



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

All Body Piercing Instruments Body Piercing Instruments


Reusable stainless-steel all body piercing instruments forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



10
PCS



10
PCS



26

PRODUCTS

\$3.00

PRICES FROM



OEM

PRIVATE LABEL

2
PCS
Global

WORLDWIDE
EXPORT

Contents

26 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

Product	Code	Price	Page
10 Pcs Ring Closing Segment Plier 6" with PVC Grip Handle	WI-70001780	\$69.99	3
10 Piercing Barbell Clamps 4.75" Bead Holder Tweezer	WI-70001798	\$49.99	3
2 Pc Ring Opening and Closing Pliers Set	WI-70002747	\$19.99	4
2 Pcs Ring Closing Segment Plier 6" with PVC Grip Handle	WI-70001782	\$17.99	4
3 Piercing Barbell Clamps 4.75" Bead Holder Tweezer	WI-70001800	\$19.99	5
5 Pcs Ring Closing Segment Plier 6" with PVC Grip Handle	WI-70001781	\$39.99	5
5 Piercing Barbell Clamps 4.75" Bead Holder Tweezer	WI-70001799	\$29.99	6
7 Pcs Body Piercing Instruments Tools Kit	WI-70002757	\$25.99	6
Body Piercing Sponge Forceps 7" Non Slotted Curved Serrated	WI-70002371	\$8.99	7
Body Piercing Sponge Forceps 7" Non Slotted Straight Serrated	WI-70002372	\$8.99	7
Body Piercing Stainless Steel Sponge Forceps Instrument Set	WI-26058427	\$3.00	8
Dermal Anchor Tube Hemostat Forceps 5" for 4mm Dermal Tops	WI-70001299	\$14.99	8
German Forester Sponge Holding Forceps Slotted Body Piercing	WI-34417582	\$8.00	9
Low MOQ Body Piercing Sponge Forceps Instrument Set	WI-27195797	\$3.00	9
Mini Forester	WI-26157271	\$3.00	10
Piercing Barbell Clamps 4.75" Bead Holder Tweezer	WI-70001801	\$9.99	10
Ring Closing Plier 5.5" Gold Jaws with Grip Handle	WI-70001303	\$10.99	11
Ring Closing Plier 5.5" with Deep Groove Tips & PVC Grip Handle	WI-70002318	\$11.24	11
Ring Closing Plier 6" Captive Bead Loop Closer with Double Groove Tips	WI-70001140	\$19.99	12
Ring Closing Segment Plier 6" with PVC Grip Handle	WI-70001783	\$10.99	12
Ring Holding & Closing Pliers 6" Gold Jaws with PVC Grip Handle	WI-70001304	\$10.99	13
Ring Opening Plier 6.25" Body Piercing Jewelry Bead Ball Capture Tool	WI-70002324	\$12.49	13
Ring Opening Plier 7" Body Piercing Jewelry Bead Ball Capture Tool	WI-70002319	\$13.74	14
Septum Forceps 6" Nasal Body Piercing Tool	WI-70001123	\$15.99	14
Stainless Steel Piercing Sponge Forceps Polish	WI-28176355	\$8.00	15
Top Notched Large Ring Closing Plier 6"	WI-70001989	\$14.99	15

**10
PCS****WI-70001780****10 Pcs Ring Closing Segment Plier 6" with PVC Grip Handle****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Stainless Steel Captive Bead Ring Closing Pliers Segment 6" with PVC Grip Handle Hoops Heavy Big Gauges Captive Bead Bar Body Piercing Tool.

The loop closing pliers is an essential tool for any jewelry maker that has to close a lot of jump rings or eye pins. This Large Captive Rings Closing Pliers have Jaws that are hollowed out with a smooth finish for closing Rings, Loops, Jump Rings, Bows, Bracelet Links, Charms Links, etc. without marring or twisting your wire. Making a repetitive job easier and quicker, loop closing pliers may be the one tool missing from your toolbox. Can be used for nose, ear, eyebrow, tongue, nipple, belly button, and lip piercings.

RING CLOSING SEGMENT PLIER FEATURES

- Loop closing pliers are primarily used to close already formed loops with a slight opening. The inner notches inside the pliers allow you to securely seed your ring or eye pin into the pliers before gently closing.
- They are autoclavable. Are used to gently close ball closure rings, bar closure rings, smooth segment rings and other similar items around their captive ball or bar. Piercing professional puncture tools.
- Open up the jaw of your loop closing pliers to seed the ring or eye pin inside the pre existing ridges. Gently place your ring into the ridges of the loop closing pliers. Gently squeeze

**10
PCS****WI-70001798****10 Piercing Barbell Clamps 4.75" Bead Holder Tweezer****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Piercing Barbell Clamps Forceps Pliers 4.75" (12cm) Surgical Steel Catch Barbell Jewelry Clamp Tweezer for Tongue Rings Labret Cartilage Earrings All Piercing.

Piercing Forceps/ Piercing Clamps Tweezer for Nose Septum Belly Nipple Tongue Eyebrow, Lips, Nose Piercing. These body piercing clamps are made of stainless steel. Useful for closing captive bead rings, these professional-grade pliers will help open rings without scratching or distorting their shape. Reusable after sterilized.

BARBELL JEWELRY CLAMP TWEEZER FEATURES

4.75" BARBELL CLAMP TWEEZER: The Piercing Barbell Clamps, used to hold

the piercing barbell jewelry. This helps to hold the barbell securely. They feature a soft easy grip for high comfort and flexibility during use.


PREMIUM QUALITY O.R. GRADE STAINLESS STEEL:

Forged from corrosion-resistant surgical stainless steel for the highest durability and product life. Smooth and tarnish-free with a high mirror finish. 100% Brand New + Autoclavable for proper sterilization. Ultrasonic cleaned + passivated.

USER FRIENDLY AND

ERGONOMICALLY DESIGNED: 100% Hygienic and easy to clean and disinfected. They are designed for precise and easy use with a high degree of flexibility during use.

DEPENDABLE BODY PIERCING TOOL: These clamps are used for a variety



**WI-70002747****2 Pc Ring Opening and Closing Pliers Set****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Set of 2 Ring Opening and Ring Closing Pliers Ball Closure Body Piercing Tattoo Jewelry Captive Bead Tools

RING OPENING AND CLOSING PLIERS FEATURES

BEAD CAPTURE RING OPENER & CLOSER SET: The Stainless Steel construction and OR quality make the Bead Capture Ring Opener & Closer Set ideal for tattoo shops and jewelry making. When you're working with sensitive new piercings, a good quality pair of opening and closing pliers is essential to cause your clients the least amount of discomfort possible. This Ring Opener and Closer Tool Set is easy to use and is built to remain dependable for body piercing professionals as well as jewelry makers

at home, thanks to the durable construction.

PREMIUM QUALITY O.R. GRADE STAINLESS STEEL: Forged from corrosion-resistant surgical stainless steel for the highest durability and product life. Smooth and tarnish free with a high mirror finish. 100% Brand New + Autoclavable for proper sterilization. Ultrasonic cleaned + passivated.

RING-OPENING PLIERS: The Ring Opening Pliers work best on rings 3/8 to 3/4 inches in diameter and up to 6 gauge in thickness. The head of the pliers has four notches to accommodate different ring diameters. It features extra grip notches in the handle for precision use Needle Nose Beak for ultimate control on big or small rings.

**WI-70001782****2 Pcs Ring Closing Segment Plier 6" with PVC Grip Handle****\$17.99** per pieceTrade price from **\$16.55**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Stainless Steel Captive Bead Ring Closing Pliers Segment 6" with PVC Grip Handle Hoops Heavy Big Gauges Captive Bead Bar Body Piercing Tool.

The loop closing pliers is an essential tool for any jewelry maker that has to close a lot of jump rings or eye pins. This Large Captive Rings Closing Pliers have Jaws that are hollowed out with a smooth finish for closing Rings, Loops, Jump Rings, Bows, Bracelet Links, Charms Links, etc. without marring or twisting your wire. Making a repetitive job easier and quicker, loop closing pliers may be the one tool missing from your toolbox. Can be used for nose, ear, eyebrow, tongue, nipple, belly button, and lip piercings.

RING CLOSING SEGMENT PLIER FEATURES

- Loop closing pliers are primarily used to close already formed loops with a slight opening. The inner notches inside the pliers allow you to securely seed your ring or eye pin into the pliers before gently closing.
- They are autoclavable. Are used to gently close ball closure rings, bar closure rings, smooth segment rings and other similar items around their captive ball or bar. Piercing professional puncture tools.
- Open up the jaw of your loop closing pliers to seed the ring or eye pin inside the pre existing ridges. Gently place your ring into the ridges of the loop closing pliers. Gently squeeze



**3
PCS****WI-70001800****3 Piercing Barbell Clamps 4.75\"/>****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Piercing Barbell Clamps
Forceps Pliers 4.75\"/>

Piercing Forceps/ Piercing Clamps
Tweezer for Nose Septum Belly Nipple
Tongue Eyebrow, Lips, Nose Piercing.
These body piercing clamps are made
of stainless steel. Useful for closing
captive bead rings, these professional-
grade pliers will help open rings without
scratching or distorting their shape.
Reusable after sterilized.

**BARBELL JEWELRY CLAMP
TWEezer FEATURES**


4.75\"/>

the piercing barbell jewelry. This helps
to hold the barbell securely. They
feature a soft easy grip for high comfort
and flexibility during use.

**PREMIUM QUALITY O.R. GRADE
STAINLESS STEEL:** Forged from
corrosion-resistant surgical stainless
steel for the highest durability and
product life. Smooth and tarnish-free
with a high mirror finish. 100% Brand
New + Autoclavable for proper
sterilization. Ultrasonic cleaned +
passivated.

**USER FRIENDLY AND
ERGONOMICALLY DESIGNED:** 100%
Hygienic and easy to clean and
disinfected. They are designed for
precise and easy use with a high
degree of flexibility during use.

DEPENDABLE BODY PIERCING TOOL:
These clamps are used for a variety

**5
PCS****WI-70001781****5 Pcs Ring Closing Segment Plier 6\"/>****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Stainless Steel Captive Bead
Ring Closing Pliers Segment 6\"/>

The loop closing pliers is an essential
tool for any jewelry maker that has to
close a lot of jump rings or eye pins.
This Large Captive Rings Closing Pliers
have Jaws that are hollowed out with a
smooth finish for closing Rings, Loops,
Jump Rings, Bows, Bracelet Links,
Charms Links, etc. without marring or
twisting your wire. Making a repetitive
job easier and quicker, loop closing
pliers may be the one tool missing from
your toolbox. Can be used for nose,
ear, eyebrow, tongue, nipple, belly
button, and lip piercings.

**RING CLOSING SEGMENT PLIER
FEATURES**

- Loop closing pliers are primarily
used to close already formed loops
with a slight opening. The inner
notches inside the pliers allow you to
securely seed your ring or eye pin
into the pliers before gently closing.
- They are autoclavable. Are used to
gently close ball closure rings, bar
closure rings, smooth segment rings
and other similar items around their
captive ball or bar. Piercing
professional puncture tools.
- Open up the jaw of your loop closing
pliers to seed the ring or eye pin
inside the pre existing ridges. Gently
place your ring into the ridges of the
loop closing pliers. Gently squeeze

LARGE RING CLOSING SEGMENT PLIER



**WI-70001799****5 Piercing Barbell Clamps 4.75\"/>****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Piercing Barbell Clamps Forceps Pliers 4.75\"/>

the piercing barbell jewelry. This helps to hold the barbell securely. They feature a soft easy grip for high comfort and flexibility during use.

PREMIUM QUALITY O.R. GRADE

STAINLESS STEEL: Forged from corrosion-resistant surgical stainless steel for the highest durability and product life. Smooth and tarnish-free with a high mirror finish. 100% Brand New + Autoclavable for proper sterilization. Ultrasonic cleaned + passivated.

USER FRIENDLY AND

ERGONOMICALLY DESIGNED: 100%

Hygienic and easy to clean and disinfected. They are designed for precise and easy use with a high degree of flexibility during use.

DEPENDABLE BODY PIERCING TOOL:

These clamps are used for a variety

Piercing Forceps/ Piercing Clamps Tweezer for Nose Septum Belly Nipple Tongue Eyebrow, Lips, Nose Piercing. These body piercing clamps are made of stainless steel. Useful for closing captive bead rings, these professional-grade pliers will help open rings without scratching or distorting their shape. Reusable after sterilized.

**BARBELL JEWELRY CLAMP
TWEezer FEATURES**

4.75\"/>

**WI-70002757****7 Pcs Body Piercing Instruments Tools Kit****\$25.99** per pieceTrade price from **\$23.91**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 7 Pcs Body Piercing Instruments Tools Kit – Premium Quality Stainless Steel Pennington Forceps, Pliers, Clamps, Ideal for Skin Jewelry, Body Piercing, Eyebrow, Belly Button, Tongue, Navel & More.


**BODY PIERCING INSTRUMENTS
KIT FEATURES****QUALITY BODY PIERCING KIT: 7**

pieces professional piercing kit is a quality value pack that comprises Pennington Forceps, Sponge Forceps, Ring Opening & Closing Plier with Kelly Forceps all you need for body piercings. This set has it all for comprehensive professional-grade piercing process instruments.

HIGH-QUALITY MATERIAL: High-quality sterilized material used for safe piercings. The presence of high-grade Surgical Stainless Steel allows for a greater level of comfort and protection against skin infections, allergic reactions, and other skin problems. All the instruments and jewelry pieces in

**7 PCS BODY PIERCING SET
INCLUDES**

- 1 Pc Pennington Forceps Unslotted 6\"/>
- 1 Pc Pennington Forceps 6\"/>
- 1 Pc Sponge Forceps 8\"/>
- 1 Pc Sponge Forceps 8\"/>
- 1 Pc Ring Closing Plier 5.5\"/>
- 1 Pc Ring-Opening Plier 6\"/>
- 1 Pc Kelly Forceps 5.5\"/>



**WI-70002371****Body Piercing Sponge Forceps 7" Non Slotted Curved Serrated****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Foerster Sponge Forceps Curved 7" (17.8cm) Serrated Jaws with the Non-Slotted end are used for a variety of body piercings. Especially for the tongue piercing.

Body Piercing Sponge Hemostat Clamps are multi-purpose forceps generally used for Body piercing, which is a form of body modification like tongue, lip, ear and naval piercings.

This one, in particular, is the non-slotted version. The closed end helps to make sure there is a firm grip on the skin and the piercing is performed safely but also accurately.


BODY PIERCING FORCEPS FEATURES.

- Sponge Forceps Non-Slotted Locking

- Oval Shape Tip.
- Serrated jaws for improved grip
- Fully autoclavable, Rust Free
- Standard non-slotted Foerster jaws
- 3 Step Adjustable ratcheting lock

PACKAGE INCLUDE:

- 1 Pc Body Piercing Sponge Forceps Non-Slotted 7" Curved Serrated Jaws

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Body Piercing Tools**APPLICATION:** Body Jewelry**TIP CONFIGURATION:** Oval Shape Fenestrated**TIP JAWS:** Serrated**CURVATURE:** Curved**HANDLE:** Ratcheted Finger Ring Handle**WI-70002372****Body Piercing Sponge Forceps 7" Non Slotted Straight Serrated****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Foerster Sponge Forceps Straight 7" (17.8cm) Serrated Jaws with the Non-Slotted end are used for a variety of body piercings. Especially for the tongue piercing.

Body Piercing Sponge Hemostat Clamps are multi-purpose forceps generally used for Body piercing, which is a form of body modification like tongue, lip, ear and naval piercings.

This one, in particular, is the non-slotted version. The closed end helps to make sure there is a firm grip on the skin and the piercing is performed safely but also accurately.



BODY PIERCING FORCEPS FEATURES.

- Sponge Forcep Non-Slotted Locking
- Oval Shape Tip.

- Serrated jaws for improved grip
- Fully autoclavable, Rust Free
- Standard non-slotted Foerster jaws
- 3 Step Adjustable ratcheting lock

PACKAGE INCLUDE:

- 1 Pc Body Piercing Sponge Forceps Non-Slotted 7" Straight Serrated Jaws

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Body Piercing Tools**APPLICATION:** Body Jewelry**TIP CONFIGURATION:** Oval Shape Fenestrated**TIP JAWS:** Serrated**CURVATURE:** Straight**HANDLE:** Ratcheted Finger Ring Handle**REUSABLE:** Yes



WI-26058427

Body Piercing Stainless Steel Sponge Forceps Instrument Set

\$3.00 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Body Piercing Stainless Steel Sponge Forceps Instrument Set, made for professional piercing studios. Instruments of this type hold and align tissue accurately during placement. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.00, falling to \$1.50 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001299

Dermal Anchor Tube Hemostat Forceps 5" for 4mm Dermal Tops

\$14.99 per piece

Trade price from \$13.79
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dermal Anchor Tube Hemostat Forceps 5" for Dermal Tops German Grade Stainless Steel Anchor Holder Body Piercing Tools for 4mm Dermal Tops.

great with this tool. Made from high quality stainless steel designed with maximum comfortable enough to use hundreds of times.

These dermal anchor forceps are perfect for holding dermal anchors and piercing needles. They can also be used for handling small beads with care. This micro dermal anchor piercing tool has been designed for easy insertion of micro dermal anchors or surface anchors. This Body piercing tool allows you to lock the dermal anchor in and easily unlock it once you have finished manipulating the insertion. Opening in these forceps is 4mm wide so you can use this with any dermal anchor top that is 4mm or larger. That means that a disc or most other dermal threaded tops will work

DERMAL ANCHOR FORCEPS FEATURES

- Dermal Anchor Tube Hemostat Forceps for Inserting Dermal Anchors.
- Size 5" Overall Length.
- For 4mm Dermal Tops.
- Holds Dermal Top Firmly for Placement.
- Ratchet Locking Mechanism.
- Finger Ring Handle
- Reusable.Perfect for Holding Needles, or Anything That is Bigger Than 3mm Wide.





Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34417582

German Forester Sponge Holding Forceps Slotted Body Piercing

\$8.00 per piece

Trade price from \$4.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of German Forester Sponge Holding Forceps Slotted Body Piercing, made for professional piercing studios. Instruments of this type hold and align tissue accurately during placement. The steel is heat treated and passivated for corrosion resistance, then each instrument is

checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-27195797

Low MOQ Body Piercing Sponge Forceps Instrument Set

\$3.00 per piece

Trade price from \$1.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Low MOQ Body Piercing Sponge Forceps Instrument Set — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type hold and align tissue accurately during placement. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$3.00, falling to \$1.50 at 1 piece — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.



WI-26157271

Mini Forester

\$3.00 per piece

Trade price from **\$1.50**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Mini Forester for professional piercing studios — made from bar stock to finished instrument in house. Instruments of this type hold and align tissue accurately during placement.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.00, falling to \$1.50 at 1 piece — about 50% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001801

Piercing Barbell Clamps 4.75" Bead Holder Tweezer

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Piercing Barbell Clamps Forceps Pliers 4.75" (12cm) Surgical Steel Catch Barbell Jewelry Clamp Tweezer for Tongue Rings Labret Cartilage Earrings All Piercing.

Piercing Forceps/ Piercing Clamps Tweezer for Nose Septum Belly Nipple Tongue Eyebrow, Lips, Nose Piercing. These body piercing clamps are made of stainless steel. Useful for closing captive bead rings, these professional-grade pliers will help open rings without scratching or distorting their shape. Reusable after sterilized.

the piercing barbell jewelry. This helps to hold the barbell securely. They feature a soft easy grip for high comfort and flexibility during use.

PREMIUM QUALITY O.R. GRADE STAINLESS STEEL:

Forged from corrosion-resistant surgical stainless steel for the highest durability and product life. Smooth and tarnish-free with a high mirror finish. 100% Brand New + Autoclavable for proper sterilization. Ultrasonic cleaned + passivated.

USER FRIENDLY AND

ERGONOMICALLY DESIGNED: 100% Hygienic and easy to clean and disinfected. They are designed for precise and easy use with a high degree of flexibility during use.

DEPENDABLE BODY PIERCING TOOL: These clamps are used for a variety



BARBELL JEWELRY CLAMP TWEEZER FEATURES

4.75" BARBELL CLAMP TWEEZER: The Piercing Barbell Clamps, used to hold

**WI-70001303**

Ring Closing Plier 5.5" Gold Jaws with Grip Handle

\$10.99 per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Ring Closing Plier 5.5 Inches Gold Jaws with Blue PVC Grip Handle Jump Ring Loop Closing Jewelry Tools German Grade Stainless Steel Body Piercing Jewelry Making Repair DIY Tools

Ring closing pliers are specially designed with grooves in the matching faces of the tips. This helps to hold the curve of a ring securely. They feature a soft easy grip for high comfort and flexibility during use. Loop closing plier is an essential tool for any jewelry maker that has to close a lot of jump rings or eye pins. Ring closing pliers have gold jaws that are hollowed out with a smooth finish for closing rings loops, Jump rings loop, Bows, Bracelet links, Charms links, etc.

GOLD RING CLOSING PLIER: When you're a professional body piercing artist working with sensitive new piercings, a good quality of plier helpful tool to ensure you cause your clients the least amount of discomfort possible. This ring closer plier is designed for piercing rings in a professional, sterile environment. They can also come in handy for jewelers looking for a dependable tool for closing jump rings of many sizes.

HIGH QUALITY STAINLESS STEEL: Forged from corrosion resistant Premium Quality Stainless Steel for the highest amount of durability and product life. Tarnish free with a high mirror finish + ultrasonic cleaned. These pliers will not bend or break on you and will last for many years of service.

RING CLOSING PLIER FEATURES

**WI-70002318**

Ring Closing Plier 5.5" with Deep Groove Tips & PVC Grip Handle

\$11.24 per pieceTrade price from **\$10.34**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Quality 5.5" Loop Closing Plier Close Jump Rings Jewelry Making Tool Beading With Red PVC Grip Handle. Ideal for Multiple Ring Closing Usage Navel Ring/Tongue Ring/Lip Ring and Many other professional uses.

with a high mirror finish. 100% Brand New + Autoclavable for proper sterilization. Ultrasonic cleaned + passivated.

USER FRIENDLY AND

ERGONOMICALLY DESIGNED: 100% Hygienic and easy to clean and disinfected. They are designed for precise and easy use with a high degree of flexibility during use.

DEPENDABLE BODY PIERCING TOOL: Loop closing plier with deeply grooved, smoothly finished jaws. Used to close loops, jump rings, and bracelet links without scratching or distorting. 5.5" total length, Box joint construction for smooth functioning. They feature a soft easy grip PVC handle for high comfort and flexibility during use.

SATISFACTION GUARANTEED: These high-quality professional grade jewelry

RING CLOSING PLIERS FEATURES

5.5" RING CLOSING PLIERS: The ring-closing pliers are specially designed with grooves in the matching faces of the tips. This helps to hold the curve of a ring securely. They feature a soft easy grip for high comfort and flexibility during use.

PREMIUM QUALITY O.R. GRADE STAINLESS STEEL: Forged from corrosion-resistant surgical stainless steel for the highest durability and product life. Smooth and tarnish-free



**WI-70001140****Ring Closing Plier 6" Captive Bead Loop Closer with Double Groove Tips****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Ring Closing Plier 6" (15.2cm) with Grip Handle – Captive Bead Loop and Ring Closer Hoop Jewelry Tools German Grade Stainless Steel Needle Nose Beak Body Piercing Jewelry Making Repair DIY Tools.

This ring closing plier feature hollowed out and smooth jaws with needle nose beak used to close jump rings, piercing rings, captive beads loops etc. Box joint style ring closing plier ensure comfortable holding with ease of open and close. Stainless steel non-slip grip handle designed for a sterile environment, this plier is fully autoclavable without fear of degradation, ensuring 100% hygiene and easy cleaning, making it a dependable tool for professional body piercing artists and jewelers.

CAPTIVE BEADS RING CLOSING PLIER FEATURES

- Captive Bead Capture & Large Ring Clothing Plier
- Overall Length 6" (15.2cm)
- Twice Hollowed Out and Smooth Jaws
- Needle Nose Beak
- Box Joint Style
- Non-slip Grip Handle for Strong Grip While in Use
- Reusable and Fully Autoclavable
- Made With German Grade Stainless Steel

PACKAGE INCLUDE

- 1 Pc Captive Bead Ring Closing Plier 6" (15.2cm) with Twice Deep Groove

**WI-70001783****Ring Closing Segment Plier 6" with PVC Grip Handle****\$10.99** per pieceTrade price from **\$10.11**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

6" Large Stainless Steel Captive Bead Ring Closing Pliers Segment with PVC Grip Handle Hoops Heavy Big Gauges Captive Bead Bar Body Piercing Tool.

The loop closing pliers is an essential tool for any jewelry maker that has to close a lot of jump rings or eye pins. This Large Captive Rings Closing Pliers have Jaws that are hollowed out with a smooth finish for closing Rings, Loops, Jump Rings, Bows, Bracelet Links, Charms Links, etc. without marring or twisting your wire. Making a repetitive job easier and quicker, loop closing pliers may be the one tool missing from your toolbox. Can be used for nose, ear, eyebrow, tongue, nipple, belly button, and lip piercings.

RING CLOSING SEGMENT PLIER FEATURES

- Loop closing pliers are primarily used to close already formed loops with a slight opening. The inner notches inside the pliers allow you to securely seed your ring or eye pin into the pliers before gently closing.
- They are autoclavable. Are used to gently close ball closure rings, bar closure rings, smooth segment rings and other similar items around their captive ball or bar. Piercing professional puncture tools.
- Open up the jaw of your loop closing pliers to seed the ring or eye pin inside the pre existing ridges. Gently place your ring into the ridges of the loop closing pliers. Gently squeeze



**WI-70001304**

Ring Holding & Closing Pliers 6" Gold Jaws with PVC Grip Handle

\$10.99 per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Ring Holding & Closing Pliers with Gold Jaws 6" Captive Bead Ring Closing Segment Pliers with Blue PVC Grip Handle Hoops Heavy Big Gauges Captive Bead Bar Body Piercing Stainless Steel Jewelry Making Repair DIY Tools

The loop closing pliers is an essential tool for any jewelry maker that has to close a lot of jump rings or eye pins. This Large Captive Rings Closing Pliers have Jaws that are hollowed out with a smooth finish for closing Rings, Loops, Jump Rings, Bows, Bracelet Links, Charms Links, etc. without marring or twisting your wire. Making a repetitive job easier and quicker, loop closing pliers may be the one tool missing from your toolbox. Can be used for nose,

ear, eyebrow, tongue, nipple, belly button, and lip piercings.

RING HOLDING & CLOSING PLIERS FEATURES

- These pliers are designed for placing and closing rings of any size or shape. Because the jaws come to a point and are nicely grooved, you can close rings without scratching or marring. These pliers also allow the user to get into tight places and hold the ring for soldering.
- Loop closing pliers are primarily used to close already formed loops with a slight opening. The inner notches inside the pliers allow you to securely seed your ring or eye pin into the pliers before gently closing.

**WI-70002324**

Ring Opening Plier 6.25" Body Piercing Jewelry Bead Ball Capture Tool

\$12.49 per pieceTrade price from **\$11.49**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Stainless steel 6.25" (15.87cm) ring opening plier for captive bead rings up to 3/4" and 6g. A must have for anyone who likes to change jewelry or even just to change the bead to another color or style. Ideal Loop opening pliers jump ring opener captive bead body piercing hobby craft jewelry making tool

PREMIUM QUALITY O.R. GRADE

STAINLESS STEEL: Forged from corrosion-resistant surgical stainless steel for the highest durability and product life. Smooth and tarnish free with a high mirror finish. 100% Brand New + Autoclavable for proper sterilization. Ultrasonic cleaned + passivated.

USER FRIENDLY AND

ERGONOMICALLY DESIGNED: 100% Hygienic and easy to clean, the Ring Opener is fully autoclavable without fear of degradation. You can sterilize them without fear. They are designed for precise and easy use with a high degree of flexibility during use.

DEPENDABLE BODY PIERCING TOOL:

The Ring Opening Plier have a needle nose beak for both big and small rings. The head of the pliers has two notches

RING OPENING PLIER FEATURES

RING-OPENING PLIERS: The Ring Opening Pliers work best on rings 3/8 to 3/4 inches in diameter and up to 6 gauge in thickness. The head of the pliers has 2 notches to accommodate different ring diameters. It features extra grip notches in the handle for precision use Needle Nose Beak for ultimate control on big or small rings.



**WI-70002319**

Ring Opening Plier 7" Body Piercing Jewelry Bead Ball Capture Tool

\$13.74 per pieceTrade price from **\$12.64**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Stainless steel 7" (17.8cm) ring opening plier for captive bead rings up to 1/2" and 6g diameter. A must have for anyone who likes to change jewelry or even just to change the bead to another color or style. Ideal Loop opening pliers jump ring opener captive bead body piercing hobby craft jewelry making tool

Needle Nose Beak for ultimate control on big or small rings.

PREMIUM QUALITY O.R. GRADE STAINLESS STEEL: Forged from corrosion-resistant surgical stainless steel for the highest durability and product life. Smooth and tarnish free with a high mirror finish. 100% Brand New + Autoclavable for proper sterilization. Ultrasonic cleaned + passivated.

USER FRIENDLY AND ERGONOMICALLY DESIGNED: 100% Hygienic and easy to clean, the Ring Opener is fully autoclavable without fear of degradation. You can sterilize them without fear. They are designed for precise and easy use with a high degree of flexibility during use.

DEPENDABLE BODY PIERCING TOOL: The Ring Opening Plier have a needle

RING OPENING PLIER FEATURES

RING-OPENING PLIERS: The Ring Opening Pliers works best with ring diameters 1/2" (13mm) up to 40mm in diameter and up to 6 gauge in thickness. The head of the pliers has 2 notches to accommodate different ring diameters. It features extra grip notches in the handle for precision use

**WI-70001123**

Septum Forceps 6" Nasal Body Piercing Tool

\$15.99 per pieceTrade price from **\$14.71**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Septum Forceps 6" (15.2cm) Stainless Steel Nose Body Piercing Septum Forceps Clamp Pliers with 5mm O.D. Tube Working Ends for Ear Lip Belly Navel Tongue. Great for tragus and other cartilage piercings.

piercing professionals. 3-step locking mechanism and ring handles conforms easy to grip and bring convenience to operate. We are using high quality stainless steel ensures maximum durability and dependability.

Septum piercing forceps allows the needle to pass through straight with ease. Septum forceps with tubes working ends that are large enough to allow a 10g needle to pass through with ease. The ergonomic design for septum piercings easy to hold and use. The precise nib can hold small objects firmly, and can also be used as a fixture for other purposes, which is hygienic and easy to grasp, such as sewing needles, soybeans, mung beans, etc. Septum forceps are delicately designed for both home DIY use and body

SEPTUM FORCEPS FEATURES.

- 5mm O.D Tube Shape Working Tips.
- 6" (15.2cm) Overall Length
- Box Joint
- 3 Step Adjustable ratcheting lock
- Fully autoclavable, Rust Free
- Finger Ring Handles
- Made with High Quality Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Septum Forceps 6" (15.2cm) with Ratchet Lock Finger Ring Handle





WI-28176355

Stainless Steel Piercing Sponge Forceps Polish

\$8.00 per piece

Trade price from \$4.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Piercing Sponge Forceps Polish is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type hold and align tissue accurately during placement.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001989

Top Notched Large Ring Closing Plier 6"

\$14.99 per piece

Trade price from \$13.79
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Top Notched Ring Closing Pliers 6" (15.2cm) with Stainless Steel Grip Handle. Ideal closing & wire bending pliers for captive beads body piercing hobby craft jewelry making tool

These large ring closing plier help to tighten captive rings with a notched area designed to prevent scratching the jewelry in the process. To close easily all kind of rings of any thickness and also helps to secure beads in captives. The large closing & wire bending pliers work best for captives from 5/8" and up. This tool is useful for Rings 20g – 6g.

TOP NOTCHED LARGE RING CLOSING PLIERS: The closing pliers are specially designed with deep grooves and captive top notch. A perfect tool to close easily all kind of rings of any thickness.

PREMIUM QUALITY O.R. GRADE STAINLESS STEEL: Forged from corrosion-resistant surgical stainless steel for the highest durability and product life. Smooth and tarnish-free with a high mirror finish. 100% Brand New + Autoclavable for proper sterilization. Ultrasonic cleaned + passivated.

USER FRIENDLY AND ERGONOMICALLY DESIGNED: 100% Hygienic and easy to clean and disinfected. They are designed for precise and easy use with a high degree of flexibility during use.

NOTCHED RING CLOSING PLIERS FEATURES





WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

Request a quotation

Send us product codes from this catalogue, quantities and your country.

We reply the same working day on WhatsApp.

Online shop	shop.waqarintl.com
Email	sales@waqarintl.com
WhatsApp /	
Phone	+92 307 6173627
Works	Roras Road, Opposite Anwar Khawaja Industries, Sialkot 51310, Pakistan

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.