



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

Elevators & Curettes Dental Instruments

Reusable stainless-steel elevators & curettes forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



95

PRODUCTS

\$3.85

PRICES FROM

OEM

PRIVATE LABEL

Global

WORLDWIDE
EXPORT

Contents

95 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

Product	Code	Price	Page
10 Piece Dental Luxating Elevator Set	WI-42900435	\$31.20	5
10 Seldin Root Elevator 34S – 4mm Tip	WI-70002147	\$62.49	5
12 Periosteal Elevator Seldin #S23 Double Ended	WI-70003904	\$76.83	6
13 Pcs Set Stainless Steel Curved Dental Elevator	WI-25639209	\$4.00	6
2 Apical Elevator #302 Left Angled 3mm Tip	WI-70002137	\$12.99	7
2 Apical Elevator 303 – Right	WI-70002136	\$12.99	7
2 Apical Root Elevator 301 Straight 2.5mm Tip	WI-70002140	\$12.99	8
2 Seldin Elevator 4R – 10mm Tip	WI-70002145	\$12.99	8
2 Seldin Root Elevator 34S – 4mm Tip	WI-70002149	\$12.99	9
3 PCS Dental Extraction Elevators Set	WI-42895925	\$13.00	9
3 PCS Dental Straight Elevator Kit	WI-42894919	\$13.00	10
3 PCS Dental Straight Extraction Elevator Set	WI-42912655	\$13.00	10
3 Piece Dental Straight Extraction Elevator Set	WI-42918470	\$13.00	11
4R Seldin Elevator – 10mm Tip	WI-70002146	\$9.99	11
5 Seldin Root Elevator 34S – 4mm Tip	WI-70002148	\$31.24	12
6 Periosteal Elevator Seldin #S23 Double Ended	WI-70003899	\$41.99	12
Apical Dental Elevator #46 Straight 3.5mm Tip	WI-70002142	\$9.99	13
Apical Elevator #302 Left Angled 3mm Tip	WI-70002139	\$9.99	13
Apical Elevator #303 – Right	WI-70002138	\$9.99	14
Apical Elevator #303HF Serrated Tip	WI-70002135	\$9.99	14
Apical Root Elevators 301 Straight 2.5mm Small Tip	WI-70003897	\$9.99	15
Berenard Curved Elevator Dental Root Extraction Tool	WI-42896583	\$7.00	15
Bernard Spade Root Tip Elevator Stainless Steel Curved Dental	WI-42892597	\$7.00	16
Cogswell A Root Elevator 6.25" Straight 4mm Tip	WI-70002141	\$9.99	16
Cryer #44 Dental Root Elevator	WI-70000008	\$9.99	17
Cryer #45 Dental Root Elevator	WI-70003885	\$9.99	17
Cryer Root Elevator #25	WI-70000005	\$9.99	18
Cryer Root Elevator #26	WI-70000004	\$9.99	18
Curved Dental Luxating Elevator Stainless Steel Oral Surgery Tool	WI-42928029	\$7.00	19
Dental Berenard Spade Elevator Concave Curved Root Extraction Tool	WI-42913338	\$7.00	19
Dental Bone Curettes Curettage	WI-70004107	\$20.53	20
Dental Crane Pick Elevator #8	WI-70003886	\$9.99	20
Dental Curved Luxating Elevator Root Extraction Instrument	WI-42910772	\$7.00	21
Dental Elevator #3	WI-70000006	\$9.99	21
Dental Elevator #81	WI-70000007	\$9.99	22
Dental Elevators	WI-32893707	\$4.00	22
Dental Extraction Spear Root Elevators	WI-01670769	\$3.85	23
Dental Heidbrink Root Elevator Stainless Steel Oral Surgery Tool	WI-42917335	\$6.00	23
Dental Luxating Elevator Curved Stainless Steel Tooth Extraction	WI-42905879	\$7.00	24
Dental Luxation Kit Surgical Luxating Tooth Extraction Root Elevators	WI-25658193	\$4.00	24
Dental Luxation Root Elevators 6-Piece	WI-15381460	\$12.00	25

Contents (continued)

Product	Code	Price	Page
Dental PDL Luxating Root Elevators Set	WI-30367248	\$25.00	26
Dental Periodontal Surgical Stainless Steel Bone Curettes (Lucas	WI-27580627	\$6.00	26
Dental Periosteal Elevator Molt M9 Blue Titanium Instrument	WI-35370689	\$7.00	27
Dental Root Elevator Set 5 Pcs	WI-32940520	\$4.00	27
Dental Root Elevator Stainless Steel Dental Instruments	WI-09976073	\$4.50	28
Dental Root Elevator Stainless Steel Surgical	WI-32963538	\$4.00	28
Dental Root Elevator Surgical Instrument Orthodontic	WI-25812497	\$4.50	29
Dental Root Elevators Set 10PCS Stainless Steel Extraction Tools	WI-42880564	\$31.00	29
Dental Root Planning Lab Stainless Steel Handle Surgical	WI-32277655	\$4.00	30
Dental Root Tip Elevator DEL 303 Stainless Steel Tooth	WI-42919101	\$7.00	30
Dental Root Tip Elevator DEL 303 Surgical Instrument	WI-42874818	\$7.00	31
Dental Scaler Stainless Steel Periodontal Curette Surgical Instrument	WI-44789087	\$4.00	31
Dental Straight Elevators Set 3PCS Tooth Extraction Surgical Tools	WI-42909751	\$13.00	32
Dental Tooth Extraction Elevator Apical Root Tip Extractor Luxating	WI-32960530	\$4.00	32
Double-Ended Dental Hygienist Sickle Scaler Curettes Machinery	WI-36929137	\$8.00	33
Heidbrink Elevator Dental Root Extraction Instrument	WI-42902709	\$6.00	33
Heidbrink Root Elevator #2 Right Angled Slightly Curved Pointed Tip	WI-70003895	\$9.99	34
Heidbrink Root Elevator 1 Straight Slightly Curved Pointed Tip	WI-70003896	\$9.99	34
Heidbrink Root Elevator 3 Left Angled Slightly Curved Pointed Tip	WI-70003894	\$9.99	35
Heidbrink Root Tip Elevator H1 Straight Hollow Handle Dental	WI-42904604	\$7.00	35
Heidbrink Root Tip Elevator Stainless Steel Dental Extraction Instrument	WI-42892741	\$6.00	36
Heidbrink Root Tip Pick Elevator H1 Straight Pointed Tip	WI-70000644	\$8.99	36
Heidbrink Root Tip Pick Elevator H2 Right Angled Pointed Tip	WI-70000643	\$8.99	37
Heidbrink Root Tip Pick Elevator H3 Left Angled Pointed Tip	WI-70000642	\$8.99	37
Heidbrink Straight Elevator Dental Oral Surgery Tool	WI-42880759	\$7.00	38
Light Weight Stainless Steel Dental Elevators Ergonomic Grip 3	WI-32966532	\$4.00	38
Molt Bone Curette #2/4 Double Ended	WI-70000540	\$9.99	39
Molt Bone Curette #2/4 Double Ended Hollow Handle	WI-70000227	\$9.99	39
Molt Bone Curette #6R Single Ended 2.9mm Spoon Tip	WI-70000541	\$10.99	40
Molt M9 Dental Periosteal Elevator Stainless Steel Surgical Tool	WI-35292882	\$7.00	40
Offset Dental Crown Splitter Remover 5.75" Elevator Handle	WI-70002024	\$12.99	41
Orthodontic Surgical Dental Bone Curettes Double Ended	WI-27587648	\$6.00	41
Periosteal Elevator Baby Molt #7A Double Ended	WI-70000705	\$7.99	42
Periosteal Elevator Baby Molt #7A Double Ended with Round Handle	WI-70003892	\$8.99	42
Prichard Periosteal Elevator PR-3 with Blue Titanium Coated Tips	WI-70000327	\$10.99	43
Reusable Dental Root Tip Elevator Stainless Steel Oral Surgery	WI-42881674	\$7.00	43
Root Elevator Stainless Steel Dental Instrument	WI-25818494	\$4.50	44
Root Tip Elevator Stainless Steel	WI-42919007	\$7.00	44
Seldin Elevator 4L – 10mm Tip	WI-70002144	\$9.99	45
Seldin Periosteal Elevator #S23 8" Double Ended	WI-70003893	\$9.99	45
Seldin Root Elevator 34 – 5mm Tip	WI-70003890	\$9.99	46
Seldin Root Elevator 34S – 4mm Tip	WI-70003240	\$9.99	46

Contents (continued)

Product	Code	Price	Page
Set 2 German Stainless Steel Dental Elevators	WI-32967470	\$4.00	47
Set 2 German Stainless Steel Hollow Handle Dental	WI-32974535	\$4.00	47
Stainless Steel DEL 303 Dental Oral Surgery Elevator	WI-42911488	\$7.00	48
Stainless Steel Dental Elevators Set 2 German Ergonomic	WI-32931581	\$4.00	48
Stainless Steel Dental Elevators Set Root	WI-32955561	\$4.00	49
Stainless Steel Dental Scaler Gracey Curette	WI-32290513	\$4.00	49
Stainless Steel Dental Scalers & Bone Curettes	WI-27606659	\$6.00	50
Stainless Steel Dental Tooth Extraction Elevator Apical Root	WI-25722160	\$4.00	50
Surgical Dental Elevators Stainless Steel Ergonomic Grip	WI-31557272	\$4.00	51
Tip Dental Bone Curettes Periodontal Surgery 2.5mm/3.5mm/4mm	WI-27573589	\$6.00	51
Volkman Bone Curettes	WI-27581634	\$6.00	52



WI-42900435

10 Piece Dental Luxating Elevator Set

\$31.20 per piece

Trade price from **\$26.50**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

10 Piece Dental Luxating Elevator Set is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$31.20, falling to \$26.50 at 1 set — about 15% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002147

10 Seldin Root Elevator 34S – 4mm Tip

\$62.49 per piece

Trade price from **\$57.49**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Pieces Seldin Elevator #34S Dental Extraction Surgery Instruments

Seldin Root Elevator 34S Straight with 4mm Tip used during extractions to loosen the tooth from the gums to ease extraction and prevent trauma to surrounding teeth and tissue. The Seldin elevator #34S is used for upper and lower anterior roots. It helps professionals in easy extraction with its multiple features.

SELDIN ROOT ELEVATOR 34S FEATURES

- 4mm sharp straight tip for upper & lower anterior roots.
- Long shaft and a very thin, rounded, slightly hollow tip.
- Ergonomic Handle with thumb hold for increased comfort.



- Superior quality Stainless Steel construction.

PACKAGE INCLUDE:

- 10 Pcs Seldin Root Elevator #34S Straight – 4mm Tip

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Root Elevators
APPLICATION: Dental Extraction
BLADE TIP SIZE: 4mm
CURVATURE: Straight
HANDLE: Octagonal hollow handle
REUSABLE: Yes
AUTOCCLAVED: Yes
RUST PREVENTION PROCEDURE: Passivated
ULTRASONIC CLEANED: Yes
TESTS PERFORMED: Boil Test,





WI-70003904

12 Periosteal Elevator Seldin #S23 Double Ended

\$76.83 per piece

Trade price from **\$70.68**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs #S23 Periosteal Elevator Seldin 8" (20.3cm) Double Ended Large Flap Retractor / Small Periosteal 10mm & 13mm wide Slightly Curved Blade with Large Flat Handle Quality Stainless Steel Dental Surgical Instruments

Seldin Periosteal Elevator #S23 are designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. Seldin S23 periosteal elevator designed with a large flat handle with double ended design approximately 10mm and 13mm wide blades. The smaller end is sharp for tissue reflecting, and the larger end is slightly curved and is a non-cutting blade, used for tissue retracting. Periosteal elevator S23 manufactured from high grade stainless steel to exacting specifications to ensure

optimal performance. Seldin elevators is autoclavable and crafted from high-quality stainless steel ensuring durability and sustained performance. 100% satisfaction guaranteed

PERIOSTEAL ELEVATOR SELDIN #S23 FEATURES

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Broad Rounded End & Broad Sharp Chisel End.
- Slight Curved 10mm and 13mm Wide Blades.
- Large Flat Handle.
- Used to Elevate the Periosteum
- Reusable.
- Made With Premium Quality Stainless Steel



WI-25639209

13 Pcs Set Stainless Steel Curved Dental Elevator

\$4.00 per piece

Trade price from **\$2.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture 13 Pcs Set Stainless Steel Curved Dental Elevator for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a curved working end for reach around structures. Supplied as a 13-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$4.00, falling to \$2.00 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70002137****2 Apical Elevator #302 Left Angled 3mm Tip****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Pcs Apical Elevator Set #302 6" (15.2cm) Long Left Angled 3mm Tip with Standard Handle Premium Quality Stainless Steel Dental Extraction Surgery Instruments



Apical Root Elevator #302 Used during extractions to loosen the tooth from the gums to ease extraction and prevent trauma to surrounding teeth and tissue. It works as a lever pushed between the gums and the tooth to loosen the surrounding tissues. Manufactured from surgical grade stainless steel provides a high degree of flexibility while conducting clinical procedures. Apical Elevator comes with highly polished finish resistance and anti-oxidation and can be autoclaved.

APICAL ELEVATOR #302 FEATURES

- Overall Length 6" (15.2cm).
- 3mm Left Angled Tip for Broken or Deep-seated Roots.
- Lightweight Stainless Steel Octagonal Handle.
- Thumb Hold for Increased Comfort.
- Fully Autoclavable & Reusable.
- Superior Quality Stainless Steel Construction.

PACKAGE INCLUDE:

- 2 Pcs Apical Root Elevator 302 Left – 3mm Tip – Standard Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 3mm**TIP ANGLE:** Left Angled**WI-70002136****2 Apical Elevator 303 – Right****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Dental Elevator Apical #303 – Right Angled 3mm Tip Dental Extraction Surgery Instruments

Apical Root Elevator #303 Used during extractions to loosen the tooth from the gums to ease extraction and prevent trauma to surrounding teeth and tissue. It works as a lever pushed between the gums and the tooth to loosen the surrounding tissues.

APICAL ELEVATOR 303 FEATURES

- 3mm right angled tip for broken or deep-seated roots
- Lightweight stainless steel handle
- Thumb hold for increased comfort
- Superior quality Stainless Steel construction

PACKAGE INCLUDE:



- 2 Pcs Apical Root Elevator 303 Right – 3mm Tip – Standard Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 3mm**TIP ANGLE:** Right Angled**HANDLE:** Octagonal hollow handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test,

Performance Test, Shape Test

QC PASSED: Yes**PACKING DETAILS:** Polybag Packed

**2
PCS****WI-70002140****2 Apical Root Elevator 301 Straight 2.5mm Tip****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Apexo Dental Elevator Apical 301 Small Tip 2.5mm Tip Surgical Denture Instruments

Apical Root Elevator 301 is designed for tooth extraction by elevating the periosteum layer. 301 Apexo Small straight elevator with a tip measuring 2.5mm is frequently used for beginning the luxation of an erupted tooth, before application of the forceps.

APICAL ELEVATOR 301 FEATURES

- 2.5mm Straight, small, smooth tip for deep-seated roots
- Lightweight stainless steel handle
- Thumb hold for increased comfort
- Superior quality Stainless Steel construction

PACKAGE INCLUDE:


- 2 Pcs Apical Root Elevator 301 Straight – 2.5mm Tip – Standard Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 3mm**TIP ANGLE:** Slight Angled**HANDLE:** Octagonal hollow handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test,

Performance Test, Shape Test

QC PASSED: Yes**PACKING DETAILS:** Polybag Packed**2
PCS****WI-70002145****2 Seldin Elevator 4R – 10mm Tip****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Root Elevator Seldin #4R Dental Extraction Surgery Instruments

Seldin Elevator 4R is used for pulling up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction. 4R elevator is a smaller version of 1R. Seldin Elevators used to remove soft tissue and periosteal membrane from the alveolar bone with pushing strokes. It has one working end, a robust shank, and a sturdy handle. Sharp tips helping in extraction procedure. These can be used for both the upper and lower molar roots.

SELDIN ELEVATOR 4R FEATURES

- 10mm sharp right angled tip
- Right upper & lower molar roots
- Smaller version of 1R


- Ergonomic standard handle with thumb hold for increased comfort.
- Superior quality Stainless Steel construction.

PACKAGE INCLUDE:

- 2 Pieces 4R Seldin Elevator 10mm Tip Standard Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 10mm**TIP ANGLE:** Right**HANDLE:** Octagonal hollow handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:**

Passivated





WI-70002149

2 Seldin Root Elevator 34S – 4mm Tip

\$12.99 per piece

Trade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Seldin Root Elevator 34S Straight with 4mm Tip used during extractions to loosen the tooth from the gums to ease extraction and prevent trauma to surrounding teeth and tissue. The Seldin elevator #34S is used for upper and lower anterior roots. It helps professionals in easy extraction with its multiple features.

PACKAGE INCLUDE:

- 2 Pcs Seldin Root Elevator #34S Straight – 4mm Tip

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Root Elevators

APPLICATION: Dental Extraction

BLADE TIP SIZE: 4mm

CURVATURE: Straight

HANDLE: Octagonal hollow handle

REUSABLE: Yes

AUTOCCLAVED: Yes



RUST PREVENTION PROCEDURE: Passivated

ULTRASONIC CLEANED: Yes

TESTS PERFORMED: Boil Test, Performance Test, Shape Test

QC PASSED: Yes

PACKING DETAILS: Polybag Packed



SELDIN ROOT ELEVATOR 34S FEATURES

- 4mm sharp straight tip for upper & lower anterior roots.
- Long shaft and a very thin, rounded, slightly hollow tip.
- Ergonomic Handle with thumb hold for increased comfort.
- Superior quality Stainless Steel construction.



WI-42895925

3 PCS Dental Extraction Elevators Set

\$13.00 per piece

Trade price from **\$6.00**

Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture 3 PCS Dental Extraction Elevators Set for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. Supplied as a 3-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$13.00, falling to \$6.00 at 1 set — about 54% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-42894919

3 PCS Dental Straight Elevator Kit

\$13.00 per piece

Trade price from \$6.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture 3 PCS Dental Straight Elevator Kit for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a straight working end for direct, in-line access.

Supplied as a set — the contents can be adjusted to your list. We work the

instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$13.00, falling to \$6.00 at 1 set — about 54% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-42912655

3 PCS Dental Straight Extraction Elevator Set

\$13.00 per piece

Trade price from \$6.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

This is our pattern of 3 PCS Dental Straight Extraction Elevator Set, made for dental, hygiene and oral surgery work. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a straight working end for direct, in-line access.

Supplied as a 3-piece set. We can add, drop or substitute items so the tray

matches your protocol rather than ours. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge.

Listed at \$13.00, falling to \$6.00 at 1 set — about 54% off the single-unit price. Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



WI-42918470

3 Piece Dental Straight Extraction Elevator Set

\$13.00 per piece

Trade price from \$6.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

3 Piece Dental Straight Extraction Elevator Set is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a straight working end for direct, in-line access.

Supplied as a set — the contents can be adjusted to your list. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$13.00, falling to \$6.00 at 1 set — about 54% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002146

4R Seldin Elevator – 10mm Tip

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Root Elevator Seldin #4R Right Angled 10mm Tip with Standard Handle Premium Quality Stainless Steel Dental Extraction Surgery Instruments.

Seldin Elevator 4R is used for pulling up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction. 4R elevator is a smaller version of 1R. Seldin Elevators used to removes soft tissue and periosteal membrane from the alveolar bone with pushing strokes. It has one working end, a robust shank, and a sturdy handle. Sharp tips helping in extraction procedure. These can be used for both the upper and lower molar roots.

SELDIN ELEVATOR 4R FEATURES

- 10mm sharp right angled tip


- Right upper & lower molar roots
- Smaller version of 1R
- Ergonomic standard handle with thumb hold for increased comfort.
- Superior quality Stainless Steel construction.

PACKAGE INCLUDE:

- 1 Piece 4R Seldin Elevator 10mm Tip Standard Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Root Elevators
APPLICATION: Dental Extraction
BLADE TIP SIZE: 10mm
TIP ANGLE: Right
HANDLE: Octagonal hollow handle
REUSABLE: Yes
AUTOCCLAVED: Yes



**WI-70002148****5 Seldin Root Elevator 34S – 4mm Tip****\$31.24** per pieceTrade price from **\$28.74**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pieces Seldin Elevator #34S
Dental Extraction Surgery Instruments

Seldin Root Elevator 34S Straight with 4mm Tip used during extractions to loosen the tooth from the gums to ease extraction and prevent trauma to surrounding teeth and tissue. The Seldin elevator #34S is used for upper and lower anterior roots. It helps professionals in easy extraction with its multiple features.



SELDIN ROOT ELEVATOR 34S FEATURES

- 4mm sharp straight tip for upper & lower anterior roots.
- Long shaft and a very thin, rounded, slightly hollow tip.
- Ergonomic Handle with thumb hold for increased comfort.

- Superior quality Stainless Steel construction.

PACKAGE INCLUDE:

- 5 Pcs Seldin Root Elevator #34S Straight – 4mm Tip

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 4mm**CURVATURE:** Straight**HANDLE:** Octagonal hollow handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test,**WI-70003899****6 Periosteal Elevator Seldin #S23 Double Ended****\$41.99** per pieceTrade price from **\$38.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs #S23 Periosteal Elevator Seldin 8" (20.3cm) Double Ended Large Flap Retractor / Small Periosteal 10mm & 13mm wide Slightly Curved Blade with Large Flat Handle Quality Stainless Steel Dental Surgical Instruments

Seldin Periosteal Elevator #S23 are designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. Seldin S23 periosteal elevator designed with a large flat handle with double ended design approximately 10mm and 13mm wide blades. The smaller end is sharp for tissue reflecting, and the larger end is slightly curved and is a non-cutting blade, used for tissue retracting.

Periosteal elevator S23 manufactured from high grade stainless steel to exacting specifications to ensure

optimal performance. Seldin elevators is autoclavable and crafted from high-quality stainless steel ensuring durability and sustained performance. 100% satisfaction guaranteed

PERIOSTEAL ELEVATOR SELDIN #S23 FEATURES

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Broad Rounded End & Broad Sharp Chisel End.
- Slight Curved 10mm and 13mm Wide Blades.
- Large Flat Handle.
- Used to Elevate the Periosteum
- Reusable.
- Made With Premium Quality Stainless Steel



**WI-70002142****Apical Dental Elevator #46 Straight 3.5mm Tip****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Apical Root Dental Elevator #46
Straight 3.5mm Tip with Standard
Handle Premium Quality Stainless Steel
Dental Extraction Surgery Instruments

46 Apical Root elevator is a dental tool designed to facilitate extractions by releasing and elevating the upper roots from the periodontal ligament, preventing trauma to surrounding teeth. With overall 6.25" (16cm) overall length, dental elevator #46 features a lightweight standard handle with a thumb hold, is made of high-quality stainless steel, is fully autoclavable, and is reusable.


APICAL DENTAL ELEVATOR 46 FEATURES

- 6.25" (16cm) Overall Length.
- 3.5mm Tip; Straight, Apical

- Lightweight Ergonomic Standard Handle
- Thumb Hold for Increased Comfort.
- Resuable.
- Superior Quality Stainless Steel Construction.

PACKAGE INCLUDE:

- 1 Pc Apical Root Dental Elevator #46 Straight 3.5mm Tip Standard Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 3.5mm**TIP ANGLE:** Straight**HANDLE:** Octagonal hollow handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes – Must Autoclave Before Use**WI-70002139****Apical Elevator #302 Left Angled 3mm Tip****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Apical Elevator #302 6" (15.2cm) Long
Left Angled 3mm Tip with Standard
Handle Premium Quality Stainless Steel
Dental Extraction Surgery Instruments



Apical Root Elevator #302 Used during extractions to loosen the tooth from the gums to ease extraction and prevent trauma to surrounding teeth and tissue. It works as a lever pushed between the gums and the tooth to loosen the surrounding tissues. Manufactured from surgical grade stainless steel provides a high degree of flexibility while conducting clinical procedures. Apical Elevator comes with highly polished finish resistance and anti-oxidation and can be autoclaved.



APICAL ELEVATOR #302 FEATURES

- Overall Length 6" (15.2cm).
- 3mm Left Angled Tip for Broken or Deep-seated Roots.
- Lightweight Stainless Steel Octagonal Handle.
- Thumb Hold for Increased Comfort.
- Fully Autoclavable & Reusable.
- Superior Quality Stainless Steel Construction.

PACKAGE INCLUDE:

- 1 Pc Apical Root Elevator 302 Left – 3mm Tip – Standard Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 3mm**TIP ANGLE:** Left Angled

**WI-70002138****Apical Elevator #303 – Right****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Elevator Apical 303 Right Angled 3mm Tip with Thumb Hold Standard Handle Premium Quality Stainless Steel Dental Extraction Surgery Instruments




Apical Root Elevator #303 Used during extractions to loosen the tooth from the gums to ease extraction and prevent trauma to surrounding teeth and tissue. It works as a lever pushed between the gums and the tooth to loosen the surrounding tissues.

APICAL ELEVATOR #303 FEATURES

- 3mm right angled tip for broken or deep-seated roots
- Lightweight stainless steel handle
- Thumb hold for increased comfort
- Superior quality Stainless Steel construction

PACKAGE INCLUDE:

- 1 Pc Apical Root Elevator 303 Right – 3mm Tip – Standard Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 3mm**TIP ANGLE:** Right Angled**HANDLE:** Octagonal hollow handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packed**WI-70002135****Apical Elevator #303HF Serrated Tip****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#303HF Apical Root Elevator Serrated Back-Action Bend with 2.5mm Blade Tip Standard Handle Premium Quality Stainless Steel Dental Extraction Surgery Instrument

Apical Root Elevator #303HF or Elevator Apexo 303 has serrated backward bend tip used to loosen the tooth from the periodontal ligament and ease extraction. It is lightweight, unbendable, and reliable for dental procedures. Manufactured with premium quality stainless steel for optimal results and precision.

APICAL ELEVATOR #303HF FEATURES



- Crafted with 6.25" (15.9cm) Overall Length.
- 2.5mm Backward Bend Tip.

- Lightweight Stainless Steel Handle.
- Thumb Hold for Increased Comfort.
- Reusable.
- Superior Quality Stainless Steel Construction

PACKAGE INCLUDE:

- 1 Pc Apical Root Elevator 303HF Serrated Back-Action Bend 2.5mm Tip

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 2.5mm**TIP:** Serrated**TIP ANGLE:** Back-Action Bend Straight**HANDLE:** Octagonal hollow handle**REUSABLE:** Yes



WI-70003897

Apical Root Elevators 301 Straight 2.5mm Small Tip

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Apical Root Elevator 301 Small 2.5mm Straight Tip with Standard Handle
Premium Quality Stainless Steel Dental Surgical Instruments.

#301 Apical Root Elevator is designed for tooth extraction by elevating the periosteum layer. 301 Apexo Small straight elevator with a tip measuring 2.5mm is frequently used for beginning the luxation of an erupted tooth, before application of the forceps. Apical Root Elevator #301 ideal for easy tooth extraction. It is lightweight, unbendable, and reliable for dental procedures. Manufactured with premium quality stainless steel for optimal results and precision.

APICAL ELEVATORS 301 FEATURES

- Crafted with 6" (15.2cm) Overall Length.
- 2.5mm Straight, Small, Smooth Tip for Deep-Seated Roots
- Lightweight Stainless Steel Handle
- Thumb Hold for Increased Comfort
- Reusable
- Superior Quality Stainless Steel Construction

PACKAGE INCLUDE:

- 1 Pc Apical Root Elevator 301 Straight – 2.5mm Tip – Standard Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Root Elevators
APPLICATION: Dental Extraction
BLADE TIP SIZE: 3mm
TIP ANGLE: Slight Angled



WI-42896583

Berenard Curved Elevator Dental Root Extraction Tool

\$7.00 per piece

Trade price from **\$3.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Berenard Curved Elevator Dental Root Extraction Tool — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a curved working end for reach around structures.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$7.00, falling to \$3.00 at 10 pieces — about 57% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42892597

Bernard Spade Root Tip Elevator Stainless Steel Curved Dental

\$7.00 per piece

Trade price from **\$3.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Bernard Spade Root Tip Elevator Stainless Steel Curved Dental — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a curved working end for reach around structures.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$7.00, falling to \$3.00 at 10 pieces — about 57% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002141

Cogswell A Root Elevator 6.25" Straight 4mm Tip

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extraction Cogswell #A Root Elevator 6.25" (15.9cm) Straight 4mm Tip with Standard Handle Premium Quality Stainless Steel Dental Surgical Instruments

Cogswell elevators are used for third molar extraction and provide a complete solution for loosening the periodontal ligament and elevating the root. The straight #A Cogswell elevator is ideal for loosening the periodontal ligament around all faces of a root. Used without leverage, a straight Cogswell dental elevator will safely sever the periodontal ligament without causing trauma to the surrounding tissue, or damage to the adjacent teeth. The handle design allows for a balance of force across the handle. All out dental elevators are handcrafted from


stainless steel to ensure excellent performance and durability.



COGSWELL A ROOT ELEVATOR FEATURES

- Crafted with 6.25" (15.9cm) Overall Length
- Straight 4mm Wide Tips
- Sharp, Thin and Narrow Spade-shaped Tip
- Lightweight Stainless Steel Standard Handle.
- Thumb Hold for Increased Comfort.
- Fully Autoclaved and Reusable.
- Superior Quality Stainless Steel Construction.

PACKAGE INCLUDE:

- 1 Pc Cogswell #A Root Elevator 6.25" Straight 4mm Tip Standard Handle



**WI-70000008****Cryer #44 Dental Root Elevator****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Cryer Root Elevator #44 6" (15.2cm)
Small Sharp Left Angle Tip with
Standard Elevator Handle. Premium
Quality Stainless Steel Tooth Extraction
Surgical Dental Instruments

#44 Cryer Elevator is a surgical dental instrument used to loosen and remove tooth roots and fragments, particularly from lower molars. Cryer Root Elevator 44 features a sharp, angled, flag-shaped tip for accessing root bifurcations or empty sockets, adaptable for right or left configurations. The ergonomic standard elevator handle provides enhanced control in confined spaces, making it vital for dental professionals. Crafted from surgical-grade stainless steel to ensure strength and corrosion resistance. The instrument is reusable

and autoclavable, meeting standard dental hygiene protocols.

44 SMALL CRYER ELEVATOR FEATURES

- 6 Inch (15.2cm) Overall Length.
- Left Angle Small Triangle Tip for Root Fragment Removal.
- Ergonomic Handle Provides Improved Grip and Control for the Dentist.
- Designed for the Extraction of Root Segments During Dental Procedures.
- Corrosion Resistance and Anti-oxidation.
- Fully Autoclavable for Reuse.
- Made With German Grade Stainless Steel

PACKAGE INCLUDE**WI-70003885****Cryer #45 Dental Root Elevator****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Cryer Root Elevator #45 6" (15.2cm)
Small Sharp Right Angle Tip with
Standard Elevator Handle. Premium
Quality Stainless Steel Tooth Extraction
Surgical Dental Instruments

Cryer Root Elevator #45 is a surgical dental instrument used to loosen and remove tooth roots and fragments, particularly from lower molars. 45 small cryer elevator features a sharp, angled, flag-shaped tip for accessing root bifurcations or empty sockets, adaptable for right or left configurations. The ergonomic standard elevator handle provides enhanced control in confined spaces, making it vital for dental professionals. Crafted from surgical-grade stainless steel to ensure strength and corrosion resistance. The instrument is reusable

and autoclavable, meeting standard dental hygiene protocols.

45 SMALL CRYER ELEVATOR FEATURES

- 6 Inch (15.2cm) Overall Length.
- Right Angle Small Triangle Tip for Root Fragment Removal.
- Ergonomic Handle Provides Improved Grip and Control for the Dentist.
- Designed for the Extraction of Root Segments During Dental Procedures.
- Corrosion Resistance and Anti-oxidation.
- Fully Autoclavable for Reuse.
- Made With German Grade Stainless Steel

PACKAGE INCLUDE

**WI-70000005****Cryer Root Elevator #25****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Cryer Root Elevator #25 Sharp Left Triangle Tip with Standard Elevator Handle. Premium Quality Stainless Steel Tooth Extraction Surgical Dental Instruments

#25 Cryer Elevator is a surgical dental instrument used to loosen and remove tooth roots and fragments, particularly from lower molars. Cryer root elevator #25 features a sharp, triangled tip for accessing root bifurcations or empty sockets, adaptable for right or left configurations. The ergonomic standard elevator handle provides enhanced control in confined spaces, making it vital for dental professionals. Crafted from surgical-grade stainless steel to ensure strength and corrosion resistance. The instrument is reusable

and autoclavable, meeting standard dental hygiene protocols.

DEL 25 CRYER ELEVATOR FEATURES

- 6 Inch (15.2cm) Overall Length.
- Medium Left Triangle Tip for Root Fragment Removal.
- Ergonomic Handle Provides Improved Grip and Control for the Dentist.
- Designed for the Extraction of Root Segments During Dental Procedures.
- Corrosion Resistance and Anti-oxidation.
- Fully Autoclavable for Reuse.
- Made With German Grade Stainless Steel

PACKAGE INCLUDE**WI-70000004****Cryer Root Elevator #26****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Cryer Root Elevator #26 Sharp Right Triangle Tip with Standard Elevator Handle. Premium Quality Stainless Steel Tooth Extraction Surgical Dental Instruments

#26 Cryer Elevator is a surgical dental instrument used to loosen and remove tooth roots and fragments, particularly from lower molars. Cryer elevator #26 features a sharp triangle-shaped tip for accessing root bifurcations or empty sockets, adaptable for right or left configurations. The ergonomic standard elevator handle provides enhanced control in confined spaces, making it vital for dental professionals. Crafted from surgical-grade stainless steel to ensure strength and corrosion resistance. The instrument is reusable

and autoclavable, meeting standard dental hygiene protocols.

DEL 25 CRYER ELEVATOR FEATURES

- 6 Inch (15.2cm) Overall Length.
- Medium Right Triangle Tip for Root Fragment Removal.
- Ergonomic Handle Provides Improved Grip and Control for the Dentist.
- Designed for the Extraction of Root Segments During Dental Procedures.
- Corrosion Resistance and Anti-oxidation.
- Fully Autoclavable for Reuse.
- Made With German Grade Stainless Steel

PACKAGE INCLUDE



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-42928029

Curved Dental Luxating Elevator Stainless Steel Oral Surgery Tool

\$7.00 per piece

Trade price from \$3.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Curved Dental Luxating Elevator Stainless Steel Oral Surgery Tool is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a curved working end for reach around structures.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$3.00 at 1 piece — about 57% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-42913338

Dental Berenard Spade Elevator Concave Curved Root Extraction Tool

\$7.00 per piece

Trade price from \$3.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Berenard Spade Elevator Concave Curved Root Extraction Tool is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a curved working end for reach around structures.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$7.00, falling to \$3.00 at 10 pieces — about 57% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



WI-70004107

Dental Bone Curettes Curettage

\$20.53 per piece

Trade price from **\$18.89**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Bone Curettes Curettage is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$20.53, falling to \$18.89 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003886

Dental Crane Pick Elevator #8

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Elevator Crane Pick #8 6" (15.2cm) Sharp Penetrating Crane Tip with Standard Elevator Handle. Premium Quality Stainless Steel Tooth Extraction Surgical Dental Instruments

Crane Pick Elevator #8 is a dental surgical instrument designed for various functions: it serves to retrieve roots and fragments by leveraging its sharp pick tip, allows for the separation of the periodontal ligament to facilitate smoother extractions, and assists in luxation and splitting of multi-rooted teeth during complex extractions.

Crafted from surgical-grade stainless steel to ensure strength and corrosion resistance. The instrument is reusable and autoclavable, meeting standard dental hygiene protocols.


CRANE PICK ELEVATOR #8 FEATURES

- 6 Inch (15.2cm) Overall Length.
- Sharp Penetrating Crane Tip.
- Ergonomic Handle Provides Improved Grip and Control for the Dentist.
- Designed for the Extraction of Root Segments During Dental Procedures.
- Corrosion Resistance and Anti-oxidation.
- Fully Autoclavable for Reuse.
- Made With German Grade Stainless Steel

PACKAGE INCLUDE

- 1 Pc #8 Crane Pick Elevator 6" Sharp Crane Tip

TECHNICAL INFORMATION





WI-42910772

Dental Curved Luxating Elevator Root Extraction Instrument

\$7.00 per piece

Trade price from \$3.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

We manufacture Dental Curved Luxating Elevator Root Extraction Instrument for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a curved working end for reach around structures.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$3.00 at 1 piece — about 57% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000006

Dental Elevator #3

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Elevator #3 Slightly Curved Sharp Tip with Standard Elevator Handle. Premium Quality Stainless Steel Tooth Extraction Surgical Dental Instruments

Dental Elevator #3 (Miltex DEL3) is a dental instrument designed for loosening teeth, expanding sockets, and assisting in the removal of root tips prior to extractions. Its main functions include tooth luxation, ligament cutting to sever the periodontal ligament, and aiding in the extraction of broken or stubborn root tips. Crafted from surgical-grade stainless steel to ensure strength and corrosion resistance. The instrument is reusable and autoclavable, meeting standard dental hygiene protocols.

- 6 Inch (15.2cm) Overall Length.
- Slightly Curved Round Edge Blade.
- Ergonomic Handle Provides Improved Grip and Control for the Dentist.
- Designed for Loosen/Removal of Root Tips Prior to Extractions.
- Corrosion Resistance and Anti-oxidation.
- Fully Autoclavable for Reuse.
- Made With German Grade Stainless Steel

PACKAGE INCLUDE


- 1 Pc Dental Elevator #3 Slightly Curved Blade

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENT TYPE: Dental Elevators

#3 DENTAL ELEVATOR FEATURES





WI-70000007

Dental Elevator #81

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Elevator #81 6" (15.2cm) Straight Flat Sharp Tip with Standard Elevator Handle. Premium Quality Stainless Steel Tooth Extraction Surgical Dental Instruments

81 Elevator is a specialized surgical instrument used in dental extractions to loosen teeth and small roots from the periodontal ligament space. It is designed to sever the periodontal ligament with its sharp tip, allowing for atraumatic extractions, especially in upper anterior teeth and small root fragments. This technique minimizes damage to surrounding bone and soft tissue by gently rocking and cutting to ease tooth removal. Crafted from surgical-grade stainless steel to ensure strength and corrosion resistance. The instrument is reusable and

autoclavable, meeting standard dental hygiene protocols.

DENTAL ELEVATOR DEL 81 FEATURES

- 6 Inch (15.2cm) Overall Length.
- Straight, Flat Blade for PDL Access and Tooth Luxation.
- Ergonomic Handle Provides Improved Grip and Control for the Dentist.
- Designed for Upper Anterior and Small Root Extraction.
- Corrosion Resistance and Anti-oxidation.
- Fully Autoclavable for Reuse.
- Made With German Grade Stainless Steel

PACKAGE INCLUDE



WI-32893707

Dental Elevators

\$4.00 per piece

Trade price from \$2.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Elevators is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$4.00, falling to \$2.00 at 25 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-01670769

Dental Extraction Spear Root Elevators

\$3.85 per piece

Trade price from \$1.98
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dental Extraction Spear Root Elevators, made for dental, hygiene and oral surgery work. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.85, falling to \$1.98 at 10 pieces — about 49% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42917335

Dental Heidbrink Root Elevator Stainless Steel Oral Surgery Tool

\$6.00 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Heidbrink Root Elevator Stainless Steel Oral Surgery Tool for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$6.00, falling to \$2.00 at 10 pieces — about 67% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42905879

Dental Luxating Elevator Curved Stainless Steel Tooth Extraction

\$7.00 per piece

Trade price from \$3.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Luxating Elevator Curved Stainless Steel Tooth Extraction is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a curved working end for reach around structures.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$7.00, falling to \$3.00 at 1 piece — about 57% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25658193

Dental Luxation Kit Surgical Luxating Tooth Extraction Root Elevators

\$4.00 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Luxation Kit Surgical Luxating Tooth Extraction Root Elevators for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.00, falling to \$2.00 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-15381460

Dental Luxation Root Elevators 6-Piece

\$12.00 per piece

Trade price from **\$8.50**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dental Luxation Root Elevators 6-Piece, made for dental, hygiene and oral surgery work. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$12.00, falling to \$8.50 at 1 set — about 29% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35359757

Dental Molt M9 Elevator Blue Titanium Coated

\$7.00 per piece

Trade price from **\$1.30**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

We manufacture Dental Molt M9 Elevator Blue Titanium Coated for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$7.00, falling to \$1.30 at 1 piece — about 81% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-30367248

Dental PDL Luxating Root Elevators Set

\$25.00 per piece

Trade price from \$19.00
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Dental PDL Luxating Root Elevators Set is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$25.00, falling to \$19.00 at 10 sets — about 24% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27580627

Dental Periodontal Surgical Stainless Steel Bone Curettes (Lucas)

\$6.00 per piece

Trade price from \$3.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Periodontal Surgical Stainless Steel Bone Curettes (Lucas is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$6.00, falling to \$3.00 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35370689

Dental Periosteal Elevator Molt M9 Blue Titanium Instrument

\$7.00 per piece

Trade price from \$1.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Periosteal Elevator Molt M9 Blue Titanium Instrument is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$7.00, falling to \$1.50 at 1 piece — about 79% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32940520

Dental Root Elevator Set 5 Pcs

\$4.00 per piece

Trade price from \$2.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Root Elevator Set 5 Pcs — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. Supplied as a 5-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$2.00 at 25 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-09976073

Dental Root Elevator Stainless Steel Dental Instruments

\$4.50 per piece

Trade price from \$1.99
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Root Elevator Stainless Steel Dental Instruments — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.50, falling to \$1.99 at 10 pieces — about 56% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32963538

Dental Root Elevator Stainless Steel Surgical

\$4.00 per piece

Trade price from \$2.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Root Elevator Stainless Steel Surgical for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$4.00, falling to \$2.00 at 25 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25812497

Dental Root Elevator Surgical Instrument Orthodontic

\$4.50 per piece

Trade price from **\$1.99**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Root Elevator Surgical Instrument Orthodontic for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$4.50, falling to \$1.99 at 10 pieces — about 56% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42880564

Dental Root Elevators Set 10PCS Stainless Steel Extraction Tools

\$31.00 per piece

Trade price from **\$27.60**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Root Elevators Set 10PCS Stainless Steel Extraction Tools is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$31.00, falling to \$27.60 at 1 set — about 11% off the single-unit price. Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32277655

Dental Root Planning Lab Stainless Steel Handle Surgical

\$4.00 per piece

Trade price from \$0.90
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Root Planning Lab Stainless Steel Handle Surgical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$0.90 at 10 pieces — about 78% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42919101

Dental Root Tip Elevator DEL 303 Stainless Steel Tooth

\$7.00 per piece

Trade price from \$3.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Root Tip Elevator DEL 303 Stainless Steel Tooth for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$3.00 at 10 pieces — about 57% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42874818

Dental Root Tip Elevator DEL 303 Surgical Instrument

\$7.00 per piece

Trade price from **\$3.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Dental Root Tip Elevator DEL 303 Surgical Instrument is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$7.00, falling to \$3.00 at 10 pieces — about 57% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44789087

Dental Scaler Stainless Steel Periodontal Curette Surgical Instrument

\$4.00 per piece

Trade price from **\$0.90**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Scaler Stainless Steel Periodontal Curette Surgical Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.00, falling to \$0.90 at 10 pieces — about 78% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-42909751

Dental Straight Elevators Set 3PCS Tooth Extraction Surgical Tools

\$13.00 per piece

Trade price from \$6.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Straight Elevators Set 3PCS Tooth Extraction Surgical Tools — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a straight working end for direct, in-line access.

Supplied as a 3-piece set. We can add, drop or substitute items so the tray

matches your protocol rather than ours. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge.

Listed at \$13.00, falling to \$6.00 at 1 set — about 54% off the single-unit price. Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-32960530

Dental Tooth Extraction Elevator Apical Root Tip Extractor Luxating

\$4.00 per piece

Trade price from \$2.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Dental Tooth Extraction Elevator Apical Root Tip Extractor Luxating is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$4.00, falling to \$2.00 at 25 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-36929137

Double-Ended Dental Hygienist Sickle Scaler Curettes Machinery

\$8.00 per piece

Trade price from \$5.00
Minimum order 5 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Double-Ended Dental Hygienist Sickle Scaler Curettes Machinery — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: two working ends on one handle, halving what you carry.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$5.00 at 5 sets — about 38% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-42902709

Heidbrink Elevator Dental Root Extraction Instrument

\$6.00 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Heidbrink Elevator Dental Root Extraction Instrument is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$6.00, falling to \$2.00 at 10 pieces — about 67% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

**WI-70003895****Heidbrink Root Elevator #2 Right Angled Slightly Curved Pointed Tip****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heidbrink Root Tip Elevator #2 Right Angled 6" (15.2cm) Slightly Curved Pointed Tip with Hexagonal Handle Premium Quality Stainless Steel Root Tips Tooth Extraction Dental Instruments.

Heidbrink elevators are a single-ended instruments designed for the effective removal of root tips and residual roots in surgical procedures. This delicate dental elevators are used to retrieve loose root tips that may break away from a tooth during extraction and are in hard-to-reach places. Heidbrink also referred to as a root tip pick elevator with a narrow sharp slightly curved design and lightweight hexagonal handle for better control during use. Heidbrink Root Tip Elevator is crafted with straight, left and right angled tips.

This heidbrink root elevators with its unique design is worthy of being a part of your dental kit. Made with premium grade stainless steel, these elevators are lightweight, rust-free, and reusable.

HEIDBRINK ROOT TIP ELEVATOR FEATURES

- Crafted with 6" (15.2cm) Overall Length.
- 2.5mm Working End.
- Narrow, Sharp and Slightly Curved Tips.
- Available in Straight, Left and Right Angled Profiles.
- Single Ended Instruments Used to Extract Root Tips.
- Hexagonal Handle Creates Less Force on Tooth With Better Gripping.

**WI-70003896****Heidbrink Root Elevator 1 Straight Slightly Curved Pointed Tip****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heidbrink Root Tip Elevator 1 Straight 6" (15.2cm) Slightly Curved Pointed Tip with Hexagonal Handle Premium Quality Stainless Steel Root Tips Tooth Extraction Dental Instruments.

Heidbrink elevators are a single-ended instruments designed for the effective removal of root tips and residual roots in surgical procedures. This delicate dental elevators are used to retrieve loose root tips that may break away from a tooth during extraction and are in hard-to-reach places. Heidbrink also referred to as a root tip pick elevator with a narrow sharp slightly curved design and lightweight hexagonal handle for better control during use. Heidbrink Root Tip Elevator is crafted with straight, left and right angled tips. This heidbrink root elevators with its

unique design is worthy of being a part of your dental kit. Made with premium grade stainless steel, these elevators are lightweight, rust-free, and reusable.

HEIDBRINK ROOT TIP ELEVATOR FEATURES

- Crafted with 6" (15.2cm) Overall Length.
- 2.5mm Working End.
- Narrow, Sharp and Slightly Curved Tips.
- Available in Straight, Left and Right Angled Profiles.
- Single Ended Instruments Used to Extract Root Tips.
- Hexagonal Handle Creates Less Force on Tooth With Better Gripping.
- Unique Shape Reducing Chance of Tooth Fractures.



**WI-70003894****Heidbrink Root Elevator 3 Left Angled Slightly Curved Pointed Tip****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heidbrink Root Tip Elevator 3 Left Angled 6" (15.2cm) Slightly Curved Pointed Tip with Hexagonal Handle Premium Quality Stainless Steel Root Tips Tooth Extraction Dental Instruments.

Heidbrink elevators are a single-ended instruments designed for the effective removal of root tips and residual roots in surgical procedures. This delicate dental elevators are used to retrieve loose root tips that may break away from a tooth during extraction and are in hard-to-reach places. Heidbrink also referred to as a root tip pick elevator with a narrow sharp slightly curved design and lightweight hexagonal handle for better control during use. Heidbrink Root Tip Elevator is crafted with straight, left and right angled tips.

This heidbrink root elevators with its unique design is worthy of being a part of your dental kit. Made with premium grade stainless steel, these elevators are lightweight, rust-free, and reusable.

HEIDBRINK ROOT TIP ELEVATOR FEATURES

- Crafted with 6" (15.2cm) Overall Length.
- 2.5mm Working End.
- Narrow, Sharp and Slightly Curved Tips.
- Available in Straight, Left and Right Angled Profiles.
- Single Ended Instruments Used to Extract Root Tips.
- Hexagonal Handle Creates Less Force on Tooth With Better Gripping.

**WI-42904604****Heidbrink Root Tip Elevator H1 Straight Hollow Handle Dental****\$7.00** per pieceTrade price from **\$3.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Heidbrink Root Tip Elevator H1 Straight Hollow Handle Dental for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a straight working end for direct, in-line access.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$3.00 at 10 pieces — about 57% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42892741

Heidbrink Root Tip Elevator Stainless Steel Dental Extraction Instrument

\$6.00 per piece

Trade price from **\$2.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Heidbrink Root Tip Elevator Stainless Steel Dental Extraction Instrument for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$6.00, falling to \$2.00 at 10 pieces — about 67% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000644

Heidbrink Root Tip Pick Elevator H1 Straight Pointed Tip

\$8.99 per piece

Trade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heidbrink Root Tip Pick Elevator H1 Straight 6.5" (16.5cm) Slightly Curved Pointed Tip with Heavy Duty Slim Hexagonal Handle Premium Quality Stainless Steel Root Tips Tooth Extraction Dental Instruments.

right angled tips. This heidbrink root elevators with its unique design is worthy of being a part of your dental kit. Made with premium grade stainless steel, these elevators are lightweight, rust-free, and reusable.

Heidbrink Root Teaser elevators are a single-ended instruments designed for the effective removal of root tips and residual roots in surgical procedures. This delicate dental elevators are used to retrieve loose root tips that may break away from a tooth during extraction and are in hard-to-reach places. Heidbrink also referred to as a root tip pick elevator with a narrow sharp slightly curved design and lightweight hexagonal handle for better control during use. Heidbrink Root Tip Elevator is crafted with straight, left and

HEIDBRINK ROOT TIP PICK ELEVATOR FEATURES

- Crafted with 6.5" (16.5cm) Overall Length.
- 2.5mm Working End.
- Narrow, Sharp and Slightly Curved Tips.
- Available in Straight, Left and Right Angled Profiles.
- Single Ended Instruments Used to Extract Root Tips.
- Heavy Duty Slim Hexagonal Long Shank Handle.



**WI-70000643**

Heidbrink Root Tip Pick Elevator H2 Right Angled Pointed Tip

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heidbrink Root Tip Pick Elevator H2 Right Angled 6.5" (16.5cm) Slightly Curved Pointed Tip with Heavy Duty Slim Hexagonal Handle Premium Quality Stainless Steel Root Tips Tooth Extraction Dental Instruments.

Heidbrink Root Teaser elevators are a single-ended instruments designed for the effective removal of root tips and residual roots in surgical procedures. This delicate dental elevators are used to retrieve loose root tips that may break away from a tooth during extraction and are in hard-to-reach places. Heidbrink also referred to as a root tip pick elevator with a narrow sharp slightly curved design and lightweight hexagonal handle for better control during use. Heidbrink Root Tip Elevator is crafted with straight, left and

right angled tips. This heidbrink root elevators with its unique design is worthy of being a part of your dental kit. Made with premium grade stainless steel, these elevators are lightweight, rust-free, and reusable.

HEIDBRINK ROOT TIP PICK ELEVATOR FEATURES

- Crafted with 6.5" (16.5cm) Overall Length.
- 2.5mm Working End.
- Narrow, Sharp and Slightly Curved Tips.
- Available in Straight, Left and Right Angled Profiles.
- Single Ended Instruments Used to Extract Root Tips.
- Heavy Duty Slim Hexagonal Long Shank Handle.

**WI-70000642**

Heidbrink Root Tip Pick Elevator H3 Left Angled Pointed Tip

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heidbrink Root Tip Pick Elevator H3 Left Angled 6.5" (16.5cm) Slightly Curved Pointed Tip with Heavy Duty Slim Hexagonal Handle Premium Quality Stainless Steel Root Tips Tooth Extraction Dental Instruments.

Heidbrink Root Teaser elevators are a single-ended instruments designed for the effective removal of root tips and residual roots in surgical procedures. This delicate dental elevators are used to retrieve loose root tips that may break away from a tooth during extraction and are in hard-to-reach places. Heidbrink also referred to as a root tip pick elevator with a narrow sharp slightly curved design and lightweight hexagonal handle for better control during use. Heidbrink Root Tip Elevator is crafted with straight, left and

right angled tips. This heidbrink root elevators with its unique design is worthy of being a part of your dental kit. Made with premium grade stainless steel, these elevators are lightweight, rust-free, and reusable.

HEIDBRINK ROOT TIP PICK ELEVATOR FEATURES

- Crafted with 6.5" (16.5cm) Overall Length.
- 2.5mm Working End.
- Narrow, Sharp and Slightly Curved Tips.
- Available in Straight, Left and Right Angled Profiles.
- Single Ended Instruments Used to Extract Root Tips.
- Heavy Duty Slim Hexagonal Long Shank Handle.





WI-42880759

Heidbrink Straight Elevator Dental Oral Surgery Tool

\$7.00 per piece

Trade price from \$3.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Heidbrink Straight Elevator Dental Oral Surgery Tool is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a straight working end for direct, in-line access.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$3.00 at 10 pieces — about 57% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32966532

Light Weight Stainless Steel Dental Elevators Ergonomic Grip 3

\$4.00 per piece

Trade price from \$2.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Light Weight Stainless Steel Dental Elevators Ergonomic Grip 3 is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: an ergonomic handle for long lists without fatigue.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$4.00, falling to \$2.00 at 25 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000540****Molt Bone Curette #2/4 Double Ended****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Molt Bone Curette 2/4 Double Ended
3.5mm & 7.0mm Working Tip Round
Pattern Handle Premium Quality
Stainless Steel Surgical Dental
Instruments.

2/4 Molt surgical curette is a double ended bone curette designed for curettage, cyst removal, and tooth socket debridement. With a 3.5 mm and 7.0 mm working ends allows the surgeon to scrape away bone with durable and sharp edges. Molt 2/4 bone curette has round pattern non-slip knurled handle provides optimum control and omits hand fatigue.

This double ended bone curette can also use in a variety of veterinary dental procedures for efficient scraping of diseased tissue from oral/dental wounds for cats and dogs. Molt Bone

curette 2-4, crafted with high-quality surgical stainless steel with a smooth surface, ensures a strong grip. Completely sterilizable and reusable

MOLT BONE CURETTE #2/4 FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- Double Ended.
- 3.5 mm – 7.0mm Spoon Sharp Blades.
- Ergonomic Design Round Pattern Grip Handle.
- Sterilizable / Reusable.
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE

- 1 Pc Molt Bone Curette #2/4 Double Ended 7"

**WI-70000227****Molt Bone Curette #2/4 Double Ended Hollow Handle****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Molt Bone Curette 2/4 Double Ended
3.5mm & 7.0mm Working Tip Rounded
Pattern Hollow Handle Premium Quality
Stainless Steel Surgical Dental
Instruments.

2/4 Molt surgical curette is a double ended bone curette designed for curettage, cyst removal, and tooth socket debridement. With a 3.5 mm and 7.0 mm working ends allows the surgeon to scrape away bone with durable and sharp edges. Molt 2/4 bone curette has round pattern hollow handle provides optimum control and omits hand fatigue.

This double ended bone curette can also use in a variety of veterinary dental procedures for efficient scraping of diseased tissue from oral/dental wounds for cats and dogs. Molt Bone

curette 2-4, crafted with high-quality surgical stainless steel with a smooth surface, ensures a strong grip. Completely sterilizable and reusable

MOLT BONE CURETTE #2/4 HOLLOW FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- Double Ended.
- 3.5mm – 7.0mm Spoon Sharp Blades.
- Ergonomic Design Lightweight Hollow Handle.
- Sterilizable / Reusable.
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE

- 1 Pc Molt Bone Curette #2/4 Double Ended – Hollow Handle

**WI-70000541****Molt Bone Curette #6R Single Ended 2.9mm Spoon Tip****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Molt Bone Curette #6R Single Ended 6" (15.2cm) Long Right Angled 2.9mm Spoon Tip Hollow Handle Premium Quality Stainless Steel Surgical Curette Dental Instruments

6R Molt surgical curette is a right-angled bone curette designed for curettage, cyst removal, and tooth socket debridement. With a 2.9 mm spoon-shaped working end Molt #6R Bone Curette allows the surgeon to scrape away bone with durable and sharp edges. Moreover, the Molt bone curette 6R hand dental instrument is designed with a hollow handle that provides optimum control and omits hand fatigue. Molt Bone curette 6R, crafted with high-quality surgical stainless steel with a smooth surface,

ensures a strong grip. Completely sterilizable and reusable

MOLT BONE CURETTE #6R FEATURES

- Crafted with 6" (15.2cm) Overall Length.
- Right Angled.
- 2.9mm Spoon Shape Blade.
- Single Ended Curette.
- Hexagonal Hollow Grip Handle.
- Sterilizable / Reusable.
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE

- 1 Pc Molt Bone Curette #6R Single Ended 2.9mm Spoon Tip

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**WI-35292882****Molt M9 Dental Periosteal Elevator Stainless Steel Surgical Tool****\$7.00** per pieceTrade price from **\$1.50**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Molt M9 Dental Periosteal Elevator Stainless Steel Surgical Tool is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$7.00, falling to \$1.50 at 1 piece — about 79% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70002024**

Offset Dental Crown Splitter Remover 5.75\"/>

\$12.99 per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Posterior Crown Spreader 5.75" (15cm) Offset Dental Crown Splitter with Elevator Handle Premium Quality Stainless Steel Restorative Dentistry Instruments

With max Posterior Crown Spreader permanent crown removing process can be quick and easy. Offset crown spreader is used to split the crown in two parts after partial occlusal cutting, when the crown is required to remove. Crown is first cut into two halves with a bur and these two parts are then pushed apart with this crown splitter. Its working end is tempered and strong enough to resist the deformity while pushing apart the crown into two halves. The unique multi-directional tips allow approach from any angle or any direction. Our tools offer a sturdy grip

and plus-shaped tip design, ensuring a comfortable experience for dentists, endodontists, and orthodontists. This Elevator-Type Posterior Crown Spreader, crafted from high-quality stainless steel, is a reliable and durable tool for dental professionals, fully autoclaved and reusable.

POSTERIOR CROWN SPREADER FEATURES

- Crafted with 5.75" (15cm) Overall Length.
- Plus-shaped Multi Directional Tip.
- Easily Pry Crown Apart.
- Elevator handle for optimal comfort and grip
- Extra long handle and stem allows better access to posterior regions.
- Reusable.

**WI-27587648**

Orthodontic Surgical Dental Bone Curettes Double Ended

\$6.00 per pieceTrade price from **\$3.00**

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Orthodontic Surgical Dental Bone Curettes Double Ended is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: two working ends on one handle, halving what you carry.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$6.00, falling to \$3.00 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000705****Periosteal Elevator Baby Molt #7A Double Ended****\$7.99** per pieceTrade price from **\$7.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Periosteal Elevator Baby Molt #7A Double Ended 6.5" (16.5cm) Hexagonal Grip Handle Premium Quality Stainless Steel Dental Instruments.

Molt-7 elevators have thin blades on both ends that are used to cut periodontal ligaments and loosen teeth. Periosteal Baby Molt #7A comprises a flat, oval tip at each end, one blunt and one sharp. Surgeons use the blunt oval tip for lifting and retracting the periosteal flap, and the sharp tip to dissect it. In surgery dentist often needs to separate a bone or tooth from the periosteum that covers it. This is done with a periosteal elevator. The dentist may also use it to gain access to retained roots and surrounding bone. Elevators are manufactured from

stainless steel to ensure optimal performance.

PERIOSTEAL ELEVATOR BABY MOLT #7A FEATURES

- Overall Length 6.5" (16.5cm).
- Double Ended.
- One Flat and One Oval Tip.
- Non-Slip Hexagonal Grip Handle.
- Use to Elevate, Retract and Dissecting Periosteal Flap.
- Corrosion Resistance.
- Fully Autoclavable & Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Periosteal Elevator Baby Molt #7A Double Ended

TECHNICAL INFORMATION**WI-70003892****Periosteal Elevator Baby Molt #7A Double Ended with Round Handle****\$8.99** per pieceTrade price from **\$8.27**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Periosteal Elevator Baby Molt #7A Double Ended 6.5" (16.5cm) Round Pattern Grip Handle Premium Quality Stainless Steel Dental Instruments.

Molt-7 elevators have thin blades on both ends that are used to cut periodontal ligaments and loosen teeth. Periosteal Baby Molt #7A comprises a flat, oval tip at each end, one blunt and one sharp. Surgeons use the blunt oval tip for lifting and retracting the periosteal flap, and the sharp tip to dissect it. In surgery dentist often needs to separate a bone or tooth from the periosteum that covers it. This is done with a periosteal elevator. The dentist may also use it to gain access to retained roots and surrounding bone. Elevators are manufactured from

stainless steel to ensure optimal performance.

PERIOSTEAL ELEVATOR BABY MOLT #7A FEATURES

- Overall Length 6.5" (16.5cm).
- Double Ended.
- One Flat and One Oval Tip.
- Non-Slip Knurled Grip Handle.
- Use to Elevate, Retract and Dissecting Periosteal Flap.
- Corrosion Resistance.
- Fully Autoclavable & Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Periosteal Elevator Baby Molt #7A Double Ended

TECHNICAL INFORMATION



WI-70000327

Prichard Periosteal Elevator PR-3 with Blue Titanium Coated Tips

\$10.99 per piece

Trade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Prichard Periosteal Elevator PC-3 Blue Titanium Coated Double Ends with Round Pattern Hollow Grip Handle Premium Quality Stainless Steel Dental Implant Surgery Instruments.

PR-3 Prichard Periosteal elevators are designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. Featuring two blades, one large flat blade and one narrow sharp blade on round tip. PR-3 prichard periosteal elevator is a specialized tool for dental implant and oral surgery, featuring durable, corrosion-resistant blue titanium-coated tips. Its lightweight, ergonomic hollow grip handle allows for optimal control and comfort. Made from premium German-grade stainless steel, this elevator provides reliable

performance and precision tips for smooth periosteal reflection with minimal tissue trauma, making it essential for flap elevation, periodontal, and implant surgery.

PRICHARD PERIOSTEAL ELEVATOR PR-3 FEATURES

- Double Ended.
- Blue Titanium Coated Tips.
- One Large Flat Blade and One Narrow Sharp Blade.
- Non-Slip Knurled Grip Hollow Handle.
- Use to Elevate, Retract and Dissecting Periosteal Flap.
- Corrosion Resistance.
- Reusable.
- Made with Premium Quality Stainless Steel.



WI-42881674

Reusable Dental Root Tip Elevator Stainless Steel Oral Surgery

\$7.00 per piece

Trade price from **\$3.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Reusable Dental Root Tip Elevator Stainless Steel Oral Surgery — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: full reusability — clean, autoclave and use again.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$7.00, falling to \$3.00 at 10 pieces — about 57% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25818494

Root Elevator Stainless Steel Dental Instrument

\$4.50 per piece

Trade price from **\$1.99**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Root Elevator Stainless Steel Dental Instrument is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.50, falling to \$1.99 at 10 pieces — about 56% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42919007

Root Tip Elevator Stainless Steel

\$7.00 per piece

Trade price from **\$3.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Root Tip Elevator Stainless Steel is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$3.00 at 10 pieces — about 57% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70002144****Seldin Elevator 4L – 10mm Tip****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


4L Seldin Root Elevator 6" (15.2cm) Left Angled 10mm Tip with Standard Handle Premium Quality Stainless Steel Dental Extraction Surgery Instruments.

Seldin Elevator 4L is used for pulling up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction. 4L elevator is a smaller version of 1L. Seldin Elevators used to removes soft tissue and periosteal membrane from the alveolar bone with pushing strokes. It has one working end, a robust shank, and a sturdy handle. Sharp tips helping in extraction procedure. These can be used for both the upper and lower molar roots.

- Crafted with 6" (15.2cm) Overall Length.
- 10mm Sharp Left Angled Tip
- Left Upper & Lower Molar Roots
- Smaller Version of 1L.
- Ergonomic Standard Handle With Thumb Hold for Increased Comfort.
- Fully Autoclaved & Reusable.
- Superior Quality Stainless Steel Construction.

PACKAGE INCLUDE:

- 1 Pc 4L Seldin Elevators 10mm Tip Standard Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 10mm**SELDIN ELEVATOR 4L FEATURES****WI-70003893****Seldin Periosteal Elevator #S23 8" Double Ended****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#S23 Periosteal Elevator Seldin 8" (20.3cm) Double Ended Large Flap Retractor / Small Periosteal 10mm & 13mm wide Slightly Curved Blade with Large Flat Handle Quality Stainless Steel Dental Surgical Instruments

Seldin Periosteal Elevator #S23 are designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. Seldin S23 periosteal elevator designed with a large flat handle with double ended design approximately 10mm and 13mm wide blades. The smaller end is sharp for tissue reflecting, and the larger end is slightly curved and is a non-cutting blade, used for tissue retracting.

Periosteal elevator S23 manufactured from high grade stainless steel to exacting specifications to ensure

optimal performance. Seldin elevators is autoclavable and crafted from high-quality stainless steel ensuring durability and sustained performance. 100% satisfaction gauranteed

PERIOSTEAL ELEVATOR SELDIN #S23 FEATURES

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Broad Rounded End & Broad Sharp Chisel End.
- Slight Curved 10mm and 13mm Wide Blades.
- Large Flat Handle.
- Used to Elevate the Periosteum
- Reusable.
- Made With Premium Quality Stainless Steel



**WI-70003890****Seldin Root Elevator 34 – 5mm Tip****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Seldin Root Elevator 34 Straight with 5mm Tip used during extractions to loosen the tooth from the gums to ease extraction and prevent trauma to surrounding teeth and tissue. The #34 Seldin is used for upper and lower anterior roots. It helps professionals in easy extraction with its multiple features.

SELDIN ELEVATOR 34 FEATURES

- 5mm sharp straight tip for upper & lower anterior roots.
- Long shaft and a very thin, rounded, slightly hollow tip.
- Ergonomic Handle with thumb hold for increased comfort.
- Superior quality Stainless Steel construction.

PACKAGE INCLUDE:

- 1 Pc Seldin Root Elevator #34 Straight – 5mm Tip

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 5mm**CURVATURE:** Straight**HANDLE:** Octagonal hollow handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packed**WI-70003240****Seldin Root Elevator 34S – 4mm Tip****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Seldin Root Elevator 34S Straight with 4mm Tip used during extractions to loosen the tooth from the gums to ease extraction and prevent trauma to surrounding teeth and tissue. The #34S Seldin is used for upper and lower anterior roots. It helps professionals in easy extraction with its multiple features.

SELDIN ELEVATOR 34S FEATURES

- 4mm sharp straight tip for upper & lower anterior roots.
- Long shaft and a very thin, rounded, slightly hollow tip.
- Ergonomic Handle with thumb hold for increased comfort.
- Superior quality Stainless Steel construction.

PACKAGE INCLUDE:

- 1 Pc Seldin Root Elevator #34S Straight – 4mm Tip

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 4mm**CURVATURE:** Straight**HANDLE:** Octagonal hollow handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packed



WI-32967470

Set 2 German Stainless Steel Dental Elevators

\$4.00 per piece

Trade price from \$2.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Set 2 German Stainless Steel Dental Elevators — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$2.00 at 25 pieces — about 50% off the single-unit price. Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32974535

Set 2 German Stainless Steel Hollow Handle Dental

\$4.00 per piece

Trade price from \$2.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Set 2 German Stainless Steel Hollow Handle Dental, made for dental, hygiene and oral surgery work. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$4.00, falling to \$2.00 at 25 pieces — about 50% off the single-unit price. Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42911488

Stainless Steel DEL 303 Dental Oral Surgery Elevator

\$7.00 per piece

Trade price from **\$3.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel DEL 303 Dental Oral Surgery Elevator — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$3.00 at 10 pieces — about 57% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32931581

Stainless Steel Dental Elevators Set 2 German Ergonomic

\$4.00 per piece

Trade price from **\$2.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Elevators Set 2 German Ergonomic is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: an ergonomic handle for long lists without fatigue.

Supplied as a set — the contents can be adjusted to your list. The steel is

heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$4.00, falling to \$2.00 at 25 pieces — about 50% off the single-unit price. Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32955561

Stainless Steel Dental Elevators Set Root

\$4.00 per piece

Trade price from **\$2.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Elevators Set Root is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$2.00 at 25 pieces — about 50% off the single-unit price. Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32290513

Stainless Steel Dental Scaler Gracey Curette

\$4.00 per piece

Trade price from **\$0.90**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Scaler Gracey Curette — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.00, falling to \$0.90 at 10 pieces — about 78% off the single-unit price. Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27606659

Stainless Steel Dental Scalers & Bone Curettes

\$6.00 per piece

Trade price from **\$3.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Scalers & Bone Curettes is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$6.00, falling to \$3.00 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25722160

Stainless Steel Dental Tooth Extraction Elevator Apical Root

\$4.00 per piece

Trade price from **\$2.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Dental Tooth Extraction Elevator Apical Root, made for dental, hygiene and oral surgery work. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.00, falling to \$2.00 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-31557272

Surgical Dental Elevators Stainless Steel Ergonomic Grip

\$4.00 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Dental Elevators Stainless Steel Ergonomic Grip is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: an ergonomic handle for long lists without fatigue.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$4.00, falling to \$2.00 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-27573589

Tip Dental Bone Curettes Periodontal Surgery 2.5mm/3.5mm/4mm

\$6.00 per piece

Trade price from \$3.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Tip Dental Bone Curettes Periodontal Surgery 2.5mm/3.5mm/4mm for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The 2.5mm size shown is our standard; other sizes in the same pattern are made to order.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$6.00, falling to \$3.00 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.



WI-27581634

Volkmann Bone Curettes

\$6.00 per piece

Trade price from \$3.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Volkmann Bone Curettes for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$6.00, falling to \$3.00 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

Request a quotation

Send us product codes from this catalogue, quantities and your country.
We reply the same working day on WhatsApp.

Online shop	shop.waqarintl.com
Email	sales@waqarintl.com
WhatsApp /	
Phone	+92 307 6173627
Works	Roras Road, Opposite Anwar Khawaja Industries, Sialkot 51310, Pakistan

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.