



WAQAR INTERNATIONAL
PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

Forceps Dental Instruments

Reusable stainless-steel forceps forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



242

PRODUCTS

\$0.70

PRICES FROM

OEM

PRIVATE LABEL

Global

WORLDWIDE
EXPORT

Contents

242 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

Product	Code	Price	Page
10 American Pattern Extracting Forceps No. 10S	WI-70003259	\$89.99	8
10 Extracting Forceps 3FS Dental Surgical Instruments	WI-70004220	\$89.99	8
10pcs Dental Surgical Instruments Set Adult	WI-29420860	\$53.00	9
12 American Pattern Extracting Forceps No. 10S	WI-70004169	\$99.99	9
12 cm Anatomic Tweezers 1x2 Teeth	WI-26507772	\$2.50	10
12 Dental Extracting Forceps #217 for 1st & 2nd Lower Molars	WI-70004132	\$99.99	10
12 Extracting Forceps 3FS Dental Surgical Instruments	WI-70004178	\$99.99	11
12 Extracting Forceps Dental Surgical 1 Hen	WI-70004185	\$99.99	11
12cm Stainless Steel Surgical Adson Toothed Forceps Dissecting	WI-26506845	\$2.50	12
16cm Stainless Steel Dental Extracting Forceps No 16	WI-36801721	\$12.50	12
2 Dental Extracting Forceps #150AS & #151AS Serrated Split Beaks	WI-70000404	\$19.99	13
2 Extracting Forceps #88R & #88L with Gold Hollow Python Grip Handle	WI-70000980	\$26.99	13
2 Extracting Forceps 18L & 18R for Upper 1st & 2nd Molars	WI-70000951	\$19.99	14
2 Extracting Forceps 3FS Dental Surgical Instruments	WI-70004126	\$19.99	14
2 Universal Extraction Forceps #150 & #151 For Lower & Upper Teeth	WI-70004137	\$19.99	15
3 American Extracting Forceps #15 Lower Molars	WI-70004149	\$29.99	15
3 American Extracting Forceps #16 Lower 1st & 2nd Molars	WI-70003315	\$29.99	16
3 Dental Extracting Forceps #10H Upper Molars Universal	WI-70004236	\$29.99	16
3 Dental Extracting Forceps #150 #151 #23	WI-70002154	\$29.99	17
3 Dental Extracting Forceps #17 Lower Molars	WI-70004252	\$29.99	17
3 Dental Extracting Forceps #18R for Upper 1st & 2nd Molars	WI-70004238	\$29.99	18
3 Dental Extracting Forceps #217 for 1st & 2nd Lower Molars	WI-70004227	\$29.99	18
3 Dental Extracting Forceps #222 Lower 3rd Molars	WI-70004250	\$29.99	19
3 Dental Extracting Forceps #88L 1st & 2nd Upper Molar	WI-70004219	\$29.99	19
3 Dental Extracting Forceps 18L for Upper 1st & 2nd Molars	WI-70004230	\$29.99	20
3 Dental Extracting Kells Forceps #99C	WI-70004232	\$29.99	20
3 Dental Extraction Forceps #150 Upper Incisors, Root Teeth	WI-70004243	\$29.99	21
3 Dental Tooth Extracting Forceps #10S with Serrated Jaws	WI-70004244	\$29.99	21
3 Extracting Forceps #103 for Lower Bicuspids, Incisors & Roots	WI-70004242	\$29.99	22
3 Extracting Forceps #150A for Upper Incisors, Cuspids and Bicuspids	WI-70004235	\$29.99	22
3 Extracting Forceps #203 Lower Incisors, Canines, Premolars & Roots	WI-70004217	\$29.99	23
3 Extracting Forceps #210S Upper 3rd Molars, Universal	WI-70004223	\$29.99	23
3 Extracting Forceps #213 for Upper Central & Canines Roots	WI-70004222	\$29.99	24
3 Extracting Forceps #286 for Upper Roots, Incisors & Premolars	WI-70004233	\$29.99	24
3 Extracting Forceps #6 Upper & Lower Teeth Splitting & Separating	WI-70004239	\$29.99	25
3 Extracting Forceps #69 Upper Lower Roots and Fragments	WI-70004229	\$29.99	25
3 Extracting Forceps #88R For Right 1st And 2nd Upper Molar	WI-70004218	\$29.99	26
3 Extracting Forceps #99A Upper Incisors, Canines & Premolar Root	WI-70004237	\$29.99	26
3 Extracting Forceps 3FS Dental Surgical Instruments	WI-70004221	\$29.99	27
3 Extracting Forceps Dental Surgical 1 Hen	WI-70004226	\$29.99	27
3 Pedo Extracting Forceps #151S Dental Surgical Instrument	WI-70004213	\$29.99	28

Contents (continued)

Product	Code	Price	Page
3 Universal Extraction Forceps #151 for Lower Teeth	WI-70004231	\$29.99	29
3 Universal Pedro Extracting Forceps #17S for Lower Molars	WI-70004253	\$29.99	29
3.5 Inch Stainless Steel Dental Backhaus Towel Forceps	WI-25943254	\$2.50	30
5 Piece Dental Extraction Forceps Set with Gold	WI-42865573	\$30.00	30
6 American Extracting Forceps #15 Lower Molars	WI-70004200	\$59.99	31
6 American Extracting Forceps #16 Lower 1st & 2nd Molars	WI-70003257	\$59.99	31
6 Dental Extracting Forceps #10H Upper Molars Universal	WI-70004194	\$59.99	32
6 Dental Extracting Forceps #150AS Upper Anteriors	WI-70004187	\$59.99	32
6 Dental Extracting Forceps #17 Lower Molars	WI-70004211	\$59.99	33
6 Dental Extracting Forceps #217 for 1st & 2nd Lower Molars	WI-70004186	\$59.99	33
6 Dental Extracting Forceps #222 Lower 3rd Molars	WI-70004209	\$59.99	34
6 Dental Extracting Forceps #88L 1st & 2nd Upper Molar	WI-70004176	\$59.99	34
6 Dental Extracting Forceps 18L for Upper 1st & 2nd Molars	WI-70004189	\$59.99	35
6 Dental Extracting Kells Forceps #99C	WI-70004191	\$59.99	35
6 Dental Extraction Forceps #150 Upper Incisors, Root Teeth	WI-70004202	\$59.99	36
6 Extracting Forceps #103 for Lower Bicuspids, Incisors & Roots	WI-70004201	\$59.99	36
6 Extracting Forceps #150A for Upper Incisors, Cuspids and Bicuspids	WI-70004195	\$59.99	37
6 Extracting Forceps #203 Lower Incisors, Canines, Premolars & Roots	WI-70004174	\$59.99	37
6 Extracting Forceps #210S Upper 3rd Molars, Universal	WI-70004181	\$59.99	38
6 Extracting Forceps #213 for Upper Central & Canines Roots	WI-70004180	\$59.99	38
6 Extracting Forceps #286 for Upper Roots, Incisors & Premolars	WI-70004192	\$59.99	39
6 Extracting Forceps #6 Upper & Lower Teeth Splitting & Separating	WI-70004197	\$59.99	39
6 Extracting Forceps #69 Upper Lower Roots and Fragments	WI-70004188	\$59.99	40
6 Extracting Forceps #88R For Right 1st And 2nd Upper Molar	WI-70004175	\$59.99	40
6 Extracting Forceps #99A Upper Incisors, Canines & Premolar Root	WI-70004196	\$59.99	41
6 Extracting Forceps 3FS Dental Surgical Instruments	WI-70004179	\$59.99	41
6 Extracting Forceps Dental Surgical 1 Hen	WI-70004184	\$59.99	42
6 Pedro Extracting Forceps #151S Dental Surgical Instrument	WI-70004170	\$59.99	42
6 Universal Dental Extracting Forceps No. 24 for Upper Molars	WI-70004206	\$59.99	43
6 Universal Extraction Forceps #151 for Lower Teeth	WI-70004190	\$59.99	43
6 Universal Pedro Extracting Forceps #17S for Lower Molars	WI-70004212	\$59.99	44
Adson Bayonet Forceps 16 cm Medical Dressing Serrated 1.7cm	WI-22083570	\$5.00	44
Adson Brown Dissecting Forceps Teeth Surgical Tissue Medical	WI-00533624	\$2.50	45
Adson Tissue Forceps 1x2 Teeth Reusable Surgical	WI-43462069	\$3.20	45
Adson Tissue Forceps 1x2 Teeth Stainless Steel Surgical Dressing	WI-43417980	\$3.20	46
Adson Tissue Forceps 1x2 Teeth Stainless Steel Surgical Instrument	WI-43465046	\$3.20	46
Adson Tissue Forceps Angled Teeth Precision Tips-Surgical	WI-33489833	\$7.00	47
Allis Atraumatic Ergonomic 3-Jaw Holding Forceps Size	WI-36617604	\$5.00	47
Allis Tissue Forceps 19cm 5x6 Teeth	WI-43452128	\$3.20	48
Allis Tissue Forceps 5x6 Teeth	WI-43427329	\$3.20	48
Allis Tissue Forceps Stainless Steel Instrument Set Dental	WI-35100238	\$5.00	49
Allis Tissue Forceps Toothed Locking Pliers Intestinal Forceps Dental	WI-25901013	\$5.00	49

Contents (continued)

Product	Code	Price	Page
American Extracting Forceps #151A	WI-70004168	\$12.99	50
American Extracting Forceps #16 Lower 1st & 2nd Molars	WI-70004167	\$11.99	50
American Extracting Forceps #16S Lower Molars	WI-70004166	\$11.99	51
American Extracting Forceps #23 with Gold Hollow Python Grip Handle	WI-70000972	\$14.99	51
American Extracting Forceps #88L 1st & 2nd Upper Molar	WI-70001099	\$14.99	52
American Extracting Forceps No. 15 for 1st & 2nd Lower Molars	WI-70003258	\$11.99	52
American Extracting Forceps No. 85A Lower Molars	WI-70004162	\$12.99	53
American Pattern Extracting Forceps #90 Upper Left Molars	WI-70004122	\$11.99	53
American Pado Extracting Forceps #23S for Primary Lower Molars	WI-70004164	\$11.99	54
Atraumatic Dressing Forceps Teeth Precision Tissue Handling Instrument Set	WI-36702823	\$3.00	54
Backhaus Stainless Steel Dental Towel Forceps	WI-34366784	\$2.50	55
Backhaus Towel Forceps Stainless Steel Dental	WI-34363801	\$2.50	55
Bone Reduction Forceps with Teeth Curved 5.5 Inch Orthopedic	WI-35539093	\$40.00	56
Children's Dental Extraction Forceps Curved Straight	WI-13330170	\$5.99	56
Dental Adson Brown Forceps Curved Thumb Dressing Tweezers	WI-43392166	\$6.00	57
Dental Basic Examination Set Mouth Mirror Handle Double-Ended Hook	WI-29087315	\$2.00	57
Dental College Tweezer Forceps-Stainless Steel Serrated Tip Multipurpose	WI-30787180	\$0.80	58
Dental Cotton Holding Tweezers Stainless Steel Clinical Instrument	WI-36801641	\$2.00	58
Dental Dressing Tweezers Precision Hand Instrument	WI-36803638	\$2.00	59
Dental Elevator Forceps Extraction 8-Piece Kit	WI-24285324	\$18.00	59
Dental Examination Kit 3	WI-29020422	\$2.00	60
Dental Examination Set Mirror Probe Tweezers Dentist College	WI-29028418	\$2.00	60
Dental Extracting Forceps #101 for Upper & Lower Premolars	WI-70003316	\$11.99	61
Dental Extracting Forceps #10H Upper Molars Universal	WI-70004144	\$11.99	61
Dental Extracting Forceps #150AS Serrated Split Beaks	WI-70004134	\$11.99	62
Dental Extracting Forceps #17 Lower Molars	WI-70004160	\$11.99	62
Dental Extracting Forceps #18R for Upper 1st & 2nd Molars	WI-70004118	\$11.99	63
Dental Extracting Forceps #1STD Gold Hollow Python Grip Handle	WI-70001370	\$14.99	63
Dental Extracting Forceps #210H Upper 3rd Molars	WI-70004142	\$11.99	64
Dental Extracting Forceps #213 for Upper Central & Canines Roots	WI-70004128	\$11.99	64
Dental Extracting Forceps #222 Lower 3rd Molars	WI-70004158	\$11.99	65
Dental Extracting Forceps #32 Upper Canines, Premolars & Molar	WI-70004159	\$11.99	65
Dental Extracting Forceps #32A Molars Root – Slim Beak	WI-70004157	\$11.99	66
Dental Extracting Forceps #349 for Upper Incisors & Root	WI-70000962	\$11.99	66
Dental Extracting Forceps #53L For Upper 1st & 2nd Left Molars	WI-70004139	\$11.99	67
Dental Extracting Forceps #53L with Gold Hollow Python Grip Handle	WI-70001236	\$14.99	67
Dental Extracting Forceps #62 Upper and Lower Anterior Teeth	WI-70003255	\$12.99	68
Dental Extracting Forceps #65 for Upper Roots Dental Tools	WI-70004156	\$11.99	68
Dental Extracting Forceps #88L 1st & 2nd Upper Molar	WI-70004125	\$11.99	69
Dental Extracting Forceps #88R For Right 1st And 2nd Upper Molar	WI-70004124	\$11.99	69
Dental Extracting Forceps #89 for Upper Right Molars	WI-70004117	\$11.99	70
Dental Extracting Forceps 18L for Upper 1st & 2nd Molars	WI-70004136	\$11.99	70

Contents (continued)

Product	Code	Price	Page
Dental Extracting Forceps No 16 16cm	WI-36864648	\$12.50	71
Dental Extracting Forceps No. 24 for Upper Molars	WI-70004155	\$11.99	71
Dental Extracting Forceps No.44 Lower Anterior	WI-70003637	\$12.99	72
Dental Extracting Forceps Set 5PCS Gold Handle Tooth Extraction	WI-42862583	\$30.00	72
Dental Extracting Forceps Set 8PCS Stainless Steel Tooth Extraction	WI-42826808	\$45.00	73
Dental Extracting Forceps Stainless Steel Oral Surgery Instrument	WI-42856850	\$10.00	73
Dental Extracting Forceps Tooth Removal Instrument for Oral Surgery	WI-46820684	\$8.00	74
Dental Extracting Kells Forceps #99C	WI-70004138	\$11.99	74
Dental Extracting Serrated Forceps	WI-42891146	\$9.00	75
Dental Extraction Forceps #150 Upper Incisors, Root Teeth	WI-70004151	\$11.99	75
Dental Extraction Forceps 10S Gold Handle Stainless Steel Upper	WI-42834901	\$9.00	76
Dental Extraction Forceps Stainless Steel Serrated Jaw Tooth Extraction	WI-42861734	\$9.50	76
Dental Extraction Forceps Stainless Steel Universal Serrated Jaw Instrument	WI-42847770	\$7.00	77
Dental Extraction Serrated Forceps Stainless Steel Tooth Extraction Instrument	WI-42887187	\$9.00	77
Dental Forceps Holder Stainless Steel Dental Instrument Rack Sterilization	WI-36729139	\$1.00	78
Dental Forceps Storage Rack Stainless Steel Sterilization Holder	WI-36755108	\$1.00	78
Dental Hygiene Kit Stainless Steel Oral Care Tools Scaler	WI-42868317	\$31.00	79
Dental Instrument Holder Stainless Steel Forceps Storage	WI-36727154	\$0.90	79
Dental Kit Steel Tweezer Set	WI-29059343	\$2.00	80
Dental Orthodontic Cheatle Forceps	WI-34949415	\$3.00	80
Dental Surgical Tweezers Steel Blade Operation	WI-27787337	\$1.00	81
Dental TC Tweezers 15cm/6" Stainless Steel Surgical Forceps	WI-25924708	\$3.00	81
Dental Tooth Extracting Forceps #10S with Serrated Jaws	WI-70004152	\$11.99	82
Dental Tooth Extraction Forceps Set 7 Gold Plated	WI-30805248	\$65.00	82
Dental Tooth Forceps 10-Piece Set Stainless Steel Silver	WI-45076284	\$42.00	83
Dental Tweezers Surgical Cotton Tweezers	WI-27812198	\$1.00	83
Dental Universal Extraction Forceps Stainless Steel Tooth Extraction Instrument	WI-42861952	\$9.00	84
Dissecting Tissue Thumb Forceps Lab Multi-Purpose Tweezers Blunt 1x2	WI-40981380	\$7.00	84
English Pattern Dental Tooth Extracting Forceps	WI-34250028	\$8.00	85
English Pattern Dental Tooth Extracting Forceps German	WI-32723692	\$8.00	85
Extracting Forceps #1 with Gold Hollow Python Grip Handle	WI-70001324	\$14.99	86
Extracting Forceps #103 for Lower Bicuspid, Incisors & Roots	WI-70004150	\$11.99	86
Extracting Forceps #150A for Upper Incisors, Cuspids and Bicuspid	WI-70004143	\$11.99	87
Extracting Forceps #150AS Diamond Dusted Python Grip Gold Hollow Handle	WI-70000411	\$14.99	87
Extracting Forceps #151 P-1 Diamond Dusted Gold Hollow Handle	WI-70000408	\$14.99	88
Extracting Forceps #151 Pattern-2 Diamond Dusted	WI-70000407	\$14.99	88
Extracting Forceps #151AS Serrated Split Beaks	WI-70004130	\$11.99	89
Extracting Forceps #203 Lower Incisors, Canines, Premolars & Roots	WI-70004123	\$11.99	89
Extracting Forceps #210S Upper 3rd Molars, Universal	WI-70004141	\$11.99	90
Extracting Forceps #217 for 1st & 2nd Lower Molars	WI-70004133	\$11.99	90
Extracting Forceps #23 Cow Horn Lower Molar	WI-70004148	\$11.99	91
Extracting Forceps #23 Lower Molars with Green Color Coated Handle	WI-70001251	\$11.99	91

Contents (continued)

Product	Code	Price	Page
Extracting Forceps #286 for Upper Roots, Incisors & Premolars	WI-70004140	\$11.99	92
Extracting Forceps #53R for Upper 1st & 2nd Right Molars	WI-70004147	\$11.99	92
Extracting Forceps #6 Upper & Lower Teeth Splitting & Separating	WI-70004146	\$11.99	93
Extracting Forceps #69 Upper Lower Roots and Fragments	WI-70004135	\$11.99	93
Extracting Forceps #88R 1st & 2nd Upper Molar with Gold Python Grip	WI-70000981	\$14.99	94
Extracting Forceps #99A Upper Incisors, Canines & Premolar Root	WI-70004145	\$11.99	94
Extracting Forceps 1 Hen Dental Surgical	WI-70004131	\$11.99	95
Extracting Forceps 3FS Dental Surgical Instrument	WI-70004127	\$11.99	95
Extracting Forceps Set Dental Instruments for Tooth Surgery	WI-28825128	\$4.00	96
Extraction Forceps #150 Upper Incisors Gold Hollow Python Grip Handle	WI-70000576	\$14.99	96
Extraction Forceps No 7 Upper Roots English Pattern	WI-31344355	\$8.00	97
Fine Tip Dressing Forceps Teeth Reusable Stainless Steel Surgical	WI-36763643	\$3.00	97
Genuine German Steel Machinery 17-Piece Oral Dental	WI-16643035	\$45.00	98
Gillies Toothed Dissecting Forceps Stainless Steel Dental Surgical Instrument	WI-35080915	\$1.20	98
HOT SALES Stainless Steel Dental Extracting Forceps Set	WI-31352471	\$8.00	99
MD1 Dental Extracting Forceps Stainless Steel Upper Tooth Extraction	WI-42867462	\$8.00	99
Michel Clip Extracting Forceps Surgical Kit for Precise Skin	WI-46955210	\$8.90	100
Micro Precision 16cm Non-Toothed Dressing Forceps Reusable Surgical	WI-36804704	\$0.90	100
Non Toothed Dissecting Forceps Medical	WI-36260860	\$1.20	101
Non Toothed Dissecting Forceps Stainless Steel Surgical Dressing Tissue	WI-36299043	\$1.20	101
Orthodontic Dental Surgical Tweezer Set 14cm Medical	WI-29050394	\$2.00	102
Pedo Extracting Forceps #151S Dental Surgical Instrument	WI-70004121	\$11.99	102
Potts Smith Tissue Forceps 1×2 Teeth	WI-35032080	\$3.00	103
Precision Dental Exam Tool Kit Stainless Steel Dental Probe	WI-29064365	\$2.00	103
Reusable Dental College Tweezers Precision Cotton Holding Instrument	WI-36827612	\$2.00	104
Reusable Dental Tooth Extraction Forceps Stainless Steel 3PCS Kit	WI-42844798	\$19.00	104
Solid Color Stainless Steel Brown	WI-27423479	\$2.50	105
Stainless Steel 16cm Dental Extracting Forceps No	WI-36811740	\$12.50	105
Stainless Steel Adson Surgical Forceps (Toothed	WI-30455551	\$1.00	106
Stainless Steel Adult Tooth Extraction Dental	WI-31374392	\$8.00	106
Stainless Steel Dental Cotton Dressing Tweezer Surgical Tissue	WI-34891604	\$2.00	107
Stainless Steel Dental Extracting Forceps	WI-31406595	\$8.00	107
Stainless Steel Dental Extracting Forceps-Professional	WI-31379371	\$8.00	108
Stainless Steel Dental Extraction Forceps Ergonomic Grip	WI-32764672	\$8.00	108
Stainless Steel Dental Extraction Forceps Oral Surgery Instrument	WI-46801758	\$9.00	109
Stainless Steel Dental Forceps Tooth Extracting	WI-32770678	\$8.00	109
Stainless Steel Dental Hygiene Tool Kit Removal	WI-29044372	\$2.00	110
Stainless Steel Dental Surgical Tweezers & Tooth Forceps	WI-26456933	\$2.50	110
Stainless Steel Dental TC Surgical Tweezers Ergonomic	WI-25875769	\$3.00	111
Stainless Steel Dental Tooth Extracting Forceps	WI-32730713	\$8.00	111
Stainless Steel Dental Tooth Upper Roots Extraction	WI-31373307	\$8.00	112
Stainless Steel Dressing Forceps Straight Surgical Thumb Tweezer	WI-34951329	\$2.00	112

Contents (continued)

Product	Code	Price	Page
Stainless Steel Gillies Dissecting Forceps Toothed Dental Surgical Instrument	WI-35092446	\$1.20	113
Stainless Steel Heaney Hysterectomy Forceps 21cm	WI-26029780	\$8.00	113
Stainless Steel Pediatric Dental Forceps Ergonomic Grip	WI-32732705	\$8.00	114
Stainless Steel Russian Toothed Dissecting Forceps	WI-40970984	\$7.00	114
Stainless Steel Serrated Curved (15mm) Dental Root	WI-25799728	\$4.50	115
Stainless Steel Surgical Tissue Forceps/Fine Teeth Grasping Claw	WI-26535751	\$2.50	115
Stainless Steel Toothless Dissecting Forceps Precision Tissue Handling Instrument	WI-36832566	\$0.70	116
Stainless-Steel Dental Tweezers Medical	WI-27811200	\$1.00	116
Steel Dental Examination Set	WI-29077349	\$2.00	117
Surgical Dental Miller Articulating Straight Steel Paper Forceps	WI-22435893	\$3.00	117
Surgical Tissue Forceps 17.5cm Stainless Steel 1×2 Teeth	WI-36917027	\$8.50	118
T/C Dental Extracting Forceps #17 Lower Molars	WI-70001254	\$14.99	118
T/C Dental Extraction Forceps #150 Upper Incisors, Root Teeth	WI-70001098	\$14.99	119
T/C Extracting Forceps #53L For Upper 1st & 2nd Left Molars	WI-70001237	\$14.99	119
T/C Extracting Forceps #53R for Upper 1st & 2nd Right Molars	WI-70001235	\$14.99	120
T/C Extracting Forceps #65 Upper Roots, Fragments & Incisors	WI-70001101	\$14.99	120
T/C Wynman Crown Remover Gripper Forceps 7" Curved Serrated Tips	WI-70002178	\$12.99	121
TC Extracting Forceps #1STD Upper Incisors and Cuspids	WI-70001371	\$14.99	121
TC Extraction Forceps #151 for Lower Teeth	WI-70000410	\$14.99	122
Thumb Dressing Forceps Teeth Surgical Tissue Pickup Forceps	WI-36762663	\$3.00	122
Tooth Extracting Forceps #1STD Upper Incisors and Cuspids	WI-70001372	\$11.99	123
Tooth Extraction Dental Forceps	WI-10243183	\$5.99	123
Tooth Extraction Forceps Dental Set Surgical Tool	WI-46793830	\$8.80	124
Toothed Dissecting Forceps Stainless Steel Surgical Tissue Forceps	WI-36275794	\$1.20	124
Toothed Dissecting Forceps Stainless Steel Tissue	WI-36277895	\$1.20	125
Toothed Forceps 6" for Hospitals-Brown Adson Surgical	WI-27296632	\$2.39	125
Universal Dental Extracting Forceps Stainless Steel Oral Surgery Tool	WI-42901035	\$9.00	126
Universal Extraction Forceps #151 for Lower Teeth	WI-70004120	\$11.99	126
Universal Extraction Forceps #151 Green Handle	WI-70000409	\$11.99	127
Universal Extraction Forceps #17 with Gold Hollow Python Grip Handle	WI-70001252	\$14.99	127
Universal Pedo Extracting Forceps #17S for Lower Molars	WI-70004161	\$11.99	128
Vulsellum Forceps 3×4 Teeth's	WI-06803073	\$3.50	128

**WI-70003259****10 American Pattern Extracting Forceps No. 10S****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Heavy Duty American Pattern Tooth Extracting Forceps Fig. 10S Upper Molar Serrated Curved Jaws German Grade Stainless Steel Dental Instrument.

Durable and sturdy, Upper molars extracting forceps #10S are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals dentists.

EXTRACTING FORCEPS 10S FEATURES

- Serrated Beak For Better Grip
- Figure 10S For Upper Molar
- Smoothly Curved Tips Making Them Suitable For Both Sides
- Non-Slip Grip Handle

- German Stainless Steel

PACKAGE INCLUDE

- 10 Pcs Extracting Forceps No. 10S Curved Serrated – Upper Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**WI-70004220****10 Extracting Forceps 3FS Dental Surgical Instruments****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Extracting Forceps #3FS Woodward's American Pattern Tooth Extraction Forceps Fig. 3FS For Lower Molars German Grade Stainless Steel Dental Surgical Instruments

Extracting forceps #3FS used for lower molars extracting teeth from the alveolar bones. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by dental professionals.

DENTAL EXTRACTING FORCEPS #3FS FEATURES

- Pointed Beak Design
- For Lower Molars
- Non-Slip Grip Handle

- German Stainless Steel

PACKAGE INCLUDE

- 10 Pcs Woodward Dental Extracting Forceps 3FS – Lower Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test

**WI-29420860****10pcs Dental Surgical Instruments Set Adult****\$53.00** per piece

Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

10pcs Dental Surgical Instruments Set Adult is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$53.00. Minimum order 10 sets, and the price moves with quantity. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70004169****12 American Pattern Extracting Forceps No. 10S****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Heavy Duty American Pattern Tooth Extracting Forceps Fig. 10S Upper Molar Serrated Curved Jaws German Grade Stainless Steel Dental Instrument.

Durable and sturdy, Upper molars extracting forceps #10S are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals dentists.

EXTRACTING FORCEPS 10S FEATURES

- Serrated Beak For Better Grip
- Figure 10S For Upper Molar
- Smoothly Curved Tips Making Them Suitable For Both Sides
- Non-Slip Grip Handle

- German Stainless Steel

PACKAGE INCLUDE

- 12 Pcs Extracting Forceps No. 10S Curved Serrated – Upper Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test



WI-26507772

12 cm Anatomic Tweezers 1x2 Teeth

\$2.50 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture 12 cm Anatomic Tweezers 1x2 Teeth for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

The 12cm size shown is our standard; other sizes in the same pattern are made to order. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.50, falling to \$2.00 at 10 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004132

12 Dental Extracting Forceps #217 for 1st & 2nd Lower Molars

\$99.99 per piece

Trade price from \$91.99
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 American Pattern Dental Extracting Forceps #217 For Lower Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#217 Extraction Forceps is used to extract lower molars. The working end of this tooth extraction instrument is designed with a tri-pointed shaped beaks, making sure it reaches deep into the socket for a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

- Tri-pointed Serrated Beak Design
- For 1st & 2nd Lower Molars
- Non-slip Grip Handle
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE

- 12 Pcs Dental Extracting Forceps #217 – Lower Molars

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Extracting Forceps
APPLICATION: Teeth Extraction
POLISHED: Yes
STERILE: Yes
REUSABLE: Yes
AUTOCCLAVED: Yes
RUST PREVENTION
PROCEDURE: Passivated



**WI-70004178****12 Extracting Forceps 3FS Dental Surgical Instruments****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Extracting Forceps #3FS Woodward's American Pattern Tooth Extraction Forceps Fig. 3FS For Lower Molars German Grade Stainless Steel Dental Surgical Instruments

Extracting forceps #3FS used for lower molars extracting teeth from the alveolar bones. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by dental professionals.

DENTAL EXTRACTING FORCEPS #3FS FEATURES

- Pointed Beak Design
- For Lower Molars
- Non-Slip Grip Handle

- German Stainless Steel

PACKAGE INCLUDE

- 12 Pcs Woodward Dental Extracting Forceps 3FS – Lower Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**WI-70004185****12 Extracting Forceps Dental Surgical 1 Hen****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Extracting Forceps 1 HEN, Upper Incisors Henahan Teeth Extraction German Grade Stainless Steel Surgical Dental Instrument

Extraction forceps are the most familiar instruments with regard to the extraction of tooth teeth. Intra-alveolar extractions are mainly accomplished using elevators and extraction forceps. They are the instruments that are used to remove the tooth out of the alveolar bone socket.

Doctor use this tool for upper incisors, #1HEN Henahan. Our tools are professional quality and used by professionals. All Dental Forceps are handcrafted from german grade stainless steel for logtivity and durability

EXTRACTING FORCEPS 1 HEN FEATURES

- #1HEN Henahan Dental Froceps
- Upper Incisors and Canines Forceps
- Non-Slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 12 Pcs Extracting Forceps 1 HEN

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated



WI-26506845

12cm Stainless Steel Surgical Adson Toothed Forceps Dissecting

\$2.50 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture 12cm Stainless Steel Surgical Adson Toothed Forceps Dissecting for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

The 12cm size shown is our standard; other sizes in the same pattern are made to order. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.50, falling to \$2.00 at 10 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36801721

16cm Stainless Steel Dental Extracting Forceps No 16

\$12.50 per piece

Trade price from \$8.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of 16cm Stainless Steel Dental Extracting Forceps No 16, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 16cm size shown is our standard; other sizes in the same pattern are made to order.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$12.50, falling to \$8.00 at 25 pieces — about 36% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000404****2 Dental Extracting Forceps #150AS & #151AS Serrated Split Beaks****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Dental Extracting Forceps #150AS & #151AS Split Beaks, Serrated for Lower & Upper Anterior German Grade Stainless Steel American Pattern Tooth Extraction Dental Surgical Instruments.

#150AS & #151AS Extraction Forceps are used for extracting teeth from the alveolar bone. Lower and Upper split serrated beaks designed to provide increased grasp beneath the crown to prevent crown or root breakage. Serrations provide increased mechanical grasp on the root without increased pressure. These dental forceps can be used left and right.

American pattern forceps have a smooth hinge so no pin is visible. The stainless steel construction ensures maximum durability and dependability.

You can reuse these #150AS and 151AS extracting forceps over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

AMERICAN EXTRACTING FORCEPS #150AS & #151AS FEATURES

- For Lower and Upper Anteriors.
- Serrated, Split Beaks Provide an Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Extracting Forceps #150AS Split Beaks Serrated – Upper Anterior
- 1 Pc Extracting Forceps #151AS Split Beaks Serrated – Lower Anterior

**WI-70000980****2 Extracting Forceps #88R & #88L with Gold Hollow Python Grip Handle****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Extracting Forceps #88R & #88L with Right and Left Pointed, Split Beak Design with Gold Python Grip Handle for 1st & 2nd Upper Molars German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.



American Pattern dental extraction forceps No. 88R and 88L used to extract 1st and 2nd molars. This extraction forceps set has pointed split beak design. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #88R & #88L FEATURES

- For Upper 1st & 2nd Molars.
- Right and Left Pointed, Split Beak Design.
- Gold Hollow Python Grip Handle.
- American Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc American Pattern Extracting Forceps #88R – Gold Python Grip Hollow Handle
- 1 Pc American Pattern Extracting Forceps #88L – Gold Python Grip Hollow Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

**WI-70000951****2 Extracting Forceps 18L & 18R for Upper 1st & 2nd Molars****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Extracting Forceps 18L & 18R for Upper 1st & 2nd Molars, Left & Right Pointed Tip Design with Hook Handle German Grade Stainless Steel Teeth Extraction Surgical Dental Instruments

American Extracting forceps No. 18L and 18R hook handle are the perfect tool to remove teeth from the alveolar bone. This 2 pcs extraction forceps set has left & right pointed beak design features an American Pattern. Bifurcated pointed beaks designed for upper first and second molars. Hooked handle provides increased leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you

#18L & #18R EXTRACTING FORCEPS FEATURES

- Left & Right Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc American Extraction Forceps #18L Left Pointed Beak with Hook Handle
- 1 Pc American Extraction Forceps #18R Right Pointed Beak with Hook Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental

**2
PCS****WI-70004126****2 Extracting Forceps 3FS Dental Surgical Instruments****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Extracting Forceps #3FS Woodward's American Pattern Tooth Extraction Forceps Fig. 3FS For Lower Molars German Grade Stainless Steel Dental Surgical Instruments

Extracting forceps #3FS used for lower molars extracting teeth from the alveolar bones. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by dental professionals.

DENTAL EXTRACTING FORCEPS #3FS FEATURES

- Pointed Beak Design
- For Lower Molars
- Non-Slip Grip Handle

- German Stainless Steel

PACKAGE INCLUDE

- 2 Pcs Woodward Dental Extracting Forceps 3FS – Lower Molars

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental
Extracting Forceps
APPLICATION: Teeth Extraction
POLISHED: Yes
STERILE: Yes
REUSABLE: Yes
AUTOCCLAVED: Yes
RUST PREVENTION
PROCEDURE: Passivated
ULTRASONIC CLEANED: Yes
TESTS PERFORMED: Boil Test, Performance Test, Shape Test



**WI-70004137**

2 Universal Extraction Forceps #150 & #151 For Lower & Upper Teeth

\$19.99 per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Pcs Universal Extraction Forceps #150 & #151 Serrated Beaks for Lower & Upper Incisors, Canines, Premolars & Roots – Premium Quality Stainless Steel American Pattern Tooth Extracting Surgical Dental Instrument.

#150 & #151 Extraction Forceps are used for extracting teeth from the alveolar bone. Lower and Upper serrated beaks designed to provide increased mechanical grasp on the root without increased pressure. These dental forceps can be used left and right. American pattern forceps have a smooth hinge so no pin is visible. The stainless steel construction ensures maximum durability and dependability. You can reuse these #150 and 151 extracting forceps over time as they will not bend or tarnish on you. High quality

stainless steel, excellent international quality standard.

UNIVERSAL EXTRACTION FORCEPS #150 & #151 FEATURES

- For Lower Incisors, Canines, Premolars & Roots.
- For Upper Incisors, Canines, Premolars & Roots.
- Thin, Tapered, and Sharp Parallel Beaks for Good Contact.
- Longitudinal Serrations for Improved Grip and Atraumatic Extraction.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Universal Extraction Forceps #150 Upper Teeth

**WI-70004149**

3 American Extracting Forceps #15 Lower Molars

\$29.99 per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Extracting Forceps #15 for 1st & 2nd Lower Molars with Hook Handle, Universal German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American pattern dental extracting forceps No. 15 used to extract 1st and 2nd lower molars. This extraction forceps has tri-pointed beak design. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #15 FEATURES

- For Lower 1st & 2nd Molars.

- Tri-Pointed Beaks Design.
- Non-Slip Grip Hook Handle.
- American Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 3 Pcs American Pattern Extracting Forceps #15 Hook Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Surgery Procedure**BEAK/JAWS:** Tri-Pointed Design**POLISHED:** Yes**HANDLE:** Hook Handle**REUSABLE:** Yes**AUTOCLAVED:** Yes**RUST PREVENTION**

**3
PCS****WI-70003315****3 American Extracting Forceps #16 Lower 1st & 2nd Molars****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Universal Extracting Forceps #16 Hook Handle for Lower 1st & 2nd Molars Cowhorn Beak Design American Pattern German Grade Stainless Steel Tooth Extraction Surgical Dental Instruments

American Extracting forceps No. 16 hook handle are the perfect tool to remove teeth from the alveolar bone. This universal extraction forceps #16 has pointed cowhorn beak design features an American Pattern. Bifurcated pointed beaks designed for lower first and second molars. Hooked handle provides increased leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over

time as they will not bend or tarnish on you

AMERICAN EXTRACTION FORCEPS #16 FEATURES

- Cow Horn Beak Design.
- For Lower 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 3 Pcs American Extraction Forceps #16 Cowhorn Beak with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**3
PCS****WI-70004236****3 Dental Extracting Forceps #10H Upper Molars Universal****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Extracting Forceps #10H Universal Upper Molars with Hook Handle Premium Quality Stainless Steel American Pattern Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 10 with hook handle is used to extract upper molars.

American extraction forceps #10H features conical beaks designed to adapt to upper molars to improve grasp. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

- Conical Beaks Design.
- Upper Molars, Universal
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 3 Pcs American Pattern Extraction Forceps #10H with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction

Forceps

BEAK/JAWS: Conical Serrated Beaks**POLISHED:** Yes**STERILE:** Yes

**WI-70002154****3 Dental Extracting Forceps #150 #151 #23****\$29.99** per pieceTrade price from **\$27.59**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty Dental Extracting Forceps #150 #151 #23, Apical Forceps, Cowhorn Forceps, Lower Molars, Atraumair Forceps 3 Pcs Set


EXTRACTION FORCEP FEATURES

- Presidential Upper Universal Forceps used to extract teeth from alveolar bone & Presidential Lower Universal Forceps used to extract teeth from alveolar bone.
- Thin, tapered and sharp parallel beaks improve root and crown contact and longitudinal serrations provide additional grip. The design of Apical Forceps facilitates an atraumatic tooth extraction.
- Surgical grade stainless steel promotes longevity and corrosion resistance

- Thin, sharp, conical beaks designed to be able to cut through surrounding bone so the beaks achieve contact with the crown and root structure on lower incisors, canines and premolars.
- Longitudinal serrations contact crown and root providing enhanced contact on the tooth reducing possible damage to crown or root

PACKAGE INCLUDE:

- 1 Pc Extraction Forceps #150
- 1 Pc Extraction Forceps #151
- 1 Pc Extraction Forceps #23

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Extraction Forceps**WI-70004252****3 Dental Extracting Forceps #17 Lower Molars****\$29.99** per pieceTrade price from **\$27.59**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Premium Quality American Pattern Dental Extracting Forceps #17 For Lower Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#17 xtraction Forceps is used to extract Lower Molars. The Working end of this tooth extraction instrument is designed with a sharp bend, making sure it reaches deep into the socket for a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.


EXTRACTING FORCEPS #17 FEATURES

- Tri-pointed Beak Design

- For Lower Molars
- Non-slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 3 Pcs Dental Extracting Forceps #17 – Lower Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test,

**WI-70004238****3 Dental Extracting Forceps #18R for Upper 1st & 2nd Molars****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Extracting Forceps #18R for Upper 1st & 2nd Molars, Right Pointed Tip Design with Hook Handle German Grade Stainless Steel Teeth Extraction Surgical Dental Instruments

American Extracting forceps No. 18R hook handle are the perfect tool to remove teeth from the alveolar bone. Extraction forceps #18R has right pointed beak design features an American Pattern. Bifurcated pointed beaks designed for upper first and second molars. Hooked handle provides increased leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you

EXTRACTION FORCEPS #18R FEATURES

- Right Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE


- 3 Pcs American Extraction Forceps #18R Right Pointed Beak with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction

Surgery Procedure

BEAK/JAWS: Right Pointed Tip Design**WI-70004227****3 Dental Extracting Forceps #217 for 1st & 2nd Lower Molars****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 American Pattern Dental Extracting Forceps #217 For Lower Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#217 Extraction Forceps is used to extract lower molars. The working end of this tooth extraction instrument is designed with a tri-pointed shaped beaks, making sure it reaches deep into the socket for a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

- Tri-pointed Serrated Beak Design
- For 1st & 2nd Lower Molars
- Non-slip Grip Handle
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE

- 3 Pcs Dental Extracting Forceps #217 – Lower Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated

**3
PCS****WI-70004250****3 Dental Extracting Forceps #222 Lower 3rd Molars****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Extracting Forceps #222 For Lower 3rd Molars Premium Quality Stainless Steel American Pattern Apical Beaks Teeth Extraction Dental Surgical Instruments

American Extraction Forceps No. 222 are designed to be used for extracting lower wisdom teeth. Extracting forceps #222 designed with a sharp apical beaks make complete root contact, cross serration give nonslip surface, quickly enter sub-gingival. It ensures no fracture or tissue damage and easy extraction. These dental forceps are perfect to be a part of your dental kit. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #222 FEATURES


- Serrated Apical Beak.
- For Lower 3rd Molars.
- Non-slip Grip Handle
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE

- 3 Pcs Dental Extracting Forceps #222 – Lower 3rd Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**3
PCS****WI-70004219****3 Dental Extracting Forceps #88L 1st & 2nd Upper Molar****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Extracting Forceps #88L with Left Pointed Split Beak Grip Handle for 1st & 2nd Upper Molars German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American pattern dental extraction forceps No. 88L used to extract 1st and 2nd molars. This extraction forceps has left pointed split beak design. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #88L FEATURES

- For Upper 1st & 2nd Molars.

- Left Pointed, Split Beak Design.
- Non-Slip Grip Handle.
- American Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 3 Pcs American Pattern Extracting Forceps #88L for 1st & 2nd Upper Molar

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Left Pointed Split Design**POLISHED:** Yes**HANDLE:** Non-Slip Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION**

**3
PCS****WI-70004230****3 Dental Extracting Forceps 18L for Upper 1st & 2nd Molars****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Extracting Forceps #18L for Upper 1st & 2nd Molars, Left Pointed Tip Design with Hook Handle German Grade Stainless Steel Teeth Extraction Surgical Dental Instruments

American Extracting forceps No. 18L hook handle are the perfect tool to remove teeth from the alveolar bone. Extraction forceps #18L has left pointed beak design features an American Pattern. Bifurcated pointed beaks designed for upper first and second molars. Hooked handle provides increased leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you

EXTRACTION FORCEPS #18L FEATURES

- Left Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 3 Pcs American Extraction Forceps #18L Left Pointed Beak with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction

Surgery Procedure

BEAK/JAWS: Left Pointed Tip Design**3
PCS****WI-70004232****3 Dental Extracting Kells Forceps #99C****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Extracting Kells Forceps #99C for Upper Incisors, Canines & Premolars German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument

Kells extraction forceps No. 99 used for extracting teeth from the alveolar bone. The Short, straight conical beaks designed to adapt to upper incisors, canines and premolars. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

KELLS EXTRACTING FORCEPS #99C FEATURES

- For Upper Incisors, Canines & Premolars.
- Short, Straight Conical Beaks for Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 3 Pcs Kells Extracting Forceps #99C – Serrated Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Short Serrated Beaks**POLISHED:** Yes**HANDLE:** Non-Slip Handle**REUSABLE:** Yes

**WI-70004243****3 Dental Extraction Forceps #150 Upper Incisors, Root Teeth****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Extracting Forceps #150 Upper Incisors, Canines, and Premolars German Grade Stainless Steel Teeth Extraction Forceps Dental Surgical Instruments.

Universal Extraction Forceps No. 150 are used for extracting teeth from the alveolar bone. Universal conical beaks designed to adapt to upper incisors, canines, premolars, roots, universal right or left. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

UNIVERSAL EXTRACTION FORCEPS #150 FEATURES

- For Upper Incisors, Canines & Premolars.
- Serrated Conical Beaks for Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 3 Pcs Universal Extraction Forceps #150 – Serrated Beaks

TECHNICAL INFORMATION

MATERIAL: High-Quality Medical Grade Stainless Steel

INSTRUMENTS TYPE: Dental Extracting Forceps

APPLICATION: Extraction Forceps

BEAK/JAWS: Serrated

POLISHED: Yes

STERILE: Yes

**WI-70004244****3 Dental Tooth Extracting Forceps #10S with Serrated Jaws****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Heavy Duty American Pattern Tooth Extracting Forceps Fig. 10S Upper Molar Serrated Curved Jaws German Grade Stainless Steel Dental Instrument.

Durable and sturdy, Upper molars extracting forceps #10S are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals dentists.

DENTAL EXTRACTING FORCEPS #10S FEATURES

- Serrated Beak For Better Grip
- Figure 10S For Upper Molar
- Smoothly Curved Tips Making Them Suitable For Both Sides
- Non-Slip Grip Handle

- German Stainless Steel

PACKAGE INCLUDE

- 3 Pcs Extracting Forceps #10S Curved Serrated – Upper Molars

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Dental Extracting Forceps

APPLICATION: Teeth Extraction

POLISHED: Yes

STERILE: Yes

REUSABLE: Yes

AUTOCCLAVED: Yes

RUST PREVENTION

PROCEDURE: Passivated

ULTRASONIC CLEANED: Yes

TESTS PERFORMED: Boil Test, Performance Test, Shape Test



**WI-70004242****3 Extracting Forceps #103 for Lower Bicuspid, Incisors & Roots****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Universal Extracting Forceps #103 for Lower Bicuspid, Incisors & Roots with Hook Handle Premium Quality Stainless Steel American Pattern Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No.103 with hook handle is used to extract teeth from the alveolar bone. Dental extracting forceps #103 features universal serrated precise beaks designed to adapt to upper & lower bicuspid, incisors & roots. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

UNIVERSAL EXTRACTION FORCEPS #103 FEATURES

- Precise Serrated Tips Meet at a Point, Ensuring Optimal Grasp.
- For Lower Bicuspid Incisor and Roots.
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 3 Pcs American Pattern Universal Extraction Forceps #103

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Serrated Beaks**WI-70004235****3 Extracting Forceps #150A for Upper Incisors, Cuspids and Bicuspid****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Universal Extracting Forceps #150A Parallel Beaks with Plain Edge for Presidential Upper Incisors, Cuspids and Bicuspid Premium Quality Stainless Steel American Pattern Teeth Extraction Dental Surgical Instruments.

#150A maxillary forceps is a universal, parallel conical beaks designed to adapt to upper incisors, canines and premolars used for extracting teeth from the alveolar bone. 150A parallel beaks allow for increased contact beneath the crown for right or left extractions. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction and easy gripping handles ensures maximum durability and dependability. You can reuse t150a

extracting forceps over time as they will not bend or tarnish on you.

UNIVERSAL EXTRACTION FORCEPS #150A FEATURES

- For Upper Incisors, Cuspids and Bicuspid
- Parallel Beaks with Plain Edge for Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 3 Pcs American Extraction Forceps #151A – Parallel Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps

**WI-70004217****3 Extracting Forceps #203 Lower Incisors, Canines, Premolars & Roots****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Extracting Forceps #203 for Lower Incisors, Canines, Premolars, and Roots Premium Quality Stainless Steel American Pattern Universal Wider Conical Beaks Teeth Extraction Surgical Dental Instruments.

American extracting forceps No. 203 used for extracting teeth from the alveolar bone. Extraction forceps #203 features a universal wider conical beaks designed to adapt to lower incisors, canines, premolars and roots, may be used left or right. American Pattern forceps have a smooth hinge so no pin is visible. All our extraction forceps can be serrated unless otherwise indicated.

AMERICAN EXTRACTING FORCEPS #203 FEATURES

- Wider Beaks for Better Grip.
- For Lower Incisors, Cuspids, Bicuspid & Root.
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 3 Pcs American Pattern Extraction Forceps #203 Wider Conical Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Conical Beaks**POLISHED:** Yes**STERILE:** Yes**WI-70004223****3 Extracting Forceps #210S Upper 3rd Molars, Universal****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Extracting Forceps #210S for Upper 3rd Molars Premium Quality Stainless Steel American Pattern Bayonet Style Serrated Beaks Teeth Extraction Surgical Dental Instruments.

American No.210S are used for extracting upper wisdom teeth and root fragments from the alveolar bone. Extraction forceps #210S features offset (Bayonet) conical beaks designed to adapt to upper third molars with longitudinal serrations to improve grasp. The unique design help in removing a tooth with a minor injury. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

NOTE: 210S is the same design as 210H without hooked handle.

AMERICAN EXTRACTION FORCEPS #210S FEATURES

- Offset (Bayonet) Conical Beaks Design.
- For Upper 3rd Molars.
- Longitudinal Serrations Improve Grasp.
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 3 Pc American Pattern Extraction Forceps #210S for Upper Molars

TECHNICAL INFORMATION

**WI-70004222****3 Extracting Forceps #213 for Upper Central & Canines Roots****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Extracting Forceps #213 for Upper Central & Canines Roots with Hook Handle Premium Quality Stainless Steel American Pattern Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 213 with hook handle is used to extract teeth from the alveolar bone. American extraction forceps #213 features conical beaks designed to adapt to upper central & canines roots with to improve grasp. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**AMERICAN PATTERN
EXTRACTION FORCEPS #213****FEATURES**

- Conical Beaks Design.
- Upper Central & Canines Roots
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 3 Pcs American Pattern Extraction Forceps #213 with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Conical Beaks**POLISHED:** Yes**WI-70004233****3 Extracting Forceps #286 for Upper Roots, Incisors & Premolars****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Extracting Forceps #286 for Upper Anterior Roots, Incisors & Bicuspid Premium Quality Stainless Steel American Pattern Bayonet Beaks Teeth Extraction Surgical Dental Instruments.

American extracting forceps No. 286 used for extracting teeth from the alveolar bone. Extraction forceps #286 features offset (Bayonet) long, narrow beaks designed to adapt to upper roots, incisors and premolars. American Pattern forceps have a smooth hinge so no pin is visible. All our extraction forceps can be serrated unless otherwise indicated.

**AMERICAN EXTRACTING
FORCEPS #286 FEATURES**

- Bayonet Pattern Long Narrow Tips

- For Upper Anterior Roots, Incisors & Bicuspid
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 3 Pcs American Pattern Extraction Forceps #286 Bayonet Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Bayonet Beaks**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes

**3
PCS****WI-70004239****3 Extracting Forceps #6 Upper & Lower Teeth Splitting & Separating****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 American Extraction Forceps #6 for Upper and Lower Teeth Splitting & Separating Premium Quality Stainless Steel Teeth Extraction Surgical Dental Instruments.

Extracting Forceps #6 are precision atraumatic extraction forceps that provide optimum tactical feel while reducing hand fatigue. The tapered tips adapt to the tooth with excellent subgingival contact. These dental forceps are perfect to be a part of your dental kit. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #6
FEATURES**

- Blade-like Beaks for Splitting & Separating.
- For Upper & Lower Teeth.
- 6.5" (16.5cm) Overall Length.
- Non-Slip Textured-grip Handle.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 3 Pcs American Extracting Forceps #6 for Upper and Lower Teeth Splitting

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**3
PCS****WI-70004229****3 Extracting Forceps #69 Upper Lower Roots and Fragments****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Extracting Forceps #69 with Long Conical Beaks for Upper and Lower Roots and Fragments – American Pattern – German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

American pattern extracting forceps 69 are designed to remove upper and lower dental fragments and roots, particularly in the posterior regions. They feature long precise beaks, a smooth hinge without a visible pin, and are made of stainless steel for durability. These forceps are fully autoclavable and reusable.

**AMERICAN EXTRACTING
FORCEPS #69 FEATURES**

- For Upper and Lower Roots and Fragments.

- Long, Precise Conical Beaks
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 3 Pcs Dental Extracting Forceps #69

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Tooth Extraction

Forceps

STERILE: Non-Sterile**REUSABLE/AUTOCCLAVABLE:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test

**3
PCS****WI-70004218****3 Extracting Forceps #88R For Right 1st And 2nd Upper Molar****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Extracting Forceps #88R For Right 1st And 2nd Upper Molar American Pattern Right Pointed, Serrated Split Beak Design German Grade Stainless Steel Tooth Extraction Dental Surgical Instrument.

Extraction Forceps No. 88R design to remove teeth from bony socket.

Suitable for maxillary right first and second Molars. Presidential forceps for 1st and Upper 2nd Molars. Bayonet design split beak for engaging palatal root. Each beak has different design to adapt to the maxillary molar roots that differ anatomically on the facial and lingual. Beak with one projection contours to bifurcation of mesial-buccal and distal-buccal roots. Beak with two projections contours to lingual root. American pattern forceps have a

smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

AMERICAN EXTRACTION FORCEPS #88R FEATURES

- For Upper 1st & 2nd Molars.
- Right Pointed, Serrated Split Beak Design
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 3 Pcs American Extraction Extracting Forceps #88R

TECHNICAL INFORMATION**3
PCS****WI-70004237****3 Extracting Forceps #99A Upper Incisors, Canines & Premolar Root****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Extracting Forceps #99A Straight Serrated Beaks for Upper Incisors, Canines, Premolar Root with Hook Handle Premium Quality American Pattern Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 99A used to extract upper incisor, cuspid and bicuspid from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #99A FEATURES

- Straight Serrated Beaks.
- For Upper Incisors, Canines, Premolar Root.

- Hook Handle.
- American Pattern.
- Made Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 3 Pcs American Pattern Extracting Forceps #99A with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Extraction Forceps**BEAK/JAWS:** Straight Serrated Beaks**POLISHED:** Yes**HANDLE:** Non-Slip Hook Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated

**WI-70004221**

3 Extracting Forceps 3FS Dental Surgical Instruments

\$29.99 per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Extracting Forceps #3FS Woodward's American Pattern Tooth Extraction Forceps Fig. 3FS For Lower Molars German Grade Stainless Steel Dental Surgical Instruments

Extracting forceps #3FS used for lower molars extracting teeth from the alveolar bones. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by dental professionals.

DENTAL EXTRACTING FORCEPS #3FS FEATURES

- Pointed Beak Design
- For Lower Molars
- Non-Slip Grip Handle

- German Stainless Steel

PACKAGE INCLUDE

- 3 Pcs Woodward Dental Extracting Forceps 3FS – Lower Molars

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**WI-70004226**

3 Extracting Forceps Dental Surgical 1 Hen

\$29.99 per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Extracting Forceps 1 HEN, Upper Incisors Henahan Teeth Extraction German Grade Stainless Steel Surgical Dental Instrument

Extraction forceps are the most familiar instruments with regard to the extraction of tooth teeth. Intra-alveolar extractions are mainly accomplished using elevators and extraction forceps. They are the instruments that are used to remove the tooth out of the alveolar bone socket.

Doctor use this tool for upper incisors, #1HEN Henahan. Our tools are professional quality and used by professionals. All Dental Forceps are handcrafted from german grade stainless steel for logtivity and durability

EXTRACTING FORCEPS 1 HEN FEATURES

- #1HEN Henahan Dental Froceps
- Upper Incisors and Canines Forceps
- Non-Slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 3 Pcs Extracting Forceps 1 HEN

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated

**WI-70004213****3 Pedo Extracting Forceps #151S Dental Surgical Instrument****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pedo Extraction Forceps #151S for Children Lower Primary Teeth and Roots Stainless Steel American Pattern Tooth Extraction Dental Surgical Instruments

Universal Pedo Forceps No. 151S used to extract lower primary teeth of children. The small, closely-situated beaks ensures maximum accessibility and a tight grasp of the tooth, while minimizing chances of trauma to the patient.

151S Pedo forceps guarantees the easiest tooth extraction procedure. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these

extracting forceps over time as they will not bend or tarnish on you.

AMERICAN EXTRACTION FORCEPS #151S FEATURES

- For Lower Primary Teeth and Roots
- Small, Conical Beaks for Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 3 Pcs American Extraction Forceps #151S

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Pedo Extraction**WI-70004247****3 Universal Dental Extracting Forceps No. 24 for Upper Molars****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Universal Extracting Forceps Set #24 for Upper Molars with Wide Rounded Beaks and Hook Handle Premium Quality American Pattern Tooth Extraction Forceps for both Human and Small Animals Dentistry Instrument.

Universal extracting forceps #24 are designed for the extraction of upper molars, including first, second, and sometimes third molars. Their beaks securely grip the teeth, allowing for controlled extraction while minimizing trauma to surrounding tissues. They feature an American pattern design with wide, rounded, serrated beaks and a hook handle for enhanced leverage, making them versatile for various upper molar extractions.

AMERICAN UNIVERSAL EXTRACTING FORCEPS #24 FEATURES

- For Upper Molars.
- Wide, Rounded Serrated Jaws.
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Autoclavable/Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 3 Pcs Universal Extraction Forceps #24 Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps

**WI-70004231****3 Universal Extraction Forceps #151 for Lower Teeth****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Universal Extraction Forceps #151 Serrated Beaks for Lower Incisors, Canines, Premolars & Roots – Premium Quality Stainless Steel American Pattern Tooth Extracting Surgical Dental Instrument.

#151 extracting forceps also known as "lower universal forceps" is a dental tool used for removing teeth in the lower jaw, including incisors, canines, premolars, and roots. American Pattern #151 forceps is designed to be versatile, enabling dentists to extract various teeth in the lower jaw. #151 dental forceps are designed with thin, tapered, and sharp parallel beaks that ensure good contact with the root and crown of the tooth. 151 dental extracting forceps feature longitudinal serrations for improved grip and atraumatic

extraction, and are typically made of surgical-grade stainless steel for durability and corrosion resistance.

AMERICAN EXTRACTION FORCEPS #151 FEATURES

- For Lower Incisors, Canines, Premolars & Roots.
- Thin, Tapered, and Sharp Parallel Beaks for Good Contact.
- Longitudinal Serrations for Improved Grip and Atraumatic Extraction.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 3 Pcs Universal Extraction Forceps #151

TECHNICAL INFORMATION**WI-70004253****3 Universal Pedo Extracting Forceps #17S for Lower Molars****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Universal Pedo Extraction Forceps #17S For Lower Molars with Tri-Pointed Beaks & Spring Handle Premium Quality German Grade Stainless Steel American Pattern Tooth Extraction Dental Surgical Instruments

Pedodontic 17S pediatric dental forceps are specifically designed for extracting lower molars, featuring sharp tri-pointed tips that ensure a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #17S FEATURES

- Universal Pedo Extraction Forceps.

- Tri-pointed Beak Design.
- For Lower Molars.
- Non-slip Grip Spring Handle.
- German Stainless Steel.

PACKAGE INCLUDE

- 3 Pcs American Pattern Pedo Extracting Forceps #17S – Lower Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Childern Teeth Extraction**STERILE:** Non-Sterile**REUSABLE/AUTOCCLAVABLE:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-25943254

3.5 Inch Stainless Steel Dental Backhaus Towel Forceps

\$2.50 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

3.5 Inch Stainless Steel Dental Backhaus Towel Forceps is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 3.5" size shown is our standard; other sizes in the same pattern are made to order.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.50, falling to \$1.50 at 10 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-42865573

5 Piece Dental Extraction Forceps Set with Gold

\$30.00 per piece

Trade price from \$16.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture 5 Piece Dental Extraction Forceps Set with Gold for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$30.00, falling to \$16.00 at 1 set — about 47% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

**6
PCS****WI-70004200****6 American Extracting Forceps #15 Lower Molars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Extracting Forceps #15 for 1st & 2nd Lower Molars with Hook Handle, Universal German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American pattern dental extracting forceps No. 15 used to extract 1st and 2nd lower molars. This extraction forceps has tri-pointed beak design. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**EXTRACTING FORCEPS #15
FEATURES**

- For Lower 1st & 2nd Molars.

- Tri-Pointed Beaks Design.
- Non-Slip Grip Hook Handle.
- American Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 6 Pcs American Pattern Extracting Forceps #15 Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction Surgery Procedure**BEAK/JAWS:** Tri-Pointed Design**POLISHED:** Yes**HANDLE:** Hook Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****6
PCS****WI-70003257****6 American Extracting Forceps #16 Lower 1st & 2nd Molars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Universal Extracting Forceps #16 Hook Handle for Lower 1st & 2nd Molars Cowhorn Beak Design American Pattern German Grade Stainless Steel Tooth Extraction Surgical Dental Instruments

American Extracting forceps No. 16 hook handle are the perfect tool to remove teeth from the alveolar bone. This universal extraction forceps #16 has pointed cowhorn beak design features an American Pattern.

Bifurcated pointed beaks designed for lower first and second molars. Hooked handle provides increased leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over

time as they will not bend or tarnish on you

**AMERICAN EXTRACTION
FORCEPS #16 FEATURES**

- Cow Horn Beak Design.
- For Lower 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 6 Pcs American Extraction Forceps #16 Cowhorn Beak with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction

**6
PCS****WI-70004194****6 Dental Extracting Forceps #10H Upper Molars Universal****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Extracting Forceps #10H Universal Upper Molars with Hook Handle Premium Quality Stainless Steel American Pattern Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 10 with hook handle is used to extract upper molars. American extraction forceps #10H features conical beaks designed to adapt to upper molars to improve grasp. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**AMERICAN PATTERN
EXTRACTION FORCEPS #10H
FEATURES**

- Conical Beaks Design.
- Upper Molars, Universal
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 6 Pcs American Pattern Extraction Forceps #10H with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Conical Serrated Beaks**POLISHED:** Yes**STERILE:** Yes**6
PCS****WI-70004187****6 Dental Extracting Forceps #150AS Upper Anteriors****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 American Pattern Dental Extracting Forceps #150AS Serrated Split Beaks For Upper Anteriors German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Extraction Forceps No. 150AS used for extracting teeth from the alveolar bone. Upper universal split serrated beaks designed to provide increased grasp beneath the crown to prevent crown or root breakage. Serrations provide increased mechanical grasp on the root without increased pressure. 150AS dental forceps can be used left and right.

American pattern forceps have a smooth hinge so no pin is visible. The stainless steel construction ensures maximum durability and dependability. You can reuse these 150AS extracting


forceps over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

**AMERICAN EXTRACTING
FORCEPS #150AS FEATURES**

- For Upper Anteriors.
- Serrated, Split Beaks Provide an Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 6 Pcs Dental Extracting Forceps #150AS Split Beaks – Serrated

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

**WI-70004211****6 Dental Extracting Forceps #17 Lower Molars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Premium Quality American Pattern Dental Extracting Forceps #17 For Lower Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#17 xtraction Forceps is used to extract Lower Molars. The Working end of this tooth extraction instrument is designed with a sharp bend, making sure it reaches deep into the socket for a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #17 FEATURES

- Tri-pointed Beak Design

- For Lower Molars
- Non-slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 6 Pcs Dental Extracting Forceps #17 – Lower Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test,**WI-70004186****6 Dental Extracting Forceps #217 for 1st & 2nd Lower Molars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 American Pattern Dental Extracting Forceps #217 For Lower Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#217 Extraction Forceps is used to extract lower molars. The working end of this tooth extraction instrument is designed with a tri-pointed shaped beaks, making sure it reaches deep into the socket for a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #217 FEATURES

- Tri-pointed Serrated Beak Design
- For 1st & 2nd Lower Molars
- Non-slip Grip Handle
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE

- 6 Pcs Dental Extracting Forceps #217 – Lower Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated

**6
PCS****WI-70004209****6 Dental Extracting Forceps #222 Lower 3rd Molars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Extracting Forceps #222 For Lower 3rd Molars Premium Quality Stainless Steel American Pattern Apical Beaks Teeth Extraction Dental Surgical Instruments

American Extraction Forceps No. 222 are designed to be used for extracting lower wisdom teeth. Extracting forceps #222 designed with a sharp apical beaks make complete root contact, cross serration give nonslip surface, quickly enter sub-gingival. It ensures no fracture or tissue damage and easy extraction. These dental forceps are perfect to be a part of your dental kit. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #222 FEATURES

- Serrated Apical Beak.
- For Lower 3rd Molars.
- Non-slip Grip Handle
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE

- 6 Pcs Dental Extracting Forceps #222 – Lower 3rd Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**6
PCS****WI-70004176****6 Dental Extracting Forceps #88L 1st & 2nd Upper Molar****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Extracting Forceps #88L with Left Pointed Split Beak Grip Handle for 1st & 2nd Upper Molars German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American pattern dental extraction forceps No. 88L used to extract 1st and 2nd molars. This extraction forceps has left pointed split beak design. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #88L FEATURES

- For Upper 1st & 2nd Molars.

- Left Pointed, Split Beak Design.
- Non-Slip Grip Handle.
- American Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 6 Pcs American Pattern Extracting Forceps #88L for 1st & 2nd Upper Molar

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Left Pointed Split Design**POLISHED:** Yes**HANDLE:** Non-Slip Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION**

**6
PCS****WI-70004189****6 Dental Extracting Forceps #18L for Upper 1st & 2nd Molars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Extracting Forceps #18L for Upper 1st & 2nd Molars, Left Pointed Tip Design with Hook Handle German Grade Stainless Steel Teeth Extraction Surgical Dental Instruments

American Extracting forceps No. 18L hook handle are the perfect tool to remove teeth from the alveolar bone. Extraction forceps #18L has left pointed beak design features an American Pattern. Bifurcated pointed beaks designed for upper first and second molars. Hooked handle provides increased leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you

EXTRACTION FORCEPS #18L FEATURES

- Left Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 6 Pcs American Extraction Forceps #18L Left Pointed Beak with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction

Surgery Procedure

BEAK/JAWS: Left Pointed Tip Design**6
PCS****WI-70004191****6 Dental Extracting Kells Forceps #99C****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Extracting Kells Forceps #99C for Upper Incisors, Canines & Premolars German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument

Kells extracting forceps No. 99 used for extracting teeth from the alveolar bone. The Short, straight conical beaks designed to adapt to upper incisors, canines and premolars. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

KELLS EXTRACTING FORCEPS #99C FEATURES

- For Upper Incisors, Canines & Premolars.
- Short, Straight Conical Beaks for Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 6 Pcs Kells Extracting Forceps #99C – Serrated Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Short Serrated Beaks**POLISHED:** Yes**HANDLE:** Non-Slip Handle**REUSABLE:** Yes

**6
PCS****WI-70004202****6 Dental Extraction Forceps #150 Upper Incisors, Root Teeth****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Extracting Forceps #150 Upper Incisors, Canines, and Premolars German Grade Stainless Steel Teeth Extraction Forceps Dental Surgical Instruments.

Universal Extraction Forceps No. 150 are used for extracting teeth from the alveolar bone. Universal conical beaks designed to adapt to upper incisors, canines, premolars, roots, universal right or left. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

UNIVERSAL EXTRACTION FORCEPS #150 FEATURES

- For Upper Incisors, Canines & Premolars.
- Serrated Conical Beaks for Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 6 Pcs Universal Extraction Forceps #150 – Serrated Beaks

TECHNICAL INFORMATION

MATERIAL: High-Quality Medical Grade Stainless Steel

INSTRUMENTS TYPE: Dental Extracting Forceps

APPLICATION: Extraction Forceps

BEAK/JAWS: Serrated

POLISHED: Yes

STERILE: Yes

**6
PCS****WI-70004201****6 Extracting Forceps #103 for Lower Bicuspid, Incisors & Roots****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Universal Extracting Forceps #103 for Lower Bicuspid, Incisors & Roots with Hook Handle Premium Quality Stainless Steel American Pattern Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 103 with hook handle is used to extract teeth from the alveolar bone. Dental extracting forceps #103 features universal serrated precise beaks designed to adapt to upper & lower bicuspid, incisors & roots. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

UNIVERSAL EXTRACTION FORCEPS #103 FEATURES

- Precise Serrated Tips Meet at a Point, Ensuring Optimal Grasp.
- For Lower Bicuspid Incisor and Roots.
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 6 Pcs American Pattern Universal Extraction Forceps #103

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Dental Extracting Forceps

APPLICATION: Teeth Extraction Forceps

BEAK/JAWS: Serrated Beaks



**WI-70004195****6 Extracting Forceps #150A for Upper Incisors, Cuspids and Bicuspid****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Universal Extracting Forceps #150A Parallel Beaks with Plain Edge for Presidential Upper Incisors, Cuspids and Bicuspid Premium Quality Stainless Steel American Pattern Teeth Extraction Dental Surgical Instruments.

#150A maxillary forceps is a universal, parallel conical beaks designed to adapt to upper incisors, canines and premolars used for extracting teeth from the alveolar bone. 150A parallel beaks allow for increased contact beneath the crown for right or left extractions. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction and easy gripping handles ensures maximum durability and dependability. You can reuse t150a

extracting forceps over time as they will not bend or tarnish on you.

UNIVERSAL EXTRACTING FORCEPS #150A FEATURES

- For Upper Incisors, Cuspids and Bicuspid
- Parallel Beaks with Plain Edge for Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 6 Pcs American Extraction Forceps #151A – Parallel Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**WI-70004174****6 Extracting Forceps #203 Lower Incisors, Canines, Premolars & Roots****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Extracting Forceps #203 for Lower Incisors, Canines, Premolars, and Roots Premium Quality Stainless Steel American Pattern Universal Wider Conical Beaks Teeth Extraction Surgical Dental Instruments.

American extracting forceps No. 203 used for extracting teeth from the alveolar bone. Extraction forceps #203 features a universal wider conical beaks designed to adapt to lower incisors, canines, premolars and roots, may be used left or right. American Pattern forceps have a smooth hinge so no pin is visible. All our extraction forceps can be serrated unless otherwise indicated.

- Wider Beaks for Better Grip.
- For Lower Incisors, Cuspids, Bicuspid & Root.
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 6 Pcs American Pattern Extraction Forceps #203 Wider Conical Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Conical Beaks**POLISHED:** Yes**STERILE:** Yes

**WI-70004181****6 Extracting Forceps #210S Upper 3rd Molars, Universal****\$59.99** per pieceTrade price from **\$55.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Extracting Forceps #210S for Upper 3rd Molars Premium Quality Stainless Steel American Pattern Bayonet Style Serrated Beaks Teeth Extraction Surgical Dental Instruments.

American No.210S are used for extracting upper wisdom teeth and root fragments from the alveolar bone. Extraction forceps #210S features offset (Bayonet) conical beaks designed to adapt to upper third molars with longitudinal serrations to improve grasp. The unique design help in removing a tooth with a minor injury.

American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

NOTE: 210S is the same design as 210H without hooked handle.

AMERICAN EXTRACTION FORCEPS #210S FEATURES

- Offset (Bayonet) Conical Beaks Design.
- For Upper 3rd Molars.
- Longitudinal Serrations Improve Grasp.
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 6 Pc American Pattern Extraction Forceps #210S for Upper Molars

TECHNICAL INFORMATION**WI-70004180****6 Extracting Forceps #213 for Upper Central & Canines Roots****\$59.99** per pieceTrade price from **\$55.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Extracting Forceps #213 for Upper Central & Canines Roots with Hook Handle Premium Quality Stainless Steel American Pattern Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 213 with hook handle is used to extract teeth from the alveolar bone. American extraction forceps #213 features conical beaks designed to adapt to upper central & canines roots with to improve grasp. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

FEATURES

- Conical Beaks Design.
- Upper Central & Canines Roots
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 6 Pcs American Pattern Extraction Forceps #213 with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Conical Beaks**POLISHED:** Yes**AMERICAN PATTERN EXTRACTION FORCEPS #213**

**WI-70004192****6 Extracting Forceps #286 for Upper Roots, Incisors & Premolars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Extracting Forceps #286 for Upper Anterior Roots, Incisors & Bicuspid Premium Quality Stainless Steel American Pattern Bayonet Beaks Teeth Extraction Surgical Dental Instruments.

American extracting forceps No. 286 used for extracting teeth from the alveolar bone. Extraction forceps #286 features offset (Bayonet) long, narrow beaks designed to adapt to upper roots, incisors and premolars. American Pattern forceps have a smooth hinge so no pin is visible. All our extraction forceps can be serrated unless otherwise indicated.

AMERICAN EXTRACTING FORCEPS #286 FEATURES

- Bayonet Pattern Long Narrow Tips

- For Upper Anterior Roots, Incisors & Bicuspid
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 6 Pcs American Pattern Extraction Forceps #286 Bayonet Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Bayonet Beaks**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**WI-70004197****6 Extracting Forceps #6 Upper & Lower Teeth Splitting & Separating****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 American Extracting Forceps #6 for Upper and Lower Teeth Splitting & Separating Premium Quality Stainless Steel Teeth Extraction Surgical Dental Instruments.

Extracting Forceps #6 are precision atraumatic extraction forceps that provide optimum tactical feel while reducing hand fatigue. The tapered tips adapt to the tooth with excellent subgingival contact. These dental forceps are perfect to be a part of your dental kit. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #6 FEATURES

- Blade-like Beaks for Splitting & Separating.
- For Upper & Lower Teeth.
- 6.5" (16.5cm) Overall Length.
- Non-Slip Textured-grip Handle.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 6 Pcs American Extracting Forceps #6 for Upper and Lower Teeth Splitting

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes

**WI-70004188****6 Extracting Forceps #69 Upper Lower Roots and Fragments****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Extracting Forceps #69 with Long Conical Beaks for Upper and Lower Roots and Fragments – American Pattern – German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

American pattern extracting forceps 69 are designed to remove upper and lower dental fragments and roots, particularly in the posterior regions. They feature long precise beaks, a smooth hinge without a visible pin, and are made of stainless steel for durability. These forceps are fully autoclavable and reusable.


AMERICAN EXTRACTING FORCEPS #69 FEATURES

- For Upper and Lower Roots and Fragments.

- Long, Precise Conical Beaks
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 6 Pcs Dental Extracting Forceps #69

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Tooth Extraction Forceps**STERILE:** Non-Sterile**REUSABLE/AUTOCCLAVABLE:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**WI-70004175****6 Extracting Forceps #88R For Right 1st And 2nd Upper Molar****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Extracting Forceps #88R For Right 1st And 2nd Upper Molar American Pattern Right Pointed, Serrated Split Beak Design German Grade Stainless Steel Tooth Extraction Dental Surgical Instrument.

Extraction Forceps No. 88R design to remove teeth from bony socket. Suitable for maxillary right first and second Molars. Presidential forceps for 1st and Upper 2nd Molars. Bayonet design split beak for engaging palatal root. Each beak has different design to adapt to the maxillary molar roots that differ anatomically on the facial and lingual. Beak with one projection contours to bifurcation of mesial-buccal and distal-buccal roots. Beak with two projections contours to lingual root. American pattern forceps have a


smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

AMERICAN EXTRACTION FORCEPS #88R FEATURES

- For Upper 1st & 2nd Molars.
- Right Pointed, Serrated Split Beak Design
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 6 Pcs American Extraction Extracting Forceps #88R

TECHNICAL INFORMATION

**WI-70004196****6 Extracting Forceps #99A Upper Incisors, Canines & Premolar Root****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Extracting Forceps #99A Straight Serrated Beaks for Upper Incisors, Canines, Premolar Root with Hook Handle Premium Quality American Pattern Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 99A used to extract upper incisor, cuspid and bicuspid from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #99A FEATURES

- Straight Serrated Beaks.
- For Upper Incisors, Canines, Premolar Root.

- Hook Handle.
- American Pattern.
- Made Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 6 Pcs American Pattern Extracting Forceps #99A with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Straight Serrated Beaks**POLISHED:** Yes**HANDLE:** Non-Slip Hook Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**WI-70004179****6 Extracting Forceps 3FS Dental Surgical Instruments****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Extracting Forceps #3FS Woodward's American Pattern Tooth Extraction Forceps Fig. 3FS For Lower Molars German Grade Stainless Steel Dental Surgical Instruments

Extracting forceps #3FS used for lower molars extracting teeth from the alveolar bones. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by dental professionals.

DENTAL EXTRACTING FORCEPS #3FS FEATURES

- Pointed Beak Design
- For Lower Molars
- Non-Slip Grip Handle

- German Stainless Steel

PACKAGE INCLUDE

- 6 Pcs Woodward Dental Extracting Forceps 3FS – Lower Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test

**WI-70004184****6 Extracting Forceps Dental Surgical 1 Hen****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Extracting Forceps 1 HEN, Upper Incisors Henahan Teeth Extraction German Grade Stainless Steel Surgical Dental Instrument

Extraction forceps are the most familiar instruments with regard to the extraction of tooth teeth. Intra-alveolar extractions are mainly accomplished using elevators and extraction forceps. They are the instruments that are used to remove the tooth out of the alveolar bone socket.

Doctor use this tool for upper incisors, #1HEN Henahan. Our tools are professional quality and used by professionals. All Dental Forceps are handcrafted from german grade stainless steel for logtivity and durability

EXTRACTING FORCEPS 1 HEN FEATURES

- #1HEN Henahan Dental Froceps
- Upper Incisors and Canines Forceps
- Non-Slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 6 Pcs Extracting Forceps 1 HEN

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**WI-70004170****6 Pedro Extracting Forceps #151S Dental Surgical Instrument****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pedro Extraction Forceps #151S for Childern Lower Primary Teeth and Roots Stainless Steel American Pattern Tooth Extraction Dental Surgical Insruments

Universal Pedro Forceps No. 151S used to extract lower primary teeth of children. The small, closely-situated beaks ensures maximum accessibility and a tight grasp of the tooth, while minimizing chances of trauma to the patient.

151S Pedro forceps guarantees the easiest tooth extraction procedure. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these

extracting forceps over time as they will not bend or tarnish on you.

AMERICAN EXTRACTION FORCEPS #151S FEATURES


- For Lower Primary Teeth and Roots
- Small, Conical Beaks for Optimal Grip.
- American Pattern.
- Non-ship Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 6 Pcs American Extraction Forceps #151S

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Pedro Extraction

**WI-70004206**

6 Universal Dental Extracting Forceps No. 24 for Upper Molars

\$59.99 per pieceTrade price from **\$55.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Universal Extracting Forceps Set #24 for Upper Molars with Wide Rounded Beaks and Hook Handle Premium Quality American Pattern Tooth Extraction Forceps for both Human and Small Animals Dentistry Instrument.

Universal extracting forceps #24 are designed for the extraction of upper molars, including first, second, and sometimes third molars. Their beaks securely grip the teeth, allowing for controlled extraction while minimizing trauma to surrounding tissues. They feature an American pattern design with wide, rounded, serrated beaks and a hook handle for enhanced leverage, making them versatile for various upper molar extractions.

AMERICAN UNIVERSAL EXTRACTING FORCEPS #24 FEATURES

- For Upper Molars.
- Wide, Rounded Serrated Jaws.
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Autoclavable/Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 6 Pcs Universal Extraction Forceps #24 Hook Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Extracting Forceps

**WI-70004190**

6 Universal Extraction Forceps #151 for Lower Teeth

\$59.99 per pieceTrade price from **\$55.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Universal Extraction Forceps #151 Serrated Beaks for Lower Incisors, Canines, Premolars & Roots – Premium Quality Stainless Steel American Pattern Tooth Extracting Surgical Dental Instrument.

#151 extracting forceps also known as "lower universal forceps" is a dental tool used for removing teeth in the lower jaw, including incisors, canines, premolars, and roots. American Pattern #151 forceps is designed to be versatile, enabling dentists to extract various teeth in the lower jaw. #151 dental forceps are designed with thin, tapered, and sharp parallel beaks that ensure good contact with the root and crown of the tooth. 151 dental extracting forceps feature longitudinal serrations for improved grip and atraumatic

extraction, and are typically made of surgical-grade stainless steel for durability and corrosion resistance.

AMERICAN EXTRACTION FORCEPS #151 FEATURES

- For Lower Incisors, Canines, Premolars & Roots.
- Thin, Tapered, and Sharp Parallel Beaks for Good Contact.
- Longitudinal Serrations for Improved Grip and Atraumatic Extraction.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 6 Pcs Universal Extraction Forceps #151

TECHNICAL INFORMATION





WI-70004212

6 Universal Pedo Extracting Forceps #17S for Lower Molars

\$59.99 per piece

Trade price from \$55.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Universal Pedo Extraction Forceps #17S For Lower Molars with Tri-Pointed Beaks & Spring Handle Premium Quality German Grade Stainless Steel American Pattern Tooth Extraction Dental Surgical Instruments

Pedodontic 17S pediatric dental forceps are specifically designed for extracting lower molars, featuring sharp tri-pointed tips that ensure a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #17S FEATURES

- Universal Pedo Extraction Forceps.

- Tri-pointed Beak Design.
- For Lower Molars.
- Non-slip Grip Spring Handle.
- German Stainless Steel.

PACKAGE INCLUDE

- 6 Pcs American Pattern Pedo Extracting Forceps #17S – Lower Molars

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Dental Extracting Forceps

APPLICATION: Children Teeth Extraction

STERILE: Non-Sterile

REUSABLE/AUTOCCLAVABLE: Yes

RUST PREVENTION

PROCEDURE: Passivated

ULTRASONIC CLEANED: Yes



WI-22083570

Adson Bayonet Forceps 16 cm Medical Dressing Serrated 1.7cm

\$5.00 per piece

Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Adson Bayonet Forceps 16 cm Medical Dressing Serrated 1.7cm, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: serrated jaws for grip that does not slip under tension.

The 16cm size shown is our standard; other sizes in the same pattern are made to order. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$5.00. Minimum order 25 pieces, and the price moves with quantity. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-00533624

Adson Brown Dissecting Forceps Teeth Surgical Tissue Medical

\$2.50 per piece

Trade price from \$1.90
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Adson Brown Dissecting Forceps Teeth Surgical Tissue Medical, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.50, falling to \$1.90 at 1 piece — about 24% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-43462069

Adson Tissue Forceps 1x2 Teeth Reusable Surgical

\$3.20 per piece

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Adson Tissue Forceps 1x2 Teeth Reusable Surgical is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: interlocking teeth for a positive hold on tough tissue, and full reusability — clean, autoclave and use again.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.20. Minimum order 1 piece, and the price moves with quantity.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-43417980

Adson Tissue Forceps 1x2 Teeth Stainless Steel Surgical Dressing

\$3.20 per piece

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Adson Tissue Forceps 1×2 Teeth Stainless Steel Surgical Dressing for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.20. Minimum order 1 piece, and the price moves with quantity.

Samples go out on request so you can check the steel and the finish before you place a real order.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-43465046

Adson Tissue Forceps 1x2 Teeth Stainless Steel Surgical Instrument

\$3.20 per piece

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Adson Tissue Forceps 1×2 Teeth Stainless Steel Surgical Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$3.20. Minimum order 1 piece, and the price moves with quantity.

Samples go out on request so you can check the steel and the finish before you place a real order.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-33489833

Adson Tissue Forceps Angled Teeth Precision Tips-Surgical

\$7.00 per piece

Trade price from \$4.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Adson Tissue Forceps Angled Teeth Precision Tips-Surgical, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: an angled shank that keeps the hand clear of the sightline, and interlocking teeth for a positive hold on tough tissue.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$7.00, falling to \$4.00 at 25 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-36617604

Allis Atraumatic Ergonomic 3-Jaw Holding Forceps Size

\$5.00 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Allis Atraumatic Ergonomic 3-Jaw Holding Forceps Size, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: an ergonomic handle for long lists without fatigue.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$5.00, falling to \$2.00 at 10 pieces — about 60% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.



WI-43452128

Allis Tissue Forceps 19cm 5x6 Teeth

\$3.20 per piece

Trade price from \$2.90
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Allis Tissue Forceps 19cm 5x6 Teeth is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

The 19cm size shown is our standard; other sizes in the same pattern are

made to order. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.20, falling to \$2.90 at 1 piece — about 9% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-43427329

Allis Tissue Forceps 5x6 Teeth

\$3.20 per piece

Trade price from \$2.90
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Allis Tissue Forceps 5x6 Teeth — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.20, falling to \$2.90 at 1 piece — about 9% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35100238

Allis Tissue Forceps Stainless Steel Instrument Set Dental

\$5.00 per piece

Trade price from \$2.90
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Allis Tissue Forceps Stainless Steel Instrument Set Dental is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$5.00, falling to \$2.90 at 10 pieces — about 42% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25901013

Allis Tissue Forceps Toothed Locking Pliers Intestinal Forceps Dental

\$5.00 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Allis Tissue Forceps Toothed Locking Pliers Intestinal Forceps Dental — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$5.00, falling to \$2.00 at 10 pieces — about 60% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70004168****American Extracting Forceps #151A****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Extracting Forceps #151A
Parallel Beaks for Lower Incisors,
Cuspids & Bicuspid Premium Quality
Stainless Steel American Pattern Tooth
Extraction Surgical Dental Instrument.


Extracting Forceps No. 151A are
specialized dental instruments
designed for the extraction of teeth
from the lower arch, particularly
incisors, canines, premolars, and roots.
They feature parallel beaks, offering a
firmer and more precise grip compared
to standard No. 151 forceps. American
pattern forceps have a smooth hinge so
no pin is visible. Made of surgical-
grade stainless steel for durability and
corrosion resistance.

**UNIVERSAL EXTRACTION
FORCEPS #151A FEATURES**

- For Lower Incisors, Cuspids &
Bicuspid
- Parallel Beaks for Firm Grip..
- American Pattern.
- Non-slip Gold Plated Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc TC Universal Extraction Forceps
#151A

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental
Extracting Forceps**APPLICATION:** Lower Teeth Extraction
Forceps**AUTOCCLAVABLE/ REUSABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION****PROCEDURE:** Passivated**WI-70004167****American Extracting Forceps #16 Lower 1st & 2nd Molars****\$11.99** per pieceTrade price from **\$11.03**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Extraction Forceps #16 Hook
Handle for Lower 1st & 2nd Molars
Cowhorn Beak Design American
Pattern German Grade Stainless Steel
Tooth Extraction Surgical Dental
Instruments

American extraction forceps No. 16
hook handle are the perfect tool to
remove teeth from the alveolar bone.
Extracting forceps #16 has pointed
cowhorn beak design features an
American Pattern. Bifurcated pointed
beaks designed for lower first and
second molars. Hooked handle
provides increased leverage when
approach is from behind the head.
American pattern forceps have a
smooth hinge so no pin is visible. You
can reuse these extracting forceps over

time as they will not bend or tarnish on
you

**UNIVERSAL COWHORN
EXTRACTION FORCEPS #16
FEATURES**

- Cow Horn Beak Design.
- For Lower 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless
Steel.

PACKAGE INCLUDE

- 1 Pc American Extraction Forceps
#16 Cowhorn Beak with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental
Extracting Forceps

**WI-70004166****American Extracting Forceps #16S Lower Molars****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pedo Extracting Forceps #16S for Lower Molars Cowhorn Beak Design American Pattern German Grade Stainless Steel Children Teeth Extraction Surgical Dental Instruments

Cowhorn Pedo Extraction Forceps No. 16S are the perfect tool to remove children teeth from the alveolar bone. American extracting forceps #16S has cowhorn pointed beaks designed for lower first and second molars. Hooked handle provides increased leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you. 100% satisfaction guarantee against defects in material and workmanship.

PEDO EXTRACTING FORCEPS #16S FEATURES

- Cow Horn Beak Design.
- For Lower 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel

PACKAGE INCLUDE

- 1 Pc American Pedo Extraction Forceps #16S Cowhorn Beak with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction

Surgery Procedure

BEAK/JAWS: Cowhorn Beaks Design**WI-70000972****American Extracting Forceps #23 with Gold Hollow Python Grip Handle****\$14.99** per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty Universal Extracting Forceps #23 Cow Horn 1st & 2nd Lower Molars with Gold Hollow Python Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#23 extraction forceps is used to extract teeth from 1st and 2nd lower molars with cow horn jaws. Forceps dimension is 7.5" (19cm) overall in length with gold hollow python grip handle. Durable and sturdy, our extraction forceps are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #23 FEATURES


- Slightly Larger Cowhorn Beak Design Than #16
- For Lower 1st & 2nd Molars
- Gold Hollow Python Grip Handle
- American Pattern
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Extracting Forceps #23 Cow Horn Jaws with Gold Hollow Python Grip Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes

**WI-70001099****American Extracting Forceps #88L 1st & 2nd Upper Molar****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps #88L with Left Pointed Split Beak Gold Python Grip Handle for 1st & 2nd Upper Molars German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American pattern dental extraction forceps No. 88L used to extract 1st and 2nd molars. This extraction forceps has left pointed split beak design.

American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #88L FEATURES

- For Upper 1st & 2nd Molars.
- Left Pointed, Split Beak Design.
- Gold Hollow Python Grip Handle.
- American Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc American Pattern Extracting Forceps #88L – Gold Python Grip Hollow Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Left Pointed Split Design**POLISHED:** Yes**HANDLE:** Gold Python Grip Hollow Handle**REUSABLE:** Yes**WI-70003258****American Extracting Forceps No. 15 for 1st & 2nd Lower Molars****\$11.99** per pieceTrade price from **\$11.03**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps No. 15 for 1st & 2nd Lower Molars with Hook Handle, Universal German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American pattern dental extracting forceps No. 15 used to extract 1st and 2nd lower molars. This extraction forceps has tri-pointed beak design. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS NO. 15 FEATURES

- For Lower 1st & 2nd Molars.

- Tri-Pointed Beaks Design.
- Non-Slip Grip Hook Handle.
- American Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc American Pattern Extracting Forceps #15 Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction

Surgery Procedure

BEAK/JAWS: Tri-Pointed Design**POLISHED:** Yes**HANDLE:** Hook Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION**

**WI-70004162****American Extracting Forceps No. 85A Lower Molars****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Extracting Forceps #85A for Lower Molars with Hook Handle
Premium Quality American Pattern
Tooth Extraction Forceps Dental
Instrument.

Extracting Forceps No. 85A are dental surgical instruments used for tooth extractions in the lower arch, specifically designed to securely grip lower molars, cuspids, and bicusps. They are known as universal forceps due to their anatomical design, featuring a hook handle for better grip and leverage during procedures. They feature an American pattern design with hook handle for enhanced leverage, making them versatile for various upper molar extractions.

**UNIVERSAL EXTRACTING
FORCEPS #85A FEATURES**

- For Lower Molars.
- Noserrated Jaws.
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Autoclavable/Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 3 Pcs Universal Extraction Forceps #85A Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental
Extracting Forceps**APPLICATION:** Teeth Extraction
Forceps**STERILE:** Non-Sterile**REUSABLE/AUTOCCLAVABLE:** Yes**WI-70004122****American Pattern Extracting Forceps #90
Upper Left Molars****\$11.99** per pieceTrade price from **\$11.03**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #90 for Upper Molars with Left Pointed Tip Design German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American pattern dental extraction forceps No. 90 used to extract upper left molars. This extraction forceps has left pointed tip design. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**EXTRACTING FORCEPS #90
FEATURES**

- For Upper Left Molars.

- Left Pointed Tip Design.
- Non-Slip Grip Handle.
- American Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc American Pattern Extracting Forceps #90 – Upper Left Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental
Extracting Forceps**APPLICATION:** Extraction Forceps**BEAK/JAWS:** Left Pointed Tip**POLISHED:** Yes**HANDLE:** Non-Slip Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated



WI-70004164

American Pedo Extracting Forceps #23S for Primary Lower Molars

\$11.99 per piece

Trade price from \$11.03
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pedo Extraction Forceps #23S
Cowhorn Beaks with Spring Handle for Primary Lower Molars American Pattern Premium Quality Stainless Steel Surgical Dental Instruments.

American pattern extracting forceps #23S are pediatric dental instruments designed for the extraction of lower primary molars. 23S forceps is pedo version of "Forceps No. 23" having a cowhorn beak design, suited for grasping lower molars due to its curved beaks. These forceps are specifically used in children's dentistry for first and second molar extractions. American pattern forceps have a smooth hinge so no pin is visible. All Our Extraction Forceps can be serrated unless otherwise indicated.

DENTAL EXTRACTION FORCEPS #23S:

- Pediatric Version of Extraction Forceps #23
- Designed Primary Lower Molars.
- Cowhorn Beaks Design.
- Single Spring Grip Handle.
- Autoclavable/ Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Pedo Extraction Forceps #23S – Spring Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Extracting Forceps
APPLICATION: Children Teeth



WI-36702823

Atraumatic Dressing Forceps Teeth Precision Tissue Handling Instrument Set

\$3.00 per piece

Trade price from \$1.30
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Atraumatic Dressing Forceps Teeth Precision Tissue Handling Instrument Set — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Supplied as a set — the contents can be adjusted to your list. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.00, falling to \$1.30 at 25 pieces — about 57% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34366784

Backhaus Stainless Steel Dental Towel Forceps

\$2.50 per piece

Trade price from \$1.50
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Backhaus Stainless Steel Dental Towel Forceps for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.50, falling to \$1.50 at 25 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34363801

Backhaus Towel Forceps Stainless Steel Dental

\$2.50 per piece

Trade price from \$1.50
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Backhaus Towel Forceps Stainless Steel Dental — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.50, falling to \$1.50 at 25 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35539093

Bone Reduction Forceps with Teeth Curved 5.5 Inch Orthopedic

\$40.00 per piece

Trade price from **\$25.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Bone Reduction Forceps with Teeth Curved 5.5 Inch Orthopedic — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: a curved working end for reach around structures, and interlocking teeth for a positive hold on tough tissue.

The 5.5" size shown is our standard; other sizes in the same pattern are made to order. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$40.00, falling to \$25.00 at 10 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-13330170

Children's Dental Extraction Forceps Curved Straight

\$5.99 per piece

Trade price from **\$2.50**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Children's Dental Extraction Forceps Curved Straight — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$5.99, falling to \$2.50 at 10 pieces — about 58% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-43392166

Dental Adson Brown Forceps Curved Thumb Dressing Tweezers

\$6.00 per piece Trade price from **\$4.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dental Adson Brown Forceps Curved Thumb Dressing Tweezers, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a curved working end for reach around structures.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$6.00, falling to \$4.00 at 1 piece — about 33% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-29087315

Dental Basic Examination Set Mouth Mirror Handle Double-Ended Hook

\$2.00 per piece Trade price from **\$1.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Basic Examination Set Mouth Mirror Handle Double-Ended Hook is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: two working ends on one handle, halving what you carry.

Supplied as a set — the contents can be adjusted to your list. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.00, falling to \$1.00 at 1 set — about 50% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.



WI-30787180

Dental College Tweezer Forceps–Stainless Steel Serrated Tip Multipurpose

\$0.80 per piece

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dental College Tweezer Forceps–Stainless Steel Serrated Tip Multipurpose, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: serrated jaws for grip that does not slip under tension.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$0.80. Minimum order 10 pieces, and the price moves with quantity.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36801641

Dental Cotton Holding Tweezers Stainless Steel Clinical Instrument

\$2.00 per piece

Trade price from \$1.40
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Cotton Holding Tweezers Stainless Steel Clinical Instrument is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.00, falling to \$1.40 at 25 pieces — about 30% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36803638

Dental Dressing Tweezers Precision Hand Instrument

\$2.00 per piece

Trade price from \$1.50
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Dressing Tweezers Precision Hand Instrument is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.00, falling to \$1.50 at 25 pieces — about 25% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-24285324

Dental Elevator Forceps Extraction 8-Piece Kit

\$18.00 per piece

Trade price from \$15.00
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Elevator Forceps Extraction 8-Piece Kit is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$18.00, falling to \$15.00 at 2 sets — about 17% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29020422

Dental Examination Kit 3

\$2.00 per piece

Trade price from \$1.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dental Examination Kit 3, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.00, falling to \$1.00 at 1 set — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29028418

Dental Examination Set Mirror Probe Tweezers Dentist College

\$2.00 per piece

Trade price from \$1.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Examination Set Mirror Probe Tweezers Dentist College for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.00, falling to \$1.00 at 1 set — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003316****Dental Extracting Forceps #101 for Upper & Lower Premolars****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Extracting Forceps #101 for Presidential Upper and Lower Premolars and Deciduous Premium Quality Stainless Steel American Pattern Teeth Extraction Surgical Dental Instruments.

101 Universal Extraction Forceps used for extracting teeth from the alveolar bone. Dental extracting forceps #101 features universal conical beaks designed to adapt to upper and lower premolars, deciduous teeth, right or left. American Pattern forceps have a smooth hinge so no pin is visible. All our extraction forceps can be serrated unless otherwise indicated.

UNIVERSAL EXTRACTING FORCEPS #101 FEATURES

- Conical Precise Tips Meet at a Point, Ensuring Optimal Grasp.

- For Upper & Lower Premolars and Deciduous Teeth.
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc American Pattern Universal Extraction Forceps #101

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Conical Beaks**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**WI-70004144****Dental Extracting Forceps #10H Upper Molars Universal****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #10H Universal Upper Molars with Hook Handle Premium Quality Stainless Steel American Pattern Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 10 with hook handle is used to extract upper molars.

American extraction forceps #10H features conical beaks designed to adapt to upper molars to improve grasp. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

AMERICAN PATTERN EXTRACTION FORCEPS #10H FEATURES

- Conical Beaks Design.
- Upper Molars, Universal
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 1 Pc American Pattern Extraction Forceps #10H with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Conical Serrated Beaks**POLISHED:** Yes**STERILE:** Yes

**WI-70004134****Dental Extracting Forceps #150AS Serrated Split Beaks****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

American Pattern Dental Extracting Forceps #150AS Split Beaks, Serrated German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Extraction Forceps No. 150AS used for extracting teeth from the alveolar bone. Upper universal split serrated beaks designed to provide increased grasp beneath the crown to prevent crown or root breakage. Serrations provide increased mechanical grasp on the root without increased pressure. This dental forceps can be used left and right.

American pattern forceps have a smooth hinge so no pin is visible. The stainless steel construction ensures maximum durability and dependability. You can reuse these 150AS extracting forceps over time as they will not bend or tarnish on you. High quality stainless


steel, excellent international quality standard.

AMERICAN EXTRACTING FORCEPS #150AS FEATURES

- For Upper Anteriors.
- Serrated, Split Beaks Provide an Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc American Pattern Dental Extracting Forceps #150AS Split Beaks – Serrated

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**WI-70004160****Dental Extracting Forceps #17 Lower Molars****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality American Pattern Dental Extracting Forceps #17 For Lower Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#17 xtraction Forceps is used to extract Lower Molars. The Working end of this tooth extraction instrument is designed with a sharp bend, making sure it reaches deep into the socket for a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.


EXTRACTING FORCEPS #17 FEATURES

- Tri-pointed Beak Design

- For Lower Molars
- Non-slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Dental Extracting Forceps #17 – Lower Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test,

**WI-70004118****Dental Extracting Forceps #18R for Upper 1st & 2nd Molars****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps #18R for Upper 1st & 2nd Molars, Right Pointed Tip Design with Hook Handle German Grade Stainless Steel Teeth Extraction Surgical Dental Instruments

American Extracting forceps No. 18R hook handle are the perfect tool to remove teeth from the alveolar bone.

Extraction forceps #18R has right pointed beak design features an American Pattern. Bifurcated pointed beaks designed for upper first and second molars. Hooked handle provides increased leverage when approach is from behind the head.

American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you

EXTRACTION FORCEPS #18R FEATURES

- Right Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE



- 1 Pc American Extraction Forceps #18R Right Pointed Beak with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction

Surgery Procedure

BEAK/JAWS: Right Pointed Tip Design**WI-70001370****Dental Extracting Forceps #1STD Gold Hollow Python Grip Handle****\$14.99** per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty Dental Extracting Forceps #1STD Upper Incisors and Cuspids, Standard, Serrated Beak with Gold Hollow Python Grip Handle German Grade Stainless Steel Dental Surgical Instruments

#1 Standard Extraction Forceps (V90-1STD) upper incisors and cuspids serrated forceps used for extracting teeth from the alveolar bones. German grade stainless steel promotes longevity and corrosion resistance. Straight, conical beaks designed to adapt upper anterior incisors and canines. Our tools are premium quality and used by professionals.

EXTRACTING FORCEPS 1STD FEATURES

- Vantage #1 Standard Extraction Forceps (V90-1STD)
- Upper Incisors and Cuspids Serrated Forceps
- Gold Hollow Python Grip Handle
- Non-Slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Extracting Forceps #1STD with Gold Hollow Python Grip Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes

**WI-70004142****Dental Extracting Forceps #210H Upper 3rd Molars****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #210H for Upper 3rd Molars with Hook Handle Premium Quality Stainless Steel American Pattern Bayonet Style Serrated Beaks Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 210H with hook handle is used to extract teeth from the alveolar bone. American extraction forceps #210H features offset (Bayonet) conical beaks designed to adapt to upper third molars with longitudinal serrations to improve grasp. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

AMERICAN PATTERN EXTRACTION FORCEPS #210H FEATURES

- Offset (Bayonet) Conical Beaks Design.
- For Upper 3rd Molars.
- Longitudinal Serrations Improve Grasp.
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 1 Pc American Pattern Extraction Forceps #210H with Hook Handle

TECHNICAL INFORMATION**WI-70004128****Dental Extracting Forceps #213 for Upper Central & Canines Roots****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #213 for Upper Central & Canines Roots with Hook Handle Premium Quality Stainless Steel American Pattern Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 213 with hook handle is used to extract teeth from the alveolar bone. American extraction forceps #213 features conical beaks designed to adapt to upper central & canines roots with to improve grasp. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

FEATURES

- Conical Beaks Design.
- Upper Central & Canines Roots
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 1 Pc American Pattern Extraction Forceps #213 with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Conical Beaks**POLISHED:** Yes

**WI-70004158****Dental Extracting Forceps #222 Lower 3rd Molars****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #222 For Lower 3rd Molars Premium Quality Stainless Steel American Pattern Apical Beaks Teeth Extraction Dental Surgical Instruments

American Extraction Forceps No. 222 are designed to be used for extracting lower wisdom teeth. Extracting forceps #222 designed with a sharp apical beaks make complete root contact, cross serration give nonslip surface, quickly enter sub-gingival. It ensures no fracture or tissue damage and easy extraction. These dental forceps are perfect to be a part of your dental kit. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #222 FEATURES

- Serrated Apical Beak.
- For Lower 3rd Molars.
- Non-slip Grip Handle
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE

- 1 Pc Dental Extracting Forceps #222 – Lower 3rd Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOC LAVED:** Yes**WI-70004159****Dental Extracting Forceps #32 Upper Canines, Premolars & Molar****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Extracting Forceps #32 for Upper Canines, Premolars & Molars with Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#32 Universal forceps used to extract teeth from Upper Cuspids, Bicuspid and Molars with sharp edged narrow beaks for breaking the alveolar process and trimming it after extraction. Forceps dimension is 7.5" (19cm) overall in length with textured grip handles. Durable and sturdy our extraction forceps are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #32 FEATURES

- Narrow beak, beveled and smoothed edges
- For Upper Canines, Premolars & Molars
- Non-slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Extracting Forceps #32 – Narrow, Bayonet Pattern Beak

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes

**WI-70004157****Dental Extracting Forceps #32A Molars Root – Slim Beak****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #32A for Upper Cuspids, Biscuspids, Molars and Roots with Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#32A Extraction forceps used to extract teeth from Upper Cuspids, Biscuspids, Molars and roots. Sharp edged Slim beaks for breaking the alveolar process and trimming it after extraction.

Forceps dimension is 7.5" (19cm) overall in length with textured grip handles. Durable and sturdy our extraction forceps are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #32A FEATURES

- Slim, bayonet pattern beaks with parallel lines for optimal grip.
- For Upper Cuspids, Biscuspids, Molars and Roots
- Non-slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Extracting Forceps #32A – Slim, Bayonet Pattern Beak

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**WI-70000962****Dental Extracting Forceps #349 for Upper Incisors & Root****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forcep #349 for Upper Incisors & Root Premium Quality Stainless Steel American Pattern Serrated Jaws with Grip Handle Teeth Extraction Forceps Surgical Dental Instruments.

Extracting Forceps No. 349 used to upper incisors & root. #349 forceps features American Pattern design with deep reaching serrated jaws. American Pattern forceps have a smooth hinge so no pin is visible. All our extraction forceps can be serrated unless otherwise indicated. They are made from premium-grade stainless steel ensures maximum durability and dependability. You can reuse these #349 forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #349 FEATURES

- Serrated Beaks for Deep & Parallel Grip.
- Upper Incisors & Root.
- Non-slip Textured-grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc American Pattern Extracting Forceps #349 – Serrated Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Tapered & Serrated Beaks

**WI-70004139****Dental Extracting Forceps #53L For Upper 1st & 2nd Left Molars****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality American Pattern Dental Extracting Forceps #53L for Upper 1st and 2nd Left Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

53L dental forceps are used to extract maxillary 1st and 2nd left molars. Left angled pointed offset (Bayonet) beaks for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

EXTRACTING FORCEPS #53L FEATURES

- Left Pointed Tip Design
- For Upper 1st & 2nd Molars
- Non-Slip Textured-Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc American Pattern Dental Extracting Forceps #53L

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**WI-70001236****Dental Extracting Forceps #53L with Gold Hollow Python Grip Handle****\$14.99** per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty American Pattern Extracting Forceps #53L for Upper 1st and 2nd Left Molars with Gold Python Grip Hollow Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Gold hollow python grip extraction forceps #53L are used to extract maxillary 1st and 2nd left molars. Left angled pointed offset (Bayonet) beaks for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

EXTRACTING FORCEPS #53L FEATURES

- Left Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Angular Beaks.
- Gold Python Grip Hollow Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Dental Extracting Forceps #53L with Gold Python Grip Hollow Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes

**WI-70003255****Dental Extracting Forceps #62 Upper and Lower Anterior Teeth****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal American Pattern Extracting Forceps No. 62 for Upper & Lower Incisors, Canines, Premolars and Deciduous Teeth German Grade Stainless Steel Tooth Extraction Dental Surgical Instrument.

American Pattern Tooth Extraction Forceps #62 used to extract Upper and Lower Incisors & Canines. The long, extra sharp and very narrow beaks allow a tight grasp of the tooth. The slightly curved working end makes this instrument capable of handling both upper and lower teeth, in addition to minimizing trauma to the patient.

American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All

our dental tools are professional quality and used by professionals.

AMERICAN EXTRACTING FORCEPS #62 FEATURES

- Precise Tips Meet at a Point, Ensuring Optimal Grip.
- Universal, Conical Beaks Design.
- For Upper & Lower Incisors, Canines, Premolars and Deciduous Teeth.
- Gently-Curved Working End Making It Suitable For Both Upper & Lower Teeth.
- Non-Slip Textured-Grip Handle.
- German Stainless Steel.

PACKAGE INCLUDE

- 1 Pc American Pattern Extracting Forceps #62

TECHNICAL INFORMATION**WI-70004156****Dental Extracting Forceps #65 for Upper Roots Dental Tools****\$11.99** per pieceTrade price from **\$11.03**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

American Pattern Dental Extraction Forceps #65 for Upper Roots, Fragments & Incisors, Slim Bayonet Pattern Jaws with Non-slip Handle German Grade Stainless Steel Teeth Extraction Dental Surgical Instrument.

American pattern dental extracting forceps No. 65 used to extract teeth from upper roots, fragments & incisors. Slim bayonet pattern beaks designed specially for hard to reach area. American pattern forceps have a smooth hinge so no pin is visible. The top quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.


- For Upper Roots, Fragments & Incisors.
- Slim, Bayonet Pattern Beaks.
- Non-Slip Textured-grip Handle.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc American Pattern Extracting Forceps #65

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Slim Bayonet Jaws**POLISHED:** Yes**HANDLE:** Non-slip Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**EXTRACTING FORCEPS #65 FEATURES**

**WI-70004125****Dental Extracting Forceps #88L 1st & 2nd Upper Molar****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps #88L with Left Pointed Split Beak Grip Handle for 1st & 2nd Upper Molars German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American pattern dental extraction forceps No. 88L used to extract 1st and 2nd molars. This extraction forceps has left pointed split beak design. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #88L FEATURES

- For Upper 1st & 2nd Molars.


- Left Pointed, Split Beak Design.
- Non-Slip Grip Handle.
- American Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc American Pattern Extracting Forceps #88L for 1st & 2nd Upper Molar

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Left Pointed Split Design**POLISHED:** Yes**HANDLE:** Non-Slip Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****WI-70004124****Dental Extracting Forceps #88R For Right 1st And 2nd Upper Molar****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #88R For Right 1st And 2nd Upper Molar American Pattern Right Pointed, Serrated Split Beak Design German Grade Stainless Steel Tooth Extraction Dental Surgical Instrument.

Extraction Forceps No. 88R design to remove teeth from bony socket. Suitable for maxillary right first and second Molars. Presidential forceps for 1st and Upper 2nd Molars. Bayonet design split beak for engaging palatal root. Each beak has different design to adapt to the maxillary molar roots that differ anatomically on the facial and lingual. Beak with one projection contours to bifurcation of mesial-buccal and distal-buccal roots. Beak with two projections contours to lingual root. American pattern forceps have a


smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

AMERICAN EXTRACTION FORCEPS #88R FEATURES

- For Upper 1st & 2nd Molars.
- Right Pointed, Serrated Split Beak Design
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc American Extraction Extracting Forceps #88R

TECHNICAL INFORMATION

**WI-70004117**

Dental Extracting Forceps #89 for Upper Right Molars

\$11.99 per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps #89 for Upper Molars with Right Pointed Tip and Textured Grip Handle Premium Quality Teeth Extraction Surgical Dental Instrument

Extracting Forceps No. 89 are used for extracting Upper Right Molars during dental procedures. Forceps have parallel grasp on the neck of the tooth and deep reach in the alveolar region which provides a secure hold for the extraction of the tooth. Durable and sturdy, our extraction forceps #89 are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #89 FEATURES

- Right Pointed Tip Design.
- For Upper Molars.
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc American Pattern Extracting Forceps #89 for Upper Molar

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Right Pointed Tip Design**POLISHED:** Yes**HANDLE:** Non-Slip Handle**WI-70004136**

Dental Extracting Forceps 18L for Upper 1st & 2nd Molars

\$11.99 per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps #18L for Upper 1st & 2nd Molars, Left Pointed Tip Design with Hook Handle German Grade Stainless Steel Teeth Extraction Surgical Dental Instruments

American Extracting forceps No. 18L hook handle are the perfect tool to remove teeth from the alveolar bone. Extraction forceps #18L has left pointed beak design features an American Pattern. Bifurcated pointed beaks designed for upper first and second molars. Hooked handle provides increased leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you

EXTRACTION FORCEPS #18L FEATURES

- Left Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc American Extraction Forceps #18L Left Pointed Beak with Hook Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction

Surgery Procedure

BEAK/JAWS: Left Pointed Tip Design



WI-36864648

Dental Extracting Forceps No 16 16cm

\$12.50 per piece

Trade price from \$8.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Extracting Forceps No 16 16cm — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 16cm size shown is our standard; other sizes in the same pattern are made to order.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$12.50, falling to \$8.00 at 25 pieces — about 36% off the single-unit price. Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004155

Dental Extracting Forceps No. 24 for Upper Molars

\$11.99 per piece

Trade price from \$11.03
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Extracting Forceps #24 for Upper Molars with Wide Rounded Beaks and Hook Handle Premium Quality American Pattern Tooth Extraction Forceps for both Human and Small Animals Dentistry Instrument.

Universal extracting forceps #24 are designed for the extraction of upper molars, including first, second, and sometimes third molars. Their beaks securely grip the teeth, allowing for controlled extraction while minimizing trauma to surrounding tissues. They feature an American pattern design with wide, rounded, serrated beaks and a hook handle for enhanced leverage, making them versatile for various upper molar extractions.

AMERICAN UNIVERSAL
EXTRACTING FORCEPS #24

FEATURES

- For Upper Molars.
- Wide, Rounded Serrated Jaws.
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Autoclavable/Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Universal Extraction Forceps #24 Hook Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Extracting Forceps
APPLICATION: Teeth Extraction Forceps





WI-70003637

Dental Extracting Forceps No.44 Lower Anterior

\$12.99 per piece

Trade price from **\$11.95**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pediatric Extracting Forceps #44
Narrow Tapered Beaks for Lower Anterior Premium Quality Stainless Steel American Pattern Child Teeth Extraction Surgical Dental Instruments.

#44 Lower anterior pediatric extracting forceps are durable tools designed for tooth removal from the alveolar bone in young patients. Their slim beaks and tips enhance maneuverability and field-of-view in small oral cavities while maintaining grip strength. These American pattern forceps feature a smooth hinge with no visible pin, optimizing precision and control in delicate dental procedures. All our extraction forceps can be serrated unless otherwise indicated.

EXTRACTING FORCEPS #44 FEATURES

- Narrow Tapered Beaks.
- For Lower Anterior.
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Pediatric Extraction Forceps #44 – Lower Anterior

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Dental

Extracting Forceps

APPLICATION: Pediatric Teeth Extraction Procedure

STERILE: Non-Sterile

AUTOCCLAVABLE/REUSABLE: Yes

RUST PREVENTION

PROCEDURE: Passivated



WI-42862583

Dental Extracting Forceps Set 5PCS Gold Handle Tooth Extraction

\$30.00 per piece

Trade price from **\$16.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Extracting Forceps Set 5PCS Gold Handle Tooth Extraction for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house.

Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Supplied as a 5-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$30.00, falling to \$16.00 at 1 set — about 47% off the single-unit price. Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42826808

Dental Extracting Forceps Set 8PCS Stainless Steel Tooth Extraction

\$45.00 per piece

Trade price from \$30.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Extracting Forceps Set 8PCS Stainless Steel Tooth Extraction is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a 8-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$45.00, falling to \$30.00 at 1 set — about 33% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42856850

Dental Extracting Forceps Stainless Steel Oral Surgery Instrument

\$10.00 per piece

Trade price from \$4.90
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Extracting Forceps Stainless Steel Oral Surgery Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$10.00, falling to \$4.90 at 1 piece — about 51% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-46820684

Dental Extracting Forceps Tooth Removal Instrument for Oral Surgery

\$8.00 per piece

Trade price from \$4.00
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Extracting Forceps Tooth Removal Instrument for Oral Surgery is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$8.00, falling to \$4.00 at 50 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



WI-70004138

Dental Extracting Kells Forceps #99C

\$11.99 per piece

Trade price from \$11.03
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Kells Forceps #99C for Upper Incisors, Canines & Premolars German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument

Kells extracting forceps No. 99 used for extracting teeth from the alveolar bone. The Short, straight conical beaks designed to adapt to upper incisors, canines and premolars. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

- For Upper Incisors, Canines & Premolars.
- Short, Straight Conical Beaks for Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Kells Extracting Forceps #99C – Serrated Beaks

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Extracting Forceps

APPLICATION: Extraction Forceps
BEAK/JAWS: Short Serrated Beaks

POLISHED: Yes
HANDLE: Non-Slip Handle
REUSABLE: Yes



WI-42891146

Dental Extracting Serrated Forceps

\$9.00 per piece

Trade price from \$4.30
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Extracting Serrated Forceps is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: serrated jaws for grip that does not slip under tension.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$9.00, falling to \$4.30 at 1 piece — about 52% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004151

Dental Extraction Forceps #150 Upper Incisors, Root Teeth

\$11.99 per piece

Trade price from \$11.03
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #150 Upper Incisors, Canines, and Premolars German Grade Stainless Steel Teeth Extraction Forceps Dental Surgical Instruments.

Universal Extraction Forceps No. 150 are used for extracting teeth from the alveolar bone. Universal conical beaks designed to adapt to upper incisors, canines, premolars, roots, universal right or left. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

- For Upper Incisors, Canines & Premolars.
- Serrated Conical Beaks for Optimal Grip.
- American Pattern.
- Non-ship Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Universal Extraction Forceps #150 – Serrated Beaks

TECHNICAL INFORMATION

MATERIAL: High-Quality Medical Grade Stainless Steel

INSTRUMENTS TYPE: Dental Extracting Forceps

APPLICATION: Extraction Forceps

BEAK/JAWS: Serrated

POLISHED: Yes

STERILE: Yes



UNIVERSAL EXTRACTION FORCEPS #150 FEATURES



WI-42834901

**Dental Extraction Forceps 10S Gold Handle
Stainless Steel Upper**

\$9.00 per piece

Trade price from **\$5.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dental Extraction Forceps 10S Gold Handle Stainless Steel Upper, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$9.00, falling to \$5.00 at 1 piece — about 44% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42861734

**Dental Extraction Forceps Stainless Steel
Serrated Jaw Tooth Extraction**

\$9.50 per piece

Trade price from **\$4.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Extraction Forceps Stainless Steel Serrated Jaw Tooth Extraction is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: serrated jaws for grip that does not slip under tension.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$9.50, falling to \$4.00 at 1 piece — about 58% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42847770

Dental Extraction Forceps Stainless Steel
Universal Serrated Jaw Instrument

\$7.00 per piece

Trade price from \$3.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dental Extraction Forceps Stainless Steel Universal Serrated Jaw Instrument, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: serrated jaws for grip that does not slip under tension.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$7.00, falling to \$3.50 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42887187

Dental Extraction Serrated Forceps Stainless
Steel Tooth Extraction Instrument

\$9.00 per piece

Trade price from \$4.30
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dental Extraction Serrated Forceps Stainless Steel Tooth Extraction Instrument, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: serrated jaws for grip that does not slip under tension.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$9.00, falling to \$4.30 at 1 piece — about 52% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36729139

Dental Forceps Holder Stainless Steel Dental Instrument Rack Sterilization

\$1.00 per piece

Trade price from **\$0.70**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Forceps Holder Stainless Steel Dental Instrument Rack Sterilization for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$1.00, falling to \$0.70 at 25 pieces — about 30% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36755108

Dental Forceps Storage Rack Stainless Steel Sterilization Holder

\$1.00 per piece

Trade price from **\$0.70**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dental Forceps Storage Rack Stainless Steel Sterilization Holder, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$1.00, falling to \$0.70 at 25 pieces — about 30% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42868317

Dental Hygiene Kit Stainless Steel Oral Care Tools Scaler

\$31.00 per piece

Trade price from **\$14.00**

Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Hygiene Kit Stainless Steel Oral Care Tools Scaler is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$31.00, falling to \$14.00 at 1 set — about 55% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36727154

Dental Instrument Holder Stainless Steel Forceps Storage

\$0.90 per piece

Trade price from **\$0.60**

Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Instrument Holder Stainless Steel Forceps Storage is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$0.90, falling to \$0.60 at 25 pieces — about 33% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29059343

Dental Kit Steel Tweezer Set

\$2.00 per piece

Trade price from \$1.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Kit Steel Tweezer Set is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.00, falling to \$1.00 at 1 set — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34949415

Dental Orthodontic Cheatle Forceps

\$3.00 per piece

Trade price from \$1.90
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Orthodontic Cheatle Forceps is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$3.00, falling to \$1.90 at 10 pieces — about 37% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27787337

Dental Surgical Tweezers Steel Blade Operation

\$1.00 per piece

Trade price from \$0.60
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

This is our pattern of Dental Surgical Tweezers Steel Blade Operation, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$1.00, falling to \$0.60 at 5 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25924708

Dental TC Tweezers 15cm/6" Stainless Steel Surgical Forceps

\$3.00 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental TC Tweezers 15cm/6" Stainless Steel Surgical Forceps — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: tungsten carbide inserts for a longer cutting life.

The 15cm size shown is our standard; other sizes in the same pattern are

made to order. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$3.00, falling to \$1.50 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70004152****Dental Tooth Extracting Forceps #10S with Serrated Jaws****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty American Pattern Tooth Extracting Forceps Fig. 10S Upper Molar Serrated Curved Jaws German Grade Stainless Steel Dental Instrument.

Durable and sturdy, Upper molars extracting forceps #10S are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals dentists.

DENTAL EXTRACTING FORCEPS 10S FEATURES

- Serrated Beak For Better Grip
- Figure 10S For Upper Molar
- Smoothly Curved Tips Making Them Suitable For Both Sides
- Non-Slip Grip Handle


- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Extracting Forceps Fig. 10S Curved Serrated – Upper Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**WI-30805248****Dental Tooth Extraction Forceps Set 7 Gold Plated****\$65.00** per piece

Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Tooth Extraction Forceps Set 7 Gold Plated for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house.

Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$65.00. Minimum order 10 sets, and the price moves with quantity.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-45076284

Dental Tooth Forceps 10-Piece Set Stainless Steel Silver

\$42.00 per piece Trade price from \$32.00 Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Tooth Forceps 10-Piece Set Stainless Steel Silver is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$42.00, falling to \$32.00 at 1 piece — about 24% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27812198

Dental Tweezers Surgical Cotton Tweezers

\$1.00 per piece Trade price from \$0.60 Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Tweezers Surgical Cotton Tweezers — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by

hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$1.00, falling to \$0.60 at 5 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42861952

Dental Universal Extraction Forceps Stainless Steel Tooth Extraction Instrument

\$9.00 per piece

Trade price from \$4.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Universal Extraction Forceps Stainless Steel Tooth Extraction Instrument for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$9.00, falling to \$4.00 at 1 piece — about 56% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-40981380

Dissecting Tissue Thumb Forceps Lab Multi-Purpose Tweezers Blunt 1x2

\$7.00 per piece

Trade price from \$4.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dissecting Tissue Thumb Forceps Lab Multi-Purpose Tweezers Blunt 1x2 is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: blunt tips for atraumatic handling.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$7.00, falling to \$4.00 at 25 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34250028

English Pattern Dental Tooth Extracting Forceps

\$8.00 per piece

Trade price from \$5.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

English Pattern Dental Tooth Extracting Forceps is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32723692

English Pattern Dental Tooth Extracting Forceps German

\$8.00 per piece

Trade price from \$5.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

English Pattern Dental Tooth Extracting Forceps German — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001324****Extracting Forceps #1 with Gold Hollow Python Grip Handle****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty Gold Python Grip
Extracting Forceps #1 Hollow Handle
American Pattern Tooth Extraction
Forceps for Upper Incisors as &
Canines German Grade Stainless Steel
Dental Instrument

#1 extraction forcep used in the
extraction of upper anterior teeth i.e
Central & Lateral Incisor. German grade
stainless steel promotes longevity and
corrosion resistance. American pattern
forceps have a smooth hinge so no pin
is visible. Our tools are professional
quality and used by professionals
dentists.

**EXTRACTING FORCEPS #1
FEATURES**

- Serrated Jaws for Firm Grip
- Figure 1 For Canines/Incisors.

- Gold Hollow Python Grip Handle
- Non-Slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Extracting Forceps #1 with Gold Hollow Python Grip Handle – Upper Incisors & Canines

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**WI-70004150****Extracting Forceps #103 for Lower Bicuspid, Incisors & Roots****\$11.99** per pieceTrade price from **\$11.03**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Extracting Forceps #103 for
Lower Bicuspid, Incisors & Roots with
Hook Handle Premium Quality Stainless
Steel American Pattern Teeth Extraction
Surgical Dental Instruments

Extracting Forceps No.103 with hook
handle is used to extract teeth from the
alveolar bone. Dental extracting
forceps #103 features universal
serrated precise beaks designed to
adapt to upper & lower bicuspid,
incisors & roots. American pattern
forceps have a smooth hinge so no pin
is visible. You can reuse these
extracting forceps over time as they will
not bend or tarnish on you.

**UNIVERSAL EXTRACTION
FORCEPS #103 FEATURES**

- Precise Serrated Tips Meet at a Point, Ensuring Optimal Grasp.
- For Lower Bicuspid Incisor and Roots.
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 1 Pc American Pattern Universal Extraction Forceps #103

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction

Forceps

BEAK/JAWS: Serrated Beaks

**WI-70004143****Extracting Forceps #150A for Upper Incisors, Cuspids and Bicuspids****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Extracting Forceps #150A Parallel Beaks with Plain Edge for Presidential Upper Incisors, Cuspids and Bicuspids Premium Quality Stainless Steel American Pattern Teeth Extraction Dental Surgical Instruments.

#150A maxillary forceps is a universal, parallel conical beaks designed to adapt to upper incisors, canines and premolars used for extracting teeth from the alveolar bone. 150A parallel beaks allow for increased contact beneath the crown for right or left extractions. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction and easy gripping handles ensures maximum durability and dependability. You can reuse t150a


extracting forceps over time as they will not bend or tarnish on you.

UNIVERSAL EXTRACTING FORCEPS #150A FEATURES

- For Upper Incisors, Cuspids and Bicuspids
- Parallel Beaks with Plain Edge for Optimal Grip.
- American Pattern.
- Non-ship Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc American Extraction Forceps #151A – Parallel Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**WI-70000411****Extracting Forceps #150AS Diamond Dusted Python Grip Gold Hollow Handle****\$14.99** per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty Diamond Dusted Extracting Forceps #150AS for Upper Anteriors American Pattern Split Serrated Beaks with Python Grip Gold Hollow Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Extraction Forceps No. 150AS used for extracting teeth from the alveolar bone. Diamond dusted upper universal split serrated beaks designed to provide increased grasp beneath the crown to prevent crown or root breakage. Serrations provide increased mechanical grasp on the root without increased pressure. #150AS dental forceps with gold color hollow python grip handles can be used left and right.

American pattern forceps have a smooth hinge so no pin is visible. The stainless steel construction ensures

maximum durability and dependability. You can reuse these 150AS extracting forceps over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

AMERICAN EXTRACTING FORCEPS #150AS FEATURES

- For Upper Anteriors.
- Diamond Dusted Tips.
- Split Serrated Beaks Provide an Optimal Grip.
- American Pattern.
- Non-ship Gold Python Grip Hollow Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Extracting Forceps #150AS with Python Grip Gold Hollow Handle



**WI-70000408****Extracting Forceps #151 P-1 Diamond Dusted Gold Hollow Handle****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Extraction Forceps #151 Diamond Dusted Serrated Beaks (Pattern #1) with Gold Hollow Python Grip Handle for Lower Incisors, Canines, Premolars & Roots Premium Quality Stainless Steel American Pattern Tooth Extracting Surgical Dental Instrument.

#151 extracting forceps, also known as lower universal forceps, are dental tools used for removing teeth in the lower jaw. #151 dental forceps feature thin, tapered, and sharp parallel beaks for good contact with tooth root and crown, Diamond dusted longitudinal serrated beaks feature an improved grip with gold-colored hollow python grip handles for durability and corrosion resistance. These forceps are meticulously crafted and finished with a

scissors hinge for easier cleaning and force distribution.

DIAMOND DUSTED #151 EXTRACTION FORCEPS FEATURES

- For Lower Incisors, Canines, Premolars & Roots.
- Diamond Dusted Tips.
- Thin, Tapered, and Sharp Parallel Beaks for Good Contact.
- Longitudinal Serrations for Improved Grip and Atraumatic Extraction.
- American Pattern.
- Gold Hollow Python Grip Handle.
- Made With Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Extracting Forceps #151 Pattern-1 Diamond Dusted

**WI-70000407****Extracting Forceps #151 Pattern-2 Diamond Dusted****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #151 Pattern-2 Diamond Dusted Serrated Grip Beaks with Gold Hollow Python Grip Handle for Lower Incisors, Canines, Premolars & Roots Premium Quality Stainless Steel American Pattern Tooth Extracting Surgical Dental Instrument.

#151 extracting forceps, also known as lower universal forceps, are dental tools used for removing teeth in the lower jaw. #151 dental forceps feature thin, tapered, and sharp parallel grip beaks for good contact with tooth root and crown, Diamond dusted longitudinal serrated with central cut beaks feature an improved grip with gold-colored hollow python grip handles for durability and corrosion resistance. These forceps are meticulously crafted

and finished with a scissors hinge for easier cleaning and force distribution.

DIAMOND DUSTED #151 P-2 EXTRACTION FORCEPS FEATURES

- For Lower Incisors, Canines, Premolars & Roots.
- Pattern-2 Diamond Dusted Tips.
- Thin, Tapered, and Sharp Parallel Beaks for Good Contact.
- Longitudinal Serrations for Improved Grip and Atraumatic Extraction.
- American Pattern.
- Gold Hollow Python Grip Handle.
- Made With Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Extracting Forceps #151 Pattern-2 Diamond Dusted



**WI-70004130****Extracting Forceps #151AS Serrated Split Beaks****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #151AS Split Beaks, Serrated for Lower Anterior – Standard – American Pattern German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Extraction Forceps No. 151AS used for extracting teeth from the alveolar bone. Lower universal split serrated beaks designed to provide increased grasp beneath the crown to prevent crown or root breakage. Serrations provide increased mechanical grasp on the root without increased pressure. This 151AS dental forceps can be used left and right.

American pattern forceps have a smooth hinge so no pin is visible. The stainless steel construction ensures maximum durability and dependability. You can reuse these 151AS extracting

forceps over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

AMERICAN EXTRACTING FORCEPS #151AS FEATURES

- For Lower Anterior.
- Serrated, Split Beaks Provide an Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Extracting Forceps #151AS Split Beaks Serrated -Lower Anterior

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental**WI-70004123****Extracting Forceps #203 Lower Incisors, Canines, Premolars & Roots****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps #203 for Lower Incisors, Canines, Premolars, and Roots Premium Quality Stainless Steel American Pattern Universal Wider Conical Beaks Teeth Extraction Surgical Dental Instruments.

American extracting forceps No. 203 used for extracting teeth from the alveolar bone. Extraction forceps #203 features a universal wider conical beaks designed to adapt to lower incisors, canines, premolars and roots, may be used left or right. American Pattern forceps have a smooth hinge so no pin is visible. All our extraction forceps can be serrated unless otherwise indicated.

AMERICAN EXTRACTING FORCEPS #203 FEATURES

- Wider Beaks for Better Grip.
- For Lower Incisors, Cuspids, Bicuspsids & Root.
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc American Pattern Extraction Forceps #203 Wider Conical Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Conical Beaks**POLISHED:** Yes**STERILE:** Yes

**WI-70004141****Extracting Forceps #210S Upper 3rd Molars, Universal****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #210S for Upper 3rd Molars Premium Quality Stainless Steel American Pattern Bayonet Style Serrated Beaks Teeth Extraction Surgical Dental Instruments.

American No.210S are used for extracting upper wisdom teeth and root fragments from the alveolar bone. Extraction forceps #210S features offset (Bayonet) conical beaks designed to adapt to upper third molars with longitudinal serrations to improve grasp. The unique design help in removing a tooth with a minor injury. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

NOTE: 210S is the same design as 210H without hooked handle.

AMERICAN EXTRACTION FORCEPS #210S FEATURES

- Offset (Bayonet) Conical Beaks Design.
- For Upper 3rd Molars.
- Longitudinal Serrations Improve Grasp.
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 1 Pc American Pattern Extraction Forceps #210S for Upper Molars

TECHNICAL INFORMATION**WI-70004133****Extracting Forceps #217 for 1st & 2nd Lower Molars****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

American Pattern Dental Extracting Forceps #217 For Lower Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#217 Extraction Forceps is used to extract lower molars. The working end of this tooth extraction instrument is designed with a tri-pointed shaped beaks, making sure it reaches deep into the socket for a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #217 FEATURES

- Tri-pointed Serrated Beak Design

- For 1st & 2nd Lower Molars
- Non-slip Grip Handle
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE

- 1 Pc Dental Extracting Forceps #217 – Lower Molars

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Dental Extracting Forceps

APPLICATION: Teeth Extraction

POLISHED: Yes

STERILE: Yes

REUSABLE: Yes

AUTOCCLAVED: Yes

RUST PREVENTION

PROCEDURE: Passivated

ULTRASONIC CLEANED: Yes



**WI-70004148****Extracting Forceps #23 Cow Horn Lower Molar****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Extracting Forceps #23 Cow Horn 1st & 2nd Lower Molars with Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#23 extraction forceps is use to extract teeth from 1st and 2nd lower molars with cow horn jaws. Forceps dimension is 7.5" (19cm) overall in length with Textured grip handles. Durable and sturdy our extraction forceps are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #23 FEATURES

- Slightly Larger Cowhorn Beak Design Than #16
- For Lower 1st & 2nd Molars
- Non-slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Extracting Forceps #23 – Cow Horn Jaws

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**WI-70001251****Extracting Forceps #23 Lower Molars with Green Color Coated Handle****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Extracting Forceps #23 Cow Horn 1st & 2nd Lower Molars with Green Color Coated Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#23 extraction forceps is use to extract teeth from 1st and 2nd lower molars with cow horn jaws. Forceps dimension is 7.5" (19cm) overall in length. Textured grip with green color coated handles. Durable and sturdy our extraction forceps are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #23 FEATURES

- Slightly Larger Cowhorn Beak Design Than #16
- For Lower 1st & 2nd Molars
- Green Color Coated Non-slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Extracting Forceps #23 – Cow Horn Jaws – Green Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCLAVED:** Yes**RUST PREVENTION**

**WI-70004140****Extracting Forceps #286 for Upper Roots, Incisors & Premolars****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps #286 for Upper Anterior Roots, Incisors & Bicuspsids Premium Quality Stainless Steel American Pattern Bayonet Beaks Teeth Extraction Surgical Dental Instruments.

American extracting forceps No. 286 used for extracting teeth from the alveolar bone. Extraction forceps #286 features offset (Bayonet) long, narrow beaks designed to adapt to upper roots, incisors and premolars. American Pattern forceps have a smooth hinge so no pin is visible. All our extraction forceps can be serrated unless otherwise indicated.

AMERICAN EXTRACTING FORCEPS #286 FEATURES

- Bayonet Pattern Long Narrow Tips

- For Upper Anterior Roots, Incisors & Bicuspsids
- Non-Slip Textured-Grip Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc American Pattern Extraction Forceps #286 Bayonet Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Forceps**BEAK/JAWS:** Bayonet Beaks**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**WI-70004147****Extracting Forceps #53R for Upper 1st & 2nd Right Molars****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality American Pattern Dental Extracting Forceps #53R for Upper 1st and 2nd Right Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

53R dental forceps are used to extract maxillary 1st and 2nd right molars. Right angled pointed offset (Bayonet) beaks for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

EXTRACTING FORCEPS #53R FEATURES

- Right Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Angular Beaks.
- Non-Slip Textured-Grip Handle.
- German Stainless Steel.

PACKAGE INCLUDE

- 1 Pc American Pattern Dental Extracting Forceps #53R

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated

**WI-70004146****Extracting Forceps #6 Upper & Lower Teeth Splitting & Separating****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

American Extraction Forceps #6 for Upper and Lower Teeth Splitting & Separating Premium Quality Stainless Steel Teeth Extraction Surgical Dental Instruments.

Extracting Forceps #6 are precision atraumatic extraction forceps that provide optimum tactical feel while reducing hand fatigue. The tapered tips adapt to the tooth with excellent subgingival contact. These dental forceps are perfect to be a part of your dental kit. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #6 FEATURES

- Blade-like Beaks for Splitting & Separating.
- For Upper & Lower Teeth.
- 6.5" (16.5cm) Overall Length.
- Non-Slip Textured-grip Handle.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc American Extracting Forceps #6 for Upper and Lower Teeth Splitting

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**WI-70004135****Extracting Forceps #69 Upper Lower Roots and Fragments****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #69 with Long Conical Beaks for Upper and Lower Roots and Fragments – American Pattern – German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

American pattern extracting forceps 69 are designed to remove upper and lower dental fragments and roots, particularly in the posterior regions. They feature long precise beaks, a smooth hinge without a visible pin, and are made of stainless steel for durability. These forceps are fully autoclavable and reusable.

AMERICAN EXTRACTING FORCEPS #69 FEATURES

- For Upper and Lower Roots and Fragments.

- Long, Precise Conical Beaks
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Dental Extracting Forceps #69

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Tooth Extraction Forceps**STERILE:** Non-Sterile**REUSABLE/AUTOCCLAVABLE:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test

**WI-70000981****Extracting Forceps #88R 1st & 2nd Upper Molar with Gold Python Grip****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps #88R with Right Pointed Split Beak Gold Python Grip Handle for 1st & 2nd Upper Molars German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American pattern dental extraction forceps No. 88R used to extract 1st and 2nd molars. This extraction forceps has right pointed split beak design. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #88R FEATURES

- For Upper 1st & 2nd Molars.
- Right Pointed, Split Beak Design.
- Gold Hollow Python Grip Handle.
- American Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc American Pattern Extracting Forceps #88R – Gold Python Grip Hollow Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Right Pointed Split Design**POLISHED:** Yes**HANDLE:** Gold Python Grip Hollow Handle**REUSABLE:** Yes**WI-70004145****Extracting Forceps #99A Upper Incisors, Canines & Premolar Root****\$11.99** per pieceTrade price from **\$11.03**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps #99A Straight Serrated Beaks for Upper Incisors, Canines, Premolar Root with Hook Handle Premium Quality American Pattern Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 99A used to extract upper incisor, cuspid and bicuspid from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

EXTRACTING FORCEPS #99A FEATURES

- Straight Serrated Beaks.
- For Upper Incisors, Canines, Premolar Root.

- Hook Handle.
- American Pattern.
- Made Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc American Pattern Extracting Forceps #99A with Hook Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps**BEAK/JAWS:** Straight Serrated Beaks**POLISHED:** Yes**HANDLE:** Non-Slip Hook Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated

**WI-70004131****Extracting Forceps 1 Hen Dental Surgical****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Henahan Dental Forceps #1HEN, Upper Incisors Teeth Extraction German Grade Stainless Steel Surgical Dental Instrument

Extraction forceps are the most familiar instruments with regard to the extraction of tooth teeth. Intra-alveolar extractions are mainly accomplished using elevators and extraction forceps. They are the instruments that are used to remove the tooth out of the alveolar bone socket.

Doctor use this tool for upper incisors, Henahan. Our tools are professional quality and used by professionals. All Dental Forceps are handcrafted from german grade stainless steel for logitivity and durability

EXTRACTING FORCEPS 1HEN FEATURES

- Upper Incisors and Canines Forceps
- Non-Slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Extracting Forceps #1HEN

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Extraction Forceps**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**WI-70004127****Extracting Forceps 3FS Dental Surgical Instrument****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #3FS Woodward's American Pattern Tooth Extraction Forceps Fig. 3FS For Lower Molars German Grade Stainless Steel Dental Surgical Instruments

Extracting forceps #3FS used for lower molars extracting teeth from the alveolar bones. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS 3FS FEATURES

- Pointed beak design
- For lower molars
- Non-Slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Woodward Dental Extracting Forceps 3FS – Lower Molars

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing



WI-28825128

Extracting Forceps Set Dental Instruments for Tooth Surgery

\$4.00 per piece

Minimum order 4 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Extracting Forceps Set Dental Instruments for Tooth Surgery — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$4.00. Minimum order 4 pieces, and the price moves with quantity.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000576

Extraction Forceps #150 Upper Incisors Gold Hollow Python Grip Handle

\$14.99 per piece

Trade price from \$13.79
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps #150 Upper Incisors, Canines, and Premolars with Gold Hollow Python Grip Handles German Grade Stainless Steel Universal Teeth Extraction Forceps Dental Surgical Instruments.

Universal Extraction Forceps No. 150 are used for extracting teeth from the alveolar bone. Universal conical beaks designed to adapt to upper incisors, canines, premolars, roots, universal right or left. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction and easy gripping gold color hollow python grip handles ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

UNIVERSAL EXTRACTION FORCEPS #150 FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- For Upper Incisors, Canines & Premolars.
- Serrated Conical Beaks for Optimal Grip.
- American Pattern.
- Non-ship Gold Hollow Python Grip Handles.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Dental Extraction Forceps #150 Serrated Beaks with Gold Hollow Python Grip Handles

TECHNICAL INFORMATION





WI-31344355

Extraction Forceps No 7 Upper Roots English Pattern

\$8.00 per piece

Trade price from **\$5.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Extraction Forceps No 7 Upper Roots English Pattern — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$8.00, falling to \$5.00 at 10 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36763643

Fine Tip Dressing Forceps Teeth Reusable Stainless Steel Surgical

\$3.00 per piece

Trade price from **\$1.20**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Fine Tip Dressing Forceps Teeth Reusable Stainless Steel Surgical is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: interlocking teeth for a positive hold on tough tissue, a fine pattern for delicate tissue, and full reusability — clean, autoclave and use again.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.00, falling to \$1.20 at 25 pieces — about 60% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-16643035

Genuine German Steel Machinery 17-Piece Oral Dental

\$45.00 per piece

Trade price from \$25.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Genuine German Steel Machinery 17-Piece Oral Dental — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$45.00, falling to \$25.00 at 1 set — about 44% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35080915

Gillies Toothed Dissecting Forceps Stainless Steel Dental Surgical Instrument

\$1.20 per piece

Trade price from \$0.70
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Gillies Toothed Dissecting Forceps Stainless Steel Dental Surgical Instrument is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.20, falling to \$0.70 at 10 pieces — about 42% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-31352471

HOT SALES Stainless Steel Dental Extracting Forceps Set

\$8.00 per piece

Trade price from **\$5.00**

Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of HOT SALES Stainless Steel Dental Extracting Forceps Set, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$8.00, falling to \$5.00 at 2 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42867462

MD1 Dental Extracting Forceps Stainless Steel Upper Tooth Extraction

\$8.00 per piece

Trade price from **\$4.00**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

MD1 Dental Extracting Forceps Stainless Steel Upper Tooth Extraction is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by

hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46955210

Michel Clip Extracting Forceps Surgical Kit for Precise Skin

\$8.90 per piece

Trade price from \$4.90
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Michel Clip Extracting Forceps Surgical Kit for Precise Skin, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$8.90, falling to \$4.90 at 25 pieces — about 45% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36804704

Micro Precision 16cm Non-Toothed Dressing Forceps Reusable Surgical

\$0.90 per piece

Trade price from \$0.45
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Micro Precision 16cm Non-Toothed Dressing Forceps Reusable Surgical is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: a micro-scale working end for fine and delicate work, interlocking teeth for a positive hold on tough tissue, and full reusability — clean, autoclave and use again.

The 16cm size shown is our standard; other sizes in the same pattern are made to order. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$0.90, falling to \$0.45 at 100 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36260860

Non Toothed Dissecting Forceps Medical

\$1.20 per piece

Trade price from \$0.60
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Non Toothed Dissecting Forceps Medical, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.20, falling to \$0.60 at 5 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36299043

Non Toothed Dissecting Forceps Stainless Steel Surgical Dressing Tissue

\$1.20 per piece

Trade price from \$0.60
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Non Toothed Dissecting Forceps Stainless Steel Surgical Dressing Tissue for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.20, falling to \$0.60 at 5 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29050394

Orthodontic Dental Surgical Tweezer Set 14cm Medical

\$2.00 per piece

Trade price from \$1.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Orthodontic Dental Surgical Tweezer Set 14cm Medical for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

The 14cm size shown is our standard; other sizes in the same pattern are

made to order. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.00, falling to \$1.00 at 1 set — about 50% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004121

Pedo Extracting Forceps #151S Dental Surgical Instrument

\$11.99 per piece

Trade price from \$11.03
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

American Pattern Pedo Extraction Forceps #151S for Children Lower Primary Teeth and Roots Stainless Steel Tooth Extraction Dental Surgical Instruments

Universal Pedo Forceps No. 151S used to extract lower primary teeth of children. The small, closely-situated beaks ensures maximum accessibility and a tight grasp of the tooth, while minimizing chances of trauma to the patient.

151S Pedo forceps guarantees the easiest tooth extraction procedure. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these

extracting forceps over time as they will not bend or tarnish on you.

AMERICAN EXTRACTION FORCEPS #151S FEATURES

- For Lower Primary Teeth and Roots
- Small, Conical Beaks for Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc American Extraction Forceps #151S

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Extracting Forceps
APPLICATION: Pedo Extraction





WI-35032080

Potts Smith Tissue Forceps 1x2 Teeth

\$3.00 per piece

Trade price from \$0.90

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Potts Smith Tissue Forceps 1x2 Teeth for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.00, falling to \$0.90 at 10 pieces — about 70% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29064365

Precision Dental Exam Tool Kit Stainless Steel Dental Probe

\$2.00 per piece

Trade price from \$1.00

Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Precision Dental Exam Tool Kit Stainless Steel Dental Probe — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.00, falling to \$1.00 at 1 set — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36827612

Reusable Dental College Tweezers Precision Cotton Holding Instrument

\$2.00 per piece

Trade price from \$1.50
Minimum order 30 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Reusable Dental College Tweezers Precision Cotton Holding Instrument is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.00, falling to \$1.50 at 30 pieces — about 25% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42844798

Reusable Dental Tooth Extraction Forceps Stainless Steel 3PCS Kit

\$19.00 per piece

Trade price from \$9.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Reusable Dental Tooth Extraction Forceps Stainless Steel 3PCS Kit — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$19.00, falling to \$9.00 at 1 set — about 53% off the single-unit price. Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Supplied as a set — the contents can be adjusted to your list. The steel is

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27423479

Solid Color Stainless Steel Brown

\$2.50 per piece

Trade price from \$1.90
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Solid Color Stainless Steel Brown — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.50, falling to \$1.90 at 1 piece — about 24% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36811740

Stainless Steel 16cm Dental Extracting Forceps No

\$12.50 per piece

Trade price from \$8.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Stainless Steel 16cm Dental Extracting Forceps No — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 16cm size shown is our standard; other sizes in the same pattern are made to order.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$12.50, falling to \$8.00 at 25 pieces — about 36% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-30455551

Stainless Steel Adson Surgical Forceps (Toothed)

\$1.00 per piece

Trade price from \$0.40

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Adson Surgical Forceps (Toothed, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.00, falling to \$0.40 at 1 piece — about 60% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-31374392

Stainless Steel Adult Tooth Extraction Dental

\$8.00 per piece

Trade price from \$5.00

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Adult Tooth Extraction Dental for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$8.00, falling to \$5.00 at 10 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34891604

Stainless Steel Dental Cotton Dressing Tweezer Surgical Tissue

\$2.00 per piece

Trade price from **\$0.90**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Dental Cotton Dressing Tweezer Surgical Tissue, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.00, falling to \$0.90 at 10 pieces — about 55% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-31406595

Stainless Steel Dental Extracting Forceps

\$8.00 per piece

Trade price from **\$5.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


This is our pattern of Stainless Steel Dental Extracting Forceps, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$8.00, falling to \$5.00 at 10 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-31379371

Stainless Steel Dental Extracting Forceps-Professional

\$8.00 per piece

Trade price from \$5.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Extracting Forceps-Professional is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$8.00, falling to \$5.00 at 10 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-32764672

Stainless Steel Dental Extraction Forceps Ergonomic Grip

\$8.00 per piece

Trade price from \$5.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Extraction Forceps Ergonomic Grip is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: an ergonomic handle for long lists without fatigue.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



WI-46801758

Stainless Steel Dental Extraction Forceps Oral Surgery Instrument

\$9.00 per piece

Trade price from **\$4.00**
Minimum order 12 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Extraction Forceps Oral Surgery Instrument is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$9.00, falling to \$4.00 at 12 pieces — about 56% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32770678

Stainless Steel Dental Forceps Tooth Extracting

\$8.00 per piece

Trade price from **\$5.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Forceps Tooth Extracting is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

We work the instrument by hand at the finishing stage — box locks, hinges and

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29044372

Stainless Steel Dental Hygiene Tool Kit Removal

\$2.00 per piece

Trade price from \$1.00

Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Hygiene Tool Kit Removal — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.00, falling to \$1.00 at 1 set — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-26456933

Stainless Steel Dental Surgical Tweezers & Tooth Forceps

\$2.50 per piece

Trade price from \$2.00

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Surgical Tweezers & Tooth Forceps is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.50, falling to \$2.00 at 10 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-25875769

Stainless Steel Dental TC Surgical Tweezers Ergonomic

\$3.00 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Dental TC Surgical Tweezers Ergonomic for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: tungsten carbide inserts for a longer cutting life, and an ergonomic handle for long lists without fatigue.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$3.00, falling to \$1.50 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-32730713

Stainless Steel Dental Tooth Extracting Forceps

\$8.00 per piece

Trade price from \$5.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Tooth Extracting Forceps — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



WI-31373307

Stainless Steel Dental Tooth Upper Roots Extraction

\$8.00 per piece

Trade price from \$5.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Tooth Upper Roots Extraction is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$5.00 at 10 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34951329

Stainless Steel Dressing Forceps Straight Surgical Thumb Tweezer

\$2.00 per piece

Trade price from \$0.90
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dressing Forceps Straight Surgical Thumb Tweezer is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a straight working end for direct, in-line access.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.00, falling to \$0.90 at 10 pieces — about 55% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-35092446

**Stainless Steel Gillies Dissecting Forceps
Toothed Dental Surgical Instrument**

\$1.20 per piece

Trade price from **\$0.70**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Gillies Dissecting Forceps Toothed Dental Surgical Instrument for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.20, falling to \$0.70 at 10 pieces — about 42% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-26029780

**Stainless Steel Heaney Hysterectomy Forceps
21cm**

\$8.00 per piece

Trade price from **\$4.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Heaney Hysterectomy Forceps 21cm, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 21cm size shown is our standard; other sizes in the same pattern are made to order.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$8.00, falling to \$4.00 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-32732705

Stainless Steel Pediatric Dental Forceps Ergonomic Grip

\$8.00 per piece

Trade price from \$5.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Pediatric Dental Forceps Ergonomic Grip for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: an ergonomic handle for long lists without fatigue.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-40970984

Stainless Steel Russian Toothed Dissecting Forceps

\$7.00 per piece

Trade price from \$4.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Russian Toothed Dissecting Forceps is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$7.00, falling to \$4.00 at 25 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-25799728

Stainless Steel Serrated Curved (15mm) Dental Root

\$4.50 per piece

Trade price from \$1.99
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Serrated Curved (15mm) Dental Root is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: a curved working end for reach around structures, and serrated jaws for grip that does not slip under tension.

The 15mm size shown is our standard; other sizes in the same pattern are

made to order. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$4.50, falling to \$1.99 at 10 pieces — about 56% off the single-unit price. Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-26535751

Stainless Steel Surgical Tissue Forceps/Fine Teeth Grasping Claw

\$2.50 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Tissue Forceps/Fine Teeth Grasping Claw is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: interlocking teeth for a positive hold on tough tissue, and a fine pattern for delicate tissue.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.50, falling to \$2.00 at 10 pieces — about 20% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-36832566

Stainless Steel Toothless Dissecting Forceps
Precision Tissue Handling Instrument

\$0.70 per piece

Trade price from \$0.60
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Toothless Dissecting Forceps Precision Tissue Handling Instrument for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance,

then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$0.70, falling to \$0.60 at 50 pieces — about 14% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-27811200

Stainless-Steel Dental Tweezers Medical

\$1.00 per piece

Trade price from \$0.60
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless-Steel Dental Tweezers Medical is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$1.00, falling to \$0.60 at 5 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Every piece is hand-finished and inspected before it leaves the bench:



WI-29077349

Steel Dental Examination Set

\$2.00 per piece

Trade price from \$1.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Steel Dental Examination Set is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.00, falling to \$1.00 at 1 set — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-22435893

Surgical Dental Miller Articulating Straight Steel Paper Forceps

\$3.00 per piece

Trade price from \$1.50
Minimum order 6 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Surgical Dental Miller Articulating Straight Steel Paper Forceps for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a straight working end for direct, in-line access.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.00, falling to \$1.50 at 6 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36917027

Surgical Tissue Forceps 17.5cm Stainless Steel 1x2 Teeth

\$8.50 per piece

Trade price from **\$4.30**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Tissue Forceps 17.5cm Stainless Steel 1x2 Teeth is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

The 17.5cm size shown is our standard; other sizes in the same pattern are

made to order. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$8.50, falling to \$4.30 at 10 pieces — about 49% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001254

T/C Dental Extracting Forceps #17 Lower Molars

\$14.99 per piece

Trade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality American Pattern T/C Universal Dental Extracting Forceps #17 For Lower 1st & 2nd Molars with Gold Plated Textured Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Tungsten Carbide Universal Extraction Forceps #17 is used to extract lower 1st & 2nd molars. The working end of this tooth extraction instrument is designed with a sharp bend, making sure it reaches deep into the socket for a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

T/C EXTRACTING FORCEPS #17 FEATURES

- Tri-pointed Beak Design
- Tungsten Carbide Inserts Jaws
- For Lower 1st & 2nd Molars
- Non-slip Grip Gold Plated Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc T/C Universal Extracting Forceps #17 – 1st & 2nd Lower Molars

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Extracting Forceps
APPLICATION: Teeth Extraction
POLISHED: Yes
STERILE: Yes
REUSABLE: Yes
AUTOCCLAVED: Yes
RUST PREVENTION

**WI-70001098****T/C Dental Extraction Forceps #150 Upper Incisors, Root Teeth****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Dental Extracting Forceps #150 Tungsten Carbide Serrated Beaks for Upper Incisors, Canines, and Premolars German Grade Stainless Steel Teeth Extraction Forceps Dental Surgical Instruments.

T/C Universal Extraction Forceps No. 150 are used for extracting teeth from the alveolar bone. Universal conical beaks designed to adapt to upper incisors, canines, premolars, roots, universal right or left. American pattern forceps have a smooth hinge so no pin is visible. The tungsten carbide stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you. High quality stainless


steel, excellent international quality standard.

T/C UNIVERSAL EXTRACTION FORCEPS #150 FEATURES

- For Upper Incisors, Canines & Premolars.
- Serrated Conical Beaks for Optimal Grip.
- Tungsten Carbide Inserts Jaws.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc T/C Universal Extraction Forceps #150 – Serrated Beaks

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental**WI-70001237****T/C Extracting Forceps #53L For Upper 1st & 2nd Left Molars****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Dental Extracting Forceps #53L for Upper 1st and 2nd Left Molars American Pattern German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Tungsten carbide inlaid 53L dental forceps are used to extract maxillary 1st and 2nd left molars. Left angled pointed offset (Bayonet) beaks for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

- Left Pointed Tip Design
- For Upper 1st & 2nd Molars
- Tungsten Carbide Inserts Angular Beaks.
- Non-Slip Textured-Grip Gold Color Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Tungsten Carbide American Pattern Dental Extracting Forceps #53L

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**T/C EXTRACTING FORCEPS #53L FEATURES**

**WI-70001235****T/C Extracting Forceps #53R for Upper 1st & 2nd Right Molars****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Tungsten Carbide (T/C) Extracting Forceps #53R for Upper 1st and 2nd Right Molars American Pattern German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Tungsten Carbide 53R dental forceps are used to extract maxillary 1st and 2nd right molars. Right angled pointed offset (Bayonet) beaks for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

T/C EXTRACTING FORCEPS #53R FEATURES

- Right Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Tungsten Carbide Inserts Angular Beaks.
- Non-Slip Textured-Grip Handle.
- German Stainless Steel.

PACKAGE INCLUDE

- 1 Pc T/C American Pattern Dental Extracting Forceps #53R

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****WI-70001101****T/C Extracting Forceps #65 Upper Roots, Fragments & Incisors****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

American Pattern T/C Dental Extracting Forceps #65 for Upper Roots, Fragments & Incisors, Tungsten Carbide Slim Bayonet Pattern Jaws with Gold Plated Handle German Grade Stainless Steel Teeth Extraction Dental Surgical Instrument.

American pattern tungsten carbide dental extraction forceps No. 65 used to extract teeth from upper roots, fragments & incisors. Slim bayonet pattern beaks designed specially for hard to reach area. The tungsten carbide stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

T/C EXTRACTING FORCEPS #65 FEATURES

- For Upper Roots, Fragments & Incisors.
- Slim, Bayonet Pattern Beaks.
- Tungsten Carbide Inserts Jaws.
- Gold Plated Textured-grip Handle.
- American Pattern.
- Made With High Quality Stainless Steel.


PACKAGE INCLUDE:

- 1 Pc T/C American Pattern Extracting Forceps #65 – Gold Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

APPLICATION: Extraction Forceps

**WI-70002178****T/C Wynman Crown Remover Gripper Forceps
7" Curved Serrated Tips****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Wynman Crown Remover Gripper Forceps 7" (17.78cm) Tungsten Carbide Tips with Gold Ring Handle Premium Quality Stainless Steel Temporary Crown Remover Surgical Dental Forceps Instruments.

Wynman Crown Remover Gripper Forceps features a serrated tip made with durable tungsten carbide tips for a long-lasting working end. Ideal to remove or place temporary crowns with ease, even in the most hard-to-reach places. Our TC Wynman crown grippers are hemostat style forceps that enable a strong grip on the temporary crown without slippage, keeping the crown intact even in the tightest areas of the mouth. Crafted from stainless steel and tungsten carbides,

this instrument offers seamless crown removal solution

**TC WYNMAN CROWN GRIPPER
FORCEPS FEATURES**

TC WYNMAN CROWN GRIPPER FORCEPS: The TC Wynman Crown Gripper Forceps are an ideal addition to your inventory. The tungsten carbide inserts on the Dental Crown Gripper Forceps will grip the temporary crown firmly without slipping, keeping the crown intact. They are useful for gripping wet glass such as in porcelain restorations and assist in more precise procedures.

O.R. GRADE T/C AND STAINLESS STEEL: Hot drop-forged construction Wynman crown gripper constructed of O.R. Grade stainless steel and tungsten carbide for the highest degree of

**WI-70001371****TC Extracting Forceps #1STD Upper Incisors
and Cuspids****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

TC Extracting Forceps #1STD Upper Incisors and Cuspids, Standard, Serrated Beak Extraction Forceps German Grade Stainless Steel Dental Surgical Instruments

#1 Standard Extraction Forceps (V90-1STD) upper incisors and cuspids serrated forceps used for extracting teeth from the alveolar bones. German grade stainless steel promotes longevity and corrosion resistance. Straight, conical beaks designed to adapt upper anterior incisors and canines. Our tools are premium quality and used by professionals.

**TC EXTRACTING FORCEPS 1STD
FEATURES**

- Vantage #1 Standard Extraction Forceps (V90-1STD)

- Upper Incisors and Cuspids Serrated Forceps
- Tungsten Carbide Inserts Jaws
- Non-Slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc TC Extracting Forceps 1STD, Standard – Gold Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Dental
Extracting Forceps

APPLICATION: Extraction Forceps

POLISHED: Yes

STERILE: Yes

REUSABLE: Yes

AUTOCCLAVED: Yes

RUST PREVENTION

PROCEDURE: Passivated





WI-70000410

TC Extraction Forceps #151 for Lower Teeth

\$14.99 per piece

Trade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

American Pattern TC Universal Extraction Forceps #151 Tungsten Carbide Serrated Beaks for Lower Incisors, Canines, Premolars & Roots Premium Quality Stainless Steel Tooth Extracting Surgical Dental Instrument.

TC #151 extracting forceps also known as "lower universal forceps" is a dental tool used for removing teeth in the lower jaw, including incisors, canines, premolars, and roots. American Pattern 151 forceps is designed to be versatile, enabling dentists to extract various teeth in the lower jaw. Tungsten Carbide #151 dental forceps are designed with thin, tapered, and sharp parallel beaks that ensure good contact with the root and crown of the tooth. 151 tungsten carbide tip extraction forcep feature longitudinal serrations

for improved grip and atraumatic extraction. Made of surgical-grade stainless steel for durability and corrosion resistance.

TC EXTRACTION FORCEPS #151 FEATURES

- For Lower Incisors, Canines, Premolars & Roots.
- Tungsten Carbide Inserts Tips.
- Thin, Tapered, and Sharp Parallel Beaks for Good Contact.
- Longitudinal Serrations for Improved Grip and Atraumatic Extraction.
- American Pattern.
- Non-ship Gold Plated Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc TC Universal Extraction Forceps #151 Serrated Beaks Gold Handle



WI-36762663

Thumb Dressing Forceps Teeth Surgical Tissue Pickup Forceps

\$3.00 per piece

Trade price from **\$1.30**

Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Thumb Dressing Forceps Teeth Surgical Tissue Pickup Forceps, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.00, falling to \$1.30 at 25 pieces — about 57% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001372

Tooth Extracting Forceps #1STD Upper Incisors and Cuspids

\$11.99 per piece

Trade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps #1STD Upper Incisors and Cuspids, Standard, Serrated Beak Extraction Forceps German Grade Stainless Steel Dental Surgical Instruments

#1 Standard Extraction Forceps (V90-1STD) upper incisors and cuspids serrated forceps used for extracting teeth from the alveolar bones. German grade stainless steel promotes longevity and corrosion resistance. Straight, conical beaks designed to adapt upper anterior incisors and canines. Our tools are premium quality and used by professionals.

EXTRACTING FORCEPS 1STD FEATURES

- Vantage #1 Standard Extraction Forceps (V90-1STD)

- Upper Incisors and Cuspids Serrated Forceps
- Non-Slip Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Extracting Forceps 1STD, Standard

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Dental Extracting Forceps
APPLICATION: Extraction Forceps
POLISHED: Yes
STERILE: Yes
REUSABLE: Yes
AUTOCCLAVED: Yes
RUST PREVENTION
PROCEDURE: Passivated
ULTRASONIC CLEANED: Yes



WI-10243183

Tooth Extraction Dental Forceps

\$5.99 per piece

Trade price from **\$2.50**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Tooth Extraction Dental Forceps — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$5.99, falling to \$2.50 at 10 pieces — about 58% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-46793830

Tooth Extraction Forceps Dental Set Surgical Tool

\$8.80 per piece

Trade price from \$3.80
Minimum order 25 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Tooth Extraction Forceps Dental Set Surgical Tool is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$8.80, falling to \$3.80 at 25 pieces — about 57% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-36275794

Toothed Dissecting Forceps Stainless Steel Surgical Tissue Forceps

\$1.20 per piece

Trade price from \$0.60
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Toothed Dissecting Forceps Stainless Steel Surgical Tissue Forceps is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$1.20, falling to \$0.60 at 5 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-36277895

Toothed Dissecting Forceps Stainless Steel Tissue

\$1.20 per piece

Trade price from \$0.60
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Toothed Dissecting Forceps Stainless Steel Tissue for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$1.20, falling to \$0.60 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-27296632

Toothed Forceps 6" for Hospitals-Brown Adson Surgical

\$2.39 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Toothed Forceps 6" for Hospitals-Brown Adson Surgical is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.39, falling to \$2.00 at 10 pieces — about 16% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



WI-42901035

Universal Dental Extracting Forceps Stainless Steel Oral Surgery Tool

\$9.00 per piece

Trade price from \$4.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Universal Dental Extracting Forceps Stainless Steel Oral Surgery Tool — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$9.00, falling to \$4.00 at 1 piece — about 56% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004120

Universal Extraction Forceps #151 for Lower Teeth

\$11.99 per piece

Trade price from \$11.03
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

American Pattern Universal Extraction Forceps #151 Serrated Beaks for Lower Incisors, Canines, Premolars & Roots – Premium Quality Stainless Steel Tooth Extracting Surgical Dental Instrument.

#151 extracting forceps also known as “lower universal forceps” is a dental tool used for removing teeth in the lower jaw, including incisors, canines, premolars, and roots. American Pattern #151 forceps is designed to be versatile, enabling dentists to extract various teeth in the lower jaw. #151 dental forceps are designed with thin, tapered, and sharp parallel beaks that ensure good contact with the root and crown of the tooth. 151 dental extracting forceps feature longitudinal serrations for improved grip and atraumatic extraction, and are typically made of

surgical-grade stainless steel for durability and corrosion resistance.

AMERICAN EXTRACTION FORCEPS #151 FEATURES

- For Lower Incisors, Canines, Premolars & Roots.
- Thin, Tapered, and Sharp Parallel Beaks for Good Contact.
- Longitudinal Serrations for Improved Grip and Atraumatic Extraction.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Universal Extraction Forceps #151

TECHNICAL INFORMATION

MATERIAL: Stainless Steel



**WI-70000409**

Universal Extraction Forceps #151 Green Handle

\$11.99 per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Extracting Forceps #151 Green Color Coated Handles with Serrated Beaks for Lower Incisors, Canines, Premolars & Roots Premium Quality Stainless Steel American Pattern Tooth Extracting Surgical Dental Instrument.

#151 extracting forceps also known as "lower universal forceps" is a dental tool used for removing teeth in the lower jaw, including incisors, canines, premolars, and roots. American Pattern 151 forceps is designed to be versatile, enabling dentists to extract various teeth in the lower jaw. Green handle #151 dental forceps are designed with thin, tapered, and sharp parallel beaks that ensure good contact with the root and crown of the tooth. 151 extraction forcep feature longitudinal serrations for improved grip and atraumatic

extraction. Made of surgical-grade stainless steel for durability and corrosion resistance.

EXTRACTION FORCEPS #151 GREEN HANDLE FEATURES

- For Lower Incisors, Canines, Premolars & Roots.
- Thin, Tapered and Sharp Parallel Beaks for Good Contact.
- Longitudinal Serrations for Improved Grip and Atraumatic Extraction.
- American Pattern.
- Non-ship Green Color Plated Grip Handle.
- Made With Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Extraction Forceps #151 Serrated Beaks Green Handle

**WI-70001252**

Universal Extraction Forceps #17 with Gold Hollow Python Grip Handle

\$14.99 per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality American Pattern Universal Extracting Forceps #17 For Lower 1st & 2nd Molars with Gold Hollow Python Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

#17 extraction forceps use for 1st and 2nd lower molars holding with serrated jaws. The working end of this tooth extraction instrument is designed with a sharp bend, making sure it reaches deep into the socket for a strong grip. Aesthetically crafted python grip hollow handles offers better grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #17 FEATURES

- Tri-pointed Beak Design
- For Lower 1st & 2nd Molars
- Gold Hollow Python Grip Handle
- German Stainless Steel

PACKAGE INCLUDE

- 1 Pc Universal Extracting Forceps #17 with Gold Hollow Python Grip Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION**



WI-70004161

Universal Pedo Extracting Forceps #17S for Lower Molars

\$11.99 per piece

Trade price from \$11.03
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Pedo Extraction Forceps #17S For Lower Molars with Tri-Pointed Beaks & Spring Handle Premium Quality German Grade Stainless Steel American Pattern Tooth Extraction Dental Surgical Instruments

Pedodontic 17S pediatric dental forceps are specifically designed for extracting lower molars, featuring sharp tri-pointed tips that ensure a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

EXTRACTING FORCEPS #17S FEATURES

- Universal Pedo Extraction Forceps.

- Tri-pointed Beak Design.
- For Lower Molars.
- Non-slip Grip Spring Handle.
- German Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Pedo Extracting Forceps #17S – Lower Molars

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Dental Extracting Forceps

APPLICATION: Childern Teeth Extraction

STERILE: Non-Sterile

REUSABLE/AUTOCCLAVABLE: Yes

RUST PREVENTION

PROCEDURE: Passivated

ULTRASONIC CLEANED: Yes

TESTS PERFORMED: Boil Test,



WI-06803073

Vulsellum Forceps 3x4 Teeth's

\$3.50 per piece

Trade price from \$1.85
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Vulsellum Forceps 3×4 Teeth's, made for dental, hygiene and oral surgery work. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$3.50, falling to \$1.85 at 10 pieces — about 47% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

Request a quotation

Send us product codes from this catalogue, quantities and your country.

We reply the same working day on WhatsApp.

Online shop	shop.waqarintl.com
Email	sales@waqarintl.com
WhatsApp /	
Phone	+92 307 6173627
Works	Roras Road, Opposite Anwar Khawaja Industries, Sialkot 51310, Pakistan

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.