



# WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

# Instrument Sets

## Dental Instruments

Reusable stainless-steel instrument sets forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



10  
PCS



10  
PCS



10  
PCS

**216**

PRODUCTS

**\$1.50**

PRICES FROM

**OEM**

PRIVATE LABEL

**Global**

WORLDWIDE  
EXPORT

# Contents

216 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

| Product                                                                   | Code        | Price    | Page |
|---------------------------------------------------------------------------|-------------|----------|------|
| 10 Extracting Forceps Set Tooth Extraction Dental Instruments             | WI-70001994 | \$89.99  | 8    |
| 10 Manual Crown Remover with 3 Tips Double Ended Set                      | WI-70002034 | \$187.49 | 8    |
| 10 Pcs Amalgam Carrier Double Ended Set Large/Medium 2.5mm/2mm            | WI-70004257 | \$69.99  | 9    |
| 10 Pcs Ball Burnisher #34 Single Ended                                    | WI-70000703 | \$39.99  | 9    |
| 10 Pcs Dental Amalgam Well Non-Slip Restorative Carrier                   | WI-70001442 | \$109.99 | 10   |
| 10 Pcs Dental Elevators Set – Periosteal Apical Molt Oral Implant         | WI-70002287 | \$35.99  | 10   |
| 10 PCS Dental Orthodontic Instrument Set                                  | WI-42944644 | \$38.00  | 11   |
| 10 PCS Dental Spoon Excavator Set Composite Filling                       | WI-35370681 | \$38.00  | 11   |
| 10 PCS German Instrument Sets                                             | WI-25126184 | \$4.00   | 12   |
| 10 PCS Hollow Handle Dental Inspection Mirror                             | WI-42929636 | \$38.00  | 12   |
| 10 PCS Replaceable Dental Mouth Mirror Set                                | WI-42919953 | \$38.00  | 13   |
| 10 Pcs Wax Carving Spatula Modeling Kit                                   | WI-70002724 | \$19.99  | 13   |
| 10 Periosteal Elevators Set                                               | WI-70000073 | \$39.99  | 14   |
| 100 Pcs Ball Burnisher #34 Single Ended                                   | WI-70000701 | \$149.99 | 14   |
| 12 P.K Thomas #2 Dental Wax Modelling Carver Set                          | WI-70003937 | \$49.99  | 15   |
| 12 Pcs Amalgam Plugger Black #1/2 Double Ended                            | WI-70003940 | \$71.99  | 15   |
| 12 Pcs American Pattern Extracting Forceps No. 17 Lower Molars            | WI-70004165 | \$99.99  | 16   |
| 12 Pcs American Pado Extracting Forceps #16S Lower Molars                 | WI-70004129 | \$99.99  | 16   |
| 12 Pcs Ball Burnisher #28/29 Set Double Ended 2.3mm – 3.6mm Football Tips | WI-70004278 | \$59.99  | 17   |
| 12 Pcs Bone File Set #12A Double Ended                                    | WI-70004104 | \$69.99  | 17   |
| 12 Pcs Castroviejo Caliper Long 0-40mm Angled 6"                          | WI-70001466 | \$155.99 | 18   |
| 12 Pcs Castroviejo Caliper Long 0-40mm Angled 7"                          | WI-70001457 | \$164.99 | 18   |
| 12 Pcs Castroviejo Caliper Long 0-40mm Straight 7"                        | WI-70001461 | \$164.99 | 19   |
| 12 Pcs Cement Spatulas 7" (17.8cm) Double Ended                           | WI-70003997 | \$69.99  | 19   |
| 12 Pcs Cement Spatulas Single Ended                                       | WI-70000731 | \$64.99  | 20   |
| 12 Pcs Dental Ball Burnisher #1 Double Ended                              | WI-70004113 | \$49.99  | 20   |
| 12 Pcs Dental Bone File #12CA Double Ended                                | WI-70004103 | \$69.99  | 21   |
| 12 Pcs Dental Bone File #33 Double Ended                                  | WI-70004102 | \$69.99  | 21   |
| 12 Pcs Dental Extracting Forceps #53L For Upper 1st & 2nd Left Molars     | WI-70004171 | \$99.99  | 22   |
| 12 Pcs Dental Extracting Forceps Tooth Extraction Instruments Kit         | WI-70000630 | \$109.99 | 22   |
| 12 Pcs Dental Hollenback Wax Carver #3                                    | WI-70003939 | \$65.99  | 23   |
| 12 PCS Dental Mouth Mirror Set                                            | WI-42941584 | \$9.00   | 23   |
| 12 Pcs Extracting Forceps Upper & Lower Dental Instruments                | WI-70000971 | \$99.99  | 24   |
| 12 Pcs Miller Bone File #45 Double Ended                                  | WI-70000685 | \$69.99  | 24   |
| 12 Pcs Molt Bone Curette Set #2/4 Double Ended                            | WI-70000537 | \$83.99  | 25   |
| 12 Pcs Periosteal Elevator Baby Molt #7A Double Ended with Round Handle   | WI-70003903 | \$59.99  | 25   |
| 12 Pcs Polokoff Bone File #1 Double Ended                                 | WI-70000682 | \$65.99  | 26   |
| 12 Pcs Spring Divider Calipers Set 4.5" Dental Lab Tool                   | WI-70003956 | \$84.99  | 26   |
| 12 Pcs Star-Like Boone Bracket Positioning Gauge Set                      | WI-70003983 | \$72.09  | 27   |
| 12 Pcs T/C Wynman Crown Remover Gripper Forceps 7" Curved Serrated Tips   | WI-70002175 | \$126.99 | 27   |
| 15 Pcs Dental Elevators Instruments Set                                   | WI-70002288 | \$49.99  | 28   |

## Contents (continued)

| Product                                                                  | Code        | Price    | Page |
|--------------------------------------------------------------------------|-------------|----------|------|
| 2 Extracting Forceps Set #53L & #53R Upper Left & Right Molars           | WI-70001233 | \$18.99  | 29   |
| 2 Extracting Forceps Set #88R & #88L for Upper 1st & 2nd Molars          | WI-70000973 | \$19.99  | 29   |
| 2 Iwanson Calipers Set 4" Sharp and Blunt Tips Improved Model            | WI-70001449 | \$15.99  | 30   |
| 2 Pcs Amalgam Carriers Double Ended + Amalgam Well Set                   | WI-70004262 | \$23.99  | 30   |
| 2 Pcs Amalgam Carriers Single Ended + Amalgam Well Set                   | WI-70004261 | \$21.99  | 31   |
| 2 Pcs Ball Burnisher #28/29 Set Double Ended 2.3mm – 3.6mm Football Tips | WI-70004260 | \$14.99  | 31   |
| 2 Pcs Christensen Crown Remover Set CRCH2 90 Degree Angled               | WI-70002028 | \$17.99  | 32   |
| 2 Pcs Cryer Root Elevator Set #25 & #26                                  | WI-70000001 | \$15.99  | 32   |
| 2 Pcs Cryer Root Elevator Set #44 & #45                                  | WI-70000002 | \$15.99  | 33   |
| 2 Pcs Dental Amalgam Well Non-Slip Restorative Carrier                   | WI-70001444 | \$26.99  | 33   |
| 2 Pcs Dental Boley Gauge & Iwanson Gauge Set Vernier Lab Instrument      | WI-70003969 | \$17.99  | 34   |
| 2 Pcs Dental Bone File #12CA Double Ended                                | WI-70004070 | \$15.99  | 34   |
| 2 Pcs Dental Bone File #33 Double Ended                                  | WI-70004075 | \$15.99  | 35   |
| 2 Pcs Dental Bone Files #12A Double Ended                                | WI-70004071 | \$15.99  | 35   |
| 2 Pcs Seldin Elevator 4L – 10mm Tip                                      | WI-70002143 | \$12.99  | 36   |
| 2 Periotome Luxating Dental Elevators Set Del #8 and #A                  | WI-70002284 | \$12.99  | 36   |
| 2 T/C Extracting Forceps Set #53L & #53R Upper Left & Right Molars       | WI-70001234 | \$24.99  | 37   |
| 20 Pcs Amalgam Restorative Set with Sterilization Cassette               | WI-70004276 | \$164.99 | 37   |
| 2023 Dental Hygiene Kit Oral Care Kit Pakistan Made                      | WI-13357194 | \$7.00   | 38   |
| 3 Bone File Kits of 5 Pcs Dental Surgical Orthopedic Veterinary Surgery  | WI-70004094 | \$109.99 | 38   |
| 3 Bone File Kits of 8 Pcs Surgical Dental Orthopedic Veterinary          | WI-70004093 | \$179.99 | 39   |
| 3 Heidbrink Root Tip Elevators Set Straight, Left, Right Angled Tips     | WI-70003888 | \$23.99  | 39   |
| 3 P.K Thomas #2 Dental Wax Modelling Carver Set                          | WI-70003930 | \$20.99  | 40   |
| 3 Pcs Amalgam Carrier Double Ended Set Large/Medium 2.5mm/2mm            | WI-70004256 | \$26.99  | 40   |
| 3 Pcs Amalgam Carrier Single Ended Set 1.5mm Mini Tip                    | WI-70004275 | \$23.99  | 41   |
| 3 Pcs Amalgam Carrier Single Ended Set 2.0mm Regular Tapered Tip         | WI-70000527 | \$23.99  | 41   |
| 3 Pcs Amalgam Carrier Single Ended Set 3.0mm Large Tapered Tip           | WI-70000524 | \$23.99  | 42   |
| 3 Pcs American Extracting Forceps #90 Upper Left Molars                  | WI-70004216 | \$29.99  | 42   |
| 3 Pcs American Pedro Extracting Forceps #16S Lower Molars                | WI-70004224 | \$29.99  | 43   |
| 3 Pcs Ball Burnisher #28/29 Set Double Ended 2.3mm – 3.6mm Football Tips | WI-70004272 | \$19.99  | 43   |
| 3 Pcs Boley Gauge 10cm Measuring 6" (15.24cm)                            | WI-70003980 | \$26.99  | 44   |
| 3 Pcs Bone File Set #12A Double Ended                                    | WI-70004097 | \$23.99  | 44   |
| 3 Pcs Castroviejo Caliper Long 0-40mm Angled 6"                          | WI-70001468 | \$47.99  | 45   |
| 3 Pcs Castroviejo Caliper Long 0-40mm Angled 7"                          | WI-70001459 | \$49.99  | 45   |
| 3 Pcs Castroviejo Caliper Long 0-40mm Straight 6"                        | WI-70003954 | \$47.99  | 46   |
| 3 Pcs Castroviejo Caliper Long 0-40mm Straight 7"                        | WI-70001463 | \$49.99  | 46   |
| 3 Pcs Cement Spatulas 7" (17.8cm) Double Ended                           | WI-70003995 | \$25.99  | 47   |
| 3 Pcs Cement Spatulas Single End                                         | WI-70003993 | \$23.99  | 47   |
| 3 Pcs Dental Bone File #12CA Double Ended                                | WI-70004096 | \$23.99  | 48   |
| 3 Pcs Dental Bone File #33 Double Ended                                  | WI-70004095 | \$23.99  | 48   |
| 3 Pcs Dental Extracting Forceps #150AS Upper Anteriors                   | WI-70004228 | \$29.99  | 49   |
| 3 Pcs Dental Extracting Forceps #210H Upper 3rd Molars                   | WI-70004234 | \$29.99  | 49   |

## Contents (continued)

| Product                                                                 | Code        | Price    | Page |
|-------------------------------------------------------------------------|-------------|----------|------|
| 3 Pcs Dental Extracting Forceps #23 Cow Horn for Lower Molar            | WI-70004241 | \$29.99  | 50   |
| 3 Pcs Dental Extracting Forceps #32                                     | WI-70004251 | \$29.99  | 50   |
| 3 Pcs Dental Extracting Forceps #53L For Upper 1st & 2nd Left Molars    | WI-70004215 | \$29.99  | 51   |
| 3 Pcs Dental Extracting Forceps #65 for Upper Roots                     | WI-70004248 | \$29.99  | 51   |
| 3 Pcs Dental Extracting Forceps 32A Molars Root – Slim Beak             | WI-70004249 | \$29.99  | 52   |
| 3 Pcs Dental Orthodontic Star-Like Boone Bracket Positioning Gauge Set  | WI-70003978 | \$21.99  | 52   |
| 3 Pcs Dental Pick Tools Set Plaque Remover Scaler                       | WI-25737862 | \$1.50   | 53   |
| 3 Pcs Extracting Forceps #53R for Upper 1st & 2nd Right Molars          | WI-70004240 | \$29.99  | 53   |
| 3 Pcs Extracting Forceps #62 Upper and Lower Anterior Teeth             | WI-70004163 | \$29.99  | 54   |
| 3 Pcs Extracting Forceps Set #151AS Serrated Split Beaks                | WI-70004225 | \$29.99  | 54   |
| 3 Pcs Heidbrink Root Tip Pick Elevator Set H1, H2 & H3 Pointed Tips     | WI-70000641 | \$23.99  | 55   |
| 3 Pcs Hollenback Wax Carver #3                                          | WI-70003932 | \$19.99  | 55   |
| 3 Pcs Miller Bone File #45 Double Ended                                 | WI-70004092 | \$23.99  | 56   |
| 3 Pcs Molt Bone Curette Set #2/4 Double Ended                           | WI-70000539 | \$23.99  | 56   |
| 3 Pcs Pado Dental Extracting Forceps 150S                               | WI-70004214 | \$29.99  | 57   |
| 3 Pcs Periosteal Elevator Baby Molt #7A Double Ended with Round Handle  | WI-70003901 | \$19.99  | 57   |
| 3 Pcs Periosteal Seldin #S23 Dental Elevator Double Ended               | WI-70003902 | \$23.99  | 58   |
| 3 Pcs Polokoff Bone File #1 Double Ended                                | WI-70000684 | \$23.99  | 58   |
| 3 Pcs Spring Divider Calipers Set 4.5" Dental Lab Tool                  | WI-70003953 | \$26.99  | 59   |
| 3 Pcs T/C Wynman Crown Remover Gripper Forceps 7" Curved Serrated Tips  | WI-70002177 | \$34.49  | 59   |
| 3 Pcs Warwick-James Elevators Set Right, Left & Straight                | WI-70000931 | \$23.99  | 60   |
| 3 Root Tip Pick Elevators Set #78 #79 & #80                             | WI-70002285 | \$18.74  | 60   |
| 3 Swivel Head Bracket Height Gauge Set 6" Adjustable Pointed Tips       | WI-70000542 | \$34.99  | 61   |
| 3 Warwick-James Elevators Set with Sterilization Cassette               | WI-70003887 | \$49.99  | 61   |
| 30 Pcs Amalgam Carrier Double Ended Set Large/Medium 2.5mm/2mm          | WI-70004270 | \$179.99 | 62   |
| 4 Dental Amalgam Carrier Double Ended Set                               | WI-70000517 | \$35.99  | 62   |
| 4 Lucas Bone Curette Set #85, #86, #87, #88 Double Ended                | WI-70000350 | \$23.99  | 63   |
| 4 Pcs Ball Burnishers Set #27/29 Double Ended                           | WI-70004281 | \$25.99  | 63   |
| 4 Pcs Dental Amalgam Well Non-Slip Restorative Carrier                  | WI-70001443 | \$51.99  | 64   |
| 4 Pcs Spatulas Wax & Clay Sculpting Tool Carver Set                     | WI-70001740 | \$9.99   | 64   |
| 4 Pcs Wax & Modelling LeCron Zahle Beale Carver Cement Spatula Set      | WI-70003919 | \$16.99  | 65   |
| 5 Lucas Bone Curette Set #84, #85, #86, #87, #88 Double Ended           | WI-70000533 | \$30.99  | 65   |
| 5 Pcs Ball Burnisher #34 Single Ended                                   | WI-70000704 | \$24.99  | 66   |
| 5 Pcs Bone Files Kit Dental Surgical Orthopedic Veterinary Surgery      | WI-70004074 | \$39.99  | 66   |
| 5 Pcs Christensen Crown Remover Set CRCH2 90 Degree Angled              | WI-70002027 | \$43.99  | 67   |
| 5 Pcs Dental Elevators #301, #34, #34S, #46 & #81 Set                   | WI-70002286 | \$24.99  | 67   |
| 5 Pcs Periosteal Dental Elevator Set                                    | WI-70001922 | \$14.99  | 68   |
| 5 Pcs Wax Carving Metal Clay Pottery Sculpting Tool Set                 | WI-70002722 | \$10.99  | 68   |
| 50 Pcs Ball Burnisher #34 Single Ended                                  | WI-70000702 | \$124.99 | 69   |
| 50 Pcs Football Burnisher #28/29 Set 2.3mm – 3.6mm Double Ended         | WI-70000603 | \$124.99 | 69   |
| 5-Piece Stainless Steel Double-Sided Cement Spatula Set Ergonomic       | WI-30694736 | \$5.00   | 70   |
| 6 Bone File Kits of 5 Pcs Dental Surgical Orthopedic Veterinary Surgery | WI-70004084 | \$199.99 | 70   |

## Contents (continued)

| Product                                                                  | Code        | Price   | Page |
|--------------------------------------------------------------------------|-------------|---------|------|
| 6 Bone Rongeurs Surgical Dental Oral Surgery Set                         | WI-70004017 | \$89.99 | 71   |
| 6 P.K Thomas #2 Dental Wax Modelling Carver Set                          | WI-70003924 | \$35.99 | 71   |
| 6 Pcs Amalgam Carrier Double Ended Set Large/Medium 2.5mm/2mm            | WI-70004258 | \$47.99 | 72   |
| 6 Pcs Amalgam Carrier Single Ended Set 1.5mm Mini Tip                    | WI-70004269 | \$41.99 | 72   |
| 6 Pcs Amalgam Carrier Single Ended Set 2.0mm Regular Tapered Tip         | WI-70000526 | \$41.99 | 73   |
| 6 Pcs Amalgam Carrier Single Ended Set 3.0mm Large Tapered Tip           | WI-70000523 | \$41.99 | 73   |
| 6 Pcs American Extracting Forceps #90 Upper Left Molars                  | WI-70004173 | \$59.99 | 74   |
| 6 Pcs American Pedro Extracting Forceps #16S Lower Molars                | WI-70004182 | \$59.99 | 74   |
| 6 Pcs Ball Burnisher #28/29 Set Double Ended 2.3mm – 3.6mm Football Tips | WI-70004266 | \$35.99 | 75   |
| 6 Pcs Bone File Set #12A Double Ended                                    | WI-70004087 | \$43.99 | 75   |
| 6 Pcs Castroviejo Caliper Long 0-40mm Angled 6"                          | WI-70001467 | \$89.99 | 76   |
| 6 Pcs Castroviejo Caliper Long 0-40mm Angled 7"                          | WI-70001458 | \$99.99 | 76   |
| 6 Pcs Castroviejo Caliper Long 0-40mm Straight 6"                        | WI-70003951 | \$89.99 | 77   |
| 6 Pcs Castroviejo Caliper Long 0-40mm Straight 7"                        | WI-70001462 | \$99.99 | 77   |
| 6 Pcs Cement Spatulas 7" (17.8cm) Double Ended                           | WI-70003996 | \$45.99 | 78   |
| 6 Pcs Cement Spatulas Single Ended                                       | WI-70003992 | \$43.99 | 78   |
| 6 Pcs Dental Bone File #12CA Double Ended                                | WI-70004086 | \$43.99 | 79   |
| 6 Pcs Dental Bone File #33 Double Ended                                  | WI-70004085 | \$43.99 | 79   |
| 6 Pcs Dental Extracting Forceps #210H Upper 3rd Molars                   | WI-70004193 | \$59.99 | 80   |
| 6 Pcs Dental Extracting Forceps #23 Cow Horn for Lower Molar             | WI-70004199 | \$59.99 | 80   |
| 6 Pcs Dental Extracting Forceps #32                                      | WI-70004210 | \$59.99 | 81   |
| 6 Pcs Dental Extracting Forceps #53L For Upper 1st & 2nd Left Molars     | WI-70004172 | \$59.99 | 81   |
| 6 Pcs Dental Extracting Forceps 32A Molars Root – Slim Beak              | WI-70004208 | \$59.99 | 82   |
| 6 Pcs Dental Extracting Forceps Set #65 for Upper Roots                  | WI-70004207 | \$59.99 | 82   |
| 6 Pcs Dental Hollenback Wax Carver #3                                    | WI-70003926 | \$35.99 | 83   |
| 6 Pcs Dental Orthodontic Boone Bracket Positioning Gauge Set             | WI-70003974 | \$41.99 | 83   |
| 6 Pcs Dental Tooth Extracting Forceps #10S with Serrated Jaws            | WI-70004203 | \$59.99 | 84   |
| 6 Pcs Extracting Forceps #53R for Upper 1st & 2nd Right Molars           | WI-70004198 | \$59.99 | 84   |
| 6 Pcs Extracting Forceps #62 Upper and Lower Anterior Teeth              | WI-70004177 | \$59.99 | 85   |
| 6 Pcs Extracting Forceps Set #151AS Serrated Split Beaks                 | WI-70004183 | \$59.99 | 85   |
| 6 Pcs Lesmann Knife Set 17cm – Large                                     | WI-70003905 | \$36.99 | 86   |
| 6 Pcs Miller Bone File #45 Double Ended                                  | WI-70004083 | \$43.99 | 86   |
| 6 Pcs Molt Bone Curette Set #2/4 Double Ended                            | WI-70000538 | \$44.99 | 87   |
| 6 Pcs Periosteal Elevator Baby Molt #7A Double Ended with Round Handle   | WI-70003898 | \$35.99 | 87   |
| 6 Pcs Polokoff Bone File #1 Double Ended                                 | WI-70000683 | \$43.99 | 88   |
| 6 Pcs Spring Divider Calipers Set 4.5" Dental Lab Tool                   | WI-70003950 | \$47.99 | 88   |
| 6 Pcs T/C Wynman Crown Remover Gripper Forceps 7" Curved Serrated Tips   | WI-70002176 | \$69.99 | 89   |
| 6pcs Kit Dental Composite Resin Filling Spatulas Double Head             | WI-34589234 | \$7.50  | 89   |
| 7 Pcs PDL Luxating Root Elevator Color Handle Set                        | WI-70002756 | \$59.99 | 90   |
| 7 Pcs Wax Carving Metal Clay Pottery Sculpting Tool Set                  | WI-70002723 | \$14.99 | 90   |
| 8 Pcs Bone File Set for Surgical Dental Orthopedic Veterinary            | WI-70004080 | \$65.99 | 91   |
| 8 Pcs Stainless Steel Color Dental Tool                                  | WI-32307453 | \$4.00  | 91   |

## Contents (continued)

| Product                                                                    | Code        | Price    | Page |
|----------------------------------------------------------------------------|-------------|----------|------|
| 8-Piece Precision Dental Instruments Endodontic Root Canal                 | WI-24322645 | \$8.00   | 92   |
| 9 Pcs Wax Carver Set – Double Ended                                        | WI-70000328 | \$16.99  | 92   |
| Adult Jennings Mouth Probe Gag Stainless Steel Oral                        | WI-36815967 | \$20.00  | 93   |
| Ainsworth Rubber Dam Punch Stainless Steel Dental Instrument               | WI-35031143 | \$4.10   | 93   |
| Amalgam Carriers Stainless Steel Single/Double Ended Syringe               | WI-30844027 | \$3.00   | 94   |
| Basic Orthodontics Dental Instrument Set                                   | WI-32307440 | \$4.00   | 94   |
| Bone Scraper Harvesters Implant Surgery Set Grafting Dental Surgical       | WI-34359269 | \$120.00 | 95   |
| Bracket Positioning Star Boone & Height Gauge Wick Type Orthodontic Set    | WI-70003964 | \$14.99  | 95   |
| Complete Instrument Sets                                                   | WI-32031427 | \$7.00   | 96   |
| Dental Autoclave Stainless Steel Cassette Tray Rack Box                    | WI-24379958 | \$15.00  | 96   |
| Dental Equipment Kit 10pcs Power Silicone Handle Composite                 | WI-24513101 | \$25.00  | 97   |
| Dental Hollowware Surgical Instruments                                     | WI-41184859 | \$15.00  | 97   |
| Dental Hygiene Kit Calculus Plaque Remover                                 | WI-09615268 | \$3.50   | 98   |
| Dental Instruments Composite Resin Filling Spatula Filler Hygiene Kit      | WI-36904798 | \$8.00   | 98   |
| Dental Micro Probe Scaler Explorer Set Stainless Steel Double              | WI-42839665 | \$30.00  | 99   |
| Dental Trephine Drills Instrument Set Stainless Steel Color                | WI-30655614 | \$40.00  | 99   |
| Double Ended Stainless Steel Dental Scaler Instrument                      | WI-32309377 | \$4.00   | 100  |
| Fork Pattern Articulator Nickel Chromium Plate                             | WI-30660169 | \$18.00  | 100  |
| German Stainless Steel Adjustable Dental Articulator Holder Instrument     | WI-30598718 | \$18.00  | 101  |
| Implant Buddy Stainless Steel Dental Implant Case Kit                      | WI-34294513 | \$120.00 | 101  |
| Maxillofacial Surgery Instruments Set Stainless Steel Oral                 | WI-36131731 | \$190.00 | 102  |
| Medical Devices Orthopedic Reusable Dental Instruments Vet Instruments Kit | WI-34083176 | \$600.00 | 102  |
| Morrel Crown Remover Set with 3 Tips                                       | WI-70002039 | \$10.99  | 103  |
| Most Selling Stainless Steel Dental Scaler Instrument                      | WI-32283577 | \$4.00   | 103  |
| Oral Dental Morrel Crown Bridges Removal Kit with Holding Forcep           | WI-70003295 | \$18.74  | 104  |
| Orthodontic Ligature Pen Dental Set                                        | WI-46731692 | \$7.00   | 104  |
| Orthopedic Dog Dental Instruments Reusable AC 110-220V Vet Instrument      | WI-34057331 | \$600.00 | 105  |
| Orthopedic Surgical Instruments Dental Instrument Set                      | WI-34382078 | \$120.00 | 105  |
| Pack 3 Stainless Steel Dental Boxes Rectangular Small                      | WI-33797689 | \$14.00  | 106  |
| Pack of 12 PCS Dental Teeth Inspection Mirrors                             | WI-42923809 | \$9.00   | 106  |
| Periodontics Dental Scaler Instrument Set                                  | WI-32303451 | \$4.00   | 107  |
| Ready Ship Micro Periodontal Oral Surgery Kit Surgical                     | WI-26215650 | \$80.00  | 107  |
| Ridge Splitting Chisel Dental Set for Oral Surgery Bone                    | WI-46803605 | \$6.90   | 108  |
| Rubber Dam Kit Dental Isolation Instrument Set                             | WI-42858541 | \$8.00   | 108  |
| Rubber Dam Punch Board Practice Tool Set Dental                            | WI-30479857 | \$28.00  | 109  |
| Set of 10 Dental Spoon Excavators Blue Titanium Coated                     | WI-35330785 | \$38.00  | 109  |
| Set of 10 PCS Hollow Handle Dental Mouth Mirrors                           | WI-42919952 | \$38.00  | 110  |
| Set of 4 Miller Dental Root Extracting East West Right Left Elevators      | WI-70003889 | \$29.99  | 110  |
| Stainless Steel Dental Articulator                                         | WI-30564806 | \$18.00  | 111  |
| Stainless Steel Dental Scaler Instrument Set Dentistry                     | WI-32279576 | \$4.00   | 111  |
| Stainless Steel Dental Scaler Set 6pc                                      | WI-32274735 | \$4.00   | 112  |
| Stainless Steel Dental Scalers                                             | WI-32309354 | \$4.00   | 112  |

## Contents (continued)

| Product                                                   | Code        | Price   | Page |
|-----------------------------------------------------------|-------------|---------|------|
| Stainless Steel Dental Set Sterilization Tray for Medical | WI-36839984 | \$12.00 | 113  |
| Stainless Steel Dental Trephine Drills Implant Kit        | WI-30583793 | \$40.00 | 113  |
| Stainless Steel Dental Trephine Drills Kit Free Bur Block | WI-30622697 | \$40.00 | 114  |
| Stainless Steel Reusable Surgical Dental Suture           | WI-27358659 | \$12.00 | 114  |
| Stainless Steel Wax Carver Set Art Supplies Kit Pottery   | WI-25115209 | \$4.00  | 115  |
| Steel Dental Tooth Cleaning Scalers                       | WI-32299493 | \$4.00  | 115  |

**WI-70001994**

## 10 Extracting Forceps Set Tooth Extraction Dental Instruments

**\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Dental Extracting Forceps  
Teeth Extraction Tools Textured-Grip  
Handle Stainless Steel Dental Surgical  
Instruments

### 10 PCS DENTAL EXTRACTING FORCEPS SET

You can make 10 extracting Forceps set of your choice. We have following forceps available and you can choose any 10 forceps;

1HEN, 1STD, 3FS, 6, 10H, 10S, 13, 15, 16, 16S, 17, 17S, 18L, 18R, 22, 23, 23S, 24, 32, 32A, 33, 53L, 53R, 62, 65, 69, 73, 74, 74N, 85A, 88L, 88R, 89, 90, 99C, 101, 103, 116, 120, 150, 150 1/2S, 150A, 150AS, 150S, 150XAS, 151, 151 1/2S, 151A, 151AS, 151S, 151XAS, 203, 210H, 210S, 217, 222, 286, A, B, E, F, H, MD1, MD2, MD3 and MD4

**NOTE:** Please specify when ordering what number of forceps you need.

### TECHNICAL INFORMATION

**MATERIAL:** Stainless Steel

**INSTRUMENTS TYPE:** Dental  
Extracting Forceps

**APPLICATION:** Extraction Forceps

**POLISHED:** Yes

**STERILE:** Yes

**REUSABLE:** Yes

**AUTOCCLAVED:** Yes

**RUST PREVENTION**

**PROCEDURE:** Passivated

**ULTRASONIC CLEANED:** Yes

**TESTS PERFORMED:** Boil Test,  
Performance Test, Shape Test

**QC PASSED:** Yes

**PACKING DETAILS:** Polybag Packing

**WI-70002034**

## 10 Manual Crown Remover with 3 Tips Double Ended Set

**\$187.49** per pieceTrade price from **\$172.49**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Manual Crown Remover  
Complete Set with 3 Attachment Tips  
Crown and Bridge Removers  
Restorative Dentistry Instruments.

Manual dental crown bridge remover kit double end is made of stainless steel and comes with 3 tips. Classic, easy-to-use crown remover is perfect for removing temporary or permanent restorations. The crown remover tips attach to the working ends of crown removers to grip the crown securely without slipping. This crown remover set is autoclavable and crafted from high-quality stainless steel ensuring durability and sustained performance. 100% satisfaction guaranteed

- For use in all areas of the mouth
- Manual, hammer mechanism
- Double Ended
- Comes with 3 tempered stainless steel tips
- Reusable
- Made with Highest quality medical grade Stainless Steel

### PACKAGE INCLUDE:

- 10 Pcs Manual Crown Remover with 3 Attachment Tips – Double Ended

### TECHNICAL INFORMATION:

**MATERIAL:** Stainless Steel

**INSTRUMENTS TYPE:** Restorative  
Dental Instruments

**APPLICATION:** Permanent and  
Temporary Crowns and Bridges  
Removal

**TIPS SHAPE:** Hook Shape / Curved



**10  
PCS****WI-70004257****10 Pcs Amalgam Carrier Double Ended Set  
Large/Medium 2.5mm/2mm****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Pcs Amalgam Carrier Double Ended Set 8" (20.3cm) Large/Medium 2.5mm – 2mm Tapered Tips with Non-Slip Textured Grip Handle Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.


Dental Amalgam Filling Carrier is a durable, high-quality instrument used to fill dental cavities with amalgam.

Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam

carrier. The textured surface helps effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

**DOUBLE ENDED AMALGAM  
CARRIER FEATURES**

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Color Bands for Easy Instruments Identification.
- Large – Regular (2.5mm – 2mm) Tapered Tips for Even Distribution of Filling Material.

**100% HIGH QUALITY  
STAINLESS STEEL****10  
PCS****WI-70000703****10 Pcs Ball Burnisher #34 Single Ended****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Pcs Ball Burnisher #34 Single Ended 6.5" (16.5cm) Long T-shape Burnisher With a Ball on One Side Premium Quality Stainless Steel Dental Amalgam Composite Filling Instruments.

Single Ended ball burnisher #34 is often used to smooth the surface of a freshly placed amalgam restoration. T-shaped working end with a small ball on one side of "T" and paddle on the other. The shape of burnisher #34 conforms to the rounded cavity surfaces and allows ease of access into the rounded corners or junctions of the cavity surfaces to condense and shape the composite against the cavity wall. All our amalgam composite filling instruments are manufactured for optimal results and precision. These

dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

**BALL BURNISHER #34 FEATURES**

- Crafted with 6.5" (16.5cm) Overall Length.
- Single Ended.
- Small Pattern.
- T-Shaped Working End.
- Grip Design Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel.
- Reusable.
- Hand Crafted From High-Grade Stainless Steel.

**PACKAGE INCLUDE:**

- 10 Pcs Ball Burnisher #34 Single Ended 6.5" (16.5cm) Long



**10  
PCS****WI-70001442****10 Pcs Dental Amalgam Well Non-Slip Restorative Carrier****\$109.99** per pieceTrade price from **\$101.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Dentist Bone Grafting Amalgam Mixing Well Pot Non-Slip Base Stainless Steel Dentist Lab Restorative Dental Surgical Instruments.

Amalgam well pot used to mix and dispense amalgam during restorative procedures. Dentist use this non-slip small bowl shaped container to holds mixed amalgam prior to its being loaded into the amalgam carrier. amalgam well used to hold amalgam until it is used in the restoration. Dental restorative amalgam well can also be used for bone grinding/crushing in dental implant procedures.

**DENTAL AMALGAM WELL FEATURES:**

**AMALGAM WELL POT:** Amalgam well have a feature is 3cm diameter, Heavy well pot with secure base with rubber friction locking ring. Machined from a solid piece of steel with a cone shaped inner well, Smooth cavity for optimal mixing and flat cone tip for easy amalgam carrier transfer.

**NON-SLIP RUBBER RING BOTTOM:** Non-slip removable rubber 'O' ring and bottom heavily weighted to avoid spillage.

**STAINLESS STEEL CONSTRUCTION:** The stainless steel construction ensures maximum durability and dependability. You can reuse these amalgam well over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

**EASY TO CLEAN:** Amalgam well is a

**WI-70002287****10 Pcs Dental Elevators Set – Periosteal Apical Molt Oral Implant****\$35.99** per pieceTrade price from **\$33.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Dental Elevators Instrument Osteotome Periotome Luxating Elevator Apical Molt Oral Implant Surgical Root Extraction Spreading – 10 Pcs Set


Elevators are used to loosen teeth prior to forceps extraction, to remove roots or impacted teeth, when teeth are compromised and susceptible to fracture or when they are malpositioned and cannot be reached with forceps. These delicate, sharp instruments cut periodontal ligament while preserving bone anatomy with thin, tapered tips to ensure the smooth widening of the alveolar. Fine blade inserts between bone and root with minimal damage to bone structure alternative to standard elevators for less traumatic extractions

**10 PCS DENTAL ELEVATORS SET INCLUDES:**

- 1 Pc Dental Elevator 4R
- 1 Pc Dental Elevator 4L
- 1 Pc Dental Elevator 303
- 1 Pc Dental Elevator 302
- 1 Pc Dental Elevator 301
- 1 Pc Dental Elevator 303HF
- 1 Pc Dental Elevator 34
- 1 Pc Dental Elevator 34S
- 1 Pc Dental Elevator 46
- 1 Pc Dental Elevator A

**DENTAL ELEVATORS FEATURES**

**DELICATE, SHARP INSTRUMENTS:** Cut periodontal ligament while preserving bone anatomy thin, tapered tips ensure smooth widening of the alveolar. Fine blade inserts between bone and root with minimal damage to bone structure





WI-42944644

**10 PCS Dental Orthodontic Instrument Set**

**\$38.00** per piece

Trade price from **\$22.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

This is our pattern of 10 PCS Dental Orthodontic Instrument Set, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. We work the instrument by

hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$38.00, falling to \$22.00 at 1 piece — about 42% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35370681

**10 PCS Dental Spoon Excavator Set Composite Filling**

**\$38.00** per piece

Trade price from **\$25.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

10 PCS Dental Spoon Excavator Set Composite Filling is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$38.00, falling to \$25.00 at 1 set — about 34% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-25126184

10 PCS German Instrument Sets

**\$4.00** per piece

Trade price from **\$3.00**  
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

10 PCS German Instrument Sets is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$3.00 at 10 sets — about 25% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-42929636

10 PCS Hollow Handle Dental Inspection Mirror

**\$38.00** per piece

Trade price from **\$22.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of 10 PCS Hollow Handle Dental Inspection Mirror, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$38.00, falling to \$22.00 at 1 set — about 42% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-42919953

### 10 PCS Replaceable Dental Mouth Mirror Set

**\$38.00** per piece

Trade price from **\$22.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

10 PCS Replaceable Dental Mouth Mirror Set is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$38.00, falling to \$22.00 at 1 set — about 42% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-70002724

### 10 Pcs Wax Carving Spatula Modeling Kit

**\$19.99** per piece

Trade price from **\$18.39**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 PCS Wax Carving Tools Dental Waxing Instruments with Zipper Case – Premium Quality Stainless Steel Wax Carvers Set – Metal Clay Sculpting Pottery Sculpture Tools Double Ended DIY Waxing Modeling Kits Steel Wax Carving Tools Dental Clay Pottery Plaster Sculpture Carvers Tools Kit Dabber Tool Set.

#### PACKAGE INCLUDE:

- 1 Pc Ward 1 + Zahle 3 Tips
- 1 Pc Lecorn 5 Wax Carver
- 1 Pc Ward 1 + Zahle Tips
- 1 Pc Pk Thomas #2
- 1 Pc Hollenback + Ward Tips
- 1 Pc Ward + Zahle Tips
- 1 Pc Excavator + Zahle Tips
- 1 Pc Vehe + Spoon Spatula Tips
- 1 Pc Lecorn 5 + PK Thomas #2 Tips

- 1 Pc Fahnestock Wax Knife with **FREE** Zipper Carrying Case

#### 10 PCS WAX CARVING TOOLS SET FEATURES:

**COMPLETE 10-PIECES KIT:** Our clay/wax sculpting tools set contains the key 10 tools you'll need, like the Probe, Scraper, Spoon, Spatula, Chisel, Spearhead, Ballpoint, Paddle, Carver, and more.

**QUALITY STEEL INSTRUMENTS:** Made of durable and long-lasting hardened stainless steel, these carving tools are heat resistant, non-stick, stain-resistant, and built to withstand your creative output.

**NON-SLIP HANDLE:** Frosted handles provide a good grip and non-slip performance. They are comfortable



**WI-70000073****10 Periosteal Elevators Set****\$39.99** per pieceTrade price from **\$36.79**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Pcs Periosteal Elevators Set – Premium Quality Stainless Steel Dental Surgical Instruments for Implants, Periodontal and Bone Grafting Procedures.

Periosteal elevators are essential tools for periodontists, implantologists, oral surgeons, and general dentists, utilized for elevating the periosteum during surgical procedures. They facilitate the separation of the periosteal layer from the bone, allowing access to underlying structures and improving surgical outcomes. This 10 pieces periosteal elevators set featuring a variety of elevators having double-ended structures, slightly curved, and thin tip to accommodate less traumatic manipulation for different clinical needs. Made with premium grade

stainless steel, these elevators offer durability, rust resistance, and precise handling to enhance the patient's comfort. These are easy cleaning and sterilization, ensuring safety and hygiene for dental practices.

**SET INCLUDES**

- 1 Pc Molt #9 Periosteal Elevator – Double Ended
- 1 Pc #OHL Periosteal Elevator – Double Ended
- 1 Pc #PR3 Prichard Periosteal Elevator – Double Ended
- 1 Pc #S23 Seldin Periosteal Elevator – Double Ended
- 1 Pc 2 West (#W2) Periosteal Elevator – Double Ended
- 1 Pc GF 14 Periosteal Elevator – Double Ended

**WI-70000701****100 Pcs Ball Burnisher #34 Single Ended****\$149.99** per pieceTrade price from **\$137.99**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 Pcs Ball Burnisher Set #34 Single Ended 6.5" (16.5cm) Long T-shape Burnisher With a Ball on One Side Premium Quality Stainless Steel Dental Amalgam Composite Filling Instruments.

Single Ended ball burnisher #34 is often used to smooth the surface of a freshly placed amalgam restoration. T-shaped working end with a small ball on one side of "T" and paddle on the other. The shape of burnisher #34 conforms to the rounded cavity surfaces and allows ease of access into the rounded corners or junctions of the cavity surfaces to condense and shape the composite against the cavity wall. All our amalgam composite filling instruments are manufactured for optimal results and precision. These

dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.


**BALL BURNISHER #34 FEATURES**

- Crafted with 6.5" (16.5cm) Overall Length.
- Single Ended.
- Small Pattern.
- T-Shaped Working End.
- Grip Design Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel.
- Reusable.
- Hand Crafted From High-Grade Stainless Steel.

**PACKAGE INCLUDE:**

- 100 Pcs Ball Burnisher #34 Single Ended 6.5" (16.5cm) Long



**12  
PCS****WI-70003937****12 P.K Thomas #2 Dental Wax Modelling Carver Set****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

PTK2 – 12 Pcs P.K Thomas #2 Dental Wax Modelling Carver with Elongated Curve on Thinner Tips Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT2 features curved on thinner tapered tips used to flow on molten wax. P.K. Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas.

Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture when working on hard-to-reach areas. Made of high



quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

**P.K THOMAS #2 WAX CARVER FEATURES**

- Elongated Curve on Thinner Tips Than P.K. Thomas #1.
- Double-ended
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

**PACKAGE INCLUDE**

- 12 Pcs Thomas Waxing Instrument – PKT2

**12  
PCS****WI-70003940****12 Pcs Amalgam Plugger Black #1/2 Double Ended****\$71.99** per pieceTrade price from **\$66.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Amalgam Plugger Black #1/2 Double Ended 7" (17.8cm) with Serrated Tips and Round Pattern Knurled Grip Handle Premium Quality Stainless Steel Surgical Dental Instruments.

An amalgam condenser is a non-cutting dental instrument used in operative dentistry and endodontics to condense silver amalgam in natural tooth cavities during dental filling procedures. Amalgam Plugger Black 1/2 is used to pack and condense amalgam into cavity preparation. With 1.5 mm – 2.5 mm double-ended cross-serrated tips, the 1/2 black amalgam plugger can also be used with other restorative materials. PLGB 1/2 Black Plugger/Condenser, handcrafted from premium quality stainless steel,



ensures durability and performance longevity. Can be reused after sterilization.

**AMALGAM PLUGGER BLACK #1/2 FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- Double Ended.
- 1.5mm – 2.5mm Cross Serrated Tips.
- Ergonomically Designed Round Pattern Grip Handle.
- Fully Autoclaved / Reusable
- Made with Premium Grade Stainless Steel.

**PACKAGE INCLUDE**

- 12 Pcs Amalgam Plugger Black #1/2 Double Ended

**TECHNICAL INFORMATION**

**WI-70004165****12 Pcs American Pattern Extracting Forceps No. 17 Lower Molars****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Premium Quality American Pattern Dental Extracting Forceps #17 For Lower Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#17 xtraction Forceps is used to extract Lower Molars. The Working end of this tooth extraction instrument is designed with a sharp bend, making sure it reaches deep into the socket for a strong grip. German grade stainless steel promotes longevity and corrosion resistance. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #17 FEATURES**

- Tri-pointed Beak Design

- For Lower Molars
- Non-slip Grip Handle
- German Stainless Steel

**PACKAGE INCLUDE**

- 12 Pcs Dental Extracting Forceps #17 – Lower Molars

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test,**WI-70004129****12 Pcs American Pado Extracting Forceps #16S Lower Molars****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pado Extracting Forceps #16S for Lower Molars Cowhorn Beak Design American Pattern German Grade Stainless Steel Childern Teeth Extraction Surgical Dental Insruments

Cowhorn Pado Extraction Forceps No. 16S are the perfect tool to remove childern teeth from the alveolar bone. American extracting forceps #16S has cowhorn pointed beaks designed for lower first and second molars. Hooked handle provides increased leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you. 100% satisfaction guarantee against defects in material and workmanship.

**PEDO EXTRACTING FORCEPS #16S FEATURES**

- Cow Horn Beak Design.
- For Lower 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel

**PACKAGE INCLUDE**

- 12 Pcs American Pado Extraction Forceps #16S Cowhorn Beak with Hook Handle

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction Surgery Procedure**BEAK/JAWS:** Cowhorn Beaks Design

**12  
PCS****WI-70004278****12 Pcs Ball Burnisher #28/29 Set Double Ended  
2.3mm – 3.6mm Football Tips****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Dental Ball Burnisher #28/29 Double Ended 2.3 mm – 3.6 mm Football Working Tips Premium Quality Stainless Steel Amalgam Restorative Composite/Plastic Filling Instruments.


28/29 football end burnisher is an operative hand instrument for dental professionals used to condense, carve and polish restorative materials such as amalgam. Ball Burnishers can also be used to form occlusal anatomy. The working end of a #28/29 ball burnisher is smooth and rounded like football (2.3 mm – 3.6 mm) made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. Fully autoclavable for reuse purposes.

**BALL BURNISHER # 28/29  
FEATURES**

- Crafted with 6.5" (16.5cm) Overall Length.
- 2.3mm Football & 3.6mm Football.
- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Reusable.
- Smooth Finish For Excellent Tactile Feel

**PACKAGE INCLUDE:**

- 12 Pcs Foot Ball Burnisher #28/29 Double Ended

**12  
PCS****WI-70004104****12 Pcs Bone File Set #12A Double Ended****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Bone File #12A Double Ended 7" (17.8cm) 5mm Cross Serrated and Plain Serrated Edge with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

12A Bone File designed to smooth and shape bone using a pull stroke. Their cross-serrated & serrated working ends make it easy to maneuver between soft tissue and bone. Using a surgical bone file to smooth the bone and narrow the canals of the interior of the tooth helps to keep the area disinfected. Round non-slip knurled handle allow strong grip while handling. Hand crafted from premium quality stainless steel, fully autoclaved and sterilized for reuse purpose.

- Crafted with 7" (17.8cm) Overall Length.
- Double-ended.
- One Cross-serrated & One Plain Serrated Ends.
- Non-slip Knurled Grip Round Handle.
- Reusable.
- Hand Crafted Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 12 Pcs Dental Bone File #12A Double Ended 7" (17.8cm)

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Bone Files**APPLICATION:** Used for Smoothing and Shaping of Bones**CURVATURE:** Curved**RESUABLE:** Yes**12A BONE FILE FEATURES**

**WI-70001466****12 Pcs Castroviejo Caliper Long 0-40mm Angled 6"****\$155.99** per pieceTrade price from **\$143.51**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Castroviejo Implant Bone Caliper 6" (15.2cm) Long Model – 30° Angled Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper Angled 15.2cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 6 inches and tapers from the measurement end to the tip. The shanks are straight with 30° Angled jaws

**LONG CASTROVIEJO CALIPER FEATURES**

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 15.2cm
- 30° Angled Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

**PACKAGE INCLUDE:**

- 12 Pcs Castroviejo Caliper 6" (15.2cm) Long Model – 30° Angled Tip Graduated 0-40mm

**TECHNICAL INFORMATION****WI-70001457****12 Pcs Castroviejo Caliper Long 0-40mm Angled 7"****\$164.99** per pieceTrade price from **\$151.79**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Castroviejo Implant Bone Caliper 7" (17.8cm) Long Model – 30° Angled Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper Angled 17.8cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 7 inches and tapers from the measurement end to the tip. The shanks are straight with angled jaws.

**LONG CASTROVIEJO CALIPER FEATURES**

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 17.8cm
- 30° Angled Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

**PACKAGE INCLUDE:**

- 12 Pcs Castroviejo Caliper 7" (17.8cm) Long Model – 30° Angled Tip Graduated 0-40mm

**TECHNICAL INFORMATION**

**WI-70001461****12 Pcs Castroviejo Caliper Long 0-40mm Straight 7"****\$164.99** per pieceTrade price from **\$151.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Castroviejo Implant Bone Caliper 7" (17.8cm) Long Model – Straight Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper Straight 17.8cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 7 inches and tapers from the measurement end to the tip. The shanks and jaws are straight.

**LONG CASTROVIEJO CALIPER FEATURES**

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 17.8cm
- Straight Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

**PACKAGE INCLUDE:**

- 12 Pcs Castroviejo Caliper 7" (17.8cm) Long Model – Straight Tip Graduated 0-40mm

**TECHNICAL INFORMATION****WI-70003997****12 Pcs Cement Spatulas 7" (17.8cm) Double Ended****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Cement Spatulas Double Ended 7" (17.8cm) Lab Mixing Modeling Premium Quality Stainless Steel Dental Restorative Instruments.


Cement spatulas used for mixing cement materials and various other dental materials like silicones, waxes, cements etc. Dental spatula comes in double ends, one tip is rounded for scooping, mixing and the other is arch headed spreading tip used to mix, spread and lift material including foods, drugs, plaster and paints.

Dental cement spatulas available in double and single ended patterns with flexible and rigid forms which are ideal for lighter, medium to heavier cements and other restorative materials. These spatulas can also be used to load materials into crowns or other

preparations. Made with premium quality stainless steel with superior craftsmanship. Provides a reliable performance for dental professionals over time. Fully Autoclavable for hygiene and reuse purpose.

**DOUBLE ENDED CEMENT SPATULAS FEATURES**

- Crafted With 7" (17.8cm) Overall Length.
- Double Ended.
- One Rounded End and Other is Arch Headed.
- Optimal Size for Handling for Optimal Maneuvering.
- Non-slip Textured Area for Better Grip
- Highly Polished Rounded Handle
- Fully Autoclave and Reusable



**WI-70000731****12 Pcs Cement Spatulas Single Ended****\$64.99** per pieceTrade price from **\$59.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 12 Pcs Dental Cement Spatulas Single Ended 7" (17.8cm) Lab Mixing Modeling Stainless Steel Restorative Instruments.

Cement spatulas are used for mixing cement materials and various other dental materials like silicones, waxes, cements etc. A spatula is a broad, flat, flexible blade used to mix, spread and lift material including foods, drugs, plaster and paints. Available in double and single ended patterns with flexible and rigid forms which are ideal for lighter, medium to heavier cements and other restorative materials. Cement spatulas can also be used to load materials into crowns or other preparations. Made with premium quality stainless steel with superior craftsmanship. Cement spatulas can be

reused over time as they will not bend or tarnish on you. We offer 100% satisfaction guaranteed.

**SINGLE ENDED CEMENT SPATULAS FEATURES**

- Crafted With 7" (17.8cm) Overall Length.
- Single Ended.
- Large Width Spatula Design.
- Optimal Size for Handling for Optimal Maneuvering
- Non-slip Textured Area for Better Grip
- Highly Polished Rounded Handle
- Fully Autoclave and Reusable
- Made With High Quality Stainless Steel

**PACKAGE INCLUDE****WI-70004113****12 Pcs Dental Ball Burnisher #1 Double Ended****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Ball Burnisher #1 Double 1.3mm & 1.8mm Ball Ends 7" (17.8cm) Long Premium Quality Stainless Steel Dental Amalgam Filling Instruments.

Dental 1 Ball Burnisher with 7" overall length used to condense, carve and polish restorative materials such as amalgam. 1 Ball Burnisher composite amalgam filling instrument is non-stick and comes with ball shaped ends. Amalgam Ball Burnisher #1 can be used to blend your composite material to the margin for final contouring. All of our composite placement instruments made with stainless steel ensuring excellent performance instruments which are autoclavable and can reuse after sterilization.


**DENTAL BALL BURNISHER #1**

- Crafted With 7" (17.8cm) Overall Length.
- 1.3mm & 1.8mm Ball Ends.
- Double-ended.
- Non-slip Knurled With Center Round Handle.
- Autoclavable & Reusable.
- Made With Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 12 Pcs Ball Burnisher #1 Double Ended 1.3mm & 1.8mm

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.



**WI-70004103****12 Pcs Dental Bone File #12CA Double Ended****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Dental Bone File #12CA Double Ended 7" (17.8cm) Cross Serrated Working Tips with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Dental Instruments.

12CA Bone File is a cross serrated working ends dental instrument used to excavate small spicules of bones from narrow spaces using a pull stroke. This bone file usually maneuvers between the soft tissues and the bones that need scraping. The knurled design allows strong grip and easy maneuvering. Made with premium quality stainless instrument for durability and long-lasting performance efficiency. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for

a variety of precise neurosurgical procedures.

**BONE FILE #12CA FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- Cross-serrated Ends.
- Double-Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 12 Pcs Bone File #12CA 7" Double Ended

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Surgical Dental Instruments**WI-70004102****12 Pcs Dental Bone File #33 Double Ended****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Dental Bone File #33 Double Ended 7" (17.8cm) Long with 4.8mm Narrow Pattern Straight Plain Serrated Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

#33 Bone File is designed to smooth and shape bone using a pull stroke. Bone File 33 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. The straight serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 4.8mm wide. Constructed from surgical-grade stainless steel for excellent performance and longevity. The round

handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

**BONE FILE #33 FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- 4.8mm Narrow Pattern.
- Double Ended.
- Straight Parallel Sides Plain Serration.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 12 Pcs Bone File #33 Straight 4.8mm Double Double Ends

**12  
PCS****12  
PCS**

**WI-70004171****12 Pcs Dental Extracting Forceps #53L For Upper 1st & 2nd Left Molars****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Dental Extracting Forceps #53L for Upper 1st and 2nd Left Molars American Pattern German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

53L dental forceps are used to extract maxillary 1st and 2nd left molars. Left angled pointed offset (Bayonet) beaks for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #53L FEATURES**

- Left Pointed Tip Design
- For Upper 1st & 2nd Molars
- Non-Slip Textured-Grip Handle
- German Stainless Steel

**PACKAGE INCLUDE**

- 12 Pcs American Pattern Dental Extracting Forceps #53L

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**12 PCS EXTRACTING FORCEPS KIT INCLUDES**

|      |       |      |       |
|------|-------|------|-------|
| #150 | #150S | #151 | #151S |
| #23  | #65   | #88L | #88R  |
| #MD1 | #MD2  | #MD3 | #MD4  |

**WI-70000630****12 Pcs Dental Extracting Forceps Tooth Extraction Instruments Kit****\$109.99** per pieceTrade price from **\$101.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Dental Extracting Forceps Teeth Extraction Kit with FREE Carrying Case Textured-Grip Handle Stainless Steel Dental Surgical Instruments.

12 pcs tooth extraction dental forceps kit assembled to cover most common dental surgery needs. Extracting forceps are used to extract tooth, designed for use in specific areas of the mouth. The beak is shaped to conform snugly to the contour of the tooth which makes it easier to reach different teeth effectively. Our american pattern dental extracting forcep's handle shaped are so well that a maximum amount of force can be applied to the beaks, while the handles are still in a comfortable position for the dental surgeon. Made with premium

grade stainless steel ensures durability and performance longevity. Dental extracting forceps can be fully autoclaved and are reusable.

**12 PCS DENTAL EXTRACTING FORCEPS KIT INCLUDES**

- 1 Pc Extracting Forceps #23
- 1 Pc Extracting Forceps #65
- 1 Pc Extracting Forceps #88L
- 1 Pc Extracting Forceps #88R
- 1 Pc Extracting Root Forceps #150
- 1 Pc Extracting Root Forceps #150S
- 1 Pc Extracting Root Forceps #151
- 1 Pc Extracting Root Forceps #150S
- 1 Pc Extracting Mead Forceps #MD1
- 1 Pc Extracting Mead Forceps #MD2
- 1 Pc Extracting Mead Forceps #MD3
- 1 Pc Extracting Mead Forceps #MD4



**WI-70003939****12 Pcs Dental Hollenback Wax Carver #3****\$65.99** per pieceTrade price from **\$60.71**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Hollenback Carver #3 Double Ended Spear Shape Sharp Blades with 6.5" Overall Length Premium Quality Stainless Steel Wax Modelling Dental Restorative Instruments.

3 hollenback carver is a dental instrument used for placing, carving, and contouring restorative materials like amalgam and composite. Hollenback #3 used for shaping restorative materials to resemble natural tooth anatomy, placing and packing materials into cavities, and trimming excess material after restoration. Double-ended design with large and sharp blades enable detailed work on both the chewing surface and sides of teeth, even in tight spaces. Stainless steel construction of hollenback carver 3 guarantees


durability and dependability, allowing for reuse without bending or tarnishing, and meets international quality standards.

**DENTAL HOLLENBACK CARVER #3 FEATURES**

- Large Pointed Spear Shape Sharp Blades.
- Double Ended.
- Non-Slip Knurled Grip Handle.
- Used for Wax and Porcelain Carving.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 6 Pcs Hollenback #3 Double Ended

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**WI-42941584****12 PCS Dental Mouth Mirror Set****\$9.00** per pieceTrade price from **\$2.50**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**


This is our pattern of 12 PCS Dental Mouth Mirror Set, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 12-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$9.00, falling to \$2.50 at 1 set — about 72% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|                  |                                                                                 |
|------------------|---------------------------------------------------------------------------------|
| <b>Material</b>  | Stainless steel — we state the grade on your quotation                          |
| <b>Reusable</b>  | Yes, autoclavable                                                               |
| <b>Finish</b>    | Satin, mirror or black — your choice at order                                   |
| <b>Marking</b>   | Your logo or part number can be applied — ask us what is possible for this item |
| <b>Lead time</b> | Typically 15–30 days depending on quantity and marking                          |

**WI-70000971****12 Pcs Extracting Forceps Upper & Lower Dental Instruments****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

12 Pcs Dental Extracting Forceps Set German Grade Stainless Steel Upper & Lower American Pattern Tooth Extraction Instruments with Textured Grip Handle.

Dental extracting forceps are used to extract tooth, designed for use in specific areas of the mouth. The beak is shaped to conform snugly to the contour of the tooth. This makes it easier to reach different teeth effectively. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**12 PCS DENTAL EXTRACTING FORCEPS SET INCLUDES**

- 1 Extracting Forceps #23
- 1 Extracting Forceps #65
- 1 Extracting Forceps #88L
- 1 Extracting Forceps #88R
- 1 Extracting Root Forceps #150
- 1 Extracting Root Forceps #151
- 1 Extracting Root Forceps #150S
- 1 Extracting Root Forceps #151S
- 1 Extracting Mead Forceps #MD1
- 1 Extracting Mead Forceps #MD2
- 1 Extracting Mead Forceps #MD3
- 1 Extracting Mead Forceps #MD4

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**WI-70000685****12 Pcs Miller Bone File #45 Double Ended****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Miller Bone File #45 Straight-Cut 7" (17.8cm) Long with 5mm Wide Straight & Curved Double Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

Bone file #45 is designed to smooth and shape bone using a pull stroke. Dental bone file #45 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. It can scrape off bony segments from soft tissues and bones. The serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features straight-cut plain serrations and is 5mm wide. Constructed from surgical-grade

stainless steel for excellent performance and longevity. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

**BONE FILE MILLER #45 FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- 5mm Straight-cut Plain Serrated Ends.
- Curved & Straight Double Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 12 Pcs Miller Bone File #45 Double Ended 5mm Straight & Cruved Ends



**WI-70000537****12 Pcs Molt Bone Curette Set #2/4 Double Ended****\$83.99** per pieceTrade price from **\$77.27**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Molt Bone Curette Set #2/4 Double Ended 3.5mm & 7.0mm Working Tip Round Pattern Handle Premium Quality Stainless Steel Surgical Dental Instruments.

2/4 Molt surgical curette is a double ended bone curette designed for curettage, cyst removal, and tooth socket debridement. With a 3.5 mm and 7.0 mm working ends allows the surgeon to scrape away bone with durable and sharp edges. Molt #2/4 bone curette has round pattern non-slip knurled handle provides optimum control and omits hand fatigue.

This double ended bone curette can also use in a variety of veterinary dental procedures for efficient scraping of diseased tissue from oral/dental wounds for cats and dogs. Molt Bone

curette 2-4, crafted with high-quality surgical stainless steel with a smooth surface, ensures a strong grip. Completely sterilizable and reusable

**MOLT BONE CURETTE #2/4 FEATURES**

- Crafted with 6" (15.2cm) Overall Length.
- Double Ended.
- 3.5 mm – 7.0mm Spoon Sharp Blades.
- Ergonomic Design Round Pattern Grip Handle.
- Sterilizable / Reusable.
- Made with Premium Quality Stainless Steel

**PACKAGE INCLUDE**

- 12 Pcs Molt Bone Curette #2/4 Double Ended

**WI-70003903****12 Pcs Periosteal Elevator Baby Molt #7A Double Ended with Round Handle****\$59.99** per pieceTrade price from **\$55.19**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Periosteal Elevator Baby Molt #7A Double Ended 6.5" (16.5cm) Round Pattern Grip Handle Premium Quality Stainless Steel Dental Instruments.

Molt-7 elevators have thin blades on both ends that are used to cut periodontal ligaments and loosen teeth. Periosteal Baby Molt #7A comprises a flat, oval tip at each end, one blunt and one sharp. Surgeons use the blunt oval tip for lifting and retracting the periosteal flap, and the sharp tip to dissect it. In surgery dentist often needs to separate a bone or tooth from the periosteum that covers it. This is done with a periosteal elevator. The dentist may also use it to gain access to retained roots and surrounding bone. Elevators are manufactured from

stainless steel to ensure optimal performance.

**PERIOSTEAL ELEVATOR BABY MOLT #7A FEATURES**

- Overall Length 6.5" (16.5cm).
- Double Ended.
- One Flat and One Oval Tip.
- Non-Slip Knurled Grip Handle.
- Use to Elevate, Retract and Dissecting Periosteal Flap.
- Corrosion Resistance.
- Fully Autoclavable & Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 12 Pcs Periosteal Elevator Baby Molt #7A Double Ended

**TECHNICAL INFORMATION**

**WI-70000682****12 Pcs Polokoff Bone File #1 Double Ended****\$65.99** per pieceTrade price from **\$60.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Polokoff Bone File #1 Double Ended 7" (17.8cm) Long with 3.4mm Narrow Pattern Straight Serrated Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

#1 Polokoff Bone File is broadly used in surgical procedures to aid in mild bone reshaping using a pull stroke. Polokoff Bone File is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. The straight serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 3.4mm wide. Constructed from surgical-grade stainless steel for excellent

performance and longevity. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

**POLOKOFF BONE FILE #1 FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- 3.4mm Narrow Pattern.
- Double Ended.
- Straight Parallel Sides Serration.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 12 Pcs Polokoff Bone File 1 Straight 3.4mm Narrow Pattern Working

**WI-70003956****12 Pcs Spring Divider Calipers Set 4.5" Dental Lab Tool****\$84.99** per pieceTrade price from **\$78.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Pack of 12 Pcs Spring Divider Calipers Set 4.5" (11.4cm) with Flat Legs Stainless Steel Precision Tool Bone for Accurate Measurements During Dental Lab Procedures.

The dental calipers measures the width of metal, wax, and crown size in dental cosmetic procedures. It has a scale or screw for the measuring of sizes in the dental lab. Spring Divider Caliper features two adjustable legs to measure hard-to-measure dimensions, made from carefully selected flat stock with heat treated points for durability, and they pivot on a roller while the bow spring ensures a precise tension setting. Accurate adjustments are made by turning the knurled nut on the threaded rod. The design of these calipers, with their fine adjustment

feature, ensures quick and positive settings when taking or transferring specific measurements. Manufactured from high quality stainless steel with a polished finish.

**SPRING DIVIDER CALIPER FEATURES**

- Crafted With 4.5" (11.4cm) Overall Length.
- Flat Legs with Sharpened Hard and Pointed Tips.
- Adjustable Lockable Steel Wing Divider.
- Rivet Design to Increase Stability.
- Finely Threaded to Achieve Fine Adjustment.
- Durable and Rustproof.
- Reusable.



**WI-70003983****12 Pcs Star-Like Boone Bracket Positioning Gauge Set****\$72.09** per pieceTrade price from **\$66.32**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Boone Bracket Positioning Gauge with Four Integrated Measurement 3.5, 4.0, 4.5 and 5.0 Millimeter – Premium Quality Stainless Steel Dental Orthodontic Instruments for Orthodontists, Dentists, and Dental Assistants involved in Orthodontic Procedures.

Boone gauge facilitates consistent bracket placement on the incisal or occlusal surfaces and promotes uniform bonding heights across the dental arch, aiding proper tooth movement during treatment. Made with premium quality stainless steel ensures sterility for hygienic orthodontic practices.

**BOONE BRACKET POSITIONING GAUGE FEATURES**

- Four Measurement Options: 3.5, 4.0, 4.5, and 5.0 Millimeters.
- Pointed Tips with Flat Base Portion for Superficial Marking
- Measure Incisal or Occlusal Surface of the Tooth to Determine the Ideal Bracket Position.
- Withstand Repeated Use and Sterilization Procedures.



Boone Bracket Positioning Gauge (Star-Like Gauge) is an orthodontic instrument used for precisely measuring and marking the vertical height for braces placement on teeth. It includes multiple tips sized at 3.5mm, 4.0mm, 4.5mm, and 5.0mm, aligning with standard distances from a tooth's biting edge, ensuring consistent bracket bonding during procedures. It also features pointed tips for accurate surface marking on the tooth crown.

**WI-70002175****12 Pcs T/C Wynman Crown Remover Gripper Forceps 7" Curved Serrated Tips****\$126.99** per pieceTrade price from **\$116.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs T/C Wynman Crown Remover Gripper Forceps Set 7" (17.78cm) Tungsten Carbide Tips with Gold Ring Handle Premium Quality Stainless Steel Temporary Crown Remover Surgical Dental Forceps Instruments.

this instrument offers seamless crown removal solution

**TC WYNMAN CROWN GRIPPER FORCEPS FEATURES**

**12 PCS TC WYNMAN CROWN GRIPPER FORCEPS:** The TC Wynman Crown Gripper Forceps are an ideal addition to your inventory. The tungsten carbide inserts on the Dental Crown Gripper Forceps will grip the temporary crown firmly without slipping, keeping the crown intact. They are useful for gripping wet glass such as in porcelain restorations and assist in more precise procedures.

**O.R. GRADE T/C AND STAINLESS STEEL:** Hot drop-forged construction Wynman crown gripper constructed of O.R. Grade stainless steel and tungsten carbide for the highest degree of



Wynman Crown Remover & Gripper Forceps features a serrated tip made from durable tungsten carbide for a long-lasting working end. Ideal to remove or place temporary crowns with ease, even in the most hard-to-reach places. Our TC Wynman crown grippers are hemostat style forceps that enable a strong grip on the temporary crown without slippage, keeping the crown intact even in the tightest areas of the mouth. Crafted from stainless steel and tungsten carbides,

**WI-70002288****15 Pcs Dental Elevators Instruments Set****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

15 Pcs Dental Elevators Set – Premium Quality Stainless Steel Various Types Dental Periotome Luxating Elevator Apical Molt Oral Implant Surgical Root Extraction Spreading Instrument Kit.

Elevators are used to loosen teeth prior to forceps extraction, to remove roots or impacted teeth, when teeth are compromised and susceptible to fracture or when they are malpositioned and cannot be reached with forceps. These delicate, sharp instruments cut periodontal ligament while preserving bone anatomy with thin, tapered tips to ensure the smooth widening of the alveolar. Fine blade inserts between bone and root with minimal damage to bone structure alternative to standard elevators for less traumatic extractions

**15 PCS DENTAL ELEVATORS SET INCLUDES:**

- 1 Pc Dental Elevator #81
- 1 Pc Dental Elevator #34
- 1 Pc Dental Elevator #46
- 1 Pc Dental Elevator #301
- 1 Pc Dental Elevator #302HF
- 1 Pc Dental Elevator #303HF
- 1 Pc Dental Elevator #302
- 1 Pc Dental Elevator #303
- 1 Pc Dental Elevator #44
- 1 Pc Dental Elevator #45
- 1 Pc Dental Elevator #A
- 1 Pc Dental Elevator #3
- 1 Pc Dental Elevator #8
- 1 Pc Dental Elevator #34S
- 1 Pc Molt Elevator #9 Double Ended

**DENTAL ELEVATORS FEATURES****WI-70003998****2 Cement Spatulas Set 7" (17.8cm) Single & Double Ended****\$16.99** per pieceTrade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Cement Spatulas 7" (17.8cm) Single & Double Ended – Premium Quality Stainless Steel Lab Mixing Wax Modeling Dental Restorative Instruments.

Cement spatulas used for mixing cement materials and various other dental materials like silicones, waxes, cements etc. Dental cement spatulas available in double and single ended patterns with flexible and rigid forms which are ideal for lighter, medium to heavier cements and other restorative materials. These spatulas can also be used to load materials into crowns or other preparations. Made with premium quality stainless steel with superior craftsmanship. Provides a reliable performance for dental professionals

over time. Fully Autoclavable for hygiene and reuse purpose.

**CEMENT SPATULAS FEATURES**

- Crafted With 7" (17.8cm) Overall Length.
- Single and Double Ended Variation.
- Optimal Size for Handling for Optimal Maneuvering.
- Non-slip Textured Area for Better Grip
- Highly Polished Rounded Handle
- Fully Autoclave and Reusable
- Made With High Quality Stainless Steel

**PACKAGE INCLUDE**

- 1 Cement Spatula 7" (17.8cm) Single Ended
- 1 Cement Spatula 7" (17.8cm) Double Ended



**WI-70001233****2 Extracting Forceps Set #53L & #53R Upper Left & Right Molars****\$18.99** per pieceTrade price from **\$17.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Pcs Extracting Forceps Set #53L & #53R for Upper Left & Right Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.



53L & 53R dental forceps are used to extract maxillary 1st and 2nd molars. With one Left and one Right angled pointed offset (Bayonet) beaks designed for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #53L & #53R FEATURES**

- Left & Right Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Angular Beaks.
- Non-Slip Textured-Grip Handle.
- German Stainless Steel.

**PACKAGE INCLUDE**

- 1 Pc American Pattern Dental Extracting Forceps #53L – Upper Left Molars
- 1 Pc American Pattern Dental Extracting Forceps #53R – Upper Right Molars

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**WI-70000973****2 Extracting Forceps Set #88R & #88L for Upper 1st & 2nd Molars****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Extracting Forceps #88R & #88L with Right and Left Pointed, Split Beak Design with Grip Handle for 1st & 2nd Upper Molars German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American Pattern dental extraction forceps No. 88R and 88L used to extract 1st and 2nd molars. This extraction forceps set has pointed split beak design. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**EXTRACTING FORCEPS #88R & #88L FEATURES**

- For Upper 1st & 2nd Molars.
- Right and Left Pointed, Split Beak Design.
- Non-Slip Grip Handle.
- American Pattern.
- Made With Stainless Steel

**PACKAGE INCLUDE:**

- 1 Pc American Pattern Extracting Forceps #88R
- 1 Pc American Pattern Extracting Forceps #88L

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Extraction Forceps**BEAK/JAWS:** Pointed Split Design**POLISHED:** Yes**HANDLE:** Non-Slip Handle

**WI-70001449****2 Iwanson Calipers Set 4" Sharp and Blunt Tips Improved Model****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack 2 Crown Gauge Iwanson Spring Calipers Set 4" (11.4cm) Blunt & Sharp Tips for Wax or Metal Improved Model Stainless Steel Orthodontics Dental Lab Instrument.

Iwanson Caliper is a dental instrument used to measure the thickness of crowns and spaces. This orthodontic lab instrument is also known as Iwanson Gauge & Iwanson Spring Caliper. Iwanson Spring Caliper features Blunt and Sharp tips mainly used for precision measuring of wax crowns and models. High quality medical standard stainless steel product with accurate readings 10 cm length. Calibrated scale capacity 0.1mm to 10.0mm. Marked In 0.1mm increments. Can be sterilized for multiple use.

**IWANSON CALIPER FEATURES**

**CROWN GAUGE CALIPER:** Iwanson Caliper mainly for precision measuring dental gauge. Iwanson Caliper has an improved caliper dial with blunt and sharp tips variations. Spring action enable smooth movement, and a serrated side for better grip. This caliper is a mandatory tool for jeweler's art and craft workers and dental labs

**HIGH-QUALITY STAINLESS STEEL:** Made with superior quality durable stainless steel. The Premium Quality Stainless Steel composition of this caliper offers a long life of strong performance, as they are built to last.

**FUNCTIONAL ERGONOMIC DESIGN:** Features an ergonomic design to complete the task at hand and make your life easy. Provides a comfortable

**WI-70004262****2 Pcs Amalgam Carriers Double Ended + Amalgam Well Set****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Dental Amalgam Filling Set Including Amalgam Carriers Double Ended Large 2.8mm/ Regular 2mm + Amalgam Well Pot Premium Quality Stainless Steel Restorative Dentistry Instruments.

Double-ended amalgam carriers and an amalgam well set designed to facilitates the transport and packing of dental amalgam into tooth cavities during restorative procedures. Amalgam well is a stable, bowl-shaped metal container that holds freshly mixed amalgam while the double-ended amalgam carrier, has two tubular tips which allows the dentist to scoop the amalgam from the well and dispense it into the cavity by pressing a lever. This amalgam filling set handcrafted with premium quality stainless steel ensures

durability and performance in dental settings.

**AMALGAM FILLING PACKAGE INCLUDE**

- 1 Amalgam Carrier Double Ended Large 2.8mm / Regular 2mm Tips
- 1 Amalgam Well with Non-Slip Base

**TECHNICAL INFORMATION**

**MATERIAL:** Stainless Steel

**INSTRUMENTS TYPE:** Dental Restorative Instruments

**APPLICATION:** Used to Fill Dental Cavities With Amalgam

**REUSABLE:** Yes

**STERILE:** Non-Sterile

**RUST PREVENTION PROCEDURE:** Passivated

**ULTRASONIC CLEANED:** Yes



**WI-70004261****2 Pcs Amalgam Carriers Single Ended + Amalgam Well Set****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Dental Amalgam Filling Set Including Amalgam Carriers Single Ended Reulgar 2mm + Amalgam Well Pot Premium Quality Stainless Steel Restorative Dentistry Instruments.

A pack consisting of two single-ended amalgam carriers and an amalgam well is utilized for the efficient handling of dental amalgam. The amalgam carrier functions similarly to a small syringe, allowing the dentist to scoop and dispense the filling material into the tooth cavity. Meanwhile, the amalgam well is a heavy metal bowl that securely holds the mixed amalgam, facilitating easy reloading during the procedure. Crafted from premium stainless steel, this amalgam filler set assures durability and optimal performance in dental settings.

**AMALGAM FILLING PACKAGE INCLUDE**

- 1 Amalgam Carrier Single Ended Regular 2mm Tip
- 1 Amalgam Well with Non-Slip Base

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Restorative Instruments**APPLICATION:** Used to Fill Dental Cavities With Amalgam**REUSABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**

Passivated

**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing**WI-70004260****2 Pcs Ball Burnisher #28/29 Set Double Ended 2.3mm – 3.6mm Football Tips****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Pcs Dental Ball Burnisher #28/29 Double Ended 2.3 mm – 3.6 mm Football Working Tips Premium Quality Stainless Steel Amalgam Restorative Composite/Plastic Filling Instruments.

28/29 football end burnisher is an operative hand instrument for dental professionals used to condense, carve and polish restorative materials such as amalgam. Ball Burnishers can also be used to form occlusal anatomy. The working end of a #28/29 ball burnisher is smooth and rounded like football (2.3 mm – 3.6 mm) made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. Fully autoclavable for reuse purposes.


**BALL BURNISHER # 28/29 FEATURES**

- Crafted with 6.5" (16.5cm) Overall Length.
- 2.3mm Football & 3.6mm Football.
- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Reusable.
- Smooth Finish For Excellent Tactile Feel

**PACKAGE INCLUDE:**

- 2 Pcs Foot Ball Burnisher #28/29 Double Ended



**2  
PCS****WI-70002028****2 Pcs Christensen Crown Remover Set CRCH2  
90 Degree Angled****\$17.99** per pieceTrade price from **\$16.55**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Pcs Christensen Crown Remover CRCH2 90 Degree Angled with Hexagonal Dental Elevator Handle Premium Quality Stainless Steel Permanent Crown Removal Restoration Dental Instruments

For permanent removal of crowns by breaking the seal between tooth and crown after sectioning with a bur. Crown Remover Christensen is ideally constructed for the disassembly method of crown removal. After sectioning the crown with a bur, the elevator handle provides a comfortable grip and maximum leverage. The notched tip then aids the single-end Christensen Crown Remover in effortless removal. Pressure against the tooth is lessened, which reduces the potential for tooth fracture.

**CHRISTENSEN CROWN  
REMOVER FEATURES**

- 90-degree angled CRCH2 Crown Remover is ideally suited for posterior teeth.
- Elevator handle for optimal comfort and grip.
- Notched tips.
- Reusable.
- Made with Highest quality medical grade Stainless Steel.

**PACKAGE INCLUDE:**

- 2 Pcs Christensen Crown Remover CRCH2 90 Degree Angled with Hexagon Handle

**TECHNICAL INFORMATION:****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Restorative Dental Instruments

CHRISTENSEN CROWN REMOVER

6"

**WI-70000001****2 Pcs Cryer Root Elevator Set #25 & #26****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 2 Cryer Root Elevator #25 (Left Angled) and #26 (Right Angled) Medium Triangled Tip with Standard Elevator Handle. Premium Quality Stainless Steel Tooth Extraction Surgical Dental Instruments

Set of 2 cryer elevators (#25 and #26) are specialized dental surgical instruments designed for luxating and loosening teeth, as well as removing stubborn root tips. They effectively detach teeth from the periodontal ligament before extraction. #25 and #26 sharp pointed tip configuration is particularly useful for fractured or retained roots which allow clean extraction of roots by engaging and rotating them. Made from surgical-grade stainless steel, they are durable,

corrosion-resistant, and fully autoclavable for sterilization.

**CRYER ELEVATORS SET 25 & 26  
FEATURES**

- 6 Inch (15.2cm) Overall Length.
- Left and Right Angled Medium Triangled Tip for Root Fragment Removal.
- Ergonomic Handle Provides Improved Grip and Control for the Dentist.
- Designed for the Extraction of Root Segments During Dental Procedures.
- Corrosion Resistance and Anti-oxidation.
- Fully Autoclavable for Reuse.
- Made With German Grade Stainless Steel

**PACKAGE INCLUDE**

**WI-70000002****2 Pcs Cryer Root Elevator Set #44 & #45****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Cryer Root Elevator #44 (Left Angled) and #45 (Right Angled) Small Triangle Tip with Standard Elevator Handle. Premium Quality Stainless Steel Tooth Extraction Surgical Dental Instruments

Set 2 Cryer elevators (#44 and #45) consists of specialized dental surgical instruments designed to luxate and extract broken or retained tooth roots. The set features a left-angled (#44) and a right-angled (#45) pointed tip, enabling access to both mesial and distal sides of a tooth socket. These instruments are particularly beneficial for lower molar extractions, where their sharp tips can engage remaining root fragments or bone using a rotational force. Made from surgical-grade stainless steel, they are durable,

corrosion-resistant, and fully autoclavable for sterilization.

**44 & 45 SMALL CRYER ELEVATORS FEATURES**

- 6 Inch (15.2cm) Overall Length.
- Left and Right Angled Small Triangle Tip for Root Fragment Removal.
- Ergonomic Handle Provides Improved Grip and Control for the Dentist.
- Designed for the Extraction of Root Segments During Dental Procedures.
- Corrosion Resistance and Anti-oxidation.
- Fully Autoclavable for Reuse.
- Made With German Grade Stainless Steel

**PACKAGE INCLUDE****WI-70001444****2 Pcs Dental Amalgam Well Non-Slip Restorative Carrier****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Dentist Bone Grafting Amalgam Mixing Well Pot Non-Slip Base Stainless Steel Dentist Lab Restorative Dental Surgical Instruments.

Amalgam well pot used to mix and dispense amalgam during restorative procedures. Dentist use this non-slip small bowl shaped container to holds mixed amalgam prior to its being loaded into the amalgam carrier. amalgam well used to hold amalgam until it is used in the restoration. Dental restorative amalgam well can also be used for bone grinding/crushing in dental implant procedures.

**DENTAL AMALGAM WELL FEATURES:**

**AMALGAM WELL POT:** Amalgam well have a feature is 3cm diameter, Heavy well pot with secure base with rubber friction locking ring. Machined from a solid piece of steel with a cone shaped inner well, Smooth cavity for optimal mixing and flat cone tip for easy amalgam carrier transfer.

**NON-SLIP RUBBER RING BOTTOM:**

Non-slip removable rubber 'O' ring and bottom heavily weighted to avoid spillage.

**STAINLESS STEEL CONSTRUCTION:**

The stainless steel construction ensures maximum durability and dependability. You can reuse these amalgam well over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

**EASY TO CLEAN:** Amalgam well is a



**WI-70003969****2 Pcs Dental Boley Gauge & Iwanson Gauge Set  
Vernier Lab Instrument****\$17.99** per pieceTrade price from **\$16.55**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Set of 2 Boley Gauge Vernier & Iwanson Wax Caliper Orthodontic Dental Lab Measuring Gauges Stainless Steel Surgical Dentist Instruments.

These dental device used to measure teeth, spaces and crowns in all types of diagnostic, restorative and orthodontic procedures. This Dental Orthodontic Gauge Set comes with a highly polished finish for aesthetics, corrosion resistance, and anti-oxidation. With ultra precision and effortless operation measuring, this dental gauge caliper set is easy to use and a valuable choice for any dental setting. All instruments are fully reusable and can be autoclaved. Made with german grade stainless steel built to last long.

**PACKAGE INCLUDE:**

- 1 Boley Gauge Sliding Vernier – Flat Ends 10mm Measuring with Overall Length 6" (15.24cm)
- 1 Iwanson For Mattel & Wax Crown Gauge, Improved Model with 2 Sided Scale – Measuring Scale: 0.1mm – 10mm

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Instruments**APPLICATION:** Dental Dimensions Measurement**SCALE:** Millimeter Measuring Scale**ADJUSTABLE:** Yes**REUSABLE:** Yes**AUTOC LAVED:** Yes**RUST PREVENTION PROCEDURE:**

Passivated

**ULTRASONIC CLEANED:** Yes**WI-70004070****2 Pcs Dental Bone File #12CA Double Ended****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Dental Bone File #12CA Double Ended 7" (17.8cm) Cross Serrated Working Tips with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Dental Instruments.

12CA Bone File is a cross serrated working ends dental instrument used to excavate small spicules of bones from narrow spaces using a pull stroke. This bone file usually maneuvers between the soft tissues and the bones that need scraping. The knurled design allows strong grip and easy maneuvering. Made with premium quality stainless instrument for durability and long-lasting performance efficiency. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for

a variety of precise neurosurgical procedures.

**BONE FILE #12CA FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- Cross-serrated Ends.
- Double-Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 2 Pcs Bone File #12CA 7" Double Ended

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Surgical Dental Instruments

**WI-70004075****2 Pcs Dental Bone File #33 Double Ended****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Pcs Dental Bone File #33 Double Ended 7" (17.8cm) Long with 4.8mm Narrow Pattern Straight Plain Serrated Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

#33 Bone File is designed to smooth and shape bone using a pull stroke. Bone File 33 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. The straight serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 4.8mm wide. Constructed from surgical-grade stainless steel for excellent performance and longevity. The round


handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

**BONE FILE #33 FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- 4.8mm Narrow Pattern.
- Double Ended.
- Straight Parallel Sides Plain Serration.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 2 Pcs Bone File #33 Straight 4.8mm Double Double Ends

**WI-70004071****2 Pcs Dental Bone Files #12A Double Ended****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 2 Pcs Dental Bone Files #12A Double Ended 7" (17.8cm) 5mm Cross Serrated and Plain Serrated Edge with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

12A Bone File designed to smooth and shape bone using a pull stroke. Dental Bone Files #12A has cross-serrated & serrated working ends make it easy to maneuver between soft tissue and bone. Using a surgical bone file to smooth the bone and narrow the canals of the interior of the tooth helps to keep the area disinfected. Round non-slip knurled handle allow strong grip while handling. Hand crafted from premium quality stainless steel, fully autoclaved and sterilized for reuse purpose.

- Crafted with 7" (17.8cm) Overall Length.
- Double-ended.
- One Cross-serrated & One Plain Serrated Ends.
- Non-slip Knurled Grip Round Handle.
- Reusable.
- Hand Crafted Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 2 Pcs Dental Bone File #12A Double Ended 7" (17.8cm)

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Bone Files**APPLICATION:** Used for Smoothing and Shaping of Bones**CURVATURE:** Curved**RESUABLE:** Yes

**WI-70002143****2 Pcs Seldin Elevator 4L – 10mm Tip****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 2 Pcs 4L Seldin Root Elevators 6" (15.2cm) Left Angled 10mm Tip with Standard Handle Premium Quality Stainless Steel Dental Extraction Surgery Instruments.

Seldin Elevator 4L is used for pulling up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction. 4L elevator is a smaller version of 1L. Seldin Elevators used to remove soft tissue and periosteal membrane from the alveolar bone with pushing strokes. It has one working end, a robust shank, and a sturdy handle. Sharp tips helping in extraction procedure. These can be used for both the upper and lower molar roots.

- Crafted with 6" (15.2cm) Overall Length.
- 10mm Sharp Left Angled Tip
- Left Upper & Lower Molar Roots
- Smaller Version of 1L.
- Ergonomic Standard Handle With Thumb Hold for Increased Comfort.
- Fully Autoclaved & Reusable.
- Superior Quality Stainless Steel Construction.

**PACKAGE INCLUDE:**

- 2 Pcs 4L Seldin Elevators 10mm Tip Standard Handle

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Root Elevators**APPLICATION:** Dental Extraction**BLADE TIP SIZE:** 10mm**SELDIN ELEVATOR 4L FEATURES****WI-70002284****2 Periotome Luxating Dental Elevators Set Del #8 and #A****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Pcs Dental Luxating Elevators Del #8 and #A with Standard Elevator Handle. Premium Quality Dental Elevators Instruments Osteotome Periotome Luxating Elevator Apical Molt Oral Implant Surgical Root Extraction & Spreading.

Elevators are used to loosen teeth prior to forceps extraction, to remove roots or impacted teeth, when teeth are compromised and susceptible to fracture or when they are malpositioned and cannot be reached with forceps. These delicate, sharp instruments cut periodontal ligament while preserving bone anatomy with thin, tapered tips to ensure the smooth widening of the alveolar. Fine blade inserts between bone and root with minimal damage to bone structure

alternative to standard elevators for less traumatic extractions


**2 PCS DENTAL ELEVATORS SET INCLUDES:**

- 1 Pc Dental Elevator DEL #8
- 1 Pc Dental Elevator #A

**DENTAL ELEVATORS FEATURES**

**DELICATE, SHARP INSTRUMENTS:** Cut periodontal ligament while preserving bone anatomy thin, tapered tips ensure smooth widening of the alveolar. Fine blade inserts between bone and root with minimal damage to bone structure alternative to standard elevators for less traumatic extractions.

**USE OF DENTAL ELEVATORS:** Elevators used to loosen teeth prior to forcep extraction, to remove roots or impacted teeth, when teeth are compromised and



**WI-70001234****2 T/C Extracting Forceps Set #53L & #53R  
Upper Left & Right Molars****\$24.99** per pieceTrade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Pcs Tungsten Carbide (T/C) Extracting Forceps Set #53L & #53R for Upper Left & Right Molars American Pattern German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Tungsten Carbide 53L & 53R dental forceps are used to extract maxillary 1st and 2nd molars. With one Left and one Right angled pointed offset (Bayonet) beaks designed for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

**T/C EXTRACTING FORCEPS #53L & #53R FEATURES**

- Left & Right Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Tungsten Carbide Inserts Angular Beaks.
- Non-Slip Textured-Grip Handle.
- German Stainless Steel.

**PACKAGE INCLUDE**

- 1 Pc T/C American Pattern Dental Extracting Forceps #53L – Upper Left Molars
- 1 Pc T/C American Pattern Dental Extracting Forceps #53R – Upper Right Molars

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental**WI-70004276****20 Pcs Amalgam Restorative Set with  
Sterilization Cassette****\$164.99** per pieceTrade price from **\$151.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 20 Pieces Amalgam Restorative Set with Sterilization Cassette – Premium Quality Dental Amalgam Composite Plastic Filling Instruments.

The 20-piece Amalgam Restorative Set with Sterilization Cassette is designed for placing and contouring restorative materials in dental procedures, blending composite and amalgam materials for final contouring. As a direct manufacturer of surgical and dental instruments, we have created a precise Amalgam Restorative Set of 20 pieces using high-grade Stainless Steel for durability and efficiency. The Sterilization Cassette is utilized to efficiently sterilize instruments, ensuring aseptic conditions during restorative procedures and allowing easy access during dental procedures.

All instruments are rust-resistant and autoclave-compatible.

**AMALGAM RESTORATIVE SET INCLUDED**

- 1 Pc Mirror Handle
- 1 Pc Mirror # 5, Front Surface
- 1 Pc Explorer # 5, Double Ended
- 1 Pc College Plier 6" Serrated
- 1 Pc Periodontal Probe X2, 3, 6, 9, 12 mm Markings
- 1 Pc Burnisher # 34, Single Ended
- 1 Pc Burnisher # 27/29, Double Ended
- 1 Pc Burnisher # 28/29, Double Ended
- 1 Pc Acorn Burnisher # 21, Double Ended
- 1 Pc Hollenback Carver # 3, Double Ended



WI-13357194

## 2023 Dental Hygiene Kit Oral Care Kit Pakistan Made

**\$7.00** per piece

Trade price from **\$3.50**

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

2023 Dental Hygiene Kit Oral Care Kit Pakistan Made is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$7.00, falling to \$3.50 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

3 Pieces



WI-70004094

## 3 Bone File Kits of 5 Pcs Dental Surgical Orthopedic Veterinary Surgery

**\$109.99** per piece

Trade price from **\$101.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

3 Set of 5 Pcs Bone Files Kit with Carrying Case Including #64, #45, #12, #1PO and #33 – Premium Quality Stainless Steel 8 Dental Bone File Double Ended for Precise and Delicate Dental Surgical Clinical Procedures.

5 piece dental bone file set comprises surgical-grade stainless steel instruments used in oral, periodontal, and orthopedic surgeries. These bone files designed to smooth and shape alveolar bone, primarily used for post-extraction smoothing, alveoloplasty and flap surgery. Features a double-ended configuration for access to hard-to-reach areas and serrated cutting edges designed for efficient bone tissue removal with minimal trauma. The instruments come in varying sizes to accommodate different surgical needs.

Made from surgical-grade stainless steel ensures outstanding durability and performance over time. These dental surgical bone files can be reused after proper sterilization.

### DENTAL BONE FILES FEATURES

- Double Ended.
- Serrated Tips.
- Round Non-Slip Handles.
- Available in Different Shapes and Diameters.
- Ideal to Smooth and Shape Rough Edges Bone Spicules During Surgical Procedure.
- Fully Autoclavable/Reusable.
- Hand Crafted from Premium Quality Stainless Steel.

### PACKAGE INCLUDE

- 3 Pcs Bone File #64 – Double Ended

**3  
Pieces****WI-70004093****3 Bone File Kits of 8 Pcs Surgical Dental Orthopedic Veterinary****\$179.99** per pieceTrade price from **\$165.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

3 Bone File Kits of 8 Pcs with Carrying Case Including #10, #33, #45, #64, #12, #12A, #12CA and #1PO – Premium Quality Stainless Steel 8 Dental Bone File Double Ended for Precise and Delicate Dental Clinical Procedures.

8-piece dental bone file set comprises surgical-grade stainless steel instruments used in oral, periodontal, and orthopedic surgeries. These bone files designed to smooth and shape alveolar bone, primarily used for post-extraction smoothing, alveoloplasty and flap surgery. Features a double-ended configuration for access to hard-to-reach areas and serrated cutting edges designed for efficient bone tissue removal with minimal trauma. The instruments come in varying sizes to accommodate different surgical needs.

Made from surgical-grade stainless steel ensures outstanding durability and performance over time. These dental surgical bone files can be reused after proper sterilization.

**DENTAL BONE FILES FEATURES**

- Double Ended.
- Serrated Tips.
- Round Non-Slip Handles.
- Available in Different Shapes and Diameters.
- Ideal to Smooth and Shape Rough Edges Bone Spicules During Surgical Procedure.
- Reusable.
- Hand Crafted from Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

3 Sets of 8 Pcs Bone File Kits Includes;

**WI-70003888****3 Heidbrink Root Tip Elevators Set Straight, Left, Right Angled Tips****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Heidbrink Root Tip Elevators Set 1, 2, 3 Straight, Left and Right Angled Slightly Curved Pointed Tip with Hexagonal Handle Premium Quality Stainless Steel Root Tips Tooth Extraction Dental Instruments.


Heidbrink elevators are a single-ended instruments designed for the effective removal of root tips and residual roots in surgical procedures. This delicate dental elevators are used to retrieve loose root tips that may break away from a tooth during extraction and are in hard-to-reach places. Heidbrink also referred to as a root tip pick elevator with a narrow sharp slightly curved design and lightweight hexagonal handle for better control during use. Heidbrink Root Tip Elevator is crafted with straight, left and right angled tips.

This heidbrink root elevators with its unique design is worthy of being a part of your dental kit. Made with premium grade stainless steel, these elevators are lightweight, rust-free, and reusable.

**HEIDBRINK ROOT TIP ELEVATOR FEATURES**

- Crafted with 6" (15.2cm) Overall Length.
- 2.5mm Working End.
- Narrow, Sharp and Slightly Curved Tips.
- Available in Straight, Left and Right Angled Profiles.
- Single Ended Instruments Used to Extract Root Tips.
- Hexagonal Handle Creates Less Force on Tooth With Better Gripping.



**3  
PCS****WI-70003930****3 P.K Thomas #2 Dental Wax Modelling Carver Set****\$20.99** per pieceTrade price from **\$19.31**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

PTK2 – 3 Pcs P.K Thomas #2 Dental Wax Modelling Carver with Elongated Curve on Thinner Tips Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PPKT2 features curved on thinner tapered tips used to flow on molten wax. P.K. Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas.

Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture when working on hard-to-reach areas. Made of high



quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

**P.K THOMAS #2 FEATURES**

- Elongated Curve on Thinner Tips Than P.K. Thomas #1.
- Double-ended
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

**PACKAGE INCLUDE**

- 3 Pcs Thomas Waxing Instrument – PTK2

**3  
PCS****WI-70004256****3 Pcs Amalgam Carrier Double Ended Set Large/Medium 2.5mm/2mm****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 3 Pcs Amalgam Carrier Double Ended Set 8" (20.3cm) Large/Medium 2.5mm – 2mm Tapered Tips with Non-Slip Textured Grip Handle Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument used to fill dental cavities with amalgam. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam

carrier. The textured surface helps effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

**DOUBLE ENDED AMALGAM CARRIER FEATURES**

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Color Bands for Easy Instruments Identification.
- Large – Regular (2.5mm – 2mm) Tapered Tips for Even Distribution of Filling Material.

**100% HIGH QUALITY  
STAINLESS STEEL**

**3  
PCS****WI-70004275****3 Pcs Amalgam Carrier Single Ended Set  
1.5mm Mini Tip****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Amalgam Carrier Single Ended 1.5mm Mini Tapered Tip 6" (15.2cm) Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument for precise and controlled amalgam filling application, available in single-ended or double-ended designs. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. The textured surface helps effectively grip the filling material, preventing it from slipping.

Amalgam filling carriers like a syringes, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental

cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

**SINGLE ENDED AMALGAM  
CARRIER FEATURES**

- 6" (15.2cm) Overall Length.
- Single Ended.
- 1.5mm Mini Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.

**100% HIGH QUALITY  
STAINLESS STEEL****3  
PCS****WI-70000527****3 Pcs Amalgam Carrier Single Ended Set  
2.0mm Regular Tapered Tip****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Amalgam Carrier Single Ended Set 2.0mm Regular Tapered Tip 6" (15.2cm) Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument for precise and controlled amalgam filling application, available in single-ended or double-ended designs. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. The textured surface helps effectively grip the filling material, preventing it from slipping.

Amalgam filling carriers like a syringes, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental

cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

**SINGLE ENDED AMALGAM  
CARRIER FEATURES**

- 6" (15.2cm) Overall Length.
- Single Ended.
- 2.0mm Regular Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.

**100% HIGH QUALITY  
STAINLESS STEEL**

**3  
PCS****WI-70000524****3 Pcs Amalgam Carrier Single Ended Set  
3.0mm Large Tapered Tip****\$23.99** per pieceTrade price from **\$22.07**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Amalgam Carrier Single Ended Set 3.0mm Large Tapered Tip 6" (15.2cm) Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument for precise and controlled amalgam filling application, available in single-ended or double-ended designs. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. The textured surface helps effectively grip the filling material, preventing it from slipping.

Amalgam filling carriers like a syringes, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental

cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

**SINGLE ENDED AMALGAM  
CARRIER FEATURES**

- 6" (15.2cm) Overall Length.
- Single Ended.
- 3.0mm Large Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.

**100% HIGH QUALITY  
STAINLESS STEEL****3  
PCS****WI-70004216****3 Pcs American Extracting Forceps #90 Upper  
Left Molars****\$29.99** per pieceTrade price from **\$27.59**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Dental Extracting Forceps #90 for Upper Molars with Left Pointed Tip Design German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American pattern dental extraction forceps No. 90 used to extract upper left molars. This extraction forceps has left pointed tip design. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**EXTRACTING FORCEPS #90  
FEATURES**

- For Upper Left Molars.

- Left Pointed Tip Design.
- Non-Slip Grip Handle.
- American Pattern.
- Made With Stainless Steel

**PACKAGE INCLUDE:**

- 3 Pcs American Pattern Extracting Forceps #90 – Upper Left Molars

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental  
Extracting Forceps**APPLICATION:** Extraction Forceps**BEAK/JAWS:** Left Pointed Tip**POLISHED:** Yes**HANDLE:** Non-Slip Handle**REUSABLE:** Yes**AUTOCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated

**WI-70004224****3 Pcs American Pedo Extracting Forceps #16S Lower Molars****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pedo Extracting Forceps #16S for Lower Molars Cowhorn Beak Design American Pattern German Grade Stainless Steel Children Teeth Extraction Surgical Dental Instruments

Cowhorn Pedo Extraction Forceps No. 16S are the perfect tool to remove children teeth from the alveolar bone. American extracting forceps #16S has cowhorn pointed beaks designed for lower first and second molars. Hooked handle provides increased leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you. 100% satisfaction guarantee against defects in material and workmanship.

**PEDO EXTRACTING FORCEPS #16S FEATURES**

- Cow Horn Beak Design.
- For Lower 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel

**PACKAGE INCLUDE**

- 3 Pcs American Pedo Extraction Forceps #16S Cowhorn Beak with Hook Handle

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

**APPLICATION:** Teeth Extraction

Surgery Procedure

**BEAK/JAWS:** Cowhorn Beaks Design**WI-70004272****3 Pcs Ball Burnisher #28/29 Set Double Ended 2.3mm – 3.6mm Football Tips****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Dental Ball Burnisher #28/29 Double Ended 2.3 mm – 3.6 mm Football Working Tips Premium Quality Stainless Steel Amalgam Restorative Composite/Plastic Filling Instruments.

28/29 football end burnisher is an operative hand instrument for dental professionals used to condense, carve and polish restorative materials such as amalgam. Ball Burnishers can also be used to form occlusal anatomy. The working end of a #28/29 ball burnisher is smooth and rounded like football (2.3 mm – 3.6 mm) made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. Fully autoclavable for reuse purposes.

**BALL BURNISHER # 28/29 FEATURES**

- Crafted with 6.5" (16.5cm) Overall Length.
- 2.3mm Football & 3.6mm Football.
- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Reusable.
- Smooth Finish For Excellent Tactile Feel

**PACKAGE INCLUDE:**

- 3 Pcs Foot Ball Burnisher #28/29 Double Ended



**WI-70003980****3 Pcs Boley Gauge 10cm Measuring 6" (15.24cm)****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Boley Gauge 10cm Measuring 6" (15.24cm) Long Orthodontic Crown Measuring Teeth Size Spaces Vernier Caliper Dental Implant Lab Orthodontics Surgical Instruments

Boley Gauge is one of the most commonly used measuring gauges in dentistry. The Boley Gauge is specialized dental Instruments helps in determining tooth dimensions of length, width, and thickness in millimeter increments. It is used for measuring teeth and spaces. The flat end is good for measuring teeth size and distance of spaces in the mouth.

This dental gauge featuring smooth operation and easy control to enhance user experience helps in a variety of


clinic, laboratory procedures and dental practices

**BOLEY GAUGE FEATURES**

- Millimeter Increments Measuring Scale Range 100mm
- Small Knob for Adjustment
- Flat End Facilitates Teeth Measurement
- A Boley Gauge Has a Fixed Scale and a Sliding Scale.
- Micro Measuring for More Accuracy
- Ideal For Space Analysis & Band Selection
- Superior Quality Stainless Steel Construction

**PACKAGE INCLUDE**

- 3 Pcs Boley Gauge 10cm Measuring with 6" (15.24cm) Overall Length

**WI-70004097****3 Pcs Bone File Set #12A Double Ended****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Bone File #12A Double Ended 7" (17.8cm) 5mm Cross Serrated and Plain Serrated Edge with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

12A Bone File designed to smooth and shape bone using a pull stroke. Their cross-serrated & serrated working ends make it easy to maneuver between soft tissue and bone. Using a surgical bone file to smooth the bone and narrow the canals of the interior of the tooth helps to keep the area disinfected. Round non-slip knurled handle allow strong grip while handling. Hand crafted from premium quality stainless steel, fully autoclaved and sterilized for reuse purpose.

- Crafted with 7" (17.8cm) Overall Length.
- Double-ended.
- One Cross-serrated & One Plain Serrated Ends.
- Non-slip Knurled Grip Round Handle.
- Reusable.
- Hand Crafted Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 3 Pcs Dental Bone File #12A Double Ended 7" (17.8cm)

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Bone Files**APPLICATION:** Used for Smoothing and Shaping of Bones**CURVATURE:** Curved**RESUABLE:** Yes**12A BONE FILE FEATURES**

**WI-70001468****3 Pcs Castroviejo Caliper Long 0-40mm Angled 6"****\$47.99** per pieceTrade price from **\$44.15**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Castroviejo Implant Bone Caliper 6" (15.2cm) Long Model – 30° Angled Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper Angled 15.2cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 6 inches and tapers from the measurement end to the tip. The shanks are straight with 30° Angled jaws

**LONG CASTROVIEJO CALIPER FEATURES**

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 15.2cm
- 30° Angled Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

**PACKAGE INCLUDE:**

- 3 Pcs Castroviejo Caliper 6" (15.2cm) Long Model – 30° Angled Tip Graduated 0-40mm

**TECHNICAL INFORMATION****WI-70001459****3 Pcs Castroviejo Caliper Long 0-40mm Angled 7"****\$49.99** per pieceTrade price from **\$45.99**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Castroviejo Implant Bone Caliper 7" (17.8cm) Long Model – 30° Angled Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper Angled 17.8cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 7 inches and tapers from the measurement end to the tip. The shanks are straight with angled jaws.

**LONG CASTROVIEJO CALIPER FEATURES**

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 17.8cm
- 30° Angled Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

**PACKAGE INCLUDE:**

- 3 Pcs Castroviejo Caliper 7" (17.8cm) Long Model – 30° Angled Tip Graduated 0-40mm

**TECHNICAL INFORMATION**

**WI-70003954****3 Pcs Castroviejo Caliper Long 0-40mm Straight 6"****\$47.99** per pieceTrade price from **\$44.15**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Castroviejo Implant Bone Caliper, Long Model 6" (15.2cm) Straight Tips, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Castroviejo Bone Caliper Straight 15.2cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 6 inches and tapers from the measurement end to the tip. The shanks and jaws are straight.

**LONG CASTROVIEJO CALIPER FEATURES**

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 15.2cm
- Pointed tips that open up for precise measurement.
- Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

**PACKAGE INCLUDE:**

- 3 Pcs Castroviejo Caliper Long 6" (15.2cm) Straight Graduated 0-40mm

**TECHNICAL INFORMATION****WI-70001463****3 Pcs Castroviejo Caliper Long 0-40mm Straight 7"****\$49.99** per pieceTrade price from **\$45.99**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Castroviejo Implant Bone Caliper 7" (17.8cm) Long Model – Straight Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper Straight 17.8cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 7 inches and tapers from the measurement end to the tip. The shanks and jaws are straight.

**LONG CASTROVIEJO CALIPER FEATURES**

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 17.8cm
- Straight Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

**PACKAGE INCLUDE:**

- 3 Pcs Castroviejo Caliper 7" (17.8cm) Long Model – Straight Tip Graduated 0-40mm

**TECHNICAL INFORMATION**

**WI-70003995****3 Pcs Cement Spatulas 7" (17.8cm) Double Ended****\$25.99** per pieceTrade price from **\$23.91**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Cement Spatulas Double Ended 7" (17.8cm) Lab Mixing Modeling Premium Quality Stainless Steel Dental Restorative Instruments.

Cement spatulas used for mixing cement materials and various other dental materials like silicones, waxes, cements etc. Dental spatula comes in double ends, one tip is rounded for scooping, mixing and the other is arch headed spreading tip used to mix, spread and lift material including foods, drugs, plaster and paints.

Dental cement spatulas available in double and single ended patterns with flexible and rigid forms which are ideal for lighter, medium to heavier cements and other restorative materials. These spatulas can also be used to load materials into crowns or other

preparations. Made with premium quality stainless steel with superior craftsmanship. Provides a reliable performance for dental professionals over time. Fully Autoclavable for hygiene and reuse purpose.

**DOUBLE ENDED CEMENT SPATULAS FEATURES**

- Crafted With 7" (17.8cm) Overall Length.
- Double Ended.
- One Rounded End and Other is Arch Headed.
- Optimal Size for Handling for Optimal Maneuvering.
- Non-slip Textured Area for Better Grip
- Highly Polished Rounded Handle
- Fully Autoclave and Reusable

**WI-70003993****3 Pcs Cement Spatulas Single End****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 3 Pcs Dental Cement Spatulas Single Ended 7" (17.8cm) Lab Mixing Modeling Stainless Steel Restorative Instruments.

Cement spatulas are used for mixing cement materials and various other dental materials like silicones, waxes, cements etc. A spatula is a broad, flat, flexible blade used to mix, spread and lift material including foods, drugs, plaster and paints. Available in double and single ended patterns with flexible and rigid forms which are ideal for lighter, medium to heavier cements and other restorative materials. Cement spatulas can also be used to load materials into crowns or other preparations. Made with premium quality stainless steel with superior craftsmanship. Cement spatulas can be

reused over time as they will not bend or tarnish on you. We offer 100% satisfaction guaranteed.

**SINGLE ENDED CEMENT SPATULAS FEATURES**

- Crafted With 7" (17.8cm) Overall Length.
- Single Ended.
- Large Width Spatula Design.
- Optimal Size for Handling for Optimal Maneuvering
- Non-slip Textured Area for Better Grip
- Highly Polished Rounded Handle
- Fully Autoclave and Reusable
- Made With High Quality Stainless Steel

**PACKAGE INCLUDE**

**WI-70004096****3 Pcs Dental Bone File #12CA Double Ended****\$23.99** per pieceTrade price from **\$22.07**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Bone File #12CA Double Ended 7" (17.8cm) Cross Serrated Working Tips with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Dental Instruments.

12CA Bone File is a cross serrated working ends dental instrument used to excavate small spicules of bones from narrow spaces using a pull stroke. This bone file usually maneuvers between the soft tissues and the bones that need scraping. The knurled design allows strong grip and easy maneuvering. Made with premium quality stainless instrument for durability and long-lasting performance efficiency. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for

a variety of precise neurosurgical procedures.

**BONE FILE #12CA FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- Cross-serrated Ends.
- Double-Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 2 Pcs Bone File #12CA 7" Double Ended

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Surgical Dental Instruments**WI-70004095****3 Pcs Dental Bone File #33 Double Ended****\$23.99** per pieceTrade price from **\$22.07**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Dental Bone File #33 Double Ended 7" (17.8cm) Long with 4.8mm Narrow Pattern Straight Plain Serrated Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

#33 Bone File is designed to smooth and shape bone using a pull stroke. Bone File 33 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. The straight serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 4.8mm wide. Constructed from surgical-grade stainless steel for excellent performance and longevity. The round


handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

**BONE FILE #33 FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- 4.8mm Narrow Pattern.
- Double Ended.
- Straight Parallel Sides Plain Serration.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 3 Pcs Bone File #33 Straight 4.8mm Double Double Ends

**3  
PCS****WI-70004228****3 Pcs Dental Extracting Forceps #150AS Upper Anteriors****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 American Pattern Dental Extracting Forceps #150AS Serrated Split Beaks For Upper Anteriors German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Extraction Forceps No. 150AS used for extracting teeth from the alveolar bone. Upper universal split serrated beaks designed to provide increased grasp beneath the crown to prevent crown or root breakage. Serrations provide increased mechanical grasp on the root without increased pressure. 150AS dental forceps can be used left and right.

American pattern forceps have a smooth hinge so no pin is visible. The stainless steel construction ensures maximum durability and dependability. You can reuse these 150AS extracting


forceps over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

**AMERICAN EXTRACTING FORCEPS #150AS FEATURES**

- For Upper Anteriors.
- Serrated, Split Beaks Provide an Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

**PACKAGE INCLUDE:**

- 3 Pcs Dental Extracting Forceps #150AS Split Beaks – Serrated

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental**3  
PCS****WI-70004234****3 Pcs Dental Extracting Forceps #210H Upper 3rd Molars****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Extracting Forceps #210H for Upper 3rd Molars with Hook Handle Premium Quality Stainless Steel American Pattern Bayonet Style Serrated Beaks Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 210H with hook handle is used to extract teeth from the alveolar bone. American extraction forceps #210H features offset (Bayonet) conical beaks designed to adapt to upper third molars with longitudinal serrations to improve grasp. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**AMERICAN PATTERN EXTRACTION FORCEPS #210H FEATURES**

- Offset (Bayonet) Conical Beaks Design.
- For Upper 3rd Molars.
- Longitudinal Serrations Improve Grasp.
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

**PACKAGE INCLUDES**

- 3 Pcs American Pattern Extraction Forceps #210H with Hook Handle

**TECHNICAL INFORMATION**

**WI-70004241****3 Pcs Dental Extracting Forceps #23 Cow Horn for Lower Molar****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Universal Extracting Forceps #23 Cow Horn 1st & 2nd Lower Molars with Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#23 extraction forceps is use to extract teeth from 1st and 2nd lower molars with cow horn jaws. Forceps dimension is 7.5" (19cm) overall in length with Textured grip handles. Durable and sturdy our extraction forceps are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #23 FEATURES**


- Slightly Larger Cowhorn Beak Design Than #16
- For Lower 1st & 2nd Molars
- Non-slip Grip Handle
- German Stainless Steel

**PACKAGE INCLUDE**

- 3 Pcs Extracting Forceps #23 – Cow Horn Jaws

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**WI-70004251****3 Pcs Dental Extracting Forceps #32****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Universal Extracting Forceps #32 for Upper Canines, Premolars & Molars with Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#32 Universal forceps used to extract teeth from Upper Cuspids, Bicuspid and Molars with sharp edged narrow beaks for breaking the alveolar process and trimming it after extraction. Forceps dimension is 7.5" (19cm) overall in length with textured grip handles. Durable and sturdy our extraction forceps are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #32 FEATURES**

- Narrow beak, beveled and smoothed edges
- For Upper Canines, Premolars & Molars
- Non-slip Grip Handle
- German Stainless Steel

**PACKAGE INCLUDE**

- 3 Pcs Extracting Forceps #32 – Marrow, Bayonet Pattern Beak

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes

**WI-70004215****3 Pcs Dental Extracting Forceps #53L For Upper 1st & 2nd Left Molars****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Extracting Forceps #53L for Upper 1st and 2nd Left Molars American Pattern German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

53L dental forceps are used to extract maxillary 1st and 2nd left molars. Left angled pointed offset (Bayonet) beaks for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #53L FEATURES**

- Left Pointed Tip Design
- For Upper 1st & 2nd Molars
- Non-Slip Textured-Grip Handle
- German Stainless Steel

**PACKAGE INCLUDE**

- 3 Pcs American Pattern Dental Extracting Forceps #53L

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**WI-70004248****3 Pcs Dental Extracting Forceps #65 for Upper Roots****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 American Pattern Dental Extraction Forceps #65 for Upper Roots, Fragments & Incisors, Slim Bayonet Pattern Jaws with Non-slip Handle German Grade Stainless Steel Teeth Extraction Dental Surgical Instrument.

American pattern dental extracting forceps No. 65 used to extract teeth from upper roots, fragments & incisors. Slim bayonet pattern beaks designed specially for hard to reach area.


American pattern forceps have a smooth hinge so no pin is visible. The top quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**EXTRACTING FORCEPS #65 FEATURES**

- For Upper Roots, Fragments & Incisors.
- Slim, Bayonet Pattern Beaks.
- Non-Slip Textured-grip Handle.
- Made With High Quality Stainless Steel.

**PACKAGE INCLUDE:**

- 3 Pcs American Pattern Extracting Forceps #65

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Extraction Forceps**BEAK/JAWS:** Slim Bayonet Jaws**POLISHED:** Yes**HANDLE:** Non-slip Handle

**WI-70004249**

### 3 Pcs Dental Extracting Forceps #32A Molars Root – Slim Beak

**\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Dental Extracting Forceps #32A for Upper Cuspids, Biscuspids, Molars and Roots with Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#32A Extraction forceps used to extract teeth from Upper Cuspids, Biscuspids, Molars and roots. Sharp edged Slim beaks for breaking the alveolar process and trimming it after extraction. Forceps dimension is 7.5" (19cm) overall in length with textured grip handles. Durable and sturdy our extraction forceps are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

#### EXTRACTING FORCEPS #32A FEATURES

- Slim, bayonet pattern beaks with parallel lines for optimal grip.
- For Upper Cuspids, Biscuspids, Molars and Roots
- Non-slip Grip Handle
- German Stainless Steel

#### PACKAGE INCLUDE

- 3 Pcs Extracting Forceps #32A – Slim, Bayonet Pattern Beak

#### TECHNICAL INFORMATION

**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**WI-70003978**

### 3 Pcs Dental Orthodontic Star-Like Boone Bracket Positioning Gauge Set

**\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 3 Pcs Boone Bracket Positioning Gauge with Four Integrated Measurement 3.5, 4.0, 4.5 and 5.0 Millimeter – Premium Quality Stainless Steel Dental Orthodontic Instruments for Orthodontists, Dentists, and Dental Assistants involved in Orthodontic Procedures.

Boone Bracket Positioning Gauge (Star-Like Gauge) is an orthodontic instrument used for precisely measuring and marking the vertical height for braces placement on teeth. It includes multiple tips sized at 3.5mm, 4.0mm, 4.5mm, and 5.0mm, aligning with standard distances from a tooth's biting edge, ensuring consistent bracket bonding during procedures. It also features pointed tips for accurate surface marking on the tooth crown.

Boone gauge facilitates consistent bracket placement on the incisal or occlusal surfaces and promotes uniform bonding heights across the dental arch, aiding proper tooth movement during treatment. Made with premium quality stainless steel ensures sterility for hygienic orthodontic practices.

#### BOONE BRACKET POSITIONING GAUGE FEATURES

- Four Measurement Options: 3.5, 4.0, 4.5, and 5.0 Millimeters.
- Pointed Tips with Flat Base Portion for Superficial Marking
- Measure Incisal or Occlusal Surface of the Tooth to Determine the Ideal Bracket Position.
- Withstand Repeated Use and Sterilization Procedures.





WI-25737862

3 Pcs Dental Pick Tools Set Plaque Remover Scaler

\$1.50 per piece Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

3 Pcs Dental Pick Tools Set Plaque Remover Scaler — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 3-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$1.50. Minimum order 10 sets, and the price moves with quantity. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-70004240

3 Pcs Extracting Forceps #53R for Upper 1st & 2nd Right Molars

\$29.99 per piece Trade price from \$27.59 Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Premium Quality American Pattern Dental Extracting Forceps #53R for Upper 1st and 2nd Right Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

53R dental forceps are used to extract maxillary 1st and 2nd right molars. Right angled pointed offset (Bayonet) beaks for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

EXTRACTING FORCEPS #53R FEATURES

- Right Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Angular Beaks.
- Non-Slip Textured-Grip Handle.
- German Stainless Steel.

PACKAGE INCLUDE

- 3 Pcs American Pattern Dental Extracting Forceps #53R

TECHNICAL INFORMATION

**MATERIAL:** Stainless Steel  
**INSTRUMENTS TYPE:** Dental Extracting Forceps  
**APPLICATION:** Teeth Extraction  
**POLISHED:** Yes  
**STERILE:** Yes  
**REUSABLE:** Yes  
**AUTOCCLAVED:** Yes  
**RUST PREVENTION PROCEDURE:** Passivated



**WI-70004163****3 Pcs Extracting Forceps #62 Upper and Lower Anterior Teeth****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Universal American Pattern Extracting Forceps #62 for Upper & Lower Incisors, Canines, Premolars and Deciduous Teeth American Pattern German Grade Stainless Steel Tooth Extraction Dental Surgical Instrument.

American Pattern Tooth Extraction Forceps No. 62 used to extract Upper and Lower Incisors & Canines. The long, extra sharp and very narrow beaks allow a tight grasp of the tooth. The slightly curved working end makes this instrument capable of handling both upper and lower teeth, in addition to minimizing trauma to the patient.

American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All

our dental tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #62 FEATURES**

- Precise Tips Meet at a Point, Ensuring Optimal Grip.
- Universal, Conical Beaks Design.
- For Upper & Lower Incisors, Canines, Premolars and Deciduous Teeth.
- Gently-Curved Working End Making It Suitable For Both Upper & Lower Teeth.
- Non-Slip Textured-Grip Handle.
- German Stainless Steel.

**PACKAGE INCLUDE**

- 3 Pcs American Pattern Extracting Forceps #62

**TECHNICAL INFORMATION****WI-70004225****3 Pcs Extracting Forceps Set #151AS Serrated Split Beaks****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Dental Extracting Forceps #151AS Split Beaks, Serrated for Lower Anterior – Standard – American Pattern German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Extraction Forceps No. 151AS used for extracting teeth from the alveolar bone. Lower universal split serrated beaks designed to provide increased grasp beneath the crown to prevent crown or root breakage. Serrations provide increased mechanical grasp on the root without increased pressure. This 151AS dental forceps can be used left and right.

American pattern forceps have a smooth hinge so no pin is visible. The stainless steel construction ensures maximum durability and dependability.

You can reuse these 151AS extracting forceps over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

**AMERICAN EXTRACTING FORCEPS #151AS FEATURES**

- For Lower Anterior.
- Serrated, Split Beaks Provide an Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

**PACKAGE INCLUDE:**

- 3 Pcs Extracting Forceps #151AS Split Beaks Serrated -Lower Anterior

**TECHNICAL INFORMATION**

**WI-70000641****3 Pcs Heidbrink Root Tip Pick Elevator Set H1, H2 & H3 Pointed Tips****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 3 Pcs Heidbrink Root Tip Pick Elevator Set H1, H2 & H3 Straight, Left and Right Angled 6.5" (16.5cm) Slightly Curved Pointed Tip with Heavy Duty Slim Hexagonal Handle Premium Quality Stainless Steel Root Tips Tooth Extraction Dental Instruments.

Heidbrink Root Teaser elevators are a single-ended instruments designed for the effective removal of root tips and residual roots in surgical procedures. This delicate dental elevators are used to retrieve loose root tips that may break away from a tooth during extraction and are in hard-to-reach places. Heidbrink also referred to as a root tip pick elevator with a narrow sharp slightly curved design and lightweight hexagonal handle for better control during use. Heidbrink Root Tip

Elevator is crafted with straight, left and right angled tips. This heidbrink root elevators with its unique design is worthy of being a part of your dental kit. Made with premium grade stainless steel, these elevators are lightweight, rust-free, and reusable.

**HEIDBRINK ROOT TIP PICK ELEVATOR FEATURES**

- Crafted with 6.5" (16.5cm) Overall Length.
- 2.5mm Working End.
- Narrow, Sharp and Slightly Curved Tips.
- Available in Straight, Left and Right Angled Profiles.
- Single Ended Instruments Used to Extract Root Tips.

**WI-70003932****3 Pcs Hollenback Wax Carver #3****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Hollenback Carver #3 Double Ended Spear Shape Sharp Blades with 6.5" Overall Length Premium Quality Stainless Steel Wax Modelling Dental Restorative Instruments.

3 hollenback carver is a dental instrument used for placing, carving, and contouring restorative materials like amalgam and composite. Hollenback #3 used for shaping restorative materials to resemble natural tooth anatomy, placing and packing materials into cavities, and trimming excess material after restoration. Double-ended design with large and sharp blades enable detailed work on both the chewing surface and sides of teeth, even in tight spaces. Stainless steel construction of hollenback carver 3 guarantees

durability and dependability, allowing for reuse without bending or tarnishing, and meets international quality standards.

**DENTAL HOLLENBACK CARVER #3 FEATURES**

- Large Pointed Spear Shape Sharp Blades.
- Double Ended.
- Non-Slip Knurled Grip Handle.
- Used for Wax and Porcelain Carving.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 3 Pcs Hollenback #3 Double Ended

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel

**WI-70004092****3 Pcs Miller Bone File #45 Double Ended****\$23.99** per pieceTrade price from **\$22.07**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Miller Bone File #45 Straight-Cut 7" (17.8cm) Long with 5mm Wide Straight & Curved Double Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

Bone file #45 is designed to smooth and shape bone using a pull stroke. Dental bone file #45 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. It can scrape off bony segments from soft tissues and bones. The serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features straight-cut plain serrations and is 5mm wide. Constructed from surgical-grade

stainless steel for excellent performance and longevity. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

**BONE FILE MILLER #45 FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- 5mm Straight-cut Plain Serrated Ends.
- Curved & Straight Double Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 3 Pcs Miller Bone File #45 Double Ended 5mm Straight & Curved Ends

**WI-70000539****3 Pcs Molt Bone Curette Set #2/4 Double Ended****\$23.99** per pieceTrade price from **\$22.07**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Molt Bone Curette Set #2/4 Double Ended 3.5mm & 7.0mm Working Tip Round Pattern Handle Premium Quality Stainless Steel Surgical Dental Instruments.

2/4 Molt surgical curette is a double ended bone curette designed for curettage, cyst removal, and tooth socket debridement. With a 3.5 mm and 7.0 mm working ends allows the surgeon to scrape away bone with durable and sharp edges. Molt #2/4 bone curette has round pattern non-slip knurled handle provides optimum control and omits hand fatigue.

This double ended bone curette can also use in a varteity of veterinary dental procdures for efficient scraping of diseased tissue from oral/dental wounds for cats and dogs. Molt Bone

curette 2-4, crafted with high-quality surgical stainless steel with a smooth surface, ensures a strong grip. Completely sterilizable and reusable

**MOLT BONE CURETTE #2/4 FEATURES**

- Crafted with 6" (15.2cm) Overall Length.
- Double Ended.
- 3.5 mm – 7.0mm Spoon Sharp Blades.
- Ergonomic Design Round Pattern Grip Handle.
- Sterilizable / Reusable.
- Made with Premium Quality Stainless Steel

**PACKAGE INCLUDE**

- 3 Pcs Molt Bone Curette #2/4 Double Ended



**WI-70004214****3 Pcs Pedo Dental Extracting Forceps #150S****\$29.99** per pieceTrade price from **\$27.59**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Universal Pedo Extracting Forceps #150S Upper Incisors & Roots German Grade Stainless Steel Teeth Extraction Dental Surgical Instruments.

Pedo extraction forceps No.150S for upper primary teeth and roots; universal forceps used for extracting teeth from the alveolar bone. The universal conical beaks designed to adapt to upper primary teeth and roots, may be used right or left. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**UNIVERSAL PEDO EXTRACTION FORCEPS #150S FEATURES**

- For Upper Incisors and Roos
- Conical Beaks for Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

**PACKAGE INCLUDE:**

- 3 Pcs Universal Pedo Extracting Forceps #150S

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Pedo Extraction Forceps**BEAK/JAWS:** Conical Beaks**POLISHED:** Yes**STERILE:** Yes**WI-70003901****3 Pcs Periosteal Elevator Baby Molt #7A Double Ended with Round Handle****\$19.99** per pieceTrade price from **\$18.39**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Periosteal Elevator Baby Molt #7A Double Ended 6.5" (16.5cm) Round Pattern Grip Handle Premium Quality Stainless Steel Dental Instruments.

Molt-7 elevators have thin blades on both ends that are used to cut periodontal ligaments and loosen teeth. Periosteal Baby Molt #7A comprises a flat, oval tip at each end, one blunt and one sharp. Surgeons use the blunt oval tip for lifting and retracting the periosteal flap, and the sharp tip to dissect it. In surgery dentist often needs to separate a bone or tooth from the periosteum that covers it. This is done with a periosteal elevator. The dentist may also use it to gain access to retained roots and surrounding bone. Elevators are manufactured from

stainless steel to ensure optimal performance.


**PERIOSTEAL ELEVATOR BABY MOLT #7A FEATURES**

- Overall Length 6.5" (16.5cm).
- Double Ended.
- One Flat and One Oval Tip.
- Non-Slip Knurled Grip Handle.
- Use to Elevate, Retract and Dissecting Periosteal Flap.
- Corrosion Resistance.
- Fully Autoclavable & Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 3 Pcs Periosteal Elevator Baby Molt #7A Double Ended

**TECHNICAL INFORMATION**

**WI-70003902****3 Pcs Periosteal Seldin #S23 Dental Elevator Double Ended****\$23.99** per pieceTrade price from **\$22.07**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs #S23 Periosteal Elevator Seldin 8" (20.3cm) Double Ended Large Flap Retractor / Small Periosteal 10mm & 13mm wide Slightly Curved Blade with Large Flat Handle Quality Stainless Steel Dental Surgical Instruments

Seldin Periosteal Elevator #S23 are designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. Seldin S23 periosteal elevator designed with a large flat handle with double ended design approximately 10mm and 13mm wide blades. The smaller end is sharp for tissue reflecting, and the larger end is slightly curved and is a non-cutting blade, used for tissue retracting. Periosteal elevator S23 manufactured from high grade stainless steel to exacting specifications to ensure

optimal performance. Seldin elevators is autoclavable and crafted from high-quality stainless steel ensuring durability and sustained performance. 100% satisfaction guaranteed

**PERIOSTEAL ELEVATOR SELDIN #S23 FEATURES**

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Broad Rounded End & Broad Sharp Chisel End.
- Slight Curved 10mm and 13mm Wide Blades.
- Large Flat Handle.
- Used to Elevate the Periosteum
- Reusable.
- Made With Premium Quality Stainless Steel

**100% HIGH QUALITY  
RUST RESISTANT STEEL****WI-70000684****3 Pcs Polokoff Bone File #1 Double Ended****\$23.99** per pieceTrade price from **\$22.07**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Polokoff Bone File #1 Double Ended 7" (17.8cm) Long with 3.4mm Narrow Pattern Straight Serrated Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

#1 Polokoff Bone File is broadly used in surgical procedures to aid in mild bone reshaping using a pull stroke. Polokoff Bone File is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. The straight serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 3.4mm wide. Constructed from surgical-grade stainless steel for excellent

performance and longevity. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

**POLOKOFF BONE FILE #1 FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- 3.4mm Narrow Pattern.
- Double Ended.
- Straight Parallel Sides Serration.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 3 Pcs Polokoff Bone File 1 Straight 3.4mm Narrow Pattern Working

DENTAL BONE FILE POLOK OFF # 1

**3  
PCS**

**WI-70003953****3 Pcs Spring Divider Calipers Set 4.5" Dental Lab Tool****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 3 Pcs Spring Divider Calipers Set 4.5" (11.4cm) with Flat Legs Stainless Steel Precision Tool Bone for Accurate Measurements During Dental Lab Procedures.

The dental calipers measures the width of metal, wax, and crown size in dental cosmetic procedures. It has a scale or screw for the measuring of sizes in the dental lab. Spring Divider Caliper features two adjustable legs to measure hard-to-measure dimensions, made from carefully selected flat stock with heat treated points for durability, and they pivot on a roller while the bow spring ensures a precise tension setting. Accurate adjustments are made by turning the knurled nut on the threaded rod. The design of these calipers, with their fine adjustment

feature, ensures quick and positive settings when taking or transferring specific measurements. Manufactured from high quality stainless steel with a polished finish.

**SPRING DIVIDER CALIPER FEATURES**

- Crafted With 4.5" (11.4cm) Overall Length.
- Flat Legs with Sharpened Hard and Pointed Tips.
- Adjustable Lockable Steel Wing Divider.
- Rivet Design to Increase Stability.
- Finely Threaded to Achieve Fine Adjustment.
- Durable and Rustproof.
- Reusable.

**WI-70002177****3 Pcs T/C Wynman Crown Remover Gripper Forceps 7" Curved Serrated Tips****\$34.49** per pieceTrade price from **\$31.73**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs T/C Wynman Crown Remover Gripper Forceps Set 7" (17.78cm) Tungsten Carbide Tips with Gold Ring Handle Premium Quality Stainless Steel Temporary Crown Remover Surgical Dental Forceps Instruments.

Wynman Crown Remover Gripper Forceps features a serrated tip made from durable tungsten carbide for a long-lasting working end. Ideal to remove or place temporary crowns with ease, even in the most hard-to-reach places. Our TC Wynman crown grippers are hemostat style forceps that enable a strong grip on the temporary crown without slippage, keeping the crown intact even in the tightest areas of the mouth. Crafted from stainless steel and tungsten carbides,

this instrument offers seamless crown removal solution

**TC WYNMAN CROWN GRIPPER FORCEPS FEATURES**

**3 PCS TC WYNMAN CROWN GRIPPER FORCEPS:** The TC Wynman Crown Gripper Forceps are an ideal addition to your inventory. The tungsten carbide inserts on the Dental Crown Gripper Forceps will grip the temporary crown firmly without slipping, keeping the crown intact. They are useful for gripping wet glass such as in porcelain restorations and assist in more precise procedures.

**O.R. GRADE T/C AND STAINLESS STEEL:** Hot drop-forged construction Wynman crown gripper constructed of O.R. Grade stainless steel and tungsten carbide for the highest degree of



**WI-70000931****3 Pcs Warwick-James Elevators Set Right, Left & Straight****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

3 Pcs Warwick-James Elevators Set Right, Left and Straight Premium Grade Stainless Steel Teeth Extraction Surgical Dental Instruments.

Warwick-James elevators used for extracting fractured roots and also used in extraction procedure for elevating, loosen or removal root fragments or teeth from periodontal ligament. This Warwick-James Root Elevators set features left, right, or straight tip with a hexagonal handle that makes it ideal for easy tooth extraction. It is lightweight, unbendable, and reliable for dental procedures. Our multifaceted elevators are widely preferred for veterinary dental procedures for tooth extraction. We can design dental extraction custom kit packs according to surgeon's specific

needs and procedure. Manufactured with premium quality stainless steel for optimal results and precision.

**WARWICK-JAMES TOOTH ELEVATORS SET INCLUDES**

- 1 Warwick-James Elevator Straight
- 1 Warwick-James Elevator Left
- 1 Warwick-James Elevator Right

**WARWICK-JAMES ELEVATORS SET FEATURES**

- Warwick-James Set Consists of Straight, Left and Right Tips.
- Commonly Utilized to Extract Fractured Roots or Loosen and 'elevate' the Teeth in Their Sockets Prior to Extraction.
- The Single Working End of the Elevator is Designed for the Function of the Instrument.

**WI-70002285****3 Root Tip Pick Elevators Set #78 #79 & #80****\$18.74** per pieceTrade price from **\$17.24**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Root Tip Pick #78 #79 & #80 with Standard Handle – Oral Surgery Dental Extraction Instruments Set

Root Tip Pick Elevators are ideal for removing root tips from narrow areas during dental procedures. These are delicate instruments for the removal of root tips in difficult-to-reach areas.

blade inserts between bone and root with minimal damage to bone structure alternative to standard elevators for less traumatic extractions.

**USE OF DENTAL ELEVATORS:** Elevators used to loosen teeth prior to forcep extraction, to remove roots or impacted teeth, when teeth are compromised and susceptible to fracture or when they are malpositioned and cannot be reached with forceps. Highly polished finish for aesthetic and corrosion resistance.


**HIGH QUALITY PRECISION:** High Degree of Precision and Flexibility while conducting the Clinical Procedure. Our Product Conforms to all quality standards. Premium Quality Stainless Steel with Superior Craftsmanship.

**3 ROOT TIP PICK ELEVATORS SET INCLUDES:**

- 1 Pc Root Tip Pick #78 Straight
- 1 Pc Root Tip Pick #79 Left Angle
- 1 Pc Root Tip Pick #80 Right Angle

**DENTAL ELEVATORS FEATURES**

**DELICATE, SHARP INSTRUMENTS:** Cut periodontal ligament while preserving bone anatomy thin, tapered tips ensure smooth widening of the alveolar. Fine



**WI-70000542**

### 3 Swivel Head Bracket Height Gauge Set 6" Adjustable Pointed Tips

**\$34.99** per pieceTrade price from **\$32.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

3 Swivel Head Bracket Height Gauge Set 6" (15.2cm) Adjustable Pointed Tips Premium Quality Stainless Steel Wick Positioning Bracket Gauge for Both .018" & .022" Brackets Dental Orthodontic Instruments.

Swivel Head Bracket Height Gauge is an orthodontic tool used to accurately measure bracket height for placement from incisal edge, featuring articulated wick type tips for posterior tooth fixation and adjustable ends for four bracket height dimensions from the occlusal surface. Bracket Height Gauge can be adjusted for easy use on anterior, cuspid, and bicuspid brackets due to its 90° angled hinge on both ends. Our swivel head bracket height gauge set combines 3 different size adjustable bracket positioner that gives

you a handy experience and ensures proper bracket placement during orthodontic procedures.

#### ADJUSTABLE BRACKET HEIGHT GAUGE FEATURES

- 6" (15.2cm) Overall Length.
- Double Ended.
- Four Bracket Height Dimensions Gauges.
- Pointed Adjustable Tips for Superficial Marking at Different Angles.
- Ideal for Use on .018" and .022" Bracket Slots.
- Non-Slip Grip Handle.
- Fully Autoclavable / Reusable.
- Made with Premium Grade Stainless Steel.

#### PACKAGE INCLUDE

**WI-70003887**

### 3 Warwick-James Elevators Set with Sterilization Cassette

**\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

3 Pcs Warwick James Elevators Set Right, Left & Straight with Sterilization Cassette Tray Premium Quality Stainless Steel Teeth Extraction Surgical Dental Instruments.




Warwick-James elevators used for extracting fractured roots and also used in extraction procedure for elevating, loosen or removal root fragments or teeth from periodontal ligament. This Warwick-James Root Elevators set features left, right, or straight tip with a hexagonal handle that makes it ideal for easy tooth extraction. It is lightweight, unbendable, and reliable for dental procedures. Our multifaceted elevators are widely preferred for veterinary dental procedures for tooth extraction. We can design dental extraction custom kit

packs according to surgeon's specific needs and procedure. This elevator set comes with stainless steel sterilization container with Rubber inserts for efficient organization. Manufactured with premium quality stainless steel for optimal results and precision.

#### PACKAGE INCLUDES

- 1 Warwick-James Elevator Straight
- 1 Warwick-James Elevator Left
- 1 Warwick-James Elevator Right
- 1 Sterilization Cassette 7" x 3.5"

#### TECHNICAL INFORMATION

**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Elevators**APPLICATION:** Teeth Extraction**HANDLE:** Hexagonal Handle**RUST PREVENTION PROCEDURE:** Passivated

**WI-70004270****30 Pcs Amalgam Carrier Double Ended Set  
Large/Medium 2.5mm/2mm****\$179.99** per pieceTrade price from **\$165.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 30 Pcs Amalgam Carrier Double Ended Set 8" (20.3cm) Large/Medium 2.5mm – 2mm Tapered Tips with Non-Slip Textured Grip Handle Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.




Dental Amalgam Filling Carrier is a durable, high-quality instrument used to fill dental cavities with amalgam.

Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam

carrier. The textured surface helps effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

**DOUBLE ENDED AMALGAM CARRIER FEATURES**

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Color Bands for Easy Instruments Identification.
- Large – Regular (2.5mm – 2mm) Tapered Tips for Even Distribution of Filling Material.

**WI-70000517****4 Dental Amalgam Carrier Double Ended Set****\$35.99** per pieceTrade price from **\$33.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 4 Amalgam Carrier Double Ended Set 8" (20.3cm) Tapered Tips with Non-Slip Textured Grip Handle Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.


This 4 Pcs dental Amalgam carrier double ended set included various tips which dental professional needs to fill dental cavities with amalgam of all levels. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back

end of the amalgam carrier. The textured surface helps effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

**DOUBLE ENDED AMALGAM CARRIER FEATURES**

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Color Bands for Easy Instruments Identification.
- Various Tapered Tips for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.



**WI-70000350****4 Lucas Bone Curette Set #85, #86, #87, #88 Double Ended****\$23.99** per pieceTrade price from **\$22.07**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 4 Pcs Lucas Bone Curette Set #84, #85, #86, #87, #88 Spoon Shaped Double Ended Premium Quality Stainless Steel Surgical Curettes Dentist Lab Periodontal Dental Intraoral Surgery Instruments.

Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° various wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium


quality stainless steel, it ensures durability and lasting performance. These curettes are reusable after sterilization.

**LUCAS CURETTE FEATURES**

- Spoon Shaped Sharp Ends.
- Double-ended.
- Non-Slip Strong Grip Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 1 Pc Lucas Bone Curette #85 Small 2.7mm Double Ended
- 1 Pc Lucas Bone Curette #86 Small 3.0mm Double Ended
- 1 Pc Lucas Bone Curette #87 Small 3.8mm Double Ended

**WI-70004281****4 Pcs Ball Burnishers Set #27/29 Double Ended****\$25.99** per pieceTrade price from **\$23.91**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 4 Dental Ball Burnisher #27/29 Double Ended 1.6 mm – 3.1 mm Working Ends Premium Quality Stainless Steel De Amalgam Restorative Composite/Plastic Filling Instruments.

27/29 ball burnisher used to blend material for final contouring, to achieve sculpting of areas like grooves, fissures or pits. Can also be used to form occlusal anatomy. The working end of a burnisher is smooth and rounded (1.6 mm – 3.1 mm). Made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. These dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

**BALL BURNISHER # 27/29 FEATURES**

- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel

**PACKAGE INCLUDE:**

- 4 Pcs Ball Burnishers #27/29 Double Ended

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Ball Burnisher Dental Instrument

**WI-70001443**

## 4 Pcs Dental Amalgam Well Non-Slip Restorative Carrier

**\$51.99** per pieceTrade price from **\$47.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 4 Dentist Bone Grafting Amalgam Mixing Well Pot Non-Slip Base Stainless Steel Dentist Lab Restorative Dental Surgical Instruments.

Amalgam well pot used to mix and dispense amalgam during restorative procedures. Dentist use this non-slip small bowl shaped container to holds mixed amalgam prior to its being loaded into the amalgam carrier. amalgam well used to hold amalgam until it is used in the restoration. Dental restorative amalgam well can also be used for bone grinding/crushing in dental implant procedures.


### DENTAL AMALGAM WELL FEATURES:

**AMALGAM WELL POT:** Amalgam well have a feature is 3cm diameter, Heavy well pot with secure base with rubber friction locking ring. Machined from a solid piece of steel with a cone shaped inner well, Smooth cavity for optimal mixing and flat cone tip for easy amalgam carrier transfer.

**NON-SLIP RUBBER RING BOTTOM:** Non-slip removable rubber 'O' ring and bottom heavily weighted to avoid spillage.

**STAINLESS STEEL CONSTRUCTION:** The stainless steel construction ensures maximum durability and dependability. You can reuse these amalgam well over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

**EASY TO CLEAN:** Amalgam well is a

**WI-70001740**

## 4 Pcs Spatulas Wax & Clay Sculpting Tool Carver Set

**\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

4 Pcs Stainless Steel Spetulas Set – Dental Cement Lab Mixing Wax Modeling Carving Spatulas Ceramic Pottery Clay Sculpting Tools Double Ended. Excellent for a Variety of Uses Within the Dentla Lab Environment.

### METAL SPETULAS SET FEATURES:

**4 PCS SPATULAS CARVER SET:** Stainless Steel Spatulas set includes 4 premium quality commonly used dental wax carving spatulas for a variety of different applications. You can cut, clean molds, engrave, make designs, repair, scrape, and shape ceramics, clay, electronics, jewelry, wood, or details of other items. Each tool included in this set is double-sided with different types of tips for creating a wide variety of shapes and projects.

**DURABLE STAINLESS STEEL**

**CONSTRUCTION:** You can depend on this spatula carver set to last you many years of constant use. Built to never rust, break or bend even with continuous use, they are a valuable addition to any medical or home setting.

**PRECISE AND EASY TO USE:** Each wax carving spatula is perfectly balanced and designed to be easy to maneuver for precise use. The stainless steel kit includes instruments with built-in finger grips for ultimate comfort and user-friendliness. The average length is 7", ideal for many uses. Great for carving, cutting, modeling, scraping, shaping, etc.

**EXCELLENT DENTAL WAX CARVERS KIT:** Perfect Dental Laboratory Equipment for wax carving and mold forming. High Quality and Durable



**WI-70003919****4 Pcs Wax & Modelling LeCron Zahle Beale Carver Cement Spatula Set****\$16.99** per pieceTrade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 4 Pcs Wax & Modelling Tools Set LeCron, Zahle & Beale Carver Cement Spatula Double Ended Dental Alignate Mixing Spatula Dental Lab Mixing Instruments.


Dental wax & modelling kit of 4 pcs Zahle, Beale, LeCron Cement Spatula double-ended wax carvers designed for efficient carving. This versatile dental tool set utilized for accurate wax modeling, sculpting, and contouring of crowns, bridges, and dentures, in addition to mixing and applying dental cements and impression materials by technicians and students in laboratories and clinics. These essential instruments enhance dental lab work, ensuring accuracy and efficiency for both professionals and students. Made of stainless steel, the carvers prevent rust

and are ergonomically shaped for comfort and balance.

**4 PCS WAX & MODELLING CARVERS SET INCLUDES**

Dental lab instruments set consists of various essential instruments.

- 1 Pc Cement Spatula – Double Ended
- 1 Pc Beale Wax Carver – Double Ended
- 1 Pc LeCron Wax Carver – Double Ended
- 1 Pc Zahle Wax Carver – Double Ended

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Lab Instruments**APPLICATION:** Alignate Mixing, Wax & Modelling, Sculpting**WI-70000533****5 Lucas Bone Curette Set #84, #85, #86, #87, #88 Double Ended****\$30.99** per pieceTrade price from **\$28.51**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Lucas Bone Curette Set #84, #85, #86, #87, #88 Spoon Shaped Double Ended Premium Quality Stainless Steel Surgical Curettes Dentist Lab Periodontal Dental Intraoral Surgery Instruments.

Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° various wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium



quality stainless steel, it ensures durability and lasting performance. These curettes are reusable after sterilization.

**LUCAS CURETTE FEATURES**

- Crafted With 6.5" (16.5cm) Overall Length.
- Spoon Shaped Sharp Ends.
- Double-ended.
- Non-Slip Strong Grip Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 1 Pc Lucas Bone Curette #84 Small 2.4mm Double Ended
- 1 Pc Lucas Bone Curette #85 Small 2.7mm Double Ended



**WI-70000704****5 Pcs Ball Burnisher #34 Single Ended****\$24.99** per pieceTrade price from **\$22.99**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Ball Burnisher #34 Single Ended 6.5" (16.5cm) Long T-shape Burnisher With a Ball on One Side Premium Quality Stainless Steel Dental Amalgam Composite Filling Instruments.

Single Ended ball burnisher #34 is often used to smooth the surface of a freshly placed amalgam restoration. T-shaped working end with a small ball on one side of "T" and paddle on the other. The shape of burnisher #34 conforms to the rounded cavity surfaces and allows ease of access into the rounded corners or junctions of the cavity surfaces to condense and shape the composite against the cavity wall. All our amalgam composite filling instruments are manufactured for optimal results and precision. These



dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

**BALL BURNISHER #34 FEATURES**

- Crafted with 6.5" (16.5cm) Overall Length.
- Single Ended.
- Small Pattern.
- T-Shaped Working End.
- Grip Design Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel.
- Reusable.
- Hand Crafted From High-Grade Stainless Steel.

**PACKAGE INCLUDE:**

- 5 Pc Ball Burnisher #34 Single Ended 6.5" (16.5cm) Long

**WI-70004074****5 Pcs Bone Files Kit Dental Surgical Orthopedic Veterinary Surgery****\$39.99** per pieceTrade price from **\$36.79**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Bone Files Set with Carrying Case Including #64, #45, #12, #1PO and #33 – Premium Quality Stainless Steel 8 Dental Bone File Double Ended for Precise and Delicate Dental Surgical Clinical Procedures.

5 piece dental bone file set comprises surgical-grade stainless steel instruments used in oral, periodontal, and orthopedic surgeries. These bone files designed to smooth and shape alveolar bone, primarily used for post-extraction smoothing, alveoloplasty and flap surgery. Features a double-ended configuration for access to hard-to-reach areas and serrated cutting edges designed for efficient bone tissue removal with minimal trauma. The instruments come in varying sizes to accommodate different surgical needs.

Made from surgical-grade stainless steel ensures outstanding durability and performance over time. These dental surgical bone files can be reused after proper sterilization.

**DENTAL BONE FILES FEATURES**

- Double Ended.
- Serrated Tips.
- Round Non-Slip Handles.
- Available in Different Shapes and Diameters.
- Ideal to Smooth and Shape Rough Edges Bone Spicules During Surgical Procedure.
- Fully Autoclavable/Reusable.
- Hand Crafted from Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 1 Pc Bone File #64 – Double Ended



**5  
PCS**

**WI-70002027**

## 5 Pcs Christensen Crown Remover Set CRCH2 90 Degree Angled

**\$43.99** per piece

Trade price from **\$40.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Christensen Crown Remover CRCH2 90 Degree Angled with Hexagonal Dental Elevator Handle Premium Quality Stainless Steel Permanent Crown Removal Restoration Dental Instruments

For permanent removal of crowns by breaking the seal between tooth and crown after sectioning with a bur. Crown Remover Christensen is ideally constructed for the disassembly method of crown removal. After sectioning the crown with a bur, the elevator handle provides a comfortable grip and maximum leverage. The notched tip then aids the single-end Christensen Crown Remover in effortless removal. Pressure against the tooth is lessened, which reduces the potential for tooth fracture.

### CHRISTENSEN CROWN REMOVER FEATURES


- 90-degree angled CRCH2 Crown Remover is ideally suited for posterior teeth.
- Elevator handle for optimal comfort and grip.
- Notched tips.
- Reusable.
- Made with Highest quality medical grade Stainless Steel.

>

### PACKAGE INCLUDE:

- 5 Pcs Christensen Crown Remover CRCH2 90 Degree Angled with Hexagon Handle

### TECHNICAL INFORMATION:



**WI-70002286**

## 5 Pcs Dental Elevators #301, #34, #34S, #46 & #81 Set

**\$24.99** per piece

Trade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Dental Elevators Set Premium Quality Stainless Steel Osteotome Periotome Luxating Elevator Apical Molt Oral Implant Surgical Root Extraction Spreading Tools for Dental Surgical Procedure

Elevators are used to loosen teeth prior to forceps extraction, to remove roots or impacted teeth, when teeth are compromised and susceptible to fracture or when they are malpositioned and cannot be reached with forceps. These delicate, sharp instruments cut periodontal ligament while preserving bone anatomy with thin, tapered tips to ensure the smooth widening of the alveolar. Fine blade inserts between bone and root with minimal damage to bone structure


alternative to standard elevators for less traumatic extractions

### 5 PCS DENTAL ELEVATORS SET INCLUDES:

- 1 Pc Dental Elevator #301
- 1 Pc Dental Elevator #34
- 1 Pc Dental Elevator #34S
- 1 Pc Dental Elevator #46
- 1 Pc Dental Elevator #81

### DENTAL ELEVATORS FEATURES

**DELICATE, SHARP INSTRUMENTS:** Cut periodontal ligament while preserving bone anatomy thin, tapered tips ensure smooth widening of the alveolar. Fine blade inserts between bone and root with minimal damage to bone structure



**WI-70001922****5 Pcs Periosteal Dental Elevator Set****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Periosteal Dental Elevator Set premium Quality German Grade Stainless Steel Dental Extraction Surgical Instruments.

Periosteal Elevator is designed to lift the full thickness of soft tissue flaps during dental practices. Likewise, it will separate the gingiva present around the tooth's neck. These elevators are reliable for tooth removal practices as they provide maximum leverage and an optimal force during dental extraction procedures. These dentist elevators come in single or double-ended structures, slightly curved, and thin tip to accommodate less traumatic manipulation. The high-quality material makes the tools suitable for health care providers and enhances the patient's comfort level.

**SET INCLUDES**

- 1 Pc Periosteal Elevator #9 Molt (#M9).
- 1 Pc Periosteal Elevator Seldin #23 Straight – Double Ended.
- 1 Pc Periosteal Elevator #OHL – Double Ended.
- 1 Pc West #W2 Periosteal Elevator – Double Ended.
- 1 Pc GF 14 Periosteal Elevator – Double Ended.

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Instruments**APPLICATION:** Teeth Extraction**RUST PREVENTION PROCEDURE:** Passivated**SURFACE:** Polished**WI-70002722****5 Pcs Wax Carving Metal Clay Pottery Sculpting Tool Set****\$10.99** per pieceTrade price from **\$10.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

5 Pcs Spatula Wax & Clay Sculpting Tool Set – Premium Quality Stainless Steel DIY Waxing Modeling Kits Stainless Steel Wax Carving Tool Dental Clay Pottery Plaster Sculpture Carvers Tools Kit Dabber Tool Set

**5 PCS WAX CARVING TOOLS SET FEATURES:**

**5 PCS WAX CARVERS KIT:** The Kit includes 5 Premium Quality commonly used wax carving instruments for a variety of different applications. You can cut, clean molds, engrave, make designs, repair, scrape, and shape ceramics, clay, electronics, jewelry, wood, or details of other items. Each tool included in this set is double-sided with different types of tips for creating a wide variety of shapes and projects.

**DURABLE STAINLESS STEEL**

**CONSTRUCTION:** The Premium Quality stainless steel construction of our Wax Carvers Kit ensures long-lasting durability over many years of use. Our wax carvers kit is built to last, corrosion-resistant with a high mirror finish, and tarnish-free.

**PRECISE AND EASY TO USE:** Each double-sided wax carving tool is perfectly balanced and designed to be easy to maneuver for precise use. The stainless steel kit includes instruments with built-in finger grips for ultimate comfort and user-friendliness. Great for carving, cutting, modeling, scraping, shaping, etc.

**EXCELLENT DENTAL WAX CARVERS KIT:** Perfect Dental Laboratory Equipment for wax carving and mold forming. High Quality and Durable Dentist Sculpture Knife Instruments

**WI-70000702****50 Pcs Ball Burnisher #34 Single Ended****\$124.99** per pieceTrade price from **\$114.99**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 50 Pcs Ball Burnisher Set #34 Single Ended 6.5" (16.5cm) Long T-shape Burnisher With a Ball on One Side Premium Quality Stainless Steel Dental Amalgam Composite Filling Instruments.

Single Ended ball burnisher #34 is often used to smooth the surface of a freshly placed amalgam restoration. T-shaped working end with a small ball on one side of "T" and paddle on the other. The shape of burnisher #34 conforms to the rounded cavity surfaces and allows ease of access into the rounded corners or junctions of the cavity surfaces to condense and shape the composite against the cavity wall. All our amalgam composite filling instruments are manufactured for optimal results and precision. These


dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

**BALL BURNISHER #34 FEATURES**

- Crafted with 6.5" (16.5cm) Overall Length.
- Single Ended.
- Small Pattern.
- T-Shaped Working End.
- Grip Design Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel.
- Reusable.
- Hand Crafted From High-Grade Stainless Steel.

**PACKAGE INCLUDE:**

- 50 Pcs Ball Burnisher #34 Single Ended 6.5" (16.5cm) Long

**WI-70000603****50 Pcs Football Burnisher #28/29 Set 2.3mm - 3.6mm Double Ended****\$124.99** per pieceTrade price from **\$114.99**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 50 Pcs Football Burnisher #28/29 Double Ended 2.3 mm - 3.6 mm Football Working Tips Premium Quality Stainless Steel Amalgam Restorative Composite/Plastic Filling Instruments.


28/29 football end burnisher is an operative hand instrument for dental professionals used to condense, carve and polish restorative materials such as amalgam. Ball Burnishers can also be used to form occlusal anatomy. The working end of a #28/29 ball burnisher is smooth and rounded like football (2.3 mm - 3.6 mm) made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. Fully autoclavable for reuse purposes.

**BALL BURNISHER # 28/29 FEATURES**

- Crafted with 6.5" (16.5cm) Overall Length.
- 2.3mm Football & 3.6mm Football.
- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Reusable.
- Smooth Finish For Excellent Tactile Feel

**PACKAGE INCLUDE:**

- 50 Pcs Football Burnisher #28/29 Double Ended





WI-30694736

### 5-Piece Stainless Steel Double-Sided Cement Spatula Set Ergonomic

**\$5.00** per piece

Trade price from **\$3.00**  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

We manufacture 5-Piece Stainless Steel Double-Sided Cement Spatula Set Ergonomic for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: an ergonomic handle for long lists without fatigue.

Supplied as a set — the contents can be adjusted to your list. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$5.00, falling to \$3.00 at 10 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-70004084

### 6 Bone File Kits of 5 Pcs Dental Surgical Orthopedic Veterinary Surgery

**\$199.99** per piece

Trade price from **\$183.99**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

6 Set of 5 Pcs Bone Files Kit with Carrying Case Including #64, #45, #12, #1PO and #33 – Premium Quality Stainless Steel 8 Dental Bone File Double Ended for Precise and Delicate Dental Surgical Clinical Procedures.

5 piece dental bone file set comprises surgical-grade stainless steel instruments used in oral, periodontal, and orthopedic surgeries. These bone files designed to smooth and shape alveolar bone, primarily used for post-extraction smoothing, alveoloplasty and flap surgery. Features a double-ended configuration for access to hard-to-reach areas and serrated cutting edges designed for efficient bone tissue removal with minimal trauma. The instruments come in varying sizes to accommodate different surgical needs.

Made from surgical-grade stainless steel ensures outstanding durability and performance over time. These dental surgical bone files can be reused after proper sterilization.

#### DENTAL BONE FILES FEATURES

- Double Ended.
- Serrated Tips.
- Round Non-Slip Handles.
- Available in Different Shapes and Diameters.
- Ideal to Smooth and Shape Rough Edges Bone Spicules During Surgical Procedure.
- Fully Autoclavable/Reusable.
- Hand Crafted from Premium Quality Stainless Steel.

#### PACKAGE INCLUDE

- 6 Pcs Bone File #64 – Double Ended

**WI-70004017****6 Bone Rongeurs Surgical Dental Oral Surgery Set****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Pack of 6 Pcs Bone Rongeurs Set Including Blumenthal, Mead, and Friedman – Premium Quality Stainless Steel Heavy-duty for Cutting and Removing Bone During Oral, Maxillofacial, Orthopedic and Neurosurgical Procedures.

This pack of 6 bone rongeurs featuring Blumenthal, Mead, and Friedman variations, is specifically designed for surgical tasks involving cutting, trimming, shaping, and scooping small bone or cartilage. They are utilized in orthopedic, spinal, neurosurgical, and dental procedures to facilitate surgical access by clearing bone. Exclusively used in remodeling bone surfaces in orthopedic surgery, removing bone fragments in spine and neurosurgery, and trimming alveolar bone in dental

surgery. The multi-pack offers various sizes and jaw configurations with spring-loaded hinged handles, allowing surgeons to adapt to different anatomical depths and angles for efficient cutting. Constructed from premium grade surgical stainless steel, these instruments ensure durability and precision.

**PACKAGE INCLUDE:**

- 1 Pc Blumenthal Bone Rongeur 6" (15.2cm) – 30 Degree
- 1 Pc Blumenthal Bone Rongeur 6" (15.2cm) – 45 Degree
- 1 Pc Blumenthal Bone Rongeur 6" (15.2cm) – 90 Degree
- 1 Pc Mead 1A Bone Rongeur 6.5" (16.5cm) – Round Jaws

**WI-70003924****6 P.K Thomas #2 Dental Wax Modelling Carver Set****\$35.99** per pieceTrade price from **\$33.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

PTK2 – 6 Pcs P.K Thomas #2 Dental Wax Modelling Carver with Elongated Curve on Thinner Tips Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PPKT2 features curved on thinner tapered tips used to flow on molten wax. P.K. Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture when working on hard-to-reach areas. Made of high



quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

**P.K THOMAS #2 WAX CARVER FEATURES**

- Elongated Curve on Thinner Tips Than P.K. Thomas #1.
- Double-ended
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

**PACKAGE INCLUDE**

- 6 Pcs Thomas Waxing Instrument – PKT2



**WI-70004258****6 Pcs Amalgam Carrier Double Ended Set  
Large/Medium 2.5mm/2mm****\$47.99** per pieceTrade price from **\$44.15**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Amalgam Carrier Double Ended Set 8" (20.3cm) Large/Medium 2.5mm – 2mm Tapered Tips with Non-Slip Textured Grip Handle Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument used to fill dental cavities with amalgam.

Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam

carrier. The textured surface helps effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

**DOUBLE ENDED AMALGAM  
CARRIER FEATURES**

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Color Bands for Easy Instruments Identification.
- Large – Regular (2.5mm – 2mm) Tapered Tips for Even Distribution of Filling Material.

**WI-70004269****6 Pcs Amalgam Carrier Single Ended Set  
1.5mm Mini Tip****\$41.99** per pieceTrade price from **\$38.63**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Amalgam Carrier Single Ended 1.5mm Mini Tapered Tip 6" (15.2cm) Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument for precise and controlled amalgam filling application, available in single-ended or double-ended designs. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. The textured surface helps effectively grip the filling material, preventing it from slipping.

Amalgam filling carriers like a syringes, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental

cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

**SINGLE ENDED AMALGAM  
CARRIER FEATURES**

- 6" (15.2cm) Overall Length.
- Single Ended.
- 1.5mm Mini Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.



**WI-70000526****6 Pcs Amalgam Carrier Single Ended Set  
2.0mm Regular Tapered Tip****\$41.99** per pieceTrade price from **\$38.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Amalgam Carrier Single Ended Set 2.0mm Regular Tapered Tip 6" (15.2cm) Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument for precise and controlled amalgam filling application, available in single-ended or double-ended designs. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. The textured surface helps effectively grip the filling material, preventing it from slipping.

Amalgam filling carriers like a syringes, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental

cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

**SINGLE ENDED AMALGAM CARRIER FEATURES**

- 6" (15.2cm) Overall Length.
- Single Ended.
- 2.0mm Regular Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.

**WI-70000523****6 Pcs Amalgam Carrier Single Ended Set  
3.0mm Large Tapered Tip****\$41.99** per pieceTrade price from **\$38.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Amalgam Carrier Single Ended Set 3.0mm Large Tapered Tip 6" (15.2cm) Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument for precise and controlled amalgam filling application, available in single-ended or double-ended designs. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. The textured surface helps effectively grip the filling material, preventing it from slipping.

Amalgam filling carriers like a syringes, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental

cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

**SINGLE ENDED AMALGAM CARRIER FEATURES**

- 6" (15.2cm) Overall Length.
- Single Ended.
- 3.0mm Large Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.



**WI-70004173****6 Pcs American Extracting Forceps #90 Upper Left Molars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Dental Extracting Forceps #90 for Upper Molars with Left Pointed Tip Design German Grade Stainless Steel American Pattern Teeth Extraction Dental Surgical Instrument.

American pattern dental extraction forceps No. 90 used to extract upper left molars. This extraction forceps has left pointed tip design. American pattern forceps have a smooth hinge so no pin is visible. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**EXTRACTING FORCEPS #90 FEATURES**

- For Upper Left Molars.

- Left Pointed Tip Design.
- Non-Slip Grip Handle.
- American Pattern.
- Made With Stainless Steel

**PACKAGE INCLUDE:**

- 6 Pcs American Pattern Extracting Forceps #90 – Upper Left Molars

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

**APPLICATION:** Extraction Forceps**BEAK/JAWS:** Left Pointed Tip**POLISHED:** Yes**HANDLE:** Non-Slip Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**WI-70004182****6 Pcs American Pedo Extracting Forceps #16S Lower Molars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pedo Extracting Forceps #16S for Lower Molars Cowhorn Beak Design American Pattern German Grade Stainless Steel Childern Teeth Extraction Surgical Dental Instruments

Cowhorn Pedo Extraction Forceps No. 16S are the perfect tool to remove children teeth from the alveolar bone. American extracting forceps #16S has cowhorn pointed beaks designed for lower first and second molars. Hooked handle provides increased leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you. 100% satisfaction guarantee against defects in material and workmanship.

**PEDO EXTRACTING FORCEPS #16S FEATURES**

- Cow Horn Beak Design.
- For Lower 1st & 2nd Molars.
- Textured-grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel

**PACKAGE INCLUDE**

- 6 Pcs American Pedo Extraction Forceps #16S Cowhorn Beak with Hook Handle

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

**APPLICATION:** Teeth Extraction

Surgery Procedure

**BEAK/JAWS:** Cowhorn Beaks Design

**6  
PCS****WI-70004266****6 Pcs Ball Burnisher #28/29 Set Double Ended  
2.3mm – 3.6mm Football Tips****\$35.99** per pieceTrade price from **\$33.11**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Dental Ball Burnisher #28/29 Double Ended 2.3 mm – 3.6 mm Football Working Tips Premium Quality Stainless Steel Amalgam Restorative Composite/Plastic Filling Instruments.


28/29 football end burnisher is an operative hand instrument for dental professionals used to condense, carve and polish restorative materials such as amalgam. Ball Burnishers can also be used to form occlusal anatomy. The working end of a #28/29 ball burnisher is smooth and rounded like football (2.3 mm – 3.6 mm) made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. Fully autoclavable for reuse purposes.

**BALL BURNISHER # 28/29  
FEATURES**

- Crafted with 6.5" (16.5cm) Overall Length.
- 2.3mm Football & 3.6mm Football.
- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Reusable.
- Smooth Finish For Excellent Tactile Feel

**PACKAGE INCLUDE:**

- 6 Pcs Foot Ball Burnisher #28/29 Double Ended

**6  
PCS****WI-70004087****6 Pcs Bone File Set #12A Double Ended****\$43.99** per pieceTrade price from **\$40.47**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Bone File #12A Double Ended 7" (17.8cm) 5mm Cross Serrated and Plain Serrated Edge with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

12A Bone File designed to smooth and shape bone using a pull stroke. Their cross-serrated & serrated working ends make it easy to maneuver between soft tissue and bone. Using a surgical bone file to smooth the bone and narrow the canals of the interior of the tooth helps to keep the area disinfected. Round non-slip knurled handle allow strong grip while handling. Hand crafted from premium quality stainless steel, fully autoclaved and sterilized for reuse purpose.

- Crafted with 7" (17.8cm) Overall Length.
- Double-ended.
- One Cross-serrated & One Plain Serrated Ends.
- Non-slip Knurled Grip Round Handle.
- Reusable.
- Hand Crafted Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 6 Pcs Dental Bone File #12A Double Ended 7" (17.8cm)

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Bone Files**APPLICATION:** Used for Smoothing and Shaping of Bones**CURVATURE:** Curved**RESUABLE:** Yes**12A BONE FILE FEATURES**

**WI-70001467****6 Pcs Castroviejo Caliper Long 0-40mm Angled 6"****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Castroviejo Implant Bone Caliper 6" (15.2cm) Long Model – 30° Angled Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper Angled 15.2cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 6 inches and tapers from the measurement end to the tip. The shanks are straight with 30° Angled jaws

**LONG CASTROVIEJO CALIPER FEATURES**

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 15.2cm
- 30° Angled Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

**PACKAGE INCLUDE:**

- 6 Pcs Castroviejo Caliper 6" (15.2cm) Long Model – 30° Angled Tip Graduated 0-40mm

**TECHNICAL INFORMATION****WI-70001458****6 Pcs Castroviejo Caliper Long 0-40mm Angled 7"****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Castroviejo Implant Bone Caliper 7" (17.8cm) Long Model – 30° Angled Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper Angled 17.8cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 7 inches and tapers from the measurement end to the tip. The shanks are straight with angled jaws.

**LONG CASTROVIEJO CALIPER FEATURES**

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 17.8cm
- 30° Angled Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

**PACKAGE INCLUDE:**

- 6 Pcs Castroviejo Caliper 7" (17.8cm) Long Model – 30° Angled Tip Graduated 0-40mm

**TECHNICAL INFORMATION**

**WI-70003951****6 Pcs Castroviejo Caliper Long 0-40mm Straight 6"****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Castroviejo Implant Bone Caliper, Long Model 6" (15.2cm) Straight Tips, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Castroviejo Bone Caliper Straight 15.2cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 6 inches and tapers from the measurement end to the tip. The shanks and jaws are straight.

**LONG CASTROVIEJO CALIPER FEATURES**

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 15.2cm
- Pointed tips that open up for precise measurement.
- Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

**PACKAGE INCLUDE:**

- 6 Pcs Castroviejo Caliper Long 6" (15.2cm) Straight Graduated 0-40mm

**TECHNICAL INFORMATION****WI-70001462****6 Pcs Castroviejo Caliper Long 0-40mm Straight 7"****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Castroviejo Implant Bone Caliper 7" (17.8cm) Long Model – Straight Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper Straight 17.8cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 7 inches and tapers from the measurement end to the tip. The shanks and jaws are straight.

**LONG CASTROVIEJO CALIPER FEATURES**

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 17.8cm
- Straight Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

**PACKAGE INCLUDE:**

- 6 Pcs Castroviejo Caliper 7" (17.8cm) Long Model – Straight Tip Graduated 0-40mm

**TECHNICAL INFORMATION**

**WI-70003996****6 Pcs Cement Spatulas 7" (17.8cm) Double Ended****\$45.99** per pieceTrade price from **\$42.31**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Cement Spatulas Double Ended 7" (17.8cm) Lab Mixing Modeling Premium Quality Stainless Steel Dental Restorative Instruments.

Cement spatulas used for mixing cement materials and various other dental materials like silicones, waxes, cements etc. Dental spatula comes in double ends, one tip is rounded for scooping, mixing and the other is arch headed spreading tip used to mix, spread and lift material including foods, drugs, plaster and paints.

Dental cement spatulas available in double and single ended patterns with flexible and rigid forms which are ideal for lighter, medium to heavier cements and other restorative materials. These spatulas can also be used to load materials into crowns or other

preparations. Made with premium quality stainless steel with superior craftsmanship. Provides a reliable performance for dental professionals over time. Fully Autoclavable for hygiene and reuse purpose.

**DOUBLE ENDED CEMENT SPATULAS FEATURES**

- Crafted With 7" (17.8cm) Overall Length.
- Double Ended.
- One Rounded End and Other is Arch Headed.
- Optimal Size for Handling for Optimal Maneuvering.
- Non-slip Textured Area for Better Grip
- Highly Polished Rounded Handle
- Fully Autoclave and Reusable

**WI-70003992****6 Pcs Cement Spatulas Single Ended****\$43.99** per pieceTrade price from **\$40.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 6 Pcs Dental Cement Spatulas Single Ended 7" (17.8cm) Lab Mixing Modeling Stainless Steel Restorative Instruments.

Cement spatulas are used for mixing cement materials and various other dental materials like silicones, waxes, cements etc. A spatula is a broad, flat, flexible blade used to mix, spread and lift material including foods, drugs, plaster and paints. Available in double and single ended patterns with flexible and rigid forms which are ideal for lighter, medium to heavier cements and other restorative materials. Cement spatulas can also be used to load materials into crowns or other preparations. Made with premium quality stainless steel with superior craftsmanship. Cement spatulas can be

reused over time as they will not bend or tarnish on you. We offer 100% satisfaction guaranteed.

**SINGLE ENDED CEMENT SPATULAS FEATURES**

- Crafted With 7" (17.8cm) Overall Length.
- Single Ended.
- Large Width Spatula Design.
- Optimal Size for Handling for Optimal Maneuvering
- Non-slip Textured Area for Better Grip
- Highly Polished Rounded Handle
- Fully Autoclave and Reusable
- Made With High Quality Stainless Steel

**PACKAGE INCLUDE**

**WI-70004086****6 Pcs Dental Bone File #12CA Double Ended****\$43.99** per pieceTrade price from **\$40.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Bone File #12CA Double Ended 7" (17.8cm) Cross Serrated Working Tips with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Dental Instruments.

12CA Bone File is a cross serrated working ends dental instrument used to excavate small spicules of bones from narrow spaces using a pull stroke. This bone file usually maneuvers between the soft tissues and the bones that need scraping. The knurled design allows strong grip and easy maneuvering. Made with premium quality stainless instrument for durability and long-lasting performance efficiency. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for

a variety of precise neurosurgical procedures.

**BONE FILE #12CA FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- Cross-serrated Ends.
- Double-Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 6 Pcs Bone File #12CA 7" Double Ended

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Surgical Dental Instruments**WI-70004085****6 Pcs Dental Bone File #33 Double Ended****\$43.99** per pieceTrade price from **\$40.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Dental Bone File #33 Double Ended 7" (17.8cm) Long with 4.8mm Narrow Pattern Straight Plain Serrated Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

#33 Bone File is designed to smooth and shape bone using a pull stroke. Bone File 33 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. The straight serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 4.8mm wide. Constructed from surgical-grade stainless steel for excellent performance and longevity. The round

handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

**BONE FILE #33 FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- 4.8mm Narrow Pattern.
- Double Ended.
- Straight Parallel Sides Plain Serration.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 6 Pcs Bone File #33 Straight 4.8mm Double Double Ends

**6  
PCS****WI-70004193****6 Pcs Dental Extracting Forceps #210H Upper 3rd Molars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Extracting Forceps #210H for Upper 3rd Molars with Hook Handle Premium Quality Stainless Steel American Pattern Bayonet Style Serrated Beaks Teeth Extraction Surgical Dental Instruments.

Extracting Forceps No. 210H with hook handle is used to extract teeth from the alveolar bone. American extraction forceps #210H features offset (Bayonet) conical beaks designed to adapt to upper third molars with longitudinal serrations to improve grasp. Hooked handle provides improved leverage when approach is from behind the head. American pattern forceps have a smooth hinge so no pin is visible. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**AMERICAN PATTERN  
EXTRACTION FORCEPS #210H  
FEATURES**

- Offset (Bayonet) Conical Beaks Design.
- For Upper 3rd Molars.
- Longitudinal Serrations Improve Grasp.
- Non-Slip Textured-Grip Hook Handle.
- American Pattern.
- Made With Premium Quality Stainless Steel.

**PACKAGE INCLUDES**

- 6 Pcs American Pattern Extraction Forceps #210H with Hook Handle

**TECHNICAL INFORMATION****6  
PCS****WI-70004199****6 Pcs Dental Extracting Forceps #23 Cow Horn for Lower Molar****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Universal Extracting Forceps #23 Cow Horn 1st & 2nd Lower Molars with Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#23 extraction forceps is used to extract teeth from 1st and 2nd lower molars with cow horn jaws. Forceps dimension is 7.5" (19cm) overall in length with Textured grip handles. Durable and sturdy our extraction forceps are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #23  
FEATURES**

- Slightly Larger Cowhorn Beak Design Than #16
- For Lower 1st & 2nd Molars
- Non-slip Grip Handle
- German Stainless Steel

**PACKAGE INCLUDE**

- 6 Pcs Extracting Forceps #23 – Cow Horn Jaws

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated

**WI-70004210****6 Pcs Dental Extracting Forceps #32****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Universal Extracting Forceps #32 for Upper Canines, Premolars & Molars with Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#32 Universal forceps used to extract teeth from Upper Cuspids, Bicuspids and Molars with sharp edged narrow beaks for breaking the alveolar process and trimming it after extraction.

Forceps dimension is 7.5" (19cm) overall in length with textured grip handles. Durable and sturdy our extraction forceps are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #32 FEATURES**

- Narrow beak, beveled and smoothed edges
- For Upper Canines, Premolars & Molars
- Non-slip Grip Handle
- German Stainless Steel

**PACKAGE INCLUDE**

- 6 Pcs Extracting Forceps #32 – Narrow, Bayonet Pattern Beak

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**WI-70004172****6 Pcs Dental Extracting Forceps #53L For Upper 1st & 2nd Left Molars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Extracting Forceps #53L for Upper 1st and 2nd Left Molars American Pattern German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

53L dental forceps are used to extract maxillary 1st and 2nd left molars. Left angled pointed offset (Bayonet) beaks for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

- Left Pointed Tip Design
- For Upper 1st & 2nd Molars
- Non-Slip Textured-Grip Handle
- German Stainless Steel

**PACKAGE INCLUDE**

- 6 Pcs American Pattern Dental Extracting Forceps #53L

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**EXTRACTING FORCEPS #53L FEATURES**

**WI-70004208****6 Pcs Dental Extracting Forceps #32A Molars Root – Slim Beak****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Dental Extracting Forceps #32A for Upper Cuspids, Biscuspids, Molars and Roots with Grip Handle German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments

#32A Extraction forceps used to extract teeth from Upper Cuspids, Biscuspids, Molars and roots. Sharp edged Slim beaks for breaking the alveolar process and trimming it after extraction. Forceps dimension is 7.5" (19cm) overall in length with textured grip handles. Durable and sturdy our extraction forceps are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #32A FEATURES**

- Slim, bayonet pattern beaks with parallel lines for optimal grip.
- For Upper Cuspids, Biscuspids, Molars and Roots
- Non-slip Grip Handle
- German Stainless Steel

**PACKAGE INCLUDE**

- 6 Pcs Extracting Forceps #32A – Slim, Bayonet Pattern Beak

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**WI-70004207****6 Pcs Dental Extracting Forceps Set #65 for Upper Roots****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 American Pattern Dental Extraction Forceps #65 for Upper Roots, Fragments & Incisors, Slim Bayonet Pattern Jaws with Non-slip Handle German Grade Stainless Steel Teeth Extraction Dental Surgical Instrument.

American pattern dental extracting forceps No. 65 used to extract teeth from upper roots, fragments & incisors. Slim bayonet pattern beaks designed specially for hard to reach area.

American pattern forceps have a smooth hinge so no pin is visible. The top quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

**EXTRACTING FORCEPS #65 FEATURES**


- For Upper Roots, Fragments & Incisors.
- Slim, Bayonet Pattern Beaks.
- Non-Slip Textured-grip Handle.
- Made With High Quality Stainless Steel.

**PACKAGE INCLUDE:**

- 6 Pcs American Pattern Extracting Forceps #65

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps

**APPLICATION:** Extraction Forceps**BEAK/JAWS:** Slim Bayonet Jaws**POLISHED:** Yes**HANDLE:** Non-slip Handle

**WI-70003926****6 Pcs Dental Hollenback Wax Carver #3****\$35.99** per pieceTrade price from **\$33.11**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Hollenback Carver #3 Double Ended Spear Shape Sharp Blades with 6.5" Overall Length Premium Quality Stainless Steel Wax Modelling Dental Restorative Instruments.

3 hollenback carver is a dental instrument used for placing, carving, and contouring restorative materials like amalgam and composite. Hollenback #3 used for shaping restorative materials to resemble natural tooth anatomy, placing and packing materials into cavities, and trimming excess material after restoration. Double-ended design with large and sharp blades enable detailed work on both the chewing surface and sides of teeth, even in tight spaces. Stainless steel construction of hollenback carver 3 guarantees


durability and dependability, allowing for reuse without bending or tarnishing, and meets international quality standards.

**DENTAL HOLLENBACK CARVER #3 FEATURES**

- Large Pointed Spear Shape Sharp Blades.
- Double Ended.
- Non-Slip Knurled Grip Handle.
- Used for Wax and Porcelain Carving.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 6 Pcs Hollenback #3 Double Ended

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**WI-70003974****6 Pcs Dental Orthodontic Boone Bracket Positioning Gauge Set****\$41.99** per pieceTrade price from **\$38.63**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Pack of 6 Pcs Boone Bracket Positioning Gauge with Four Integrated Measurement 3.5, 4.0, 4.5 and 5.0 Millimeter – Premium Quality Stainless Steel Dental Orthodontic Instruments for Orthodontists, Dentists, and Dental Assistants involved in Orthodontic Procedures.

Boone Bracket Positioning Gauge (Star-Like Gauge) is an orthodontic instrument used for precisely measuring and marking the vertical height for braces placement on teeth. It includes multiple tips sized at 3.5mm, 4.0mm, 4.5mm, and 5.0mm, aligning with standard distances from a tooth's biting edge, ensuring consistent bracket bonding during procedures. It also features pointed tips for accurate surface marking on the tooth crown.

Boone gauge facilitates consistent bracket placement on the incisal or occlusal surfaces and promotes uniform bonding heights across the dental arch, aiding proper tooth movement during treatment. Made with premium quality stainless steel ensures sterility for hygienic orthodontic practices.

**BOONE BRACKET POSITIONING GAUGE FEATURES**

- Four Measurement Options: 3.5, 4.0, 4.5, and 5.0 Millimeters.
- Pointed Tips with Flat Base Portion for Superficial Marking
- Measure Incisal or Occlusal Surface of the Tooth to Determine the Ideal Bracket Position.
- Withstand Repeated Use and Sterilization Procedures.



**WI-70004203****6 Pcs Dental Tooth Extracting Forceps #10S with Serrated Jaws****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Heavy Duty American Pattern Tooth Extracting Forceps Fig. 10S Upper Molar Serrated Curved Jaws German Grade Stainless Steel Dental Instrument.

Durable and sturdy, Upper molars extracting forceps #10S are the perfect tool to remove teeth from the alveolar bone. American pattern forceps have a smooth hinge so no pin is visible. Our tools are professional quality and used by professionals dentists.

**EXTRACTING FORCEPS #10S FEATURES**

- Serrated Beak For Better Grip
- Figure 10S For Upper Molar
- Smoothly Curved Tips Making Them Suitable For Both Sides
- Non-Slip Grip Handle

- German Stainless Steel

**PACKAGE INCLUDE**

- 6 Pcs Extracting Forceps #10S Curved Serrated – Upper Molars

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**WI-70004198****6 Pcs Extracting Forceps #53R for Upper 1st & 2nd Right Molars****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Premium Quality American Pattern Dental Extracting Forceps #53R for Upper 1st and 2nd Right Molars German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.


53R dental forceps are used to extract maxillary 1st and 2nd right molars. Right angled pointed offset (Bayonet) beaks for engaging bifurcated buccal root. Deep cavity in beaks prevents crown from being crushed. American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All our dental tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #53R FEATURES**

- Right Pointed Tip Design.
- For Upper 1st & 2nd Molars.
- Angular Beaks.
- Non-Slip Textured-Grip Handle.
- German Stainless Steel.

**PACKAGE INCLUDE**

- 6 Pcs American Pattern Dental Extracting Forceps #53R

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Teeth Extraction**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION****PROCEDURE:** Passivated

**WI-70004177****6 Pcs Extracting Forceps #62 Upper and Lower Anterior Teeth****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Universal American Pattern Extracting Forceps #62 for Upper & Lower Incisors, Canines, Premolars and Deciduous Teeth American Pattern German Grade Stainless Steel Tooth Extraction Dental Surgical Instrument.

American Pattern Tooth Extraction Forceps No. 62 used to extract Upper and Lower Incisors & Canines. The long, extra sharp and very narrow beaks allow a tight grasp of the tooth. The slightly curved working end makes this instrument capable of handling both upper and lower teeth, in addition to minimizing trauma to the patient.

American pattern forceps have a smooth hinge so no pin is visible. German grade stainless steel promotes longevity and corrosion resistance. All

our dental tools are professional quality and used by professionals.

**EXTRACTING FORCEPS #62 FEATURES**

- Precise Tips Meet at a Point, Ensuring Optimal Grip.
- Universal, Conical Beaks Design.
- For Upper & Lower Incisors, Canines, Premolars and Deciduous Teeth.
- Gently-Curved Working End Making It Suitable For Both Upper & Lower Teeth.
- Non-Slip Textured-Grip Handle.
- German Stainless Steel.

**PACKAGE INCLUDE**

- 6 Pcs American Pattern Extracting Forceps #62

**TECHNICAL INFORMATION****WI-70004183****6 Pcs Extracting Forceps Set #151AS Serrated Split Beaks****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Dental Extracting Forceps #151AS Split Beaks, Serrated for Lower Anterior – Standard – American Pattern German Grade Stainless Steel Tooth Extraction Dental Surgical Instruments.

Extraction Forceps No. 151AS used for extracting teeth from the alveolar bone. Lower universal split serrated beaks designed to provide increased grasp beneath the crown to prevent crown or root breakage. Serrations provide increased mechanical grasp on the root without increased pressure. This 151AS dental forceps can be used left and right.

American pattern forceps have a smooth hinge so no pin is visible. The stainless steel construction ensures maximum durability and dependability.


You can reuse these 151AS extracting forceps over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

**AMERICAN EXTRACTING FORCEPS #151AS FEATURES**

- For Lower Anterior.
- Serrated, Split Beaks Provide an Optimal Grip.
- American Pattern.
- Non-slip Grip Handle.
- Made With Stainless Steel

**PACKAGE INCLUDE:**

- 6 Pcs Extracting Forceps #151AS Split Beaks Serrated -Lower Anterior

**TECHNICAL INFORMATION**

**WI-70003905****6 Pcs Lesmann Knife Set 17cm – Large****\$36.99** per pieceTrade price from **\$34.03**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Lesmann Knives Large 17cm Double Ended Blade with Wooden Handle Premium Quality Stainless Steel Dental Technician Laboratory Instrument.

Lessmann knife (Lessmann Knife) is a specialized dental laboratory tool used for carving, shaping, trimming, and modeling dental wax patterns for crowns, bridges, and dentures. Its common applications include sculpting wax-ups for dental work, refining contours of dentures, adding detailing to biting surfaces of model teeth, and using heat to manipulate wax.

Lessmann Knife features a double ended stainless steel blade ensuring sharpness for precise cuts with lightweight wooden handle for balanced handling and reduced hand

fatigue. Lessmann wax knife designed to withstand regular sterilization protocols.

**LESSMANN WAX KNIFE FEATURES:**

- Wax & Porcelain Carvers for Dental Lab Technicians.
- Overall Length 6.7" (17cm) – Large.
- Double-ended, One Knife End & Other is Spoon End.
- Smooth Edges Minimize the Risk of Hand Injury.
- Lightweight Wooden Handle for Easy Grip.
- Reusable.
- High Quality Stainless Steel Blades.

**PACKAGE INCLUDE**

- 6 Pcs Lesmann Knife 17cm Large – Double Ended

**WI-70004083****6 Pcs Miller Bone File #45 Double Ended****\$43.99** per pieceTrade price from **\$40.47**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Miller Bone File #45 Straight-Cut 7" (17.8cm) Long with 5mm Wide Straight & Curved Double Working Ends Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

Bone file #45 is designed to smooth and shape bone using a pull stroke. Dental bone file #45 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. It can scrape off bony segments from soft tissues and bones. The serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features straight-cut plain serrations and is 5mm wide. Constructed from surgical-grade

stainless steel for excellent performance and longevity. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

**BONE FILE MILLER #45 FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- 5mm Straight-cut Plain Serrated Ends.
- Curved & Straight Double Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 6 Pcs Miller Bone File #45 Double Ended 5mm Straight & Cruved Ends



**WI-70000538****6 Pcs Molt Bone Curette Set #2/4 Double Ended****\$44.99** per pieceTrade price from **\$41.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Molt Bone Curette Set #2/4 Double Ended 3.5mm & 7.0mm Working Tip Round Pattern Handle Premium Quality Stainless Steel Surgical Dental Instruments.

2/4 Molt surgical curette is a double ended bone curette designed for curettage, cyst removal, and tooth socket debridement. With a 3.5 mm and 7.0 mm working ends allows the surgeon to scrape away bone with durable and sharp edges. Molt #2/4 bone curette has round pattern non-slip knurled handle provides optimum control and omits hand fatigue.

This double ended bone curette can also use in a variety of veterinary dental procedures for efficient scraping of diseased tissue from oral/dental wounds for cats and dogs. Molt Bone

curette 2-4, crafted with high-quality surgical stainless steel with a smooth surface, ensures a strong grip. Completely sterilizable and reusable

**MOLT BONE CURETTE #2/4 FEATURES**

- Crafted with 6" (15.2cm) Overall Length.
- Double Ended.
- 3.5 mm – 7.0mm Spoon Sharp Blades.
- Ergonomic Design Round Pattern Grip Handle.
- Sterilizable / Reusable.
- Made with Premium Quality Stainless Steel

**PACKAGE INCLUDE**

- 6 Pcs Molt Bone Curette #2/4 Double Ended

**WI-70003898****6 Pcs Periosteal Elevator Baby Molt #7A Double Ended with Round Handle****\$35.99** per pieceTrade price from **\$33.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Periosteal Elevator Baby Molt #7A Double Ended 6.5" (16.5cm) Round Pattern Grip Handle Premium Quality Stainless Steel Dental Instruments.

Molt-7 elevators have thin blades on both ends that are used to cut periodontal ligaments and loosen teeth. Periosteal Baby Molt #7A comprises a flat, oval tip at each end, one blunt and one sharp. Surgeons use the blunt oval tip for lifting and retracting the periosteal flap, and the sharp tip to dissect it. In surgery dentist often needs to separate a bone or tooth from the periosteum that covers it. This is done with a periosteal elevator. The dentist may also use it to gain access to retained roots and surrounding bone. Elevators are manufactured from

stainless steel to ensure optimal performance.

**PERIOSTEAL ELEVATOR BABY MOLT #7A FEATURES**

- Overall Length 6.5" (16.5cm).
- Double Ended.
- One Flat and One Oval Tip.
- Non-Slip Knurled Grip Handle.
- Use to Elevate, Retract and Dissecting Periosteal Flap.
- Corrosion Resistance.
- Fully Autoclavable & Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 6 Pcs Periosteal Elevator Baby Molt #7A Double Ended

**TECHNICAL INFORMATION**

**WI-70000683****6 Pcs Polokoff Bone File #1 Double Ended****\$43.99** per pieceTrade price from **\$40.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Polokoff Bone File #1 Double Ended 7" (17.8cm) Long with 3.4mm Narrow Pattern Straight Serrated Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

#1 Polokoff Bone File is broadly used in surgical procedures to aid in mild bone reshaping using a pull stroke. Polokoff Bone File is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. The straight serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 3.4mm wide. Constructed from surgical-grade stainless steel for excellent



performance and longevity. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

**POLOKOFF BONE FILE #1 FEATURES**

- Crafted with 7" (17.8cm) Overall Length.
- 3.4mm Narrow Pattern.
- Double Ended.
- Straight Parallel Sides Serration.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

**PACKAGE INCLUDE**

- 6 Pcs Polokoff Bone File 1 Straight 3.4mm Narrow Pattern Working

**WI-70003950****6 Pcs Spring Divider Calipers Set 4.5" Dental Lab Tool****\$47.99** per pieceTrade price from **\$44.15**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Pack of 6 Pcs Spring Divider Calipers Set 4.5" (11.4cm) with Flat Legs Stainless Steel Precision Tool Bone for Accurate Measurements During Dental Lab Procedures.

The dental calipers measures the width of metal, wax, and crown size in dental cosmetic procedures. It has a scale or screw for the measuring of sizes in the dental lab. Spring Divider Caliper features two adjustable legs to measure hard-to-measure dimensions, made from carefully selected flat stock with heat treated points for durability, and they pivot on a roller while the bow spring ensures a precise tension setting. Accurate adjustments are made by turning the knurled nut on the threaded rod. The design of these calipers, with their fine adjustment

feature, ensures quick and positive settings when taking or transferring specific measurements. Manufactured from high quality stainless steel with a polished finish.

**SPRING DIVIDER CALIPER FEATURES**

- Crafted With 4.5" (11.4cm) Overall Length.
- Flat Legs with Sharpened Hard and Pointed Tips.
- Adjustable Lockable Steel Wing Divider.
- Rivet Design to Increase Stability.
- Finely Threaded to Achieve Fine Adjustment.
- Durable and Rustproof.
- Reusable.





WI-70002176

6 Pcs T/C Wynman Crown Remover Gripper Forceps 7" Curved Serrated Tips

\$69.99 per piece

Trade price from \$64.39

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs T/C Wynman Crown Remover Gripper Forceps Set 7" (17.78cm) Tungsten Carbide Tips with Gold Ring Handle Premium Quality Stainless Steel Temporary Crown Remover Surgical Dental Forceps Instruments.

Wynman Crown Remover Gripper Forceps features a serrated tip made from durable tungsten carbide for a long-lasting working end. Ideal to remove or place temporary crowns with ease, even in the most hard-to-reach places. Our TC Wynman crown grippers are hemostat style forceps that enable a strong grip on the temporary crown without slipe, keeping the crown intact even in the tightest areas of the mouth. Crafted from stainless steel and tungsten carbides,

this instrument offers seamless crown removal solution


TC WYNMAN CROWN GRIPPER FORCEPS FEATURES

6 PCS TC WYNMAN CROWN GRIPPER

**FORCEPS:** The TC Wynman Crown Gripper Forceps are an ideal addition to your inventory. The tungsten carbide inserts on the Dental Crown Gripper Forceps will grip the temporary crown firmly without slipping, keeping the crown intact. They are useful for gripping wet glass such as in porcelain restorations and assist in more precise procedures.

O.R. GRADE T/C AND STAINLESS

**STEELL:** Hot drop-forged construction Wynman crown gripper constructed of O.R. Grade stainless steel and tungsten carbide for the highest degree of



WI-34589234

6pcs Kit Dental Composite Resin Filling Spatulas Double Head

\$7.50 per piece

Trade price from \$6.50

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

6pcs Kit Dental Composite Resin Filling Spatulas Double Head is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.50, falling to \$6.50 at 1 piece — about 13% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

**WI-70002756****7 Pcs PDL Luxating Root Elevator Color Handle Set****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Luxating Elevator Set of 7 Pcs Implant PDL Luxating Root Elevator Proximators Tooth Extracting Surgical Dental Set Stainless Steel Color Handle German Dental Instruments.

**PDL LUXATING ELEVATOR SET FEATURES**

**DENTAL PROXIMATORS LUXATING ELEVATOR SET:** The 7pc Proximotor Luxating Elevators Set is used for effective dental extractions while reducing the amount of potential damage during a procedure. Choosing the right instrument is critical for the success of a dental procedure. With the number of tooth extractions increasing every year, it is important to have a dental tool that can provide a precise outcome.

**HIGH GRADE GERMAN STAINLESS STEEL:** Forged from Premium High-Grade O.R. Quality German Stainless Steel. Corrosion Resistant, Tarnish Free with a Smooth Finish. Ultrasonic Cleaned – Fully Autoclavable For Proper Disinfection

**ERGONOMICALLY DESIGNED HANDLE:**

The color coated + ergonomic handle design encourages proper directional force while allowing for full control during the apical progression of the sharp working edge. The solid, one-piece design provides a firm and balanced true feel during the procedure. The high grip surface provides minimal rotation for precise applications.

**PRECISE AND EASY TO USE:** The 7PC Dental Luxating Elevators Proximators Kit gives the best possible access to

**WI-70002723****7 Pcs Wax Carving Metal Clay Pottery Sculpting Tool Set****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

7 PCS Wax Carving Tools Dental Waxing Instruments, Stainless Steel Wax Carvers Set – Metal Clay Sculpting Pottery Sculpture Tools – Double Ended DIY Waxing Modeling Kits Steel Wax Carving Tool Dental Clay Pottery Plaster Sculpture Carvers Tools Kit Dabber Tool Set.

**7 PCS WAX CARVING TOOLS SET FEATURES:**

**WAX CARVERS KIT:** The Kit includes 7 Premium Quality commonly used wax carving instruments for a variety of different applications. You can cut, clean molds, engrave, make designs, repair, scrape, and shape ceramics, clay, electronics, jewelry, wood, or details of other items. Each tool included in this set is double-sided with different types of tips for creating a

wide variety of shapes and projects.

**DURABLE STAINLESS STEEL CONSTRUCTION:**

The Premium Quality stainless steel construction of our Wax Carvers Kit ensures long-lasting durability over many years of use. Our wax carvers kit is built to last, corrosion-resistant with a high mirror finish, and tarnish-free.

**PRECISE AND EASY TO USE:** Each double-sided wax carving tool is perfectly balanced and designed to be easy to maneuver for precise use. The stainless steel kit includes instruments with built-in finger grips for ultimate comfort and user-friendliness. The average length is 6", ideal for many uses. Great for carving, cutting, modeling, scraping, shaping, etc.

**EXCELLENT DENTAL WAX CARVERS KIT:** Perfect Dental Laboratory





WI-70004080

## 8 Pcs Bone File Set for Surgical Dental Orthopedic Veterinary

**\$65.99** per piece

Trade price from **\$60.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 8 Pcs Bone Files Set with Carrying Case Including #10, #33, #45, #64, #12, #12A, #12CA and #1PO – Premium Quality Stainless Steel 8 Dental Bone File Double Ended for Precise and Delicate Dental Surgical Clinical Procedures.

8-piece dental bone file set comprises surgical-grade stainless steel instruments used in oral, periodontal, and orthopedic surgeries. These bone files designed to smooth and shape alveolar bone, primarily used for post-extraction smoothing, alveoloplasty and flap surgery. Features a double-ended configuration for access to hard-to-reach areas and serrated cutting edges designed for efficient bone tissue removal with minimal trauma. The instruments come in varying sizes to

accommodate different surgical needs. Made from surgical-grade stainless steel ensures outstanding durability and performance over time. These dental surgical bone files can be reused after proper sterilization.

### DENTAL BONE FILES FEATURES

- Double Ended.
- Serrated Tips.
- Round Non-Slip Handles.
- Available in Different Shapes and Diameters.
- Ideal to Smooth and Shape Rough Edges Bone Spicules During Surgical Procedure.
- Fully Autoclavable/Reusable.
- Hand Crafted from Premium Quality Stainless Steel.

### PACKAGE INCLUDE

WI-32307453

## 8 Pcs Stainless Steel Color Dental Tool

**\$4.00** per piece

Trade price from **\$0.90**

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

8 Pcs Stainless Steel Color Dental Tool — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$0.90 at 10 pieces — about 78% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-24322645

## 8-Piece Precision Dental Instruments Endodontic Root Canal

**\$8.00** per piece

Trade price from **\$5.00**

Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

We manufacture 8-Piece Precision Dental Instruments Endodontic Root Canal for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house.

Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$8.00, falling to \$5.00 at 10 sets — about 38% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-70000328

## 9 Pcs Wax Carver Set – Double Ended

**\$16.99** per piece

Trade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

9 Pcs Wax Carver Set – Double Ended – Premium Quality Stainless Steel Wax Carving & Modelling Dental Lab Mixing Instruments Designed for Carving, Cutting, Modeling, and Shaping Wax With Precision.

Waqar International's offer Carver kit includes tools that are each carefully selected for their versatility across a wide range of applications. This wax carver kit is ideal for clay or wax and can also be used for dabber, jewelry, sculpting, carving, cutting, modeling, scraping, grinding, shaping and whipping, marking, carpentry, engraving, and styling. The double-headed design comes with different cutting tips, making them ideal for dental, crafts and DIY projects. All instruments are crafted from durable

stainless steel, ensuring their long-lasting use. Frosted handles provide a good grip and non-slip performance. They are comfortable and easy to use. Can be reuse after sterilization.

### 9 PCS WAX CARVER SET INCLUDES

- 1 Pc Ward 1 + Zahle 3 Tips
- 1 Pc Lecorn 5 Wax Carver
- 1 Pc Ward 1 + Zahle Tips
- 1 Pc Pk Thomas #2
- 1 Pc Hollenback + Ward Tips
- 1 Pc Ward + Zahle Tips
- 1 Pc Excavator + Zahle Tips
- 1 Pc Vehe + Spoon Spatula Tips
- 1 Pc Lecorn 5 + PK Thomas #2 Tips

### TECHNICAL INFORMATION





WI-36815967

### Adult Jennings Mouth Probe Gag Stainless Steel Oral

**\$20.00** per piece

Trade price from **\$12.00**  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Adult Jennings Mouth Probe Gag Stainless Steel Oral is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$20.00, falling to \$12.00 at 10 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35031143

### Ainsworth Rubber Dam Punch Stainless Steel Dental Instrument

**\$4.10** per piece

Trade price from **\$3.80**  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Ainsworth Rubber Dam Punch Stainless Steel Dental Instrument is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.10, falling to \$3.80 at 10 pieces — about 7% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Every piece is hand-finished and inspected before it leaves the bench:

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-30844027

## Amalgam Carriers Stainless Steel Single/Double Ended Syringe

**\$3.00** per piece

Trade price from **\$2.00**  
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

Amalgam Carriers Stainless Steel Single/Double Ended Syringe — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: two working ends on one handle, halving what you carry.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.00, falling to \$2.00 at 10 sets — about 33% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-32307440

## Basic Orthodontics Dental Instrument Set

**\$4.00** per piece

Trade price from **\$0.90**  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

Basic Orthodontics Dental Instrument Set is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$0.90 at 10 pieces — about 78% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-34359269

## Bone Scraper Harvesters Implant Surgery Set Grafting Dental Surgical

**\$120.00** per piece

Trade price from **\$70.00**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

This is our pattern of Bone Scraper Harvesters Implant Surgery Set Grafting Dental Surgical, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$120.00, falling to \$70.00 at 1 piece — about 42% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-70003964

## Bracket Positioning Star Boone & Height Gauge Wick Type Orthodontic Set

**\$14.99** per piece

Trade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Bracket Height Gauge Wick Type & Boone Positioning Star-Like Gauge Premium Quality Dental Tools for Professional Orthodontists and Clinical Use.

Wick Type Bracket Height Gauge and Star Boone Bracket Positioning Gauge Set enables orthodontists to accurately measure and position braces brackets at the correct vertical height on teeth. The Boone Gauge features a star-shaped design with calibrated measurement points (3.5mm, 4.0mm, 4.5mm, and 5.0mm) for gauging distance from the incisal or occlusal edge. The Wick Type Height Gauge utilizes a rod tip to verify the precise vertical slot position for .022" compatible brackets, suitable for both posterior and anterior teeth. This

combined purpose eliminates guesswork during bonding, ensuring uniform bracket placement on the clinical crown for effective tooth movement.

### PACKAGE INCLUDE

- 1 Pc Bracket Height Gauge 3.5mm – 5.0mm for 022 Brackets 6" (15.2cm)
- 1 Pc 1 Pc Star-Like Boone Bracket Positioning Gauge 3.5mm – 5.0mm

### TECHNICAL INFORMATION

**MATERIAL:** Stainless Steel, Aluminum  
**INSTRUMENTS TYPE:** Orthodontic Instruments

**APPLICATION:** Bracket Height Positioning Measurement

**AUTOCCLAVABLE/REUSABLE:** Yes

**STERILE:** Non-Sterile

**RUST PREVENTION PROCEDURE:**



WI-32031427

Complete Instrument Sets

**\$7.00** per piece

Trade price from **\$4.00**

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Complete Instrument Sets is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$4.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-24379958

Dental Autoclave Stainless Steel Cassette Tray Rack Box

**\$15.00** per piece

Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Autoclave Stainless Steel Cassette Tray Rack Box is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$15.00. Minimum order 2 pieces, and the price moves with quantity.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-24513101

### Dental Equipment Kit 10pcs Power Silicone Handle Composite

**\$25.00** per piece

Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Dental Equipment Kit 10pcs Power Silicone Handle Composite is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$25.00. Minimum order 2 sets, and the price moves with quantity.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-41184859

### Dental Hollowware Surgical Instruments

**\$15.00** per piece

Trade price from **\$12.00**  
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Dental Hollowware Surgical Instruments — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$15.00, falling to \$12.00 at 25 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-09615268

Dental Hygiene Kit Calculus Plaque Remover

\$3.50 per piece

Trade price from \$1.65  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Hygiene Kit Calculus Plaque Remover is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$3.50, falling to \$1.65 at 10 pieces — about 53% off the single-unit price. Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36904798

Dental Instruments Composite Resin Filling Spatula Filler Hygiene Kit

\$8.00 per piece

Trade price from \$5.00  
Minimum order 5 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Instruments Composite Resin Filling Spatula Filler Hygiene Kit is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$5.00 at 5 sets — about 38% off the single-unit price. Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-42839665

### Dental Micro Probe Scaler Explorer Set Stainless Steel Double

**\$30.00** per piece

Trade price from **\$17.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Dental Micro Probe Scaler Explorer Set Stainless Steel Double is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: a micro-scale working end for fine and delicate work.

Supplied as a set — the contents can be adjusted to your list. The steel is

heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$30.00, falling to \$17.00 at 1 set — about 43% off the single-unit price. Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-30655614

### Dental Trephine Drills Instrument Set Stainless Steel Color

**\$40.00** per piece

Trade price from **\$33.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Dental Trephine Drills Instrument Set Stainless Steel Color is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$40.00, falling to \$33.00 at 1 piece — about 18% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



**WI-32309377**

## Double Ended Stainless Steel Dental Scaler Instrument

**\$4.00** per piece

Trade price from **\$0.90**  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

We manufacture Double Ended Stainless Steel Dental Scaler Instrument for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: two working ends on one handle, halving what you carry.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$4.00, falling to \$0.90 at 10 pieces — about 78% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



**WI-30660169**

## Fork Pattern Articulator Nickel Chromium Plate

**\$18.00** per piece

Trade price from **\$12.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

Fork Pattern Articulator Nickel Chromium Plate is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$18.00, falling to \$12.00 at 1 piece — about 33% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-30598718

### German Stainless Steel Adjustable Dental Articulator Holder Instrument

**\$18.00** per piece

Trade price from **\$12.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

This is our pattern of German Stainless Steel Adjustable Dental Articulator Holder Instrument, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$18.00, falling to \$12.00 at 1 piece — about 33% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-34294513

### Implant Buddy Stainless Steel Dental Implant Case Kit

**\$120.00** per piece

Trade price from **\$70.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

This is our pattern of Implant Buddy Stainless Steel Dental Implant Case Kit, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$120.00, falling to \$70.00 at 1 piece — about 42% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36131731

## Maxillofacial Surgery Instruments Set Stainless Steel Oral

**\$190.00** per piece

Trade price from **\$170.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

This is our pattern of Maxillofacial Surgery Instruments Set Stainless Steel Oral, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$190.00, falling to \$170.00 at 1 set — about 11% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-34083176

## Medical Devices Orthopedic Reusable Dental Instruments Vet Instruments Kit

**\$600.00** per piece

Trade price from **\$597.00**  
Minimum order 5 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

This is our pattern of Medical Devices Orthopedic Reusable Dental Instruments Vet Instruments Kit, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Supplied as a set — the contents can be adjusted to your list. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$600.00, falling to \$597.00 at 5 sets — about 1% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-70002039

Morrell Crown Remover Set with 3 Tips

\$10.99 per piece

Trade price from \$10.11  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Morrell Crown Remover with 3 Point Attachment Tips Premium Quality Stainless Steel Dental Crowns and Bridge Removers Restorative Dentistry Instruments.

The Morrell Crown Remover is made of stainless steel and comes with 3 tip points. Classic, easy-to-use crown remover is perfect for removing temporary or permanent restorations. The crown remover tips attach to the working end of crown removers to grip the crown securely without slipping.

MORRELL CROWN REMOVER SET FEATURES

- For use in all areas of the mouth
- Manual, hammer mechanism
- Single Ended



- Comes with 3 tempered stainless steel tips
- Made with Highest quality medical grade Stainless Steel

PACKAGE INCLUDE:

- 1 Morrell Crown Remover Handle
- 3 Attachment Tips

TECHNICAL INFORMATION:

**MATERIAL:** Stainless Steel  
**INSTRUMENTS TYPE:** Restorative Dental Instruments  
**APPLICATION:** Permanent and Temporary Crowns and Bridges Removal  
**TIPS SHAPE:** Hook Shape / Curved /Angled  
**HANDLE:** Hand Held  
**REUSABLE:** Yes  
**AUTOCCLAVABLE:** Yes, Must Autoclave



WI-32283577

Most Selling Stainless Steel Dental Scaler Instrument

\$4.00 per piece

Trade price from \$0.90  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Most Selling Stainless Steel Dental Scaler Instrument, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$0.90 at 10 pieces — about 78% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

**WI-70003295****Oral Dental Morrel Crown Bridges Removal Kit with Holding Forcep****\$18.74** per pieceTrade price from **\$17.24**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Orthodontic Oral Dental Morrel Crown Bridges Removal Kit with 5 Attachment Tips & 1 Holding Forcep Dentistry Lab Restorative Dentistry Instruments.

7 Pcs Morrel Crown Bridges Removal Kit. The Morrel Crown Remover is made of stainless steel and comes with 5 tip points. Classic, easy-to-use crown remover is perfect for removing temporary or permanent restorations. The crown remover tips attach to the working end of crown removers to grip the crown securely without slipping.

- Comes with 5 tempered stainless steel tips
- Made with Highest quality medical grade Stainless Steel

**PACKAGE INCLUDE:**

- 1 Morrel Crown Remover
- 5 Attachment Tips
- 1 Hemostat Forcep

**TECHNICAL INFORMATION:****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Restorative Dental Instruments**APPLICATION:** Permanent and Temporary Crowns and Bridges Removal**TIPS SHAPE:** Hook Shape / Curved / Serrated**HANDLE:** Hand Held**REUSABLE:** Yes**MORREL CROWN REMOVER SET FEATURES**

- For use in all areas of the mouth
- Manual, hammer mechanism

**WI-46731692****Orthodontic Ligature Pen Dental Set****\$7.00** per pieceTrade price from **\$3.50**

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

This is our pattern of Orthodontic Ligature Pen Dental Set, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$3.50 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



|                  |                                                                                 |
|------------------|---------------------------------------------------------------------------------|
| <b>Material</b>  | Stainless steel — we state the grade on your quotation                          |
| <b>Reusable</b>  | Yes, autoclavable                                                               |
| <b>Finish</b>    | Satin, mirror or black — your choice at order                                   |
| <b>Marking</b>   | Your logo or part number can be applied — ask us what is possible for this item |
| <b>Lead time</b> | Typically 15–30 days depending on quantity and marking                          |



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-34057331

Orthopedic Dog Dental Instruments Reusable  
AC 110–220V Vet Instrument

**\$600.00** per piece

Trade price from **\$597.00**  
Minimum order 5 sets


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Orthopedic Dog Dental Instruments Reusable AC 110–220V Vet Instrument is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: full reusability — clean, autoclave and use again.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$600.00, falling to \$597.00 at 5 sets — about 1% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-34382078

Orthopedic Surgical Instruments Dental  
Instrument Set

**\$120.00** per piece

Trade price from **\$70.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Orthopedic Surgical Instruments Dental Instrument Set is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$120.00, falling to \$70.00 at 1 piece — about 42% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



WI-33797689

### Pack 3 Stainless Steel Dental Boxes Rectangular Small

**\$14.00** per piece

Trade price from **\$8.00**  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

This is our pattern of Pack 3 Stainless Steel Dental Boxes Rectangular Small, made for dental, hygiene and oral surgery work. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$14.00, falling to \$8.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-42923809

### Pack of 12 PCS Dental Teeth Inspection Mirrors

**\$9.00** per piece

Trade price from **\$2.50**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Pack of 12 PCS Dental Teeth Inspection Mirrors is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$9.00, falling to \$2.50 at 1 set — about 72% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-32303451

Periodontics Dental Scaler Instrument Set

\$4.00 per piece

Trade price from \$0.90  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Periodontics Dental Scaler Instrument Set — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.00, falling to \$0.90 at 1 piece — about 78% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-26215650

Ready Ship Micro Periodontal Oral Surgery Kit Surgical

\$80.00 per piece

Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Ready Ship Micro Periodontal Oral Surgery Kit Surgical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: a micro-scale working end for fine and delicate work.

Supplied as a set — the contents can be adjusted to your list. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$80.00.

Minimum order 1 set, and the price moves with quantity. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-46803605

## Ridge Splitting Chisel Dental Set for Oral Surgery Bone

**\$6.90** per piece

Trade price from **\$3.50**  
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

This is our pattern of Ridge Splitting Chisel Dental Set for Oral Surgery Bone, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$6.90, falling to \$3.50 at 50 pieces — about 49% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-42858541

## Rubber Dam Kit Dental Isolation Instrument Set

**\$8.00** per piece

Trade price from **\$4.10**  
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

Rubber Dam Kit Dental Isolation Instrument Set — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$4.10 at 10 sets — about 49% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-30479857

Rubber Dam Punch Board Practice Tool Set  
Dental

**\$28.00** per piece

Trade price from **\$25.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Rubber Dam Punch Board Practice Tool Set Dental is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$28.00, falling to \$25.00 at 1 piece — about 11% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35330785

Set of 10 Dental Spoon Excavators Blue  
Titanium Coated

**\$38.00** per piece

Trade price from **\$25.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Set of 10 Dental Spoon Excavators Blue Titanium Coated is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$38.00, falling to \$25.00 at 1 set — about 34% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-42919952

## Set of 10 PCS Hollow Handle Dental Mouth Mirrors

**\$38.00** per piece

Trade price from **\$22.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

Set of 10 PCS Hollow Handle Dental Mouth Mirrors is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. We work the instrument by

hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$38.00, falling to \$22.00 at 1 set — about 42% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-70003889

## Set of 4 Miller Dental Root Extracting East West Right Left Elevators

**\$29.99** per piece

Trade price from **\$27.59**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


4 Miller Dental Tooth Root Elevator #71+72+73+74 Angled & Curved Left Right with Rounded Baldes Dental Surgical Instruments. Perfect for pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligaments for easy extraction. The lightweight standard handle with no seams, which prevents leaking or tip from caving into the handle.

an elevator consist of three parts; Handle, Shank and a Blade which is the working end of the instrument. These instruments are used to exert immense force on the dental roots and the bones around them. The thin and sharp working end is used to sever the periodontal ligament so the tooth can be easily and swiftly extracted.

### MILLER DENTAL ROOT ELEVATORS SET INCLUDES

- 1 Miller Elevator #71 Angled, Left, Flat Rounded Blade
- 1 Miller Elevator #72 Angled, Right, Flat Rounded Blade
- 1 Miller Elevator #73 Curved, Left, Flat Rounded Blade
- 1 Miller Elevator #74 Curved, Right, Flat Rounded Blade

Elevators (also known as Luxators) are Oral Surgical Instruments used in dental extractions. They may be used to loosen teeth prior to forceps extraction, to remove roots or impacted teeth, when teeth are compromised and susceptible to fracture, or when they are malpositioned and cannot be reached with forceps. Most commonly,





WI-30564806

### Stainless Steel Dental Articulator

**\$18.00** per piece

Trade price from **\$12.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Dental Articulator for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$18.00, falling to \$12.00 at 1 piece — about 33% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-32279576

### Stainless Steel Dental Scaler Instrument Set Dentistry

**\$4.00** per piece

Trade price from **\$0.90**  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Stainless Steel Dental Scaler Instrument Set Dentistry — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.00, falling to \$0.90 at 10 pieces — about 78% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-32274735

## Stainless Steel Dental Scaler Set 6pc

**\$4.00** per piece

Trade price from **\$0.90**  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Dental Scaler Set 6pc for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 6-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$4.00, falling to \$0.90 at 10 pieces — about 78% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-32309354

## Stainless Steel Dental Scalers

**\$4.00** per piece

Trade price from **\$0.90**  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

Stainless Steel Dental Scalers is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$4.00, falling to \$0.90 at 10 pieces — about 78% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36839984

### Stainless Steel Dental Set Sterilization Tray for Medical

**\$12.00** per piece

Trade price from **\$6.00**  
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Stainless Steel Dental Set Sterilization Tray for Medical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$12.00, falling to \$6.00 at 25 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-30583793

### Stainless Steel Dental Trephine Drills Implant Kit

**\$40.00** per piece

Trade price from **\$33.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Stainless Steel Dental Trephine Drills Implant Kit — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$40.00, falling to \$33.00 at 1 piece — about 18% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-30622697

## Stainless Steel Dental Trephine Drills Kit Free Bur Block

**\$40.00** per piece

Trade price from **\$33.00**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Dental Trephine Drills Kit Free Bur Block, made for dental, hygiene and oral surgery work. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$40.00, falling to \$33.00 at 1 piece — about 18% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-27358659

## Stainless Steel Reusable Surgical Dental Suture

**\$12.00** per piece

Trade price from **\$10.00**

Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


### ABOUT THIS INSTRUMENT

Stainless Steel Reusable Surgical Dental Suture is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$12.00, falling to \$10.00 at 1 set — about 17% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-25115209

## Stainless Steel Wax Carver Set Art Supplies Kit Pottery

**\$4.00** per piece

Trade price from **\$2.80**

Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

Stainless Steel Wax Carver Set Art Supplies Kit Pottery is part of our dental instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$4.00, falling to \$2.80 at 10 sets — about 30% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-32299493

## Steel Dental Tooth Cleaning Scalers

**\$4.00** per piece

Trade price from **\$0.90**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

Steel Dental Tooth Cleaning Scalers — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.00, falling to \$0.90 at 1 piece — about 78% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



# WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

## Request a quotation

Send us product codes from this catalogue, quantities and your country.  
We reply the same working day on WhatsApp.

|                    |                                                                           |
|--------------------|---------------------------------------------------------------------------|
| <b>Online shop</b> | shop.waqarintl.com                                                        |
| <b>Email</b>       | sales@waqarintl.com                                                       |
| <b>WhatsApp /</b>  |                                                                           |
| <b>Phone</b>       | +92 307 6173627                                                           |
| <b>Works</b>       | Roras Road, Opposite Anwar Khawaja<br>Industries, Sialkot 51310, Pakistan |

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.