



WAQAR INTERNATIONAL


PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

Specialty Instruments

Dental Instruments

Reusable stainless-steel specialty instruments forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



230

PRODUCTS

\$0.04

PRICES FROM

OEM

PRIVATE LABEL

Global

WORLDWIDE
EXPORT

Contents

230 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

Product	Code	Price	Page
10 Acorn Burnisher 21 Double Ended with Hollow Handle	WI-70002091	\$59.99	8
12 Amalgam Carrier Single Ended 1.5mm Mini Surgical Dental Instrument	WI-70004280	\$71.99	8
12 Ball Burnisher #27/29 De Dental Amalgam Instruments	WI-70004279	\$59.99	9
12 Ball Burnisher #34 Single Ended Dental Amalgam Instrument	WI-70004277	\$59.99	9
12 Beerendonk Caliper 0-80mm Measuring Curved Pointed Tips	WI-70003958	\$119.99	10
12 Boley Gauge Dental Implant Lab Orthodontic Instruments	WI-70003985	\$89.99	10
12 Bone File #12 Surgical Dental Dentist	WI-70004105	\$69.99	11
12 Bracket Height Gauge 3.5mm-5mm for .022" Brackets	WI-70003967	\$89.99	11
12 Castroviejo Caliper Long 0-40mm Straight 6"	WI-70003957	\$155.99	12
12 Crown Gauge Iwanson Spring Caliper 4" for Metal Improved Model	WI-70003984	\$59.99	12
12 Crown Gauge Iwanson Spring Caliper 4" for Wax Improved Model	WI-70001450	\$59.99	13
12 Dental Bone File #10 Double Ended 7" Curved Up & Straight	WI-70004106	\$69.99	13
12 Hollenback Wax Carver #3 Double Ended Hollow Handle	WI-70001082	\$89.99	14
12 Lecron Inlay Carver 6.5" Wax Carving Dental Tool	WI-70003938	\$61.99	14
12 Miller Bone File #64 Double Ended	WI-70004101	\$69.99	15
12 P.K Thomas #3 Dental Wax Modelling Carver	WI-70003936	\$49.99	15
12 P.K Thomas #4 Dental Wax Modelling Carver	WI-70003935	\$49.99	16
12 P.K Thomas #5 Dental Wax Modelling Carver	WI-70003934	\$49.99	16
12 Palti Bone Graft Scoop Dental Implant	WI-70004100	\$92.71	17
12 Ridge Mapping Caliper Surgical Dental	WI-70003982	\$154.18	17
2 Acorn Burnisher 21 Double Ended with Hollow Handle	WI-70002094	\$13.99	18
2 Automatic Crown Remover Gun Type with 5 Attachments	WI-70002032	\$36.99	18
2 Castroviejo Caliper Long 6" Straight + Angled Graduated 0-40mm	WI-70001465	\$32.99	19
2 Castroviejo Caliper Long 7" Straight + Angled Graduated 0-40mm	WI-70001456	\$35.99	19
2 Crown Remover Manual with 3 Tips Double Ended	WI-70002037	\$43.74	20
2 Krekeler Sliding Calipers Round Gauge 5" 0-80mm Straight Curved	WI-70000864	\$35.99	20
3 Acorn Burnisher 21 Double Ended with Hollow Handle	WI-70002093	\$19.99	21
3 Amalgam Carriers Single Ended 1.5mm, 2.0mm, 3.0mm Tip	WI-70000522	\$24.99	21
3 Ball Burnisher #27/29 Dental Amalgam Instruments	WI-70004273	\$19.99	22
3 Ball Burnisher #1 Double Ended Dental Amalgam Instruments	WI-70004274	\$15.99	22
3 Ball Burnisher #34 Single Ended Dental Amalgam Instrument	WI-70004271	\$20.99	23
3 Beerendonk Caliper 0-80mm Measuring Curved Pointed Tips	WI-70003955	\$35.99	23
3 Bracket Positioning Height Gauge 3.5mm-5mm for .022" Brackets	WI-70003966	\$29.99	24
3 Crown Gauge Iwanson Spring Caliper 4" for Metal Improved Model	WI-70003979	\$19.99	24
3 Crown Gauge Iwanson Spring Caliper 4" for Wax Improved Model	WI-70001452	\$19.99	25
3 Crown Remover Manual with 3 Tips Double Ended	WI-70002036	\$62.49	25
3 Dental Bone File #10 Double Ended 7" Curved Up & Straight	WI-70004099	\$23.99	26
3 Hollenback Wax Carver #3 Double Ended Hollow Handle	WI-70001084	\$25.99	26
3 Howard Bone File #12 Double Ended	WI-70004098	\$23.99	27
3 Lecron Carver Wax 6.5" Double Ended Carving Dental Tool	WI-70003931	\$20.99	27
3 Miller Bone File #64 Double Ended	WI-70004091	\$25.99	28

Contents (continued)

Product	Code	Price	Page
3 P.K Thomas #4 Dental Wax Modelling Carver	WI-70003928	\$20.99	29
3 P.K Thomas #5 Dental Wax Modelling Carver	WI-70003927	\$20.99	29
3 Palti Bone Graft Scoop Dental Implant	WI-70004090	\$30.84	30
3 Ridge Mapping Caliper Surgical Dental	WI-70003977	\$39.99	30
5 Acorn Burnisher 21 Double Ended with Hollow Handle	WI-70002092	\$29.99	31
5 Automatic Crown Remover Gun Type with 5 Attachments	WI-70002031	\$79.99	31
5 Bone Ridge Mapping Calipers 6.5" Ring Style Measures Upto 0-35mm	WI-70003946	\$89.99	32
5 Crown Remover Manual with 3 Tips Double Ended	WI-70002035	\$99.99	32
50 Ball Burnisher #27/29 De Dental Amalgam Instruments	WI-70002331	\$124.99	33
50 Lecron Inlay Wax & Porcelain Carver 6.5" Double Ended	WI-70003917	\$208.99	33
6 Ball Burnisher #1 Double Ended Dental Amalgam Instruments	WI-70004268	\$29.99	34
6 Ball Burnisher #27/29 De Dental Amalgam Instruments	WI-70004267	\$35.99	34
6 Ball Burnisher #34 Single Ended Dental Amalgam Instrument	WI-70004265	\$35.99	35
6 Beerendonk Caliper 0-80mm Measuring Curved Pointed Tips	WI-70003952	\$65.99	35
6 Boley Gauge Dental Implant Lab Orthodontic Instruments	WI-70003976	\$47.99	36
6 Bracket Positioning Height Gauge 3.5mm-5mm for .022" Brackets	WI-70003965	\$49.99	36
6 Crown Gauge Iwanson Spring Caliper 4" for Metal Improved Model	WI-70003975	\$35.99	37
6 Crown Gauge Iwanson Spring Caliper 4" for Wax Improved Model	WI-70001451	\$35.99	37
6 Dental Bone File #10 Double Ended 7" Curved Up & Straight	WI-70004089	\$43.99	38
6 Hollenback Wax Carver #3 Double Ended Hollow Handle	WI-70001083	\$46.99	38
6 Howard Bone File #12 Double Ended	WI-70004088	\$43.99	39
6 Lecron Inlay Wax & Porcelain Carver 6.5" Double Ended	WI-70003925	\$35.99	39
6 Miller Bone File #64 Double Ended	WI-70004082	\$46.30	40
6 P.K Thomas #3 Dental Wax Modelling Carver	WI-70003923	\$35.99	40
6 P.K Thomas #4 Dental Wax Modelling Carver	WI-70003922	\$35.99	41
6 P.K Thomas #5 Dental Wax Modelling Carver	WI-70003921	\$35.99	41
6 Palti Bone Graft Scoop Dental Implant	WI-70004081	\$51.46	42
6 Ridge Mapping Caliper Surgical Dental	WI-70003973	\$87.14	42
7 Inch Curved Surgical Dental Instrument	WI-35080040	\$3.50	43
Acorn Burnisher #21B Double Ended Hollow Handle	WI-70000060	\$9.99	43
Acorn Burnisher 21 Double Ended with Hollow Handle	WI-70002095	\$9.99	44
Amalgam Carrier Double Ended 8" (20.3cm) Jumbo/Jumbo Tips	WI-70000521	\$10.99	44
Amalgam Carrier Double Ended 8" (20.3cm) Large/Jumbo Tips	WI-70000518	\$10.99	45
Amalgam Carrier Double Ended 8" (20.3cm) Medium/Jumbo Tips	WI-70000519	\$10.99	45
Amalgam Carrier Double Ended 8" (20.3cm) Medium/Large Tips	WI-70000520	\$10.99	46
Amalgam Carrier Double Ended Large 2.8mm /Regular 2mm	WI-70004264	\$10.99	46
Amalgam Carrier Double Ended Large/Medium 2.5mm/2mm	WI-70004259	\$10.99	47
Amalgam Carrier Double Ended Mini 1.5mm / Regular 2mm	WI-70004254	\$10.99	47
Amalgam Carrier Single Ended 1.5mm Mini Tip	WI-70004263	\$8.99	48
Amalgam Carrier Single Ended 2.0mm Regular Tapered Tip	WI-70000528	\$8.99	48
Amalgam Carrier Single Ended 3.0mm Large Tapered Tip	WI-70000525	\$8.99	49
Automatic Crown Remover Gun Type with 5 Attachments	WI-70002033	\$21.99	49

Contents (continued)

Product	Code	Price	Page
Ball Burnisher #1 Double Ended Dental Amalgam Instruments	WI-70004282	\$7.99	50
Ball Burnisher #27/29 DE Dental Amalgam Instrument	WI-70004112	\$7.99	50
Ball Burnisher #28/29 Double Ended 2.3 mm – 3.6 mm Football Tips	WI-70004111	\$7.99	51
Ball Burnisher #34 Single Ended Dental Amalgam Instrument	WI-70004110	\$7.99	51
Ball Burnisher Double Ended 1mm-2mm Tips Hollow Handle	WI-70000330	\$9.99	52
Beerendonk Caliper 0-80mm Measuring Double Sided Scale	WI-70003949	\$12.99	52
Boley Gauge 10cm Measuring 6" (15.24cm)	WI-70003972	\$9.99	53
Bone File #12A Double Ended	WI-70004077	\$8.99	53
Bone Ridge Mapping Caliper 6.5" Ring Style with Gold Handle	WI-70003947	\$19.99	54
Bracket Height Gauge 2.0mm-3.5mm Adjustable Swivel Head	WI-70000545	\$12.99	54
Bracket Height Gauge 2.5mm-5.0mm Adjustable Swivel Head	WI-70000544	\$12.99	55
Bracket Height Gauge Moveable Head 4.0mm-6.0mm Pointed Tips	WI-70000543	\$12.99	55
Bracket Positioning Height Gauge 3.5mm-5mm for .022" Brackets	WI-70003963	\$10.99	56
Burnisher #21 Double Ended Solid Handle	WI-70000231	\$7.99	56
Castroviejo Caliper 3.5" (8.89cm) Angled Graduated 0-20mm	WI-70003945	\$10.99	57
Castroviejo Caliper 3.5" Graduated 0-40mm Angled	WI-70003942	\$12.99	57
Castroviejo Caliper 3.5" Graduated 0-40mm Straight	WI-70003941	\$12.99	58
Castroviejo Caliper 6" Round Tips Graduated 0-40mm	WI-70001455	\$14.99	58
Castroviejo Caliper Long 0-40mm Angled 6"	WI-70001469	\$16.99	59
Castroviejo Caliper Long 0-40mm Angled 7"	WI-70001460	\$18.99	59
Castroviejo Caliper Long 0-40mm Straight 6"	WI-70003948	\$16.99	60
Castroviejo Caliper Long 0-40mm Straight 7"	WI-70001464	\$18.99	60
Cement Spatulas 7" (17.8cm) Double Ended	WI-70003994	\$9.99	61
Cement Spatulas Single Ended	WI-70003991	\$8.99	61
Christensen Crown Remover CRCH2 90 Degree Angled	WI-70002029	\$9.99	62
Christensen Crown Remover Straight – CRCH1	WI-70002030	\$9.99	62
Cone Burnisher #21B Double Ended Solid Handle	WI-70004109	\$7.99	63
Crown Gauge Iwanson Spring Caliper 4" for Metal Improved Model	WI-70003971	\$8.99	63
Crown Gauge Iwanson Spring Caliper 4" for Wax Improved Model	WI-70001453	\$8.99	64
Dental Amalgam Well with Non-Slip Base	WI-70004255	\$14.99	64
Dental Autoclave Sterilization Cassette Tray	WI-25693507	\$14.00	65
Dental Bone File #10 Double Ended 7" Curved Up & Straight	WI-70004079	\$8.99	65
Dental Bone File #12 Double-Ended Medical Surgical Scaler	WI-21890483	\$1.50	66
Dental Bone File #12CA Cross Serrated Double Ended	WI-70004076	\$8.99	66
Dental Bone File #33 Double Ended	WI-70004069	\$8.99	67
Dental Bone Files	WI-70004072	\$15.36	67
Dental Bone Mill Bone Crusher Grinder Silver/Gold Implants Grafting Instrument	WI-70000692	\$39.99	68
Dental Cement Mixing Spatula Restorative Laboratory Instrument	WI-43373329	\$2.59	68
Dental Cement Spatula Double Ended	WI-43373326	\$2.59	69
Dental Cement Spatula Double Ended Stainless Steel Mixing Modeling	WI-43347847	\$2.59	69
Dental Edentulous Stainless Steel Perforated High Precision Denture	WI-29044198	\$4.00	70
Dental Flat Plastic Filling Instruments	WI-36933034	\$8.00	70

Contents (continued)


Product	Code	Price	Page
Dental Flex Periotomes Root Extraction Screw Extracting	WI-34335291	\$120.00	71
Dental Hollenback Wax Carver #3	WI-70003918	\$7.99	71
Dental Impression Tray	WI-06303870	\$4.50	72
Dental Impression Trays Autoclave	WI-29063159	\$4.00	72
Dental Impression Tray-Stainless Steel Perforated	WI-29022207	\$4.00	73
Dental Inspection Mirror Hollow Handle Oral Examination Tool	WI-42924596	\$5.00	73
Dental Inspection Mirror Hollow Handle Replaceable Mirror Instrument	WI-42931489	\$6.00	74
Dental Instruments	WI-35452239	\$9.00	74
Dental Laboratory Articulator Plain Line Ball Bearing Mortimer	WI-30574792	\$18.00	75
Dental Magnetic Articulator Adjustable Denture	WI-30655901	\$10.00	75
Dental Mouth Inspection Mirror Stainless Steel Oral Examination	WI-36338525	\$3.00	76
Dental Mouth Mirror with Private Labeling Light Weight	WI-13320359	\$2.50	76
Dental Mouth Teeth Inspection Mirror Reusable Dental	WI-42932414	\$3.50	77
Dental Scaler Tray Stainless Steel Autoclavable Instrument Sterilization Organizer	WI-36853954	\$1.00	77
Dental Scaler Tray Stainless Steel Sterilization Instrument Tray	WI-36877851	\$2.00	78
Dental Spoon Excavators Double Ended Blue Titanium	WI-35293845	\$38.00	78
Design Cassette IMS 1/2 DIN IME12DIN Orthopedic Dog	WI-34049315	\$600.00	79
Design Stainless Steel Dental Scaler	WI-10236086	\$2.50	79
Double-Sided Dental Mouth Examination Mirror	WI-10256215	\$2.50	80
First Class Dental Handpiece Crown Removal Tool	WI-35450235	\$9.00	80
Full Medical Stainless Steel Dental Implant	WI-34318383	\$120.00	81
German Stainless Steel Medical Instrument Tray Model	WI-28829372	\$140.00	81
German Stainless Steel NON-PERFORATED	WI-29036187	\$4.00	82
Graft Packer Large 6mm-8mm Serrated Tip Hollow Handle	WI-70000055	\$9.99	82
Graft Packer Small 3mm-4mm Serrated Tip Hollow Handle	WI-70000056	\$9.99	83
Grisdale Bone Condenser Plugger 3.3mm-4mm Double Ended	WI-70000324	\$9.99	83
High Efficient Steel/Metal Dental Mirror	WI-06348016	\$3.80	84
High Precision Dental Autoclave Impression Trays	WI-29011233	\$4.00	84
High Precision Dental Impression Perforated Large Trays Upper	WI-29081153	\$4.00	85
High Precision Dental Impression Trays Upper Lower Teeth	WI-29044192	\$4.00	85
High Precision Stainless Steel	WI-29057170	\$4.00	86
High Precision Stainless Steel Autoclave Impression Trays	WI-29047203	\$4.00	86
High Precision Stainless Steel Dental	WI-01090818	\$4.00	87
High Precision Stainless Steel Dental Impression Trays	WI-29057171	\$4.00	87
High Precision Stainless Steel Reusable Dental Orthodontic	WI-29035205	\$4.00	88
Hollenback Wax Carver #3 Blue Tips with Blue & Silver Hollow Handle	WI-70001043	\$10.99	88
Hollenback Wax Carver #3 Double Ended Hollow Handle	WI-70001085	\$9.99	89
Hollow Handle Dental Mouth Inspection Mirror Stainless Steel Hygiene	WI-42948395	\$5.00	89
Hollow Handle Dental Mouth Mirror with Replaceable Mirror	WI-42963158	\$4.00	90
Howard Bone File #12 Double Ended	WI-70004078	\$8.99	90
Krekeler Sliding Caliper Round Gauge 5" 0-80mm Curved	WI-70000866	\$19.99	91
Krekeler Sliding Caliper Round Gauge 5" 0-80mm Straight	WI-70000865	\$19.99	91

Contents (continued)

Product	Code	Price	Page
Lore Enterprises A+ Stainless Steel Adjustable Articulator	WI-30624645	\$18.00	92
Manual Crown Remover with 3 Tips Double Ended	WI-70002038	\$24.99	92
Marquette 0/1 Amalgam Plugger Double Ended	WI-70000678	\$7.99	93
Medical Dental Surgical Storage Tray for Lab Instruments	WI-32715105	\$14.00	93
Medical Surgical Dental Use Stainless Steel Instrument Storage Boxes	WI-25690521	\$14.00	94
Medical Surgical Tray Organizer Holder Stainless Steel Dental Instruments	WI-22428104	\$8.00	94
Medical Surgical Trays 201 Stainless Steel Dental Instruments	WI-33871516	\$14.00	95
Miller Bone File #45 Double Ended	WI-70004073	\$8.99	95
Miller Bone File #64 Double Ended	WI-70004068	\$9.65	96
Mini Handle Pakistan Made Dentist Mouth Mirror Bright	WI-13329342	\$2.50	96
Orthodontic Ligature Applicator Dental Instrument for Accurate Elastic Ligature	WI-46733664	\$6.80	97
P.K Thomas #1 Dental Wax Modeling Carver	WI-70003916	\$7.99	97
P.K Thomas #2 Dental Wax Modelling Carver	WI-70003915	\$7.99	98
P.K Thomas #3 Dental Wax Modelling Carver	WI-70003914	\$7.99	98
P.K Thomas #4 Dental Wax Modelling Carver	WI-70003913	\$7.99	99
P.K Thomas #5 Dental Wax Modelling Carver	WI-70003912	\$7.99	99
P.K. Thomas #3 Dental Wax Carver Gold Tip Hollow Handle	WI-70000059	\$10.99	100
P.K. Thomas #3 Dental Wax Carver Hollow Handle	WI-70000057	\$9.99	100
Pakistan Manufacture	WI-29068196	\$4.00	101
Palti Bone Graft Scoop Dental Implant	WI-70004067	\$13.10	101
Palti Bone Graft Scoop Dental Implant Instrument	WI-70000477	\$12.99	102
Periodontal Montana Jack Stainless Steel Dental Hygiene Instrument	WI-35310804	\$9.00	102
Polokoff Bone File #1 Double Ended	WI-70004066	\$8.99	103
Private Label Dental Crown Remover Attachments Medical	WI-35438360	\$9.00	103
Resin Dental Cement Spatula Stainless Steel Laboratory	WI-44644850	\$2.50	104
Reusable Hollow Handle Dental Mouth Mirror Stainless Steel Diagnostic	WI-42944448	\$5.00	104
Ridge Mapping Caliper Surgical Dental	WI-70003968	\$14.99	105
Semi-Adjustable Denture Magnetic Anatomic	WI-30639652	\$18.00	105
Single-Sided Dental Mouth Mirror 2023	WI-10239292	\$2.50	106
Size Kidney Shaped Emesis Basin Stainless Steel Surgical	WI-20375053	\$4.00	106
Spring Divider Caliper 4.5" Dental Lab Tool	WI-70003959	\$9.99	107
Stainless Steel Bone Scraper Surgical Instrument for Grafting	WI-36382864	\$4.00	107
Stainless Steel Dental Anterior Posterior Scaler Montana	WI-35323822	\$9.00	108
Stainless Steel Dental Cement Spatula Filling Instrument	WI-44664727	\$2.50	108
Stainless Steel Dental Examination Mirror	WI-36335510	\$3.00	109
Stainless Steel Dental Examination Mirror Diagnostic Octagon Handle	WI-25861757	\$1.50	109
Stainless Steel Dental Hygiene Scaler Tray Reusable	WI-36824967	\$2.00	110
Stainless Steel Dental Implant Mixing Bowl Size	WI-35102909	\$1.50	110
Stainless Steel Dental Impression Trays	WI-10245159	\$2.25	111
Stainless Steel Dental Instrument Sterilization Tray	WI-34971261	\$6.00	111
Stainless Steel Dental Instrument Tray Organizer for Clinic Hospital	WI-36867910	\$12.00	112
Stainless Steel Dental Instrument Trays Surgical Treatment Boxes	WI-27818100	\$8.00	112

Contents (continued)

Product	Code	Price	Page
Stainless Steel Dental Mouth Mirror Single Side	WI-34105560	\$0.04	113
Stainless Steel Dental Mouth Mirror Teeth	WI-42929441	\$3.00	113
Stainless Steel Dental Scaler-High-quality Surgical Instrument for General	WI-32309392	\$4.00	114
Stainless Steel Dental Surgical Storage Tray for Clinic	WI-40089398	\$14.00	114
Stainless Steel Instrument Tray for Medical Surgical Dental	WI-28839302	\$8.75	115
Stainless Steel Jacquette Scaler Dental Hygiene Periodontal	WI-43364645	\$1.15	115
Stainless Steel Medical Surgical Dental Lab Instruments Holder	WI-33820635	\$14.00	116
Stainless Steel Perforated Dental Surgical Instrument Box	WI-25649494	\$14.00	116
Stainless Steel Perforated Surgical Instrument Box for Dental	WI-27780137	\$8.00	117
Stainless Steel Single Ended Beauty Mixing Spatula	WI-36993453	\$1.00	117
Stainless Steel Surgical Dental Instruments Sterilization Box Cover	WI-27763216	\$8.00	118
Stainless Steel Surgical Instrument Box Perforated Dental Tool Box	WI-27742545	\$8.00	118
Stainless Steel Surgical Sterilization Trays Dental Instrument Boxes	WI-25638480	\$14.00	119
Stainless Steel Surgical Tray for Dental Lab Instruments Holder	WI-28798338	\$8.75	119
Star-Like Boone Bracket Positioning Gauge	WI-70003970	\$9.99	120
Surgical Tooth Crown Remover 3-Tip Bridge Dental Instrument	WI-29080566	\$9.00	120
Titanium Coated Hollenback Wax Carver #2 Double Ended	WI-70001081	\$10.99	121
VB Stainless Steel Dental Alloy VeraBond with Be/Dental Alloy	WI-34636056	\$2.00	121
Weiss Modified Castroviejo Caliper 6.75" Curved 0-40mm	WI-70001454	\$15.99	122
Weston #3 Wax Spatula – Double Ended	WI-70000329	\$9.99	122

**WI-70002091****10 Acorn Burnisher #21 Double Ended with Hollow Handle****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Acorn Burnisher #21 with Hollow Handle Dental Composite Restoration Instruments

Acorn burnisher #21 is a double-ended (2.7 mm – 1.7 mm) dental instrument used to condense, smooth, carve and polish amalgam. Acorn-shaped ends are excellent for carving occlusal anatomy. Acorn Burnishers used to assure that amalgam restorations are free of irregularities and are smooth

ACORN BURNISHER #21 FEATURES

- 1.7mm & 2.7mm acorn-shaped ends.
- Double-Ended Silver Color
- Hand crafted from High Quality stainless steel.
- Available in hollow, ergonomic handle.

- 10 Pcs Acorn Burnisher #21 Double Ended with Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Burnisher**APPLICATION:** Composite Placement and Contouring**COLOR:** Silver**TIP:** Double Ended (2.7 mm – 1.7 mm)**TIP TYPE:** Acorn**HANDLE:** Hollow Handle**FINISH:** Polished**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**PACKAGE INCLUDE:**

ACORN BURNISHER #21

**WI-70004280****12 Amalgam Carrier Single Ended 1.5mm Mini Surgical Dental Instrument****\$71.99** per pieceTrade price from **\$66.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Amalgam Carrier Single Ended 1.5mm Mini Tapered Tip 6" (15.2cm) Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.


Dental Amalgam Filling Carrier is a durable, high-quality instrument for precise and controlled amalgam filling application, available in single-ended or double-ended designs. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. The textured surface helps effectively grip the filling material, preventing it from slipping.

Amalgam filling carriers like a syringes, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental

cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

SINGLE ENDED AMALGAM CARRIER FEATURES

- 6" (15.2cm) Overall Length.
- Single Ended.
- 1.5mm Mini Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.

**12
PCS****100% HIGH QUALITY
STAINLESS STEEL**

**WI-70004279****12 Ball Burnisher #27/29 De Dental Amalgam Instruments****\$59.99** per pieceTrade price from **\$55.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Dental Ball Burnisher #27/29 Double Ended 1.6 mm – 3.1 mm Working Ends Premium Quality Stainless Steel De Amalgam Restorative Composite/Plastic Filling Instruments.

27/29 ball burnisher used to blend material for final contouring, to achieve sculpting of areas like grooves, fissures or pits. Can also be used to form occlusal anatomy. The working end of a burnisher is smooth and rounded (1.6 mm – 3.1 mm). Made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. These dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

BALL BURNISHER # 27/29 FEATURES

- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel

PACKAGE INCLUDE:

- 12 Pcs Ball Burnishers #27/29 Double Ended

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Ball Burnisher Dental Instrument**WI-70004277****12 Ball Burnisher #34 Single Ended Dental Amalgam Instrument****\$59.99** per pieceTrade price from **\$55.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Ball Burnisher #34 Single Ended 6.5" (16.5cm) Long T-shape Round Pattern Burnisher With a Ball on One Side Premium Quality Stainless Steel Dental Amalgam Composite Filling Instruments.

Single Ended ball burnisher #34 is often used to smooth the surface of a freshly placed amalgam restoration. T-shaped working end with a small ball on one side of "T" and paddle on the other. The shape of burnisher #34 conforms to the rounded cavity surfaces and allows ease of access into the rounded corners or junctions of the cavity surfaces to condense and shape the composite against the cavity wall. All our amalgam composite filling instruments are manufactured for optimal results and precision. These

dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

BALL BURNISHER #34 FEATURES

- Crafted with 6.5" (16.5cm) Overall Length.
- Single Ended.
- Small Pattern.
- T-Shaped Working End.
- Round Pattern Grip Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel.
- Reusable.
- Hand Crafted From High-Grade Stainless Steel.

PACKAGE INCLUDE:

- 12 Pcs Ball Burnisher #34 Single Ended 6.5" (16.5cm) – Round Pattern



WI-70003958

12 Beerendonk Caliper 0-80mm Measuring Curved Pointed Tips

\$119.99 per piece

Trade price from **\$110.39**

Minimum order 21 pieces



- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Beerendonk Caliper 5" (12.7cm) Sliding Gauge 0-80mm Measuring Dental Lab Instruments Premium Grade Stainless Steel Precision Measurement Tool with Adjustable Calibration for Dental Technicians.

Beerendonk Caliper is a high precision dental-calipers with double sided points curved tips for internal & external measurement directly into the mouth. This caliper have intra oral measurements, scale on the back for measurements of sagittal steps and depths, measuring range 0-80mm. Beerendonk caliper double-sided sliding gauge used during various types of procedures, especially in dentistry, orthodontics and prosthetics. The basic purpose of the measure is to

make precise intraoral measurements, e.g. when fitting a fixed appliance. For this purpose, the caliper has long, thin tips, bent at the ends, which allow for accurate measurement in hard-to-reach places. It features an adjustable calibration mechanism, empowering users to fine-tune measurements with precision and reliability. The intuitive layout and clear, easy-to-read displays enable users to obtain accurate measurements swiftly and efficiently.

Whether you're engaged in metalworking, machining intricate components, or overseeing engineering projects, the Beerendonk Caliper 5" 0-80mm stands ready to meet your exacting demands. Moreover, Premium Grade Stainless Steel durable



WI-70003985

12 Boley Gauge Dental Implant Lab Orthodontic Instruments

\$89.99 per piece

Trade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Dental Boley Gauge 10cm Measuring 6" (15.24cm) Long Orthodontic Crown Measuring Teeth Size Spaces Vernier Caliper Dental Implant Lab Orthodontics Surgical Instruments

Boley Gauge is one of the most commonly used measuring gauges in dentistry. The Boley Gauge is specialized dental Instruments helps in determining tooth dimensions of length, width, and thickness in millimeter increments. It is used for measuring teeth and spaces. The flat end is good for measuring teeth size and distance of spaces in the mouth.

This dental gauge featuring smooth operation and easy control to enhance user experience helps in a variety of


clinic, laboratory procedures and dental practices

BOLET GAUGE FEATURES

- Millimeter Increments Measuring Scale Range 100mm
- Small Knob for Adjustment
- Flat End Facilitates Teeth Measurement
- A Boley Gauge Has a Fixed Scale and a Sliding Scale.
- Micro Measuring for More Accuracy
- Ideal For Space Analysis & Band Selection
- Superior Quality Stainless Steel Construction

PACKAGE INCLUDE

- 12 Pcs Boley Gauge 10cm Measuring with 6" (15.24cm) Overall Length



**WI-70004105****12 Bone File #12 Surgical Dental Dentist****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Howard Bone File #12 Double Ended 7" (17.8cm) Straight-Cut 6.4mm Wide Working Tips with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Dental Instruments.

Double Ended Bone File #12 is a straight-cut 6.4mm wide tip dental instrument used for the final smoothing and reshaping of bones before sewing the wound. Straight cut bone files reduce and smooths bone with only a pull stroke due to the parallel position of the cutting edges. The knurled design allows strong grip and easy maneuvering. Made with premium quality stainless instrument for durability and long-lasting performance efficiency. The round handle design helps to reduce the risk of slippage,


making this instrument ideally suited for a variety of precise neurosurgical procedures.

HOWARD BONE FILE #12 FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- 6.4mm Wide Straight-Cut Serrated Ends.
- Double-Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 12 Pcs Howard Bone File #12 Double Ended

TECHNICAL INFORMATION**WI-70003967****12 Bracket Height Gauge 3.5mm-5mm for .022" Brackets****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 12 Bracket Positioning Height Gauge Set 3.5mm – 5.0mm for .022" Brackets Double Ended Premium Quality Aluminum Height Gauge Bracket Positioner Height Placement Dental Orthodontic Tool.

Aluminum Height Gauge featuring double ends with four measurement level tips at 3.5mm, 4mm, 4.5mm, 5mm slotted precisely for .022" brackets. Every Orthodontist must have this gauge during bracket placement procedures. Used to accurately measure height for bracket placement on the posterior and anterior teeth from the insical edge. Made from high-quality, durable Aluminum material, this gauge features a clear, easy-to-read scale that facilitates accurate measurement of bracket height on

each tooth. With precision engineering, our bracket height gauge guarantees optimal bracket positioning for effective treatment outcomes for 022 brackets.

BRACKET POSITIONING HEIGHT GAUGE FEATURES

- 6" (15.2cm) Overall Length.
- Double Ended.
- Four Bracket Height Dimensions 3.5, 4.0, 4.5 and 5.0 Millimeter Gauges.
- .022" Brackets Slot Compatible.
- Pointed Tips for Superficial Marking.
- Long Flat Handle.
- Lightweight.
- Reusable.
- Made with Top Quality Aluminum Material.

PACKAGE INCLUDE

**WI-70003957****12 Castroviejo Caliper Long 0-40mm Straight 6"****\$155.99** per pieceTrade price from **\$143.51**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Castroviejo Implant Bone Caliper, Long Model 6" (15.2cm) Straight Tips, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.


This Castroviejo Bone Caliper Straight 15.2cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 6 inches and tapers from the measurement end to the tip. The shanks and jaws are straight.

LONG CASTROVIEJO CALIPER FEATURES

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 15.2cm
- Pointed tips that open up for precise measurement.
- Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 12 Pcs Castroviejo Caliper Long 6" (15.2cm) Straight Graduated 0-40mm

TECHNICAL INFORMATION**WI-70003984****12 Crown Gauge Iwanson Spring Caliper 4" for Metal Improved Model****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Crown Gauge Iwanson Spring Caliper 4" (11.4cm) Sharp Pointed Tips for Metal Crowns Improved Model Stainless Steel Orthodontics Dental Lab Instrument.

Iwanson caliper is used by dentists for measuring the thickness of a metal so that the caliper has a pointed tip. Iwanson style calipers were originally used by dentists to measure metal and wax thickness. However, it can be used just as effectively for measuring smaller items and model engineering in arts, crafts, jewelry, and DIY projects. Made from stainless steel with sharp tip design gives accurate readings. Easy and convenient to use, soft in holding.

IWANSON CALIPER FEATURES

CROWN GAUGE CALIPER: This product is a dental caliper mainly for precision measuring. This Iwanson Caliper has an improved caliber dial, sharp pointed tips and spring action for smooth movement, and a serrated side for better grip.

HIGH-QUALITY STAINLESS STEEL: Made with superior quality durable stainless steel. The Premium Quality Stainless Steel composition of this caliper offers a long life of strong performance, as they are built to last.

FUNCTIONAL ERGONOMIC DESIGN: Features an ergonomic design to complete the task at hand and make your life easy. Provides a comfortable experience and offers more control for precise and easy operation. The Iwanson Crown Gauge measure range is 0-10mm, which can be used for



**12
PCS****WI-70001450****12 Crown Gauge Iwanson Spring Caliper 4" for Wax Improved Model****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Crown Gauge Iwanson Spring Caliper 4" (11.4cm) Blunt Blunt Tips For Wax Improved Model Stainless Steel Orthodontics Dental Lab Instrument.

Iwanson Caliper is a dental instrument used to measure the thickness of crowns and spaces. This orthodontic lab instrument is also known as Iwanson Gauge & Iwanson Spring Caliper. Iwanson Spring Caliper for Wax features Blunt tips that won't mark wax specially used for measuring wax crowns and models. High quality medical standard stainless steel product with accurate readings 10 cm length. Calibrated scale capacity 0.1mm to 10.0mm. Marked In 0.1mm increments. Can be sterilized for multiple use.

IWANSON CALIPER FEATURES


CROWN GAUGE CALIPER: Iwanson Caliper mainly for precision measuring. Iwanson Caliper has an improved caliber dial, blunt tips and spring action for smooth movement, and a serrated side for better grip.

HIGH-QUALITY STAINLESS STEEL:

Made with superior quality durable stainless steel. The Premium Quality Stainless Steel composition of this caliper offers a long life of strong performance, as they are built to last.

FUNCTIONAL ERGONOMIC DESIGN:

Features an ergonomic design to complete the task at hand and make your life easy. Provides a comfortable experience and offers more control for precise and easy operation. The Dental Orthodontic Gauge measure range is 0-

**12
PCS****WI-70004106****12 Dental Bone File #10 Double Ended 7" Curved Up & Straight****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Dental Bone File #10 Double Ended 7" (17.8cm) Long with 7mm & 6.5mm Curved Up and Straight Serrated Tips and No-Slip Round Handle Premium Quality Stainless Steel Surgical Dentist Instruments.

Bone file #10 is designed to smooth and shape bone using a pull stroke. Dental bone file #10 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. It can scrape off bony segments from soft tissues and bones. The double-ended end makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 7.0mm and 6.5mm wide. Constructed from surgical-grade stainless steel for excellent performance and longevity.


The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

BONE FILE DOUBE ENDED #10 FEATURES

- Bone File #10.
- 7" (17.8cm) Overall Length.
- 7mm & 6.5mm Serrated Ends.
- Curved Up and Straight Double-Ended.
- Non-Slip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.


PACKAGE INCLUDE

- 12 Pcs Bone File Double Ended #10 Serrated Curved Up/Straight Ends





**12
PCS**



WI-70001082

12 Hollenback Wax Carver #3 Double Ended Hollow Handle

\$89.99 per piece

Trade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Hollenback Carver #3 Double Ended Hollow Handle Premium Quality Stainless Steel Dental Wax Modeling Carver Amalgam Restorative Instruments.

Hollenback wax carver used to place, carve and contour plastic dental materials such as amalgam during restorative procedures. Having universal adaptability this hollow round handle hollenback carver 3 is ideal for placing, carving and contouring amalgam. Mainly used for the contour and carve occlusal surfaces & interproximal surfaces. With flat profile hollenback wax carver has spear shape blade that is sharp around the entire periphery. Our hollenback carver known for high degree of precision and flexibility while conducting the clinical

procedure. Made with premium quality stainless steel ensure maxium durability and dependibility, last long and can be reused overtime.

HOLLENBACK CARVER #3 FEATURES

- Overall 7.5" (19.1cm).
- Double Ended.
- Large Pointed Spear Shape Sharp Blades.
- Round Shaped Hollow Handle.
- Reusable.
- For Wax and Porcelain Carving
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 12 Pcs Hollenback Carver 3 – Double Ended – Hollow Handle



**12
PCS**



WI-70003938

12 Lecron Inlay Carver 6.5" Wax Carving Dental Tool

\$61.99 per piece

Trade price from **\$57.03**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 LeCron Wax Carver Set 6.5" (16.5cm) Double Ended One Sharp Knife End & Small Scoop End Premium Quality Stainless Steel Wax & Porcelain Carvers for Dental Technicians Dental Laboratory Instruments

Lecron Inlay Wax & Porcelain Carvers #5 is a dental laboratory tool used for shaping and carving wax patterns during the fabrication of dental prosthesis. Double-ended Lecron wax carvers are used to carve and cut wax during impressions and various other dental applications that involve dental wax. Dental technicians can use this to easily shape the wax of fixed and removable prosthodontics into the desired shape. Lecron Inlay Carver #5 features an ergonomically designed lightweight handle with a comfortable

non-slip grip and is made from premium quality stainless steel, making it sterilizable and reusable.

LECRON INLAY CARVER #5 FEATURES

- Crafted with 6.5" (16.5cm) Overall Length.
- Double Ended.
- One Sharp Knife End & One Small Scoop End.
- Ergonomic, Round & Lightweight Non-Slip Grip Handle.
- Maximum Comfort & Minimizing The Risk Of Carpal Tunnel Syndrome.
- Ideal for Carve Wax Blocks and Dentures.
- Sterilizable and Reusable.
- Made with Premium Grade Stainless Steel.



WI-70004101

12 Miller Bone File #64 Double Ended

\$69.99 per piece

Trade price from \$64.39
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

12 Miller Bone File #64 Double Ended is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$69.99, falling to \$64.39 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003936

12 P.K Thomas #3 Dental Wax Modelling Carver

\$49.99 per piece

Trade price from \$45.99
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

PTK3 – 12 Pcs P.K Thomas #3 Dental Wax Modelling Carver Set with Sharp, Conicular Ponited End Premium Quality Dental Restorative P.K. Thomas Waxing Instrument



PKT3 features a pointed burnisher used to perfect and enhance the supplemental and developmental grooves. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture

when working on hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #3 WAX CARVER FEATURES

- Sharp, Conicular End to Precisely Carve and Shape Anatomical Features.
- Ponited Double-ended
- Non-slip Knurled Handlee for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE



**WI-70003935****12 P.K Thomas #4 Dental Wax Modelling Carver****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking




PTK4 – 12 Pcs P.K Thomas #4 Dental Wax Modelling Carver with Precise, Spatula-like Ends Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT4 is a modified carver used to perfect the external contours and remove excess wax at the cavo-surface margins. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture

when working on hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #4 WAX CARVER FEATURES

- Use to Shape External Contours and Remove Excess Wax.
- Precise, Spatula Shaped Double Ends.
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE**WI-70003934****12 P.K Thomas #5 Dental Wax Modelling Carver****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



PTK5 – 12 Pcs P.K Thomas #5 Dental Wax Modelling Carver with Unique Design Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT5 is a special carver used to remove excess wax as cusp ridges are developed; its contour maintains desired convexity at these ridges. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote

neutral wrist posture when working on hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #5 WAX CARVER FEATURES

- Use to Contour Tooth Anatomy and Remove Excess Wax.
- Unique Carver Tips Design.
- Double Ended.
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.





12 Pieces

WI-70004100

12 Palti Bone Graft Scoop Dental Implant

\$92.71 per piece

Trade price from \$85.29
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

12 Palti Bone Graft Scoop Dental Implant is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$92.71, falling to \$85.29 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



12 Pieces

WI-70003982

12 Ridge Mapping Caliper Surgical Dental

\$154.18 per piece

Trade price from \$141.85
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

12 Ridge Mapping Caliper Surgical Dental is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$154.18, falling to \$141.85 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70002094**

2 Acorn Burnisher #21 Double Ended with Hollow Handle

\$13.99 per pieceTrade price from **\$12.87**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Acorn Burnisher #21 with Hollow Handle Dental Composite Restoration Instruments

Acorn burnisher 21 is a double-ended (2.7 mm – 1.7 mm) dental instrument used to condense, smooth, carve and polish amalgam. Acorn-shaped ends are excellent for carving occlusal anatomy. Acorn Burnishers used to assure that amalgam restorations are free of irregularities and are smooth

ACORN BURNISHER 21 FEATURES

- 1.7mm & 2.7mm acorn-shaped ends.
- Double-Ended Silver Color
- Hand crafted from High Quality stainless steel.
- Available in hollow, ergonomic handle.

PACKAGE INCLUDE:

- 2 Pcs Acorn Burnisher 21 Double Ended with Hollow Handle

TECHNICAL SPECIFICATIONS


MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Dental Burnisher**APPLICATION:** Composite Placement and Contouring**COLOR:** Silver**TIP:** Double Ended (2.7 mm – 1.7 mm)**TIP TYPE:** Acorn**HANDLE:** Hollow Handle**FINISH:** Polished**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test,

Performance Test, Shape Test

ACORN BURNISHER 21

**2
PCS****WI-70002032**

2 Automatic Crown Remover Gun Type with 5 Attachments

\$36.99 per pieceTrade price from **\$34.03**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Dental Automatic Crown & Bridge Remover Gun Spring Loaded with 3 Tips 1 Wrench and Assemble Key Restoration Dental Instruments

Automatic crown and bridge remover is used for removal of temporary and fixed restorations. The automatic crown remover is a spring loaded mechanical device to remove the bonded crowns and bridges by an intensive impulse of the impact stress caused by the stored energy in the spring located inside the instrument. Cement layer breaks with the impulse and the crown or bridge can be removed easily. Packet contain spring loaded action sold with 3 working tips of different styles, a wrench to open it for cleaning if required. It also contain 1 assembly key.

AUTOMATIC CROWN REMOVER SET FEATURES

- Ideally designed for all areas of the mouth
- Spring-loaded mechanism delivers a quick force that releases the crown
- Comes with three tempered stainless steel tips, one wrench and assemble key
- Made with Highest quality medical grade Stainless Steel

PACKAGE INCLUDE:

- 2 Pcs Automatic Crowns Bridges Remover Gun Type with 5 Attachments

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel**INSTRUMENTS TYPE:** Restorative Dental Instruments

**WI-70001465****2 Castroviejo Caliper Long 6" Straight + Angled Graduated 0-40mm****\$32.99** per pieceTrade price from **\$30.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Castroviejo Implant Bone Caliper 6" (15.2cm) Long Model Straight & Angled Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper 15.2cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 6 inches and tapers from the measurement end to the tip. The shanks are straight with 30° Angled jaws

LONG CASTROVIEJO CALIPER FEATURES

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 15.2cm
- Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 1 Pc Castroviejo Caliper 6" (15.2cm) Long Model – 30° Angled Tip Graduated 0-40mm
- 1 Pc Castroviejo Caliper 6" (15.2cm) Long Model – Straight Tip Graduated 0-40mm

**WI-70001456****2 Castroviejo Caliper Long 7" Straight + Angled Graduated 0-40mm****\$35.99** per pieceTrade price from **\$33.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Castroviejo Implant Bone Caliper 7" (17.8cm) Long Model Straight & Angled Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper 17.8cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 7 inches and tapers from the measurement end to the tip. The shanks are straight with Straight & 30° Angled jaws.

LONG CASTROVIEJO CALIPER FEATURES

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 17.8cm
- Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 1 Pc Castroviejo Caliper 7" (17.8cm) Long Model – 30° Angled Tip Graduated 0-40mm
- 1 Pc Castroviejo Caliper 7" (17.8cm) Long Model – Straight Tip Graduated 0-40mm

**2
PCS****WI-70002037**

2 Crown Remover Manual with 3 Tips Double Ended

\$43.74 per pieceTrade price from **\$40.24**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Manual Crown Remover Complete Set with 3 Attachment Tips Crown and Bridge Removers Restorative Dentistry Instruments.

Manual dental crown bridge remover kit double end is made of stainless steel and comes with 3 tips. Classic, easy-to-use crown remover is perfect for removing temporary or permanent restorations. The crown remover tips attach to the working ends of crown removers to grip the crown securely without slipping.

- Comes with 3 tempered stainless steel tips
- Made with Highest quality medical grade Stainless Steel

PACKAGE INCLUDE:

- 2 Pcs Manual Crown Remover with 3 Attachment Tips – Double Ended

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Restorative Dental Instruments**APPLICATION:** Permanent and Temporary Crowns and Bridges Removal**TIPS SHAPE:** Hook Shape / Curved / Serrated**HANDLE:** Hand Held – Double Ended**REUSABLE:** Yes**AUTOCCLAVED:** Yes**MANUAL CROWN REMOVER SET FEATURES**

- For use in all areas of the mouth
- Manual, hammer mechanism
- Double Ended

**WI-70000864**

2 Krekeler Sliding Calipers Round Gauge 5" 0-80mm Straight Curved

\$35.99 per pieceTrade price from **\$33.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Krekeler Sliding Calipers Round Gauge 5" (12.7cm) Measures up to 0-80mm Straight & Curved Pointed Tips with Screw Hold Mechanism Premium Grade Stainless Steel Dental Instruments for Implantology, Surgical, Restorative or Orthodontic Procedures.


Krekeler sliding implant caliper gauge to measure interval and length of implants. This multi purpose krekeler sliding caliper round gauge is specially designed to accurately measure bone and teeth or spaces between teeth and implants. Accurately measure the width and depth of the osteotomy. Facilitate positioning of implants. Krekeler sliding caliper measures 0-80mm range with pointed tips can be used to make markings. The narrow tips ideal for implants. With its versatile functionality,

this Krekeler caliper is suitable for a wide range of applications from metal fabrication to woodworking and engineering tasks as well. Crafted from stainless steel, this caliper is built to last. Its robust construction ensures longevity and reliability.

KREKELER SLIDING CALIPER ROUND GAUGE FEATURES

- Krekeler Sliding Caliper Round Gauge is 5" (12.7cm) Long
- 0-80mm Measuring Scale.
- Equipped With Easy-to-Read Markings
- Pointed Tips for Accurate Measurement.
- Screw Hold Mechanism
- Narrow Tips Ideal for Implants



**WI-70002093**

3 Acorn Burnisher #21 Double Ended with Hollow Handle

\$19.99 per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Acorn Burnisher #21 with Hollow Handle Dental Composite Restoration Instruments

Acorn burnisher 21 is a double-ended (2.7 mm – 1.7 mm) dental instrument used to condense, smooth, carve and polish amalgam. Acorn-shaped ends are excellent for carving occlusal anatomy. Acorn Burnishers used to assure that amalgam restorations are free of irregularities and are smooth

ACORN BURNISHER 21 FEATURES

- 1.7mm & 2.7mm acorn-shaped ends.
- Double-Ended Silver Color
- Hand crafted from High Quality stainless steel.
- Available in hollow, ergonomic handle.

PACKAGE INCLUDE:

- 3 Pcs Acorn Burnisher 21 Double Ended with Hollow Handle

TECHNICAL SPECIFICATIONS

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Dental Burnisher**APPLICATION:** Composite Placement and Contouring**COLOR:** Silver**TIP:** Double Ended (2.7 mm – 1.7 mm)**TIP TYPE:** Acorn**HANDLE:** Hollow Handle**FINISH:** Polished**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test

ACORN BURNISHER 21

**WI-70000522**

3 Amalgam Carriers Single Ended 1.5mm, 2.0mm, 3.0mm Tip

\$24.99 per pieceTrade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Amalgam Carriers Single Ended 1.5mm, 2.0mm & 3.0mm Tapered Tip 6" (15.2cm) Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument for precise and controlled amalgam filling application, available in single-ended or double-ended designs. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. The textured surface helps effectively grip the filling material, preventing it from slipping.

Amalgam filling carriers like a syringes, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental

cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

SINGLE ENDED AMALGAM CARRIER FEATURES

- 6" (15.2cm) Overall Length.
- Single Ended.
- Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.



**3
PCS****WI-70004273**

3 Ball Burnisher #27/29 Dental Amalgam Instruments

\$19.99 per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Ball Burnisher #27/29 Double Ended 1.6 mm – 3.1 mm Working Ends Premium Quality Stainless Steel De Amalgam Restorative Composite/Plastic Filling Instruments.

27/29 ball burnisher used to blend material for final contouring, to achieve sculpting of areas like grooves, fissures or pits. Can also be used to form occlusal anatomy. The working end of a burnisher is smooth and rounded (1.6 mm – 3.1 mm). Made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. These dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.


BALL BURNISHER # 27/29 FEATURES

- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel

PACKAGE INCLUDE:

- 3 Pcs Ball Burnishers #27/29 Double Ended

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENTS TYPE:** Ball Burnisher Dental Instrument**3
PCS****WI-70004274**

3 Ball Burnisher #1 Double Ended Dental Amalgam Instruments

\$15.99 per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Ball Burnisher #1 Double 1.3mm & 1.8mm Ball Ends 7" (17.8cm) Long Premium Quality Stainless Steel Dental Amalgam Filling Instruments.

Dental 1 Ball Burnisher with 7" overall length used to condense, carve and polish restorative materials such as amalgam. 1 Ball Burnisher composite amalgam filling instrument is non-stick and comes with ball shaped ends. Amalgam Ball Burnisher #1 can be used to blend your composite material to the margin for final contouring. All of our composite placement instruments made with stainless steel ensuring excellent performance instruments which are autoclavable and can reuse after sterilization.

- Crafted With 7" (17.8cm) Overall Length.
- 1.3mm & 1.8mm Ball Ends.
- Double-ended.
- Non-slip Knurled With Center Round Handle.
- Autoclavable & Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 3 Pcs Ball Burnisher #1 Double Ended 1.3mm & 1.8mm

Discover more **Amalgam Filling Ball Burnisher**

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.

DENTAL BALL BURNISHER #1

**WI-70004271****3 Ball Burnisher #34 Single Ended Dental Amalgam Instrument****\$20.99** per pieceTrade price from **\$19.31**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Ball Burnisher #34 Single Ended 6.5" (16.5cm) Long T-shape Round Pattern Burnisher With a Ball on One Side Premium Quality Stainless Steel Dental Amalgam Composite Filling Instruments.

Single Ended ball burnisher #34 is often used to smooth the surface of a freshly placed amalgam restoration. T-shaped working end with a small ball on one side of "T" and paddle on the other. The shape of burnisher #34 conforms to the rounded cavity surfaces and allows ease of access into the rounded corners or junctions of the cavity surfaces to condense and shape the composite against the cavity wall. All our amalgam composite filling instruments are manufactured for optimal results and precision. These

dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

BALL BURNISHER #34 FEATURES

- Crafted with 6.5" (16.5cm) Overall Length.
- Single Ended.
- Small Pattern.
- T-Shaped Working End.
- Round Pattern Grip Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel.
- Reusable.
- Hand Crafted From High-Grade Stainless Steel.

PACKAGE INCLUDE:

- 3 Pcs Ball Burnisher #34 Single Ended 6.5" (16.5cm) – Round Pattern

**WI-70003955****3 Beerendonk Caliper 0-80mm Measuring Curved Pointed Tips****\$35.99** per pieceTrade price from **\$33.11**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Beerendonk Caliper 5" (12.7cm) Sliding Gauge 0-80mm Measuring Dental Lab Instruments Premium Grade Stainless Steel Precision Measurement Tool with Adjustable Calibration for Dental Technicians.

Beerendonk Caliper is a high precision dental-calipers with double sided points curved tips for internal & external measurement directly into the mouth. This caliper have intra oral measurements, scale on the back for measurements of sagittal steps and depths, measuring range 0-80mm. Beerendonk caliper double-sided sliding gauge used during various types of procedures, especially in dentistry, orthodontics and prosthetics. The basic purpose of the measure is to

make precise intraoral measurements, e.g. when fitting a fixed appliance. For this purpose, the caliper has long, thin tips, bent at the ends, which allow for accurate measurement in hard-to-reach places. It features an adjustable calibration mechanism, empowering users to fine-tune measurements with precision and reliability. The intuitive layout and clear, easy-to-read displays enable users to obtain accurate measurements swiftly and efficiently.

Whether you're engaged in metalworking, machining intricate components, or overseeing engineering projects, the Beerendonk Caliper 5" 0-80mm stands ready to meet your exacting demands. Moreover, Premium Grade Stainless Steel durable



**WI-70003966****3 Bracket Positioning Height Gauge 3.5mm-5mm for .022\"/>****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 3 Bracket Positioning Height Gauge Set 3.5mm – 5.0mm for .022\"/>

Aluminum Height Gauge featuring double ends with four measurement level tips at 3.5mm, 4mm, 4.5mm, 5mm slotted precisely for .022\"/>

each tooth. With precision engineering, our bracket height gauge guarantees optimal bracket positioning for effective treatment outcomes for 022 brackets.

BRACKET POSITIONING HEIGHT GAUGE FEATURES

- 6\"/>
- Double Ended.
- Four Bracket Height Dimensions 3.5, 4.0, 4.5 and 5.0 Millimeter Gauges.
- .022\"/>
- Pointed Tips for Superficial Marking.
- Long Flat Handle.
- Lightweight.
- Reusable.
- Made with Top Quality Aluminum Material.

PACKAGE INCLUDE**WI-70003979****3 Crown Gauge Iwanson Spring Caliper 4\"/>****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Crown Gauge Iwanson Spring Caliper 4\"/>

Iwanson caliper is used by dentists for measuring the thickness of a metal so that the caliper has a pointed tip. Iwanson style calipers were originally used by dentists to measure metal and wax thickness. However, it can be used just as effectively for measuring smaller items and model engineering in arts, crafts, jewelry, and DIY projects. Made from stainless steel with sharp tip design gives accurate readings. Easy and convenient to use, soft in holding.

CROWN GAUGE CALIPER: This product is a dental caliper mainly for precision measuring. This Iwanson Caliper has an improved caliber dial, sharp pointed tips and spring action for smooth movement, and a serrated side for better grip.

HIGH-QUALITY STAINLESS STEEL: Made with superior quality durable stainless steel. The Premium Quality Stainless Steel composition of this caliper offers a long life of strong performance, as they are built to last.

FUNCTIONAL ERGONOMIC DESIGN: Features an ergonomic design to complete the task at hand and make your life easy. Provides a comfortable experience and offers more control for precise and easy operation. The Iwanson Crown Gauge measure range is 0-10mm, which can be used for

**IWANSON CALIPER FEATURES**

**3
PCS****WI-70001452****3 Crown Gauge Iwanson Spring Caliper 4" for Wax Improved Model****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Crown Gauge Iwanson Spring Caliper 4" (11.4cm) Blunt Blunt Tips For Wax Improved Model Stainless Steel Orthodontics Dental Lab Instrument.

Iwanson Caliper is a dental instrument used to measure the thickness of crowns and spaces. This orthodontic lab instrument is also known as Iwanson Gauge & Iwanson Spring Caliper. Iwanson Spring Caliper for Wax features Blunt tips that won't mark wax product with accurate readings 10 cm length. Calibrated scale capacity 0.1mm to 10.0mm. Marked In 0.1mm increments. Can be sterilized for multiple use.

CROWN GAUGE CALIPER: Iwanson Caliper mainly for precision measuring. Iwanson Caliper has an improved caliber dial, blunt tips and spring action for smooth movement, and a serrated side for better grip.

HIGH-QUALITY STAINLESS STEEL: Made with superior quality durable stainless steel. The Premium Quality Stainless Steel composition of this caliper offers a long life of strong performance, as they are built to last.

FUNCTIONAL ERGONOMIC DESIGN: Features an ergonomic design to complete the task at hand and make your life easy. Provides a comfortable experience and offers more control for precise and easy operation. The Dental Orthodontic Gauge measure range is 0-10mm, which can be used for precise measurement.

IWANSON CALIPER FEATURES**WI-70002036****3 Crown Remover Manual with 3 Tips Double Ended****\$62.49** per pieceTrade price from **\$57.49**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Manual Crown Remover Complete Set with 3 Attachment Tips Crown and Bridge Removers Restorative Dentistry Instruments.

Manual dental crown bridge remover kit double end is made of stainless steel and comes with 3 tips. Classic, easy-to-use crown remover is perfect for removing temporary or permanent restorations. The crown remover tips attach to the working ends of crown removers to grip the crown securely without slipping.


- Comes with 3 tempered stainless steel tips
- Made with Highest quality medical grade Stainless Steel

PACKAGE INCLUDE:

- 3 Pcs Manual Crown Remover with 3 Attachment Tips – Double Ended

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Restorative Dental Instruments**APPLICATION:** Permanent and Temporary Crowns and Bridges Removal**TIPS SHAPE:** Hook Shape / Curved /Serrated**HANDLE:** Hand Held – Double Ended**REUSABLE:** Yes**AUTOC LAVED:** Yes**MANUAL CROWN REMOVER SET FEATURES**

- For use in all areas of the mouth
- Manual, hammer mechanism
- Double Ended

**3
PCS**

**3
PCS**

DENTAL BONE FILE #10 DOUBLE ENDED

**WI-70004099****3 Dental Bone File #10 Double Ended 7" Curved Up & Straight****\$23.99** per pieceTrade price from **\$22.07**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Bone File #10 Double Ended 7" (17.8cm) Long with 7mm & 6.5mm Curved Up and Straight Serrated Tips and No-Slip Round Handle Premium Quality Stainless Steel Surgical Dentist Instruments.

Bone file #10 is designed to smooth and shape bone using a pull stroke. Dental bone file #10 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. It can scrape off bony segments from soft tissues and bones. The double-ended end makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 7.0mm and 6.5mm wide. Constructed from surgical-grade stainless steel for excellent performance and longevity.

The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

BONE FILE DOUBE ENDED #10 FEATURES

- Bone File #10.
- 7" (17.8cm) Overall Length.
- 7mm & 6.5mm Serrated Ends.
- Curved Up and Straight Double-Ended.
- Non-Slip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 3 Pcs Bone File Double Ended #10 Serrated Curved Up/Straight Ends

**3
PCS****WI-70001084****3 Hollenback Wax Carver #3 Double Ended Hollow Handle****\$25.99** per pieceTrade price from **\$23.91**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Hollenback Carver #3 Double Ended Hollow Handle Premium Quality Stainless Steel Dental Wax Modeling Carver Amalgam Restorative Instruments.

Hollenback wax carver used to place, carve and contour plastic dental materials such as amalgam during restorative procedures. Having universal adaptability this hollow round handle hollenback carver 3 is ideal for placing, carving and contouring amalgam. Mainly used for the contour and carve occlusal surfaces & interproximal surfaces. With flat profile hollenback wax carver has spear shape blade that is sharp around the entire periphery. Our hollenback carver known for high degree of precision and flexibility while conducting the clinical

procedure. Made with premium quality stainless steel ensure maxium durability and dependability, last long and can be reused overtime.

HOLLENBACK CARVER #3 FEATURES

- Overall 7.5" (19.1cm).
- Double Ended.
- Large Pointed Spear Shape Sharp Blades.
- Round Shaped Hollow Handle.
- Reusable.
- For Wax and Porcelain Carving
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 3 Pcs Hollenback Carver 3 – Double Ended – Hollow Handle

**WI-70004098****3 Howard Bone File #12 Double Ended****\$23.99** per pieceTrade price from **\$22.07**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Howard Bone File #12 Double Ended 7" (17.8cm) Straight-Cut 6.4mm Wide Working Tips with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Dental Instruments.

Double Ended Bone File #12 is a straight-cut 6.4mm wide tip dental instrument used for the final smoothing and reshaping of bones before sewing the wound. Straight cut bone files reduce and smooths bone with only a pull stroke due to the parallel position of the cutting edges. The knurled design allows strong grip and easy maneuvering. Made with premium quality stainless instrument for durability and long-lasting performance efficiency. The round handle design helps to reduce the risk of slippage,

making this instrument ideally suited for a variety of precise neurosurgical procedures.

HOWARD BONE FILE #12 FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- 6.4mm Wide Straight-Cut Serrated Ends.
- Double-Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 3 Pcs Howard Bone File #12 Double Ended

TECHNICAL INFORMATION**WI-70003931****3 Lecron Carver Wax 6.5" Double Ended Carving Dental Tool****\$20.99** per pieceTrade price from **\$19.31**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 LeCron Wax Carver 6.5" (16.5cm) Double Ended One Sharp Knife End & Small Scoop End Premium Quality Stainless Steel Wax & Porcelain Carvers for Dental Technicians Dental Laboratory Instruments

Lecron Inlay Wax & Porcelain Carvers #5 is a dental laboratory tool used for shaping and carving wax patterns during the fabrication of dental prosthesis. Double-ended Lecron wax carvers are used to carve and cut wax during impressions and various other dental applications that involve dental wax. Dental technicians can use this to easily shape the wax of fixed and removable prosthodontics into the desired shape. Lecron Inlay Carver #5 features an ergonomically designed lightweight handle with a comfortable

non-slip grip and is made from premium quality stainless steel, making it sterilizable and reusable.

LECRON INLAY CARVER #5 FEATURES

- Crafted with 6.5" (16.5cm) Overall Length.
- Double Ended.
- One Sharp Knife End & One Small Scoop End.
- Ergonomic, Round & Lightweight Non-Slip Grip Handle.
- Maximum Comfort & Minimizing The Risk Of Carpal Tunnel Syndrome.
- Ideal for Carve Wax Blocks and Dentures.
- Sterilizable and Reusable.
- Made with Premium Grade Stainless Steel.

**WI-70004091****3 Miller Bone File #64 Double Ended****\$25.99** per pieceTrade price from **\$23.91**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

3 Miller Bone File #64 Double Ended is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$25.99, falling to \$23.91 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003929****3 P.K Thomas #3 Dental Wax Modelling Carver****\$20.99** per pieceTrade price from **\$19.31**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



PTK3 – 3 Pcs P.K Thomas #3 Dental Wax Modelling Carver Set with Sharp, Conicular Ponited End Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT3 features a pointed burnisher used to perfect and enhance the supplemental and developmental grooves. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture

when working on hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #3 WAX CARVER FEATURES

- Sharp, Conicular End to Precisely Carve and Shape Anatomical Features.
- Ponited Double-ended
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE

**WI-70003928****3 P.K Thomas #4 Dental Wax Modelling Carver****\$20.99** per pieceTrade price from **\$19.31**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



PTK4 – 3 Pcs P.K Thomas #4 Dental Wax Modelling Carver with Precise, Spatula-like Ends Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT4 is a modified carver used to perfect the external contours and remove excess wax at the cavo-surface margins. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture

when working on hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #4 WAX CARVER FEATURES

- Use to Shape External Contours and Remove Excess Wax.
- Precise, Spatula Shaped Double Ends.
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE**WI-70003927****3 P.K Thomas #5 Dental Wax Modelling Carver****\$20.99** per pieceTrade price from **\$19.31**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


PTK5 – 3 Pcs P.K Thomas #5 Dental Wax Modelling Carver with Unique Design Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT5 is a special carver used to remove excess wax as cusp ridges are developed; its contour maintains desired convexity at these ridges. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote

neutral wrist posture when working on hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #5 WAX CARVER FEATURES

- Use to Contour Tooth Anatomy and Remove Excess Wax.
- Unique Carver Tips Design.
- Double Ended.
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.



**WI-70004090****3 Palti Bone Graft Scoop Dental Implant****\$30.84** per pieceTrade price from **\$28.37**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

3 Palti Bone Graft Scoop Dental Implant is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$30.84, falling to \$28.37 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003977****3 Ridge Mapping Caliper Surgical Dental****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

WAQAR INTERNATIONAL 3 PCS RIDGE MAPPING CALIPER SET 6" (15.24CM) IMPLANT SURGERY INSTRUMENT

Ridge Mapping Caliper Measures Upto 0-25mm Range 6" (15.24cm) Long Dental Implants Surgery Bone Measuring Instrument

Ridge Mapping Caliper used to measure bone and teeth or spaces between teeth. It can be used in dental surgical, restorative or orthodontic procedures. During Implantology, this caliper provides the surgeon with information which can be required during dental treatment. It is used as measuring instruments in preparation for an implant procedure required for various dental treatments.

RIDGE MAPPING CALIPER FEATURES

- Measures up to 25mm
- 6" length
- Pointed tips for measuring
- Flat finger ring design
- Made with High Quality Stainless Steel

PACKAGE INCLUDES:

- 3 Pcs Ridge Mapping Caliper Set Measures upto 0.25mm

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Instruments**APPLICATION:** Implantology**TIPS:** Pointed Measuring Tups**SCALE:** upto -0-25mm (Millimeter) Measuring Scale

**WI-70002092****5 Acorn Burnisher 21 Double Ended with Hollow Handle****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Acorn Burnisher #21 with Hollow Handle Dental Composite Restoration Instruments

Acorn burnisher 21 is a double-ended (2.7 mm – 1.7 mm) dental instrument used to condense, smooth, carve and polish amalgam. Acorn-shaped ends are excellent for carving occlusal anatomy. Acorn Burnishers used to assure that amalgam restorations are free of irregularities and are smooth

ACORN BURNISHER 21 FEATURES

- 1.7mm & 2.7mm acorn-shaped ends.
- Double-Ended Silver Color
- Hand crafted from High Quality stainless steel.
- Available in hollow, ergonomic handle.

PACKAGE INCLUDE:

- 5 Pcs Acorn Burnisher 21 Double Ended with Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Burnisher**APPLICATION:** Composite Placement and Contouring**COLOR:** Silver**TIP:** Double Ended (2.7 mm – 1.7 mm)**TIP TYPE:** Acorn**HANDLE:** Hollow Handle**FINISH:** Polished**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**ACORN BURNISHER 21****WI-70002031****5 Automatic Crown Remover Gun Type with 5 Attachments****\$79.99** per pieceTrade price from **\$73.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Dental Automatic Crown & Bridge Remover Gun Spring Loaded with 3 Tips 1 Wrench and Assemble Key Restoration Dental Instruments

Automatic crown and bridge remover is used for removal of temporary and fixed restorations. The automatic crown remover is a spring loaded mechanical device to remove the bonded crowns and bridges by an intensive impulse of the impact stress caused by the stored energy in the spring located inside the instrument. Cement layer breaks with the impulse and the crown or bridge can be removed easily. Packet contain spring loaded action sold with 3 working tips of different styles, a wrench to open it for cleaning if required. It also contain 1 assembly key.


AUTOMATIC CROWN REMOVER SET FEATURES

- Ideally designed for all areas of the mouth
- Spring-loaded mechanism delivers a quick force that releases the crown
- Comes with three tempered stainless steel tips, one wrench and assemble key
- Made with Highest quality medical grade Stainless Steel

PACKAGE INCLUDE:

- 5 Pcs Automatic Crowns Bridges Remover Gun Type with 5 Attachments

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Restorative Dental Instruments**5
PCS**

**WI-70003946****5 Bone Ridge Mapping Calipers 6.5" Ring Style Measures Upto 0-35mm****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Dental Implant Bone Ridge Mapping Calipers 6.5" (16.5cm) Ring Style Measures Upto 0-35mm With Gold Handle Precision Stainless Steel Measurement Tool for Mapping Bone Width Implantology, Dental Restorative Orthodontic Procedures.

Bone Ridge Mapping Caliper is expertly crafted with a ring style that allows for easy handling and accurate measurements. This caliper is designed for dental implant procedures, specifically for mapping bone ridges with precision. The innovative ring-style design, coupled with adjustable arms, offers flexibility and precision, enabling you to tailor measurements to your exact specifications with ease and accuracy. The gold-plated handles not

only provide a luxurious touch but also ensure durability and longevity.

Dental Implant Bone Ridge Mapping Caliper is made of high-quality surgical instruments stainless steel, guaranteeing reliability and accuracy. The caliper is lightweight and easy to use, making it a valuable tool for any dental professional. Whether you are a seasoned dentist or a dental student, this is a must-have tool for ensuring successful dental implant procedures.

BONE RIDGE MAPPING CALIPER FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- Measures Upto 0-35mm.
- Clear and Easy-to-read Markings.

**WI-70002035****5 Crown Remover Manual with 3 Tips Double Ended****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Manual Crown Remover Complete Set with 3 Attachment Tips Crown and Bridge Removers Restorative Dentistry Instruments.

Manual dental crown bridge remover kit double end is made of stainless steel and comes with 3 tips. Classic, easy-to-use crown remover is perfect for removing temporary or permanent restorations. The crown remover tips attach to the working ends of crown removers to grip the crown securely without slipping.

- Comes with 3 tempered stainless steel tips
- Made with Highest quality medical grade Stainless Steel

PACKAGE INCLUDE:

- 5 Pcs Manual Crown Remover with 3 Attachment Tips – Double Ended

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Restorative Dental Instruments**APPLICATION:** Permanent and Temporary Crowns and Bridges Removal**TIPS SHAPE:** Hook Shape / Curved / Serrated**HANDLE:** Hand Held – Double Ended**REUSABLE:** Yes**AUTOCCLAVED:** Yes**MANUAL CROWN REMOVER SET FEATURES**

- For use in all areas of the mouth
- Manual, hammer mechanism
- Double Ended



**WI-70002331****50 Ball Burnisher #27/29 De Dental Amalgam Instruments****\$124.99** per pieceTrade price from **\$114.99**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ OEM / private-label marking
- ✓ Reusable, autoclavable

Pack of 50 Dental Ball Burnisher #27/29 Double Ended 1.6 mm – 3.1 mm Working Ends Premium Quality Stainless Steel De Amalgam Restorative Composite/Plastic Filling Instruments.

27/29 ball burnisher used to blend material for final contouring, to achieve sculpting of areas like grooves, fissures or pits. Can also be used to form occlusal anatomy. The working end of a burnisher is smooth and rounded (1.6 mm – 3.1 mm). Made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. These dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

BALL BURNISHER # 27/29 FEATURES

- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel

PACKAGE INCLUDE:

- 50 Pcs Ball Burnishers #27/29 Double Ended

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Ball Burnisher Dental Instrument**WI-70003917****50 Lecron Inlay Wax & Porcelain Carver 6.5\"/>****\$208.99** per pieceTrade price from **\$192.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ OEM / private-label marking
- ✓ Reusable, autoclavable

Pack of 50 LeCron Wax Carver Set 6.5\"/>

Lecron Inlay Wax & Porcelain Carvers #5 is a dental laboratory tool used for shaping and carving wax patterns during the fabrication of dental prosthesis. Double-ended Lecron wax carvers are used to carve and cut wax during impressions and various other dental applications that involve dental wax. Dental technicians can use this to easily shape the wax of fixed and removable prosthodontics into the desired shape. Lecron Inlay Carver #5 features an ergonomically designed lightweight handle with a comfortable

non-slip grip and is made from premium quality stainless steel, making it sterilizable and reusable.

LECRON INLAY CARVER #5 FEATURES

- Crafted with 6.5\"/>
- Double Ended.
- One Sharp Knife End & One Small Scoop End.
- Ergonomic, Round & Lightweight Non-Slip Grip Handle.
- Maximum Comfort & Minimizing The Risk Of Carpal Tunnel Syndrome.
- Ideal for Carve Wax Blocks and Dentures.
- Sterilizable and Reusable.
- Made with Premium Grade Stainless Steel.



**6
PCS****WI-70004268****6 Ball Burnisher #1 Double Ended Dental Amalgam Instruments****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Ball Burnisher #1 Double 1.3mm & 1.8mm Ball Ends 7" (17.8cm) Long Premium Quality Stainless Steel Dental Amalgam Filling Instruments.

Dental 1 Ball Burnisher with 7" overall length used to condense, carve and polish restorative materials such as amalgam. 1 Ball Burnisher composite amalgam filling instrument is non-stick and comes with ball shaped ends.

Amalgam Ball Burnisher #1 can be used to blend your composite material to the margin for final contouring. All of our composite placement instruments made with stainless steel ensuring excellent performance instruments which are autoclavable and can reuse after sterilization.

DENTAL BALL BURNISHER #1

- Crafted With 7" (17.8cm) Overall Length.
- 1.3mm & 1.8mm Ball Ends.
- Double-ended.
- Non-slip Knurled With Center Round Handle.
- Autoclavable & Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 6 Pcs Ball Burnisher #1 Double Ended 1.3mm & 1.8mm

Discover more **Amalgam Filling Ball Burnisher**

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.

**6
PCS****WI-70004267****6 Ball Burnisher #27/29 De Dental Amalgam Instruments****\$35.99** per pieceTrade price from **\$33.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Ball Burnisher #27/29 Double Ended 1.6 mm – 3.1 mm Working Ends Premium Quality Stainless Steel De Amalgam Restorative Composite/Plastic Filling Instruments.

27/29 ball burnisher used to blend material for final contouring, to achieve sculpting of areas like grooves, fissures or pits. Can also be used to form occlusal anatomy. The working end of a burnisher is smooth and rounded (1.6 mm – 3.1 mm). Made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. These dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

BALL BURNISHER # 27/29 FEATURES

- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel

PACKAGE INCLUDE:

- 6 Pcs Ball Burnishers #27/29 Double Ended

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Ball Burnisher Dental Instrument

**WI-70004265**

6 Ball Burnisher #34 Single Ended Dental Amalgam Instrument

\$35.99 per pieceTrade price from **\$33.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Ball Burnisher #34 Single Ended 6.5" (16.5cm) Long T-shape Round Pattern Burnisher With a Ball on One Side Premium Quality Stainless Steel Dental Amalgam Composite Filling Instruments.

Single Ended ball burnisher #34 is often used to smooth the surface of a freshly placed amalgam restoration. T-shaped working end with a small ball on one side of "T" and paddle on the other. The shape of burnisher #34 conforms to the rounded cavity surfaces and allows ease of access into the rounded corners or junctions of the cavity surfaces to condense and shape the composite against the cavity wall. All our amalgam composite filling instruments are manufactured for optimal results and precision. These

dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

BALL BURNISHER #34 FEATURES

- Crafted with 6.5" (16.5cm) Overall Length.
- Single Ended.
- Small Pattern.
- T-Shaped Working End.
- Round Pattern Grip Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel.
- Reusable.
- Hand Crafted From High-Grade Stainless Steel.

PACKAGE INCLUDE:

- 6 Pcs Ball Burnisher #34 Single Ended 6.5" (16.5cm) – Round Pattern

**WI-70003952**

6 Beerendonk Caliper 0-80mm Measuring Curved Pointed Tips

\$65.99 per pieceTrade price from **\$60.71**
Minimum order 21 pieces



- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Beerendonk Caliper 5" (12.7cm) Sliding Gauge 0-80mm Measuring Dental Lab Instruments Premium Grade Stainless Steel Precision Measurement Tool with Adjustable Calibration for Dental Technicians.

Beerendonk Caliper is a high precision dental-calipers with double sided points curved tips for internal & external measurement directly into the mouth. This caliper have intra oral measurements, scale on the back for measurements of sagittal steps and depths, measuring range 0-80mm. Beerendonk caliper double-sided sliding guage used during various types of procedures, especially in dentistry, orthodontics and prosthetics. The basic purpose of the measure is to

make precise intraoral measurements, e.g. when fitting a fixed appliance. For this purpose, the caliper has long, thin tips, bent at the ends, which allow for accurate measurement in hard-to-reach places. It features an adjustable calibration mechanism, empowering users to fine-tune measurements with precision and reliability. The intuitive layout and clear, easy-to-read displays enable users to obtain accurate measurements swiftly and efficiently.

Whether you're engaged in metalworking, machining intricate components, or overseeing engineering projects, the Beerendonk Caliper 5" 0-80mm stands ready to meet your exacting demands. Moreover, Premium Grade Stainless Steel durable



**6
PCS****WI-70003976**

6 Boley Gauge Dental Implant Lab Orthodontic Instruments

\$47.99 per pieceTrade price from **\$44.15**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Boley Gauge 10cm Measuring 6" (15.24cm) Long Orthodontic Crown Measuring Teeth Size Spaces Vernier Caliper Dental Implant Lab Orthodontics Surgical Instruments

Boley Gauge is one of the most commonly used measuring gauges in dentistry. The Boley Gauge is specialized dental Instruments helps in determining tooth dimensions of length, width, and thickness in millimeter increments. It is used for measuring teeth and spaces. The flat end is good for measuring teeth size and distance of spaces in the mouth.

This dental gauge featuring smooth operation and easy control to enhance user experience helps in a variety of



clinic, laboratory procedures and dental practices

BOLEY GAUGE FEATURES

- Millimeter Increments Measuring Scale Range 100mm
- Small Knob for Adjustment
- Flat End Facilitates Teeth Measurement
- A Boley Gauge Has a Fixed Scale and a Sliding Scale.
- Micro Measuring for More Accuracy
- Ideal For Space Analysis & Band Selection
- Superior Quality Stainless Steel Construction

PACKAGE INCLUDE

- 6 Pcs Boley Gauge 10cm Measuring with 6" (15.24cm) Overall Length

**WI-70003965**

6 Bracket Positioning Height Gauge 3.5mm-5mm for .022" Brackets

\$49.99 per pieceTrade price from **\$45.99**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Bracket Positioning Height Gauge Set 3.5mm – 5.0mm for .022" Brackets Double Ended Premium Quality Aluminum Height Gauge Bracket Positioner Height Placement Dental Orthodontic Tool.

Aluminum Height Gauge featuring double ends with four measurement level tips at 3.5mm, 4mm, 4.5mm, 5mm slotted precisely for .022" brackets. Every Orthodontist must have this gauge during bracket placement procedures. Used to accurately measure height for bracket placement on the posterior and anterior teeth from the insical edge. Made from high-quality, durable Aluminum material, this gauge features a clear, easy-to-read scale that facilitates accurate measurement of bracket height on

each tooth. With precision engineering, our bracket height gauge guarantees optimal bracket positioning for effective treatment outcomes for 022 brackets.

BRACKET POSITIONING HEIGHT GAUGE FEATURES

- 6" (15.2cm) Overall Length.
- Double Ended.
- Four Bracket Height Dimensions 3.5, 4.0, 4.5 and 5.0 Millimeter Gauges.
- .022" Brackets Slot Compatible.
- Pointed Tips for Superficial Marking.
- Long Flat Handle.
- Lightweight.
- Reusable.
- Made with Top Quality Aluminum Metrial.

PACKAGE INCLUDE



**6
PCS****WI-70003975****6 Crown Gauge Iwanson Spring Caliper 4" for Metal Improved Model****\$35.99** per pieceTrade price from **\$33.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Crown Gauge Iwanson Spring Caliper 4" (11.4cm) Sharp Pointed Tips for Metal Crowns Improved Model Stainless Steel Orthodontics Dental Lab Instrument.

Iwanson caliper is used by dentists for measuring the thickness of a metal so that the caliper has a pointed tip. Iwanson style calipers were originally used by dentists to measure metal and wax thickness. However, it can be used just as effectively for measuring smaller items and model engineering in arts, crafts, jewelry, and DIY projects. Made from stainless steel with sharp tip design gives accurate readings. Easy and convenient to use, soft in holding.

IWANSON CALIPER FEATURES

CROWN GAUGE CALIPER: This product is a dental caliper mainly for precision measuring. This Iwanson Caliper has an improved caliber dial, sharp pointed tips and spring action for smooth movement, and a serrated side for better grip.

HIGH-QUALITY STAINLESS STEEL: Made with superior quality durable stainless steel. The Premium Quality Stainless Steel composition of this caliper offers a long life of strong performance, as they are built to last.

FUNCTIONAL ERGONOMIC DESIGN: Features an ergonomic design to complete the task at hand and make your life easy. Provides a comfortable experience and offers more control for precise and easy operation. The Iwanson Crown Gauge measure range is 0-10mm, which can be used for

**6
PCS****WI-70001451****6 Crown Gauge Iwanson Spring Caliper 4" for Wax Improved Model****\$35.99** per pieceTrade price from **\$33.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Crown Gauge Iwanson Spring Caliper 4" (11.4cm) Blunt Blunt Tips For Wax Improved Model Stainless Steel Orthodontics Dental Lab Instrument.

Iwanson Caliper is a dental instrument used to measure the thickness of crowns and spaces. This orthodontic lab instrument is also known as Iwanson Gauge & Iwanson Spring Caliper. Iwanson Spring Caliper for Wax features Blunt tips that won't mark wax specially used for measuring wax crowns and models. High quality medical standard stainless steel product with accurate readings 10 cm length. Calibrated scale capacity 0.1mm to 10.0mm. Marked in 0.1mm increments. Can be sterilized for multiple use.

IWANSON CALIPER FEATURES

CROWN GAUGE CALIPER: Iwanson Caliper mainly for precision measuring. Iwanson Caliper has an improved caliber dial, blunt tips and spring action for smooth movement, and a serrated side for better grip.

HIGH-QUALITY STAINLESS STEEL: Made with superior quality durable stainless steel. The Premium Quality Stainless Steel composition of this caliper offers a long life of strong performance, as they are built to last.

FUNCTIONAL ERGONOMIC DESIGN: Features an ergonomic design to complete the task at hand and make your life easy. Provides a comfortable experience and offers more control for precise and easy operation. The Dental Orthodontic Gauge measure range is 0-10mm, which can be used for precise measurement.

**6
PCS****WI-70004089****6 Dental Bone File #10 Double Ended 7\"/>****\$43.99** per pieceTrade price from **\$40.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Dental Bone File #10 Double Ended 7" (17.8cm) Long with 7mm & 6.5mm Curved Up and Straight Serrated Tips and No-Slip Round Handle Premium Quality Stainless Steel Surgical Dentist Instruments.

Bone file #10 is designed to smooth and shape bone using a pull stroke. Dental bone file #10 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. It can scrape off bony segments from soft tissues and bones. The double-ended end makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 7.0mm and 6.5mm wide. Constructed from surgical-grade stainless steel for excellent performance and longevity.

The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

BONE FILE DOUBE ENDED #10 FEATURES

- Bone File #10.
- 7" (17.8cm) Overall Length.
- 7mm & 6.5mm Serrated Ends.
- Curved Up and Straight Double-Ended.
- Non-Slip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 6 Pcs Bone File Double Ended #10 Serrated Curved Up/Straight Ends

**6
PCS****WI-70001083****6 Hollenback Wax Carver #3 Double Ended Hollow Handle****\$46.99** per pieceTrade price from **\$43.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Hollenback Carver #3 Double Ended Hollow Handle Premium Quality Stainless Steel Dental Wax Modeling Carver Amalgam Restorative Instruments.

Hollenback wax carver used to place, carve and contour plastic dental materials such as amalgam during restorative procedures. Having universal adaptability this hollow round handle hollenback carver 3 is ideal for placing, carving and contouring amalgam. Mainly used for the contour and carve occlusal surfaces & interproximal surfaces. With flat profile hollenback wax carver has spear shape blade that is sharp around the entire periphery. Our hollenback carver known for high degree of precision and flexibility while conducting the clinical

procedure. Made with premium quality stainless steel ensure maxium durability and dependability, last long and can be reused overtime.

HOLLENBACK CARVER #3 FEATURES

- Overall 7.5" (19.1cm).
- Double Ended.
- Large Pointed Spear Shape Sharp Blades.
- Round Shaped Hollow Handle.
- Reusable.
- For Wax and Porcelain Carving
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 6 Pcs Hollenback Carver 3 – Double Ended – Hollow Handle



**WI-70004088****6 Howard Bone File #12 Double Ended****\$43.99** per pieceTrade price from **\$40.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Howard Bone File #12 Double Ended 7" (17.8cm) Straight-Cut 6.4mm Wide Working Tips with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Dental Instruments.

Double Ended Bone File #12 is a straight-cut 6.4mm wide tip dental instrument used for the final smoothing and reshaping of bones before sewing the wound. Straight cut bone files reduce and smooths bone with only a pull stroke due to the parallel position of the cutting edges. The knurled design allows strong grip and easy maneuvering. Made with premium quality stainless instrument for durability and long-lasting performance efficiency. The round handle design helps to reduce the risk of slippage,

making this instrument ideally suited for a variety of precise neurosurgical procedures.

HOWARD BONE FILE #12 FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- 6.4mm Wide Straight-Cut Serrated Ends.
- Double-Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 6 Pcs Howard Bone File #12 Double Ended

TECHNICAL INFORMATION**WI-70003925****6 Lecron Inlay Wax & Porcelain Carver 6.5" Double Ended****\$35.99** per pieceTrade price from **\$33.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 6 LeCron Wax Carver Set 6.5" (16.5cm) Double Ended One Sharp Knife End & Small Scoop End Premium Quality Stainless Steel Wax & Porcelain Carvers for Dental Technicians Dental Laboratory Instruments

Lecron Inlay Wax & Porcelain Carvers #5 is a dental laboratory tool used for shaping and carving wax patterns during the fabrication of dental prosthesis. Double-ended Lecron wax carvers are used to carve and cut wax during impressions and various other dental applications that involve dental wax. Dental technicians can use this to easily shape the wax of fixed and removable prosthodontics into the desired shape. Lecron Inlay Carver #5 features an ergonomically designed lightweight handle with a comfortable

non-slip grip and is made from premium quality stainless steel, making it sterilizable and reusable.

LECRON INLAY CARVER #5 FEATURES

- Crafted with 6.5" (16.5cm) Overall Length.
- Double Ended.
- One Sharp Knife End & One Small Scoop End.
- Ergonomic, Round & Lightweight Non-Slip Grip Handle.
- Maximum Comfort & Minimizing The Risk Of Carpal Tunnel Syndrome.
- Ideal for Carve Wax Blocks and Dentures.
- Sterilizable and Reusable.
- Made with Premium Grade Stainless Steel.



**6
Pieces****WI-70004082****6 Miller Bone File #64 Double Ended****\$46.30** per pieceTrade price from **\$42.60**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

6 Miller Bone File #64 Double Ended is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$46.30, falling to \$42.60 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**6
PCS****WI-70003923****6 P.K Thomas #3 Dental Wax Modelling Carver****\$35.99** per pieceTrade price from **\$33.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



PTK3 – 6 Pcs P.K Thomas #3 Dental Wax Modelling Carver Set with Sharp, Conicular Ponited End Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT3 features a pointed burnisher used to perfect and enhance the supplemental and developmental grooves. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture

when working on hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #3 WAX CARVER FEATURES

- Sharp, Conicular End to Precisely Carve and Shape Anatomical Features.
- Ponited Double-ended
- Non-slip Knurled Handlee for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE

**WI-70003922****6 P.K Thomas #4 Dental Wax Modelling Carver****\$35.99** per pieceTrade price from **\$33.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



PTK4 – 6 Pcs P.K Thomas #4 Dental Wax Modelling Carver with Precise, Spatula-like Ends Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT4 is a modified carver used to perfect the external contours and remove excess wax at the cavo-surface margins. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture

when working on hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #4 WAX CARVER FEATURES

- Use to Shape External Contours and Remove Excess Wax.
- Precise, Spatula Shaped Double Ends.
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE**WI-70003921****6 P.K Thomas #5 Dental Wax Modelling Carver****\$35.99** per pieceTrade price from **\$33.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


PTK5 – 6 Pcs P.K Thomas #5 Dental Wax Modelling Carver with Unique Design Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT5 is a special carver used to remove excess wax as cusp ridges are developed; its contour maintains desired convexity at these ridges. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote

neutral wrist posture when working on hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #5 WAX CARVER FEATURES

- Use to Contour Tooth Anatomy and Remove Excess Wax.
- Unique Carver Tips Design.
- Double Ended.
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.





6 Pieces

WI-70004081

6 Palti Bone Graft Scoop Dental Implant

\$51.46 per piece

Trade price from \$47.34
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

6 Palti Bone Graft Scoop Dental Implant is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$51.46, falling to \$47.34 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



6 Pieces

WI-70003973

6 Ridge Mapping Caliper Surgical Dental

\$87.14 per piece

Trade price from \$80.17
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

6 Ridge Mapping Caliper Surgical Dental is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$87.14, falling to \$80.17 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35080040

7 Inch Curved Surgical Dental Instrument

\$3.50 per piece

Minimum order 250 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture 7 Inch Curved Surgical Dental Instrument for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. The name tells you what you are getting: a curved working end for reach around structures. The 7" size shown is our standard; other sizes in the same pattern are made to order.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.50. Minimum order 250 pieces, and the price moves with quantity.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000060

Acorn Burnisher #21B Double Ended Hollow Handle

\$9.99 per piece

Trade price from \$9.19

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Burnisher #21B Acorn Shaped Double Ended 2.7mm – 3.0mm with Non-Slip Hollow Handle Dental Composite Restoration Instruments.

Burnisher #21B (Acorn Style) is a double-ended dental instrument designed for condensing, contouring, and smoothing amalgam or composite materials. 21B burinsher features acorn-shaped tips are particularly effective for shaping occlusal anatomy in posterior restorations, refining margins, and polishing surfaces for a smooth finish. Non-slip grip hollow handle allow less hand fatigue with, comfortable grip during long procedures. Made from high-grade stainless steel ensure optimal strength, corrosion resistance, and

precision. They are fully autoclavable for reuse.

ACORN BURNISHER #21B HOLLOW FEATURES

- Double-Ended Design.
- 2.7mm & 3.0mm Acorn-shaped ends.
- Non-Slip Grip Hollow Handle.
- Reusable.
- Hand crafted from High Quality stainless steel.

PACKAGE INCLUDE:

- 1 Pc Acorn Burnisher #21B Double Ended – Hollow Handle

TECHNICAL SPECIFICATIONS

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Dental Burnisher

APPLICATION: Composite Placement and Contouring





ACORN BURNISHER 21



WI-70002095

Acorn Burnisher 21 Double Ended with Hollow Handle

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ OEM / private-label marking
- ✓ Reusable, autoclavable

Acorn Burnisher 21 with Hollow Handle Dental Composite Restoration Instruments

Acorn burnisher 21 is a double-ended (2.7 mm – 1.7 mm) dental instrument used to condense, smooth, carve and polish amalgam. Acorn-shaped ends are excellent for carving occlusal anatomy. Acorn Burnishers used to assure that amalgam restorations are free of irregularities and are smooth

ACORN BURNISHER #21 FEATURES

- 1.7mm & 2.7mm acorn-shaped ends.
- Double-Ended Silver Color
- Hand crafted from High Quality stainless steel.
- Available in hollow, ergonomic handle.

PACKAGE INCLUDE:

- 1 Piece Acorn Burnisher 21 Double Ended with Hollow Handle

TECHNICAL SPECIFICATIONS

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Dental Burnisher

APPLICATION: Composite Placement and Contouring

COLOR: Silver

TIP: Double Ended (2.7 mm – 1.7 mm)

TIP TYPE: Acorn

HANDLE: Hollow Handle

FINISH: Polished

REUSABLE: Yes

AUTOCLAVED: Yes

RUST PREVENTION PROCEDURE:

Passivated

ULTRASONIC CLEANED: Yes

TESTS PERFORMED: Boil Test,

Performance Test, Shape Test

WI-70000521

Amalgam Carrier Double Ended 8" (20.3cm) Jumbo/Jumbo Tips

\$10.99 per piece

Trade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ OEM / private-label marking
- ✓ Reusable, autoclavable

Amalgam Carrier Double Ended 8"
(20.3cm) Jumbo/Jumbo Tapered Tips
with Non-Slip Textured Grip Handle
Premium Quality Stainless Steel
Amalgam Filling Dental Restorative
Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument used to fill dental cavities with amalgam. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. The textured surface helps

effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

DOUBLE ENDED AMALGAM CARRIER FEATURES

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Color Bands for Easy Instruments Identification.
- Jumbo – Jumbo (3mm – 3mm) Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.



**WI-70000518****Amalgam Carrier Double Ended 8\"/>****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Amalgam Carrier Double Ended 8\"/>

Dental Amalgam Filling Carrier is a durable, high-quality instrument used to fill dental cavities with amalgam. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. The textured surface helps

effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

DOUBLE ENDED AMALGAM CARRIER FEATURES

- Crafted with 8\"/>
- Double Ended.
- Color Bands for Easy Instruments Identification.
- Large – Jumbo (2.5mm – 3mm) Tapered Tips for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.

**WI-70000519****Amalgam Carrier Double Ended 8\"/>****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Amalgam Carrier Double Ended 8\"/>

Dental Amalgam Filling Carrier is a durable, high-quality instrument used to fill dental cavities with amalgam. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. The textured surface helps

effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

DOUBLE ENDED AMALGAM CARRIER FEATURES

- Crafted with 8\"/>
- Double Ended.
- Color Bands for Easy Instruments Identification.
- Medium – Jumbo (2mm – 3mm) Tapered Tips for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.



**WI-70000520****Amalgam Carrier Double Ended 8" (20.3cm) Medium/Large Tips****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Amalgam Carrier Double Ended 8" (20.3cm) Medium/Large Tapered Tips with Non-Slip Textured Grip Handle Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument used to fill dental cavities with amalgam. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. The textured surface helps

effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

DOUBLE ENDED AMALGAM CARRIER FEATURES

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Color Bands for Easy Instruments Identification.
- Medium – Large (2mm – 2.5mm) Tapered Tips for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.

**100% HIGH QUALITY STAINLESS STEEL****WI-70004264****Amalgam Carrier Double Ended Large 2.8mm /Regular 2mm****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Amalgam Carrier Double Ended 8" (20.3cm) Large/Regular (2.8mm – 2.0mm) Tapered Tips with Non-Slip Textured Grip Handle Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument used to fill dental cavities with amalgam. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. The textured surface helps

effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

DOUBLE ENDED AMALGAM CARRIER FEATURES

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Large – Regular (2.8mm – 2.0mm) Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.

**100% HIGH QUALITY STAINLESS STEEL**

**WI-70004259****Amalgam Carrier Double Ended Large/Medium 2.5mm/2mm****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Amalgam Carrier Double Ended 8" (20.3cm) Large/Medium 2.5mm – 2mm Tapered Tips with Non-Slip Textured Grip Handle Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument used to fill dental cavities with amalgam. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. The textured surface helps

effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

DOUBLE ENDED AMALGAM CARRIER FEATURES

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Color Bands for Easy Instruments Identification.
- Large – Regular (2.5mm – 2mm) Tapered Tips for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.

**WI-70004254****Amalgam Carrier Double Ended Mini 1.5mm / Regular 2mm****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Amalgam Carrier Double Ended 8" (20.3cm) Mini/Regular 1.5mm – 2mm Tapered Tips with Non-Slip Textured Grip Handle Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument used to fill dental cavities with amalgam. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. Amalgam filling carriers act like a syringe, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. The textured surface helps

effectively grip the filling material, preventing it from slipping. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

DOUBLE ENDED AMALGAM CARRIER FEATURES

- Crafted with 8" (20.3cm) Overall Length.
- Double Ended.
- Mini 1.5mm / Regular 2mm Tapered Tips for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.



**WI-70004263****Amalgam Carrier Single Ended 1.5mm Mini Tip****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Amalgam Carrier Single Ended 1.5mm Mini Tapered Tip 6" (15.2cm) Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument for precise and controlled amalgam filling application, available in single-ended or double-ended designs. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. The textured surface helps effectively grip the filling material, preventing it from slipping.

Amalgam filling carriers like a syringes, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental

cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

SINGLE ENDED AMALGAM CARRIER FEATURES

- 6" (15.2cm) Overall Length.
- Single Ended.
- 1.5mm Mini Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.

**WI-70000528****Amalgam Carrier Single Ended 2.0mm Regular Tapered Tip****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Amalgam Carrier Single Ended 2.0mm Regular Tapered Tip 6" (15.2cm) Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument for precise and controlled amalgam filling application, available in single-ended or double-ended designs. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. The textured surface helps effectively grip the filling material, preventing it from slipping.

Amalgam filling carriers like a syringes, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental

cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

SINGLE ENDED AMALGAM CARRIER FEATURES

- 6" (15.2cm) Overall Length.
- Single Ended.
- 2.0mm Regular Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.

**WI-70000525****Amalgam Carrier Single Ended 3.0mm Large Tapered Tip****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Amalgam Carrier Single Ended 3.0mm Large Tapered Tip 6" (15.2cm) Premium Quality Stainless Steel Amalgam Filling Dental Restorative Endodontics Instruments.

Dental Amalgam Filling Carrier is a durable, high-quality instrument for precise and controlled amalgam filling application, available in single-ended or double-ended designs. Amalgam carrier features a tapered tip to distribute the filling material evenly and a curved body for spill prevention. The textured surface helps effectively grip the filling material, preventing it from slipping.

Amalgam filling carriers like a syringes, use to carry and dispense amalgam material in their tubular tip which is inserted and pressed into the dental

cavity by pressing on the lever which activates the piston inside the tube or the back end of the amalgam carrier. We are using high quality stainless steel ensures maximum durability and dependability. You can reuse these amalgam carriers over time as they will not bend or tarnish on you.

SINGLE ENDED AMALGAM CARRIER FEATURES

- 6" (15.2cm) Overall Length.
- Single Ended.
- 3.0mm Large Tapered Tip for Even Distribution of Filling Material.
- Hollow Curved Body to Prevent Material From Spilling Out.
- Non-Slip Textured Surface for Better Grip.

**100% HIGH QUALITY STAINLESS STEEL****WI-70002033****Automatic Crown Remover Gun Type with 5 Attachments****\$21.99** per pieceTrade price from **\$20.23**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

6 Pcs Complete Set of Dental Automatic Crown & Bridge Remover Gun Spring Loaded with 3 Tips 1 Wrench and Assemble Key Restoration Dental Instruments.

Automatic crown and bridge remover is used for removal of temporary and fixed restorations. The automatic crown remover is a spring loaded mechanical device to remove the bonded crowns and bridges by an intensive impulse of the impact stress caused by the stored energy in the spring located inside the instrument. Cement layer breaks with the impulse and the crown or bridge can be removed easily. Packet contain spring loaded action sold with 3 working tips of different styles, a wrench to open it for cleaning if required. It also contain 1 assembly key.


AUTOMATIC CROWN REMOVER SET FEATURES

- Ideally designed for all areas of the mouth.
- Spring-loaded mechanism delivers a quick force that releases the crown
- Comes with three tempered stainless steel tips, one wrench and assemble key
- Made with Highest quality medical grade Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Automatic Crowns Bridges Remover Gun Type with 5 Attachments

TECHNICAL INFORMATION:

**WI-70004282****Ball Burnisher #1 Double Ended Dental Amalgam Instruments****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Ball Burnisher #1 Double 1.3mm & 1.8mm Ball Ends 7" (17.8cm) Long Premium Quality Stainless Steel Dental Amalgam Filling Instruments.

Dental 1 Ball Burnisher with 7" overall length used to condense, carve and polish restorative materials such as amalgam. 1 Ball Burnisher composite amalgam filling instrument is non-stick and comes with ball shaped ends.

Amalgam Ball Burnisher #1 can be used to blend your composite material to the margin for final contouring. All of our composite placement instruments made with stainless steel ensuring excellent performance tools which are autoclavable and can be reused after sterilization.

DENTAL BALL BURNISHER #1

- Crafted With 7" (17.8cm) Overall Length.
- 1.3mm & 1.8mm Ball Ends.
- Double-ended.
- Non-slip Knurled With Center Round Handle.
- Autoclavable & Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Ball Burnisher #1 Double Ended 1.3mm & 1.8mm

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.

**WI-70004112****Ball Burnisher #27/29 DE Dental Amalgam Instrument****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Ball Burnisher #27/29 Double Ended 1.6 mm – 3.1 mm Working Ends Premium Quality Stainless Steel De Amalgam Restorative Composite/Plastic Filling Instruments.

27/29 ball burnisher used to blend material for final contouring, to achieve sculpting of areas like grooves, fissures or pits. Can also be used to form occlusal anatomy. The working end of a burnisher is smooth and rounded (1.6 mm – 3.1 mm). Made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. These dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

BALL BURNISHER # 27/29 FEATURES

- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel

PACKAGE INCLUDE:

- 1 Pc Ball Burnishers #27/29 Double Ended

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Ball Burnisher Dental Instrument

**WI-70004111****Ball Burnisher #28/29 Double Ended 2.3 mm – 3.6 mm Football Tips****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Ball Burnisher #28/29 Double Ended 2.3 mm – 3.6 mm Football Working Tips Premium Quality Stainless Steel Amalgam Restorative Composite/Plastic Filling Instruments.


28/29 football end burnisher is an operative hand instrument for dental professionals used to condense, carve and polish restorative materials such as amalgam. Ball Burnishers can also be used to form occlusal anatomy. The working end of a #28/29 ball burnisher is smooth and rounded like football (2.3 mm – 3.6 mm) made from high grade stainless steel for optimal strength and corrosion resistance. All our amalgam composite filling instruments are manufactured for optimal results and precision. Fully autoclavable for reuse purposes.

BALL BURNISHER # 28/29 FEATURES

- Crafted with 6.5" (16.5cm) Overall Length.
- 2.3mm Football & 3.6mm Football.
- Double Ended.
- Hand Crafted From High-Grade Stainless Steel.
- Long Light Weight Designed.
- Less Hand and Wrist Fatigue.
- Grip Design Handle for Better Handling.
- Reusable.
- Smooth Finish For Excellent Tactile Feel

PACKAGE INCLUDE:

- 1 Pc Foot Ball Burnisher #28/29 Double Ended

**WI-70004110****Ball Burnisher #34 Single Ended Dental Amalgam Instrument****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Ball Burnisher #34 Single Ended 6.5" (16.5cm) Long T-shape Round Pattern Burnisher With a Ball on One Side Premium Quality Stainless Steel Dental Amalgam Composite Filling Instruments.

Single Ended ball burnisher #34 is often used to smooth the surface of a freshly placed amalgam restoration. T-shaped working end with a small ball on one side of "T" and paddle on the other. The shape of burnisher #34 conforms to the rounded cavity surfaces and allows ease of access into the rounded corners or junctions of the cavity surfaces to condense and shape the composite against the cavity wall. All our amalgam composite filling instruments are manufactured for optimal results and precision. These

dental hand instruments are used for long-term duration. Fully autoclavable for reuse purposes.

BALL BURNISHER #34 FEATURES

- Crafted with 6.5" (16.5cm) Overall Length.
- Single Ended.
- Small Pattern.
- T-Shaped Working End.
- Round Pattern Grip Handle for Better Handling.
- Smooth Finish For Excellent Tactile Feel.
- Reusable.
- Hand Crafted From High-Grade Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Ball Burnisher #34 Single Ended 6.5" (16.5cm) – Round Pattern



**WI-70000330****Ball Burnisher Double Ended 1mm-2mm Tips Hollow Handle****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Ball Burnisher Double Ended 1mm – 2mm Ball-Shaped Ends 7" (17.8cm) Long Premium Quality Stainless Steel Dental Restorative Amalgam and Composite Filling Instruments.

Ball Burnisher 1mm-2mm ball ends with 7" overall length used to condense, carve and polish restorative materials such as amalgam & composite. This ball-shaped hand instrument offers precise manipulation for anatomically correct surfaces. Its double-ended 1mm-2mm version provides greater control and accuracy. Burnisher 1mm-2mm specific design allows dentists to sculpt the material, smooth surfaces, and contour the restoration to match the tooth's anatomy, ensuring a precise and aesthetically pleasing finish.


All of our composite placement instruments made with stainless steel ensuring excellent performance tools which are autoclavable and can be reused after sterilization.

DENTAL BALL BURNISHER 1MM – 2MM FEATURES

- Crafted With 7" (17.8cm) Overall Length.
- 1mm & 2mm Ball Ends.
- Double-ended.
- Non-slip Hollow Handle.
- Autoclavable & Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Ball Burnisher Double Ended 1mm-2mm Hollow Handle

**WI-70003949****Beerendonk Caliper 0-80mm Measuring Double Sided Scale****\$12.99** per pieceTrade price from **\$11.95**
Minimum order 21 pieces



- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Beerendonk Caliper 5" (12.7cm) Sliding Gauge 0-80mm Measuring Dental Lab Instruments Premium Grade Stainless Steel Precision Measurement Tool with Adjustable Calibration for Dental Technicians.

Beerendonk Caliper is a high precision dental-calipers with double sided points curved tips for internal & external measurement directly into the mouth. This caliper have intra oral measurements, scale on the back for measurements of sagittal steps and depths, measuring range 0-80mm. Beerendonk caliper double-sided sliding guage used during various types of procedures, especially in dentistry, orthodontics and prosthetics. The basic purpose of the measure is to make precise intraoral measurements,

e.g. when fitting a fixed appliance. For this purpose, the caliper has long, thin tips, bent at the ends, which allow for accurate measurement in hard-to-reach places. It features an adjustable calibration mechanism, empowering users to fine-tune measurements with precision and reliability. The intuitive layout and clear, easy-to-read displays enable users to obtain accurate measurements swiftly and efficiently.

Whether you're engaged in metalworking, machining intricate components, or overseeing engineering projects, the Beerendonk Caliper 5" 0-80mm stands ready to meet your exacting demands. Moreover, Premium Grade Stainless Steel durable construction of the caliper ensures consistent performance over time.



**WI-70003972****Boley Gauge 10cm Measuring 6" (15.24cm)****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Boley Gauge 10cm Measuring 6" (15.24cm) Long Orthodontic Crown Measuring Teeth Size Spaces Vernier Caliper Dental Implant Lab Orthodontics Surgical Instruments

Boley Gauge is one of the most commonly used measuring gauges in dentistry. The Boley Gauge is specialized dental Instruments helps in determining tooth dimensions of length, width, and thickness in millimeter increments. It is used for measuring teeth and spaces. The flat end is good for measuring teeth size and distance of spaces in the mouth.

This dental gauge featuring smooth operation and easy control to enhance user experience helps in a variety of clinic, laboratory procedures and dental practices

BOLET GAUGE FEATURES


- Millimeter Increments Measuring Scale Range 100mm
- Small Knob for Adjustment
- Flat End Facilitates Teeth Measurement
- A Boley Gauge Has a Fixed Scale and a Sliding Scale.
- Micro Measuring for More Accuracy
- Ideal For Space Analysis & Band Selection
- Superior Quality Stainless Steel Construction

PACKAGE INCLUDE


- 1 Pc Boley Gauge 10cm Measuring Scale 6" (15.24cm)

TECHNICAL INFORMATION

BOLEY GAUGE



BOLEY GAUGE

**WI-70004077****Bone File #12A Double Ended****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Bone File #12A Double Ended 7" (17.8cm) 5mm Cross Serrated and Plain Serrated Edge with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

12A Bone File designed to smooth and shape bone using a pull stroke. Their cross-serrated & serrated working ends make it easy to maneuver between soft tissue and bone. Using a surgical bone file to smooth the bone and narrow the canals of the interior of the tooth helps to keep the area disinfected. Round non-slip knurled handle allow strong grip while handling. Hand crafted from premium quality stainless steel, fully autoclaved and sterilized for reuse purpose.

12A BONE FILE FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- Double-ended.
- 1 Cross-serrated & 1 Plain Serrated Ends.
- Non-slip Knurled Grip Round Handle.
- Reusable.
- Hand Crafted Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Dental Bone File #12A Double Ended 7" (17.8cm)

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Bone Files**APPLICATION:** Used for Smoothing and Shaping of Bones**CURVATURE:** Curved**RESUABLE:** Yes

DENTAL BONE FILE HOWARD #12A



**WI-70003947****Bone Ridge Mapping Caliper 6.5" Ring Style with Gold Handle****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Implant Bone Ridge Mapping Caliper 6.5" (16.5cm) Ring Style Measures Upto 0-35mm With Gold Handle Precision Stainless Steel Measurement Tool for Mapping Bone Width Implantology, Dental Restorative Orthodontic Procedures.


Bone Ridge Mapping Caliper is expertly crafted with a ring style that allows for easy handling and accurate measurements. This caliper is designed for dental implant procedures, specifically for mapping bone ridges with precision. The innovative ring-style design, coupled with adjustable arms, offers flexibility and precision, enabling you to tailor measurements to your exact specifications with ease and accuracy. The gold-plated handles not

only provide a luxurious touch but also ensure durability and longevity.

Dental Implant Bone Ridge Mapping Caliper is made of high-quality surgical instruments stainless steel, guaranteeing reliability and accuracy. The caliper is lightweight and easy to use, making it a valuable tool for any dental professional. Whether you are a seasoned dentist or a dental student, this is a must-have tool for ensuring successful dental implant procedures.

BONE RIDGE MAPPING CALIPER FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- Measures Upto 0-35mm.
- Clear and Easy-to-read Markings.

**WI-70000545****Bracket Height Gauge 2.0mm-3.5mm Adjustable Swivel Head****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Bracket Height Gauge 2.0, 2.5, 3.0, 3.5 mm Adjustable Swivel Head 6" (15.2cm) Premium Quality Stainless Steel Wick Positioning Bracket Gauge for Both .018" & .022" Brackets Dental Orthodontic Instruments.

Swivel Head Bracket Height Gauge is an orthodontic tool used to accurately measure bracket height for placement from incisal edge, featuring articulated wick type tips for posterior tooth fixation and adjustable ends for four bracket height dimensions from the occlusal surface. Bracket Height Gauge can be adjusted for easy use on anterior, cuspid, and bicuspid brackets due to its 90° angled hinge on both ends. This adjustable bracket height gauge is an ergonomic double-sided device that ensures proper bracket

placement during orthodontic procedures.

ADJUSTABLE BRACKET HEIGHT GAUGE FEATURES

- 6" (15.2cm) Overall Length.
- Double Ended.
- Four Bracket Height Dimensions 2.0, 2.5, 3.0, and 3.5 Millimeter Gauges.
- Pointed Adjustable Tips for Superficial Marking at Different Angles.
- Ideal for Use on .018" and .022" Bracket Slots.
- Non-Slip Grip Handle.
- Fully Autoclavable / Reusable.
- Made with Premium Grade Stainless Steel.

PACKAGE INCLUDE

**WI-70000544****Bracket Height Gauge 2.5mm-5.0mm
Adjustable Swivel Head****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Bracket Height Gauge 2.5mm-5.0mm Adjustable Swivel Head 6" (15.2cm) Premium Quality Stainless Steel Wick Positioning Bracket Gauge for Both .018" & .022" Brackets Dental Orthodontic Instruments.

Swivel Head Bracket Height Gauge is an orthodontic tool used to accurately measure bracket height for placement from incisal edge, featuring articulated wick type tips for posterior tooth fixation and adjustable ends for four bracket height dimensions from the occlusal surface. Bracket Height Gauge can be adjusted for easy use on anterior, cuspid, and bicuspid brackets due to its 90° angled hinge on both ends. This adjustable bracket height gauge is an ergonomic double-sided device that ensures proper bracket

placement during orthodontic procedures.

**ADJUSTABLE BRACKET HEIGHT
GAUGE FEATURES**

- 6" (15.2cm) Overall Length.
- Double Ended.
- Four Bracket Height Dimensions 2.5, 4.0, 4.5, and 5.0 Millimeter Gauges.
- Pointed Adjustable Tips for Superficial Marking at Different Angles.
- Ideal for Use on .018" and .022" Bracket Slots.
- Non-Slip Grip Handle.
- Fully Autoclavable / Reusable.
- Made with Premium Grade Stainless Steel.

PACKAGE INCLUDE**WI-70000543****Bracket Height Gauge Moveable Head 4.0mm-
6.0mm Pointed Tips****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Bracket Height Gauge Moveable Head 4.0mm-6.0mm Adjustable Pointed Tips 6" (15.2cm) Long Premium Quality Stainless Steel Wick Positioning Bracket Gauge for Both .018" & .022" Brackets Dental Orthodontic Instruments.

Swivel Head Bracket Height Gauge is an orthodontic tool used to accurately measure bracket height for placement from incisal edge, featuring articulated wick type tips for posterior tooth fixation and adjustable ends for four bracket height dimensions from the occlusal surface. Bracket Height Gauge can be adjusted for easy use on anterior, cuspid, and bicuspid brackets due to its 90° angled hinge on both ends. This adjustable bracket height gauge is an ergonomic double-sided

device that ensures proper bracket placement during orthodontic procedures.

**ADJUSTABLE BRACKET HEIGHT
GAUGE FEATURES**

- 6" (15.2cm) Overall Length.
- Double Ended.
- Four Bracket Height Dimensions 4.0, 4.5, 5.0 and 6.0 Millimeter Gauges.
- Pointed Adjustable Tips for Superficial Marking at Different Angles.
- Ideal for Use on .018" and .022" Bracket Slots.
- Non-Slip Grip Handle.
- Fully Autoclavable / Reusable.
- Made with Premium Grade Stainless Steel.

PACKAGE INCLUDE

**WI-70003963****Bracket Positioning Height Gauge 3.5mm-5mm for .022" Brackets****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Bracket Positioning Height Gauge 3.5mm – 5.0mm for .022" Brackets
Double Ended Premium Quality Aluminum Height Gauge Bracket Positioner Height Placement Dental Orthodontic Tool.

Aluminum Height Gauge featuring double ends with four measurement level tips at 3.5mm, 4mm, 4.5mm, 5mm slotted precisely for .022" brackets. Every Orthodontist must have this gauge during bracket placement procedures. Used to accurately measure height for bracket placement on the posterior and anterior teeth from the incisal edge. Made from high-quality, durable Aluminum material, this gauge features a clear, easy-to-read scale that facilitates accurate measurement of bracket height on

each tooth. With precision engineering, our bracket height gauge guarantees optimal bracket positioning for effective treatment outcomes for 022 brackets.

BRACKET POSITIONING HEIGHT GAUGE FEATURES

- 6" (15.2cm) Overall Length.
- Double Ended.
- Four Bracket Height Dimensions 3.5, 4.0, 4.5 and 5.0 Millimeter Gauges.
- .022" Brackets Slot Compatible.
- Pointed Tips for Superficial Marking.
- Long Flat Handle.
- Lightweight.
- Reusable.
- Made with Top Quality Aluminum Material.

PACKAGE INCLUDE**WI-70000231****Burnisher #21 Double Ended Solid Handle****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Burnisher #21 Acorn Shaped Double Ended 1.7mm – 2.7mm with Non-Slip Solid Handle Dental Composite Restoration Instruments.


Double ended Burnisher #21 is a premium dental instrument made of durable, corrosion-resistant stainless steel. Its acorn-shaped tips, measuring 2.7 mm and 1.7 mm, are designed for condensing, smoothing, and shaping composite or amalgam fillings in restorative dental procedures. Non-slip knurled grip solid handle ideal for a firm, comfortable grip. Made from high-grade stainless steel, our amalgam composite filling instruments ensure optimal strength, corrosion resistance, and precision. They are fully autoclavable for reuse.

ACORN BURNISHER #21 FEATURES

- Double-Ended Design.
- 1.7mm & 2.7mm Acorn-shaped ends.
- Non-Slip Knurled Grip Handle.
- Reusable.
- Hand crafted from High Quality stainless steel.

PACKAGE INCLUDE:

- 1 Pc Acorn Burnisher #21 Double Ended

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Burnisher**APPLICATION:** Composite Placement and Contouring**AUTOCCLAVED/REUSABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**

**WI-70003945****Castroviejo Caliper 3.5" (8.89cm) Angled Graduated 0-20mm****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Angled Castroviejo Implant Bone Caliper 3.5" (8.89cm), a Combination Measuring and Marking Caliper Graduated 0-20mm Stainless Steel Surgical Dental Implant Instrument.


This Castroviejo Bone Caliper 8.89cm Angled tips uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-20mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 3-1/2 inches and tapers from the measurement end to the tip. The straight shanks with angled jaws.

CASTROVIEJO CALIPER FEATURES

- 0 to 20mm marking scale for small & wide measurements.
- Overall length 8.89cm
- Pointed tips that open up for precise measurement.
- Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 1 Pc Castroviejo Bone Caliper 3.5" (8.89cm) Angled Graduated 0-20mm

TECHNICAL INFORMATION**WI-70003942****Castroviejo Caliper 3.5" Graduated 0-40mm Angled****\$12.99** per pieceTrade price from **\$11.95**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Angled Castroviejo Implant Bone Caliper 3.5" (8.89cm), a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.


This Castroviejo Bone Caliper 8.89cm Angled tips uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 3-1/2 inches and tapers from the measurement end to the tip. The straight shanks with angled jaws.

CASTROVIEJO CALIPER FEATURES

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 8.89cm
- Pointed tips that open up for precise measurement.
- Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 1 Pc Castroviejo Bone Caliper 3.5" (8.89cm) Angled Graduated 0-40mm

TECHNICAL INFORMATION

**WI-70003941****Castroviejo Caliper 3.5\"/>****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Castroviejo Implant Bone Caliper 3.5\"/>

This Castroviejo Bone Caliper Straight 8.89Cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-20mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 3-1/2 inches and tapers from the measurement end to the tip. The shanks and jaws are straight.

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 8.89cm
- Pointed tips that open up for precise measurement.
- Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 1 Pc Castroviejo Caliper 3.5\"/>

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Implant Instruments**CASTROVIEJO CALIPER FEATURES****WI-70001455****Castroviejo Caliper 6\"/>****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Castroviejo Implant Caliper 6\"/>

This Long Castroviejo Caliper Round Tips 15.2cm used to measure bone and teeth or spaces between teeth. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 6 inches and tapers from the measurement end to the tip.

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 15.2cm
- Round Pointed tips that open up for precise measurement.
- Ridge Mapping V Shape Round Tips.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 1 Pc Long Castroviejo Caliper 6\"/>

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Implant Instruments**APPLICATION:** Dental Implant,**LONG CASTROVIEJO CALIPER FEATURES**

**WI-70001469****Castroviejo Caliper Long 0-40mm Angled 6"****\$16.99** per pieceTrade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Castroviejo Implant Bone Caliper 6" (15.2cm) Long Model – 30° Angled Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper Angled 15.2cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 6 inches and tapers from the measurement end to the tip. The shanks are straight with 30° Angled jaws

LONG CASTROVIEJO CALIPER FEATURES

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 15.2cm
- 30° Angled Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 1 Pc Castroviejo Caliper 6" (15.2cm) Long Model – 30° Angled Tip Graduated 0-40mm

TECHNICAL INFORMATION**WI-70001460****Castroviejo Caliper Long 0-40mm Angled 7"****\$18.99** per pieceTrade price from **\$17.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Castroviejo Implant Bone Caliper 7" (17.8cm) Long Model – 30° Angled Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.

This Long Castroviejo Bone Caliper Angled 17.8cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 7 inches and tapers from the measurement end to the tip. The shanks are straight with angled jaws.


LONG CASTROVIEJO CALIPER FEATURES

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 17.8cm
- 30° Angled Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 1 Pc Castroviejo Caliper 7" (17.8cm) Long Model – 30° Angled Tip Graduated 0-40mm

TECHNICAL INFORMATION

**WI-70003948****Castroviejo Caliper Long 0-40mm Straight 6"****\$16.99** per pieceTrade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Castroviejo Implant Bone Caliper, Long Model 6" (15.2cm) Straight Tips, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.


This Castroviejo Bone Caliper Straight 15.2cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 6 inches and tapers from the measurement end to the tip. The shanks and jaws are straight.

LONG CASTROVIEJO CALIPER FEATURES

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 15.2cm
- Straight Pointed tips that open up for precise measurement.
- Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 1 Pc Castroviejo Caliper Long 6" (15.2cm) Straight Graduated 0-40mm

TECHNICAL INFORMATION**WI-70001464****Castroviejo Caliper Long 0-40mm Straight 7"****\$18.99** per pieceTrade price from **\$17.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Castroviejo Implant Bone Caliper 7" (17.8cm) Long Model – Straight Tip, a Combination Measuring and Marking Caliper Graduated 0-40mm Stainless Steel Surgical Dental Implant Instrument.


This Long Castroviejo Bone Caliper Straight 17.8cm uses accurate calibrations to measure bone and other oral structures precisely. Can be used in surgical, restorative or orthodontic procedures. This caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 7 inches and tapers from the measurement end to the tip. The shanks and jaws are straight.

LONG CASTROVIEJO CALIPER FEATURES

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 17.8cm
- Straight Pointed tips that open up for precise measurement.
- Available in Straight & Angled pattern to suit all needs.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 1 Pc Castroviejo Caliper 7" (17.8cm) Long Model – Straight Tip Graduated 0-40mm

TECHNICAL INFORMATION

**WI-70003994****Cement Spatulas 7" (17.8cm) Double Ended****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Cement Spatula Double Ended 7" (17.8cm) Lab Mixing Modeling Premium Quality Stainless Steel Dental Restorative Instruments.

Cement spatulas used for mixing cement materials and various other dental materials like silicones, waxes, cements etc. Dental spatula comes in double ends, one tip is rounded for scooping, mixing and the other is arch headed spreading tip used to mix, spread and lift material including foods, drugs, plaster and paints.

Dental cement spatulas available in double and single ended patterns with flexible and rigid forms which are ideal for lighter, medium to heavier cements and other restorative materials. These spatulas can also be used to load materials into crowns or other

preparations. Made with premium quality stainless steel with superior craftsmanship. Provides a reliable performance for dental professionals over time. Fully Autoclavable for hygiene and reuse purpose.

DOUBLE ENDED CEMENT SPATULAS FEATURES

- Crafted With 7" (17.8cm) Overall Length.
- Double Ended.
- One Rounded End and Other is Arch Headed.
- Optimal Size for Handling for Optimal Maneuvering.
- Non-slip Textured Area for Better Grip
- Highly Polished Rounded Handle
- Fully Autoclave and Reusable

WI-70003991**Cement Spatulas Single Ended****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Cement Spatula Single Ended 7" (17.8cm) Lab Mixing Modeling Stainless Steel Restorative Instruments.

Cement spatulas are used for mixing cement materials and various other dental materials like silicones, waxes, cements etc. A spatula is a broad, flat, flexible blade used to mix, spread and lift material including foods, drugs, plaster and paints. Available in double and single ended patterns with flexible and rigid forms which are ideal for lighter, medium to heavier cements and other restorative materials. Cement spatulas can also be used to load materials into crowns or other preparations. Made with premium quality stainless steel with superior craftsmanship. Cement spatulas can be reused over time as they will not bend

or tarnish on you. We offer 100% satisfaction guaranteed.

SINGLE ENDED CEMENT SPATULAS FEATURES

- Crafted With 7" (17.8cm) Overall Length.
- Single Ended.
- Large Width Spatula Design.
- Optimal Size for Handling for Optimal Maneuvering
- Non-slip Textured Area for Better Grip
- Highly Polished Rounded Handle
- Fully Autoclave and Reusable
- Made With High Quality Stainless Steel

PACKAGE INCLUDE

**WI-70002029****Christensen Crown Remover CRCH2 90 Degree Angled****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Christensen Crown Remover CRCH2 90 Degree Angled with Hexagonal Dental Elevator Handle Permanent Crown Removal Restoration Dental Instruments


For permanent removal of crowns by breaking the seal between tooth and crown after sectioning with a bur. Crown Remover Christensen is ideally constructed for the disassembly method of crown removal. After sectioning the crown with a bur, the elevator handle provides a comfortable grip and maximum leverage. The notched tip then aids the single-end Christensen Crown Remover in effortless removal. Pressure against the tooth is lessened, which reduces the potential for tooth fracture.

CHRISTENSEN CROWN REMOVER FEATURES

- The 90-degree angled CRCH2 Crown Remover is ideally suited for posterior teeth.
- Elevator handle for optimal comfort and grip
- Notched tips
- Made with Highest quality medical grade Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Christensen Crown Remover CRCH2 90 Degree Angled with Hexagon Handle

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Restorative Dental Instruments**APPLICATION:** Permanent Crown**WI-70002030****Christensen Crown Remover Straight – CRCH1****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Christensen Crown Remover Straight (CRCH1) with Hexagonal Dental Elevator Handle Crown and Bridge Removal Restoration Dental Instruments

For permanent removal of crowns by breaking the seal between tooth and crown after sectioning with a bur. Crown Remover Christensen is ideally constructed for the disassembly method of crown removal. After sectioning the crown with a bur, the elevator handle provides a comfortable grip and maximum leverage. The notched tip then aids the single-end Christensen Crown Remover in effortless removal. Pressure against the tooth is lessened, which reduces the potential for tooth fracture.

CHRISTENSEN CROWN REMOVER FEATURES

- The Straight tip Crown Remover CRCH1 ideal for anterior teeth.
- Elevator handle for optimal comfort and grip
- Notched tips
- Made with Highest quality medical grade Stainless Steel

PACKAGE INCLUDE:

- 1 Pc CRCH1 Christensen Crown Remover – Straight

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Restorative Dental Instruments**APPLICATION:** Permanent Crown Removal**TIPS SHAPE:** Straight – Notched Tips

**WI-70004109**

Cone Burnisher #21B Double Ended Solid Handle

\$7.99 per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Burnisher #21B Acorn Shaped Double Ended 2.7mm – 3.0mm with Non-Slip Solid Handle Dental Composite Restoration Instruments.

Burnisher #21B (Acorn Style) is a double-ended dental instrument designed for condensing, contouring, and smoothing amalgam or composite materials. 21B burinsher features acorn-shaped tips are particularly effective for shaping occlusal anatomy in posterior restorations, refining margins, and polishing surfaces for a smooth finish. Non-slip knurled grip solid handle ideal for a firm, comfortable grip. Made from high-grade stainless steel, our amalgam composite filling instruments ensure optimal strength, corrosion resistance, and precision. They are fully autoclavable for reuse.

ACORN BURNISHER #21B FEATURES

- Double-Ended Design.
- 2.7mm & 3.0mm Acorn-shaped ends.
- Non-Slip Knurled Grip Handle.
- Reusable.
- Hand crafted from High Quality stainless steel.

PACKAGE INCLUDE:

- 1 Pc Acorn Burnisher #21B Double Ended

TECHNICAL SPECIFICATIONS

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Dental Burnisher**APPLICATION:** Composite Placement and Contouring**AUTOCCLAVED/REUSABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:****WI-70003971**

Crown Gauge Iwanson Spring Caliper 4" for Metal Improved Model

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Crown Gauge Iwanson Spring Caliper 4" (11.4cm) Sharp Pointed Tips for Metal Crowns Improved Model Stainless Steel Orthodontics Dental Lab Instrument.

Iwanson caliper is used by dentists for measuring the thickness of a metal so that the caliper has a pointed tip. Iwanson style calipers were originally used by dentists to measure metal and wax thickness. However, it can be used just as effectively for measuring smaller items and model engineering in arts, crafts, jewelry, and DIY projects. Made from stainless steel with sharp tip design gives accurate readings. Easy and convenient to use with its 1-10mm calibrated scale and soft in holding.

IWANSON CALIPER FEATURES

CROWN GAUGE CALIPER: This product is a dental caliper mainly for precision measuring. This Iwanson Caliper has an improved caliber dial, sharp pointed tips and spring action for smooth movement, and a serrated side for better grip.

HIGH-QUALITY STAINLESS STEEL: Made with superior quality durable stainless steel. The Premium Quality Stainless Steel composition of this caliper offers a long life of strong performance, as they are built to last.

FUNCTIONAL ERGONOMIC DESIGN: Features an ergonomic design to complete the task at hand and make your life easy. Provides a comfortable experience and offers more control for precise and easy operation. The Iwanson Crown Gauge measure range is 0-10mm, which can be used for



**WI-70001453****Crown Gauge Iwanson Spring Caliper 4" for Wax Improved Model****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Crown Gauge Iwanson Spring Caliper 4" (11.4cm) Blunt Blunt Tips For Wax Improved Model Stainless Steel Orthodontics Dental Lab Instrument.

Iwanson Caliper is a dental instrument used to measure the thickness of crowns and spaces. This orthodontic lab instrument is also known as Iwanson Gauge & Iwanson Spring Caliper. Iwanson Spring Caliper for Wax features Blunt tips that won't mark wax product with accurate readings 10 cm length. Calibrated scale capacity 0.1mm to 10.0mm. Marked In 0.1mm increments. Can be sterilized for multiple use.

CROWN GAUGE CALIPER: Iwanson Caliper mainly for precision measuring. Iwanson Caliper has an improved caliber dial, blunt tips and spring action for smooth movement, and a serrated side for better grip.

HIGH-QUALITY STAINLESS STEEL: Made with superior quality durable stainless steel. The Premium Quality Stainless Steel composition of this caliper offers a long life of strong performance, as they are built to last.

FUNCTIONAL ERGONOMIC DESIGN: Features an ergonomic design to complete the task at hand and make your life easy. Provides a comfortable experience and offers more control for precise and easy operation. The Dental Orthodontic Gauge measure range is 0-10mm, which can be used for precise measurement.

IWANSON CALIPER FEATURES**WI-70004255****Dental Amalgam Well with Non-Slip Base****\$14.99** per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dentist Bone Grafting Amalgam Mixing Well Pot Non-Slip Base Stainless Steel Dentist Lab Restorative Dental Surgical Instruments.

Amalgam well pot used to mix and dispense amalgam during restorative procedures. Dentist use this non-slip small bowl shaped container to holds mixed amalgam prior to its being loaded into the amalgam carrier. amalgam well used to hold amalgam until it is used in the restoration. Dental restorative amalgam well can also be used for bone grinding/crushing in dental implant procedures.

DENTAL AMALGAM WELL FEATURES:

AMALGAM WELL POT: Amalgam well have a feature is 3cm diameter, Heavy

well pot with secure base with rubber friction locking ring. Machined from a solid piece of steel with a cone shaped inner well, Smooth cavity for optimal mixing and flat cone tip for easy amalgam carrier transfer.

NON-SLIP RUBBER RING BOTTOM: Non-slip removable rubber 'O' ring and bottom heavily weighted to avoid spillage.

STAINLESS STEEL CONSTRUCTION: The stainless steel construction ensures maximum durability and dependability. You can reuse these amalgam well over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

EASY TO CLEAN: Amalgam well is a user friendly and easy to clean dental tool. Can be wipe it even with a damp

**EASY TO CLEAN**
WIPE IT WITH A DAMP CLOTH**ANTI-SLIP**
ANTI-SLIP RUBBER BASE



WI-25693507

Dental Autoclave Sterilization Cassette Tray

\$14.00 per piece

Trade price from \$8.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Autoclave Sterilization Cassette Tray for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$14.00, falling to \$8.00 at 10 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004079

Dental Bone File #10 Double Ended 7" Curved Up & Straight

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Bone File #10 Double Ended 7" (17.8cm) Long with 7mm & 6.5mm Curved Up and Straight Serrated Tips and No-Slip Round Handle Premium Quality Stainless Steel Surgical Dentist Instruments.

Bone file #10 is designed to smooth and shape bone using a pull stroke. Dental bone file #10 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. It can scrape off bony segments from soft tissues and bones. The double-ended end makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 7.0mm and 6.5mm wide. Constructed from surgical-grade stainless steel for excellent performance and longevity.


The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

BONE FILE DOUBE ENDED #10 FEATURES

- Bone File #10.
- 7" (17.8cm) Overall Length.
- 7mm & 6.5mm Serrated Ends.
- Curved Up and Straight Double-Ended.
- Non-Slip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Bone File Double Ended #10 Serrated Curved Up/Straight Ends





WI-21890483

Dental Bone File #12 Double-Ended Medical Surgical Scaler

\$1.50 per piece

Trade price from \$0.90
Minimum order 8 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Bone File #12 Double-Ended Medical Surgical Scaler — one of 4,770 instruments we manufacture ourselves rather than buy in. The name tells you what you are getting: two working ends on one handle, halving what you carry.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.50, falling to \$0.90 at 8 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004076

Dental Bone File #12CA Cross Serrated Double Ended

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Bone File #12CA Double Ended 7" (17.8cm) Cross Serrated Working Tips with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Dental Instruments.

12CA Bone File is a cross serrated working ends dental instrument used to excavate small spicules of bones from narrow spaces using a pull stroke. This bone file usually maneuvers between the soft tissues and the bones that need scraping. The knurled design allows strong grip and easy maneuvering. Made with premium quality stainless instrument for durability and long-lasting performance efficiency. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for

a variety of precise neurosurgical procedures.

BONE FILE #12CA FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- Cross-serrated Ends.
- Double-Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Bone File #12CA 7" Cross Serrated Double Ended

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Surgical Dental Instruments





WI-70004069

Dental Bone File #33 Double Ended

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Bone File #33 Double Ended 7" (17.8cm) Long with 4.8mm Narrow Pattern Straight Plain Serrated Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

#33 Bone File is designed to smooth and shape bone using a pull stroke. Bone File 33 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. The straight serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 4.8mm wide. Constructed from surgical-grade stainless steel for excellent performance and longevity. The round handle design helps to reduce the risk

of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

BONE FILE #33 FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- 4.8mm Narrow Pattern.
- Double Ended.
- Straight Parallel Sides Plain Serration.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Bone File #33 Straight 4.8mm Double Double Ends

TECHNICAL INFORMATION



WI-70004072

Dental Bone Files

\$15.36 per piece

Trade price from \$14.13
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Dental Bone Files is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$15.36, falling to \$14.13 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000692

Dental Bone Mill Bone Crusher Grinder
Silver/Gold Implants Grafting Instrument

\$39.99 per piece

Trade price from \$36.79
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Bone Mill Crusher Grinder Heavy Duty Stainless Steel Dental Implant Surgery Bone Morselizers Graft Augmentation Instruments – Silver/Gold

Bone Mill or Dental Bone Morselizers (also known as Dental Bone Crushers and Bone Smashers) are used to prepare materials for guided tissue regeneration. Dental Bone Crusher Mill Grinder primarily used to grind blocks of bone into smaller particle sizes needed for grafting procedures for small practices of implant surgery. Dental bone mill is essential tool for any dental settings that frequently perform grafting procedures. The stainless steel construction ensures maximum durability and dependability. You can reuse this bone crusher over

time as they will not bend or tarnish on you.

DENTAL MILL BONE CRUSHER GRINDER FEATURES

- Size: 43mm x 30mm.
- Silver and Gold Color.
- Cylindrical Cutting Surface.
- Durable Steel Drum.
- Non-slip Tension.
- Rotating Knobe.
- Reusable.
- Fully Autoclaved.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 1 Pc Dental Bone Morselizers (Bone Crusher) – Silver/Gold

TECHNICAL INFORMATION



WI-43373329

Dental Cement Mixing Spatula Restorative
Laboratory Instrument

\$2.59 per piece

Trade price from \$2.15
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Cement Mixing Spatula Restorative Laboratory Instrument for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what

separates a good instrument from a cheap one.

Listed at \$2.59, falling to \$2.15 at 1 piece — about 17% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-43373326

Dental Cement Spatula Double Ended

\$2.59 per piece

Trade price from \$2.15
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Cement Spatula Double Ended — one of 4,770 instruments we manufacture ourselves rather than buy in. The name tells you what you are getting: two working ends on one handle, halving what you carry. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$2.59, falling to \$2.15 at 1 piece — about 17% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-43347847

Dental Cement Spatula Double Ended
Stainless Steel Mixing Modeling

\$2.59 per piece

Trade price from \$2.15
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Cement Spatula Double Ended Stainless Steel Mixing Modeling is a dental instrument we forge, grind and finish in our own Sialkot workshop. The name tells you what you are getting: two working ends on one handle, halving what you carry.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.59, falling to \$2.15 at 1 piece — about 17% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29044198

Dental Edentulous Stainless Steel Perforated High Precision Denture

\$4.00 per piece

Trade price from \$2.00

Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dental Edentulous Stainless Steel Perforated High Precision Denture, made for dental, hygiene and oral surgery work. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36933034

Dental Flat Plastic Filling Instruments

\$8.00 per piece

Trade price from \$5.00

Minimum order 5 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Flat Plastic Filling Instruments is part of our dental instruments range, produced start to finish at our Sialkot facility. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge.

Listed at \$8.00, falling to \$5.00 at 5 sets — about 38% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34335291

Dental Flex Periotomes Root Extraction Screw Extracting

\$120.00 per piece

Trade price from **\$70.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

We manufacture Dental Flex Periotomes Root Extraction Screw Extracting for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$120.00, falling to \$70.00 at 1 piece — about 42% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003918

Dental Hollenback Wax Carver #3

\$7.99 per piece

Trade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hollenback Carver #3 Double Ended Spear Shape Sharp Blades with 6.5" Overall Length Premium Quality Stainless Steel Wax Modelling Dental Restorative Instruments.

3 hollenback carver is a dental instrument used for placing, carving, and contouring restorative materials like amalgam and composite. Hollenback #3 used for shaping restorative materials to resemble natural tooth anatomy, placing and packing materials into cavities, and trimming excess material after restoration. Double-ended design with large and sharp blades enable detailed work on both the chewing surface and sides of teeth, even in tight spaces. Stainless steel construction of hollenback carver 3 guarantees

durability and dependability, allowing for reuse without bending or tarnishing, and meets international quality standards.

DENTAL HOLLENBACK CARVER #3 FEATURES

- Large Pointed Spear Shape Sharp Blades.
- Double Ended.
- Non-Slip Knurled Grip Handle.
- Used for Wax and Porcelain Carving.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Hollenback #3 Double Ended

TECHNICAL INFORMATION

MATERIAL: Stainless Steel





WI-06303870

Dental Impression Tray

\$4.50 per piece

Trade price from \$2.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Impression Tray — one of 4,770 instruments we manufacture ourselves rather than buy in. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$4.50, falling to \$2.50 at 1 piece — about 44% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29063159

Dental Impression Trays Autoclave

\$4.00 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Impression Trays Autoclave for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29022207

Dental Impression Tray-Stainless Steel Perforated

\$4.00 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Impression Tray-Stainless Steel Perforated for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip

alignment and surface finish are checked on each unit, not on a sample.

Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42924596

Dental Inspection Mirror Hollow Handle Oral Examination Tool

\$5.00 per piece

Trade price from \$0.80
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dental Inspection Mirror Hollow Handle Oral Examination Tool, made for dental, hygiene and oral surgery work. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a

dimension or a finish without a tooling charge.

Listed at \$5.00, falling to \$0.80 at 1 piece — about 84% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42931489

Dental Inspection Mirror Hollow Handle
Replaceable Mirror Instrument

\$6.00 per piece

Trade price from \$2.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Inspection Mirror Hollow Handle Replaceable Mirror Instrument is a dental instrument we forge, grind and finish in our own Sialkot workshop. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$6.00, falling to \$2.00 at 1 piece — about 67% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35452239

Dental Instruments

\$9.00 per piece

Trade price from \$6.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Instruments — one of 4,770 instruments we manufacture ourselves rather than buy in. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge.

Listed at \$9.00, falling to \$6.00 at 1 piece — about 33% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-30574792

Dental Laboratory Articulator Plain Line Ball Bearing Mortimer

\$18.00 per piece Trade price from \$12.00 Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Laboratory Articulator Plain Line Ball Bearing Mortimer — one of 4,770 instruments we manufacture ourselves rather than buy in. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$18.00, falling to \$12.00 at 1 piece — about 33% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-30655901

Dental Magnetic Articulator Adjustable Denture

\$10.00 per piece Trade price from \$7.00 Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Magnetic Articulator Adjustable Denture — one of 4,770 instruments we manufacture ourselves rather than buy in. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$10.00, falling to \$7.00 at 2 pieces — about 30% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36338525

Dental Mouth Inspection Mirror Stainless Steel Oral Examination

\$3.00 per piece

Trade price from \$0.80
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Mouth Inspection Mirror Stainless Steel Oral Examination is part of our dental instruments range, produced start to finish at our Sialkot facility. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and

surface finish are checked on each unit, not on a sample.

Listed at \$3.00, falling to \$0.80 at 25 pieces — about 73% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-13320359

Dental Mouth Mirror with Private Labeling Light Weight

\$2.50 per piece

Trade price from \$0.99
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Mouth Mirror with Private Labeling Light Weight for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a

dimension or a finish without a tooling charge.

Listed at \$2.50, falling to \$0.99 at 10 pieces — about 60% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42932414

Dental Mouth Teeth Inspection Mirror
Reusable Dental

\$3.50 per piece

Trade price from \$0.90
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Mouth Teeth Inspection Mirror Reusable Dental is part of our dental instruments range, produced start to finish at our Sialkot facility. From the pattern: interlocking teeth for a positive hold on tough tissue, and full reusability — clean, autoclave and use again. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$3.50, falling to \$0.90 at 1 piece — about 74% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36853954

Dental Scaler Tray Stainless Steel
Autoclavable Instrument Sterilization
Organizer

\$1.00 per piece

Trade price from \$0.80
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Dental Scaler Tray Stainless Steel Autoclavable Instrument Sterilization Organizer for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by

hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$1.00, falling to \$0.80 at 50 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36877851

Dental Scaler Tray Stainless Steel Sterilization Instrument Tray

\$2.00 per piece

Trade price from \$1.40
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Dental Scaler Tray Stainless Steel Sterilization Instrument Tray — one of 4,770 instruments we manufacture ourselves rather than buy in. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what

separates a good instrument from a cheap one.

Listed at \$2.00, falling to \$1.40 at 25 pieces — about 30% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35293845

Dental Spoon Excavators Double Ended Blue Titanium

\$38.00 per piece

Trade price from \$23.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Dental Spoon Excavators Double Ended Blue Titanium — one of 4,770 instruments we manufacture ourselves rather than buy in. The name tells you what you are getting: two working ends on one handle, halving what you carry.

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$38.00, falling to \$23.00 at 1 set — about 39% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Forging, milling, heat treatment, hand finishing and passivation all happen

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34049315

Design Cassette IMS 1/2 DIN IME12DIN
Orthopedic Dog

\$600.00 per piece

Trade price from **\$597.00**
Minimum order 5 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Design Cassette IMS 1/2 DIN IME12DIN Orthopedic Dog — one of 4,770 instruments we manufacture ourselves rather than buy in. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$600.00, falling to \$597.00 at 5 sets — about 1% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-10236086

Design Stainless Steel Dental Scaler

\$2.50 per piece

Trade price from **\$0.99**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Design Stainless Steel Dental Scaler is a dental instrument we forge, grind and finish in our own Sialkot workshop. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one.

Listed at \$2.50, falling to \$0.99 at 10 pieces — about 60% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-10256215

Double-Sided Dental Mouth Examination Mirror

\$2.50 per piece

Trade price from \$0.99
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Double-Sided Dental Mouth Examination Mirror, made for dental, hygiene and oral surgery work. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one.

Listed at \$2.50, falling to \$0.99 at 10 pieces — about 60% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35450235

First Class Dental Handpiece Crown Removal Tool

\$9.00 per piece

Trade price from \$6.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

First Class Dental Handpiece Crown Removal Tool is part of our dental instruments range, produced start to finish at our Sialkot facility. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$9.00, falling to \$6.00 at 1 piece — about 33% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34318383

Full Medical Stainless Steel Dental Implant

\$120.00 per piece

Trade price from **\$70.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Full Medical Stainless Steel Dental Implant for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why

we can change a dimension or a finish without a tooling charge.

Listed at \$120.00, falling to \$70.00 at 1 piece — about 42% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28829372

German Stainless Steel Medical Instrument Tray Model

\$140.00 per piece

Trade price from **\$70.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of German Stainless Steel Medical Instrument Tray Model, made for dental, hygiene and oral surgery work. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$140.00, falling to \$70.00 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29036187

German Stainless Steel NON-PERFORATED

\$4.00 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture German Stainless Steel NON-PERFORATED for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000055

Graft Packer Large 6mm-8mm Serrated Tip Hollow Handle

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Graft Packer Large 6mm-8mm Serrated Tip with Silver Hollow Handle – Premium Quality Dental Implant Bone Grafting Instruments.


Bone graft packer with serrated tips used in dental and oral surgery to compact and condense bone graft material into precise sites, like extraction sockets. Their serrated ends prevent slippage, enable tight packing to eliminate air pockets, and facilitate access to narrow spaces without harming surrounding tissue. Additionally, they assist in shaping bone during ridge augmentation procedures. The dual bone grafting instrument features 6.0mm and 8.0mm serrated ends designed for effective bone packing and compacting. It has a lightweight hollow handle, is made from

100% premium quality stainless steel for enhanced durability, and is fully autoclavable for repeated clinical use.

BONE GRAFT PACKER FEATURES

- Double Ended
- 6mm-6mm Diameter Serrated Tips.
- Designed for Optimal Access to Hard to Reach Oral Regions.
- Ergonomic, Lightweight Handle Designed to Reduce Hand Fatigue.
- Used for Precise Placement and Packing of Graft Material.
- Ideal for Bone Grafting, Dental Implantology, Periodontal Surgery and Educational Training.
- Fully Compatible with Standard Sterilization Procedures.
- Made with Premium Quality Stainless Steel.



**WI-70000056**

Graft Packer Small 3mm-4mm Serrated Tip Hollow Handle

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Graft Packer Small 3mm-4mm Serrated Tip with Silver Hollow Handle – Premium Quality Dental Implant Bone Grafting Instruments.

Small bone graft packers with serrated tips are used in dental and oral surgery to compact and condense bone graft material into precise sites, like extraction sockets. Their serrated ends prevent slippage, enable tight packing to eliminate air pockets, and facilitate access to narrow spaces without harming surrounding tissue.

Additionally, they assist in shaping bone during ridge augmentation procedures. The dual bone grafting instrument features 3.0mm and 4.0mm serrated ends designed for effective bone packing and compacting. It has a lightweight hollow handle, is made from

100% premium quality stainless steel for enhanced durability, and is fully autoclavable for repeated clinical use.

BONE GRAFT PACKER FEATURES

- Double Ended
- 3mm-4mm Diameter Serrated Tips.
- Designed for Optimal Access to Hard to Reach Oral Regions.
- Ergonomic, Lightweight Handle Designed to Reduce Hand Fatigue.
- Used for Precise Placement and Packing of Graft Material.
- Ideal for Bone Grafting, Dental Implantology, Periodontal Surgery and Educational Training.
- Fully Compatible with Standard Sterilization Procedures.
- Made with Premium Quality Stainless Steel.

**WI-70000324**

Grisdale Bone Condenser Plugger 3.3mm-4mm Double Ended

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Bone Graft Plugger Grisdale Condenser 3.3mm – 4mm Flat Double-Ended Extra Large Bone Graft Packer with 8-11-13-15-18 mm Markings – Premium Quality Dental Implant Bone Grafting Instruments.

Grisdale bone condenser is a dental surgical instrument, specifically a bone plugger used during bone grafting procedures, such as for dental implants. Its primary function is to compact and place bone graft material accurately and efficiently into a specific area. The double-ended design with tips of various sizes have markings to measure graft placement (8-11-13-15-18 mm) to condense bone graft material and help achieve primary stability for implants. Crafted from high-quality surgical stainless steel, it offers


durability, ergonomic handling, and optimal control for implant procedures.

GRISDALE BONE CONDENSER PLUGGER FEATURES

- Crafter with 7.25" Overall Length.
- 3.3mm-4mm Double Ended Tips.
- 8-11-13-15-18mm Depth Markings Measure for Graft Placement.
- Non-Slip Grip Hollow Handle.
- Used for Precise Placement and Packing of Graft Material.
- Corrosion Resistance.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Grisdale Bone Condenser Plugger 3.3/4mm Double Ended with





WI-06348016

High Efficient Steel/Metal Dental Mirror

\$3.80 per piece

Trade price from **\$1.60**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

High Efficient Steel/Metal Dental Mirror — one of 4,770 instruments we manufacture ourselves rather than buy in. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$3.80, falling to \$1.60 at 10 pieces — about 58% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29011233

High Precision Dental Autoclave Impression Trays

\$4.00 per piece

Trade price from **\$2.00**
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of High Precision Dental Autoclave Impression Trays, made for dental, hygiene and oral surgery work. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one.

Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29081153

High Precision Dental Impression Perforated Large Trays Upper

\$4.00 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

High Precision Dental Impression Perforated Large Trays Upper — one of 4,770 instruments we manufacture ourselves rather than buy in. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29044192

High Precision Dental Impression Trays Upper Lower Teeth

\$4.00 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

High Precision Dental Impression Trays Upper Lower Teeth — one of 4,770 instruments we manufacture ourselves rather than buy in. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29057170

High Precision Stainless Steel

\$4.00 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

High Precision Stainless Steel is a dental instrument we forge, grind and finish in our own Sialkot workshop. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29047203

High Precision Stainless Steel Autoclave Impression Trays

\$4.00 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

High Precision Stainless Steel Autoclave Impression Trays — one of 4,770 instruments we manufacture ourselves rather than buy in. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what

separates a good instrument from a cheap one.

Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-01090818

High Precision Stainless Steel Dental

\$4.00 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

High Precision Stainless Steel Dental is a dental instrument we forge, grind and finish in our own Sialkot workshop. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29057171

High Precision Stainless Steel Dental Impression Trays

\$4.00 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

High Precision Stainless Steel Dental Impression Trays — one of 4,770 instruments we manufacture ourselves rather than buy in. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29035205

High Precision Stainless Steel Reusable Dental Orthodontic

\$4.00 per piece

Trade price from **\$2.00**

Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


High Precision Stainless Steel Reusable Dental Orthodontic — one of 4,770 instruments we manufacture ourselves rather than buy in. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001043

Hollenback Wax Carver #3 Blue Tips with Blue & Silver Hollow Handle

\$10.99 per piece

Trade price from **\$10.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hollenback Wax Carver #3 Double Ended Titanium Coated Blue Tips with Blue & Silver Hollow Handle Premium Quality Stainless Steel Dental Wax Modeling Carver Amalgam Restorative Instruments.

Wax Carver "Hollenback #3" used to place, carve and contour dental plastic materials such as amalgam during restorative procedures. Hollenback #3 feature double ended design combines a sharp stiff metal blades pointed ends are protrude at differen tangles carves other restoratvie materials.

Blue titanium coated tips gauranteed not to rust for the life of the instrument. The ergonomically designed rounded hollow handles makes this instrument super light weight and ensures the most comfortable no-slip grip with


minimum fatigue. All our amalgam restorative instruments are highly polished stainless steel dental instruments is sterilizable and can be reused.

HOLLENBACK CARVER #3 FEATURES

- Crafted 7" (17.8cm) Overall Length.
- Large Pointed Spear Shape Sharp Blades.
- Double Ended Titanium Coated Blue Tips.
- Round Shaped Blue & Silver Hollow Handle.
- Reusable.
- Used for Wax and Porcelain Carving.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE



**WI-70001085**

Hollenback Wax Carver #3 Double Ended Hollow Handle

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Hollenback Carver #3 Double Ended Hollow Handle Stainless Steel Dental Wax Modeling Carver Amalgam Restorative Instruments.

Hollenback wax carver used to place, carve and contour plastic dental materials such as amalgam during restorative procedures. Having universal adaptability this hollow round handle hollenback carver 3 is ideal for placing, carving and contouring amalgam. Mainly used for the contour and carve occlusal surfaces & interproximal surfaces. With flat profile hollenback wax carver has spear shape blade that is sharp around the entire periphery. Our hollenback carver known for high degree of precision and flexibility while conducting the clinical procedure. Made with premium quality

stainless steel ensure maxium durability and dependibility, last long and can be reused overtime.


HOLLENBACK CARVER #3 FEATURES

- Overall 7.5" (19.1cm).
- Double Ended.
- Large Pointed Spear Shape Sharp Blades.
- Round Shaped Hollow Handle.
- Reusable.
- For Wax and Porcelain Carving
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Hollenback Carver 3 – Double Ended – Hollow Handle

TECHNICAL INFORMATION

**WI-42948395**

Hollow Handle Dental Mouth Inspection Mirror Stainless Steel Hygiene

\$5.00 per pieceTrade price from **\$0.90**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hollow Handle Dental Mouth Inspection Mirror Stainless Steel Hygiene is part of our dental instruments range, produced start to finish at our Sialkot facility. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface

finish are checked on each unit, not on a sample.

Listed at \$5.00, falling to \$0.90 at 1 piece — about 82% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42963158

Hollow Handle Dental Mouth Mirror with Replaceable Mirror

\$4.00 per piece

Trade price from \$2.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Hollow Handle Dental Mouth Mirror with Replaceable Mirror, made for dental, hygiene and oral surgery work. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a

dimension or a finish without a tooling charge.

Listed at \$4.00, falling to \$2.00 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004078

Howard Bone File #12 Double Ended

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Howard Bone File #12 Double Ended 7" (17.8cm) Straight-Cut 6.4mm Wide Working Tips with Non-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Dental Instruments.

Double Ended Bone File #12 is a straight-cut 6.4mm wide tip dental instrument used for the final smoothing and reshaping of bones before sewing the wound. Straight cut bone files reduce and smooths bone with only a pull stroke due to the parallel position of the cutting edges. The knurled design allows strong grip and easy maneuvering. Made with premium quality stainless instrument for durability and long-lasting performance efficiency. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for

a variety of precise neurosurgical procedures.

HOWARD BONE FILE #12 FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- 6.4mm Wide Straight-Cut Serrated Ends.
- Double-Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Howard Bone File #12 Double Ended

TECHNICAL INFORMATION



**WI-70000866****Krekeler Sliding Caliper Round Gauge 5" 0-80mm Curved****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces



- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Krekeler Sliding Caliper Round Gauge 5" (12.7cm) Measures up to 0-80mm Curved Pointed Tips with Screw Hold Mechanism Premium Grade Stainless Steel Dental Instruments for Implantology, Surgical, Restorative or Orthodontic Procedures.

this Krekeler caliper is suitable for a wide range of applications from metal fabrication to woodworking and engineering tasks as well. Crafted from stainless steel, this caliper is built to last. Its robust construction ensures longevity and reliability.

KREKELER SLIDING CALIPER ROUND GAUGE FEATURES

- Krekeler Sliding Caliper Round Gauge is 5" (12.7cm) Long
- 0-80mm Measuring Scale.
- Equipped With Easy-to-Read Markings
- Pointed Tips for Accurate Measurement.
- Screw Hold Mechanism
- Narrow Tips Ideal for Implants



Krekeler sliding implant caliper gauge to measure interval and length of implants. This multi purpose krekeler sliding caliper round gauge is specially designed to accurately measure bone and teeth or spaces between teeth and implants. Accurately measure the width and depth of the osteotomy. Facilitate positioning of implants. Krekeler sliding caliper measures 0-80mm range with pointed tips can be used to make markings. The narrow tips ideal for implants. With its versatile functionality,

WI-70000865**Krekeler Sliding Caliper Round Gauge 5" 0-80mm Straight****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Krekeler Sliding Caliper Round Gauge 5" (12.7cm) Measures up to 0-80mm Straight Pointed Tips with Screw Hold Mechanism Premium Grade Stainless Steel Dental Instruments for Implantology, Surgical, Restorative or Orthodontic Procedures.

this Krekeler caliper is suitable for a wide range of applications from metal fabrication to woodworking and engineering tasks as well. Crafted from stainless steel, this caliper is built to last. Its robust construction ensures longevity and reliability.

KREKELER SLIDING CALIPER ROUND GAUGE FEATURES

- Krekeler Sliding Caliper Round Gauge is 5" (12.7cm) Long
- 0-80mm Measuring Scale.
- Equipped With Easy-to-Read Markings
- Pointed Tips for Accurate Measurement.
- Screw Hold Mechanism
- Narrow Tips Ideal for Implants



Krekeler sliding implant caliper gauge to measure interval and length of implants. This multi purpose krekeler sliding caliper round gauge is specially designed to accurately measure bone and teeth or spaces between teeth and implants. Accurately measure the width and depth of the osteotomy. Facilitate positioning of implants. Krekeler sliding caliper measures 0-80mm range with pointed tips can be used to make markings. The narrow tips ideal for implants. With its versatile functionality,



WI-30624645

Lore Enterprises A+ Stainless Steel Adjustable Articulator

\$18.00 per piece

Trade price from \$12.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Lore Enterprises A+ Stainless Steel Adjustable Articulator, made for dental, hygiene and oral surgery work. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$18.00, falling to \$12.00 at 1 piece — about 33% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002038

Manual Crown Remover with 3 Tips Double Ended

\$24.99 per piece

Trade price from \$22.99
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Manual Crown Remover Complete Set with 3 Attachment Tips Crown and Bridge Removers Restorative Dentistry Instruments.

- Comes with 3 tempered stainless steel tips
- Made with Highest quality medical grade Stainless Steel

PACKAGE INCLUDE:

- 1 Crown Remover Handle – Double Ended
- 3 Attachment Tips

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Restorative Dental Instruments

APPLICATION: Permanent and Temporary Crowns and Bridges Removal


TIPS SHAPE: Hook Shape / Curved /Serrated

HANDLE: Hand Held – Double Ended
REUSABLE: Yes

MANUAL CROWN REMOVER SET FEATURES

- For use in all areas of the mouth
- Manual, hammer mechanism
- Double Ended





WI-70000678

Marquette 0/1 Amalgam Plugger Double Ended

\$7.99 per piece

Trade price from \$7.35
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Marquette 0/1 Amalgam Plugger 7" (17.8cm) Serrated Double Ended Premium Quality Stainless Steel Dental Filling Instruments.

#0/1 Marquette is a double ended condenser used to pack restorative materials into cavity preparations with flat and serrated ends. Plugger also know as "Condenser"Is a dental instrument with serrated or flat working ends designed to pack amalgam or other filling materials into cavity preparations during restorative dentistry procedures. These amalgam pluggers/condensers are available in a wide range of combinations with plain or serrated tips. Dental condenser have light weight grip handle for maximum comfort, reduced hand fatigue and increased control. Made of high quality

stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency in restorative dentistry.

MARQUETTE 0/1 AMALGAM PLUGGER FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- Double Ended.
- 1mm – 1.5mm Serrated/Flat Ends.
- Knurled Grip with Center Plain Handle.
- Resuable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Marquette #0/1 Amalgam Plugger – Double Ended



WI-32715105

Medical Dental Surgical Storage Tray for Lab Instruments

\$14.00 per piece

Trade price from \$8.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Medical Dental Surgical Storage Tray for Lab Instruments, made for dental, hygiene and oral surgery work. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$14.00, falling to \$8.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25690521

Medical Surgical Dental Use Stainless Steel Instrument Storage Boxes

\$14.00 per piece

Trade price from **\$8.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Medical Surgical Dental Use Stainless Steel Instrument Storage Boxes, made for dental, hygiene and oral surgery work. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what

separates a good instrument from a cheap one.

Listed at \$14.00, falling to \$8.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-22428104

Medical Surgical Tray Organizer Holder Stainless Steel Dental Instruments

\$8.00 per piece

Trade price from **\$4.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Medical Surgical Tray Organizer Holder Stainless Steel Dental Instruments, made for dental, hygiene and oral surgery work. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-33871516

Medical Surgical Trays 201 Stainless Steel Dental Instruments

\$14.00 per piece

Trade price from **\$8.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Medical Surgical Trays 201 Stainless Steel Dental Instruments, made for dental, hygiene and oral surgery work. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$14.00, falling to \$8.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004073

Miller Bone File #45 Double Ended

\$8.99 per piece

Trade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Miller Bone File #45 Straight-Cut 7" (17.8cm) Long with 5mm Wide Straight & Curved Double Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

Bone file #45 is designed to smooth and shape bone using a pull stroke. Dental bone file #45 is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. It can scrape off bony segments from soft tissues and bones. The serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features straight-cut plain serrations and is 5mm wide. Constructed from surgical-grade stainless steel for excellent


performance and longevity. The round handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

BONE FILE MILLER #45 FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- 5mm Straight-cut Plain Serrated Ends.
- Curved & Straight Double Ended.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Miller Bone File #45 Double Ended 5mm Straight & Cruved Ends





WI-70004068

Miller Bone File #64 Double Ended

\$9.65 per piece

Trade price from **\$8.88**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Miller Bone File #64 Double Ended is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.65, falling to \$8.88 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-13329342

Mini Handle Pakistan Made Dentist Mouth Mirror Bright

\$2.50 per piece

Trade price from **\$0.99**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Mini Handle Pakistan Made Dentist Mouth Mirror Bright is a dental instrument we forge, grind and finish in our own Sialkot workshop. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what

separates a good instrument from a cheap one.

Listed at \$2.50, falling to \$0.99 at 10 pieces — about 60% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46733664

Orthodontic Ligature Applicator Dental Instrument for Accurate Elastic Ligature

\$6.80 per piece

Trade price from \$3.50
Minimum order 25 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Orthodontic Ligature Applicator Dental Instrument for Accurate Elastic Ligature — one of 4,770 instruments we manufacture ourselves rather than buy in. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why

we can change a dimension or a finish without a tooling charge.
Listed at \$6.80, falling to \$3.50 at 25 pieces — about 49% off the single-unit price. Larger quantities are priced on request.
Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003916

P.K Thomas #1 Dental Wax Modeling Carver

\$7.99 per piece

Trade price from \$7.35
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

PTK1 – P.K Thomas #1 Dental Wax Modeling Carver with Elongated Curve on Tapered Tips Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT1 features curved tapered tips used to flow on molten wax. P.K. Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture when working on hard-to-reach areas. Made of high

quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.



P.K THOMAS #1 FEATURES



- Elongated Curve on Tapered Tips.
- Double-ended
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Thomas Waxing Instrument – PKT1

TECHNICAL INFORMATION



**WI-70003915****P.K Thomas #2 Dental Wax Modelling Carver****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

PTK2 – P.K Thomas #2 Dental Wax Modelling Carver with Elongated Curve on Thinner Tips Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT2 features curved on thinner tapered tips used to flow on molten wax. P.K. Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture when working on hard-to-reach areas. Made of high




quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #2 WAX CARVER FEATURES

- Elongated Curve on Thinner Tips Than P.K. Thomas #1.
- Double-ended
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Thomas Waxing Instrument – PKT2

**WI-70003914****P.K Thomas #3 Dental Wax Modelling Carver****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

PTK3 – P.K Thomas #3 Dental Wax Modelling Carver with Sharp, Conicular Ponited End Premium Quality Dental Restorative P.K. Thomas Waxing Instrument



PKT3 features a pointed burnisher used to perfect and enhance the supplemental and developmental grooves. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture

when working on hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #3 WAX CARVER FEATURES

- Sharp, Conicular End to Precisely Carve and Shape Anatomical Features.
- Ponited Double-ended
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE

**WI-70003913****P.K Thomas #4 Dental Wax Modelling Carver****\$7.99** per pieceTrade price from **\$7.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking




PTK4 – P.K Thomas #4 Dental Wax Modelling Carver with Precise, Spatula-like Ends Premium Quality Dental Restorative P.K. Thomas Waxing Instrument

PKT4 is a modified carver used to perfect the external contours and remove excess wax at the cavo-surface margins. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture

when working on hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #4 WAX CARVER FEATURES

- Use to Shape External Contours and Remove Excess Wax.
- Precise, Spatula Shaped Double Ends.
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE**WI-70003912****P.K Thomas #5 Dental Wax Modelling Carver****\$7.99** per pieceTrade price from **\$7.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

PTK5 – P.K Thomas #5 Dental Wax Modelling Carver with Unique Design Premium Quality Dental Restorative P.K. Thomas Waxing Instrument


PKT5 is a special carver used to remove excess wax as cusp ridges are developed; its contour maintains desired convexity at these ridges. P.K Thomas waxing instrument series are perfect for all kinds of waxing procedures. They are designed to use as part of restorative procedures to carve anatomical features and remove access materials. PK Thomas wax and modelling technician carving lab instruments giving you more control and access to tight areas. Ergonomically designed knurled handle provide firm grip and helps to promote neutral wrist posture when working on

hard-to-reach areas. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #5 WAX CARVER FEATURES

- Use to Contour Tooth Anatomy and Remove Excess Wax.
- Unique Carver Tips Design.
- Double Ended.
- Non-slip Knurled Handle for a Firm Grip.
- Polished Finish.
- Easy to Sterilize, Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE

**WI-70000059****P.K. Thomas #3 Dental Wax Carver Gold Tip Hollow Handle****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



P.K. Thomas #3 Dental Wax Modelling Carver with Sharp, Conicular Ponited End Gold Color Coated Tips with Silver & Gold Hollow Handle Premium Quality Dental Restorative P.K. Thomas Waxing Instruments.

P.K. Thomas #3 (PKT-3) is a specialized double-ended dental instrument used for wax modeling and restoration techniques. It features gold-coated cone-shaped tips for sculpting grooves in teeth, and blunt and sharp tips for smoothing and polishing biting surfaces in crowns and prosthetics. The ergonomic silver and color hollow handle reduces hand and wrist fatigue during during meticulous lab work. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with

precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #3 GOLD TIPS WAX CARVER FEATURES

- Double Ended Design.
- Gold Color Coated Sharp, Conicular Ponited Ends.
- Sharp Ends to Precisely Carve and Shape Anatomical Features.
- Ponited Burnisher for Refining Supplemental and Developmental Grooves.
- Non-slip Silver & Gold Color Coated Handle for a Firm Grip.
- PKT3 Ideal for Waxing Procedures and Techniques.
- Easy to Sterilize, Reusable.
- Made with High Quality Stainless Steel.

**WI-70000057****P.K. Thomas #3 Dental Wax Carver Hollow Handle****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

P.K. Thomas #3 Dental Wax Modelling Carver with Sharp, Conicular Ponited End Tips with Silver Hollow Handle Premium Quality Stainless Steel Dental Restorative Waxing Instruments.

P.K. Thomas #3 (PKT-3) is a specialized double-ended dental instrument used for wax modeling and restoration techniques. It features silver-coated cone-shaped tips for sculpting grooves in teeth, and blunt and sharp tips for smoothing and polishing biting surfaces in crowns and prosthetics. The ergonomic hollow handle reduces hand and wrist fatigue during during meticulous lab work. Made of high quality stainless steel for optimum utility, easy to sterilize, and fully autoclavable. Crafted with

precision to maximize efficiency for practicing waxing skills.

P.K THOMAS #3 HOLLOW HANDLE WAX CARVER FEATURES

- Double Ended Design.
- Sharp, Conicular Ponited Ends.
- Sharp Ends to Precisely Carve and Shape Anatomical Features.
- Ponited Burnisher for Refining Supplemental and Developmental Grooves.
- Non-slip Silver Color Coated Handle for a Firm Grip.
- PKT3 Ideal for Waxing Procedures and Techniques.
- Easy to Sterilize, Reusable.
- Made with High Quality Stainless Steel.

PACKAGE INCLUDE



WI-29068196

Pakistan Manufacture

\$4.00 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Pakistan Manufacture — one of 4,770 instruments we manufacture ourselves rather than buy in. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.00, falling to \$2.00 at 2 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004067

Palti Bone Graft Scoop Dental Implant

\$13.10 per piece

Trade price from \$12.05
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Palti Bone Graft Scoop Dental Implant is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$13.10, falling to \$12.05 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000477

Palti Bone Graft Scoop Dental Implant Instrument

\$12.99 per piece

Trade price from \$11.95
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Palti Bone Graft Scoop Double Ended 10×25 mm & 5×21 mm with Non-Slip Grip Handle Premium Quality Dental Implants Surgical Instrument.

This is a surgical tool utilized in bone grafting procedures to transport and apply the graft material to the surgical site during dental surgeries. With Two different sized working ends (10×25 mm & 5×21 mm), our instrument is very versatile and can be used in different areas of the mouth with ease.

Our Palti Bone Graft Scoop features ergonomically designed, lightweight handles for comfortable, no-slip grip. With minimal fatigue and premium quality stainless steel construction, dental implant instruments play a pivotal role in delivering exceptional oral care, enabling dental experts to

achieve optimal dental hygiene and patient satisfaction.

PALTI BONE GRAFT SCOOP FEATURES

- Double Ended 10×25mm & 5×21mm.
- Diamond Pattern Knurled Handle for an Optimum Grip.
- Ergonomic, Round, Hollow & Light Weight Handles
- Maximum Comfort & Minimizing the Risk of Carpal Tunnel Syndrome.
- Different Sized Working Ends for Versatility & Accessibility.
- Versatile Design for Exceptional Oral Care.
- Enable Optimal Dental Hygiene and Patient Satisfaction.
- Sterilizable and Re-usable.
- Durable Stainless Steel Construction.



WI-35310804

Periodontal Montana Jack Stainless Steel Dental Hygiene Instrument

\$9.00 per piece

Trade price from \$1.30
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Periodontal Montana Jack Stainless Steel Dental Hygiene Instrument for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a

dimension or a finish without a tooling charge.

Listed at \$9.00, falling to \$1.30 at 1 piece — about 86% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004066

Polokoff Bone File #1 Double Ended

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Polokoff Bone File #1 Double Ended 7" (17.8cm) Long with 3.4mm Narrow Pattern Straight Serrated Working Ends No-Slip Knurled Grip Round Handle Premium Quality Stainless Steel Surgical Dental Instruments.

#1 Polokoff Bone File is broadly used in surgical procedures to aid in mild bone reshaping using a pull stroke. Polokoff Bone File is 7" long and used for reshaping, smoothing, and removing a small bony object from narrow areas. The straight serrated ends makes it quick and easy to switch to another unused file without having to switch instrumentation. Each end features plain serrations and is 3.4mm wide. Constructed from surgical-grade stainless steel for excellent performance and longevity. The round

handle design helps to reduce the risk of slippage, making this instrument ideally suited for a variety of precise neurosurgical procedures.

POLOKOFF BONE FILE #1 FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- 3.4mm Narrow Pattern.
- Double Ended.
- Straight Parallel Sides Serration.
- Non-Slip Knurled Grip Round Handle.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Polokoff Bone File 1 Straight 3.4mm Narrow Pattern Working Double Ends



WI-35438360

Private Label Dental Crown Remover Attachments Medical

\$9.00 per piece

Trade price from \$6.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Private Label Dental Crown Remover Attachments Medical is part of our dental instruments range, produced start to finish at our Sialkot facility. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can

change a dimension or a finish without a tooling charge.

Listed at \$9.00, falling to \$6.00 at 1 piece — about 33% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44644850

Resin Dental Cement Spatula Stainless Steel Laboratory

\$2.50 per piece

Trade price from **\$2.30**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Resin Dental Cement Spatula Stainless Steel Laboratory is a dental instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$2.50, falling to \$2.30 at 1 piece — about 8% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42944448

Reusable Hollow Handle Dental Mouth Mirror Stainless Steel Diagnostic

\$5.00 per piece

Trade price from **\$0.80**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Reusable Hollow Handle Dental Mouth Mirror Stainless Steel Diagnostic is part of our dental instruments range, produced start to finish at our Sialkot facility. The name tells you what you are getting: full reusability — clean, autoclave and use again. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$5.00, falling to \$0.80 at 1 piece — about 84% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003968****Ridge Mapping Caliper Surgical Dental****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

WAQAR INTERNATIONAL BONE RIDGE MAPPING CALIPER 6" (15.24CM) IMPLANT SURGERY INSTRUMENT

Ridge Mapping Caliper Measures Upto 0-25mm Range 6" (15.24cm) Long Dental Implants Surgery Bone Measuring Instrument

Ridge Mapping Calipers used to measure bone and teeth or spaces between teeth. It can be used in dental surgical, restorative or orthodontic procedures. During Implantology, this caliper provides the surgeon with information which can be required during dental treatment. It is used as measuring instruments in preparation for an implant procedure required for various dental treatments.

BONE RIDGE MAPPING CALIPER FEATURES

- Measures up to 25mm
- 6" length
- Pointed tips for measuring
- Flat finger ring design
- Made with High Quality Stainless Steel

PACKAGE INCLUDES:

- 1 Pc Ridge Mapping Caliper Set Measures upto 0.25mm

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Instruments**APPLICATION:** Implantology**TIPS:** Pointed Measuring Tups**SCALE:** upto -0-25mm (Millimeter) Measuring Scale**WI-30639652****Semi-Adjustable Denture Magnetic Anatomic****\$18.00** per pieceTrade price from **\$12.00**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Semi-Adjustable Denture Magnetic Anatomic is part of our dental instruments range, produced start to finish at our Sialkot facility. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$18.00, falling to \$12.00 at 1 piece — about 33% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-10239292

Single-Sided Dental Mouth Mirror 2023

\$2.50 per piece

Trade price from \$0.99
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Single-Sided Dental Mouth Mirror 2023 is a dental instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$2.50, falling to \$0.99 at 10 pieces — about 60% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-20375053

Size Kidney Shaped Emesis Basin Stainless Steel Surgical

\$4.00 per piece

Trade price from \$2.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Size Kidney Shaped Emesis Basin Stainless Steel Surgical, made for dental, hygiene and oral surgery work. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.00, falling to \$2.50 at 1 piece — about 38% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003959

Spring Divider Caliper 4.5" Dental Lab Tool

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Spring Divider Caliper 4.5" with Flat Legs Stainless Steel Precision Tool Bone for Accurate Measurements During Dental Lab Procedures.

The dental caliper measures the width of metal, wax, and crown size in dental cosmetic procedures. It has a scale or screw for the measuring of sizes in the dental lab. Spring Divider Caliper features two adjustable legs to measure hard-to-measure dimensions, made from carefully selected flat stock with heat treated points for durability, and they pivot on a roller while the bow spring ensures a precise tension setting. Accurate adjustments are made by turning the knurled nut on the threaded rod. The design of these calipers, with their fine adjustment feature, ensures quick and positive

settings when taking or transferring specific measurements. Manufactured from high quality stainless steel with a polished finish.

SPRING DIVIDER CALIPER FEATURES

- Crafted With 4.5" (11.4cm) Overall Length.
- Flat Legs with Sharpened Hard and Pointed Tips.
- Adjustable Lockable Steel Wing Divider.
- Rivet Design to Increase Stability.
- Finely Threaded to Achieve Fine Adjustment.
- Durable and Rustproof.
- Reusable.
- Made with Premium Grade Stainless Steel.



WI-36382864

Stainless Steel Bone Scraper Surgical Instrument for Grafting

\$4.00 per piece

Trade price from \$3.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Bone Scraper Surgical Instrument for Grafting — one of 4,770 instruments we manufacture ourselves rather than buy in. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$4.00, falling to \$3.00 at 1 piece — about 25% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35323822

Stainless Steel Dental Anterior Posterior Scaler Montana

\$9.00 per piece

Trade price from \$1.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Anterior Posterior Scaler Montana is a dental instrument we forge, grind and finish in our own Sialkot workshop. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what

separates a good instrument from a cheap one.

Listed at \$9.00, falling to \$1.50 at 1 piece — about 83% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44664727

Stainless Steel Dental Cement Spatula Filling Instrument

\$2.50 per piece

Trade price from \$2.30
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Cement Spatula Filling Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$2.50, falling to \$2.30 at 1 piece — about 8% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36335510

Stainless Steel Dental Examination Mirror

\$3.00 per piece

Trade price from \$0.80
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Dental Examination Mirror for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$3.00, falling to \$0.80 at 25 pieces — about 73% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25861757

Stainless Steel Dental Examination Mirror Diagnostic Octagon Handle

\$1.50 per piece

Trade price from \$0.70
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Examination Mirror Diagnostic Octagon Handle — one of 4,770 instruments we manufacture ourselves rather than buy in. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is

what separates a good instrument from a cheap one.

Listed at \$1.50, falling to \$0.70 at 1 piece — about 53% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36824967

Stainless Steel Dental Hygiene Scaler Tray Reusable

\$2.00 per piece

Trade price from \$1.30
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Hygiene Scaler Tray Reusable is a dental instrument we forge, grind and finish in our own Sialkot workshop. The name tells you what you are getting: full reusability — clean, autoclave and use again.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.00, falling to \$1.30 at 25 pieces — about 35% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35102909

Stainless Steel Dental Implant Mixing Bowl Size

\$1.50 per piece

Trade price from \$0.90
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Dental Implant Mixing Bowl Size, made for dental, hygiene and oral surgery work. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$1.50, falling to \$0.90 at 10 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-10245159

Stainless Steel Dental Impression Trays

\$2.25 per piece

Trade price from \$0.99
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Impression Trays — one of 4,770 instruments we manufacture ourselves rather than buy in. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge.

Listed at \$2.25, falling to \$0.99 at 10 sets — about 56% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34971261

Stainless Steel Dental Instrument Sterilization Tray

\$6.00 per piece

Trade price from \$4.70
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Instrument Sterilization Tray — one of 4,770 instruments we manufacture ourselves rather than buy in. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$6.00, falling to \$4.70 at 10 pieces — about 22% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36867910

Stainless Steel Dental Instrument Tray Organizer for Clinic Hospital

\$12.00 per piece

Trade price from **\$6.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Instrument Tray Organizer for Clinic Hospital is part of our dental instruments range, produced start to finish at our Sialkot facility. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$12.00, falling to \$6.00 at 25 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27818100

Stainless Steel Dental Instrument Trays Surgical Treatment Boxes

\$8.00 per piece

Trade price from **\$5.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Dental Instrument Trays Surgical Treatment Boxes for dental, hygiene and oral surgery work — made from bar stock to finished instrument in house. Every piece is hand-finished and inspected before it leaves the bench: joint action,

tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$8.00, falling to \$5.00 at 1 piece — about 38% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34105560

Stainless Steel Dental Mouth Mirror Single Side

\$0.04 per piece

Minimum order 8 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Mouth Mirror Single Side is part of our dental instruments range, produced start to finish at our Sialkot facility. Each unit is inspected by hand for alignment and action before packing, and we can

photograph your actual pieces before they ship.

Listed at \$0.04. Minimum order 8 pieces, and the price moves with quantity.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42929441

Stainless Steel Dental Mouth Mirror Teeth

\$3.00 per piece

Trade price from \$0.80
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Dental Mouth Mirror Teeth, made for dental, hygiene and oral surgery work. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.00, falling to \$0.80 at 1 piece — about 73% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32309392

Stainless Steel Dental Scaler-High-quality Surgical Instrument for General

\$4.00 per piece

Trade price from **\$0.90**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Dental Scaler-High-quality Surgical Instrument for General is a dental instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.00, falling to \$0.90 at 10 pieces — about 78% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-40089398

Stainless Steel Dental Surgical Storage Tray for Clinic

\$14.00 per piece

Trade price from **\$8.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Dental Surgical Storage Tray for Clinic, made for dental, hygiene and oral surgery work. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$14.00, falling to \$8.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-28839302

Stainless Steel Instrument Tray for Medical Surgical Dental

\$8.75 per piece

Trade price from **\$6.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Instrument Tray for Medical Surgical Dental, made for dental, hygiene and oral surgery work. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what

separates a good instrument from a cheap one.

Listed at \$8.75, falling to \$6.00 at 1 piece — about 31% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-43364645

Stainless Steel Jacquette Scaler Dental Hygiene Periodontal

\$1.15 per piece

Trade price from **\$0.70**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Jacquette Scaler Dental Hygiene Periodontal is a dental instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$1.15, falling to \$0.70 at 1 piece — about 39% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.



WI-33820635

Stainless Steel Medical Surgical Dental Lab Instruments Holder

\$14.00 per piece

Trade price from \$8.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Medical Surgical Dental Lab Instruments Holder, made for dental, hygiene and oral surgery work. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what

separates a good instrument from a cheap one.

Listed at \$14.00, falling to \$8.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25649494

Stainless Steel Perforated Dental Surgical Instrument Box

\$14.00 per piece

Trade price from \$8.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Perforated Dental Surgical Instrument Box — one of 4,770 instruments we manufacture ourselves rather than buy in. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$14.00, falling to \$8.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27780137

Stainless Steel Perforated Surgical Instrument Box for Dental

\$8.00 per piece

Trade price from \$5.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Perforated Surgical Instrument Box for Dental, made for dental, hygiene and oral surgery work. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample.

Listed at \$8.00, falling to \$5.00 at 1 piece — about 38% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36993453

Stainless Steel Single Ended Beauty Mixing Spatula

\$1.00 per piece

Trade price from \$0.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Single Ended Beauty Mixing Spatula is a dental instrument we forge, grind and finish in our own Sialkot workshop. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$1.00, falling to \$0.50 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27763216

Stainless Steel Surgical Dental Instruments Sterilization Box Cover

\$8.00 per piece

Trade price from \$5.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Stainless Steel Surgical Dental Instruments Sterilization Box Cover is part of our dental instruments range, produced start to finish at our Sialkot facility. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and

surface finish are checked on each unit, not on a sample.

Listed at \$8.00, falling to \$5.00 at 1 piece — about 38% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27742545

Stainless Steel Surgical Instrument Box Perforated Dental Tool Box

\$8.00 per piece

Trade price from \$5.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Surgical Instrument Box Perforated Dental Tool Box, made for dental, hygiene and oral surgery work. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$8.00, falling to \$5.00 at 1 piece — about 38% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25638480

**Stainless Steel Surgical Sterilization Trays
Dental Instrument Boxes**

\$14.00 per piece

Trade price from **\$8.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Sterilization Trays Dental Instrument Boxes is a dental instrument we forge, grind and finish in our own Sialkot workshop. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can

change a dimension or a finish without a tooling charge.

Listed at \$14.00, falling to \$8.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28798338

**Stainless Steel Surgical Tray for Dental Lab
Instruments Holder**

\$8.75 per piece

Trade price from **\$6.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Tray for Dental Lab Instruments Holder is a dental instrument we forge, grind and finish in our own Sialkot workshop. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what

separates a good instrument from a cheap one.

Listed at \$8.75, falling to \$6.00 at 1 piece — about 31% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003970

Star-Like Boone Bracket Positioning Gauge

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Boone Bracket Positioning Gauge With Four Integrated Measurement 3.5, 4.0, 4.5 and 5.0 Millimeter – Premium Quality Stainless Steel Dental Orthodontic Instruments for Orthodontists, Dentists, and Dental Assistants involved in Orthodontic Procedures.

Boone Bracket Positioning Gauge (Star-Like Gauge) is an orthodontic instrument used for precisely measuring and marking the vertical height for braces placement on teeth. It includes multiple tips sized at 3.5mm, 4.0mm, 4.5mm, and 5.0mm, aligning with standard distances from a tooth's biting edge, ensuring consistent bracket bonding during procedures. It also features pointed tips for accurate surface marking on the tooth crown.

Boone gauge facilitates consistent bracket placement on the incisal or occlusal surfaces and promotes uniform bonding heights across the dental arch, aiding proper tooth movement during treatment. Made with premium quality stainless steel ensures sterility for hygienic orthodontic practices.

BOONE BRACKET POSITIONING GAUGE FEATURES

- Four Measurement Options: 3.5, 4.0, 4.5, and 5.0 Millimeters.
- Pointed Tips with Flat Base Portion for Superficial Marking
- Measure Incisal or Occlusal Surface of the Tooth to Determine the Ideal Bracket Position.
- Withstand Repeated Use and Sterilization Procedures.



WI-29080566

Surgical Tooth Crown Remover 3-Tip Bridge Dental Instrument

\$9.00 per piece

Trade price from \$6.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Tooth Crown Remover 3-Tip Bridge Dental Instrument is part of our dental instruments range, produced start to finish at our Sialkot facility. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface

finish are checked on each unit, not on a sample.

Listed at \$9.00, falling to \$6.00 at 1 piece — about 33% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001081

Titanium Coated Hollenback Wax Carver #2 Double Ended

\$10.99 per piece

Trade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Titanium Coated Hollenback Carver #2 Double Ended Hollow Handle Premium Quality Stainless Steel Dental Wax Modeling Carver Amalgam Restorative Instruments.

Hollenback wax carver used to place, carve and contour plastic dental materials such as amalgam during restorative procedures. Dental Amalgam Carvers are commonly used to carve anatomical features into amalgam of the teeth and trim excess materials. To Countour and crave occlusal and interproximal anatomy in amlgam restorations characteristics. Double ended sharp stiff metal blade sharp point ends are protrude at differen tangles carves other restoratvie materials.

Our hollenback carver known for high degree of precision and flexibility while conducting the clinical procedure. Made with premium quality stainless steel ensure maxium durability and dependibility, last long and can be reused overtime.

HOLLENBACK CARVER #2 FEATURES

- Double Ended.
- Blue Color Pointed Spear Shape Sharp Blades.
- Round Shaped Hollow Handle.
- 7.5" (19.1cm) Overall Length.
- Reusable.
- For Wax and Porcelain Carving.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE



WI-34636056

VB Stainless Steel Dental Alloy VeraBond with Be/Dental Alloy

\$2.00 per piece

Trade price from **\$1.50**
Minimum order 15 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

VB Stainless Steel Dental Alloy VeraBond with Be/Dental Alloy is a dental instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$2.00, falling to \$1.50 at 15 pieces — about 25% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001454****Weiss Modified Castroviejo Caliper 6.75\"/>****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Weiss Modified Castroviejo Caliper 6.75" (17cm) Curved Graduated 0-40 Stainless Steel Surgical Dental Instrument.


Castroviejo ridge mapping caliper (Weiss Modified Caliper, Curved) is an ideal implant caliper with curved pointed tips penetrate through soft tissue and measure the original thickness of the ridge. This Ridge Mapping Castroviejo Caliper is available with graduated 0-40mm markings and a screw lock that will lock the instrument in place when the measurement is obtained. The overall length of this caliper is 6.75 inches and tapers from the measurement end to the tip. With curved, box tip design tips this caliper used to measure bone or other oral structures.

WEISS MODIFIED CURVED BONE CALIPER FEATURES

- 0 to 40mm marking scale for small & wide measurements.
- Overall length 17cm
- Curved box tip design pointed tips that open up for precise measurement.
- Ridge mapping V shape.
- Screw hold mechanism.
- Reusable.
- Crafted with German grade stainless steel.

PACKAGE INCLUDE:

- 1 Pc Castroviejo Ridge Mapping Caliper (Weiss Modified Curved) 6.75" (17cm) Curved Tips Graduated 0-40mm

TECHNICAL INFORMATION**WI-70000329****Weston #3 Wax Spatula – Double Ended****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Weston #3 Wax/Cement Spatula 6.75" Double Ended – Premium Quality Stainless Steel Dental Mixing Instrument Suitable for Dentists, Prosthodontists, and Dental Technicians for Restorative and Prosthetic Dental Work.

Cement Spatula Weston #3 is a high-quality dental tool designed for precise manipulation of dental materials, such as temporary cements, filling materials, and impression compounds. Commonly used in both clinical and laboratory settings to support efficient cement mixing and spreading. Weston #3 features versatile use during dental procedures to facilitate seamless transitions between mixing and application tasks. The wide blunt-edge blade feature for efficient blending,

while the sharp chisel edge ensures precision in material application, shaping, or removal. The ergonomically shaped handle ensures a firm grip and controlled application.

Crafted from German-grade stainless steel, is a crucial tool for dentists, prosthodontists, and dental technicians, offering precision, longevity, and ease of sterilization

WESTON #3 DENTAL WAX MIXING CEMENT SPATULA

- Crafted with 6.75" (17.1cm) Overall Length.
- Double-ended Design.
- Ideal Wax & Modelling Carvers Cement Spatula for Mixing and Applications





WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

Request a quotation

Send us product codes from this catalogue, quantities and your country.
We reply the same working day on WhatsApp.

Online shop	shop.waqarintl.com
Email	sales@waqarintl.com
WhatsApp /	
Phone	+92 307 6173627
Works	Roras Road, Opposite Anwar Khawaja Industries, Sialkot 51310, Pakistan

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.