



# WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

# Instrument Sets

## Orthopedic Instruments

Reusable stainless-steel instrument sets forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



10  
PCS



**102**

PRODUCTS

**\$5.90**

PRICES FROM

**OEM**

PRIVATE LABEL

**Global**

WORLDWIDE  
EXPORT

# Contents

102 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

| Product                                                                   | Code        | Price    | Page |
|---------------------------------------------------------------------------|-------------|----------|------|
| 10 Pcs Bone Reduction Forceps Set 5.5" Curved Pointed Tip Gold Handle     | WI-70001501 | \$159.99 | 5    |
| 10PCS Orthopedic Surgery Set Bone Surgery Instruments                     | WI-36028338 | \$180.00 | 5    |
| 10PCS Orthopedic Surgery Set Stainless Steel Surgical Instruments         | WI-36010323 | \$180.00 | 6    |
| 10PCS Orthopedic Surgical Instrument Set Stainless Steel Orthopedic Tools | WI-36005345 | \$180.00 | 6    |
| 12PCS Orthopedic Surgery Set Stainless Steel Orthopedic                   | WI-36067241 | \$260.00 | 7    |
| 12PCS Orthopedic Surgical Instrument Set Stainless Steel Bone             | WI-36093014 | \$260.00 | 7    |
| 2 Lane Bone Holding Forceps Set 13" & 17" with Ratchet                    | WI-70000632 | \$127.99 | 8    |
| 2 Pcs Bone Reduction Forceps Set 5.5" Curved Pointed Tip Gold Handle      | WI-70001503 | \$49.99  | 8    |
| 2 Pcs Gigli Saw Handle Set                                                | WI-70001813 | \$9.99   | 9    |
| 20PCS Orthopedic Surgery Set                                              | WI-36056304 | \$160.00 | 9    |
| 20PCS Orthopedic Surgery Set Bone Surgery Instruments                     | WI-36088265 | \$160.00 | 10   |
| 20PCS Orthopedic Surgery Set Surgical Instruments                         | WI-35998304 | \$160.00 | 10   |
| 22 Pcs Mueller Cement Removal Instrument Set                              | WI-45575604 | \$550.00 | 11   |
| 22PCS Shoulder Surgery Instrument Set Orthopedic Arthroscopy Surgical     | WI-45916812 | \$600.00 | 11   |
| 25 PCS Basic Orthopedic Surgery Set                                       | WI-35528201 | \$300.00 | 12   |
| 25 Piece Orthopedic Surgery Instrument Set                                | WI-35257162 | \$190.00 | 12   |
| 25PCS Orthopedic Surgery Set Stainless Steel Orthopedic                   | WI-36061249 | \$150.00 | 13   |
| 25PCS Orthopedic Surgical Set Stainless Steel Orthopedic Instruments      | WI-36049244 | \$150.00 | 13   |
| 25-Piece Steel Major Basic Orthopedic Surgery Set Power                   | WI-30873075 | \$190.00 | 14   |
| 3 Pcs Dingman Cartilage Bone Clamp 7½" Serrated 2x2 Teeth Jaws            | WI-70001250 | \$49.99  | 14   |
| 3 Pcs Gigli Saw Wire 12", 20" & 30"                                       | WI-70001814 | \$15.99  | 15   |
| 3 Plate Bone Holding Reduction Forceps Set 5.5" Gold Handle               | WI-70000649 | \$79.99  | 15   |
| 36-Piece Orthopedic Surgical Set Precision Instruments 2-Year             | WI-27305681 | \$250.00 | 16   |
| 5 Kern Bone Holding Forceps Set 5.5", 6", 9", 13" & 17" with Ratchet      | WI-70001790 | \$149.99 | 16   |
| 5 Pcs Bone Reduction Forceps Set 5.5" Curved Pointed Tip Gold Handle      | WI-70001502 | \$99.99  | 17   |
| 5 Pcs Cobb Spinal Elevators Set 10mm – 32mm Wide Blade                    | WI-70000149 | \$89.99  | 17   |
| 5 Pcs Gigli Saw Wire Set 12", 16", 20", 25" & 30"                         | WI-70000179 | \$24.99  | 18   |
| 5 Pcs Plaster & Cast Removal Set                                          | WI-70001137 | \$99.99  | 18   |
| 6 Pcs Dingman Cartilage Bone Clamp 7½" Serrated 2x2 Teeth Jaws            | WI-70001249 | \$79.99  | 19   |
| 7 Pcs Hohmann Retractor Set 6, 8, 10, 18, 25, 43 & 70mm Blades            | WI-70001746 | \$39.99  | 19   |
| 7 Pcs Smith Peterson Chisel Set 8" Straight 0.64cm-3.2cm Blades           | WI-70001122 | \$99.99  | 20   |
| Advanced Orthopedic Biology Lab Dissecting Instrument Set Medical Student | WI-35367738 | \$9.00   | 20   |
| Amputation Instrument Set-03                                              | WI-46492870 | \$248.00 | 21   |
| Amputation Surgery Instrument Set 23 Pieces                               | WI-27858712 | \$150.00 | 21   |
| Amputation Surgical Instrument Set Large Fragment Medical Tools           | WI-27331774 | \$270.00 | 22   |
| Basic for Lumbar Instrument Set Precision Spine Surgery Kit               | WI-45890366 | \$240.00 | 22   |
| Basic for Lumbar Instruments Set Minimally Invasive Spine                 | WI-45870590 | \$240.00 | 23   |
| Basic for Lumbar Surgical Instrument Set                                  | WI-45889397 | \$240.00 | 23   |
| Basic Major Orthopedic Surgical Instruments Set 75 Pieces                 | WI-27324175 | \$180.00 | 24   |
| Basic Major Orthopedic Surgical Instruments Set Stainless Steel Medical   | WI-35467832 | \$30.00  | 24   |
| Basic Orthopedic Surgery Set 25 PCS                                       | WI-24339638 | \$90.00  | 25   |

## Contents (continued)

| Product                                                                       | Code        | Price    | Page |
|-------------------------------------------------------------------------------|-------------|----------|------|
| Cardiac Vascular Thoracic Surgery Instrument Kit Precision Heart Operation    | WI-46197258 | \$80.00  | 26   |
| Complete Basic for Lumbar Instrument Set Minimally Invasive Spine             | WI-45564529 | \$240.00 | 26   |
| Complete Laminectomy Instrument Set 14PCS Precision Spine Surgery Kit         | WI-45882337 | \$300.00 | 27   |
| Complete Mueller Type Cement Removal Set 22PCS Orthopedic Revision            | WI-45886250 | \$550.00 | 27   |
| Complete Orthopedic Surgical Set                                              | WI-35535185 | \$400.00 | 28   |
| Complete Periacetabular Osteotomy Orthopedic Instrument Set 10PCS Hip Surgery | WI-45595260 | \$300.00 | 28   |
| Complete Small Fragment Orthopedic Instrument Set 30PCS Bone Trauma           | WI-45875269 | \$499.00 | 29   |
| Complete Stainless Steel Shoulder Surgery Set 22Pcs Orthopedic Arthroscopy    | WI-45955039 | \$600.00 | 29   |
| DHS Box Medium Orthopedic Trauma Instrument Set Dynamic                       | WI-45858783 | \$150.00 | 30   |
| Excellent Broken Screw Removal Instrument Set Orthopedic                      | WI-28844453 | \$424.00 | 30   |
| Genecology Instruments                                                        | WI-26509143 | \$110.00 | 31   |
| German Stainless Steel Surgical Orthopedic Instrument                         | WI-25752471 | \$50.00  | 31   |
| Gigli Wire Saw Set 25" with 2 Handle                                          | WI-70001819 | \$10.99  | 32   |
| Large Fragment Instrument Set Stainless Steel Orthopedic Surgical             | WI-34729565 | \$200.00 | 32   |
| Large Stainless Steel Pad 14 Pre-Cut Suture Wounds Practice                   | WI-32595825 | \$10.00  | 33   |
| Large Stainless Steel Pad 14 Pre-Cut Wounds Suture Practice                   | WI-32659810 | \$10.00  | 33   |
| Lumbar Spine Surgery Set Stainless Steel Orthopedic Instruments               | WI-32883708 | \$900.00 | 34   |
| Lumbar Spine Surgical Instrument Kit for Spinal Operation Decompression       | WI-46296869 | \$850.00 | 34   |
| Major Orthopedic Surgery Instrument Set Hospital Surgical Tools               | WI-35485713 | \$30.00  | 35   |
| Mueller Type Cement Removal Set of 22 PCS                                     | WI-45872364 | \$550.00 | 35   |
| Orthopedic Cerclage Wire Passing Surgical Kit 3-Piece                         | WI-45862402 | \$120.00 | 36   |
| Orthopedic Hip Knee Spine Tibia Femur                                         | WI-28866024 | \$420.00 | 36   |
| Orthopedic Instrument Set Steel 56PCS Broken Screw Removal                    | WI-34585765 | \$415.00 | 37   |
| Orthopedic Intramedullary Nailing Instrument Set Stainless Steel Surgical     | WI-36932285 | \$699.00 | 37   |
| Orthopedic Small Fragment Set Stainless Steel Locking Plate                   | WI-34735529 | \$200.00 | 38   |
| Orthopedic Stainless Steel Wire Passer Fiber Handle                           | WI-45880190 | \$120.00 | 38   |
| Orthopedic Surgery Instruments Kit Stainless Steel Medical Tools              | WI-35540173 | \$400.00 | 39   |
| Orthopedic Surgical Instrument Set                                            | WI-35395938 | \$30.00  | 39   |
| Orthopedic Surgical Instrument Set 25 PCS Medical                             | WI-35525197 | \$300.00 | 40   |
| Orthopedic Surgical Instrument Set Stainless Steel Bone Surgery               | WI-36088232 | \$128.00 | 40   |
| Orthopedic Surgical Instruments Basic Stainless Steel Re-usable               | WI-25837418 | \$50.00  | 41   |
| Orthopedic Surgical Instruments Solid Handle Stainless Steel Hudson           | WI-27352365 | \$42.00  | 41   |
| Orthopedic Trauma Instrument Set German                                       | WI-36969261 | \$125.00 | 42   |
| Orthopedic Trauma Surgical Instrument Set Bone Fixation Surgery               | WI-36913377 | \$125.00 | 42   |
| Periacetabular Osteotomy Orthopedic Instruments Set 10 Pcs Hip Preservation   | WI-45601048 | \$300.00 | 43   |
| Periacetabular Osteotomy Orthopedic Instruments Set of 10 Pcs                 | WI-45886322 | \$300.00 | 43   |
| Periacetabular Osteotomy Surgical Instrument Set 10 Piece Orthopedic          | WI-45589374 | \$300.00 | 44   |
| Reusable Stainless Steel Orthopedic Trauma Surgery Set Fracture               | WI-36935339 | \$125.00 | 44   |
| Shoulder Joint Arthroscopy Instrument Set Orthopedic Surgical Kit             | WI-35493270 | \$25.00  | 45   |
| Shoulder Surgery Set 22Pcs Orthopedic Shoulder Arthroscopy Surgical           | WI-45937620 | \$600.00 | 45   |
| Small Fragment 3.5 4.0mm Stainless Steel Instrument Set                       | WI-34308005 | \$200.00 | 46   |
| Small Fragment Instrument Set Stainless Steel Orthopedic Surgical             | WI-42812777 | \$500.00 | 46   |

## Contents (continued)

| Product                                                                     | Code        | Price    | Page |
|-----------------------------------------------------------------------------|-------------|----------|------|
| Small Fragment Instrument Set Stainless Steel Orthopedic Trauma             | WI-36177725 | \$500.00 | 47   |
| Small Fragment Instruments Set 30 Pcs Orthopedic Trauma                     | WI-45881225 | \$499.00 | 47   |
| Small Fragment Orthopedic Instrument Set Stainless Steel Bone Fixation      | WI-36131764 | \$500.00 | 48   |
| Small Fragment Surgical Instrument Set Orthopedic Bone Plating              | WI-36168695 | \$500.00 | 48   |
| Spine Fixation Box Spinal Surgery Instrument                                | WI-45963958 | \$199.00 | 49   |
| Spine Fixation Instrument Set Spinal Fusion Surgery                         | WI-45996496 | \$199.00 | 49   |
| Stainless Steel Amputation Surgery Instrument Set for Limb                  | WI-46521683 | \$249.90 | 50   |
| Stainless Steel Amputation Surgical Instrument Kit Orthopedic Surgery Tools | WI-46517704 | \$250.00 | 50   |
| Stainless Steel Bone Mallet Hammer Surgical Orthopedic                      | WI-09914253 | \$11.35  | 51   |
| Stainless Steel Orthopedic Shoulder Arthroscopy Instrument Set              | WI-35554136 | \$25.00  | 51   |
| Stainless Titanium Cement Removal Instruments Set Orthopedic                | WI-45575603 | \$550.00 | 52   |
| Steel Bone Chisel Orthopedic Surgical Kit for Bone Cutting                  | WI-46819581 | \$5.90   | 52   |
| Suture Practice Kit Stainless Steel Medical Instrument Small                | WI-32640778 | \$10.00  | 53   |
| Thoracotomy Surgical Instrument Set Chest Surgery Kit Thoracic Operation    | WI-46524106 | \$399.00 | 53   |
| Trauma Instrument Set Stainless Steel Orthopedic Surgical                   | WI-28792056 | \$280.00 | 54   |
| Wire Passer Fiber Handle Set 3 Pcs Orthopedic                               | WI-45592253 | \$120.00 | 54   |
| Wire Passer Fiber Handle Set 3 Precision Orthopedic                         | WI-45850861 | \$120.00 | 55   |
| Wood Handle Wire Passer Set 3-Piece Orthopedic Surgical Instrument          | WI-45560591 | \$120.00 | 55   |



WI-70001501

10 Pcs Bone Reduction Forceps Set 5.5\"/>

\$159.99 per piece

Trade price from \$147.19  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Bone Holding Reduction Forceps 5.5\"/>

The bone reduction forceps with curved tips allow the gripping of denser and deeper bones. Reduction forceps assist in grasping dense bones located in deeper regions. These Bone Reduction Forceps also helpful for the soft tissue attachment intended to reduce distal radius fractures.

The instrument aids in the alignment and setting of bone fractures while firmly holding them. The ergonomic design of the forceps allows precision and reduces the chances of fracture during the surgical procedure. These forceps are available in different

variations. The long ratchet mechanism gives a precise performance during bone manipulation.

BONE REDUCTION FORCEPS FEATURES:

- Minimal damage during the surgical procedure
- Curved Light Serrated Jaws
- Precise Pointed Tips
- Size and tip variations are available
- German Grade Stainless Steel Finish



PACKAGE INCLUDE:

- 10 Pcs Bone Reduction Forceps 5.5\"/>

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Bone Reduction Forceps



WI-36028338

10PCS Orthopedic Surgery Set Bone Surgery Instruments

\$180.00 per piece

Trade price from \$150.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

10PCS Orthopedic Surgery Set Bone Surgery Instruments is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. The steel is heat treated and

passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$180.00, falling to \$150.00 at 1 set — about 17% off the single-unit price. Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

**WI-36010323****10PCS Orthopedic Surgery Set Stainless Steel Surgical Instruments****\$180.00** per pieceTrade price from **\$150.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

10PCS Orthopedic Surgery Set Stainless Steel Surgical Instruments is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$180.00, falling to \$150.00 at 1 set — about 17% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

**WI-36005345****10PCS Orthopedic Surgical Instrument Set Stainless Steel Orthopedic Tools****\$180.00** per pieceTrade price from **\$150.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

10PCS Orthopedic Surgical Instrument Set Stainless Steel Orthopedic Tools is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$180.00, falling to \$150.00 at 1 set — about 17% off the single-unit price. Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36067241

**12PCS Orthopedic Surgery Set Stainless Steel Orthopedic**

**\$260.00** per piece

Trade price from **\$220.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

12PCS Orthopedic Surgery Set Stainless Steel Orthopedic is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 12-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$260.00, falling to \$220.00 at 1 set — about 15% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36093014

**12PCS Orthopedic Surgical Instrument Set Stainless Steel Bone**

**\$260.00** per piece

Trade price from **\$220.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

This is our pattern of 12PCS Orthopedic Surgical Instrument Set Stainless Steel Bone, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 12-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. We work the instrument by

hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$260.00, falling to \$220.00 at 1 set — about 15% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

**WI-70000632****2 Lane Bone Holding Forceps Set 13" & 17" with Ratchet****\$127.99** per pieceTrade price from **\$117.75**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Pcs Lane Bone Holding Forceps Set 13" & 17" Heavy Serrated 2x2 Teeth Jaws with Bar Ratchet Mechanism Premium Quality Stainless Steel Surgical Orthopedic Instruments


Lane Bone Holding Forceps is a specialized orthopedic Instrument that surgeons commonly use for grasping and mobilizing fractured bone pieces, in order to align them into the anatomical stabilization spot during surgical procedures.

This bone holding forceps comes with a pair of 2x2 toothed jaws with a square-shaped design. Moreover, the jaws have horizontal heavy serration on their inner surfaces for fixating onto the bone surfaces with ease. As a result, surgeons can reposition the desired bone fragments into their original

anatomical place. Moreover Lane Bone Holding Forceps have a plier style handle with a bar ratchet system at the bottom which allows the surgeons to statically lock the instrument at the desired width for controlling the clamping pressure.

**LANE BONE HOLDING FORCEPS FEATURES**

- Available in 13" and 17" Overall Length.
- 2x2 Toothed Tips for a Secure Grasps.
- Atraumatic Jaws Serrations for Avoiding Local Injury.
- Single Spring Action Handle.
- Bar Ratchet Mechanism for Statically Locking the Jaws.

**WI-70001503****2 Pcs Bone Reduction Forceps Set 5.5" Curved Pointed Tip Gold Handle****\$49.99** per pieceTrade price from **\$45.99**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Bone Holding Reduction Forceps 5.5" Curved Serrated Jaws Pointed Tips with Ratcheted Finger Ring Gold Handles Stainless Steel Surgical Orthopedic Instrument.

The bone reduction forceps with curved tips allow the gripping of denser and deeper bones. Reduction forceps assist in grasping dense bones located in deeper regions. These Bone Reduction Forceps also helpful for the soft tissue attachment intended to reduce distal radius fractures.

The instrument aids in the alignment and setting of bone fractures while firmly holding them. The ergonomic design of the forceps allows precision and reduces the chances of fracture during the surgical procedure. These forceps are available in different


variations. The long ratchet mechanism gives a precise performance during bone manipulation.

**BONE REDUCTION FORCEPS FEATURES:**

- Minimal damage during the surgical procedure
- Curved Light Serrated Jaws
- Precise Pointed Tips
- Size and tip variations are available
- German Grade Stainless Steel Finish

**PACKAGE INCLUDE:**

- 2 Pcs Bone Reduction Forceps 5.5" (14cm) Long Ratchet Gold Handles

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Bone Reduction Forceps



WI-70001813

## 2 Pcs Gigli Saw Handle Set

**\$9.99** per piece

Trade price from **\$9.19**  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 GIGLI Saw Handles Surgical Veterinary Neurology Orthopedic Instruments

Gigli saw handle with T-shape pattern is made for amputation practices that require continuous and seamless bone cutting.. Gigli Wire Saw is a specialized surgical device that orthopedic and neurosurgeons can use for cutting bone surgeries during a wide range of procedures, especially during amputation. The Gigli Wire Saw feature small rings to attach the handles. Surgeons can insert T-shaped or loop-shaped handles before starting the bone cutting.

### GIGLI SAW HANDLE FEATURES

**ORTHOPEDIC GIGLI SAW SET:** The premium quality construction of Gigli

Saw Handle offers a valuable and durable addition to any hospital or clinic. It is a popular orthopedic instrument used for a wide variety of different procedures comprised of flexible wires twisted in order to create multiple cutting surfaces.

**OR GRADE STAINLESS STEEL:** The Gigli Saw Handle is made of high quality O.R. Grade stainless steel that will not bend or break on you. It is fully autoclavable so that it can be used safely again and again without degradation.

**PRECISE AND DURABLE:** The handles are is corrosion resistant, features a high mirror finish, and is tarnish free. Ultrasonic cleaned for the highest possible quality delivered.

**HANDLE DESIGN:** Gigli saw handle is utilized in several amputation



WI-36056304

## 20PCS Orthopedic Surgery Set

**\$160.00** per piece

Trade price from **\$135.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

20PCS Orthopedic Surgery Set is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 20-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. The steel is heat treated and

passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$160.00, falling to \$135.00 at 1 set — about 16% off the single-unit price. Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36088265

## 20PCS Orthopedic Surgery Set Bone Surgery Instruments

**\$160.00** per piece

Trade price from **\$135.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT


20PCS Orthopedic Surgery Set Bone Surgery Instruments — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 20-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. Forging, milling, heat

treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$160.00, falling to \$135.00 at 1 set — about 16% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35998304

## 20PCS Orthopedic Surgery Set Surgical Instruments

**\$160.00** per piece

Trade price from **\$135.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

We manufacture 20PCS Orthopedic Surgery Set Surgical Instruments for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 20-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$160.00, falling to \$135.00 at 1 set — about 16% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45575604

**22 Pcs Mueller Cement Removal Instrument Set**

**\$550.00** per piece

Trade price from **\$480.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

22 Pcs Mueller Cement Removal Instrument Set — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 22-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. Forging, milling, heat

treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$550.00, falling to \$480.00 at 1 set — about 13% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45916812

**22PCS Shoulder Surgery Instrument Set Orthopedic Arthroscopy Surgical**

**\$600.00** per piece

Trade price from **\$350.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

22PCS Shoulder Surgery Instrument Set Orthopedic Arthroscopy Surgical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 22-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$600.00, falling to \$350.00 at 1 set — about 42% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35528201

25 PCS Basic Orthopedic Surgery Set

\$300.00 per piece

Trade price from \$210.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of 25 PCS Basic Orthopedic Surgery Set, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 25-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. The steel is heat treated and

passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$300.00, falling to \$210.00 at 1 set — about 30% off the single-unit price. Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35257162

25 Piece Orthopedic Surgery Instrument Set

\$190.00 per piece

Trade price from \$100.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of 25 Piece Orthopedic Surgery Instrument Set, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$190.00, falling to \$100.00 at 1 set — about 47% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

**WI-36061249**

**25PCS Orthopedic Surgery Set Stainless Steel Orthopedic**

**\$150.00** per piece

Trade price from **\$120.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


**ABOUT THIS INSTRUMENT**

This is our pattern of 25PCS Orthopedic Surgery Set Stainless Steel Orthopedic, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 25-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. Forging, milling, heat

treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$150.00, falling to \$120.00 at 1 set — about 20% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

**WI-36049244**

**25PCS Orthopedic Surgical Set Stainless Steel Orthopedic Instruments**

**\$150.00** per piece

Trade price from **\$120.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

25PCS Orthopedic Surgical Set Stainless Steel Orthopedic Instruments — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 25-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$150.00, falling to \$120.00 at 1 set — about 20% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

**WI-30873075****25-Piece Steel Major Basic Orthopedic Surgery Set Power****\$190.00** per pieceTrade price from **\$160.00**

Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

We manufacture 25-Piece Steel Major Basic Orthopedic Surgery Set Power for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$190.00, falling to \$160.00 at 1 set — about 16% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|                  |                                                                                 |
|------------------|---------------------------------------------------------------------------------|
| <b>Material</b>  | Stainless steel — we state the grade on your quotation                          |
| <b>Reusable</b>  | Yes, autoclavable                                                               |
| <b>Finish</b>    | Satin, mirror or black — your choice at order                                   |
| <b>Marking</b>   | Your logo or part number can be applied — ask us what is possible for this item |
| <b>Lead time</b> | Typically 15–30 days depending on quantity and marking                          |

**WI-70001250****3 Pcs Dingman Cartilage Bone Clamp 7½" Serrated 2×2 Teeth Jaws****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Premium Quality DINGMAN Small Bone Holding & Cartilage Clamp Forceps 7.5" (19cm) Serrated 2×2 Teeth, 4mm Wide Jaws German Grade Stainless Steel Orthopedic Surgical Instruments


Dingmann Cartilage Seizing Forceps is a specialized surgical tool that orthopedic surgeons frequently use to grasp, retract and hold back cartilages and other tissues, especially during bone fracture repairing and joint surgery.

Dingman bone clamp comes with curved sideways shanks for engaging in difficult approaches. This orthopedic instrument has atraumatic 2×2 toothed jaws. Moreover, the jaws have fine and deeper serrations along their inner surfaces. As a result, surgeons can

manipulate cartilages and other tissues with a strong grasp that prevents slippage. The finger ring ratchet lock helps in performing static clamping.

**DINGMAN BONE CLAMP FEATURES**

- 2×2 Toothed Tips Granting Secure Grasps.
- 4mm Wide Jaws
- 7.5" (19cm) Overall Length
- Atraumatic Delicate Serrated Jaws for Avoiding Local Damage.
- Sideways Curved Shanks
- Perfect for Holding Small Bone Fragments and Cartilage
- Ergonomic Finger Rings Ratchet Lock Ensuring Superior Handling
- Made With German Grade Stainless Steel



**WI-70001814****3 Pcs Gigli Saw Wire 12", 20" & 30"****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 GIGLI Saw Wire 12", 20" & 30" Surgical Veterinary Neurology Orthopedic Instruments

Gigli Wire Saw is a specialized surgical device that orthopedic and neurosurgeons can use for cutting bone surgeries during a wide range of procedures, especially during amputation. The Gigli Wire Saw feature small rings to attach the handles. Surgeons can insert T-shaped or loop-shaped handles before starting the bone cutting. The saw consists of a bundle of wires twisted in a cutting fashion. As a result, surgeons can slide the wire repeatedly for obtaining precise cuts.

**GIGLI SAW WIRE FEATURES**

**ORTHOPEDIC GIGLI SAW SET:** The premium quality construction of Gigli Saw Wire offers a valuable and durable addition to any hospital or clinic. It is a popular orthopedic instrument used for a wide variety of different procedures comprised of flexible wires twisted in order to create multiple cutting surfaces.

**OR GRADE STAINLESS STEEL:** The Gigli Saw wire is made of high quality O.R. Grade stainless steel that will not bend or break on you. It is fully autoclavable so that it can be used safely again and again without degradation.

**PRECISE AND DURABLE:** The Gigli Saw wire is corrosion resistant, features a high mirror finish, and is tarnish free. Ultrasonic cleaned for the highest possible quality delivered.

**WI-70000649****3 Plate Bone Holding Reduction Forceps Set 5.5" Gold Handle****\$79.99** per pieceTrade price from **\$73.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Plate Bone Holding Reduction Forceps Set 5.5" (14cm) with Ratcheted Finger Ring Gold Handles Premium Quality Stainless Steel Surgical Orthopedic Instruments.

Bone forceps are widely used for orthopedic surgical procedures. This plate bone holding reduction forceps set consists of three different types of forceps used for bone holding, bone reduction and bone plate holding. These forceps are ideal to grip and manipulate bones and bone fragments during orthopedic procedures. The ergonomic Jaws are engineered to enhance their grip on separate bone surfaces, which minimizes the risk of slippage. The finger ring handles reduce the fatigue during extended procedures. Made with premium grade

stainless steel ensuring the highest quality and durability standards in your operating rooms. Fully autoclaved and can be reuse.

**BONE HOLDING FORCEPS FEATURES**

- Crafted with 5.5" (14cm) Overall Length.
- Curved Tips Profile
- Variety of Tips for Strong Grasping.
- Ratchet Lock for Controlling the Clamping Degree.
- Ergonomic Finger Ring Gold Plated Handle Ensuring Superior Control.
- Ideal for Grasping Bone and Tissue Fragments.
- Fully Autoclaved and Reusable.
- Made with German Grade Stainless Steel.





|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-27305681

36-Piece Orthopedic Surgical Set Precision Instruments 2-Year

\$250.00 per piece

Trade price from \$230.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

36-Piece Orthopedic Surgical Set Precision Instruments 2-Year is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$250.00, falling to \$230.00 at 1 set — about 8% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.



WI-70001790

5 Kern Bone Holding Forceps Set 5.5", 6", 9", 13" & 17" with Ratchet

\$149.99 per piece

Trade price from \$137.99  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 KERN Bone Holding Forceps 5.5", 6", 9", 13" and 17" Long 2×2 Teeth Tips Serrated Jaws with Bar Ratchet Mechanism – Orthopedic Surgical Instruments

The Kern Bone Holding Forceps is a specialized orthopedic device that surgeons commonly use for grasping and mobilizing fractured bone pieces, in order to align them into the anatomical stabilization spot during surgical procedures.

This instrument comes with a pair of 2×2 toothed jaws. Moreover, the jaws have fine horizontal serration on their inner surfaces for fixating onto the bone surfaces with ease. As a result, surgeons can reposition the desired bone fragments into their original anatomical place. Moreover Kern Bone


Holding Forceps have a plier style handle with a bar ratchet system at the bottom which allows the surgeons to statically lock the instrument at the desired width for controlling the clamping pressure.

KERN BONE HOLDING FORCEPS FEATURES

- Secure Grasps with 2×2 Toothed Tips.
- Heavy Serrated Jaws.
- Single Spring Action Handle.
- Bar Ratchet Mechanism for Statically Locking the Jaws.
- Ergonomic Plier Style Handle Ensuring Maximum Control.
- Made with German Grade Stainless Steel

PACKAGE INCLUDE



**5  
PCS****WI-70001502****5 Pcs Bone Reduction Forceps Set 5.5" Curved Pointed Tip Gold Handle****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Bone Holding Reduction Forceps 5.5" Curved Serrated Jaws Pointed Tips with Ratcheted Finger Ring Gold Handles Stainless Steel Surgical Orthopedic Instrument.

The bone reduction forceps with curved tips allow the gripping of denser and deeper bones. Reduction forceps assist in grasping dense bones located in deeper regions. These Bone Reduction Forceps also helpful for the soft tissue attachment intended to reduce distal radius fractures.

The instrument aids in the alignment and setting of bone fractures while firmly holding them. The ergonomic design of the forceps allows precision and reduces the chances of fracture during the surgical procedure. These forceps are available in different

variations. The long ratchet mechanism gives a precise performance during bone manipulation.

**BONE REDUCTION FORCEPS FEATURES:**

- Minimal damage during the surgical procedure
- Curved Light Serrated Jaws
- Precise Pointed Tips
- Size and tip variations are available
- German Grade Stainless Steel Finish

**PACKAGE INCLUDE:**

- 5 Pcs Bone Reduction Forceps 5.5" (14cm) Long Ratchet Gold Handles

**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Bone Reduction Forceps**WI-70000149****5 Pcs Cobb Spinal Elevators Set 10mm – 32mm Wide Blade****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Cobb Spinal Elevators Set 11" (27.9cm) Slightly Curved 10mm, 13mm, 19mm, 26mm, 32mm Wide Oval Tip Blades with Non-Slip Hollow Round Handle – Premium Quality Stainless Steel Orthopedic Spinal Surgical Instruments.

Cobb elevators are specialized surgical instruments essential in orthopedic and spinal surgeries for elevating, dissecting, and retracting the periosteum, muscles, and soft tissues from the bone. These durable stainless steel tools feature wide, slightly curved blades, providing a clear view of the surgical site during procedures such as spinal fusions or fracture repairs. Cobb spinal elevators have a round hollow handle for a secure grip and prevent slipping. The blades are slightly curved,

with round cutting edge widths ranging from 10mm to 32mm for various surgical applications. Made from stainless steel, they are autoclavable and reusable.

**COBB SPINAL ELEVATORS FEATURES:**

- Size: 11 Inches (27.9 cm) Overall Length.
- Slightly Curved Multiple Size Blade.
- Oval Tip Round Cutting Edge for Opening Incisions.
- Knurled Grip Non-slip Hollow Round Handle for Enhanced Gripping.
- Long Straight Shank Allow Deep Hard-to-reach Surgical Access.
- Multiple Blade Sizes for Reaching Narrow Spaces.



**WI-70000179****5 Pcs Gigli Saw Wire Set 12", 16", 20", 25" & 30"****\$24.99** per pieceTrade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 GIGLI Saw Wires Set 12", 16", 20", 25" & 30" – Premium Quality Stainless Steel Surgical Veterinary Neurology Orthopedic Instruments.

Gigli Wire Saw is a specialized surgical device that orthopedic and neurosurgeons can use for cutting bone surgeries during a wide range of procedures, especially during amputation. The Gigli Wire Saw feature small rings to attach the handles. Surgeons can insert T-shaped or loop-shaped handles before starting the bone cutting. The saw consists of a bundle of wires twisted in a cutting fashion. As a result, surgeons can slide the wire repeatedly for obtaining precise cuts.

**GIGLI SAW WIRE FEATURES**

**ORTHOPEDIC GIGLI SAW SET:** The premium quality construction of Gigli Saw Wire offers a valuable and durable addition to any hospital or clinic. It is a popular orthopedic instrument used for a wide variety of different procedures comprised of flexible wires twisted in order to create multiple cutting surfaces.

**OR GRADE STAINLESS STEEL:** The Gigli Saw wire is made of high quality O.R. Grade stainless steel that will not bend or break on you. It is fully autoclavable so that it can be used safely again and again without degradation.

**PRECISE AND DURABLE:** The Gigli Saw wire is corrosion resistant, features a high mirror finish, and is tarnish free. Ultrasonic cleaned for the highest possible quality delivered.

**WI-70001137****5 Pcs Plaster & Cast Removal Set****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Plaster Cast Removal Instruments Set Premium Quality Cast Cutting Shears, Plaster Remover Cutter, Spreader and Breaker Orthopedic Surgical Instruments.

This plaster cast removal set specially designed for efficient removal of plaster casts during orthopedic surgeries. Made with high quality stainless steel ensures maximum durability and dependability. You can reuse this cast plaster instrument set over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

**PLASTER CAST REMOVAL SET INCLUDES:**

- 1 Henning Cast Spreader Plaster Breaker 11"
- 1 Bruns Esmarch Plaster Cast Cutting Shears 9.5"
- 1 Cast & Plaster Utility Shears 8"
- 1 Wolfe-Boehler Plaster Breaker 9"
- 1 Three Prong Cast Spreader 9"

**TECHNICAL INFORMATION**

**MATERIAL:** Stainless Steel

**INSTRUMENT TYPE:** Orthopedic Surgical Instrument

**APPLICATION:** Cut and Remove the Plaster Cast


**REUSABLE:** Yes

**AUTOCCLAVED:** Yes

**RUST PREVENTION PROCEDURE:** Passivated

**ULTRASONIC CLEANED:** Yes

**TESTS PERFORMED:** Boil Test,

**TOOLS SIZE**

**WI-70001249****6 Pcs Dingman Cartilage Bone Clamp 7½" Serrated 2×2 Teeth Jaws****\$79.99** per pieceTrade price from **\$73.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Premium Quality DINGMAN Small Bone Holding & Cartilage Clamp Forceps 7.5" (19cm) Serrated 2×2 Teeth, 4mm Wide Jaws German Grade Stainless Steel Orthopedic Surgical Instruments


Dingmann Cartilage Seizing Forceps is a specialized surgical tool that orthopedic surgeons frequently use to grasp, retract and hold back cartilages and other tissues, especially during bone fracture repairing and joint surgery.

Dingman bone clamp comes with curved sideways shanks for engaging in difficult approaches. This orthopedic instrument has atraumatic 2×2 toothed jaws. Moreover, the jaws have fine and deeper serrations along their inner surfaces. As a result, surgeons can

manipulate cartilages and other tissues with a strong grasp that prevents slippage. The finger ring ratchet lock helps in performing static clamping.

**DINGMAN BONE CLAMP FEATURES**

- 2×2 Toothed Tips Granting Secure Grasps.
- 4mm Wide Jaws
- 7.5" (19cm) Overall Length
- Atraumatic Delicate Serrated Jaws for Avoiding Local Damage.
- Sideways Curved Shanks
- Perfect for Holding Small Bone Fragments and Cartilage
- Ergonomic Finger Rings Ratchet Lock Ensuring Superior Handling
- Made With German Grade Stainless Steel

**WI-70001746****7 Pcs Hohmann Retractor Set 6, 8, 10, 18, 25, 43 & 70mm Blades****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 7 Hohmann Retractor Set (6mm, 8mm, 10mm, 18mm, 25mm, 43mm, 70mm) Blades Width Premium Quality Veterinary Surgical Orthopedic Instrument

Hohmann retractors, a crucial handheld orthopedic surgical instruments are designed to minimise tissue damage in arthroplasty surgery, and are widely used in minimally invasive hip and knee surgeries, where they retract the tissues around the margins of the joint. It Hohmann Bone Retractor features a curved, flat blade with a rounded tip useful for retracting tissues and organs gently.

These instruments can protect the soft tissue of the surgical site and increase exposure for the surgeons. The Hohmann Bone Liver Retractor

Veterinary Surgical features a slim body and a flat handle. The tip is leaf-shaped, round, square, etc. They are made of high-grade German stainless steel and multiple sizes and variations of this retractor available to needs of the patient.

**HOHMANN RETRACTOR FEATURES**

- 6mm, 8mm, 10mm, 18mm, 25mm, 43mm, 70mm Blades Width.
- Rectangular, Curved, Flattened Blade Tips.
- Ergonomic Flat Handle Design For Convenient Handling.
- Round Edges Minimizes Damage to the Soft Tissues
- Crafted with High-Grade German Stainless Steel





WI-70001122

7 Pcs Smith Peterson Chisel Set 8" Straight 0.64cm-3.2cm Blades

\$99.99 per piece

Trade price from \$91.99  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

7 Pcs Smith Peterson Chisel Set 8" (20.3cm) Straight 0.64cm – 3.2cm wide Blades Round Handle with FREE Carrying Pouch Stainless Steel Surgical Veterinary Orthopedic Instruments.

Smith Peterson Chisel designed to scrape and make fine cuts without damaging the surrounding structure. Smith Peterson Bone Chisel has a one-side beveled blade featured with a sharp cutting edge available in different sizes, ranging from 0.64cm to 3.2cm. As a result, surgeons can pick the adequate size for treating small or large bones, also contemplating the deep access to narrow spaces.

Moreover, Smith-Peterson Bone Chisel features a solid, bulky handle supplied with a smooth surface. Also, the handle includes a bottom circular impacting

platform. Surgeons can push the blade against the bone and hammer the platform with a mallet to cut cortical bone tissue. Specialty Used in Neurosurgery / Spine, Orthopedic, Post Mortem.

SMITH PETERSON CHISEL FEATURES

- Atraumatic Sharp One-side Beveled Cutting Blade for Smooth Cutting Action.
- From 0.64cm to 3.2cm Different Working End Sizes for Deep Surgical Access.
- Ergonomic Solid Handle for Superior Control.
- Ideal for Neurosurgery / Spine, Orthopedic, Post Mortem



WI-35367738

Advanced Orthopedic Biology Lab Dissecting Instrument Set Medical Student

\$9.00 per piece

Trade price from \$6.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Advanced Orthopedic Biology Lab Dissecting Instrument Set Medical Student is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$9.00, falling to \$6.00 at 1 set — about 33% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15-30 days depending on quantity and marking                          |



WI-46492870

Amputation Instrument Set-03

\$248.00 per piece

Trade price from \$188.00  
Minimum order 5 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Amputation Instrument Set-03 is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$248.00, falling to \$188.00 at 5 sets — about 24% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-27858712

Amputation Surgery Instrument Set 23 Pieces

\$150.00 per piece

Trade price from \$127.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Amputation Surgery Instrument Set 23 Pieces — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 23-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. The steel is heat treated and

passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$150.00, falling to \$127.00 at 1 set — about 15% off the single-unit price. Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-27331774

## Amputation Surgical Instrument Set Large Fragment Medical Tools

**\$270.00** per piece

Trade price from **\$255.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


### ABOUT THIS INSTRUMENT

We manufacture Amputation Surgical Instrument Set Large Fragment Medical Tools for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$270.00, falling to \$255.00 at 1 set — about 6% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45890366

## Basic for Lumbar Instrument Set Precision Spine Surgery Kit

**\$240.00** per piece

Trade price from **\$210.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


### ABOUT THIS INSTRUMENT

Basic for Lumbar Instrument Set Precision Spine Surgery Kit is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$240.00, falling to \$210.00 at 1 set — about 12% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-45870590

Basic for Lumbar Instruments Set Minimally Invasive Spine


**\$240.00** per piece Trade price from **\$210.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Basic for Lumbar Instruments Set Minimally Invasive Spine, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$240.00, falling to \$210.00 at 1 set — about 12% off the single-unit price. Larger quantities are priced on request.  
  
Order it on its own or as part of a mixed shipment; we do not insist on a full container.



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-45889397

Basic for Lumbar Surgical Instrument Set


**\$240.00** per piece Trade price from **\$210.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Basic for Lumbar Surgical Instrument Set, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$240.00, falling to \$210.00 at 1 set — about 12% off the single-unit price. Larger quantities are priced on request.  
  
Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



WI-27324175

### Basic Major Orthopedic Surgical Instruments Set 75 Pieces

**\$180.00** per piece

Trade price from **\$173.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Basic Major Orthopedic Surgical Instruments Set 75 Pieces is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 75-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$180.00, falling to \$173.00 at 1 set — about 4% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35467832

### Basic Major Orthopedic Surgical Instruments Set Stainless Steel Medical

**\$30.00** per piece

Trade price from **\$17.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Basic Major Orthopedic Surgical Instruments Set Stainless Steel Medical is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$30.00, falling to \$17.00 at 1 set — about 43% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-24339638

Basic Orthopedic Surgery Set 25 PCS

\$90.00 per piece

Trade price from \$70.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Basic Orthopedic Surgery Set 25 PCS is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 25-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. Forging, milling, heat

treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$90.00, falling to \$70.00 at 1 set — about 22% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-28705390

Basic Orthopedic Surgery Set 25 Pcs Complete Surgical Kit

\$200.00 per piece

Trade price from \$150.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Basic Orthopedic Surgery Set 25 Pcs Complete Surgical Kit is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 25-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. The steel is heat treated and

passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$200.00, falling to \$150.00 at 1 set — about 25% off the single-unit price. Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-46197258

Cardiac Vascular Thoracic Surgery Instrument Kit Precision Heart Operation

\$80.00 per piece

Trade price from \$50.00  
Minimum order 5 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

We manufacture Cardiac Vascular Thoracic Surgery Instrument Kit Precision Heart Operation for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$80.00, falling to \$50.00 at 5 sets — about 38% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45564529

Complete Basic for Lumbar Instrument Set Minimally Invasive Spine

\$240.00 per piece

Trade price from \$210.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Complete Basic for Lumbar Instrument Set Minimally Invasive Spine — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$240.00, falling to \$210.00 at 1 set — about 12% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45882337

### Complete Laminectomy Instrument Set 14PCS Precision Spine Surgery Kit

**\$300.00** per piece

Trade price from **\$250.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Complete Laminectomy Instrument Set 14PCS Precision Spine Surgery Kit — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 14-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$300.00, falling to \$250.00 at 1 set — about 17% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45886250

### Complete Mueller Type Cement Removal Set 22PCS Orthopedic Revision

**\$550.00** per piece

Trade price from **\$480.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

We manufacture Complete Mueller Type Cement Removal Set 22PCS Orthopedic Revision for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 22-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$550.00, falling to \$480.00 at 1 set — about 13% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35535185

Complete Orthopedic Surgical Set

\$400.00 per piece

Trade price from \$303.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Complete Orthopedic Surgical Set is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$400.00, falling to \$303.00 at 1 set — about 24% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45595260

Complete Periacetabular Osteotomy  
Orthopedic Instrument Set 10PCS Hip Surgery

\$300.00 per piece

Trade price from \$270.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Complete Periacetabular Osteotomy Orthopedic Instrument Set 10PCS Hip Surgery — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$300.00, falling to \$270.00 at 1 set — about 10% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45875269

Complete Small Fragment Orthopedic Instrument Set 30PCS Bone Trauma

\$499.00 per piece

Trade price from \$450.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Complete Small Fragment Orthopedic Instrument Set 30PCS Bone Trauma is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 30-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$499.00, falling to \$450.00 at 1 set — about 10% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45955039

Complete Stainless Steel Shoulder Surgery Set 22Pcs Orthopedic Arthroscopy

\$600.00 per piece

Trade price from \$350.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Complete Stainless Steel Shoulder Surgery Set 22Pcs Orthopedic Arthroscopy for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 22-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$600.00, falling to \$350.00 at 1 set — about 42% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45858783

## DHS Box Medium Orthopedic Trauma Instrument Set Dynamic

**\$150.00** per piece

Trade price from **\$130.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


### ABOUT THIS INSTRUMENT

DHS Box Medium Orthopedic Trauma Instrument Set Dynamic is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$150.00, falling to \$130.00 at 1 piece — about 13% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-28844453

## Excellent Broken Screw Removal Instrument Set Orthopedic

**\$424.00** per piece

Trade price from **\$400.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

### ABOUT THIS INSTRUMENT

This is our pattern of Excellent Broken Screw Removal Instrument Set Orthopedic, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$424.00, falling to \$400.00 at 1 set — about 6% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-26509143

Genecology Instruments

\$110.00 per piece

Trade price from \$100.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Genecology Instruments — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$110.00, falling to \$100.00 at 1 set — about 9% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-25752471

German Stainless Steel Surgical Orthopedic Instrument

\$50.00 per piece

Trade price from \$30.00  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

German Stainless Steel Surgical Orthopedic Instrument is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$50.00, falling to \$30.00 at 1 piece — about 40% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-70001819

Gigli Wire Saw Set 25" with 2 Handle

\$10.99 per piece

Trade price from \$10.11  
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

GIGLI Saw Set (2 Handles & 1 Wire)  
Premium Quality Surgical Neurology  
Orthopedic Instruments

Gigli Wire Saw is a specialized surgical device that orthopedic and neurosurgeons can use for cutting bone surgeries during a wide range of procedures, especially during amputation. The Gigli Wire Saw feature small rings to attach the handles. Surgeons can insert T-shaped or loop-shaped handles before starting the bone cutting. The saw consists of a bundle of wires twisted in a cutting fashion. As a result, surgeons can slide the wire repeatedly for obtaining precise cuts.

GIGLI SAW SET FEATURES

**ORTHOPEDIC GIGLI SAW SET:** The premium quality construction of Gigli Saw Set (2 Handles 1 Wire) offers a valuable and durable addition to any hospital or clinic. It is a popular orthopedic instrument used for a wide variety of different procedures comprised of flexible wires twisted in order to create multiple cutting surfaces.

**OR GRADE STAINLESS STEEL:** The Gigli Saw Set is made of high quality O.R. Grade stainless steel that will not bend or break on you. It is fully autoclavable so that it can be used safely again and again without degradation.

**PRECISE AND DURABLE:** The Gigli Saw set is corrosion resistant, features a high mirror finish, and is tarnish free. Ultrasonic cleaned for the highest



WI-34729565

Large Fragment Instrument Set Stainless Steel  
Orthopedic Surgical

\$200.00 per piece

Trade price from \$120.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Large Fragment Instrument Set Stainless Steel Orthopedic Surgical, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$200.00, falling to \$120.00 at 1 set — about 40% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-32595825

Large Stainless Steel Pad 14 Pre-Cut Suture Wounds Practice

\$10.00 per piece

Trade price from \$7.00  
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Large Stainless Steel Pad 14 Pre-Cut Suture Wounds Practice is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$10.00, falling to \$7.00 at 25 pieces — about 30% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-32659810

Large Stainless Steel Pad 14 Pre-Cut Wounds Suture Practice

\$10.00 per piece

Trade price from \$7.00  
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Large Stainless Steel Pad 14 Pre-Cut Wounds Suture Practice is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$10.00, falling to \$7.00 at 25 pieces — about 30% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-32883708

Lumbar Spine Surgery Set Stainless Steel Orthopedic Instruments

\$900.00 per piece

Trade price from \$650.00  
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Lumbar Spine Surgery Set Stainless Steel Orthopedic Instruments for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$900.00, falling to \$650.00 at 25 sets — about 28% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-46296869

Lumbar Spine Surgical Instrument Kit for Spinal Operation Decompression

\$850.00 per piece

Trade price from \$750.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Lumbar Spine Surgical Instrument Kit for Spinal Operation Decompression is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$850.00, falling to \$750.00 at 1 set — about 12% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-35485713

Major Orthopedic Surgery Instrument Set  
Hospital Surgical Tools

\$30.00 per piece

Trade price from \$16.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Major Orthopedic Surgery Instrument Set Hospital Surgical Tools for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$30.00, falling to \$16.00 at 1 set — about 47% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-45872364

Mueller Type Cement Removal Set of 22 PCS

\$550.00 per piece

Trade price from \$480.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Mueller Type Cement Removal Set of 22 PCS — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 22-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. We work the instrument by

hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$550.00, falling to \$480.00 at 1 set — about 13% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.



WI-45862402

### Orthopedic Cerclage Wire Passing Surgical Kit 3-Piece

**\$120.00** per piece

Trade price from **\$80.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

This is our pattern of Orthopedic Cerclage Wire Passing Surgical Kit 3-Piece, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$120.00, falling to \$80.00 at 1 set — about 33% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-28866024

### Orthopedic Hip Knee Spine Tibia Femur

**\$420.00** per piece

Trade price from **\$405.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

We manufacture Orthopedic Hip Knee Spine Tibia Femur for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$420.00, falling to \$405.00 at 1 set — about 4% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-34585765

### Orthopedic Instrument Set Steel 56PCS Broken Screw Removal

**\$415.00** per piece

Trade price from **\$275.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Orthopedic Instrument Set Steel 56PCS Broken Screw Removal is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 56-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. Forging, milling, heat

treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$415.00, falling to \$275.00 at 1 piece — about 34% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36932285

### Orthopedic Intramedullary Nailing Instrument Set Stainless Steel Surgical

**\$699.00** per piece

Trade price from **\$685.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


#### ABOUT THIS INSTRUMENT

This is our pattern of Orthopedic Intramedullary Nailing Instrument Set Stainless Steel Surgical, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$699.00, falling to \$685.00 at 1 piece — about 2% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-34735529

Orthopedic Small Fragment Set Stainless Steel Locking Plate

\$200.00 per piece

Trade price from \$120.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Orthopedic Small Fragment Set Stainless Steel Locking Plate, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$200.00, falling to \$120.00 at 1 set — about 40% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45880190

Orthopedic Stainless Steel Wire Passer Fiber Handle

\$120.00 per piece

Trade price from \$80.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Orthopedic Stainless Steel Wire Passer Fiber Handle for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$120.00, falling to \$80.00 at 1 set — about 33% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35540173

Orthopedic Surgery Instruments Kit Stainless Steel Medical Tools

\$400.00 per piece

Trade price from \$303.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Orthopedic Surgery Instruments Kit Stainless Steel Medical Tools is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$400.00, falling to \$303.00 at 1 set — about 24% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35395938

Orthopedic Surgical Instrument Set

\$30.00 per piece

Trade price from \$14.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Orthopedic Surgical Instrument Set, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$30.00, falling to \$14.00 at 1 set — about 53% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35525197

Orthopedic Surgical Instrument Set 25 PCS Medical

\$300.00 per piece Trade price from \$210.00 Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Orthopedic Surgical Instrument Set 25 PCS Medical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 25-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. The steel is heat treated and

passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$300.00, falling to \$210.00 at 1 set — about 30% off the single-unit price. Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36088232

Orthopedic Surgical Instrument Set Stainless Steel Bone Surgery

\$128.00 per piece Trade price from \$100.00 Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Orthopedic Surgical Instrument Set Stainless Steel Bone Surgery for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$128.00, falling to \$100.00 at 1 set — about 22% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-25837418

### Orthopedic Surgical Instruments Basic Stainless Steel Re-usable

**\$50.00** per piece

Trade price from **\$30.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Orthopedic Surgical Instruments Basic Stainless Steel Re-usable is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$50.00, falling to \$30.00 at 1 piece — about 40% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-27352365

### Orthopedic Surgical Instruments Solid Handle Stainless Steel Hudson

**\$42.00** per piece

Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

Orthopedic Surgical Instruments Solid Handle Stainless Steel Hudson is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$42.00. Minimum order 1 set, and the price moves with quantity.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36969261

Orthopedic Trauma Instrument Set German

**\$125.00** per piece

Trade price from **\$90.00**  
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Orthopedic Trauma Instrument Set German is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$125.00, falling to \$90.00 at 2 pieces — about 28% off the single-unit price. Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36913377

Orthopedic Trauma Surgical Instrument Set Bone Fixation Surgery

**\$125.00** per piece

Trade price from **\$100.00**  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Orthopedic Trauma Surgical Instrument Set Bone Fixation Surgery for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$125.00, falling to \$100.00 at 1 piece — about 20% off the single-unit price. Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-45601048

Periacetabular Osteotomy Orthopedic Instruments Set 10 Pcs Hip Preservation

\$300.00 per piece

Trade price from \$270.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Periacetabular Osteotomy Orthopedic Instruments Set 10 Pcs Hip Preservation for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$300.00, falling to \$270.00 at 1 set — about 10% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

WI-45886322

Periacetabular Osteotomy Orthopedic Instruments Set of 10 Pcs

\$300.00 per piece

Trade price from \$270.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Periacetabular Osteotomy Orthopedic Instruments Set of 10 Pcs — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 10-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$300.00, falling to \$270.00 at 1 set — about 10% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



WI-45589374

Periacetabular Osteotomy Surgical Instrument Set 10 Piece Orthopedic

\$300.00 per piece

Trade price from \$270.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Periacetabular Osteotomy Surgical Instrument Set 10 Piece Orthopedic is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$300.00, falling to \$270.00 at 1 set — about 10% off the single-unit price. Larger quantities are priced on request.  
  
Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36935339

Reusable Stainless Steel Orthopedic Trauma Surgery Set Fracture

\$125.00 per piece

Trade price from \$90.00  
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Reusable Stainless Steel Orthopedic Trauma Surgery Set Fracture is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Supplied as a set — the contents can be adjusted to your list. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$125.00, falling to \$90.00 at 2 pieces — about 28% off the single-unit price.  
  
Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35493270

Shoulder Joint Arthroscopy Instrument Set  
Orthopedic Surgical Kit

\$25.00 per piece

Trade price from \$12.00  
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Shoulder Joint Arthroscopy Instrument Set Orthopedic Surgical Kit, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$25.00, falling to \$12.00 at 10 sets — about 52% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45937620

Shoulder Surgery Set 22Pcs Orthopedic  
Shoulder Arthroscopy Surgical

\$600.00 per piece

Trade price from \$350.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Shoulder Surgery Set 22Pcs Orthopedic Shoulder Arthroscopy Surgical is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 22-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. The steel is heat treated and

passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$600.00, falling to \$350.00 at 1 set — about 42% off the single-unit price. Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-34308005

Small Fragment 3.5 4.0mm Stainless Steel Instrument Set

\$200.00 per piece

Trade price from \$120.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Small Fragment 3.5 4.0mm Stainless Steel Instrument Set — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The 4.0mm size shown is our standard; other sizes in the same pattern are

made to order. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$200.00, falling to \$120.00 at 1 set — about 40% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-42812777

Small Fragment Instrument Set Stainless Steel Orthopedic Surgical

\$500.00 per piece

Trade price from \$300.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Small Fragment Instrument Set Stainless Steel Orthopedic Surgical for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$500.00, falling to \$300.00 at 1 set — about 40% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36177725

Small Fragment Instrument Set Stainless Steel  
Orthopedic Trauma

**\$500.00** per piece

Trade price from **\$300.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Small Fragment Instrument Set Stainless Steel Orthopedic Trauma — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$500.00, falling to \$300.00 at 1 set — about 40% off the single-unit price. Larger quantities are priced on request.  
Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45881225

Small Fragment Instruments Set 30 Pcs  
Orthopedic Trauma

**\$499.00** per piece

Trade price from **\$450.00**  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Small Fragment Instruments Set 30 Pcs Orthopedic Trauma is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 30-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. The steel is heat treated and

passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.  
Listed at \$499.00, falling to \$450.00 at 1 set — about 10% off the single-unit price. Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36131764

Small Fragment Orthopedic Instrument Set  
Stainless Steel Bone Fixation

\$500.00 per piece

Trade price from \$300.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Small Fragment Orthopedic Instrument Set Stainless Steel Bone Fixation is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$500.00, falling to \$300.00 at 1 set — about 40% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-36168695

Small Fragment Surgical Instrument Set  
Orthopedic Bone Plating

\$500.00 per piece

Trade price from \$300.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Small Fragment Surgical Instrument Set Orthopedic Bone Plating is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$500.00, falling to \$300.00 at 1 set — about 40% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45963958

Spine Fixation Box Spinal Surgery Instrument

\$199.00 per piece

Trade price from \$150.00  
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Spine Fixation Box Spinal Surgery Instrument is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$199.00, falling to \$150.00 at 1 piece — about 25% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45996496

Spine Fixation Instrument Set Spinal Fusion Surgery

\$199.00 per piece

Trade price from \$150.00  
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Spine Fixation Instrument Set Spinal Fusion Surgery for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$199.00, falling to \$150.00 at 2 pieces — about 25% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |

**WI-46521683****Stainless Steel Amputation Surgery Instrument Set for Limb****\$249.90** per pieceTrade price from **\$189.90**

Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

We manufacture . Stainless Steel Amputation Surgery Instrument Set for Limb for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Supplied as a set — the contents can be adjusted to your list. Every piece is

hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$249.90, falling to \$189.90 at 2 sets — about 24% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

|                  |                                                                                 |
|------------------|---------------------------------------------------------------------------------|
| <b>Material</b>  | Stainless steel — we state the grade on your quotation                          |
| <b>Reusable</b>  | Yes, autoclavable                                                               |
| <b>Finish</b>    | Satin, mirror or black — your choice at order                                   |
| <b>Marking</b>   | Your logo or part number can be applied — ask us what is possible for this item |
| <b>Lead time</b> | Typically 15–30 days depending on quantity and marking                          |

**WI-46517704****Stainless Steel Amputation Surgical Instrument Kit Orthopedic Surgery Tools****\$250.00** per pieceTrade price from **\$190.00**

Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

**ABOUT THIS INSTRUMENT**

This is our pattern of Stainless Steel Amputation Surgical Instrument Kit Orthopedic Surgery Tools, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$250.00, falling to \$190.00 at 1 set — about 24% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|                  |                                                                                 |
|------------------|---------------------------------------------------------------------------------|
| <b>Material</b>  | Stainless steel — we state the grade on your quotation                          |
| <b>Reusable</b>  | Yes, autoclavable                                                               |
| <b>Finish</b>    | Satin, mirror or black — your choice at order                                   |
| <b>Marking</b>   | Your logo or part number can be applied — ask us what is possible for this item |
| <b>Lead time</b> | Typically 15–30 days depending on quantity and marking                          |



WI-09914253

### Stainless Steel Bone Mallet Hammer Surgical Orthopedic

**\$11.35** per piece

Trade price from **\$5.68**  
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Bone Mallet Hammer Surgical Orthopedic for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$11.35, falling to \$5.68 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-35554136

### Stainless Steel Orthopedic Shoulder Arthroscopy Instrument Set

**\$25.00** per piece

Trade price from **\$12.00**  
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

#### ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Orthopedic Shoulder Arthroscopy Instrument Set for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$25.00, falling to \$12.00 at 10 sets — about 52% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45575603

Stainless Titanium Cement Removal  
Instruments Set Orthopedic

\$550.00 per piece

Trade price from \$480.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Titanium Cement Removal Instruments Set Orthopedic is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$550.00, falling to \$480.00 at 1 set — about 13% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-46819581

Steel Bone Chisel Orthopedic Surgical Kit for  
Bone Cutting

\$5.90 per piece

Trade price from \$3.50  
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Steel Bone Chisel Orthopedic Surgical Kit for Bone Cutting — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$5.90, falling to \$3.50 at 50 pieces — about 41% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-32640778

Suture Practice Kit Stainless Steel Medical Instrument Small

\$10.00 per piece

Trade price from \$7.00  
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

We manufacture Suture Practice Kit Stainless Steel Medical Instrument Small for bone, joint and trauma surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$10.00, falling to \$7.00 at 25 pieces — about 30% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-46524106

Thoracotomy Surgical Instrument Set Chest Surgery Kit Thoracic Operation

\$399.00 per piece

Trade price from \$299.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Thoracotomy Surgical Instrument Set Chest Surgery Kit Thoracic Operation is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$399.00, falling to \$299.00 at 1 set — about 25% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-28792056

Trauma Instrument Set Stainless Steel  
Orthopedic Surgical

\$280.00 per piece

Trade price from \$250.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Trauma Instrument Set Stainless Steel Orthopedic Surgical is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$280.00, falling to \$250.00 at 1 set — about 11% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45592253

Wire Passer Fiber Handle Set 3 Pcs Orthopedic

\$120.00 per piece

Trade price from \$75.00  
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Wire Passer Fiber Handle Set 3 Pcs Orthopedic is part of our orthopedic instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 3-piece set.

than ours. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$120.00, falling to \$75.00 at 2 sets — about 38% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45850861

Wire Passer Fiber Handle Set 3 Precision Orthopedic

\$120.00 per piece

Trade price from \$80.00  
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Wire Passer Fiber Handle Set 3 Precision Orthopedic, made for bone, joint and trauma surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$120.00, falling to \$80.00 at 1 set — about 33% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



WI-45560591

Wood Handle Wire Passer Set 3-Piece Orthopedic Surgical Instrument

\$120.00 per piece

Trade price from \$75.00  
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Wood Handle Wire Passer Set 3-Piece Orthopedic Surgical Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$120.00, falling to \$75.00 at 2 sets — about 38% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

|           |                                                                                 |
|-----------|---------------------------------------------------------------------------------|
| Material  | Stainless steel — we state the grade on your quotation                          |
| Reusable  | Yes, autoclavable                                                               |
| Finish    | Satin, mirror or black — your choice at order                                   |
| Marking   | Your logo or part number can be applied — ask us what is possible for this item |
| Lead time | Typically 15–30 days depending on quantity and marking                          |



# WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

## Request a quotation

Send us product codes from this catalogue, quantities and your country.  
We reply the same working day on WhatsApp.

|                    |                                                                           |
|--------------------|---------------------------------------------------------------------------|
| <b>Online shop</b> | shop.waqarintl.com                                                        |
| <b>Email</b>       | sales@waqarintl.com                                                       |
| <b>WhatsApp /</b>  |                                                                           |
| <b>Phone</b>       | +92 307 6173627                                                           |
| <b>Works</b>       | Roras Road, Opposite Anwar Khawaja<br>Industries, Sialkot 51310, Pakistan |

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.