



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

Clamps Surgical Instruments

Reusable stainless-steel clamps forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



10
PCS



10
PCS



5
PCS



261

PRODUCTS

\$0.80

PRICES FROM

OEM

PRIVATE LABEL

Global

WORLDWIDE
EXPORT

Contents

261 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

Product	Code	Price	Page
10 Jones Towel Clamps 2.5"	WI-70000876	\$39.99	9
10 Jones Towel Clamps 2"	WI-70000883	\$39.99	9
10 Jones Towel Clamps 3.5"	WI-70000874	\$41.99	10
105mm Pediatric Vascular Clamps German Stainless Steel Surgical Instrument	WI-35040585	\$1.50	10
12.5cm Curved Mosquito Hemostatic Forceps Stainless Steel Locking Surgical	WI-34861976	\$1.50	11
12cm Curved Mosquito Hemostatic Forceps Stainless Steel Reusable Surgical	WI-34846965	\$1.50	11
2 Dieffenbach Serrefine Clamps 1.5" (3.8cm) Straight Curved	WI-70000848	\$32.99	12
2 Dieffenbach Serrefine Clamps 2.25" (5.7cm) Straight Curved	WI-70000847	\$34.99	12
21cm Foster Sponge Forceps Atraumatic Dressing Clamp Surgical Holding	WI-36791585	\$6.20	13
5 Inch Curved Mosquito Hemostatic Forceps Stainless Steel Locking	WI-34860963	\$1.50	13
5 Jones Towel Clamps 2.5"	WI-70000877	\$24.99	14
5 Jones Towel Clamps 2"	WI-70000884	\$24.99	14
5.5 Inch Curved Steel Multi-Use Locking Mosquito Hemostat Surgical	WI-30416011	\$1.50	15
9cm Curved Mosquito Hemostatic Forceps Stainless Steel Surgical Artery	WI-34901070	\$1.50	15
Acland Vessel Clamp Set 6-Piece Precision Kit Microsurgery	WI-24305786	\$8.00	16
Acrylic Nail Clipper 5.5" Stainless Steel False Nails Cutter	WI-70001247	\$7.99	16
Adson Hemostatic Surgical Forceps Instrument Set	WI-46788893	\$7.80	17
Adson Steel Hemostatic Forceps Surgical Instrument for Precise	WI-46828006	\$7.00	17
ADSON Steel Reusable Hemostatic Forceps for Controlled Tissue Handling	WI-46809738	\$8.00	18
Allis Tissue Forceps 19cm Reusable Surgical Clamp Medical	WI-43454129	\$3.20	18
Artery Forceps Mosquito Hemostat Stainless Steel Surgical Clamp Curved	WI-36659791	\$1.00	19
Backhaus Towel Clamp and General Surgical Instrument Set-Manual Medical	WI-29266362	\$4.00	19
Backhaus Towel Forceps Ball Socket 11" Instrument Set	WI-29264363	\$4.00	20
Burlisher Clamp Reusable Surgical Instrument	WI-47144772	\$5.80	20
Burlisher Clamp Surgical Instrument for Precise Tissue Handling Vessel	WI-47166528	\$6.00	21
Burlisher Clamp-02 Surgical Clamp Precision Instrument for Tissue	WI-47145714	\$5.90	21
Carmalt Angiotribe Stainless Steel Forceps Surgical Hemostat Kit	WI-46957115	\$2.99	22
Carmalt Angiotribe Titanium Surgical Forceps	WI-46965076	\$2.90	22
Circumcision Clamp 0.8cm Urology Instrument	WI-70002243	\$14.99	23
Circumcision Clamp 1.1cm Urology Instrument	WI-70002242	\$14.99	23
Circumcision Clamp 1.2cm Urology Instrument	WI-70002241	\$14.99	24
Circumcision Clamp 1.3cm Urology Instrument	WI-70002240	\$14.99	24
Circumcision Clamp 1.4cm Urology Instrument	WI-70002239	\$14.99	25
Circumcision Clamp 1.6cm Urology Instrument	WI-70002238	\$20.99	25
Circumcision Clamp 2.1cm Urology Instrument	WI-70002237	\$24.49	26
Circumcision Clamp 2.6cm Urology Instrument	WI-70002236	\$25.49	26
Circumcision Clamp 2.9cm Urology Instrument	WI-70002235	\$25.99	27
Circumcision Clamp 3.5cm Urology Instrument	WI-70002234	\$27.49	27
Cooley Derra Pediatric Vascular Clamp 6.5" For Anastomosis – Large	WI-70001034	\$21.99	28
Cooley Derra Pediatric Vascular Clamp 6.5" For Anastomosis – Medium	WI-70001035	\$19.99	28
Cooley Derra Pediatric Vascular Clamp 6.5" For Anastomosis – Small	WI-70001036	\$19.99	29

Contents (continued)

Product	Code	Price	Page
Crile Artery Forceps Straight Stainless Steel Reusable Surgical Hemostatic	WI-34901642	\$1.70	30
Crile Forceps Stainless Steel Hemostatic Artery Forceps Surgical Clamp	WI-36628867	\$0.80	30
Crile Hemostat Forceps Straight Curved Sponge Forceps	WI-27261492	\$2.55	31
Crile Hemostatic Forceps 5.5 Inch Straight Stainless Steel Surgical	WI-35084646	\$1.30	31
Crile Surgical Forceps Stainless Steel Hemostatic Clamp	WI-36622951	\$0.90	32
Debakey Aortic Aneurysm Clamp 10" Curved 95mm Atraumatic 1x2 Teeth Jaws	WI-70000753	\$24.99	32
Debakey Aortic Aneurysm Clamp 12" Curved 110mm Atraumatic 1x2 Teeth Jaws	WI-70000752	\$26.99	33
Debakey Atraumatic Aortic Aneurysm Clamps 10.75" Full Curved 111mm Jaws	WI-70000713	\$26.99	33
DeBakey Atraumatic Clamp Vascular Surgery Instrument Non-Crushing Vessel	WI-47174075	\$17.00	34
DeBakey Atraumatic Clamp Vascular Surgical Kit	WI-47167608	\$18.00	34
Debakey Atraumatic Multi-Purpose Clamps 6.75" Curved Downwards 78mm Jaws	WI-70000747	\$14.99	35
DeBakey Atraumatic Ring Handle Bulldog Clamp 5" Curved	WI-70001002	\$19.99	35
DeBakey Atraumatic Vascular Clamp Surgical Instrument for Controlled	WI-47142823	\$17.90	36
Debakey Atraumatic Vascular Clamps 9" Angled 80mm Jaws	WI-70000748	\$17.99	36
Debakey Atraumatic Vascular Tissue Forceps Clamp Stainless Steel Instrument	WI-34885189	\$17.00	37
Debakey Beck Atraumatic Vena-Cava Clamp 8" Double Angled 65mm Jaws	WI-70000716	\$24.99	37
DeBakey Bulldog Clamp Atraumatic Forceps 1.5 Straight	WI-34944365	\$7.00	38
Debakey Forceps Atraumatic Vascular Tissue Forceps Clamp The Basis	WI-34849412	\$17.00	38
Debakey Multi-Purpose Clamp 6.5" Angled 30° 3.5cm Long Atraumatic Jaws	WI-70001033	\$24.99	39
Debakey Peripheral Vascular Clamps 6.5" Curved 56mm Atraumatic Jaws	WI-70000751	\$14.99	39
Debakey Peripheral Vascular Clamps 7.75" Curved 90mm Atraumatic Jaws	WI-70000750	\$17.99	40
Debakey Peripheral Vascular Clamps 8" Curved 93mm Atraumatic Jaws	WI-70000749	\$18.99	40
DeBAKEY Ring Handled Bulldog Clamps 45° Angled 4 7/8" Curved Shanks	WI-70001003	\$19.99	41
DeBakey Tangential Occlusion Clamps 11" Double Angled 75mm Long x 17mm Deep Jaws	WI-70000746	\$28.99	41
DeBakey-Satinsky Vena Cava Clamp 8.25" Double Angled	WI-70001106	\$24.99	42
DeBakey-Satinsky Vena Cava Clamp 9.5" Double Angled	WI-70001105	\$26.99	42
Dieffenbach Serrefine Bulldog Clamp 2 1/4" (5.7cm) Straight	WI-70000851	\$18.99	43
Dieffenbach Serrefine Clamp 1 1/2" (3.8cm) Curved	WI-70000850	\$17.99	43
Dieffenbach Serrefine Clamp 1 1/2" (3.8cm) Straight	WI-70000852	\$17.99	44
Dieffenbach Serrefine Clamp 2 1/4" (5.7cm) Curved	WI-70000849	\$18.99	44
Duval Clamp Atraumatic Surgical Instrument for Gentle Tissue	WI-47078580	\$7.90	45
Duval Lung Clamp Atraumatic Surgical Tool	WI-47091453	\$7.50	45
Endodontic Rubber Dam Clamp #0 – Premolar	WI-70001597	\$8.99	46
Endodontic Rubber Dam Clamp #00 Winged	WI-70003674	\$8.99	46
Endodontic Rubber Dam Clamp #1 – Premolar	WI-70001603	\$8.99	47
Endodontic Rubber Dam Clamp #10 – General Molars	WI-70001592	\$8.99	47
Endodontic Rubber Dam Clamp #11 – General Molars	WI-70001591	\$8.99	48
Endodontic Rubber Dam Clamp #12 – Lower Molars	WI-70001590	\$8.99	48
Endodontic Rubber Dam Clamp #12A for Lower Right Molars	WI-70003673	\$8.99	49
Endodontic Rubber Dam Clamp #13 – Lower Left Molars	WI-70001589	\$8.99	49
Endodontic Rubber Dam Clamp #138 for Left-Side Third Molars	WI-70000314	\$8.99	50
Endodontic Rubber Dam Clamp #139 for Right-Side Third Molars	WI-70003672	\$8.99	50

Contents (continued)

Product	Code	Price	Page
Endodontic Rubber Dam Clamp #13A for Lower Left Molars	WI-70003671	\$8.99	51
Endodontic Rubber Dam Clamp #14 for General Molars	WI-70003669	\$8.99	51
Endodontic Rubber Dam Clamp #14A for General Molars	WI-70003670	\$8.99	52
Endodontic Rubber Dam Clamp #15 for Big Lower Molars	WI-70003668	\$8.99	52
Endodontic Rubber Dam Clamp #18 for Upper & Lower Molars	WI-70003667	\$8.99	53
Endodontic Rubber Dam Clamp #1A – Canines and Premolars	WI-70001595	\$8.99	53
Endodontic Rubber Dam Clamp #1T – Edged Premolars	WI-70001594	\$8.99	54
Endodontic Rubber Dam Clamp #2 – Premolar	WI-70001602	\$8.99	54
Endodontic Rubber Dam Clamp #200 for Lower Molars	WI-70003666	\$8.99	55
Endodontic Rubber Dam Clamp #201 for Upper Molars	WI-70003665	\$8.99	55
Endodontic Rubber Dam Clamp #202 for Large Lower Molars	WI-70000313	\$8.99	56
Endodontic Rubber Dam Clamp #203 for Small Lower Molars	WI-70003664	\$8.99	56
Endodontic Rubber Dam Clamp #205	WI-70003662	\$8.99	57
Endodontic Rubber Dam Clamp #206	WI-70003661	\$8.99	57
Endodontic Rubber Dam Clamp #208	WI-70003660	\$8.99	58
Endodontic Rubber Dam Clamp #209 for Premolars	WI-70003647	\$8.99	58
Endodontic Rubber Dam Clamp #210 for Upper Centrals and Incisors	WI-70000312	\$8.99	59
Endodontic Rubber Dam Clamp #211 for Adult Lower Anterior	WI-70000311	\$8.99	59
Endodontic Rubber Dam Clamp #212 – Anterior	WI-70001586	\$8.99	60
Endodontic Rubber Dam Clamp #212SA for Anteriors and Premolars	WI-70000310	\$8.99	60
Endodontic Rubber Dam Clamp #22 for Upper & Lower Bicuspid	WI-70000316	\$8.99	61
Endodontic Rubber Dam Clamp #26N – Upper Molars	WI-70001585	\$8.99	61
Endodontic Rubber Dam Clamp #27 for Large Molars	WI-70000315	\$8.99	62
Endodontic Rubber Dam Clamp #27N – Premolars	WI-70001584	\$8.99	62
Endodontic Rubber Dam Clamp #28 for Lower Molars	WI-70003655	\$8.99	63
Endodontic Rubber Dam Clamp #29 for Lower Premolars	WI-70003654	\$8.99	63
Endodontic Rubber Dam Clamp #2A – Larger Premolar	WI-70001601	\$8.99	64
Endodontic Rubber Dam Clamp #3 for Small Lower Molars	WI-70003653	\$8.99	64
Endodontic Rubber Dam Clamp #30 for Upper Right Molars	WI-70003652	\$8.99	65
Endodontic Rubber Dam Clamp #31 for Upper Left Molars	WI-70003651	\$8.99	65
Endodontic Rubber Dam Clamp #4	WI-70003650	\$8.99	66
Endodontic Rubber Dam Clamp #5	WI-70003649	\$8.99	66
Endodontic Rubber Dam Clamp #51 for Upper & Lower Molars	WI-70003648	\$8.99	67
Endodontic Rubber Dam Clamp #7 – Lower Molars	WI-70001599	\$8.99	67
Endodontic Rubber Dam Clamp #7A – Larger Lower Molars	WI-70001598	\$8.99	68
Endodontic Rubber Dam Clamp #8A	WI-70003646	\$8.99	68
Endodontic Rubber Dam Clamp #8AD – Molar	WI-70001587	\$8.99	69
Endodontic Rubber Dam Clamp #9 – Anterior	WI-70001593	\$8.99	69
Endodontic Rubber Dam Clamp #W0	WI-70001607	\$8.99	70
Endodontic Rubber Dam Clamp #W14 – Wingless	WI-70003645	\$8.99	70
Endodontic Rubber Dam Clamp #W3	WI-70001606	\$8.99	71
Endodontic Rubber Dam Clamp #W4	WI-70001605	\$8.99	71

Contents (continued)

Product	Code	Price	Page
Endodontic Rubber Dam Clamp #W7	WI-70001604	\$8.99	72
Endodontic Rubber Dam Clamp #W8	WI-70003644	\$8.99	72
Endodontic Rubber Dam Clamp Labial #6	WI-70001600	\$8.99	73
Endodontic Rubber Dam Clamps #W1	WI-70001608	\$8.99	73
Ferguson Angiotribe Forceps Vascular Surgical Kit Steel Components	WI-46417415	\$13.00	74
Fogarty Clamp Surgical Kit Vascular Surgery Instrument	WI-47073655	\$13.00	74
Fogarty Clamp Vascular Surgical Instrument Set for Precise Vessel	WI-47112217	\$13.90	75
Fogarty Vascular Clamp (Model Fogarty Clamp-01) Steel Surgical Instrument	WI-47117134	\$14.00	75
Foster Sponge Forceps 21cm Locking Dressing Clamp Precision Surgical	WI-36701790	\$6.20	76
German Debakey Bulldog Clamp Set Atraumatic Curved Vascular Forceps	WI-34986137	\$7.00	76
German Stainless Steel Satinsky DeBakey Clamp 34	WI-34858395	\$17.00	77
Glover Atraumatic Bulldog Clamps 5.5	WI-34907469	\$7.00	77
Glover Bulldog Clamp with Spring Adjusting Screw 9cm Curved 4cm Jaws	WI-70000391	\$43.99	78
Halsted Mosquito Artery Forceps Straight Curved Stainless Steel Hemostatic	WI-34952003	\$1.50	78
Halsted Mosquito Forceps 12.5cm Straight Curved Stainless Steel Hemostatic	WI-34851876	\$1.50	79
Harken Vena-Cava Clamp 9.5" Curved #1 Serrated Jaws	WI-70000854	\$24.99	79
Hartman Hemostat Forceps Curved Surgical Clamping Tool	WI-46377602	\$3.90	80
Hartman Hemostatic Forceps Surgical Hemostat Precise Artery Clamping Bleeding	WI-46421374	\$4.00	80
Heavy Duty Steel Straight Forceps Tissue Grasping Clamp Surgical	WI-45271703	\$8.00	81
Heavy Pattern Rochester Ochsner Curved Clamp Precision Surgical	WI-36752650	\$8.00	81
Hemostatic Forceps Surgical Instrument	WI-36938012	\$5.00	82
Hemostatic Mixer Artery Forceps Curved Stainless Steel Surgical Instruments	WI-28116498	\$8.00	82
Hemostatic Mixer Artery Straight Stainless Steel Surgical Forceps	WI-28135348	\$8.00	83
Ingrown Toenail Clippers & Cutter 4.5" (11.4cm) For Thick Nails	WI-70000936	\$9.99	83
Javid Carotid Clamp Vascular Surgery Instrument for Secure Vessel	WI-47080544	\$17.00	84
Javid Carotid Clamp Vascular Surgical Instrument for Precise Carotid	WI-47073640	\$18.00	84
Javid Carotid Surgical Clamp Precision Vascular Instrument for Controlled	WI-47091422	\$17.90	85
Jones Towel Clamp 2.5"	WI-70000887	\$8.99	85
Jones Towel Clamp 2"	WI-70002921	\$8.99	86
Jones Towel Clamp 3.5"	WI-70000886	\$8.99	86
Kelly Artery Forceps Surgical Hemostat Clamp	WI-43456093	\$1.00	87
Kelly Curved Clamp Stainless Steel Hemostatic Artery Forceps	WI-36657806	\$1.30	87
Kelly Curved Hemostatic Forceps Medical Stainless Steel Artery	WI-36692761	\$1.30	88
Kelly Forceps Clamp Stainless Steel 3-Year Shelf	WI-27089727	\$3.00	88
Kelly Forceps Stainless Steel Hemostatic Artery Forceps Surgical Instrument	WI-44570126	\$1.30	89
Kelly Hemostatic Forceps 5.50" Stainless Steel Straight A+	WI-27135104	\$3.00	89
Kelly Hemostatic Forceps 5.50" Straight Surgical Instruments	WI-27673962	\$1.00	90
Kelly Hemostatic Forceps Stainless Steel Surgical Artery	WI-44557721	\$1.30	90
Kelly Hemostatic Forceps Straight & Curved Stainless Steel Surgical	WI-43448128	\$1.00	91
Kocher Ochsner Forceps Surgical Clamps Straight Curved	WI-01168086	\$3.00	91
Kocher Tissue Forceps-02 Stainless Steel Hemostatic	WI-35079891	\$2.20	92
Lahey Artery Forceps Stainless Steel Hemostatic Forceps	WI-44579068	\$3.20	92

Contents (continued)

Product	Code	Price	Page
Lahey Artery Forceps Stainless Steel Hemostatic Surgical Clamp	WI-44578087	\$3.20	93
Lahey Artery Forceps Surgical Hemostatic Clamp	WI-35106242	\$3.20	93
Lorna Non-Perforating Towel Clamp 4.5" Inward Curved Jaws	WI-70000958	\$9.99	94
Lorna Non-Perforating Towel Clamp 4" Inward Curved Jaws	WI-70000959	\$8.99	94
Lorna Non-Perforating Towel Clamp 5.25" Inward Curved Jaws	WI-70000957	\$10.99	95
Lorna Towel Clamp 3.5" Inward Curved Jaws	WI-70002920	\$8.99	95
Metal Angiotribe Forceps Hemostatic Surgical Instrument for Vessel Clamping	WI-46929399	\$3.00	96
Mikulicz Hemostat Forceps Model 03 Reusable Stainless Steel Surgical	WI-46597023	\$4.00	96
Mikulicz Hemostatic Forceps Surgical Hemostat for Tissue Clamping Bleeding	WI-46578175	\$9.00	97
Mikulicz Surgical Forceps Hemostatic Clamp for General Surgery Tissue	WI-46587072	\$8.90	97
Mosquito Artery Forceps Stainless Steel Surgical	WI-25889597	\$0.90	98
Mosquito Forceps Powder Coated	WI-25865740	\$1.00	98
Mosquito Forceps-19 Hemostatic Mosquito Artery Clamp	WI-36626810	\$1.00	99
Moynihan Artery Forceps Stainless Steel Surgical Hemostatic Instrument	WI-35110212	\$3.20	99
Multi Color Rainbow Kelly Hemostat Forceps Instrument Set 5	WI-34927695	\$3.50	100
Ochsner Kocher Hemostatic Forceps Curved 20cm Stainless Steel Surgical	WI-34903836	\$1.50	100
Payr Pylorus Clamps 11.5" with Pin 4.75" (120mm) Serrated Jaws	WI-70000689	\$45.99	101
Pean Forceps Locking Hemostat Clamp	WI-45169351	\$1.00	101
Pean Forceps Straight Curved Hemostatic Artery	WI-46148546	\$0.80	102
Pean Forceps Surgical Hemostat Clamp Reusable German	WI-46137586	\$0.90	102
Pean Rochester Forceps 6.5" Straight Curved	WI-27670944	\$1.00	103
Pean Rochester Straight Curved Forceps Hemostat Surgical	WI-27276391	\$2.55	103
Professional Cantilever Nail Clipper 6" Black Coated Double Action	WI-70003230	\$21.99	104
Rainbow Cantilever Nail Clipper 6" Double Action	WI-70000857	\$19.99	104
Reusable Stainless Steel Mosquito Artery Forceps Straight Surgical Hemostatic	WI-34893156	\$1.50	105
Robert Artery Forceps Stainless Steel Hemostatic Forceps	WI-44576095	\$2.50	105
Robert Artery Forceps Stainless Steel Hemostatic Surgical Clamp	WI-44584009	\$2.50	106
Rochester Ochsner Curved Forceps Heavy Duty Tissue Grasping Clamp	WI-36705771	\$8.00	106
Rubber Dam Clamp #204 Pedo Primary Tooth & Adult Premolar	WI-70003663	\$8.99	107
Rubber Dam Clamp #212R for Right Anterior	WI-70003659	\$8.99	107
Rubber Dam Clamp #24 for Upper Right & Lower Left Molars	WI-70003658	\$8.99	108
Rubber Dam Clamp #25 for Upper Left & Lower Right Molars	WI-70003657	\$8.99	108
Rubber Dam Clamp #26 for Small Lower Molars	WI-70003656	\$8.99	109
Rubber Dam Clamp #W14A for Partially Erupted or Irregularly Shaped Molars	WI-70000309	\$8.99	109
Satinsky Vascular Clamp Surgical 9" Forceps DeBakey	WI-24277762	\$17.50	110
Seller Roeder Towel Clamp Ball Jaws 11"	WI-29256369	\$4.00	110
Single Use Stainless Steel Artery Mosquito Forceps Surgical	WI-06914644	\$2.50	111
Spencer Wells Artery Forceps 16cm Stainless Steel Surgical Hemostatic	WI-34984273	\$1.70	111
Spencer Wells Artery Forceps Straight 12.5cm Stainless Steel Hemostatic	WI-34904633	\$2.50	112
Spencer Wells Artery Hemostatic Forceps Stainless Steel Surgical Instrument	WI-34982225	\$2.50	112
Spencer Wells Forceps 5 Inch Straight Stainless Steel Hemostatic	WI-35047678	\$1.80	113
Spencer Wells Forceps 8 Inch Curved Stainless Steel Hemostatic	WI-34935540	\$2.50	113

Contents (continued)

Product	Code	Price	Page
Spencer Wells Hemostatic Forceps Stainless Steel Surgical	WI-44544624	\$1.80	114
Spencer Wells Stainless Steel Hemostatic Artery Forceps Surgical Instrument	WI-44553598	\$1.80	114
Sponge Forceps Stainless Steel Surgical Dressing Clamp	WI-35059671	\$2.90	115
Stainless Pean Hemostatic Forceps Surgical Artery Clamp	WI-46164525	\$1.00	115
Stainless Steel 4-Piece Hair Extension Pliers	WI-37714333	\$20.00	116
Stainless Steel Babcock Surgical Forceps Tissue Holding Clamp	WI-35098433	\$2.20	116
Stainless Steel Babcock Tissue Forceps Atraumatic Surgical Clamp Set	WI-34990043	\$4.00	117
Stainless Steel Backhaus Towel Clamp Forceps	WI-29202455	\$4.00	117
Stainless Steel Benguela Hemostatic Forceps Kocher Straight/Curved	WI-26092304	\$3.00	118
Stainless Steel Cheatle Forceps Surgical Clamp	WI-34926496	\$3.00	118
Stainless Steel Crile Hemostat Forceps Surgical Locking	WI-36678845	\$0.80	119
Stainless Steel Curved Mosquito Forceps Locking Hemostatic Surgical Instrument	WI-34917508	\$2.00	119
Stainless Steel Diethrich Bulldog Clamps Use Surgery	WI-34940385	\$7.00	120
Stainless Steel Hair Extension Pliers 4-Piece	WI-37726309	\$15.00	120
Stainless Steel Hair Extension Pliers Kit	WI-39821326	\$15.00	121
Stainless Steel Hair Extension Pliers Kit 2 Clamp	WI-39815363	\$15.00	121
Stainless Steel Hair Extension Pliers Kit with 2	WI-39816352	\$15.00	122
Stainless Steel Hartman Hemostatic Surgical Forceps Precision Artery Clamp	WI-46420370	\$3.50	122
Stainless Steel Hemostat Kocher Forceps 7"	WI-01169093	\$3.00	123
Stainless Steel Hemostatic Forceps Surgical Instrument	WI-35076339	\$3.50	123
Stainless Steel Kelly Forceps	WI-33358800	\$3.00	124
Stainless Steel Mosquito Forceps Model 1	WI-06206494	\$2.50	124
Stainless Steel Mosquito Forceps Powder Coated	WI-06247465	\$2.50	125
Stainless Steel Mosquito Hemostat Locking Forceps	WI-06958100	\$1.50	125
Stainless Steel Moynihan Artery Forceps Surgical Hemostatic Clamp Medical	WI-35027754	\$3.20	126
Stainless Steel Spencer Wells Artery Forceps Curved Surgical	WI-34971373	\$1.70	126
Stainless Steel Sponge Forceps Surgical Holding Clamp Reusable Medical	WI-35014734	\$2.90	127
Stainless Steel Surgical Bulldog Clamps & Clips-Precision	WI-34958310	\$7.00	127
Stainless Steel Surgical Instrument Set	WI-00889157	\$2.00	128
Stainless Steel Surgical Micro Straight Bulldog Clamps	WI-34970224	\$7.00	128
Stainless Steel Surgical Towel Clip Clamp	WI-36646844	\$0.90	129
Stainless Steel Surgical Towel Clip Clamp Forceps	WI-44539597	\$1.50	129
Stainless Steel Towel Clamp	WI-29303341	\$4.00	130
Stainless Steel Towel Clip Surgical Draping Clamp Reusable Medical	WI-36649828	\$1.00	130
Stainless Steel Weft Extensions Tape Sealing	WI-33067066	\$15.00	131
Steel Curved Hemostatic Surgical Instruments	WI-28123441	\$8.00	131
Stokes Rubber Dam Clamp 7" with Lock Handle	WI-70003312	\$10.99	132
Straight Curved Mosquito Hemostatic Forceps Stainless Steel Artery Pean	WI-34855926	\$1.50	132
Straight Spencer Wells Hemostatic Forceps Stainless Steel Surgical Instrument	WI-35043709	\$1.80	133
Surgical Hemostatic Forceps With Blue Clip	WI-36822994	\$17.00	133
Surgical Instruments Micro Straight Bulldog Clamps Micro	WI-34890616	\$7.00	134
Surgical Instruments Set for Micro	WI-29308399	\$15.00	134

Contents (continued)

Product	Code	Price	Page
Surgical Instruments Stainless Steel	WI-01174096	\$3.00	135
Surgical Pean Hemostatic Forceps Straight Curved Artery Clamp	WI-45181285	\$1.00	135
Surgical Rochester Ochsner Straight Clamp Precision Tissue Holding Forceps	WI-45269768	\$8.00	136
Surgical Stainless Steel Metal Artery Forceps	WI-06214445	\$2.50	136
TC Babcock Forceps Stainless Steel Surgical Tissue	WI-34945610	\$4.00	137
Towel Clamp Stainless Steel Surgical Draping Clip	WI-36660872	\$1.00	137
Vascular Surgery Set-07 Vascular Surgical Kit Clamp Precision Artery	WI-46234797	\$220.00	138
Wertheim Atraumatic Parametrium Clamp 10" Straight 2×3 Teeth Jaws	WI-70000715	\$14.99	138
Wertheim Atraumatic Parametrium Clamp 9.75" Curved 2×3 Teeth Jaws	WI-70000714	\$14.99	139

**10
PCS****WI-70000876****10 Jones Towel Clamps 2.5"****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Jones Towel Clamps 2.5" (6.4cm) with Sharp Pointed Inwardly Curved Jaws and Cross-Action Spring Handle Premium Quality Stainless Steel Surgical Veterinary Instruments.

Jones Towel Clamps have sharp pointed tips which are curved and inward provides a secure and firm grip that doesn't let the object slip away. The tips overlap each other and are used to secure drapes and towels during surgical procedures. These instruments are very lightweight and are widely used in veterinary surgical procedures. The jones towel clamp forceps features a cross-action spring handle to open the jaws with comfort. Surgeons can loosen the pressure to slowly close the jaws and fixate the drapes without damaging the skin.

Manufactured in premium grade stainless steel to ensure the highest quality instruments in your operating rooms.

**JONES TOWEL CLAMPS
FEATURES**

- Crafted with 2.5" (6.4cm) Overall Length.
- Inward Curved Jaws.
- 1x1 Prong.
- Overlapping Sharp Pointed Tips.
- Cross Action Spring Handle.
- Lightweight and Comfortable to Hold.
- Widely Used in All Sorts of Surgical Procedures.
- Fully Autoclavable and Reusable.
- Made With Premium Grade Stainless Steel.

**10
PCS****WI-70000883****10 Jones Towel Clamps 2"****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Jones Towel Clamp 2" (5.1cm) with Sharp Pointed Inwardly Curved Jaws and Cross-Action Spring Handle Premium Quality Stainless Steel Surgical Veterinary Instruments.

Jones Towel Clamps have sharp pointed tips which are curved and inward provides a secure and firm grip that doesn't let the object slip away. The tips overlap each other and are used to secure drapes and towels during surgical procedures. These instruments are very lightweight and are widely used in veterinary surgical procedures. The jones towel clamp forceps features a cross-action spring handle to open the jaws with comfort. Surgeons can loosen the pressure to slowly close the jaws and fixate the drapes without damaging the skin.

Manufactured in premium grade stainless steel to ensure the highest quality instruments in your operating rooms.

**JONES TOWEL CLAMPS
FEATURES**

- Crafted with 2" (5.1cm) Overall Length.
- Inward Curved Jaws.
- 1x1 Prong.
- Overlapping Sharp Pointed Tips.
- Cross Action Spring Handle.
- Lightweight and Comfortable to Hold.
- Widely Used in All Sorts of Surgical Procedures.
- Fully Autoclavable and Reusable.
- Made With Premium Grade Stainless Steel.





WI-70000874

10 Jones Towel Clamps 3.5"

\$41.99 per piece

Trade price from \$38.63
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Jones Towel Clamps 3.5" (8.9cm) with Sharp Pointed Inwardly Curved Jaws and Cross-Action Spring Handle Premium Quality Stainless Steel Surgical Veterinary Instruments.

Jones Towel Clamps have sharp pointed tips which are curved and inward provides a secure and firm grip that doesn't let the object slip away. The tips overlap each other and are used to secure drapes and towels during surgical procedures. These instruments are very lightweight and are widely used in veterinary surgical procedures. The jones towel clamp forceps features a cross-action spring handle to open the jaws with comfort. Surgeons can loosen the pressure to slowly close the jaws and fixate the drapes without damaging the skin.

Manufactured in premium grade stainless steel to ensure the highest quality instruments in your operating rooms.

JONES TOWEL CLAMPS FEATURES

- Crafted with 3.5" (8.9cm) Overall Length.
- Inward Curved Jaws.
- 1x1 Prong.
- Overlapping Sharp Pointed Tips.
- Cross Action Spring Handle.
- Lightweight and Comfortable to Hold.
- Widely Used in All Sorts of Surgical Procedures.
- Fully Autoclavable and Reusable.
- Made With Premium Grade Stainless Steel.



WI-35040585

105mm Pediatric Vascular Clamps German Stainless Steel Surgical Instrument

\$1.50 per piece

Trade price from \$0.70
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

105mm Pediatric Vascular Clamps German Stainless Steel Surgical Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The 105mm size shown is our standard; other sizes in the same pattern are made to order.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.50, falling to \$0.70 at 10 pieces — about 53% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34861976

12.5cm Curved Mosquito Hemostatic Forceps Stainless Steel Locking Surgical

\$1.50 per piece

Trade price from **\$1.20**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

12.5cm Curved Mosquito Hemostatic Forceps Stainless Steel Locking Surgical is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

The 12.5cm size shown is our standard; other sizes in the same pattern are

made to order. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.50, falling to \$1.20 at 10 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34846965

12cm Curved Mosquito Hemostatic Forceps Stainless Steel Reusable Surgical

\$1.50 per piece

Trade price from **\$1.20**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

12cm Curved Mosquito Hemostatic Forceps Stainless Steel Reusable Surgical is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and full reusability — clean, autoclave and use again.

The 12cm size shown is our standard; other sizes in the same pattern are

made to order. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.50, falling to \$1.20 at 10 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

**WI-70000848****2 Dieffenbach Serrefine Clamps 1.5" (3.8cm)
Straight Curved****\$32.99** per pieceTrade price from **\$30.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Dieffenbach Serrefine Clamps Set 1.5" (3.8cm) Straight & Curved Serrated Round Tip Jaws with Spring Action Self Locking Mechanism Stainless Steel Coronary Artery Cardiovascular Surgery Instruments.

Dieffenbach Serrefine Clamp also known as "bulldog clamps" are a versatile surgical instrument used for precise and reliable occlusion (blocking) of blood vessels or tissue to control bleeding during surgical procedures. They used in oculoplastics to secure suture and prevent entangling with other suture as is done for the rectus muscle sutures in an enucleation. With spring-loaded mechanism and a series of sharp serrated teeth on the clamp's jaws, these dieffenbach serrefine forcep

clamp grip the blood vessel or tissue firmly, while the spring-loaded mechanism keeps the clamp closed.

Bulldog Dieffenbach Serrefine clamp commonly used in vascular surgeries, including open-heart surgery, where they are used to clamp off the coronary artery during a bypass procedure. They are also used in neurosurgery and orthopedic surgery. These clamps comes with different variations to help the professionals in multiple surgical procedures. Dieffenbach Serrefine available in straight and curved profiles while straight pattern provides easy manipulation. Made with premium quality stainless steel to ensure reliable and excellent performance

**DIEFFENBACH SERREFINE CLAMP
FEATURES****WI-70000847****2 Dieffenbach Serrefine Clamps 2.25" (5.7cm)
Straight Curved****\$34.99** per pieceTrade price from **\$32.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 2 Dieffenbach Serrefine Clamps Set 2.25" (5.7cm) Straight & Curved Serrated Round Tip Jaws with Spring Action Self Locking Mechanism Stainless Steel Coronary Artery Cardiovascular Surgery Instruments.

Dieffenbach Serrefine Clamp also known as "bulldog clamps" are a versatile surgical instrument used for precise and reliable occlusion (blocking) of blood vessels or tissue to control bleeding during surgical procedures. They used in oculoplastics to secure suture and prevent entangling with other suture as is done for the rectus muscle sutures in an enucleation. With spring-loaded mechanism and a series of sharp serrated teeth on the clamp's jaws, these dieffenbach serrefine forcep

clamp grip the blood vessel or tissue firmly, while the spring-loaded mechanism keeps the clamp closed.

Bulldog Dieffenbach Serrefine clamp commonly used in vascular surgeries, including open-heart surgery, where they are used to clamp off the coronary artery during a bypass procedure. They are also used in neurosurgery and orthopedic surgery. These clamps comes with different variations to help the professionals in multiple surgical procedures. Dieffenbach Serrefine available in straight and curved profiles while straight pattern provides easy manipulation. Made with premium quality stainless steel to ensure reliable and excellent performance

**DIEFFENBACH SERREFINE CLAMP
FEATURES**



WI-36791585

21cm Foster Sponge Forceps Atraumatic Dressing Clamp Surgical Holding

\$6.20 per piece

Trade price from \$4.50
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

21cm Foster Sponge Forceps Atraumatic Dressing Clamp Surgical Holding is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The 21cm size shown is our standard; other sizes in the same pattern are made to order.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$6.20, falling to \$4.50 at 25 pieces — about 27% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34860963

5 Inch Curved Mosquito Hemostatic Forceps Stainless Steel Locking

\$1.50 per piece

Trade price from \$1.20
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of 5 Inch Curved Mosquito Hemostatic Forceps Stainless Steel Locking, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

The 5" size shown is our standard; other sizes in the same pattern are

made to order. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$1.50, falling to \$1.20 at 10 pieces — about 20% off the single-unit price. Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**5
PCS****WI-70000877****5 Jones Towel Clamps 2.5"****\$24.99** per pieceTrade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Jones Towel Clamps 2.5" (6.4cm) with Sharp Pointed Inwardly Curved Jaws and Cross-Action Spring Handle Premium Quality Stainless Steel Surgical Veterinary Instruments.

Jones Towel Clamps have sharp pointed tips which are curved and inward provides a secure and firm grip that doesn't let the object slip away. The tips overlap each other and are used to secure drapes and towels during surgical procedures. These instruments are very lightweight and are widely used in veterinary surgical procedures. The jones towel clamp forceps features a cross-action spring handle to open the jaws with comfort. Surgeons can loosen the pressure to slowly close the jaws and fixate the drapes without damaging the skin.

Manufactured in premium grade stainless steel to ensure the highest quality instruments in your operating rooms.

**JONES TOWEL CLAMPS
FEATURES**

- Crafted with 2.5" (6.4cm) Overall Length.
- Inward Curved Jaws.
- 1x1 Prong.
- Overlapping Sharp Pointed Tips.
- Cross Action Spring Handle.
- Lightweight and Comfortable to Hold.
- Widely Used in All Sorts of Surgical Procedures.
- Fully Autoclavable and Reusable.
- Made With Premium Grade Stainless Steel.

**5
PCS****WI-70000884****5 Jones Towel Clamps 2"****\$24.99** per pieceTrade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Jones Towel Clamp 2" (5.1cm) with Sharp Pointed Inwardly Curved Jaws and Cross-Action Spring Handle Premium Quality Stainless Steel Surgical Veterinary Instruments.

Jones Towel Clamps have sharp pointed tips which are curved and inward provides a secure and firm grip that doesn't let the object slip away. The tips overlap each other and are used to secure drapes and towels during surgical procedures. These instruments are very lightweight and are widely used in veterinary surgical procedures. The jones towel clamp forceps features a cross-action spring handle to open the jaws with comfort. Surgeons can loosen the pressure to slowly close the jaws and fixate the drapes without damaging the skin.

Manufactured in premium grade stainless steel to ensure the highest quality instruments in your operating rooms.

**JONES TOWEL CLAMPS
FEATURES**

- Crafted with 2" (5.1cm) Overall Length.
- Inward Curved Jaws.
- 1x1 Prong.
- Overlapping Sharp Pointed Tips.
- Cross Action Spring Handle.
- Lightweight and Comfortable to Hold.
- Widely Used in All Sorts of Surgical Procedures.
- Fully Autoclavable and Reusable.
- Made With Premium Grade Stainless Steel.





Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-30416011

5.5 Inch Curved Steel Multi-Use Locking Mosquito Hemostat Surgical

\$1.50 per piece

Trade price from \$0.85
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

5.5 Inch Curved Steel Multi-Use Locking Mosquito Hemostat Surgical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

The 5.5" size shown is our standard; other sizes in the same pattern are

made to order. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.50, falling to \$0.85 at 1 piece — about 43% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34901070

9cm Curved Mosquito Hemostatic Forceps Stainless Steel Surgical Artery

\$1.50 per piece

Trade price from \$1.20
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of 9cm Curved Mosquito Hemostatic Forceps Stainless Steel Surgical Artery, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

The 9cm size shown is our standard; other sizes in the same pattern are made to order. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$1.50, falling to \$1.20 at 10 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



WI-24305786

Acland Vessel Clamp Set 6-Piece Precision Kit Microsurgery

\$8.00 per piece

Trade price from \$5.00
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Acland Vessel Clamp Set 6-Piece Precision Kit Microsurgery for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$8.00, falling to \$5.00 at 10 sets — about 38% off the single-unit price. Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001247

Acrylic Nail Clipper 5.5" Stainless Steel False Nails Cutter

\$7.99 per piece

Trade price from \$7.35
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Professional Acrylic Nail Clipper 5.5" (14cm) Stainless Steel False Nails Cutter Fake Nail Clippers Nail Tip Trimmer for Artificial Nail Art, Manicure Clip Tool Perfect for Nails Salon, Manicure Shop or Personal Home Diy Nail Art

tip to the nail at a 45-degree angle. Easy to use and efficient tool, portable and easy to clean. This nai cutter has comfortable, anti-slip handle and are safe and accurate to cut nails

ACRYLIC NAIL CLIPPER FEATURES

- Sharp & Long-lasting Blades
- Double-bladed Guillotine Style Acrylic False Nails Trimmer
- Spring Action with Anti-slip Handle
- Suitable for Acrylic Nails, Gel Nails, Fake Nails, Artificial Nails, False & Natural Nails
- Safe, Accurate and Easy to Use
- Made With High Quality Stainless Steel

PACKAGE INCLUDE

- 1 Pc Acrylic Nail Clipper 5.5" (14cm)





WI-46788893

Adson Hemostatic Surgical Forceps Instrument Set

\$7.80 per piece

Trade price from **\$3.80**
Minimum order 12 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Adson Hemostatic Surgical Forceps Instrument Set is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$7.80, falling to \$3.80 at 12 pieces — about 51% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46828006

Adson Steel Hemostatic Forceps Surgical Instrument for Precise

\$7.00 per piece

Trade price from **\$3.50**
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Adson Steel Hemostatic Forceps Surgical Instrument for Precise, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$7.00, falling to \$3.50 at 50 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46809738

ADSON Steel Reusable Hemostatic Forceps for Controlled Tissue Handling

\$8.00 per piece

Trade price from **\$4.00**
Minimum order 20 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of ADSON Steel Reusable Hemostatic Forceps for Controlled Tissue Handling, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$4.00 at 20 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-43454129

Allis Tissue Forceps 19cm Reusable Surgical Clamp Medical

\$3.20 per piece

Trade price from **\$2.90**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Allis Tissue Forceps 19cm Reusable Surgical Clamp Medical for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: full reusability — clean, autoclave and use again.

The 19cm size shown is our standard; other sizes in the same pattern are made to order. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.20, falling to \$2.90 at 1 piece — about 9% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36659791

Artery Forceps Mosquito Hemostat Stainless Steel Surgical Clamp Curved

\$1.00 per piece

Trade price from \$0.60
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Artery Forceps Mosquito Hemostat Stainless Steel Surgical Clamp Curved for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$1.00, falling to \$0.60 at 100 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29266362

Backhaus Towel Clamp and General Surgical Instrument Set-Manual Medical

\$4.00 per piece

Trade price from \$1.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Backhaus Towel Clamp and General Surgical Instrument Set-Manual Medical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$4.00, falling to \$1.50 at 1 piece — about 62% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29264363

Backhaus Towel Forceps Ball Socket 11\"/>

\$4.00 per piece

Trade price from \$1.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Backhaus Towel Forceps Ball Socket 11\"/>

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.00, falling to \$1.50 at 1 piece — about 62% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47144772

Burlisher Clamp Reusable Surgical Instrument

\$5.80 per piece

Trade price from \$1.70
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Burlisher Clamp Reusable Surgical Instrument, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: full reusability — clean, autoclave and use again.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$5.80, falling to \$1.70 at 50 pieces — about 71% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47166528

Burlisher Clamp Surgical Instrument for Precise Tissue Handling Vessel

\$6.00 per piece

Trade price from **\$2.00**
Minimum order 12 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Burlisher Clamp Surgical Instrument for Precise Tissue Handling Vessel — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$6.00, falling to \$2.00 at 12 pieces — about 67% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47145714

Burlisher Clamp-02 Surgical Clamp Precision Instrument for Tissue

\$5.90 per piece

Trade price from **\$1.90**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Burlisher Clamp-02 Surgical Clamp Precision Instrument for Tissue is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$5.90, falling to \$1.90 at 25 pieces — about 68% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46957115

Carmalt Angiotribe Stainless Steel Forceps
Surgical Hemostat Kit

\$2.99 per piece

Trade price from \$1.49
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Carmalt Angiotribe Stainless Steel Forceps Surgical Hemostat Kit is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.99, falling to \$1.49 at 50 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46965076

Carmalt Angiotribe Titanium Surgical Forceps

\$2.90 per piece

Trade price from \$1.50
Minimum order 30 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Carmalt Angiotribe Titanium Surgical Forceps — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.90, falling to \$1.50 at 30 pieces — about 48% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Forging, milling, heat treatment, hand finishing and passivation all happen

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70002243****Circumcision Clamp 0.8cm Urology Instrument****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Newborn Nursery Circumcision Clamp 0.8cm (8mm) Bell Shaped Premium Quality Chrome Plated Urology Surgical Instruments.

Circumcision is the oldest known surgical procedure but it is the most frequently performed surgery around the world in males. Circumcision Clamps is a specialized device that urology surgeons use to remove the glans "foreskin" from neonatal patients, as a preventive method to minimize the risk of urinary tract infection, HIV acquisition and other benefits.

CIRCUMCISION CLAMPS FEATURES

CIRCUMCISION CLAMP: With dimensions of: 4"(10.2cm) x 1"(2.54cm) x 1.5"(3.81cm) this clamp has a

circumcision bell shaped bevel hole base plate. Rocker arm with yoke to help the bell holding with the nut screw. We are using high-quality chrome plated metal.

CHROME PLATED CONSTRUCTION:

Made with high quality chrome-plated brass for a smooth finish. Brass is stronger and stiffer than pure copper, it is more susceptible to developing stress cracks. Compared to bronze, brass is more malleable to your operating rooms. The instrument can be sterilized and reused.

WHERE WE USE:: Its principal use is to provide adequate haemostasis throughout the neonatal circumcision procedure. Circumcision performed using the this clamp is usually quick and effective, and results in very little bleeding.

WI-70002242**Circumcision Clamp 1.1cm Urology Instrument****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Newborn Nursery Circumcision Clamp 1.1cm (11mm) Bell Shaped Premium Quality Chrome Plated Urology Surgical Instruments.

Circumcision is the oldest known surgical procedure but it is the most frequently performed surgery around the world in males. Circumcision Clamps is a specialized device that urology surgeons use to remove the glans "foreskin" from neonatal patients, as a preventive method to minimize the risk of urinary tract infection, HIV acquisition and other benefits.

CIRCUMCISION CLAMPS FEATURES

CIRCUMCISION CLAMP: With dimensions of: 4"(10.2cm) x 1"(2.54cm) x 1.5"(3.81cm) this clamp has a

circumcision bell shaped bevel hole base plate. Rocker arm with yoke to help the bell holding with the nut screw. We are using high-quality chrome plated metal.

CHROME PLATED CONSTRUCTION:

Made with high quality chrome-plated brass for a smooth finish. Brass is stronger and stiffer than pure copper, it is more susceptible to developing stress cracks. Compared to bronze, brass is more malleable to your operating rooms. The instrument can be sterilized and reused.

WHERE WE USE:: Its principal use is to provide adequate haemostasis throughout the neonatal circumcision procedure. Circumcision performed using the this clamp is usually quick and effective, and results in very little bleeding.

**WI-70002241****Circumcision Clamp 1.2cm Urology Instrument****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Newborn Nursery Circumcision Clamp 1.2cm (12mm) Bell Shaped Premium Quality Chrome Plated Urology Surgical Instruments.

Circumcision is the oldest known surgical procedure but it is the most frequently performed surgery around the world in males. Circumcision Clamps is a specialized device that urology surgeons use to remove the glans "foreskin" from neonatal patients, as a preventive method to minimize the risk of urinary tract infection, HIV acquisition and other benefits.

CIRCUMCISION CLAMPS FEATURES


CIRCUMCISION CLAMP: With dimensions of: 4"(10.2cm) x 1"(2.54cm) x 1.5"(3.81cm) this clamp has a

circumcision bell shaped bevel hole base plate. Rocker arm with yoke to help the bell holding with the nut screw. We are using high-quality chrome plated metal.

CHROME PLATED CONSTRUCTION:

Made with high quality chrome-plated brass for a smooth finish. Brass is stronger and stiffer than pure copper, it is more susceptible to developing stress cracks. Compared to bronze, brass is more malleable to your operating rooms. The instrument can be sterilized and reused.

WHERE WE USE:: Its principal use is to provide adequate haemostasis throughout the neonatal circumcision procedure. Circumcision performed using the this clamp is usually quick and effective, and results in very little bleeding.

**WI-70002240****Circumcision Clamp 1.3cm Urology Instrument****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Newborn Nursery Circumcision Clamp 1.3cm (13mm) Bell Shaped Premium Quality Chrome Plated Urology Surgical Instruments.

Circumcision is the oldest known surgical procedure but it is the most frequently performed surgery around the world in males. Circumcision Clamps is a specialized device that urology surgeons use to remove the glans "foreskin" from neonatal patients, as a preventive method to minimize the risk of urinary tract infection, HIV acquisition and other benefits.

CIRCUMCISION CLAMPS FEATURES


CIRCUMCISION CLAMP: With dimensions of: 4"(10.2cm) x 1"(2.54cm) x 1.5"(3.81cm) this clamp has a

circumcision bell shaped bevel hole base plate. Rocker arm with yoke to help the bell holding with the nut screw. We are using high-quality chrome plated metal.

CHROME PLATED CONSTRUCTION:

Made with high quality chrome-plated brass for a smooth finish. Brass is stronger and stiffer than pure copper, it is more susceptible to developing stress cracks. Compared to bronze, brass is more malleable to your operating rooms. The instrument can be sterilized and reused.

WHERE WE USE:: Its principal use is to provide adequate haemostasis throughout the neonatal circumcision procedure. Circumcision performed using the this clamp is usually quick and effective, and results in very little bleeding.



**WI-70002239****Circumcision Clamp 1.4cm Urology Instrument****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Newborn Nursery Circumcision Clamp 1.4cm (14mm) Bell Shaped Premium Quality Chrome Plated Urology Surgical Instruments.

Circumcision is the oldest known surgical procedure but it is the most frequently performed surgery around the world in males. Circumcision Clamps is a specialized device that urology surgeons use to remove the glans "foreskin" from neonatal patients, as a preventive method to minimize the risk of urinary tract infection, HIV acquisition and other benefits.

CIRCUMCISION CLAMPS FEATURES

CIRCUMCISION CLAMP: With dimensions of: 4"(10.2cm) x 1"(2.54cm) x 1.5"(3.81cm) this clamp has a

circumcision bell shaped bevel hole base plate. Rocker arm with yoke to help the bell holding with the nut screw. We are using high-quality chrome plated metal.

CHROME PLATED CONSTRUCTION:

Made with high quality chrome-plated brass for a smooth finish. Brass is stronger and stiffer than pure copper, it is more susceptible to developing stress cracks. Compared to bronze, brass is more malleable to your operating rooms. The instrument can be sterilized and reused.

WHERE WE USE:: Its principal use is to provide adequate haemostasis throughout the neonatal circumcision procedure. Circumcision performed using the this clamp is usually quick and effective, and results in very little bleeding.

**WI-70002238****Circumcision Clamp 1.6cm Urology Instrument****\$20.99** per pieceTrade price from **\$19.31**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Newborn Nursery Circumcision Clamp 1.6cm (16mm) Bell Shaped Premium Quality Chrome Plated Urology Surgical Instruments.

Circumcision is the oldest known surgical procedure but it is the most frequently performed surgery around the world in males. Circumcision Clamps is a specialized device that urology surgeons use to remove the glans "foreskin" from neonatal patients, as a preventive method to minimize the risk of urinary tract infection, HIV acquisition and other benefits.

CIRCUMCISION CLAMPS FEATURES


CIRCUMCISION CLAMP: With dimensions of: 4"(10.2cm) x 1"(2.54cm) x 1.5"(3.81cm) this clamp has a

circumcision bell shaped bevel hole base plate. Rocker arm with yoke to help the bell holding with the nut screw. We are using high-quality chrome plated metal.

CHROME PLATED CONSTRUCTION:

Made with high quality chrome-plated brass for a smooth finish. Brass is stronger and stiffer than pure copper, it is more susceptible to developing stress cracks. Compared to bronze, brass is more malleable to your operating rooms. The instrument can be sterilized and reused.

WHERE WE USE:: Its principal use is to provide adequate haemostasis throughout the neonatal circumcision procedure. Circumcision performed using the this clamp is usually quick and effective, and results in very little bleeding.



**WI-70002237****Circumcision Clamp 2.1cm Urology Instrument****\$24.49** per pieceTrade price from **\$22.53**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Newborn Nursery Circumcision Clamp 2.1cm (21mm) Bell Shaped Premium Quality Chrome Plated Urology Surgical Instruments.

Circumcision is the oldest known surgical procedure but it is the most frequently performed surgery around the world in males. Circumcision Clamps is a specialized device that urology surgeons use to remove the glans "foreskin" from neonatal patients, as a preventive method to minimize the risk of urinary tract infection, HIV acquisition and other benefits.

CIRCUMCISION CLAMPS FEATURES

CIRCUMCISION CLAMP: With dimensions of: 4"(10.2cm) x 1"(2.54cm) x 1.5"(3.81cm) this clamp has a

circumcision bell shaped bevel hole base plate. Rocker arm with yoke to help the bell holding with the nut screw. We are using high-quality chrome plated metal.

CHROME PLATED CONSTRUCTION:

Made with high quality chrome-plated brass for a smooth finish. Brass is stronger and stiffer than pure copper, it is more susceptible to developing stress cracks. Compared to bronze, brass is more malleable to your operating rooms. The instrument can be sterilized and reused.

WHERE WE USE:: Its principal use is to provide adequate haemostasis throughout the neonatal circumcision procedure. Circumcision performed using the this clamp is usually quick and effective, and results in very little bleeding.

WI-70002236**Circumcision Clamp 2.6cm Urology Instrument****\$25.49** per pieceTrade price from **\$23.45**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Newborn Nursery Circumcision Clamp 2.6cm (26mm) Bell Shaped Premium Quality Chrome Plated Urology Surgical Instruments.

Circumcision is the oldest known surgical procedure but it is the most frequently performed surgery around the world in males. Circumcision Clamps is a specialized device that urology surgeons use to remove the glans "foreskin" from neonatal patients, as a preventive method to minimize the risk of urinary tract infection, HIV acquisition and other benefits.

CIRCUMCISION CLAMPS FEATURES

CIRCUMCISION CLAMP: With dimensions of: 4"(10.2cm) x 1"(2.54cm) x 1.5"(3.81cm) this clamp has a

circumcision bell shaped bevel hole base plate. Rocker arm with yoke to help the bell holding with the nut screw. We are using high-quality chrome plated metal.

CHROME PLATED CONSTRUCTION:

Made with high quality chrome-plated brass for a smooth finish. Brass is stronger and stiffer than pure copper, it is more susceptible to developing stress cracks. Compared to bronze, brass is more malleable to your operating rooms. The instrument can be sterilized and reused.

WHERE WE USE:: Its principal use is to provide adequate haemostasis throughout the neonatal circumcision procedure. Circumcision performed using the this clamp is usually quick and effective, and results in very little bleeding.



**WI-70002235**

Circumcision Clamp 2.9cm Urology Instrument

\$25.99 per pieceTrade price from **\$23.91**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Newborn Nursery Circumcision Clamp 2.9cm (29mm) Bell Shaped Premium Quality Chrome Plated Urology Surgical Instruments.

Circumcision is the oldest known surgical procedure but it is the most frequently performed surgery around the world in males. Circumcision Clamps is a specialized device that urology surgeons use to remove the glans "foreskin" from neonatal patients, as a preventive method to minimize the risk of urinary tract infection, HIV acquisition and other benefits.

CIRCUMCISION CLAMPS FEATURES

CIRCUMCISION CLAMP: With dimensions of: 4"(10.2cm) x 1"(2.54cm) x 1.5"(3.81cm) this clamp has a

circumcision bell shaped bevel hole base plate. Rocker arm with yoke to help the bell holding with the nut screw. We are using high-quality chrome plated metal.

CHROME PLATED CONSTRUCTION: Made with high quality chrome-plated brass for a smooth finish. Brass is stronger and stiffer than pure copper, it is more susceptible to developing stress cracks. Compared to bronze, brass is more malleable to your operating rooms. The instrument can be sterilized and reused.

WHERE WE USE:: Its principal use is to provide adequate haemostasis throughout the neonatal circumcision procedure. Circumcision performed using the this clamp is usually quick and effective, and results in very little bleeding.

**WI-70002234**

Circumcision Clamp 3.5cm Urology Instrument

\$27.49 per pieceTrade price from **\$25.29**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Newborn Nursery Circumcision Clamp 3.5cm (35mm) Bell Shaped Premium Quality Chrome Plated Urology Surgical Instruments.

Circumcision is the oldest known surgical procedure but it is the most frequently performed surgery around the world in males. Circumcision Clamps is a specialized device that urology surgeons use to remove the glans "foreskin" from neonatal patients, as a preventive method to minimize the risk of urinary tract infection, HIV acquisition and other benefits.

CIRCUMCISION CLAMPS FEATURES

CIRCUMCISION CLAMP: With dimensions of: 4"(10.2cm) x 1"(2.54cm) x 1.5"(3.81cm) this clamp has a

circumcision bell shaped bevel hole base plate. Rocker arm with yoke to help the bell holding with the nut screw. We are using high-quality chrome plated metal.

CHROME PLATED CONSTRUCTION: Made with high quality chrome-plated brass for a smooth finish. Brass is stronger and stiffer than pure copper, it is more susceptible to developing stress cracks. Compared to bronze, brass is more malleable to your operating rooms. The instrument can be sterilized and reused.

WHERE WE USE:: Its principal use is to provide adequate haemostasis throughout the neonatal circumcision procedure. Circumcision performed using the this clamp is usually quick and effective, and results in very little bleeding.



**WI-70001034****Cooley Derra Pediatric Vascular Clamp 6.5"
For Anastomosis – Large****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Cooley Derra Pediatric Vascular Clamp For Anastomosis 6 1/2" (16.5cm) Large 27mm Wide Angled Jaws Curved Shank Ratched Lock Finger Ring Handle Premium Quality Stainless Steel Surgical Cardiovascular Instruments.

Cooley Derra Pediatric Vascular Clamp is a specially designed instrument to assist cardiothoracic surgeons in manipulating specific blood vessels in pediatric patients. The serrated jaws of the instrument ensure a firm but gentle grip on the vessels. Cooley-Derra Atrauma Anastomosis Clamps features a pair of jaws supplied with two angles to hold the artery. As a result, the jaws form a vaulted, empty space in which surgeons can isolate the obstruction. Also, their inner surfaces include atraumatic serrations to minimize

injuries. This clamp has a ratchet lock to achieve self-retaining property with boxlock joint permits smooth movements. Made with premium quality stainless steel ensures maximum durability and can reuse these cooley derra clamp over time after sterilization.

COOLEY-DERRA PEDIATRIC VASCULAR CLAMP FEATURES

- 6.5" (16.5cm) Overall Length
- Cooley Serration Large 27mm Wide Double Angled Jaws.
- Inner Microserrations to Reduce Endothelial Trauma.
- Curved Shank.
- Multiple Working End Sizes to Suit All Applications.
- Ratchet Ratchet Mechanism to Achieve Self-retaining Clamping

WI-70001035**Cooley Derra Pediatric Vascular Clamp 6.5"
For Anastomosis – Medium****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Cooley Derra Pediatric Vascular Clamp For Anastomosis 6 1/2" (16.5cm) Medium 21mm Wide Angled Jaws Curved Shank Ratched Lock Finger Ring Handle Premium Quality Stainless Steel Surgical Cardiovascular Instruments.

Cooley Derra Pediatric Vascular Clamp is a specially designed instrument to assist cardiothoracic surgeons in manipulating specific blood vessels in pediatric patients. The serrated jaws of the instrument ensure a firm but gentle grip on the vessels. Cooley-Derra Atrauma Anastomosis Clamps features a pair of jaws supplied with two angles to hold the artery. As a result, the jaws form a vaulted, empty space in which surgeons can isolate the obstruction. Also, their inner surfaces include

atraumatic serrations to minimize injuries. This clamp has a ratchet lock to achieve self-retaining property with boxlock joint permits smooth movements. Made with premium quality stainless steel ensures maximum durability and can reuse these cooley derra clamp over time after sterilization.

COOLEY-DERRA PEDIATRIC VASCULAR CLAMP FEATURES

- 6.5" (16.5cm) Overall Length
- Cooley Serration Medium 21mm Wide Double Angled Jaws.
- Inner Microserrations to Reduce Endothelial Trauma.
- Curved Shank.
- Multiple Working End Sizes to Suit All Applications.



WI-70001036

Cooley Derra Pediatric Vascular Clamp 6.5\"/>

\$19.99 per piece

Trade price from \$18.39
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Cooley Derra Pediatric Vascular Clamp For Anastomosis 6 1/2" (16.5cm) Small 16mm Wide Angled Jaws Curved Shank Ratched Lock Finger Ring Handle Premium Quality Stainless Steel Surgical Cardiovascular Instruments.

Cooley Derra Pediatric Vascular Clamp is a specially designed instrument to assist cardiothoracic surgeons in manipulating specific blood vessels in pediatric patients. The serrated jaws of the instrument ensure a firm but gentle grip on the vessels. Cooley-Derra Atrauma Anastomosis Clamps features a pair of jaws supplied with two angles to hold the artery. As a result, the jaws form a vaulted, empty space in which surgeons can isolate the obstruction. Also, their inner surfaces include atraumatic serrations to minimize

injuries. This clamp has a ratchet lock to achieve self-retaining property with boxlock joint permits smooth movements. Made with premium quality stainless steel ensures maximum durability and can reuse these cooley derra clamp over time after sterilization.

COOLEY-DERRA PEDIATRIC VASCULAR CLAMP FEATURES

- 6.5" (16.5cm) Overall Length
- Cooley Serration Small 16mm Wide Double Angled Jaws.
- Inner Microserrations to Reduce Endothelial Trauma.
- Curved Shank.
- Multiple Working End Sizes to Suit All Applications.
- Ratchet Ratchet Mechanism to Achieve Self-retaining Clamping



WI-35019757

Crile Artery Forceps Fine Jaw Hemostatic Clamp

\$1.30 per piece

Trade price from \$0.70
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Crile Artery Forceps Fine Jaw Hemostatic Clamp is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a fine pattern for delicate tissue.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.30, falling to \$0.70 at 10 pieces — about 46% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34901642

Crile Artery Forceps Straight Stainless Steel Reusable Surgical Hemostatic

\$1.70 per piece Trade price from **\$0.90** Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Crile Artery Forceps Straight Stainless Steel Reusable Surgical Hemostatic for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a straight working end for direct, in-line access, and full reusability — clean, autoclave and use again.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$1.70, falling to \$0.90 at 10 pieces — about 47% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-36628867

Crile Forceps Stainless Steel Hemostatic Artery Forceps Surgical Clamp

\$0.80 per piece Trade price from **\$0.70** Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Crile Forceps Stainless Steel Hemostatic Artery Forceps Surgical Clamp is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$0.80, falling to \$0.70 at 100 pieces — about 13% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.



WI-27261492

Crile Hemostat Forceps Straight Curved
Sponge Forceps

\$2.55 per piece

Trade price from **\$1.55**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Crile Hemostat Forceps Straight Curved Sponge Forceps for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.55, falling to \$1.55 at 1 piece — about 39% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35084646

Crile Hemostatic Forceps 5.5 Inch Straight
Stainless Steel Surgical

\$1.30 per piece

Trade price from **\$0.70**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Crile Hemostatic Forceps 5.5 Inch Straight Stainless Steel Surgical for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a straight working end for direct, in-line access.

The 5.5" size shown is our standard; other sizes in the same pattern are made to order. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.30, falling to \$0.70 at 10 pieces — about 46% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36622951

Crile Surgical Forceps Stainless Steel Hemostatic Clamp

\$0.90 per piece

Trade price from **\$0.70**
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Crile Surgical Forceps Stainless Steel Hemostatic Clamp for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$0.90, falling to \$0.70 at 100 pieces — about 22% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000753

DeBakey Aortic Aneurysm Clamp 10" Curved 95mm Atraumatic 1x2 Teeth Jaws

\$24.99 per piece

Trade price from **\$22.99**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



DeBakey Aortic Aneurysm Clamps 10" (25cm) 95mm Curved DeBakey Atraumatic 1x2 Teeth Jaws with Reversed Shank Curvature Premium Quality Stainless Steel Ratcheted Handle Hemostatic Cardiovascular Surgical Instrument.


DeBakey Atrauma Aortic Aneurysm Clamps are specialized hemostatic devices that allow surgeons to temporarily occlude the aortic lumen during vascular resection and anastomosis procedures. Aortic aneurysm clamp with 9.5cm (95mm) atraumatic debakey jaws designed to provide the perfect balance between control and reach, ensuring comfortable and precise tissue handling during a wide range of surgical procedures.

With its specially designed non-traumatic slight curve tips debakey aortic aneurysm forceps are tailored for delicate tissue handling. This clamp feature horizontal serrations on their inner surfaces, which help to avoid vessel injury when clamping. With elongated curve shanks with boxlock and a ratchet finger ring handle ensures self-retaining occlusions. DeBakey Atrauma Aortic Aneurysm Clamps are available in different working end sizes, ranging from 95 to 140mm. Made with premium quality stainless steel ensures exceptional precision and accuracy. Fully autoclaved and can be reuse.

DEBAKEY AORTIC ANEURYSM CLAMPS FEATURES:

- 10" (25cm) Overall Length.



**WI-70000752****DeBakey Aortic Aneurysm Clamp 12" Curved 110mm Atraumatic 1x2 Teeth Jaws****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



DeBakey Aortic Aneurysm Clamps 12" (30cm) Curved 110mm Atraumatic 1x2 Teeth Jaws with Reversed Shank Curvature Premium Quality Stainless Steel Ratcheted Handle Hemostatic Cardiovascular Surgical Instrument.

DeBakey Atrauma Aortic Aneurysm Clamps are specialized hemostatic devices that allow surgeons to temporarily occlude the aortic lumen during vascular resection and anastomosis procedures. Aortic aneurysm clamp with 11cm (110mm) atraumatic debakey jaws designed to provide the perfect balance between control and reach, ensuring comfortable and precise tissue handling during a wide range of surgical procedures.

With its specially designed non-traumatic slight curve tips debakey aortic aneurysm forceps are tailored for delicate tissue handling. This clamp feature horizontal serrations on their inner surfaces, which help to avoid vessel injury when clamping. With elongated curve shanks with boxlock and a ratchet finger ring handle ensures self-retaining occlusions. DeBakey Atrauma Aortic Aneurysm Clamps are available in different working end sizes, ranging from 95 to 140mm. Made with premium quality stainless steel ensures exceptional precision and accuracy. Fully autoclaved and can be reuse.

DEBAKEY AORTIC ANEURYSM CLAMPS FEATURES:

- 12" (30cm) Overall Length.

**WI-70000713****DeBakey Atraumatic Aortic Aneurysm Clamps 10.75" Full Curved 111mm Jaws****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

DeBakey Atraumatic Aortic Aneurysm Clamps 10.75" (27.3cm) Full Curved 111mm Atraumatic 1x2 Teeth Jaws with Curved Shanks Premium Quality Stainless Steel Ratcheted Handle Hemostatic Cardiovascular Surgical Instrument.



DeBakey Aortic Aneurysm Clamp is a specialized surgical tool that cardiothoracic surgeons commonly use for diverting the blood flow in the aorta artery in order to perform heart surgeries. The clamp has a full curved atraumatic DeBakey Jaws with curved shanks that easily adapts to the aorta artery. The inner surfaces of the jaws have fine serrations for minimizing endothelial damage.

DeBakey Aortic Aneurysm Clamp has ratchet system located between the

stems ensures self-retaining clamping with ergonomically designed finger rings which permit a firm grip. DeBakey Aortic Clamps manufactured from premium grade stainless steel to provide the highest quality and product performance during cardiovascular surgeries. Fully autoclaved and can be reuse.

FULLY CURVED DEBAKEY AORTIC ANEURYSM CLAMPS FEATURES:

- 10.75" (27.3cm) Overall Length.
- 4.37" (11.1cm) Long 1X2 Atraumatic Fine Teeth for Secure Grip.
- Fully Curved Long Jaws.
- Elongated Curved Shanks for Deep Access.
- Ratchet Step Locking for Self-retaining Properties.





Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-47174075

DeBakey Atraumatic Clamp Vascular Surgery Instrument Non-Crushing Vessel

\$17.00 per piece

Trade price from \$13.00
Minimum order 35 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of DeBakey Atraumatic Clamp Vascular Surgery Instrument Non-Crushing Vessel, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$17.00, falling to \$13.00 at 35 pieces — about 24% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-47167608

DeBakey Atraumatic Clamp Vascular Surgical Kit

\$18.00 per piece

Trade price from \$14.00
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

DeBakey Atraumatic Clamp Vascular Surgical Kit is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$18.00, falling to \$14.00 at 5 pieces — about 22% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

**WI-70000747****DeBakey Atraumatic Multi-Purpose Clamps
6.75" Curved Downwards 78mm Jaws****\$14.99** per pieceTrade price from **\$13.79**


Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

DeBakey Atraumatic Multi-Purpose Vascular Clamps 6.75" (171mm) Curved Downwards 3.1" (78mm) Long 1x2 Teeth Jaws with Curved Shanks Premium Quality Stainless Steel Ratcheted Finger Ring Handle Hemostatic Cardiovascular Surgical Instrument.

DeBakey Atraumatic Multi-Purpose Curved Clamps is a specialized instrument used in a variety of surgical procedures to temporary clamp blood vessels for hemostasis. It is also often used for clamping tissues as well as ducts during general, gynecological and cardiovascular surgeries. Atraumatic Multipurpose Vascular Clamps especially includes curved down angled jaws with curved shanks for deep access. This vascular clamp

has fine serrations on the internal side for a secure hold and smooth outer edges to protect adjacent soft tissue from damage. Moreover, it comes with a secure ratchet lock which gives it its self-retaining property. The instrument also includes ergonomic finger ring handles for comfortable handling. The DeBakey Atraumatic Multipurpose Vascular Clamps is available in multiple lengths for different surgical preferences. These clamps manufactured with high-quality German surgical stainless ensure corrosion resistance and durability. All our cardiothoracic instruments are fully autoclavable and reusable

**DEBAKEY ATRAUMATIC
MULTIPURPOSE VASCULAR
CLAMPS FEATURES:****WI-70001002****DeBakey Atraumatic Ring Handle Bulldog
Clamp 5" Curved****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

DeBakey Atraumatic Ring Handle Bulldog Clamp 5" (12.7cm) Curved 3" (76.2mm) Jaws Straight Shanks Ratcheted Handle Cardiovascular Surgical Instrument.

DeBakey Ring Handle Bulldog Curved Clamp are used to cease blood flow during cardiothoracic procedures. The ring-handles and a locking mechanism allow the instrument to clamp the vessels without any damage. The 7.62cm long atraumatic curved jaws and overall instrument length is 5" (12.7cm). DeBakey bulldog clamps are reusable and lightweight instruments which are available in various sizes to meet the surgery's requirements.

- 5" (12.7cm) Overall Length.
- 7.62cm Long 1X2 Rows of Nontraumatic Fine Teeth for Secure Grip.
- Curved Jaws.
- Straight Shanks for Deep Access.
- Ratchet Step Locking for Self-retaining Properties.
- Ergonomic Finger Ring Handles for Superior Control & Comfortable Handle.
- Rust-free and Autoclavable.
- Reusable After Sterilization.
- Made With High-quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc DeBakey Atraumatic Ring Handle Bulldog Clamp 5" Curved Jaws – Straight Shanks

**DEBAKEY RING HANDLE
BULLDOG CLAMP CURVED
FEATURES:**



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-47142823

DeBakey Atraumatic Vascular Clamp Surgical Instrument for Controlled

\$17.90 per piece

Trade price from \$13.90
Minimum order 12 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

DeBakey Atraumatic Vascular Clamp Surgical Instrument for Controlled — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$17.90, falling to \$13.90 at 12 pieces — about 22% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



WI-70000748

Debakey Atraumatic Vascular Clamps 9" Angled 80mm Jaws

\$17.99 per piece

Trade price from \$16.55
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



DeBakey Atrauma Vascular Clamps 9" (230mm) Angled 3.1" (80mm) Long 1x2 Teeth Jaws with Curved Shanks Premium Quality Stainless Steel Ratcheted Finger Ring Handle Hemostatic Cardiovascular Surgical Instrument.


DeBakey Atraumatic Vascular Clamps is a specialized cardiovascular surgical tool that commonly use for occluding the blood flow in the major cardiothoracic blood vessels, in order to prevent haemorrhage. Atraumatic Vascular Clamps have 80mm long working ends features an angled jaws that adapt to the mediastinum's anatomy. Moreover jaw's inner surface has 1x2 atraumatic serrations designed to prevent the vessel trauma. DeBakey Atrauma Vascular Clamp forceps

features curved shanks with a ratchet self-lock system provides comfort while the finger rings offer stability and reduce fatigue during lengthy procedures. These clamps manufactured with high-quality German surgical stainless ensure corrosion resistance and durability. All our cardiothoracic instruments are fully autoclavable and reusable

DEBAKEY ATRAUMATIC VASCULAR CLAMPS FEATURES:

- Crafted with 9" (230mm) Overall Length.
- 3.1" (80mm) Long 1X2 Atraumatic Fine Teeth for Secure Grip.
- Angled Jaws.
- Curved Shanks for Deep Access.





WI-34885189

Debakey Atraumatic Vascular Tissue Forceps Clamp Stainless Steel Instrument

\$17.00 per piece

Trade price from **\$11.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Debakey Atraumatic Vascular Tissue Forceps Clamp Stainless Steel Instrument is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$17.00, falling to \$11.00 at 25 pieces — about 35% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000716

Debakey Beck Atraumatic Vena-Cava Clamp 8" Double Angled 65mm Jaws

\$24.99 per piece

Trade price from **\$22.99**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

DeBaKey Beck Atraumatic Vena-Cava Clamp 8" (20.3cm) Double Angled 1x2 Atraumatic Jaws with Straight Shanks Premium Quality Stainless Steel Ratcheted Handle Hemostatic Cardiovascular Surgical Instrument.


DeBaKey-Beck atrauma anastomosis clamps are surgical tools that vascular and thoracic surgeons frequently use to grasp, hold, and manipulate blood vessels, in order to control hemostasis while performing anastomosis procedures. DeBaKey-Beck features double-angled jaws for improved sightlines and great access with fine serrations designed to achieve a strong but atraumatic grip.

DeBaKey Beck Atraumatic Vena-Cava Clamp has straight shafts to reach deep spaces. The boxlock joint that attaches

the shafts permits smooth closure and opening jaw movements. Manufactured in premium grade stainless steel to ensure the highest quality instruments in your operating rooms. Available in multiple sizes to suit different applications. Fully autoclaved and can be reuse.

DEBAKEY BECK ATRAUMATIC VENA-CAVA CLAMP FEATURES:

- 8" (20.3cm) Overall Length.
- 65mm Long 1X2 Atraumatic Fine Teeth for Secure Grip.
- Double-Angled Jaws for Optimal Clamping.
- Straight Shanks for Deep Access.
- Ratchet System for Self-Retaining Occlusions.





WI-34944365

DeBakey Bulldog Clamp Atraumatic Forceps
1.5 Straight

\$7.00 per piece

Trade price from **\$3.99**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of DeBakey Bulldog Clamp Atraumatic Forceps 1.5 Straight, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a straight working end for direct, in-line access.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$3.99 at 25 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34849412

Debakey Forceps Atraumatic Vascular Tissue
Forceps Clamp The Basis

\$17.00 per piece

Trade price from **\$11.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Debakey Forceps Atraumatic Vascular Tissue Forceps Clamp The Basis is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$17.00, falling to \$11.00 at 25 pieces — about 35% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001033****DeBakey Multi-Purpose Clamp 6.5\"/>****\$24.99** per pieceTrade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

DeBakey Atraumatic Multi-Purpose Clamp 6 1/2" (16.5cm) Angled 30°, 3.5cm Long Atraumatic Jaws Straight Shank Ratched Lock Finger Ring Handle Premium Quality Stainless Steel Surgical Cardiovascular Instruments.

The DeBakey Atraumatic Multi-Purpose Vascular Clamps are multipurpose instruments used to perform a variety of tasks during cardiothoracic surgical procedures. DeBakey is xcellent instrument commonly used during cardiovascular and thoracic surgeries. The DeBakey atraumatic clamp stunning design consists of finger-ring handles and a ratchet lock mechanism, and angled jaws aid in maneuvering some smaller spaces where other clamps cannot serve. Multiple variations of these clamps are available

for the ease of the customers' preferences. It is made from premium-grade German surgical stainless material. So, it is reusable after sterilization.

DEBAKEY MULTI-PURPOSE CLAMP FEATURES

- 6.5" (16.5cm) Overall Length
- 3.5cm DeBakey Atraumatic Serration 30° Angled Jaws.
- 1x2 DeBakey Serration to Reduce Endothelial Trauma.
- Straight Shank.
- Multiple Working End Sizes to Suit All Applications.
- Ratchet Ratchet Mechanism to Achieve Self-retaining Clamping
- Ergonomic Finger Rings Ensuring Excellent Handling.

**WI-70000751****DeBakey Peripheral Vascular Clamps 6.5\"/>****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

DeBakey Peripheral Vascular Clamps 6.5" (165mm) Curved 2.2" (56mm) Long 1x2 Atraumatic Serration Jaws Premium Quality Stainless Steel Ratcheted Finger Ring Handle Hemostatic Cardiovascular Surgical Instrument.


DeBakey Peripheral Vascular Clamps Forcep is designed for use in cardiothoracic surgical procedures where it is necessary to block or occlude blood flow. Although it is a multi-purpose instrument is is best suited for clamping of peripheral vessels. The curved jaws as well as the 1x2 jaw teeth allow the clamp to remain securely fastened at the site where it is affixed through the procedure.


Peripheral Vascular Clamps have multiple length sizes from 16.5cm to

20.3cm and with different working ends ranging from 56mm to 93mm. As a result, the surgical staff can treat both peripheral and major vessels. The ratchet lock mechanism of the instrument keeps it in its place during longer surgeries. Made with premium quality stainless steel ensures exceptional precision and accuracy. Fully autoclaved and can be reuse.

DEBAKEY PERIPHERAL VASCULAR CLAMPS FEATURES:

- 6.5" (165mm) Overall Length.
- 2.2" (56mm) Long 1X2 Atraumatic Fine Teeth for Secure Grip.
- Curved Jaws.
- Straight Shanks for Deep Access.
- Ratchet Step Locking for Self-retaining Properties.



**WI-70000750****DeBakey Peripheral Vascular Clamps 7.75\"/>****\$17.99** per pieceTrade price from **\$16.55**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

DeBakey Peripheral Vascular Clamps 7.75" (197mm) Curved 3.5" (90mm) Long 1x2 Atraumatic Serration Jaws Premium Quality Stainless Steel S-Shaped Ratcheted Finger Ring Handle Hemostatic Cardiovascular Surgical Instrument.



DeBakey Peripheral Vascular Clamps Forcep is designed for use in cardiothoracic surgical procedures where it is necessary to block or occlude blood flow. Although it is a multi-purpose instrument is best suited for clamping of peripheral vessels. The curved jaws as well as the 1x2 jaw teeth allow the clamp to remain securely fastened at the site where it is affixed through the procedure.

Peripheral Vascular Clamps have multiple length sizes from 16.5cm to

20.3cm and with different working ends ranging from 56mm to 93mm. As a result, the surgical staff can treat both peripheral and major vessels. The versatile s-shaped outline for reaching narrow spaces and ratchet lock mechanism of the instrument keeps it in its place during longer surgeries. Made with premium quality stainless steel ensures exceptional precision and accuracy. Fully autoclaved and can be reuse.

DEBAKEY PERIPHERAL VASCULAR CLAMPS FEATURES:

- 7.75" (197mm) Overall Length.
- 3.5" (90mm) Long 1X2 Atraumatic Fine Teeth for Secure Grip.
- Long Curved Jaws.
- Curved Shanks for Deep Access.

**WI-70000749****DeBakey Peripheral Vascular Clamps 8\"/>****\$18.99** per pieceTrade price from **\$17.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

DeBakey Peripheral Vascular Clamps 8" (203mm) Curved 3.6" (93mm) Long 1x2 Atraumatic Serration Jaws Premium Quality Stainless Steel S-Shaped Ratcheted Finger Ring Handle Hemostatic Cardiovascular Surgical Instrument.


DeBakey Peripheral Vascular Clamps Forcep is designed for use in cardiothoracic surgical procedures where it is necessary to block or occlude blood flow. Although it is a multi-purpose instrument is best suited for clamping of peripheral vessels. The curved jaws as well as the 1x2 jaw teeth allow the clamp to remain securely fastened at the site where it is affixed through the procedure.

Peripheral Vascular Clamps have multiple length sizes from 16.5cm to

20.3cm and with different working ends ranging from 56mm to 93mm. As a result, the surgical staff can treat both peripheral and major vessels. The versatile s-shaped outline for reaching narrow spaces and ratchet lock mechanism of the instrument keeps it in its place during longer surgeries. Made with premium quality stainless steel ensures exceptional precision and accuracy. Fully autoclaved and can be reuse.

DEBAKEY PERIPHERAL VASCULAR CLAMPS FEATURES:

- 8" (203mm) Overall Length.
- 3.6" (93mm) Long 1X2 Atraumatic Fine Teeth for Secure Grip.
- Long Curved Jaws.
- Curved Shanks for Deep Access.



**WI-70001003****DeBAKEY Ring Handled Bulldog Clamps 45° Angled 4 7/8\"/>****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

DEBAKEY Ring Handled Bulldog Clamps 4 7/8\" (12.4cm) Angled 45°, 4cm Long Atraumatic Jaws Curved Shanks Ratcheted Handle Cardiovascular Surgical Instrument.


DeBAKEY Ring Handle Bulldog Clamp are used to cease blood flow during cardiothoracic procedures. The ring-handles and a locking mechanism allow the instrument to clamp the vessels without any damage. The 4cm long atraumatic jaws make an angle of 45°. Overall instrument length is 4 7/8\". DeBAKEY bulldog clamps are reusable and lightweight instruments which are available in various sizes to meet the surgery's requirements.

DEBAKEY RING HANDLED BULLDOG CLAMPS FEATURES:

- 4 7/8\" (12.4cm) Overall Length.
- 4.0cm Long 1X2 Rows of Nontraumatic Fine Teeth for Secure Grip.
- 45 Degree Angled Jaws.
- Curved Shanks for Deep Access.
- Ratchet Step Locking for Self-retaining Properties.
- Ergonomic Finger Ring Handles for Superior Control & Comfortable Handle.
- Rust-free and Autoclavable.
- Reusable After Sterilization.
- Made With High-quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc DeBAKEY Ring Handled Bulldog Clamps Forcep 45 Degree Angle – 4 7/8\" (12.4cm) Curved Shanks

**WI-70000746****DeBAKEY Tangential Occlusion Clamps 11\" Double Angled 75mm Long x 17mm Deep Jaws****\$28.99** per pieceTrade price from **\$26.67**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


DeBAKEY Tangential Occlusion Clamps 11\" (280mm) Double Angled 75mm Long x 17mm Deep Atraumatic Jaws Slightly Curved Shanks Premium Quality Stainless Steel Ratcheted Finger Ring Handle Hemostatic Cardiovascular Surgical Instrument.

DeBAKEY Tangential Occlusion Clamp is used for partial tangential occlusion during cardiothoracic surgical procedures. Tangential Occlusion Clamp especially includes double angled blunt tips to prevent piercing injuries to tissues. The smooth outer surface of slender jaws prevent inadvertent trauma to surrounding tissue. DeBAKEY Atraumatic Tangential Clamp have ergonomic ratchet lock finger ring handles for superior control and also available in multiple sizes as

well as different working end lengths to suit different surgical needs. The high-quality German surgical stainless material is used to make these cardiothoracic instruments. So, all our instruments are reusable.

DEBAKEY TANGENTIAL OCCLUSION CLAMP FEATURES:

- Crafted with 11\" (280mm) Overall Length.
- 75mm Long x 17mm Deep Atraumatic Fine 1x2 Teeth Jaws for Secure Grip.
- Double Angled Tips For Precise Action.
- Slightly Curved Long Shanks for Deep Access.
- Ratchet Step Locking for Self-retaining Properties.



**WI-70001106****DeBakey-Satinsky Vena Cava Clamp 8.25\"/>****\$24.99** per pieceTrade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

DEBAKEY-SATINSKY Vena Cava Clamp 8.25" (21cm) Vascular Atraumatic Forceps Double Angled Curved Shank Ratcheted Handle Cardiovascular Surgical Instrument.


The Satinsky Debakey Vena Cava Clamp designed to clamping of the partial or total occlusion of blood vessels, especially the vena cava. The Debakey style toothed tips of the vascular clamp ensure a secure but gentle grip on the blood vessels, preventing any trauma or complication. The double angled jaw provides adequate grip on the vessel without causing injury to the vessel walls.

The vena cava clamp holds a section of the vessel wall while blood continues to flow through the open section. The ratcheted handle allows the surgeon to

adjust and maintain blood flow as necessary for the presenting conditions and ensure comfortable handling. The stainless steel construction ensures maximum durability and dependability, So you can reuse these debakey vena cava clamp forceps over time as they will not bend or tarnish on you.

DEBAKEY-SATINSKY VENA CAVA CLAMP FEATURES:

- 8.25" (21cm) Overall Length.
- Double Angled Narrow Tips for Precise Action.
- 1x2 Rows of Nontraumatic Debakey Styl Fine Serrated Jaws for Secure Grip.
- Slim Curved Shanks for Deep Access.

**WI-70001105****DeBakey-Satinsky Vena Cava Clamp 9.5\"/>****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

DEBAKEY-SATINSKY Vena Cava Clamp 9.5" (24.1cm) Vascular Atraumatic Forceps Double Angled Curved Shank Ratcheted Handle Cardiovascular Surgical Instrument.


The Satinsky Debakey Vena Cava Clamp designed to clamping of the partial or total occlusion of blood vessels, especially the vena cava. The Debakey style toothed tips of the vascular clamp ensure a secure but gentle grip on the blood vessels, preventing any trauma or complication. The double angled jaw provides adequate grip on the vessel without causing injury to the vessel walls.

The vena cava clamp holds a section of the vessel wall while blood continues to flow through the open section. The ratcheted handle allows the surgeon to

adjust and maintain blood flow as necessary for the presenting conditions and ensure comfortable handling. The stainless steel construction ensures maximum durability and dependability, So you can reuse these debakey vena cava clamp forceps over time as they will not bend or tarnish on you.

DEBAKEY SATINSKY VENA CAVA CLAMP FEATURES:

- 9.5" (24.1cm) Overall Length.
- Double Angled Narrow Tips for Precise Action.
- 1x2 Rows of Nontraumatic Debakey Style Fine Serrated Jaws for Secure Grip.
- Slim Curved Shanks for Deep Access.



**WI-70000851****Dieffenbach Serrefine Bulldog Clamp 2 1/4\"**
(5.7cm) Straight**\$18.99** per pieceTrade price from **\$17.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dieffenbach Serrefine Clamp 2 1/4\" (5.7cm) Straight Serrated Round Tip Jaws with Spring Action Self Locking Mechanism Stainless Steel Coronary Artery Cardiovascular Surgery Instruments.

Dieffenbach Serrefine Clamp also known as "bulldog clamps" are a versatile surgical instrument used for precise and reliable occlusion (blocking) of blood vessels or tissue to control bleeding during surgical procedures. They used in oculoplastics to secure suture and prevent entangling with other suture as is done for the rectus muscle sutures in an enucleation. With spring-loaded mechanism and a series of sharp serrated teeth on the clamp's jaws, these dieffenbach serrefine forcep

clamp grip the blood vessel or tissue firmly, while the spring-loaded mechanism keeps the clamp closed.

Bulldog Dieffenbach Serrefine clamp commonly used in vascular surgeries, including open-heart surgery, where they are used to clamp off the coronary artery during a bypass procedure. They are also used in neurosurgery and orthopedic surgery. These clamps comes with different variations to help the professionals in multiple surgical procedures. Dieffenbach Serrefine available in straight and curved profiles while straight pattern provides easy manipulation. Made with premium quality stainless steel to ensure reliable and excellent performance

DIEFFENBACH SERREFINE CLAMP FEATURES**WI-70000850****Dieffenbach Serrefine Clamp 1 1/2\"**
Curved**\$17.99** per pieceTrade price from **\$16.55**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dieffenbach Serrefine Clamp 1 1/2\" (3.8cm) Curved Serrated Round Tip Jaws with Spring Action Self Locking Mechanism Stainless Steel Coronary Artery Cardiovascular Surgery Instruments.

Dieffenbach Serrefine Clamp also known as "bulldog clamps" are a versatile surgical instrument used for precise and reliable occlusion (blocking) of blood vessels or tissue to control bleeding during surgical procedures. They used in oculoplastics to secure suture and prevent entangling with other suture as is done for the rectus muscle sutures in an enucleation. With spring-loaded mechanism and a series of sharp serrated teeth on the clamp's jaws, these dieffenbach serrefine forcep

clamp grip the blood vessel or tissue firmly, while the spring-loaded mechanism keeps the clamp closed.

Bulldog Dieffenbach Serrefine clamp commonly used in vascular surgeries, including open-heart surgery, where they are used to clamp off the coronary artery during a bypass procedure. They are also used in neurosurgery and orthopedic surgery. These clamps comes with different variations to help the professionals in multiple surgical procedures. Dieffenbach Serrefine available in straight and curved profiles while curved pattern provides easy manipulation. Made with premium quality stainless steel to ensure reliable and excellent performance

DIEFFENBACH SERREFINE CLAMP FEATURES

**WI-70000852****Dieffenbach Serrefine Clamp 1 1/2" (3.8cm) Straight****\$17.99** per pieceTrade price from **\$16.55**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dieffenbach Serrefine Clamp 1 1/2" (3.8cm) Straight Serrated Round Tip Jaws with Spring Action Self Locking Mechanism Stainless Steel Coronary Artery Cardiovascular Surgery Instruments.

Dieffenbach Serrefine Clamp also known as "bulldog clamps" are a versatile surgical instrument used for precise and reliable occlusion (blocking) of blood vessels or tissue to control bleeding during surgical procedures. They used in oculoplastics to secure suture and prevent entangling with other suture as is done for the rectus muscle sutures in an enucleation. With spring-loaded mechanism and a series of sharp serrated teeth on the clamp's jaws, these dieffenbach serrefine forcep

clamp grip the blood vessel or tissue firmly, while the spring-loaded mechanism keeps the clamp closed.

Bulldog Dieffenbach Serrefine clamp commonly used in vascular surgeries, including open-heart surgery, where they are used to clamp off the coronary artery during a bypass procedure. They are also used in neurosurgery and orthopedic surgery. These clamps comes with different variations to help the professionals in multiple surgical procedures. Dieffenbach Serrefine available in straight and curved profiles while straight pattern provides easy manipulation. Made with premium quality stainless steel to ensure reliable and excellent performance

DIEFFENBACH SERREFINE CLAMP FEATURES**WI-70000849****Dieffenbach Serrefine Clamp 2 1/4" (5.7cm) Curved****\$18.99** per pieceTrade price from **\$17.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dieffenbach Serrefine Clamp 2 1/4" (5.7cm) Curved Serrated Round Tip Jaws with Spring Action Self Locking Mechanism Stainless Steel Coronary Artery Cardiovascular Surgery Instruments.

Dieffenbach Serrefine Clamp also known as "bulldog clamps" are a versatile surgical instrument used for precise and reliable occlusion (blocking) of blood vessels or tissue to control bleeding during surgical procedures. They used in oculoplastics to secure suture and prevent entangling with other suture as is done for the rectus muscle sutures in an enucleation. With spring-loaded mechanism and a series of sharp serrated teeth on the clamp's jaws, these dieffenbach serrefine forcep

clamp grip the blood vessel or tissue firmly, while the spring-loaded mechanism keeps the clamp closed.

Bulldog Dieffenbach Serrefine clamp commonly used in vascular surgeries, including open-heart surgery, where they are used to clamp off the coronary artery during a bypass procedure. They are also used in neurosurgery and orthopedic surgery. These clamps comes with different variations to help the professionals in multiple surgical procedures. Dieffenbach Serrefine available in straight and curved profiles while curved pattern provides easy manipulation. Made with premium quality stainless steel to ensure reliable and excellent performance

DIEFFENBACH SERREFINE CLAMP FEATURES



WI-47078580

Duval Clamp Atraumatic Surgical Instrument for Gentle Tissue

\$7.90 per piece

Trade price from \$3.90
Minimum order 30 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture . Duval Clamp Atraumatic Surgical Instrument for Gentle Tissue for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$7.90, falling to \$3.90 at 30 pieces — about 51% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47091453

Duval Lung Clamp Atraumatic Surgical Tool

\$7.50 per piece

Trade price from \$3.50
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Duval Lung Clamp Atraumatic Surgical Tool — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$7.50, falling to \$3.50 at 50 pieces — about 53% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Forging, milling, heat treatment, hand finishing and passivation all happen

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001597****Endodontic Rubber Dam Clamp #0 – Premolar****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 0
Winged for Longer Premolars Teeth
Stainless Steel Dental Endodontic
Instruments

the flexibility to stretch around
specified teeth and the strength to
retain their shape while holding down
rubber dam.

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help in
preventing saliva and blood from
interfering with dental work, keep the
patient's mouth open to improve
visibility and ensure instruments are not
inhaled and do not damage the mouth.
Rubber dam clamps also help to
provide patient comfort, since a suction
in the mouth is not needed.

Our stainless steel clamps are crafted
with precise tension to give each clamp

**RUBBER DAM CLAMP #0
FEATURES:**

- For Bicuspid (Premolar).
- Flat Jaws and High Bow for Longer Premolars.
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**WI-70003674****Endodontic Rubber Dam Clamp #00 Winged****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #00 Winged Flat
Beaks High Bow for Small Lower
Canines and Premolars Teeth Premium
Quality Stainless Steel Dental
Endodontic Instruments

#00 rubber dam clamp is a dental tool
designed for isolating small lower
canines and premolars, ensuring a dry,
clear field for restorative, endodontic,
and cosmetic procedures, while
preventing debris and saliva from
interfering, thereby enhancing visibility
for the dentist. Our stainless steel
clamps are crafted with precise tension
to give each clamp the flexibility to
stretch around specified teeth and the
strength to retain their shape while
holding down rubber dam.

**RUBBER DAM CLAMP
#00 FEATURES:**

- For Small Lower Canines and Premolars
- Flat Jaws, High Bow.
- Winged, Smooth.
- Tempered Stainless Steel.
- Reusable.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**
Passivated**ULTRASONIC CLEANED:** Yes

**WI-70001603****Endodontic Rubber Dam Clamp #1 – Premolar****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 1
Winged for Upper Premolar Tooth 4-5-
12-13 – Stainless Steel Dental
Endodontic Instruments

the flexibility to stretch around
specified teeth and the strength to
retain their shape while holding down
rubber dam.

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help in
preventing saliva and blood from
interfering with dental work, keep the
patient's mouth open to improve
visibility and ensure instruments are not
inhaled and do not damage the mouth.
Rubber dam clamps also help to
provide patient comfort, since a suction
in the mouth is not needed.

Our stainless steel clamps are crafted
with precise tension to give each clamp

**RUBBER DAM CLAMP #1
FEATURES:**

- Slightly Festooned Jaws for Upper
Premolars.
- Winged, Smooth.
- Tempered stainless steel
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**TOOTH:** Upper Premolars 4-5-12-13**WI-70001592****Endodontic Rubber Dam Clamp #10 – General
Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 10 –
Winged – for Upper Left & Lower Right
Molars Teeth Stainless Steel Dental
Endodontic Instruments

the flexibility to stretch around
specified teeth and the strength to
retain their shape while holding down
rubber dam.

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help in
preventing saliva and blood from
interfering with dental work, keep the
patient's mouth open to improve
visibility and ensure instruments are not
inhaled and do not damage the mouth.
Rubber dam clamps also help to
provide patient comfort, since a suction
in the mouth is not needed.

Our stainless steel clamps are crafted
with precise tension to give each clamp

**RUBBER DAM CLAMP #10
FEATURES:**

- For Upper Left & Lower Right Molars.
- Flat Jaws Designed for Small Upper
Left and Lower Right Molars.
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures

**WI-70001591****Endodontic Rubber Dam Clamp #11 – General Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 11 – Winged – for Small Upper Right & Lower Left Molars Teeth Stainless Steel Dental Endodontic Instruments

the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

Rubber Dam Clamps are used to isolate the restorative site from the rest of the mouth when performing restorative, cosmetic, or endodontic procedures such as root canals. They help in preventing saliva and blood from interfering with dental work, keep the patient's mouth open to improve visibility and ensure instruments are not inhaled and do not damage the mouth. Rubber dam clamps also help to provide patient comfort, since a suction in the mouth is not needed.

Our stainless steel clamps are crafted with precise tension to give each clamp

RUBBER DAM CLAMP #11 FEATURES:

- For Small Upper Right & Lower Left Molars.
- Flat Jaws Designed for Small Upper Right and Lower Left Molars.
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic,**WI-70001590****Endodontic Rubber Dam Clamp #12 – Lower Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 12 – Winged – for Lower Molars Teeth Stainless Steel Dental Endodontic Instruments

the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

Rubber Dam Clamps are used to isolate the restorative site from the rest of the mouth when performing restorative, cosmetic, or endodontic procedures such as root canals. They help in preventing saliva and blood from interfering with dental work, keep the patient's mouth open to improve visibility and ensure instruments are not inhaled and do not damage the mouth. Rubber dam clamps also help to provide patient comfort, since a suction in the mouth is not needed.

Our stainless steel clamps are crafted with precise tension to give each clamp

RUBBER DAM CLAMP #12 FEATURES:

- For Lower Molars
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** Lower Molars**SIZE:** 12

**WI-70003673****Endodontic Rubber Dam Clamp #12A for Lower Right Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #12A Winged Serrated Jaws for Lower Right Molars Teeth Stainless Steel Dental Endodontic Instruments.

Rubber Dam Clamp #12A is used for isolating lower right molars during dental procedures, such as restorative, endodontic, and prosthodontic treatments. Its specifically designed shape helps to keep the surgical site clean and dry by holding the rubber dam in place and separating it from the rest of the mouth. The clamp has serrated jaws and winged jaws, making it suitable for lower molars and third molars.

Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around

specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #12A FEATURES:

- For Lower Right Molars
- Winged, Serrated.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** Lower Right Molars**SIZE:** 12A**REUSABLE/AUTOCCLAVABLE:** Yes**WI-70001589****Endodontic Rubber Dam Clamp #13 – Lower Left Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 13 – Winged – for Lower Left Molars Teeth Stainless Steel Dental Endodontic Instruments

Rubber Dam Clamps are used to isolate the restorative site from the rest of the mouth when performing restorative, cosmetic, or endodontic procedures such as root canals. They help in preventing saliva and blood from interfering with dental work, keep the patient's mouth open to improve visibility and ensure instruments are not inhaled and do not damage the mouth. Rubber dam clamps also help to provide patient comfort, since a suction in the mouth is not needed.

Our stainless steel clamps are crafted with precise tension to give each clamp

the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #13 FEATURES:

- For Lower Left Molars
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** Lower Left Molars**SIZE:** 13

**WI-70000314**

Endodontic Rubber Dam Clamp #138 for Left-Side Third Molars

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #138 Wingless with Stiff Spring and Serrated Beaks for Left Side Wisdom Teeth (Third Molars) Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber dam clamp #138 is designed for use during dental procedures involving left-sided wisdom teeth (Third Molars). It secures a rubber dam to isolate these teeth, creating a dry and clean work environment by preventing saliva and blood from affecting the treatment. clamp #138 features a wingless design tailored specifically for left-sided third molars. Rubber dam clamps made from stainless steel offer strength and flexibility, resisting corrosion. They are available in a wide variety of sizes for various applications.

RUBBER DAM CLAMP #138 FEATURES:

- Used for Left-Side Third Molars.
- Stiff Spring and Serrated Beaks.
- Wingless, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH:** Left Side Wisdom Teeth (Third Molars)**SIZE:** 138**REUSABLE:** Yes**AUTOCLAVED:** Yes**WI-70003672**

Endodontic Rubber Dam Clamp #139 for Right-Side Third Molars

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #138 Wingless with Stiff Spring and Serrated Beaks for Right Side Wisdom Teeth (Third Molars) Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber dam clamp #139 is designed for use during dental procedures involving right-sided wisdom teeth (Third Molars). It secures a rubber dam to isolate these teeth, creating a dry and clean work environment by preventing saliva and blood from affecting the treatment. Clamp #139 features a wingless design tailored specifically for right-sided third molars. Rubber dam clamps made from stainless steel offer strength and flexibility, resisting corrosion. They are available in a wide variety of sizes for various applications.

RUBBER DAM CLAMP #139 FEATURES:

- Used for Right-Side Third Molars.
- Stiff Spring and Serrated Beaks.
- Wingless, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH:** Right Side Wisdom Teeth (Third Molars)**SIZE:** 139**REUSABLE:** Yes**AUTOCLAVED:** Yes

**WI-70003671****Endodontic Rubber Dam Clamp #13A for Lower Left Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #13A Winged Serrated Jaws for Lower Left Molars Teeth Stainless Steel Dental Endodontic Instruments.

Clamp #13A is used to hold a rubber dam in place for posterior teeth, particularly the lower left molar and upper right molar. It is a winged, serrated clamp with an offset design, meaning it has a wider jaw for the buccal (cheek side) surface and a narrower jaw for the lingual (tongue side) surface. The serrations provide a more secure grip, and its features are ideal for procedures like restorative work and endodontic therapy on these specific teeth

Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around

specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #13A FEATURES:

- Used for the Lower Left Molar and Upper Right Molar
- Winged, Serrated Jaws for Improved Grip.
- Can Be Paired With Rubber Dam Clamp #12A.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**WI-70003669****Endodontic Rubber Dam Clamp #14 for General Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #14 Deeply Festooned Jaws for Partially Erupted or Irregularly Shaped Molars Stainless Steel Dental Endodontic Instruments

Rubber Dam Clamp #14 is used for general molar applications, particularly for isolating partially erupted or irregularly shaped molars during restorative, cosmetic, or endodontic (root canal) procedures. Its deeply festooned jaws and strong, flexible design ensure a secure fit while protecting the tooth and soft tissues

Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #14 FEATURES:

- Used for Partially Erupted or Irregular Molars
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** General Molars**SIZE:** 14**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**

**WI-70003670**

Endodontic Rubber Dam Clamp #14A for General Molars

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #14A Deeply Festooned Jaws for Larger Partially Erupted or Irregularly Shaped Molars Stainless Steel Dental Endodontic Instruments

Clamp #14A is used for isolating larger, partially erupted or irregularly shaped molars during dental procedures like endodontic therapy or extensive restorations. Its deeply festooned jaws are specifically designed to securely grip these challenging teeth, while other features like padded or darkly coated jaws help protect surrounding tissue and improve visibility during the procedure.

Our stainless steel rubber dam clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength

to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #14A FEATURES:

- Used for Larger Partially Erupted or Irregular Molars
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** General Molars**SIZE:** 14A**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**WI-70003668**

Endodontic Rubber Dam Clamp #15 for Big Lower Molars

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #15 Wingless Smooth Wide Beaks for Larger Lower Molars Premium Quality Stainless Steel Dental Endodontic Instruments

Clamp #15 is a wingless clamp designed for isolating large lower molars during dental procedures such as endodontic treatments and restorations. Its wide beaks ensure a stable and secure fit, effectively creating a dry and isolated field by preventing saliva and debris from contaminating the working area. Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #15 FEATURES:

- Used for Larger Lower Molars.
- Wingless Smooth, Wide Beak.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** Big Lower Molars**SIZE:** 15**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test,

Performance Test, Shape Test

QC PASSED: Yes

**WI-70003667****Endodontic Rubber Dam Clamp #18 for Upper & Lower Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Endodontic Rubber Dam Clamp #18 Wingless Flat Jaws and Stiff Bow for Upper and Lower Molars Premium Quality Stainless Steel Dental Instruments

Rubber Dam Clamp #18 is designed to secure a dental dam around a molar tooth, providing isolation for various dental procedures such as restorative, endodontic, and fixed prosthetic treatments. Its flat jaws and stiff bow design enhance stability and suitability for molars. Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

- Used for Upper & Lower Molars.
- Wingless Smooth.
- Flat Jaws and Stiff Bow.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** For Molars**SIZE:** 18**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**RUBBER DAM CLAMPS #18
FEATURES:****WI-70001595****Endodontic Rubber Dam Clamp #1A – Canines and Premolars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 1A – Winged – for Canines and Premolars Teeth Stainless Steel Dental Endodontic Instruments

Rubber Dam Clamps are used to isolate the restorative site from the rest of the mouth when performing restorative, cosmetic, or endodontic procedures such as root canals. They help in preventing saliva and blood from interfering with dental work, keep the patient's mouth open to improve visibility and ensure instruments are not inhaled and do not damage the mouth. Rubber dam clamps also help to provide patient comfort, since a suction in the mouth is not needed.

Our stainless steel clamps are crafted with precise tension to give each clamp

the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

**RUBBER DAM CLAMP #1A
FEATURES:**

- For Canines and Premolars
- Deeply Festooned Jaws for Canines and Premolars
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures

**WI-70001594****Endodontic Rubber Dam Clamp #1T – Edged Premolars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 1T – Winged – for Edged Premolars Teeth
Stainless Steel Dental Endodontic Instruments

the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

Rubber Dam Clamps are used to isolate the restorative site from the rest of the mouth when performing restorative, cosmetic, or endodontic procedures such as root canals. They help in preventing saliva and blood from interfering with dental work, keep the patient's mouth open to improve visibility and ensure instruments are not inhaled and do not damage the mouth. Rubber dam clamps also help to provide patient comfort, since a suction in the mouth is not needed.

Our stainless steel clamps are crafted with precise tension to give each clamp

RUBBER DAM CLAMP #1T FEATURES:

- For Edged Premolars
- Serrated Rubber Dam Clamp for general use upper with adjustable jaws according to gingival line shapes.
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**WI-70001602****Endodontic Rubber Dam Clamp #2 – Premolar****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 2 Winged for Lower Premolar Tooth
Stainless Steel Dental Endodontic Instruments

the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

Rubber Dam Clamps are used to isolate the restorative site from the rest of the mouth when performing restorative, cosmetic, or endodontic procedures such as root canals. They help in preventing saliva and blood from interfering with dental work, keep the patient's mouth open to improve visibility and ensure instruments are not inhaled and do not damage the mouth. Rubber dam clamps also help to provide patient comfort, since a suction in the mouth is not needed.

Our stainless steel clamps are crafted with precise tension to give each clamp

RUBBER DAM CLAMP #2 FEATURES:

- Slightly Festooned Jaws for Lower Premolars.
- Winged, Smooth.
- Tempered stainless steel
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TOOTH:** Lower Premolars

**WI-70003666****Endodontic Rubber Dam Clamp #200 for Lower Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #200 Winged Smooth Beaks for Lower Molars Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber dam clamp 200 is used to isolate lower molars in endodontic or restorative dental procedures. Its winged design secures the rubber dam sheet, ensuring a dry and clean operating field by isolating the tooth from the mouth. Our clamps made from stainless steel offer strength and flexibility, resisting corrosion. They are available in a wide variety of sizes for various applications.

RUBBER DAM CLAMP #200 FEATURES:

- Used for Lower Molars.

- Wingless, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH:** Lower Molars)**SIZE:** 200**REUSABLE:** Yes**AUTOCLAVED:** Yes**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing**WI-70003665****Endodontic Rubber Dam Clamp #201 for Upper Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #201 Winged Slightly Festooned Jaws for Upper Molars Premium Quality Stainless Steel Dental Endodontic Instruments.

Endodontic Rubber dam clamp #201 secures a rubber dam to isolate an upper molar during dental procedures, featuring festooned jaws for a secure fit, enabling a clean and dry environment for treatments such as fillings and root canals. Our tempered stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam. We offer a wide range of rubber dam clamps suitable for anterior, premolar, and molar applications, made from tempered stainless steel for superior

performance, durability, and corrosion resistance.

RUBBER DAM CLAMPS #201 FEATURES:

- Used for Upper Molars.
- Festooned Jaws.
- Winged Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH:** Upper Molars**SIZE:** 201**REUSABLE/AUTOCLAVABLE:** Yes**STERILE:** Non-Sterile

**WI-70000313****Endodontic Rubber Dam Clamp #202 for Large Lower Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #202 Winged with Large Flat Jaws for Large Adult Lower Molars Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber Dam Clamp #202 is designed for anchoring the rubber dam during dental procedures on large lower molars, specifically for teeth 17-18-19 and 30-31-32. Its large, flat jaws ensure a secure grip, allowing for a dry and isolated working field for the dentist. We offer a wide range of rubber dam clamps suitable for anterior, premolar, and molar applications, made from tempered stainless steel for superior performance, durability, and corrosion resistance.

RUBBER DAM CLAMP #202 FEATURES:

- Used for Large Lower Molars.
- Large, Flat Jaws.
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH:** Large Lower Molars)**SIZE:** 202**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:** Passivated**WI-70003664****Endodontic Rubber Dam Clamp #203 for Small Lower Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #203 Winged for Small Lower Molars Including Posterior for Children Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber Dam Clamp #203 is designed to secure a rubber dam to small lower molars, including those of children or adult premolars, effectively isolating the tooth during dental procedures. Clamp #203 is a winged, smooth clamp used for the lower right and upper left sides, and functions as a general posterior clamp for children. We offer a wide range of rubber dam clamps suitable for anterior, premolar, and molar applications, made from tempered stainless steel for superior performance, durability, and corrosion resistance.

RUBBER DAM CLAMPS #203 FEATURES:

- Used for Small Lower Molars.
- Winged Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH:** Small Lower Molars**SIZE:** 203**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes

**WI-70003662****Endodontic Rubber Dam Clamp #205****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Endodontic Rubber Dam Clamp #205
Winged Festooned Jaws for Large
Upper Molars Teeth Premium Quality
Stainless Steel Dental Instruments.

Rubber Dam Clamp 205 is designed for
isolating the cervical area of large
upper molars during dental procedures,
featuring slightly curved, festooned
jaws that securely and comfortably grip
these specific teeth. Our stainless steel
clamps are crafted with precise tension
to give each clamp the flexibility to
stretch around specified teeth and the
strength to retain their shape while
holding down rubber dam.

**RUBBER DAM CLAMP #205
FEATURES:**

- For Large Upper Molars.
- Winged, Smooth.

- Beveled, Festooned Jaws
- Autoclavable.
- Tempered Stainless Steel.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**
Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test,
Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag PackingExplore more **Rubber Dam Clamps****WI-70003661****Endodontic Rubber Dam Clamp #206****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Endodontic Rubber Dam Clamp #206
Winged for Upper and Lower Bicuspid
(Premolars) Teeth Premium Quality
Stainless Steel Dental Instruments.

Rubber Dam Clamp #206 is designed
primarily for use on premolars in both
upper and lower arches. Its winged
design facilitates easier placement and
removal, securely holding the dam
material. Our stainless steel clamps are
crafted with precise tension to give
each clamp the flexibility to stretch
around specified teeth and the strength
to retain their shape while holding down
rubber dam.

**RUBBER DAM CLAMP #206
FEATURES:**

- For Upper and Lower Premolars
(Bicuspid).

- Winged, Smooth.
- Autoclavable.
- Tempered Stainless Steel.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**
Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test,
Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag PackingExplore more **Rubber Dam Clamps**

**WI-70003660****Endodontic Rubber Dam Clamp #208****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Endodontic Rubber Dam Clamp #208 Winged for Big Bicuspid & Primary Molars Teeth Premium Quality Stainless Steel Dental Instruments.

Rubber Dam Clamp 208 is utilized in dentistry to anchor a rubber dam around the neck of large premolars (bicuspid) or primary molars, ensuring a clean and dry field for procedures such as fillings and root canals, while preventing saliva and blood interference and enhancing visibility. Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #208 FEATURES:

- For Big Bicuspid.
- Winged, Smooth.
- Autoclavable.
- Tempered Stainless Steel.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing**WI-70003647****Endodontic Rubber Dam Clamp #209 for Premolars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #209 Winged Flat Beaks for Premolars Teeth Premium Quality Stainless Steel Dental Endodontic Instruments.

#209 rubber dam clamp is a dental instrument designed to isolate premolars or sometimes small pediatric anterior teeth (Front Teeth). It features small, flat jaws that securely grip these tooth structures during procedures such as fillings or root canals, ensuring the area remains dry and clean. Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #209 FEATURES:

- For Anterior, Premolars and Molars.
- Small, Flat Jaws
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test

**WI-70000312****Endodontic Rubber Dam Clamp #210 for Upper Centrals and Incisors****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #210 Winged Anterior Upper Centrals and Incisors Teeth Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber dam clamp #210 is designed to anchor a rubber dam to the anterior teeth, specifically the upper centrals and incisors, providing a stable and secure grip in the cervical region to facilitate dental procedures such as fillings or endodontic treatments. We offer a wide range of rubber dam clamps suitable for anterior, premolar, and molar applications, made from tempered stainless steel for superior performance, durability, and corrosion resistance.

RUBBER DAM CLAMP #210 FEATURES:

- Used for Anterior Upper Centrals and Incisors.
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH:** Anterior Upper Centrals**SIZE:** 210**REUSABLE:** Yes**AUTOCLAVED:** Yes**RUST PREVENTION PROCEDURE:** Passivated**WI-70000311****Endodontic Rubber Dam Clamp #211 for Adult Lower Anterior****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #211 Winged for Adult Anterior Lower Jaw Premium Quality Stainless Steel Dental Endodontic Instruments.

#211 rubber dam clamp is utilized to isolate adult anterior teeth, primarily on the lower jaw, during dental procedures such as endodontic and restorative work. It secures the rubber dam at the tooth's cervical area, ensuring a dry and isolated operating field. We offer a wide range of rubber dam clamps suitable for anterior, premolar, and molar applications, made from tempered stainless steel for superior performance, durability, and corrosion resistance.

RUBBER DAM CLAMP #211 FEATURES:

- Used for Adult Anterior Lower Jaw.
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH:** Universal for Anterior**SIZE:** 211**REUSABLE/AUTOCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes

**WI-70001586****Endodontic Rubber Dam Clamp #212 – Anterior****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 212 – Wingless – for Labial Caries on Anterior Teeth Stainless Steel Dental Endodontic Instruments

Rubber Dam Clamps are used to isolate the restorative site from the rest of the mouth when performing restorative, cosmetic, or endodontic procedures such as root canals. They help in preventing saliva and blood from interfering with dental work, keep the patient's mouth open to improve visibility and ensure instruments are not inhaled and do not damage the mouth. Rubber dam clamps also help to provide patient comfort, since a suction in the mouth is not needed.

Our stainless steel clamps are crafted with precise tension to give each clamp

the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #212 FEATURES:

- For Cervical Caries on Anterior Teeth
- Flat jaws for labial caries on anterior teeth
- Wingless, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**WI-70000310****Endodontic Rubber Dam Clamp #212SA for Anteriors and Premolars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #212SA Winged, Butterfly-shaped for Anteriors and Premolars Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber Dam Clamp #212SA is a specialized winged, butterfly-shaped clamp designed to hold dental dams. It is primarily used for isolating and retracting anterior teeth, facilitating access to conditions such as labial caries (cavities on the front of the tooth). The "SA" designation typically indicates a stronger pressure on the gingiva (gums) than the standard #212 clamp. We offer a wide range of rubber dam clamps suitable for anterior, premolar, and molar applications, made from tempered stainless steel for superior performance, durability, and corrosion resistance.

RUBBER DAM CLAMP #212SA FEATURES:

- Used for Anteriors and Premolars.
- Winged, Butterfly-shaped.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH:** Dor Anteriors**SIZE:** 212SA**REUSABLE:** Yes**AUTOCLAVED:** Yes**RUST PREVENTION PROCEDURE:** Passivated

**WI-70000316****Endodontic Rubber Dam Clamp #22 for Upper & Lower Bicuspid****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #22 Wingless For Upper & Lower Bicuspid Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber Dam Clamp #22 is used for isolating premolars, as it is a wingless clamp specifically designed for upper and lower bicuspid. Its purpose is to hold a rubber dam sheet in place, creating a dry and isolated field for dental procedures like fillings, crowns, and root canals. This isolation prevents saliva, blood, and debris from interfering with the work, improves visibility for the dentist, and protects the patient from accidentally swallowing instruments. Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified

teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMPS #22 FEATURES:

- Used for Upper & Lower Bicuspid (Premolars)
- Wingless Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** For Premolars**SIZE:** 22**REUSABLE/AUTOCALVABLE:** Yes**STERILE:** Non-Sterile**WI-70001585****Endodontic Rubber Dam Clamp #26N – Upper Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp 26N Wingless Flat Jaws for Upper Molar Teeth Premium Quality Stainless Steel Dental Endodontic Instruments

Endodontic Rubber Dam Clamp #26N is designed for isolating upper molars during dental procedures, such as restorations and endodontics. It maintains a dry working area, enhances visibility for the dentist, and promotes patient comfort. Clamp #26N features flat jaws that securely grip the surfaces of upper molars. Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #26N FEATURES:

- For Upper Molars
- Flat Jaws to Secure the Rubber Dam to the Upper Molars
- Wingless, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** Upper Molars**SIZE:** 26N**REUSABLE:** Yes**AUTOCALVED:** Yes**RUST PREVENTION PROCEDURE:**

**WI-70000315****Endodontic Rubber Dam Clamp #27 for Large Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp 27 Wingless Flat Jaws for Large Molars Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber Dam Clamp #27 is designed to anchor the rubber dam for isolating teeth during dental procedures, primarily restorations and endodontics. It is suitable for distal access on erupted molars and can be used in conjunction with another clamp in a "piggybacking" technique. Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #27 FEATURES:

- Used for Large Molars.
- Flat Jaws.
- Wingless, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** Large Molars**SIZE:** 27**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:** Passivated**WI-70001584****Endodontic Rubber Dam Clamp #27N – Premolars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 27N – Wingless – for Premolar Teeth Stainless Steel Dental Endodontic Instruments

Rubber Dam Clamps are used to isolate the restorative site from the rest of the mouth when performing restorative, cosmetic, or endodontic procedures such as root canals. They help in preventing saliva and blood from interfering with dental work, keep the patient's mouth open to improve visibility and ensure instruments are not inhaled and do not damage the mouth. Rubber dam clamps also help to provide patient comfort, since a suction in the mouth is not needed.

Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around

specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #27N FEATURES:

- For Upper Premolars
- Slightly Festooned Jaws for Premolars
- Wingless, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** Premolars

**WI-70003655**

Endodontic Rubber Dam Clamp #28 for Lower Molars

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #28 Wingless Flat Jaws for Lower Molars Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber dam clamp #28 is a wingless dental tool used to stabilize and hold a rubber dam on lower molars during procedures. Its flat jaws conform to the tooth's shape, ensuring smooth operation and reducing the risk of pinching soft tissue. Rubber dam clamps are typically crafted from tempered stainless steel, ensuring both durability and corrosion resistance.

RUBBER DAM CLAMP #28 FEATURES:

- Used for Lower Molars.
- Flat Jaws.
- Wingless, Smooth.

- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** Lower Molars**SIZE:** 28**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing**WI-70003654**

Endodontic Rubber Dam Clamp #29 for Lower Premolars

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #29 Wingless Flat Jaws for Lower Premolars Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber Dam Clamp #29 is designed for anchoring a dental rubber dam to lower premolars during endodontic procedures. It ensures precise tooth isolation, maintaining a dry and clean working area, enhances visibility for the dentist, and prevents the swallowing of instruments. Rubber dam clamps are typically crafted from tempered stainless steel, ensuring both durability and corrosion resistance.

RUBBER DAM CLAMP #29 FEATURES:

- Used for Lower Premolars.
- Flat Jaws.

- Wingless, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** Lower Premolars**SIZE:** 29**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing

**WI-70001601****Endodontic Rubber Dam Clamp #2A – Larger Premolar****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 2A
Winged for Larger Premolar Tooth
Stainless Steel Dental Endodontic
Instruments

the flexibility to stretch around
specified teeth and the strength to
retain their shape while holding down
rubber dam.

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help in
preventing saliva and blood from
interfering with dental work, keep the
patient's mouth open to improve
visibility and ensure instruments are not
inhaled and do not damage the mouth.
Rubber dam clamps also help to
provide patient comfort, since a suction
in the mouth is not needed.

Our stainless steel clamps are crafted
with precise tension to give each clamp

**RUBBER DAM CLAMP #2A
FEATURES:**

- For Larger Premolars.
- Winged, Smooth.
- Tempered stainless steel
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**TOOTH:** Larger Premolars**SIZE:** 2A**WI-70003653****Endodontic Rubber Dam Clamp #3 for Small Lower Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #3 Flat Narrow
Jaws for Small Lower Molars Stainless
Steel Dental Endodontic Instruments

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help in
preventing saliva and blood from
interfering with dental work, keep the
patient's mouth open to improve
visibility and ensure instruments are not
inhaled and do not damage the mouth.
Rubber dam clamps also help to
provide patient comfort, since a suction
in the mouth is not needed.

Rubber dam clamp #3 is primarily used
to isolate and secure the rubber dam to
small lower molars during dental
procedures. Its flat, narrow jaws are

specifically designed to fit these teeth,
such as the #17, #18, #19, #30, #31, and
#32 molars. This clamp helps provide a
dry and stable working field by keeping
saliva and debris out of the treatment
area. Made with stainless steel crafted
with precise tension to give each clamp
the flexibility to stretch around
specified teeth and the strength to
retain their shape while holding down
rubber dam.

**RUBBER DAM CLAMP #3
FEATURES:**

- For Small Lower Molars
- Flat Narrow Jaws.
- Tempered Stainless Steel.
- Autoclavable.
- Rust Free.

TECHNICAL INFORMATION:

**WI-70003652****Endodontic Rubber Dam Clamp #30 for Upper Right Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #30 Offset Jaws for Upper Right Molars Premium Quality Stainless Steel Dental Endodontic Instruments

Rubber dam clamp #30 is designed to secure a dental dam around an upper right molar or a lower molar in its wingless design. Its main function is to isolate the targeted tooth or teeth during dental procedures, such as endodontics and restorative work, by separating the surgical site from the rest of the mouth.. Made with stainless steel crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

- For Upper Right Molars
- Offset Jaws.
- Wingless Design.
- Tempered Stainless Steel.
- Autoclavable.
- Rust Free.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**SIZE:** 30**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**RUBBER DAM CLAMP #30
FEATURES:****WI-70003651****Endodontic Rubber Dam Clamp #31 for Upper Left Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #31 Offset Jaws for Upper Left Molars Premium Quality Stainless Steel Dental Endodontic Instruments

Clamp #31 is designed to secure a dental dam around an upper left molar or a lower molar in its wingless design. Its main function is to isolate the targeted tooth or teeth during dental procedures, such as endodontics and restorative work, by separating the surgical site from the rest of the mouth.. Made with stainless steel crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

- For Upper Left Molars
- Offset Jaws.
- Wingless Design.
- Tempered Stainless Steel.
- Autoclavable.
- Rust Free.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**SIZE:** 31**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**RUBBER DAM CLAMP #31
FEATURES:**



WI-70003650

Endodontic Rubber Dam Clamp #4

\$8.99 per piece

Trade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Endodontic Rubber Dam Clamp #4 is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$8.99, falling to \$8.27 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003649

Endodontic Rubber Dam Clamp #5

\$8.99 per piece

Trade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Endodontic Rubber Dam Clamp #5 is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$8.99, falling to \$8.27 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003648****Endodontic Rubber Dam Clamp #51 for Upper & Lower Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #51 Wingless Jaws for Upper & Lower Molars Premium Quality Stainless Steel Dental Endodontic Instruments.

Endodontic Rubber Dam Clamp 51 is a stainless steel, wingless clamp utilized to secure a dental dam, effectively isolating and stabilizing upper and lower molars during restorative dental procedures. It is designed for durability, cleanliness, and can be autoclaved. Made with stainless steel crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #51 FEATURES:

- For for Upper & Lower Molars.

- Wingless Design.
- Tempered Stainless Steel.
- Autoclavable.
- Rust Free.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**SIZE:** 51**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing**WI-70001599****Endodontic Rubber Dam Clamp #7 – Lower Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 7 Winged for Lower Molars Teeth Stainless Steel Dental Endodontic Instruments

Rubber Dam Clamps are used to isolate the restorative site from the rest of the mouth when performing restorative, cosmetic, or endodontic procedures such as root canals. They help in preventing saliva and blood from interfering with dental work, keep the patient's mouth open to improve visibility and ensure instruments are not inhaled and do not damage the mouth. Rubber dam clamps also help to provide patient comfort, since a suction in the mouth is not needed.

Our stainless steel clamps are crafted with precise tension to give each clamp

the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #7 FEATURES:

- For Lower Molars
- Winged, Smooth.
- Tempered stainless steel
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TOOTH** Lower Molars**SIZE:** 7

**WI-70001598****Endodontic Rubber Dam Clamp #7A – Larger Lower Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 7A
Winged for Larger Lower Molars Teeth
Stainless Steel Dental Endodontic
Instruments

the flexibility to stretch around
specified teeth and the strength to
retain their shape while holding down
rubber dam.

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help in
preventing saliva and blood from
interfering with dental work, keep the
patient's mouth open to improve
visibility and ensure instruments are not
inhaled and do not damage the mouth.
Rubber dam clamps also help to
provide patient comfort, since a suction
in the mouth is not needed.

Our stainless steel clamps are crafted
with precise tension to give each clamp

**RUBBER DAM CLAMP #7
FEATURES:**

- For Larger Lower Molars
- Winged, Smooth.
- Tempered stainless steel
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**TEETH** Larger Lower Molars**SIZE:** 7A**WI-70003646****Endodontic Rubber Dam Clamp #8A****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #8A Winged for
Upper & Lower Molars Teeth Premium
Quality Stainless Steel Dental
Endodontic Instruments.

8A rubber dam clamp is utilized for
upper and lower molar roots,
particularly on smaller or partially
erupted molars with compromised
crowns. It aims to maintain a dry, stable
working area for dental procedures
such as fillings, root canals, and
restorations, improving visibility and
reducing contamination. Our stainless
steel clamps are crafted with precise
tension to give each clamp the
flexibility to stretch around specified
teeth and the strength to retain their
shape while holding down rubber dam.

**RUBBER DAM CLAMP #8A
FEATURES:**

- For Upper & Lower Molars.
- Deeply Festooned Four-point Jaws
for Molars
- Winged, Smooth.
- Autoclavable.
- Tempered Stainless Steel.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**
Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test,
Performance Test, Shape Test

**WI-70001587****Endodontic Rubber Dam Clamp #8AD – Molar****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 8AD – Winged – for Molars teeth Stainless Steel Dental Endodontic Instruments


specified teeth and the strength to retain their shape while holding down rubber dam.

Rubber Dam Clamps are used to isolate the restorative site from the rest of the mouth when performing restorative, cosmetic, or endodontic procedures such as root canals. They help in preventing saliva and blood from interfering with dental work, keep the patient's mouth open to improve visibility and ensure instruments are not inhaled and do not damage the mouth. Rubber dam clamps also help to provide patient comfort, since a suction in the mouth is not needed.

Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around

RUBBER DAM CLAMP #8AD FEATURES:

- For General Molars
- Deeply Festooned Jaws and Extended Bow for Distal Access to Molars.
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**WI-70001593****Endodontic Rubber Dam Clamp #9 – Anterior****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 9 – Winged – for Labial Caries on Anterior Teeth Stainless Steel Dental Endodontic Instruments

the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

Rubber Dam Clamps are used to isolate the restorative site from the rest of the mouth when performing restorative, cosmetic, or endodontic procedures such as root canals. They help in preventing saliva and blood from interfering with dental work, keep the patient's mouth open to improve visibility and ensure instruments are not inhaled and do not damage the mouth. Rubber dam clamps also help to provide patient comfort, since a suction in the mouth is not needed.

Our stainless steel clamps are crafted with precise tension to give each clamp

RUBBER DAM CLAMP #9 FEATURES:

- For Labial Caries on Anterior Teeth
- Offset Flat Jaws and Stiff Bow for Upper and Lower Anteriors Teeth.
- Winged, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures



WI-70001607

Endodontic Rubber Dam Clamp #W0

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Ivory Rubber Dam Clamp
W0 Wingless Bicuspid Stainless Steel
Dental Endodontic Instruments

specified teeth and the strength to
retain their shape while holding down
rubber dam.

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help in
preventing saliva and blood from
interfering with dental work, keep the
patient's mouth open to improve
visibility and ensure instruments are not
inhaled and do not damage the mouth.
Rubber dam clamps also help to
provide patient comfort, since a suction
in the mouth is not needed.

Our stainless steel clamps are crafted
with precise tension to give each clamp
the flexibility to stretch around

RUBBER DAM CLAMP #W0
FEATURES:

- For upper & lower Bicuspids.
- Wingless, smooth.
- Tempered stainless steel uniquely
designed to retain shape.
- Autoclavable.

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Endodontic
Instruments

APPLICATION: Restorative, Cosmetic,
or Endodontic Procedures

TOOTH: Bicuspids

SIZE: W0



WI-70003645

Endodontic Rubber Dam Clamp #W14 –
Wingless

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Endodontic Rubber Dam Clamp #W14 –
Wingless is an instrument we forge,
grind and finish in our own Sialkot
workshop. The steel is heat treated and
passivated for corrosion resistance,
then each instrument is checked by
hand for action and finish.

Listed at \$8.99, falling to \$8.27 at 21
pieces - about 8% off the single-unit
price. Larger quantities are priced on
request.

Small trial orders are welcome — most
buyers start with a sample or a single
carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001606****Endodontic Rubber Dam Clamp #W3****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp W3
Wingless For Small Lower Molars
Stainless Steel Dental Endodontic
Instruments

the flexibility to stretch around
specified teeth and the strength to
retain their shape while holding down
rubber dam.

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help in
preventing saliva and blood from
interfering with dental work, keep the
patient's mouth open to improve
visibility and ensure instruments are not
inhaled and do not damage the mouth.
Rubber dam clamps also help to
provide patient comfort, since a suction
in the mouth is not needed.

Our stainless steel clamps are crafted
with precise tension to give each clamp

**RUBBER DAM CLAMP #W3
FEATURES:**

- For small lower molars.
- Wingless, smooth.
- Tempered stainless steel
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**TOOTH:** Small Lower Molars**SIZE:** W3**WI-70001605****Endodontic Rubber Dam Clamp #W4****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp W4
Wingless For Small Upper Molars
Stainless Steel Dental Endodontic
Instruments

the flexibility to stretch around
specified teeth and the strength to
retain their shape while holding down
rubber dam.

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help in
preventing saliva and blood from
interfering with dental work, keep the
patient's mouth open to improve
visibility and ensure instruments are not
inhaled and do not damage the mouth.
Rubber dam clamps also help to
provide patient comfort, since a suction
in the mouth is not needed.

Our stainless steel clamps are crafted
with precise tension to give each clamp

**RUBBER DAM CLAMP #W4
FEATURES:**

- Flat Narrow Jaws For Small Upper
Molars.
- Wingless, smooth.
- Tempered stainless steel
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**TOOTH:** Small Upper Molars



WI-70001604

Endodontic Rubber Dam Clamp #W7

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp W7
Wingless For Lower Molars Stainless
Steel Dental Endodontic Instruments

specified teeth and the strength to
retain their shape while holding down
rubber dam.

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help in
preventing saliva and blood from
interfering with dental work, keep the
patient's mouth open to improve
visibility and ensure instruments are not
inhaled and do not damage the mouth.
Rubber dam clamps also help to
provide patient comfort, since a suction
in the mouth is not needed.

Our stainless steel clamps are crafted
with precise tension to give each clamp
the flexibility to stretch around

RUBBER DAM CLAMP #W7
FEATURES:

- For Lower molars.
- Wingless, smooth.
- Tempered stainless steel
- Autoclavable.

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel
INSTRUMENT TYPE: Endodontic
Instruments
APPLICATION: Restorative, Cosmetic,
or Endodontic Procedures
TOOTH: Lower Molars
SIZE: W7
REUSABLE: Yes



WI-70003644

Endodontic Rubber Dam Clamp #W8

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Endodontic Rubber Dam Clamp #W8 is
an instrument we forge, grind and finish
in our own Sialkot workshop. The steel
is heat treated and passivated for
corrosion resistance, then each
instrument is checked by hand for
action and finish.

Listed at \$8.99, falling to \$8.27 at 21
pieces - about 8% off the single-unit
price. Larger quantities are priced on
request.

Small trial orders are welcome — most
buyers start with a sample or a single
carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001600****Endodontic Rubber Dam Clamp Labial #6****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Rubber Dam Clamp 6
Winged for Labial Caries on Anterior
Teeth Stainless Steel Dental Endodontic
Instruments

the flexibility to stretch around
specified teeth and the strength to
retain their shape while holding down
rubber dam.

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help in
preventing saliva and blood from
interfering with dental work, keep the
patient's mouth open to improve
visibility and ensure instruments are not
inhaled and do not damage the mouth.
Rubber dam clamps also help to
provide patient comfort, since a suction
in the mouth is not needed.

Our stainless steel clamps are crafted
with precise tension to give each clamp

**RUBBER DAM CLAMP #6
FEATURES:**

- For Labial Caries on Anterior Teeth
- Winged, Smooth.
- Tempered stainless steel
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**TOOTH** Labial Caries on Anterior Teeth
SIZE: 6**WI-70001608****Endodontic Rubber Dam Clamps #W1****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

1 Pieces of Ivory Rubber Dam Clamp
W1 Wingless Bicuspid Stainless Steel
Dental Endodontic Instruments

around specified teeth and the strength
to retain their shape while holding down
rubber dam.

Rubber Dam Clamps are used to isolate
the restorative site from the rest of the
mouth when performing restorative,
cosmetic, or endodontic procedures
such as root canals. They help prevent
saliva and blood from interfering with
dental work, keep the patient's mouth
open to improve visibility and ensure
instruments are not inhaled and do not
damage the mouth. Rubber dam
clamps also help provide patient
comfort, since a suction in the mouth is
not needed.

Our stainless steel rubber dam clamps
are crafted with precise tension to give
each clamp the flexibility to stretch

**RUBBER DAM CLAMP #W1
FEATURES:**

- For upper & lower Bicuspids.
- Wingless, smooth.
- Tempered stainless steel uniquely
designed to retain shape.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic
Instruments**APPLICATION:** Restorative, Cosmetic,
or Endodontic Procedures**TOOTH:** Bicuspids**SIZE:** W1



WI-46417415

Ferguson Angiotribe Forceps Vascular Surgical Kit Steel Components

\$13.00 per piece

Trade price from **\$7.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Ferguson Angiotribe Forceps Vascular Surgical Kit Steel Components is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$13.00, falling to \$7.00 at 10 pieces — about 46% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47073655

Fogarty Clamp Surgical Kit Vascular Surgery Instrument

\$13.00 per piece

Trade price from **\$7.00**
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Fogarty Clamp Surgical Kit Vascular Surgery Instrument is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$13.00, falling to \$7.00 at 50 pieces — about 46% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47112217

Fogarty Clamp Vascular Surgical Instrument Set for Precise Vessel

\$13.90 per piece

Trade price from **\$8.90**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Fogarty Clamp Vascular Surgical Instrument Set for Precise Vessel for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$13.90, falling to \$8.90 at 25 pieces — about 36% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47117134

Fogarty Vascular Clamp (Model Fogarty Clamp-01) Steel Surgical Instrument

\$14.00 per piece

Trade price from **\$9.00**
Minimum order 12 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Fogarty Vascular Clamp (Model Fogarty Clamp-01) Steel Surgical Instrument is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by

hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$14.00, falling to \$9.00 at 12 pieces — about 36% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36701790

Foster Sponge Forceps 21cm Locking Dressing Clamp Precision Surgical

\$6.20 per piece

Trade price from \$4.50
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Foster Sponge Forceps 21cm Locking Dressing Clamp Precision Surgical is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The 21cm size shown is our standard; other sizes in the same pattern are made to order.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$6.20, falling to \$4.50 at 25 pieces — about 27% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34986137

German Debakey Bulldog Clamp Set Atraumatic Curved Vascular Forceps

\$7.00 per piece

Trade price from \$3.99
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

German Debakey Bulldog Clamp Set Atraumatic Curved Vascular Forceps is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$7.00, falling to \$3.99 at 25 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Supplied as a set — the contents can be adjusted to your list. We work the

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34858395

German Stainless Steel Satinsky DeBakey Clamp 34

\$17.00 per piece

Trade price from \$11.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

German Stainless Steel Satinsky DeBakey Clamp 34 is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$17.00, falling to \$11.00 at 25 pieces — about 35% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34907469

Glover Atraumatic Bulldog Clamps 5.5

\$7.00 per piece

Trade price from \$3.99
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Glover Atraumatic Bulldog Clamps 5.5 is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$3.99 at 25 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000391

Glover Bulldog Clamp with Spring Adjusting Screw 9cm Curved 4cm Jaws

\$43.99 per piece Trade price from \$40.47 Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Glover Bulldog Clamp with Spring Adjusting Screw 3.5" (9cm) Size Curved 4cm DeBakey Atraumatic Serrated Jaws Premium Quality Stainless Steel Coronary Artery Cardiovascular Surgery Instruments.

The Glover bulldog clamp is a surgical instrument primarily used to temporarily occlude (stop) blood flow in vascular and cardiovascular surgeries. Glover DeBakey bulldog clamp designed to minimize trauma to the blood vessels while providing controlled pressure during procedures like coronary artery bypass grafting (CABG) and other cardiac procedures. This clamp used to temporarily blocks blood flow to a specific area, enabling surgeons to work on blood vessels or surrounding tissue without

interference, commonly used in cardiovascular surgeries. Glover DeBakey bulldog clamp equipped with 1×2 atraumatic serrated jaws, features a tension-adjusting mechanism that allows surgeons to regulate blood vessel pressure, preventing excessive clamping. The jaws, either straight or slightly curved, feature atraumatic serrations, enabling surgeons to safely clamp small and medium-sized vessels without damaging their walls. Made with premium quality stainless steel to ensure reliable and excellent performance

GLOVER BULLDOG CLAMP FEATURES



WI-34952003

Halsted Mosquito Artery Forceps Straight Curved Stainless Steel Hemostatic

\$1.50 per piece Trade price from \$1.20 Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Halsted Mosquito Artery Forceps Straight Curved Stainless Steel Hemostatic is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$1.50, falling to \$1.20 at 10 pieces — about 20% off the single-unit price. Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34851876

Halsted Mosquito Forceps 12.5cm Straight Curved Stainless Steel Hemostatic

\$1.50 per piece

Trade price from **\$1.20**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Halsted Mosquito Forceps 12.5cm Straight Curved Stainless Steel Hemostatic for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

The 12.5cm size shown is our standard; other sizes in the same pattern are made to order. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.50, falling to \$1.20 at 10 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000854

Harken Vena-Cava Clamp 9.5" Curved #1 Serrated Jaws

\$24.99 per piece

Trade price from **\$22.99**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Harken Vana-Cava Clamp 9.5" (24cm) Curved #1 Serrated Jaws Slightly Curved Shank with Ring Handle Premium Quality Stainless Steel Hemostatic Vascular Clamp Surgical Instruments.

Harken Vena Cava Clamp is a specialized 9.5-inch surgical instrument used in cardiovascular, thoracic, and transplant surgeries to temporarily occlude major blood vessels. Features a curved, finely serrated jaws allow for secure handling of delicate vessels without causing damage or crushing tissue. With finger ring ratcheted handle, the surgeons can lock the instrument in place and perform static blood vessel occlusions. You can reuse these Harken Vana-Cava Clamp over time as they will not bend or tarnish on

you. Made with High quality stainless steel, meet all international quality standards.

HARKEN VANA-CAVA CLAMP FEATURES:

- Crafted with 9.5" (24cm) Overall Length.
- Curved #1 Profile.
- Curved Jaws with Fine Serration.
- Slightly Curved Shank For Reaching Deep Vessels.
- Finger Ring Ratcheted Handle for Convenient Holding
- Specially Designed to Avoid Vascular Tissue Damage
- Ideal for Cardiovascular Procedures and Atraumatic Tissue Manipulation.
- Fully Autoclaved, Rust Free and Reusable





WI-46377602

Hartman Hemostat Forceps Curved Surgical Clamping Tool

\$3.90 per piece

Trade price from \$2.00
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hartman Hemostat Forceps Curved Surgical Clamping Tool is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.90, falling to \$2.00 at 50 pieces — about 49% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46421374

Hartman Hemostatic Forceps Surgical Hemostat Precise Artery Clamping Bleeding

\$4.00 per piece

Trade price from \$2.00
Minimum order 20 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hartman Hemostatic Forceps Surgical Hemostat Precise Artery Clamping Bleeding — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$2.00 at 20 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-45271703

Heavy Duty Steel Straight Forceps Tissue Grasping Clamp Surgical

\$8.00 per piece

Trade price from \$5.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Heavy Duty Steel Straight Forceps Tissue Grasping Clamp Surgical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a straight working end for direct, in-line access, and a heavy pattern built for load.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36752650

Heavy Pattern Rochester Ochsner Curved Clamp Precision Surgical

\$8.00 per piece

Trade price from \$5.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Heavy Pattern Rochester Ochsner Curved Clamp Precision Surgical is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a heavy pattern built for load.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36938012

Hemostatic Forceps Surgical Instrument

\$5.00 per piece

Trade price from \$4.20
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hemostatic Forceps Surgical Instrument is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$5.00, falling to \$4.20 at 2 pieces — about 16% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28116498

Hemostatic Mixer Artery Forceps Curved Stainless Steel Surgical Instruments

\$8.00 per piece

Trade price from \$4.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hemostatic Mixer Artery Forceps Curved Stainless Steel Surgical Instruments is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28135348

Hemostatic Mixer Artery Straight Stainless Steel Surgical Forceps

\$8.00 per piece

Trade price from \$4.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hemostatic Mixer Artery Straight Stainless Steel Surgical Forceps — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a straight working end for direct, in-line access.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000936

Ingrown Toenail Clippers & Cutter 4.5" (11.4cm) For Thick Nails

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Ingrown Toenail Clippers & Cutter 4.5" (11.4cm) Spring Action Handle Heavy Duty Thick Toe Nail Nipper German Grade Stainless Steel Dermatology Podiatry Instruments.

Toenail cutter clippers are specifically designed for working on hard, brittle, out of reach toenails. Ingrown toenail clippers are a high quality professional tools. The stainless-steel cutters are hand-filed for precision to make even the most awkward of jobs easy. They are perfectly aligned so that they close evenly, no matter the amount of pressure applied. Extremely sharp, hand-filed, pointed blades allow little pressure needed when cutting thick, hard and stubborn nails. High quality stainless steel used to manufacture this professional podiatrist tool which is

easy to clean and wipe down after use. An essential tool for adult personal hygiene and can also be used in nail salons or at home by seniors.

INGROWN TOENAIL CLIPPER CUTTER FEATURES

- 4.5" (11.4cm) Overall Length.
- Sharp Curved Blades.
- Single Spring Action Design.
- Adjustable Screw.
- Non-Slip Plier Grip Handle.
- Fully Autoclavable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Ingrown Toenail Clipper 4.5" (11.4cm) Spring Action Handle

TECHNICAL INFORMATION





WI-47080544

Javid Carotid Clamp Vascular Surgery Instrument for Secure Vessel

\$17.00 per piece

Trade price from \$9.00
Minimum order 30 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Javid Carotid Clamp Vascular Surgery Instrument for Secure Vessel is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$17.00, falling to \$9.00 at 30 pieces — about 47% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47073640

Javid Carotid Clamp Vascular Surgical Instrument for Precise Carotid

\$18.00 per piece

Trade price from \$11.90
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Javid Carotid Clamp Vascular Surgical Instrument for Precise Carotid is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$18.00, falling to \$11.90 at 10 pieces — about 34% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47091422

Javid Carotid Surgical Clamp Precision Vascular Instrument for Controlled

\$17.90 per piece

Trade price from \$11.00
Minimum order 12 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Javid Carotid Surgical Clamp Precision Vascular Instrument for Controlled, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$17.90, falling to \$11.00 at 12 pieces — about 39% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000887

Jones Towel Clamp 2.5"

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Jones Towel Clamp 2.5" (6.4cm) with Sharp Pointed Inwardly Curved Jaws and Cross-Action Spring Handle Premium Quality Stainless Steel Surgical Veterinary Instruments.

Jones Towel Clamps have sharp pointed tips which are curved and inward provides a secure and firm grip that doesn't let the object slip away. The tips overlap each other and are used to secure drapes and towels during surgical procedures. These instruments are very lightweight and are widely used in veterinary surgical procedures. The jones towel clamp forceps features a cross-action spring handle to open the jaws with comfort. Surgeons can loosen the pressure to slowly close the jaws and fixate the drapes without damaging the skin.

Manufactured in premium grade stainless steel to ensure the highest quality instruments in your operating rooms.

JONES TOWEL CLAMPS FEATURES

- Crafted with 2.5" (6.4cm) Overall Length.
- Inward Curved Jaws.
- 1x1 Prong.
- Overlapping Sharp Pointed Tips.
- Cross Action Spring Handle.
- Lightweight and Comfortable to Hold.
- Widely Used in All Sorts of Surgical Procedures.
- Fully Autoclavable and Reusable.
- Made With Premium Grade Stainless Steel.



**WI-70002921****Jones Towel Clamp 2"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Jones Towel Clamp 2" (5.1cm) with Sharp Pointed Inwardly Curved Jaws and Cross-Action Spring Handle Premium Quality Stainless Steel Surgical Veterinary Instruments.

Jones Towel Clamps have sharp pointed tips which are curved and inward provides a secure and firm grip that doesn't let the object slip away. The tips overlap each other and are used to secure drapes and towels during surgical procedures. These instruments are very lightweight and are widely used in veterinary surgical procedures. The jones towel clamp forceps features a cross-action spring handle to open the jaws with comfort. Surgeons can loosen the pressure to slowly close the jaws and fixate the drapes without damaging the skin.

Manufactured in premium grade stainless steel to ensure the highest quality instruments in your operating rooms.

JONES TOWEL CLAMPS FEATURES

- Crafted with 2" (5.1cm) Overall Length.
- Inward Curved Jaws.
- 1x1 Prong.
- Overlapping Sharp Pointed Tips.
- Cross Action Spring Handle.
- Lightweight and Comfortable to Hold.
- Widely Used in All Sorts of Surgical Procedures.
- Fully Autoclavable and Reusable.
- Made With Premium Grade Stainless Steel.

**WI-70000886****Jones Towel Clamp 3.5"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Jones Towel Clamp 3.5" (8.9cm) with Sharp Pointed Inwardly Curved Jaws and Cross-Action Spring Handle Premium Quality Stainless Steel Surgical Veterinary Instruments.

Jones Towel Clamps have sharp pointed tips which are curved and inward provides a secure and firm grip that doesn't let the object slip away. The tips overlap each other and are used to secure drapes and towels during surgical procedures. These instruments are very lightweight and are widely used in veterinary surgical procedures. The jones towel clamp forceps features a cross-action spring handle to open the jaws with comfort. Surgeons can loosen the pressure to slowly close the jaws and fixate the drapes without damaging the skin.

Manufactured in premium grade stainless steel to ensure the highest quality instruments in your operating rooms.

JONES TOWEL CLAMPS FEATURES

- Crafted with 3.5" (8.9cm) Overall Length.
- Inward Curved Jaws.
- 1x1 Prong.
- Overlapping Sharp Pointed Tips.
- Cross Action Spring Handle.
- Lightweight and Comfortable to Hold.
- Widely Used in All Sorts of Surgical Procedures.
- Fully Autoclavable and Reusable.
- Made With Premium Grade Stainless Steel.





WI-43456093

Kelly Artery Forceps Surgical Hemostat Clamp

\$1.00 per piece

Trade price from \$0.90
Minimum order 11 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kelly Artery Forceps Surgical Hemostat Clamp is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.00, falling to \$0.90 at 11 pieces — about 10% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36657806

Kelly Curved Clamp Stainless Steel Hemostatic Artery Forceps

\$1.30 per piece

Trade price from \$0.60
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kelly Curved Clamp Stainless Steel Hemostatic Artery Forceps is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.30, falling to \$0.60 at 25 pieces — about 54% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-36692761

Kelly Curved Hemostatic Forceps Medical Stainless Steel Artery

\$1.30 per piece

Trade price from \$0.60
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kelly Curved Hemostatic Forceps Medical Stainless Steel Artery is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$1.30, falling to \$0.60 at 25 pieces — about 54% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-27089727

Kelly Forceps Clamp Stainless Steel 3-Year Shelf

\$3.00 per piece

Trade price from \$1.90
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kelly Forceps Clamp Stainless Steel 3-Year Shelf is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$3.00, falling to \$1.90 at 1 piece — about 37% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.



WI-44570126

Kelly Forceps Stainless Steel Hemostatic Artery Forceps Surgical Instrument

\$1.30 per piece

Trade price from **\$0.60**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Kelly Forceps Stainless Steel Hemostatic Artery Forceps Surgical Instrument, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.
We state the grade we used on your quotation. Listed at \$1.30, falling to \$0.60 at 10 pieces — about 54% off the single-unit price. Larger quantities are priced on request.
Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27135104

Kelly Hemostatic Forceps 5.50" Stainless Steel Straight A+

\$3.00 per piece

Trade price from **\$1.90**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Kelly Hemostatic Forceps 5.50" Stainless Steel Straight A+, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a straight working end for direct, in-line access.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.00, falling to \$1.90 at 10 pieces — about 37% off the single-unit price. Larger quantities are priced on request.
Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27673962

Kelly Hemostatic Forceps 5.50" Straight Surgical Instruments

\$1.00 per piece

Trade price from \$0.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Kelly Hemostatic Forceps 5.50" Straight Surgical Instruments, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a straight working end for direct, in-line access.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$1.00, falling to \$0.50 at 1 piece — about 50% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44557721

Kelly Hemostatic Forceps Stainless Steel Surgical Artery

\$1.30 per piece

Trade price from \$0.60
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Kelly Hemostatic Forceps Stainless Steel Surgical Artery, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by

hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$1.30, falling to \$0.60 at 10 pieces — about 54% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-43448128

Kelly Hemostatic Forceps Straight & Curved Stainless Steel Surgical

\$1.00 per piece

Minimum order 11 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kelly Hemostatic Forceps Straight & Curved Stainless Steel Surgical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$1.00.

Minimum order 11 pieces, and the price moves with quantity. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-01168086

Kocher Ochsner Forceps Surgical Clamps Straight Curved

\$3.00 per piece

Trade price from **\$1.50**

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kocher Ochsner Forceps Surgical Clamps Straight Curved — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$3.00, falling to \$1.50 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35079891

Kocher Tissue Forceps-02 Stainless Steel Hemostatic

\$2.20 per piece

Trade price from \$1.55
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kocher Tissue Forceps-02 Stainless Steel Hemostatic is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.20, falling to \$1.55 at 10 pieces — about 30% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44579068

Lahey Artery Forceps Stainless Steel Hemostatic Forceps

\$3.20 per piece

Trade price from \$2.10
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Lahey Artery Forceps Stainless Steel Hemostatic Forceps, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$3.20, falling to \$2.10 at 10 pieces — about 34% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-44578087

Lahey Artery Forceps Stainless Steel Hemostatic Surgical Clamp

\$3.20 per piece

Trade price from \$2.10
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Lahey Artery Forceps Stainless Steel Hemostatic Surgical Clamp — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.20, falling to \$2.10 at 10 pieces — about 34% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-35106242

Lahey Artery Forceps Surgical Hemostatic Clamp

\$3.20 per piece

Trade price from \$1.90
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Lahey Artery Forceps Surgical Hemostatic Clamp for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$3.20, falling to \$1.90 at 10 pieces — about 41% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

**WI-70000958****Lorna Non-Perforating Towel Clamp 4.5" Inward Curved Jaws****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lorna Towel Clamp 4.5" (11.4cm) Curved Inward Serrated Jaws with Finger Ring Ratcheted Lock Premium Quality Stainless Steel Surgical Instrument.

Lorna Non-Perforating Towel Clamps are used for securing drapes or towels. Edna towel clamp is a specialized surgical tool that surgeons in multiple fields can use hold surgical drapes around the operative site, in order to enhance the sightlines and minimize the risk of infection. This towel clamp is suitable for many other electrosurgical applications, attaching suction and camera cords to drapes. This non-perforating instrument feature inward curved, serrated jaws. The fine serrated tips teeth are blunted, so the instrument gives a tight grip and does

not penetrate the surgical drape. The finger ring ratcheted handles assure firm holding. Made with premium quality stainless steel keeps it lightweight and rustproof and can be reused after sterilization.

LORNA NON PERFORATING TOWEL CLAMP FEATURES

- 4.5" (11.4cm) Overall Length.
- Wide Inward Curved Jaws.
- Fine Serrated Jaw Edges for Atraumatic Fixation.
- Non-perforating Blunt Teeth Which Not Penetrate the Surgical Drapes.
- Ratchet System for Self-retaining the Jaws.
- Ergonomic Finger Rings Ensuring Superior Handling.
- Reusable.

**WI-70000959****Lorna Non-Perforating Towel Clamp 4" Inward Curved Jaws****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lorna Towel Clamp 4" (10.2cm) Curved Inward Serrated Jaws with Finger Ring Ratcheted Lock Premium Quality Stainless Steel Surgical Instrument.

Lorna Non-Perforating Towel Clamps are used for securing drapes or towels. Edna towel clamp is a specialized surgical tool that surgeons in multiple fields can use hold surgical drapes around the operative site, in order to enhance the sightlines and minimize the risk of infection. This towel clamp is suitable for many other electrosurgical applications, attaching suction and camera cords to drapes. This non-perforating instrument feature inward curved, serrated jaws. The fine serrated tips teeth are blunted, so the instrument gives a tight grip and does not penetrate the surgical drape. The

finger ring ratcheted handles assure firm holding. Made with premium quality stainless steel keeps it lightweight and rustproof and can be reused after sterilization.

LORNA NON PERFORATING TOWEL CLAMP FEATURES

- 4" (10.2cm) Overall Length.
- Wide Inward Curved Jaws.
- Fine Serrated Jaw Edges for Atraumatic Fixation.
- Non-perforating Blunt Teeth Which Not Penetrate the Surgical Drapes.
- Ratchet System for Self-retaining the Jaws.
- Ergonomic Finger Rings Ensuring Superior Handling.
- Reusable.



**WI-70000957****Lorna Non-Perforating Towel Clamp 5.25\"/>****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lorna Towel Clamp 5.25" (13.3cm) Curved Inward Serrated Jaws with Finger Ring Ratcheted Lock Premium Quality Stainless Steel Surgical Instrument.

Lorna Non-Perforating Towel Clamps are used for securing drapes or towels. Edna towel clamp is a specialized surgical tool that surgeons in multiple fields can use hold surgical drapes around the operative site, in order to enhance the sightlines and minimize the risk of infection. This towel clamp is suitable for many other electrosurgical applications, attaching suction and camera cords to drapes. This non-perforating instrument feature inward curved, serrated jaws. The fine serrated tips teeth are blunted, so the instrument gives a tight grip and does

not penetrate the surgical drape. The finger ring ratcheted handles assure firm holding. Made with premium quality stainless steel keeps it lightweight and rustproof and can be reused after sterilization.

LORNA NON PERFORATING TOWEL CLAMP FEATURES

- 5.25" (13.3cm) Overall Length.
- Wide Inward Curved Jaws.
- Fine Serrated Jaw Edges for Atraumatic Fixation.
- Non-perforating Blunt Teeth Which Not Penetrate the Surgical Drapes.
- Ratchet System for Self-retaining the Jaws.
- Ergonomic Finger Rings Ensuring Superior Handling.
- Reusable.

**WI-70002920****Lorna Towel Clamp 3.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lorna Towel Clamp 3.5" (8.9cm) Curved Inward Serrated Jaws with Finger Ring Ratcheted Lock Premium Quality Stainless Steel Surgical Instrument.

Lorna Non-Perforating Towel Clamps are used for securing drapes or towels. Edna towel clamp is a specialized surgical tool that surgeons in multiple fields can use hold surgical drapes around the operative site, in order to enhance the sightlines and minimize the risk of infection. This towel clamp is suitable for many other electrosurgical applications, attaching suction and camera cords to drapes. This non-perforating instrument feature inward curved, serrated jaws. The fine serrated tips teeth are blunted, so the instrument gives a tight grip and does

not penetrate the surgical drape. The finger ring ratcheted handles assure firm holding. Made with premium quality stainless steel keeps it lightweight and rustproof and can be reused after sterilization.

LORNA NON PERFORATING TOWEL CLAMP FEATURES

- 3.5" (8.9cm) Overall Length.
- Wide Inward Curved Jaws.
- Fine Serrated Jaw Edges for Atraumatic Fixation.
- Non-perforating Blunt Teeth Which Not Penetrate the Surgical Drapes.
- Ratchet System for Self-retaining the Jaws.
- Ergonomic Finger Rings Ensuring Superior Handling.
- Reusable.





WI-46929399

Metal Angiotribe Forceps Hemostatic Surgical Instrument for Vessel Clamping

\$3.00 per piece

Trade price from \$1.50
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Metal Angiotribe Forceps Hemostatic Surgical Instrument for Vessel Clamping is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$3.00, falling to \$1.50 at 25 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46597023

Mikulicz Hemostat Forceps Model 03 Reusable Stainless Steel Surgical

\$4.00 per piece

Trade price from \$2.50
Minimum order 500 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Mikulicz Hemostat Forceps Model 03 Reusable Stainless Steel Surgical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$2.50 at 500 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46578175

Mikulicz Hemostatic Forceps Surgical Hemostat for Tissue Clamping Bleeding

\$9.00 per piece

Trade price from \$5.00
Minimum order 14 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Mikulicz Hemostatic Forceps Surgical Hemostat for Tissue Clamping Bleeding is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by

hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$9.00, falling to \$5.00 at 14 pieces — about 44% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46587072

Mikulicz Surgical Forceps Hemostatic Clamp for General Surgery Tissue

\$8.90 per piece

Trade price from \$4.90
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Mikulicz Surgical Forceps Hemostatic Clamp for General Surgery Tissue is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$8.90, falling to \$4.90 at 50 pieces — about 45% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Every piece is hand-finished and inspected before it leaves the bench:

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25889597

Mosquito Artery Forceps Stainless Steel Surgical

\$0.90 per piece

Trade price from **\$0.80**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Mosquito Artery Forceps Stainless Steel Surgical, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$0.90, falling to \$0.80 at 10 pieces — about 11% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25865740

Mosquito Forceps Powder Coated

\$1.00 per piece

Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Mosquito Forceps Powder Coated is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$1.00.

Minimum order 2 pieces, and the price moves with quantity. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36626810

Mosquito Forceps-19 Hemostatic Mosquito Artery Clamp

\$1.00 per piece

Trade price from \$0.60
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Mosquito Forceps-19 Hemostatic Mosquito Artery Clamp — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$1.00, falling to \$0.60 at 100 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35110212

Moynihan Artery Forceps Stainless Steel Surgical Hemostatic Instrument

\$3.20 per piece

Trade price from \$2.20
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Moynihan Artery Forceps Stainless Steel Surgical Hemostatic Instrument is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$3.20, falling to \$2.20 at 10 pieces — about 31% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34927695

Multi Color Rainbow Kelly Hemostat Forceps Instrument Set 5

\$3.50 per piece

Trade price from \$2.00

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Multi Color Rainbow Kelly Hemostat Forceps Instrument Set 5 — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.50, falling to \$2.00 at 1 piece — about 43% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34903836

Ochsner Kocher Hemostatic Forceps Curved 20cm Stainless Steel Surgical

\$1.50 per piece

Trade price from \$0.90

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Ochsner Kocher Hemostatic Forceps Curved 20cm Stainless Steel Surgical, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

The 20cm size shown is our standard; other sizes in the same pattern are

made to order. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$1.50, falling to \$0.90 at 10 pieces — about 40% off the single-unit price. Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000689

Payr Pylorus Clamps 11.5" with Pin 4.75" (120mm) Serrated Jaws

\$45.99 per piece

Trade price from \$42.31
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


PAYR Pylorus Clamps 11.5" (29.2cm) with Pin 4.75" (120mm) Long Longitudinal Serrated Angled Jaws with Special Hinge Handles Premium Quality Stainless Steel Intestinal Surgical Instruments.

Payr Intestinal Clamps is an incredible surgical instrument intended to use compressing or clamping bowel to remove during surgical procedures. Payr Pylorus Clamps commonly used during gut surgeries for secure closure of bowl loops and pylorus. These clamps often used in cardiothoracic surgeries as a valuable hemostat. Payr Pylorus clamp includes long sender blades having longitudinal serrations that minimize trauma to the tissues and pin allows a proper and secure closure of jaws with maximum force. Payr

Pylorus Intestinal Clamp is available in multiple lengths for various surgical needs. Made with premium quality stainless steel ensuring precision and durability. 100% customer satisfaction guaranteed.

PAYR PYLORUS INTESTINAL CLAMPS FEATURES

- Crafted With 11.5" (29.2cm) Overall Length.
- 4.75" (120mm) Long Sender Blades.
- Pin in Jaw for Secure Lock.
- Longitudinal Serration Blades to Prevent Soft Tissue Slipage and Damage.
- Special Hinge Mechanism for Secure Grip on Tissues.
- Smooth Without Lock Handle to Reduce User Fatigue.



WI-45169351

Pean Forceps Locking Hemostat Clamp

\$1.00 per piece

Trade price from \$0.70
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Pean Forceps Locking Hemostat Clamp — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.00, falling to \$0.70 at 25 pieces — about 30% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46148546

Pean Forceps Straight Curved Hemostatic Artery

\$0.80 per piece

Minimum order 500 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Pean Forceps Straight Curved Hemostatic Artery — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$0.80. Minimum order 500 pieces, and the price moves with quantity.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46137586

Pean Forceps Surgical Hemostat Clamp Reusable German

\$0.90 per piece

Minimum order 500 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Pean Forceps Surgical Hemostat Clamp Reusable German is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: full reusability — clean, autoclave and use again.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$0.90. Minimum order 500 pieces, and the price moves with quantity.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27670944

Pean Rochester Forceps 6.5" Straight Curved

\$1.00 per piece

Trade price from \$0.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Pean Rochester Forceps 6.5" Straight Curved, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.00, falling to \$0.50 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27276391

Pean Rochester Straight Curved Forceps Hemostat Surgical

\$2.55 per piece

Trade price from \$1.55
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Pean Rochester Straight Curved Forceps Hemostat Surgical for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.55, falling to \$1.55 at 1 piece — about 39% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003230****Professional Cantilever Nail Clipper 6" Black Coated Double Action****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Professional Cantilever Nail Clipper 6" Sharp Concave Jaws Double Action Black Color with Gold Barrel Spring German Grade Stainless Steel Chiropody Podiatry Tool for Thick and Tough Toenail Fungus Nails.

Besides the pedicurists podiatrists elderly care professionals, Seniors, home personal care, this heavy duty cantilever toe nail clipper is especially designed for the who are suffering from fungus infections. Regular toenail clippers are proven to be ineffective painful and cause even more damage. Fungus infections nails requires a right tool that can do the job and prevent injury. Ideal for cutting of problematic nails.

CANTILEVER TOE NAIL NIPPER FEATURES

PROFESSIONAL CANTILEVER NAIL CLIPPER – BLACK: Podiatrist double action mycotic nail nipper is the most durable pedicure nipper designed to work on extremely thick hard damaged and fungus infected toenails. The heavy duty cantilever nail cutter toenail clippers are ideal for both professionals and home use – for men, women and seniors. Effortlessly and safely cut thick toenails, ingrown nails and hangnails with the razor sharp edge and double barrel spring clippers.

PREMIUM QUALITY STAINLESS STEEL: Forged from O.R. Grade Stainless Steel, the Podiatrist Toenail Clippers for the best product life and performance. The curved razor sharp edge ensures a precise smooth cut through thick nails with ease. Corrosion resistant, tarnish free with a locking

**WI-70000857****Rainbow Cantilever Nail Clipper 6" Double Action****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty Rainbow Cantilever Nail Clipper 6" (15.2cm) Sharp Concave Jaws with Double Spring Action German Grade Stainless Steel Chiropody Podiatry Tool for Thick and Tough Toenail Fungus Nails.


Besides the pedicurists podiatrists elderly care professionals, Seniors, home personal care, this heavy duty multi titanium rainbow color cantilever toenail clipper is especially designed for the who are suffering from fungus infections. Regular toenail clippers are proven to be ineffective painful and cause even more damage. Fungus infections nails requires a right tool that can do the job and prevent injury. Ideal for cutting of problematic nails.

Podiatrist double action mycotic nail nipper is the most durable pedicure

nipper designed to work on extremely thick hard damaged and fungus infected toenails. Our rainbow cantilever toenail nippera are ideal for both professionals and home use – for men, women and seniors. Effortlessly and safely cut thick toenails, ingrown nails and hangnails with the razor sharp edge and double spring clippers.

RAINBOW CANTILEVER NAIL CLIPPER FEATURES

- Crafted with 6" (15.2cm) Overall Length.
- Multi Titanium Coated Rainbow Color.
- Sharp Concave Jaws.
- Double Spring Action.
- Locking Mechanism for Safety.
- Gripped Handle for Stronger Hold.





WI-34893156

Reusable Stainless Steel Mosquito Artery Forceps Straight Surgical Hemostatic

\$1.50 per piece

Trade price from \$1.20
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Reusable Stainless Steel Mosquito Artery Forceps Straight Surgical Hemostatic for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a straight working end for direct, in-line access, and full reusability — clean, autoclave and use again.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.50, falling to \$1.20 at 10 pieces — about 20% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44576095

Robert Artery Forceps Stainless Steel Hemostatic Forceps

\$2.50 per piece

Trade price from \$1.40
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Robert Artery Forceps Stainless Steel Hemostatic Forceps — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$2.50, falling to \$1.40 at 10 pieces — about 44% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44584009

Robert Artery Forceps Stainless Steel Hemostatic Surgical Clamp

\$2.50 per piece

Trade price from **\$1.40**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Robert Artery Forceps Stainless Steel Hemostatic Surgical Clamp is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.50, falling to \$1.40 at 10 pieces — about 44% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36705771

Rochester Ochsner Curved Forceps Heavy Duty Tissue Grasping Clamp

\$8.00 per piece

Trade price from **\$5.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Rochester Ochsner Curved Forceps Heavy Duty Tissue Grasping Clamp is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a heavy pattern built for load.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003663****Rubber Dam Clamp #204 Pedo Primary Tooth & Adult Premolar****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #204 Winged for Small Lower Molars, Pedo Primary Tooth & Adult Premolar Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber Dam Clamp #204 is specifically designed for holding a rubber dam on small lower molars, children's primary molars, and adult premolars. This winged clamp is utilized in dental procedures, such as root canals, to effectively isolate the tooth and maintain a clean working area free from saliva and blood. We offer a wide range of rubber dam clamps suitable for anterior, premolar, and molar applications, made from tempered stainless steel for superior performance, durability, and corrosion resistance.

RUBBER DAM CLAMPS #204 FEATURES:

- Used for Small Lower Molars.
- Winged Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH:** Pedo Primary Tooth, Adult Premolar**SIZE:** 204**REUSABLE/AUTOCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:** Passivated**WI-70003659****Rubber Dam Clamp #212R for Right Anterior****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #212R Wingless for Right Anterior Teeth Premium Quality Stainless Steel Dental Endodontic Instruments.

Endodontic Rubber Dam Clamp 212R is a specialized instrument for anterior teeth and premolars, providing retraction and access for labial caries. It secures a rubber dam, ensuring the treatment area remains dry and visible. The "R" signifies it's for the right side and is a wingless clamp designed to securely hold a dental dam around the tooth to create a dry and clean field, preventing interference from saliva and blood. We offer a wide range of rubber dam clamps suitable for anterior, premolar, and molar applications, made from tempered stainless steel for

superior performance, durability, and corrosion resistance.

RUBBER DAM CLAMP #212R FEATURES:

- Used for Right Anterior.
- Springs Back to Original Position.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH:** Anterior Teeth and Premolars**SIZE:** 212R**REUSABLE/AUTOCLAVABLE:** Yes**STERILE:** Non-Sterile

**WI-70003658****Rubber Dam Clamp #24 for Upper Right & Lower Left Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #24 Wingless Flat Jaws for Upper Right and Lower Left Molars Premium Quality Stainless Steel Dental Endodontic Instruments.

Endodontic rubber dam clamp #24 stabilizes a rubber dam on upper right and lower left molars. With flat jaws, clamp #24 is commonly used in endodontic procedures such as root canals and general dental restorations to isolate teeth and maintain a dry, stable operating field. Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMPS #24 FEATURES:

- Used for Upper Right and Lower Left Molars
- Wingless Smooth.
- Flat Jaws.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** For Molars**SIZE:** 24**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**WI-70003657****Rubber Dam Clamp #25 for Upper Left & Lower Right Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #25 Wingless Flat Jaws for Upper Left & Lower Right Molars Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber Dam Clamp #25 is utilized in dentistry to isolate upper left and lower right molars, particularly during endodontic treatments and restorations. Its flat jaws ensure a stable grip on the dental dam, enhancing visibility, maintaining dryness, and safeguarding the patient from the risk of aspirating instruments or debris. Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMPS #25 FEATURES:

- Used for Upper Left & Lower Right Molars.
- Wingless Smooth.
- Flat Jaws.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** For Molars**SIZE:** 25**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile

**WI-70003656****Rubber Dam Clamp #26 for Small Lower Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp 26 Wingless Slightly Reversed Festooned Jaws for Small Lower Molars Premium Quality Stainless Steel Dental Endodontic Instruments

Rubber Dam Clamp #26 is designed to secure a rubber dam to small lower molars, with variations for upper and pediatric molars. It isolates the tooth during restorative, endodontic, and prosthodontic procedures. Clamp #26 featuring beaks that snugly fit the contours of small lower molars to ensure secure retention without excessive pressure. Our stainless steel clamps are crafted with precise tension to give each clamp the flexibility to stretch around specified teeth and the strength to retain their shape while holding down rubber dam.

RUBBER DAM CLAMP #26 FEATURES:

- For Small Lower Molars
- Slightly Reversed Festooned Jaws.
- Wingless, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**TEETH** Small Lower Molars**SIZE:** 26**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:****WI-70000309****Rubber Dam Clamp #W14A for Partially Erupted or Irregularly Shaped Molars****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp #W14A Wingless for Partially Erupted or Irregularly Shaped Molars Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber dam clamp W14A isolates teeth, especially partially erupted or irregularly shaped molars, in dental procedures. Its wingless design features deeply festooned jaws for a secure grip, making it suitable for endodontic, fixed prosthetic, and general restorative treatments. We offer a wide range of rubber dam clamps suitable for anterior, premolar, and molar applications, made from tempered stainless steel for superior performance, durability, and corrosion resistance.

RUBBER DAM CLAMP #W14A FEATURES:

- Used for Partially Erupted or Irregularly Shaped Molars.
- Deeply Festooned Jaws.
- Wingless, Smooth.
- Tempered Stainless Steel.
- Autoclavable.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Endodontic Instruments**APPLICATION:** Restorative, Cosmetic, or Endodontic Procedures**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes



WI-24277762

Satinsky Vascular Clamp Surgical 9" Forceps DeBakey

\$17.50 per piece

Trade price from \$8.85
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Satinsky Vascular Clamp Surgical 9" Forceps DeBakey is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$17.50, falling to \$8.85 at 10 pieces — about 49% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29256369

Seller Roeder Towel Clamp Ball Jaws 11"

\$4.00 per piece

Trade price from \$1.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Seller Roeder Towel Clamp Ball Jaws 11" is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.00, falling to \$1.50 at 1 piece — about 62% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-06914644

Single Use Stainless Steel Artery Mosquito Forceps Surgical

\$2.50 per piece

Trade price from **\$1.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Single Use Stainless Steel Artery Mosquito Forceps Surgical, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.50, falling to \$1.00 at 10 pieces — about 60% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34984273

Spencer Wells Artery Forceps 16cm Stainless Steel Surgical Hemostatic

\$1.70 per piece

Trade price from **\$0.90**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Spencer Wells Artery Forceps 16cm Stainless Steel Surgical Hemostatic is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The 16cm size shown is our standard; other sizes in the same pattern are made to order.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$1.70, falling to \$0.90 at 10 pieces — about 47% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34904633

Spencer Wells Artery Forceps Straight 12.5cm
Stainless Steel Hemostatic

\$2.50 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Spencer Wells Artery Forceps Straight 12.5cm Stainless Steel Hemostatic, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a straight working end for direct, in-line access.

The 12.5cm size shown is our standard; other sizes in the same pattern are

made to order. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.50, falling to \$1.50 at 10 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34982225

Spencer Wells Artery Hemostatic Forceps
Stainless Steel Surgical Instrument

\$2.50 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Spencer Wells Artery Hemostatic Forceps Stainless Steel Surgical Instrument, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.50, falling to \$1.50 at 10 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.



WI-35047678

Spencer Wells Forceps 5 Inch Straight
Stainless Steel Hemostatic

\$1.80 per piece

Trade price from **\$0.90**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Spencer Wells Forceps 5 Inch Straight Stainless Steel Hemostatic for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a straight working end for direct, in-line access.

The 5" size shown is our standard; other sizes in the same pattern are made to order. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$1.80, falling to \$0.90 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34935540

Spencer Wells Forceps 8 Inch Curved Stainless
Steel Hemostatic

\$2.50 per piece

Trade price from **\$1.50**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Spencer Wells Forceps 8 Inch Curved Stainless Steel Hemostatic for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

The 8" size shown is our standard; other sizes in the same pattern are made to order. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.50, falling to \$1.50 at 10 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44544624

Spencer Wells Hemostatic Forceps Stainless Steel Surgical

\$1.80 per piece

Trade price from \$0.80
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Spencer Wells Hemostatic Forceps Stainless Steel Surgical is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$1.80, falling to \$0.80 at 10 pieces — about 56% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44553598

Spencer Wells Stainless Steel Hemostatic Artery Forceps Surgical Instrument

\$1.80 per piece

Trade price from \$0.70
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Spencer Wells Stainless Steel Hemostatic Artery Forceps Surgical Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.80, falling to \$0.70 at 10 pieces — about 61% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35059671

Sponge Forceps Stainless Steel Surgical Dressing Clamp

\$2.90 per piece

Trade price from \$1.70
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Sponge Forceps Stainless Steel Surgical Dressing Clamp is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.90, falling to \$1.70 at 10 pieces — about 41% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46164525

Stainless Pean Hemostatic Forceps Surgical Artery Clamp

\$1.00 per piece

Trade price from \$0.90
Minimum order 500 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Pean Hemostatic Forceps Surgical Artery Clamp is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.00, falling to \$0.90 at 500 pieces — about 10% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-37714333

Stainless Steel 4-Piece Hair Extension Pliers

\$20.00 per piece

Trade price from \$12.00
Minimum order 20 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel 4-Piece Hair Extension Pliers is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$20.00, falling to \$12.00 at 20 sets — about 40% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35098433

Stainless Steel Babcock Surgical Forceps Tissue Holding Clamp

\$2.20 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Babcock Surgical Forceps Tissue Holding Clamp — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$2.20, falling to \$1.50 at 10 pieces — about 32% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34990043

Stainless Steel Babcock Tissue Forceps Atraumatic Surgical Clamp Set

\$4.00 per piece

Trade price from **\$2.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Babcock Tissue Forceps Atraumatic Surgical Clamp Set for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$2.00 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29202455

Stainless Steel Backhaus Towel Clamp Forceps

\$4.00 per piece

Trade price from **\$1.50**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Backhaus Towel Clamp Forceps, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$1.50 at 1 piece — about 62% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-26092304

Stainless Steel Benguela Hemostatic Forceps Kocher Straight/Curved

\$3.00 per piece

Trade price from **\$1.50**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Benguela Hemostatic Forceps Kocher Straight/Curved, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.00, falling to \$1.50 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34926496

Stainless Steel Cheatle Forceps Surgical Clamp

\$3.00 per piece

Trade price from **\$1.80**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Cheatle Forceps Surgical Clamp — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.00, falling to \$1.80 at 10 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36678845

Stainless Steel Crile Hemostat Forceps Surgical Locking

\$0.80 per piece

Trade price from **\$0.70**
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Crile Hemostat Forceps Surgical Locking — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$0.80, falling to \$0.70 at 100 pieces — about 13% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34917508

Stainless Steel Curved Mosquito Forceps Locking Hemostatic Surgical Instrument

\$2.00 per piece

Trade price from **\$0.90**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Curved Mosquito Forceps Locking Hemostatic Surgical Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.00, falling to \$0.90 at 10 pieces — about 55% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34940385

Stainless Steel Diethrich Bulldog Clamps Use Surgery

\$7.00 per piece

Trade price from \$3.99
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Diethrich Bulldog Clamps Use Surgery is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$7.00, falling to \$3.99 at 25 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-37726309

Stainless Steel Hair Extension Pliers 4-Piece

\$15.00 per piece

Trade price from \$12.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Hair Extension Pliers 4-Piece for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$15.00, falling to \$12.00 at 1 set — about 20% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.



WI-39821326

Stainless Steel Hair Extension Pliers Kit

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Hair Extension Pliers Kit, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-39815363

Stainless Steel Hair Extension Pliers Kit 2 Clamp

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Hair Extension Pliers Kit 2 Clamp for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-39816352

Stainless Steel Hair Extension Pliers Kit with 2

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Hair Extension Pliers Kit with 2 is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46420370

Stainless Steel Hartman Hemostatic Surgical Forceps Precision Artery Clamp

\$3.50 per piece

Trade price from \$2.50
Minimum order 500 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Hartman Hemostatic Surgical Forceps Precision Artery Clamp — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$3.50, falling to \$2.50 at 500 pieces — about 29% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-01169093

Stainless Steel Hemostat Kocher Forceps 7"

\$3.00 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Hemostat Kocher Forceps 7" is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by

hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$3.00, falling to \$1.50 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35076339

Stainless Steel Hemostatic Forceps Surgical Instrument

\$3.50 per piece

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Hemostatic Forceps Surgical Instrument for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for

corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$3.50. Minimum order 1 piece, and the price moves with quantity.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-33358800

Stainless Steel Kelly Forceps

\$3.00 per piece

Trade price from **\$1.90**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Kelly Forceps is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.00, falling to \$1.90 at 25 pieces — about 37% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-06206494

Stainless Steel Mosquito Forceps Model 1

\$2.50 per piece

Trade price from **\$1.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Mosquito Forceps Model 1 — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.50, falling to \$1.00 at 10 pieces — about 60% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-06247465

Stainless Steel Mosquito Forceps Powder Coated

\$2.50 per piece

Trade price from **\$1.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Mosquito Forceps Powder Coated — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.50, falling to \$1.00 at 10 pieces — about 60% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-06958100

Stainless Steel Mosquito Hemostat Locking Forceps

\$1.50 per piece

Trade price from **\$0.70**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Mosquito Hemostat Locking Forceps is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$1.50, falling to \$0.70 at 10 pieces — about 53% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35027754

Stainless Steel Moynihan Artery Forceps
Surgical Hemostatic Clamp Medical

\$3.20 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Moynihan Artery Forceps Surgical Hemostatic Clamp Medical, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.20, falling to \$2.00 at 10 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34971373

Stainless Steel Spencer Wells Artery Forceps
Curved Surgical

\$1.70 per piece

Trade price from \$0.90
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Spencer Wells Artery Forceps Curved Surgical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.70, falling to \$0.90 at 10 pieces — about 47% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35014734

Stainless Steel Sponge Forceps Surgical Holding Clamp Reusable Medical

\$2.90 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Sponge Forceps Surgical Holding Clamp Reusable Medical is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.90, falling to \$2.00 at 10 pieces — about 31% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34958310

Stainless Steel Surgical Bulldog Clamps & Clips-Precision

\$7.00 per piece

Trade price from \$3.99
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Bulldog Clamps & Clips-Precision is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by

hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$7.00, falling to \$3.99 at 25 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-00889157

Stainless Steel Surgical Instrument Set

\$2.00 per piece

Trade price from \$0.90
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Instrument Set is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.00, falling to \$0.90 at 10 pieces — about 55% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34970224

Stainless Steel Surgical Micro Straight Bulldog Clamps

\$7.00 per piece

Trade price from \$3.99
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Surgical Micro Straight Bulldog Clamps for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a straight working end for direct, in-line access, and a micro-scale working end for fine and delicate work.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$7.00, falling to \$3.99 at 25 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36646844

Stainless Steel Surgical Towel Clip Clamp

\$0.90 per piece

Trade price from \$0.80
Minimum order 200 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Surgical Towel Clip Clamp for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$0.90, falling to \$0.80 at 200 pieces — about 11% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44539597

Stainless Steel Surgical Towel Clip Clamp Forceps

\$1.50 per piece

Trade price from \$0.70
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Stainless Steel Surgical Towel Clip Clamp Forceps is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$1.50, falling to \$0.70 at 10 pieces — about 53% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29303341

Stainless Steel Towel Clamp

\$4.00 per piece

Trade price from \$1.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Towel Clamp, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.00, falling to \$1.50 at 1 piece — about 62% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36649828

Stainless Steel Towel Clip Surgical Draping Clamp Reusable Medical

\$1.00 per piece

Trade price from \$0.80
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Towel Clip Surgical Draping Clamp Reusable Medical, made for general and specialty operating theatre use. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$1.00, falling to \$0.80 at 100 pieces — about 20% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-33067066

Stainless Steel Weft Extensions Tape Sealing

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Weft Extensions Tape Sealing — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28123441

Steel Curved Hemostatic Surgical Instruments

\$8.00 per piece

Trade price from \$4.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Steel Curved Hemostatic Surgical Instruments is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a curved working end for reach around structures.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003312****Stokes Rubber Dam Clamp 7" with Lock Handle****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Stokes Rubber Dam Clamp Forceps 7" (17.8cm) Small Projection Curved Ends with Slide Ring Lock Mechanism Premium Quality Stainless Steel Dental Endodontic Instruments.

Rubber Dam Clamp Stokes Forceps are dental instruments used for placing and removing rubber dam clamps to isolate the working area during restorative and endodontic procedures. Rubber Dam Stokes forcep feature a sliding ring locking mechanism that holds the forceps at any desired position, allowing for secure and accurate clamp placement. The small beak projections at the curved ends engage securely into the clamp holes for a firm, non-slip grip allow easier clamp manipulation. Made from premium stainless steel, these forceps are fully autoclavable for

maintaining hygiene and safety standards.

STOKES RUBBER DAM CLAMP FORCEPS FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- Design for Efficient and Secure Rubber Dam Placement & Removal.
- Smooth Spring Makes the Placement of Clamps Easy.
- Curved Tips with Small Projections Fit Into Holes on the Rubber Dam Clamp.
- Sliding Ring Locks Handles in Position During Tooth Placement.
- Fully Autoclaved / Reusable.
- Made of Surgical Stainless Steel.

PACKAGE INCLUDES**WI-34855926****Straight Curved Mosquito Hemostatic Forceps Stainless Steel Artery Pean****\$1.50** per pieceTrade price from **\$1.20**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Straight Curved Mosquito Hemostatic Forceps Stainless Steel Artery Pean — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.50, falling to \$1.20 at 10 pieces — about 20% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35043709

Straight Spencer Wells Hemostatic Forceps
Stainless Steel Surgical Instrument

\$1.80 per piece

Trade price from \$0.90
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Straight Spencer Wells Hemostatic Forceps Stainless Steel Surgical Instrument for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a straight working end for direct, in-line access.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.80, falling to \$0.90 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36822994

Surgical Hemostatic Forceps With Blue Clip

\$17.00 per piece

Trade price from \$8.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Hemostatic Forceps With Blue Clip is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$17.00, falling to \$8.00 at 1 piece — about 53% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

We work the instrument by hand at the finishing stage — box locks, hinges and

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34890616

Surgical Instruments Micro Straight Bulldog Clamps Micro

\$7.00 per piece

Trade price from \$3.99
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Instruments Micro Straight Bulldog Clamps Micro is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a straight working end for direct, in-line access, and a micro-scale working end for fine and delicate work.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$7.00, falling to \$3.99 at 25 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-29308399

Surgical Instruments Set for Micro

\$15.00 per piece

Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Instruments Set for Micro is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a micro-scale working end for fine and delicate work.

Supplied as a set — the contents can be adjusted to your list. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$15.00.

Minimum order 10 sets, and the price moves with quantity. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



WI-01174096

Surgical Instruments Stainless Steel

\$3.00 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Instruments Stainless Steel is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$3.00, falling to \$1.50 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-45181285

Surgical Pean Hemostatic Forceps Straight Curved Artery Clamp

\$1.00 per piece

Trade price from \$0.70
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Pean Hemostatic Forceps Straight Curved Artery Clamp — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$1.00, falling to \$0.70 at 25 pieces — about 30% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-45269768

Surgical Rochester Ochsner Straight Clamp Precision Tissue Holding Forceps

\$8.00 per piece

Trade price from \$5.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Surgical Rochester Ochsner Straight Clamp Precision Tissue Holding Forceps for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: a straight working end for direct, in-line access.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-06214445

Surgical Stainless Steel Metal Artery Forceps

\$2.50 per piece

Trade price from \$1.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Stainless Steel Metal Artery Forceps is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.50, falling to \$1.50 at 1 piece — about 40% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Every piece is hand-finished and inspected before it leaves the bench:

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34945610

TC Babcock Forceps Stainless Steel Surgical Tissue

\$4.00 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

TC Babcock Forceps Stainless Steel Surgical Tissue is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. The name tells you what you are getting: tungsten carbide inserts for a longer cutting life.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$4.00, falling to \$2.00 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-36660872

Towel Clamp Stainless Steel Surgical Draping Clip

\$1.00 per piece

Trade price from \$0.80
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Towel Clamp Stainless Steel Surgical Draping Clip is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place.

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$1.00, falling to \$0.80 at 100 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Every piece is hand-finished and inspected before it leaves the bench:



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-46234797

Vascular Surgery Set-07 Vascular Surgical Kit Clamp Precision Artery

\$220.00 per piece

Trade price from \$200.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Vascular Surgery Set-07 Vascular Surgical Kit Clamp Precision Artery is a clamp we forge, grind and finish in our own Sialkot workshop. Instruments of this type occlude a vessel or duct, or hold tissue and drapes securely in place. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$220.00, falling to \$200.00 at 1 set — about 9% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



WI-70000715

Wertheim Atraumatic Parametrium Clamp 10" Straight 2x3 Teeth Jaws

\$14.99 per piece

Trade price from \$13.79
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Wertheim Atraumatic Parametrium Clamp 10" (25.4cm) Straight 2x3 Teeth Debakey Serration Premium Quality Hysterectomy and Parametric Gynecology and Obstetric Surgical Instruments

Wertheim parametrium forceps is a specialized device that gynecology's and obstetric surgeons frequently use to manipulate uterus ligaments and blood vessels, in order to permit their dissection. This forceps can also used as a hemostat to control bleeding blood vessels during general procedures as well. Parametrium clamp comes with an atraumatic toothting (Debakey style 2x3 jaws) for superior grip on parametrial tissues, ensuring secure manipulation and dissection during surgical procedures.

Wertheim Atraumatic Parametrium Clamp features a pair atraumatic 2x3 fine serration jaws with elongated shafts for deep surgical reaching conform to the anatomy of the pelvis. Also, the box lock system helps to face the jaws against each other with smooth joint movements. Wertheim parametrium forceps manufactured with premium grade stainless steel to ensure the highest quality and product performance during surgical procedures.

WERTHEIM ATRAUMATIC PARAMETRIUM CLAMP FEATURES

- Crafted with 10" (25.4cm) Overall Length.
- Debakey Atraumatic 2x3 Fine Serration.



**WI-70000714****Wertheim Atraumatic Parametrium Clamp
9.75" Curved 2×3 Teeth Jaws****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Wertheim Atraumatic Parametrium Clamp 9.75" (24.8cm) Straight 2×3 Teeth Debaquey Serration Premium Quality Hysterectomy and Parametric Gynecology and Obstetric Surgical Instruments

Wertheim parametrium forceps is a specialized device that gynecology's and obstetric surgeons frequently use to manipulate uterus ligaments and blood vessels, in order to permit their dissection. This forceps can also used as a hemostat to control bleeding blood vessels during general procedures as well. Parametrium clamp comes with an atraumatic toothing (Debaquey style 2×3 jaws) for superior grip on parametrial tissues, ensuring secure manipulation and dissection during surgical procedures.

Wertheim Atraumatic Parametrium Clamp features a pair atraumatic 2×3 fine serration jaws with elongated shafts for deep surgical reaching conform to the anatomy of the pelvis. Also, the box lock system helps to face the jaws against each other with smooth joint movements. Wertheim parametrium forceps manufactured with premium grade stainless steel to ensure the highest quality and product performance during surgical procedures.

**WERTHEIM ATRAUMATIC
PARAMETRIUM CLAMP
FEATURES**

- Crafted with 9.75" (24.8cm) Overall Length.
- Debaquey Atraumatic 2×3 Fine Serration.





WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

Request a quotation

Send us product codes from this catalogue, quantities and your country.

We reply the same working day on WhatsApp.

Online shop	shop.waqarintl.com
Email	sales@waqarintl.com
WhatsApp /	
Phone	+92 307 6173627
Works	Roras Road, Opposite Anwar Khawaja Industries, Sialkot 51310, Pakistan

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.