




WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

Elevators & Currettes Surgical Instruments

Reusable stainless-steel elevators & currettes forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



123

PRODUCTS

\$4.00

PRICES FROM

OEM

PRIVATE LABEL

Global

WORLDWIDE
EXPORT

Contents

123 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

Product	Code	Price	Page
10 Gracey Curettes 11/12 Blue Titanium Coated Tips Hollow Handle	WI-70002077	\$59.99	5
10 Gracey Curettes 13/14 Gold Coated Double Ended	WI-70002072	\$59.99	5
10 Periosteal Dental Elevator Molt No. 9	WI-70003852	\$39.99	6
11.34 Inch Stainless Steel Rectus Curette No	WI-32397761	\$5.00	6
12 Lucas Bone Curette #84 Small 2.4mm Spoon Double Ended	WI-70003879	\$59.99	7
12 Lucas Curette #86 Dentist Lab Dental	WI-70003877	\$71.46	7
12 Molt Periosteal M9 Dental Surgery Tooth Elevator	WI-70003876	\$44.99	8
12 Sayre Periosteal Elevator 6.75" Double Ended	WI-70003875	\$131.99	8
2 Gracey Curettes 13/14 Gold Coated Double Ended	WI-70002075	\$13.99	9
3 Gracey Curettes 11/12 Blue Titanium Coated Tips Hollow Handle	WI-70002079	\$19.99	9
3 Gracey Curettes 13/14 Gold Coated Double Ended	WI-70002074	\$19.99	10
3 Lucas Bone Curette #84 Small 2.4mm Spoon Double Ended	WI-70003871	\$19.99	10
3 Novak Endometrial Suction Biopsy Curettes 9.75" 2, 3 & 4mm Dia	WI-70002395	\$24.99	11
3 Sayre Periosteal Elevators 6.75" Double Ended	WI-70003867	\$38.99	11
5 Buck Ear Curettes 6.5" Straight Sharp Size #00, #0, #1, #2, #3	WI-70000466	\$34.99	12
5 Gracey Curettes 11/12 Blue Titanium Coated Tips Hollow Handle	WI-70002078	\$29.99	12
5 Gracey Curettes 13/14 Gold Coated Double Ended	WI-70002073	\$29.99	13
5 Kevorkian Younge Biopsy Curette 12" with Basket Loop Size 3x12mm	WI-70002390	\$49.99	13
5 Kevorkian-Younge Biopsy Curette 12" Loop 3X12mm with Basket Half Gold Handle	WI-70002248	\$49.99	14
5 Kevorkian-Younge Biopsy Curette 12" Loop 3X12mm Without Basket Half Gold Handle	WI-70002249	\$49.99	14
5 Offset Sinus Osteotomes with Sterilization Cassette 7.25" X 5"	WI-70003124	\$69.99	15
5 Pc Molt Periosteal M9 Dental Elevator	WI-70003865	\$27.49	15
6 Lucas Bone Curette #84 Small 2.4mm Spoon Double Ended	WI-70003861	\$35.99	16
6 Molt Periosteal M9 Dental Surgery Tooth Elevator	WI-70003858	\$29.99	16
6 Sayre Periosteal Elevator 6.75" Double Ended	WI-70003857	\$71.99	17
Adson Periosteal Elevator 7.5" Semi Sharp Curved 10mm Wide	WI-70000611	\$24.99	17
Adson Periosteal Elevator 7.5" Semi Sharp Curved 7mm Wide	WI-70000612	\$24.99	18
Alexander Chisel	WI-70003278	\$33.09	18
Billeau Ear Loop Curette Size #1 – Small	WI-70000178	\$9.99	19
Billeau Ear Loop Curette Size #2 – Medium	WI-70000177	\$9.99	19
Billeau Ear Loop Curette Size #3 – Large	WI-70000176	\$9.99	20
Cage Curette Surgical	WI-46975459	\$15.00	20
Cobb Osteotome	WI-70003276	\$38.19	21
Cone Ring Curette 9 Inch Angled	WI-35248110	\$8.00	21
Cone Ring Curette 9" Straight 6mm Size #2	WI-70000627	\$12.99	22
Coupland Root Elevators 2 0 mm	WI-70003254	\$9.99	22
Coupland Root Elevators 2.5 mm	WI-70003253	\$9.99	23
Coupland Root Elevators 3 0 mm	WI-70003252	\$9.99	23
Curette With Cage Reusable Surgical Tool for Precise Curettage	WI-46993252	\$14.80	24
Curette with Cage Surgical Kit	WI-47001203	\$14.90	24
Dental Columbia Universal Curette 2R/2L	WI-70003227	\$7.99	25

Contents (continued)

Product	Code	Price	Page
Dental Columbia Universal Curette Scalars 13/14	WI-70003228	\$7.99	26
Dental Elevators Autoclave Sterilization Cassette Tray 7" x 8" for 8 Instruments	WI-70000441	\$32.99	26
Dental Gracey Curettes 11/12 Double Ended	WI-70003623	\$9.99	27
Double Ended Micro Curettes Medical	WI-35495799	\$4.00	27
Double Ended Micro Curettes Precision Surgical Tools	WI-35506825	\$4.00	28
Double Ended Micro Curettes Stainless Steel Surgical Instrument Set	WI-35538762	\$4.00	28
Fox Dermal Curette 1mm Round Tip – 5.5" (14cm) Straight	WI-70002649	\$8.99	29
Fox Dermal Curette 2mm Round Tip – 5.5" (14cm) Straight	WI-70002648	\$8.99	29
Fox Dermal Curette 3mm Round Tip – 5.5" (14cm) Straight	WI-70002647	\$8.99	30
Fox Dermal Curette 4mm Round Tip – 5.5" (14cm) Straight	WI-70002646	\$8.99	30
Fox Dermal Curette 5mm Round Tip – 5.5" (14cm) Straight	WI-70002645	\$8.99	31
Fox Dermal Curette 6mm Round Tip – 5.5" (14cm) Straight	WI-70002644	\$8.99	31
Gracey Curette 1/2 Blue Tips Double Ended	WI-70000306	\$10.99	32
Gracey Curette 1/2 Gold Tips Double Ended	WI-70000319	\$10.99	32
Gracey Curette 1/2 Titanium Coated Double Ended	WI-70002090	\$10.99	33
Gracey Curette 11/12 Blue Titanium Coated Tips Double Ended	WI-70002081	\$10.99	33
Gracey Curette 13/14 Gold Coated Double Ended	WI-70002076	\$10.99	34
Gracey Curette 3/4 Blue Tips Double Ended	WI-70000307	\$10.99	34
Gracey Curette 5/6 Blue Tips Double Ended	WI-70000295	\$10.99	35
Gracey Curette 7/8 Blue Tips Double Ended	WI-70000294	\$10.99	35
Gracey Curette 9/10 Blue Coated Tips Double Ended	WI-70000305	\$10.99	36
Gracey Curette 9/10 Gold Coated Double Ended	WI-70002085	\$10.99	36
Heaney Uterine Biopsy Curette 9" Serrated Loop	WI-70002394	\$10.99	37
Kevorkian-Younger Biopsy Curette 12" Loop 3X12mm with Basket Half Gold Handle	WI-70002252	\$12.99	37
Kevorkian-Younger Biopsy Curette 12" Loop 3X12mm Without Basket Half Gold Handle	WI-70002251	\$12.99	38
Kevorkian-Younger Endocervical Biopsy Curette 12" Loop 3X12mm with Basket	WI-70002393	\$11.99	38
Kevorkian-Younger Endocervical Biopsy Curette 12" Loop 3X12mm Without Basket	WI-70002392	\$11.99	39
Lambotte Osteotome 9" Straight 13mm Tip	WI-70001915	\$20.99	39
Lambotte Osteotome 9" Straight 19mm Tip	WI-70001914	\$21.99	40
Lambotte Osteotome 9" Straight 25mm Tip	WI-70001913	\$22.99	40
Lambotte Osteotome 9" Straight 32mm Tip	WI-70001912	\$23.99	41
Lambotte Osteotome 9" Straight 38mm Tip	WI-70000202	\$24.99	41
Lambotte Osteotome 9" Straight 6mm Tip	WI-70001916	\$19.99	42
Lexer Mini Chisel	WI-70003273	\$33.59	42
Lucas Bone Curette #84 Small 2.4mm Spoon Double Ended	WI-70003838	\$7.99	43
Lucas Bone Curette #85 Double Ended 2.7mm Smooth Spoon Tips	WI-70003837	\$7.99	43
Lucas Bone Curette #85 Double Ended 2.7mm Tips Hexagonal Handle	WI-70000609	\$7.99	44
Lucas Bone Curette #86 Double Ended 3mm Smooth Spoon Tips	WI-70003836	\$7.99	44
Lucas Bone Curette #86 Double Ended 3mm Tips Hexagonal Handle	WI-70000608	\$7.99	45
Lucas Bone Curette #86 Double Ended Spoon Tips – Hollow	WI-70000072	\$9.99	45
Lucas Bone Curette #87 Double Ended 3.8mm Smooth Spoon Tips	WI-70003840	\$7.99	46
Lucas Bone Curette #88 Double Ended 5mm Smooth Spoon Tips	WI-70003839	\$7.99	46

Contents (continued)

Product	Code	Price	Page
McCalls Curette 13S/14S Blue Coated Tips with Hollow Handle	WI-70000301	\$10.99	47
Meyhoefer Chalazion Curette 1.5mm Size #0 Flat Handle 5"	WI-70002416	\$9.99	47
Meyhoefer Chalazion Curette 1.7mm Size #1 Flat Handle 5"	WI-70002415	\$9.99	48
Meyhoefer Chalazion Curette 1mm Size #00 Flat Handle 5"	WI-70002417	\$9.99	48
Meyhoefer Chalazion Curette 2.2mm Size #2 Flat Handle 5"	WI-70002414	\$9.99	49
Meyhoefer Chalazion Curette 2.5mm Size #3 Flat Handle 5"	WI-70002413	\$9.99	49
Meyhoefer Chalazion Curette 3.5mm Size #4 Flat Handle 5"	WI-70002412	\$9.99	50
Molt Periosteal #M9 Elevator Blue Tips Hollow Handle	WI-70000228	\$10.99	50
Molt Periosteal M9 Dental Surgery Tooth Elevator	WI-70003835	\$7.99	51
Novak Endometrial Suction Biopsy Curette 2mm 9 3/4" with Luer Hub	WI-70002398	\$8.99	51
Novak Endometrial Suction Biopsy Curette 3mm 9 3/4" with Luer Hub	WI-70002397	\$8.99	52
Novak Endometrial Suction Biopsy Curette 4mm 9 3/4" with Luer Hub	WI-70002396	\$8.99	52
Rectus Curette 11-Inch Stainless Steel Surgical	WI-32394847	\$800.00	53
Rectus Flat Back Spinal Fusion Curette 11" Straight Size	WI-35858438	\$5.00	53
Root Elevator German Stainless Steel Ergonomic Handle	WI-25002864	\$8.00	54
Sayre Periosteal Elevator 6.75" Double Ended	WI-70003891	\$15.99	54
Spratt Brun Bone Curette 6.5" Oval Cup 2.2×3.4mm Size #5/0	WI-70000667	\$9.99	55
Spratt Brun Bone Curette 6.5" Oval Cup 2.5×4.0mm Size #4/0	WI-70000666	\$9.99	55
Spratt Brun Bone Curette 6.5" Oval Cup 2.8×4.5mm Size #3/0	WI-70000665	\$9.99	56
Spratt Brun Bone Curette 6.5" Oval Cup 3.3×5.1mm Size #2/0	WI-70000664	\$9.99	56
Spratt Brun Bone Curette 6.5" Oval Cup 3.7×5.8mm Size #0	WI-70000663	\$9.99	57
Spratt Brun Bone Curette 6.5" Oval Cup 4.3×6.6mm Size #1	WI-70000662	\$9.99	57
Spratt Brun Bone Curette 6.5" Oval Cup 4.8×7.5mm Size #2	WI-70000661	\$9.99	58
Spratt Brun Bone Curette 6.5" Oval Cup 5.6×9.0mm Size #3	WI-70000660	\$9.99	58
Spratt Brun Bone Curette 6.5" Oval Cup 6.1×9.5mm Size #4	WI-70000659	\$9.99	59
Spratt Brun Bone Curette 6.5" Oval Cup 6.7×10.4mm Size #5	WI-70000658	\$9.99	59
Spratt Brun Bone Curette 6.5" Oval Cup 8.7×12.6mm Size #6	WI-70000657	\$9.99	60
Stainless Steel Cone Ring Curette Hollow Handle	WI-35255073	\$8.00	60
Stainless Steel Heidbrink H1 Straight Elevator Hollow Handle	WI-42909585	\$7.00	61
Stille Type Chisel	WI-70003271	\$76.48	61
Stille Type Osteotome	WI-70003269	\$37.68	62
Trending Elegant Looking Rectus Curette Medical	WI-01193293	\$5.00	62
U.S Army Pattern Osteotome 10mm Straight	WI-70000364	\$10.99	63
U.S Army Pattern Osteotome 12mm Straight	WI-70003268	\$10.99	63
U.S Army Pattern Osteotome 14mm Straight	WI-70000363	\$11.99	64
U.S Army Pattern Osteotome 16mm Straight	WI-70000362	\$13.99	64
U.S Army Pattern Osteotome 18mm Straight	WI-70000361	\$14.99	65
U.S Army Pattern Osteotome 6mm Straight	WI-70000366	\$9.99	65
U.S Army Pattern Osteotome 8mm Straight	WI-70000365	\$9.99	66

**WI-70002077****10 Gracey Curettes 11/12 Blue Titanium Coated Tips Hollow Handle****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Pcs Gracey Curettes 11/12 Titanium Blue Coated Tips Double Ended with Blue & Silver Color Hollow Handle Premium Quality Stainless Steel Periodontal Dental Instruments



Gracey Curette 11/12 is sharp and precise instrument used to remove both supragingival and subgingival calculus in specific areas of the mouth. Gracey 11/12 used for on mesial surfaces. The small contra angle and return angle are complemented by a third angle to offset the blade, which forms a smooth bladed toe. The curette tips are notably long, permitting access to the posterior teeth while keeping the clinician's hand comfortably out of the mouth.

11/12 BLUE TIPS GRACEY CURETTES FEATURES

- Titanium Blue Coated Contra-angle Tips.
- Double-ended Design.
- Ergonomic Blue & Silver Color Hollow Handle.
- Reusable.
- Hand Crafted from High Quality Stainless Steel.

PACKAGE INCLUDE:

- 10 Pcs Gracey Curette 11/12 Titanium Coated Blue Tips – Double Ended – Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curette**APPLICATION:** For Premolars' and Molars' Mesial Surfaces**COLOR:** Titanium Blue Coated Color Tips**WI-70002072****10 Gracey Curettes 13/14 Gold Coated Double Ended****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Pcs Gracey Curettes 13/14 Gold Coated Tips Double Ended Titanium Blue & Silver Color Hollow Handle Periodontal Dental Instruments



The Gracey 13/14 Curette is an area-specific curette that is great for moderate to heavy calculus removal on the distal surfaces of the posterior teeth. Gracey curettes 13/14 are sharp and precise instrument used to remove both supragingival and subgingival calculus in specific areas of the mouth. The small contra angle and return angle are complemented by a third angle to offset the blade, which forms a smooth bladed toe. The curette tips are notably long, permitting access to the posterior teeth while keeping the clinician's hand comfortably out of the mouth.

13/14 GRACEY CURETTES FEATURES

- Goldcoated contra-angle tips
- Double-ended
- Hand crafted from High Quality stainless steel.
- Ergonomic Titanium blue & silver color Hollow Handle.

PACKAGE INCLUDE:

- 10 Pcs Gracey Curettes 13/14 Gold Coated Double Ended Titanium Blue & Silver Color Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curette**APPLICATION:** For Premolars' and Molars' Distal surfaces**COLOR:** Gold Color Coated Tips**TIP:** Double Ended



WI-70003852

10 Periosteal Dental Elevator Molt No. 9

\$39.99 per piece

Trade price from **\$36.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Periosteal Molt 9 (M9) Elevator German Grade Stainless Steel Extraction Surgery Surgical Veterinary Dental Instruments

Molt Periosteal Elevator is an oral instrument commonly used to retract soft tissue in dental procedures. It has a sharp, pointed, narrow end that is used to begin the process of lifting the mucoperiosteum and a slightly wider, more rounded end that completes the tissue removal enough for bone exposure. The rounded edge of the Molt 9 periosteal elevator minimise the tearing of tissue when lifting the gum flap. The sharp end makes the initial lifting of the flap easier.

PERIOSTEAL MOLT #9 FEATURES

- Double end design for optimal lifting & dissecting
- 1 Wider round sharp blade
- 1 sharp, pointed and narrow blade
- Non-reflective surface
- Ergonomic lightweight handle
- Handcrafted with German grade stainless steel
- Corrosion resistance.

PACKAGE INCLUDE:

- 10 Pcs Periosteal Elevator Molt 9 Double Ended 7.5" (19cm) Long

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENT TYPE: Dental Instruments

APPLICATION: Teeth Extraction
RUST PREVENTION PROCEDURE: Passivated



WI-32397761

11.34 Inch Stainless Steel Rectus Curette No

\$5.00 per piece

Trade price from **\$2.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

11.34 Inch Stainless Steel Rectus Curette No is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The 11.34" size shown is our standard; other sizes in the same pattern are made to order.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$5.00, falling to \$2.00 at 25 pieces — about 60% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003879****12 Lucas Bone Curette #84 Small 2.4mm Spoon Double Ended****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Lucas Bone Curette #84 Small 2.4mm Spoon Shaped Double Ended Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 2.4 mm wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium quality stainless steel, it ensures

durability and lasting performance. These curettes are reusable after sterilization.

LUCAS CURETTE #84 FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 2.4mm Spoon Shaped Sharp Ends.
- Double-ended.
- Non-Slip Strong Grip Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 12 Pcs Lucas Bone Curette #84 Small 2.4mm Double Ended

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**WI-70003877****12 Lucas Curette #86 Dentist Lab Dental****\$71.46** per pieceTrade price from **\$65.74**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Lucas Bone Curette #86 Large 3mm Double Ended Smooth Spoon Shaped Tips with Round Pattern Knurled Grip Handle Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Dental Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 3 mm wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium


quality stainless steel, it ensures durability and lasting performance. These curettes are reusable after sterilization.

LUCAS CURETTE #86 FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 3mm Spoon Shaped Sharp Ends.
- Double-ended.
- Smooth Pattern.
- Non-Slip Knurled Grip Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 12 Pcs Lucas Bone Curette #86 Large 3mm Double Ended



**12
PCS****WI-70003876****12 Molt Periosteal M9 Dental Surgery Tooth Elevator****\$44.99** per pieceTrade price from **\$41.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Periosteal Molt 9 (M9) Elevator German Grade Stainless Steel Extraction Surgery Surgical Veterinary Dental Instruments



Molt Periosteal Elevator is an oral instrument commonly used to retract soft tissue in dental procedures. It has a sharp, pointed, narrow end that is used to begin the process of lifting the mucoperiosteum and a slightly wider, more rounded end that completes the tissue removal enough for bone exposure. The rounded edge of the Molt 9 periosteal elevator minimise the tearing of tissue when lifting the gum flap. The sharp end makes the initial lifting of the flap easier.

PERIOSTEAL MOLT #9 FEATURES

- Double end design for optimal lifting & dissecting
- 1 Wider round sharp blade
- 1 sharp, pointed and narrow blade
- Non-reflective surface
- Ergonomic lightweight handle
- Handcrafted with German grade stainless steel
- Corrosion resistance.

PACKAGE INCLUDE:

- 12 Pcs Periosteal Elevator Molt 9 Double Ended 7.5" (19cm) Long

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Instruments**APPLICATION:** Teeth Extraction**RUST PREVENTION PROCEDURE:** Passivated**12
PCS****WI-70003875****12 Sayre Periosteal Elevator 6.75" Double Ended****\$131.99** per pieceTrade price from **\$121.43**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Sayre Periosteal Elevator Set 6.75" (17.1cm) Double-ended (Curved Sharp and Straight Blunt Blades) with Flat Handle – Premium Quality Stainless Steel Surgical Orthopedic Veterinary Instrument.

Sayre Periosteal Elevator is a double-ended surgical instrument designed to elevate, separate, and reflect the periosteum and soft tissues from bones. Its primary applications include orthopedic surgery for fracture repairs and osteotomies, facial and Neuro/ENT surgeries for tissue isolation, and veterinary procedures for soft-tissue separation. Sayre elevator features one straight blunt end and a sharp curved end, allowing versatility in tasks. Made from surgical-grade stainless steel, it is

sterilizable and suitable for various surgical environments.

SAYRE PERIOSTEAL ELEVATOR FEATURES

- Crafted with 6.75" (17.14cm) Overall Length.
- Double Ended Profile
- One Sharp Curved End & One Straight Blunt End.
- Ergonomic Solid Handle for Maximum Control.
- Used to Elevate, Separate and Reflect the Periosteum and Soft Tissue from Bones.
- Sterilizable and Reusable.
- Made With Premium Quality Stainless Steel

PACKAGE INCLUDE:

**WI-70002075****2 Gracey Curettes 13/14 Gold Coated Double Ended****\$13.99** per pieceTrade price from **\$12.87**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Pcs Gracey Curettes 13/14 Gold Coated Tips Double Ended Titanium Blue & Silver Color Hollow Handle Periodontal Dental Instruments

The Gracey 13/14 Curette is an area-specific curette that is great for moderate to heavy calculus removal on the distal surfaces of the posterior teeth. Gracey curettes 13/14 are sharp and precise instrument used to remove both supragingival and subgingival calculus in specific areas of the mouth. The small contra angle and return angle are complemented by a third angle to offset the blade, which forms a smooth bladed toe. The curette tips are notably long, permitting access to the posterior teeth while keeping the clinician's hand comfortably out of the mouth.

13/14 GRACEY CURETTES FEATURES


- Goldcoated contra-angle tips
- Double-ended
- Hand crafted from High Quality stainless steel.
- Ergonomic Titanium blue & silver color Hollow Handle.

PACKAGE INCLUDE:

- 2 Pcs Gracey Curettes 13/14 Gold Coated Double Ended Titanium Blue & Silver Color Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curette**APPLICATION:** For Premolars' and Molars' Distal surfaces**COLOR:** Gold Color Coated Tips**TIP:** Double Ended

GRACEY CURETTES 13/14

**3 PCS****WI-70002079****3 Gracey Curettes 11/12 Blue Titanium Coated Tips Hollow Handle****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Gracey Curettes 11/12 Titanium Blue Coated Tips Double Ended with Blue & Silver Color Hollow Handle Premium Quality Stainless Steel Periodontal Dental Instruments

Gracey Curette 11/12 is sharp and precise instrument used to remove both supragingival and subgingival calculus in specific areas of the mouth. Gracey 11/12 used for on mesial surfaces. The small contra angle and return angle are complemented by a third angle to offset the blade, which forms a smooth bladed toe. The curette tips are notably long, permitting access to the posterior teeth while keeping the clinician's hand comfortably out of the mouth.

- Titanium Blue Coated Contra-angle Tips.
- Double-ended Design.
- Ergonomic Blue & Silver Color Hollow Handle.
- Reusable.
- Hand Crafted from High Quality Stainless Steel.

PACKAGE INCLUDE:

- 3 Pcs Gracey Curette 11/12 Titanium Coated Blue Tips – Double Ended – Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curette**APPLICATION:** For Premolars' and Molars' Mesial Surfaces**COLOR:** Titanium Blue Coated Color Tips**11/12 BLUE TIPS GRACEY CURETTES FEATURES**

**WI-70002074****3 Gracey Curettes 13/14 Gold Coated Double Ended****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Gracey Curettes 13/14 Gold Coated Tips Double Ended Titanium Blue & Silver Color Hollow Handle Periodontal Dental Instruments

The Gracey 13/14 Curette is an area-specific curette that is great for moderate to heavy calculus removal on the distal surfaces of the posterior teeth. Gracey curettes 13/14 are sharp and precise instrument used to remove both supragingival and subgingival calculus in specific areas of the mouth. The small contra angle and return angle are complemented by a third angle to offset the blade, which forms a smooth bladed toe. The curette tips are notably long, permitting access to the posterior teeth while keeping the clinician's hand comfortably out of the mouth.

13/14 GRACEY CURETTES FEATURES


- Goldcoated contra-angle tips
- Double-ended
- Hand crafted from High Quality stainless steel.
- Ergonomic Titanium blue & silver color Hollow Handle.

PACKAGE INCLUDE:

- 3 Pcs Gracey Curettes 13/14 Gold Coated Double Ended Titanium Blue & Silver Color Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curette**APPLICATION:** For Premolars' and Molars' Distal surfaces**COLOR:** Gold Color Coated Tips**TIP:** Double Ended

GRACEY CURETTES 13/14

**WI-70003871****3 Lucas Bone Curette #84 Small 2.4mm Spoon Double Ended****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Lucas Bone Curette #84 Small 2.4mm Spoon Shaped Double Ended Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 2.4 mm wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium quality stainless steel, it ensures

durability and lasting performance. These curettes are reusable after sterilization.

LUCAS CURETTE #84 FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 2.4mm Spoon Shaped Sharp Ends.
- Double-ended.
- Non-Slip Strong Grip Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 3 Pcs Lucas Bone Curette #84 Small 2.4mm Double Ended

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**3
PCS**

**WI-70002395****3 Novak Endometrial Suction Biopsy Curettes
9.75" 2, 3 & 4mm Dia****\$24.99** per pieceTrade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Set of 3 Endometrial Biopsy Curette NOVAK 9 3/4" (24.8cm) Length Single Ended, 2mm, 3mm & 4mm Tip Curved Round Tip with Sharp Serrated Teeth OB/Gynecology Surgical Instrument

Novak Suction Biopsy Curette is primarily used for scraping and removing a portion of the lining of the uterus. Novak Endocervical Suction Biopsy Curette features an angled, serrated tip to aid in the removal and collection of biopsy tissue from the endocervix. The tip is at the end of a long, slender handle. The Novak curette is used with a suction device that contains the tissue in a syringe after it has been scraped and removed. It comes in multiple variations such as a curette with Leur Hub in 2mm, 3mm and 4mm diameter

**NOVAK BIOPSY CURETTE
FEATURES:**

- Long and thin metallic tube narrow built for deep access.
- Broad suction hole for easy suctioning.
- Blunt tip to prevent piercing injuries.
- Ergonomic handle for better control
- Lightweight for comfortable handling

PACKAGE INCLUDE:

- Novak Endometrial Suction Biopsy Curette Size 9 3/4" (24.8cm) with Luer Hub 2mm Diameter.
- Novak Endometrial Suction Biopsy Curette Size 9 3/4" (24.8cm) with Luer Hub 3mm Diameter.
- Novak Endometrial Suction Biopsy Curette Size 9 3/4" (24.8cm) with Luer Hub 4mm Diameter.

**WI-70003867****3 Sayre Periosteal Elevators 6.75" Double Ended****\$38.99** per pieceTrade price from **\$35.87**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Sayre Periosteal Elevator Set 6.75" (17.1cm) Double-ended (Curved Sharp and Straight Blunt Blades) with Flat Handle – Premium Quality Stainless Steel Surgical Orthopedic Veterinary Instrument.

Sayre Periosteal Elevator is a double-ended surgical instrument designed to elevate, separate, and reflect the periosteum and soft tissues from bones. Its primary applications include orthopedic surgery for fracture repairs and osteotomies, facial and Neuro/ENT surgeries for tissue isolation, and veterinary procedures for soft-tissue separation. Sayre elevator features one straight blunt end and a sharp curved end, allowing versatility in tasks. Made from surgical-grade stainless steel, it is

sterilizable and suitable for various surgical environments.

**SAYRE PERIOSTEAL ELEVATOR
FEATURES**

- Crafted with 6.75" (17.14cm) Overall Length.
- Double Ended Profile
- One Sharp Curved End & One Straight Blunt End.
- Ergonomic Solid Handle for Maximum Control.
- Used to Elevate, Separate and Reflect the Periosteum and Soft Tissue from Bones.
- Sterilizable and Reusable.
- Made With Premium Quality Stainless Steel

PACKAGE INCLUDE:



WI-70000466

5 Buck Ear Curettes 6.5" Straight Sharp Size #00, #0, #1, #2, #3

\$34.99 per piece

Trade price from **\$32.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Buck Ear Curettes 6.5" (16.5cm) Straight Sharp Ring Blade Size #00, #0, #1, #2, #3 with Round Pattern Knurled Grip Handle Stainless Steel ENT Surgical Veterinary Instruments.

Buck Ear Curettes are specialized surgical tools used in ear surgeries to safely remove tumors, clean infections, treat pain and prevent hearing loss through biopsies. This Ear Curette offer surgical options with blunt ring blades for sampling and sharp blades for deep cleaning, available with straight or angled necks.

Buck Ear 6.5" (16.5cm) curettes, made from premium grade stainless steel which making it reusable and autoclavable, are designed to effectively treat deep diseases and have slender shafts for clear working

tips during delicate samples. The veterinary practitioners also use Buck Ear Curettes for ENT procedures like vertical cleaning cealing to remove ear wax and other foreign objects. The long handle to give precise functioning during multiple ENT practices.

BUCK EAR CURETTE FEATURES

- 6.5" (16.5cn) Long.
- Straight Sharp Profile.
- Slender Shafts For Deep Surgical Access.
- Non-Slip Handle For Controlled Grips.
- Multiple Model and Sizes For Diverse Indications.
- Designed for Ear Canal Cleaning & Removes Cerumen.



WI-70002078

5 Gracey Curettes 11/12 Blue Titanium Coated Tips Hollow Handle

\$29.99 per piece

Trade price from **\$27.59**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Gracey Curettes 11/12 Titanium Blue Coated Tips Double Ended with Blue & Silver Color Hollow Handle Premium Quality Stainless Steel Periodontal Dental Instruments

Gracey Curette 11/12 is sharp and precise instrument used to remove both supragingival and subgingival calculus in specific areas of the mouth. Gracey 11/12 used for on mesial surfaces. The small contra angle and return angle are complemented by a third angle to offset the blade, which forms a smooth bladed toe. The curette tips are notably long, permitting access to the posterior teeth while keeping the clinician's hand comfortably out of the mouth.

- Titanium Blue Coated Contra-angle Tips.
- Double-ended Design.
- Ergonomic Blue & Silver Color Hollow Handle.
- Reusable.
- Hand Crafted from High Quality Stainless Steel.

PACKAGE INCLUDE:

- 5 Pcs Gracey Curette 11/12 Titanium Coated Blue Tips – Double Ended – Hollow Handle


TECHNICAL SPECIFICATIONS

MATERIAL: Stainless Steel


INSTRUMENT TYPE: Dental Curette

APPLICATION: For Premolars' and Molars' Mesial Surfaces

COLOR: Titanium Blue Coated Color Tips



11/12 BLUE TIPS GRACEY CURETTES FEATURES

**WI-70002073****5 Gracey Curettes 13/14 Gold Coated Double Ended****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Gracey Curettes 13/14 Gold Coated Tips Double Ended Titanium Blue & Silver Color Hollow Handle Periodontal Dental Instruments

The Gracey 13/14 Curette is an area-specific curette that is great for moderate to heavy calculus removal on the distal surfaces of the posterior teeth. Gracey curettes 13/14 are sharp and precise instrument used to remove both supragingival and subgingival calculus in specific areas of the mouth. The small contra angle and return angle are complemented by a third angle to offset the blade, which forms a smooth bladed toe. The curette tips are notably long, permitting access to the posterior teeth while keeping the clinician's hand comfortably out of the mouth.

13/14 GRACEY CURETTES FEATURES



- Goldcoated contra-angle tips
- Double-ended
- Hand crafted from High Quality stainless steel.
- Ergonomic Titanium blue & silver color Hollow Handle.

PACKAGE INCLUDE:

- 5 Pcs Gracey Curettes 13/14 Gold Coated Double Ended Titanium Blue & Silver Color Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curette**APPLICATION:** For Premolars' and Molars' Distal surfaces**COLOR:** Gold Color Coated Tips**TIP:** Double Ended

GRACEY CURETTES 13/14

**WI-70002390****5 Kevorkian Youngs Biopsy Curette 12" with Basket Loop Size 3x12mm****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Kevorkian-Youngs Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, with Basket, Angular Premium Quality Stainless Steel Obstetrical/Gynecology Surgical Instruments

Kevorkian Biopsy Curette is a specialized gynecological instrument used for dissecting and removing part of the uterine wall for diagnostic purposes. Kevorkian-Youngs Endocervical Biopsy Curette is designed for the excision of tissue from the endocervix for biopsy. The curette is available with or without a basket for use in excision or biopsy. The Kevorkian Biopsy curette also comes with sharp edges that effectively scrape off the uterine tissue and the shape allows the surgeon to easily

scoop up and collect the dissected tissue for secure removal.


KEVORKIAN YOUNG BIOPSY CURETTE FEATURES

- Long shank for deep tissue access.
- Ergonomic solid handle for superior control.
- With Or Without basket for optimum dissection.
- Smooth outer curves to prevent inadvertent soft tissue damage.
- Sharp edges for effective scraping off the uterine tissues.

PACKAGE INCLUDE:

- 5 Pcs Kevorkian-Youngs Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, with Basket

KEVORKIAN BIOPSY CURETTE WITH BASKET



**WI-70002248****5 Kevorkian-Younger Biopsy Curette 12" Loop 3X12mm with Basket Half Gold Handle****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pieces Kevorkian-Younger Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, with Basket Half Gold Handle – Angular Premium Quality Premium Quality Stainless Steel Obstetrical/Gynecology Surgical Instruments

Kevorkian Biopsy Curette is a specialized gynecological instrument used for dissecting and removing part of the uterine wall for diagnostic purposes. Kevorkian-Younger Endocervical Biopsy Curette is designed for the excision of tissue from the endocervix for biopsy. The curette is available with or without a basket for use in excision or biopsy. The Kevorkian Biopsy curette also comes with sharp edges that effectively scrape off the uterine tissue and the

shape allows the surgeon to easily scoop up and collect the dissected tissue for secure removal.

KEVORKIAN-YOUNGER ENDOCERVICAL BIOPSY CURETTE FEATURES

- Long shank for deep tissue access.
- Ergonomic solid half gold handle for superior control.
- With Or Without basket for optimum dissection.
- Smooth outer curves to prevent inadvertent soft tissue damage.
- Sharp edges for effective scraping off the uterine tissues.

PACKAGE INCLUDE:

- 5 Pcs Kevorkian-Younger Endocervical Biopsy Curette 12"

WI-70002249**5 Kevorkian-Younger Biopsy Curette 12" Loop 3X12mm Without Basket Half Gold Handle****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Kevorkian-Younger Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, without Basket – Half Gold Handle – Angular Premium Quality Stainless Steel Obstetrical/Gynecology Surgical Instruments

Kevorkian Biopsy Curette is a specialized gynecological instrument used for dissecting and removing part of the uterine wall for diagnostic purposes. Kevorkian-Younger Endocervical Biopsy Curette is designed for the excision of tissue from the endocervix for biopsy. The curette is available with or without a basket for use in excision or biopsy. The Kevorkian Biopsy curette also comes with sharp edges that effectively scrape off the uterine tissue and the

shape allows the surgeon to easily scoop up and collect the dissected tissue for secure removal.

KEVORKIAN-YOUNGER ENDOCERVICAL BIOPSY CURETTE FEATURES

- Long shank for deep tissue access.
- Ergonomic solid Half Gold Handle for superior control.
- With Or Without basket for optimum dissection.
- Smooth outer curves to prevent inadvertent soft tissue damage.
- Sharp edges for effective scraping off the uterine tissues.

PACKAGE INCLUDE:

- 5 Pcs Kevorkian-Younger Endocervical Biopsy Curette 12"

**5 PCS****5 PCS**

**WI-70003124****5 Offset Sinus Osteotomes with Sterilization Cassette 7.25" X 5"****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Sterilization Cassette with 5 Offset Sinus Lift Osteotomes Set with L-Key Adjustment Stops Concave Tips Diameters: 2.0, 3.0, 3.3, 4.0, 5.0mm Depths Marked at 8, 10, 13, 15, 17 & 20mm Surgical Dental Implant Instruments.

Sterilization cassette & Osteotome set includes 5 Offset concave tips sinus osteotomes. The tips of the osteotomes are tapered and concave, making them able to compress bone graft into the sinus window. The osteotomes provided have L-key adjustable stoppers, preventing sinus perforation. The stainless steel construction ensures maximum durability and dependability. You can reuse these sinus lift osteotomes over time as they will not bend or tarnish on you.

Sterilization cassette ensuring the safety and efficiency of dental and medical procedures. In the fast-evolving world of healthcare, maintaining impeccable hygiene standards is paramount, and the sterilization cassette emerges as an unsung hero in this endeavor. These sleek, often stainless steel with removable plastic inserts, not only organize and protect instruments during sterilization but also ensure that each tool is ready for use in its optimal condition.

STERILIZATION CASSETTE WITH SINUS LIFT OSTEOTOMES SET INCLUDES

- 1 x Offset Sinus Osteotomes L-Key Adjustment Stop (Concave) Diameter 2.0

**WI-70003865****5 Pc Molt Periosteal M9 Dental Elevator****\$27.49** per pieceTrade price from **\$25.29**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Periosteal Molt 9 (M9) Elevator German Grade Stainless Steel Extraction Surgery Surgical Veterinary Dental Instruments

Molt Periosteal Elevator is an oral instrument commonly used to retract soft tissue in dental procedures. It has a sharp, pointed, narrow end that is used to begin the process of lifting the mucoperiosteum and a slightly wider, more rounded end that completes the tissue removal enough for bone exposure. The rounded edge of the Molt 9 periosteal elevator minimise the tearing of tissue when lifting the gum flap. The sharp end makes the initial lifting of the flap easier.

PERIOSTEAL MOLT #9 FEATURES

- Double end design for optimal lifting & dissecting
- 1 Wider round sharp blade
- 1 sharp, pointed and narrow blade
- Non-reflective surface
- Ergonomic lightweight handle
- Handcrafted with German grade stainless steel
- Corrosion resistance.


PACKAGE INCLUDE:

- 5 Pcs Periosteal Elevator Molt 9 Double Ended 7.5" (19cm) Long

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Instruments**APPLICATION:** Teeth Extraction**RUST PREVENTION PROCEDURE:**

Passivated



**WI-70003861****6 Lucas Bone Curette #84 Small 2.4mm Spoon Double Ended****\$35.99** per pieceTrade price from **\$33.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Lucas Bone Curette #84 Small 2.4mm Spoon Shaped Double Ended Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 2.4 mm wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium quality stainless steel, it ensures

durability and lasting performance. These curettes are reusable after sterilization.

LUCAS CURETTE #84 FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 2.4mm Spoon Shaped Sharp Ends.
- Double-ended.
- Non-Slip Strong Grip Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 6 Pcs Lucas Bone Curette #84 Small 2.4mm Double Ended

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**WI-70003858****6 Molt Periosteal M9 Dental Surgery Tooth Elevator****\$29.99** per pieceTrade price from **\$27.59**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Periosteal Molt 9 (M9) Elevator German Grade Stainless Steel Extraction Surgery Surgical Veterinary Dental Instruments

Molt Periosteal Elevator is an oral instrument commonly used to retract soft tissue in dental procedures. It has a sharp, pointed, narrow end that is used to begin the process of lifting the mucoperiosteum and a slightly wider, more rounded end that completes the tissue removal enough for bone exposure. The rounded edge of the Molt 9 periosteal elevator minimise the tearing of tissue when lifting the gum flap. The sharp end makes the initial lifting of the flap easier.

- Double end design for optimal lifting & dissecting
- 1 Wider round sharp blade
- 1 sharp, pointed and narrow blade
- Non-reflective surface
- Ergonomic lightweight handle
- Handcrafted with German grade stainless steel
- Corrosion resistance.

PACKAGE INCLUDE:

- 6 Pcs Periosteal Elevator Molt 9 Double Ended 7.5" (19cm) Long

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Instruments**APPLICATION:** Teeth Extraction**RUST PREVENTION PROCEDURE:** Passivated**PERIOSTEAL MOLT #9 FEATURES**

**WI-70003857****6 Sayre Periosteal Elevator 6.75\"/>****\$71.99** per pieceTrade price from **\$66.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Sayre Periosteal Elevator Set 6.75\" (17.1cm) Double-ended (Curved Sharp and Straight Blunt Blades) with Flat Handle – Premium Quality Stainless Steel Surgical Orthopedic Veterinary Instrument.

Sayre Periosteal Elevator is a double-ended surgical instrument designed to elevate, separate, and reflect the periosteum and soft tissues from bones. Its primary applications include orthopedic surgery for fracture repairs and osteotomies, facial and Neuro/ENT surgeries for tissue isolation, and veterinary procedures for soft-tissue separation. Sayre elevator features one straight blunt end and a sharp curved end, allowing versatility in tasks. Made from surgical-grade stainless steel, it is

sterilizable and suitable for various surgical environments.

SAYRE PERIOSTEAL ELEVATOR FEATURES

- Crafted with 6.75\" (17.14cm) Overall Length.
- Double Ended Profile
- One Sharp Curved End & One Straight Blunt End.
- Ergonomic Solid Handle for Maximum Control.
- Used to Elevate, Separate and Reflect the Periosteum and Soft Tissue from Bones.
- Sterilizable and Reusable.
- Made With Premium Quality Stainless Steel

PACKAGE INCLUDE:**WI-70000611****Adson Periosteal Elevator 7.5\"/>****\$24.99** per pieceTrade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Adson Periosteal Elevator 7.5\" (19.1cm) Semi-Sharp Curved 10mm Wide Blade with Hollow Handle Premium Grade Stainless Steel Orthopedic Neurosurgical Instrument.

Adson periosteal elevator is a specialized instrument used for lifting and cutting periosteum tissue to scrape flaps and reveal diseased or fractured bone structures during neurosurgery and orthopedic surgical procedures. The Adson periosteal elevator has a cylindrical, straight shank with a 10mm-wide, semi-sharp, round-ended, curved blade, which allows surgeons to dissect periosteum from multiple types of bone. The Adson Elevator features a hollow handle with three vertical grooves to secure a firm grip. Additionally, the handle transforms into

a cylindrical, straight shank, which enhances accuracy. Made with premium-grade stainless steel, ensuring a long-lasting experience and well-suited for multiple medical approaches. Can be sterilized for repeated use.

ADSON PERIOSTEAL ELEVATOR FEATURES

- Crafted with 7.5\" (19.1cm) Overall Length.
- Round End Semi-Sharp Curved Profile.
- 10mm Wide Blade.
- One-Sided Beveled Blade for Smooth Tissue Dissection.
- Cylindrical Straight Shank Easy and Deep Access.





WI-70000612

Adson Periosteal Elevator 7.5" Semi Sharp Curved 7mm Wide

\$24.99 per piece

Trade price from **\$22.99**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Adson Periosteal Elevator 7.5" (19.1cm) Semi-Sharp Curved 7mm Wide Blade with Hollow Handle Premium Grade Stainless Steel Orthopedic Neurosurgical Instrument.

Adson periosteal elevator is a specialized instrument used for lifting and cutting periosteum tissue to scrape flaps and reveal diseased or fractured bone structures during neurosurgery and orthopedic surgical procedures. The Adson periosteal elevator has a cylindrical, straight shank with a 7mm-wide, semi-sharp, round-ended, curved blade, which allows surgeons to dissect periosteum from multiple types of bone. The Adson Elevator features a hollow handle with three vertical grooves to secure a firm grip. Additionally, the handle transforms into

a cylindrical, straight shank, which enhances accuracy. Made with premium-grade stainless steel, ensuring a long-lasting experience and well-suited for multiple medical approaches. Can be sterilized for repeated use.

ADSON PERIOSTEAL ELEVATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Round End Semi-Sharp Curved Profile.
- 7mm Wide Blade.
- One-Sided Beveled Blade for Smooth Tissue Dissection.
- Cylindrical Straight Shank Easy and Deep Access.



WI-70003278

Alexander Chisel

\$33.09 per piece

Trade price from **\$30.44**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Alexander Chisel is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$33.09, falling to \$30.44 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000178****Billeau Ear Loop Curette Size #1 – Small****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Billeau Ear Loop Curette 6.5" (16.5cm)
Size #1 Small Loop Shaped Blade
Premium Quality Stainless Steel Ear
Cleaning ENT Otology Surgery
Instruments.

Billeau Ear Loop is a surgical instrument used to remove cerumen (ear-wax) and foreign bodies from the external auditory meatus during various otology procedures particularly for reconstructive surgery and tympanoplasty (eardrum repair). Billeau Ear Curette is a large, straight, 6-1/2 inch straight curette designed for cleaning the ear canal and removing cerumen, with a sharp, teardrop, fenestrated wire tip. This Ear Loop Curette handles feature a round profile for comfort and dexterity, and outer serrations for confidence and accuracy


in device movement. Made from premium grade stainless steel, ensuring years of reliable and high-quality service.

BILLEAU EAR LOOP FEATURES:

- 6.5" (16.5cm) Long.
- Various Size Blade Tips: #1, #2, #3.
- Loop Shaped Ring Blades for Cleaning or Sampling Procedures.
- Rounded Knurled Handle for Balanced Movements.
- Slender Shafts for Easy Access to the Deep Surgical Sites.
- Easily Sterilizable and Reusable.
- Made with Premium Grade Stainless Steel.

PACKAGE INCLUDES:

- Billeau Ear Loop Curette 6.5" (16.5cm) Size #1 – Small

**WI-70000177****Billeau Ear Loop Curette Size #2 – Medium****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Billeau Ear Loop Curette 6.5" (16.5cm)
Size #2 Medium Loop Shaped Blade
Premium Quality Stainless Steel Ear
Cleaning ENT Otology Surgery
Instruments.

Billeau Ear Loop is a surgical instrument used to remove cerumen (ear-wax) and foreign bodies from the external auditory meatus during various otology procedures particularly for reconstructive surgery and tympanoplasty (eardrum repair). Billeau Ear Curette is a large, straight, 6-1/2 inch straight curette designed for cleaning the ear canal and removing cerumen, with a sharp, teardrop, fenestrated wire tip. This Ear Loop Curette handles feature a round profile for comfort and dexterity, and outer serrations for confidence and accuracy


in device movement. Made from premium grade stainless steel, ensuring years of reliable and high-quality service.

BILLEAU EAR LOOP FEATURES:

- 6.5" (16.5cm) Long.
- Various Size Blade Tips: #1, #2, #3.
- Loop Shaped Ring Blades for Cleaning or Sampling Procedures.
- Rounded Knurled Handle for Balanced Movements.
- Slender Shafts for Easy Access to the Deep Surgical Sites.
- Easily Sterilizable and Reusable.
- Made with Premium Grade Stainless Steel.

PACKAGE INCLUDES:

- Billeau Ear Loop Curette 6.5" (16.5cm) Size #2 – Medium





WI-70000176

Billeau Ear Loop Curette Size #3 – Large

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Billeau Ear Loop Curette 6.5" (16.5cm)
Size #3 Large Loop Shaped Blade
Premium Quality Stainless Steel Ear
Cleaning ENT Otolaryngology Surgery
Instruments.

Billeau Ear Loop is a surgical instrument used to remove cerumen (ear-wax) and foreign bodies from the external auditory meatus during various otology procedures particularly for reconstructive surgery and tympanoplasty (eardrum repair). Billeau Ear Curette is a large, straight, 6-1/2 inch straight curette designed for cleaning the ear canal and removing cerumen, with a sharp, teardrop, fenestrated wire tip. This Ear Loop Curette handles feature a round profile for comfort and dexterity, and outer serrations for confidence and accuracy

in device movement. Made from premium grade stainless steel, ensuring years of reliable and high-quality service.

BILLEAU EAR LOOP FEATURES:

- 6.5" (16.5cm) Long.
- Various Size Blade Tips: #1, #2, #3.
- Loop Shaped Ring Blades for Cleaning or Sampling Procedures.
- Rounded Knurled Handle for Balanced Movements.
- Slender Shafts for Easy Access to the Deep Surgical Sites.
- Easily Sterilizable and Reusable.
- Made with Premium Grade Stainless Steel.

PACKAGE INCLUDES:

- Billeau Ear Loop Curette 6.5" (16.5cm) Size #3 – Large



WI-46975459

Cage Curette Surgical

\$15.00 per piece

Trade price from \$12.00
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Cage Curette Surgical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$15.00, falling to \$12.00 at 5 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003276

Cobb Osteotome

\$38.19 per piece

Trade price from \$35.13
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Cobb Osteotome is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$38.19, falling to \$35.13 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35248110

Cone Ring Curette 9 Inch Angled

\$8.00 per piece

Trade price from \$6.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Cone Ring Curette 9 Inch Angled for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: an angled shank that keeps the hand clear of the sightline.

The 9" size shown is our standard; other sizes in the same pattern are

made to order. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$8.00, falling to \$6.00 at 10 pieces — about 25% off the single-unit price. Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000627

Cone Ring Curette 9" Straight 6mm Size #2

\$12.99 per piece

Trade price from \$11.95
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Cone Ring Curette 9" (22.9cm) Straight 6mm Round Fenestrated Cup Size #2 with Knurled Handle Premium Quality Stainless Steel Neurosurgery Orthopedic Instruments.

Cone ring curettes straight remove the unwanted tissue growth in the neurosurgeries. These curettes are designed to effectively remove excess tissue for sampling or to remove growths from regions of the body encountered in neurosurgery. Cone ring curette features a ring at the end of the shaft that enables it to precisely and accurately remove tissue for collection or to increase access to the needed surgical field. Additionally, this curette is available in different diameters as well as either angled or straight, making it suited for many

neurosurgical/orthopedic procedures and surgeon preferences. This instrument has a knurled handle to provide precise performance. Made with German-grade stainless steel material, ensuring performance durability and making it sterilizable/easy to clean.

CONE RING CURETTES FEATURES

- Crafted with 9" (22.9cm) Overall Length.
- Round Fenestrated Straight Profile.
- 6mm Ring Diameter Size #2.
- Various Pattern Cup Sizes.
- Knurled Grip Round Handle Ensure Better Control.
- Rustproof and Reusable.
- Made With Premium Quality Stainless Steel.



WI-70003254

Coupland Root Elevators 2 0 mm

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Coupland Root Elevators 2 0 mm is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003253

Coupland Root Elevators 2.5 mm

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Coupland Root Elevators 2.5 mm is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003252

Coupland Root Elevators 3 0 mm

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Coupland Root Elevators 3 0 mm is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46993252

Curette With Cage Reusable Surgical Tool for Precise Curettage

\$14.80 per piece

Trade price from \$11.80
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Curette With Cage Reusable Surgical Tool for Precise Curettage for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$14.80, falling to \$11.80 at 5 pieces — about 20% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47001203

Curette with Cage Surgical Kit

\$14.90 per piece

Trade price from \$11.90
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Curette with Cage Surgical Kit — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. Supplied as a set — the contents can be adjusted to your list.

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$14.90, falling to \$11.90 at 5 pieces — about 20% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Forging, milling, heat treatment, hand finishing and passivation all happen

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003227****Dental Columbia Universal Curette 2R/2L****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Anterior Columbia Universal Curette 2R/2L Double Ended 6.5" (16.5cm) Premium Quality Stainless Steel Heavy Calculus Removal Dental Instruments.

Columbia Curette 2R-2L is a sharp and precise dental instrument used to remove both supragingival and subgingival calculus. 2R/2L Columbia University Curette feature a longer shank and shorter thicker blade for use on heavier anterior deposits. 2L-2R Columbia Curette is a dentist choice scaler for heavy calculus removal. The sharp mirror ends are elongated with round toe able to reach deep into the mouth and easily cut away heavy and challenging buildup. 2R/2L Columbia University Curette manufactured with premium quality stainless steel to last long and guaranteed to satisfy standard

practices of dental settings with exceptional performance. Fully autoclaved and reusable.

COLUMBIA UNIVERSAL CURETTE 2R/2L FEATURES

- Crafted With 6.5" Overall Length.
- Double Ended.
- Short Thick Blade.
- 14mm Long Terminal Shank.
- Rounded Tips Sharp on Both Sides.
- Designed for Heavier Anterior Deposits.
- Round Pattern Non-Slip Handle for Better Handling.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE:**WI-70003226****Dental Columbia Universal Curette 4R/4L Double Ended****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Posterior Columbia Universal Curette 4R/4L Double Ended 6.5" (16.5cm) Premium Quality Stainless Steel Heavy Calculus Removal Dental Instruments.


Columbia 4R/4L Curette is the universal curette made for scaling posterior teeth. With super sharp rounded tips 4R/4L columbia universal curette is designed for removing heavy and moderate calculus from sub gingival and supragingival tooth surface. It is considered an ideal instrument for clinicians who are looking for low-maintenance dental instruments that provide great precision and durable sharpness. Columbia Curette 4R-4L manufactured with premium quality stainless steel to last long and guaranteed to satisfy standard practices of dental settings with

exceptional performance. Fully autoclaved and reusable.

COLUMBIA UNIVERSAL CURETTE 4L/4R FEATURES

- Crafted With 6.5" Overall Length.
- Double Ended.
- 12mm Long Terminal Shank.
- Two Sharp Cutting Edges With Rounded Toe.
- Designed for Heavier Posterior Deposits.
- Round Pattern Non-Slip Handle for Better Handling.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE:

**WI-70003228****Dental Columbia Universal Curette Scalars 13/14****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Columbia Universal Curette 13/14 Double Ended 6.5" (16.5cm) Premium Quality Stainless Steel Anterior and Posterior Moderate to Heavy Calculus Removal Dental Instruments.

Columbia 13/14 Universal Curette used for scaling, root planing and periodontal debridement of supragingival and subgingival tooth surfaces. 13/14 Columbia Curette means you can use it for both posterior and anterior on all surfaces. These instruments remove moderate to heavy calculus from supra- and subgingival tooth surfaces. The sharp, mirror-ended tips of 13/14 dental universal curette allow you to achieve every angle with ease. The rounded toe gives you room to navigate tight spaces and remove every last bit of calculus. Columbia University

Curette 13-14 manufactured with premium quality stainless steel to last long and guaranteed to satisfy standard practices of dental settings with exceptional performance. Fully autoclaved and reusable.

COLUMBIA UNIVERSAL CURETTE 13/14 FEATURES

- Crafted With 6.5" Overall Length.
- Double Ended.
- Rounded Tips Sharp on Both Sides.
- For Anterior and Posterior Calculus Removal
- Designed for Moderate to Heavy Deposits.
- Round Pattern Non-Slip Handle for Better Handling.
- Reusable.

**WI-70000441****Dental Elevators Autoclave Sterilization Cassette Tray 7" x 8" for 8 Instruments****\$32.99** per pieceTrade price from **\$30.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Elevators Sterilization Cassette Tray 7" X 8" X 1.5" for 8 Instruments Instruments – Premium Quality Stainless Steel Container With Removable Rubber/ Silicone Inserts – Heavy Duty Lockable Sterilization Tray Rack Box Dental Tools – Efficient Dental Surgical Lab Organizer.


The dental sterilization tray measures 7" x 8" x 1.5" and features removable rubber inserts and a secure hinged lid for safe instrument enclosure during sterilization processes. The Dental Sterilization cassette is a compact, versatile tool that optimizes space utilization in dental clinics by accommodating various dental instruments, promoting efficient organization and workflow management.

DENTAL ELEVATORS STERILIZATION CASSETTE RACK BOX FEATURES

- Stainless Steel Cassettes, 7" x 8" x 1.5"
- Hold Up to 8 Instruments.
- Protective Silicone Rubber Inserts Scratches and Impact Damage.
- Secure Hinged Lid with Locking Mechanism.
- Designed for Autoclave Sterilization for Easy Maintenance and Safety.
- Can Withstand the Rigors of Repeated Use.
- Compatible With Diverse Dental and Surgical Instruments.

PACKAGE INCLUDE

- 1 x Elevators Sterilization Cassette Tray 7" x 8" x 1.5" for 8 Instrument

**MAIN FEATURES**

- Corrosion & rust resistant
- Premium quality construction
- Easy carry
- Durable material
- Comes with removable silicon inserts



WI-70003623

Dental Gracey Curettes 11/12 Double Ended

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Gracey Curettes 11/12 Double Ended is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35495799

Double Ended Micro Curettes Medical

\$4.00 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Double Ended Micro Curettes Medical — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. From the pattern: a micro-scale working end for fine and delicate work, and two working ends on one handle, halving what you carry.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$4.00, falling to \$2.00 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-35506825

Double Ended Micro Curettes Precision Surgical Tools

\$4.00 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Double Ended Micro Curettes Precision Surgical Tools — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. From the pattern: a micro-scale working end for fine and delicate work, and two working ends on one handle, halving what you carry.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$2.00 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-35538762

Double Ended Micro Curettes Stainless Steel Surgical Instrument Set

\$4.00 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Double Ended Micro Curettes Stainless Steel Surgical Instrument Set for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. From the pattern: a micro-scale working end for fine and delicate work, and two working ends on one handle, halving what you carry.

Supplied as a set — the contents can be adjusted to your list. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$4.00, falling to \$2.00 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

**WI-70002649****Fox Dermal Curette 1mm Round Tip – 5.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Fox Dermal Curettes Dermatology 1mm Round Tip – 5.5\"/>

Fox Dermal Curette is commonly used in dermatology to scrape the skin growths such as excess skin surface, warts, and cancer. It features a ring-shaped cutting edge and a flat handle. It comes with multiple variations, so you can pick the best one for your surgical kit. It is crafted from high-grade surgical stainless steel.

Fox Dermal Curette is used for scraping skin lesions and growths such as warts. It is a single-ended instrument available with an oval or round loop tip. The loop is sharp on the inside and blunt on the outside so the lesion can be removed,

but the surrounding skin is not damaged. A variety of tip sizes are available to accommodate varying lesion sizes. The handle of the Fox curette is flat and lightweight.

We offer an extensive range of surgical curettes with multiple variations. Our versatile instruments are perfect to be a part of your amazing surgical kit for years if cared for properly. These instruments are widely preferred to perform a broad range of surgical operations. We can also design our surgeons' custom kit packs according to their needs. The instruments are all fine, durable, and sharp.

FOX DERMAL CURETTE FEATURES**WI-70002648****Fox Dermal Curette 2mm Round Tip – 5.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Fox Dermal Curettes Dermatology 2mm Round Tip – 5.5\"/>

Fox Dermal Curette is commonly used in dermatology to scrape the skin growths such as excess skin surface, warts, and cancer. It features a ring-shaped cutting edge and a flat handle. It comes with multiple variations, so you can pick the best one for your surgical kit. It is crafted from high-grade surgical stainless steel.

Fox Dermal Curette is used for scraping skin lesions and growths such as warts. It is a single-ended instrument available with an oval or round loop tip. The loop is sharp on the inside and blunt on the outside so the lesion can be removed,

but the surrounding skin is not damaged. A variety of tip sizes are available to accommodate varying lesion sizes. The handle of the Fox curette is flat and lightweight.

We offer an extensive range of surgical curettes with multiple variations. Our versatile instruments are perfect to be a part of your amazing surgical kit for years if cared for properly. These instruments are widely preferred to perform a broad range of surgical operations. We can also design our surgeons' custom kit packs according to their needs. The instruments are all fine, durable, and sharp.

FOX DERMAL CURETTE FEATURES

**WI-70002647****Fox Dermal Curette 3mm Round Tip – 5.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Fox Dermal Curettes Dermatology 3mm Round Tip – 5.5\"/>

Fox Dermal Curette is commonly used in dermatology to scrape the skin growths such as excess skin surface, warts, and cancer. It features a ring-shaped cutting edge and a flat handle. It comes with multiple variations, so you can pick the best one for your surgical kit. It is crafted from high-grade surgical stainless steel.

Fox Dermal Curette is used for scraping skin lesions and growths such as warts. It is a single-ended instrument available with an oval or round loop tip. The loop is sharp on the inside and blunt on the outside so the lesion can be removed,

but the surrounding skin is not damaged. A variety of tip sizes are available to accommodate varying lesion sizes. The handle of the Fox curette is flat and lightweight.

We offer an extensive range of surgical curettes with multiple variations. Our versatile instruments are perfect to be a part of your amazing surgical kit for years if cared for properly. These instruments are widely preferred to perform a broad range of surgical operations. We can also design our surgeons' custom kit packs according to their needs. The instruments are all fine, durable, and sharp.

FOX DERMAL CURETTE FEATURES**WI-70002646****Fox Dermal Curette 4mm Round Tip – 5.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Fox Dermal Curettes Dermatology 4mm Round Tip – 5.5\"/>

Fox Dermal Curette is commonly used in dermatology to scrape the skin growths such as excess skin surface, warts, and cancer. It features a ring-shaped cutting edge and a flat handle. It comes with multiple variations, so you can pick the best one for your surgical kit. It is crafted from high-grade surgical stainless steel.

Fox Dermal Curette is used for scraping skin lesions and growths such as warts. It is a single-ended instrument available with an oval or round loop tip. The loop is sharp on the inside and blunt on the outside so the lesion can be removed,

but the surrounding skin is not damaged. A variety of tip sizes are available to accommodate varying lesion sizes. The handle of the Fox curette is flat and lightweight.

We offer an extensive range of surgical curettes with multiple variations. Our versatile instruments are perfect to be a part of your amazing surgical kit for years if cared for properly. These instruments are widely preferred to perform a broad range of surgical operations. We can also design our surgeons' custom kit packs according to their needs. The instruments are all fine, durable, and sharp.

FOX DERMAL CURETTE FEATURES

**WI-70002645****Fox Dermal Curette 5mm Round Tip – 5.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Fox Dermal Curettes Dermatology 5mm Round Tip – 5.5\"/>

Fox Dermal Curette is commonly used in dermatology to scrape the skin growths such as excess skin surface, warts, and cancer. It features a ring-shaped cutting edge and a flat handle. It comes with multiple variations, so you can pick the best one for your surgical kit. It is crafted from high-grade surgical stainless steel.

Fox Dermal Curette is used for scraping skin lesions and growths such as warts. It is a single-ended instrument available with an oval or round loop tip. The loop is sharp on the inside and blunt on the outside so the lesion can be removed,

but the surrounding skin is not damaged. A variety of tip sizes are available to accommodate varying lesion sizes. The handle of the Fox curette is flat and lightweight.

We offer an extensive range of surgical curettes with multiple variations. Our versatile instruments are perfect to be a part of your amazing surgical kit for years if cared for properly. These instruments are widely preferred to perform a broad range of surgical operations. We can also design our surgeons' custom kit packs according to their needs. The instruments are all fine, durable, and sharp.

FOX DERMAL CURETTE FEATURES**WI-70002644****Fox Dermal Curette 6mm Round Tip – 5.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Fox Dermal Curettes Dermatology 6mm Round Tip – 5.5\"/>


Fox Dermal Curette is commonly used in dermatology to scrape the skin growths such as excess skin surface, warts, and cancer. It features a ring-shaped cutting edge and a flat handle. It comes with multiple variations, so you can pick the best one for your surgical kit. It is crafted from high-grade surgical stainless steel.

Fox Dermal Curette is used for scraping skin lesions and growths such as warts. It is a single-ended instrument available with an oval or round loop tip. The loop is sharp on the inside and blunt on the outside so the lesion can be removed,

but the surrounding skin is not damaged. A variety of tip sizes are available to accommodate varying lesion sizes. The handle of the Fox curette is flat and lightweight.

We offer an extensive range of surgical curettes with multiple variations. Our versatile instruments are perfect to be a part of your amazing surgical kit for years if cared for properly. These instruments are widely preferred to perform a broad range of surgical operations. We can also design our surgeons' custom kit packs according to their needs. The instruments are all fine, durable, and sharp.

FOX DERMAL CURETTE FEATURES

**WI-70000306****Gracey Curette 1/2 Blue Tips Double Ended****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Gracey Curette 1/2 Blue Coated Tips Double Ended with Blue and Silver Color Hollow Handle Stainless Steel Periodontal Dental Instruments.


Sharp and precise, dental gracey curettes are used to remove both supragingival and subgingival calculus in specific areas of the mouth. The Gracey 1/2 is designed for use on all aspects of anterior teeth (including canines). The 1/2 gracey curette has a small contra-angle and a return angle resulting in an instrument toe which is close to the mid-line of the instrument handle for better stability. The shank is particularly long on Gracey 1/2s to permit access to lower anteriors with the patient in a supine position.

1/2 GRACEY CURETTE BLUE TIPS FEATURES

- Used for Anterior Incisors & Canines
- Blue Color Coated Tips.
- Double-ended.
- Ergonomic Blue & Silver Color Hollow Handle.
- Hand Crafted From High Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Gracey Curette 1/2 Blue Tips with Blue & Silver Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curettes**APPLICATION:** Anterior Incisors and Canines**TIP COLOR:** Blue Color Coated**FINISH:** Polished**WI-70000319****Gracey Curette 1/2 Gold Tips Double Ended****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Gracey Curette 1/2 Gold Coated Tips Double Ended with Blue and Silver Hollow Handle Premium Quality Stainless Steel Periodontal Dental Instruments.

Sharp and precise, dental gracey curettes are used to remove both supragingival and subgingival calculus in specific areas of the mouth. The Gracey 1/2 Gold Tips is designed for use on all aspects of anterior teeth (including canines). The 1/2 gracey curette has a small contra-angle and a return angle resulting in an instrument toe which is close to the mid-line of the instrument handle for better stability. The shank is particularly long on Gracey 1/2s to permit access to lower anteriors with the patient in a supine position.

1/2 GRACEY CURETTE GOLD TIPS FEATURES

- Gold Color Coated Contra-angle Tips Double-ended
- Ergonomic Blue & Silver Color Hollow Handle.
- Hand Crafted From High Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Gracey Curette 1/2 Gold Tips with Blue & Silver Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curette**APPLICATION:** Anterior Incisors and Canines**TIP TYPE:** Contra-Angle**FINISH:** Polished**REUSABLE/AUTOCCLAVABLE:** Yes

**WI-70002090****Gracey Curette 1/2 Titanium Coated Double Ended****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Gracey Curette 1/2 Titanium Blue Coated Tip Double Ended Gold and Silver Color Hollow Handle Periodontal Dental Instruments

Sharp and precise, Gracey Curettes are used to remove both supragingival and subgingival calculus in specific areas of the mouth. The Gracey 1/2 is designed for use on all aspects of anterior teeth (including canines). The Gracey Curette has a small contra-angle and a return angle resulting in an instrument toe which is close to the mid-line of the instrument handle for better stability. The shank is particularly long on Gracey 1/2s to permit access to lower anteriors with the patient in a supine position.

1/2 GRACEY CURETTE FEATURES

- Titanium Blue coated contra-angle tips
- Double-ended
- Hand crafted from High Quality stainless steel.
- Ergonomic gold & silver color hollow handle.

PACKAGE INCLUDE:

- 1 Pc Gracey Curette 1/2 Titanium Blue Coated Double Ended with Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curette**APPLICATION:** Anterior Incisors and Canines**COLOR:** Titanium Color Coated Tip**TIP:** Double Ended**TIP TYPE:** Contra-Angle

GRACEY CURETTES 1/2

**WI-70002081****Gracey Curette 11/12 Blue Titanium Coated Tips Double Ended****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Gracey Curette 11/12 Titanium Blue Coated Tips Double Ended with Blue & Silver Color Hollow Handle Premium Quality Stainless Steel Periodontal Dental Instruments

Gracey Curette 11/12 is sharp and precise instrument used to remove both supragingival and subgingival calculus in specific areas of the mouth. Gracey 11/12 used for on mesial surfaces. The small contra angle and return angle are complemented by a third angle to offset the blade, which forms a smooth bladed toe. The curette tips are notably long, permitting access to the posterior teeth while keeping the clinician's hand comfortably out of the mouth.

11/12 BLUE TIPS GRACEY CURETTES FEATURES

- Titanium Blue Coated Contra-angle Tips.
- Double-ended Design.
- Ergonomic Blue & Silver Color Hollow Handle.
- Reusable.
- Hand Crafted from High Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Gracey Curette 11/12 Titanium Coated Blue Tips – Double Ended – Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curette**APPLICATION:** For Premolars' and Molars' Mesial Surfaces**COLOR:** Titanium Blue Coated Color Tips**MIRROR FINISH**

100% HIGH GRADE STAINLESS STEEL

**WI-70002076****Gracey Curette 13/14 Gold Coated Double Ended****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Gracey Curette 13/14 Gold Coated Tips Double Ended Titanium Blue & Silver Color Hollow Handle Periodontal Dental Instruments

The Gracey 13/14 Curette is an area-specific curette that is great for moderate to heavy calculus removal on the distal surfaces of the posterior teeth. Gracey curettes 13/14 are sharp and precise instrument used to remove both supragingival and subgingival calculus in specific areas of the mouth. The small contra angle and return angle are complemented by a third angle to offset the blade, which forms a smooth bladed toe. The curette tips are notably long, permitting access to the posterior teeth while keeping the clinician's hand comfortably out of the mouth.

13/14 GRACEY CURETTES FEATURES

- Goldcoated contra-angle tips
- Double-ended
- Hand crafted from High Quality stainless steel.
- Ergonomic Titanium blue & silver color Hollow Handle.

PACKAGE INCLUDE:

- 1 Pc Gracey Curette 13/14 Gold Coated Double Ended Titanium Blue & Silver Color Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curette**APPLICATION:** For Premolars' and Molars' Distal surfaces**COLOR:** Gold Color Coated Tips**TIP:** Double Ended**WI-70000307****Gracey Curette 3/4 Blue Tips Double Ended****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Gracey Curette 3/4 Blue Tips Double Ended with Blue and Silver Hollow Handle Premium Quality Stainless Steel Periodontal Dental Instruments.





Sharp and precise, dental gracey curettes are used to remove both supragingival and subgingival calculus in specific areas of the mouth. Gracey 3/4 blue tips is a dental instrument specifically designed for scaling and root planing procedures on anterior teeth, including incisors, canines, and premolars. The 3/4 gracey curette has slight contra-angled tips for anterior incisors and canines with mid-line easy grip hollow handle for better stability. All our composite/plastic filling instruments are highly polished stainless steel dental instruments is sterilizable and can be reused.


3/4 BLUE COATED TIPS GRACEY CURETTE FEATURES

- Used for Anterior Teeth (Incisors and Canines) and Premolars
- Blue Color Coated Slight Contra-angle Tips.
- Double-ended.
- Ergonomic Blue & Silver Color Hollow Handle.
- Reusable.
- Hand Crafted From High Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Gracey Curette 3/4 Blue Coated Tips with Blue & Silver Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Curette

**WI-70000295****Gracey Curette 5/6 Blue Tips Double Ended****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Gracey Curette 5/6 Blue Coated Tips Double Ended with Blue and Silver Color Hollow Handle Stainless Steel Periodontal Dental Instruments.

Sharp and precise blue tips Gracey Curette 5/6 is a specialized periodontal tool designed for the removal of subgingival calculus on the anterior teeth and premolars, particularly effective on their proximal surfaces. Gracey 5/6 enabling precise access to difficult areas for early periodontal disease treatment and root planing. Features a slight contra-angle and a unique blade curvature for precise adaptation to specific tooth surfaces, minimizing tissue distention. Its ergonomically designed lightweight rounded blue and silver hollow handle



reduces hand fatigue during dental procedures.

5/6 GRACEY CURETTE BLUE TIPS FEATURES

- Used for Anterior Teeth and Premolars.
- Blue Color Coated Tips.
- Double-ended.
- Ergonomic Blue & Silver Color Hollow Handle.
- Reusable.
- Hand Crafted From High Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Gracey Curette 5/6 Blue Tips with Blue & Silver Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**WI-70000294****Gracey Curette 7/8 Blue Tips Double Ended****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Gracey Curette 7/8 Blue Coated Tips Double Ended with Blue and Silver Color Hollow Handle Stainless Steel Periodontal Dental Instruments.

Gracey curette 7/8 is a specialized dental instrument designed for hygienists and dentists, effective on the buccal and lingual surfaces of premolars and molars, where standard curettes may be less effective. Its offset blade adapts accurately to tooth surfaces, with a cutting edge angled towards the root. The instrument's uniquely bent shank enhances access to curved surfaces, while its rounded toe and back minimize gum tissue trauma during use. The ergonomically designed lightweight rounded blue and

silver hollow handle reduces hand fatigue during dental procedures.


TITANIUM BLUE TIPS 7/8 GRACEY CURETTE FEATURES

- Used for Buccal and Lingual (Inner) Surfaces of Premolars and Molars.
- Blue Color Coated Tips.
- Double-ended.
- Ergonomic Blue & Silver Color Hollow Handle.
- Reusable.
- Hand Crafted From High Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Gracey Curette 7/8 Blue Tips with Blue & Silver Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel

**WI-70000305**

Gracey Curette 9/10 Blue Coated Tips Double Ended

\$10.99 per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Gracey Curette 9/10 Blue Color Coated Tips Double Ended with Blue and Silver Hollow Handle Periodontal Dental Instruments

Gracey curettes 9/10 are Sharp and precise. The 9/10 Gracey is predominantly used on the buccal and lingual molar aspects, especially around the bifurcation, where access can be difficult. Its wide contra-angle, wide return angle make it distinctive instrument. The toe is very nearly in line with the instrument handle, making this pattern especially suited to the least accessible root surfaces, while still allowing good control.

9/10 GRACEY CURETTE FEATURES

- Used for Buccal and Lingual Molars

- Blue Coated Wide Contra-angle Tips.
- Double-ended.
- Hand Crafted from High Quality Stainless Steel.
- Ergonomic Blue & Silver Color Hollow Handle.

PACKAGE INCLUDE:

- 1 Pc Gracey Curette 9/10 Blue Tips with Blue & Silver Color Hollow Handle

TECHNICAL SPECIFICATIONS

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Dental Curettes**APPLICATION:** For Molars' Buccal and Lingual Surfaces**TIPS COLOR:** Blue Color Coated Tip**FINISH:** Polished**REUSABLE/AUTOCCLAVABLE:** Yes**STERILE:** Non-Sterile**WI-70002085**

Gracey Curette 9/10 Gold Coated Double Ended

\$10.99 per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Gracey Curette 9/10 Gold Color Coated Tips Double Ended Blue and Silver Hollow Handle Periodontal Dental Instruments

Gracey curettes 9/10 are Sharp and precise. The 9/10 Gracey is predominantly used on the buccal and lingual molar aspects, especially around the bifurcation, where access can be difficult. Its wide contra-angle, wide return angle make it distinctive instrument. The toe is very nearly in line with the instrument handle, making this pattern especially suited to the least accessible root surfaces, while still allowing good control.

9/10 GRACEY CURETTE FEATURES


- Gold coated wide contra-angle tips

- Double-ended
- Hand crafted from High Quality stainless steel.
- Ergonomic blue & silver color hollow handle.

PACKAGE INCLUDE:

- 1 Piece Gracey Curette 9/10/2 Gold Coated Double Ended Blue & Silver Color Hollow Handle

TECHNICAL SPECIFICATIONS

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Dental Curette**APPLICATION:** For Molars' Buccal and Lingual Surfaces**COLOR:** Gold Color Coated Tip**TIP:** Double Ended**TIP TYPE:** Wide Contra-Angle**HANDLE:** Blue & Silver Color Hollow Handle

**WI-70002394**

Heaney Uterine Biopsy Curette 9" Serrated Loop

\$10.99 per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heaney Uterine Biopsy Tissue Curette 9" (22.9 cm) Open Oval Loop Serrated Tip 5.0mm Premium Quality Stainless Steel OB/Gynecology Surgical Instruments

Heaney Uterine Biopsy Curette is a gynecological tool used to scrape and collect uterine lining tissue samples. Its key uses include endometrial biopsies for detecting abnormal cells and diagnosing causes of abnormal uterine bleeding and infertility. Additionally, it is utilized in dilation and curettage (D&C) procedures. Heaney curette features a serrated loop tip for efficient tissue cutting, a malleable shaft for adaptation to uterine angles, and a flat handle for a secure grip during use. Made with high-quality Stainless Steel and thus can be reused after sterilization.



HEANEY UTERINE BIOPSY CURETTES FEATURES:

- 9 Inch (22.9cm) Overall Length.
- Lightweight long handle with serrated tip built for deep access.
- 5.0mm angled open tip serrated loop.
- Blunt tip to prevent piercing injuries.
- Solid flat handle for comfortable handling.
- Reusable.
- Made with premium quality stainless steel.

PACKAGE INCLUDE:

- 1 Pc Heaney Uterine Biopsy Curette 9" (22.9 cm) Serrated Loop.

TECHNICAL INFORMATION

**WI-70002252**

Kevorkian-Youngue Biopsy Curette 12" Loop 3X12mm with Basket Half Gold Handle

\$12.99 per pieceTrade price from **\$11.95**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kevorkian-Youngue Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, with Basket Half Gold Handle, Angular Premium Quality Stainless Steel Obstetrical/Gynecology Surgical Instruments

Kevorkian Biopsy Curette is a specialized gynecological instrument used for dissecting and removing part of the uterine wall for diagnostic purposes. Kevorkian-Youngue Endocervical Biopsy Curette is designed for the excision of tissue from the endocervix for biopsy. The curette is available with or without a basket for use in excision or biopsy. The Kevorkian Biopsy curette also comes with sharp edges that effectively scrape off the uterine tissue and the shape allows the surgeon to easily

scoop up and collect the dissected tissue for secure removal.

KEVORKIAN-YOUNGUE ENDOCERVICAL BIOPSY CURETTE FEATURES

- Long shank for deep tissue access.
- Ergonomic solid half gold handle for superior control.
- With Or Without basket for optimum dissection.
- Smooth outer curves to prevent inadvertent soft tissue damage.
- Sharp edges for effective scraping off the uterine tissues.

PACKAGE INCLUDE:

- 1 Pc Kevorkian-Youngue Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, with Basket Half Gold Handle



**WI-70002251****Kevorkian-Youngue Biopsy Curette 12" Loop 3X12mm Without Basket Half Gold Handle****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kevorkian-Youngue Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, without Basket – Half Gold Handle – Angular Premium Quality Stainless Steel Obstetrical/Gynecology Surgical Instruments

Kevorkian Biopsy Curette is a specialized gynecological instrument used for dissecting and removing part of the uterine wall for diagnostic purposes. Kevorkian-Youngue Endocervical Biopsy Curette is designed for the excision of tissue from the endocervix for biopsy. The curette is available with or without a basket for use in excision or biopsy. The Kevorkian Biopsy curette also comes with sharp edges that effectively scrape off the uterine tissue and the shape allows the surgeon to easily



scoop up and collect the dissected tissue for secure removal.

**KEVORKIAN-YOUNGUE
ENDOCERVICAL BIOPSY
CURETTE FEATURES**

- Long shank for deep tissue access.
- Ergonomic solid Half Gold Handle for superior control.
- With Or Without basket for optimum dissection.
- Smooth outer curves to prevent inadvertent soft tissue damage.
- Sharp edges for effective scraping off the uterine tissues.

PACKAGE INCLUDE:

- 1 Pc Kevorkian-Youngue Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, without Basket – Half Gold Handle

**WI-70002393****Kevorkian-Youngue Endocervical Biopsy Curette 12" Loop 3X12mm with Basket****\$11.99** per pieceTrade price from **\$11.03**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kevorkian-Youngue Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, with Basket, Angular Premium Quality Stainless Steel Obstetrical/Gynecology Surgical Instruments

Kevorkian Biopsy Curette is a specialized gynecological instrument used for dissecting and removing part of the uterine wall for diagnostic purposes. Kevorkian-Youngue Endocervical Biopsy Curette is designed for the excision of tissue from the endocervix for biopsy. The curette is available with or without a basket for use in excision or biopsy. The Kevorkian Biopsy curette also comes with sharp edges that effectively scrape off the uterine tissue and the shape allows the surgeon to easily


scoop up and collect the dissected tissue for secure removal.


**KEVORKIAN YOUNG BIOPSY
CURETTE FEATURES**

- Long shank for deep tissue access.
- Ergonomic solid handle for superior control.
- With Or Without basket for optimum dissection.
- Smooth outer curves to prevent inadvertent soft tissue damage.
- Sharp edges for effective scraping off the uterine tissues.

PACKAGE INCLUDE:

- 1 Pc Kevorkian-Youngue Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, with Basket

TECHNICAL INFORMATION

**WI-70002392****Kevorkian-Young Endocervical Biopsy Curette 12" Loop 3X12mm Without Basket****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kevorkian-Young Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, without Basket, Angular Premium Quality Stainless Steel Obstetrical/Gynecology Surgical Instruments

Kevorkian Biopsy Curette is a specialized gynecological instrument used for dissecting and removing part of the uterine wall for diagnostic purposes. Kevorkian-Young Endocervical Biopsy Curette is designed for the excision of tissue from the endocervix for biopsy. The curette is available with or without a basket for use in excision or biopsy. The Kevorkian Biopsy curette also comes with sharp edges that effectively scrape off the uterine tissue and the shape allows the surgeon to easily



scoop up and collect the dissected tissue for secure removal.

KEVORKIAN YOUNG BIOPSY CURETTE FEATURES

- Long shank for deep tissue access.
- Ergonomic solid handle for superior control.
- With Or Without basket for optimum dissection.
- Smooth outer curves to prevent inadvertent soft tissue damage.
- Sharp edges for effective scraping off the uterine tissues.

PACKAGE INCLUDE:

- 1 Pc Kevorkian-Young Endocervical Biopsy Curette 12" (30.5 cm), Loop Size 3 X 12 mm, without Basket

TECHNICAL INFORMATION**WI-70001915****Lambotte Osteotome 9" Straight 13mm Tip****\$20.99** per pieceTrade price from **\$19.31**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lambotte Osteotome 9" (22.9cm) Straight 1/2" (13mm) Tip Flat Solid Handle Premium Quality Stainless Steel Orthopedic Surgical Instruments.


Lambotte Osteotome is a specialized surgical tool that orthopedic surgeons can use to cut, sculpture and reshaped hard compact bone, especially during osteotomy and fracture repairing. Lambotte Bone Osteotomes have a beveled tip equipped with a sharp cutting edge. Moreover, the working end features straight and curved patterns. As a result, surgeons can reshape flat and round bone surfaces. This handheld orthopedic tool includes a solid, flat handle supplied with a smooth working surface. This particular design ensures enhanced gripping and prevents finger strain. All our Lambotte

Bone Osteotomes are made with premium grade stainless steel, granting the highest quality and performance.

LAMBOTTE OSTEOTOME FEATURES

- Size: 9 Inch (22.9cm) Long.
- Atraumatic Beveled Sharp Working End for Minimizing Local Injury.
- Straight Tip Designs for Accessing Narrow Spaces.
- Ergonomic Solid Handle For Superior Control.
- Reusable, Long-Lasting, Durable
- Suitable For Various Osteotomy Procedures.
- Made with German grade Stainless Steel

PACKAGE INCLUDE:

**WI-70001914****Lambotte Osteotome 9" Straight 19mm Tip****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Lambotte Osteotome 9" (22.9cm) Straight 3/4" (19mm) Tip Flat Solid Handle Premium Quality Stainless Steel Orthopedic Surgical Instruments.

Lambotte Osteotome is a specialized surgical tool that orthopedic surgeons can use to cut, sculpture and reshaped hard compact bone, especially during osteotomy and fracture repairing. Lambotte Bone Osteotomes have a beveled tip equipped with a sharp cutting edge. Moreover, the working end features straight and curved patterns. As a result, surgeons can reshape flat and round bone surfaces. This handheld orthopedic tool includes a solid, flat handle supplied with a smooth working surface. This particular design ensures enhanced gripping and prevents finger strain. All our Lambotte

Bone Osteotomes are made with premium grade stainless steel, granting the highest quality and performance.

LAMBOTTE OSTEOTOME FEATURES

- Size: 9 Inch (22.9cm) Long.
- Atraumatic Beveled Sharp Working End for Minimizing Local Injury.
- Straight Tip Designs for Accessing Narrow Spaces.
- Ergonomic Solid Handle For Superior Control.
- Reusable, Long-Lasting, Durable
- Suitable For Various Osteotomy Procedures.
- Made with German grade Stainless Steel

PACKAGE INCLUDE:**WI-70001913****Lambotte Osteotome 9" Straight 25mm Tip****\$22.99** per pieceTrade price from **\$21.15**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Lambotte Osteotome 9" (22.9cm) Straight 1" (25mm) Tip Flat Solid Handle Premium Quality Stainless Steel Orthopedic Surgical Instruments.


Lambotte Osteotome is a specialized surgical tool that orthopedic surgeons can use to cut, sculpture and reshaped hard compact bone, especially during osteotomy and fracture repairing. Lambotte Bone Osteotomes have a beveled tip equipped with a sharp cutting edge. Moreover, the working end features straight and curved patterns. As a result, surgeons can reshape flat and round bone surfaces. This handheld orthopedic tool includes a solid, flat handle supplied with a smooth working surface. This particular design ensures enhanced gripping and prevents finger strain. All our Lambotte

Bone Osteotomes are made with premium grade stainless steel, granting the highest quality and performance.

LAMBOTTE OSTEOTOME FEATURES

- Size: 9 Inch (22.9cm) Long.
- Atraumatic Beveled Sharp Working End for Minimizing Local Injury.
- Straight Tip Designs for Accessing Narrow Spaces.
- Ergonomic Solid Handle For Superior Control.
- Reusable, Long-Lasting, Durable
- Suitable For Various Osteotomy Procedures.
- Made with German grade Stainless Steel

PACKAGE INCLUDE:

**WI-70001912****Lambotte Osteotome 9" Straight 32mm Tip****\$23.99** per pieceTrade price from **\$22.07**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Lambotte Osteotome 9" (22.9cm) Straight 1 1/4" (32mm) Tip Flat Solid Handle Premium Quality Stainless Steel Orthopedic Surgical Instruments.

Lambotte Osteotome is a specialized surgical tool that orthopedic surgeons can use to cut, sculpture and reshaped hard compact bone, especially during osteotomy and fracture repairing. Lambotte Bone Osteotomes have a beveled tip equipped with a sharp cutting edge. Moreover, the working end features straight and curved patterns. As a result, surgeons can reshape flat and round bone surfaces. This handheld orthopedic tool includes a solid, flat handle supplied with a smooth working surface. This particular design ensures enhanced gripping and prevents finger strain. All our Lambotte

Bone Osteotomes are made with premium grade stainless steel, granting the highest quality and performance.

LAMBOTTE OSTEOTOME FEATURES

- Size: 9 Inch (22.9cm) Long.
- Atraumatic Beveled Sharp Working End for Minimizing Local Injury.
- Straight Tip Designs for Accessing Narrow Spaces.
- Ergonomic Solid Handle For Superior Control.
- Reusable, Long-Lasting, Durable
- Suitable For Various Osteotomy Procedures.
- Made with German grade Stainless Steel

PACKAGE INCLUDE:**WI-70000202****Lambotte Osteotome 9" Straight 38mm Tip****\$24.99** per pieceTrade price from **\$22.99**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lambotte Osteotome 9" (22.9cm) Straight 1 1/2" (38mm) Tip Flat Solid Handle Premium Quality Stainless Steel Orthopedic Surgical Instruments.


Lambotte Osteotome is a specialized surgical tool that orthopedic surgeons can use to cut, sculpture and reshaped hard compact bone, especially during osteotomy and fracture repairing. Lambotte Bone Osteotomes have a beveled tip equipped with a sharp cutting edge. Moreover, the working end features straight and curved patterns. As a result, surgeons can reshape flat and round bone surfaces. This handheld orthopedic tool includes a solid, flat handle supplied with a smooth working surface. This particular design ensures enhanced gripping and prevents finger strain. All our Lambotte

Bone Osteotomes are made with premium grade stainless steel, granting the highest quality and performance.

LAMBOTTE OSTEOTOME FEATURES

- Size: 9 Inch (22.9cm) Long.
- Atraumatic Beveled Sharp Working End for Minimizing Local Injury.
- Straight Tip Designs for Accessing Narrow Spaces.
- Ergonomic Solid Handle For Superior Control.
- Reusable, Long-Lasting, Durable
- Suitable For Various Osteotomy Procedures.
- Made with German grade Stainless Steel

PACKAGE INCLUDE:



WI-70001916

Lambotte Osteotome 9" Straight 6mm Tip

\$19.99 per piece

Trade price from \$18.39
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lambotte Osteotome 9" (22.9cm) Straight 1/4" (6mm) Tip Flat Solid Handle Premium Quality Stainless Steel Orthopedic Surgical Instruments.

Lambotte Osteotome is a specialized surgical tool that orthopedic surgeons can use to cut, sculpture and reshaped hard compact bone, especially during osteotomy and fracture repairing. Lambotte Bone Osteotomes have a beveled tip equipped with a sharp cutting edge. Moreover, the working end features straight and curved patterns. As a result, surgeons can reshape flat and round bone surfaces. This handheld orthopedic tool includes a solid, flat handle supplied with a smooth working surface. This particular design ensures enhanced gripping and prevents finger strain. All our Lambotte

Bone Osteotomes are made with premium grade stainless steel, granting the highest quality and performance.

LAMBOTTE OSTEOTOME FEATURES

- Size: 9 Inch (22.9cm) Long.
- Atraumatic Beveled Sharp Working End for Minimizing Local Injury.
- Straight Tip Designs for Accessing Narrow Spaces.
- Ergonomic Solid Handle For Superior Control.
- Reusable, Long-Lasting, Durable
- Suitable For Various Osteotomy Procedures.
- Made with German grade Stainless Steel

PACKAGE INCLUDE:



WI-70003273

Lexer Mini Chisel

\$33.59 per piece

Trade price from \$30.90
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Lexer Mini Chisel is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$33.59, falling to \$30.90 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003838****Lucas Bone Curette #84 Small 2.4mm Spoon Double Ended****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lucas Bone Curette #84 small 2.4mm Double Ended Smooth Spoon Shaped Tips with Round Pattern Grip Handle Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 2.4 mm wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium


quality stainless steel, it ensures durability and lasting performance. These curettes are reusable after sterilization.

LUCAS CURETTE #84 FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 2.4mm Spoon Shaped Sharp Ends.
- Double-ended.
- Smooth Pattern.
- Non-Slip Strong Grip Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Lucas Bone Curette #84 Small 2.4mm Double Ended

**WI-70003837****Lucas Bone Curette #85 Double Ended 2.7mm Smooth Spoon Tips****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lucas Bone Curette #85 Medium 2.7mm Double Ended Smooth Spoon Shaped Tips with Round Pattern Grip Handle Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Dental Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 2.7 mm wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium



quality stainless steel, it ensures durability and lasting performance. These curettes are reusable after sterilization.

LUCAS CURETTE #85 FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 2.7mm Spoon Shaped Sharp Ends.
- Double-ended.
- Smooth Pattern.
- Non-Slip Strong Grip Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Lucas Bone Curette #85 Medium 2.7mm Double Ended

**WI-70000609****Lucas Bone Curette #85 Double Ended 2.7mm Tips Hexagonal Handle****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lucas Bone Curette #85 Medium 2.7mm Double Ended Smooth Spoon Shaped Tips with Hexagonal Handle Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Dental Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 2.7 mm wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium



quality stainless steel, it ensures durability and lasting performance. These curettes are reusable after sterilization.

LUCAS CURETTE #85 FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 2.7mm Spoon Shaped Sharp Ends.
- Double-ended.
- Smooth Pattern.
- Non-Slip Strong Hexagonal Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Lucas Bone Curette #85 Medium 2.7mm Double Ended

**WI-70003836****Lucas Bone Curette #86 Double Ended 3mm Smooth Spoon Tips****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lucas Bone Curette #86 Large 3mm Double Ended Smooth Spoon Shaped Tips with Round Pattern Knurled Grip Handle Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Dental Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 3 mm wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium


quality stainless steel, it ensures durability and lasting performance. These curettes are reusable after sterilization.

LUCAS CURETTE #86 FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 3mm Spoon Shaped Sharp Ends.
- Double-ended.
- Smooth Pattern.
- Non-Slip Knurled Grip Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Lucas Bone Curette #86 Large 3mm Double Ended

**WI-70000608****Lucas Bone Curette #86 Double Ended 3mm Tips Hexagonal Handle****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lucas Bone Curette #86 Large 3mm Double Ended Smooth Spoon Shaped Tips with Hexagonal Handle Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Dental Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 3 mm wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium



quality stainless steel, it ensures durability and lasting performance. These curettes are reusable after sterilization.

LUCAS CURETTE #86 FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 3mm Spoon Shaped Sharp Ends.
- Double-ended.
- Smooth Pattern.
- Non-Slip Strong Hexagonal Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Lucas Bone Curette #86 Large 3mm Double Ended

**WI-70000072****Lucas Bone Curette #86 Double Ended Spoon Tips - Hollow****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lucas Bone Curette #86 Large 3mm Double Ended Smooth Spoon Shaped Tips with Round Pattern Hollow Handle Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Dental Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 3 mm wide spoon-shaped blades that have an elongated radius. The ergonomic hollow handle provides comfort for the dentist with an easy hold non-slip grip allow less hand

fatigue during long procedures. Made with premium quality stainless steel, it ensures durability and lasting performance. These curettes are reusable after sterilization.

LUCAS CURETTE #86 FEATURES

- 3mm Spoon Shaped Sharp Ends.
- Double-ended.
- Smooth Pattern.
- Lightweight Non-Slip Hollow Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Lucas Bone Curette #86 Hollow Handle

TECHNICAL INFORMATION

**WI-70003840****Lucas Bone Curette #87 Double Ended 3.8mm Smooth Spoon Tips****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lucas Bone Curette #87 Double Ended 3.8mm Smooth Jumbo Spoon Shaped Tips with Round Pattern Knurled Grip Handle Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Dental Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 3.8 mm wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium


quality stainless steel, it ensures durability and lasting performance. These curettes are reusable after sterilization.

LUCAS CURETTE #87 FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 3.8mm Spoon Shaped Sharp Ends.
- Double-ended.
- Smooth Pattern.
- Non-Slip Knurled Grip Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Lucas Bone Curette #87 Jumbo 3.8mm Double Ended

**WI-70003839****Lucas Bone Curette #88 Double Ended 5mm Smooth Spoon Tips****\$7.99** per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lucas Bone Curette #88 Double Ended 5mm Smooth Spoon Shaped Tips with Round Pattern Knurled Grip Handle Premium Quality Stainless Steel Surgical Curettes Dentist Lab Dental Teeth Instruments.

Dental Lucas curette used to remove unwanted tissue, cysts or other debris from the tooth socket. Lucas surgical curettes are designed for curettage, cyst removal, tooth socket debridement and removing soft tissue from bony defects. Double-ended Lucas bone curettes have spoon-shaped mirror image ends. The terminal shank is angled at 50° with 5 mm wide spoon-shaped blades that have an elongated radius. The ergonomic handle provides comfort for the dentist with an easy hold non-slip grip. Made with premium


quality stainless steel, it ensures durability and lasting performance. These curettes are reusable after sterilization.


LUCAS CURETTE #88 FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 5mm Spoon Shaped Sharp Ends.
- Double-ended.
- Smooth Pattern.
- Non-Slip Knurled Grip Handle.
- Easy to Hold.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Lucas Bone Curette #88 Double Ended 5mm Spoon Tips



**WI-70000301****McCall's Curette 13S/14S Blue Coated Tips with Hollow Handle****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

McCall's Curette #13S/14S Sharp Blue Coated Tips Double Ended with Blue & Silver Hollow Handle Premium Quality Stainless Steel Periodontal Dental Scalers.

Universal McCall's #13S/14S curette is a dental tool for removing moderate to heavy calculus from supra- and subgingival tooth surfaces. McCall's 13S-14S has sharpened point which allows effective cleaning of deep root pockets, particularly in posterior teeth, making it ideal for tight, narrow areas. 13S/14S pointed McCall's Curette can also be used for removing excess filling material or cement. The blue titanium-coated tips are rust-resistant, and the ergonomically rounded hollow handles provide a lightweight, comfortable grip. Made of highly polished stainless steel,

these instruments are sterilizable and reusable.

MCCALLS 13S-14S CURETTE FEATURES:

- Sharp Pointed Cutting Edges.
- Double Ended.
- Titanium Coated Blue Tips.
- Slightly Thinner and Shorten Blade than 13/14 McCall.
- Ergonomic Non-slip Blue & Silver Hollow Handle.
- Reusable.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc McCall's Curette #13S/14S Blue Tips – Hollow Handle

**WI-70002416****Meyhoefer Chalazion Curette 1.5mm Size #0 Flat Handle 5"****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Meyhoefer Chalazion Curettes 5" (12.7cm) Size #0 Cupped 1.5mm Solid Handle Premium Quality Stainless Steel Surgical Medical Ophthalmic ENT Instruments.

Meyhoefer Chalazion Curette is an ophthalmic instrument used for removing chalazion. A small cyst formed by inflammation of a gland in the margin of the eyelid, a chalazion is similar to a sty but often requires an incision and drainage to be properly treated. Meyhoefer curette features a sharp, rounded scoop tip in various sizes, often used with Chalazion Forceps. Made from high-quality stainless steel, these instruments resist corrosion and rust, allowing for reuse after sterilization.

MEYHOEFER CHALAZION CURETTE FEARES

- Size: 5 Inch (12.7cm) Overall Length.
- Multiple sizes #00 (1mm) to #4(3.5mm)
- Cup-shaped scoop Tip with Semi-sharp Edges for Enhanced Scraping and Scooping.
- Flat Solid Handle Improves Finger Grip and Reduce Hand Fatigue.
- Primarily for Extracting Cysts, Lid Masses and Tissue Samples.
- Useful for Pushing Small Tumors Down.
- Fully Autoclavable and Safe for Repeated Clinical Use
- Made From High Quality Stainless Steel

PACKAGE INCLUDE:

**WI-70002415****Meyhoefer Chalazion Curette 1.7mm Size #1 Flat Handle 5"****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Meyhoefer Chalazion Curettes 5" (12.7cm) Size #1 Cupped 1.7mm Solid Handle Premium Quality Stainless Steel Surgical Medical Ophthalmic ENT Instruments.

Meyhoefer Chalazion Curette is an ophthalmic instrument used for removing chalazion. A small cyst formed by inflammation of a gland in the margin of the eyelid, a chalazion is similar to a sty but often requires an incision and drainage to be properly treated. Meyhoefer curette features a sharp, rounded scoop tip in various sizes, often used with Chalazion Forceps. Made from high-quality stainless steel, these instruments resist corrosion and rust, allowing for reuse after sterilization.

MEYHOEFER CHALAZION CURETTE FEARES

- Size: 5 Inch (12.7cm) Overall Length.
- Multiple sizes #00 (1mm) to #4(3.5mm)
- Cup-shaped scoop Tip with Semi-sharp Edges for Enhanced Scraping and Scooping.
- Flat Solid Handle Improves Finger Grip and Reduce Hand Fatigue.
- Primarily for Extracting Cysts, Lid Masses and Tissue Samples.
- Useful for Pushing Small Tumors Down.
- Fully Autoclavable and Safe for Repeated Clinical Use
- Made From High Quality Stainless Steel

PACKAGE INCLUDE:**WI-70002417****Meyhoefer Chalazion Curette 1mm Size #00 Flat Handle 5"****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Meyhoefer Chalazion Curettes 5" (12.7cm) Size #00 Cupped 1mm Solid Handle Premium Quality Stainless Steel Surgical Medical Ophthalmic ENT Instruments.

Meyhoefer Chalazion Curette is an ophthalmic instrument used for removing chalazion. A small cyst formed by inflammation of a gland in the margin of the eyelid, a chalazion is similar to a sty but often requires an incision and drainage to be properly treated. Meyhoefer curette features a sharp, rounded scoop tip in various sizes, often used with Chalazion Forceps. Made from high-quality stainless steel, these instruments resist corrosion and rust, allowing for reuse after sterilization.

MEYHOEFER CHALAZION CURETTE FEARES

- Size: 5 Inch (12.7cm) Overall Length.
- Multiple sizes #00 (1mm) to #4(3.5mm)
- Cup-shaped scoop Tip with Semi-sharp Edges for Enhanced Scraping and Scooping.
- Flat Solid Handle Improves Finger Grip and Reduce Hand Fatigue.
- Primarily for Extracting Cysts, Lid Masses and Tissue Samples.
- Useful for Pushing Small Tumors Down.
- Fully Autoclavable and Safe for Repeated Clinical Use
- Made From High Quality Stainless Steel

PACKAGE INCLUDE:

**WI-70002414****Meyhoefer Chalazion Curette 2.2mm Size #2 Flat Handle 5"****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Meyhoefer Chalazion Curettes 5" (12.7cm) Size #2 Cupped 2.2mm Solid Handle Premium Quality Stainless Steel Surgical Medical Ophthalmic ENT Instruments.

Meyhoefer Chalazion Curette is an ophthalmic instrument used for removing chalazion. A small cyst formed by inflammation of a gland in the margin of the eyelid, a chalazion is similar to a sty but often requires an incision and drainage to be properly treated. Meyhoefer curette features a sharp, rounded scoop tip in various sizes, often used with Chalazion Forceps. Made from high-quality stainless steel, these instruments resist corrosion and rust, allowing for reuse after sterilization.

MEYHOEFER CHALAZION CURETTE FEARES

- Size: 5 Inch (12.7cm) Overall Length.
- Multiple sizes #00 (1mm) to #4(3.5mm)
- Cup-shaped scoop Tip with Semi-sharp Edges for Enhanced Scraping and Scooping.
- Flat Solid Handle Improves Finger Grip and Reduce Hand Fatigue.
- Primarily for Extracting Cysts, Lid Masses and Tissue Samples.
- Useful for Pushing Small Tumors Down.
- Fully Autoclavable and Safe for Repeated Clinical Use
- Made From High Quality Stainless Steel

PACKAGE INCLUDE:**WI-70002413****Meyhoefer Chalazion Curette 2.5mm Size #3 Flat Handle 5"****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Meyhoefer Chalazion Curettes 5" (12.7cm) Size #3 Cupped 2.5mm Solid Handle Premium Quality Stainless Steel Surgical Medical Ophthalmic ENT Instruments.

Meyhoefer Chalazion Curette is an ophthalmic instrument used for removing chalazion. A small cyst formed by inflammation of a gland in the margin of the eyelid, a chalazion is similar to a sty but often requires an incision and drainage to be properly treated. Meyhoefer curette features a sharp, rounded scoop tip in various sizes, often used with Chalazion Forceps. Made from high-quality stainless steel, these instruments resist corrosion and rust, allowing for reuse after sterilization.

MEYHOEFER CHALAZION CURETTE FEARES

- Size: 5 Inch (12.7cm) Overall Length.
- Multiple sizes #00 (1mm) to #4(3.5mm)
- Cup-shaped scoop Tip with Semi-sharp Edges for Enhanced Scraping and Scooping.
- Flat Solid Handle Improves Finger Grip and Reduce Hand Fatigue.
- Primarily for Extracting Cysts, Lid Masses and Tissue Samples.
- Useful for Pushing Small Tumors Down.
- Fully Autoclavable and Safe for Repeated Clinical Use
- Made From High Quality Stainless Steel

PACKAGE INCLUDE:

**WI-70002412****Meyhoefer Chalazion Curette 3.5mm Size #4 Flat Handle 5"****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Meyhoefer Chalazion Curettes 5" (12.7cm) Size #4 Cupped 3.5mm Solid Handle Premium Quality Stainless Steel Surgical Medical Ophthalmic ENT Instruments.

Meyhoefer Chalazion Curette is an ophthalmic instrument used for removing chalazion. A small cyst formed by inflammation of a gland in the margin of the eyelid, a chalazion is similar to a sty but often requires an incision and drainage to be properly treated. Meyhoefer curette features a sharp, rounded scoop tip in various sizes, often used with Chalazion Forceps. Made from high-quality stainless steel, these instruments resist corrosion and rust, allowing for reuse after sterilization.

MEYHOEFER CHALAZION CURETTE FEARES

- Size: 5 Inch (12.7cm) Overall Length.
- Multiple sizes #00 (1mm) to #4(3.5mm)
- Cup-shaped scoop Tip with Semi-sharp Edges for Enhanced Scraping and Scooping.
- Flat Solid Handle Improves Finger Grip and Reduce Hand Fatigue.
- Primarily for Extracting Cysts, Lid Masses and Tissue Samples.
- Useful for Pushing Small Tumors Down.
- Fully Autoclavable and Safe for Repeated Clinical Use
- Made From High Quality Stainless Steel

PACKAGE INCLUDE:**WI-70000228****Molt Periosteal #M9 Elevator Blue Tips Hollow Handle****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Molt 9 Periosteal Elevator (M9) Blue Coated Tips Double Ended with Blue & Silver Color Hollow Handle German Grade Stainless Steel Extractions and Impant Surgery Dental Surgical Veterinary Procedures

Molt Periosteal Elevator is an oral instrument commonly used to retract soft tissue in dental procedures. It has a sharp, pointed, narrow end that is used to begin the process of lifting the mucoperiosteum and a slightly wider, more rounded end that completes the tissue removal enough for bone exposure. The rounded edge of the Molt 9 periosteal elevator minimise the tearing of tissue when lifting the gum flap. The sharp end makes the initial lifting of the flap easier. The blue, titanium-coated tips improve visibility

and durability, while the hollow handle offers an ergonomic, non-slip grip to reduce hand fatigue.

BLUE TIPS PERIOSTEAL MOLT #9 (M9) FEATURES

- Double-ended Design.
- Titanium Blue Coated Tips.
- 1 Broad, Rounded Wider Blade, 1 Sharp, Pointed and Narrow Blade.
- Ergonomic Lightweight Blue & Silver Color Hollow Handle
- Corrosion Resistance.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Molt 9 Periosteal Elevator Blue Double Ended – Hollow Handle

**WI-70003835**

Molt Periosteal M9 Dental Surgery Tooth Elevator

\$7.99 per pieceTrade price from **\$7.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Periosteal Molt #9 (M9) Elevator
German Grade Stainless Steel
Extraction Surgery Surgical Veterinary
Dental Instruments

Molt Periosteal Elevator is an oral instrument commonly used to retract soft tissue in dental procedures. It has a sharp, pointed, narrow end that is used to begin the process of lifting the mucoperiosteum and a slightly wider, more rounded end that completes the tissue removal enough for bone exposure. The rounded edge of the Molt 9 periosteal elevator minimise the tearing of tissue when lifting the gum flap. The sharp end makes the initial lifting of the flap easier.

PERIOSTEAL MOLT #9 FEATURES

- Double-ended

- 1 Wide sharp blade
- 1 sharp, pointed and narrow blade
- Non-reflective surface
- Ergonomic lightweight handle
- Handcrafted with German grade stainless steel
- Corrosion resistance.

PACKAGE INCLUDE:

- 1 Periosteal Elevator Molt 9 Double Ended 7.5" (19cm) Long

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Dental Instruments**APPLICATION:** Teeth Extraction**RUST PREVENTION PROCEDURE:**

Passivated

SURFACE: Polished**AUTOCLAVED:** Yes**WI-70002398**

Novak Endometrial Suction Biopsy Curette 2mm 9 3/4" with Luer Hub

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Endometrial Biopsy Curette NOVAK 9 3/4" (24.8cm) Length Single Ended, 2mm Tip Curved Round Tip with Sharp Serrated Teeth OB/Gynecology Surgical Instrument

Novak Biopsy Curette is primarily used for scraping and removing a portion of the lining of the uterus. Novak Endocervical Suction Biopsy Curette features an angled, serrated tip to aid in the removal and collection of biopsy tissue from the endocervix. The tip is at the end of a long, slender handle. The Novak curette is used with a suction device that contains the tissue in a syringe after it has been scraped and removed. It comes in multiple variations such as a curette with Leur Hub in 2mm, 3mm and 4mm diameter

NOVAK BIOPSY CURETTE FEATURES:


- Long and thin metallic tube narrow built for deep access.
- Broad suction hole for easy suctioning.
- Blunt tip to prevent piercing injuries.
- Ergonomic handle for better control
- Lightweight for comfortable handling

PACKAGE INCLUDE:

- Novak Endometrial Suction Biopsy Curette Size 9 3/4" (24.8cm) with Luer Hub 2mm Diameter.

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENTS TYPE:** Uterine Curettes Gynecology & Obstetrics Instruments**APPLICATION:** Endometrial Biopsy Curette

**NOVAK ENDOMETRIAL
SUCTION BIOPSY CURETTE****NOVAK ENDOMETRIAL
SUCTION BIOPSY CURETTE****WI-70002397****Novak Endometrial Suction Biopsy Curette
3mm 9 3/4" with Luer Hub****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Endometrial Biopsy Curette NOVAK 9 3/4" (24.8cm) Length Single Ended, 3mm Tip Curved Round Tip with Sharp Serrated Teeth OB/Gynecology Surgical Instrument


Novak Biopsy Curette is primarily used for scraping and removing a portion of the lining of the uterus. Novak Endocervical Suction Biopsy Curette features an angled, serrated tip to aid in the removal and collection of biopsy tissue from the endocervix. The tip is at the end of a long, slender handle. The Novak curette is used with a suction device that contains the tissue in a syringe after it has been scraped and removed. It comes in multiple variations such as a curette with Leur Hub in 2mm, 3mm and 4mm diameter

**NOVAK BIOPSY CURETTE
FEATURES:**

- Long and thin metallic tube narrow built for deep access.
- Broad suction hole for easy suctioning.
- Blunt tip to prevent piercing injuries.
- Ergonomic handle for better control
- Lightweight for comfortable handling

PACKAGE INCLUDE:

- Novak Endometrial Suction Biopsy Curette Size 9 3/4" (24.8cm) with Luer Hub 3mm Diameter.

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Uterine Curettes
Gynecology & Obstetrics Instruments**APPLICATION:** Endometrial Biopsy Curette**NOVAK ENDOMETRIAL
SUCTION BIOPSY CURETTE****NOVAK ENDOMETRIAL
SUCTION BIOPSY CURETTE****WI-70002396****Novak Endometrial Suction Biopsy Curette
4mm 9 3/4" with Luer Hub****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Endometrial Biopsy Curette NOVAK 9 3/4" (24.8cm) Length Single Ended, 4mm Tip Curved Round Tip with Sharp Serrated Teeth OB/Gynecology Surgical Instrument

Novak Biopsy Curette is primarily used for scraping and removing a portion of the lining of the uterus. Novak Endocervical Suction Biopsy Curette features an angled, serrated tip to aid in the removal and collection of biopsy tissue from the endocervix. The tip is at the end of a long, slender handle. The Novak curette is used with a suction device that contains the tissue in a syringe after it has been scraped and removed. It comes in multiple variations such as a curette with Leur Hub in 2mm, 3mm and 4mm diameter

**NOVAK BIOPSY CURETTE
FEATURES:**

- Long and thin metallic tube narrow built for deep access.
- Broad suction hole for easy suctioning.
- Blunt tip to prevent piercing injuries.
- Ergonomic handle for better control
- Lightweight for comfortable handling

PACKAGE INCLUDE:

- Novak Endometrial Suction Biopsy Curette Size 9 3/4" (24.8cm) with Luer Hub 4mm Diameter.

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Uterine Curettes
Gynecology & Obstetrics Instruments**APPLICATION:** Endometrial Biopsy Curette



WI-32394847

Rectus Curette 11-Inch Stainless Steel Surgical

\$800.00 per piece

Trade price from **\$300.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Rectus Curette 11-Inch Stainless Steel Surgical is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The 11" size shown is our standard; other sizes in the same pattern are made to order.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$800.00, falling to \$300.00 at 1 piece — about 62% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-35858438

Rectus Flat Back Spinal Fusion Curette 11" Straight Size

\$5.00 per piece

Trade price from **\$3.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Rectus Flat Back Spinal Fusion Curette 11" Straight Size is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a straight working end for direct, in-line access.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$5.00, falling to \$3.00 at 1 piece — about 40% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25002864

Root Elevator German Stainless Steel Ergonomic Handle

\$8.00 per piece

Trade price from \$4.00
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Root Elevator German Stainless Steel Ergonomic Handle, made for general and specialty operating theatre use. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: an ergonomic handle for long lists without fatigue.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$8.00, falling to \$4.00 at 10 sets — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003891

Sayre Periosteal Elevator 6.75" Double Ended

\$15.99 per piece

Trade price from \$14.71
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Sayre Periosteal Elevator 6.75" (17.1cm) Double-ended (Curved Sharp and Straigh Blunt Blades) with Flat Handle – Premium Quality Stainless Steel Surgical Orthopedic Veterinary Instrument.

Sayre Periosteal Elevator is a double-ended surgical instrument designed to elevate, separate, and reflect the periosteum and soft tissues from bones. Its primary applications include orthopedic surgery for fracture repairs and osteotomies, facial and Neuro/ENT surgeries for tissue isolation, and veterinary procedures for soft-tissue separation. Sayre elevator features one straight blunt end and a sharp curved end, allowing versatility in tasks. Made from surgical-grade stainless steel, it is


sterilizable and suitable for various surgical environments.

SAYRE PERIOSTEAL ELEVATOR FEATURES

- Crafted with 6.75" (17.14cm) Overall Length.
- Double Ended Profile
- One Sharp Curved End & One Straight Blunt End.
- Ergonomic Solid Handle for Maximum Control.
- Used to Elevate, Separate and Reflect the Periosteum and Soft Tissue from Bones.
- Sterilizable and Reusable.
- Made With Premium Quality Stainless Steel

PACKAGE INCLUDE:



**WI-70000667****Spratt Brun Bone Curette 6.5" Oval Cup
2.2×3.4mm Size #5/0****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking




Spratt Brun Bone Curette 6.5" (16.5cm) Oval Cup Size 5/0 (00000) 2.2mm × 3.4mm Single Ended Hollow Handle Premium Quality Stainless Steel Neurosurgical, Orthopedic, ENT Surgical Instrument.

These curettes with a small scope shape on the tip is used for scraping hard tissues from bone surfaces to expose narrow surgical sites. Spratt Bone Curette is mostly used in biopsy procedures and generally for cleaning debris or clearing excess tissue from bone or other structures to increase visibility and access to a particular area. This Bone Curette is available in a wide range of sizes from 5/0 to 6. The straight, tapered shaft to easily traverse tight passages. Moreover, the hollow handle with a thumb depression to

minimize finger strain. Spratt Curettes are crafted from surgical grade stainless steel, ideal for hospitals and specialty clinics. These are sterilizable and you can reuse it.

**SPRATT BRUN BONE CURETTE
FEATURES**

- Crafted With 6.5" (16.5cm) Overall Length.
- 2.2mm × 3.4mm Size #5/0.
- Oval Shaped Cup Blade Design.
- Sharp Blads.
- Straight Tapered Shaft.
- Very Small Tip for Scooping in Tight Spaces.
- Single Working End.
- Hollow Handle With Grooves.
- Reusable.

**WI-70000666****Spratt Brun Bone Curette 6.5" Oval Cup
2.5×4.0mm Size #4/0****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Spratt Brun Bone Curette 6.5" (16.5cm) Oval Cup Size 4/0 (0000) 2.5mm × 4.0mm Single Ended Hollow Handle Premium Quality Stainless Steel Neurosurgical, Orthopedic, ENT Surgical Instrument.

These curettes with a small scope shape on the tip is used for scraping hard tissues from bone surfaces to expose narrow surgical sites. Spratt Bone Curette is mostly used in biopsy procedures and generally for cleaning debris or clearing excess tissue from bone or other structures to increase visibility and access to a particular area. This Bone Curette is available in a wide range of sizes from 5/0 to 6. The straight, tapered shaft to easily traverse tight passages. Moreover, the hollow handle with a thumb depression to

minimize finger strain. Spratt Curettes are crafted from surgical grade stainless steel, ideal for hospitals and specialty clinics. These are sterilizable and you can reuse it.

**SPRATT BRUN BONE CURETTE
FEATURES**

- Crafted With 6.5" (16.5cm) Overall Length.
- 2.5mm × 4.0mm Size #4/0.
- Oval Shaped Cup Blade Design.
- Sharp Blads.
- Straight Tapered Shaft.
- Very Small Tip for Scooping in Tight Spaces.
- Single Working End.
- Hollow Handle With Grooves.
- Reusable.



**WI-70000665****Spratt Brun Bone Curette 6.5" Oval Cup
2.8×4.5mm Size #3/0****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Spratt Brun Bone Curette 6.5" (16.5cm) Oval Cup Size 3/0 (000) 2.8mm × 4.5mm Single Ended Hollow Handle Premium Quality Stainless Steel Neurosurgical, Orthopedic, ENT Surgical Instrument.

These curettes with a small scope shape on the tip is used for scraping hard tissues from bone surfaces to expose narrow surgical sites. Spratt Bone Curette is mostly used in biopsy procedures and generally for cleaning debris or clearing excess tissue from bone or other structures to increase visibility and access to a particular area. This Bone Curette is available in a wide range of sizes from 5/0 to 6. The straight, tapered shaft to easily traverse tight passages. Moreover, the hollow handle with a thumb depression to

minimize finger strain. Spratt Curettes are crafted from surgical grade stainless steel, ideal for hospitals and specialty clinics. These are sterilizable and you can reuse it.

**SPRATT BRUN BONE CURETTE
FEATURES**

- Crafted With 6.5" (16.5cm) Overall Length.
- 2.8mm × 4.5mm Size #3/0.
- Oval Shaped Cup Blade Design.
- Sharp Blads.
- Straight Tapered Shaft.
- Very Small Tip for Scooping in Tight Spaces.
- Single Working End.
- Hollow Handle With Grooves.
- Reusable.

**WI-70000664****Spratt Brun Bone Curette 6.5" Oval Cup
3.3×5.1mm Size #2/0****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Spratt Brun Bone Curette 6.5" (16.5cm) Oval Cup Size 2/0 (00) 3.3mm × 5.1mm Single Ended Hollow Handle Premium Quality Stainless Steel Neurosurgical, Orthopedic, ENT Surgical Instrument.


These curettes with a small scope shape on the tip is used for scraping hard tissues from bone surfaces to expose narrow surgical sites. Spratt Bone Curette is mostly used in biopsy procedures and generally for cleaning debris or clearing excess tissue from bone or other structures to increase visibility and access to a particular area. This Bone Curette is available in a wide range of sizes from 5/0 to 6. The straight, tapered shaft to easily traverse tight passages. Moreover, the hollow handle with a thumb depression to minimize finger strain. Spratt Curettes

are crafted from surgical grade stainless steel, ideal for hospitals and specialty clinics. These are sterilizable and you can reuse it.

**SPRATT BRUN BONE CURETTE
FEATURES**

- Crafted With 6.5" (16.5cm) Overall Length.
- 3.3mm × 5.1mm Size #2/0.
- Oval Shaped Cup Blade Design.
- Sharp Blads.
- Straight Tapered Shaft.
- Very Small Tip for Scooping in Tight Spaces.
- Single Working End.
- Hollow Handle With Grooves.
- Reusable.
- Made With Premium Grade Stainless Steel.



**SIZE
0****WI-70000663****Spratt Brun Bone Curette 6.5" Oval Cup
3.7×5.8mm Size #0****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking




Spratt Brun Bone Curette 6.5" (16.5cm) Oval Cup Size #0 (3.7mm × 5.8mm) Single Ended Hollow Handle Premium Quality Stainless Steel Neurosurgical, Orthopedic, ENT Surgical Instrument.

These curettes with a small scope shape on the tip is used for scraping hard tissues from bone surfaces to expose narrow surgical sites. Spratt Bone Curette is mostly used in biopsy procedures and generally for cleaning debris or clearing excess tissue from bone or other structures to increase visibility and access to a particular area. This Bone Curette is available in a wide range of sizes from 5/0 to 6. The straight, tapered shaft to easily traverse tight passages. Moreover, the hollow handle with a thumb depression to minimize finger strain. Spratt Curettes

are crafted from surgical grade stainless steel, ideal for hospitals and specialty clinics. These are sterilizable and you can reuse it.

**SPRATT BRUN BONE CURETTE
FEATURES**

- Crafted With 6.5" (16.5cm) Overall Length.
- 3.7mm × 5.8mm Size #0.
- Oval Shaped Cup Blade Design.
- Sharp Blads.
- Straight Tapered Shaft.
- Very Small Tip for Scooping in Tight Spaces.
- Single Working End.
- Hollow Handle With Grooves.
- Reusable.
- Made With Premium Grade Stainless Steel.

**SIZE
1****WI-70000662****Spratt Brun Bone Curette 6.5" Oval Cup
4.3×6.6mm Size #1****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Spratt Brun Bone Curette 6.5" (16.5cm) Oval Cup Size #1 (4.3mm × 6.6mm) Single Ended Hollow Handle Premium Quality Stainless Steel Neurosurgical, Orthopedic, ENT Surgical Instrument.


These curettes with a small scope shape on the tip is used for scraping hard tissues from bone surfaces to expose narrow surgical sites. Spratt Bone Curette is mostly used in biopsy procedures and generally for cleaning debris or clearing excess tissue from bone or other structures to increase visibility and access to a particular area. This Bone Curette is available in a wide range of sizes from 5/0 to 6. The straight, tapered shaft to easily traverse tight passages. Moreover, the hollow handle with a thumb depression to minimize finger strain. Spratt Curettes

are crafted from surgical grade stainless steel, ideal for hospitals and specialty clinics. These are sterilizable and you can reuse it.

**SPRATT BRUN BONE CURETTE
FEATURES**

- Crafted With 6.5" (16.5cm) Overall Length.
- 4.3mm × 6.6mm Size #1.
- Oval Shaped Cup Blade Design.
- Sharp Blads.
- Straight Tapered Shaft.
- Very Small Tip for Scooping in Tight Spaces.
- Single Working End.
- Hollow Handle With Grooves.
- Reusable.
- Made With Premium Grade Stainless Steel.



**SIZE
2****WI-70000661****Spratt Brun Bone Curette 6.5" Oval Cup
4.8x7.5mm Size #2****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking




Spratt Brun Bone Curette 6.5" (16.5cm) Oval Cup Size #2 (4.8mm x 7.5mm) Single Ended Hollow Handle Premium Quality Stainless Steel Neurosurgical, Orthopedic, ENT Surgical Instrument.

These curettes with a small scope shape on the tip is used for scraping hard tissues from bone surfaces to expose narrow surgical sites. Spratt Bone Curette is mostly used in biopsy procedures and generally for cleaning debris or clearing excess tissue from bone or other structures to increase visibility and access to a particular area. This Bone Curette is available in a wide range of sizes from 5/0 to 6. The straight, tapered shaft to easily traverse tight passages. Moreover, the hollow handle with a thumb depression to minimize finger strain. Spratt Curettes

are crafted from surgical grade stainless steel, ideal for hospitals and specialty clinics. These are sterilizable and you can reuse it.

**SPRATT BRUN BONE CURETTE
FEATURES**

- Crafted With 6.5" (16.5cm) Overall Length.
- 4.8mm x 7.5mm Size #2.
- Oval Shaped Cup Blade Design.
- Sharp Blads.
- Straight Tapered Shaft.
- Very Small Tip for Scooping in Tight Spaces.
- Single Working End.
- Hollow Handle With Grooves.
- Reusable.
- Made With Premium Grade Stainless Steel.

**SIZE
3****WI-70000660****Spratt Brun Bone Curette 6.5" Oval Cup
5.6x9.0mm Size #3****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Spratt Brun Bone Curette 6.5" (16.5cm) Oval Cup Size #3 (5.6mm x 9.0mm) Single Ended Hollow Handle Premium Quality Stainless Steel Neurosurgical, Orthopedic, ENT Surgical Instrument.


These curettes with a small scope shape on the tip is used for scraping hard tissues from bone surfaces to expose narrow surgical sites. Spratt Bone Curette is mostly used in biopsy procedures and generally for cleaning debris or clearing excess tissue from bone or other structures to increase visibility and access to a particular area. This Bone Curette is available in a wide range of sizes from 5/0 to 6. The straight, tapered shaft to easily traverse tight passages. Moreover, the hollow handle with a thumb depression to minimize finger strain. Spratt Curettes

are crafted from surgical grade stainless steel, ideal for hospitals and specialty clinics. These are sterilizable and you can reuse it.

**SPRATT BRUN BONE CURETTE
FEATURES**

- Crafted With 6.5" (16.5cm) Overall Length.
- 5.6mm x 9.0mm Size #3.
- Oval Shaped Cup Blade Design.
- Sharp Blads.
- Straight Tapered Shaft.
- Very Small Tip for Scooping in Tight Spaces.
- Single Working End.
- Hollow Handle With Grooves.
- Reusable.



**SIZE
4****WI-70000659****Spratt Brun Bone Curette 6.5" Oval Cup
6.1×9.5mm Size #4****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking




Spratt Brun Bone Curette 6.5" (16.5cm) Oval Cup Size #4 (6.1mm × 9.5mm) Single Ended Hollow Handle Premium Quality Stainless Steel Neurosurgical, Orthopedic, ENT Surgical Instrument.

These curettes with a small scope shape on the tip is used for scraping hard tissues from bone surfaces to expose narrow surgical sites. Spratt Bone Curette is mostly used in biopsy procedures and generally for cleaning debris or clearing excess tissue from bone or other structures to increase visibility and access to a particular area. This Bone Curette is available in a wide range of sizes from 5/0 to 6. The straight, tapered shaft to easily traverse tight passages. Moreover, the hollow handle with a thumb depression to minimize finger strain. Spratt Curettes

are crafted from surgical grade stainless steel, ideal for hospitals and specialty clinics. These are sterilizable and you can reuse it.

**SPRATT BRUN BONE CURETTE
FEATURES**

- Crafted With 6.5" (16.5cm) Overall Length.
- 6.1mm × 9.5mm Size #4.
- Oval Shaped Cup Blade Design.
- Sharp Blads.
- Straight Tapered Shaft.
- Very Small Tip for Scooping in Tight Spaces.
- Single Working End.
- Hollow Handle With Grooves.
- Reusable.
- Made With Premium Grade Stainless Steel.

**SIZE
5****WI-70000658****Spratt Brun Bone Curette 6.5" Oval Cup
6.7×10.4mm Size #5****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Spratt Brun Bone Curette 6.5" (16.5cm) Oval Cup Size #5 (6.7mm × 10.4mm) Single Ended Hollow Handle Premium Quality Stainless Steel Neurosurgical, Orthopedic, ENT Surgical Instrument.


These curettes with a small scope shape on the tip is used for scraping hard tissues from bone surfaces to expose narrow surgical sites. Spratt Bone Curette is mostly used in biopsy procedures and generally for cleaning debris or clearing excess tissue from bone or other structures to increase visibility and access to a particular area. This Bone Curette is available in a wide range of sizes from 5/0 to 6. The straight, tapered shaft to easily traverse tight passages. Moreover, the hollow handle with a thumb depression to minimize finger strain. Spratt Curettes

are crafted from surgical grade stainless steel, ideal for hospitals and specialty clinics. These are sterilizable and you can reuse it.

**SPRATT BRUN BONE CURETTE
FEATURES**

- Crafted With 6.5" (16.5cm) Overall Length.
- 6.7mm × 10.4mm Size #5.
- Oval Shaped Cup Blade Design.
- Sharp Blads.
- Straight Tapered Shaft.
- Very Small Tip for Scooping in Tight Spaces.
- Single Working End.
- Hollow Handle With Grooves.
- Reusable.
- Made With Premium Grade Stainless Steel.





WI-70000657

Spratt Brun Bone Curette 6.5" Oval Cup 8.7×12.6mm Size #6

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Spratt Brun Bone Curette 6.5" (16.5cm) Oval Cup Size #6 (8.7mm×12.6mm) Single Ended Hollow Handle Premium Quality Stainless Steel Neurosurgical, Orthopedic, ENT Surgical Instrument.

These curettes with a small scope shape on the tip is used for scraping hard tissues from bone surfaces to expose narrow surgical sites. Spratt Bone Curette is mostly used in biopsy procedures and generally for cleaning debris or clearing excess tissue from bone or other structures to increase visibility and access to a particular area. This Bone Curette is available in a wide range of sizes from 5/0 to 6. The straight, tapered shaft to easily traverse tight passages. Moreover, the hollow handle with a thumb depression to minimize finger strain. Spratt Curettes

are crafted from surgical grade stainless steel, ideal for hospitals and specialty clinics. These are sterilizable and you can reuse it.

SPRATT BRUN BONE CURETTE FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- 8.7mm × 12.6mm Size #6.
- Oval Shaped Cup Blade Design.
- Sharp Blads.
- Straight Tapered Shaft.
- Very Small Tip for Scooping in Tight Spaces.
- Single Working End.
- Hollow Handle With Grooves.
- Reusable.
- Made With Premium Grade Stainless Steel.



WI-35255073

Stainless Steel Cone Ring Curette Hollow Handle

\$8.00 per piece

Trade price from **\$6.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Cone Ring Curette Hollow Handle is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$6.00 at 10 pieces — about 25% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42909585

Stainless Steel Heidbrink H1 Straight Elevator Hollow Handle

\$7.00 per piece

Trade price from \$3.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Heidbrink H1 Straight Elevator Hollow Handle is an elevators & curette we forge, grind and finish in our own Sialkot workshop. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets. The name tells you what you are getting: a straight working end for direct, in-line access.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$7.00, falling to \$3.00 at 10 pieces — about 57% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003271

Stille Type Chisel

\$76.48 per piece

Trade price from \$70.36
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stille Type Chisel is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$76.48, falling to \$70.36 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003269

Stille Type Osteotome

\$37.68 per piece

Trade price from \$34.67
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stille Type Osteotome is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$37.68, falling to \$34.67 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-01193293

Trending Elegant Looking Rectus Curette Medical

\$5.00 per piece

Trade price from \$3.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Trending Elegant Looking Rectus Curette Medical, made for general and specialty operating theatre use. Instruments of this type lift periosteum or soft tissue, and scrape out cavities and sockets.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$5.00, falling to \$3.00 at 1 piece — about 40% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000364****U.S Army Pattern Osteotome 10mm Straight****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality U.S. Army Pattern Osteotome 6.25" (16cm) Straight 10mm Wide Blade Stainless Steel Orthopedic Surgery Instrument.

U.S. Army Pattern Osteotome is a surgical instrument used in orthopedic surgery to cut, scrape, and chisel bones, often in difficult-to-reach areas. This Osteotome is a straight, beveled instrument designed to remove small bone portions-to-reach difficult areas, using a mallet. The handle is typically flat and square-shaped, with a cap-like base for easier force application with a mallet. U.S. Army Pattern Osteotome, with a chisel blade, is versatile for surgical procedures and comes in different sizes to meet various patient needs.

**U.S. ARMY PATTERN
OSTEOTOMES FEATURES**

- Length: 6.25" (16cm) Overall Length.
- Blade: Chisel-type 10mm Wide Straight Blade.
- Functionality: Used for Cutting and Sculpting Bone, Allowing Surgeons to Reshape and Create Cavities.
- Beveling: Beveled on Both Sides for Cutting and Scraping.
- Handle: Flat, Square-shaped With a Cap-like Base, Designed for Increased Force.
- Durability: Made of High-tensile Materials Like German Stainless Steel.
- Slip-Resistant Design: Provides Secure Grip During Surgery.
- Variations: Multiple Tip Sizes.

**WI-70003268****U.S Army Pattern Osteotome 12mm Straight****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Premium Quality U.S. Army Pattern Osteotome 6.25" (16cm) Straight 12mm Wide Blade Stainless Steel Orthopedic Surgery Instrument

U.S. Army Pattern Osteotome is a surgical instrument used in orthopedic surgery to cut, scrape, and chisel bones, often in difficult-to-reach areas. This Osteotome is a straight, beveled instrument designed to remove small bone portions-to-reach difficult areas, using a mallet. The handle is typically flat and square-shaped, with a cap-like base for easier force application with a mallet. U.S. Army Pattern Osteotome, with a chisel blade, is versatile for surgical procedures and comes in different sizes to meet various patient needs.

**U.S. ARMY PATTERN
OSTEOTOMES FEATURES**

- Length: 6.25" (16cm) Overall Length.
- Blade: Chisel-type 12mm Wide Straight Blade.
- Functionality: Used for Cutting and Sculpting Bone, Allowing Surgeons to Reshape and Create Cavities.
- Beveling: Beveled on Both Sides for Cutting and Scraping.
- Handle: Flat, Square-shaped With a Cap-like Base, Designed for Increased Force.
- Durability: Made of High-tensile Materials Like German Stainless Steel.
- Slip-Resistant Design: Provides Secure Grip During Surgery.
- Variations: Multiple Tip Sizes.



**WI-70000363****U.S Army Pattern Osteotome 14mm Straight****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces



- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality U.S. Army Pattern Osteotome 6.25" (16cm) Straight 14mm Wide Blade Stainless Steel Orthopedic Surgery Instrument.

U.S. Army Pattern Osteotome is a surgical instrument used in orthopedic surgery to cut, scrape, and chisel bones, often in difficult-to-reach areas. This Osteotome is a straight, beveled instrument designed to remove small bone portions-to-reach difficult areas, using a mallet. The handle is typically flat and square-shaped, with a cap-like base for easier force application with a mallet. U.S. Army Pattern Osteotome, with a chisel blade, is versatile for surgical procedures and comes in different sizes to meet various patient needs.

**U.S. ARMY PATTERN
OSTEOTOMES FEATURES**

- Length: 6.25" (16cm) Overall Length.
- Blade: Chisel-type 14mm Wide Straight Blade.
- Functionality: Used for Cutting and Sculpting Bone, Allowing Surgeons to Reshape and Create Cavities.
- Beveling: Beveled on Both Sides for Cutting and Scraping.
- Handle: Flat, Square-shaped With a Cap-like Base, Designed for Increased Force.
- Durability: Made of High-tensile Materials Like German Stainless Steel.
- Slip-Resistant Design: Provides Secure Grip During Surgery.
- Variations: Multiple Tip Sizes.

**WI-70000362****U.S Army Pattern Osteotome 16mm Straight****\$13.99** per pieceTrade price from **\$12.87**
Minimum order 21 pieces



- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking




Premium Quality U.S. Army Pattern Osteotome 6.25" (16cm) Straight 16mm Wide Blade Stainless Steel Orthopedic Surgery Instrument.

U.S. Army Pattern Osteotome is a surgical instrument used in orthopedic surgery to cut, scrape, and chisel bones, often in difficult-to-reach areas. This Osteotome is a straight, beveled instrument designed to remove small bone portions-to-reach difficult areas, using a mallet. The handle is typically flat and square-shaped, with a cap-like base for easier force application with a mallet. U.S. Army Pattern Osteotome, with a chisel blade, is versatile for surgical procedures and comes in different sizes to meet various patient needs.

**U.S. ARMY PATTERN
OSTEOTOMES FEATURES**

- Length: 6.25" (16cm) Overall Length.
- Blade: Chisel-type 16mm Wide Straight Blade.
- Functionality: Used for Cutting and Sculpting Bone, Allowing Surgeons to Reshape and Create Cavities.
- Beveling: Beveled on Both Sides for Cutting and Scraping.
- Handle: Flat, Square-shaped With a Cap-like Base, Designed for Increased Force.
- Durability: Made of High-tensile Materials Like German Stainless Steel.
- Slip-Resistant Design: Provides Secure Grip During Surgery.
- Variations: Multiple Tip Sizes.



**WI-70000361****U.S Army Pattern Osteotome 18mm Straight****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality U.S. Army Pattern Osteotome 6.25" (16cm) Straight 18mm Wide Blade Stainless Steel Orthopedic Surgery Instrument.

U.S. Army Pattern Osteotome is a surgical instrument used in orthopedic surgery to cut, scrape, and chisel bones, often in difficult-to-reach areas. This Osteotome is a straight, beveled instrument designed to remove small bone portions-to-reach difficult areas, using a mallet. The handle is typically flat and square-shaped, with a cap-like base for easier force application with a mallet. U.S. Army Pattern Osteotome, with a chisel blade, is versatile for surgical procedures and comes in different sizes to meet various patient needs.

U.S. ARMY PATTERN OSTEOTOMES FEATURES

- Length: 6.25" (16cm) Overall Length.
- Blade: Chisel-type 18mm Wide Straight Blade.
- Functionality: Used for Cutting and Sculpting Bone, Allowing Surgeons to Reshape and Create Cavities.
- Beveling: Beveled on Both Sides for Cutting and Scraping.
- Handle: Flat, Square-shaped With a Cap-like Base, Designed for Increased Force.
- Durability: Made of High-tensile Materials Like German Stainless Steel.
- Slip-Resistant Design: Provides Secure Grip During Surgery.
- Variations: Multiple Tip Sizes.

**WI-70000366****U.S Army Pattern Osteotome 6mm Straight****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality U.S. Army Pattern Osteotome 6" (15.2cm) Straight 6mm Wide Blade Stainless Steel Orthopedic Surgery Instrument

U.S. Army Pattern Osteotome is a surgical instrument used in orthopedic surgery to cut, scrape, and chisel bones, often in difficult-to-reach areas. This Osteotome is a straight, beveled instrument designed to remove small bone portions-to-reach difficult areas, using a mallet. The handle is typically flat and square-shaped, with a cap-like base for easier force application with a mallet. U.S. Army Pattern Osteotome, with a chisel blade, is versatile for surgical procedures and comes in different sizes to meet various patient needs.

U.S. ARMY PATTERN OSTEOTOMES FEATURES

- Length: 6" (15.2cm) Overall Length.
- Blade: Chisel-type 6mm Wide Straight Blade.
- Functionality: Used for Cutting and Sculpting Bone, Allowing Surgeons to Reshape and Create Cavities.
- Beveling: Beveled on Both Sides for Cutting and Scraping.
- Handle: Flat, Square-shaped With a Cap-like Base, Designed for Increased Force.
- Durability: Made of High-tensile Materials Like German Stainless Steel.
- Slip-Resistant Design: Provides Secure Grip During Surgery.
- Variations: Multiple Tip Sizes.

**WI-70000365****U.S Army Pattern Osteotome 8mm Straight****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality U.S. Army Pattern Osteotome 6.25" (16cm) Straight 8mm Wide Blade Stainless Steel Orthopedic Surgery Instrument.

U.S. Army Pattern Osteotome is a surgical instrument used in orthopedic surgery to cut, scrape, and chisel bones, often in difficult-to-reach areas. This Osteotome is a straight, beveled instrument designed to remove small bone portions-to-reach difficult areas, using a mallet. The handle is typically flat and square-shaped, with a cap-like base for easier force application with a mallet. U.S. Army Pattern Osteotome, with a chisel blade, is versatile for surgical procedures and comes in different sizes to meet various patient needs.

**U.S. ARMY PATTERN
OSTEOTOMES FEATURES**

- Length: 6.25" (16cm) Overall Length.
- Blade: Chisel-type 8mm Wide Straight Blade.
- Functionality: Used for Cutting and Sculpting Bone, Allowing Surgeons to Reshape and Create Cavities.
- Beveling: Beveled on Both Sides for Cutting and Scraping.
- Handle: Flat, Square-shaped With a Cap-like Base, Designed for Increased Force.
- Durability: Made of High-tensile Materials Like German Stainless Steel.
- Slip-Resistant Design: Provides Secure Grip During Surgery.
- Variations: Multiple Tip Sizes.





WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

Request a quotation

Send us product codes from this catalogue, quantities and your country.
We reply the same working day on WhatsApp.

Online shop	shop.waqarintl.com
Email	sales@waqarintl.com
WhatsApp /	
Phone	+92 307 6173627
Works	Roras Road, Opposite Anwar Khawaja Industries, Sialkot 51310, Pakistan

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.