



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026 · PART 2 OF 2

Forceps Surgical Instruments

Reusable stainless-steel forceps forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



331

PRODUCTS

\$1.45

PRICES FROM

OEM

PRIVATE LABEL

Global

WORLDWIDE
EXPORT

Contents

331 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

Product	Code	Price	Page
Hartmann Mosquito Forceps 3.5" (8.9cm) Straight Gold Plated	WI-70000386	\$8.99	10
Hartmann Mosquito Forceps 3.5" Curved	WI-70002625	\$8.99	10
Hartmann Mosquito Forceps 3.5" Straight	WI-70003528	\$8.99	11
Hartmann-Wullstein Ear Forceps 8.5cm	WI-70003260	\$46.86	11
HEANEY FORCEPS 20.0cm Stainless Steel Surgical Instruments	WI-26032810	\$8.00	12
Heaney Hysterectomy Forceps 203mm Stainless Steel Surgical Instruments	WI-26019784	\$8.00	12
Heaney Hysterectomy Forceps 8.25" Double Tooth Curved Jaws	WI-70001292	\$14.99	13
Heaney Hysterectomy Forceps 8.25" Double Tooth Straight Jaws	WI-70000068	\$14.99	13
Heaney Hysterectomy Forceps 8.25" Single Tooth Curved Jaws	WI-70001293	\$12.99	14
Heaney Hysterectomy Forceps 8.25" Single Tooth Straight Jaws	WI-70000069	\$12.99	14
Heaney-Ballentine Hysterectomy Forcep 8.5" Curved Longitudinal Serrated Double Tooth	WI-70000066	\$14.99	15
Heaney-Ballentine Hysterectomy Forcep 8.5" Curved Longitudinal Serrated Single Tooth	WI-70000067	\$12.99	15
Heaney-Ballentine Hysterectomy Forcep 8.5" Straight Longitudinal Serrated Double Tooth	WI-70000064	\$14.99	16
Heaney-Ballentine Hysterectomy Forcep 8.5" Straight Longitudinal Serrated Single Tooth	WI-70000065	\$12.99	16
Hemostat Forceps Straight & Curved with Grip Handle	WI-70002510	\$8.99	17
High Quality Anti Static Stain Less Steel Satin Eyelash Tweezer Eyebrow Electric Tweezers	WI-70003599	\$14.75	17
Hoff Towel Clamp Forceps 5.25"	WI-70000382	\$9.99	18
Hoff Towel Forceps 3.5"	WI-70002922	\$8.99	18
Hot Sales Debakey Forceps Straight Debakey Vascular Tissue Forceps	WI-34850411	\$17.00	19
Hot Sales Stainless Steel Debakey Forceps Straight Vascular Tissue	WI-34796578	\$17.00	19
Hulka Controlling Tenaculum Uterine Forceps 11"	WI-70002014	\$11.99	20
Hysterectomy Forceps Surgical Stainless Steel Medical	WI-42586281	\$6.50	20
Iris Tissue Forceps 4" 1x2 Teeth Fully Curved	WI-70000913	\$9.99	21
Iris Tissue Forceps 4" 1x2 Teeth Half Curved	WI-70000914	\$9.99	21
Iris Tissue Forceps 4" 1x2 Teeth Straight	WI-70000915	\$9.99	22
Ivory Light Rubber Dam Clamp Forceps 6.7"	WI-70001700	\$8.99	22
Jackson Endoscopic Forceps 24" Angled Cup Jaw 4mm Bite	WI-70000199	\$109.99	23
Jackson Endoscopic Forceps 24" Straight Cup Jaw 4mm Bite	WI-70000200	\$109.99	23
Jansen Ear Dressing Forceps 6.25"	WI-70003636	\$8.99	24
Jar Forceps Stainless Steel Medical Instrument for Precise Handling	WI-47140991	\$7.90	24
Jar Forceps Surgical Instrument for Secure Grasping Handling	WI-47170605	\$8.00	25
Jar Surgical Forceps Reusable Medical Instrument for Controlled Grasping	WI-47135996	\$7.50	25
Jewelers Forceps Fine Tip Stainless Steel Precision Handling	WI-43520264	\$6.00	26
Jewelers Forceps Reusable Stainless Steel Precision Instrument	WI-43522144	\$6.00	26
Kelly Forceps 5.5" Curved Full Black Coated	WI-70000270	\$9.99	27
Kelly Forceps 5.5" Straight Full Black Coated	WI-70003520	\$9.99	27
Kelly Forceps Stainless Steel Curved Straight Full-Body	WI-27077854	\$3.00	28
Kelly Hemostat Forceps 5.5" Curved Fishing Pliers	WI-70002600	\$8.99	28

Contents (continued)

Product	Code	Price	Page
Kelly Hemostat Forceps 6.25" Straight Fishing Plier	WI-70002605	\$9.99	31
Kelly Hemostatic Forceps Curved 5.5" Blue PVC Handle	WI-70000266	\$9.99	31
Kelly Hemostatic Forceps Curved 5.5" Red PVC Handle	WI-70000259	\$9.99	32
Kelly Hemostatic Forceps Curved 6.25" Blue PVC Handle	WI-70000265	\$10.99	32
Kelly Hemostatic Forceps Curved 6.25" Red PVC Handle	WI-70000258	\$10.99	33
Kelly Hemostatic Forceps Straight 5.5" Blue PVC Handle	WI-70000267	\$9.99	33
Kelly Hemostatic Forceps Straight 5.5" Red PVC Handle	WI-70000261	\$9.99	34
Kelly Hemostatic Forceps Straight 6.25" Blue PVC Handle	WI-70000264	\$10.99	34
Kelly Hemostatic Forceps Straight 6.25" Red PVC Handle	WI-70000260	\$10.99	35
Kelly Placenta Forceps 12.5" Curved Serrated Oval Tip Jaws	WI-70000398	\$13.99	35
Kelly Placenta Forceps 12.5" Straight Serrated Oval Tip Jaws	WI-70000397	\$13.99	36
Kelly Placenta Forceps Set 4-Piece OB/GYN Tools	WI-24316734	\$20.00	36
Keramo Grip Forceps 14.5 Cm TC Root Curved	WI-45053571	\$20.00	37
Kitchen Cooking Tongs 10" Round Tip BBQ Grill Tweezers – Gold	WI-70001532	\$15.99	37
Kitchen Cooking Tongs 10" Straight Tip BBQ Grill Tweezers – Gold	WI-70001534	\$14.99	38
Kitchen Cooking Tongs 12" Round Tip BBQ Grill Tweezers – Gold	WI-70001531	\$18.99	38
Kitchen Cooking Tongs 12" Straight Tip BBQ Grill Tweezers – Gold	WI-70001533	\$17.99	39
Kitchen Tweezer 10" Round Tip BBQ Grill Plating Long Food Tongs – Silver	WI-70001538	\$14.99	39
Kitchen Tweezer 10" Straight Tip BBQ Grill Plating Long Food Tongs – Silver	WI-70001540	\$13.99	40
Kitchen Tweezer 12" Round Tip BBQ Grill Plating Long Food Tongs – Silver	WI-70001537	\$17.99	40
Kitchen Tweezer 12" Straight Tip BBQ Grill Plating Long Food Tongs – Silver	WI-70001539	\$16.99	41
Kocher Forceps 14cm	WI-70003526	\$12.99	41
Kocher Tissue Forceps Stainless Steel Surgical Grasping Instrument	WI-35038993	\$2.30	42
Kocher Tissue Forceps Stainless Steel Surgical Instrument Set	WI-44538586	\$2.30	42
Laborde Tracheal Dilator Forceps 5.5" 3 Prong with Gold Handle	WI-70001544	\$14.99	43
Laborde Tracheal Dilator Forceps 5.5" with 3 Blades	WI-70001549	\$8.99	43
Lahey Cholecystectomy Forceps Surgical Instrument	WI-47040507	\$8.00	44
Lahey Cholecystectomy Forceps Surgical Tool	WI-47060312	\$7.00	44
Lahey Gall Duct Forceps 7.5" Curved Longitudinal Serrated Jaws	WI-70000867	\$9.99	45
Lahey Gall Duct Forceps 9" Curved Longitudinal Serrated Jaws	WI-70000869	\$10.99	45
Lahey Surgical Forceps Cholecystectomy Instrument	WI-47054374	\$7.90	46
Lane Tissue Forceps Atraumatic Surgical Tissue Holding Forceps	WI-36254879	\$2.50	46
Lane Tissue Forceps Reusable Surgical Instrument	WI-35090434	\$2.30	47
Lane Tissue Forceps Stainless Steel Surgical Tissue Manipulation Instrument	WI-35076893	\$2.30	47
Light Silver Eyebrow Tweezers Eye Lash Extension Tweezer	WI-70003596	\$9.99	48
Light Silver Eyebrow Tweezers Eyelash Extension Tweezer	WI-70003595	\$17.43	48
London College Tweezer 6" (15.2cm) Dental Surgical	WI-70003756	\$8.99	49
Lower Molar Extraction Forceps Serrated Handles Ergonomic	WI-34249008	\$8.00	49
Lower Molar Extraction Forceps Serrated Handles Firm	WI-32732715	\$8.00	50
Lucae Ear Forceps 5.5" Bayonet Shape with Serrated Tips	WI-70001812	\$8.99	50
Magill Catheter Forceps 10" (25.4cm) – Adult	WI-70002132	\$11.99	51
Magill Catheter Forceps 8" (20.3cm) – Child	WI-70002133	\$10.99	51

Contents (continued)

Product	Code	Price	Page
Magill Forceps 6" (15.2cm) – Infant	WI-70004116	\$8.99	52
Magnifying Splinter Forcep Tweezer 3.5" Fine Pointed Tips	WI-70000532	\$8.99	52
Maier Uterine Polypus Dressing Forceps 10" Serrated 5X30mm Straight Jaws	WI-70000723	\$14.99	53
Mashing Tweezer 8" Glass Work Bead Masher	WI-70001738	\$9.99	53
Maxillofacial Surgical Instrument Set Stainless Steel Oral Surgery	WI-36155701	\$190.00	54
McClintock Ovum Forceps Gynecological Surgical Instrument	WI-47075025	\$10.00	54
McClintock Ovum Forceps Gynecology Surgical Instrument	WI-47060663	\$9.50	55
McClintock Steel Ovum Forceps Reusable Gynecological Instrument	WI-47072473	\$9.90	55
Mcindoe Tissue Forceps 6" Straight 1.5mm Wide 1x2 Teeth	WI-70003735	\$8.99	56
Mclean-Tucker Obstetrical Forceps	WI-70003374	\$45.99	56
McPherson Suture Tying Forceps Precision Eye Surgery Tool	WI-43629101	\$9.00	57
Medic Instrument Foerster Stainless Steel	WI-28242334	\$8.00	57
Medical Stainless Steel Sponge Forceps Curved	WI-26119318	\$3.00	58
Medical Surgical Stainless Steel Hospital Forceps	WI-27738897	\$5.50	58
Meriam Tweezer 16cm – Serrated Angled Tip	WI-70002134	\$8.99	59
Michel Clip Forceps Surgical Instrument for Precise Skin Clip	WI-46972100	\$9.90	59
Michel Clip Removal Forceps	WI-46951208	\$9.00	60
Michel Clip Surgical Forceps Precision Instrument for Skin Closure	WI-46970135	\$9.50	60
Michel Skin Clip Forceps Stainless Steel Surgical Instrument	WI-46933488	\$10.00	61
Micro Alligator Forceps 3.3" Serrated Jaws	WI-70001651	\$8.99	61
Micro Alligator Forceps 3.5" Serrated Jaws	WI-70001650	\$12.99	62
Micro Alligator Forceps 3"	WI-70002749	\$8.99	62
Micro Alligator Forceps Black Coated 3.5"	WI-70004283	\$12.99	63
Micro Extracting Forceps #133L Left Angled Tapered Tips	WI-70000278	\$14.99	63
Micro Extracting Forceps #133R Right Angled Tapered Tips	WI-70000277	\$14.99	64
Microsurgery Forceps 7.25" TC Micro Tying Forcep	WI-70003734	\$18.74	64
Miller Articulating Paper Forceps 6" Straight Surgical Dental Tweezer	WI-70003762	\$8.99	65
Millin Vulsellum Forceps Gynecological Surgical Instrument for Secure Tissue	WI-47039525	\$12.00	65
Millin Vulsellum Forceps Reusable Gynecological Surgical Tool for Controlled	WI-47056422	\$11.90	66
Millin Vulsellum Surgical Forceps with Clip for Gynecology-Precise	WI-47062362	\$12.00	66
Mixer Forceps 11" (27.9cm) Right Angled Longitudinal Serration Jaws	WI-70000160	\$15.99	67
Mixer Forceps 7.25" (18.4cm) Fully Curved Serrated Jaws	WI-70000147	\$9.99	67
Mixer Forceps 9" (22.9cm) Fully Curved Serrated Jaws	WI-70000146	\$12.99	68
Mixer Forceps 9" (22.9cm) Right Angled Longitudinal Serration Jaws	WI-70000161	\$13.99	68
Mixer Petit Surgical Forceps Curved Instrument for Controlled Vessel	WI-47249272	\$4.90	69
Multi Rainbow Color Kelly Hemostat Forceps 5.5" Curved	WI-70001725	\$10.99	69
Multi Rainbow Color Kelly Hemostat Forceps 5.5" Straight	WI-70001726	\$10.99	70
Multi Rainbow Mosquito Forceps 5" Curved	WI-70001730	\$9.99	70
Multi Rainbow Mosquito Forceps 5" Straight	WI-70001731	\$9.99	71
Museux Tenaculum Forcep 6mm Wide 2x2 Teeth 9.5" Curved	WI-70000977	\$14.99	71
Museux Tenaculum Forceps 6mm Wide 2x2 Teeth 9.5" Straight	WI-70000976	\$14.99	72
Naegele Obstetrical Forceps	WI-70003373	\$127.51	72

Contents (continued)

Product	Code	Price	Page
New Eye Lash Straight Fine Point Color Tweezers	WI-70003593	\$9.99	73
New Fine Point Eyelash Extension Xtreme Lashes Tweezers	WI-70003592	\$23.20	73
Pean Hemostat 12" Curved Fishing Forceps Locking Clamps	WI-70002500	\$15.99	74
Pean Hemostat 12" Straight Fishing Forceps Locking Clamps	WI-70002501	\$15.99	74
Pean Hemostat 14" Curved Fishing Forceps Locking Clamps	WI-70002497	\$18.99	75
Pean Hemostat 14" Straight Fishing Forceps Locking Clamps	WI-70002498	\$18.99	75
Pean Hemostat 16" Curved Fishing Forceps Locking Clamps	WI-70002495	\$21.99	76
Pean Hemostat 16" Straight Fishing Forceps Locking Clamps	WI-70002496	\$21.99	76
Pedo Extracting Forceps #150XAS Split Beak Serrated	WI-70000405	\$11.99	77
Pedo Extracting Forceps #151XAS Diamond Dusted Gold Hollow Handle	WI-70000140	\$14.99	77
Pedo Extracting Forceps #151XAS for Lower Anterior	WI-70000406	\$11.99	78
Pedo Extracting Forceps #74N Narrow Beaks	WI-70000132	\$12.99	78
Pedo Extracting Forceps #A Dental Surgical Instrument	WI-70001405	\$11.99	79
Pedo Extracting Forceps #B Lower Roots Universal Mini	WI-70000415	\$11.99	79
Pedo Extracting Forceps #C Lower Molars Universal Mini	WI-70000414	\$11.99	80
Pedo Extracting Forceps #D Lower Incisors Klein Pedodontic	WI-70000412	\$11.99	80
Pedo Extracting Forceps #E Upper Front Universal Mini	WI-70000413	\$11.99	81
Pedo Extracting Forceps #F Dental Surgical Instrument	WI-70001404	\$11.99	81
Pedo Extracting Forceps #H Upper & Lower Universal Mini	WI-70000418	\$11.99	82
Peers-Bertram Towel Forceps 5.75" Serrated Tips with Grooves	WI-70002919	\$17.99	82
Piper Obstetrical Forceps 17 1/2" (44.5cm)	WI-70003372	\$49.99	83
Plate and Bone Holding Forceps 5.5" with Footplate	WI-70004064	\$11.99	83
Plate and Bone Holding Forceps 5.5" with Footplate & Gold Handle	WI-70000650	\$39.99	84
Platypus Nail Pulling Forceps 5.5" Narrow Jaws	WI-70000432	\$9.99	84
Platypus Nail Pulling Forceps 5.5" Standard Wide Jaws	WI-70000999	\$10.99	85
Pozzi Tenaculum Forceps 25cm	WI-00849593	\$2.98	85
Pozzi Tenaculum Forceps Gynecology Surgical Instrument for Secure Cervical	WI-47256124	\$4.90	86
Puller Serrated Hemostat Forceps 2 1/2" Endodontic Instrument with Pulling Cap	WI-	\$12.99	86

Contents (continued)

Product	Code	Price	Page
Reusable Medical Equipment Surgical Stainless Steel Forceps Sponge	WI-28193341	\$8.00	94
Reusable Surgical Dressing Forceps Stainless Steel Medical Tissue	WI-45080905	\$2.40	94
ROBA English Pattern Extracting Forceps #35A Upper Premolars	WI-70001104	\$14.99	95
ROBA Extracting Forceps #34M Upper Incisors and Cuspids Serrated Jaws	WI-70001241	\$14.99	95
ROBA Extracting Forceps #34N Upper Incisors and Eye Teeth	WI-70001240	\$14.99	96
ROBA Extracting Forceps #35M Upper Premolars	WI-70001103	\$14.99	96
Robin Anchoring Forceps Reusable Stainless Steel Surgical Tool	WI-46953081	\$35.00	97
Robin Anchoring Forceps Surgical Instrument Set	WI-46958051	\$40.00	97
Robin Anchoring Surgical Forceps Precision Instrument for Tissue Fixation	WI-46952103	\$32.00	98
Rochester Carmalt Forceps 6.25" Curved	WI-70001072	\$9.99	98
Rochester Carmalt Forceps 6.25" Straight	WI-70001073	\$9.99	99
Rochester Carmalt Forceps 8" Curved	WI-70001070	\$9.99	99
Rochester Carmalt Forceps 8" Straight	WI-70001071	\$9.99	100
Rochester Pean / Spencer Well Artery Forceps Curved/Straight 6/7/8/9/10	WI-27261542	\$2.55	100
Rochester Pean Forceps 14cm Curved	WI-70003524	\$8.99	101
Rochester Pean Hemostat Forceps 14cm Straight	WI-70002515	\$8.99	101
Rochester Pean Hemostat Forceps 18" Curved Locking Clamps	WI-70002493	\$23.99	102
Rochester Pean Hemostat Forceps 18" Straight Locking Clamps	WI-70002494	\$23.99	102
Rochester Pean Hemostat Forceps 20" Curved Locking Clamps	WI-70002491	\$26.99	103
Rochester Pean Hemostat Forceps 20" Straight Locking Clamps	WI-70002492	\$26.99	103
Rochester Pean Hemostat Forceps 24" (61cm) Curved	WI-70002489	\$31.99	104
Rochester Pean Hemostat Forceps 24" Straight Locking Clamps	WI-70002490	\$31.99	104
Rochester Pean Hemostat Forceps 6.25" Curved	WI-70002513	\$8.99	105
Rochester Pean Hemostat Forceps 6.25" Straight	WI-70002514	\$8.99	105
Rochester-Ochsner Forceps 5.5" Curved 1x2 Teeth Serrated Jaws	WI-70000083	\$10.99	106
Rochester-Ochsner Forceps 5.5" Straight 1x2 Teeth Serrated Jaws	WI-70000084	\$10.99	106
Rochester-Ochsner Forceps 6.25" Curved 1x2 Teeth Serrated Jaws	WI-70000081	\$10.99	107
Rochester-Ochsner Forceps 6.25" Straight 1x2 Teeth Serrated Jaws	WI-70000082	\$11.99	107
Rochester-Ochsner Forceps 7.25" Curved 1x2 Teeth Serrated Jaws	WI-70000079	\$11.99	108
Rochester-Ochsner Forceps 7.25" Straight 1x2 Teeth Serrated Jaws	WI-70000080	\$11.99	108
Rochester-Ochsner Forceps 8" Curved 1x2 Teeth Serrated Jaws	WI-70000077	\$12.99	109
Rochester-Ochsner Forceps 8" Straight 1x2 Teeth Serrated Jaws	WI-70000078	\$12.99	109
Ruskin Kleinkert Kutz Bone Splitting Forceps 6" Straight – Double Action	WI-70001305	\$49.99	110
Russian Thumb Forceps 6" (15.2cm) Round Serrated Tips	WI-70003733	\$8.99	110
Russian Tissue Forceps 8" Round Serrated Tip	WI-70002218	\$9.99	111
Satinsky DeBakey Dissecting Forceps 6" (15cm)	WI-34791613	\$17.00	111
Schaedel Towel Forceps	WI-70002918	\$30.21	112
Schroeder Braun Uterine Tenaculum Forceps 9.5" Straight	WI-70002335	\$12.49	112
Semken Dressing Forceps Stainless Steel Surgical Instrument Controlled Tissue	WI-47228647	\$2.20	113
Semken Steel Dressing Forceps Surgical Instrument	WI-47258231	\$2.50	113
Semken Thump Forceps	WI-70003732	\$8.99	114
Simpson Obstetrical Forceps 14" (35.6cm) – Fenestrated Blades	WI-70003371	\$39.99	114

Contents (continued)

Product	Code	Price	Page
Simpson-Braun Obstetrical Forceps	WI-70003370	\$127.51	115
Simpson-Luikart Obstetrical Forceps 14" Solid Blade	WI-70003369	\$43.99	115
Singley Surgical Tissue Forceps	WI-47222924	\$3.90	116
Singley Thump Forceps	WI-70003731	\$12.28	116
Singley Tissue Forceps Stainless Steel Surgical Instrument	WI-47235367	\$4.00	117
Singley Tissue Grasping Forceps Reusable Surgical Instrument Secure	WI-47258123	\$3.80	117
Skin Staple Remover Forceps Surgical Kit for Precise Staple	WI-47016212	\$15.00	118
Skin Staple Remover Forceps-Reusable	WI-47002321	\$14.90	118
Slanted Eyebrow Tweezers 4" Angled Tips	WI-70000144	\$9.99	119
Splinter Forceps 3.5" Fine Point Serrated Tips	WI-70000818	\$8.99	119
Splinter Forceps 4.5" Fine Point Curved Serrated Tips	WI-70000814	\$8.99	120
Splinter Forceps 4.5" Fine Point Serrated Tips	WI-70000819	\$8.99	120
Splinter Forceps 5" Fine Point Curved Serrated Tips	WI-70000813	\$8.99	121
Splinter Forceps 5" Fine Point Straight Serrated Tips	WI-70000816	\$8.99	121
Sponge Forceps 12 Inch Straight	WI-28223336	\$8.00	122
Sponge Forceps 9.5" Straight Serrated Surgical	WI-28180345	\$8.00	122
Sponge Forceps Blue Coated 9.5" Straight Serrated Jaws	WI-70002228	\$14.99	123
Sponge Forceps Stainless Steel Medical Surgical Dressing Forceps	WI-44567176	\$2.90	123
Stainless Forceps	WI-25867723	\$5.00	124
Stainless Steel 6.4" Babcock Tissue Grasping Forceps Power	WI-30847995	\$3.00	124
Stainless Steel 8 Inch Babcock Forceps Atraumatic Tissue	WI-36329018	\$2.30	125
Stainless Steel Artery Forceps Straight Tip Surgical Instrument Set	WI-34445744	\$5.50	125
Stainless Steel Bayonet Dressing Forceps Surgical Tissue Forceps	WI-36115438	\$5.00	126
Stainless Steel Brown Tissue Forceps Surgical Instrument Medical Forceps	WI-36176181	\$3.00	126
Stainless Steel Cheatle Forceps (10	WI-39628463	\$3.00	127
Stainless Steel Debakey	WI-34844439	\$17.00	127
Stainless Steel Debakey Forceps Instrument	WI-34795638	\$17.00	128
Stainless Steel Debakey Forceps Reusable Surgical	WI-35031400	\$3.27	128
Stainless Steel Entomology Dissection Kit Teaser Needles Adson	WI-24100124	\$30.00	129
Stainless Steel Foerster Sponge Forceps 14cm	WI-34972290	\$3.50	129
Stainless Steel Foerster Sponge Holding Forceps 10.5	WI-34904583	\$3.50	130
Stainless Steel Foldable Portable Surgical Cheatle Forceps	WI-34963423	\$3.00	130
Stainless Steel Hair Transplant Forceps Surgical	WI-33270313	\$12.50	131
Stainless Steel Kocher Tissue Forceps Surgical Instrument Set	WI-44547556	\$2.30	131
Stainless Steel Magill Forceps for Medical Intubation	WI-34952384	\$6.00	132
Stainless Steel Pediatric Mixer Dissecting Forceps	WI-25991828	\$8.00	132
Stainless Steel Skin Staple Removal Forceps	WI-46980504	\$14.80	133
Stainless Steel Sponge Dressing Forceps Curved Slims Forceps	WI-26095341	\$3.00	133
Stainless Steel Sponge Holding Forceps	WI-34936434	\$3.50	134
Stainless Steel Surgical Dressing Forceps	WI-10247516	\$2.99	134
Stainless Steel Surgical Forceps Sponge Holding	WI-27187853	\$3.00	135
Stainless Steel Surgical Instruments Set Foerster Sponge	WI-27191808	\$3.00	135

Contents (continued)

Product	Code	Price	Page
Stainless Steel TC Babcock Forceps for Atraumatic	WI-30909034	\$3.00	136
Stainless Steel Utrata Capsulorrhesis Forceps Eye Surgery	WI-43480659	\$12.00	136
Standard Articulating Paper Forcep 6" Cross Serrated Tips	WI-70000552	\$8.99	137
Steel Babcock Tissue Intestinal Holding Forceps Surgical	WI-30931001	\$3.00	137
Steel Magill Catheter Forceps Surgical Instrument for Precise Airway	WI-47223591	\$8.00	138
Steel Surgical Sponge Holding Dressing Forceps	WI-30785999	\$3.00	138
Stille Thumb Forceps	WI-70003730	\$8.99	139
Stille Thump Forceps	WI-70003729	\$9.90	139
Straight ESD 12 Fine Pointed Eyelash Tweezer	WI-70003583	\$9.99	140
Straight Polish Silver Straight Eyelash Tweezers Custom Eyelash Extension Straight Tweezer	WI-70003587	\$9.99	140
Straight Sponge Holding Forceps Stainless Steel Surgical	WI-44565366	\$2.90	141
Straight Surgical Instrument Dressing Forceps	WI-08601496	\$3.00	141
Straight Surgical Tissue Forceps Brown Pattern	WI-36183192	\$3.00	142
Straight Tweezers for Eyelash Extension	WI-70003586	\$9.99	142
Strong Curved Volume Eye Care Curved Tweezers for Volume Eyelash Extension	WI-70003585	\$9.99	143
Surgical Babcock Tissue Forceps Straight Size	WI-30866109	\$3.00	143
Surgical Dressing Tissue Forceps	WI-27282717	\$2.39	144
Surgical Kelly Locking Forceps Stainless Steel Reusable Different	WI-33339752	\$3.00	144
Surgical Stainless Steel Gross Forceps 20cm	WI-30819959	\$3.00	145
Surgical Use Reusable Steel Sponge Holding Forceps	WI-26066344	\$3.00	145
T/C Adson Brown Forceps 4.75" 9×9 Teeth Gold Handle	WI-70001477	\$12.99	146
T/C Adson Brown Tissue Forceps 7×7 Teeth 4.75"	WI-70001480	\$12.99	146
T/C Adson Dressing Forceps 4.75" Serrated Gold Handle	WI-70001479	\$12.99	147
T/C Adson Dressing Forceps 6" Serrated Gold Handle	WI-70001478	\$13.99	147
T/C Adson Micro Tissue Forceps 1×2 Teeth 4.75"	WI-70003743	\$12.99	148
T/C Adson Tissue Forceps 1×2 Teeth 4.75" Gold Handle	WI-70001476	\$12.99	148
T/C Adson Tissue Forceps 1×2 Teeth 6" Gold Handle	WI-70001475	\$12.99	149
T/C Extracting Forceps #859 for Lower Root & Molars, Delicate	WI-70003314	\$14.99	149
T/C Micro Suture Tying Forceps 5.5" with Platform Gold Plated Handle	WI-70002212	\$12.99	150
TC Potts Smith Forceps Gold Handle Stainless Steel Surgical	WI-34992064	\$3.00	150
TC Potts Smith Forceps Straight Tungsten Carbide Tips	WI-35015094	\$3.00	151
Tenaculum Pozzi Forceps Gynecological Surgical Instrument for Precise Cervical	WI-47263071	\$5.00	151
Tenaculum Pozzi Surgical Forceps Precision Gynecological Instrument Cervical	WI-47243244	\$4.88	152
Thomas-Gaylor Uterine Biopsy Punch Forceps 8.5" Curved Jaws 5mm Bite	WI-70000694	\$24.99	152
Thomas-Gaylor Uterine Biopsy Punch Forceps 9.5" Curved Jaws 5mm Bite	WI-70000693	\$26.99	153
Three Prong Sterilizer Forceps 11"	WI-70000335	\$19.99	153
Three Prong Sterilizer Forceps 8"	WI-70000336	\$15.99	154
Thumb Tissue Forceps 1×2 Teeth 4.5"	WI-70000794	\$8.99	154
Thumb Tissue Forceps 1×2 Teeth 5.5"	WI-70000792	\$8.99	155

Contents (continued)

Product	Code	Price	Page
Thumb Tissue Forceps 2×3 Teeth 5.5"	WI-70000787	\$8.99	157
Thumb Tissue Forceps 2×3 Teeth 5"	WI-70000788	\$8.99	157
Thumb Tissue Forceps 2×3 Teeth 6"	WI-70000786	\$8.99	158
Thumb Tissue Forceps 3×4 Teeth 5.5"	WI-70000783	\$8.99	158
Thumb Tissue Forceps 3×4 Teeth 5"	WI-70000784	\$8.99	159
Thumb Tissue Forceps 3×4 Teeth 6"	WI-70000782	\$9.99	159
Thumb Tissue Forceps 4×5 Teeth 5.5"	WI-70000780	\$8.99	160
Tischler Kevorkian Biopsy Punch Forceps 3mm x 9.5mm Bite	WI-70002654	\$59.99	160
Tischler-Morgan Biopsy Forceps Oval Shape 3mmx7mm Bite	WI-70002653	\$59.99	161
Tissue Forceps Brown Style Stainless Steel Surgical	WI-36089992	\$3.00	161
Tohoku Towel Forceps	WI-70002917	\$30.21	162
Towel Forceps	WI-70002916	\$8.99	162
Troeltsch Wilde Ear Forceps 5" Angular 1×2 Teeth	WI-70001053	\$9.99	163
Trousseau Tracheal Dilator Forceps 5.5" ENT Tracheostomy Instrument	WI-70001824	\$9.99	163
Ulrich Dressing Forceps	WI-70003699	\$30.84	164
Universal Pado Extracting Forceps #150 1/2S for Upper Teeth	WI-70004154	\$11.99	164
Universal Pado Extracting Forceps #151 1/2S For Lower Teeth	WI-70004119	\$11.99	165
University of Washington Rubber Dam Clamp Forceps 7"	WI-70001703	\$10.99	165
Vasectomy Meatotomy Instrument Set Stainless Steel Forceps	WI-36123694	\$150.00	166
Verbrugge Bone Holding Forceps 10" Self Centering with Speed Lock Premium Quality	WI-70004063	\$72.70	166
Verbrugge Bone Holding Forceps 6" Self Centering with Speed Lock Premium Quality	WI-70004062	\$41.76	167
Virtus Splinter Forceps 6" (15.2cm) Straight Serrated Tips	WI-70000486	\$9.99	167
Virtus Splinter Forceps 6" Straight Serrated Triangle Tips	WI-70000810	\$10.99	168
Volume Lashes Tweezers Best Volum Lashes Tweezers for Professional Eyelash Extension	WI-70003581	\$9.99	168
Whole Sale Eyelash Tweezers Strong Curved and Straight Eye Lash Extension Tweezers	WI-70003579	\$15.99	169
Wilde Blakesley Ethmoid Forceps 5" Shaft 45° Up 4×12mm Pointed Fenestrated Cup	WI-70000624	\$29.99	169
Wilde Blakesley Ethmoid Forceps 5" Shaft 4×12mm Straight Pointed Fenestrated Cup	WI-70000625	\$29.99	170
Wilde Blakesley Ethmoid Forceps 5" Shaft 45° Up 4×12mm Pointed Fenestrated Cup	WI-	\$29.99	170

**WI-70000386****Hartmann Mosquito Forceps 3.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hemostat Mosquito Forceps 3.5\"/>

Full gold hartmann mosquito forceps feature fine tips, short, and fully serrated jaws used as a clamping instrument to halt blood flow so that cauterization or ligation can be performed. "Hartmann" hemostat mosquito forceps are used as both a hemostat and for fine tissue dissection in shallow wounds and procedures. They are also used for holding small, fine sutures. Hartmann Mosquito Forceps are shorter and heavier than halsted mosquito forceps. Ratcheted handles can be locked in multiple positions in order to maintain variable

levels of constant pressure. Forged from durable stainless steel and have large finger loops with serrated jaws that are ideal for grasping and holding small objects securely. Easy to clean and store in any tackle box, backpack, belt, pocket + other outdoor gear.

FULL GOLD HARTMANN MOSQUITO FORCEPS FEATURES

- 3.5\"/>
- Fully Gold Plated.
- Fine Serrated Jaws for Firm Grip and Minimal Tissue Trauma.
- Ratchet Lock Mechanism for Secure Locking.
- Finger Ring Handle for Comfortable, Secure Grip.
- Available in Straight and Curved Profiles

**WI-70002625****Hartmann Mosquito Forceps 3.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Mosquito Hemostat Locking Forceps 3.5\"/>

Hartmann Mosquito forceps feature fine tips, short, and fully serrated jaws, these Hartmann Mosquito Forceps can be used as a clamping instrument to halt blood flow so that cauterization or ligation can be performed. Hartmann Mosquito Forceps are used as both a hemostat and for fine tissue dissection in shallow wounds and procedures. They are also used for holding small, fine sutures. Hartmann Mosquito Forceps are shorter and heavier than Halsted Mosquito Forceps.

HARTMANN MOSQUITO FORCEPS FEATURES

3 1/2\"/>

THREE LOCKING POSITIONS: The high grip and locking ratchet mechanism reduce the possibility of the vessels or collection of tissues sliding or moving during use. This hemostat has handles that can be held in place by a locking mechanism. The locking mechanism consists of a series of 3 interlocking teeth that allow the user to adjust the clamping force of the pliers. Once





WI-70003528

Hartmann Mosquito Forceps 3.5" Straight

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Mosquito Hemostat Locking Forceps 3.5" (8.9cm) Stainless Steel – Multipurpose Ideal Hemostats for Nurses, Fishing Forceps, Crafts and Hobby


Hartmann Mosquito forceps feature fine tips, short, and fully serrated jaws, these Hartmann Mosquito Forceps can be used as a clamping instrument to halt blood flow so that cauterization or ligation can be performed. Hartmann Mosquito Forceps are used as both a hemostat and for fine tissue dissection in shallow wounds and procedures. They are also used for holding small, fine sutures. Hartmann Mosquito Forceps are shorter and heavier than Halsted Mosquito Forceps.

HARTMANN MOSQUITO FORCEPS FEATURES

3 1/2" STRAIGHT MOSQUITO

FORCEPS: Hartman Mosquito Forceps are used as both a hemostat and for blunt tissue dissection. These straight forceps are only 3.5" long, making them easy to work with in tight spaces. Ratcheted handles can be locked in multiple positions in order to maintain variable levels of constant pressure. Delicate hemostats have thin jaws with fine serrations

THREE LOCKING POSITIONS: The high grip and locking ratchet mechanism reduce the possibility of the vessels or collection of tissues sliding or moving during use. This hemostat has handles that can be held in place by a locking mechanism. The locking mechanism consists of a series of 3 interlocking teeth that allow the user to adjust the clamping force of the pliers. Once



WI-70003260

Hartmann-Wullstein Ear Forceps 8.5cm

\$46.86 per piece

Trade price from \$43.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hartmann-Wullstein Ear Forceps 8.5cm is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$46.86, falling to \$43.11 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-26032810

HEANEY FORCEPS 20.0cm Stainless Steel Surgical Instruments

\$8.00 per piece

Trade price from **\$4.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

HEANEY FORCEPS 20.0cm Stainless Steel Surgical Instruments is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 20.0cm size shown is our standard; other sizes in the same pattern are made to order.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-26019784

Heaney Hysterectomy Forceps 203mm Stainless Steel Surgical Instruments

\$8.00 per piece

Trade price from **\$4.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

This is our pattern of Heaney Hysterectomy Forceps 203mm Stainless Steel Surgical Instruments, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 203mm size shown is our standard; other sizes in the same pattern are made to order.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$4.00 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001292**

Heaney Hysterectomy Forceps 8.25\"/>

\$14.99 per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Heaney Hysterectomy Forceps 8.25\"/>




Heaney Hysterectomy Forceps are widely used to clamp off a ligament on both sides of the uterus during a hysterectomy procedure. The gynaecology and obstetrics specialists use Heaney Hysterectomy Forceps to grasp, mobilize and clamp the blood vessels and ligaments that attach the uterus to the pelvic cavity.

Heaney Forceps features curved jaws that conform to the anatomy of the pelvic cavity. These jaws include a cross-serration with single and double tooth patterns that provides a strong gripping, minimizes trauma and

enhance the clamping action. The finger ring handle with ratchet lock system allows self-locking the jaws so surgeons can apply different levels of pressure to the clamped tissues with maximum control.

HEANEY HYSTERECTOMY FORCEPS FEATURES

- 8.25 Inch (21cm) Long.
- Curved Fully Serrated Double Tooth Jaws.
- Cross-serrated Pattern for Enhanced Gripping.
- Boxlock System Promoting Smooth Jaw Movements.
- Ratchet Locking Mechanism to Grasp the Ligament Firmly.
- Ergonomic Finger Rings Ensuring Maximum Control.

**WI-70000068**

Heaney Hysterectomy Forceps 8.25\"/>

\$14.99 per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Heaney Hysterectomy Forceps 8.25\"/>


Heaney Hysterectomy Forceps are widely used to clamp off a ligament on both sides of the uterus during a hysterectomy procedure. The gynaecology and obstetrics specialists use Heaney Hysterectomy Forceps to grasp, mobilize and clamp the blood vessels and ligaments that attach the uterus to the pelvic cavity.

Straight Jaws Heaney Hysterectomy forceps provide excellent access and clear sightlines, with controlled clamping depths beneficial for various surgical techniques. They are available

in double-toothed and single-toothed variations, ensuring a secure grip on tissues and effectively interlocking to hold thick ligaments, preventing slippage during manipulation or cutting. The finger ring handle with ratchet lock system allows self-locking the jaws so surgeons can apply different levels of pressure to the clamped tissues with maximum control.

HEANEY HYSTERECTOMY FORCEPS FEATURES

- 8.25 Inch (21cm) Long.
- Straight Fully Serrated Double-Toothed Jaws.
- Cross-serrated Pattern for Enhanced Gripping.
- Boxlock System Promoting Smooth Jaw Movements.



**WI-70001293**

Heaney Hysterectomy Forceps 8.25" Single Tooth Curved Jaws

\$12.99 per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Heaney Hysterectomy Forceps 8.25" (21cm) Curved Heavy Serrated Single Tooth Jaws German Grade Stainless Steel Gynaecology and Obstetrics Instrument.

Heaney Hysterectomy Forceps are widely used to clamp off a ligament on both sides of the uterus during a hysterectomy procedure. The gynaecology and obstetrics specialists use Heaney Hysterectomy Forceps to grasp, mobilize and clamp the blood vessels and ligaments that attach the uterus to the pelvic cavity.

Heaney Forceps features curved jaws that conform to the anatomy of the pelvic cavity. These jaws include a cross-serration with single and double tooth patterns that provides a strong gripping, minimizes trauma and

enhance the clamping action. The finger ring handle with ratchet lock system allows self-locking the jaws so surgeons can apply different levels of pressure to the clamped tissues with maximum control.

HEANEY HYSTERECTOMY FORCEPS FEATURES

- 8.25 Inch (21cm) Long.
- Curved Fully Serrated Single-Toothed Jaws.
- Cross-serrated Pattern for Enhanced Gripping.
- Boxlock System Promoting Smooth Jaw Movements.
- Ratchet Locking Mechanism to Grasp the Ligament Firmly.
- Ergonomic Finger Rings Ensuring Maximum Control.

WI-70000069

Heaney Hysterectomy Forceps 8.25" Single Tooth Straight Jaws

\$12.99 per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Heaney Hysterectomy Forceps 8.25" (21cm) Straight Serrated Single Tooth Jaws with Ratcheted Finger Ring Handle – German Grade Stainless Steel Gynaecology and Obstetrics Instrument.


Heaney Hysterectomy Forceps are widely used to clamp off a ligament on both sides of the uterus during a hysterectomy procedure. The gynaecology and obstetrics specialists use Heaney Hysterectomy Forceps to grasp, mobilize and clamp the blood vessels and ligaments that attach the uterus to the pelvic cavity.

Straight Jaws Heaney Hysterectomy forceps provide excellent access and clear sightlines, with controlled clamping depths beneficial for various surgical techniques. They are available

in double-toothed and single-toothed variations, ensuring a secure grip on tissues and effectively interlocking to hold thick ligaments, preventing slippage during manipulation or cutting. The finger ring handle with ratchet lock system allows self-locking the jaws so surgeons can apply different levels of pressure to the clamped tissues with maximum control.

HEANEY HYSTERECTOMY FORCEPS FEATURES

- 8.25 Inch (21cm) Long.
- Straight Fully Serrated Single-Toothed Jaws.
- Cross-serrated Pattern for Enhanced Gripping.
- Boxlock System Promoting Smooth Jaw Movements.

**WI-70000066****Heaney-Ballentine Hysterectomy Forcep 8.5\"/>****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Heaney-Ballentine Hysterectomy Forceps 8.5" (21.6cm) Curved Longitudinal Serration Double Tooth/Groove Jaws with Ratcheted Finger Ring Handle – Premium Quality German Grade Stainless Steel Gynaecology and Obstetrics Instrument for Hysterectomy Procedures.

Heaney-Ballentine hysterectomy forceps are specialized gynecological instruments used to clamp ligaments at both ends of the uterus, facilitating secure excision while minimizing bleeding. They feature longitudinal serrations along the jaws, accompanied single tooth and double toothed and come in curved and straight profiles, which aid in maneuverability within the pelvic cavity. Varying lengths and groove configurations allow surgeons

to maneuver efficiently and reach deep into the pelvic cavity. The ratcheted ring handle locks the clamp securely to prevent hemorrhaging. Heaney-Ballentine forceps are made from premium quality stainless steel, offering a fully autoclavable and reliable option for OB/GYN clinics and professional settings, ensuring high quality and product performance during surgical procedures.

HEANEY-BALLENTINE HYSTERECTOMY FORCEPS FEATURES

- 8.5 Inch (21.6cm) Long.
- Longitudinal Serrated Double-Toothed Jaws Improve Grip & Hemostatic Control.
- Curved Pattern Enhances Maneuverability Around Delicate

**WI-70000067****Heaney-Ballentine Hysterectomy Forcep 8.5\"/>****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Heaney-Ballentine Hysterectomy Forceps 8.5" (21.6cm) Curved Longitudinal Serration Single Tooth/Groove Jaws with Ratcheted Finger Ring Handle – Premium Quality German Grade Stainless Steel Gynaecology and Obstetrics Instrument for Hysterectomy Procedures.

Heaney-Ballentine hysterectomy forceps are specialized gynecological instruments used to clamp ligaments at both ends of the uterus, facilitating secure excision while minimizing bleeding. They feature longitudinal serrations along the jaws, accompanied single tooth and double toothed and come in curved and straight profiles, which aid in maneuverability within the pelvic cavity. Varying lengths and groove configurations allow surgeons

to maneuver efficiently and reach deep into the pelvic cavity. The ratcheted ring handle locks the clamp securely to prevent hemorrhaging. Heaney-Ballentine forceps are made from premium quality stainless steel, offering a fully autoclavable and reliable option for OB/GYN clinics and professional settings, ensuring high quality and product performance during surgical procedures.

HEANEY-BALLENTINE HYSTERECTOMY FORCEPS FEATURES

- 8.5 Inch (21.6cm) Long.
- Longitudinal Serrated Single-Toothed Jaws Improve Grip & Hemostatic Control.
- Curved Pattern Enhances Maneuverability Around Delicate



**WI-70000064****Heaney-Ballentine Hysterectomy Forcep 8.5\"/>****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Heaney-Ballentine Hysterectomy Forceps 8.5" (21.6cm) Straight Longitudinal Serration Double Tooth/Groove Jaws with Ratcheted Finger Ring Handle – Premium Quality German Grade Stainless Steel Gynaecology and Obstetrics Instrument for Hysterectomy Procedures.

Heaney-Ballentine hysterectomy forceps are specialized gynecological instruments used to clamp ligaments at both ends of the uterus, facilitating secure excision while minimizing bleeding. They feature longitudinal serrations along the jaws, accompanied single tooth and double toothed and come in curved and straight profiles, which aid in maneuverability within the pelvic cavity. Varying lengths and groove configurations allow surgeons

to maneuver efficiently and reach deep into the pelvic cavity. The ratcheted ring handle locks the clamp securely to prevent hemorrhaging. Heaney-Ballentine forceps are made from premium quality stainless steel, offering a fully autoclavable and reliable option for OB/GYN clinics and professional settings, ensuring high quality and product performance during surgical procedures.

HEANEY-BALLENTINE HYSTERECTOMY FORCEPS FEATURES

- 8.5 Inch (21.6cm) Long.
- Longitudinal Serrated Double-Toothed Jaws Improve Grip & Hemostatic Control.
- Straight Pattern Enhances Maneuverability Around Delicate

**WI-70000065****Heaney-Ballentine Hysterectomy Forcep 8.5\"/>****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Heaney-Ballentine Hysterectomy Forceps 8.5" (21.6cm) Straight Longitudinal Serration Single Tooth/Groove Jaws with Ratcheted Finger Ring Handle – Premium Quality German Grade Stainless Steel Gynaecology and Obstetrics Instrument for Hysterectomy Procedures.

Heaney-Ballentine hysterectomy forceps are specialized gynecological instruments used to clamp ligaments at both ends of the uterus, facilitating secure excision while minimizing bleeding. They feature longitudinal serrations along the jaws, accompanied single tooth and double toothed and come in curved and straight profiles, which aid in maneuverability within the pelvic cavity. Varying lengths and groove configurations allow surgeons

to maneuver efficiently and reach deep into the pelvic cavity. The ratcheted ring handle locks the clamp securely to prevent hemorrhaging. Heaney-Ballentine forceps are made from premium quality stainless steel, offering a fully autoclavable and reliable option for OB/GYN clinics and professional settings, ensuring high quality and product performance during surgical procedures.

HEANEY-BALLENTINE HYSTERECTOMY FORCEPS FEATURES

- 8.5 Inch (21.6cm) Long.
- Longitudinal Serrated Single-Toothed Jaws Improve Grip & Hemostatic Control.
- Straight Pattern Enhances Maneuverability Around Delicate





WI-70002510

Hemostat Forceps Straight & Curved with Grip Handle

\$8.99 per piece

Trade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Set of 2 Hemostat Forceps Straight & Curved with Red Grip Handle Serrated Multipurpose Locking Tweezer Clamps – Stainless Steel – Hemostats for Nurses, Fishing Forceps, Crafts and Hobby

HEMOSTAT FORCEPS FEATURES

SET OF 2 HEMOSTAT FORCEPS: This package includes one straight + one curved Serrated Hemostat Forceps. Hemostat is a multipurpose instrument used to clamp off blood vessels, grasp loose objects, and remove small root tips. Handy to have in an office, hospital, or practice – they are versatile instruments.

PREMIUM QUALITY HEMOSTATIC FORCEPS: These are Premium Grade quality with an ultra-smooth finish and fine craftsmanship.

NO SLIP GRIP HANDLE DESIGN: The ergonomic handle design provides extra comfort and reduces user fatigue while performing many different procedures.

RUST AND CORROSION RESISTANT: These Hemostat Forceps will not rust or tarnish on you, they are ultrasonic cleaned and hygienic. Built to last and serve you for many years.

SERRATED JAWS FOR PRECISION USE: The serrated jaws create an even more secure grip, allowing for easier operation and precision use. 100% Satisfaction Guaranteed.

MULTI-PURPOSE FISHING/ MEDICAL/ HOBBY FORCEPS: Waqar International is the leading source for hemostats straight and curved in many sizes. A versatile tool that can be used for so many activities, from efficient and safe



WI-70003599

High Quality Anti Static Stain Less Steel Satin Eyelash Tweezer Eyebrow Electric Tweezers

\$14.75 per piece

Trade price from **\$13.57**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

High Quality Anti Static Stain Less Steel Satin Eyelash Tweezer Eyebrow Electric Tweezers is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$14.75, falling to \$13.57 at 21 pieces – about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000382****Hoff Towel Clamp Forceps 5.25"****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hoff Towel Clamp Forceps 5.25" (13.3cm) Curved Jaws with Rounded Cross-Serration Tips Flat Finger Ring Ratcht Handle – Premium Quality Stainless Steel Non-Perforated Surgical Towel Clamps.

Hoff Non-Perforated Towel Clamps, made from premium-grade German surgical stainless material, are versatile tools for tying down towels to wrap around the surgical site. This Towel Clamp forceps features bent tips and jaws, flat handle, and cross-serrations rounded tips for secure hold and easy towel movement. The tightened finger ring handles assured firm grasping of towel or drapes. These forceps should never be used on the patient's skin. Made from premium-grade German

surgical stainless, making it lightweight and rustproof.

HOFF NON PERFORATED TOWEL CLAMP FEATURES

- 5.25" (13.3cm) Overall Length.
- Curved Jaws with Rounded Cross-Serration Tps
- Straight Shanks.
- Flat Finger Ring Handles.
- Ratcht Lock Ensure Firm Grasping.
- Not for Use on Patient's Skin.
- Reusable.
- Made with High-quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Hoff Towel Clamp Forceps 5.25" (13.3cm)

TECHNICAL INFORMATION**WI-70002922****Hoff Towel Forceps 3.5"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hoff Towel Clamp Forceps 3.5" (8.9cm) Curved Jaws with Rounded Cross-Serration Tips Flat Finger Ring Ratcht Handle – Premium Quality Stainless Steel Non-Perforated Surgical Towel Clamps.

Hoff Non-Perforated Towel Clamps, made from premium-grade German surgical stainless material, are versatile tools for tying down towels to wrap around the surgical site. This Towel Clamp forceps features bent tips and jaws, flat handle, and cross-serrations rounded tips for secure hold and easy towel movement. The tightened finger ring handles assured firm grasping of towel or drapes. These forceps should never be used on the patient's skin. Made from premium-grade German

surgical stainless, making it lightweight and rustproof.


HOFF NON PERFORATED TOWEL CLAMP FEATURES

- 3.5" (8.9cm) Overall Length.
- Curved Jaws with Rounded Cross-Serration Tps
- Straight Shanks.
- Flat Finger Ring Handles.
- Ratcht Lock Ensure Firm Grasping.
- Not for Use on Patient's Skin.
- Reusable.
- Made with High-quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Hoff Towel Clamp Forceps 3.5" (8.9cm)

TECHNICAL INFORMATION



WI-34850411

Hot Sales Debakey Forceps Straight Debakey Vascular Tissue Forceps

\$17.00 per piece

Trade price from **\$11.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hot Sales Debakey Forceps Straight Debakey Vascular Tissue Forceps — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a straight working end for direct, in-line access.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$17.00, falling to \$11.00 at 25 pieces — about 35% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34796578

Hot Sales Stainless Steel Debakey Forceps Straight Vascular Tissue

\$17.00 per piece

Trade price from **\$11.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hot Sales Stainless Steel Debakey Forceps Straight Vascular Tissue — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a straight working end for direct, in-line access.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$17.00, falling to \$11.00 at 25 pieces — about 35% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002014

Hulka Controlling Tenaculum Uterine Forceps 11"

\$11.99 per piece

Trade price from \$11.03
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



ABOUT THIS INSTRUMENT

Hulka Controlling Tenaculum Uterine Forceps 11" is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$11.99, falling to \$11.03 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.


Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

HULKA CONTROLLING TENACULUM



Material Stainless steel — we state the grade on your quotation

Reusable Yes, autoclavable



WI-42586281

Hysterectomy Forceps Surgical Stainless Steel Medical

\$6.50 per piece

Trade price from \$5.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hysterectomy Forceps Surgical Stainless Steel Medical is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$6.50, falling to \$5.00 at 10 pieces — about 23% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000913****Iris Tissue Forceps 4" 1x2 Teeth Fully Curved****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Iris Tissue Forceps 4" (10.2cm) Fully Curved Standard 0.8mm 1x2 Teeth with Grip Handle Premium Quality Stainless Steel Eye Surgery Surgical Ophthalmic Instruments.

Iris Tissue Forceps are specially designed to deal with iris tissue during ophthalmologic procedures. The serration aid in gripping iris tissue for manipulation or extraction. The narrow jaws of iris thumb forceps allow for fine manipulation of structures like the iris of the eye. These forceps are 10cm long with 0.8mm tips, and they come with straight or curved tips. These Iris Forceps have teeth (1x2). One blade has one tooth and the other has two. The blades align so that the teeth fit together. If the tissue to be grasped is

slippery or hard to hold, surgical forceps with teeth offer a firmer grasp.

Iris Forceps also used in veterinary ophthalmic procedures which are delicate and critical. The eye organ is considered fragile, requires extreme care during surgery. These corrective eye surgeries turn out to be successful when high-quality instruments play their vital role. Multiple half curved, fully curved, and straight are available to meet the surgical specialties for the veterinarians/ophthalmologists.

IRIS TISSUE FORCEPS FEATURES

- Crafted with 4" (10.2cm) Overall Length.
- Fully Curved Standard Pattern Profile.
- 1x2 Teeth for Safe Tissue Holding.

WI-70000914**Iris Tissue Forceps 4" 1x2 Teeth Half Curved****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Iris Tissue Forceps 4" (10.2cm) Half Curved Standard 0.8mm 1x2 Teeth with Grip Handle Premium Quality Stainless Steel Eye Surgery Surgical Ophthalmic Instruments.

Iris Tissue Forceps are specially designed to deal with iris tissue during ophthalmologic procedures. The serration aid in gripping iris tissue for manipulation or extraction. The narrow jaws of iris thumb forceps allow for fine manipulation of structures like the iris of the eye. These forceps are 10cm long with 0.8mm tips, and they come with straight or curved tips. These Iris Forceps have teeth (1x2). One blade has one tooth and the other has two. The blades align so that the teeth fit together. If the tissue to be grasped is

slippery or hard to hold, surgical forceps with teeth offer a firmer grasp.

Iris Forceps also used in veterinary ophthalmic procedures which are delicate and critical. The eye organ is considered fragile, requires extreme care during surgery. These corrective eye surgeries turn out to be successful when high-quality instruments play their vital role. Multiple half curved, fully curved, and straight are available to meet the surgical specialties for the veterinarians/ophthalmologists.

IRIS TISSUE FORCEPS FEATURES

- Crafted with 4" (10.2cm) Overall Length.
- Half Curved Standard Pattern Profile.
- 1x2 Teeth for Safe Tissue Holding.
- 0.8mm Wide Tips.

**WI-70000915****Iris Tissue Forceps 4" 1x2 Teeth Straight****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Iris Tissue Forceps 4" (10.2cm) Straight Standard 0.8mm 1x2 Teeth with Grip Handle Premium Quality Stainless Steel Eye Surgery Surgical Ophthalmic Instruments.

Iris Tissue Forceps are specially designed to deal with iris tissue during ophthalmologic procedures. The serration aid in gripping iris tissue for manipulation or extraction. The narrow jaws of iris thumb forceps allow for fine manipulation of structures like the iris of the eye. These forceps are 10cm long with 0.8mm tips, and they come with straight or curved tips. These Iris Forceps have teeth (1x2). One blade has one tooth and the other has two. The blades align so that the teeth fit together. If the tissue to be grasped is

slippery or hard to hold, surgical forceps with teeth offer a firmer grasp.

Iris Forceps also used in veterinary ophthalmic procedures which are delicate and critical. The eye organ is considered fragile, requires extreme care during surgery. These corrective eye surgeries turn out to be successful when high-quality instruments play their vital role. Multiple half curved, fully curved, and straight are available to meet the surgical specialties for the veterinarians/ophthalmologists.

IRIS TISSUE FORCEPS FEATURES

- Crafted with 4" (10.2cm) Overall Length.
- Straight Standard Pattern Profile.
- 1x2 Teeth for Safe Tissue Holding.
- 0.8mm Wide Tips.

**WI-70001700****Ivory Light Rubber Dam Clamp Forceps 6.7"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Ivory Light Weight Rubber Dam Clamp Forceps Plier 6.7" (17cm) Restorative Endodontic tools

The Ivory rubber dam forceps are able to reach around the bow of any clamp without tilting, which is the most important factor when reversing a clamp. The straight head design of Ivory rubber dam clamp forceps allows easy gripping and placement of any clamp size. The ivory pattern has stabilizers that prevent the clamp from rotating on the beaks. It limits the use of Ivory forceps to teeth that are within a range of normal angulation..

IVORY DAM CLAMP FORCEPS FEATURES

- Ivory Straight Head Design.

- Smooth Spring Enables Quick and Easy Placement and Removal of Clamps.
- Locking Ring Mechanism to Hold Forceps in Any Position.
- Long Shanks with Locking Mechanism.
- Completely Sterilizable.
- Made with High Quality Stainless Steel

PACKAGE INCLUDES

- 1 Pc Ivory Light Weight Rubber Dam Clamp Forceps 6.7" (17cm)

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Endodontic Dental Instruments**APPLICATION:** Place and Remove Rubber Dam Clamps

**WI-70000199****Jackson Endoscopic Forceps 24\"/>****\$109.99** per pieceTrade price from **\$101.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Jackson Endoscopic Forceps 24\"/>

Jackson Biopsy Forceps facilitate precise tissue sampling in ENT and endoscopic procedures, featuring an extended length for access to deep areas. Jackson Endoscopic Forceps are specialized for grasping and manipulating tissues or retrieving foreign bodies during ENT (ear, nose, and throat) and endoscopic procedures, with long, thin shafts designed for reaching deep anatomical regions. They are commonly used in laryngoscopy, bronchoscopy, and bronchoscopic biopsies.

Features circular cup jaws used to take samples in areas like the larynx, trachea, and bronchi providing controlled handling and dependable performance. Jackson biopsy forceps are designed with large finger ring handles for better grip and reduced finger strain. They come in both straight and angled profiles, with tip shapes including cupped or alligator jaws. Made from high-quality stainless steel, these forceps are sterilizable and reusable, ensuring long-lasting service.

JACKSON ENDOSCOPIC FORCEPS FEATURES

- Size: 24 Inch (61cm) Long Shaft
- 4mm Diameter Circular Cup-shaped Jaws For Reliable Biopsies.
- Straight and Angled Profiles.

**WI-70000200****Jackson Endoscopic Forceps 24\"/>****\$109.99** per pieceTrade price from **\$101.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Jackson Endoscopic Forceps 24\"/>

Jackson Biopsy Forceps facilitate precise tissue sampling in ENT and endoscopic procedures, featuring an extended length for access to deep areas. Jackson Endoscopic Forceps are specialized for grasping and manipulating tissues or retrieving foreign bodies during ENT (ear, nose, and throat) and endoscopic procedures, with long, thin shafts designed for reaching deep anatomical regions. They are commonly used in laryngoscopy, bronchoscopy, and bronchoscopic biopsies.

Features circular cup jaws used to take samples in areas like the larynx, trachea, and bronchi providing controlled handling and dependable performance. Jackson biopsy forceps are designed with large finger ring handles for better grip and reduced finger strain. They come in both straight and angled profiles, with tip shapes including cupped or alligator jaws. Made from high-quality stainless steel, these forceps are sterilizable and reusable, ensuring long-lasting service.

JACKSON ENDOSCOPIC FORCEPS FEATURES

- Size: 24 Inch (61cm) Long Shaft
- 4mm Diameter Circular Cup-shaped Jaws For Reliable Biopsies.
- Straight and Angled Profiles.





WI-70003636

Jansen Ear Dressing Forceps 6.25"

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Jansen Ear Dressing Forceps 6.25" (15.9cm) Bayonet Style Serrated Tips ENT Medical Otology Surgical Instruments

Jansen ear forceps are bayonet grasping devices that ear doctors use for aural procedures, like inserting gauze packs and removing foreign bodies. Jansen ear forceps assist doctors to treat ear diseases. Doctors use these forceps to remove foreign bodies that close the ear canal.

The shaft is long in order to improve access into deep spaces. Moreover, the bayonet design gives a more direct vision of the jaws for accuracy. The serrations on the handles offer a steadfast grip. On the other hand, the serrations in the jaws hold a wide range of materials securely. We manufacture

ear forceps from the finest class of stainless steel that you can depend on for years of efficient performance.

JANSEN EAR FORCEPS FEATURES

- Bayonet Design For Optimal Accuracy.
- Serrated Handles For Non-slip Grips.
- Long Slender Jaws For Deep Access.
- Lightweight
- Made with High Quality Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Jansen Ear Forceps 5.5" (14cm) Bayonet Style Serrated Tips

TECHNICAL INFORMATION

MATERIAL: Stainless Steel



WI-47140991

Jar Forceps Stainless Steel Medical Instrument for Precise Handling

\$7.90 per piece

Trade price from \$4.90
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Jar Forceps Stainless Steel Medical Instrument for Precise Handling is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$7.90, falling to \$4.90 at 25 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47170605

Jar Forceps Surgical Instrument for Secure Grasping Handling

\$8.00 per piece

Trade price from **\$5.00**
Minimum order 12 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Jar Forceps Surgical Instrument for Secure Grasping Handling for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$8.00, falling to \$5.00 at 12 pieces — about 38% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47135996

Jar Surgical Forceps Reusable Medical Instrument for Controlled Grasping

\$7.50 per piece

Trade price from **\$4.50**
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Jar Surgical Forceps Reusable Medical Instrument for Controlled Grasping — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$7.50, falling to \$4.50 at 50 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-43520264

Jewelers Forceps Fine Tip Stainless Steel Precision Handling

\$6.00 per piece

Trade price from **\$4.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Jewelers Forceps Fine Tip Stainless Steel Precision Handling for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a fine pattern for delicate tissue.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$6.00, falling to \$4.00 at 10 pieces — about 33% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-43522144

Jewelers Forceps Reusable Stainless Steel Precision Instrument

\$6.00 per piece

Trade price from **\$4.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Jewelers Forceps Reusable Stainless Steel Precision Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$6.00, falling to \$4.00 at 10 pieces — about 33% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000270****Kelly Forceps 5.5\"/>****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Kelly Hemostat 5.5" (14cm) Curved Full Black Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Serrated Jaws, Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobbies.

Full black Kelly hemostat is a 5.5-inch stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. The matte black coating minimizes glare and enhances corrosion resistance.

These hemostat clamps feature a locking forceps mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

BLACK KELLY HEMOSTAT FORCEPS FEATURES

- Size: 5.5" (14cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- Matte Full Black Finish Reduces Glare and Enhances Corrosion Resistance.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Fully Autoclavable and Reusable.

**WI-70003520****Kelly Forceps 5.5\"/>****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostat 5.5" (14cm) Straight Full Black Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Serrated Jaws, Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobbies.

Full black Kelly hemostat is a 5.5-inch stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. The matte black coating minimizes glare and enhances corrosion resistance.

These hemostat clamps feature a locking forceps mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

BLACK KELLY HEMOSTAT FORCEPS FEATURES

- Size: 5.5" (14cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- Matte Full Black Finish Reduces Glare and Enhances Corrosion Resistance.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Fully Autoclavable and Reusable.





WI-27077854

Kelly Forceps Stainless Steel Curved Straight Full-Body

\$3.00 per piece

Trade price from \$1.90
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kelly Forceps Stainless Steel Curved Straight Full-Body is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.00, falling to \$1.90 at 1 piece — about 37% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002609

Kelly Hemostat Forceps 5.5" Curved Fishing Plier

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostat 5.5" (14cm) Curved Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Serrated Jaws, Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobby


KELLY FORCEPS FEATURES

KELLY HEMOSTAT FORCEPS: Waqar International's Kelly Forceps 14cm Curved hemostat clamps are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely. **HIGH-QUALITY STAINLESS STEEL CONSTRUCTION:** Hot-drop forged from medical-grade stainless steel,

these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

THREE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once locked in place, even the smallest objects can be held securely in the serrated jaws with great force.

A MUST-HAVE TOOL FOR REMOVING HOOKS: These Fishing clamps are one of the most essential pieces of equipment that you can have with you on a day on the water. They are ideal



**WI-70000268****Kelly Hemostat Forceps 5.5\"/>****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostat 5.5\" (14cm) Curved Full Gold Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Serrated Jaws, Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobbies.

Full Gold Kelly hemostat is a 5.5-inch stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. Full gold color coating minimizes glare and enhances corrosion resistance.

These hemostat clamps feature a locking forceps mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

FULL GOLD KELLY FORCEPS FEATURES

- Size: 5.5\" (14cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- High-quality Full Gold Color Coating Reduces Glare and Enhances Corrosion Resistance.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Fully Autoclavable and Reusable.

**WI-70002608****Kelly Hemostat Forceps 5.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostat 5.5\" (14cm) Straight Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Serrated Jaws, Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobby

KELLY FORCEPS FEATURES

KELLY HEMOSTAT FORCEPS: Waqar International's Kelly Forceps 14cm Straight hemostat clamps are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel,

these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

THREE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once locked in place, even the smallest objects can be held securely in the serrated jaws with great force.

A MUST-HAVE TOOL FOR REMOVING HOOKS: These Fishing clamps are one of the most essential pieces of equipment that you can have with you on a day on the water. They are ideal



**WI-70000269****Kelly Hemostat Forceps 5.5\"/>****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostat 5.5\" (14cm) Straight Full Gold Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Serrated Jaws, Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobbies.

Full Gold Kelly hemostat is a 5.5-inch stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. Full gold color coating minimizes glare and enhances corrosion resistance.

These hemostat clamps feature a locking forceps mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

FULL GOLD KELLY FORCEPS FEATURES

- Size: 5.5\" (14cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- High-quality Full Gold Color Coating Reduces Glare and Enhances Corrosion Resistance.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Fully Autoclavable and Reusable.

**WI-70002607****Kelly Hemostat Forceps 6.25\"/>****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostat 6.25\" (15.9cm) Curved Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Serrated Jaws, Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobby

KELLY FORCEPS FEATURES


KELLY HEMOSTAT FORCEPS: Waqar International's Kelly Forceps 15.9cm Curved hemostat clamps are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel,

these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

THREE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once locked in place, even the smallest objects can be held securely in the serrated jaws with great force.

A MUST-HAVE TOOL FOR REMOVING HOOKS: These Fishing clamps are one of the most essential pieces of equipment that you can have with you on a day on the water. They are ideal



**WI-70002605**

Kelly Hemostat Forceps 6.25" Straight Fishing Plier

\$9.99 per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostat 6.25" (15.9cm) Straight Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Serrated Jaws, Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobby

KELLY FORCEPS FEATURES



KELLY HEMOSTAT FORCEPS: Waqar International's Kelly Forceps 15.9cm Straight hemostat clamps are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel,

these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

THREE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once locked in place, even the smallest objects can be held securely in the serrated jaws with great force.

A MUST-HAVE TOOL FOR REMOVING HOOKS: These Fishing clamps are one of the most essential pieces of equipment that you can have with you on a day on the water. They are ideal

**WI-70000266**

Kelly Hemostatic Forceps Curved 5.5" Blue PVC Handle

\$9.99 per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostatic Forceps 5.5" (14cm) Curved Serrated Jaws with Blue PVC Ratcheted Handle Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Stainless Steel Premium Quality Hemostat for Nurses, Fishing Forceps, Crafts and Hobbies.


5.5 Inch Kelly forceps is a stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. The non-slip PVC grip handles enhance comfort


and reduce hand fatigue during procedures.

These hemostat clamps feature a locking mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

PVC HANDLE KELLY FORCEPS FEATURES

- Size: 5.5" (14cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- Blue PVC Dipped Handle Reduces Hand Fatigue.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Reusable.



**WI-70000259****Kelly Hemostatic Forceps Curved 5.5" Red PVC Handle****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostatic Forceps 5.5" (14cm) Curved Serrated Jaws with Red PVC Ratcheted Handle Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Stainless Steel Premium Quality Hemostat for Nurses, Fishing Forceps, Crafts and Hobbies.



5.5 Inch Kelly forceps is a stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. The non-slip PVC grip handles enhance comfort

and reduce hand fatigue during procedures.

These hemostat clamps feature a locking mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

PVC HANDLE KELLY FORCEPS FEATURES

- Size: 5.5" (14cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- Red PVC Dipped Handle Reduces Hand Fatigue.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Reusable.

**WI-70000265****Kelly Hemostatic Forceps Curved 6.25" Blue PVC Handle****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rakin Kelly Forceps 6.25" (16cm) Curved Serrated Jaws with Blue PVC Ratcheted Handle Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Stainless Steel Premium Quality Hemostat for Nurses, Fishing Forceps, Crafts and Hobbies.


6.25 Inch Kelly Rankin forceps is a stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. The non-slip PVC grip handles enhance comfort

and reduce hand fatigue during procedures.

These hemostat clamps feature a locking mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

PVC HANDLE RANKIN KELLY FORCEPS FEATURES

- Size: 6.25" (16cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- Blue PVC Dipped Handle Reduces Hand Fatigue.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Reusable.



**WI-70000258****Kelly Hemostatic Forceps Curved 6.25\"/>****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Kelly Hemostatic Forceps 6.25" (16cm) Curved Serrated Jaws with Red PVC Ratcheted Handle Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Stainless Steel Premium Quality Hemostat for Nurses, Fishing Forceps, Crafts and Hobbies.

and reduce hand fatigue during procedures.

These hemostat clamps feature a locking mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

PVC HANDLE KELLY FORCEPS FEATURES

- Size: 6.25" (16cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- Red PVC Dipped Handle Reduces Hand Fatigue.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Reusable.



6.25 Inch Kelly forceps is a stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. The non-slip PVC grip handles enhance comfort

**WI-70000267****Kelly Hemostatic Forceps Straight 5.5\"/>****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Kelly Hemostatic Forceps 5.5" (14cm) Straight Serrated Jaws with Blue PVC Handle Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobbies.

and reduce hand fatigue during procedures.

These hemostat clamps feature a locking mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

PVC HANDLE KELLY FORCEPS FEATURES

- Size: 5.5" (14cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- Blue PVC Dipped Handle Reduces Hand Fatigue.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Reusable.



5.5 Inch Kelly forceps is a stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. The non-slip PVC grip handles enhance comfort

**WI-70000261****Kelly Hemostatic Forceps Straight 5.5" Red PVC Handle****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostatic Forceps 5.5" (14cm) Straight Serrated Jaws with Red PVC Handle Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobbies.


5.5 Inch Kelly forceps is a stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. The non-slip PVC grip handles enhance comfort

and reduce hand fatigue during procedures.

These hemostat clamps feature a locking mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

PVC HANDLE KELLY FORCEPS FEATURES

- Size: 5.5" (14cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- Red PVC Dipped Handle Reduces Hand Fatigue.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Reusable.

**WI-70000264****Kelly Hemostatic Forceps Straight 6.25" Blue PVC Handle****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostatic Forceps 6.25" (16cm) Straight Serrated Jaws with Blue PVC Handle Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobbies.

6.25 Inch Kelly forceps is a stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. The non-slip PVC grip handles enhance comfort

and reduce hand fatigue during procedures.

These hemostat clamps feature a locking mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

PVC HANDLE KELLY FORCEPS FEATURES

- Size: 6.25" (16cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- Blue PVC Dipped Handle Reduces Hand Fatigue.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Reusable.



**WI-70000260****Kelly Hemostatic Forceps Straight 6.25" Red PVC Handle****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Hemostatic Forceps 6.25" (16cm) Straight Serrated Jaws with Red PVC Handle Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobbies.

and reduce hand fatigue during procedures.

These hemostat clamps feature a locking mechanism with three interlocking teeth, ensuring a secure hold. They are easy to clean, store, and can be reused.

PVC HANDLE KELLY FORCEPS FEATURES

- Size: 6.25" (16cm) Overall Length.
- Available in Straight & Curved Profile.
- Serrated Jaws for Secure Grip.
- Ratcheted Locking Mechanism for Adjustable Pressure.
- Red PVC Dipped Handle Reduces Hand Fatigue.
- Suitable for Nurses, Fishing, Crafts, and Hobby Enthusiasts.
- Reusable.



6.25 Inch Kelly forceps is a stainless steel locking clamp featuring serrated jaws, utilized for clamping blood vessels, controlling bleeding, and grasping tissues in various medical, dental, and surgical procedures. Kelly forceps jaws design offers a secure grip on sutures and small instruments. Beyond medical uses, it is also popular in hobbies, fishing, electronics repair, and for removing pet ticks. The non-slip PVC grip handles enhance comfort

**WI-70000398****Kelly Placenta Forceps 12.5" Curved Serrated Oval Tip Jaws****\$13.99** per pieceTrade price from **\$12.87**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Placenta Forceps 12.5" (31.8cm) Curved Serrated Fenestrated Oval Tip Jaws with Straight Shank Non-locking OB/GYN Ovum Forceps – Premium Quality Stainless Steel Obstetrics and Gynecology Surgical Instruments.

Forceps, made of stainless steel, are non-locking finger ring handle forceps which requiring sterilization before reuse.

KELLY PLACENTA FORCEPS FEATURES

- 12.5" (31.8cm) Overall Length.
- Slightly Curved Serrated Jaws with Fenestrated Oval Tips.
- Long Straight Shank for Deep Access.
- Comfortable Ring Handles for Precise Control and Manipulation.
- Highly Durable.
- Sterilization-friendly
- High-quality Stainless Steel Construction

PACKAGE INCLUDE

- 1 Pc Kelly Placenta Forceps 12.5" (31.8cm) Curved Serrated Jaws



Kelly Placenta Forceps are commonly utilized post-childbirth to assist the mother in pushing out the placenta. These forceps are primarily used for placenta removal but can also be used in other OB/GYN procedures and for handling or manipulating tissues in other surgical settings. Kelly Placenta OB/GYN forceps have ring handles, slightly curved serrated jaws and fenestrated oval tips with straight shank, accommodating a wide range of cases with their straight and curved profiles. OB GYN Placenta Ovum



WI-70000397

Kelly Placenta Forceps 12.5" Straight Serrated Oval Tip Jaws

\$13.99 per piece

Trade price from **\$12.87**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Kelly Placenta Forceps 12.5" (31.8cm) Straight Serrated Fenestrated Oval Tip Jaws with Non-locking OB/GYN Ovum Forceps – Premium Quality Stainless Steel Obstetrics and Gynecology Surgical Instruments.

Kelly Placenta Forceps are commonly utilized post-childbirth to assist the mother in pushing out the placenta. These forceps are primarily used for placenta removal but can also be used in other OB/GYN procedures and for handling or manipulating tissues in other surgical settings. Kelly Placenta OB/GYN forceps have ring handles, straight serrated jaws and fenestrated oval tips with straight shank, accommodating a wide range of cases with their straight and curved profiles. OB GYN Placenta Ovum Forceps, made

of stainless steel, are non-locking finger ring handle forceps which requiring sterilization before reuse.

KELLY PLACENTA FORCEPS FEATURES

- 12.5" (31.8cm) Overall Length.
- Straight Serrated Jaws with Fenestrated Oval Tips.
- Long Straight Shank for Deep Access.
- Comfortable Ring Handles for Precise Control and Manipulation.
- Highly Durable.
- Sterilization-friendly
- High-quality Stainless Steel Construction

PACKAGE INCLUDE

- 1 Pc Kelly Placenta Forceps 12.5" (31.8cm) Straight Serrated Jaws



WI-24316734

Kelly Placenta Forceps Set 4-Piece OB/GYN Tools

\$20.00 per piece

Trade price from **\$8.00**
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kelly Placenta Forceps Set 4-Piece OB/GYN Tools is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$20.00, falling to \$8.00 at 10 sets — about 60% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-45053571

Keramo Grip Forceps 14.5 Cm TC Root Curved

\$20.00 per piece

Trade price from \$17.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Keramo Grip Forceps 14.5 Cm TC Root Curved for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: a curved working end for reach around structures, and tungsten carbide inserts for a longer cutting life.

The 14.5cm size shown is our standard; other sizes in the same pattern are made to order. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$20.00, falling to \$17.00 at 2 pieces — about 15% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001532

Kitchen Cooking Tongs 10" Round Tip BBQ Grill Tweezers – Gold

\$15.99 per piece

Trade price from \$14.71
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Inch Tweezer Tongs Gold Color
Stainless Steel Extra-Long Food Tongs, Kitchen Cooking Tweezers Tongs with Precision Serrated Rounded Tips for BBQ Grill Baking Decorating & Sea Food Multi-use

Tweezer tongs used for garnishing pastries, main courses, and mixed drinks without cutting or damaging your smallest or most delicate ingredients. Designed for detail work, these Long Food Tongs let you make sure the position of every component of your finished creation is perfect. They are ideal for precise plating applications, so you can add detail to any dish or drink at your hotel, event, bakery, bar, or upscale restaurant.

KITCHEN TWEEZER TONGS FEATURES:

TWEEZER TONGS FOR COOKING: You can use these 10 inch tweezers tongs to pick up meat, fruits, vegetables and so on, in addition, not only suitable for kitchen grilling, but also for taking aquarium weeding tweezers and more.
STUDY TO USE: These BBQ tweezer tongs are made of stainless steel, not easy to get rust, break or deform, can contact with food directly, safe to use and no harm to body, can serve for a long time use if keep well.
WHERE TO USE: You can use these tweezer tongs for cooking, roasting, serving, grilling and more usages, must-have items for kitchen and bring you more convenience when stand by oven, compare to shovel, it's easier to



**WI-70001534**

Kitchen Cooking Tongs 10" Straight Tip BBQ Grill Tweezers – Gold

\$14.99 per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Inch Tweezer Tongs Gold Color Stainless Steel Extra-Long Food Tongs, Kitchen Cooking Tweezer Tongs with Precision Serrated Tips for BBQ Grill Baking Decorating & Sea Food Multi-use

Tweezer tongs used for garnishing pastries, main courses, and mixed drinks without cutting or damaging your smallest or most delicate ingredients. Designed for detail work, these Long Food Tongs let you make sure the position of every component of your finished creation is perfect. They're ideal for precise plating applications, so you can add detail to any dish or drink at your hotel, event, bakery, bar, or upscale restaurant.

KITCHEN TWEEZER TONGS FEATURES:

PREMIUM MATERIAL, EASY ON YOUR HAND: Our cooking tweezers are made of high quality 100% food grade stainless steel, rust -proof, anti-corrosive, sturdy, durable and built to last. Precision serrated tooth is enough to be reliable to touch and even press, which prevents food from dropping, easy to grip and clamp fine. While remaining lightweight and avoids slipping and sliding.

MULTIFUNCTIONAL USE: Perfect for outdoor barbecue or grilling in the kitchen. Works well for turning and lifting appetizer cheese, fish, crab, chicken, turkey, burger, asparagus, vegetable salad, ice block, escargot, pastry, cake, and taco.

**WI-70001531**

Kitchen Cooking Tongs 12" Round Tip BBQ Grill Tweezers – Gold

\$18.99 per pieceTrade price from **\$17.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

12 Inch Tweezer Tongs Gold Color Stainless Steel Extra-Long Food Tongs, Kitchen Cooking Tweezers Tongs with Precision Serrated Rounded Tips for BBQ Grill Baking Decorating & Sea Food Multi-use

Tweezer tongs used for garnishing pastries, main courses, and mixed drinks without cutting or damaging your smallest or most delicate ingredients. Designed for detail work, these Long Food Tongs let you make sure the position of every component of your finished creation is perfect. They are ideal for precise plating applications, so you can add detail to any dish or drink at your hotel, event, bakery, bar, or upscale restaurant.

KITCHEN TWEEZER TONGS FEATURES:

TWEEZER TONGS FOR COOKING: You can use these 12 inch tweezers tongs to pick up meat, fruits, vegetables and so on, in addition, not only suitable for kitchen grilling, but also for taking aquarium weeding tweezers and more.

STUDY TO USE: These BBQ tweezers tongs are made of stainless steel, not easy to get rust, break or deform, can contact with food directly, safe to use and no harm to body, can serve for a long time use if kept well.

WHERE TO USE: You can use these tweezers tongs for cooking, roasting, serving, grilling and more usages, must-have items for kitchen and bring you more convenience when standing by oven, compare to shovel, it's easier to



**WI-70001533**

Kitchen Cooking Tongs 12" Straight Tip BBQ Grill Tweezers – Gold

\$17.99 per pieceTrade price from **\$16.55**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


12 Inch Tweezer Tongs Gold Color Stainless Steel Extra-Long Food Tongs, Kitchen Cooking Tweezer Tongs with Precision Serrated Tips for BBQ Grill Baking Decorating & Sea Food Multi-use

Tweezer tongs used for garnishing pastries, main courses, and mixed drinks without cutting or damaging your smallest or most delicate ingredients. Designed for detail work, these Long Food Tongs let you make sure the position of every component of your finished creation is perfect. They're ideal for precise plating applications, so you can add detail to any dish or drink at your hotel, event, bakery, bar, or upscale restaurant.

KITCHEN TWEEZER TONGS FEATURES:

PREMIUM MATERIAL, EASY ON YOUR HAND: Our cooking tweezer are made of high quality 100% food grade stainless steel, rust -proof, anti-corrosive, sturdy, durable and built to last. Precision serrated tooth is enough to be reliable to touch and even press, which is prevent food to drop, easy to grip and clamp fine. While remaining lightweight and avoid to slipping and sliding.

MULTIFUNCTIONAL USE: Perfect for outdoor barbecue or grilling in the kitchen. Works well for turning and lifting appetizer cheese, fish, crab, chicken, turkey, burger, asparagus, vegetable salad, ice block, escargot, pastry, cake, and taco.

**WI-70001538**

Kitchen Tweezer 10" Round Tip BBQ Grill Plating Long Food Tongs – Silver

\$14.99 per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Inch Tweezer Tongs, Stainless Steel Extra-Long Tweezers Tongs, Kitchen Cooking Tweezer Tongs with Precision Serrated Rounded Tips for BBQ Grill Baking Decorating & Sea Food Multi-use

Tweezer tongs used for garnishing pastries, main courses, and mixed drinks without cutting or damaging your smallest or most delicate ingredients. Designed for detail work, these Long Food Tongs let you make sure the position of every component of your finished creation is perfect. They are ideal for precise plating applications, so you can add detail to any dish or drink at your hotel, event, bakery, bar, or upscale restaurant.

KITCHEN TWEEZER TONGS FEATURES:

TWEEZER TONGS FOR COOKING: You can use these 10 inch tweezers tongs to pick up meat, fruits, vegetables and so on, in addition, not only suitable for kitchen grilling, but also for taking aquarium weeding tweezers and more.

STUDY TO USE: These BBQ tweezer tongs are made of stainless steel, not easy to get rust, break or deform, can contact with food directly, safe to use and no harm to body, can serve for a long time use if keep well.

WHERE TO USE: You can use these tweezer tongs for cooking, roasting, serving, grilling and more usages, must-have items for kitchen and bring you more convenience when stand by oven, compare to shovel, it's easier to



**WI-70001540**

Kitchen Tweezer 10" Straight Tip BBQ Grill Plating Long Food Tongs – Silver

\$13.99 per pieceTrade price from **\$12.87**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Inch Tweezer Tongs, Stainless Steel Extra-Long Tweezers Tongs, Kitchen Cooking Tweezer Tongs with Precision Serrated Tips for BBQ Grill Baking Decorating & Sea Food Multi-use

Tweezer tongs used for garnishing pastries, main courses, and mixed drinks without cutting or damaging your smallest or most delicate ingredients. Designed for detail work, these Long Food Tongs let you make sure the position of every component of your finished creation is perfect. They're ideal for precise plating applications, so you can add detail to any dish or drink at your hotel, event, bakery, bar, or upscale restaurant.

PREMIUM MATERIAL, EASY ON YOUR HAND: Our cooking tweezer are made of high quality 100% food grade stainless steel, rust -proof, anti-corrosive, sturdy, durable and built to last. Precision serrated tooth is enough to be reliable to touch and even press, which is prevent food to drop, easy to grip and clamp fine. While remaining lightweight and avoid to slipping and sliding.

MULTIFUNCTIONAL USE: Perfect for outdoor barbecue or grilling in the kitchen. Works well for turning and lifting appetizer cheese, fish, crab, chicken, turkey, burger, asparagus, vegetable salad, ice block, escargot, pastry, cake, and taco.

ERGONOMIC DESIGN: These tweezers design with comfortable ergonomic shape, precision tweezer food tongs to

KITCHEN TWEEZER TONGS FEATURES:

**WI-70001537**

Kitchen Tweezer 12" Round Tip BBQ Grill Plating Long Food Tongs – Silver

\$17.99 per pieceTrade price from **\$16.55**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

12 Inch Tweezer Tongs, Stainless Steel Extra-Long Tweezers Tongs, Kitchen Cooking Tweezer Tongs with Precision Serrated Rounded Tips for BBQ Grill Baking Decorating & Sea Food Multi-use

Tweezer tongs used for garnishing pastries, main courses, and mixed drinks without cutting or damaging your smallest or most delicate ingredients. Designed for detail work, these Long Food Tongs let you make sure the position of every component of your finished creation is perfect. They are ideal for precise plating applications, so you can add detail to any dish or drink at your hotel, event, bakery, bar, or upscale restaurant.

KITCHEN TWEEZER TONGS FEATURES:

TWEEZER TONGS FOR COOKING: You can use these 12 inch tweezers tongs to pick up meat, fruits, vegetables and so on, in addition, not only suitable for kitchen grilling, but also for taking aquarium weeding tweezers and more.

STUDY TO USE: These BBQ tweezer tongs are made of stainless steel, not easy to get rust, break or deform, can contact with food directly, safe to use and no harm to body, can serve for a long time use if keep well.

WHERE TO USE: You can use these tweezer tongs for cooking, roasting, serving, grilling and more usages, must-have items for kitchen and bring you more convenience when stand by oven, compare to shovel, it's easier to





WI-70001539

Kitchen Tweezer 12" Straight Tip BBQ Grill Plating Long Food Tongs – Silver

\$16.99 per piece

Trade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

12 Inch Tweezer Tongs, Stainless Steel Extra-Long Tweezers Tongs, Kitchen Cooking Tweezer Tongs with Precision Serrated Tips for BBQ Grill Baking Decorating & Sea Food Multi-use

Tweezer tongs used for garnishing pastries, main courses, and mixed drinks without cutting or damaging your smallest or most delicate ingredients. Designed for detail work, these Long Food Tongs let you make sure the position of every component of your finished creation is perfect. They're ideal for precise plating applications, so you can add detail to any dish or drink at your hotel, event, bakery, bar, or upscale restaurant.

PREMIUM MATERIAL, EASY ON YOUR HAND:

Our cooking tweezer are made of high quality 100% food grade stainless steel, rust -proof, anti-corrosive, sturdy, durable and built to last. Precision serrated tooth is enough to be reliable to touch and even press, which is prevent food to drop, easy to grip and clamp fine. While remaining lightweight and avoid to slipping and sliding.

MULTIFUNCTIONAL USE:

Perfect for outdoor barbecue or grilling in the kitchen. Works well for turning and lifting appetizer cheese, fish, crab, chicken, turkey, burger, asparagus, vegetable salad, ice block, escargot, pastry, cake, and taco.

ERGONOMIC DESIGN: These tweezers design with comfortable ergonomic shape, precision tweezer food tongs to



KITCHEN TWEEZER TONGS FEATURES:



WI-70003526

Kocher Forceps 14cm

\$12.99 per piece

Trade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kocher Forceps 14cm is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$12.99, falling to \$11.95 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35038993

Kocher Tissue Forceps Stainless Steel Surgical Grasping Instrument

\$2.30 per piece

Trade price from **\$1.40**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kocher Tissue Forceps Stainless Steel Surgical Grasping Instrument is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.30, falling to \$1.40 at 10 pieces — about 39% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44538586

Kocher Tissue Forceps Stainless Steel Surgical Instrument Set

\$2.30 per piece

Trade price from **\$0.80**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kocher Tissue Forceps Stainless Steel Surgical Instrument Set is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.30, falling to \$0.80 at 10 pieces — about 65% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001544****Laborde Tracheal Dilator Forceps 5.5\"/>****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Laborde Tracheal Dilator Forceps 5.5\"/>

Laborde Tracheal Dilator Forceps is an incredible surgical instrument assists surgeons to make critical tracheotomy incisions in order to treat upper airway obstructions. The device helps dilate or widen the trachea and facilitate for easy tube insertion and placement. In tracheotomy, surgeons make an opening in the anterior wall of the trachea to bypass obstructions in the upper airways or to manage chronic breathing problems.

Laborde Tracheal Dilator helps to widen the trachea and facilitate tube insertion with ease. This prevents

damaging the back of the airways with sharp force. The large finger rings prevent accidental slips and injuries. The unique three prongs design with round tips to control the tissue separation and make a uniform opening. The blunt tips also avoid tracheal tears and injury. Moreover, the jaws have narrow angled necks so you can view the working tips clearly

LABORDE TRACHEAL DILATOR FEATURES

- Finger Ring Gold Handles for Controlled Movements.
- Three prongs to ensure uniform tissue separation.
- Angled Jaws for Easy Surgical Access.
- Blunt Prongs for Safe Manipulations

**WI-70001549****Laborde Tracheal Dilator Forceps 5.5\"/>****\$8.99** per pieceTrade price from **\$8.27**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Laborde Tracheal Dilator Forceps 5.5\"/>


Laborde Tracheal Dilator Forceps is an incredible surgical instrument assists surgeons to make critical tracheotomy incisions in order to treat upper airway obstructions. The device helps dilate or widen the trachea and facilitate for easy tube insertion and placement. In tracheotomy, surgeons make an opening in the anterior wall of the trachea to bypass obstructions in the upper airways or to manage chronic breathing problems.

Laborde Tracheal Dilator helps to widen the trachea and facilitate tube insertion with ease. This prevents

damaging the back of the airways with sharp force. The large finger rings prevent accidental slips and injuries. The unique three prongs design with round tips to control the tissue separation and make a uniform opening. The blunt tips also avoid tracheal tears and injury. Moreover, the jaws have narrow angled necks so you can view the working tips clearly

LABORDE TRACHEAL DILATOR FEATURES

- Finger Ring Handles for Controlled Movements.
- Three prongs to ensure uniform tissue separation.
- Angled Jaws for Easy Surgical Access.
- Blunt Prongs for Safe Manipulations





WI-47040507

Lahey Cholecystectomy Forceps Surgical Instrument

\$8.00 per piece

Trade price from \$4.00
Minimum order 14 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Lahey Cholecystectomy Forceps Surgical Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each instrument is

checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$8.00, falling to \$4.00 at 14 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47060312

Lahey Cholecystectomy Forceps Surgical Tool

\$7.00 per piece

Trade price from \$3.80
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Lahey Cholecystectomy Forceps Surgical Tool for general and specialty operating theatre use — made from bar stock to finished instrument in house.

Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$7.00, falling to \$3.80 at 50 pieces — about 46% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000867****Lahey Gall Duct Forceps 7.5" Curved Longitudinal Serrated Jaws****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lahey Gall Duct Forceps 7.5" (19.1cm) Curved Longitudinal Serrated Jaws with Ratched Lock Ring Handle German Grade Stainless Steel Surgical Medical Gall Bladder Instrument.

Lahey Gall Duct Forceps 7.5" are used during a wide range of gastroenteritis surgical procedures. Hepatobiliary surgeons use these forceps to grasp, mobilize and occlude the common bile duct, especially during cholecystectomy procedures. They're a top-tier necessity for anyone performing gallbladder or biliary tract surgeries. Lahey Gall Duct Clamp When it comes to clamping down on bile ducts, you want something that's reliable, durable, and efficient. That's where these forceps come in.

Lahey Clamp Forceps has a pair of partially longitudinally serrated jaws that ensure optimal vessel grasping and minimize local injury. The self-retaining ratchet lock and boxlock joint that attaches the shafts ensures smooth closing movements.

LAHEY GALL DUCT FORCEPS FEATURES

- Crafted With 7.5" (19.1cm) Overall Length.
- Longitudinal Serrations for Enhanced Occlusion.
- Curved Tips Jaw.
- Long, Slim Profile for Improved Reach and Maneuverability.
- Ergonomic Finger Rings Ensuring Maximum Control.

WI-70000869**Lahey Gall Duct Forceps 9" Curved Longitudinal Serrated Jaws****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lahey Gall Duct Forceps 9" (22.9cm) Curved Longitudinal Serrated Jaws with Ratched Lock Ring Handle German Grade Stainless Steel Surgical Medical Gall Bladder Instrument.

Lahey Gall Duct Forceps 9" are used during a wide range of gastroenteritis surgical procedures. Hepatobiliary surgeons use these forceps to grasp, mobilize and occlude the common bile duct, especially during cholecystectomy procedures. They're a top-tier necessity for anyone performing gallbladder or biliary tract surgeries. Lahey Gall Duct Clamp When it comes to clamping down on bile ducts, you want something that's reliable, durable, and efficient. That's where these forceps come in.

Lahey Clamp Forceps has a pair of partially longitudinally serrated jaws that ensure optimal vessel grasping and minimize local injury. The self-retaining ratchet lock and boxlock joint that attaches the shafts ensures smooth closing movements.

LAHEY GALL DUCT FORCEPS FEATURES

- Crafted With 9" (22.9cm) Overall Length.
- Longitudinal Serrations for Enhanced Occlusion.
- Curved Tips Jaw.
- Long, Slim Profile for Improved Reach and Maneuverability.
- Ergonomic Finger Rings Ensuring Maximum Control.



WI-47054374

Lahey Surgical Forceps Cholecystectomy Instrument

\$7.90 per piece

Trade price from \$3.90
Minimum order 8 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Lahey Surgical Forceps Cholecystectomy Instrument, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$7.90, falling to \$3.90 at 8 pieces — about 51% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36254879

Lane Tissue Forceps Atraumatic Surgical Tissue Holding Forceps

\$2.50 per piece

Trade price from \$1.50
Minimum order 8 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Lane Tissue Forceps Atraumatic Surgical Tissue Holding Forceps is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.50, falling to \$1.50 at 8 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-35090434

Lane Tissue Forceps Reusable Surgical Instrument

\$2.30 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Lane Tissue Forceps Reusable Surgical Instrument for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.30, falling to \$1.50 at 10 pieces — about 35% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-35076893

Lane Tissue Forceps Stainless Steel Surgical Tissue Manipulation Instrument

\$2.30 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Lane Tissue Forceps Stainless Steel Surgical Tissue Manipulation Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.30, falling to \$1.50 at 10 pieces — about 35% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



WI-70003596

Light Silver Eyebrow Tweezers Eye Lash Extension Tweezer

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Light Silver Eyebrow Tweezers Eye Lash Extension Tweezer is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003595

Light Silver Eyebrow Tweezers Eyelash Extension Tweezer

\$17.43 per piece

Trade price from **\$16.04**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Light Silver Eyebrow Tweezers Eyelash Extension Tweezer is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$17.43, falling to \$16.04 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003756

London College Tweezer 6" (15.2cm) Dental Surgical

\$8.99 per piece

Trade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

London College Cotton Dressing Plier Tweezer 6" (15.2cm) Angled Serrated Tips Surgical Dental Endodontic Instruments

The college tweezers are primarily used by dentists moving cotton pellets, wedges, and articulation paper to and from the surgical tray to the mouth. The knurled handles provide a firm, consistent grip even with wet hands.


The sharp tips are serrated to ensure a firm grip on most small items, including jewelry, small parts, dissection, splinters, and more!

LONDON COLLEGE TWEEZER FEATURES:

- Cotton Plier is primarily used to grasp or transfer materials in and out of the oral cavity. Made of high

quality and safe stainless steel, which is very strong, rustproof and has excellent flexibility.

- This professional dental cotton plier is not only used for personal care, medical dressings, fake teeth cleaning, freshening breath, and whitening teeth it also can be used to keep your pet's teeth clean and refresh and avoid those annoying.
- These tooth tweezers are ergonomic design. It is Priority Recommended by lots of dentists. Millions of customers are loving in these useful dental tweezers and also share them with their family and friends.
- Suitable for household, craft, and DIY. Can be used for cleaning out the fish tank, electrical repair, crafting,



WI-34249008

Lower Molar Extraction Forceps Serrated Handles Ergonomic

\$8.00 per piece

Trade price from **\$5.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Lower Molar Extraction Forceps Serrated Handles Ergonomic is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: serrated jaws for grip that does not slip under tension, and an ergonomic handle for long lists without fatigue.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32732715

Lower Molar Extraction Forceps Serrated Handles Firm

\$8.00 per piece

Trade price from **\$5.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Lower Molar Extraction Forceps Serrated Handles Firm for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: serrated jaws for grip that does not slip under tension.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$5.00 at 25 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001812

Lucae Ear Forceps 5.5" Bayonet Shape with Serrated Tips

\$8.99 per piece

Trade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lucae Ear Forceps 5.5" (14cm) Bayonet Shape with Serrated Tips Stainless Steel ENT Otolaryngology Surgical Instruments

Lucae-Gruenwald Bayonet Ear Dressing Forceps assist professionals to treat ear diseases. In fine ear surgeries, doctors use them to place medication packs and remove foreign bodies in the ear canal. This is very versatile forcep that can be used in wound dressing and general tissue grasping as well.

The serrations on the handles give a strong grip for efficient and swift packing. In addition, the serrations prevent slips during delicate movements. Moreover, the bayonet design gives a more convenient vision of the jaws for accuracy. These ear forceps made from the German grade

stainless steel that you can depend on for years of efficient performance.

LUCAE EAR FORCEPS FEATURES

- Bayonet Shaft Design For Optimal Visualization.
- Long Serrated Jaws For Better Access & Grip
- Flat Handles For Enhanced Dexterity.
- Lightweight
- Made with High Quality Stainless Steel


PACKAGE INCLUDE:

- 1 Pc Lucae Ear Forceps 5.5" (14cm) – Serrated Tips

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Ear Forceps – Otolaryngology Surgical Instrument



**WI-70002132****Magill Catheter Forceps 10" (25.4cm) – Adult****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Magill Catheter Forceps 10" (25.4cm)
Adult Size – Premium Quality Stainless
Steel ENT Anesthesia Instruments
EMT/EMS Supplies Intubation Forcep.

Magill catheter forceps are specialized, angled surgical instruments used in emergency medicine and airway management. They assist in nasotracheal intubation, foreign body retrieval, catheter and tube placement, and packing placement. The design allows the user's hand to remain out of the line of sight. They are available in various sizes for different patient needs and are made from durable, high-quality stainless steel, making them autoclavable and easy to clean.

**MAGILL INTUBATION FORCEPS
FEATURES**

- 10 Inch (25.4cm) Overall Length.
- Oblique Angles Designed to Maintain Visibility of the Larynx During Use.
- Atraumatic Serration Fenestrated Jaws.
- Rounded Ends Provide Secure Gripping Without Damaging Tissue or Equipment.
- Ring Handles Enable Effective Gripping and Manipulation by the User.
- Available in Sizes for Infants, Children and Adults for Versatile Application.
- Utilized for Guiding Endotracheal Tubes and Removing Airway Obstructions.
- Suitable for Sterilization Processes.
- Constructed From Corrosion-resistant Premium Grade Stainless

**WI-70002133****Magill Catheter Forceps 8" (20.3cm) – Child****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Magill Catheter Forceps 8" (20.3cm)
Child Size – Premium Quality Stainless
Steel ENT Anesthesia Instruments
EMT/EMS Supplies Intubation Forcep.

Magill catheter forceps are specialized, angled surgical instruments used in emergency medicine and airway management. They assist in nasotracheal intubation, foreign body retrieval, catheter and tube placement, and packing placement. The design allows the user's hand to remain out of the line of sight. They are available in various sizes for different patient needs and are made from durable, high-quality stainless steel, making them autoclavable and easy to clean.

**MAGILL INTUBATION FORCEPS
FEATURES**

- 8 Inch (20.3cm) Overall Length.
- Oblique Angles Designed to Maintain Visibility of the Larynx During Use.
- Atraumatic Serration Fenestrated Jaws.
- Rounded Ends Provide Secure Gripping Without Damaging Tissue or Equipment.
- Ring Handles Enable Effective Gripping and Manipulation by the User.
- Available in Sizes for Infants, Children and Adults for Versatile Application.
- Utilized for Guiding Endotracheal Tubes and Removing Airway Obstructions.
- Suitable for Sterilization Processes.
- Constructed From Corrosion-resistant Premium Grade Stainless



**WI-70004116****Magill Forceps 6" (15.2cm) – Infant****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Magill Catheter Forceps 6" (15.2cm)
Infant Size – Premium Quality Stainless
Steel ENT Anesthesia Instruments
EMT/EMS Supplies Intubation Forcep.

Magill catheter forceps are specialized, angled surgical instruments used in emergency medicine and airway management. They assist in nasotracheal intubation, foreign body retrieval, catheter and tube placement, and packing placement. The design allows the user's hand to remain out of the line of sight. They are available in various sizes for different patient needs and are made from durable, high-quality stainless steel, making them autoclavable and easy to clean.

**MAGILL INTUBATION FORCEPS
FEATURES**

- 6 Inch (15.2cm) Overall Length.
- Oblique Angles Designed to Maintain Visibility of the Larynx During Use.
- Atraumatic Serration Fenestrated Jaws.
- Rounded Ends Provide Secure Gripping Without Damaging Tissue or Equipment.
- Ring Handles Enable Effective Gripping and Manipulation by the User.
- Available in Sizes for Infants, Children and Adults for Versatile Application.
- Utilized for Guiding Endotracheal Tubes and Removing Airway Obstructions.
- Suitable for Sterilization Processes.
- Constructed From Corrosion-resistant Premium Grade Stainless

**WI-70000532****Magnifying Splinter Forcep Tweezer 3.5" Fine Pointed Tips****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Splinter Tweezer Forceps 3.5" with Magnifying Glass – Premium Quality Stainless Steel Precision-Tipped Multipurpose Magnifying Tweezer for Removal of Splinters and Foreign Object, Crafts, Laboratory, EMS, EMT First Aid Kit Instrument.

Magnifying Splinter Forcep is a versatile, precision-tipped tweezer with a magnifying glass, ideal for home or workshop use. It's perfect for removing splinters, foreign objects, and fine detail work, with a compact hinged design for easy storage. Splinter Magnifying Tweezer is high-quality stainless steel tool with a folding pocket design, precision fine point tips, serrated non-slip handle, and a magnifying glass tip making this


tweezer a must have tool for any home first aid kit.

This handheld magnifying glass tweezer is an essential beauty tool for shaping eyebrows, removing ingrown hairs along with removing wood splinters due to its secure grip and hand-sharpened tip. Also suitable for industrial electronics, field investigation, reading, laboratory use, watch repair, jewelry identification, coins, stamps, and antiques.

**MAGNIFYING SPLINTER FORCEP
TWEEZERS FEATURES**

- Crafted With 3.5" (8.9cm) Overall Length.
- Fine Pointed Tips to Hold and Remove Extremely Small Particles.



**WI-70000723****Maier Uterine Polypus Dressing Forceps 10" Serrated 5X30mm Straight Jaws****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



MAIER Uterine Polypus Dressing Forceps 10" (25.4cm) Serrated Jaws 5X30mm Straight with Ratchet Lock System Medical Grade Stainless Steel Gynecological Surgical Instruments.

Maier Polypus Forceps is a specialized instrument commonly used during gynecological surgeries to remove tissues from the uterus, handle internal organs as well as manipulate gauze pieces during surgery. Maier forceps is also frequently used during general surgery as well. Maier uterine dressing forceps feature a ratcheted mechanism towards the ring handles and serrated tips with internal hollow groove. This forcep is ideal for tasks such as cervical disinfection and IUD removal, as well as general dressing and gripping in clinical settings.

Maier forceps comes in straight and curved profiles with and without ratchet according to the surgical approach. MAIER uterine polypus dressing forceps made with premium grade stainless steel, are fully autoclavable ensures lasting performance with 100% guaranteed satisfaction.

MAIER UTERINE POLYPUS FORCEPS FEATURES

- 10" (25.1cm) Overall Length.
- 5×30mm Serrated Jaws.
- Straight Oval Shapped with Inner Grooved Tips.
- Elongated Straight Shafts With Boxjoint Grating Deep Surgical Access.
- Ratchet Lock System.

**WI-70001738****Mashing Tweezer 8" Glass Work Bead Masher****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

8" Stainless Steel Mashing Tweezers 3/4" Square Pad Glass Work Bead Masher Tweezer for Glassblowing Boro Lampworking.

Mashing Tweezers with 3/4" square pads are specialized tools for glass lampworking, flameworking, beadmaking, lamp-working in laboratory, electronics and many other applications designed to flatten, shape and compress hot, soft glass. Their parallel flat square pads are perfect for creating flower petals and other floral details by applying even pressure. They help square off beads or incorporate flat facets in small sculptures without leaving deep tool marks. Crafted from premium-grade stainless steel, these tweezers ensure gentle handling and

controlled pressure on delicate glass pieces.

MASHING TWEEZER FEATURES

- 3/4" X 3/4" Parallel Flat Square Pads.
- Overall length: 8" (20.3cm).
- Flat Evenly Close Tips Apply Broad, Controlled Pressure During Operation.
- Great for Mashing Flower Pedals Leaf Shapes.
- Made with High Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Masher Tweezers 8 " with 3/4" Flat Square Pads

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Glasswork Tools

Check our entire product range by visiting this



WI-36155701

Maxillofacial Surgical Instrument Set Stainless Steel Oral Surgery

\$190.00 per piece

Trade price from **\$170.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Maxillofacial Surgical Instrument Set Stainless Steel Oral Surgery is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$190.00, falling to \$170.00 at 1 set — about 11% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47075025

McClintock Ovum Forceps Gynecological Surgical Instrument

\$10.00 per piece

Trade price from **\$6.00**
Minimum order 12 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

McClintock Ovum Forceps Gynecological Surgical Instrument is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by

hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$10.00, falling to \$6.00 at 12 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47060663

McClintock Ovum Forceps Gynecology Surgical Instrument

\$9.50 per piece

Trade price from **\$5.50**
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of McClintock Ovum Forceps Gynecology Surgical Instrument, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$9.50, falling to \$5.50 at 50 pieces — about 42% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47072473

McClintock Steel Ovum Forceps Reusable Gynecological Instrument

\$9.90 per piece

Trade price from **\$5.90**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture McClintock Steel Ovum Forceps Reusable Gynecological Instrument for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$9.90, falling to \$5.90 at 25 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003735

McIndoe Tissue Forceps 6" Straight 1.5mm Wide 1x2 Teeth

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


McIndoe Tissue Forceps 6" (15.2cm) Straight Shaft 1.5mm Wide 1x2 Teeth With Flat Handle Premium Quality Stainless Steel Delicate Tissue Dissection/ Manipulation Ophthalmic Surgery or Plastic Surgery Instrument.

Delicate mcindoe tissue forceps are precision instruments designed for handling fine tissues with accuracy. They feature a narrow 1.5mm tip with 1x2 teeth offers a secure grip, preventing slippage on slippery tissues for a secure grip that minimizes trauma. Mcindoe forceps helps surgeons in holding tissue for stitching, placing or removing dressings, and clearing debris to ensure precision in small areas. Commonly used in ophthalmic, plastic and reconstructive surgeries, wound care and microdissection where

precision is essential. Their lightweight and balanced construction reduces hand fatigue during prolonged use, making them easy to sterilize and reusable.

MCINDOE DELICATE TISSUE FORCEPS FEATURES

- Overall Length 6 Inches (15.2cm).
- Straight Shaft.
- 1.5 mm Wide 1x2 Teeth Delicate Tips.
- Spring Action Thumb Grip Flat Handle .
- Perfect for Delicate Tissue Dissection/manipulation.
- Suitable for Ophthalmic Surgery or Plastic Surgery and Micro Dissection.
- Lightweight.
- Autoclavable/Reusable.



WI-70003374

Mclean-Tucker Obstetrical Forceps

\$45.99 per piece

Trade price from \$42.31
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Mclean-Tucker Obstetrical Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$45.99, falling to \$42.31 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-43629101

McPherson Suture Tying Forceps Precision Eye Surgery Tool

\$9.00 per piece

Trade price from \$5.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

McPherson Suture Tying Forceps Precision Eye Surgery Tool — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$9.00, falling to \$5.00 at 10 pieces — about 44% off the single-unit price. Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-28242334

Medic Instrument Foerster Stainless Steel

\$8.00 per piece

Trade price from \$4.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Medic Instrument Foerster Stainless Steel is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



WI-26119318

Medical Stainless Steel Sponge Forceps Curved

\$3.00 per piece

Trade price from **\$1.50**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Medical Stainless Steel Sponge Forceps Curved, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a curved working end for reach around structures.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.00, falling to \$1.50 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27738897

Medical Surgical Stainless Steel Hospital Forceps

\$5.50 per piece

Trade price from **\$2.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Medical Surgical Stainless Steel Hospital Forceps — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$5.50, falling to \$2.00 at 1 piece — about 64% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002134

Meriam Tweezer 16cm – Serrated Angled Tip

\$8.99 per piece

Trade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Meriam Foil Dental Cotton Tweezer
Serrated Angled Beak Dressing Pliers
Forceps Stainless Steel Surgical Dental
Instruments

MERIAM FOIL TWEEZER FEATURES:

PREMIUM QUALITY MERIAM FOIL TWEEZERS: Meriam Foil Tweezers are a high quality and valuable addition to dental clinics and practices. Meriam Tweezers are used for grasping or transferring materials in and out of the oral cavity during clinical procedures. They have a strong hold with serrated beaks to firmly grasp gauze, cotton and other objects with outstanding performance. The slender body and angled tips are ideal for hard to reach and tight areas.

DURABLE O.R. GRADE STAINLESS

STEEL CONSTRUCTION: Beautifully forged from corrosion-resistant Premium Quality Stainless Steel for long product life and supreme durability. High value for the price and a great addition to your inventory. Strong high quality steel that does not tarnish, degrade or bend over time. Perfectly aligned sturdy tips for even and dependable results.

FULLY AUTOCLAVABLE + HIGH MIRROR FINISH: Highly polished finish for wonderful aesthetic + ultrasonic cleaned. Easy to clean and fully autoclavable for proper sterilization and sanitation procedures.

PRECISE AND COMFORTABLE TO USE: The comfortable design features a grip panel for ultimate control and precision. Lightweight design offers outstanding



WI-46972100

Michel Clip Forceps Surgical Instrument for Precise Skin Clip

\$9.90 per piece

Trade price from **\$4.90**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Michel Clip Forceps Surgical Instrument for Precise Skin Clip, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$9.90, falling to \$4.90 at 25 pieces — about 51% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46951208

Michel Clip Removal Forceps

\$9.00 per piece

Trade price from **\$5.00**
Minimum order 12 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Michel Clip Removal Forceps, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$9.00, falling to \$5.00 at 12 pieces — about 44% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46970135

Michel Clip Surgical Forceps Precision Instrument for Skin Closure

\$9.50 per piece

Trade price from **\$5.00**
Minimum order 35 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

We manufacture Michel Clip Surgical Forceps Precision Instrument for Skin Closure for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$9.50, falling to \$5.00 at 35 pieces — about 47% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46933488

Michel Skin Clip Forceps Stainless Steel Surgical Instrument

\$10.00 per piece

Trade price from **\$5.00**
Minimum order 12 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Michel Skin Clip Forceps Stainless Steel Surgical Instrument, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$10.00, falling to \$5.00 at 12 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001651

Micro Alligator Forceps 3.3" Serrated Jaws

\$8.99 per piece

Trade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Micro Alligator Forceps 3.3" (8.3cm) Ophthalmic Lab Instruments, ENT Working Serrated. These Alligator Forceps is used for cerumen Removal. Pointed, Very Small Serrated Jaws.

Crocodile Micro Alligator Forceps is a multipurpose ENT device that surgeons commonly use to grasp, hold and manipulate the most delicate tissues within the middle ear. The forceps are also useful for grasping tubes and positioning them precisely inside the ear cavity.

MICRO ALLIGATOR FORCEPS FEATURES:


MICRO ALLIGATOR FORCEPS – The Ophthalmic Lab Micro Alligator Forceps are used for cerumen removal. They are great for grasping small objects or

removing items from small cavities in the body. They have a Serrated jaw. **PREMIUM QUALITY STAINLESS STEEL**– Forged from O.R. Grade Stainless Steel with Superior Craftsmanship. Tarnish free, corrosion resistant polished finish.

PRECISE & EASY TO USE – The Micro Alligator Forceps are pointed and very small which allows them a high degree of precision and flexibility while conducting the clinical procedure.

EASY TO CLEAN – 100% brand new, reusable, and fully autoclavable Micro Alligator Ophthalmic Lab Surgical Forceps.

QUALITY AND SAFETY STANDARDS: All our instruments are compliant with all National and International quality and safety standards making this item safe and world class.



**WI-70001650****Micro Alligator Forceps 3.5" Serrated Jaws****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Micro Alligator Forceps 3.5" (8.9cm) Ophthalmic Lab Instruments, ENT Working Serrated. These Alligator Forceps is used for cerumen Removal. Pointed, Very Small Serrated Jaws.

Crocodile Micro Alligator Forceps is a multipurpose ENT device that surgeons commonly use to grasp, hold and manipulate the most delicate tissues within the middle ear. The forceps are also useful for grasping tubes and positioning them precisely inside the ear cavity.

MICRO ALLIGATOR FORCEPS FEATURES:

MICRO ALLIGATOR FORCEPS – The Ophthalmic Lab Micro Alligator Forceps are used for cerumen removal. They are great for grasping small objects or

removing items from small cavities in the body. They have a Serrated jaw.

PREMIUM QUALITY STAINLESS STEEL– Forged from O.R. Grade Stainless Steel with Superior Craftsmanship. Tarnish free, corrosion resistant polished finish.

PRECISE & EASY TO USE – The Micro Alligator Forceps are pointed and very small which allows them a high degree of precision and flexibility while conducting the clinical procedure.

EASY TO CLEAN – 100% brand new, reusable, and fully autoclavable Micro Alligator Ophthalmic Lab Surgical Forceps.

QUALITY AND SAFETY STANDARDS: All our instruments are compliant with all National and International quality and safety standards making this item safe and world class.

**WI-70002749****Micro Alligator Forceps 3"****\$8.99** per pieceTrade price from **\$8.27**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Micro Alligator Forceps 3" (7.6cm) Ophthalmic Lab Instruments, ENT Working Serrated. These Alligator Forceps is used for cerumen Removal. Pointed, Very Small Serrated Jaw (0.8mm X 4.5mm)

Crocodile Micro Alligator Forceps is a multipurpose ENT device that surgeons commonly use to grasp, hold and manipulate the most delicate tissues within the middle ear. The forceps are also useful for grasping tubes and positioning them precisely inside the ear cavity.

MICRO ALLIGATOR FORCEPS FEATURES:

MICRO ALLIGATOR FORCEPS – The Ophthalmic Lab Micro Alligator Forceps are used for cerumen removal. They


are great for grasping small objects or removing items from small cavities in the body. They have a Serrated jaw.

PREMIUM QUALITY STAINLESS STEEL– Forged from O.R. Grade Stainless Steel with Superior Craftsmanship. Tarnish free, corrosion resistant polished finish.

PRECISE & EASY TO USE – The Micro Alligator Forceps are pointed and very small which allows them a high degree of precision and flexibility while conducting the clinical procedure.

EASY TO CLEAN – 100% brand new, reusable, and fully autoclavable Micro Alligator Ophthalmic Lab Surgical Forceps.

QUALITY AND SAFETY STANDARDS: All our instruments are compliant with all National and International quality and safety standards making this item safe



**WI-70004283****Micro Alligator Forceps Black Coated 3.5"****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hartmann Alligator Ear Forceps 3.5" (8.89cm) Curved Serrated Jaw Design, Non-Slip Grip Handle, Fine-Tip, Reusable and Easy to Clean (Black Coated)

MICRO ALLIGATOR FORCEPS FEATURES

PREMIUM QUALITY: ENT Black Coated Ear Alligator Forceps made of high-grade stainless steel, ensuring top quality and durability. Ideal for professionals and medical students.

CURVED JAW DESIGN: Unique curved jaw design allows easy access to hard-to-reach areas, perfect for precise ENT procedures involving the ear, nose, and throat. Suitable for various ENT procedures, a must-have tool for medical professionals, clinics, and hospitals.

FINE-TIP: Precision forceps with fine tips ideal for delicate procedures like myringotomy, foreign object removal, and ear tube placement.

NON-SLIP GRIP HANDLE: Comfortable, secure grip with a non-slip handle, enabling better control during procedures and improved user experience.

REUSABLE & EASY TO CLEAN: 100% brand new, fully autoclavable forceps for easy cleaning and reuse, ensuring hygiene and safety.

PACKAGE INCLUDE:

- 1 Pc Black Coated Hartman Alligator Micro Tip Forceps 3.5" (8.89cm) Serrated Jaws

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

**WI-70000278****Micro Extracting Forceps #133L Left Angled Tapered Tips****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Micro Extracting Forceps #133L Left Angled Tapered Tips for Incisor & Root Fragments with Non-Slip Finger Grip Handle – German Grade Stainless Steel English Pattern Teeth Extraction for Veterinary Dentists and General Practitioners.

Micro Extracting Forceps #133L are precision dental tools designed for delicate extractions, offering excellent control for removing small roots and fragments. Their finely tapered tips facilitate access in tight areas, making them ideal for endodontic, periodontal, and surgical dentistry. These fully autoclavable forceps are durable, crafted from high-quality German-grade stainless steel, ensuring strength and corrosion resistance, while their

ergonomic design enhances comfort during extended use.

English pattern forceps are designed so that the fulcrum is closer to the working end, providing greater leverage, easing extraction. These extracting forceps provide reliable performance for precise micro-extractions in veterinary dentistry and general atraumatic surgical extraction practice.

MICRO EXTRACTING FORCEPS #133L FEATURES

- Size: 7.25" (18.4cm)
- For Upper Incisor & Root Fragments.
- Left Angled Finally Tapered Tips.
- Non-Slip Finger Grip Handle.
- English Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:

**WI-70000277****Micro Extracting Forceps #133R Right Angled Tapered Tips****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Micro Extracting Forceps #133R Right Angled Tapered Tips for Incisor & Root Fragments with Non-Slip Finger Grip Handle – German Grade Stainless Steel English Pattern Teeth Extraction for Veterinary Dentists and General Practitioners.

Micro Extracting Forceps #133R are precision dental tools designed for delicate extractions, offering excellent control for removing small roots and fragments. Their finely tapered tips facilitate access in tight areas, making them ideal for endodontic, periodontal, and surgical dentistry. These fully autoclavable forceps are durable, crafted from high-quality German-grade stainless steel, ensuring strength and corrosion resistance, while their

ergonomic design enhances comfort during extended use.

English pattern forceps are designed so that the fulcrum is closer to the working end, providing greater leverage, easing extraction. These forceps provide reliable performance for precise micro-extractions in veterinary dentistry and general atraumatic surgical extraction practice.

MICRO EXTRACTING FORCEPS #133R FEATURES

- Size: 7.25" (18.4cm)
- For Upper Incisor & Root Fragments.
- Right Angled Finally Tapered Tips.
- Non-Slip Finger Grip Handle.
- English Pattern.
- Made With Stainless Steel

PACKAGE INCLUDE:**WI-70003734****Microsurgery Forceps 7.25" TC Micro Tying Forcep****\$18.74** per pieceTrade price from **\$17.24**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Tungsten Carbide Inserts Micro Suture Tying Forceps 7.25" with Tying Platform Straight Round Grip Gold Plated Handle Microsurgery Surgical Instruments

TUNGSTEN CARBIDE MICRO SUTURE TYING FORCEPS FEATURES

TC MICRO SUTURE TYING FORCEPS: T/C Micro Suture Tying Forceps 7.25" with Platform – Premium O.R. Grade Quality

MADE OF LONG LASTING MATERIALS: Forged from corrosion-resistant surgical stainless steel and tough tungsten carbide, these popular t/c micro suture tying forceps are built to last. The added tungsten carbide tip ensures that this instruments will not break or degrade on you. T/C is 10x

stronger than steel and is the toughest alloy on earth.

PLATFORM FOR PRECISE GRIP: These micro suture forceps have a platform near the tip of the forceps that aids in suturing. The elevated surface holds the suture securely to help make a knot. Unlike regular forceps that meet at their tips when they are closed, forceps with suture tying platforms close flat along along surface from the tip. This feature allows for a reliable, easy and gentle grip, which may be useful in applications other than suture work, too.

EASY TO USE + CLEAN: These forceps are user friendly and easy to clean as well as fully autoclavable. They can be counted on for many years in many successful applications.

**WI-70003762****Miller Articulating Paper Forceps 6" Straight Surgical Dental Tweezer****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Miller Articulating Paper Forceps 6" (15.2cm) Straight Reverse-Action Serrated Tips Stainless Steel Dental Restorative Instruments.


Miller articulating paper forceps is a reverse-action tweezer like dental instruments used to hold and place articulating paper to check occlusal surfaces of teeth before and after restorative procedures. Holding the paper with an instrument rather than fingers assures better and complete coverage of the restoration or teeth to be verified.

Articulating paper forceps is an essential tool for handling articulating paper. Facilitates the placement paper in the area of posterior teeth and helps with the detection of points of contact between opposing teeth, making it

possible to discover those that need to be adjusted and corrected. The articulating paper plier is suitable for dentists, art & craftwork, and other laboratory works. You can reuse these articulating paper forceps over time as they will not bend or tarnish on you. High quality stainless steel construction meet all international quality standard.

ARTICULATING PAPER FORCEPS FEATURES

- 6" (15.2cm) Overall Length
- Straight Tips
- Serrated Blades to Ensure Strong Grip
- Reverse-action Opening & Closing Mechanism for Firm Grip
- Long-lasting, Rust Free
- Reusable

**WI-47039525****Millin Vulsellum Forceps Gynecological Surgical Instrument for Secure Tissue****\$12.00** per pieceTrade price from **\$10.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Millin Vulsellum Forceps Gynecological Surgical Instrument for Secure Tissue is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$12.00, falling to \$10.00 at 10 pieces — about 17% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47056422

**Millin Vulsellum Forceps Reusable
Gynecological Surgical Tool for Controlled**

\$11.90 per piece

Trade price from **\$7.90**
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Millin Vulsellum Forceps Reusable Gynecological Surgical Tool for Controlled is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$11.90, falling to \$7.90 at 50 pieces — about 34% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47062362

**Millin Vulsellum Surgical Forceps with Clip for
Gynecology-Precise**

\$12.00 per piece

Trade price from **\$8.50**
Minimum order 8 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Millin Vulsellum Surgical Forceps with Clip for Gynecology-Precise for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$12.00, falling to \$8.50 at 8 pieces — about 29% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000160****Mixer Forceps 11" (27.9cm) Right Angled Longitudinal Serration Jaws****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Mixer Forceps 11" (27.9cm) Right Angled with Longitudinal Serration Premium Quality Stainless Steel Dissecting and Ligature Forcep for Cardiovascular, Thoracic, and Gynecological Surgeries.

Mixer Forceps with longitudinal serrations are surgical instruments used for deep tissue dissection, clamping blood vessels, and handling both heavy and delicate tissues. Their right-angle jaw design facilitates access to difficult areas, making them essential in thoracic, cardiovascular, and gynecological procedures for holding, isolating, or ligating structures. Longitudinally serrated mixte forcep is designed with a ratchet mechanism that allows for self-locking jaws, enabling surgeons to exert four distinct

pressure levels on clamped tissues. Crafted from premium grade stainless steel to guarantee superior quality and performance in surgical operations.

MIXTER THORACIC HEMOSTAT FORCEPS FEATURES

- 11 Inches (27.9 cm) in Length.
- Right Angle Design With Longitudinal Serration.
- Ratchet Mechanism for Self-locking Jaws
- Ergonomic Finger Rings Ensuring Maximum Control.
- Used for Dissecting and Ligature During Surgeries.
- Suitable for Cardiovascular, Thoracic, and Gynecological Procedures.
- Fully Autoclavable/Reusable.

**WI-70000147****Mixer Forceps 7.25" (18.4cm) Fully Curved Serrated Jaws****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Mixer Forceps 7.25" (18.4cm) Right Angled Vertical Serration Jaws with Finger Ring Handle – Premium Quality Stainless Steel Blood Vessels Clamping and Delicate Tissues Dissecting Forceps for Gallbladder, General Surgery, Orthopedic and Gynecological Surgeries.

Mixer hemostatic forceps are versatile surgical instruments designed for effective hemostasis used to clamp blood vessels, control bleeding, and dissect or manipulate delicate tissues. They feature right angled, fully curved, vertical serrated jaws with ring handles. The design enhances maneuverability and provides precise control in delicate procedures, while the serrated jaws ensure a secure grip and minimize slippage. The ratcheted lock ring

handle allows for easy handling and accurate pressure control. Made from high-quality stainless steel, mixer forceps are durable and reliable for repeated use.

MIXTER FORCEPS FEATURES

- 7.25 Inches (18.4 cm) in Length.
- Right Angle Design With Vertical Serration for a Secure Grip.
- Elongated Straight Shafts for Hard-to-reach Surgical Areas.
- Boxlock System for Smooth Joint Movements.
- Ratcheted Mechanism for Self-locking Jaws.
- Ergonomic Finger Rings Ensuring Maximum Control.
- Designed to Clamp Blood Vessels and Control Bleeding.



**WI-70000146****Mixer Forceps 9" (22.9cm) Fully Curved Serrated Jaws****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Mixer Forceps 9" (22.9cm) Right Angled Vertical Serration Jaws with Finger Ring Handle – Premium Quality Stainless Steel Blood Vessels Clamping and Delicate Tissues Dissecting Forceps for Gallbladder, General Surgery, Orthopedic and Gynecological Surgeries.

Mixer hemostatic forceps are versatile surgical instruments designed for effective hemostasis used to clamp blood vessels, control bleeding, and dissect or manipulate delicate tissues. They feature right angled, fully curved, vertical serrated jaws with ring handles. The design enhances maneuverability and provides precise control in delicate procedures, while the serrated jaws ensure a secure grip and minimize slippage. The ratcheted lock ring

handle allows for easy handling and accurate pressure control. Made from high-quality stainless steel, mixer forceps are durable and reliable for repeated use.

MIXTER FORCEPS FEATURES

- 9 Inches (22.9 cm) in Length.
- Right Angle Design With Vertical Serration for a Secure Grip.
- Elongated Straight Shafts for Hard-to-reach Surgical Areas.
- Boxlock System for Smooth Joint Movements.
- Ratcheted Mechanism for Self-locking Jaws.
- Ergonomic Finger Rings Ensuring Maximum Control.
- Designed to Clamp Blood Vessels and Control Bleeding.

**MIXTER HEMOSTATIC FORCEPS
HALF SERRATED**

90°

**WI-70000161****Mixer Forceps 9" (22.9cm) Right Angled Longitudinal Serration Jaws****\$13.99** per pieceTrade price from **\$12.87**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Mixer Forceps 9" (22.9cm) Right Angled with Longitudinal Serration Premium Quality Stainless Steel Dissecting and Ligature Forcep for Cardiovascular, Thoracic, and Gynecological Surgeries.

Mixer Forceps with longitudinal serrations are surgical instruments used for deep tissue dissection, clamping blood vessels, and handling both heavy and delicate tissues. Their right-angle jaw design facilitates access to difficult areas, making them essential in thoracic, cardiovascular, and gynecological procedures for holding, isolating, or ligating structures. Longitudinally serrated mixte forcep is designed with a ratchet mechanism that allows for self-locking jaws, enabling surgeons to exert four distinct

pressure levels on clamped tissues. Crafted from premium grade stainless steel to guarantee superior quality and performance in surgical operations.

**MIXTER THORACIC HEMOSTAT
FORCEPS FEATURES**

- 9 Inches (22.9 cm) in Length.
- Right Angle Design With Longitudinal Serration.
- Ratchet Mechanism for Self-locking Jaws
- Ergonomic Finger Rings Ensuring Maximum Control.
- Used for Dissecting and Ligature During Surgeries.
- Suitable for Cardiovascular, Thoracic, and Gynecological Procedures.
- Fully Autoclavable/Reusable.

**MIXTER HEMOSTATIC FORCEPS**

90°



WI-47249272

Mixer Petit Surgical Forceps Curved Instrument for Controlled Vessel

\$4.90 per piece

Trade price from \$3.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Mixer Petit Surgical Forceps Curved Instrument for Controlled Vessel is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a curved working end for reach around structures.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$4.90, falling to \$3.00 at 25 pieces — about 39% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001725

Multi Rainbow Color Kelly Hemostat Forceps 5.5" Curved

\$10.99 per piece

Trade price from \$10.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Multi Titanium Rainbow Color Kelly Hemostat 5.5" (14cm) Curved Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Serrated Jaws, Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobby


KELLY FORCEPS FEATURES

RAINBOW COLOR COATED KELLY HEMOSTAT FORCEPS: Waqar International's Rainbow Color Kelly Forceps 14cm Curved hemostat clamps are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

THREE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once locked in place, even the smallest objects can be held securely in the serrated jaws with great force.

A MUST-HAVE TOOL FOR REMOVING HOOKS: These Fishing clamps are one



**WI-70001726****Multi Rainbow Color Kelly Hemostat Forceps 5.5\"/>****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Multi Titanium Rainbow Color Kelly Hemostat 5.5\" (14cm) Straight Multipurpose Hemostat Clamping Plier Locking Tweezer Clamps Serrated Jaws, Stainless Steel Premium Quality Hemostats for Nurses, Fishing Forceps, Crafts and Hobby

KELLY FORCEPS FEATURES**RAINBOW COLOR COATED KELLY HEMOSTAT FORCEPS:**

Waqar International's Rainbow Color Kelly Forceps 14cm Straight hemostat clamps are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL

CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

THREE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once locked in place, even the smallest objects can be held securely in the serrated jaws with great force.

A MUST-HAVE TOOL FOR REMOVING HOOKS: These Fishing clamps are one

**WI-70001730****Multi Rainbow Mosquito Forceps 5\"/>****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Multi Titanium Rainbow Color Halsted Mosquito Micro Forcep Curved 5\" (12.7cm) Serrated Jaws Hemostat Clamps Locking Forceps High-Quality Stainless Steel Surgical Dental Instruments Tool

Halsted Mosquito Hemostatic Forceps are versatile, ratcheted, finger ring forceps used to clamp off smaller vessels that control the flow of blood and other clamping tasks in general surgical procedures. The ratcheted finger ring handle of Halstead Forceps provides a secure, locking grip to halt blood flow. A toothed variation is available for grasping delicate tissue in skin grafting, biopsies, or ophthalmologic procedures. They are not appropriate for use with large vessels or arteries.

Mosquito forceps are known for their fine tips, and short, fully serrated jaws. They are available curved or straight. The standard length for this instrument is 5 inches, although does offer some with a length of 5-1/2 inches. Halsted Mosquito Forceps are very similar to Hartmann Mosquito Forceps, however, the Hartmann style is finer and more delicate.

HALSTED MOSQUITO FORCEPS FEATURES**RAINBOW HEMOSTAT MOSQUITO**

FORCEPS: Our rainbow color coated 5\" Curved Halsted Mosquito are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips,



**WI-70001731****Multi Rainbow Mosquito Forceps 5" Straight****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Multi Titanium Rainbow Color Halsted Mosquito Micro Forcep Straight 5" (12.7cm) Serrated Jaws Hemostat Clamps Locking Forceps High-Quality Stainless Steel Surgical Dental Instruments Tool

Halsted Mosquito Hemostatic Forceps are versatile, ratcheted, finger ring forceps used to clamp off smaller vessels that control the flow of blood and other clamping tasks in general surgical procedures. The ratcheted finger ring handle of Halstead Forceps provides a secure, locking grip to halt blood flow. A toothed variation is available for grasping delicate tissue in skin grafting, biopsies, or ophthalmologic procedures. They are not appropriate for use with large vessels or arteries.

Mosquito forceps are known for their fine tips, and short, fully serrated jaws. They are available curved or straight. The standard length for this instrument is 5 inches, although does offer some with a length of 5-1/2 inches. Halsted Mosquito Forceps are very similar to Hartmann Mosquito Forceps, however, the Hartmann style is finer and more delicate.

HALSTED MOSQUITO FORCEPS FEATURES

RAINBOW HEMOSTAT MOSQUITO FORCEPS: Our rainbow color coated 5" Straight Halsted Mosquito are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips,

**WI-70000977****Museux Tenaculum Forcep 6mm Wide 2x2 Teeth 9.5" Curved****\$14.99** per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Museux Tenaculum Forceps 9.5" (24.1cm) with 6mm Wide 2x2 Overlapping Teeth Sideway Curved Jaws Finger Ring Ratcheted Handle Premium Quality Stainless Steel Obstetrics & Gynecology Instruments.

Museux uterine vulsellum forceps is specialized OB/GYN Instruments with finger rings & ratcheting mechanism use for grasping the cervix, fallopian tubes, ovaries and other internal female organs. Museux Tenaculum Forceps can also use for manipulating and mobilizing intrapelvic soft tissues and expose the surgical field predominantly in surgery to seize blood vessels or to steady the cervix or uterus during hysterectomies, oophorectomies, amniotomies and other procedures.

Museux Tenaculum Forceps with 2x2 overlapping Teeth is a dependable solution for grasping and manipulating the uterus, available in 6, 8 and 10mm izes for accommodating a variety of tissues. With boxlock joint that permits smooth jaw closure movements, this tenaculum forceps are made from premium quality stainless steel and are manufactured to the highest standards of quality and precision.

MUSEUX TENACULUM FORCEPS FEATURES

- 9.5 Inches (24.1cm) Overall Length.
- 2 X 2 Overlapping Teeth for Firm Grip.
- 6mm Wide Teeth.
- Available in Straight and Curved Profile.





WI-70000976

Museux Tenaculum Forceps 6mm Wide 2x2 Teeth 9.5" Straight

\$14.99 per piece

Trade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Museux Tenaculum Forceps 9.5" (24.1cm) Straight Jaws 6mm Wide 2x2 Overlapping Teeth Finger Ring Ratchted Handle Premium Quality Stainless Steel Obstetrics & Gynecology Instruments.


Museux uterine vulsellum forceps is specialized OB/GYN Instruments with finger rings & ratcheting mechanism use for grasping the cervix, fallopian tubes, ovaries and other internal female organs. Museux Forceps can also use for manipulating and mobilizing intrapelvic soft tissues and expose the surgical field predominantly in surgery to seize blood vessels or to steady the cervix or uterus during hysterectomies, oophorectomies, amniotomies and other procedures.

With 2x2 overlapping Teeth, museux tenaculum is a dependable solution for

grasping and manipulating the uterus, available in 6mm, 8mm and 10mm sizes for accommodating a variety of tissues. With boxlock joint that permits smooth jaw closure movements, this tenaculum forceps are made from premium quality stainless steel and are manufactured to the highest standards of quality and precision.

MUSEUX UTERINE VULSELLUM FORCEPS FEATURES

- 9.5 Inches (24.1cm) Overall Length.
- 2 X 2 Overlapping Teeth for Firm Grip.
- 6mm Wide Teeth.
- Available in Straight and Curved Profile.
- Boxlock Joint for Smooth Jaw's Moment.



WI-70003373

Naegele Obstetrical Forceps

\$127.51 per piece

Trade price from **\$117.31**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Naegele Obstetrical Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$127.51, falling to \$117.31 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003593

New Eye Lash Straight Fine Point Color Tweezers

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

New Eye Lash Straight Fine Point Color Tweezers is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003592

New Fine Point Eyelash Extension Xtreme Lashes Tweezers

\$23.20 per piece

Trade price from \$21.34
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

New Fine Point Eyelash Extension Xtreme Lashes Tweezers is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$23.20, falling to \$21.34 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70002500**

Pean Hemostat 12\"/>

\$15.99 per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 12\"/>


FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other

tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once

**WI-70002501**

Pean Hemostat 12\"/>

\$15.99 per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 12\"/>



FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other

tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once



**WI-70002497**

Pean Hemostat 14\"/>

\$18.99 per pieceTrade price from **\$17.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 14\"/>



FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other

tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once

**WI-70002498**

Pean Hemostat 14\"/>

\$18.99 per pieceTrade price from **\$17.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 12\"/>


FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other

tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users



**WI-70002495**

Pean Hemostat 16\"/>

\$21.99 per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 16\"/>



FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other

tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once

**WI-70002496**

Pean Hemostat 16\"/>

\$21.99 per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 16\"/>


FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other

tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once



**WI-70000405****Pedo Extracting Forceps #150XAS Split Beak Serrated****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Pedo Extracting Forceps #150XAS Serrated Split Beaks for Upper Anterior – Premium Quality Stainless Steel American Pattern Child Tooth Extraction Dental Surgical Instrument.

Pedo extraction forceps #150XAS are specifically designed for extracting upper anteriors. Dental forceps no. 150XAS are universal tools for extracting teeth from alveolar bone, with an split beak serrated design with no spring design for optimal access and engagement. Our pediatric extraction forceps are designed to fit comfortably in the palm, reducing patient anxiety and American pattern feature a smooth hinge without a visible pin. 150XAS extracting forceps are made of high-quality stainless steel,

ensuring they can withstand daily use without causing any damage or trauma to the alveolar bone.

PEDO EXTRACTING FORCEPS #150XAS FEATURES

- 6" Overall Length.
- Pediatrics Forceps For Upper Anterior.
- Serrated Split Beaks Design.
- American Pattern Design.
- Textured Grip Textured Grip Handle.
- Fully Autoclavable and Reusable.
- Made of High Grade Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Pedo Extracting Forceps #150XAS Upper Anterior Split Beak Serrated

TECHNICAL INFORMATION**WI-70000140****Pedo Extracting Forceps #151XAS Diamond Dusted Gold Hollow Handle****\$14.99** per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Pedo Extracting Forceps #151XAS Diamond Dusted Serrated Split Beaks for Lower Anterior with Gold Hollow Python Grip Handle – Premium Quality Stainless Steel American Pattern Child Tooth Extraction Dental Surgical Instrument.

Universal pedo tooth extraction forceps #151XAS are specifically designed for extracting lower anteriors. Dental forceps no. 151XAS are universal tools for extracting teeth from alveolar bone, with an split beak serrated design with no spring design for optimal access and engagement. Our pediatric extraction forceps are designed to fit comfortably in the palm, reducing patient anxiety and American pattern feature a smooth hinge without a visible pin. Diamond dusted tips of 151XAS

extracting forceps ensuring they can withstand daily use without causing any damage or trauma to the alveolar bone. The gold hollow python grip handle enhances grip, reduces fatigue, and improves control, while stainless steel construction ensures durability and reliability, allowing for reuse over time.

PEDO EXTRACTING FORCEPS #151XAS FEATURES

- 6" Overall Length.
- Pediatrics Forceps For Lower Incisor, Premolars.
- Diamond Dusted Serrated Split Beaks Design.
- American Pattern Design.
- Gold Hollow Python Grip Handle.
- Fully Autoclavable and Reusable.
- Made of High Grade Stainless Steel.



**WI-70000406****Pedo Extracting Forceps #151XAS for Lower Anterior****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Pedo Extracting Forceps #151XAS Serrated Split Beaks for Lower Anterior – Premium Quality Stainless Steel American Pattern Child Tooth Extraction Dental Surgical Instrument.

Universal pedo tooth extraction forceps #151XAS are specifically designed for extracting lower anteriors. Dental forceps no. 151XAS are universal tools for extracting teeth from alveolar bone, with an split beak serrated design with no spring design for optimal access and engagement. Our pediatric extraction forceps are designed to fit comfortably in the palm, reducing patient anxiety and American pattern feature a smooth hinge without a visible pin. 151XAS extracting forceps are made of high-quality stainless steel, ensuring they can withstand daily use

without causing any damage or trauma to the alveolar bone.

PEDO EXTRACTING FORCEPS #151XAS FEATURES

- 6" Overall Length.
- Pediatrics Forceps For Lower Incisor, Premolars.
- Serrated Split Beaks Design.
- American Pattern Design.
- Textured Grip Textured Grip Handle.
- Fully Autoclavable and Reusable.
- Made of High Grade Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Universal Pedo Extracting Forceps #151XAS Serrated Split Beak

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**WI-70000132****Pedo Extracting Forceps #74N Narrow Beaks****\$12.99** per pieceTrade price from **\$11.95**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pedo Extracting Forceps #74N Extra Narrow Beaks for Lower Interiors and Roots with Textured Grip Handle Premium Quality Stainless Steel English Pattern Tooth Extraction Surgical Dental Instruments.

English pattern extraction forceps #74N designed for tooth extraction with long, sharp, serrated beaks for effective gum expansion and secure tooth grasp. Their extra narrow beaks are suitable for small, deeply situated roots. English extraction forceps use a visible screw at the hinge that brings the fulcrum closer to the working end. Pedo Extracting Forceps No. 74N is often used in minor oral surgeries to extract fractured roots or teeth. The ergonomically designed, knurled grip ensures comfortable handling. Made

from premium quality stainless steel, these forceps are sterilizable and reusable.

ENGLISH EXTRACTION FORCEPS #74N FEATURES

- English Pattern.
- Used for Lower Incisors and Roots.
- Narrow Thin Beaks Then #74.
- Serrated Jaws.
- Non-slip Textured-grip Handle.
- Tooth Extraction Tool for Minor Oral Surgeries.
- Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc English Extracting Forceps #74N – Extra Narrow Beaks



**WI-70001405****Pedo Extracting Forceps #A Dental Surgical Instrument****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pedo Upper Molar #A Universal, Mini Pedodontic-English Pattern Stainless Steel Dental Surgical Instruments.

Upper Primary English Pediatric Forceps used for extracting teeth from the alveolar bone. Tooth extracting forceps are used to extract tooth, designed for use in specific areas of the mouth. The beak is shaped to conform snugly to the contour of the tooth. This makes it easier to reach different teeth effectively. Dental Forceps handles are shaped so that a maximum amount of force can be applied to the beaks, while the handles are still in a comfortable position for the oral surgeon.

- Pedo Upper Molar.
- Universal Forceps used for Extracting Teeth From the Alveolar Bone.
- English Pattern Design.
- Serrated Handles – Better Grip.
- Reusable.
- Made of High Grade Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc English Pattern Pedo Extraction Forceps #A

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps**APPLICATION:** Extraction Forceps**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes**PEDO EXTRACTING FORCEPS #A FEATURES****WI-70000415****Pedo Extracting Forceps #B Lower Roots Universal Mini****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

English Pattern Pedo Extracting Forceps #B Lower Roots with Textured Grip Single Spring Handle – Premium Quality Stainless Steel Universal Mini Pedodontic Dental Surgical Childern Tooth Extraction Instrument.

Pedo lower root universal mini #B pedodontic english pattern extraction forceps is used to remove teeth from the alveolar bone.

Pedodontic extracting forceps #B are durable and sturdy, specifically designed for pediatric dentistry, ensuring precise and comfortable extractions. The pedo extracting forceps number B features a serrated beak for precise extraction, a practical spring action device for faster procedure times, and comfortable handles for the oral surgeon.

PEDO EXTRACTING FORCEPS #B FEATURES

- 4.5" Overall Length.
- For Lower Roots.
- Universal Mini Pedodontic.
- Serrated Beak for Strong Graps.
- English Pattern Design.
- Textured Grip Single Spring Handle.
- Fully Autoclavable and Reusable.
- Made of High Grade Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc English Pattern Pedo Extraction Forceps #B

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps – Pediatric**APPLICATION:** Lower Roots Extraction

**WI-70000414****Pedo Extracting Forceps #C Lower Molars Universal Mini****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

English Pattern Pedo Extracting Forceps #C Lower Molars with Textured Grip Single Spring Handle – Premium Quality Stainless Steel Universal Mini Pedodontic Dental Surgical Children Tooth Extraction Instrument.

Universal extracting forceps #C are dental instruments for precise and efficient extraction of lower molars in pediatric patients, designed for versatility and durability. Pedodontic extracting forceps #C are durable and sturdy, specifically designed for pediatric dentistry, ensuring precise and comfortable extractions. The pedo extracting forceps C features a Tri-pointed tip design on top beak for precise extraction, a practical spring action device for faster procedure

times, and comfortable handles for the oral surgeon.

PEDO EXTRACTING FORCEPS #C FEATURES

- 4.5" Overall Length.
- For Lower Molars.
- Universal Mini Pedodontic.
- Tri-pointed Tip Design on Top Beak.
- English Pattern Design.
- Textured Grip Single Spring Handle.
- Fully Autoclavable and Reusable.
- Made of High Grade Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc English Pattern Pedo Extraction Forceps #C

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental**WI-70000412****Pedo Extracting Forceps #D Lower Incisors Klein Pedodontic****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pedo Extracting Forcep #D Lower Incisors with Textured Grip Single Spring Handle – Premium Quality Stainless Steel Baby Extracting Forceps English Pattern Klein Fig. 5 Pedodontic Dental Instrument

Fig.5 Klein Pattern Tooth Extraction Pedo Forceps #D is a precision and efficient surgical instrument for atraumatic extractions of lower incisors in pediatric patients. Pedodontic extracting forceps D, is designed for children's smallest teeth, with closely situated beaks and small working end. Its spring minimizes pressure force, increases joint lifetime, and features ergonomic grip handles for comfort. English pattern forceps are designed with a fulcrum closer to the working

end, enhancing leverage and facilitating extraction.

PEDO EXTRACTING FORCEPS #D FEATURES

- 4.5" Overall Length.
- For Lower Incisors & Anterior.
- Klein Pattern Tooth Extraction Forceps Fig. 5.
- Small, Closely-situated, Serrated Beaks for Maximum Control.
- English Pattern Design.
- Textured Grip Single Spring Handle.
- Fully Autoclavable and Reusable.
- Made of High Grade Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc English Pattern Pedo Extraction Forceps #D

TECHNICAL INFORMATION

**WI-70000413****Pedo Extracting Forceps #E Upper Front Universal Mini****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

English Pattern Pedo Extracting Forcep #E Upper Front with Textured Grip Single Spring Handle – Premium Quality Stainless Steel Universal Mini Pedodontic Dental Surgical Children Tooth Extraction Instrument.

#E Pedo extracting forceps are durable and versatile tools designed for pediatric dentistry, ensuring precise and comfortable extractions of teeth from alveolar bone. Pedodontic extracting forceps E, designed for comfort and convenience, feature a standard closed beak, a practical spring action, and comfortable handles for faster extraction procedure. English pattern forceps are designed with a fulcrum closer to the working end, enhancing leverage and facilitating extraction.

PEDO EXTRACTING FORCEPS #E FEATURES

- 4.5" Overall Length.
- For Upper Front .
- Universal Mini Pedodontic.
- Standard Closed Beak.
- English Pattern Design.
- Textured Grip Single Spring Handle.
- Fully Autoclavable and Reusable.
- Made of High Grade Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc English Pattern Pedo Extraction Forceps #E

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps – Pediatric

APPLICATION: Upper Front Extraction**WI-70001404****Pedo Extracting Forceps #F Dental Surgical Instrument****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Extracting Forceps Pedo Upper Front, Universal Mini Pedodontic-English Pattern Stainless Steel Dental Surgical Instruments.

This Pedo upper front #F Pedodontic English Extraction Forceps is used to remove teeth from upper primary molars. Tooth extracting forceps are used to extract tooth, designed for use in specific areas of the mouth. The beak is shaped to conform snugly to the contour of the tooth. This makes it easier to reach different teeth effectively. Dental Forceps handles are shaped so that a maximum amount of force can be applied to the beaks, while the handles are still in a comfortable position for the oral surgeon.

PEDO EXTRACTING FORCEPS #F FEATURES

- Pedo Upper Front
- Used to Extract Upper Primary Molars
- Open Tip English Pattern Design
- Serrated Handles – Better Grip
- Made of High Grade Stainless Steel

PACKAGE INCLUDE:

- 1 Pc English Pattern Pedo Extraction Forceps #F

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps – Pediatric

APPLICATION: Extraction Forceps**POLISHED:** Yes**STERILE:** Yes**REUSABLE:** Yes

**WI-70000418****Pedo Extracting Forceps #H Upper & Lower Universal Mini****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

English Pattern Pedo Extracting Forceps #H for Upper and Lower with Textured Grip Single Spring Handle – Premium Quality Stainless Steel Universal Mini Pedodontic Dental Surgical Children Tooth Extraction Instrument.

Pedo upper and lower #H pedodontic english pattern universal mini extraction forceps is used to remove teeth from the alveolar bone. The #H pediatric extraction forceps are durable and sturdy, specifically designed for pediatric dentistry, ensuring precise and comfortable extractions. The pedo extracting forceps number H features a serrated beak for precise extraction, a practical spring action device for faster procedure times, and comfortable handles for the oral surgeon.

PEDO EXTRACTING FORCEPS #H FEATURES

- 4.5" Overall Length.
- For Upper and Lower.
- Universal Mini Pedodontic.
- Serrated Beak for Strong Graps.
- English Pattern Design.
- Textured Grip Single Spring Handle.
- Fully Autoclavable and Reusable.
- Made of High Grade Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc English Pattern Pedo Extraction Forceps #H

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

Extracting Forceps – Pediatric

APPLICATION: Extraction Forceps**WI-70002919****Peers-Bertram Towel Forceps 5.75" Serrated Tips with Grooves****\$17.99** per pieceTrade price from **\$16.55**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Peers-Bertram Towel Forceps 5.75" (15cm) Inward Curved Serrated Jaws with Groove for Tube Fixation Premium Quality Stainless Steel Surgical Medical Obstetrics & Gynecology Instruments.

Peers Towel Clamp is a specialized tool utilized by surgeons to secure sterile drapes and suction tubes, enhancing the operative site's visibility. The Peers Non-Perforating Towel Clamp is a high-quality instrument with inward curved jaws, serrations, and a central groove for tube handling. This towel clamp has ergonomically designed finger rings to prevent slippage and fatigue during procedures.


PEERS TOWEL CLAMP FEATURES

- 5.75" (15cm) Overall Length.
- Inward Curved Jaws

- Serrated Central Groove Tip for Tube Fixation.
- Box Joint.
- Ratchet Locking Mechanism.
- Round Finger Ring Handles.
- Fully Autoclavable and Reusable.
- Made From High-quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Peers-Bertram Towel Clamp 5.75"

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Towel Clamps**APPLICATION:** Use For Surgical Drapes**TIPS:** Serrated Central Groove Tip
HANDLE: Ratchet Lock Finger Ring Handle**REUSABLE:** Yes

**WI-70003372****Piper Obstetrical Forceps 17 1/2" (44.5cm)****\$49.99** per pieceTrade price from **\$45.99**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Piper Obstetrical Forceps 17 1/2" (44.5cm) Stainless Steel Surgical Gynecology & Obstetrics Instruments Forceps for Delivery of Aftercoming Head of the Breech presentation

align the fetal head with the birth canal. Also, the handle features a lock system for self-retaining the jaws, which grants excellent control to rotate the baby if necessary.


Piper Obstetrical Forcep is a specialized instrument that obstetric specialists used during complicated vaginal delivery. Its principal use is to delicately clamp the fetal head to permit the baby's delivery. Piper Obstetrical Forceps are used to grasp, rotate, and pull the fetal head out of the vaginal canal, especially during forceps delivery and breech delivery procedures.

PIPER OBSTETRICAL FORCEPS FEATURES

- Long Shank, Curved Blades Promoting Deep Surgical Reach.
- Fenestrated Perineal Curve Jaws to Delicately Align the Fetal Head
- Lock System to Achieve a Self-retaining Action.
- Ergonomic Grip Handle Ensuring Maximum Control.

PACKAGE INCLUDE:

- 1 Piece Piper Obstetrical Forceps 17 1/2" (44.5cm) Gynecology &



Piper forceps have long shanks that feature a pair of large fenestrated blades including a perineal curve to

**WI-70004064****Plate and Bone Holding Forceps 5.5" with Footplate****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Plate and Bone Holding Forceps 5.5" (14cm) One Pointed and One Footplate Tip with Long Ratchet Finger Ring Handle Premium Quality Stainless Steel Surgical Orthopedic Instruments.


for optimal maneuvering. Made with premium grade stainless steel ensuring the highest quality and durability standards in your operating rooms. Fully autoclaved and can be reuse.

Plate bone holding forceps is a specialized orthopedic tool that surgeons commonly use to hold and guide plates and bone fragments into an optimal position, in order to align the bone structure and perform transfixion and other procedures.

PLATE BONE HOLDING FORCEPS FEATURES

- Crafted with 5.5" (14cm) Overall Length.
- Inward Curved Tips.
- One Pointed Tip and Footplate Serrated Tip for Strong Grasping.
- Long Ratchet for Controlling the Clamping Degree.
- Ergonomic Finger Rings Ensuring Superior Control.
- Ideal for Grasping Bone and Tissue Fragments.
- Fully Autoclaved and Reusable.

These forceps comes with a pair of long jaws that curve inward at the tips. The one tip has a pointed profile and the opposite has a hammer-like shape with serrations. As a result, surgeons can hold stabilizing plates and align long bone fragments. Long Ratcheted Finger rings ensure a comfortable grip



**WI-70000650****Plate and Bone Holding Forceps 5.5" with Footplate & Gold Handle****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Plate and Bone Holding Forceps 5.5" (14cm) One Pointed and One Footplate Tip with Long Ratchet Finger Ring Gold Handle Premium Quality Stainless Steel Surgical Orthopedic Instruments.

Plate bone holding forceps is a specialized orthopedic tool that surgeons commonly use to hold and guide plates and bone fragments into an optimal position, in order to align the bone structure and perform transfixion and other procedures. These bone holding forceps comes with a pair of long jaws that curve inward at the tips. The one tip has a pointed profile and the opposite has a hammer-like shape with serrations. As a result, surgeons can hold stabilizing plates and align long bone fragments. Long ratcheted finger rings gold handles ensure a

comfortable grip for optimal maneuvering. Made with premium grade stainless steel ensuring the highest quality and durability standards in your operating rooms. Fully autoclaved and can be reuse.

PLATE BONE HOLDING FORCEPS FEATURES

- Crafted with 5.5" (14cm) Overall Length.
- Inward Curved Tips.
- One Pointed Tip and Footplate Serrated Tip for Strong Grasping.
- Long Ratchet for Controlling the Clamping Degree.
- Ergonomic Finger Rings Gold Plated Handle Ensuring Superior Control.
- Ideal for Grasping Bone and Tissue Fragments.

**WI-70000432****Platypus Nail Pulling Forceps 5.5" Narrow Jaws****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Platypus Nail Pulling Forcep 5½" (14cm) Narrow Flat Jaws with Pulling Teeth Premium Quality Stainless Steel Toe's Ingrown Nails Puller/Remover Surgical Podiatry Instruments.

Platypus Nail Puller are highly versatile instruments used during many nail surgeries. It is widely used to split and cut the toe's ingrown nails. Platypus nail pulling forceps feature narrow jaws used to split ingrown toenails and then remove the severed portion of the nail from the soft tissue. This platypus nail forceps has straight flat jaw provides twice as many pulling teeth as the narrower designed to help users to feel ease and less fatigue during surgical procedures. The plier style smooth handle ensure easy and better gripping. This podiatry nail tool made of stainless

steel to ensures maximum durability and dependability. You can reuse these platypus nail forceps over time as they will not bend or tarnish on you.


PLATYPUS NAIL PULLING FORCEP FEATURES:

- 5.5" (14cm) Overal Length.
- Straight Narrow Jaws.
- Pulling Inner Teeth for Strong Grip.
- Smooth Plier Style Handle.
- Adjustable Screw.
- Reusable and Rust-Proof.
- Autoclavable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Platypus Nail Pulling Forceps 5½" (14cm) Narrow Jaws





WI-70000999

Platypus Nail Pulling Forceps 5.5" Standard Wide Jaws

\$10.99 per piece

Trade price from \$10.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Platypus Nail Pulling Forceps 5½" (14cm) Standard 9mm Wide Flat Jaws with Pulling Teeth Premium Quality Stainless Steel Toe's Ingrown Nails Puller/Remover Surgical Podiatry Instruments.

Platypus Nail Puller are highly versatile instruments used during many nail surgeries. It is widely used to split and cut the toe's ingrown nails. Platypus nail pulling forceps feature standard or narrow jaws used to split ingrown toenails and then remove the severed portion of the nail from the soft tissue. This platypus nail forceps has slightly wider 9mm straight flat jaw provides twice as many pulling teeth as the narrower designed to help users to feel ease and less fatigue during surgical procedures. The plier style smooth

handle ensure easy and better gripping. This podiatry nail tool made of stainless steel to ensures maximum durability and dependability. You can reuse these platypus nail forceps over time as they will not bend or tarnish on you.

PLATYPUS NAIL PULLING FORCEPS FEATURES:

- 5.5" (14cm) Overall Length.
- Straight 9mm Standard Wide Jaws.
- Pulling Inner Teeth for Strong Grip.
- Smooth Plier Style Handle.
- Adjustable Screw.
- Reusable and Rust-Proof
- Autoclavable
- Made with Premium Quality Stainless Steel

TECHNICAL INFORMATION



WI-00849593

Pozzi Tenaculum Forceps 25cm

\$2.98 per piece

Trade price from \$1.80
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Pozzi Tenaculum Forceps 25cm is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 25cm size shown is our standard; other sizes in the same pattern are made to order.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$2.98, falling to \$1.80 at 10 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47256124

Pozzi Tenaculum Forceps Gynecology Surgical Instrument for Secure Cervical

\$4.90 per piece

Trade price from \$2.50
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Pozzi Tenaculum Forceps Gynecology Surgical Instrument for Secure Cervical is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$4.90, falling to \$2.50 at 25 pieces — about 49% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001428

Pratt Smith Hemostatic Forceps 6" T-shaped Jaws with Delicate Serration

\$13.99 per piece

Trade price from \$12.87
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pratt Smith Hemostatic Forceps 6" (15.24cm) T-shaped Concave Jaws with Delicate Serration Stainless Steel Hemorrhoidal, Rectal, Surgical Gynecology Instruments

Pratt T-Shaped Clamp, also known as Pratt-Smith Grasping Forceps, is a ratcheted, finger ring, grasping instrument typically used in vaginal or perineal procedures for grasping mucosa or hemorrhoidal tissue. The Pratt Artery Forceps is used to hold and steady the cervix, uterus or blood vessels during gynecological procedures including myomectomy, oophorectomy, hysterectomy, uterine sling operation and cervical tear repair. These forceps can be sterilized for reuse in surgical procedures. Pratt Artery Forceps can be sterilized for

reuse in surgical procedures. The Pratt Haemostatic Forceps features finger rings handle that avoid accidental slippage. also has a ratchet mechanism to self-lock the jaws at different levels of pressure.

PRATT SMITH HAEMOSTATIC FORCEPS FEATURES

- T-Shaped Concave Jaws Straight Hemostatic Clamp.
- Fine Edge Serrations to Ensure Strong Hold and Maximum Tissue Fixation.
- Ergonomic Finger Rings for Superior Handling
- Ratcheted Lock Mechanism.
- Autoclaved
- Made with German Grade Stainless Steel





WI-70000932

Precision Gripper Splinter Tweezer 2" (5.1cm) Sharp Tip

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Precision Gripper Splinter Tweezer 2" (5.1cm) Sharp Pointed Tips Spring-Tempered with Tube Holder Stainless Steel Splinter Remover Tweezers First Aid Surgical Instrument

This multipurpose precision splinter tweezer can be used in numerous ways, including removing the most stubborn and tiniest of splinters, slivers, ticks, threads, hair, stingers, pricklers, and more. This Precision Splinter Tweezer ideal for camping, boating, gardening, sewing, woodworking/carpentry, machining, first aid, crafts, and more. A must for any parent or child caregiver. It comes with a convenient and reclosable plastic tube with a cap. This Sliver Gripper Precision Splinter Tweezer with tube packing is recommended for kits,

pockets, packs, etc. The sharp-pointed splinter tweezer is made from high-quality stainless steel with machined pin-point precision, which makes this rugged, durable tweezer a much-needed piece of equipment for adventurers and survivalists.

PRECISION SPLINTER TWEEZER FEATURES:

- 2" (5.1cm) Long Overall.
- Precision Sharp Pointed Tips for Pin-point Gripping.
- Ideal to Remove the Tiniest of Splinters, Threads, Hairs, Stingers, and More.
- Single Piece Spring-tempered Stainless Steel Tweezer to Keeps Perfect Alignment.



WI-46195509

Precision Vasectomy Meatotomy Instrument Set Forceps Male

\$150.00 per piece

Trade price from \$100.00
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Precision Vasectomy Meatotomy Instrument Set Forceps Male — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$150.00, falling to \$100.00 at 2 sets — about 33% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003589

Professional Eyelash Extension Tweezers
Perfect Extra Fine Point Curved Best Suppliers
Eyelash Extension Tweezers

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Professional Eyelash Extension Tweezers Perfect Extra Fine Point Curved Best Suppliers Eyelash Extension Tweezers is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion

resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000822

Rainbow Iris Forceps 4" Serrated Straight

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Rainbow Iris Forceps 4" (10.2cm) Serrated Straight Standard Pattern 0.8mm Tips Premium Quality Stainless Steel Slip-Resistant Grip Eye Dressing Ophthalmological Instrument


Iris Dressing Forceps are used for holding dressing materials such as cotton and gauze during surgical procedures, changing dressings, or packing wounds. These iris eye dressing serrated forceps are uniquely designed for ophthalmological and wound care procedures. The rainbow multicolor finish brings a touch of vibrancy to your tools kit. Featuring fine tips made of stainless steel, these forceps ensure meticulous handling. Surgeons widely use them while removing bandages. Its delicate design assures easy manipulation. The tools

are available in multiple variations. You can pick the best one of your choice according to the procedure. Crafted from high quality German surgical stainless steel ensure corrosion resistance and superior craftsmanship.

RAINBOW IRIS FORCEPS
FEATURES

- Crafted with 4" (10.2cm) Overall Length.
- Rainbow Multicolor Finish.
- Standard Pattern Straight Tips.
- Serrated Jaws for Secure Holding.
- 0.8mm Wide Tips.
- Long-Reaching Jaws for Deep Surgical Access.
- Ergonomic Spring Handle Ensuring Superior Control.



**WI-70000821****Rainbow Splinter Forceps 4.5" Pointed Tips****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Rainbow Splinter Forceps 4.5" (11.4cm)
Needle Point Serrated Straight Tips
Premium Quality Stainless Steel Slip-Resistant Grip Medical Surgical Instrument.

Splinter forceps use for removing splinters and debris on the skin. Splinter Tweezers are 4.5 inches long. Great forcep choice for first aid boxes. With fine-pointed tips and a precision gripping mechanism, these splinter forceps excel in delicate tasks, making them perfect for DIY enthusiasts in various crafting projects. Splinter forceps features tapered precision tips can be used to remove ingrown hairs and particles embedded under the skin including splinters, thorns and glass. Amazingly precise hand-finished points.

The durable stainless steel construction ensures longevity. The rainbow design adds a touch of aesthetics to these tweezers, while the non-slip handle provides a secure grip for precise maneuvers. Crafted from high quality German surgical stainless steel ensure corrosion resistance and superior craftsmanship.

RAINBOW SPLINTER FORCEPS FEATURES

- Crafted with 4.5" (11.4cm) Overall Length.
- Rainbow Multicolor Finish.
- Pointed Straight Tips.
- Serrated Jaws for Secure Holding.
- Long-Reaching Jaws for Deep Surgical Access.

**WI-70000553****Ramsey Sponge Holding Forcep 7" Straight Serrated Jaws****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Ramsey Sponge Forceps 7" (17.8cm)
Straight Serrated Loop-Shaped Oval Fenestration Jaws Premium Grade Stainless Steel Dressing & Sponge Holding Medical, Surgical Instruments.

Ramsey Sponge Holding Forcep is a versatile surgical tool used by surgeons to grab and hold swabs and sponges or to clean blood and other fluids from the operating site. This sponge holding forceps features a long loop-shaped jaws with oval fenestration and horizontal serrations for efficient sponge handling. Ramsey Sponge Forceps has elongated shanks connect with the boxlock and a ratchet mechanism allows a smooth clamping action for preventing inadvertent injury. The Stainless Steel construction ensures maximum durability and

dependability. You can reuse these Hemostat Sponge Ramsey Forceps over time as they will not bend or tarnish on you. Easy to clean and can be reused after sterilization.

RAMSEY SPONGE HOLDING FORCEPS FEATURES

- 7" (17.8cm) Overall Length.
- Loop-shaped Oval Fenestrated Jaws.
- Straight Profile.
- Atraumatic Serrated Jaws for Securing a Strong Grasp.
- Versatile Ratchet Mechanism for Self-retaining Clamping.
- Ergonomic Finger Rings for Enhanced Handling.
- Elongated Shanks Connected with Box Joint.



**WI-70002019****Randall Kidney Stone Forceps 7.5\"/>****\$13.99** per pieceTrade price from **\$12.87**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Randall Kidney Stone Forceps 7.5\"/>

Randall Kidney Stone Forceps is a specialized instrument that urology and general surgeons commonly use to grasp, mobilize and extract stones from the kidney. Its principal use is to assist surgeons in stone removal during percutaneous nephrostomy and nephrolithotomy.

These kidney stone forceps mainly hold the deep tissues into the body cavity during surgical procedures. The jaws have an oval fenestration that improves sightlines and is ideal to visualize the grabbed stone and finger rings help to minimize the risk of slippage and fatigue. These surgical forceps

standard pattern have ability that's why is highly suitable for removal of kidney stones, gall stones, or polyps. The jaws are available with four varying degrees of curvature.

RANDALL KIDNEY STONE FORCEPS FEATURES

- Serrated Jaws for enhanced gripping.
- Slightly & strongly curved jaws to suit all needs.
- Full & Slightly inward curved shank with varying degrees of curvature.
- Ergonomic non-ratcheted finger rings ensuring maximum control.
- Made with High-Quality Stainless Steel

PACKAGE INCLUDES:**WI-70002017****Randall Kidney Stone Forceps 7.75\"/>****\$13.99** per pieceTrade price from **\$12.87**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Randall Kidney Stone Forceps 7.75\"/>

Randall Kidney Stone Forceps is a specialized instrument that urology and general surgeons commonly use to grasp, mobilize and extract stones from the kidney. Its principal use is to assist surgeons in stone removal during percutaneous nephrostomy and nephrolithotomy.

These kidney stone forceps mainly hold the deep tissues into the body cavity during surgical procedures. The jaws have an oval fenestration that improves sightlines and is ideal to visualize the grabbed stone and finger rings help to minimize the risk of slippage and fatigue. These surgical forceps

standard pattern have ability that's why is highly suitable for removal of kidney stones, gall stones, or polyps. The jaws are available with four varying degrees of curvature.

RANDALL KIDNEY STONE FORCEPS FEATURES

- Serrated Jaws for enhanced gripping.
- Slightly & strongly curved jaws to suit all needs.
- Full & Slightly inward curved shanks with varying degrees of curvature.
- Ergonomic non-ratcheted finger rings ensuring maximum control.
- Made with High-Quality Stainless Steel

PACKAGE INCLUDES:

**WI-70002020****Randall Kidney Stone Forceps 8.25\" Half Curved****\$13.99** per pieceTrade price from **\$12.87**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Randall Kidney Stone Forceps 8.25\" (21cm) Half Curved Serrated Jaws Stainless Steel Urology Surgical Instruments.

Randall Kidney Stone Forceps is a specialized device that urology and general surgeons commonly use to grasp, mobilize and extract stones from the kidney. Its principal use is to assist surgeons in stone removal during percutaneous nephrostomy and nephrolithotomy.

These kidney stone forceps mainly hold the deep tissues into the body cavity during surgical procedures. The jaws have an oval fenestration that improves sightlines and is ideal to visualize the grabbed stone and finger rings help to minimize the risk of slippage and fatigue. These surgical forceps

standard pattern have ability that's why is highly suitable for removal of kidney stones, gall stones, or polyps. The jaws are available with four varying degrees of curvature.

RANDALL KIDNEY STONE FORCEPS FEATURES

- Serrated Jaws for enhanced gripping.
- Slightly & strongly curved jaws to suit all needs.
- Full & Slightly inward curved shanks with varying degrees of curvature.
- Ergonomic non-ratcheted finger rings ensuring maximum control.
- Made with High-Quality Stainless Steel

PACKAGE INCLUDES:**WI-70002018****Randall Kidney Stone Forceps 9.25\" Quarter Curved****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Randall Kidney Stone Forceps 9.25\" (23.5cm) Quarter Curved Fenestration Serrated Jaws Stainless Steel Urology Surgical Instruments.

Randall Kidney Stone Forceps is a specialized instrument that urology and general surgeons commonly use to grasp, mobilize and extract stones from the kidney. Its principal use is to assist surgeons in stone removal during percutaneous nephrostomy and nephrolithotomy.


These kidney stone forceps mainly hold the deep tissues into the body cavity during surgical procedures. The jaws have an oval fenestration that improves sightlines and is ideal to visualize the grabbed stone and finger rings help to minimize the risk of slippage and fatigue. These surgical forceps

standard pattern have ability that's why is highly suitable for removal of kidney stones, gall stones, or polyps. The jaws are available with four varying degrees of curvature.

RANDALL KIDNEY STONE FORCEPS FEATURES

- Serrated Jaws for enhanced gripping.
- Slightly & strongly curved jaws to suit all needs.
- Full & Slightly inward curved shanks with varying degrees of curvature.
- Ergonomic non-ratcheted finger rings ensuring maximum control.
- Made with High-Quality Stainless Steel

PACKAGE INCLUDES:



WI-70003525

Rankin Kelly Forceps

\$13.99 per piece

Trade price from \$12.87
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Rankin Kelly Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$13.99, falling to \$12.87 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36263969

Reusable 8 Inch Babcock Forceps General Surgery Instrument

\$2.30 per piece

Trade price from \$1.50
Minimum order 8 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Reusable 8 Inch Babcock Forceps General Surgery Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

The 8" size shown is our standard; other sizes in the same pattern are made to order. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.30, falling to \$1.50 at 8 pieces — about 35% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36900116

Reusable Axis Traction Forceps for Obstetrical Surgery

\$80.00 per piece

Trade price from **\$60.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Reusable Axis Traction Forceps for Obstetrical Surgery is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$80.00, falling to \$60.00 at 10 pieces — about 25% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35071878

Reusable Babcock Tissue Grasping Forceps Stainless Steel Surgical Tool

\$1.45 per piece

Trade price from **\$0.76**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Reusable Babcock Tissue Grasping Forceps Stainless Steel Surgical Tool — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.45, falling to \$0.76 at 10 pieces — about 48% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28193341

Reusable Medical Equipment Surgical Stainless Steel Forceps Sponge

\$8.00 per piece

Trade price from **\$4.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Reusable Medical Equipment Surgical Stainless Steel Forceps Sponge is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-45080905

Reusable Surgical Dressing Forceps Stainless Steel Medical Tissue

\$2.40 per piece

Trade price from **\$2.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Reusable Surgical Dressing Forceps Stainless Steel Medical Tissue, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.40, falling to \$2.00 at 25 pieces — about 17% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001104****ROBA English Pattern Extracting Forceps #35A Upper Premolars****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty ROBA Extracting Forceps 35A Upper Premolars Engling Pattern Special Beaks Far-reaching German Grade Stainless Steel Teeth Extraction Forceps with Gold Hollow Python Grip Handle.

English pattern dental extracting forceps #35A for Upper Premolars features far-reaching jaws with special beaks with easy gripping gold color hollow python grip handles..

RoBa extracting forceps have been designed to allow subgingival access to the root while removing a tooth. The beaks are anatomically shaped, tapered, serrated and notched to allow the forceps to penetrate deeper into the tooth socket and to maximize the tooth-to- forceps contact. The notched beaks provide an enhanced grip that assists

atraumatic extractions and reduces the risk of fracturing the root structure or causing trauma to the surrounding hard and soft tissues. Surgical grade stainless steel promotes longevity and corrosion resistance.

ROBA EXTRACTING FORCEPS #35A FEATURES

- For Upper Premolars
- Special Beaks
- Far-reaching Model
- English Pattern
- Gold Hollow Python Grip Handle
- Made with High Quality German Stainless Steel

PACKAGE INCLUDE

- 1 Pc English Pattern Extraction Forceps #35A – Serrated – Gold Hollow Python Grip Handle

**WI-70001241****ROBA Extracting Forceps #34M Upper Incisors and Cuspids Serrated Jaws****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty ROBA Extraction Forceps Fig.34M Modif. English Pattern, ROBA Edition by Dr. Beck, Upper Incisors and Eye Teeth German Grade Stainless Steel Dental Extracting Forceps with Gold Hollow Python Grip Handle.

Roba Edition Extracting Forceps 34M use for upper jaw, incisors & cuspids with wider eye teeth tips. The special beaks provide deep & parallel grip. It has gold color hollow python grip handles.

RoBa extracting forceps have been designed to allow subgingival access to the root while removing a tooth. The beaks are anatomically shaped, tapered, serrated and notched to allow the forceps to penetrate deeper into the tooth socket and to maximize the tooth-to- forceps contact. The notched

beaks provide an enhanced grip that assists atraumatic extractions and reduces the risk of fracturing the root structure or causing trauma to the surrounding hard and soft tissues. Surgical grade stainless steel promotes longevity and corrosion resistance.

ROBA EXTRACTING FORCEPS #34M FEATURES

- Upper Incisors and Canines
- Wider Eye Teeth Tips
- Special Serrated Beaks for Deep & Parallel Grip
- English Pattern
- Gold Hollow Python Grip Handle
- Made with High Quality German Stainless Steel

PACKAGE INCLUDE

**WI-70001240****ROBA Extracting Forceps #34N Upper Incisors and Eye Teeth****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty ROBA Extraction Forceps Fig.34N Modif. English Pattern, ROBA Edition by Dr. Beck, Upper Incisors and Eye Teeth German Grade Stainless Steel Dental Extracting Forceps with Gold Hollow Python Grip Handle.

Roba Edition Extracting Forceps 34N use for upper jaw, incisors & cuspids with narrow eye teeth tips. The special beaks provide deep & parallel grip. It has gold color hollow python grip handles.

RoBa extracting forceps have been designed to allow subgingival access to the root while removing a tooth. The beaks are anatomically shaped, tapered, serrated and notched to allow the forceps to penetrate deeper into the tooth socket and to maximize the tooth-to- forceps contact. The notched

beaks provide an enhanced grip that assists atraumatic extractions and reduces the risk of fracturing the root structure or causing trauma to the surrounding hard and soft tissues. Surgical grade stainless steel promotes longevity and corrosion resistance.

ROBA EXTRACTING FORCEPS #34N FEATURES

- Upper Incisors and Canine.
- Narrow Eye Teeth Tips.
- Special Serrated Beaks for Deep & Parallel Grip
- English Pattern
- Gold Hollow Python Grip Handle
- Made with High Quality German Stainless Steel

PACKAGE INCLUDE**WI-70001103****ROBA Extracting Forceps #35M Upper Premolars****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

English Pattern Dental Extracting Forceps, 35M Modif. Roba-Edition by Dr. Beck, Upper Premolars, Deep Serrated Angled Jaws German Grade Stainless Steel Teeth Extraction Forceps Dental Surgical Instruments.

ROBA extrating forceps #35M for Upper Premolars features deep serrated jaws with easy gripping gold color hollow python grip handles.

RoBa edition extraction forceps have been designed to allow subgingival access to the root while removing a tooth. The beaks are anatomically shaped, tapered, serrated and notched to allow the forceps to penetrate deeper into the tooth socket and to maximize the tooth-to- forceps contact. The notched beaks provide an enhanced grip that assists atraumatic

extractions and reduces the risk of fracturing the root structure or causing trauma to the surrounding hard and soft tissues. Surgical grade stainless steel promotes longevity and corrosion resistance.

ROBA EXTRACTING FORCEPS #35M FEATURES

- 35M Modif. Model
- For Upper Premolars
- Deep Serrated Jaws
- English Pattern
- Gold Hollow Python Grip Handle
- Made with High Quality German Stainless Steel

PACKAGE INCLUDE

- 1 Pc ROBA Extraction Forceps #35M – Deep Serrated Jaws





WI-46953081

Robin Anchoring Forceps Reusable Stainless Steel Surgical Tool

\$35.00 per piece

Trade price from \$21.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Robin Anchoring Forceps Reusable Stainless Steel Surgical Tool — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$35.00, falling to \$21.00 at 25 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46958051

Robin Anchoring Forceps Surgical Instrument Set

\$40.00 per piece

Trade price from \$25.00
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Robin Anchoring Forceps Surgical Instrument Set — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$40.00, falling to \$25.00 at 5 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46952103

Robin Anchoring Surgical Forceps Precision Instrument for Tissue Fixation

\$32.00 per piece

Trade price from \$18.00
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Robin Anchoring Surgical Forceps Precision Instrument for Tissue Fixation — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$32.00, falling to \$18.00 at 50 pieces — about 44% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001072

Rochester Carmalt Forceps 6.25" Curved

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Rochester Carmalt Forceps 6.25" (15.9cm) Curved Hemostatic Artery Forcep with Ratcheted Finger Ring Handle Premium Quality Stainless Steel Surgical Instrument.

Rochester Carmalt Forceps is a specialized hemostatic artery forcep that surgeons used to manipulate blood vessels, clamp their lumen and control the blood flow to prevent blood loss. It comes in straight and curved patterns. Carmalt Forceps features a longitudinal serrations that run the length of the blade with grid style pattern at the tip to provide the strongest grip. These large, crushing hemostatic forceps are a choice instrument for clamping blood vessels and large tissues or ligating pedicles.

The textured blades along with ratchet handles ensure a strong grip of the clamped tissue with minimal chance of cutting the tissue during a procedure. One application for the rochester carmalt forcep is an ovariectomy, where this instrument is used for crushing the vessels of the ovarian pedicle. Rochester Carmalt Forceps manufactured in premium grade stainless steel to ensure longevity with highest quality performance during surgical procedures.

ROCHESTER CARMALT FORCEPS FEATURES

- 6.25" (15.9cm) Overall Length.
- Straight or Curved Tip Patterns
- Longitudinal Serrated Blades with Grid style pattern on Tips



**WI-70001073****Rochester Carmalt Forceps 6.25" Straight****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Carmalt Forceps 6.25" (15.9cm) Straight Hemostatic Artery Forcep with Ratcheted Finger Ring Handle Premium Quality Stainless Steel Surgical Instrument.

Rochester Carmalt Forceps is a specialized hemostatic artery forcep that surgeons used to manipulate blood vessels, clamp their lumen and control the blood flow to prevent blood loss. It comes in straight and curved patterns. Carmalt Forceps features a longitudinal serrations that run the length of the blade with grid style pattern at the tip to provide the strongest grip. These large, crushing hemostatic forceps are a choice instrument for clamping blood vessels and large tissues or ligating pedicles.

The textured blades along with ratchet handles ensure a strong grip of the clamped tissue with minimal chance of cutting the tissue during a procedure. One application for the rochester carmalt forcep is an ovariectomy, where this instrument is used for crushing the vessels of the ovarian pedicle. Rochester Carmalt Forceps manufactured in premium grade stainless steel to ensure longevity with highest quality performance during surgical procedures.

ROCHESTER CARMALT FORCEPS FEATURES

- 6.25" (15.9cm) Overall Length.
- Straight or Curved Tip Patterns
- Longitudinal Serrated Blades with Grid style pattern on Tips

WI-70001070**Rochester Carmalt Forceps 8" Curved****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Rochester Carmalt Forceps 8" (20.3cm) Curved Hemostatic Artery Forcep with Ratcheted Finger Ring Handle Premium Quality Stainless Steel Surgical Instrument.

Rochester Carmalt Forceps is a specialized hemostatic artery forcep that surgeons used to manipulate blood vessels, clamp their lumen and control the blood flow to prevent blood loss. It comes in straight and curved patterns. Carmalt Forceps features a longitudinal serrations that run the length of the blade with grid style pattern at the tip to provide the strongest grip. These large, crushing hemostatic forceps are a choice instrument for clamping blood vessels and large tissues or ligating pedicles.

The textured blades along with ratchet handles ensure a strong grip of the clamped tissue with minimal chance of cutting the tissue during a procedure. One application for the rochester carmalt forcep is an ovariectomy, where this instrument is used for crushing the vessels of the ovarian pedicle. Rochester Carmalt Forceps manufactured in premium grade stainless steel to ensure longevity with highest quality performance during surgical procedures.

ROCHESTER CARMALT FORCEPS FEATURES

- 8" (20.3cm) Overall Length.
- Straight or Curved Tip Patterns
- Longitudinal Serrated Blades with Grid style pattern on Tips



WI-70001071

Rochester Carmalt Forceps 8" Straight

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Carmalt Forceps 8" (20.3cm) Straight Hemostatic Artery Forcep with Ratcheted Finger Ring Handle Premium Quality Stainless Steel Surgical Instrument.

Rochester Carmalt Forceps is a specialized hemostatic artery forcep that surgeons used to manipulate blood vessels, clamp their lumen and control the blood flow to prevent blood loss. It comes in straight and curved patterns. Carmalt Forceps features a longitudinal serrations that run the length of the blade with grid style pattern at the tip to provide the strongest grip. These large, crushing hemostatic forceps are a choice instrument for clamping blood vessels and large tissues or ligating pedicles.

The textured blades along with ratchet handles ensure a strong grip of the clamped tissue with minimal chance of cutting the tissue during a procedure. One application for the rochester carmalt forcep is an ovariectomy, where this instrument is used for crushing the vessels of the ovarian pedicle. Rochester Carmalt Forceps manufactured in premium grade stainless steel to ensure longevity with highest quality performance during surgical procedures.

ROCHESTER CARMALT FORCEPS FEATURES

- 8" (20.3cm) Overall Length.
- Straight or Curved Tip Patterns
- Longitudinal Serrated Blades with Grid style pattern on Tips



WI-27261542

Rochester Pean / Spencer Well Artery Forceps Curved/Straight 6/7/8/9/10

\$2.55 per piece

Trade price from \$1.55
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Rochester Pean / Spencer Well Artery Forceps Curved/Straight 6/7/8/9/10 is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: a curved working end for reach around structures, and a straight working end for direct, in-line access.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.55, falling to \$1.55 at 1 piece — about 39% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003524****Rochester Pean Forceps 14cm Curved****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-Pean Hemostat Forceps Curved 5.5" (14cm) Multipurpose Locking Tweezer Clamps Serrated, Stainless Steel – Hemostats for Nurses, Fishing Forceps, Crafts and Hobby.

ROCHESTER PEAN FORCEPS FEATURES



PEAN HEMOSTAT FORCEPS: Waqar International's Rochester Pean 14cm Curved hemostat clamps are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel,

these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

THREE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once locked in place, even the smallest objects can be held securely in the serrated jaws with great force.

A MUST-HAVE TOOL FOR REMOVING HOOKS: These Fishing clamps are one of the most essential pieces of equipment that you can have with you on a day on the water. They are ideal

**WI-70002515****Rochester Pean Hemostat Forceps 14cm Straight****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-Pean Hemostat Forceps Straight 5.5" (14cm) Multipurpose Locking Tweezer Clamps Serrated, Stainless Steel – Hemostats for Nurses, Fishing Forceps, Crafts and Hobby.

ROCHESTER PEAN FORCEPS FEATURES


PEAN HEMOSTAT FORCEPS: Waqar International's Rochester Pean 14cm Straight hemostat clamps are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel,

these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

THREE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once locked in place, even the smallest objects can be held securely in the serrated jaws with great force.

A MUST-HAVE TOOL FOR REMOVING HOOKS: These Fishing clamps are one of the most essential pieces of equipment that you can have with you on a day on the water. They are ideal



**WI-70002493****Rochester Pean Hemostat Forceps 18\"/>****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 18\"/>



tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once

FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other

**WI-70002494****Rochester Pean Hemostat Forceps 18\"/>****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 18\"/>


tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once

FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other



**WI-70002491**

Rochester Pean Hemostat Forceps 20\"/>

\$26.99 per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 20\"/>


tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once

FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other

**WI-70002492**

Rochester Pean Hemostat Forceps 20\"/>

\$26.99 per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 20\"/>


tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once

FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other



**WI-70002489****Rochester Pean Hemostat Forceps 24\"/>****\$31.99** per pieceTrade price from **\$29.43**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 24\"/>



tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once

FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other

**WI-70002490****Rochester Pean Hemostat Forceps 24\"/>****\$31.99** per pieceTrade price from **\$29.43**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester Pean Hemostat Forceps 24\"/>


tubular structures to obstruct the flow of blood or fluids. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel, these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

FIVE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once

FISHING HEMOSTAT FORCEPS FEATURES

PEAN HEMOSTAT FORCEPS: The Waqar International's Hemostat Forceps Locking Clamps are an essential in the clinical, medical, lab or veterinary environment. They are made of Premium Quality materials with superior craftsmanship. Ready to sterilize and use surgical-grade, stainless steel instrument that is used to clamp blood vessels or tag sutures. Typically used to compress blood vessels or other



**WI-70002513****Rochester Pean Hemostat Forceps 6.25\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-Pean Hemostat Forceps Curved 6.25" (16cm) Multipurpose Locking Tweezer Clamps Serrated, Stainless Steel – Hemostats for Nurses, Fishing Forceps, Crafts and Hobby.

ROCHESTER PEAN FORCEPS FEATURES



PEAN HEMOSTAT FORCEPS: Waqar International's Rochester Pean 16cm Curved hemostat clamps are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel,

these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

THREE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once locked in place, even the smallest objects can be held securely in the serrated jaws with great force.

A MUST-HAVE TOOL FOR REMOVING HOOKS: These Fishing clamps are one of the most essential pieces of equipment that you can have with you on a day on the water. They are ideal

**WI-70002514****Rochester Pean Hemostat Forceps 6.25\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-Pean Hemostat Forceps Straight 6.25" (16cm) Multipurpose Locking Tweezer Clamps Serrated, Stainless Steel – Hemostats for Nurses, Fishing Forceps, Crafts and Hobby.

ROCHESTER PEAN FORCEPS FEATURES


PEAN HEMOSTAT FORCEPS: Waqar International's Rochester Pean 16cm Straight hemostat clamps are ideal for personal or professional use as fishing forceps or in electronic repair, crafting, medical, model building, or any operation that requires a delicate touch. The large finger loops, precision tips, and serrated self-locking jaws are ideal for grasping small objects securely.

HIGH-QUALITY STAINLESS STEEL CONSTRUCTION: Hot-drop forged from medical-grade stainless steel,

these hemostats are autoclavable and won't rust, break, or even bend under pressure. Our production process passes the national and international quality and safety standards, these hemostat clamps are ideal for medical and veterinary settings as well as a key part of any fishing tackle kit.

THREE LOCKING POSITIONS: These hemostat clamps can be held in place by a locking forceps mechanism with three interlocking teeth that allow users to adjust the clamping force, once locked in place, even the smallest objects can be held securely in the serrated jaws with great force.

A MUST-HAVE TOOL FOR REMOVING HOOKS: These Fishing clamps are one of the most essential pieces of equipment that you can have with you on a day on the water. They are ideal



**WI-70000083****Rochester-Ochsner Forceps 5.5" Curved 1×2 Teeth Serrated Jaws****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-ochsner Hemostatic Artery Forceps 5.5" (14cm) Curved 1×2 Teeth Serrated Jaws with Ratcheted Ring Handle – Premium Quality Kocher Forceps for Clamping of Arteries during General, Orthopedic, and Gynecological Surgeries

Rochester-Ochsner forceps, also known as Kocher forceps, are robust surgical tools designed for grasping thick tissue and clamping heavy blood vessels. Commonly used in general, orthopedic, and gynecological surgeries, they prevent slippage while controlling bleeding. The forceps feature entire length horizontal serrations and a 1 × 2 tooth tip for secure gripping. With a ratchet locking mechanism, their ring handles allow hands-free operation. Available in

straight and curved patterns and various lengths, these stainless steel forceps are durable and fully sterilizable for repeated use in medical environments.

ROCHESTER OCHSNER FORCEPS FEATURES

- 5.5" (14cm) Overall Length.
- 1×2 "Rat-tooth" Tips for Secure Grip.
- Curved Profile.
- Horizontal Serrations Jaws on Entire Length for a Firm, Non-slip Grip..
- Boxlock Joint to Promote Smooth Closure Movements.
- Ratchet Lock Provides Self-retaining System.
- Ergonomic Ring Handles for Comfortable Handling.
- Anti Rust and Reusable.

WI-70000084**Rochester-Ochsner Forceps 5.5" Straight 1×2 Teeth Serrated Jaws****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-ochsner Hemostatic Artery Forceps 5.5" (14cm) Straight 1×2 Teeth Serrated Jaws with Ratcheted Ring Handle – Premium Quality Kocher Forceps for Clamping of Arteries during General, Orthopedic, and Gynecological Surgeries

Rochester-Ochsner forceps, also known as Kocher forceps, are robust surgical tools designed for grasping thick tissue and clamping heavy blood vessels. Commonly used in general, orthopedic, and gynecological surgeries, they prevent slippage while controlling bleeding. The forceps feature entire length horizontal serrations and a 1 × 2 tooth tip for secure gripping. With a ratchet locking mechanism, their ring handles allow hands-free operation. Available in

straight and curved patterns and various lengths, these stainless steel forceps are durable and fully sterilizable for repeated use in medical environments.

ROCHESTER OCHSNER FORCEPS FEATURES

- 5.5" (14cm) Overall Length.
- 1×2 "Rat-tooth" Tips for Secure Grip.
- Horizontal Serrations Jaws on Entire Length for a Firm, Non-slip Grip..
- Boxlock Joint to Promote Smooth Closure Movements.
- Ratchet Lock Provides Self-retaining System.
- Ergonomic Ring Handles for Comfortable Handling.
- Anti Rust and Reusable.

**WI-70000081****Rochester-Ochsner Forceps 6.25\"/>****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces



- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-ochsner Hemostatic Artery Forceps 6.25\" (15.9cm) Curved 1x2 Teeth Serrated Jaws with Ratcheted Ring Handle – Premium Quality Kocher Forceps for Clamping of Arteries during General, Orthopedic, and Gynecological Surgeries

straight and curved patterns and various lengths, these stainless steel forceps are durable and fully sterilizable for repeated use in medical environments.

ROCHESTER OCHSNER FORCEPS FEATURES

- 6.25\" (14cm) Overall Length.
- 1x2 "Rat-tooth" Tips for Secure Grip.
- Curved Profile.
- Horizontal Serrations Jaws on Entire Length for a Firm, Non-slip Grip..
- Boxlock Joint to Promote Smooth Closure Movements.
- Ratchet Lock Provides Self-retaining System.
- Ergonomic Ring Handles for Comfortable Handling.
- Anti Rust and Reusable.



Rochester-Ochsner forceps, also known as Kocher forceps, are robust surgical tools designed for grasping thick tissue and clamping heavy blood vessels. Commonly used in general, orthopedic, and gynecological surgeries, they prevent slippage while controlling bleeding. The forceps feature entire length horizontal serrations and a 1 x 2 tooth tip for secure gripping. With a ratchet locking mechanism, their ring handles allow hands-free operation. Available in

**WI-70000082****Rochester-Ochsner Forceps 6.25\"/>****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-ochsner Hemostatic Artery Forceps 6.25\" (15.9cm) Straight 1x2 Teeth Serrated Jaws with Ratcheted Ring Handle – Premium Quality Kocher Forceps for Clamping of Arteries during General, Orthopedic, and Gynecological Surgeries

straight and curved patterns and various lengths, these stainless steel forceps are durable and fully sterilizable for repeated use in medical environments.

ROCHESTER OCHSNER FORCEPS FEATURES

- 6.25\" (15.9cm) Overall Length.
- 1x2 "Rat-tooth" Tips for Secure Grip.
- Horizontal Serrations Jaws on Entire Length for a Firm, Non-slip Grip..
- Boxlock Joint to Promote Smooth Closure Movements.
- Ratchet Lock Provides Self-retaining System.
- Ergonomic Ring Handles for Comfortable Handling.
- Anti Rust and Reusable.



Rochester-Ochsner forceps, also known as Kocher forceps, are robust surgical tools designed for grasping thick tissue and clamping heavy blood vessels. Commonly used in general, orthopedic, and gynecological surgeries, they prevent slippage while controlling bleeding. The forceps feature entire length horizontal serrations and a 1 x 2 tooth tip for secure gripping. With a ratchet locking mechanism, their ring handles allow hands-free operation. Available in

**WI-70000079****Rochester-Ochsner Forceps 7.25\"/>****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-ochsner Hemostatic Artery Forceps 7.25\" (18.4cm) Curved 1x2 Teeth Serrated Jaws with Ratcheted Ring Handle – Premium Quality Kocher Forceps for Clamping of Arteries during General, Orthopedic, and Gynecological Surgeries

straight and curved patterns and various lengths, these stainless steel forceps are durable and fully sterilizable for repeated use in medical environments.

ROCHESTER OCHSNER FORCEPS FEATURES

- 7.25\" (18.4cm) Overall Length.
- 1x2 "Rat-tooth" Tips for Secure Grip.
- Curved Profile.
- Horizontal Serrations Jaws on Entire Length for a Firm, Non-slip Grip..
- Boxlock Joint to Promote Smooth Closure Movements.
- Ratchet Lock Provides Self-retaining System.
- Ergonomic Ring Handles for Comfortable Handling.
- Anti Rust and Reusable.



Rochester-Ochsner forceps, also known as Kocher forceps, are robust surgical tools designed for grasping thick tissue and clamping heavy blood vessels. Commonly used in general, orthopedic, and gynecological surgeries, they prevent slippage while controlling bleeding. The forceps feature entire length horizontal serrations and a 1 x 2 tooth tip for secure gripping. With a ratchet locking mechanism, their ring handles allow hands-free operation. Available in

**WI-70000080****Rochester-Ochsner Forceps 7.25\"/>****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-ochsner Hemostatic Artery Forceps 7.25\" (18.4cm) Straight 1x2 Teeth Serrated Jaws with Ratcheted Ring Handle – Premium Quality Kocher Forceps for Clamping of Arteries during General, Orthopedic, and Gynecological Surgeries


straight and curved patterns and various lengths, these stainless steel forceps are durable and fully sterilizable for repeated use in medical environments.

ROCHESTER OCHSNER FORCEPS FEATURES

- 7.25\" (18.4cm) Overall Length.
- 1x2 "Rat-tooth" Tips for Secure Grip.
- Straight Profile.
- Horizontal Serrations Jaws on Entire Length for a Firm, Non-slip Grip..
- Boxlock Joint to Promote Smooth Closure Movements.
- Ratchet Lock Provides Self-retaining System.
- Ergonomic Ring Handles for Comfortable Handling.
- Anti Rust and Reusable.



Rochester-Ochsner forceps, also known as Kocher forceps, are robust surgical tools designed for grasping thick tissue and clamping heavy blood vessels. Commonly used in general, orthopedic, and gynecological surgeries, they prevent slippage while controlling bleeding. The forceps feature entire length horizontal serrations and a 1 x 2 tooth tip for secure gripping. With a ratchet locking mechanism, their ring handles allow hands-free operation. Available in

**WI-70000077****Rochester-Ochsner Forceps 8\"/>****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-ochsner Hemostatic Artery Forceps 8" (20.3cm) Curved 1×2 Teeth Serrated Jaws with Ratcheted Ring Handle – Premium Quality Kocher Forceps for Clamping of Arteries during General, Orthopedic, and Gynecological Surgeries

straight and curved patterns and various lengths, these stainless steel forceps are durable and fully sterilizable for repeated use in medical environments.

ROCHESTER OCHSNER FORCEPS FEATURES

- 8" (20.3cm) Overall Length.
- 1×2 "Rat-tooth" Tips for Secure Grip.
- Curved Profile.
- Horizontal Serrations Jaws on Entire Length for a Firm, Non-slip Grip..
- Boxlock Joint to Promote Smooth Closure Movements.
- Ratchet Lock Provides Self-retaining System.
- Ergonomic Ring Handles for Comfortable Handling.
- Anti Rust and Reusable.



Rochester-Ochsner forceps, also known as Kocher forceps, are robust surgical tools designed for grasping thick tissue and clamping heavy blood vessels. Commonly used in general, orthopedic, and gynecological surgeries, they prevent slippage while controlling bleeding. The forceps feature entire length horizontal serrations and a 1 × 2 tooth tip for secure gripping. With a ratchet locking mechanism, their ring handles allow hands-free operation. Available in

**WI-70000078****Rochester-Ochsner Forceps 8\"/>****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rochester-ochsner Hemostatic Artery Forceps 8" (20.3cm) Straight 1×2 Teeth Serrated Jaws with Ratcheted Ring Handle – Premium Quality Kocher Forceps for Clamping of Arteries during General, Orthopedic, and Gynecological Surgeries

straight and curved patterns and various lengths, these stainless steel forceps are durable and fully sterilizable for repeated use in medical environments.

ROCHESTER OCHSNER FORCEPS FEATURES

- 8" (20.3cm) Overall Length.
- 1×2 "Rat-tooth" Tips for Secure Grip.
- Straight Profile.
- Horizontal Serrations Jaws on Entire Length for a Firm, Non-slip Grip..
- Boxlock Joint to Promote Smooth Closure Movements.
- Ratchet Lock Provides Self-retaining System.
- Ergonomic Ring Handles for Comfortable Handling.
- Anti Rust and Reusable.



Rochester-Ochsner forceps, also known as Kocher forceps, are robust surgical tools designed for grasping thick tissue and clamping heavy blood vessels. Commonly used in general, orthopedic, and gynecological surgeries, they prevent slippage while controlling bleeding. The forceps feature entire length horizontal serrations and a 1 × 2 tooth tip for secure gripping. With a ratchet locking mechanism, their ring handles allow hands-free operation. Available in

**WI-70001305****Ruskin Kleinkert Kutz Bone Splitting Forceps
6" Straight – Double Action****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Ruskin Kleinert-Kutz Bone Splitting Forceps 6 Inch Length Straight Delicate Pattern NonLocking Plier Handle with Spring Action German Grade Stainless Steel Surgical Podiatry Bone Cutting Forceps

During multiple surgical procedures, the Ruskin (Kleinert Kutz) bone splitting forceps are used to grasp, split, and extract bony structures. The Ruskin (Kleinert Kutz) bone-splitting forceps come with two types of head orientation. The straight tip forceps are ideal for smaller and narrow bones, whereas the angled tip Ruskin forceps are good for deeper bones manipulation. In addition, the ergonomically designed comes with anti-slippage technology ensuring that

the instrument does not slip or drop during a procedure.

RUSKIN KLEINKERT KUTZ BONE SPLITTING FORCEPS FEATURES

- Ruskin Style Bone-Splitting Forceps.
- Delicate Blade Straight Jaws.
- Double Action, Spring-loaded Handle.
- Two Chisel-shaped Blades Approximating Each Other for Use as a Standard Bone Forcep.
- Exact Blade Alignment and Precision Sharpening.
- Plier-style Non Slip Handles.
- Corrosion Resistance.
- Reusable, Autoclavable Design.
- Premium German Grade Stainless Steel Construction.

PACKAGE INCLUDE:**WI-70003733****Russian Thumb Forceps 6" (15.2cm) Round Serrated Tips****\$8.99** per pieceTrade price from **\$8.27**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Russian Thumb Forcep 6" (15.2cm) Multi-Tooth Circular Serrated Tips Tissue Forceps for OB/GYN, Dental Surgical

Russian Tissue Forceps are commonly used to hold heavy or thick tissues during wound debridement and closure. Additionally, they may be used during dental procedures. The large, rounded, scoop ends have deep teeth that provide surface area and a secure grip which is useful for grasping tooth fragments and placing or removing gauze. This multi-purpose tissue forceps are commonly used to hold or grasp the tissues during gynecological procedures. The round broad teeth and a head around the rim, make it an ideal tool for performing surgeries without

complications and the serrated thumb provide a secure grip.

RUSSIAN TISSUE FORCEPS FEATURES

- Rounded Head Tip With Teeth Around the Rim
- Circular Serrated Tips Allow a Secure Grip on Hard and Slippery Tissues.
- Cupped Shaped Tips to Protect Soft Tissues from Damage
- Non-Slip Thumb for Control Grip
- Ideal for Atraumatic Tissue Grasping During Dissection

PACKAGE INCLUDE:

- 1 Pc Russian Tissue Forceps 6" (15.2cm) Rounded Serrated Tip

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel



WI-70002218

Russian Tissue Forceps 8" Round Serrated Tip

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Russian Thumb Forcep 8" (20.3cm)
Multi-Tooth Circular Serrated Tips
Premium Quality Stainless Steel Tissue
Forceps for OB/GYN, Dental Surgical
Procedures.

Russian Tissue Forceps are commonly used to hold heavy or thick tissues during wound debridement and closure. Additionally, they may be used during dental procedures. The large, rounded, scoop ends have deep teeth that provide surface area and a secure grip which is useful for grasping tooth fragments and placing or removing gauze. This multi-purpose tissue forceps are commonly used to hold or grasp the tissues during gynecological procedures. The round broad teeth and a head around the rim, make it an ideal tool for performing surgeries without

complications and the serrated thumb provide a secure grip.

RUSSIAN TISSUE FORCEPS
FEATURES

- 8" (20.3cm) Overall Length.
- Rounded Head Tip With Teeth Around the Rim.
- Circular Serrated Tips Allow a Secure Grip on Hard and Slippery Tissues.
- Cupped Shaped Tips to Protect Soft Tissues from Damage.
- Non-Slip Thumb for Control Grip.
- Ideal for Atraumatic Tissue Grasping During Dissection.
- Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE:



WI-34791613

Satinsky DeBakey Dissecting Forceps 6" (15cm)

\$17.00 per piece

Trade price from \$11.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Satinsky DeBakey Dissecting Forceps 6" (15cm) for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 15cm size shown is our standard; other sizes in the same pattern are made to order.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$17.00, falling to \$11.00 at 25 pieces — about 35% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002918

Schaedel Towel Forceps

\$30.21 per piece

Trade price from \$27.79

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Schaedel Towel Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$30.21, falling to \$27.79 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002335

Schroeder Braun Uterine Tenaculum Forceps
9.5" Straight

\$12.49 per piece

Trade price from \$11.49

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Schroeder Braun Uterine Tenaculum Forceps 9.5" Straight Economy Stainless Steel OB/GYN Surgical Instruments

The unique design ensures less traumatic application.

SCHROEDER UTERINE
TENACULUM FORCEPS FEATURES

- Long Reaching Straight Shaft Ensure Deep Surgical Access.
- Sharp Pointed Single-Toothed Jaws for Firm Grip.
- Ratchet System Allows For Self-Retaining Clamping.
- Ideal for Grasping Delicate and Heavy Tissues During Gynecology Procedures.

PACKAGE INCLUDES:

- 1 Pc SCHROEDER-BRAUN Uterine Tenaculum Forceps 9 1/2" (24cm) Straight Economy

TECHNICAL INFORMATION

Schroeder Braun Tenaculum Forceps come with finger rings & ratcheting mechanism commonly used in OB/GYN for performing tubal ligation, hysterectomies and other procedures requiring manipulating the uterus or cervix in abdominal or vaginal hysterectomy. These instruments are commonly used to hold, grasp, and move arteries and vessels during gynecology surgical procedures. Schroeder Uterine Tenaculum Forceps features two flattened overlap prongs that provide a stronghold on the tissue. These sharply pointed hooks/prongs have an inward angle of 90 Degrees.





WI-47228647

Semken Dressing Forceps Stainless Steel Surgical Instrument Controlled Tissue

\$2.20 per piece

Trade price from **\$1.00**
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Semken Dressing Forceps Stainless Steel Surgical Instrument Controlled Tissue is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.20, falling to \$1.00 at 50 pieces — about 55% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47258231

Semken Steel Dressing Forceps Surgical Instrument

\$2.50 per piece

Trade price from **\$1.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Semken Steel Dressing Forceps Surgical Instrument, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$2.50, falling to \$1.00 at 25 pieces — about 60% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003732

Semken Thump Forceps

\$8.99 per piece

Trade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Semken Thump Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$8.99, falling to \$8.27 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003371

Simpson Obstetrical Forceps 14" (35.6cm) – Fenestrated Blades

\$39.99 per piece

Trade price from **\$36.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Simpson Obstetrical Forceps Long Model 14" (35.6cm) Fenestrated Blades Stainless Steel Surgical Gynecology & Obstetrics Instruments Caesarean Section Forceps.

Simpson Obstetrical Forceps are commonly used to aid the mother in childbirth. Simpson forceps have an elongated cephalic curve. These are used when there is substantial molding, that is, temporary elongation of the fetal head as it moves through the birth canal.

The forceps feature elongated, curved cephalic blades to grasp the morphed fetal head in its migration through the birth canal. They are the most common type of forceps utilized when the use of forceps is necessary. They feature an inline handle and long, fenestrated,

smooth, teardrop-shaped cephalic curve blades. The overall length of this caesarean section obstetrical forceps is 256mm.

SIMPSON OBSTETRICAL FORCEPS FEATURES

- Elongated Cephalic Curve Obstetrical Forceps
- Fenestrated Blade for Enhanced Surgical View.
- Teardrop-Shaped Blades Accommodate the Elongation of the Birth Canal.
- Lock System Ensuring Adequate Fetal Clamping.
- Long Plier Grip Style Handle Ensuring Maximum Control.

PACKAGE INCLUDE:



WI-70003370

Simpson-Braun Obstetrical Forceps

\$127.51 per piece

Trade price from \$117.31
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Simpson-Braun Obstetrical Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$127.51, falling to \$117.31 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003369

Simpson-Luikart Obstetrical Forceps 14" Solid Blade

\$43.99 per piece

Trade price from \$40.47
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Simpson-Luikart Obstetrical Forceps 14" (36cm) Solid Blades with Grip Handle Premium Quality Stainless Steel Surgical Gynecology & Obstetrics Instruments Caesarean Section Forceps.

Simpson-Luikart Obstetrical Forceps is used to expedite delivery during prolonged labor. That forceps is also used during shoulder dystocia, fetal distress and meconium staining of the liquor. Simpson-Luikart have an elongated cephalic curved solid blades to grasp the morphed fetal head in its migration through the birth canal. They are the most common type of forceps utilized when the use of forceps is necessary. They feature an inline handle and long non-fenestrated teardrop-shaped blades. The overall

length of this caesarean section obstetrical forceps is 14inch.

SIMPSON-LUIKART OBSTETRICAL FORCEPS FEATURES

- 14 Inch (36cm) Overall Length.
- Elongated Cephalic Curve Ideal Occiput Anterior Position
- Solid Blade for a Safer Withdraw of the Baby's Head.
- Teardrop-Shaped Blades Accommodate the Elongation of the Birth Canal.
- Lock System Ensuring Adequate Fetal Clamping.
- Grooved Grip Handle Ensuring Maximum Control.
- Reusable.





WI-47222924

Singley Surgical Tissue Forceps

\$3.90 per piece

Trade price from \$1.90
Minimum order 35 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Singley Surgical Tissue Forceps is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by

hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$3.90, falling to \$1.90 at 35 pieces — about 51% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003731

Singley Thump Forceps

\$12.28 per piece

Trade price from \$11.30
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Singley Thump Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$12.28, falling to \$11.30 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47235367

Single Tissue Forceps Stainless Steel Surgical Instrument

\$4.00 per piece

Trade price from **\$1.94**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Single Tissue Forceps Stainless Steel Surgical Instrument is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each instrument is

checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.00, falling to \$1.94 at 25 pieces — about 52% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47258123

Single Tissue Grasping Forceps Reusable Surgical Instrument Secure

\$3.80 per piece

Trade price from **\$1.89**
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Single Tissue Grasping Forceps Reusable Surgical Instrument Secure is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.80, falling to \$1.89 at 50 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47016212

Skin Staple Remover Forceps Surgical Kit for Precise Staple

\$15.00 per piece

Trade price from **\$13.00**
Minimum order 12 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Skin Staple Remover Forceps Surgical Kit for Precise Staple for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$15.00, falling to \$13.00 at 12 pieces — about 13% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47002321

Skin Staple Remover Forceps-Reusable

\$14.90 per piece

Trade price from **\$12.90**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



ABOUT THIS INSTRUMENT

We manufacture Skin Staple Remover Forceps-Reusable for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.


Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$14.90, falling to \$12.90 at 25 pieces — about 13% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**EYEBROW SLANT TWEEZERS****EYEBROW SLANT TWEEZERS****FINE HAIR UPROOTING**

Designed for versatile use, grip individual hairs easily at the root, making them a must-have tool for precise brow definition and effortless hair removal.

**WI-70000144****Slanted Eyebrow Tweezers 4" Angled Tips****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Slanted Eyebrow Tweezers 4" (10.2cm) Classic, Angled Slant Tip Stainless Steel Beauty Tool Designed for Shaping and Grooming Eyebrows.

These eyebrow tweezers are versatile tools for precise eyebrow shaping and hair removal, featuring an angled tip for optimal control. Slanted Eyebrow Tweezers designed to pull hair in the direction of growth, minimizing pain and ensuring effective root-level removal. Made of premium stainless steel, featuring sturdy, sharp, and durable tips that are easy to clean and reuse, ensuring lasting precision.

SLANT TIP EYEBROW TWEEZERS FEATURES

- Size: 4 Inches (10.2cm).
- 25 Degree Angled Classic Slant Tip.

- Perfectly Aligned Tips to Close Evenly and Firmly to Grab Hair.
- Ergonomic Handle for Better Handling.
- Designed for Shaping and Grooming Eyebrows.
- Suitable for Both Thick and Fine, Stubborn, or Ingrown Hairs.
- Corrosion Resistance.
- Reusable.
- Made with High-quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Slanted Eyebrow Tweezers 4" Angled Tips

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Beauty & Personal Care Instruments**WI-70000818****Splinter Forceps 3.5" Fine Point Serrated Tips****\$8.99** per pieceTrade price from **\$8.27**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Splinter Forceps 3.5" (8.9cm) Fine Point Serrated Straight Tips Premium Quality Stainless Steel Slip-Resistant Grip Medical Surgical Instrument.

Splinter forceps are tweezer like surgical instrument use for removing splinters and debris on the skin. Splinter forceps are 3.5 inches long and a great choice for first aid boxes. With fine-pointed tips and a precision gripping mechanism, splinter forceps excel in delicate tasks, making them perfect for DIY enthusiasts in various crafting projects. Splinter forceps features tapered precision tips can be used to remove ingrown hairs and particles embedded under the skin including splinters, thorns and glass. Amazingly precise hand-finished points.

The durable stainless steel construction ensures longevity. The finest finishing design adds a touch of aesthetics to these tweezers, while the non-slip handle provides a secure grip for precise maneuvers. Easy to clean and can be reused after sterilization.

SPLINTER FORCEPS FEATURES

- Crafted with 3.5" (8.9cm) Overall Length.
- Fine Point Straight Tips.
- Serrated Jaws for Secure Holding.
- Long-Reaching Jaws for Deep Surgical Access.
- Ergonomic Spring Handle Ensuring Superior Control.
- Ideal for Delicately Surgical Procedure.
- Lightweight and Reusable.

**WI-70000814****Splinter Forceps 4.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Splinter Forceps 4.5\"/>

Splinter forceps are tweezer like surgical instrument use for removing splinters and debris on the skin. Splinter forceps are 4.5 inches long and a great choice for first aid boxes. With fine-pointed tips and a precision gripping mechanism, splinter forceps excel in delicate tasks, making them perfect for DIY enthusiasts in various crafting projects. Splinter forceps features tapered precision tips can be used to remove ingrown hairs and particles embedded under the skin including splinters, thorns and glass. Amazingly precise hand-finished points.

The durable stainless steel construction ensures longevity. The finest finishing design adds a touch of aesthetics to these tweezers, while the non-slip handle provides a secure grip for precise maneuvers. Available in straight and curved patterns. Easy to clean and can be reused after sterilization.

SPLINTER FORCEPS FEATURES

- Crafted with 4.5\"/>
- Fine Point Curved Tips.
- Serrated Jaws for Secure Holding.
- Long-Reaching Jaws for Deep Surgical Access.
- Ergonomic Spring Handle Ensuring Superior Control.
- Ideal for Delicate Surgical Procedure.

**WI-70000819****Splinter Forceps 4.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Splinter Forceps 4.5\"/>


Splinter forceps are tweezer like surgical instrument use for removing splinters and debris on the skin. Splinter forceps are 4.5 inches long and a great choice for first aid boxes. With fine-pointed tips and a precision gripping mechanism, splinter forceps excel in delicate tasks, making them perfect for DIY enthusiasts in various crafting projects. Splinter forceps features tapered precision tips can be used to remove ingrown hairs and particles embedded under the skin including splinters, thorns and glass. Amazingly precise hand-finished points.

The durable stainless steel construction ensures longevity. The finest finishing design adds a touch of aesthetics to these tweezers, while the non-slip handle provides a secure grip for precise maneuvers. Easy to clean and can be reused after sterilization.

SPLINTER FORCEPS FEATURES

- Crafted with 4.5\"/>
- Fine Point Straight Tips.
- Serrated Jaws for Secure Holding.
- Long-Reaching Jaws for Deep Surgical Access.
- Ergonomic Spring Handle Ensuring Superior Control.
- Ideal for Delicate Surgical Procedure.
- Lightweight and Reusable.



**WI-70000813**

Splinter Forceps 5" Fine Point Curved Serrated Tips

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Splinter Forceps 5" (12.7cm) Fine Point Curved Serrated Tips Premium Quality Stainless Steel Slip-Resistant Grip Medical Surgical Instrument.

Splinter forceps are tweezer like surgical instrument use for removing splinters and debris on the skin. Splinter forceps are 4.5 inches long and a great choice for first aid boxes. With fine-pointed tips and a precision gripping mechanism, splinter forceps excel in delicate tasks, making them perfect for DIY enthusiasts in various crafting projects. Splinter forceps features tapered precision tips can be used to remove ingrown hairs and particles embedded under the skin including splinters, thorns and glass. Amazingly precise hand-finished points.

The durable stainless steel construction ensures longevity. The finest finishing design adds a touch of aesthetics to these tweezers, while the non-slip handle provides a secure grip for precise maneuvers. Available in straight and curved patterns. Easy to clean and can be reused after sterilization.

SPLINTER FORCEPS FEATURES

- Crafted with 5" (12.7cm) Overall Length.
- Fine Point Curved Tips.
- Serrated Jaws for Secure Holding.
- Long-Reaching Jaws for Deep Surgical Access.
- Ergonomic Spring Handle Ensuring Superior Control.
- Ideal for Delicate Surgical Procedure.

**WI-70000816**

Splinter Forceps 5" Fine Point Straight Serrated Tips

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Splinter Forceps 5" (12.7cm) Fine Point Serrated Straight Tips Premium Quality Stainless Steel Slip-Resistant Grip Medical Surgical Instrument.

Splinter forceps are tweezer like surgical instrument use for removing splinters and debris on the skin. Splinter forceps are 5 inches long and a great choice for first aid boxes. With fine-pointed tips and a precision gripping mechanism, splinter forceps excel in delicate tasks, making them perfect for DIY enthusiasts in various crafting projects. Splinter forceps features tapered precision tips can be used to remove ingrown hairs and particles embedded under the skin including splinters, thorns and glass. Amazingly precise hand-finished points.

The durable stainless steel construction ensures longevity. The finest finishing design adds a touch of aesthetics to these tweezers, while the non-slip handle provides a secure grip for precise maneuvers. Easy to clean and can be reused after sterilization.

SPLINTER FORCEPS FEATURES

- Crafted with 5" (12.7cm) Overall Length.
- Fine Point Straight Tips.
- Serrated Jaws for Secure Holding.
- Long-Reaching Jaws for Deep Surgical Access.
- Ergonomic Spring Handle Ensuring Superior Control.
- Ideal for Delicate Surgical Procedure.
- Lightweight and Reusable.



WI-28223336

Sponge Forceps 12 Inch Straight

\$8.00 per piece

Trade price from \$4.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Sponge Forceps 12 Inch Straight for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a straight working end for direct, in-line access.

The 12" size shown is our standard; other sizes in the same pattern are

made to order. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28180345

Sponge Forceps 9.5" Straight Serrated Surgical

\$8.00 per piece

Trade price from \$4.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Sponge Forceps 9.5" Straight Serrated Surgical is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: a straight working end for direct, in-line access, and serrated jaws for grip that does not slip under tension.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$8.00, falling to \$4.00 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002228

Sponge Forceps Blue Coated 9.5 inch Straight Serrated Jaws

\$14.99 per piece

Trade price from \$13.79
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Blue Coated Foerster Sponge Holding Forceps Straight 9.5 inch (24.13cm) Serrated Jaws Medical Grade Stainless Steel Gynecology Electrosurgical instruments

Sponge forceps or sponge holding forceps also known as Foerster Sponge Forceps used by doctors to pick up swabs, cotton or sponges during gynecology surgical procedures. The ratchet handle helps to hold the sponges firmly. These forceps are favorites by doctors because of their easy control. Foerster Sponge Dressing & Holding Forceps are essential in many different types of procedures. The ratchet mechanism towards the handles and have fenestrated looped jaws can have either smooth or serrated inner surfaces to prevent local

damage. Also comes in straight and curved profiles according to the surgical approach.

Sponge Forceps Blue Coated Features

Sponge Forceps Straight 9.5 inch: Sponge Forceps feature as blue coated Oval tip scissor like shape straight serrated jaws, 3-step locking ratchet for better control with holding rounded rings. It helps to hold the forceps firmly and place it in the right position. They can be used to manipulate the tissue while the doctor performs the procedure.

Stainless Steel Construction: The Stainless Steel construction ensures maximum durability and dependability along with blue non-conductive coating for added safety.



WI-44567176

Sponge Forceps Stainless Steel Medical Surgical Dressing Forceps

\$2.90 per piece

Trade price from \$0.90
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

About This Instrument

We manufacture Sponge Forceps Stainless Steel Medical Surgical Dressing Forceps for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.90, falling to \$0.90 at 10 pieces — about 69% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25867723

Stainless Forceps

\$5.00 per piece

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Forceps for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each instrument is

checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$5.00. Minimum order 1 piece, and the price moves with quantity.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-30847995

Stainless Steel 6.4" Babcock Tissue Grasping Forceps Power

\$3.00 per piece

Trade price from **\$1.90**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel 6.4" Babcock Tissue Grasping Forceps Power, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.00, falling to \$1.90 at 1 piece — about 37% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Forging, milling, heat treatment, hand finishing and passivation all happen

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36329018

**Stainless Steel 8 Inch Babcock Forceps
Atraumatic Tissue**

\$2.30 per piece

Trade price from **\$1.50**
Minimum order 8 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel 8 Inch Babcock Forceps Atraumatic Tissue, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 8" size shown is our standard; other sizes in the same pattern are made to order.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.30, falling to \$1.50 at 8 pieces — about 35% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34445744

**Stainless Steel Artery Forceps Straight Tip
Surgical Instrument Set**

\$5.50 per piece

Trade price from **\$2.00**
Minimum order 3 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Artery Forceps Straight Tip Surgical Instrument Set is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a straight working end for direct, in-line access.

Supplied as a set — the contents can be adjusted to your list. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$5.50, falling to \$2.00 at 3 pieces — about 64% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36115438

Stainless Steel Bayonet Dressing Forceps Surgical Tissue Forceps

\$5.00 per piece

Trade price from **\$2.00**
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Bayonet Dressing Forceps Surgical Tissue Forceps — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$5.00, falling to \$2.00 at 100 pieces — about 60% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36176181

Stainless Steel Brown Tissue Forceps Surgical Instrument Medical Forceps

\$3.00 per piece

Trade price from **\$0.90**
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Brown Tissue Forceps Surgical Instrument Medical Forceps, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.00, falling to \$0.90 at 5 pieces — about 70% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-39628463

Stainless Steel Cheatle Forceps (10)

\$3.00 per piece

Trade price from \$1.90
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Cheatle Forceps (10) is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$3.00, falling to \$1.90 at 25 pieces — about 37% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34844439

Stainless Steel Debakey

\$17.00 per piece

Trade price from \$11.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Debakey is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$17.00, falling to \$11.00 at 25 pieces — about 35% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34795638

Stainless Steel Debakey Forceps Instrument

\$17.00 per piece

Trade price from \$11.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Debakey Forceps Instrument is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$17.00, falling to \$11.00 at 25 pieces — about 35% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35031400

Stainless Steel Debakey Forceps Reusable Surgical

\$3.27 per piece

Trade price from \$2.93
Minimum order 4 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Debakey Forceps Reusable Surgical is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.27, falling to \$2.93 at 4 pieces — about 10% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-24100124

Stainless Steel Entomology Dissection Kit Teaser Needles Adson

\$30.00 per piece

Trade price from **\$19.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Entomology Dissection Kit Teaser Needles Adson is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$30.00, falling to \$19.00 at 1 set — about 37% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34972290

Stainless Steel Foerster Sponge Forceps 14cm

\$3.50 per piece

Trade price from **\$2.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Foerster Sponge Forceps 14cm for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 14cm size shown is our standard; other sizes in the same pattern are made to order.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$3.50, falling to \$2.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34904583

Stainless Steel Foerster Sponge Holding Forceps 10.5

\$3.50 per piece

Trade price from **\$2.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Foerster Sponge Holding Forceps 10.5, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$3.50, falling to \$2.00 at 10 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34963423

Stainless Steel Foldable Portable Surgical Cheatle Forceps

\$3.00 per piece

Trade price from **\$1.90**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Foldable Portable Surgical Cheatle Forceps is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.00, falling to \$1.90 at 10 pieces — about 37% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-33270313

Stainless Steel Hair Transplant Forceps Surgical

\$12.50 per piece

Trade price from **\$10.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Hair Transplant Forceps Surgical is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$12.50, falling to \$10.00 at 1 piece — about 20% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44547556

Stainless Steel Kocher Tissue Forceps Surgical Instrument Set

\$2.30 per piece

Trade price from **\$0.80**
Minimum order 8 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Kocher Tissue Forceps Surgical Instrument Set is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.30, falling to \$0.80 at 8 pieces — about 65% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34952384

Stainless Steel Magill Forceps for Medical Intubation

\$6.00 per piece

Trade price from **\$4.60**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Magill Forceps for Medical Intubation, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$6.00, falling to \$4.60 at 10 pieces — about 23% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25991828

Stainless Steel Pediatric Mixer Dissecting Forceps

\$8.00 per piece

Trade price from **\$4.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Pediatric Mixer Dissecting Forceps is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$8.00, falling to \$4.00 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46980504

Stainless Steel Skin Staple Removal Forceps

\$14.80 per piece

Trade price from \$9.00
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Skin Staple Removal Forceps is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$14.80, falling to \$9.00 at 50 pieces — about 39% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-26095341

Stainless Steel Sponge Dressing Forceps Curved Slims Forceps

\$3.00 per piece

Trade price from \$1.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Sponge Dressing Forceps Curved Slims Forceps, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a curved working end for reach around structures.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$3.00, falling to \$1.50 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34936434

Stainless Steel Sponge Holding Forceps

\$3.50 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Sponge Holding Forceps is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.50, falling to \$2.00 at 10 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-10247516

Stainless Steel Surgical Dressing Forceps

\$2.99 per piece

Trade price from \$1.99
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Dressing Forceps is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.99, falling to \$1.99 at 10 pieces — about 33% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27187853

Stainless Steel Surgical Forceps Sponge Holding

\$3.00 per piece

Trade price from **\$1.50**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Forceps Sponge Holding is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.00, falling to \$1.50 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27191808

Stainless Steel Surgical Instruments Set Foerster Sponge

\$3.00 per piece

Trade price from **\$1.50**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Instruments Set Foerster Sponge is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$3.00, falling to \$1.50 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-30909034

Stainless Steel TC Babcock Forceps for Atraumatic

\$3.00 per piece

Trade price from **\$1.90**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel TC Babcock Forceps for Atraumatic is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: tungsten carbide inserts for a longer cutting life.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$3.00, falling to \$1.90 at 1 piece — about 37% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-43480659

Stainless Steel Utrata Capsulorrhesis Forceps Eye Surgery

\$12.00 per piece

Trade price from **\$9.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Utrata Capsulorrhesis Forceps Eye Surgery is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$12.00, falling to \$9.00 at 10 pieces — about 25% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

We work the instrument by hand at the finishing stage — box locks, hinges and

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000552

Standard Articulating Paper Forcep 6" Cross Serrated Tips

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Articulating Paper Forcep 6" (15.2cm) Cross & Longitudinal Serrated Straight Jaws Premium Quality Stainless Steel Surgical Dental Instruments



Articulating paper forceps is an essential dentist tool for handling articulating paper. Articulating paper forceps is a reverse-action tweezer like plier used to hold and place articulating paper to check occlusal surfaces of teeth before and after restorative procedures. Facilitates the placement paper in the area of posterior teeth and helps with the detection of points of contact between opposing teeth, making it possible to discover those that need to be adjusted and corrected. The articulating paper plier is suitable for dentists, art & craftwork, and other

laboratory works. You can reuse standard articulating paper forceps over time as they will not bend or tarnish on you. High quality stainless steel construction meet all international quality standards.

STANDARD ARTICULATING PAPER FORCEPS FEATURES

- 6" (15.2cm) Overall Length
- Cross & Longitudinal Serration Jaws
- Serrated Blades to Ensure Strong Grip
- Reverse-action Opening & Closing Mechanism for Firm Grip
- Long-lasting, Rust Free
- Reusable
- Made With Premium Quality Stainless Steel

PACKAGE INCLUDE:



WI-30931001

Steel Babcock Tissue Intestinal Holding Forceps Surgical

\$3.00 per piece

Trade price from \$1.90
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Steel Babcock Tissue Intestinal Holding Forceps Surgical is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$3.00, falling to \$1.90 at 1 piece — about 37% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47223591

Steel Magill Catheter Forceps Surgical Instrument for Precise Airway

\$8.00 per piece

Trade price from \$4.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Steel Magill Catheter Forceps Surgical Instrument for Precise Airway, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$8.00, falling to \$4.00 at 10 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-30785999

Steel Surgical Sponge Holding Dressing Forceps

\$3.00 per piece

Trade price from \$1.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Steel Surgical Sponge Holding Dressing Forceps — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.00, falling to \$1.50 at 1 piece — about 50% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003730

Stille Thumb Forceps

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stille Thumb Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$8.99, falling to \$8.27 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003729

Stille Thump Forceps

\$9.90 per piece

Trade price from \$9.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stille Thump Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.90, falling to \$9.11 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003583

Straight ESD 12 Fine Pointed Eyelash Tweezer

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Straight ESD 12 Fine Pointed Eyelash Tweezer is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003587

Straight Polish Silver Straight Eyelash Tweezers Custom Eyelash Extension Straight Tweezer

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Straight Polish Silver Straight Eyelash Tweezers Custom Eyelash Extension Straight Tweezer is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-44565366

Straight Sponge Holding Forceps Stainless Steel Surgical

\$2.90 per piece

Trade price from **\$0.90**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Straight Sponge Holding Forceps Stainless Steel Surgical is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a straight working end for direct, in-line access.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.90, falling to \$0.90 at 10 pieces — about 69% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-08601496

Straight Surgical Instrument Dressing Forceps

\$3.00 per piece

Trade price from **\$1.50**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Straight Surgical Instrument Dressing Forceps for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a straight working end for direct, in-line access.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$3.00, falling to \$1.50 at 10 pieces — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36183192

Straight Surgical Tissue Forceps Brown Pattern

\$3.00 per piece

Trade price from **\$0.90**

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Straight Surgical Tissue Forceps Brown Pattern for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a straight working end for direct, in-line access.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.00, falling to \$0.90 at 10 pieces — about 70% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003586

Straight Tweezers for Eyelash Extension

\$9.99 per piece

Trade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Fine Point Black Coated Straight Tweezers for Eyelash Extension specially designed for professional precision task in grafting eyelash, designed to meet the requirement of the eyelash artists.

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Eyelash Tweezers – Beauty & Personal Care
APPLICATION: Eyelash Extension, Eyebrow, Beauty, Makeup, Salon
COLOR : Black Coated
TIP: Fine Pointed Tip

RUST PREVENTION PROCEDURE:

Passivated
ULTRASONIC CLEANED: Yes
TESTS PERFORMED: Boil Test, Performance Test, Shape Test
QC PASSED: Yes
PACKING DETAILS: Individually Packed

Explore more **Eyelash Extension Tweezers**

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.



WI-70003585

Strong Curved Volume Eye Care Curved Tweezers for Volume Eyelash Extension

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Strong Curved Volume Eye Care Curved Tweezers for Volume Eyelash Extension is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-30866109

Surgical Babcock Tissue Forceps Straight Size

\$3.00 per piece

Trade price from \$1.90
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Surgical Babcock Tissue Forceps Straight Size, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: a straight working end for direct, in-line access.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$3.00, falling to \$1.90 at 1 piece — about 37% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27282717

Surgical Dressing Tissue Forceps

\$2.39 per piece

Trade price from \$2.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Surgical Dressing Tissue Forceps is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.39, falling to \$2.00 at 10 pieces — about 16% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-33339752

Surgical Kelly Locking Forceps Stainless Steel Reusable Different

\$3.00 per piece

Trade price from \$1.90
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Kelly Locking Forceps Stainless Steel Reusable Different is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$3.00, falling to \$1.90 at 25 pieces — about 37% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-30819959

Surgical Stainless Steel Gross Forceps 20cm

\$3.00 per piece

Trade price from **\$1.50**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Stainless Steel Gross Forceps 20cm is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The 20cm size shown is our standard; other sizes in the same pattern are made to order.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.00, falling to \$1.50 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-26066344

Surgical Use Reusable Steel Sponge Holding Forceps

\$3.00 per piece

Trade price from **\$1.50**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Use Reusable Steel Sponge Holding Forceps is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.00, falling to \$1.50 at 1 piece — about 50% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001477****T/C Adson Brown Forceps 4.75\"/>****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Adson Brown Tissue Forceps 9×9 Teeth Plastic Surgery Instrument 4.75\"/>

The Adson Brown Tissue Forceps are used for holding and manipulating delicate tissues and light materials during suturing or other surgical procedures. These tissue forceps have a wide, flat thumb grasp area with horizontal ridges for a better grip. The jaws are short and the tips are narrow and they come in 4-3/4\"/>

material can be autoclaved and sterilized easily for reuse purposes.

T/C ADSON BROWN TISSUE FORCEPS FEATURES**ADSON BROWN TISSUE FORCEPS 9X9**

TEETH: The TC Adson Brown Tissue Forceps are the ideal addition to any hospital or clinical environment. They feature a gold TC handle with a wide flat thumb grasp area and are designed to allow the user to expertly hold the tissue in place during delicate procedures. The premium quality combined with affordable price sets the value high for these T/C Adson Brown Tissue Forceps with straight tip having 9×9 Teeth.

O.R. GRADE T/C AND STAINLESS STEEL: Constructed of O.R. Grade high quality stainless steel and tungsten

**WI-70001480****T/C Adson Brown Tissue Forceps 7×7 Teeth 4.75\"/>****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Adson Brown Tissue Forceps 7×7 Teeth Plastic Surgery Instrument 4.75\"/>

The Adson Brown Tissue Forceps are used for holding and manipulating delicate tissues and light materials during suturing or other surgical procedures. These tissue forceps have a wide, flat thumb grasp area with horizontal ridges for a better grip. The jaws are short and the tips are narrow and they come in 4-3/4\"/>

material can be autoclaved and sterilized easily for reuse purposes.

T/C ADSON BROWN TISSUE FORCEPS FEATURES**ADSON BROWN TISSUE FORCEPS 7X7**

TEETH: The TC Adson Brown Tissue Forceps are the ideal addition to any hospital or clinical environment. They feature a gold TC handle with a wide flat thumb grasp area and are designed to allow the user to expertly hold the tissue in place during delicate procedures. The premium quality combined with affordable price sets the value high for these T/C Adson Brown Tissue Forceps with straight tip having 7×7 Teeth.

O.R. GRADE T/C AND STAINLESS STEEL: Constructed of O.R. Grade high quality stainless steel and tungsten



**WI-70001479****T/C Adson Dressing Forceps 4.75" Serrated Gold Handle****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Adson Dressing Forceps 4.75 Inch (12cm) Straight Serrated Tips with Tungsten Carbide Inserts Plastic Surgery Instruments with Gold Handle.

Adson dressing forceps tweezer used for grasping, holding and mobilizing different types of dressings, in order to absorb blood and other fluids. Its main role during surgery is to pack wounds and surgical incisions. This dressing forcep comes with a pair of tapering jaws for grasping smaller dressings with tungsten carbide inserts for an impeccable grip experience. The gold handle has flat stems with horizontal ridges that ensure a secure grip.

T/C ADSON DRESSING FORCEPS FEATURES

TC ADSON DRESSING FORCEPS: The TC Adson Dressing Forceps are the ideal addition to any hospital or clinical environment. They feature a fine serrated point, gold TC handle and are designed to allow the user to expertly move and hold tissue/dressings in place during delicate procedures. The premium quality combined with affordable price sets the value high for these widely used T/C Adson Dressing Serrated Forceps.

O.R. GRADE T/C AND STAINLESS STEEL: Constructed of O.R. Grade high quality stainless steel and tungsten carbide for the highest degree of durability and product lifespan. The Adson Dressing Forceps won't rust, tarnish break or bend on you thanks to corrosion resistant properties. Tungsten Carbide is more dense than

**WI-70001478****T/C Adson Dressing Forceps 6" Serrated Gold Handle****\$13.99** per pieceTrade price from **\$12.87**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Adson Dressing Forceps 6 Inch (15.24cm) Straight Serrated Tips with Tungsten Carbide Inserts Plastic Surgery Instruments with Gold Handle.

Adson dressing forceps tweezer used for grasping, holding and mobilizing different types of dressings, in order to absorb blood and other fluids. Its main role during surgery is to pack wounds and surgical incisions. This dressing forcep comes with a pair of tapering jaws for grasping smaller dressings with tungsten carbide inserts for an impeccable grip experience. The gold handle has flat stems with horizontal ridges that ensure a secure grip.

T/C ADSON DRESSING FORCEPS FEATURES

TC ADSON DRESSING FORCEPS: The TC Adson Dressing Forceps are the ideal addition to any hospital or clinical environment. They feature a fine serrated point, gold TC handle and are designed to allow the user to expertly move and hold tissue/dressings in place during delicate procedures. The premium quality combined with affordable price sets the value high for these widely used T/C Adson Dressing Serrated Forceps.

O.R. GRADE T/C AND STAINLESS STEEL: Constructed of O.R. Grade high quality stainless steel and tungsten carbide for the highest degree of durability and product lifespan. The Adson Dressing Forceps won't rust, tarnish break or bend on you thanks to corrosion resistant properties. Tungsten Carbide is more dense than

**WI-70003743****T/C Adson Micro Tissue Forceps 1x2 Teeth 4.75"****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Adson Tissue Micro Kocher Forceps 1X2 Teeth Plastic Surgery Forceps 4.75" (12cm) Straight Fine Point with Tungsten Carbide Inserts Surgical Veterinary Instruments with Gold Handle.

Adson tissue forceps are specialized instrument used to grasp and stabilize soft tissue. These toothed Adson forceps used on tissues that require a firm grip such as fascia, skin or the abdominal wall. 1x2 Teeth Adson Micro Tissue Forceps are helpful in across a range of surgical procedures including biopsy, dermatology, minor surgery, vasectomy, and carpal tunnel. Available in a variety of forms with a cross-serrated 1x2 teeth and straight or angled tips. These forceps made with

german grade stainless steel for lasting performance.

T/C ADSON MICRO TISSUE FORCEPS FEATURES

ADSON TISSUE FORCEPS 1X2 TEETH: The TC Micro Adson Kocher Forceps are the ideal addition to any hospital or clinical environment. They feature grip thumb grasp area which designed to allow the user to expertly hold tissue in place during delicate procedures. The premium quality combined with affordable price sets the value high for these Adson Kocher Forceps 4.75". Perfect for bending, separating and moving small items.

O.R. GRADE T/C AND STAINLESS STEEL: Constructed of O.R. Grade high quality stainless steel for the highest degree of durability and product

**WI-70001476****T/C Adson Tissue Forceps 1x2 Teeth 4.75" Gold Handle****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Adson Tissue Kocher Forceps 1X2 Teeth Plastic Surgery Forceps 4.75" (12cm) Straight Fine Point with Tungsten Carbide Inserts Surgical Veterinary Instruments with Gold Handle

Adson tissue forceps are specialized instrument used to grasp and stabilize soft tissue. These toothed Adson forceps used on tissues that require a firm grip such as fascia, skin or the abdominal wall. The Adson Kocher Forceps are helpful in across a range of surgical procedures including biopsy, dermatology, minor surgery, vasectomy, and carpal tunnel. Available in a variety of forms with a cross-serrated 1x2 teeth and straight or angled tips. These forceps made with

german grade stainless steel for lasting performance.

T/C ADSON TISSUE KOCHER FORCEPS FEATURES

ADSON TISSUE FORCEPS 1X2 TEETH: The TC Micro Adson Kocher Forceps are the ideal addition to any hospital or clinical environment. They feature a gold handle with a high grip thumb grasp area which designed to allow the user to expertly hold tissue in place during delicate procedures. The premium quality combined with affordable price sets the value high for these T/C Micro Adson Kocher Forceps 4.75". Perfect for bending, separating and moving small items.

O.R. GRADE T/C AND STAINLESS STEEL: Constructed of O.R. Grade high quality stainless steel and tungsten



**WI-70001475****T/C Adson Tissue Forceps 1x2 Teeth 6" Gold Handle****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Adson Tissue Kocher Forceps 1X2 Teeth Plastic Surgery Forceps 6" (15.24cm) Straight Fine Point with Tungsten Carbide Inserts Surgical Veterinary Instruments with Gold Handle

Adson tissue forceps are specialized instrument used to grasp and stabilize soft tissue. These toothed Adson forceps used on tissues that require a firm grip such as fascia, skin or the abdominal wall. The Adson Kocher Forceps are helpful in across a range of surgical procedures including biopsy, dermatology, minor surgery, vasectomy, and carpal tunnel. Available in a variety of forms with a cross-serrated 1x2 teeth and straight or angled tips. These forceps made with

german grade stainless steel for lasting performance.

T/C ADSON TISSUE KOCHER FORCEPS FEATURES**ADSON TISSUE FORCEPS 1X2 TEETH:**

The TC Micro Adson Kocher Forceps are the ideal addition to any hospital or clinical environment. They feature a gold handle with a high grip thumb grasp area which designed to allow the user to expertly hold tissue in place during delicate procedures. The premium quality combined with affordable price sets the value high for these T/C Micro Adson Kocher Forceps 6". Perfect for bending, separating and moving small items.

O.R. GRADE T/C AND STAINLESS STEEL: Constructed of O.R. Grade high quality stainless steel and tungsten

**WI-70003314****T/C Extracting Forceps #859 for Lower Root & Molars, Delicate****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

English Pattern T/C Extracting Forceps #859 Lower Root & Molars, Delicate – Tungsten Carbide Inserts Serrated Jaws with Gold Plated Angled Grip Handle Premium Quality Stainless Steel Surgical Dental Instruments.

Extraction Forceps No. 859 designed with narrow, precise beaks meet at a point ensuring a firm grip. Durable and sturdy, our extraction forceps are the perfect tool to remove teeth from the alveolar bone. English pattern forceps are designed so that the fulcrum is closer to the working end, providing greater leverage, easing extraction. The high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

T/C ENGLISH EXTRACTING FORCEPS #859 FEATURES

- Right Angled, Serrated Tips.
- Delicate Jaws.
- Tungsten Carbide Inserts.
- For Lower Root Tips & Molars.
- English Pattern.
- Gold Plated Angled Grip Handle.
- Made with High Quality German Stainless Steel

PACKAGE INCLUDE

- 1 Pc T/C English Pattern Extraction Forceps #859 – Serrated Jaws – Gold Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Dental

Extracting Forceps





WI-70002212

T/C Micro Suture Tying Forceps 5.5" with Platform Gold Plated Handle

\$12.99 per piece

Trade price from \$11.95

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


T/C Micro Suture Tying Forceps 5.5" (14cm) Tungsten Carbide Inserts Tips with Tying Platform Straight Round Grip Gold Plated Handle Premium Quality Stainless Steel Microsurgery Surgical Instruments .

Micro suture tying forceps are used to hold, manipulate and tie delicate suture threads during microsurgery, ophthalmology, plastic surgery, and neurovascular procedures. Suture tying forceps comes with built-in platform near the tips for secure, precise and controlled tying of fine-gauge sutures. Featuring atraumatic tips, ergonomic gold plated ribbed-grip handle for comfortable control. The tungsten carbide tips ensure durability, making them 10x stronger than steel. Microsurgery suture tying forceps are

made from premium stainless steel, provide a reliable, easy, and gentle grip for suture work. They are user-friendly, easy to clean, and fully autoclavable, ensuring years of reliable use.

T/C MICRO SUTURE TYING FORCEPS FEATURES

- Length: 5.5 Inches (14 cm).
- Straight Pattern.
- Tying Platform for Secure, Precise and Controlled Tying of Fine Sutures.
- Tungsten Carbide Inserts Tips.
- Round Grip Gold Plated Handle for Superior Ergonomic Control.
- Suitable for Ophthalmology, Cardiovascular, Neurosurgery, and Plastic Surgery.
- Fully Autoclavable/Reusable.



WI-34992064

TC Potts Smith Forceps Gold Handle Stainless Steel Surgical

\$3.00 per piece

Trade price from \$1.00

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

TC Potts Smith Forceps Gold Handle Stainless Steel Surgical is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. The name tells you what you are getting: tungsten carbide inserts for a longer cutting life.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$3.00, falling to \$1.00 at 10 pieces — about 67% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35015094

TC Potts Smith Forceps Straight Tungsten Carbide Tips

\$3.00 per piece

Trade price from **\$1.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture TC Potts Smith Forceps Straight Tungsten Carbide Tips for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. From the pattern: a straight working end for direct, in-line access, and tungsten carbide inserts for a longer cutting life.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.00, falling to \$1.00 at 10 pieces — about 67% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47263071

Tenaculum Pozzi Forceps Gynecological Surgical Instrument for Precise Cervical

\$5.00 per piece

Trade price from **\$2.50**
Minimum order 14 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Tenaculum Pozzi Forceps Gynecological Surgical Instrument for Precise Cervical, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$5.00, falling to \$2.50 at 14 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47243244

Tenaculum Pozzi Surgical Forceps Precision Gynecological Instrument Cervical

\$4.88 per piece

Trade price from **\$2.95**
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Tenaculum Pozzi Surgical Forceps Precision Gynecological Instrument Cervical is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.88, falling to \$2.95 at 50 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000694

Thomas-Gaylor Uterine Biopsy Punch Forceps 8.5" Curved Jaws 5mm Bite

\$24.99 per piece

Trade price from **\$22.99**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thomas-Gaylor Uterine Biopsy Punch Forceps 8.5" (21.6cm) Angled Shank, 5mm Bite Curved Interfitted Cups Jaws with Ratchet Lock Finger Ring Handle Premium Quality Stainless Steel Obstetrics & Gynecology Surgical Instruments.

Thoms-Gaylor uterine biopsy forceps is a specialized OB/GYN instrument used by surgeon during uterine biopsy surgical procedures. It used to take precise biopsies from the difficult to reach anatomical areas including the cervix during cervical biopsy. Thoms Gaylor biopsy punch exclusively used for excising, grasping and obtaining uterine and cervix tissue samples. This biopsy forceps features a circular punch tip that is ideal to excise the tissue sample with precision. The upper

jaw has a pin that efficiently holds the specimen in place during removal. The long-reaching angled shafts designed for easy accessing the cervix that connect with boxlock to the finger rings for a stable grip. The boxlock permits to close the interfitted cup tips against each other with a smooth movement. Made with German Grade Stainless Steel and can be sterilized for repeat use.

THOMAS GAYLOR BIOPSY PUNCH FORCEPS FEATURES

- Crafted With 8.5" (21.6cm) Overall Length.
- 5mm Bite Interfitted Circular Cups.
- Curved Jaws.
- Long-reaching Angled Shafts for Accessing the Cervix.



**WI-70000693****Thomas-Gaylor Uterine Biopsy Punch Forceps
9.5" Curved Jaws 5mm Bite****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thomas-Gaylor Uterine Biopsy Punch Forceps 9.5" (24.1cm) Angled Shank, 5mm Bite Curved Interfitted Cups Jaws with Ratchet Lock Finger Ring Handle Premium Quality Stainless Steel Obstetrics & Gynecology Surgical Instruments.

Thoms-Gaylor uterine biopsy forceps is a specialized OB/GYN instrument used by surgeon during uterine biopsy surgical procedures. It used to take precise biopsies from the difficult to reach anatomical areas including the cervix during cervical biopsy. Thoms Gaylor biopsy punch exclusively used for excising, grasping and obtaining uterine and cervix tissue samples. This biopsy forceps features a circular punch tip that is ideal to excise the tissue sample with precision. The upper

jaw has a pin that efficiently holds the specimen in place during removal. The long-reaching angled shafts designed for easy accessing the cervix that connect with boxlock to the finger rings for a stable grip. The boxlock permits to close the interfitted cup tips against each other with a smooth movement. Made with German Grade Stainless Steel and can be sterilized for repeat use.

**THOMAS GAYLOR BIOPSY
PUNCH FORCEPS FEATURES**

- Crafted With 9.5" (24.1cm) Overall Length.
- 5mm Bite Interfitted Circular Cups.
- Curved Jaws.
- Long-reaching Angled Shafts for Accessing the Cervix.

**WI-70000335****Three Prong Sterilizer Forceps 11"****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Three Prong Sterilizer Forceps 11" (27.9cm) Premium Quality Stainless Steel Surgical Medical Sterilizing Aid.


3 Prong Sterilizer Forceps is essential surgical tools for safely handling sterile items, primarily used to retrieve sterilized instruments from autoclaves or hot solutions without contamination. These utility sterilization forceps featuring feature three prongs that curve inward, resembling a claw, providing a secure grip on specific items. Their three prong design ensures a secure grip on small or delicate items, while the 11" length provides comfortable reach and control. These high-quality stainless steel tools are durable, heat-resistant, and ideal for maintaining sterile conditions in medical and laboratory settings,

making them a must-have for ensuring aseptic techniques.

**3 PRONGS STERILIZER FORCEPS
FEATURES**

- Sterilizer Utility Forceps 11" (20.3cm) Long.
- Slightly Curved Jaws Inward Towards Each Other.
- Three Prongs Resembling Claw for Secure Grip.
- Boxlock Joint Promoting Smooth Jaws Movements.
- Long Curved Shanks to Hold Instruments From a Safe Distance.
- Ergonomic Finger Rings Without Lock Handle Ensuring Maximum Control.
- User-friendly, Easy to Clean, Autoclavable.



**WI-70000336****Three Prong Sterilizer Forceps 8"****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Three Prong Sterilizer Forceps 8" (20.3cm) Premium Quality Stainless Steel Surgical Medical Sterilizing Aid.

3 Prong Sterilizer Forceps is essential surgical tools for safely handling sterile items, primarily used to retrieve sterilized instruments from autoclaves or hot solutions without contamination. These utility sterilization forceps featuring feature three prongs that curve inward, resembling a claw, providing a secure grip on specific items. Their three prong design ensures a secure grip on small or delicate items, while the 8" length provides comfortable reach and control. These high-quality stainless steel tools are durable, heat-resistant, and ideal for maintaining sterile conditions in medical and laboratory settings,

making them a must-have for ensuring aseptic techniques.

3 PRONGS STERILIZER FORCEPS FEATURES

- Sterilizer Utility Forceps 8" (20.3cm) Long.
- Slightly Curved Jaws Inward Towards Each Other.
- Three Prongs Resembling Claw for Secure Grip.
- Boxlock Joint Promoting Smooth Jaws Movements.
- Long Curved Shanks to Hold Instruments From a Safe Distance.
- Ergonomic Finger Rings Without Lock Handle Ensuring Maximum Control.
- User-friendly, Easy to Clean, Autoclavable.

**WI-70000794****Thumb Tissue Forceps 1x2 Teeth 4.5"****\$8.99** per pieceTrade price from **\$8.27**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thumb Tissue Forceps 1x2 Teeth 4.5" (11.4cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Rat Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 1x2 rat tooth blunt tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 4.5" (11.4cm) Overall Length.
- 1x2 Rat Tooth Blunt Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.

**WI-70000792****Thumb Tissue Forceps 1x2 Teeth 5.5"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thumb Tissue Forceps 1x2 Teeth 5.5" (14cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Rat Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 1x2 rat tooth blunt tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 5.5" (14cm) Overall Length.
- 1x2 Rat Tooth Blunt Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.

WI-70000793**Thumb Tissue Forceps 1x2 Teeth 5"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thumb Tissue Forceps 1x2 Teeth 5" (12.7cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Rat Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 1x2 rat tooth blunt tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 5" (12.7cm) Overall Length.
- 1x2 Rat Tooth Blunt Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.

**WI-70000791****Thumb Tissue Forceps 1x2 Teeth 6"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Thumb Tissue Forceps 1x2 Teeth 6" (15.2cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Rat Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 1x2 rat tooth blunt tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 6" (15.2cm) Overall Length.
- 1x2 Rat Tooth Blunt Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.

**WI-70000790****Thumb Tissue Forceps 1x2 Teeth 8"****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thumb Tissue Forceps 1x2 Teeth 8" (20.3cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Rat Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 1x2 rat tooth blunt tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 8" (20.3cm) Overall Length.
- 1x2 Rat Tooth Blunt Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.



**WI-70000787****Thumb Tissue Forceps 2×3 Teeth 5.5"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thumb Tissue Forceps 2×3 Teeth 5.5" (14cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 2×3 blunt tooth tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 5.5" (14cm) Overall Length.
- 2×3 Blunt Teeth Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.

WI-70000788**Thumb Tissue Forceps 2×3 Teeth 5"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thumb Tissue Forceps 2×3 Teeth 5" (12.7cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 2×3 blunt tooth tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 5" (12.7cm) Overall Length.
- 2×3 Blunt Teeth Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.

**WI-70000786****Thumb Tissue Forceps 2×3 Teeth 6"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thumb Tissue Forceps 2×3 Teeth 6" (15.2cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 2×3 blunt teeth tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 6" (15.2cm) Overall Length.
- 2×3 Blunt Teeth Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.

WI-70000783**Thumb Tissue Forceps 3×4 Teeth 5.5"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thumb Tissue Forceps 3×4 Teeth 5.5" (14cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 3×4 blunt teeth tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 5.5" (14cm) Overall Length.
- 3×4 Blunt Teeth Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.

**WI-70000784****Thumb Tissue Forceps 3×4 Teeth 5"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thumb Tissue Forceps 3×4 Teeth 5" (12.7cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 3×4 blunt teeth tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 5" (12.7cm) Overall Length.
- 3×4 Blunt Teeth Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.

WI-70000782**Thumb Tissue Forceps 3×4 Teeth 6"****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Thumb Tissue Forceps 3×4 Teeth 6" (15.2cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 3×4 blunt teeth tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 6" (15.2cm) Overall Length.
- 3×4 Blunt Teeth Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.

**WI-70000780****Thumb Tissue Forceps 4×5 Teeth 5.5"****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Thumb Tissue Forceps 4×5 Teeth 5.5" (14cm) Multipurpose Dissection Tweezer Forceps with Precision Blunt Tooth Straight Tips for Medical Surgical and Lab Settings.

Thumb tissue forceps is a multipurpose surgical tweezer extensively used by professionals to manipulate tissues, organs, and other anatomical structures. In laboratory research, forensic science and cell culture work this can be used for handling small specimens, transferring samples, or manipulating fragile materials such as collecting small fibers, hairs, or other trace materials from crime scenes or bodies. Thumb tissue features a 4×5 blunt teeth tips with non-slip grip knurled handle are widely within the lab and medical settings due to their

premium quality materials and long lasting properties. We have wide range of teeth and length variations to choose from according to each of your procedure's requirements. All instruments are independently tested to ensure the performance you expect, all at affordable prices.

THUMB TISSUE FORCEPS FEATURES:

- Crafted with 5.5" (14cm) Overall Length.
- 4×5 Blunt Teeth Tips.
- Straight Jaws Profile for Enhanced Sightlines.
- Ergonomic Spring Handle Promoting Superior Control.
- Non-Slip Surface.
- Convenient to Hold.

**WI-70002654****Tischler Kevorkian Biopsy Punch Forceps 3mm x 9.5mm Bite****\$59.99** per pieceTrade price from **\$55.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Tischler-Kevorkian Biopsy Forceps Gold Handle 3.0 mm x 9.5 mm Large Rectangular Bite, 9.75" (25cm) Long Shaft with Locking Tip B/GYN Medical Surgical Instrument.

Cervical biopsy is a surgical procedure in which a small amount of tissue is removed from the cervix. During the colposcopy, the clinician looks at the cervix and biopsies any areas that look abnormal. Tischler Biopsy Punch is an OB/GYN instrument used in colposcopy procedures to obtain cervical tissue samples for biopsy.

The Tischler-Kevorkian Biopsy Forceps is reusable, with a shaft length of 9.75" (250mm) designed with a wider bit other than standard punches. The 3.0mm x 9.5mm bite opened with 60

degree will maximize tissue sampling while being less traumatic at the same time, hence reducing the chances of post-procedure bleeding. The Tischler-Kevorkian Forceps is equipped with a Tip-Locking system and feature a single upper and lower teeth on the rectangular bite for added traction of the tissue for biopsy sample. The full finger spring handle allows for precise, controlled cutting of the tissue sample.

TISCHLER-KEVORKIAN BIOPSY PUNCH FEATURES

TISCHLER BIOPSY PUNCH: Tischler Cervical Punch Biopsy Forceps is a Non-ratcheted, Reusable double spring OB/GYN instrument used in colposcopy procedures to obtain cervical tissue samples for biopsy. With a built-in





WI-70002653

Tischler-Morgan Biopsy Forceps Oval Shape 3mmx7mm Bite

\$59.99 per piece

Trade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Tischler-Morgan Biopsy Forceps 9.75" (25cm) Length Oval Shape 3.0mm x 7.0mm Bite with Gold Handle Premium Quality Stainless Steel OB/GYN Medical Surgical Instrument.


Biopsy Forceps are also called Biopsy Punch and commonly used in Gynecology Surgical Procedures. Tischler-Morgan Forceps are Grasping Biopsy Forceps with Inter fitted Cups. The length of the shank is 25cm. Tischler-Morgan Biopsy Forceps are used to grasp and remove cervical tissue during a cervical biopsy.

Morgan Tischler Punch forceps have standard handle and one large ring for multiple fingers, the handle has a thumb rest for precise control. Forceps have a straight shank and an oval tip with double wire. The upper and lower

jaw has single teeth which prevents slipping of the tissue. Forceps have double leaf spring for easy working. The forceps cups can fit into each other, it reduces the risk of losing the tissue sample as the instrument is withdrawn.

TISCHLER MORGAN BIOPSY FORCEPS FEATURES

TISCHLER-MORGAN BIOPSY PUNCH: Tischler Cervical Punch Biopsy Forceps is a Non-ratcheted, Reusable double spring OB/GYN instrument used in colposcopy procedures to obtain cervical tissue samples for biopsy. With a built-in thumb lock mechanism that secures and protects the cutting tip, our Tischler Biopsies meet and exceed the industry standard for cervical biopsy instruments.



WI-36089992

Tissue Forceps Brown Style Stainless Steel Surgical

\$3.00 per piece

Trade price from **\$0.90**

Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Tissue Forceps Brown Style Stainless Steel Surgical is a forcep we forge, grind and finish in our own Sialkot workshop. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.00, falling to \$0.90 at 5 pieces — about 70% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002917

Tohoku Towel Forceps

\$30.21 per piece

Trade price from \$27.79
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Tohoku Towel Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$30.21, falling to \$27.79 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002916

Towel Forceps

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Towel Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$8.99, falling to \$8.27 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001053****Troeltsch Wilde Ear Forceps 5" Angular 1x2 Teeth****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Troeltsch Ear Forceps 1x2 Teeth 5 Inches (12.7cm) Angular German Grade Stainless Steel Surgical ENT Rhinology Instruments.

Troeltsch ear forceps or Wilde-Troeltsch forceps are thumb forceps commonly used in ENT procedures. Troeltsch (Wilde) Ear Forceps assist doctors to pack gauzes in ear procedures and deliver medications into the external ear canal. This ear forceps feature long angled jaws and added outer serrations on the handles to prevent slips and slides. This forcep designed with a wide angle between the handles and the jaws to give a clear field of vision. Made with high quality stainless steel with multiple variation suit to various surgical indications.

TROELTSCH WILDE EAR FORCEPS FEATURES

- 5" (12.7cm) Overall Length.
- Long Angular Shaped Jaws.
- 1x2 Teeth Tips.
- Outer Serrated Handles.
- Multiple Jaws For Diverse Indications.
- Reusable.
- Made with Premium Quality Stainless Steel

PACJAGE INCLUDE

- 1 Pc Troeltsch Wilde Ear Forceps 5" Angular 1x2 Teeth

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Ear Forceps**APPLICATION:** Surgical ENT**WI-70001824****Trousseau Tracheal Dilator Forceps 5.5" ENT Tracheostomy Instrument****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Trousseau Bowlby Tracheal (Tracheostomy) Dilator Forceps 5.5" (14cm) Blunt, Smooth Curved Jaws with Single Spring Non-Ratcheted Style Premium Quality Stainless Steel ENT Tracheostomy Surgical Instruments.

The Trousseau Bowlby tracheal dilating forceps used in ENT to dilate tracheostomy tracks. This is highly versatile single spring, non-ratcheted, and finger ring instrument widely used to widen the tracheal incision for easy tube insertion placement. Trousseau Tracheal Dilators help surgeons to place tubes in the trachea in cases of airway obstruction. This is when doctors make an opening in the trachea to bypass a block in the upper throat. They are crafted from high-grade German grade surgical stainless


material, and thus can be reused after sterilization.

TROUSSEAU TRACHEAL DILATOR FEATURES

- Crafted with 5.5" (14cm) Overall Length.
- Finger Ring Handles for Controlled Movements.
- Angled Jaws for Easy Surgical Access.
- Blunt Prongs for Safe Manipulations.
- Single Spring Non-Ratcheted Style Handle.
- Reusable.
- Made with German Grade Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Trousseau Tracheal Dilator 14cm





WI-70003699

Ulrich Dressing Forceps

\$30.84 per piece

Trade price from \$28.37
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Ulrich Dressing Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$30.84, falling to \$28.37 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004154

Universal Pado Extracting Forceps #150 1/2S for Upper Teeth

\$11.99 per piece

Trade price from \$11.03
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Pado Extracting Forceps #150 1/2S for Upper Teeth (Incisors, Bicuspid and Roots) with Spring Handle Stainless Steel American Pattern Tooth Extraction Dental Surgical Instruments for Children.

American pattern child tooth extraction forceps No. 150 1/2S (#150SK) used to extract upper primary teeth of children. The small, closely-situated beaks ensures maximum accessibility and a tight grasp of the tooth, while minimizing chances of trauma to the patient. #150 1/2S universal pado extraction forceps guarantees the easiest tooth extraction procedure designed to fit in the palm of your hands for patient anxiety ease.

American pattern forceps have a smooth hinge so no pin is visible. The

high quality stainless steel construction ensures maximum durability and dependability. You can reuse these extracting forceps over time as they will not bend or tarnish on you.

UNIVERSAL PADO EXTRACTING FORCEPS #150 1/2S FEATURES

- Universal Tooth Extraction Forceps for Children.
- For Upper Incisors, Bicuspid and Roots.
- Small, Conical Beaks for Optimal Grip.
- American Pattern.
- Non-slip Grip Spring Handle.
- Made With Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Pado Extracting Forceps #150 1/2S – Universal

**WI-70004119**

Universal Pedo Extracting Forceps #151 1/2S For Lower Teeth

\$11.99 per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Universal Pedo Extracting Forceps #151-1/2S With Spring Handle for Lower Incisors, Bicuspid and Root – Premium Quality Stainless Steel American Pattern Child Tooth Extraction Dental Surgical Instrument

Pedo extracting forceps #151-1/2S are specifically designed for extracting lower incisors, canines, premolars, and roots. Dental forceps no. 151 1/2S are universal tools for extracting teeth from alveolar bone, with an enhanced beak geometry for optimal access and engagement. Our pediatric extraction forceps are designed to fit comfortably in the palm, reducing patient anxiety and American pattern feature a smooth hinge without a visible pin.

151-1/2S PEDO EXTRACTING FORCEPS FEATURES

- Pediatrics Forceps For Lower Incisor, Bicuspid and Root.
- 1/2 Size Smaller Than Adult #151 Forceps.
- American Pattern Design.
- Textured Grip Single Spring Handle.
- Fully Autoclavable and Reusable.
- Made of High Grade Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Universal Pedo Extracting Forceps #151 1/2S With Spring Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENTS TYPE:** Dental Extracting Forceps – Pediatric**APPLICATION:** Lower Teeth Extraction Forceps**REUSABLE/AUTOCCLAVABLE:** Yes**WI-70001703**

University of Washington Rubber Dam Clamp Forceps 7"

\$10.99 per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Rubber Dam Clamp Forceps – University of Washington Pattern 7" (17.8m) Stainless Steel Endodontic Dental Instruments


University of Washington rubber dam forceps is a great addition to your endodontic or dental practice.

University of Washington type rubber dam forceps have a spring to allow for quick and efficient placement and removal of clamps. A locking ring design can hold the forceps in any position, leaving you room to work. The length of these rubber dam clamp forceps allows you to reach any tooth, while the slim design gives you plenty of maneuverability and increased visibility. Made from german grade stainless steel, you can use these rubber dam forceps to place and hold

rubber dam clamps around teeth. These forceps are designed for comfort and performance in mind. These forceps are also autoclavable.

WASHINGTON DAM CLAMP FORCEPS FEATURES

- Perfect for quick and efficient clamp placement and removal.
- Slim pattern for easy access and visibility.
- Smooth spring enables quick and easy placement and removal of clamps.
- Locking ring mechanism to hold forceps in any position.
- Have a Curved Handle for Comfort and Efficiency to Reach Any Tooth.
- Autoclavable.





WI-36123694

Vasectomy Meatotomy Instrument Set Stainless Steel Forceps

\$150.00 per piece

Trade price from **\$105.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Vasectomy Meatotomy Instrument Set Stainless Steel Forceps, made for general and specialty operating theatre use. Instruments of this type grasp, hold and transfer tissue, dressings or small components without crushing them. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$150.00, falling to \$105.00 at 1 set — about 30% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004063

Verbrugge Bone Holding Forceps 10" Self Centering with Speed Lock Premium Quality

\$72.70 per piece

Trade price from **\$66.88**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Verbrugge Bone Holding Forceps 10" Self Centering with Speed Lock Premium Quality is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$72.70, falling to \$66.88 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70004062

Verbrugge Bone Holding Forceps 6" Self Centering with Speed Lock Premium Quality

\$41.76 per piece

Trade price from \$38.42
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Verbrugge Bone Holding Forceps 6" Self Centering with Speed Lock Premium Quality is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$41.76, falling to \$38.42 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000486

Virtus Splinter Forceps 6" (15.2cm) Straight Serrated Tips

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Waqar International Virtus Splinter Forceps 6" Straight Serrated for Medical and Veterinary Field to Remove Larger Foreign Objects.

Virtus splinter forceps measure 6" and are commonly used in the medical surgical and veterinary field to remove larger foreign objects. This Splinter forceps have a spring action pliers with a durable handle for precise operation. They are straight with serrated tip and sharp triangle point. The fine point jaws ensure access into deep cavities and wounds. The transverse serrations and vertical grooves on the jaws ensure a firm grip. Made with premium quality stainless steel that is durable and corrosion resistant. These forceps can be easily cleaned and autoclaved for many uses.

VIRTUS SPLINTER FORCEPS FEATURES

- 6" (15.2cm) Forceps Used in Medical, Surgical, and Veterinary Fields.
- Straight With Serrated Tip and Sharp Triangle Point.
- Fine Point Jaws for Deep Cavity Access.
- Transverse Serrations and Vertical Grooves for Firm Grip.
- Features a Spring Action Pliers With Durable Handle.
- Easy to Clean and Autoclave for Multiple Uses.
- Made From Durable, Corrosion-resistant Stainless Steel.

PACKAGE INCLUDES:

- 1 Pc Virtus Splinter Forceps 6" (15.2cm) Plier Handle





WI-70000810

Virtus Splinter Forceps 6" Straight Serrated Triangle Tips

\$10.99 per piece

Trade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Virtus Splinter Forceps 6" (15.2cm) Fine Point Straight Serrated Triangle Tips with Spring Action Pliers Grip Handle Premium Quality Stainless Steel Medical, Veterinary and Surgical Instruments.

Virtus splinter forceps measure 6" and are commonly used in the medical surgical and veterinary field to remove larger foreign objects. This Splinter forceps have a spring action pliers with a durable handle for precise operation. They are straight with serrated tip and sharp triangle point. The fine point jaws ensure access into deep cavities and wounds. The transverse serrations and vertical grooves on the jaws ensure a firm grip. Made with premium quality stainless steel that is durable and corrosion resistant. These forceps can

be easily cleaned and autoclaved for many uses.

VIRTUS SPLINTER FORCEPS FEATURES

- Crafted with 6" (15.2cm) Overall Length.
- Fine Point Tips for for Deep Surgical Access.
- Triangle Shaped Straight Serrated Jaws for Secure Holding.
- Ergonomic Single Spring Action Plier Style Handle Ensuring Superior Control.
- Ideal for Removal of Larger Foreign Objects.
- Full Autoclaved and Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDES:



WI-70003581

Volume Lashes Tweezers Best Volum Lashes Tweezers for Professional Eyelash Extension

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Volume Lashes Tweezers Best Volum Lashes Tweezers for Professional Eyelash Extension is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003579**

Whole Sale Eyelash Tweezers Strong Curved and Straight Eye Lash Extension Tweezers

\$15.99 per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Tweezers for Eyelash Extension – Straight and Curved Pointed Tweezers – Professional Stainless Steel Precision Tweezers Set – 2 Pcs – (Silver)

Perfect for professional eyelash extensions, these tweezers allow most precise application of designs and lashes strips, help to prevent lash extensions from sticking together during application. These eyelash tweezers are great for times when you want perfect isolation or great grip for picking up lashes. Great for Beauty Salon and individual eyelash extensions.

EYELASH EXTENSION TWEEZERS FEATURES:



EYELASH EXPERT TWEEZER: Tools that a lash artist/beginner must have.

Durable enough for beauty salons or home use. Perfect gift for your loved ones a beautiful pair of modern, sleek tweezers.

PRECISE: The edge of every pointed tweezer is aligned to ensure a great tweeze every time. protective tip covers included so tweezers never lose alignment.

PROFESSIONAL MAKING: Perfect for selecting, isolating and picking up eyelashes. The fine pointed tips make it easier to maneuver in any direction or angle. Light weight so does not cause hand fatigue.

EASY PICKUP: The straight tips are specifically designed for the easy pick up of multiple hairs off the tray. The angled tip is perfect for selecting an individual lash or to separate a client's natural lashes in order to isolate one

**WI-70000624**

Wilde Blakesley Ethmoid Forceps 5" Shaft 45° Up 4x12mm Pointed Fenestrated Cup

\$29.99 per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Wilde Blakesley Ethmoid Forceps 5" (12.7cm) Shaft 4 × 12mm Angled 45 Degree Up Pointed Fenestrated Cup Bite With Ring Handle Premium Quality Stainless Steel ENT Rhinology Instruments.


Wilde-Blakesley Ethmoid Forceps is a specialized ENT surgical tool used for grasping, biopsying, and excising delicate tissues, nasal polyps, and bone fragments in the ethmoid sinus cavity. It is essential during ENT surgical procedures such as Functional Endoscopic Sinus Surgery (FESS), Ethmoidectomy, and Nasal Polypectomy, allowing surgeons to navigate narrow passages and remove obstructions with minimal trauma.

Weil-Blakesley forceps feature 45-degree angled tip enables surgeons to

access tight lateral sinus recesses that are unreachable with straight. The fenestration tip design allows for better visualization and maneuverability in confined spaces. Wide finger ring handles provide exceptional control and comfort, reducing hand fatigue during prolonged procedures. Made with premium-grade stainless steel, these forceps are sterilized and can be reused.

WILDE BLAKESLEY NASAL FORCEPS FEATURES

- 5" (12.7cm) Working Shaft for Easy Surgical Access.
- Angled 45 Degree Profile.
- 4x12mm Pointed Fenestrated Cup Remove Soft Tissues or Bones in Narrow Sinuses.



**WI-70000625****Wilde Blakesley Ethmoid Forceps 5" Shaft
4x12mm Straight Pointed Fenestrated Cup****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Wilde Blakesley Ethmoid Forceps 5" (12.7cm) Shaft 4 × 12mm Bite Straight Pointed Fenestrated Cup With Ring Handle Premium Quality Stainless Steel ENT Rhinology Instruments.

Wilde-Blakesley Ethmoid Forceps is a specialized ENT surgical tool used for grasping, biopsying, and excising delicate tissues, nasal polyps, and bone fragments in the ethmoid sinus cavity.


It is essential during ENT surgical procedures such as Functional Endoscopic Sinus Surgery (FESS), Ethmoidectomy, and Nasal Polypectomy, allowing surgeons to navigate narrow passages and remove obstructions with minimal trauma.

Wilde-Blakesley features a 4×12mm straight pointed fenestrated cup providing excellent reach and access to

deeper nasal structures. The fenestration tip design allows for better visualization and maneuverability in confined spaces. Wide finger ring handles provide exceptional control and comfort, reducing hand fatigue during prolonged procedures. Made with premium-grade stainless steel, these forceps are sterilized and can be reused.

WILDE BLAKESLEY NASAL FORCEPS FEATURES

- 5.0" (12.7cm) Working Shaft for Easy Surgical Access.
- Straight Slender Profile to Reach Deep.
- 4×12mm Pointed Fenestrated Cup Bite Remove Soft Tissues or Bones in Narrow Sinuses.

**WI-70000623****Wilde Blakesley Ethmoid Forceps 5" Shaft 90°
Up 4x12mm Pointed Fenestrated Cup****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Wilde Blakesley Ethmoid Forceps 5" (12.7cm) Shaft 4 × 12mm Angled 90 Degree Up Pointed Fenestrated Cup Bite With Ring Handle Premium Quality Stainless Steel ENT Rhinology Instruments.

Wilde-Blakesley Ethmoid Forceps is a specialized ENT surgical tool used for grasping, biopsying, and excising delicate tissues, nasal polyps, and bone fragments in the ethmoid sinus cavity. It is essential during ENT surgical procedures such as Functional Endoscopic Sinus Surgery (FESS), Ethmoidectomy, and Nasal Polypectomy, allowing surgeons to navigate narrow passages and remove obstructions with minimal trauma.

Weil-Blakesley forceps feature 90-degree angled tip enables surgeons to

access superior anatomical structures like high ethmoid cells and frontal recess without obstructing their line of sight.. The fenestration tip design allows for better visualization and maneuverability in confined spaces. Wide finger ring handles provide exceptional control and comfort, reducing hand fatigue during prolonged procedures. Made with premium-grade stainless steel, these forceps are sterilized and can be reused.

WILDE BLAKESLEY NASAL FORCEPS FEATURES

- 5" (12.7cm) Working Shaft for Easy Surgical Access.
- Angled 90 Degree Profile.
- 4×12mm Pointed Fenestrated Cup Remove Soft Tissues or Bones in Narrow Sinuses.



**WI-70001278****William Splinter Forceps 5.5" (14cm) 45° Angle Serrated Tips****\$10.49** per pieceTrade price from **\$9.65**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

William Splinter Forceps 5.5" (14cm) 45° Angle Serrated Tips is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$10.49, falling to \$9.65 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable

**WI-70001277****William Splinter Forceps 5.5" (14cm) 90° Right Angle Pointed Serrated Tips****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Stieglitz William Splinter Forceps 5.5" (14cm) 90° Right Angle Pointed Serrated Tips with Grooves with Ratchet Ring Handle German Grade Stainless Steel Surgical Instrument.

Stieglitz (William) Splinter Forceps used in multiple areas to grasp, mobilize and remove foreign objects and other materials from the surgical field. This forcep provides a secure way to hold tissues, debris and foreign bodies for their atraumatic removal. Stieglitz Splinter Forceps are frequently used in dental procedures for removing root tips, holding pins and gutta-percha points.

William Splinter Forceps features a pointed jaws that ensure access into deep cavities and wounds. The transverse serrations and vertical

grooves on the jaws ensure a firm grip. Williams splinter forceps consists of straight, curved, and right-angle pattern so that the surgeons can reach the surgical target from multiple approaches.

STIEGLITZ WILLIAM SPLINTER FORCEPS FEATURES

- 90° Right Angle Pattern.
- Overall Length 5 1/2" (14cm).
- Serrated Tips with Vertical Grooves to Grasp Tissues.
- Pointed Jaws To Ensure Deep Access Reach.
- Ratchet System for Self-Retaining Clamping.
- Ergonomic Finger Rings Ensuring Optimal Handling.
- Reusable and Sterilizable.





WI-70001279

William Splinter Forceps 5.5" (14cm) Straight Serrated Tips

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

William Splinter Forceps 5.5" (14cm) Straight Serrated Tips is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable



WI-70003368

Wrigley Obstetrical Forceps

\$79.99 per piece

Trade price from \$73.59
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Wrigley Obstetrical Forceps is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$79.99, falling to \$73.59 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003578

X Type Curved Tweezers Eyelash Extension Angle Cross Tweezer

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

X Type Curved Tweezers Eyelash Extension Angle Cross Tweezer is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001223

XTRM CRAFT Tissue Forceps | Lung Holding Intestinal Grasping Forceps – Stainless Steel – Ergonomic, Precision-Engineered, Available in 6" and 9"

\$12.99 per piece

Trade price from \$11.95
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

XTRM CRAFT Tissue Forceps | Lung Holding Intestinal Grasping Forceps – Stainless Steel – Ergonomic, Precision-Engineered, Available in 6" and 9"

XTRM CRAFT introduces the Tissue Forceps | Lung Holding Intestinal Grasping Forceps Kit, an elite instrument set for medical professionals. Ergonomically designed for a firm, non-slip grip, these forceps are crucial for accurate tissue holding and manipulation. Crafted from high-quality stainless steel, they assure longevity and effective sterilization, crucial for maintaining cleanliness in medical settings. The unique triangular serrated jaw design enables precise control of tissue. Ideal for various

in 6" and 9" sizes, they meet procedural requirements.

LUNG GRASPING FORCEPS FEATURES

- **Ergonomic Handle:** Designed for medical professionals, XTRM CRAFT's Intestinal Grasping Forceps feature an ergonomic handle for a comfortable, non-slip grip, ensuring precise tissue manipulation.
- **Superior Stainless Steel Build:** Made with top-grade stainless steel, these forceps deliver exceptional strength and cleanliness, vital for sterile medical environments.
- **Triangular Serrated Jaws:** The precision-engineered, triangular serrated jaws offer excellent tissue



**WI-70003427****Young Tongue Forceps 16.5cm****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Young Tongue Forceps Surgical Oral Medical Clamp 6.5" (16.5cm) Tongue Seizing Forceps with Rubber Jaws Stainless Steel Dental Surgical Oral Instrument

Young Tongue Seizing Forceps is a specialized device that ENT and maxillofacial use for grasping, retracting and exposing the ventral surface of the tongue, especially during surgical procedures that require its excision or tissue harvesting. The forceps have rubber inserts to reduce trauma with curved shanks, finger rings and locking ratchets

YOUNG TONGUE FORCEPS FEATURES

QUALITY MATERIALS: These Young Tongue Holding Forceps are made from


premium quality aluminum steel, ensuring superior durability and corrosion resistance.

USER FRIENDLY: Designed to be intuitive and easy to handle, allowing for precision during oral procedures. Its user-friendly design makes it a must-have for medical professionals.

EASY TO CLEAN: These forceps feature a high mirror finish that is not only aesthetically pleasing, but also facilitates easy cleaning and sterilization, ensuring safe usage every time.

RUBBER INSERTS: Young forceps feature soft, grooved, rubber jaws. The rubber inserts provide a comfortable grip, allowing for more secure handling and control during surgical procedures, enhancing accuracy and safety.

VERSATILE USAGE: Perfect for various



SOFT, GROOVED RUBBER INSERTS TO REDUCE TONGUE TRAUMA | 100% MIRROR FINISH, STAINLESS STEEL | CURVED SHANKS FOR BETTER GRASP | EXCELLENT LOCKING MECHANISM

**WI-70000485****Young Tongue Forceps 7" (17.8cm) with Silicone Grip Tips****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Young Tongue Forceps 7 (17.8cm) Silicone Grip Tongue Seizing Fenestrated Jaws Stainless Steel Oral Dental Surgical Oral Instrument

Young Tongue Seizing Forceps are specialized tools used by ENT and maxillofacial surgeons to grasp, retract, and expose the ventral surface of the tongue during surgical procedures. Tongue holding forceps have rubber inserts to reduce trauma with curved shanks, finger rings and locking ratchets. Young Tongue Forceps offer unparalleled precision with an ergonomic design, ensuring optimal control during dental, ENT, and medical procedures. These Forceps are made from premium quality aluminum steel, ensuring superior durability and corrosion resistance facilitates easy

cleaning and sterilization, ensuring safe usage every time.

PACKAGE INCLUDE

- 1 Pc Young Tongue Forceps 7" (17.8cm)

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Dental Oral Instruments

APPLICATION: Dental, Oral & Maxillofacial, Hemiglossectomy

JAWS: Fenestrated Jaws With Rubber Inserts

SHANK: Curved

LOCK: Ratchet Lock Mechanism

HANDLE: Finger Ring

AUTOCCLAVABLE: Yes

REUSABLE: Yes

RUST PREVENTION PROCEDURE:



SOFT, GROOVED RUBBER INSERTS TO REDUCE TONGUE TRAUMA | 100% MIRROR FINISH, STAINLESS STEEL | CURVED SHANKS FOR BETTER GRASP | EXCELLENT LOCKING MECHANISM



WI-70003576

Your Best Choice High Quality Multicolor Eyelash Extension Tweezer

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Your Best Choice High Quality Multicolor Eyelash Extension Tweezer is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$9.99, falling to \$9.19 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

Request a quotation

Send us product codes from this catalogue, quantities and your country.

We reply the same working day on WhatsApp.

Online shop	shop.waqarintl.com
Email	sales@waqarintl.com
WhatsApp /	
Phone	+92 307 6173627
Works	Roras Road, Opposite Anwar Khawaja Industries, Sialkot 51310, Pakistan

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.