



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

Hooks & Needles Surgical Instruments

Reusable stainless-steel hooks & needles forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



87

PRODUCTS

\$0.50

PRICES FROM

OEM

PRIVATE LABEL

Global

WORLDWIDE
EXPORT

Contents

87 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

Product	Code	Price	Page
10 Metal Teasing Needles 6" Straight Pointed Tips	WI-70000778	\$19.99	5
10 Periodontal Probe CP15 UNC Color Coded Markings Single Ended	WI-70002062	\$59.99	5
10 Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended	WI-70002067	\$59.99	6
100 Dental Probe Explorer #5 – Double Ended	WI-70002155	\$116.13	6
12 Dental Probe Explorer #5 – Double Ended	WI-70002158	\$16.24	7
12 Rainbow Color Dental Explorer Probe #5 – Double Ended	WI-70002162	\$47.99	7
2 Periodontal Probe CP15 UNC Color Coded Markings Single Ended	WI-70002065	\$13.99	8
2 Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended	WI-70002070	\$13.99	8
23/CP11 Expro (Explorer Probe)	WI-70003223	\$16.49	9
25 Dental Probe Explorer #5 – Double Ended	WI-70002157	\$28.74	9
3 Dental Probe Explorer #5 – Double Ended	WI-70002160	\$9.99	10
3 IUD Extractor Hook 10 1/4" (26cm)	WI-70002387	\$19.99	10
3 Periodontal Probe CP15 UNC Color Coded Markings Single Ended	WI-70002064	\$19.99	11
3 Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended	WI-70002069	\$19.99	11
4 Joseph Skin Hooks 6.25" Double Sharp Prong 2, 5, 7 & 10mm	WI-70001020	\$29.99	12
5 Dental Probe Explorer #5 – Double Ended	WI-70002159	\$12.49	12
5 IUD Extractor Hook 10 1/4" (26cm)	WI-70002386	\$37.99	13
5 Periodontal Probe CP15 UNC Color Coded Markings Single Ended	WI-70002063	\$29.99	13
5 Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended	WI-70002068	\$29.99	14
50 Dental Probe Explorer #5 – Double Ended	WI-70002156	\$57.49	14
50 Teasing Needles 5.5" Bent Tips with Plastic Handle	WI-70000808	\$29.99	15
6 Rainbow Color Dental Explorer Probe #5 – Double Ended	WI-70002163	\$29.99	15
Autopsy Hammer with Hook Stainless Steel Post Mortem Examination	WI-36793611	\$28.00	16
Babinski Percussion Reflex Hammer 9" with Needle	WI-70000904	\$9.99	16
Buck Percussion Reflex Hammer with Needle and Brush	WI-70000903	\$14.99	17
Color Hair Extension Pliers Pulling Hook Tool	WI-28935429	\$15.00	17
Dejerine Percussion Hammer 8" without Needle	WI-70000889	\$9.99	18
Dejerine Percussion Hammer 9" with Needle	WI-70000890	\$12.99	18
Dental Probe Explorer #5 – Double Ended	WI-70002161	\$7.99	19
Double Ended Pratt Uterine Dilator 13/15 FR Curved	WI-70003555	\$10.30	19
Double Ended Pratt Uterine Dilator 17/19 FR Curved	WI-70003554	\$10.30	20
Double Ended Pratt Uterine Dilator 21/23 FR Curved	WI-70003553	\$15.46	20
Double Ended Pratt Uterine Dilator 25/27 FR Curved	WI-70003552	\$14.43	21
Double Ended Pratt Uterine Dilator 29/31 FR Curved	WI-70003551	\$14.43	21
Double Ended Pratt Uterine Dilator 33/35 FR Curved	WI-70003550	\$14.43	22
Double Ended Pratt Uterine Dilator 37/39 FR Curved	WI-70003549	\$14.43	22
Double Ended Pratt Uterine Dilator 41/43 FR Curved	WI-70003548	\$14.43	23
Explorer Probe #23 Single Ended Hollow Handle	WI-70000297	\$9.99	23
Explorer Probe #5 Double Ended Hollow Handle	WI-70000296	\$9.99	24

Contents (continued)

Product	Code	Price	Page
Frazier Dura Hook 5" (12.7cm) Sharp Point	WI-70000023	\$9.99	26
Gillies Skin Hook 7" Large 4mm Dia Sharp Round Tip	WI-70000763	\$9.99	26
Gillies Skin Hook 7" Small 3mm Dia Sharp Round Tip	WI-70000764	\$9.99	27
Graham Nerve Hook Blunt 6mm Deep 6.5" Neurosurgical Instruments	WI-70001407	\$9.99	27
Hegar Uterine Dilator 3mm-4mm Double Ended 7.5" (19.1cm)	WI-70003542	\$8.99	28
Hegar Uterine Dilator 5mm-6mm Double Ended 7.5" (19.1cm)	WI-70003541	\$8.99	28
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 11mm-12mm – Hollow	WI-70003546	\$8.99	29
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 13mm-14mm – Hollow	WI-70003545	\$8.99	29
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 15mm-16mm – Hollow	WI-70003544	\$8.99	30
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 17mm-18mm – Hollow	WI-70003543	\$8.99	30
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 19mm-20mm – Hollow	WI-70000773	\$9.99	31
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 21mm-22mm – Hollow	WI-70000772	\$9.99	31
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 23mm-24mm – Hollow	WI-70000771	\$9.99	32
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 25mm-26mm – Hollow	WI-70000770	\$9.99	32
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 27mm-28mm – Hollow	WI-70000769	\$9.99	33
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 29mm-30mm – Hollow	WI-70000472	\$10.99	33
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 7mm-8mm – Hollow	WI-70003540	\$8.99	34
Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 9mm-10mm – Hollow	WI-70003539	\$8.99	34
IUD Extractor Hook 10 1/4" (26cm)	WI-70002388	\$9.99	35
Jameson Large Muscle Hook Precision Eye Surgery Tool	WI-43635731	\$15.00	35
Jameson Muscle Hook Large Eye Surgery Instrument	WI-36409088	\$15.00	36
Joseph Skin Hook 6.25" Single Sharp Prong	WI-70001026	\$8.99	36
Joseph Skin Hook 6.25" Two Sharp Prong 10mm Wide	WI-70001022	\$8.99	37
Joseph Skin Hook 6.25" Two Sharp Prong 12mm Wide	WI-70001021	\$8.99	37
Joseph Skin Hook 6.25" Two Sharp Prong 2mm Wide	WI-70001025	\$8.99	38
Joseph Skin Hook 6.25" Two Sharp Prong 5mm Wide	WI-70001024	\$8.99	38
Joseph Skin Hook 6.25" Two Sharp Prong 7mm Wide	WI-70001023	\$8.99	39
Nabers Probe #Q-2N Color-Coated Marking 3,6,9,12 Hollow Handle	WI-70000224	\$10.99	39
News Tracheal Hook 6" Curved Sharp Point	WI-70000394	\$9.99	40
Periodontal Probe CP15 UNC Color Coded Markings Single Ended	WI-70002066	\$10.99	40
Periodontal Probe UNC15/23 Double Ended Hollow Handle	WI-70000142	\$10.99	41
Probe Eye-01 Stainless Steel Surgical Exploration Probe for Clinic	WI-36666994	\$0.50	41
Rainbow Color Dental Explorer Probe #5 – Double Ended	WI-70002164	\$10.99	42
Rentsch Boat Hook with Guard Precision Eye Surgery	WI-43668376	\$15.00	42
Reusable Probe with Eye Surgical Examination Instrument Medical	WI-45125898	\$0.59	43
Sinskey II Angled Lens Hook Precision Eye Surgery	WI-43644620	\$15.00	43
Sinskey II Lens Hook Angled-06 Ophthalmology Surgical Instrument	WI-46285630	\$15.00	44
Sinskey Lens Hook Straight Precision Eye Surgery Tool	WI-43638650	\$15.00	44
Stainless Steel Comfort Grip	WI-27945324	\$20.00	45
Stainless Steel Mini Hair Extension Plier	WI-25801265	\$3.50	45
Surgical	WI-28938462	\$15.00	46
Surgical Probe with Eye Reusable Stainless Steel Medical Diagnostic	WI-45137869	\$0.50	46

Contents (continued)

Product	Code	Price	Page
Tracheotomy Hook 15cm Sharp Point Stainless Steel Surgical Instrument	WI-43358584	\$7.00	47
Vera Hook Reusable Stainless Steel Eye Surgery	WI-43650631	\$15.00	47
Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended	WI-70002071	\$10.99	48

**WI-70000778****10 Metal Teasing Needles 6" Straight Pointed Tips****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Pcs Metal Teasing Needles Set Straight Pointed Tips – Medical Lab Dissecting Needle Stainless Steel Single Ended Micro Fine Precision Probe Teasing Needles Lab Biology Experiment Equipment for Medical Students.

Teasing needles are used to manipulate small samples. These are long stainless steel dissection needle available with straight sharp tips. Dissecting teasing needle is a multi-functional tool that can be used for dissecting, sampling, cutting, puncturing and other operations. With its sharp tip, it can easily penetrate tissues made of different materials and can also be used for sampling and cutting tissues. With very fine sharp tips which can precisely penetrate cells and tissues for fine operations teasing needles come with comfortable to hold solid non-slip handle

design which reduces hand fatigue and improve operation precision.

METAL DISSECTING TEASING NEEDLES FEATURES

- 6" (15.2cm) Overall Length.
- Sharp-pointed, Straight Needle.
- Single Ended Stainless Steel Tips.
- Non-Slip Central Textured Grip Handle
- Reusable.
- Easy to Clean.
- Perfect for Use in Dissection, Biology Lab and Tissue Manipulation.
- Made with Stainless Steel Tips.

PACKAGE INCLUDE

- 10 Pcs Metal Teasing Needle Set Sharp Pointed Straight Tips

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**WI-70002062****10 Periodontal Probe CP15 UNC Color Coded Markings Single Ended****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Dental Probe UNC15 Color Coded Marking 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 Single Ended with Hollow Handle Periodontal Diagnostic Instruments

Dental Probe CP-15 Periodontal Probe UNC15 color coded markings are used to assess pocket depths, periodontal attachment levels, and gingival bleeding. It is used to check the depth of periodontal pocket vertically or horizontally. Color markings help determine the depth easily. Accurate & smooth black markings finish, which will not chip, flake or fade.

UNC15 PROBE FEATURES

- Single Ended.
- CP-15 UNC probe markings: 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 on one end.


- Designed to assess periodontal pocket depths.
- Working end has calibrated markings in millimeters
- Made of high-quality stainless steel.
- Ergonomic Hollow Handle.

PACKAGE INCLUDE:

- 10 Pcs Periodontal Probe UNC15 Color Coded Marking 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 Single Ended with Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Probe**APPLICATION:** Use to Measure Periodontal Pocket Depths & Attachment Levels**COLOR:** Silver Color**TIP:** Single Ended

PERIODONTAL PROBE CP15



**WI-70002067****10 Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Dental Williams Probe Color Coded 1-2-3-5-7-8-9-10 Marking Root Depth Measurement Periodontal Examination Instruments

Williams Probe is also called the periodontal probe it is the first probe to be used to assess probing depth of a pocket with color codes. Williams Probe used to assess periodontal pocket depths, attachment levels, anatomy configurations and gingival bleeding. It has the following markings – 1,2,3,5,7,8,9,10 mm. It measures the sulcus or pocket depth of the periodontium of each tooth. The working end has calibrated markings in millimeters.


WILLIAMS PROBE FEATURES

- Single Ended.

- Williams probe markings: 1-2-3-5-7-8-9-10 on one end.
- Designed to assess periodontal pocket depths.
- Working end has calibrated markings in millimeters
- Made of high-quality stainless steel.
- Ergonomic Hollow Handle.

PACKAGE INCLUDE:

- 10 Pcs Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended with Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Probe**APPLICATION:** Use to Measure Periodontal Pocket Depths & Attachment Levels**COLOR:** Silver Color**100
PCS****WI-70002155****100 Dental Probe Explorer #5 – Double Ended****\$116.13** per pieceTrade price from **\$106.84**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 Dental Probe 17/23 Double Ended Diagnostic Examination Dental Instrument

The Explorer 5 is one of the most popular explorers for caries detection. The explore #5 designed for caries and calculus detection, exploration of pocket characteristics, furcations, and restorations. This dual-sided explorer features a Shepards hook #23 end for caries detection and a #17 end for anterior interproximal caries detection. Featuring the most ergonomic handle with knurling to the instrument shank to provide increased comfort with less hand fatigue.

DENTAL EXPLORER PROBE 17/23 FEATURES

THE FIRST QUALITY CARE TOOL: Our dental tarter scraper is made of 100%


high-purity medical stainless steel metal, which allows rust resistance, durability and safety for the tool. It also has double-ended design, cambered and curving heads that can meet your different needs.

EASY TO USE: Wipe the plaque off the blade with a tissue, then repeat until there is no more plaque on the blade. With this amazing stainless steel dental scarper, you can easily remove those food residues located in the teeth without causing pain and damaging tooth enamel.

CONVENIENT TO BRING AWAY: Thanks to the small size, you can bring it anywhere as personal belongings, no matter where you travel to make you always keep a fresh tone.

100% SATISFACTION GUARANTEE: Customer satisfaction is our utmost priority.



**WI-70002158****12 Dental Probe Explorer #5 – Double Ended****\$16.24** per pieceTrade price from **\$14.94**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Dental Probe 17/23 Double Ended Diagnostic Examination Dental Instrument

The Explorer 5 is one of the most popular explorers for caries detection. The explore #5 designed for caries and calculus detection, exploration of pocket characteristics, furcations, and restorations. This dual-sided explorer features a Shepards hook #23 end for caries detection and a #17 end for anterior interproximal caries detection. Featuring the most ergonomic handle with knurling to the instrument shank to provide increased comfort with less hand fatigue.

DENTAL EXPLORER PROBE 17/23 FEATURES



THE FIRST QUALITY CARE TOOL: Our dental tarter scraper is made of 100%

high-purity medical stainless steel metal, which allows rust resistance, durability and safety for the tool. It also has double-ended design, cambered and curving heads that can meet your different needs.

EASY TO USE: Wipe the plaque off the blade with a tissue, then repeat until there is no more plaque on the blade. With this amazing stainless steel dental scarper, you can easily remove those food residues located in the teeth without causing pain and damaging tooth enamel.

CONVENIENT TO BRING AWAY: Thanks to the small size, you can bring it anywhere as personal belongings, no matter where you travel to make you always keep a fresh tone.

100% SATISFACTION GUARANTEE: Customer satisfaction is our utmost priority.

**WI-70002162****12 Rainbow Color Dental Explorer Probe #5 – Double Ended****\$47.99** per pieceTrade price from **\$44.15**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pieces Multi Rainbow Color Dental Probe 17/23 Double Ended Diagnostic Examination Dental Instrument

Dental Explorer Probe is a professional dental scraper tool helps you to remove stains, tartar, and debris from teeth and gums. Regular care of teeth can reduce the generation of bad breath. This hygiene tool is also used for keeping your pet clean and oral fresh. You can be comfortable to be closer to your pet (cat, dog) without any fears about your health problems.

DENTAL EXPLORER PROBE 17/23 FEATURES

THE FIRST QUALITY CARE TOOL: Our dental tarter scraper is made of 100% high-purity medical stainless steel metal, which allows rust resistance, durability and

safety for the tool. It also has double-ended design, cambered and curving heads that can meet your different needs.

EASY TO USE: Wipe the plaque off the blade with a tissue, then repeat until there is no more plaque on the blade. With this amazing stainless steel dental scarper, you can easily remove those food residues located in the teeth without causing pain and damaging tooth enamel.

CONVENIENT TO BRING AWAY: Thanks to the small size, you can bring it anywhere as personal belongings, no matter where you travel to make you always keep a fresh tone.

100% SATISFACTION GUARANTEE: Customer satisfaction is our utmost priority.

PACKAGE INCLUDE:

**WI-70002065**

2 Periodontal Probe CP15 UNC Color Coded Markings Single Ended

\$13.99 per pieceTrade price from **\$12.87**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Dental Probe UNC15 Color Coded Marking 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 Single Ended with Hollow Handle Periodontal Diagnostic Instruments

Dental Probe CP-15 Periodontal Probe UNC15 color coded markings are used to assess pocket depths, periodontal attachment levels, and gingival bleeding. It is used to check the depth of periodontal pocket vertically or horizontally. Color markings help determine the depth easily. Accurate & smooth black markings finish, which will not chip, flake or fade.

UNC15 PROBE FEATURES

- Single Ended.
- CP-15 UNC probe markings: 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 on one end.

- Designed to assess periodontal pocket depths.
- Working end has calibrated markings in millimeters
- Made of high-quality stainless steel.
- Ergonomic Hollow Handle.


PACKAGE INCLUDE:

- 2 Pcs Periodontal Probe UNC15 Color Coded Marking 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 Single Ended with Hollow Handle

TECHNICAL SPECIFICATIONS

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Dental Probe**APPLICATION:** Use to Measure Periodontal Pocket Depths & Attachment Levels**COLOR:** Silver Color**TIP:** Single Ended

PERIODONTAL PROBE CP15

**WI-70002070**

2 Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended

\$13.99 per pieceTrade price from **\$12.87**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Dental Williams Probe Color Coded 1-2-3-5-7-8-9-10 Marking Root Depth Measurement Periodontal Examination Instruments

Williams Probe is also called the periodontal probe it is the first probe to be used to assess probing depth of a pocket with color codes. Williams Probe used to assess periodontal pocket depths, attachment levels, anatomy configurations and gingival bleeding. It has the following markings – 1,2,3,5,7,8,9,10 mm. It measures the sulcus or pocket depth of the periodontium of each tooth. The working end has calibrated markings in millimeters.

WILLIAMS PROBE FEATURES

- Single Ended.

- Williams probe markings: 1-2-3-5-7-8-9-10 on one end.
- Designed to assess periodontal pocket depths.
- Working end has calibrated markings in millimeters
- Made of high-quality stainless steel.
- Ergonomic Hollow Handle.


PACKAGE INCLUDE:

- 2 Pcs Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended with Hollow Handle

TECHNICAL SPECIFICATIONS

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Dental Probe**APPLICATION:** Use to Measure Periodontal Pocket Depths & Attachment Levels**COLOR:** Silver Color

PERIODONTAL PROBE WILLIAMS



**WI-70003223****23/CP11 Expro (Explorer Probe)****\$16.49** per pieceTrade price from **\$15.17**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

23/CP11 Expro (Explorer Probe) is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$16.49, falling to \$15.17 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70002157****25 Dental Probe Explorer #5 – Double Ended****\$28.74** per pieceTrade price from **\$26.44**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 25 Dental Probe 17/23 Double Ended Diagnostic Examination Dental Instrument

The Explorer 5 is one of the most popular explorers for caries detection. The explore #5 designed for caries and calculus detection, exploration of pocket characteristics, furcations, and restorations. This dual-sided explorer features a Shepards hook #23 end for caries detection and a #17 end for anterior interproximal caries detection. Featuring the most ergonomic handle with knurling to the instrument shank to provide increased comfort with less hand fatigue.

high-purity medical stainless steel metal, which allows rust resistance, durability and safety for the tool. It also has double-ended design, cambered and curving heads that can meet your different needs.

EASY TO USE: Wipe the plaque off the blade with a tissue, then repeat until there is no more plaque on the blade. With this amazing stainless steel dental scarper, you can easily remove those food residues located in the teeth without causing pain and damaging tooth enamel.

CONVENIENT TO BRING AWAY: Thanks to the small size, you can bring it anywhere as personal belongings, no matter where you travel to make you always keep a fresh tone.

100% SATISFACTION GUARANTEE: Customer satisfaction is our utmost priority.

DENTAL EXPLORER PROBE 17/23 FEATURES

THE FIRST QUALITY CARE TOOL: Our dental tarter scraper is made of 100%

**WI-70002160****3 Dental Probe Explorer #5 – Double Ended****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Probe 17/23 Double Ended Diagnostic Examination Dental Instrument

The Explorer 5 is one of the most popular explorers for caries detection. The explore #5 designed for caries and calculus detection, exploration of pocket characteristics, furcations, and restorations. This dual-sided explorer features a Shepards hook #23 end for caries detection and a #17 end for anterior interproximal caries detection. Featuring the most ergonomic handle with knurling to the instrument shank to provide increased comfort with less hand fatigue.

DENTAL EXPLORER PROBE 17/23 FEATURES


THE FIRST QUALITY CARE TOOL: Our dental tarter scraper is made of 100%

high-purity medical stainless steel metal, which allows rust resistance, durability and safety for the tool. It also has double-ended design, cambered and curving heads that can meet your different needs.

EASY TO USE: Wipe the plaque off the blade with a tissue, then repeat until there is no more plaque on the blade. With this amazing stainless steel dental scarper, you can easily remove those food residues located in the teeth without causing pain and damaging tooth enamel.

CONVENIENT TO BRING AWAY: Thanks to the small size, you can bring it anywhere as personal belongings, no matter where you travel to make you always keep a fresh tone.

100% SATISFACTION GUARANTEE: Customer satisfaction is our utmost priority.

**WI-70002387****3 IUD Extractor Hook 10 1/4" (26cm)****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pieces IUD Removal Hook with Double Extractor 10-1/4" (26.0cm) Stainless Steel Obstetrical & Gynecology Surgical Instruments.

This is a specialized instrument used for the removal of an intrauterine device (I.U.D). It is a narrow instrument with a slight angle in the shaft and two crochet style hooks on each side of the tip used to grasp the IUD strings. The IUD hook is inserted into the cervical canal, and an effort is made to hook the strings and pull them into the vagina. The long, slender handle allows for easy access to the vaginal canal.

IUD EXTRACTOR HOOK FEATURES

- Dual grasping hooks for better I.U.D. retrieval.
- Slightly angle design tip.

- Long shaft for deep access.
- Rounded tip for easy cervical insertion.
- Ergonomic design for superior control.

PACKAGE INCLUDE:

- 3 Pcs IUD Extractor Hook 10 1/4" (26cm)

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Extractor – Gynecological & Obstetrical Instruments

APPLICATION: IUD Removal

TIP CONFIGURATION: Hook Point

TIP STYLE: Dual Grasping Hooks

TIP: Blunt



CURVATURE: Slightly Angled

HANDLE: Slender Handle

REUSABLE: Yes

AUTOCLOAVED: Yes

RUST PREVENTION PROCEDURE:



**WI-70002064**

3 Periodontal Probe CP15 UNC Color Coded Markings Single Ended

\$19.99 per pieceTrade price from **\$18.39**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Probe UNC15 Color Coded Marking 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 Single Ended with Hollow Handle Periodontal Diagnostic Instruments

Dental Probe CP-15 Periodontal Probe UNC15 color coded markings are used to assess pocket depths, periodontal attachment levels, and gingival bleeding. It is used to check the depth of periodontal pocket vertically or horizontally. Color markings help determine the depth easily. Accurate & smooth black markings finish, which will not chip, flake or fade.

UNC15 PROBE FEATURES

- Single Ended.
- CP-15 UNC probe markings: 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 on one end.

- Designed to assess periodontal pocket depths.
- Working end has calibrated markings in millimeters
- Made of high-quality stainless steel.
- Ergonomic Hollow Handle.


PACKAGE INCLUDE:

- 3 Pcs Periodontal Probe UNC15 Color Coded Marking 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 Single Ended with Hollow Handle

TECHNICAL SPECIFICATIONS

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Dental Probe**APPLICATION:** Use to Measure Periodontal Pocket Depths & Attachment Levels**COLOR:** Silver Color**TIP:** Single Ended

PERIODONTAL PROBE CP15

**WI-70002069**

3 Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended

\$19.99 per pieceTrade price from **\$18.39**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Dental Williams Probe Color Coded 1-2-3-5-7-8-9-10 Marking Root Depth Measurement Periodontal Examination Instruments

Williams Probe is also called the periodontal probe it is the first probe to be used to assess probing depth of a pocket with color codes. Williams Probe used to assess periodontal pocket depths, attachment levels, anatomy configurations and gingival bleeding. It has the following markings – 1,2,3,5,7,8,9,10 mm. It measures the sulcus or pocket depth of the periodontium of each tooth. The working end has calibrated markings in millimeters.

WILLIAMS PROBE FEATURES

- Single Ended.

- Williams probe markings: 1-2-3-5-7-8-9-10 on one end.
- Designed to assess periodontal pocket depths.
- Working end has calibrated markings in millimeters
- Made of high-quality stainless steel.
- Ergonomic Hollow Handle.


PACKAGE INCLUDE:

- 3 Pcs Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended with Hollow Handle

TECHNICAL SPECIFICATIONS

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Dental Probe**APPLICATION:** Use to Measure Periodontal Pocket Depths & Attachment Levels**COLOR:** Silver Color

PERIODONTAL PROBE WILLIAMS



**WI-70001020****4 Joseph Skin Hooks 6.25" Double Sharp Prong 2, 5, 7 & 10mm****\$29.99** per pieceTrade price from **\$27.59**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 4 Joseph Skin Hooks 6.25" (15.9cm) 2mm, 5mm, 7mm and 10mm Wide Double Sharp Prong Retractor Premium Quality Stainless Steel Plastic Dermal ENT Surgical Instrument.

Joseph skin hook is a handheld retractor commonly used for holding back the edges of skin during intranasal and pharyngeal procedures. This is also known as "Joseph Mucosa Hook" which surgeons can use to retract and hold flaps of tissue, in order to access deeper cavities within the surgical field during ENT, plastic and general surgery procedures.

Joseph double hook feature two sharp hook-shaped tips that surgeons can embed in the desired tissue and easily retract a flap to reveal hidden structures. The handle is often round and knurled and it tapers to the hook ends that permits to

reach deep cavities. Joseph skin hook manufactured with premium grade stainless steel to ensure durability and product performance during surgical procedures.

JOSEPH DOUBLE PRONG SKIN HOOK FEATURES

- JOSEPH Hook 6.25" (15.9cm) Overall Length.
- 2mm – 10mm Wide Working Ends
- Double Sharp Hook Shaped Tips.
- Tapered Shaft Ensuring Deep Surgical Access.
- Non-Slip Knurled Handle for Maximum Control.
- Reusable.
- High Degree of Aesthetic and Corrosion Resistance.

**WI-70002159****5 Dental Probe Explorer #5 – Double Ended****\$12.49** per pieceTrade price from **\$11.49**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Dental Probe 17/23 Double Ended Diagnostic Examination Dental Instrument

The Explorer 5 is one of the most popular explorers for caries detection. The Explorer #5 designed for caries and calculus detection, exploration of pocket characteristics, furcations, and restorations. This dual-sided explorer features a Shepards hook #23 end for caries detection and a #17 end for anterior interproximal caries detection. Featuring the most ergonomic handle with knurling to the instrument shank to provide increased comfort with less hand fatigue.

DENTAL EXPLORER PROBE 17/23 FEATURES



THE FIRST QUALITY CARE TOOL: Our dental tartar scraper is made of 100%

high-purity medical stainless steel metal, which allows rust resistance, durability and safety for the tool. It also has double-ended design, cambered and curving heads that can meet your different needs.

EASY TO USE: Wipe the plaque off the blade with a tissue, then repeat until there is no more plaque on the blade. With this amazing stainless steel dental scraper, you can easily remove those food residues located in the teeth without causing pain and damaging tooth enamel.

CONVENIENT TO BRING AWAY: Thanks to the small size, you can bring it anywhere as personal belongings, no matter where you travel to make you always keep a fresh tone.

100% SATISFACTION GUARANTEE: Customer satisfaction is our utmost priority.

**WI-70002386****5 IUD Extractor Hook 10 1/4" (26cm)****\$37.99** per pieceTrade price from **\$34.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pieces IUD Removal Hook with Double Extractor 10-1/4" (26.0cm)
Stainless Steel Obstetrical & Gynecology Surgical Instruments.

This is a specialized instrument used for the removal of an intrauterine device (I.U.D.). It is a narrow instrument with a slight angle in the shaft and two crochet style hooks on each side of the tip used to grasp the IUD strings. The IUD hook is inserted into the cervical canal, and an effort is made to hook the strings and pull them into the vagina. The long, slender handle allows for easy access to the vaginal canal.

IUD EXTRACTOR HOOK FEATURES

- Dual grasping hooks for better I.U.D. retrieval.
- Slightly angle design tip.


- Long shaft for deep access.
- Rounded tip for easy cervical insertion.
- Ergonomic design for superior control.

PACKAGE INCLUDE:

- 5 Pcs IUD Extractor Hook 10 1/4" (26cm)

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Extractor –

Gynecological & Obstetrical Instruments

APPLICATION: IUD Removal**TIP CONFIGURATION:** Hook Point**TIP STYLE:** Dual Grasping Hooks**TIP:** Blunt**CURVATURE:** Slightly Angled**HANDLE:** Slender Handle**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:****WI-70002063****5 Periodontal Probe CP15 UNC Color Coded Markings Single Ended****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Dental Probe UNC15 Color Coded Marking 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 Single Ended with Hollow Handle Periodontal Diagnostic Instruments

Dental Probe CP-15 Periodontal Probe UNC15 color coded markings are used to assess pocket depths, periodontal attachment levels, and gingival bleeding. It is used to check the depth of periodontal pocket vertically or horizontally. Color markings help determine the depth easily. Accurate & smooth black markings finish, which will not chip, flake or fade.

UNC15 PROBE FEATURES

- Single Ended.
- CP-15 UNC probe markings: 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 on one end.

- Designed to assess periodontal pocket depths.
- Working end has calibrated markings in millimeters
- Made of high-quality stainless steel.
- Ergonomic Hollow Handle.

PACKAGE INCLUDE:

- 5 Pcs Periodontal Probe UNC15 Color Coded Marking 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 Single Ended with Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Probe**APPLICATION:** Use to Measure

Periodontal Pocket Depths & Attachment Levels

COLOR: Silver Color**TIP:** Single Ended

**WI-70002068****5 Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended****\$29.99** per pieceTrade price from **\$27.59**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Dental Williams Probe Color Coded 1-2-3-5-7-8-9-10 Marking Root Depth Measurement Periodontal Examination Instruments

Williams Probe is also called the periodontal probe it is the first probe to be used to assess probing depth of a pocket with color codes. Williams Probe used to assess periodontal pocket depths, attachment levels, anatomy configurations and gingival bleeding. It has the following markings – 1,2,3,5,7,8,9,10 mm. It measures the sulcus or pocket depth of the periodontium of each tooth. The working end has calibrated markings in millimeters.

WILLIAMS PROBE FEATURES

- Single Ended.


- Williams probe markings: 1-2-3-5-7-8-9-10 on one end.
- Designed to assess periodontal pocket depths.
- Working end has calibrated markings in millimeters
- Made of high-quality stainless steel.
- Ergonomic Hollow Handle.

PACKAGE INCLUDE:

- 5 Pcs Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended with Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Probe**APPLICATION:** Use to Measure Periodontal Pocket Depths & Attachment Levels**COLOR:** Silver Color

PERIODONTAL PROBE WILLIAMS

**50
PCS****WI-70002156****50 Dental Probe Explorer #5 – Double Ended****\$57.49** per pieceTrade price from **\$52.89**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 50 Dental Explorer Probe #5 (17/23) Double Ended Premium Quality Stainless Steel Diagnostic Examination Dental Instruments

The Explorer 5 is one of the most popular explorers for caries detection. The explore #5 designed for caries and calculus detection, exploration of pocket characteristics, furcations, and restorations. This dual-sided explorer features a Shepards hook #23 end for caries detection and a #17 end for anterior interproximal caries detection. Featuring the most ergonomic handle with knurling to the instrument shank to provide increased comfort with less hand fatigue.


DENTAL EXPLORER PROBE 17/23 FEATURES

THE FIRST QUALITY CARE TOOL: Our dental tarter scraper is made of 100% high-purity medical stainless steel metal, which allows rust resistance, durability and safety for the tool. It also has double-ended design, cambered and curving heads that can meet your different needs.

EASY TO USE: Wipe the plaque off the blade with a tissue, then repeat until there is no more plaque on the blade. With this amazing stainless steel dental scarper, you can easily remove those food residues located in the teeth without causing pain and damaging tooth enamel.

CONVENIENT TO BRING AWAY: Thanks to the small size, you can bring it anywhere as personal belongings, no matter where you travel to make you always keep a fresh tone.

100% SATISFACTION GUARANTEE:

**WI-70000808****50 Teasing Needles 5.5" Bent Tips with Plastic Handle****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

50 Pcs Bent Teasing Needles Set 5.5" (14cm) with Plastic Handle – Medical Lab Dissecting Needle Stainless Steel Single Ended Micro Fine Precision Probe Teasing Needles Lab Biology Experiment Equipment for Medical Students.

Teasing needles are used to manipulate small samples. These are long dissection needle with plastic handles and are available with straight and bent tips. Dissecting teasing needle is a multi-functional tool that can be used for dissecting, sampling, cutting, puncturing and other operations. With its sharp tip, it can easily penetrate tissues made of different materials and can also be used for sampling and cutting tissues. With very fine sharp tips which can precisely penetrate cells and tissues for fine operations teasing needles come with




comfortable to hold plastic handle design which reduces hand fatigue and improve operation precision.

DISSECTING TEASING NEEDLES FEATURES

- 5.5" (14cm) Overall Length.
- Sharp-pointed, Bent Needles.
- Single Ended Stainless Steel Tips.
- Plastic Handle.
- Reusable.
- Easy to Clean.
- Perfect for Use in Dissection, Biology Lab and Tissue Manipulation.
- Made with Stainless Steel Tips.

PACKAGE INCLUDE

- 50 Pcs Dissection Teasing Needles Set Bent Tips with Plastic Handle

TECHNICAL INFORMATION**WI-70002163****6 Rainbow Color Dental Explorer Probe #5 – Double Ended****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pieces Multi Rainbow Color Dental Probe 17/23 Double Ended Diagnostic Examination Dental Instrument

Dental Explorer Probe is a professional dental scraper tool helps you to remove stains, tartar, and debris from teeth and gums. Regular care of teeth can reduce the generation of bad breath. This hygiene tool is also used for keeping your pet clean and oral fresh. You can be comfortable to be closer to your pet (cat, dog) without any fears about your health problems.

DENTAL EXPLORER PROBE 17/23 FEATURES

THE FIRST QUALITY CARE TOOL: Our dental tarter scraper is made of 100% high-purity medical stainless steel metal, which allows rust resistance, durability and

safety for the tool. It also has double-ended design, cambered and curving heads that can meet your different needs.

EASY TO USE: Wipe the plaque off the blade with a tissue, then repeat until there is no more plaque on the blade. With this amazing stainless steel dental scarper, you can easily remove those food residues located in the teeth without causing pain and damaging tooth enamel.

CONVENIENT TO BRING AWAY: Thanks to the small size, you can bring it anywhere as personal belongings, no matter where you travel to make you always keep a fresh tone.

100% SATISFACTION GUARANTEE: Customer satisfaction is our utmost priority.

PACKAGE INCLUDE:

**WI-36793611**

Autopsy Hammer with Hook Stainless Steel Post Mortem Examination

\$28.00 per pieceTrade price from **\$20.00**
Minimum order 14 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Autopsy Hammer with Hook Stainless Steel Post Mortem Examination is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type retract, hold or pass suture in fine and delicate work.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip

alignment and surface finish are checked on each unit, not on a sample. Listed at \$28.00, falling to \$20.00 at 14 pieces — about 29% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000904**

Babinski Percussion Reflex Hammer 9" with Needle

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Babinski Reflex Hammer 9" (22.9cm) Short Metal Handle with Needle Diagnostic Reflex Percussion for Superficial Responses in Neurological Examinations

The Babinski Reflex Hammer features an adjustable head that can be pivoted in two directions and a pointed bottom. Ideal for hospitals and clinics, the Babinski Percussion Hammer is used for not only for diagnosis but for tracking the progress of therapy. The Babinsky reflex hammer has concealed needle which allow for additional reflex and neurological testing. The needles check for pain sensitivity. The hammer induces a reflex which can then be assessed to come up with a diagnosis. Babinski hammer is made of high-quality stainless steel to ensure durability and longevity.

BABINSKY PERCUSSION HAMMER FEATURES

- Stainless Steel Short Metal Handle.
- 9" (22.9cm) Long with Pointed Base.
- Comes with Needle to Check for Pain Sensitivity.
- Rubber Head Can Be Used on Horizontal and Vertical Postions.
- Reusable.

PACKAGE INCLUDE

- 1 Pc Babinski Reflex Hammer 9" with Needle

TECHNICAL INFORMATION

MATERIAL: Brass with Chrome Plated
INSTRUMENTS TYPE: Percussion Hammers

APPLICATION: Neurological, Medical, Diagnostic

HEAD: Rubber Head



**WI-70000903**

Buck Percussion Reflex Hammer with Needle and Brush

\$14.99 per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Buck Reflex Hammer with Needle and Brush, Dual Rubber Head with Stainless Steel Handle for Neurologic Reflex Diagnosis and Testing.

Buck Percussion Reflex Hammer is a versatile diagnostic instrument designed to evaluate myotatic and cutaneous responses in patients. In addition to testing deep tendon reflexes the Neurological Buck Reflex Hammer can assist in the detection of abnormalities of the central or peripheral nervous system. It can also be used for chest percussion.

The double head of the Buck Percussion Reflex Hammer consists of two silicone grade rubber tips of different sizes. This particular 3 in 1 Buck Reflex Hammer model with Needle and Brush, can spin out the needle at the top, a brush can spin out at the tail. Our reflex hammer features

a sturdy stainless steel handle; it can be sterilized easily for physical examinations. It is designed in a way that allows you to hold it with a firm and comfortable grip.

BUCK PERCUSSION REFLEX HAMMER FEATURES

- 3 in 1 Reflex Hammer Paired With Pin & Brush.
- Dual Rubber Head of Different Sizes.
- Ergonomic & Lightweight Stainless Steel Handle.
- Can Be Sterilized and Reusable.
- Firm and Comfortable Grip Handle.
- For Neurologiv Reflex Diagnosis and Testing.
- Deliver a Soft Bag Storage and Protect the Reflex Hammer.

PACKAGE INCLUDE

**WI-28935429**

Color Hair Extension Pliers Pulling Hook Tool

\$15.00 per pieceTrade price from **\$12.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Color Hair Extension Pliers Pulling Hook Tool is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type retract, hold or pass suture in fine and delicate work.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$15.00, falling to \$12.00 at 1 set — about 20% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000889****Dejerine Percussion Hammer 8" without Needle****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dejerine Percussion Hammer without Needle 8" (20.3cm) Dual Rubber Head with Stainless Steel Handle for Neurological Reflex Diagnosis and Testing.

Dejerine hammer is ideally used for testing of deep tendon reflexes in a diagnostic evaluation. The dejerine percussion reflex hammer features a double-sided head composed of two rounded rubber surfaces for increased patient comfort. This hammer includes a slim neck that flares into a wider handle for improved gripping. Our reflex hammer made with premium quality materials ideal for clinics and hospitals. Dejerine hammer can be sterilized easily for physical examinations.

DEJERINE PERCUSSION HAMMER FEATURES

- Reflex Hammer without Pin Model.
- Dual Rubber Head of Different Sizes.
- Ergonomic Designed Flare Handle.
- Can Be Sterilized and Reusable.
- Firm and Comfortable Grip Wide Handle.
- For Neurologiv Reflex Diagnosis and Testing.

PACKAGE INCLUDE

- 1 Pc Dejerine Percussion Hammer 8" without Needle

TECHNICAL INFORMATION

MATERIAL: Brass with Chrome Plated
INSTRUMENTS TYPE: Percussion Hammers

APPLICATION: Neurological, Medical, Diagnostic

HEAD: Dual Rubber Head

HANDLE: Solid Handle

**WI-70000890****Dejerine Percussion Hammer 9" with Needle****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dejerine Percussion Hammer with Needle, Dual Rubber Head with Stainless Steel Flare Handle for Neurologic Reflex Diagnosis and Testing.

Dejerine hammer is ideally used for testing of deep tendon reflexes in a diagnostic evaluation. The dejerine percussion reflex hammer features a double-sided head composed of two rounded rubber surfaces for increased patient comfort. This hammer includes a slim neck that flares into a wider handle for improved gripping. Our reflex hammer made with premium quality materials ideal for clinics and hospitals. Dejerine hammer can be sterilized easily for physical examinations.

DEJERINE PERCUSSION HAMMER FEATURES

- Reflex Hammer Paired With Pin.
- Dual Rubber Head of Different Sizes.
- Ergonomic Designed Flare Handle.
- Can Be Sterilized and Reusable.
- Firm and Comfortable Grip Handle.
- For Neurologiv Reflex Diagnosis and Testing.

PACKAGE INCLUDE

- 1 Pc Dejerine Percussion Reflex Hammer 9" with Needle

TECHNICAL INFORMATION


MATERIAL: Brass with Chrome Plated
INSTRUMENTS TYPE: Percussion Hammers

APPLICATION: Neurological, Medical, Diagnostic

HEAD: Dual Rubber Head

HANDLE: Solid Handle

REUSABLE: Yes

**WI-70002161****Dental Probe Explorer #5 – Double Ended****\$7.99** per pieceTrade price from **\$7.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Probe 17/23 Double Ended
Diagnostic Examination Dental Instrument

The Explorer 5 is also known as Shepherd's Hook. One of the most popular explorers for caries detection. The explore #5 designed for caries and calculus detection, exploration of pocket characteristics, furcations, and restorations. This dual-sided explorer features a Shepards hook #23 end for caries detection and a #17 end for anterior interproximal caries detection. Featuring the most ergonomic handle with knurling to the instrument shank to provide increased comfort with less hand fatigue.

DENTAL EXPLORER PROBE 17/23 FEATURES




THE FIRST QUALITY CARE TOOL: Our dental tarter scraper is made of 100%

high-purity medical stainless steel metal, which allows rust resistance, durability and safety for the tool. It also has double-ended design, cambered and curving heads that can meet your different needs.

EASY TO USE: Wipe the plaque off the blade with a tissue, then repeat until there is no more plaque on the blade. With this amazing stainless steel dental scarper, you can easily remove those food residues located in the teeth without causing pain and damaging tooth enamel.

CONVENIENT TO BRING AWAY: Thanks to the small size, you can bring it anywhere as personal belongings, no matter where you travel to make you always keep a fresh tone.

100% SATISFACTION GUARANTEE: Customer satisfaction is our utmost priority.

**WI-70003555****Double Ended Pratt Uterine Dilator 13/15 FR Curved****\$10.30** per pieceTrade price from **\$9.48**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pratt Uterine Dilator Curved 11.5" (29.2cm)
Double Ended 13/15 FR with Ultra Smooth Finish Premium Quality Stainless Steel
OB/Gynecology Surgical Instrument.

Pratt Uterine Dilator is an incredible OB/GYN instrument and is widely used to measure dilation or facilitate further uterus dilation. The Double-ended Pratt Uterine Dilator is Straight in the middle, slightly curved at each end, with a smooth, rounded tip. The size of these dilators is determined by the circumference of the tip and measured in French units, with one French unit equaling 1/3 mm. Suitable for D&C (Dilation and Curettage), D&E (Dilation and Evacuation), and Other Intrauterine Interventions. These are made with high-grade surgical stainless steel and can be reused after sterilization.

PRATT UTERINE DILATOR FEATURES


- 11.5 Inches (29.2cm) Overall Length.
- Double-ended Metal Rod With Two Graduated Sizes for Easy Switching by Surgeons.
- Smooth, Rounded Tips That Gently Taper to Reduce Patient Discomfort and Tissue Trauma.
- Sizes Are in French (Fr) Units, Where 1 Fr Equals 1/3 of a Millimeter.
- Solid Central Gripping Platform Offers a Strong and Secure Hold.
- Fully Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

1 Pc Pratt Uterine Dilator Curved 11.5"
Double Ended 13/15 FR – Solid

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

**WI-70003554****Double Ended Pratt Uterine Dilator 17/19 FR Curved****\$10.30** per pieceTrade price from **\$9.48**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking




Pratt Uterine Dilator Curved 11.5" (29.2cm) Double Ended 17/19 FR with Ultra Smooth Finish Premium Quality Stainless Steel OB/Gynecology Surgical Instrument.

Pratt Uterine Dilator is an incredible OB/GYN instrument and is widely used to measure dilation or facilitate further uterus dilation. The Double-ended Pratt Uterine Dilator is Straight in the middle, slightly curved at each end, with a smooth, rounded tip. The size of these dilators is determined by the circumference of the tip and measured in French units, with one French unit equaling 1/3 mm. Suitable for D&C (Dilation and Curettage), D&E (Dilation and Evacuation), and Other Intrauterine Interventions. These are made with high-grade surgical stainless steel and can be reused after sterilization.

- 11.5 Inches (29.2cm) Overall Length.
- Double-ended Metal Rod With Two Graduated Sizes for Easy Switching by Surgeons.
- Smooth, Rounded Tips That Gently Taper to Reduce Patient Discomfort and Tissue Trauma.
- Sizes Are in French (Fr) Units, Where 1 Fr Equals 1/3 of a Millimeter.
- Solid Central Gripping Platform Offers a Strong and Secure Hold.
- Fully Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

1 Pc Pratt Uterine Dilator Curved 11.5" Double Ended 17/19 FR – Solid

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**PRATT UTERINE DILATOR FEATURES****WI-70003553****Double Ended Pratt Uterine Dilator 21/23 FR Curved****\$15.46** per pieceTrade price from **\$14.22**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Pratt Uterine Dilator Curved 11.5" (29.2cm) Double Ended 21/23 FR with Ultra Smooth Finish Premium Quality Stainless Steel OB/Gynecology Surgical Instrument.


Pratt Uterine Dilator is an incredible OB/GYN instrument and is widely used to measure dilation or facilitate further uterus dilation. The Double-ended Pratt Uterine Dilator is Straight in the middle, slightly curved at each end, with a smooth, rounded tip. The size of these dilators is determined by the circumference of the tip and measured in French units, with one French unit equaling 1/3 mm. Suitable for D&C (Dilation and Curettage), D&E (Dilation and Evacuation), and Other Intrauterine Interventions. These are made with high-grade surgical stainless steel and can be reused after sterilization.

- 11.5 Inches (29.2cm) Overall Length.
- Double-ended Metal Rod With Two Graduated Sizes for Easy Switching by Surgeons.
- Smooth, Rounded Tips That Gently Taper to Reduce Patient Discomfort and Tissue Trauma.
- Sizes Are in French (Fr) Units, Where 1 Fr Equals 1/3 of a Millimeter.
- Solid Central Gripping Platform Offers a Strong and Secure Hold.
- Fully Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

1 Pc Pratt Uterine Dilator Curved 11.5" Double Ended 21/23 FR – Solid

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**PRATT UTERINE DILATOR FEATURES**

**WI-70003552****Double Ended Pratt Uterine Dilator 25/27 FR Curved****\$14.43** per pieceTrade price from **\$13.28**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Pratt Uterine Dilator Curved 11.5" (29.2cm) Double Ended 25/27 FR with Ultra Smooth Finish Premium Quality Stainless Steel OB/Gynecology Surgical Instrument.

Pratt Uterine Dilator is an incredible OB/GYN instrument and is widely used to measure dilation or facilitate further uterus dilation. The Double-ended Pratt Uterine Dilator is Straight in the middle, slightly curved at each end, with a smooth, rounded tip. The size of these dilators is determined by the circumference of the tip and measured in French units, with one French unit equaling 1/3 mm. Suitable for D&C (Dilation and Curettage), D&E (Dilation and Evacuation), and Other Intrauterine Interventions. These are made with high-grade surgical stainless steel and can be reused after sterilization.

- 11.5 Inches (29.2cm) Overall Length.
- Double-ended Metal Rod With Two Graduated Sizes for Easy Switching by Surgeons.
- Smooth, Rounded Tips That Gently Taper to Reduce Patient Discomfort and Tissue Trauma.
- Sizes Are in French (Fr) Units, Where 1 Fr Equals 1/3 of a Millimeter.
- Solid Central Gripping Platform Offers a Strong and Secure Hold.
- Fully Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

1 Pc Pratt Uterine Dilator Curved 11.5" Double Ended 25/27 FR – Solid

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**PRATT UTERINE DILATOR FEATURES****WI-70003551****Double Ended Pratt Uterine Dilator 29/31 FR Curved****\$14.43** per pieceTrade price from **\$13.28**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pratt Uterine Dilator Curved 11.5" (29.2cm) Double Ended 29/31 FR with Ultra Smooth Finish Premium Quality Stainless Steel OB/Gynecology Surgical Instrument.


Pratt Uterine Dilator is an incredible OB/GYN instrument and is widely used to measure dilation or facilitate further uterus dilation. The Double-ended Pratt Uterine Dilator is Straight in the middle, slightly curved at each end, with a smooth, rounded tip. The size of these dilators is determined by the circumference of the tip and measured in French units, with one French unit equaling 1/3 mm. Suitable for D&C (Dilation and Curettage), D&E (Dilation and Evacuation), and Other Intrauterine Interventions. These are made with high-grade surgical stainless steel and can be reused after sterilization.

- 11.5 Inches (29.2cm) Overall Length.
- Double-ended Metal Rod With Two Graduated Sizes for Easy Switching by Surgeons.
- Smooth, Rounded Tips That Gently Taper to Reduce Patient Discomfort and Tissue Trauma.
- Sizes Are in French (Fr) Units, Where 1 Fr Equals 1/3 of a Millimeter.
- Solid Central Gripping Platform Offers a Strong and Secure Hold.
- Fully Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

1 Pc Pratt Uterine Dilator Curved 11.5" Double Ended 29/31 FR – Solid

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**PRATT UTERINE DILATOR FEATURES**

**WI-70003550****Double Ended Pratt Uterine Dilator 33/35 FR Curved****\$14.43** per pieceTrade price from **\$13.28**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Pratt Uterine Dilator Curved 11.5" (29.2cm) Double Ended 33/35 FR with Ultra Smooth Finish Premium Quality Stainless Steel OB/Gynecology Surgical Instrument.

Pratt Uterine Dilator is an incredible OB/GYN instrument and is widely used to measure dilation or facilitate further uterus dilation. The Double-ended Pratt Uterine Dilator is Straight in the middle, slightly curved at each end, with a smooth, rounded tip. The size of these dilators is determined by the circumference of the tip and measured in French units, with one French unit equaling 1/3 mm. Suitable for D&C (Dilation and Curettage), D&E (Dilation and Evacuation), and Other Intrauterine Interventions. These are made with high-grade surgical stainless steel and can be reused after sterilization.

- 11.5 Inches (29.2cm) Overall Length.
- Double-ended Metal Rod With Two Graduated Sizes for Easy Switching by Surgeons.
- Smooth, Rounded Tips That Gently Taper to Reduce Patient Discomfort and Tissue Trauma.
- Sizes Are in French (Fr) Units, Where 1 Fr Equals 1/3 of a Millimeter.
- Solid Central Gripping Platform Offers a Strong and Secure Hold.
- Fully Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

1 Pc Pratt Uterine Dilator Curved 11.5" Double Ended 33/35 FR – Solid

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**PRATT UTERINE DILATOR FEATURES****WI-70003549****Double Ended Pratt Uterine Dilator 37/39 FR Curved****\$14.43** per pieceTrade price from **\$13.28**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pratt Uterine Dilator Curved 11.5" (29.2cm) Double Ended 37/39 FR with Ultra Smooth Finish Premium Quality Stainless Steel OB/Gynecology Surgical Instrument.

Pratt Uterine Dilator is an incredible OB/GYN instrument and is widely used to measure dilation or facilitate further uterus dilation. The Double-ended Pratt Uterine Dilator is Straight in the middle, slightly curved at each end, with a smooth, rounded tip. The size of these dilators is determined by the circumference of the tip and measured in French units, with one French unit equaling 1/3 mm. Suitable for D&C (Dilation and Curettage), D&E (Dilation and Evacuation), and Other Intrauterine Interventions. These are made with high-grade surgical stainless steel and can be reused after sterilization.

- 11.5 Inches (29.2cm) Overall Length.
- Double-ended Metal Rod With Two Graduated Sizes for Easy Switching by Surgeons.
- Smooth, Rounded Tips That Gently Taper to Reduce Patient Discomfort and Tissue Trauma.
- Sizes Are in French (Fr) Units, Where 1 Fr Equals 1/3 of a Millimeter.
- Solid Central Gripping Platform Offers a Strong and Secure Hold.
- Fully Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

1 Pc Pratt Uterine Dilator Curved 11.5" Double Ended 37/39 FR – Solid

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**PRATT UTERINE DILATOR FEATURES**

**WI-70003548****Double Ended Pratt Uterine Dilator 41/43 FR Curved****\$14.43** per pieceTrade price from **\$13.28**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Pratt Uterine Dilator Curved 11.5" (29.2cm) Double Ended 41/43 FR with Ultra Smooth Finish Premium Quality Stainless Steel OB/Gynecology Surgical Instrument.

Pratt Uterine Dilator is an incredible OB/GYN instrument and is widely used to measure dilation or facilitate further uterus dilation. The Double-ended Pratt Uterine Dilator is straight in the middle, slightly curved at each end, with a smooth, rounded tip. The size of these dilators is determined by the circumference of the tip and measured in French units, with one French unit equaling 1/3 mm. Suitable for D&C (Dilation and Curettage), D&E (Dilation and Evacuation), and Other Intrauterine Interventions. These are made with high-grade surgical stainless steel and can be reused after sterilization.

- 11.5 Inches (29.2cm) Overall Length.
- Double-ended Metal Rod With Two Graduated Sizes for Easy Switching by Surgeons.
- Smooth, Rounded Tips That Gently Taper to Reduce Patient Discomfort and Tissue Trauma.
- Sizes Are in French (Fr) Units, Where 1 Fr Equals 1/3 of a Millimeter.
- Solid Central Gripping Platform Offers a Strong and Secure Hold.
- Fully Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

1 Pc Pratt Uterine Dilator Curved 11.5" Double Ended 41/43 FR – Solid

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**PRATT UTERINE DILATOR FEATURES****WI-70000297****Explorer Probe #23 Single Ended Hollow Handle****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Explorer Probe #23 Shepherd Single Ended Hollow Handle Dental Diagnostic and Examination Tool for Minimally Invasive Procedures.

Shepherd single ended explorer #23 is a dental diagnostic tool designed to detect decay and calculus, explore tooth surfaces and examine restorations. Dental explorer probe #23 features 0.2mm fine point, single-ended design with uniquely shepherds hook which enhances access in small, hard-to-reach areas for tactile assessment. The ergonomic hollow handle provides increased comfort with less hand fatigue. Made of 100% high-quality stainless steel, ensuring rust resistance, durability, and safety for repeated clinical use.

- Single Ended.
- 0.2mm Fine Point Precise Tip.
- Angled Shepherds Hook.
- Ergonomic Non-slip Hollow Handle.
- Reusable.
- Made With Premium Quality Stainless Steel.


PACKAGE INCLUDE:

- 1 Pc Dental Explorer #23 Single Ended – Hollow Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Instruments**APPLICATION:** Dental Examination, Restorations and Furcations**STERILE:** Non-Sterile**REUSABLE:** Yes**RUST PREVENTION PROCEDURE:**

Passivated

EXPLORER PROBE #23 FEATURES

**WI-70000296****Explorer Probe #5 Double Ended Hollow Handle****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Explorer Probe #5 (17/23) Double Ended Hollow Handle Premium Quality Stainless Steel Dental Diagnostic and Examination Dentist Instrument.


Dental explorer #5 is one of the most popular explorers for caries detection. The explore #5 designed for caries and calculus detection, exploration of pocket characteristics, furcations, and restorations. This dual-sided explorer features a Shepards hook #23 end for caries detection and a #17 end for anterior interproximal caries detection. The rounded hollow handle design makes the instrument lighter, helping to minimize operator fatigue during long procedures. Our diagnostic explorer #5 made of 100% high-quality stainless steel, ensuring rust resistance, durability, and safety for repeated clinical use.

EXPLORER #5 HOLLOW HANDLE FEATURES

- Double Ended.
- Combination of Shepherds Hook and #17 Explorer.
- Ergonomic Non-slip Rounded Hollow Handle.
- Ideal for Most Hard to Reach Areas.
- Lightweight and Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Dental Explorer #5 Double Ended Hollow Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Instruments**APPLICATION:** Dental Examination**WI-70000051****Explorer Probe Moon Shaped – Double Ended****\$7.99** per pieceTrade price from **\$7.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Explorer Probe Moon Shaped Double Ended Sharp Tips with Solid Grip Handle – Premium Quality Calculus Removal Dental Deep Cleaning Instrument.

This stainless steel dental explorer probe has double-ended moon-shaped tips for detecting calculus, caries, and irregularities. It features an ergonomic, non-slip knurled handle for precise control during examinations. The instrument is autoclavable and reusable, making it suitable for repeated clinical use, and is an essential addition to any teeth cleaning kit.

MOON-SHAPED EXPLORER PROBE FEATURES

- Double Ended Design.
- Moon-Shaped Sharp Tips.

- Ergonomic, Non-slip Knurled Handle for Precise Control.
- Used to Detect Calculus, Caries, and Irregularities.
- Fully Autoclavable and Reusable for Clinical Use.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Explorer Probe Moon-Shaped Double Ended – Solid Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Instruments**APPLICATION:** Dental Examination to

Detect Calculus, Caries, and Irregularities

AUTOCLEAVE/REUSABLE: Yes**STERILE:** Non-Sterile**RUST PREVENTION PROCEDURE:**



WI-36754729

Forensic Autopsy Hammer Hook Stainless Steel Medical Pathology Dissection

\$28.00 per piece

Trade price from **\$20.00**
Minimum order 14 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Forensic Autopsy Hammer Hook Stainless Steel Medical Pathology Dissection, made for general and specialty operating theatre use. Instruments of this type retract, hold or pass suture in fine and delicate work.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip

alignment and surface finish are checked on each unit, not on a sample. Listed at \$28.00, falling to \$20.00 at 14 pieces — about 29% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000022

Frazier Dura Hook 5" (12.7cm) Blunt Point

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Frazier Dura Hook 5" (12.7cm) Blunt Point – Premium Quality Stainless Steel Single Prong Skin Hook Retractor Surgical Instrument.


Frazier Dura Skin Hook is a specialized tool used by neurosurgeons to grasp and mobilize the cranial dura during procedures like craniotomy. It provides an atraumatic method to lift the dura, minimizing the risk of hemorrhage and injuries. FRAZIER design hook featuring blunt tip ensures a smooth and safe edge for delicate tissues. Its slim shaft enables precise surgical access in narrow cavities, and the ergonomic handle provides optimal control for surgeons during microsurgeries. Made from high-quality stainless steel, it ensures optimal performance in surgical settings.

FRAZIER DURA SKIN HOOK FEATURES

- 5 Inch (12.7cm) Overall Length.
- Blunt Point Single Prong Tip Designed for Smooth and Safe Retracting.
- Extra-rounded Hook Tip for Secure Grasping of Delicate Tissue.
- Slim Straight Shaft Facilitates Deep Surgical Access.
- Ergonomic Solid Handle Promotes Maximum Control.
- Suitable for Precise Neurosurgical and Plastic Surgery Applications.
- Sterilizable and Reusable.
- Made With High-quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Frazier Dura Hook 5" (12.7cm) Blunt Point

**WI-70000023****Frazier Dura Hook 5" (12.7cm) Sharp Point****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Frazier Dura Hook 5" (12.7cm) Sharp Point – Premium Quality Stainless Steel Single Prong Skin Hook Retractor Surgical Instrument.

Frazier Dura Skin Hook is a specialized tool used by neurosurgeons to grasp and mobilize the cranial dura during procedures like craniotomy. It provides an atraumatic method to lift the dura, minimizing the risk of hemorrhage and injuries. FRAZIER design hook featuring a sharp tip, specifically engineered for the precise holding and retraction of the dura mater and skin during neurosurgery and plastic surgery procedures. Ergonomically designed solid handle enhances grip and reduces fatigue. Made from high-quality stainless steel, it ensures optimal performance in surgical settings.

FRAZIER DURA SKIN HOOK FEATURES

- 5 Inch (12.7cm) Overall Length.
- Sharp Point Single Prong Tip Designed for Precision Work.
- Extra-rounded Hook Tip for Secure Grasping of Delicate Tissue.
- Slim Straight Shaft Facilitates Deep Surgical Access.
- Ergonomic SOLID Handle Promotes Maximum Control.
- Suitable for Precise Neurosurgical and Plastic Surgery Applications.
- Sterilizable and Reusable.
- Made With High-quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Frazier Dura Hook 5" (12.7cm) Sharp Point

**WI-70000763****Gillies Skin Hook 7" Large 4mm Dia Sharp Round Tip****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Gillies Skin Hook 7" (17.8cm) 4mm Diameter Large Sharp Round Single Prong Premium Quality Stainless Steel Decimal Retractor Surgical Instrument.

Gillies Skin Hook is an extremely useful, multipurpose surgical instrument which allows grasping, holding, and elevating flaps of the skin and other soft tissues, especially during wound exploration, biopsy harvesting, and reconstructive surgery procedures. This hand-held retractor used to retract the skin during surgical procedures such as the dissection of lesions.

Gillies Hook features a sharp rounded single prong hook, which is available in small and large sizes. The sharp-tipped design is ideal to atraumatically traverse the desired tissue. This skin hook has a long, cylindrical shaft supplied with a

tapered design to reach deep cavities with ease. Manufactured in premium grade stainless steel to ensure the highest quality and product performance during surgical procedures.

GILLIES SKIN HOOK FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- 4mm Diameter Large Sharp Hook Tip.
- Sharp Rounded Single Prong.
- Long Cylindrical Shaft Allows Deeper Surgical Access.
- Ergonomic Flat Handle Ensuring Maximum Control.
- Reusable.
- Made with Premium Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Gillies Skin Hook 7" (17.8cm) 4mm Dia Large Size

**WI-70000764**

Gillies Skin Hook 7" Small 3mm Dia Sharp Round Tip

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Gillies Skin Hook 7" (17.8cm) 3mm Diameter Sharp Round Single Prong Premium Quality Stainless Steel Retractor Surgical Instrument.

Gillies Skin Hook is an extremely useful, multipurpose surgical instrument which allows grasping, holding, and elevating flaps of the skin and other soft tissues, especially during wound exploration, biopsy harvesting, and reconstructive surgery procedures. This hand-held retractor used to retract the skin during surgical procedures such as the dissection of lesions.

Gillies Hook features a sharp rounded single prong hook, which is available in small and large sizes. The sharp-tipped design is ideal to atraumatically traverse the desired tissue. This skin hook has a long, cylindrical shaft supplied with a


tapered design to reach deep cavities with ease. Manufactured in premium grade stainless steel to ensure the highest quality and product performance during surgical procedures.

GILLIES SKIN HOOK FEATURES

- Crafted with 7" (17.8cm) Overall Length.
- 3mm Diameter Small Sharp Hook Tip.
- Sharp Rounded Single Prong.
- Long Cylindrical Shaft Allows Deeper Surgical Access.
- Ergonomic Flat Handle Ensuring Maximum Control.
- Reusable.
- Made with Premium Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Gillies Skin Hook 7" (17.8cm) 3mm Dia Small Size

**WI-70001407**

Graham Nerve Hook Blunt 6mm Deep 6.5" Neurosurgical Instruments

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Graham Nerve Hook 6.5 Inches Blunt 6mm Deep Nerve/Vessel Retractor German Grade Stainless Steel Neurosurgical Instrument.

Graham Nerve Hook is used in neurosurgical and spinal procedures that need retraction of delicate nerve tissues. It features a blunt tip that allows easy retraction and more exploration with the reduced injury. Graham nerve retractor hook is a small, handheld instrument with a long, thin, curved hook at one end features a straight shaft that tapers to the tip, ensuring optimal access into narrow cavities. Furthermore, the ergonomically designed solid handle prevents slippage and finger strain.


- 6mm Deep Tip
- Round Blunt Hook Promoting Atraumatic Retraction.
- Straight Tapering Shaft for Deep Surgical Acces.
- Ergonomic Solid Handle Ensuring Superior Control.
- Resuable.
- Made with German Grade Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Graham Nerve Hook 6.5" (16.5cm) Blunt , 6mm Deep Tip

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENT TYPE:** Nerve Hooks**APPLICATION:** Neurosurgery / Spine**TIP:** Blunt – 6mm Deep**HANDLE:** Solid Handle

**WI-70003542**

Hegar Uterine Dilator 3mm-4mm Double Ended 7.5" (19.1cm)

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 3mm-4mm Double Ended 7.5" (19.1cm) Solid Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.


Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3.0 to 30.0 mm.
- Round Solid, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70003541**

Hegar Uterine Dilator 5mm-6mm Double Ended 7.5" (19.1cm)

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 5mm-6mm Double Ended 7.5" (19.1cm) Solid Premium Quality Urethral Surgical Gynecology Diagnostic Instrument

Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3.0 to 30.0 mm.
- Round Solid, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70003546****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 11mm-12mm – Hollow****\$8.99** per pieceTrade price from **\$8.27**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 11mm-12mm Double Ended 7.5" (19.1cm) Hollow Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.


Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70003545****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 13mm-14mm – Hollow****\$8.99** per pieceTrade price from **\$8.27**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 13mm-14mm Double Ended 7.5" (19.1cm) Hollow Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.



Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available



individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.



**WI-70003544****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 15mm-16mm – Hollow****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 15mm-16mm Double Ended 7.5" (19.1cm) Hollow Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.


Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70003543****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 17mm-18mm – Hollow****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 17mm-18mm Double Ended 7.5" (19.1cm) Hollow Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.



Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70000773****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 19mm-20mm – Hollow****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 19mm-20mm Double Ended 7.5" (19.1cm) Hollow Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.



Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70000772****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 21mm-22mm – Hollow****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 21mm-22mm Double Ended 7.5" (19.1cm) Hollow Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.



Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70000771****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 23mm-24mm – Hollow****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 23mm-24mm Double Ended 7.5" (19.1cm) Hollow Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.



Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70000770****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 25mm-26mm – Hollow****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 25mm-26mm Double Ended 7.5" (19.1cm) Hollow – Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.



Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70000769****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 27mm-28mm – Hollow****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 27mm-28mm Double Ended 7.5" (19.1cm) Hollow – Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.



Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70000472****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 29mm-30mm – Hollow****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 29mm-30mm Double Ended 7.5" (19.1cm) Hollow – Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.


Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70003540****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 7mm-8mm – Hollow****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 7mm-8mm Double Ended 7.5" (19.1cm) Hollow Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.


Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.

**WI-70003539****Hegar Uterine Dilator 7.5" (19.1cm) Double Ended 9mm-10mm – Hollow****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hegar Uterine Sounds Dilator 9mm-10mm Double Ended 7.5" (19.1cm) Hollow Premium Quality Urethral Surgical Gynecology Diagnostic Instruments.



Hegar sounds are frequently used for gently dilating the mouth of the cervix during numerous gynaecological and obstetric procedures. This uterine dilator used in outpatient procedures such as intrauterine contraceptive device insertion as well as surgeries such as dilatation and evacuation, dilation and curettage and uterine biopsy. Hegar double-ended uterine sound dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. Uterine dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. The double ended hegar dilators available

individually as well as in 8 pcs and 14 pcs sounding kit.

We manufacture hegar dilators with premium quality surgical grade stainless steel with varying thickness ranging from 3-4mm to 29-30mm . Hegar uterine sound varying can be sterilised and reused.

HEGAR SOUNDS DILATOR FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Available in Single/ Double Ended Options.
- Thickness Ranging from 3mm to 30mm.
- Round Hollow, Slightly Curved & Non-Flexible.
- Mirror Polished Finished.
- Autoclaved and Reusable.



WI-70002388

IUD Extractor Hook 10 1/4" (26cm)

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

IUD Removal Hook with Double Extractor 10-1/4" (26.0cm) Stainless Steel Obstetrical & Gynecology Surgical Instruments.

This is a specialized instrument used for the removal of an intrauterine device (I.U.D). It is a narrow instrument with a slight angle in the shaft and two crochet style hooks on each side of the tip used to grasp the IUD strings. The IUD hook is inserted into the cervical canal, and an effort is made to hook the strings and pull them into the vagina. The long, slender handle allows for easy access to the vaginal canal.

IUD EXTRACTOR HOOK FEATURES

- Dual grasping hooks for better I.U.D. retrieval.
- Slightly angle design tip.
- Long shaft for deep access.

- Rounded tip for easy cervical insertion.
- Ergonomic design for superior control.

PACKAGE INCLUDE:

- IUD Extractor Hook 10 1/4" (26cm)

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Extractor – Gynecological & Obstetrical Instruments
APPLICATION: IUD Removal
TIP CONFIGURATION: Hook Point
TIP STYLE: Dual Grasping Hooks
TIP: Blunt
CURVATURE: Slightly Angled
HANDLE: Slender Handle
REUSABLE: Yes
AUTOCCLAVED: Yes
RUST PREVENTION PROCEDURE: Passivated
ULTRASONIC CLEANED: Yes



WI-43635731

Jameson Large Muscle Hook Precision Eye Surgery Tool

\$15.00 per piece

Trade price from **\$9.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Jameson Large Muscle Hook Precision Eye Surgery Tool is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type retract, hold or pass suture in fine and delicate work.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$15.00, falling to \$9.00 at 10 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-36409088****Jameson Muscle Hook Large Eye Surgery Instrument****\$15.00** per pieceTrade price from **\$9.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Jameson Muscle Hook Large Eye Surgery Instrument for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type retract, hold or pass suture in fine and delicate work.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$15.00, falling to \$9.00 at 10 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001026****Joseph Skin Hook 6.25" Single Sharp Prong****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Joseph Skin Hook 6.25" (15.9cm) Single Sharp Prong Retractor Premium Quality Stainless Steel Plastic Dermal ENT Surgical Instrument.

Joseph skin hook is a handheld retractor commonly used for holding back the edges of skin during intranasal and pharyngeal procedures. This is also known as "Joseph Mucosa Hook" which surgeons can use to retract and hold flaps of tissue, in order to access deeper cavities within the surgical field during ENT, plastic and general surgery procedures.

Joseph single hook feature a sharp hook-shaped tip that surgeons can embed in the desired tissue and easily retract a flap to reveal hidden structures. The handle is often round and knurled and it tapers to the hook ends that permits to reach deep cavities. Joseph skin hook manufactured


with premium grade stainless steel to ensure durability and product performance during surgical procedures.

JOSEPH SKIN HOOK FEATURES

- JOSEPH Hook 6.25" (15.9cm) Overall Length.
- Single Sharp Hook Shaped Tip.
- Tapered Shaft Ensuring Deep Surgical Access.
- Non-Slip Knurled Handle for Maximum Control.
- Reusable.
- High Degree of Aesthetic and Corrosion Resistance.
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDES

- 1 Pc JOSEPH Single Hook 6.25" (15.9cm) One Sharp Prong

**WI-70001022****Joseph Skin Hook 6.25" Two Sharp Prong 10mm Wide****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Joseph Double Hook 6.25" (15.9cm) Two Sharp Prong 10mm Wide Retractor Premium Quality Stainless Steel Plastic Dermal ENT Surgical Instrument.


Joseph skin hook is a handheld retractor commonly used for holding back the edges of skin during intranasal and pharyngeal procedures. This is also known as "Joseph Mucosa Hook" which surgeons can use to retract and hold flaps of tissue, in order to access deeper cavities within the surgical field during ENT, plastic and general surgery procedures.

Joseph double hook features two sharp hook-shaped tips that surgeons can embed in the desired tissue and easily retract a flap to reveal hidden structures. The handle is often round and knurled and it tapers to the hook ends that permits to reach deep cavities. Joseph skin hook

manufactured with premium grade stainless steel to ensure durability and product performance during surgical procedures.

JOSEPH DOUBLE PRONG SKIN HOOK FEATURES

- JOSEPH Hook 6.25" (15.9cm) Overall Length.
- 10mm Wide Working Ends
- Double Sharp Hook Shaped Tips.
- Tapered Shaft Ensuring Deep Surgical Access.
- Non-Slip Knurled Handle for Maximum Control.
- Reusable.
- High Degree of Aesthetic and Corrosion Resistance.
- Made with Premium Quality Stainless Steel

**WI-70001021****Joseph Skin Hook 6.25" Two Sharp Prong 12mm Wide****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Joseph Double Hook 6.25" (15.9cm) Two Sharp Prong 12mm Wide Retractor Premium Quality Stainless Steel Plastic Dermal ENT Surgical Instrument.

Joseph skin hook is a handheld retractor commonly used for holding back the edges of skin during intranasal and pharyngeal procedures. This is also known as "Joseph Mucosa Hook" which surgeons can use to retract and hold flaps of tissue, in order to access deeper cavities within the surgical field during ENT, plastic and general surgery procedures.

Joseph double hook features two sharp hook-shaped tips that surgeons can embed in the desired tissue and easily retract a flap to reveal hidden structures. The handle is often round and knurled and it tapers to the hook ends that permits to reach deep cavities. Joseph skin hook

manufactured with premium grade stainless steel to ensure durability and product performance during surgical procedures.

JOSEPH DOUBLE PRONG SKIN HOOK FEATURES

- JOSEPH Hook 6.25" (15.9cm) Overall Length.
- 12mm Wide Working Ends.
- Double Sharp Hook Shaped Tips.
- Tapered Shaft Ensuring Deep Surgical Access.
- Non-Slip Knurled Handle for Maximum Control.
- Reusable.
- High Degree of Aesthetic and Corrosion Resistance.
- Made with Premium Quality Stainless Steel

**WI-70001025****Joseph Skin Hook 6.25" Two Sharp Prong 2mm Wide****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Joseph Double Hook 6.25" (15.9cm) Two Sharp Prong 2mm Wide Retractor
Premium Quality Stainless Steel Plastic
Dermal ENT Surgical Instrument.


Joseph skin hook is a handheld retractor commonly used for holding back the edges of skin during intranasal and pharyngeal procedures. This is also known as "Joseph Mucosa Hook" which surgeons can use to retract and hold flaps of tissue, in order to access deeper cavities within the surgical field during ENT, plastic and general surgery procedures.

Joseph double hook features two sharp hook-shaped tips that surgeons can embed in the desired tissue and easily retract a flap to reveal hidden structures. The handle is often round and knurled and it tapers to the hook ends that permits to reach deep cavities. Joseph skin hook

manufactured with premium grade stainless steel to ensure durability and product performance during surgical procedures.

JOSEPH DOUBLE PRONG SKIN HOOK FEATURES

- JOSEPH Hook 6.25" (15.9cm) Overall Length.
- 2mm Wide Working Ends
- Double Sharp Hook Shaped Tips.
- Tapered Shaft Ensuring Deep Surgical Access.
- Non-Slip Knurled Handle for Maximum Control.
- Reusable.
- High Degree of Aesthetic and Corrosion Resistance.
- Made with Premium Quality Stainless Steel

**WI-70001024****Joseph Skin Hook 6.25" Two Sharp Prong 5mm Wide****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Joseph Double Hook 6.25" (15.9cm) Two Sharp Prong 5mm Wide Retractor
Premium Quality Stainless Steel Plastic
Dermal ENT Surgical Instrument.


Joseph skin hook is a handheld retractor commonly used for holding back the edges of skin during intranasal and pharyngeal procedures. This is also known as "Joseph Mucosa Hook" which surgeons can use to retract and hold flaps of tissue, in order to access deeper cavities within the surgical field during ENT, plastic and general surgery procedures.

Joseph double hook features two sharp hook-shaped tips that surgeons can embed in the desired tissue and easily retract a flap to reveal hidden structures. The handle is often round and knurled and it tapers to the hook ends that permits to reach deep cavities. Joseph skin hook

manufactured with premium grade stainless steel to ensure durability and product performance during surgical procedures.

JOSEPH DOUBLE PRONG SKIN HOOK FEATURES

- JOSEPH Hook 6.25" (15.9cm) Overall Length.
- 5mm Wide Working Ends
- Double Sharp Hook Shaped Tips.
- Tapered Shaft Ensuring Deep Surgical Access.
- Non-Slip Knurled Handle for Maximum Control.
- Reusable.
- High Degree of Aesthetic and Corrosion Resistance.
- Made with Premium Quality Stainless Steel

**WI-70001023****Joseph Skin Hook 6.25" Two Sharp Prong 7mm Wide****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Joseph Double Hook 6.25" (15.9cm) Two Sharp Prong 7mm Wide Retractor Premium Quality Stainless Steel Plastic Dermal ENT Surgical Instrument.

Joseph skin hook is a handheld retractor commonly used for holding back the edges of skin during intranasal and pharyngeal procedures. This is also known as "Joseph Mucosa Hook" which surgeons can use to retract and hold flaps of tissue, in order to access deeper cavities within the surgical field during ENT, plastic and general surgery procedures.

Joseph double hook features two sharp hook-shaped tips that surgeons can embed in the desired tissue and easily retract a flap to reveal hidden structures. The handle is often round and knurled and it tapers to the hook ends that permits to reach deep cavities. Joseph skin hook

manufactured with premium grade stainless steel to ensure durability and product performance during surgical procedures.

JOSEPH DOUBLE PRONG SKIN HOOK FEATURES

- JOSEPH Hook 6.25" (15.9cm) Overall Length.
- 7mm Wide Working Ends
- Double Sharp Hook Shaped Tips.
- Tapered Shaft Ensuring Deep Surgical Access.
- Non-Slip Knurled Handle for Maximum Control.
- Reusable.
- High Degree of Aesthetic and Corrosion Resistance.
- Made with Premium Quality Stainless Steel

**WI-70000224****Nabers Probe #Q-2N Color-Coated Marking 3,6,9,12 Hollow Handle****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Nabers Probe #Q-2N Color Coated Marking 3, 6, 9, 12 Double Ended Hollow Handle Premium Quality Stainless Steel Furcation Probe Dental Periodontal Procedures.


Nabers Probe Q-2N is a double-ended dental instrument designed for diagnosing and measuring furcation involvement in multi-rooted teeth. Its curved, blunt tip ensures safe access and enables accurate assessment of periodontal disease, with depth markings at 3, 6, 9, and 12 mm. Q-2N Nabers Probe is crucial for detecting Class I, II, III, and IV furcation involvements in periodontal treatment planning. Its hollow handle offers a lightweight, non-slip grip that minimizes hand fatigue, and it is made from surgical-grade stainless steel, ensuring durability and resistance to repeated sterilization for clinical use.

NABERS PROB Q-2N FEATURES

- Double-ended Design.
- Color Coated Markings 3, 6, 9, and 12mm for Measurement.
- Curved, Blunt Tip for Safe Access to Ensure Accurate Assessment of Periodontal Disease.
- Non-slip Grip Hollow Handle.
- Use to Diagnose and Measure Furcation Involvement (Bone Loss) in Multi-rooted Teeth.
- Lightweight and Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES

- 1 Pc Nabers Probe #Q-2N Color Coated Marking 3, 6, 9, 12 Double Ended – Hollow Handle

TECHNICAL SPECIFICATIONS

**WI-70000394****News Tracheal Hook 6" Curved Sharp Point****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

NEWS Tracheal Hook 6" (15.2cm) Curved Sharp Point Single Prong with Thick Hook Premium Quality Stainless Steel ENT Surgical Cricothyroidotomy Procedure Instrument.

This is a specialized surgical tool used to hold cartilage during procedures like tracheostomy and cricothyroidotomy, facilitating tube insertion and securing the airway. News tracheal hook feature a thin, curved instrument with a sharp tip for retraction and stabilization of the trachea improving visibility and facilitating the insertion of a tracheostomy tube or other airway devices. During cricothyroidotomy procedures, the hook secures the airway after incision, preventing tissue damage and ensuring proper tube placement.

Made of top quality stainless steel and is reusable after cleaning and sterilization.

This hook is a trusted choice for surgeons seeking precision and excellence in tracheal surgeries.

NEWS HOOK FEATURES

- Crafted with 6" (15.2cm) Overall Length.
- Curved Profile.
- 1 Sharp Pointed Single Prong with Thick Hook End.
- Delicate Hook Design for Precise Tracheal Manipulation
- Ergonomic Handle for Surgeon Comfort
- Fully Autoclaved / Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc NEWS Tracheal Hook 6" (15.2cm)

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**WI-70002066****Periodontal Probe CP15 UNC Color Coded Markings Single Ended****\$10.99** per pieceTrade price from **\$10.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Probe UNC15 Color Coded Marking 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 Single Ended with Hollow Handle Periodontal Diagnostic Instruments

Dental Probe CP-15 Periodontal Probe UNC15 color coded markings are used to assess pocket depths, periodontal attachment levels, and gingival bleeding. It is used to check the depth of periodontal pocket vertically or horizontally. Color markings help determine the depth easily. Accurate & smooth black markings finish, which will not chip, flake or fade.


UNC15 PROBE FEATURES


- Single Ended.
- CP-15 UNC probe markings: 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 on one end.

- Designed to assess periodontal pocket depths.
- Working end has calibrated markings in millimeters
- Made of high-quality stainless steel.
- Ergonomic Hollow Handle.

PACKAGE INCLUDE:

- 1 Pc Periodontal Probe UNC15 Color Coded Marking 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 Single Ended with Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Probe**APPLICATION:** Use to Measure Periodontal Pocket Depths & Attachment Levels**COLOR:** Silver Color**TIP:** Single Ended



WI-70000142

Periodontal Probe UNC15/23 Double Ended Hollow Handle

\$10.99 per piece

Trade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Periodontal UNC15/23 Probe Explorer with Color Coded Marking Double Ended with Hollow Handle Premium Quality Stainless Steel Periodontal Diagnostic Instruments.

CP15/23 probe, also known as UNC-15/Explorer 23 is a double-ended dental instrument used for periodontal diagnostics and caries detection. One end measures periodontal pocket depths (up to 15mm) with precise 1mm markings, while the other end, the #23 explorer tip, detects cavities, calculus, and restoration margins. UNC15/23 probe evaluates periodontal disease severity and assesses furcation lesions and attachment loss. The lightweight, ergonomic hollow handle improves grip and minimizes hand fatigue. Made from high-quality, autoclavable stainless steel, it is designed for daily clinical tasks in diagnostics, restoration,

periodontics surgery, ensuring efficient use over time.

PERIODONTAL UNC15/23 PROBE EXPLORER FEATURES

- Approx 7" (17.8cm) Overall Length.
- Double-ended Handheld Dental Instrument.
- Unc15 Color-coded Marking (1-15mm) for Periodontal Pocket Depths Assessment.
- Explorer #23 Shepherd's Hook for Detecting Tooth Decay at Tooth Enamel.
- Non-slip Hollow Handle for Better Grip.
- Suitable for Periodontal Diagnostic Purposes.
- Lightweight, Easy to Clean, Reusable.
- Made From Premium Quality Stainless Steel

PACKAGE INCLUDE:

WI-36666994

Probe Eye-01 Stainless Steel Surgical Exploration Probe for Clinic

\$0.50 per piece

Trade price from **\$0.20**
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking




ABOUT THIS INSTRUMENT

This is our pattern of Probe Eye-01 Stainless Steel Surgical Exploration Probe for Clinic, made for general and specialty operating theatre use. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one.

Listed at \$0.50, falling to \$0.20 at 100 pieces — about 60% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002164

Rainbow Color Dental Explorer Probe #5 – Double Ended

\$10.99 per piece

Trade price from \$10.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Multi Rainbow Color Dental Probe 17/23
Double Ended Diagnostic Examination
Dental Instrument

Dental Explorer Probe is a professional dental scraper tool helps you to remove stains, tartar, and debris from teeth and gums. Regular care of teeth can reduce the generation of bad breath. This hygiene tool is also used for keeping your pet clean and oral fresh. You can be comfortable to be closer to your pet (cat, dog) without any fears about your health problems.

DENTAL EXPLORER PROBE 17/23 FEATURES

THE FIRST QUALITY CARE TOOL: Our dental tarter scraper is made of 100% high-purity medical stainless steel metal, which allows rust resistance, durability and

safety for the tool. It also has double-ended design, cambered and curving heads that can meet your different needs.

EASY TO USE: Wipe the plaque off the blade with a tissue, then repeat until there is no more plaque on the blade. With this amazing stainless steel dental scarper, you can easily remove those food residues located in the teeth without causing pain and damaging tooth enamel.

CONVENIENT TO BRING AWAY: Thanks to the small size, you can bring it anywhere as personal belongings, no matter where you travel to make you always keep a fresh tone.

100% SATISFACTION GUARANTEE: Customer satisfaction is our utmost priority.

PACKAGE INCLUDE:



WI-43668376

Rentsch Boat Hook with Guard Precision Eye Surgery

\$15.00 per piece

Trade price from \$9.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Rentsch Boat Hook with Guard Precision Eye Surgery — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type retract, hold or pass suture in fine and delicate work. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$15.00, falling to \$9.00 at 10 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-45125898**

Reusable Probe with Eye Surgical Examination Instrument Medical

\$0.59 per pieceTrade price from **\$0.29**

Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Reusable Probe with Eye Surgical Examination Instrument Medical is a surgical instrument we forge, grind and finish in our own Sialkot workshop. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Forging, milling, heat treatment, hand finishing and passivation all happen under

one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$0.59, falling to \$0.29 at 100 pieces — about 51% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-43644620**

Sinskey II Angled Lens Hook Precision Eye Surgery

\$15.00 per pieceTrade price from **\$9.00**

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Sinskey II Angled Lens Hook Precision Eye Surgery for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type retract, hold or pass suture in fine and delicate work. The name tells you what you are getting: an angled shank that keeps the hand clear of the sightline.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$15.00, falling to \$9.00 at 10 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-46285630

Sinskey II Lens Hook Angled-06 Ophthalmology Surgical Instrument

\$15.00 per piece

Trade price from **\$9.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Sinskey II Lens Hook Angled-06 Ophthalmology Surgical Instrument, made for general and specialty operating theatre use. Instruments of this type retract, hold or pass suture in fine and delicate work. The name tells you what you are getting: an angled shank that keeps the hand clear of the sightline.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$15.00, falling to \$9.00 at 10 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-43638650

Sinskey Lens Hook Straight Precision Eye Surgery Tool

\$15.00 per piece

Trade price from **\$9.00**
Minimum order 10 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Sinskey Lens Hook Straight Precision Eye Surgery Tool is a hooks & needle we forge, grind and finish in our own Sialkot workshop. Instruments of this type retract, hold or pass suture in fine and delicate work. The name tells you what you are getting: a straight working end for direct, in-line access.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$15.00, falling to \$9.00 at 10 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

**WI-27945324****Stainless Steel Comfort Grip****\$20.00** per pieceTrade price from **\$12.00**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Comfort Grip — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type retract, hold or pass suture in fine and delicate work.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$20.00, falling to \$12.00 at 1 piece — about 40% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-25801265**Stainless Steel Mini Hair Extension Plier****\$3.50** per pieceTrade price from **\$2.00**

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Mini Hair Extension Plier is a hooks & needle we forge, grind and finish in our own Sialkot workshop. Instruments of this type retract, hold or pass suture in fine and delicate work.

Each unit is inspected by hand for alignment and action before packing, and

we can photograph your actual pieces before they ship. Listed at \$3.50, falling to \$2.00 at 10 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28938462

Surgical

\$15.00 per piece

Trade price from \$12.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type retract, hold or pass suture in fine and delicate work.

Each unit is inspected by hand for alignment and action before packing, and

we can photograph your actual pieces before they ship. Listed at \$15.00, falling to \$12.00 at 1 set — about 20% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-45137869

Surgical Probe with Eye Reusable Stainless Steel Medical Diagnostic

\$0.50 per piece

Trade price from \$0.20
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Probe with Eye Reusable Stainless Steel Medical Diagnostic is part of our surgical instruments range, produced start to finish at our Sialkot facility. The name tells you what you are getting: full reusability — clean, autoclave and use again.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$0.50, falling to \$0.20 at 100 pieces — about 60% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-43358584**

Tracheotomy Hook 15cm Sharp Point Stainless Steel Surgical Instrument

\$7.00 per pieceTrade price from **\$3.00**
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Tracheotomy Hook 15cm Sharp Point Stainless Steel Surgical Instrument is a hooks & needle we forge, grind and finish in our own Sialkot workshop. Instruments of this type retract, hold or pass suture in fine and delicate work. The name tells you what you are getting: sharp points for precise penetration.

The 15cm size shown is our standard; other sizes in the same pattern are made

to order. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$7.00, falling to \$3.00 at 50 pieces — about 57% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-43650631**

Vera Hook Reusable Stainless Steel Eye Surgery

\$15.00 per pieceTrade price from **\$9.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Vera Hook Reusable Stainless Steel Eye Surgery for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type retract, hold or pass suture in fine and delicate work. The name tells you what you are getting: full reusability — clean, autoclave and use again.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$15.00, falling to \$9.00 at 10 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70002071****Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Williams Probe Color Coded 1-2-3-5-7-8-9-10 Marking Root Depth Measurement Periodontal Examination Instruments

Williams Probe is also called the periodontal probe it is the first probe to be used to assess probing depth of a pocket with color codes. Williams Probe used to assess periodontal pocket depths, attachment levels, anatomy configurations and gingival bleeding. It has the following markings – 1,2,3,5,7,8,9,10 mm. It measures the sulcus or pocket depth of the periodontium of each tooth. The working end has calibrated markings in millimeters.


WILLIAMS PROBE FEATURES

- Single Ended.

- Williams probe markings: 1-2-3-5-7-8-9-10 on one end.
- Designed to assess periodontal pocket depths.
- Working end has calibrated markings in millimeters
- Made of high-quality stainless steel.
- Ergonomic Hollow Handle.

PACKAGE INCLUDE:

- 1 Williams Probe Color Coded Marking 1-2-3-5-7-8-9-10 Single Ended with Hollow Handle

TECHNICAL SPECIFICATIONS**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Probe**APPLICATION:** Use to Measure Periodontal Pocket Depths & Attachment Levels**COLOR:** Silver Color



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

Request a quotation

Send us product codes from this catalogue, quantities and your country.

We reply the same working day on WhatsApp.

Online shop

shop.waqarintl.com

Email

sales@waqarintl.com

WhatsApp /

Phone

+92 307 6173627

Works

Roras Road, Opposite Anwar Khawaja Industries,
Sialkot 51310, Pakistan

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.