



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

Pliers & Cutters Surgical Instruments

Reusable stainless-steel pliers & cutters forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



10
PCS



10
PCS



10
PCS



10
PCS

194

PRODUCTS

\$1.50

PRICES FROM

OEM

PRIVATE LABEL

Global

WORLDWIDE
EXPORT

Contents

194 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

Product	Code	Price	Page
10 College Cotton Plier 6" Titanium Rainbow Color	WI-70002169	\$69.99	7
10 College Cotton Plier with Lock 6" Dental Endodontic Instruments	WI-70002167	\$49.99	7
10 Hard Wire Cutter 5" (12.7cm) Straight Tips Orthodontic Pliers	WI-70001692	\$64.99	8
10 T/C Universal Cut & Hold Distal End Cutter Orthodontic Pliers	WI-70001771	\$89.99	8
10-Piece Stainless Steel Ergonomic Orthodontic Pliers Kit	WI-42966263	\$38.00	9
12 Cow Horn Plaster Nipper 8" (20.3cm) Surgical Dental	WI-70003882	\$124.99	9
12 Mathieu Ligature Elastic Placing Pliers 5.5" (14cm)	WI-70002369	\$99.99	10
2 TC Flush Cut Distal End Cutter 4.5" Orthodontic Plier	WI-70001029	\$19.99	10
3 Cow Horn Plaster Nipper 8" (20.3cm) Surgical Dental	WI-70003874	\$31.24	11
3 Elastic Placing Mathieu Needel Holder Plier 5.5"	WI-70003392	\$19.88	11
3 Light Wire Bird Beak Pliers with Groove 5.25" (13.3cm)	WI-70001950	\$19.99	12
3 Mycotie Toe Nail Nippers 6" Double Action	WI-70001992	\$39.99	12
3 Orthodontic Band Seating Placing Aderer Plier Bird Beak Plier Wire Bending Tool	WI-70003291	\$21.99	13
3 Prong Dental Aderer Plier 5.25" Orthodontic Instrument	WI-70003139	\$10.99	13
3 Wire Bending Ortho Pliers of Orthodontic Weingart – Bird Beak – Loop forming	WI-70003279	\$23.99	14
5 College Cotton Plier 6" Titanium Rainbow Color	WI-70002170	\$35.99	14
5 College Cotton Plier with Lock 6" Dental Endodontic Instruments	WI-70002168	\$29.99	15
5 London-College Cotton Pliers 6" Angled Serrated Tips	WI-70003758	\$32.99	15
5 Mycotie Toe Nail Nippers 6" Double Action	WI-70001991	\$64.99	16
5 T/C Universal Cut & Hold Distal End Cutter Orthodontic Pliers	WI-70001772	\$49.99	16
6 Cow Horn Plaster Nipper 8" (20.3cm) Surgical Dental	WI-70003864	\$62.49	17
6 Elastic Placing Mathieu Pliers 5.5" (14cm)	WI-70002370	\$49.99	17
6 Light Wire Bird Beak Pliers with Groove 5.25" (13.3cm)	WI-70001949	\$35.99	18
Abdominal Depressor Stainless Steel General Surgery	WI-36260911	\$35.00	18
Adams Prosthetics Orthodontic Plier 5" Flat Pyramid Beaks	WI-70003367	\$10.99	19
Adhesive Removing Pliers 5.5" Dental Orthodontic Instrument	WI-70001314	\$9.99	19
Angulated Bracket Remover Orthodontic Plier 5.5"	WI-70003364	\$9.99	20
Appendix Removal Surgical Instrument Set Abdominal Operation Tools	WI-46243645	\$150.00	20
Baade Band & Shell Removing Pliers 6" with Pointed Beaks	WI-70003363	\$9.99	21
Ball Hook Crimping Plier 5" Angled Stops Tip	WI-70003362	\$10.99	21
Balloon Expansion Plier 8.26" DIY Home Party Decor Tool	WI-70002802	\$15.18	22
Barrel Spring Toenail Nipper 4.5" (11.43cm)	WI-70001966	\$8.99	22
Barrel Spring Toenail Nipper 5.5" (14cm)	WI-70001965	\$9.99	23
Bird Beak Pliers 5.25" (13.3cm) Orthodontic Instrument	WI-70003361	\$10.99	23
Black Stainless Steel Pliers Clips	WI-25753349	\$20.00	24
Bracket Remover Plier Straight 5.5" Orthodontic Instrument	WI-70003360	\$9.99	24
Cannulated Pin and Wire Cutter 7.5" End Cutting Upto 1.6mm (1/16")	WI-70001201	\$26.99	25
Carbon Steel Stainless Steel Flat Nose Jewelry	WI-29079318	\$15.00	25
Carbon Steel Stainless Steel Multi-Function Jewelry Pliers	WI-28998429	\$15.00	26
Chain Nose Jewelry Plier 4.5" Long Jaws with Blue Grip Handle	WI-70001139	\$10.99	26
College Cotton Plier 6" Titanium Rainbow Color	WI-70002171	\$10.99	27

Contents (continued)

Product	Code	Price	Page
Color Stainless Steel Carbon Steel Mini Pliers	WI-29068379	\$30.00	28
Color Stainless Steel Hair	WI-39799555	\$15.00	28
Color Stainless Steel Hair Extension	WI-37699467	\$15.00	29
Complete Surgical Suture Practice Kit Medical Student Training Model	WI-02950495	\$4.50	29
Coon Style Ligature Tying Plier 6.5" for Wires Up to .015" (.38 mm)	WI-70000826	\$10.99	30
Coral Frag Cutter 7.5" SPS Hard Reef Propagation Tool	WI-70000982	\$14.99	30
Cow Horn Plaster Nipper 8" with Handle Lock	WI-70003853	\$12.49	31
Crown and Shell Crimping Plier 5" Dental Orthodontics	WI-70003333	\$10.99	31
Cuticle Nipper 4.5" Double Spring with Back Safety Lock	WI-70000169	\$9.99	32
Dela Rosa Pliers Dental Orthodontic	WI-70003359	\$12.76	32
Dental Crown Gripper Remover Pliers 6" Angled w/Replaceable Silicone Tips	WI-70003358	\$10.99	33
Dental Orthodontic Band Seating Plier	WI-70003357	\$23.39	33
Dental Orthodontic Pliers Stand Holder	WI-70000450	\$12.99	34
Design Plier Hair Extension Tools Silver	WI-25812292	\$10.00	34
Design Stainless Steel Hair Extension	WI-39803457	\$15.00	35
Design Stainless Steel Medical Trays	WI-40089331	\$14.00	35
Donnenfeld Soft IOL Cutter Precision Eye Surgery Surgical	WI-43493851	\$70.00	36
Double Action Pin Cutter 15" Cuts Up To 4.8mm	WI-70001928	\$79.99	36
Double Spring Toenail Nipper 4.5" for Thick Ingrown Nails	WI-70000935	\$9.99	37
Extra Large Stainless Steel Wire Mesh Surgical Tray Medical	WI-35030741	\$85.00	37
Fine Dental TC Distal End & Wire Cutter Plier Orthodontic	WI-70003298	\$44.81	38
Fischer Plier 13 cm Dental Orthodontics	WI-70003354	\$19.91	38
Fisherman Fishing Plier 6" Serrated Jaws with Black Grip Handle	WI-70001282	\$9.99	39
Fisherman Fishing Plier 8.5" Fish Hook Remover Split Ring Plier	WI-70001286	\$14.99	39
Force Module Separating Pliers 6.5"	WI-70001678	\$10.99	40
Good Carbon Steel Stainless Steel Jewelry	WI-29018395	\$15.00	40
Gordon Prosthetics Contouring Pliers 5.5" Orthodontia	WI-70003351	\$9.99	41
Hair Extension Pliers Hair Extension Installation Tool	WI-36021075	\$5.00	41
Hair Extension Pliers Hold Tool Keratin Hair Extensions Microlink	WI-14915623	\$10.00	42
Hair Extension Pliers Kit Color	WI-28979253	\$15.00	42
Hair Extension Pliers Kit Surgical	WI-39797573	\$15.00	43
Hair Extension Pliers Stainless Steel Hair	WI-36006074	\$4.00	43
Hammerhead NiTi Tip Back Plier 13.5cm Serrated	WI-70000426	\$9.99	44
Hard Wire Cutter Straight Tips Orthodontic Plier	WI-70001696	\$9.99	44
Hawley Retainer Plier Orthodontic Dental	WI-70003350	\$19.81	45
Heavy Duty Ingrown Toenail Nipper 4.5" Spring Action Handle	WI-70000937	\$9.99	45
High Demand Stainless Steel Color	WI-39831159	\$15.00	46
High Demand Stainless Steel Hair Extension	WI-39828285	\$15.00	46
Hollow Chop Contouring Plier 5.5" for Wires Up to .030"	WI-70000434	\$10.99	47
Hollow Chop Plier 5.5" with Reduced Radius	WI-70003349	\$10.99	47
How Crown Plier Curved Dental Orthodontics Surgical	WI-70003348	\$14.19	48
How Crown Plier Straight Dental Orthodontics	WI-70003347	\$10.99	48

Contents (continued)


Product	Code	Price	Page
How Crown Pliers 5.5" Curved Serrated Tips	WI-70001294	\$10.99	49
How Crown Pliers 5.5" Straight Serrated Tips	WI-70001302	\$10.99	49
J2 Japanese Stainless Steel Surgical Instrument Set	WI-25817410	\$12.00	50
Jarabak Light Wire Bending Plier 5.5" for Wire Up to .020"	WI-70003346	\$9.99	50
Jewelry Nose Plier Beading Wire Shaping	WI-34568261	\$2.73	51
Jewelry Nose Plier Fine Precision Jewelry	WI-34511497	\$2.73	51
Jewelry Nose Plier Fine Precision Jewelry Making Tool	WI-34490662	\$2.73	52
Jewelry Nose Plier Fine Tip Precision Jewelry Making Pliers	WI-34519489	\$2.73	52
Jewelry Nose Plier Fine Tip Precision Making Pliers	WI-34557391	\$2.73	53
Jewelry Nose Plier Jewelry Repair and Crafting Tool	WI-34504562	\$2.73	53
Jewelry Nose Plier Mini Precision Tool for Jewelry Making	WI-34492621	\$2.73	54
Jewelry Nose Plier Stainless Steel Precision Jewelry Making Tool	WI-34529523	\$2.73	54
Jewelry Nose Plier Stainless Steel Precision Pliers	WI-34541440	\$2.73	55
Jewelry Pliers Kit Silver	WI-29053358	\$15.00	55
Jewelry Pliers Kit with Carbon	WI-28994443	\$15.00	56
Johnson Contouring Plier 5.5" Dental Orthodontics Instrument	WI-70003345	\$9.99	56
Kim Plier Orthodontic Dental	WI-70003343	\$12.05	57
Lab & Office Plier #121 Serrated Dental Orthodontics	WI-70003342	\$10.99	57
Large Stainless Steel Wire Mesh Tray Reusable Surgical	WI-35035702	\$72.00	58
Leaf Masher Plier 6" with PVC Grip for Glass Work	WI-70001739	\$9.99	58
Ligature forming Plier Orthodontic Orthodontic	WI-70003341	\$14.81	59
Light Wire Bird Beak Plier with Groove 5.25" (13.3cm)	WI-70001951	\$10.99	59
Lingual Arch Forming Plier 5" for Lingual Arch Wires Bending Up to .036"	WI-70000431	\$10.99	60
Long Nose Fisherman Fishing Plier 8" Serrated Jaws with Black Grip Handle	WI-70001281	\$10.99	60
Luxury Stainless Steel Clip Pliers	WI-25767353	\$3.50	61
Major Bowl Stainless Steel Surgical Bowl Medical Utility	WI-36335006	\$45.00	61
Male Torquing Plier 5" (12.7cm)	WI-70000429	\$10.99	62
Manufacturing Hair Extension Pliers Kit Surgical	WI-37716307	\$15.00	62
Marcotte Looping Plier 5.25" Orthodontic Dental Instruments	WI-70003338	\$9.99	63
Mathieu Ligature Elastic Placing Plier 5.5" Curved Serrated Tip	WI-70003296	\$9.99	63
Mathieu Ligature Elastic Placing Plier 5.5" Slim Narrow Tip	WI-70003378	\$9.99	64
Mathieu Ligature Elastic Placing Plier 5.5" Wide Tip	WI-70003337	\$10.61	64
Mini Bracket Remover Plier 5.5" Orthodontic Instrument	WI-70000421	\$9.99	65
Mycotic Toe Nail Nipper 6" Double Action	WI-70003238	\$14.99	65
Mycotic Toe Nail Nipper 6" Double Action Black with Gold Spring	WI-70001313	\$15.99	66
Nance Loop Forming Plier 5" for Wires Up to .022"	WI-70003336	\$10.99	66
O'Brien Loop Forming Orthodontic Pliers	WI-70000430	\$9.99	67
Omega Loop Forming Plier 5.5" Tweed Style Helical Tip	WI-70003335	\$9.99	67
Optical Pliers 5.25" for Helical Closing Loops & Arch Adjustments	WI-70001308	\$9.99	68
Orthodontic Tweed Rect. Arch forming Plier	WI-70003334	\$15.30	68
Orthopedic Pliers Orthodontic	WI-70003282	\$138.75	69
Perry Dressing Plier Surgical Instrument (Model Perry Dressing Plier-01)	WI-46769520	\$3.00	69

Contents (continued)

Product	Code	Price	Page
Perry Dressing Plier Surgical Kit for Fine Tissue	WI-46758578	\$1.50	70
Pin Puller Pliers 5.5" (14cm) with Groove Serrated Jaws	WI-70000370	\$10.99	70
Podiatry Tool Kit 5 PCS Set Curved	WI-26032655	\$25.00	71
Posterior Band Removing Pliers 5.5" Long	WI-70000938	\$10.99	71
Posterior Band Removing Pliers 5.5" Short	WI-70000939	\$10.99	72
Primary Crown Crimping Pliers 5" Slim Profile Tips	WI-70001307	\$9.99	72
Pro-Quality Stainless Steel Hair Extension Plier Top-Rated Strong Grip	WI-25800288	\$3.50	73
Rainbow 3 Prong Dental Aderer Plier 5"	WI-70000860	\$10.99	73
Rainbow Bracket Removing Plier 5.5" Straight	WI-70000858	\$10.99	74
Rainbow Hard Wire Cutter 15 Degree Curved Orthodontic Plier	WI-70000862	\$10.99	74
Rainbow Multi-Color Hair Extension Pliers Kit Surgical	WI-28898462	\$15.00	75
Reusable Stainless Steel Wire Mesh Surgical Tray Medical Sterilization	WI-35085619	\$55.00	75
Reynolds Contouring Pliers 5.5" Curved Tips	WI-70001306	\$9.99	76
Salon Tools Color	WI-25778338	\$3.50	76
Side Cutter Parallel Action Plier 4.5" (11.4cm) Serrated Jaws	WI-70000422	\$59.99	77
Slim Weingart Pliers 5.5" Orthodontic Dental Braces Instrument	WI-70003330	\$10.99	77
Solid Color Stainless Steel Hair	WI-39828267	\$15.00	78
Stainless Steel Carbon Steel Jewelry Making Pliers Kit	WI-29038420	\$15.00	78
Stainless Steel Carbon Steel Jewelry Pliers Set	WI-29040382	\$15.00	79
Stainless Steel Double Drawn Remy Extension Pliers	WI-36033077	\$5.00	79
Stainless Steel Hair Extension	WI-25748314	\$6.00	80
Stainless Steel Hair Extension Clip Pliers	WI-25779279	\$3.50	80
Stainless Steel Hair Extension Plier Clips	WI-25807282	\$7.50	81
Stainless Steel Hair Extension Plier for Hair	WI-25760361	\$10.00	81
Stainless Steel Hair Extension Pliers	WI-39791653	\$15.00	82
Stainless Steel Hair Extension Pliers and Tools	WI-25841232	\$3.50	82
Stainless Steel Hair Extension Pliers Clip Non-Slip Gripped	WI-25769320	\$3.50	83
Stainless Steel Hair Extension Pliers Hair	WI-25792255	\$35.00	83
Stainless Steel Keratin Tape Sealing Pliers	WI-39817363	\$15.00	84
Stainless Steel Keratin Tape Sealing Pliers for Hair Extensions	WI-39810413	\$15.00	84
Stainless Steel Micro Bead	WI-37706370	\$15.00	85
Stainless Steel Multi-Function Waterproof Hair Extension Pliers	WI-25827269	\$3.50	85
Stainless Steel Wire Mesh Sterilization Basket	WI-34929439	\$6.00	86
Stainless Steel Wire Mesh Tray Reusable Surgical Instrument Basket	WI-35023709	\$55.00	86
Standard Pin and Ligature Wire Cutter Pliers Box Joint Dental Orthodontics Tools	WI-70003138	\$14.33	87
Steiner Ligature Tying Plier 6" with Tapered Grooved Tips	WI-70003340	\$9.99	87
T/C Angulated Bracket Remover Plier 5.25" Orthodontic Instrument	WI-70001143	\$12.99	88
T/C Bird Beak Pliers with Cutter 5.5" Orthodontic Instrument	WI-70001400	\$12.99	88
T/C Cannulated Pin and Wire Cutter 7.5" End Cutting Upto 1.6mm (1/16")	WI-70001203	\$43.99	89
T/C Cannulated Pin and Wire Cutter 9" End Cutting Upto 1.6mm (1/16")	WI-70001202	\$49.99	89
T/C End Cutter Double Action 8.5" Max Cutting Capacity 2.4mm (3/32")	WI-70001288	\$49.99	90
T/C Male Torquing Plier 5" (12.7cm) Gold Handle	WI-70001248	\$12.99	90

Contents (continued)

Product	Code	Price	Page
T/C Mathieu Ligature Elastic Placing Plier 5.5" Elastomeric Orthodontic	WI-70003387	\$12.99	91
T/C Pin & Wire Cutter Double Action 7" Angled Max Cutting Capacity 1.6mm (1/16")	WI-70001290	\$39.99	91
T/C Pin & Wire Cutter Double Action 9" Angled Max Cutting Capacity 2.4mm (3/32")	WI-70001289	\$49.99	92
T/C Side Cutter Double Action 9" Max Cutting Capacity 2.4mm 3/32"	WI-70001287	\$49.99	92
T/C Torquing Plier Female/ Outer 5" (12.7cm)	WI-70001067	\$12.99	93
T/C Universal Cut & Hold Distal End Cutter 4.5"	WI-70003140	\$12.99	93
T/C Wire & Pin Cutter 18.5" Side Cutting With Adjustable Bolt	WI-70001927	\$149.99	94
TC Adams Pliers 5" Orthodontic Surgical Dental	WI-70003329	\$12.99	94
TC Bracket Remover Plier 5.5" Straight	WI-70003328	\$12.99	95
TC Flush Cut Distal End Cutter 4.5" Orthodontic Plier	WI-70001030	\$12.99	95
TC Hard Wire Cutter 15 Degree Curved Orthodontic Pliers	WI-70003326	\$12.99	96
TC How Crown Plier Curved Dental Orthodontics Surgical	WI-70003325	\$30.53	96
TC How Crown Plier Straight Dental Orthodontics Surgical	WI-70003324	\$27.46	97
TC Pin and Ligature Hardwire Cutter Straight Tip	WI-70003313	\$12.99	97
TC Pin Puller 5.5" Orthodontic Veterinary Plier	WI-70003323	\$12.99	98
TC Weingart Pliers 5.5" Utility Dental Orthodontic Instrument	WI-70003322	\$12.49	98
Three Prong Plier Orthodontic Dental Instrument	WI-70003321	\$9.99	99
Toe Nail Cutter 5.5" (14cm) Moon Shaped with Back Lock	WI-70003489	\$9.99	99
Torquing Plier Female/ Outer 5" (12.7cm)	WI-70000433	\$10.99	100
Ultimate 6-Piece Stainless Steel Pliers Set 8"	WI-14719250	\$12.00	100
Universal Cut and Hold Distal End Cutter 4.5" Brace Wire Cutting Plier	WI-70003327	\$10.99	101
Universal Pliers Orthodontic Surgical Dental	WI-70003320	\$22.78	101
V-Stop Pliers 5.5" with V Notch Tip	WI-70003319	\$9.99	102
Weingart Plier 5.5" Wire Bending Orthodontic Instrument	WI-70002447	\$10.99	102
Wire Mesh Surgical Tools Sterilization Tray Medical	WI-44552379	\$65.00	103
Young Plier 13.5cm Loop Forming Bending Orthodontic Instrument	WI-70003317	\$10.99	103

**10
PCS****WI-70002169****10 College Cotton Plier 6" Titanium Rainbow Color****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Pieces London College Cotton Dressing Plier Tweezer Serrated Tip Multi Titanium Rainbow Color Surgical Dental Endodontic Instruments 6" Without Lock


The college tweezers are primarily used by dentists moving cotton pellets, wedges, and articulation paper to and from the surgical tray to the mouth. The knurled handles provide a firm, consistent grip even with wet hands. The sharp tips are serrated to ensure a firm grip on most small items, including jewelry, small parts, dissection, splinters, and more!

**LONDON COLLEGE TWEEZERS
FEATURES:**

- Cotton Plier is primarily used to grasp or transfer materials in and out

of the oral cavity. Made of high quality and safe stainless steel, which is very strong, rustproof and has excellent flexibility.

- This professional dental cotton plier is not only used for personal care, medical dressings, fake teeth cleaning, freshening breath, and whitening teeth it also can be used to keep your pet's teeth clean and refresh and avoid those annoying.
- These tooth tweezers are ergonomic design. It is Priority Recommended by lots of dentists. Millions of customers are loving in these useful dental tweezers and also share them with their family and friends.
- Suitable for household, craft, and DIY. Can be used for cleaning out the fish tank, electrical repair, crafting,



RAINBOW DENTAL SURGICAL TWEEZER

6"
(152mm)**10
PCS****WI-70002167****10 College Cotton Plier with Lock 6" Dental Endodontic Instruments****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Pieces London College Cotton Dressing Plier Tweezer 6" (15.2cm) Angled Grooved Tips Self-locking Surgical Dental Endodontic Instruments

The college tweezers are primarily used by dentists moving cotton pellets, wedges, and articulation paper to and from the surgical tray to the mouth. The knurled handles provide a firm, consistent grip even with wet hands. The sharp tips are serrated to ensure a firm grip on most small items, including jewelry, small parts, dissection, splinters, and more!

**COLLEGE COTTON PLIER
FEATURES:**

- Cotton Plier is primarily used to grasp or transfer materials in and out

of the oral cavity. Made of high quality and safe stainless steel, which is very strong, rustproof and has excellent flexibility.

- This professional dental cotton plier is not only used for personal care, medical dressings, fake teeth cleaning, freshening breath, and whitening teeth it also can be used to keep your pet's teeth clean and refresh and avoid those annoying.
- These tooth tweezers are ergonomic design. It is Priority Recommended by lots of dentists. Millions of customers are loving in these useful dental tweezers and also share them with their family and friends.
- Suitable for household, craft, and DIY. Can be used for cleaning out the fish tank, electrical repair, crafting,



DENTAL SELF LOCKING TWEEZERS

6"
(152mm)

**WI-70001692****10 Hard Wire Cutter 5" (12.7cm) Straight Tips Orthodontic Pliers****\$64.99** per pieceTrade price from **\$59.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Hard Wire Cutter 5"
Designed to Cut Hard Wires up to .022"
x .028" with Straight Tip #151
Orthodontic Plier Premium Quality
German Grade Dental Orthodontics
Instruments.

STRAIGHT HARD WIRE CUTTER FEATURES:

ORTHODONTIC HARD WIRE CUTTER:
Easily trim hard orthodontic wires to the proper length with Waqar International's Hard Wire Cutter pliers. They are sharp and have the capability of cutting ligatures at the tip. Straight tips design allows for direct visibility and access to both front and back teeth during adjustments.

PREMIUM QUALITY STAINLESS STEEL: Made with german grade anti-corrosive stainless steel that is built to

withstand daily use of cutting hard wires. Fully reusable, easy to clean, and autoclavable.

HIGH DEGREE OF FLEXIBILITY: The premium quality construction and fine craftsmanship allows the Hard Wire Cutter to cut hard wires up to .022" x .028". Capability of cutting ligatures at the tip + great for making and cutting different fine loops and spring loops.

PRECISE AND USER FRIENDLY DESIGN: The ergonomic, non-slip handles enhance comfort and reduce user fatigue during long procedures.

100% SATISFACTION GUARANTEED:
We take pride in the quality of our products, each item is inspected in the USA. 100% satisfaction guaranteed.

PACKAGE INCLUDES:**WI-70001771****10 T/C Universal Cut & Hold Distal End Cutter Orthodontic Pliers****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 T/C Tip Distal End Cutter
4.5" (11.4cm) with Tungsten Carbide
Inserts designed to cut soft wire pins
and ligature wire up to 0.12"/.30mm.
Slim profile is a plus to access hard-to-
reach areas.

TC DISTAL END CUTTER PLIER FEATURES


TC UNIVERSAL CUT & HOLD DISTAL END CUTTER: The Universal Cut And Hold Distal End Cutters are a sturdy handheld orthodontic tool designed to cut soft wire pins and ligature wire from 0.012"/.30mm to 0.021" x 0.025" (0.53mm x 0.64mm)

DURABLE STAINLESS STEEL CONSTRUCTION: The Premium Quality Stainless Steel and Tungsten Carbide composition of the Universal Cut & Hold Distal End Cutter offers a long life

of strong performance. It is built to last and reliable for long term use and will not bend or break on you. Made with anti-corrosion properties + tarnish free.

FULLY REUSABLE + AUTOCLAVABLE:
A hygienic and valuable choice, the TC Distal End Cutters can be fully autoclaved without fear, for long term use without degradation. They are ultrasonic cleaned and have a beautiful high mirror finish.

DESIGNED FOR PRECISE + EASY OPERATION: The TC Distal End Cutters are designed for easy accessibility, with a comfortable grip and less chance of getting loose at the joints even after long usage. The precision joint with adjustable screw makes it possible to disassemble the plier/cutter and adjust it for wear and tear. These features make this product a valuable





Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-42966263

10-Piece Stainless Steel Ergonomic Orthodontic Pliers Kit

\$38.00 per piece

Trade price from \$22.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of 10-Piece Stainless Steel Ergonomic Orthodontic Pliers Kit, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material. The name tells you what you are getting: an ergonomic handle for long lists without fatigue.

Supplied as a set — the contents can be adjusted to your list. Every piece is

hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$38.00, falling to \$22.00 at 1 set — about 42% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.



WI-70003882

12 Cow Horn Plaster Nipper 8" (20.3cm) Surgical Dental

\$124.99 per piece

Trade price from \$114.99
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pieces Cow Horn Jaws
Plaster Nipper Dental Pliers Plaster
Cutting Nipper Shear 8" (20.3cm)
Gypsum Plaster Scissors Heavy Duty
Stainless Steel Surgical Dental
Instruments

Cow Horn Plaster Cutting Nipper dental instrument designed for accurate divesting of castings. It usually comes with double spring strap lock closure. Hand-held, manually operated heavy cutting trimming scissor applications instrument with a short, curved, toothed blade that is used to remove plaster casts. This dentist's choice instrument is lightweight and easy to grip in dental procedures.

COW HORN PLASTER NIPPER FEATURES

PLASTER CUTTING NIPPER: Cow horn plaster cutting nipper dimension is 8" (20.32cm) in length. The plaster nipper has 2 cow horn jaws one side is serrated and the other is smooth. The double spring ensures a smooth, consistent action and a gripe locking handle for safe storage. We are using high-quality stainless steel.

STAINLESS STEEL CONSTRUCTION:

The stainless steel construction ensures maximum durability and dependability. You can reuse these nipper shear over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.



**WI-70002369****12 Mathieu Ligature Elastic Placing Pliers 5.5" (14cm)****\$99.99** per pieceTrade price from **\$91.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Mathieu Ligature Elastic Placing Plier 5.5" (14cm) Needle Holder Driver, Holding Suture Needles Ligature Separator Placer, Elastomeric Orthodontic Dental Instruments.

Dental Mathieu Ligature plier with serrated tips makes this ideal for gripping and placing elastic auxiliaries. Mathieu Ligature Placing orthodontic plier is designed for placing, holding, and tying ligatures and elastics. These dental pliers are also suitable for guiding, grasping a needle when suturing. Mathieu Orthodontic Pliers are available in various tips for specific size ligatures and elastics and feature extra fine serrations for superior grip and holding of ligatures and elastics.

ELASTIC PLACING MATHIEU PLIERS FEATURES

PREMIUM QUALITY MATHIEU PLIERS 5.5": These Ligature Elastic Placing Pliers and Needle Holder is a high-value Orthodontic Instruments used for various applications such as elastic placing and wire holding.

STAINLESS STEEL CONSTRUCTION:

The Stainless Steel construction ensures maximum durability and dependability. You can reuse these pliers over time as they will not bend or tarnish on you.

FULLY AUTOCLAVABLE: A hygienic choice for the orthodontic practice, these Mathieu Ligature Elastic Placing Pliers are fully autoclavable and corrosion-resistant. They are easy to clean and use over and over for the

WI-70001029**2 TC Flush Cut Distal End Cutter 4.5" Orthodontic Plier****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Pack of 2 T/C Tip Distal End Cutter Flush Cut 4.5" (11.4cm) with Tungsten Carbide Inserts Universal Cutter for All Types of Wires Up to Maximum of 0.021" x 0.025 (.53mm x .64mm) Premium Grade Stainless Steel Dental Orthodontic Instrument.

Flush cut distal end cutter is designed to cut wires as small as 0.30 mm, and as large as 0.53x0.64. The extra thin beaks can cut even the smallest amount of jutting wire. Ensuring no trauma to the patients tissues. The tungsten carbide inserts features gives the beaks maximum sharpness and ensures a long lifetime of operation. The ergonomically designed handles of our instrument, along with the unique L-key screw design, ensure the easiest cutting with the smallest amount of

pressure. Manufactured with premium quality stainless steel to ensure durability. Flush cut distal end cutter is sterilizable and can be reused.

TC FLUSH CUT DISTAL END CUTTER FEATURES

- Universal Cutter with Ultra Thin Beaks Leaving No Trace Of The Wire.
- Maximum Cutting from 0.30mm to 0.53x0.64mm.
- Tungsten Carbide Inserts Cutting Edge For Maximum Efficiency & Durability.
- L-Key Screw Design.
- Ergonomic Grip To Decrease Effort.
- Resuable.
- Made with Premium Quality Stainless Steel.

**12
PCS****2
PCS**

**WI-70003874****3 Cow Horn Plaster Nipper 8" (20.3cm)
Surgical Dental****\$31.24** per pieceTrade price from **\$28.74**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pieces Cow Horn Jaws
Plaster Nipper Dental Pliers Plaster
Cutting Nipper Shear 8" (20.3cm)
Gypsum Plaster Scissors Heavy Duty
Stainless Steel Surgical Dental
Instruments

Cow Horn Plaster Cutting Nipper dental instrument designed for accurate divesting of castings. It usually comes with double spring strap lock closure. Hand-held, manually operated heavy cutting trimming scissor applications instrument with a short, curved, toothed blade that is used to remove plaster casts. This dentist's choice instrument is lightweight and easy to grip in dental procedures.

**COW HORN PLASTER NIPPER
FEATURES**

PLASTER CUTTING NIPPER: Cow horn plaster cutting nipper dimension is 8" (20.32cm) in length. The plaster nipper has 2 cow horn jaws one side is serrated and the other is smooth. The double spring ensures a smooth, consistent action and a gripe locking handle for safe storage. We are using high-quality stainless steel.

STAINLESS STEEL CONSTRUCTION: The stainless steel construction ensures maximum durability and dependability. You can reuse these nipper shear over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

**WI-70003392****3 Elastic Placing Mathieu Needel Holder Plier
5.5"****\$19.88** per pieceTrade price from **\$18.29**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Mathieu Ligature Elastic
Placing Plier 5.5" (14cm) Needle Holder
Driver, Holding Suture Needles Ligature
Separator Placer, Elastomeric
Orthodontic Instruments.


Dental Mathieu Ligature plier with serrated tips makes this ideal for gripping and placing elastic auxiliaries.. Mathieu Ligature Placing orthodontic plier is designed for placing, holding, and tying ligatures and elastics. These dental pliers are also suitable for guiding, grasping a needle when suturing. Mathieu Orthodontic Pliers are available in various tips for specific size ligatures and elastics and feature extra fine serrations for superior grip and holding of ligatures and elastics.

**MATHIEU LIGATURE PLIER
FEATURES**

PREMIUM QUALITY MATHIEU PLIERS 5.5": These Ligature Elastic Placing Pliers and Needle Holder is a high-value Orthodontic Instruments used for various applications such as elastic placing and wire holding.

STAINLESS STEEL CONSTRUCTION: The Stainless Steel construction ensures maximum durability and dependability. You can reuse these pliers over time as they will not bend or tarnish on you.

FULLY AUTOCLAVABLE: A hygienic choice for the orthodontic practice, these Mathieu Ligature Elastic Placing Pliers are fully autoclavable and corrosion-resistant. They are easy to clean and use over and over for the



**3
PCS****WI-70001950****3 Light Wire Bird Beak Pliers with Groove 5.25\"/>****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Light Wire Bird Beak Pliers with Groove. A longer cone than a standard bird beak is useful for forming precise loops in light wires up to 0.20\"/>

Bird Beak Plier is the most commonly used orthodontic instrument used to bend and shape arch wires in various thickness. It can do a variety of wire handling and bending, as well as creating steps and loops. It can also be used in forming retentive clasps for removable devices. These dental orthodontic pliers are available in different patterns and sizes.

BIRD BEAK PLIERS FEATURES

LIGHT WIRE BIRD BEAK PLIER: Light Wire Bird Beak Plier with Groove is

used to give sharp bends, it has 1 flat ended beaks with grooves in the inner side of the beaks which help in holding the wire in place. The Long Bird Beak Dental Pliers are designed to grip the wire up to 0.20\"/>

DURABLE STAINLESS

STEEL CONSTRUCTION: These pliers will not bend or break. Manufactured with anti-corrosion properties and extra hard steel for optimum performance.

FULLY REUSABLE + AUTOCLAVABLE: A hygienic and valuable option, these orthodontic loop forming pliers can be fully sterilized without fear, for long-term use without degradation.

TOUGH TAPERED NOSE FOR PRECISE LOOPS: The hard tapered nose of these dental pliers offers maximum precision

**3
PCS****WI-70001992****3 Mycotic Toe Nail Nippers 6\"/>****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Mycotic Toe Nail Clippers Double Barrel Spring Heavy-duty Podiatry Professional Clippers for Thick and Tough Fingernails and Toenails

Mycotic Toenail Nippers Podiatrist Pedicure Nippers for Toenail Fungus Surgical Stainless Steel Ingrown Toenail Clipper Cutter Trimmer. Besides the pedicurists podiatrists elderly care professionals, Seniors, home personal care, this nail nipper is especially designed for the who are suffering from fungus infections. Regular toenail clippers are proven to be ineffective painful and cause even more damage. Fungus infections nails requires a right tool that can do the job and prevent injury. Ideal for cutting of problematic nails.

MYCOTIC TOENAIL NIPPERS FEATURES

MYCOTIC TOE NAIL CLIPPERS: The Mycotic Nipper (also known as "Cantilever Toe Nail Clipper") is the most durable pedicure nipper designed to work on extremely thick hard damaged and fungus infected toenails. Mycotic Double Action heavy-duty toenail Nippers are designed for thick and tough fingernails and toenails. The podiatrist toenail clippers are ideal for both professionals and home use – for Men, Women and Seniors. Effortlessly and safely cut thick toenails, ingrown nails and hangnails with the razor sharp edge and double barrel spring clippers.

PREMIUM QUALITY STAINLESS

STEEL: Forged from O.R. Grade Stainless Steel, the Podiatrist Toenail



**WI-70003291****3 Orthodontic Band Seating Placing Aderer Plier Bird Beak Plier Wire Bending Tool****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

3 Pcs Orthodontic Band Seating Placing and Wire Bending Tool Set Including 3 Prong Aderer Plier, Bird Beak Plier and Band Pusher Scaler Premium Quality Dental Orthodontics Instrument.

The set of pliers with band pusher is ergonomically designed for precise wire bending in orthodontic procedures, ensuring accuracy in loop and spring creation and bracket placement. This 3 pcs set offers high-quality orthodontic instruments designed for precise wire bending, bracket placement, and band adjustments, ensuring superior results. These high-quality instruments are designed for long-term use and consistent performance, are fully autoclavable and reusable.

ORTHODONTIC WIRE BENDING TOOL SET INCLUDES

- 1 Pc 3 Prong Aderer Plier
- 1 Pc Bird Beak Plier
- 1 Pc Band Pusher Scaler

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Orthodontic Instruments**APPLICATION:** Wire Bending Procedure**REUSABLE:** Yes**AUTOCCLAVABLE:** Yes, Must Autoclave Before Use**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**WI-70003139****3 Prong Dental Aderer Plier 5.25" Orthodontic Instrument****\$10.99** per pieceTrade price from **\$10.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Aderer Plier 5.25" (13.3cm) Three Prong Orthodontic Wire Bending Triple Beak Plier for Light Wires upto 0.030" (0.76 mm) Stainless Steel Dental Supply Instruments

THREE PRONG ORTHODONTIC PLIER

Aderer pliers are also called three-jaw or three-prong pliers and are designed for consistent bends of wires, precision aligned tips, and tapered prongs for fine wire bending. Tips are rounded for ease of contouring light wires up to 0.030" (0.76 mm).

3 Prong Orthodontic Pliers useful for grasping wires and clasps for intraoral arch wire adjustments as well as for lingual arch adjustments. Three prong

orthodontic pliers can have a wide range of beak characteristics, and may be long and tapered for hard to reach quadrants or short and stubby for firmer control. They also can be offset, double sided, mini, Aderer style, round pronged or a combination. Pliers Made Of High Quality stainless steel alloy with exactly machined tool tips and precisely guided sheaths. Ergonomic design with rounded edges for fatigue-free working.

DENTAL ADERER PLIER FEATURES

DURABLE HIGH QUALITY STAINLESS STEEL: Forged from medical stainless steel, you can depend on this set of hemostats to last you many years, even with constant use. Built to never rust, break or even bend under pressure,



**WI-70003279****3 Wire Bending Ortho Pliers of Orthodontic Weingart – Bird Beak – Loop forming****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

3 Pcs Wire Bending Ortho Pliers Set for Orthodontics Weingart, Bird Beak and Hollow Chop Countouring Plier for Bending, Shaping and Loop Forming Orthodontic Archwires – Premium Quality Dental Instruments.



The set of 3 wire bending ortho pliers is a robust handheld tool for contouring and bending stainless steel wires up to 0.030". Dental pliers are essential tools for bending and shaping orthodontic wires during procedures, allowing orthodontists to customize wire shapes for individual patients, ensuring optimal tooth movement and alignment. Bird Beak wire bending plier is ideal for precise loops forming archwires. Weingart pliers with tapered

and serrated tips useful for placement and removal of archwires while hollow chop contouring plier standard size instrument for contouring and forming of archwires without use of scoring and/or torque.

These 3 Wire Bending Ortho Pliers are made of Stainless Steel offer precision in bending and adjusting light and soft wires, with anti-corrosion properties for long-lasting performance.

3 WIRE BENDING ORTHO PLIERS SET INCLUDES

- 1 Pc Weingart Plier
- 1 Pc Bird Beak
- 1 Pc Hollow Chop Countouring Plier

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Orthodontic Pliers**WI-70002170****5 College Cotton Plier 6" Titanium Rainbow Color****\$35.99** per pieceTrade price from **\$33.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pieces London College Cotton Dressing Plier Tweezer Serrated Tip Multi Titanium Rainbow Color Surgical Dental Endodontic Instruments 6" Without Lock


The college tweezers are primarily used by dentists moving cotton pellets, wedges, and articulation paper to and from the surgical tray to the mouth. The knurled handles provide a firm, consistent grip even with wet hands. The sharp tips are serrated to ensure a firm grip on most small items, including jewelry, small parts, dissection, splinters, and more!

LONDON COLLEGE TWEEZERS FEATURES:

- Cotton Plier is primarily used to grasp or transfer materials in and out

of the oral cavity. Made of high quality and safe stainless steel, which is very strong, rustproof and has excellent flexibility.

- This professional dental cotton plier is not only used for personal care, medical dressings, fake teeth cleaning, freshening breath, and whitening teeth it also can be used to keep your pet's teeth clean and refresh and avoid those annoying.
- These tooth tweezers are ergonomic design. It is Priority Recommended by lots of dentists. Millions of customers are loving in these useful dental tweezers and also share them with their family and friends.
- Suitable for household, craft, and DIY. Can be used for cleaning out the fish tank, electrical repair, crafting,



**WI-70002168****5 College Cotton Plier with Lock 6" Dental Endodontic Instruments****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pieces London College Cotton Dressing Plier Tweezer 6" (15.2cm) Angled Grooved Tips Self-locking Surgical Dental Endodontic Instruments


The college tweezers are primarily used by dentists moving cotton pellets, wedges, and articulation paper to and from the surgical tray to the mouth. The knurled handles provide a firm, consistent grip even with wet hands. The sharp tips are serrated to ensure a firm grip on most small items, including jewelry, small parts, dissection, splinters, and more!

COLLEGE COTTON PLIER FEATURES:

- Cotton Plier is primarily used to grasp or transfer materials in and out

of the oral cavity. Made of high quality and safe stainless steel, which is very strong, rustproof and has excellent flexibility.

- This professional dental cotton plier is not only used for personal care, medical dressings, fake teeth cleaning, freshening breath, and whitening teeth it also can be used to keep your pet's teeth clean and refresh and avoid those annoying.
- These tooth tweezers are ergonomic design. It is Priority Recommended by lots of dentists. Millions of customers are loving in these useful dental tweezers and also share them with their family and friends.
- Suitable for household, craft, and DIY. Can be used for cleaning out the fish tank, electrical repair, crafting,

**WI-70003758****5 London-College Cotton Pliers 6" Angled Serrated Tips****\$32.99** per pieceTrade price from **\$30.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs London College Cotton Pliers Tweezer 6" (15.2cm) Angled Serrated Tips Premium Quality Stainless Steel Surgical Dental Endodontic Instruments


London College Cotton Pliers are primarily used by dentists moving cotton pellets, wedges, and articulation paper to and from the surgical tray to the mouth. The knurled handles provide a firm, consistent grip even with wet hands. Cotton pliers has sharp angled serrated tips to ensure a firm grip on most small items, including jewelry, small parts, dissection, splinters, and more. Made with premium quality stainless steel ensures durability and lasting performance. Fully autoclaved and reusable.

COLLEGE COTTON PLIER FEATURES:

- 6" (15.2cm) Overall Length.
- Angled Tips.
- Serrated Tapered Jaws for Deep Access.
- Non-slip Spring Action Thumb Grip Handle.
- Fully Autoclavable
- Reusable
- Made with Premium Quality Stainless Steel


PACKAGE INCLUDE:

- 5 Pcs London College Cotton Plier 6" (15.2cm) Angled – Serrated Tips

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental

**5
PCS**

DOUBLE ACTION NIPPER

**WI-70001991****5 Mycotic Toe Nail Nippers 6" Double Action****\$64.99** per pieceTrade price from **\$59.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Mycotic Toe Nail Clippers
Double Barrel Spring Heavy-duty
Podiatry Professional Clippers for Thick
and Tough Fingernails and Toenails

Mycotic Toenail Nippers Podiatrist
Pedicure Nippers for Toenail Fungus
Surgical Stainless Steel Ingrown Toenail
Clippers Thick Toenails Fingernail
Clipper Cutter Trimmer. Besides the
pedicurists podiatrists elderly care
professionals, Seniors, home personal
care, this nail nipper is especially
designed for the who are suffering from
fungus infections. Regular toenail
clippers are proven to be ineffective
painful and cause even more damage.
Fungus infections nails requires a right
tool that can do the job and prevent
injury. Ideal for cutting of problematic
nails.

**MYCOTIC TOENAIL NIPPERS
FEATURES**


MYCOTIC TOE NAIL CLIPPERS: The
Mycotic Nipper (also known as
"Cantilever Toe Nail Clipper") is the
most durable pedicure nipper designed
to work on extremely thick hard
damaged and fungus infected toenails.
Mycotic Double Action heavy-duty
toenail Nippers are designed for thick
and tough fingernails and toenails. The
podiatrist toenail clippers are ideal for
both professionals and home use – for
Men, Women and Seniors. Effortlessly
and safely cut thick toenails, ingrown
nails and hangnails with the razor sharp
edge and double barrel spring clippers.

PREMIUM QUALITY STAINLESS

STEEL: Forged from O.R. Grade
Stainless Steel, the Podiatrist Toenail

**5
PCS**

T/C DISTAL END CUTTER

**WI-70001772****5 T/C Universal Cut & Hold Distal End Cutter
Orthodontic Pliers****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 T/C Tip Distal End Cutter 4.5"
(11.4cm) with Tungsten Carbide Inserts
designed to cut soft wire pins and
ligature wire up to 0.12"/.30mm. Slim
profile is a plus to access hard-to-
reach areas.

**TC DISTAL END CUTTER PLIER
FEATURES**

**TC UNIVERSAL CUT & HOLD DISTAL
END CUTTER:** The Universal Cut And
Hold Distal End Cutters are a sturdy
handheld orthodontic tool designed to
cut soft wire pins and ligature wire from
0.012"/.30mm to 0.021" x 0.025"
(0.53mm x 0.64mm)

DURABLE STAINLESS STEEL

CONSTRUCTION: The Premium Quality
Stainless Steel and Tungsten Carbide
composition of the Universal Cut &
Hold Distal End Cutter offers a long life

of strong performance. It is built to last
and reliable for long term use and will
not bend or break on you. Made with
anti-corrosion properties + tarnish free.

FULLY REUSABLE + AUTOCLAVABLE:

A hygienic and valuable choice, the TC
Distal End Cutters can be fully
autoclaved without fear, for long term
use without degradation. They are
ultrasonic cleaned and have a beautiful
high mirror finish.

DESIGNED FOR PRECISE + EASY

OPERATION: The TC Distal End Cutters
are designed for easy accessibility,
with a comfortable grip and less
chance of getting loose at the joints
even after long usage. The precision
joint with adjustable screw makes it
possible to disassemble the plier/cutter
and adjust it for wear and tear. These
features make this product a valuable

**WI-70003864****6 Cow Horn Plaster Nipper 8" (20.3cm)
Surgical Dental****\$62.49** per pieceTrade price from **\$57.49**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pieces Cow Horn Jaws
Plaster Nipper Dental Pliers Plaster
Cutting Nipper Shear 8" (20.3cm)
Gypsum Plaster Scissors Heavy Duty
Stainless Steel Surgical Dental
Instruments

Cow Horn Plaster Cutting Nipper dental instrument designed for accurate divesting of castings. It usually comes with double spring strap lock closure. Hand-held, manually operated heavy cutting trimming scissor applications instrument with a short, curved, toothed blade that is used to remove plaster casts. This dentist's choice instrument is lightweight and easy to grip in dental procedures.

**COW HORN PLASTER NIPPER
FEATURES**

PLASTER CUTTING NIPPER: Cow horn plaster cutting nipper dimension is 8" (20.32cm) in length. The plaster nipper has 2 cow horn jaws one side is serrated and the other is smooth. The double spring ensures a smooth, consistent action and a gripe locking handle for safe storage. We are using high-quality stainless steel.

STAINLESS STEEL CONSTRUCTION: The stainless steel construction ensures maximum durability and dependability. You can reuse these nipper shear over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

**WI-70002370****6 Elastic Placing Mathieu Pliers 5.5" (14cm)****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Mathieu Ligature Elastic
Placing Plier 5.5" (14cm) Needle Holder
Driver, Holding Suture Needles Ligature
Separator Placer, Elastomeric
Orthodontic Dental Instruments.


Dental Mathieu Ligature plier with serrated tips makes this ideal for gripping and placing elastic auxiliaries. Mathieu Ligature Placing orthodontic plier is designed for placing, holding, and tying ligatures and elastics. These dental pliers are also suitable for guiding, grasping a needle when suturing. Mathieu Orthodontic Pliers are available in various tips for specific size ligatures and elastics and feature extra fine serrations for superior grip and holding of ligatures and elastics.

**ELASTIC PLACING MATHIEU
PLIERS FEATURES**

PREMIUM QUALITY MATHIEU PLIERS 5.5": These Ligature Elastic Placing Pliers and Needle Holder is a high-value Orthodontic Instruments used for various applications such as elastic placing and wire holding.

STAINLESS STEEL CONSTRUCTION: The Stainless Steel construction ensures maximum durability and dependability. You can reuse these pliers over time as they will not bend or tarnish on you.

FULLY AUTOCLAVABLE: A hygienic choice for the orthodontic practice, these Mathieu Ligature Elastic Placing Pliers are fully autoclavable and corrosion-resistant. They are easy to clean and use over and over for the



**6
PCS****WI-70001949****6 Light Wire Bird Beak Pliers with Groove 5.25\"/>****\$35.99** per pieceTrade price from **\$33.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Light Wire Bird Beak Pliers with Groove. A longer cone than a standard bird beak is useful for forming precise loops in light wires up to 0.20\"/>

Bird Beak Plier is the most commonly used orthodontic instrument used to bend and shape arch wires in various thickness. It can do a variety of wire handling and bending, as well as creating steps and loops. It can also be used in forming retentive clasps for removable devices. These dental orthodontic pliers are available in different patterns and sizes.

BIRD BEAK PLIERS FEATURES

LIGHT WIRE BIRD BEAK PLIER: Light Wire Bird Beak Plier with Groove is

used to give sharp bends, it has 1 flat ended beaks with grooves in the inner side of the beaks which help in holding the wire in place. The Long Bird Beak Dental Pliers are designed to grip the wire up to 0.20\"/>

DURABLE STAINLESS

STEEL CONSTRUCTION: These pliers will not bend or break. Manufactured with anti-corrosion properties and extra hard steel for optimum performance.

FULLY REUSABLE + AUTOCLAVABLE:

A hygienic and valuable option, these orthodontic loop forming pliers can be fully sterilized without fear, for long-term use without degradation.

TOUGH TAPERED NOSE FOR PRECISE LOOPS: The hard tapered nose of these dental pliers offers maximum precision

**WI-36260911****Abdominal Depressor Stainless Steel General Surgery****\$35.00** per pieceTrade price from **\$26.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Abdominal Depressor Stainless Steel General Surgery is a pliers & cutter we forge, grind and finish in our own Sialkot workshop. Instruments of this type grip, bend or cut wire, pins and hard material. The steel is heat treated and passivated for corrosion resistance, then each instrument is

checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$35.00, falling to \$26.00 at 10 pieces — about 26% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003367****Adams Prosthetics Orthodontic Plier 5" Flat Pyramid Beaks****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Adams Plier 5" (13cm) Wire and Prosthetics Plier For Orthodontic Appliance Construction and Adjustment on Stainless Steel Wires Up to .040" (1.02 mm) – Premium Quality Stainless Steel Dental Orthodontic Instruments

Adam Pliers are precision tools designed for right angle bends and appliance adaptation on stainless steel wires up to .040" (1.02 mm). The Adam plier is a versatile orthodontic instrument that enables the construction of Adam Clasps, precise wire adjustment, and efficient tooth alignment for patient care. Adams wire bending pliers offer precision and control with a flat, pyramid-shaped tip, allowing easy adjustment of orthodontic appliances and bending wires up to .040" at a 90° angle. This

premium quality surgical grade stainless steel product is meticulously crafted for aesthetic and corrosion resistance, making it a popular choice among professionals.

ADAMS PLIER FEATURES

- 5" (13cm) Overall Length.
- 2 Smooth Pyramidal Shaped Beaks.
- Max. Wire Thickness Up to .040" (1.02 mm).
- Designed for Right Angle Bends and Appliance Adaptation.
- Enables Construction of Adam Clasps, Precise Wire Adjustment, and Efficient Tooth Alignment.
- Build for Aesthetic and Corrosion Resistance.
- Fully Autoclavable/ Reusable.

**WI-70001314****Adhesive Removing Pliers 5.5" Dental Orthodontic Instrument****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Adhesive Removing Plier 5.5" (14cm) with Replaceable Occlusal Pad Tip German Grade Stainless Steel Excess Adhesive Remover After Debonding Dental Orthodontic Instruments

Adhesive Removing Plier features an occlusal pad on one side and has a scraper on the other allowing you to remove any adhesive remaining after debonding, without damaging the tooth enamel. Reversible Tip is used to remove the bulk of composite material after debonding attachments. The Ergonomically designed handles of our instrument, along with the box joint design, ensure the easiest movement and minimum fatigue. Our adhesive removing plier is sterilizable and reusable. With anti-oxidation and anti-


corrosion properties, these plier ensure durability.

ADHESIVE REMOVING PLIERS FEATURES

- Reversible Occlusal Pad Tip.
- 1 Sharp Metal Blade Removes the Composite With Ease.
- Ergonomically Designed Handles To Decrease Fatigue.
- Light Weight For Maximum Comfort.
- Fully Autoclaved.
- Made with anti-corrosion German Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Adhesive Removing Plier 5.5" (14cm)

TECHNICAL INFORMATION



WI-70003364

Angulated Bracket Remover Orthodontic Plier 5.5"

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Angulated Bracket Remover Plier 5.5" (14cm) Angled Jaws with Long Handle – Premium Quality Dental Orthodontic Braces Removal Tool for Ceramic + Plastic Brackets.

Angulated bracket removal pliers are a sturdy orthodontics handheld instrument designed to quickly and efficiently remove brackets. The high angle + long handle design is effective for removal of all bracket styles, both anterior and posterior. Suitable for steel ceramic and plastic brackets of all types. This versatile tool for placement and removal of archwires, pins, bands, adhesive, brackets, crowns, and other auxiliaries. Cut time in half and get the braces off faster than ever.

Our angular bracket remover plier provides an extra level of resilience and

strength for ultimate reliability. These Angular Bracket Removing Pliers are made of premium quality stainless steel, offering long-lasting performance and being fully sterilizable and reusable.

ANGULATED BRACKET REMOVER PLIER FEATURES:

- Crafted with 5.5" (14cm) Overall Length.
- High Angled Tips for Precise Operation.
- Suitable for Steel, Ceramic or Plastic Brackets of All Types.
- Long Smooth Handles.
- Perfect for Posterior and Anterior Brackets Removal.
- Fully Reusable + Autoclavable



WI-46243645

Appendix Removal Surgical Instrument Set Abdominal Operation Tools

\$150.00 per piece

Trade price from **\$90.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Appendix Removal Surgical Instrument Set Abdominal Operation Tools is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$150.00, falling to \$90.00 at 1 set — about 40% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003363****Baade Band & Shell Removing Pliers 6" with Pointed Beaks****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Baade Band and Shell Removing Pliers 6" (15.2cm) Pointed Jaws with Setting Screw German Grade Stainless Steel Permanent Band & Crown Removing Dental Orthodontic Instrument.

BAADE band and shell removing plier designed for removing permanent crowns and orthodontic bands. BAADE Pliers offer excellent grip and accuracy when removing cemented appliances from teeth. The pointed beaks provide concentrated force to the tip, while the setting screw aids in grasping the appliance firmly. These pliers apply just the right amount of force to disengage the crown or band without damaging the gingiva or abutment tooth. The Dental Plier Improved No. 255, 74-80 baade pliers are ideally constructed to crimp the gingival margin of the crown

without damaging the restoration or the natural tooth underneath. The high quality stainless steel construction ensures maximum durability.


BAADE BAND & SHELL REMOVING PLIERS FEATURES:

- 6" (15.2cm) Overall Length.
- Pointed Beaks Provide Concentrated Force at the Tip.
- Setting Screw Aids in Grasping.
- Used for Removing Cemented Appliances From Teeth.
- Smooth Grip Handle.
- Reusable.
- Made With Premium Quality Stainless Steel.


PACKAGE INCLUDE:

- 1 Pc BAADE Band and Shell Removing Pliers 6" (15.2cm) with

BAADE BAND & SHELL REMOVING PLIERS



BAADE BAND & SHELL REMOVING PLIERS

**WI-70003362****Ball Hook Crimping Plier 5" Angled Stops Tip****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Ball Hook Crimping Plier 5" (12.7cm) Angled Stops Tips with Smooth Handle Premium Quality Stainless Steel Dental Orthodontic Instruments.


Ball Hook Crimping Pliers are used to create exact hooks and stops in archwire. Angled tip aids in intraoral wire bending and the slotted design fits over the base of the hook while holding it tight during crimping. The slotted design fits over the base of the hook while holding it tight during crimping. These precisely hand crafted dental pliers allow orthodontists to effortlessly crimp hooks, auxiliary stops and posts to arch wires. Made with premium quality stainless steel ensures precise long lasting performance. Fully autoclavable and reusable.

BALL HOOK CRIMPING PLIER FEATURES

- Crafted with 5" (12.7cm) Overall Length
- Angled Stop Tip.
- Slotted Design.
- L-key Design Construction.
- Smooth Plier Grip Handle.
- Fully Autoclavable & Reusable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc all Hook Crimping Plier 5" (12.7cm) Angled

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Orthodontic Plier**APPLICATION:** Placement of Hook on**100% HIGH QUALITY RUST RESISTANT STEEL**
Imported premium grade for maximum and precise resistance

**WI-70002802****Balloon Expansion Plier 8.26\"/>****\$15.18** per pieceTrade price from **\$13.97**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Balloon Opening Tool Balloon Expansion Tool Balloon Opener for Balloon Petals Filling Sequins Multi-Function Pliers, Balloon Expander, Snap Pliers for Party & Wedding Decoration.

Balloon expansion plier is a multifunctional tool that can also be used as balloon expander flared pliers latex balloon pop ball filling sequins tools for home party decoration. Made with high-quality stainless steel, the grapple design is easy to hold and use with less effort. this serves as a time-saving balloon stuffing tool that makes your party decor easy and convenient. Just hold the balloon expansion plier and stretch the neck of the plastic balloon to make it wider to add various things into the balloon.

BALLOON EXPANSION PLIER FEATURES

- Ergonomic design
- Non-slip handle
- Comfortable grip
- Integral metal extended tongs head
- Strong and Durable

TECHNICAL INFORMATION:**LENGTH:** 8.26"**OPENING:** 1.6", 1.3"**MATERIAL:** Stainless Steel Metal Frame with Plastic Handle**INSTRUMENTS TYPE:** Home Decor Instruments**COLOR:** Blue Plastic Handle**REUSABLE:** Yes**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**WI-70001966****Barrel Spring Toenail Nipper 4.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Barrel Spring Toe Nail Nipper 4.5\"/>

BARREL SPRING TOE NAIL NIPPER FEATURES

HEAVY DUTY TOENAIL CLIPPERS: 4.5\"/>

PREMIUM QUALITY STAINLESS STEEL: Forged from O.R. Grade Stainless Steel, the Podiatrist Toenail

Clippers for the best product life and performance. The curved razor sharp edge ensures a precise smooth cut through thick nails with ease. Corrosion resistant, tarnish free with a locking safety clasp for ultimate portability and convenience. High quality curved blade design stays sharper longer for amazing results over many years!

USER FRIENDLY AND EASY TO CLEAN: Ergonomically designed for easy use with an elegantly polished mirror finish and non-slip grip handles. Use it on your own nails with no fuss, thanks to the long handle and blade construction. Slip it into your pocket, bag, backpack suitcase + more for on the go use. Sterilize it with alcohol, or fully autoclave it in a professional clinical setting if needed.

GREAT FOR THICK AND INGROWN

**WI-70001965****Barrel Spring Toenail Nipper 5.5" (14cm)****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Barrel Spring Toe Nail Nipper 5.5" (14cm) Heavy Duty Thick & Tough Toenails & Fingernails Clipper Cutter Stainless Steel manicure Pedicure Chiropody Podiatry Instruments

BARREL SPRING TOE NAIL NIPPER FEATURES**HEAVY DUTY TOENAIL CLIPPERS: 5.5"**

Double Barrel Spring heavy-duty toe nail clippers are designed for thick and tough fingernails and toenails. The podiatrist toenail clippers are ideal for both professionals and home use – for Men, Women and Seniors. Effortlessly and safely cut thick toenails, ingrown nails and hangnails with the razor sharp edge and double barrel spring clippers.

PREMIUM QUALITY STAINLESS STEEL: Forged from O.R. Grade Stainless Steel, the Podiatrist Toenail

Clippers for the best product life and performance. The curved razor sharp edge ensures a precise smooth cut through thick nails with ease. Corrosion resistant, tarnish free with a locking safety clasp for ultimate portability and convenience. High quality curved blade design stays sharper longer for amazing results over many years!

USER FRIENDLY AND EASY TO CLEAN:

Ergonomically designed for easy use with an elegantly polished mirror finish and non-slip grip handles. Use it on your own nails with no fuss, thanks to the long handle and blade construction. Slip it into your pocket, bag, backpack suitcase + more for on the go use. Sterilize it with alcohol, or fully autoclave it in a professional clinical setting if needed.

GREAT FOR THICK AND INGROWN**WI-70003361****Bird Beak Pliers 5.25" (13.3cm) Orthodontic Instrument****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Bird Beak Plier 5.25" (13.3cm) Dental Orthodontic Standard Loop Forming Plier for Round Wires up to 0.30" (0.76mm) with Round and pyramid beaks are smooth to prevent wire scoring.

Bird Beak Plier is the most commonly used orthodontic instrument used to bend and shape arch wires in various thickness. It can do a variety of wire handling and bending, as well as creating steps and loops. It can also be used in forming retentive clasps for removable devices. These dental orthodontic pliers are available in different patterns and sizes.

BIRD BEAK PLIERS FEATURES

BIRD BEAK PLIER: Bird Beak Dental Pliers are two handed instruments

designed to grip, bend or hold wires and metal securely. They feature precisely aligned hard tips make them the perfect loop forming pliers for round wires up to 0.30" (0.76mm). Designed with a tip to offer a uniform pressure on the wire while looping.

DURABLE STAINLESS

STEEL CONSTRUCTION: These pliers will not bend or break. Manufactured with anti-corrosion properties and extra hard steel for optimum performance.

FULLY REUSABLE + AUTOCLAVABLE:

A hygienic and valuable option, these orthodontic loop forming pliers can be fully sterilized without fear, for long-term use without degradation.

TOUGH TAPERED NOSE FOR PRECISE LOOPS: The hard tapered nose of these dental pliers offers maximum precision and control for easy operation. Create





WI-25753349

Black Stainless Steel Pliers Clips

\$20.00 per piece

Trade price from \$15.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Black Stainless Steel Pliers Clips, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$20.00, falling to \$15.00 at 2 pieces — about 25% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003360

Bracket Remover Plier Straight 5.5" Orthodontic Instrument

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Bracket Remover Plier 5.5" (14cm) Straight Beaks Braces Removing Premium Quality Stainless Steel Dental Orthodontic Instruments

DIRECT BOND BRACKET REMOVER PLIER FEATURES

DB BRACKET REMOVER (STRAIGHT): Bracket Removal Plier is a handheld instrument designed to quickly and efficiently remove brackets. Cut time in half and get the braces off faster than ever with this handy instrument that is suitable for steel ceramic and plastic brackets of all types.


DURABLE STAINLESS STEEL CONSTRUCTION: The Premium Quality Stainless Steel composition of direct bond bracket removing plier offers a long life of strong performance, as they are built to last. These pliers will not

bend or break on you. Made with anti-corrosion properties.

FULLY REUSABLE + AUTOCLAVABLE: A hygienic and valuable choice, these orthodontic pliers can be fully autoclaved without fear, for long term use without degradation.

DESIGNED FOR PRECISE + EASY OPERATION: Bracket Removing Ortho Pliers are designed for easy accessibility, with a comfortable grip and working tip to engage the bracket under the wings. They offer smooth maneuvering, giving you more control for precise operation. They will help you to remove braces quickly and efficiently. These features make this product a valuable addition to your inventory.

SAFE AND VALUABLE INSTRUMENT: Our production process has attained





WI-70001201

Cannulated Pin and Wire Cutter 7.5" End Cutting Upto 1.6mm (1/16")

\$26.99 per piece

Trade price from \$24.83

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Cannulated Pin & Wire Cutter 7.5" (19.1cm), End Cutting Upto 3.2mm (1/8") Soft Wire or 1.6mm (1/16") Hard Wire German Grade Stainless Steel Surgical Veterinary Orthopedic Instruments.

Cannulated Pin & wire cutter widely used in fracture repairs and other standard surgical procedures. It is used to cut the inserted pins from the bone. The unique design and sharp jaws of cannulated allow surgeons for effective cutting and removal of pins. Cannulated Pin Cutter is very useful tool for cutting surgical pins. The entire portion of the pinhead is inserted in the cutter and then removed with force. The maximum cutting capacity of this Cannulated pin cutter is around 0.062" in hard wire and 0.125" in soft wire.

CANNULATED PIN & WIRE CUTTER FEATURES

- Sharp Cutting Edges.
- Overall Length 7.5" (19.1cm).
- Maximum Cutting Capacity 3.2mm (1/8") Soft Wire or 1.6mm (1/16") Hard Wire.
- Smooth Plier Grip Handle.
- Reusable.
- Corrosion-resistant.
- German Grade Stainless Steel Construction.

PACKAGE INCLUDE

- 1 Pc Cannulated Pin & Wire Cutter 7.5" Maximum Cutting Capacity 1.6mm (1/16")

TECHNICAL INFORMATION

MATERIAL: Stainless Steel



WI-29079318

Carbon Steel Stainless Steel Flat Nose Jewelry

\$15.00 per piece

Trade price from \$12.00

Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Carbon Steel Stainless Steel Flat Nose Jewelry, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$15.00, falling to \$12.00 at 2 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28998429

Carbon Steel Stainless Steel Multi-Function Jewelry Pliers

\$15.00 per piece

Trade price from \$12.00
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Carbon Steel Stainless Steel Multi-Function Jewelry Pliers is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$15.00, falling to \$12.00 at 2 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001139

Chain Nose Jewelry Plier 4.5" Long Jaws with Blue Grip Handle

\$10.99 per piece

Trade price from \$10.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Chain Nose Jewelry Plier 4.5" (11.4 cm), Long Jaws with Blue PVC Grip Handle Premium Quality Stainless Steel Hobby, Wire Work, DIY Craft & Vinyl Jewelry Making and Repair Tools.


Chain nose jewelry pliers are another must for any jewelry artist or craftsman. They are ideal for closing bead tips, getting into tight spaces, closing and opening jump rings. The flat, tapered jaws are perfect for gripping and bending items such as wire. The slim long jaws design used to grip, hold, and pick small objects without scratch. Optimum size and design is perfect for bending, crimping, opening and closing jump rings, placing objects, and shaping wire. Commonly used for forming loops in wire, these pliers are especially handy when creating curves


and geometric designs. With non-slip PVC grip handle makes it comfortable to grip and reduce manual accidents when crafts.

CHAIN NOSE JEWELRY PLIER FEATURES

- Chain Nose Plier.
- 4.5" (11.4cm) Overall in Length.
- Slim Long Jaws.
- Box Joint Style for Easy Handling.
- Soft Blue PVC Finger Grip Handles.
- Rust Proof for Long Term Use.
- Perfect for All Kinds of Diy Jewelry Making and Repair.
- Made With High Quality Stainless Steel.

PACKAGE INCLUDE



**WI-70002171**

College Cotton Plier 6" Titanium Rainbow Color

\$10.99 per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

London College Cotton Dressing Plier
Tweezer Serrated Tip Multi Titanium
Rainbow Color Surgical Dental
Endodontic Instruments 6" Without
Lock


The college tweezers are primarily used by dentists moving cotton pellets, wedges, and articulation paper to and from the surgical tray to the mouth. The knurled handles provide a firm, consistent grip even with wet hands. The sharp tips are serrated to ensure a firm grip on most small items, including jewelry, small parts, dissection, splinters, and more!

LONDON COLLEGE TWEEZERS FEATURES:

- Cotton Plier is primarily used to grasp or transfer materials in and out

of the oral cavity. Made of high quality and safe stainless steel, which is very strong, rustproof and has excellent flexibility.

- This professional dental cotton plier is not only used for personal care, medical dressings, fake teeth cleaning, freshening breath, and whitening teeth it also can be used to keep your pet's teeth clean and refresh and avoid those annoying.
- These tooth tweezers are ergonomic design. It is Priority Recommended by lots of dentists. Millions of customers are loving in these useful dental tweezers and also share them with their family and friends.
- Suitable for household, craft, and DIY. Can be used for cleaning out the fish tank, electrical repair, crafting,

**WI-70003759**

College Cotton Plier with Lock 6"

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

London College Cotton Dressing Plier
Tweezer 6" (15.2cm) Angled Serrated
Tips with Self-locking Mechanism
Surgical Dental Endodontic Instruments

College tweezers are primarily used by dentists moving cotton pellets, wedges, and articulation paper to and from the surgical tray to the mouth. Knurled handles of this cotton plier provides a firm, consistent grip even with wet hands. The sharp tips are serrated to ensure a firm grip on most small items, including jewelry, small parts, dissection, splinters, and more. Made with premium quality stainless steel to ensure durability and excellent craftsmanship. Fully Autoclaved and reusable.

COLLEGE COTTON PLIER FEATURES:

- London College Cotton Plier is primarily used to grasp or transfer materials in and out of the oral cavity. Made of high quality and safe stainless steel, which is very strong, rustproof and has excellent flexibility.
- This professional dental cotton plier tweezer is not only used for personal care, medical dressings, fake teeth cleaning, freshening breath, and whitening teeth it also can be used to keep your pet's teeth clean and refresh and avoid those annoying.
- These dental tweezers are ergonomic design and priority recommended by lots of dentists. Millions of customers are loving in these useful dental tweezers and also share them with their family and friends.





WI-29068379

Color Stainless Steel Carbon Steel Mini Pliers

\$30.00 per piece

Trade price from \$16.00
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Color Stainless Steel Carbon Steel Mini Pliers — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$30.00, falling to \$16.00 at 2 sets — about 47% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-39799555

Color Stainless Steel Hair

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Color Stainless Steel Hair is a pliers & cutter we forge, grind and finish in our own Sialkot workshop. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-37699467

Color Stainless Steel Hair Extension

\$15.00 per piece

Trade price from \$12.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Color Stainless Steel Hair Extension is a pliers & cutter we forge, grind and finish in our own Sialkot workshop. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$15.00, falling to \$12.00 at 1 set — about 20% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-02950495

Complete Surgical Suture Practice Kit Medical Student Training Model

\$4.50 per piece

Trade price from \$4.20
Minimum order 50 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Complete Surgical Suture Practice Kit Medical Student Training Model for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.50, falling to \$4.20 at 50 sets — about 7% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000826****Coon Style Ligature Tying Plier 6.5" for Wires Up to .015" (.38 mm)****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Coon Style Ligature Tying Pliers 6.5" (16.5cm) For Tying Stainless Steel Ligature Wires Up to .015"(.38 Mm) Premium Quality Stainless Steel Dental Orthodontic Instruments.

Coon style ligature tying plier is designed with the precise thin gripping beaks ensuring no slip. This is the most popular orthodontic pliers for tying stainless steel ligature wires up to .015"(.38 mm). Coon style ligature tying plier could tie ligature wire ties anywhere in the mouth, thanks to the long slender beaks. Blunt ends minimize the chances for patient trauma. The handle of instrument is designed along with the unique rivet screw design, ensure the easiest movement and minimum fatigue.

The premium quality stainless steel construction of this standard utility plier offers a long life of strong performance and made it a reliable orthodontic tool for ligature tying procedures.

COON STYLE LIGATURE TYING PLIER FEATURES

- Crafted With 6.5" (16.5cm) Overall Length.
- For Ligature Wires Up to .015"(.38 Mm).
- Long Slender Beaks.
- Blunt Tips.
- Reverse-action Handles.
- Precise Function.
- Longevity of Functional Life.
- Reusable.
- Made With Premium Grade Stainless Steel

**WI-70000982****Coral Frag Cutter 7.5" SPS Hard Reef Propagation Tool****\$14.99** per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty Coral Frag Cutter 7.5" (19cm) Stainless Steel SPS Hard Reef Frag Cutter Propagation Tool for Fragging Hard and Semi-Hard Stony Corals.

7.5" heavy duty coral cutters are designed for fragging hard and semi-hard hard stony corals. This coral cutter is made with fine quality stainless steel and has 1.2" jaws for cutting large and thicker corals. Precision ground jaws help to cut corals precisely with little to no effort. Coral frag cutter crafted with special rust resistant coated stainless steel for everlasting look and performance..

- Ground Jaws for Precise Coral Cutting.
- Sharp Cutting Tips.
- Saltwater Corrosion Resistant.
- Lightweight.
- Double Spring Action Plier Style Smooth Handles.
- Made with Stainless Steel.

PACKAGE INCLUDES:

- 1 Pc Coral Frag Cutter 7.5" (19cm)

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Fish & Aquatic Pets Instruments**APPLICATION:** Coral Frag Cutting**TIP:** Sharp Ground Jaws**HANDLE:** Spring Action Plier Style Handle**REUSABLE:** Yes**CORAL FRAG CUTTER FEATURES**

- 7.5" (19cm) Long.



**WI-70003853****Cow Horn Plaster Nipper 8" with Handle Lock****\$12.49** per pieceTrade price from **\$11.49**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Cow Horn Jaws Plaster Nipper Dental Pliers Plaster Cutting Nipper Shear 8" (20.3cm) Gypsum Plaster (Plaster of Paris) Scissors Heavy Duty Stainless Steel Surgical Dental Instruments.


Cow Horn Plaster Cutting Nipper dental instrument designed for accurate divesting of castings. It usually comes with double spring strap lock closure. Hand-held, manually operated heavy cutting trimming scissor applications instrument with a short, curved, toothed blade that is used to remove plaster casts. This dentist's choice instrument is lightweight and easy to grip in dental procedures.

COW HORN PLASTER NIPPER FEATURES

PLASTER CUTTING NIPPER: Cow horn plaster cutting nipper is 8" (20.32cm) in length. Plaster nipper has 2 cow horn jaws one side is serrated and the other is smooth. The double spring ensures a smooth, consistent action and a gripe locking handle for safe plaster cutting. We are using high-quality stainless steel.

STAINLESS STEEL CONSTRUCTION: The stainless steel construction ensures maximum durability and dependability. You can reuse these nipper shear over time as they will not bend or tarnish on you. High quality stainless steel, excellent international quality standard.

WHERE WE USE PLASTER CUTTING NIPPER: Plaster Cutting Nippers are used to cut plaster and trim models. Dentists mainly use plaster nippers

**WI-70003333****Crown and Shell Crimping Plier 5" Dental Orthodontics****\$10.99** per pieceTrade price from **\$10.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Crown and Shell Crimping Plier 5" (12.7cm) Regular For Crimping of the Gingival Margin of Stainless Steel and Temporary Crowns Dental Plier Orthodontics Instruments.

Crown & Shell Crimping Pliers #417 are specially designed to crimp the gingival margin of stainless steel and temporary crowns. These crown & band crimping pliers allow orthodontists to effortlessly crimp bands and crowns to ensure a tight fit. With its versatile tip design, these expertly constructed orthodontic pliers can also be used to place stops and tubing on arch wires. These crown and shell plier have a wider tip for the best intraoral access. Made of high-grade stainless steel for its durability and is designed for easy handling.

CROWN AND SHELL CRIMPING PLIER FEATURES

- 5" (12.7cm) Long.
- Crimping Stainless Steel and Temporary Crowns.
- Tips Sit Together Perfectly.
- Creating Accurate and Precise Crimps.
- For Lab and Office Use.
- Ergonomically Designed Handles To Decrease Fatigue.
- Light Weight For Maximum Comfort.
- Fully Autoclaved.
- Made with Anti-corrosion German Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Crown and Shell Crimping Plier 5" (12.7cm)



WI-70000169

Cuticle Nipper 4.5" Double Spring with Back Safety Lock

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Cuticle Nipper 4.5" (11.4cm) Spring Action Plier Handle with Back Lock – Premium Quality Stainless Steel Nail Care Cuticle Trimmer for for Spas and Salons.

Cuticle nippers feature ultra-precise 5mm blades for accurate removal of dead cuticles and hangnails. Their ergonomic design enhances comfort and control, while a built-in safety clasp protects the blades when not in use. With 5mm jaws and an overall length of 4.5", this cuticle nipper efficiently trim nails and remove various dead materials. The back lock feature keeps the nipper jaws closed when not in use, ensuring safety and protecting the tips from damage. Made from lightweight, rust-free German stainless steel, these nippers are also reusable.


CUTICLE NIPPER WITH BACK LOCK FEATURES

- Size: 4.5 Inches (11.4 cm)
- 5mm Sharp Blade for Precise Trimming.
- Spring Action Plier Handle for Comfortable Handling.
- Built-in Safety Clasp Ensures Protection of the Blades When Not in Use.
- Ideal for Removing Dead Cuticle and Hangnails from Sidewalls.
- Autoclavable/Reusable
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Cuticle Nipper 4.5" Spring Action with Handle Lock

TECHNICAL INFORMATION



WI-70003359

Dela Rosa Pliers Dental Orthodontic

\$12.76 per piece

Trade price from **\$11.74**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dela Rosa Pliers Dental Orthodontic is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$12.76, falling to \$11.74 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003358

Dental Crown Gripper Remover Pliers 6" Angled w/Replaceable Silicone Tips

\$10.99 per piece

Trade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Crown Gripper Remover Pliers Angled with Replaceable Silicone Tips & Double Spring Action Handle – Premium Quality Stainless Steel Dental Restoration Instruments.

Temporary Crown Remover Angled Pliers with Replaceable Silicone Tips are specialized dental tools designed for the safe and gentle removal of temporary crowns or provisional bridges. The silicone pads ensure even pressure is applied, preventing damage to porcelain restorations. Made from premium medical stainless steel, these pliers are anti-corrosive and durable, ensuring strong and long-lasting performance without bending or breaking.

- 6" (15.2cm) Overall Length.
- Long Beak With Angled Tips for Hard to Reach Area.
- Interchangeable Silicone/Rubber Pads Tips to Grip Crowns Without Scratching or Chipping.
- Double Spring Opening & Closing Mechanism.
- Extensively Used in Dental Restoration and Implant Procedures.
- Reusable and Full Autoclavability.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Dental Crown Gripper Remover Pliers Angled w/Replacable Silicone Tips

CROWN & BRIDGE REMOVING PLIER FEATURES

TECHNICAL INFORMATION



WI-70003357

Dental Orthodontic Band Seating Plier

\$23.39 per piece

Trade price from **\$21.52**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking




ABOUT THIS INSTRUMENT

Dental Orthodontic Band Seating Plier is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$23.39, falling to \$21.52 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000450

Dental Orthodontic Pliers Stand Holder

\$12.99 per piece

Trade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dental Orthodontic Pliers Stand Holder – Premium Quality Stainless Steel Pliers Rack Shelf Forceps Scissors Organizer for Dental Orthodontic Instrument Equipment.

Stainless steel dental stand holder is a useful tool for organizing and storing various dental orthodontic or jewelry pliers, enhancing safety and efficiency in your work. Dental Plier stand is easy to assemble, lightweight, and portable, featuring two brackets, three crossbars, and six screws, making it a versatile tool for any workspace.

This dental orthodontic pliers stand holder is crafted from durable stainless steel and corrosion-resistant hard alloy material, ensuring maximum durability and stability. We are dedicated to delivering top-notch products and

exceptional customer service that guarantee satisfaction and value for money.

DENTAL ORTHODONTIC PLIERS STAND HOLDER FEATURES

- Simple Assembly and Easy-to-carry Components.
- Allows Reasonable Storage and Work.
- Space Saving and Functional.
- Ensures Clean and Tidy Working Environment.
- Ideal for Orthodontic Craft Jewelry Plier Forceps Scissors.
- Reusable.
- High-quality Stainless Steel Construction.

TECHNICAL INFORMATION

MATERIAL: Stainless Steel



WI-25812292

Design Plier Hair Extension Tools Silver

\$10.00 per piece

Trade price from **\$8.00**

Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Design Plier Hair Extension Tools Silver, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$10.00, falling to \$8.00 at 2 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-39803457

Design Stainless Steel Hair Extension

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Design Stainless Steel Hair Extension is a pliers & cutter we forge, grind and finish in our own Sialkot workshop. Instruments of this type grip, bend or cut wire, pins and hard material.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-40089331

Design Stainless Steel Medical Trays

\$14.00 per piece

Trade price from \$8.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Design Stainless Steel Medical Trays is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$14.00, falling to \$8.00 at 10 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-43493851

Donnenfeld Soft IOL Cutter Precision Eye Surgery Surgical

\$70.00 per piece

Trade price from **\$59.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Donnenfeld Soft IOL Cutter Precision Eye Surgery Surgical for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$70.00, falling to \$59.00 at 1 piece — about 16% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001928

Double Action Pin Cutter 15" Cuts Up To 4.8mm

\$79.99 per piece

Trade price from **\$73.59**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Double Action Pin Cutter 15" (38.1cm) will Cut Up To 4.8mm Diameter with Max Opening 7.8mm Premium Quality Stainless Steel Orthopedic Surgical Instruments.

Bone Pin & Wire Cutter 15", end cutting maximum capacity 4.8mm is suitable for cutting wires in bone fractures in several surgical procedures. This instrument has sterilizable material due to high quality stainless steel crafting. It is quite resilient to last for many years and reused after sterilization.

PIN CUTTER FEATURES

- 15" (38.1cm) Overall Length
- Double Action Pin & Wire Cutter
- End Cutting Maximum Capacity Up To 3/16" (4.8mm) Diameter.
- Max Jaws Opening Upto 8mm.

- Grip Handle for Convenient Handling.
- Made with High Quality Stainless Steel

PACKAGE INCLUDE:

- Pin Cutter Double Action 15" (38.1cm) Cuts Up To 4.8mm Diameter

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENT TYPE: Orthopedic Instruments

APPLICATION: Bone Pin Cutting

CURVATURE: Straight

CUTTING END: Upto 4.8mm Diameter


HANDLE: Grip Handle

REUSABLE: Yes

STERILIZED: Yes

AUTOCCLAVED: Yes

RUST PREVENTION PROCEDURE: Passivated





WI-70000935

Double Spring Toenail Nipper 4.5" for Thick Ingrown Nails

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Double Spring Ingrown Toenail Clippers & Cutter 4.5" (11.4cm) Heavy Duty Professional Toe Nail Nipper for Thick Nails German Grade Stainless Steel Dermatology Podiatry Instruments.

Double Spring Toenail Cutter Clippers are specifically designed for working on hard, brittle, out of reach toenails. Ingrown toenail clippers are a high quality professional tools. The stainless-steel cutters are hand-filed for precision to make even the most awkward of jobs easy. They are perfectly aligned so that they close evenly, no matter the amount of pressure applied. Extremely sharp, hand-filed, pointed blades allow little pressure needed when cutting thick, hard and stubborn nails. High quality stainless steel used to manufacture this


professional podiatrist tool which is easy to clean and wipe down after use. An essential tool for adult personal hygiene and can also be used in nail salons or at home by seniors.

INGROWN TOENAIL CLIPPER CUTTER FEATURES

- 4.5" (11.4cm) Overall Length.
- Sharp Curved Blades.
- Double Spring Action Design.
- Adjustable Screw.
- Smooth Plier Grip Handle.
- Fully Autoclavable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Ingrown Toenail Clipper 4.5" (11.4cm) Double Spring Action Handle



WI-35030741

Extra Large Stainless Steel Wire Mesh Surgical Tray Medical

\$85.00 per piece

Trade price from **\$76.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Extra Large Stainless Steel Wire Mesh Surgical Tray Medical — one of 4,770 instruments we manufacture ourselves rather than buy in. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$85.00, falling to \$76.00 at 10 pieces — about 11% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003298

Fine Dental TC Distal End & Wire Cutter Plier Orthodontic

\$44.81 per piece

Trade price from **\$41.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Fine Dental TC Distal End & Wire Cutter Plier Orthodontic is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$44.81, falling to \$41.23 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003354

Fischer Plier 13 cm Dental Orthodontics

\$19.91 per piece

Trade price from **\$18.32**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Fischer Plier 13 cm Dental Orthodontics is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$19.91, falling to \$18.32 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001282****Fisherman Fishing Plier 6" Serrated Jaws with Black Grip Handle****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Fisherman Fishing Pliers 6" (15.2cm) – Premium Quality Stainless Steel Fishing Plier Long Serrated Jaws Deep-Seated Hooks Remover with Build in Line Cutter Corrosion Resistant Fishing Gear

The needle nose fisherman fishing plier 6" comes with long serrated jaws to help remove deep-seated hooks. Made of stainless steel, built in line cutter great for fishing, making jewelry and wiring. Stainless steel needle nose pliers are integrated with a wire cutter for cutting wires and are designed to aid in wrapping or looping in confined spaces.

This pliers can help you save time and energy, allowing you to better enjoy the fun of fishing and outdoor travel. The top high-quality stainless steel construction with PVC grip handle

ensure good hardness and corrosion resistance. These pliers for fishing measures 6 inch long, and have PVC gripped handles with finger rest for a comfortable use.

FISHING PLIER FEATURES

- 6" Long Needle Fishing Plier
- Serrated Jaws for Deep-Seated Hooks Remover
- Extra Sharp Line Cutters
- Saltwater Corrosion Resistant
- Lightweight
- Black PVC Grip Handle with Finger Rest
- Made with Stainless Steel

PACKAGE INCLUDES:

- 1 Pc Fisherman Fishing Pliers 6" Serrated Jaws with Line Cutter – Black Grip Handle

**WI-70001286****Fisherman Fishing Plier 8.5" Fish Hook Remover Split Ring Plier****\$14.99** per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Fisherman Fishing Pliers 8.5" (21.6cm) Stainless Steel Multitool Fish Hook Remover Line Cutter Split Ring Opener Lure Tuning Tool Saltwater Corrosion Resistant Fishing Gear

This fishing pliers is the ideal choice for every fishing enthusiast. It can not only cut lines, but also open rings, remove hooks, and adjust hook angles. Therefore, this pliers is not only a fishing tool, but also a practical tool for outdoor travel. It can help you save time and energy, allowing you to better enjoy the fun of fishing and outdoor travel. Add this multi-purpose fishing plier to your tackle box for tightening or removing hooks and more.


can be used in freshwater and saltwater environments. These pliers for fishing measures 8.5 inch long, and have pVC gripped handles for a comfortable use.


FISHING PLIER FEATURES

- 8.5"x3" Long Nose Fishing Plier
- Easier Split Rings & Hook Remover
- Strong Gripping Serrated Jaws
- Grip and Roation Holes
- Extra Sharp Line Cutters
- Saltwater Corrosion Resistant
- Lightweight
- Spring Action with PVC Grip Handle
- Made with Stainless Steel

PACKAGE INCLUDES:

- 1 Pc Fisherman Fishing Pliers 8.5" Fish Hook Remover Split Ring Plier





WI-70001678

Force Module Separating Pliers 6.5"

\$10.99 per piece

Trade price from \$10.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Force Module Separating Placing Pliers 6.5" (16.5cm) Stainless Steel Orthodontic Braces Dental Instruments

The Force Module Separating Plier is a blunt ended instrument used to place orthodontic elastics between teeth to create space. Separating Plier is designed for both gripping and stretching elastic separators. The Blunt Ends minimizes the chances for patient trauma.

FORCE MODULE SEPARATING PLIERS FEATURES

PERFECT SEPARATING PLIERS: Separating orthodontic pliers are designed to spread and hold elastics during placement. These pliers feature a long barrel-shaped, tapered tip for

prime accessibility and placement in any quadrant. The beaks may be angled, and the tip is grooved to prevent slippage until the elastic is ready to be placed.

EASY TO USE: This Separating Pliers is user friendly and easy to clean as it is autoclavable.

BEAUTIFULLY FORGED: High mirror finish Separating Pliers.

MADE OF SURGICAL QUALITY STAINLESS STEEL: Forged from corrosion-resistant surgical stainless steel. These popular floor grade Separating Pliers are built to last.

PACKAGE INCLUDES

- 1 Pc Force Module Separating Plier 6.5" (16.5cm)

TECHNICAL INFORMATION



WI-29018395

Good Carbon Steel Stainless Steel Jewelry

\$15.00 per piece

Trade price from \$12.00
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

We manufacture Good Carbon Steel Stainless Steel Jewelry for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$15.00, falling to \$12.00 at 2 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003351

Gordon Prosthetics Contouring Pliers 5.5\"/>

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Gordon Contouring Pliers (Gordon Prosthetics Orthodontia) 5.5 Inches (14cm) Curved Tips German Grade Stainless Steel Dental Orthodontic Instruments

Gordon Pliers are used to shape and contour wires, bands, crowns and prosthetics. Curved beaks with flat working edges. For lab and office use.

Contouring orthodontic pliers are ideal for forming and contouring arch wires, reshaping molar bands and enhancing crown form contours. They also are especially useful for placing accentuated curves in wires. These orthodontic pliers are perfect for the general dentist and cosmetic dentist who require instruments to make minor adjustments to removable prosthetics such as partials and retainers. From

wires to crowns, these pliers ensure expert contouring, enhancing patient comfort and treatment outcomes.

GORDON CONTOURING PLIERS FEATURES

- Curved Beaks With Flat Working Edges.
- Designed for Shaping and Contouring Wires, Bands, Crowns, and Prosthetics.
- Box Joint Style Construction
- Ergonomically Designed Handles To Decrease Fatigue.
- Fully Autoclaved.
- Made with anti-corrosion German Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Gordon Contouring Pliers 5.5\"/>



WI-36021075

Hair Extension Pliers Hair Extension Installation Tool

\$5.00 per piece

Trade price from \$1.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Hair Extension Pliers Hair Extension Installation Tool, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$5.00, falling to \$1.00 at 10 pieces — about 80% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-14915623

Hair Extension Pliers Hold Tool Keratin Hair Extensions Microlink

\$10.00 per piece

Trade price from \$8.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hair Extension Pliers Hold Tool Keratin Hair Extensions Microlink is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$10.00, falling to \$8.00 at 10 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28979253

Hair Extension Pliers Kit Color

\$15.00 per piece

Trade price from \$12.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hair Extension Pliers Kit Color is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$15.00, falling to \$12.00 at 1 set — about 20% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-39797573

Hair Extension Pliers Kit Surgical

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Hair Extension Pliers Kit Surgical for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36006074

Hair Extension Pliers Stainless Steel Hair

\$4.00 per piece

Trade price from \$1.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hair Extension Pliers Stainless Steel Hair — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$1.00 at 10 pieces — about 75% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000426****Hammerhead NiTi Tip Back Plier 13.5cm Serrated****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Hammerhead Niti Tie Back Pliers 13.5cm Serrated Jaws NITI Cinch Back Wire Forming Plier German Grade Stainless Steel Dental Orthodontic Instruments.

Hammerhead NiTi Tip Back Pliers are versatile pliers for bending NiTi wire up to .025" (.64 mm). These wire-forming pliers simplify a multi-step process into one step, simply squeeze and you are done; NO heat required. Dental Hammer Head Plier Serrated Tip NITI Cinch Back Wire Forming Pliers are a robust tool for forming the distal end of Niti orthodontic arch wires. Hammer Head Niti Wire Pliers are ergonomically designed for easy accessibility, providing precise control and smooth operation, making them a valuable addition to your orthodontics toolkit.

Hammerhead serrated tip NITI pliers are made of premium quality stainless steel, providing long-lasting performance and anti-corrosion properties.

HAMMERHEAD NITI TIE BACK PLIER FEATURES

- 5.3" (13.5cm) Overall Length.
- Hammerhead Serrated Tip Design.
- For Wires Up to .025".
- For NiTi Wire Only.
- Multi-use Tie-back NiTi Bender with No Heat.
- Perfect for Bending Nickel-Titanium Arches/wires in the Molar Region
- Durable, Strong, and Corrosion Resistant.
- Fully Autoclavable/Reusable.

**WI-70001696****Hard Wire Cutter Straight Tips Orthodontic Plier****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality German Grade Hard Wire Cutter Designed to Cut Hard Wires up to .022" x .028" with Straight Tip #151 Orthodontic Plier Dental Orthodontics Instruments.

STRAIGHT HARD WIRE CUTTER FEATURES:

ORTHODONTIC HARD WIRE CUTTER: Easily trim hard orthodontic wires to the proper length with Waqar International's Hard Wire Cutter pliers. They are sharp and have the capability of cutting ligatures at the tip. Straight tips design allows for direct visibility and access to both front and back teeth during adjustments.

PREMIUM QUALITY STAINLESS STEEL: Made with german grade anti-corrosive stainless steel that is built to withstand daily use of cutting hard

wires. Fully reusable, easy to clean, and autoclavable.


HIGH DEGREE OF FLEXIBILITY: The premium quality construction and fine craftsmanship allows the Hard Wire Cutter to cut hard wires up to .022" x .028". Capability of cutting ligatures at the tip + great for making and cutting different fine loops and spring loops.

PRECISE AND USER FRIENDLY DESIGN: The ergonomic, non-slip handles enhance comfort and reduce user fatigue during long procedures.

100% SATISFACTION GUARANTEED: We take pride in the quality of our products, each item is inspected in the USA. 100% satisfaction guaranteed.

PACKAGE INCLUDES:

- 1 Pc Hard Wire Cutter 5" (12.7cm) – Straight Tips





WI-70003350

Hawley Retainer Plier Orthodontic Dental

\$19.81 per piece

Trade price from \$18.23
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Hawley Retainer Plier Orthodontic Dental is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$19.81, falling to \$18.23 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000937

Heavy Duty Ingrown Toenail Nipper 4.5" Spring Action Handle

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Ingrown Toenail Nipper 4.5" (11.4cm) Spring Action Heavy Duty Thick Toe Nail Clippers & Cutters German Grade Stainless Steel Dermatology Podiatry Instruments.

Ingrown toenail nippers designed to combat even the toughest nails. Perfect for extra-thick and irregular nails misshapen by fungus, age, or other impediments. Heavy duty stainless steel toe nail nippers perfect tool for men, women, seniors & adults. Toe nail cutter features sharp curved blades, adjustable screw, spring action, & non-slip handles. High quality stainless steel used to manufacture this professional podiatrist tool which is easy to clean and wipe down after use. An essential tool for adult personal hygiene and can

also be used in nail salons or at home by seniors.

HEAVY DUTY INGROWN TOENAIL NIPPER FEATURES

- 4.5" (11.4cm) Overall Length.
- Extra Sharp Curved Blades.
- Single Spring Action Design.
- Adjustable Screw
- Non-Slip Plier Grip Handle.
- Fully Autoclavable
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDE

- 1 Pc Heavy Duty Ingrown Toenail Nipper 4.5" (11.4cm) Spring Action Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel





WI-39831159

High Demand Stainless Steel Color

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of High Demand Stainless Steel Color, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-39828285

High Demand Stainless Steel Hair Extension

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of High Demand Stainless Steel Hair Extension, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000434**

Hollow Chop Contouring Plier 5.5" for Wires Up to .030"

\$10.99 per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Hollow Chop Contouring Plier 5.5" (14cm) with Standard Radius for Contouring of Stainless Steel Arch Wires up to .030" (.76 mm) Premium Quality Stainless Steel Orthodontic Dental Instruments.

Hollow Chop Contouring Plier is designed for easy contouring and shaping of stainless steel arch wires. This orthodontic contouring plier feature a concave and convex tip design, similar to De La Rosa Pliers, primarily used for dealing with arch wires. The heavy beak of hollow chop contouring pliers enabling the shaping and contouring of arch wires up to 0.030" in diameter. Made from high quality Stainless Steel, are meticulously crafted for optimal efficiency in

orthodontics, making them ideal for long-term use.

HOLLOW CHOP CONTOURING PLIER FEATURES

- Crafted with 5.5" (14cm) Overall Length.
- Pliers Has a Standard Radius.
- Perfect For Stainless Steel Arch Wires Upto .030" (.76 Mm).
- Heavy Concave & Convex Beak Design.
- 8.5mm Concave & Convex Tip for Intricate Wire Contouring.
- Balanced Joint Assembly for Precise Tip Alignment.
- Ergonomic Handle for Balance and Grip.
- Highly Corrosion-resistant. Honed.
- Reusable.

**WI-70003349**

Hollow Chop Plier 5.5" with Reduced Radius

\$10.99 per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Hollow Chop Plier 5.5" (14cm) with Recued Radius for Contouring of Stainless Steel Arch Wires up to .030" (.76 mm) Premium Quality Stainless Steel Orthodontic Dental Instruments.

Hollow chop contouring pliers are specialized orthodontic tools used to shape or contour stainless steel arch wires, with a smaller radius than standard sizes. They are highly effective in making subtle adjustments and forming palatal arches in dental applications. The hollow chop arch forming plier features a slightly long concave and convex tip design with a smooth inner surface for subtle contouring and forming archwires without scoring or torquing, making it ideal for accentuated curves.

- Crafted with 5.5" (14cm) Overall Length.
- Pliers Has a Reduced Radius.
- For Stainless Steel Arch Wires Upto .030" (.76 Mm).
- Slightly Slim & Long Concave & Convex Beak Design.
- 6mm Concave & Convex Tip for Intricate Wire Contouring.
- Balanced Joint Assembly for Precise Tip Alignment.
- Ergonomic Handle for Balance and Grip.
- Highly Corrosion-resistant. Honed.
- Reusable.
- Made With Premium Grade Stainless Steel.

PACKAGE INCLUDE:

HOLLOW CHOP PLIERS FEATURES





WI-70003348

How Crown Plier Curved Dental Orthodontics Surgical

\$14.19 per piece

Trade price from **\$13.05**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

How Crown Plier Curved Dental Orthodontics Surgical is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$14.19, falling to \$13.05 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003347

How Crown Plier Straight Dental Orthodontics

\$10.99 per piece

Trade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

How Crown Pliers Straight 5.5" (14cm) with Serrated Tips German Grade Stainless Steel 110 Dental Orthodontic Plier for Wire Bending, Arch Forming, Crown Gripper Multipurpose Dental Oral Orthodontic Instruments

HOW orthodontic pliers are considered general utility pliers and are used for a variety of tasks, they are especially designed for grasping auxiliary. These multipurpose pliers are used to place and remove pins, archwires, and other orthodontic materials. How Plier's feature a fine serrated tip for a firm grip with smoothened edges to minimize patient discomfort. These dental pliers are available in Straight and Curved variations.

HOW CROWN PLIERS FEATURES

110 STRAIGHT HOW PLIERS: The Dental Orthodontic Pliers 110 How Crown Straight is hand held instrument designed to grip, bend or hold wires and metal securely. These multipurpose dental pliers great to use to effectively place and remove pins, archwires, and other auxiliaries.

SERRATED TIPS FOR ULTIMATE CONTROL: The round pattern serrated tips on these dental pliers offer more control for precise and easy operation.

DURABLE STAINLESS STEEL CONSTRUCTION: The Premium Quality Stainless Steel composition of these How Crown Straight Dental Pliers offers a long life of strong performance, as they are built to last. These pliers will not bend or break on you. Made with anti-corrosion properties.

FULLY REUSABLE + AUTOCLAVABLE:

**WI-70001294****How Crown Pliers 5.5\"/>****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

How Crown Pliers Curved 5.5\" (14cm) with Serrated Tips German Grade Stainless Steel 111 Dental Orthodontic Plier for Wire Bending, Arch Forming, Crown Gripper Multipurpose Dental Oral Orthodontic Instruments

HOW orthodontic pliers are considered general utility pliers and are used for a variety of tasks, they are especially designed for grasping auxiliary. These multipurpose pliers are used to place and remove pins, archwires, and other orthodontic materials. How Plier's feature a fine serrated tip for a firm grip with smoothened edges to minimize patient discomfort. These dental pliers are available in Straight and Curved variations.

HOW CROWN PLIERS FEATURES

111 CURVED HOW PLIERS: The Dental Orthodontic Pliers 111 How Crown Curved is hand held instrument designed to grip, bend or hold wires and metal securely. These multipurpose dental pliers great to use to effectively place and remove pins, archwires, and other auxiliaries.

SERRATED TIPS FOR ULTIMATE CONTROL: The round pattern serrated tips on these dental pliers offer more control for precise and easy operation.

DURABLE STAINLESS STEEL CONSTRUCTION: The Premium Quality Stainless Steel composition of these How Crown Curved Dental Pliers offers a long life of strong performance, as they are built to last. These pliers will not bend or break on you. Made with anti-corrosion properties.

FULLY REUSABLE + AUTOCLAVABLE:**WI-70001302****How Crown Pliers 5.5\"/>****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

How Crown Pliers Straight 5.5\" (14cm) with Serrated Tips German Grade Stainless Steel 110 Dental Orthodontic Plier for Wire Bending, Arch Forming, Crown Gripper Multipurpose Dental Oral Orthodontic Instruments


HOW orthodontic pliers are considered general utility pliers and are used for a variety of tasks, they are especially designed for grasping auxiliary. These multipurpose pliers are used to place and remove pins, archwires, and other orthodontic materials. How Plier's feature a fine serrated tip for a firm grip with smoothened edges to minimize patient discomfort. These dental pliers are available in Straight and Curved variations.

HOW CROWN PLIERS FEATURES

110 STRAIGHT HOW PLIERS: The Dental Orthodontic Pliers 110 How Crown Straight is hand held instrument designed to grip, bend or hold wires and metal securely. These multipurpose dental pliers great to use to effectively place and remove pins, archwires, and other auxiliaries.

SERRATED TIPS FOR ULTIMATE CONTROL: The round pattern serrated tips on these dental pliers offer more control for precise and easy operation.

DURABLE STAINLESS STEEL CONSTRUCTION: The Premium Quality Stainless Steel composition of these How Crown Straight Dental Pliers offers a long life of strong performance, as they are built to last. These pliers will not bend or break on you. Made with anti-corrosion properties.

FULLY REUSABLE + AUTOCLAVABLE:



WI-25817410

J2 Japanese Stainless Steel Surgical Instrument Set

\$12.00 per piece

Trade price from **\$9.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture J2 Japanese Stainless Steel Surgical Instrument Set for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$12.00, falling to \$9.00 at 10 pieces — about 25% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003346

Jarabak Light Wire Bending Plier 5.5" for Wire Up to .020"

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Jarabak Light Wire Bending Plier 5.5" Serrated Grip With Grooves to Assists in Precise Wire Bending and Loop Closing on Stainless Steel Wire Up to .020" (.51mm) Premium Quality Orthodontic Instrument.

designed and crafted from premium quality stainless steel, ensuring durability, strength, and resistance to corrosion.

JARABAK LIGHT WIRE BENDING PLIER FEATURES

- 5.5" (14cm) Overall Length.
- Light Wire Bending and Loop Forming Up to 0.020" Wire.
- Can Cuts Lightweight Wires Up to 0.025".
- Fine Serrations Tips Provide Firm Grip.
- Grooves Ensure Precise Bending and Closing of Loops.
- Tapered Beak Provide Precise, Multipurpose Functionality.
- Ideal for Light Wire Bending, Loop Closing, and Cutting.

Jarabak Pliers with cutter allow you to have a non-slip grip on the wire to form more intricate bends and loops. Precise and versatile, these two-in-one orthodontic pliers are used to make accurate bends, close loops and cut light wires. Precision grooves ensure uniform loop bending and closing, while full tool steel inserts prevent splay and fine serrations provide a firm grip. Jarabak Plier designed forming wire up to 0.020" and the cutting of wire up to 0.025" with the same instrument. Jarabak pliers are expertly





WI-34568261

Jewelry Nose Plier Beading Wire Shaping

\$2.73 per piece

Trade price from \$1.50
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Jewelry Nose Plier Beading Wire Shaping — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.73, falling to \$1.50 at 5 pieces — about 45% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34511497

Jewelry Nose Plier Fine Precision Jewelry

\$2.73 per piece

Trade price from \$1.50
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Jewelry Nose Plier Fine Precision Jewelry is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material. The name tells you what you are getting: a fine pattern for delicate tissue.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.73, falling to \$1.50 at 5 pieces — about 45% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34490662

Jewelry Nose Plier Fine Precision Jewelry Making Tool

\$2.73 per piece

Trade price from **\$1.50**
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Jewelry Nose Plier Fine Precision Jewelry Making Tool is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material. The name tells you what you are getting: a fine pattern for delicate tissue.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.73, falling to \$1.50 at 5 pieces — about 45% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34519489

Jewelry Nose Plier Fine Tip Precision Jewelry Making Pliers

\$2.73 per piece

Trade price from **\$1.50**
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Jewelry Nose Plier Fine Tip Precision Jewelry Making Pliers — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material. The name tells you what you are getting: a fine pattern for delicate tissue.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.73, falling to \$1.50 at 5 pieces — about 45% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34557391

Jewelry Nose Plier Fine Tip Precision Making Pliers

\$2.73 per piece

Trade price from **\$1.50**
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Jewelry Nose Plier Fine Tip Precision Making Pliers, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material. The name tells you what you are getting: a fine pattern for delicate tissue.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.73, falling to \$1.50 at 5 pieces — about 45% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34504562

Jewelry Nose Plier Jewelry Repair and Crafting Tool

\$2.73 per piece

Trade price from **\$1.50**
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Jewelry Nose Plier Jewelry Repair and Crafting Tool for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material.

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.73, falling to \$1.50 at 5 pieces — about 45% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Forging, milling, heat treatment, hand finishing and passivation all happen

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34492621

Jewelry Nose Plier Mini Precision Tool for Jewelry Making

\$2.73 per piece

Trade price from **\$1.50**
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Jewelry Nose Plier Mini Precision Tool for Jewelry Making for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.73, falling to \$1.50 at 5 pieces — about 45% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34529523

Jewelry Nose Plier Stainless Steel Precision Jewelry Making Tool

\$2.73 per piece

Trade price from **\$1.50**
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Jewelry Nose Plier Stainless Steel Precision Jewelry Making Tool, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$2.73, falling to \$1.50 at 5 pieces — about 45% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34541440

Jewelry Nose Plier Stainless Steel Precision Pliers

\$2.73 per piece

Trade price from \$1.50
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Jewelry Nose Plier Stainless Steel Precision Pliers — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.73, falling to \$1.50 at 5 pieces — about 45% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29053358

Jewelry Pliers Kit Silver

\$15.00 per piece

Trade price from \$12.00
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Jewelry Pliers Kit Silver is a pliers & cutter we forge, grind and finish in our own Sialkot workshop. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$15.00, falling to \$12.00 at 2 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28994443

Jewelry Pliers Kit with Carbon

\$15.00 per piece

Trade price from \$12.00
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Jewelry Pliers Kit with Carbon, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$15.00, falling to \$12.00 at 2 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003345

Johnson Contouring Plier 5.5" Dental Orthodontics Instrument

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Crown & Band Contouring Plier or Johnson Contouring Pliers 5.5" with Concave & Convex Tips German Grade Stainless Steel Dental Orthodontics Instruments

Crown & Band Contouring Plier is also called as Johnson Contouring Pliers used to contour and shape crowns and bands and dealing with fitting Stainless Steel crown and bands. Can also be used to make minor adjustment in partials and retainers. The Crown Contour Pliers is designed to add contour or curvature when molding the crowns during the fit test. This contouring and molding ensures the adjustment of the crown on the patient's tooth. The stainless steel construction ensures maximum durability and dependability. You can

reuse these Johnson contouring plier over time as they will not bend or tarnish on you.



JOHNSON CONTOURING PLIER FEATURES

- 5.5" (14cm) Overall in Length.
- One Concave and One Convex Tips.
- For Contouring & Adapting Stainless Steel Band and Crown.
- Box Joint Style Pliers.
- Ergonomic Smooth Handles.
- Made With High Quality Stainless Steel

PACKAGE INCLUDE

- 1 Pc Crown & Band Contouring Plier 5.5"

TECHNICAL INFORMATION





WI-70003343

Kim Plier Orthodontic Dental

\$12.05 per piece

Trade price from \$11.09

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Kim Plier Orthodontic Dental is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$12.05, falling to \$11.09 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003342

Lab & Office Plier #121 Serrated Dental Orthodontics

\$10.99 per piece

Trade price from \$10.11

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lab & Office Plier #121 Straight Serrated Tips Wire Bending Utility Plier Premium Quality Stainless Steel Dental Orthodontic Instruments.

The #121 Lab and Office Plier is a stainless steel dental utility plier designed for versatile use in dental lab and office settings, allowing gripping, bending, cutting, or holding wires and metal. This lab & office orthodontic plier features a box lock design with straight serrated beaks, preventing hand fatigue and ensuring a secure grip for general use. Lab & office pliers #121 are capable of being autoclaved for use in sterile environments.

LAB & OFFICE PLIER #121

- Overall Length: 5.5".
- Straight Serrated Tips.


- Box Lock Design Prevents Hand Fatigue.
- Smooth Handle for Easy Gripping.
- Ideal for Gripping, Bending, Cutting, and Holding Wires/Metal.
- Made From High-quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Lab & Office Plier #121 Serrated

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENT TYPE: Orthodontic Plier
APPLICATION: Wire Bending Plier
REUSABLE: Yes
AUTOCCLAVABLE: Yes, Must Autoclave Before Use
RUST PREVENTION PROCEDURE: Passivated
ULTRASONIC CLEANED: Yes





WI-35035702

Large Stainless Steel Wire Mesh Tray Reusable Surgical

\$72.00 per piece

Trade price from **\$66.30**

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Large Stainless Steel Wire Mesh Tray Reusable Surgical is part of our surgical instruments range, produced start to finish at our Sialkot facility. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$72.00, falling to \$66.30 at 10 pieces — about 8% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001739

Leaf Masher Plier 6" with PVC Grip for Glass Work

\$9.99 per piece

Trade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Glassblowing Boro Lampwork Bead Leaf Masher Flat Plier 5.5" (14cm) Stainless Steel with PVC Grip – Great for Bead Work, Gemstones and Jewelry Design, Flameworking Tools.

Leaf Masher is great for special applications such as the techniques of glass-blowing, lamp-working, flame-burning, torch-working. Angled & shaped leaf-mashing pliers are great for creating textured leaves with respect to glass creation. Rigid Stainless Steel Construction for Durability

LEAF MASHER PLIER FEATURES

- Angled Leaf Pattern Mashing Plier
- Anti-slip PVC Grip Handle – Black
- Great For Making Leaves with Serious Details

- Perfect for Glassblowing Boro Lampwork, Bead Work, Gemstones and Jewelry Design.
- Made with High Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Leaf Mashers Plier 5.5" (14cm) with Black PVC Grip Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Flameworking Tools


APPLICATION: Glassblowing Boro Lampwork Bead Making

TIPS: Angled & Shaped Leaf Details

HANDLE: PVC Grip Handle

REUSABLE: Yes

RUST PREVENTION PROCEDURE: Passivated





WI-70003341

Ligature forming Plier Orthodontic

\$14.81 per piece

Trade price from \$13.63
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Ligature forming Plier Orthodontic is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$14.81, falling to \$13.63 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001951

Light Wire Bird Beak Plier with Groove 5.25" (13.3cm)

\$10.99 per piece

Trade price from \$10.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Light Wire Bird Beak Plier with Groove. A longer cone than a standard bird beak is useful for forming precise loops in light wires up to 0.20" (0.51mm). The groove allows for locating and duplicating loops.

Bird Beak Plier is the most commonly used orthodontic instrument used to bend and shape arch wires in various thickness. It can do a variety of wire handling and bending, as well as creating steps and loops. It can also be used in forming retentive clasps for removable devices. These dental orthodontic pliers are available in different patterns and sizes.

BIRD BEAK PLIERS FEATURES

LIGHT WIRE BIRD BEAK PLIER: Light Wire Bird Beak Plier with Groove is

used to give sharp bends, it has 1 flat ended beaks with grooves in the inner side if the beaks which help in holding the wire in place. The Long Bird Beak Dental Pliers are designed to grip the wire upto 0.20" (.51mm) can easily bend or hold wires on position and metal securely.

DURABLE STAINLESS

STEEL CONSTRUCTION: These pliers will not bend or break. Manufactured with anti-corrosion properties and extra hard steel for optimum performance.

FULLY REUSABLE + AUTOCLAVABLE: A hygienic and valuable option, these orthodontic loop forming pliers can be fully sterilized without fear, for long-term use without degradation.

TOUGH TAPERED NOSE FOR PRECISE LOOPS: The hard tapered nose of these dental pliers offers maximum precision



**WI-70000431****Lingual Arch Forming Plier 5" for Lingual Arch Wires Bending Up to .036"****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lingual Arch Forming Plier 5" (12.7cm) With 4 Precision Grooves to Create Double Back and Triple Back Bends in Lingual Arch Wires Up to 0.9mm in Diameter Dental Orthodontic Instrument.

Lingual arch forming pliers are orthodontic tools used to create double-back and triple-back bends in lingual arch wires, crucial for fitting and adjusting lingual braces that are placed on the interior surfaces of teeth. This orthodontic plier is designed to precisely and efficiently shape arch wires with wire thicknesses ranging from 0.030" to 036". The orthodontic lingual arch forming plier, measuring 5" in length, features 4 precision grooves in tips for double and triple back bends. Box joint and ergonomic smooth

handles allow better gripping and ease of handling. Made with high quality stainless steel are fully sterilizable for repeated use.

LINGUAL ARCH BENDING PLIER FEATURES:

- 5" Overall Length.
- 4 Precision Groove Tips for Double & Triple Back Bends Forming.
- Box Joint Enables Easy Opening and Closing.
- Ergonomic Smooth Handles for Better Grip.
- Fully Sterilizable for Repeated Use.
- Made From High-quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Lingual Arch Forming Plier 5" (12.7cm)

**WI-70001281****Long Nose Fisherman Fishing Plier 8" Serrated Jaws with Black Grip Handle****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Fisherman Fishing Pliers 8" (20.3cm) – Premium Quality Stainless Steel Long Nose Fishing Plier, Deep-Seated Hooks Remover Long Serrated Jaws with Build in Line Cutter Corrosion Resistant Fishing Gear

The Long nose fisherman fishing plier 8" comes with long serrated jaws to help remove deep-seated hooks. Made of stainless steel, built in line cutter great for fishing, making jewelry and wiring. Stainless steel needle nose pliers are integrated with a wire cutter for cutting wires and are designed to aid in wrapping or looping in confined spaces.

This pliers can help you save time and energy, allowing you to better enjoy the fun of fishing and outdoor travel. The top high-quality stainless steel

construction with PVC grip handle ensure good hardness and corrosion resistance. These pliers for fishing measures 8 inch long, and have PVC gripped handles with finger rest for a comfortable use.

FISHING PLIER FEATURES

- 8" Long Nose Fishing Plier
- Serrated Jaws for Deep-Seated Hooks Remover
- Extra Sharp Line Cutters
- Saltwater Corrosion Resistant
- Lightweight
- Black PVC Grip Handle with Finger Rest
- Made with Stainless Steel

PACKAGE INCLUDES:

- 1 Pc Fisherman Fishing Pliers 8" Serrated Jaws with Line Cutter –



WI-25767353

Luxury Stainless Steel Clip Pliers

\$3.50 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Luxury Stainless Steel Clip Pliers — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material. Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.50, falling to \$2.00 at 2 pieces — about 43% off the single-unit price. Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36335006

Major Bowl Stainless Steel Surgical Bowl Medical Utility

\$45.00 per piece

Trade price from \$32.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Major Bowl Stainless Steel Surgical Bowl Medical Utility is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material. Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$45.00, falling to \$32.00 at 10 pieces — about 29% off the single-unit price. Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000429

Male Torquing Plier 5" (12.7cm)

\$10.99 per piece

Trade price from \$10.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Male Torquing Plier 5" (12.7cm) Dental Torque Adjustment Arch Wire Torquing Pliers German Grade Stainless Steel Dentist Orthodontic Instruments

Male Torquing orthodontic pliers are most often used to place torque into an archwire for single tooth movement. These pliers come in handy when your patient has a row of teeth that needs an extra effort to straight or align. To achieve desired torque, torquing orthodontic pliers used in conjunction with the pliers for torquing or a male and female torquing pliers set. Male torquing pliers are able to apply necessary torque without nicking or fraying the wire.

Our stainless steel torquing plier male / inner feature a streamlined design, sturdy construction and improved

system that you'd expect in an orthodontic torquing tool made for professional orthodontists.


TORQUING PLIER MALE/INNER FEATURES

- Torquing Plier Male / Inner.
- 5" (12.7cm) Overall Length.
- Perfect for Applying Torque Without Nicking or Fraying the Wire.
- Durable, Strong, and Corrosion Resistant.
- Fully Autoclavable/Reusable.
- Made with German Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Male Torquing Plier 5" (12.7cm)

TECHNICAL INFORMATION



WI-37716307

Manufacturing Hair Extension Pliers Kit Surgical

\$15.00 per piece

Trade price from \$12.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Manufacturing Hair Extension Pliers Kit Surgical is a pliers & cutter we forge, grind and finish in our own Sialkot workshop. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$15.00, falling to \$12.00 at 1 set — about 20% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003338****Marcotte Looping Plier 5.25\"/>****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Marcotte Looping Pliers 5.25" German Grade Stainless Orthodontic Dental Instruments for Precise Forming Of Helical Loops And T Loops On Wires Up To .022" X .025" (.053mm X .064mm). No Right Angles. Prevents Stress Fracture Of Wires.


amount of fatigue for the operator due to its ergonomically designed handles and box joint design which makes the opening and closing movements extremely easy.

MARCOTTE LOOPING FORMING PLIERS FEATURES

- 3 Different Cone Sizes For Versatility
- Ergonomically Designed Handles To Decrease Fatigue
- Light Weight For Maximum Comfort
- Fully Autoclaved
- Made with anti-corrosion German Grade Stainless Steel

PACKAGE INCLUDE

- 1 Pc Marcotte Looping Forming Pliers 5.25", Wire sizes up to .022" x .025" (.56 mm x .64 mm)

**WI-70003296****Mathieu Ligature Elastic Placing Plier 5.5\"/>****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Curved Serrated Tip Mathieu Ligature Elastic Placing Plier 5.5" (14cm) Needle Holder Driver, Ligature Separator Placer Premium Quality Stainless Steel Elastomeric Orthodontic Instruments

superior grip and holding of ligatures and elastics.


MATHIEU LIGATURE PLIER FEATURES

PREMIUM QUALITY MATHIEU PLIERS 5.5": These Ligature Elastic Placing Pliers and Needle Holder is a high-value Orthodontic Instruments used for various applications such as elastic placing and wire holding.

STAINLESS STEEL CONSTRUCTION: The Stainless Steel construction ensures maximum durability and dependability. You can reuse these pliers over time as they will not bend or tarnish on you.

FULLY AUTOCLAVABLE: A hygienic choice for the orthodontic practice, these Mathieu Ligature Elastic Placing



**WI-70003378**

Mathieu Ligature Elastic Placing Plier 5.5" Slim Narrow Tip

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Slim Narrow Tip Mathieu Ligature Elastic Placing Plier 5.5" (14cm) Needle Holder Driver, Holding Suture Needles Ligature Separator Placer, Elastomeric Orthodontic Instruments.

Dental Serrated Slim Narrow Tip Mathieu Ligature plier specially used to holding the extra fine beak also can be used for ligature ties on tooth braces and holding the light wires moreover placing the elastics comfortably. Mathieu Ligature Placing orthodontic plier is designed for placing, holding, and tying ligatures and elastics. These dental pliers are also suitable for guiding, grasping a needle when suturing. Mathieu Orthodontic Pliers are available in various tips for specific size ligatures and elastics and feature extra


fine serrations for superior grip and holding of ligatures and elastics.

MATHIEU LIGATURE PLIER FEATURES

PREMIUM QUALITY MATHIEU PLIERS 5.5": These Ligature Elastic Placing Pliers and Needle Holder is a high-value Orthodontic Instruments used for various applications such as elastic placing and wire holding.

STAINLESS STEEL CONSTRUCTION: The Stainless Steel construction ensures maximum durability and dependability. You can reuse these pliers over time as they will not bend or tarnish on you.

FULLY AUTOCLAVABLE: A hygienic choice for the orthodontic practice, these Mathieu Ligature Elastic Placing

**WI-70003337**

Mathieu Ligature Elastic Placing Plier 5.5" Wide Tip

\$10.61 per pieceTrade price from **\$9.76**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Wide Tip Mathieu Ligature Elastic Placing Plier 5.5" (14cm) Needle Holder Driver, Holding Suture Needles Ligature Separator Placer, Elastomeric Orthodontic Instruments.


Dental Elastic Placing Mathieu Plier with wide serrated tips makes this ideal for gripping and placing elastic auxiliaries. Mathieu Ligature Placing orthodontic plier is designed for placing, holding, and tying ligatures and elastics. These dental pliers are also suitable for guiding, grasping a needle when suturing. Mathieu Orthodontic Pliers are available in various tips for specific size ligatures and elastics and feature extra fine serrations for superior grip and holding of ligatures and elastics.

MATHIEU LIGATURE PLIER FEATURES

PREMIUM QUALITY MATHIEU PLIERS 5.5": These Ligature Elastic Placing Pliers and Needle Holder is a high-value Orthodontic Instruments used for various applications such as elastic placing and wire holding.

STAINLESS STEEL CONSTRUCTION: The Stainless Steel construction ensures maximum durability and dependability. You can reuse these pliers over time as they will not bend or tarnish on you.

FULLY AUTOCLAVABLE: A hygienic choice for the orthodontic practice, these Mathieu Ligature Elastic Placing Pliers are fully autoclavable and corrosion-resistant. They are easy to clean and use over and over for the

**WI-70000421**

Mini Bracket Remover Plier 5.5" Orthodontic Instrument

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Mini Bracket Remover Plier 5.5" (14cm) with Precision Gripping Beak for Braces Removing Premium Quality Stainless Steel Dental Orthodontic Instruments

Micro Mini bracket remover plier is a precision-engineered instrument designed for efficient bracket removal on anterior teeth. The streamlined design allows for visibility and patient comfort, and a scraping edge removes bonding adhesive for clean surfaces. Mini pliers are used to gently remove brackets by gripping the upper and lower wings between their jaws and pressing the handle. Narrow sharp tips allow for safe removal of direct bond brackets between the bracket base and tooth surface edges. Made with anti-corrosion premium quality stainless steel. This plier is autoclavable, enabling

easy sterilization between uses, a crucial aspect for maintaining a hygienic practice environment.

PACKAGE INCLUDE

- 1 Pc Mini Bracket Remover Plier 5.5" with Precision Gripping Beak

TECHNICAL INFORMATION

MATERIAL: Stainless Steel**INSTRUMENTS TYPE:** Orthodontic Pliers**APPLICATION:** Bracket Removing/Braces Removal**JAWS:** Bracket Gripper**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**SURFACE TYPE:** Polished**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**WI-70003238**

Mycotic Toe Nail Nipper 6" Double Action

\$14.99 per pieceTrade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty Mycotic Toe Nail Clipper 6" (15.2cm) Double Action Barrel Spring Model Toenail Cutter German Grade Stainless Steel Dermatology Podiatry Instruments

Mycotic Toenail Nippers Podiatrist Pedicure Nippers for Toenail Fungus Surgical Stainless Steel Ingrown Toenail Clippers Thick Toenails Fingernail Clipper Cutter Trimmer. Besides the pedicurists podiatrists elderly care professionals, Seniors, home personal care, this nail nipper is especially designed for the who are suffering from fungus infections. Regular toenail clippers are proven to be ineffective painful and cause even more damage. Fungus infections nails requires a right tool that can do the job and prevent

injury. Ideal for cutting of problematic nails.

MYCOTIC TOE NAIL NIPPER FEATURES

MYCOTIC TOE NAIL CLIPPERS: The Mycotic Nipper (also known as "Cantilever Toe Nail Clipper") is the most durable pedicure nipper designed to work on extremely thick hard damaged and fungus infected toenails. Mycotic Double Action heavy-duty toenail Nippers are designed for thick and tough fingernails and toenails. The podiatrist toenail clippers are ideal for both professionals and home use – for Men, Women and Seniors. Effortlessly and safely cut thick toenails, ingrown nails and hangnails with the razor sharp edge and double barrel spring clippers.

PREMIUM QUALITY STAINLESS

**WI-70001313**

Mycotic Toe Nail Nipper 6" Double Action Black with Gold Spring

\$15.99 per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Mycotic Nail Nipper 6" Sharp Concave Jaws Double Action Black Color Cantilever Nail Clipper with Gold Spring Barrel Spring German Grade Stainless Steel Toenail Nipper Clipper Podiatry Pedicure Tool for Thick and Tough Toenail Fungus Nails

Besides the pedicurists podiatrists elderly care professionals, Seniors, home personal care, this heavy duty cantilever toe nail clipper is especially designed for the who are suffering from fungus infections. Regular toenail clippers are proven to be ineffective painful and cause even more damage. Fungus infections nails requires a right tool that can do the job and prevent injury. Ideal for cutting of problematic nails.

MYCOTIC TOE NAIL NIPPER FEATURES

PROFESSIONAL CANTILEVER NAIL CLIPPER – BLACK: The Mycotic Nipper (also known as "Double Action Cantilever Nail Clipper") is the most durable pedicure nipper designed to work on extremely thick hard damaged and fungus infected toenails. Mycotic Double Action heavy-duty toenail Nippers are designed for thick and tough fingernails and toenails. The podiatrist toenail clippers are ideal for both professionals and home use – for Men, Women and Seniors. Effortlessly and safely cut thick toenails, ingrown nails and hangnails with the razor sharp edge and double barrel spring clippers.

PREMIUM QUALITY STAINLESS STEEL: Forged from O.R. Grade

**WI-70003336**

Nance Loop Forming Plier 5" for Wires Up to .022"

\$10.99 per pieceTrade price from **\$10.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Nance Loop Forming Plier 5" (12.7cm) Long Wire Bending Orthodontic Pliers For Precise Forming of Loops in Four Lengths: 3, 4, 5 and 6mm on Stainless Steel Wire Up to .022" (.56 mm) – Premium Quality Stainless Steel Dental Orthodontic Instrument.

Nance Loop Forming Pliers are precision dental tools for shaping arch wires with a four-step beak design that creates loops in sizes of 3mm, 4mm, 5mm, and 6mm with a groove at the base to aid in closing loops. This plier is capable of handling stainless steel wires up to .022 inches (0.56 mm) thick. Nance pliers have beveled edges to prevent wire scoring, making them ideal for orthodontic procedures. Nance loop forming wire bending plier with ergonomically designed handles

and a unique rivet screw design, make opening and closing movements easy, ensuring optimal wire bending and minimal operator fatigue. Our Nance Loop Forming Plier, crafted from high-quality stainless steel, is both sterilizable and reusable.

NANCE LOOP FORMING WIRE BENDING PLIER FEATURES

- Crafted with 5" (12.7cm) Overall Length.
- Four-step Beak Design Creates Loops in Sizes 3.0 – 6.0 Mm.
- Groove at Base Aids in Closing Loops.
- Handles Stainless Steel Wires Up to .022" (.56 mm)
- Beveled Edges Prevent Wire Scoring.



**WI-70000430****O'Brien Loop Forming Orthodontic Pliers****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

O'Brien Pliers 5" (13cm) for Precise Forming of Small Helical Loops and Intraoral Tip Back Bends on Stainless Steel Wires Up to .020" (.51 mm). Premium Quality Orthodontic Instrument.


O'Brien plier used for creating precise, small helical loops and tip-back bends in orthodontic wires up to 0.51 mm in diameter. O'Brien loop forming plier features tapered and conical tips for precise wire manipulation, micro-serrated beaks to prevent slippage during bending, and a graduated cone for precise omega and closing loops. These pliers are specialized for creating tip-back bends, forming small helical loops, and manipulating and contouring archwires to adjust tooth

position, angle, wire shape, and deliver corrective forces.

This plier also suitable for bending square or rectangular archwires, utility archwires, detorquing ends, and bending 90° angles. Made with premium grade stainless steel, features a rivet joint for smooth movements and an ergonomic handle for excellent performance during surgical procedures.

O'BRIEN ORTHODONTIC PLIERS FEATURES

- 5" (13cm) Overall Length.
- Max. Wires Up to .020" (.51 mm) in Diameter.
- Tapered & Conical Tips.
- 3 Different Cone Sizes to Suit All Applications.

**WI-70003335****Omega Loop Forming Plier 5.5" Tweed Style Helical Tip****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Omega Loop Forming Plier 5.5" (14cm) Tweed Style, Serrated & Helical Tip with L-Key Joint for Wires Up to 0.56x0.64mm in diameter. German Grade Stainless Steel Dental Orthodontic Plier Dentist Oral Braces Instruments


Tweed Omega Loop Forming Pliers are a sturdy handheld orthodontic tool designed for making bull and helical closing loops in square or rectangular wire. Helical beak has three step sections. Cylindrical beak is electro-etched to prevent wire slippage. Ingroove beak has slight parallel serrations to hold wire at 90° angle while forming loops.

Omega loop forming plier comes with 3 different-sized barrel shaped endings. This enables the operator to create

omega loops of graduated cones .045", .060", and .075" (1.14 mm, 1.52 mm, and 1.91 mm) diameters are designed to form precise omega loops on wires up to .022" x .025" (.56 mm x .64 mm) diameters as desired.

OMEGA LOOP FORMING PLIER FEATURES

- Overall Length 5.5" (14cm)
- For Wires Upto .022" x .025" (.56 mm x .64 mm) Diameters
- 3 Different Cone Sizes For Versatility
- Serrated Opposing Beak For Maximum Wire Control
- Ergonomically Designed Handles To Decrease Fatigue
- Light Weight For Maximum Comfort
- Fully Autoclaved





WI-70001308

Optical Pliers 5.25" for Helical Closing Loops & Arch Adjustments

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Optical Pliers For Forming Helical Closing Loops & Arch Adjustments
Designed for All Wires Up to 0.028"
German Grade Stainless Steel Dental Orthodontic Instrument

Orthodontic optical pliers features one concave and one convex beak to create loops in wires. These precisely handcrafted optical pliers have finely finished smooth edges allowing you effortlessly form accurate loops and bends in wires up to 0.028" without the risk of nicking the wire. Built for versatility, our optical pliers can also be used to create loop ties, gable bends and form stops in wires. Convex-concave design adjusts all types of loops, providing your patient with maximum comfort. No more impingement.

ORTHODONTIC OPTICAL PLIERS FEATURES


- One Concave and One Convex Beak
- Designed for All Wires Up to 0.028"
- Lap-joint style construction
- Ergonomically Designed Handles To Decrease Fatigue.
- Light Weight For Maximum Comfort.
- Fully Autoclaved.
- Made with anti-corrosion German Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Optical Plier 5.25" (13.3cm)

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENT TYPE: Orthodontic Plier
APPLICATION: For Forming Helical Closing Loops & Arch Adjustments



WI-70003334

Orthodontic Tweed Rect. Arch forming Plier

\$15.30 per piece

Trade price from \$14.08
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Orthodontic Tweed Rect. Arch forming Plier is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$15.30, falling to \$14.08 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003282

Orthopedic Pliers Orthodontic

\$138.75 per piece

Trade price from **\$127.65**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Orthopedic Pliers Orthodontic is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$138.75, falling to \$127.65 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46769520

Perry Dressing Plier Surgical Instrument (Model Perry Dressing Plier-01)

\$3.00 per piece

Trade price from **\$1.50**
Minimum order 30 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

We manufacture Perry Dressing Plier Surgical Instrument (Model Perry Dressing Plier-01) for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.00, falling to \$1.50 at 30 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46758578

Perry Dressing Plier Surgical Kit for Fine Tissue

\$1.50 per piece

Trade price from \$0.90
Minimum order 500 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Perry Dressing Plier Surgical Kit for Fine Tissue for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material. The name tells you what you are getting: a fine pattern for delicate tissue.

Supplied as a set — the contents can be adjusted to your list. We work the

instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.50, falling to \$0.90 at 500 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000370

Pin Puller Pliers 5.5" (14cm) with Groove Serrated Jaws

\$10.99 per piece

Trade price from \$10.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pin Puller Pliers 5 1/2" (14cm) with Groove Serrated Jaws Smooth Grip Handle Premium Quality Stainless Steel Surgical Orthopedic Instruments

The pin puller is a sturdy instrument used in bone fracture and stabilization procedures, used to remove surgical pins after proper bone treatment. Pin pulling pliers used in many orthopedic fracture repair procedures. These pliers have grooved jaws that provide a firm fitting grip that easily pulls the pins from bones. Smooth handle offers users superior comfort while in use and the serrated jaw allows for better grip. Made with premium quality stainless steel making it an efficient and time-saving tool that last for many years and reused after sterilization.

PIN PULLER PLIERS FEATURES

- 5.5" (14cm) Overall Length.
- Grooved Serrated Jaws for Secure Grip.
- Smooth Grip Handles.
- Performance Tested Surgical Pin Removal Tool.
- Reusable.
- Corrosion-resistant.
- Premium Grade Stainless Steel Construction.


PACKAGE INCLUDE

- 1 Pc Pin Puller Pliers 5 1/2" (14cm) with Groove Serrated Jaws

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENT TYPE: Orthopedic Instruments

APPLICATION: Bone Pin & Wire Cutting
REUSABLE/AUTOCCLAVABLE: Yes





WI-26032655

Podiatry Tool Kit 5 PCS Set Curved

\$25.00 per piece

Trade price from **\$15.00**
Minimum order 5 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Podiatry Tool Kit 5 PCS Set Curved, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material. The name tells you what you are getting: a curved working end for reach around structures. Supplied as a 5-piece set.

We can add, drop or substitute items so the tray matches your protocol rather

than ours. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$25.00, falling to \$15.00 at 5 sets — about 40% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000938

Posterior Band Removing Pliers 5.5" Long

\$10.99 per piece

Trade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Posterior Band Removing Pliers Long 5.5" (14cm) with Nylon Padded Tip Premium Quality Stainless Steel Orthodontic Pliers Dental Surgical Instruments.

Posterior band removing plier designed to efficiently remove posterior bands with minimum patient discomfort. This posterior band removing plier featuring a sleek long tip design which provide greater leverage during usage. The long tips allowing for a more easy access to molar region and greater visibility during operation, enabling orthodontists to remove bands with minimal effort while their plastic padded tip ensures maximum ease and comfort for patients. These orthodontic pliers can be fully autoclaved without

fear, for long term use without degradation.

POSTERIOR BAND REMOVING PLIERS FEATURES


- 5.5" (14cm) Overall Length.
- Long Tip Design.
- Nylon Padded Tip.
- Lap-Joint Style Construction.
- Smooth Grip Handle.
- Fully Autoclavable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Posterior Band Removing Plier 5.5" (14cm) Long Tip

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Orthodontic



**WI-70000939****Posterior Band Removing Pliers 5.5\"/>****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Posterior Band Removing Pliers Short 5.5\" (14cm) with Nylon Nylon Paded Tip Premium Quality Stainless Steel Orthodontic Pliers Dental Surgical Instruments.

Posterior band removing plier designed to efficiently remove posterior bands with minimum patient discomfort. This posterior band removing plier featuring a short tip design which provide greater leverage during usage. The short tips allowing for a more easy access to molar region and greater visibility during operation, enabling orthodontists to remove bands with minimal effort while their plastic padded tip ensures maximum ease and comfort for patients. These orthodontic pliers can be fully autoclaved without

fear, for long term use without degradation.

POSTERIOR BAND REMOVING PLIERS FEATURES

- 5.5\" (14cm) Overall Length.
- Short Tip Design.
- Nylon Padded Tip.
- Lap-Joint Style Construction.
- Smooth Grip Handle.
- Fully Autoclavable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE


- 1 Pc Posterior Band Removing Plier 5.5\" (14cm) Short Tip

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Orthodontic

POSTERIOR BAND REMOVING PLIERS



POSTERIOR BAND REMOVING PLIERS

**WI-70001307****Primary Crown Crimping Pliers 5\"/>****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Primary Crown Crimping Pliers 5\" Slim-Profile Tips with Spring Action Handle Stainless Steel German Grade Dental Plier Orthodontics Instruments

Primary Crown Crimping Pliers used to crimp bands and crowns to ensure a tight fit. Crimping orthodontic pliers have a variety of uses including placing stops and tubing on arch wires, crimping bands and gingival margins of crowns to insure a tight fit, and placing and crimping hooks. These multi-purpose orthodontic plier is precisely handcrafted a must-have for every professional orthodontist.

Featuring a slim-profile tip for preciseness, these crown & band crimping pliers allow orthodontists to effortlessly crimp bands and crowns to ensure a tight fit. With its versatile tip

design, these expertly constructed orthodontic pliers can also be used to place stops and tubing on arch wires.

PRIMARY CROWN CRIMPING PLIERS FEATURES

- Slim-Profile Tips
- Tips Sit Together Perfectly.
- Creating Accurate and Precise Crimps
- For Lab and Office Use
- Ergonomically Designed Handles To Decrease Fatigue.
- Spring Action Handle.
- Fully Autoclaved.
- Made with anti-corrosion German Grade Stainless Steel.

PACKAGE INCLUDE

PRIMARY CROWN CRIMPING PLIERS





WI-25800288

Pro-Quality Stainless Steel Hair Extension Plier Top-Rated Strong Grip

\$3.50 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Pro-Quality Stainless Steel Hair Extension Plier Top-Rated Strong Grip — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$3.50, falling to \$2.00 at 2 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000860

Rainbow 3 Prong Dental Aderer Plier 5"

\$10.99 per piece

Trade price from \$10.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Rainbow Dental Aderer Plier 5" (12.7cm) Three Prong Orthodontic Wire Bending Triple Beak Pliers for Light Wires upto 0.030" (0.76 mm) Stainless Steel Dental Orthodontic Instruments.


Aderer pliers are also called three-jaw or three-prong pliers and are designed for consistent bends of wires, precision aligned tips, and tapered prongs for fine wire bending. Tips are rounded for ease of contouring light wires up to 0.030" (0.76 mm). 3 prong orthodontic Aderer plier useful for grasping wires and clasps for intraoral arch wire adjustments as well as for lingual arch adjustments. Three prong orthodontic

pliers can have a wide range of beak characteristics, and may be long and tapered for hard to reach quadrants or short and stubby for firmer control. They also can be offset, double sided, mini, Aderer style, round pronged or a combination.

Our Rainbow Aderer Pliers made Of high quality stainless steel alloy with exactly machined tool tips and precisely guided sheaths. Ergonomic design with rounded edges for fatigue-free working.

RAINBOW 3 PRONG DENTAL ADERER PLIER FEATURES



**WI-70000858****Rainbow Bracket Removing Plier 5.5" Straight****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Multi Titanium Rainbow Color Bracket Remover Plier 5.5" (14cm) Straight Premium Quality Stainless Steel Braces Remover Bracket Gripper Plier for Dental Orthodontic Procedures. Efficient DIY Tool for DIY Enthusiasts and Crafters.

Our Bracket Remover Plier is a handheld orthodontic instrument designed to quickly and efficiently remove brackets. Cut time in half and get the braces off faster than ever with this handy instrument that is suitable for steel ceramic and plastic brackets of all types. Featuring a vibrant multicolor and rainbow design, these bracket removing pliers not only stand out visually but also symbolize their versatility across different crafting materials. The eye-catching aesthetics

make them a distinctive addition to any crafting project. Bracket Removing Pliers Rainbow serves as a comprehensive DIY tool kit designed for DIY enthusiasts with guarantee long-lasting performance without the risk of bending or breaking.

RAINBOW BRACKET REMOVING PLIER FEATURES

- Crafted with 5.5" (14cm) Overall Length.
- Multi Titanium Coated Rainbow Color.
- Remove Braces Quickly and Efficiently.
- Straight Bracket Gripper Jaws for Cutting, Gripping, Twisting.
- L-Key Screw Adjustment.
- Fully Autoclavable and Reusable.

**WI-70000862****Rainbow Hard Wire Cutter 15 Degree Curved Orthodontic Plier****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Multi Titanium Rainbow Hard Wire Cutter for Wires Up to .022" X .028" (.56mm X .71mm) 15 Degree Offset Curve Tip Premium Quality Dental Orthodontic Pliers.


Hard wire cutter is a specialized orthodontic instrument designed with heavy duty cutting surface which can easily cut wires and pins. Strong enough to bend archwire ends after distal cutting. The Hard wire cutter designed to cut hard wires up to .022" x .028" (.56 mm x .71 mm) with 15 degree-offset tip which enhances access in difficult areas. Our rainbow color hard wire cutter made of high-quality stainless steel. It is highly durable, anti-corrosion, will not bend, and can withstand the daily use of archwires cutting.

RAINBOW HARD WIRE CUTTER FEATURES

- Crafted with 5" (12.7cm) Overall Length.
- Multi Titanium Coated Rainbow Color.
- Cuts Wires Up to .022" X .028" (.56mm X .71mm).
- 15 Degree Offset Curve Tips.
- Rivet Screw Design for Allow Ease of Opening and Closing
- Reusable
- Made with Premium Quality Stainless Steel

PACKAGE INCLUDES:

- 1 Pc Rainbow Color Hard Wire Cutter 15° offset Tip 5" (12.7cm)

TECHNICAL INFORMATION



WI-28898462

Rainbow Multi-Color Hair Extension Pliers Kit Surgical

\$15.00 per piece

Trade price from **\$12.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Rainbow Multi-Color Hair Extension Pliers Kit Surgical, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$15.00, falling to \$12.00 at 1 set — about 20% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35085619

Reusable Stainless Steel Wire Mesh Surgical Tray Medical Sterilization

\$55.00 per piece

Trade price from **\$48.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

We manufacture Reusable Stainless Steel Wire Mesh Surgical Tray Medical Sterilization for general and specialty operating theatre use — made from bar stock to finished instrument in house. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$55.00, falling to \$48.00 at 10 pieces — about 13% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001306

Reynolds Contouring Pliers 5.5" Curved Tips

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Reynolds Contouring Pliers 5.5" (14cm)
Curved Falcon Tips Orthodontic Wire
Bending German Grade Stainless Steel
Dental Pliers 115 Orthodontic
Instrument

The Reynolds contouring dental pliers are two handed instruments designed to grip, bend or hold wires and metal securely. They are used in the orthodontic field, to make minor adjustments to removable prosthetics such as partials and retainers.

Contouring orthodontic pliers are ideal for forming and contouring arch wires, reshaping molar bands and enhancing crown form contours. They also are especially useful for placing accentuated curves in wires.

REYNOLDS CONTOURING PLIERS FEATURES

- Curved Falcon Tips.
- Designed to Contour and Shape Crowns and Bands
- Box Joint Style Construction
- Ergonomically Designed Handles To Decrease Fatigue.
- Fully Autoclaved.
- Made with anti-corrosion German Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Reynolds Contouring Pliers 5.5" (14cm)

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Orthodontic Plier

APPLICATION: To Contour and Shape



WI-25778338

Salon Tools Color

\$3.50 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Salon Tools Color is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$3.50, falling to \$2.00 at 2 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000422****Side Cutter Parallel Action Plier 4.5" (11.4cm) Serrated Jaws****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Parallel Action Wire Cutter 4.5" (11.4cm) with V-Slot Serrated Jaws and Built-in Side Wire Cutter Premium Quality Nickel Plated Steel Wire Bending, Holding, and Cutting Wires Veterinary Orthopedic Surgical Instrument.

Double action parallel wire cutter plier is a hand tool for bending, holding, and cutting small pins and wires, offering a more precise and clean cut. Nickel-plated cutters are used in veterinary orthopedics, electronics, jewelry making, and other wire cutting applications due to their corrosion resistance and durability.

Parallel Action Wire Cutter Pliers feature a compound lever action for a consistent grip, ideal for holding and gripping wire, with wide, flat nose serrated jaws featuring a V-Slot. Our

Side Wire Cutter Parallel Plier features an ergonomically designed beak and plier-style handle, enhancing efficiency and preventing errors in the wire-cutting procedure. Parallel action wire cutter is a versatile tool that combines gripping pliers and built-in side cutters, making it an efficient and time-saving tool. It is quite resilient to last for many years and reused after sterilization.

PARALLEL ACTION WIRE CUTTER FEATURES

- 4.5" (11.4cm) Overall Length.
- Wide, Flat Nose Serrated Jaws With V-slot.
- Compound Lever Action for Consistent Grip.
- Cut Both Hard and Soft Wire With Ease.

**WI-70003330****Slim Weingart Pliers 5.5" Orthodontic Dental Braces Instrument****\$10.99** per pieceTrade price from **\$10.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Slim Weingart Plier with L-Key Dental Wire Bending Orthodontics Braces Placement Stainless Steel Ortho Dental Instrument

Versatile utility pliers with serrated tips hold wires firmly at a convenient working angle. The slim Weingart plier provides easy access to lingual and posterior areas. Tapered beaks fit easily between brackets and hard-to-reach areas.

ADVANTAGES OF PLIERS WITH A HEXAGONAL SCREW/L-KEY

- Precise function.
- Angular Design is convenient for the placement and removal of arch wires.
- Longevity of life

- Less chance of getting loose at the joints even after long usage.
- When found loose, can be tightened with the aid of an Allen-key easily.
- The precision joint with a screw makes it adjustable.

SLIM WEINGART PLIERS FEATURES

PERFECT WEINGART PLIER: Weingart utility orthodontic plier is a versatile serrated tool. The tapered angular design is convenient for placement and removal of arch wires.

WEINGART L-KEY PLIERS WITH HEXAGONAL SCREW: Used to shorten arch wire. Provides an accurate 1mm stop.

EASY TO USE: This Weingart Plier is user-friendly and easy to clean and sterilize.





WI-39828267

Solid Color Stainless Steel Hair

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Solid Color Stainless Steel Hair for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29038420

Stainless Steel Carbon Steel Jewelry Making Pliers Kit

\$15.00 per piece

Trade price from \$12.00
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Stainless Steel Carbon Steel Jewelry Making Pliers Kit — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$15.00, falling to \$12.00 at 2 sets — about 20% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29040382

Stainless Steel Carbon Steel Jewelry Pliers Set

\$15.00 per piece

Trade price from **\$12.00**
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Carbon Steel Jewelry Pliers Set is a pliers & cutter we forge, grind and finish in our own Sialkot workshop. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$15.00, falling to \$12.00 at 2 sets — about 20% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36033077

Stainless Steel Double Drawn Remy Extension Pliers

\$5.00 per piece

Trade price from **\$1.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Double Drawn Remy Extension Pliers for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$5.00, falling to \$1.00 at 10 pieces — about 80% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25748314

Stainless Steel Hair Extension

\$6.00 per piece

Trade price from \$4.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Stainless Steel Hair Extension — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$6.00, falling to \$4.00 at 10 pieces — about 33% off the single-unit price. Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25779279

Stainless Steel Hair Extension Clip Pliers

\$3.50 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Hair Extension Clip Pliers is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$3.50, falling to \$2.00 at 2 pieces — about 43% off the single-unit price. Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25807282

Stainless Steel Hair Extension Plier Clips

\$7.50 per piece

Trade price from \$4.00
Minimum order 20 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Hair Extension Plier Clips — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$7.50, falling to \$4.00 at 20 pieces — about 47% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25760361

Stainless Steel Hair Extension Plier for Hair

\$10.00 per piece

Trade price from \$8.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Hair Extension Plier for Hair is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$10.00, falling to \$8.00 at 2 pieces — about 20% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-39791653

Stainless Steel Hair Extension Pliers

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Hair Extension Pliers is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25841232

Stainless Steel Hair Extension Pliers and Tools

\$3.50 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Stainless Steel Hair Extension Pliers and Tools, made for general and specialty operating theatre use. Instruments of this type grip, bend or cut wire, pins and hard material. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by

hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$3.50, falling to \$2.00 at 2 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25769320

Stainless Steel Hair Extension Pliers Clip Non-Slip Grippd

\$3.50 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Hair Extension Pliers Clip Non-Slip Grippd is a pliers & cutter we forge, grind and finish in our own Sialkot workshop. Instruments of this type grip, bend or cut wire, pins and hard material.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$3.50, falling to \$2.00 at 2 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25792255

Stainless Steel Hair Extension Pliers Hair

\$35.00 per piece

Trade price from \$12.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Hair Extension Pliers Hair is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type grip, bend or cut wire, pins and hard material.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$35.00, falling to \$12.00 at 2 pieces — about 66% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-39817363

Stainless Steel Keratin Tape Sealing Pliers

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Keratin Tape Sealing Pliers for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type grip, bend or cut wire, pins and hard material.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-39810413

Stainless Steel Keratin Tape Sealing Pliers for Hair Extensions

\$15.00 per piece

Trade price from \$12.00
Minimum order 25 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Stainless Steel Keratin Tape Sealing Pliers for Hair Extensions — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$15.00, falling to \$12.00 at 25 sets — about 20% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-37706370

Stainless Steel Micro Bead

\$15.00 per piece

Trade price from \$12.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Micro Bead — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type grip, bend or cut wire, pins and hard material. The name tells you what you are getting: a micro-scale working end for fine and delicate work.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$15.00, falling to \$12.00 at 1 set — about 20% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25827269

Stainless Steel Multi-Function Waterproof Hair Extension Pliers

\$3.50 per piece

Trade price from \$2.00
Minimum order 2 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Multi-Function Waterproof Hair Extension Pliers is a pliers & cutter we forge, grind and finish in our own Sialkot workshop. Instruments of this type grip, bend or cut wire, pins and hard material.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$3.50, falling to \$2.00 at 2 pieces — about 43% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34929439

Stainless Steel Wire Mesh Sterilization Basket

\$6.00 per piece

Trade price from **\$4.60**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Wire Mesh Sterilization Basket — one of 4,770 instruments we manufacture ourselves rather than buy in. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one.

Listed at \$6.00, falling to \$4.60 at 10 pieces — about 23% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35023709

Stainless Steel Wire Mesh Tray Reusable Surgical Instrument Basket

\$55.00 per piece

Trade price from **\$49.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Wire Mesh Tray Reusable Surgical Instrument Basket is part of our surgical instruments range, produced start to finish at our Sialkot facility. The name tells you what you are getting: full reusability — clean, autoclave and use again. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$55.00, falling to \$49.00 at 10 pieces — about 11% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003138

Standard Pin and Ligature Wire Cutter Pliers
Box Joint Dental Orhtodontics Tools

\$14.33 per piece

Trade price from \$13.18
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Standard Pin and Ligature Wire Cutter Pliers Box Joint Dental Orhtodontics Tools is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$14.33, falling to \$13.18 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003340

Steiner Ligature Tying Plier 6" with Tapered
Grooved Tips

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Steiner Ligature Tying Pliers 6" (15.2cm) For Tying Stainless Steel Ligature Wires Up to 0.015" (.38 mm) Premium Quality Stainless Steel Dental Orthodontic Instruments.

Steiner tying plier is a specialized orthodontic tool designed to grip, wrap, twist, and tie fine stainless steel ligature wires up to 0.015" (.38 mm) around orthodontic brackets, securing the main archwire in place. This ligature tying pliers designed with tapered, grooved smooth tips for gripping thin wire effectively without distortion or slippage. They feature a lightweight design that is slimmer than traditional Coon-style tying pliers and are made from surgical-grade, autoclavable stainless steel for easy sterilization and enhanced durability.

STEINER LIGATURE TYING PLIER
FEATURES



- Crafted With 6" (16.5cm) Overall Length.
- For Ligature Wires Up to .015"(.38 Mm).
- Long Slender Beaks with Grooved Blunt Tips.
- Central Spring Reverse-action Plier Handles.
- Reusable.
- Made With Premium Grade Stainless Steel

PACKAGE INCLUDE

- 1 Pc Steiner Ligature Tying Plier

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENT TYPE: Orthodontic Plier



**WI-70001143****T/C Angulated Bracket Remover Plier 5.25" Orthodontic Instrument****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Angulated Bracket Remover Plier 5.25" (13.3cm) Angled Tungsten Carbide Inserts Jaws with L-Key Adjustment, Long Gold Plated Handle Dental Orthodontic Braces Removal Tool for Ceramic + Plastic Brackets.

Tungsten Carbide Angled Bracket Removal Pliers are a sturdy orthodontics handheld instrument designed to quickly and efficiently remove brackets. The high angle + long handle design is effective for removal of all bracket styles, both anterior and posterior. Suitable for steel ceramic and plastic brackets of all types. Cut time in half and get the braces off faster than ever. Our angulated bracket remover plier comes with Tungsten Carbide Tips which provide an extra level of resilience and strength for ultimate

reliability. The premium quality stainless steel and tungsten carbide composition of these angulated bracket removing plier offers a long life of dependable performance.

T/C ANGULATED BRACKET REMOVER PLIER FEATURES:

- Angled for Easy Posterior Reach.
- Tungsten Carbide Inserts Jaws.
- Suitable for Steel, Ceramic or Plastic Brackets of All Types.
- Hexagonal Screw (L-Key) Adjustment.
- Gold Plated Smooth Handle.
- Perfect for Posterior and Anterior Brackets Removal.
- Fully Reusable + Autoclavable
- Made With Premium Quality Stainless Steel.

WI-70001400**T/C Bird Beak Pliers with Cutter 5.5" Orthodontic Instrument****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Bird Beak Plier 5.5" (14cm) with Tungsten Carbide Inserts Cutter Dental Orthodontic Standard Loop Forming Plier for Round Wires up to 0.30" (0.76mm) with Round and pyramid beaks are smooth to prevent wire scoring.

Bird Beak Plier is the most commonly used orthodontic instrument used to bend and shape arch wires in various thickness. It can do a variety of wire handling and bending, as well as creating steps and loops. It can also be used in forming retentive clasps for removable devices. These dental orthodontic pliers are available in different patterns and sizes.



BIRD BEAK PLIER WITH CUTTER: Bird Beak Orthodontic Pliers with T/C Cutter has a smooth bending surfaces for which a range of loop sizes while providing a sharp blade for cutting. They feature precisely aligned hard tips make them the perfect loop forming pliers for round wires up to 0.30" (0.76mm). Designed with a tip to offer a uniform pressure on the wire while looping.

DURABLE STAINLESS

STEEL CONSTRUCTION: These pliers will not bend or break. Manufactured with anti-corrosion properties and extra hard steel for optimum performance. The T/C Tps designed to stay sharp for a significantly longer period of time than standard instruments.

FULLY REUSABLE + AUTOCLAVABLE:

A hygienic and valuable option, these

T/C BIRD BEAK PLIER WITH CUTTER FEATURES

**WI-70001203****T/C Cannulated Pin and Wire Cutter 7.5" End Cutting Upto 1.6mm (1/16")****\$43.99** per pieceTrade price from **\$40.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty T/C Cannulated Pin & Wire Cutter 7.5" (19.1cm), End Cutting Upto 3.2mm (1/8") Soft Wire or 1.6mm (1/16") Hard Wire Gold Color Handle German Grade Stainless Steel Surgical Veterinary Orthopedic Instruments.

Cannulated Pin & wire cutter widely used in fracture repairs and other standard surgical procedures. It is used to cut the inserted pins from the bone.

The unique design and sharp jaws of cannulated allow surgeons for effective cutting and removal of pins. This comes with tungsten carbide inserts blades. Cannulated Pin Cutter is very useful tool for cutting surgical pins. The entire portion of the pinhead is inserted in the cutter and then removed with force. The maximum cutting capacity of this Cannulated pin cutter is around

0.062" in hard wire and 0.125" in soft wire.

TC CANNULATED PIN & WIRE CUTTER FEATURES

- Sharp Cutting Edges.
- Overall Length 7.5" (19.1cm).
- Maximum Cutting Capacity 3.2mm (1/8") Soft Wire or 1.6mm (1/16") Hard Wire.
- T/c (Tungsten Carbide Inserts) Cutting Jaws for Effective Cutting Action.
- Gold Color Plier Grip Handle.
- Reusable.
- Corrosion-resistant.
- German Grade Stainless Steel Construction.

PACKAGE INCLUDE**WI-70001202****T/C Cannulated Pin and Wire Cutter 9" End Cutting Upto 1.6mm (1/16")****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty T/C Cannulated Pin & Wire Cutter 9" (22.9cm), End Cutting Upto 3.2mm (1/8") Soft Wire or 1.6mm (1/16") Hard Wire Gold Color Handle German Grade Stainless Steel Surgical Veterinary Orthopedic Instruments.

Cannulated Pin & wire cutter widely used in fracture repairs and other standard surgical procedures. It is used to cut the inserted pins from the bone.

The unique design and sharp jaws of cannulated allow surgeons for effective cutting and removal of pins. This comes with tungsten carbide inserts blades. Cannulated Pin Cutter is very useful tool for cutting surgical pins. The entire portion of the pinhead is inserted in the cutter and then removed with force. The maximum cutting capacity of this Cannulated pin cutter is around

0.062" in hard wire and 0.125" in soft wire.

TC CANNULATED PIN & WIRE CUTTER FEATURES

- Sharp Cutting Edges.
- Overall Length 9" (22.9cm).
- Maximum Cutting Capacity 3.2mm (1/8") Soft Wire or 1.6mm (1/16") Hard Wire.
- T/c (Tungsten Carbide Inserts) Cutting Jaws for Effective Cutting Action.
- Gold Color Plier Grip Handle.
- Reusable.
- Corrosion-resistant.
- German Grade Stainless Steel Construction.

PACKAGE INCLUDE

**WI-70001288****T/C End Cutter Double Action 8.5" Max Cutting Capacity 2.4mm (3/32")****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Front End Tungsten Carbide Inserted Wire Pin Cutter 8.5" (22cm) Cut Up To 2.4mm (3/32") Spring Action Gold Color Handle German Grade Stainless Steel Surgical Veterinary Orthopedic Instruments

Front End cutter double-action tungsten carbide inserted orthopedic plier features a sharp cutting edge and scratch-resistant body. This sharp instrument has been assembled with German stainless to become corrosion resistant and robust for many years. T/C End Cutter is broadly used in spinal surgeries for wire placement in grafted areas. As well, it assists in cutting and manipulating the wires.

The self-opening spring action system and non-slip plier grip handles provides better safety an ergonomic grip to the

surgeon because of their perfectly balanced handle.

TC END CUTTER FEATURES

- Sharp Cutting Edges.
- Overall Length 8.5 Inches.
- Maximum Cutting Capacity 2.4mm (0.094").
- T/C (Tungsten Carbide Inserts) Cutting Jaws for Effective Cutting Action.
- Double Action
- Double Spring Action Plier Grip Handle.
- Reusable.
- Corrosion-resistant.
- German Grade Stainless Steel Construction.

PACKAGE INCLUDES:**WI-70001248****T/C Male Torquing Plier 5" (12.7cm) Gold Handle****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Male Torquing Plier 5" Tungsten Carbide Dental Torque Adjustment Arch Wire Torquing Pliers German Grade Stainless Steel Dentist Orthodontic Instruments

Male Torquing orthodontic pliers are most often used to place torque into an archwire for single tooth movement. These pliers come in handy when your patient has a row of teeth that needs an extra effort to straight or align. To achieve desired torque, torquing orthodontic pliers used in conjunction with the pliers for torquing or a male and female torquing pliers set. Our torquing plier male / inner has tungsten carbide inserts tips for superior grip and longevity. The pliers should be able to apply necessary torque without nicking or fraying the wire.

Our stainless steel Tungsten carbide torquing plier Male / Inner feature a streamlined design, sturdy construction and improved system that you'd expect in an orthodontic torquing tool made for professional orthodontists.

MALE TORQUING PLIER FEATURES

- Torquing Plier Male / Inner
- Tungsten Carbide Inserts For Better Grip & Longevity.
- 5" (12.7cm) Overall Length
- Perfect for Applying Torque Without Nicking or Fraying the Wire.
- L-Key Screw
- Durable, Strong, and Corrosion Resistant.
- Fully Autoclavable/Reusable



**WI-70003387****T/C Mathieu Ligature Elastic Placing Plier 5.5" Elastomeric Orthodontic****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality T/C Mathieu Ligature Elastic Placing Plier 5.5" (14cm) Tungsten Carbide Inserts Mathieu Needle Holder for Elastic Placing and Wire Holding Elastomeric Orthodontic Dental Instruments

T/C Mathieu Plier 5.5" (14cm) Elastomeric Tungsten Carbide Orthodontic Dental Needle Holder Braces Instruments. This versatile and affordable high quality orthodontic pliers instrument is a valuable addition to your orthodontic kit. All our products passes the National and International quality and safety standards to deliver best results.




TC MATHIEU LIGATURE ELASTIC PLACING PLIERS FEATURES**PREMIUM QUALITY TC MATHIEU PLIERS 5.5":**

These Ligature Elastic Placing Pliers and Needle Holder is a high value Orthodontic Instrument used for various applications such as elastic placing and wire holding.

STAINLESS STEEL CONSTRUCTION:

The Stainless Steel construction ensures maximum durability and dependability. You can re use these pliers over time as they will not bend or tarnish on you.

FULLY AUTOCLAVABLE: A hygienic choice for the orthodontic practice, these Mathieu Ligature Elastic Placing Pliers are fully autoclavable and corrosion resistant. They are easy to clean and use over and over for the highest possible value. The corrosion resistant properties allow them to stay in good condition for a long time.

**WI-70001290****T/C Pin & Wire Cutter Double Action 7" Angled Max Cutting Capacity 1.6mm (1/16")****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Pin and Wire Cutter 7" (17.8cm) Double Action CARB-N-SERT Cutting Blades Angled with Spring Action Gold Color Handle Maximum Cutting Capacity 1.6mm (1/16") Premium Quality Stainless Steel Surgical Veterinary Orthopedic Instruments

TC Pin and wire cutters (C-N-S pin and wire cutter) are used for cutting and removing of surgical pins & wires (Steinmann Pins and Kirschner Wires) during bone repair and fracture surgeries. It features tungsten carbide inserts for durability and precision, capable of cutting soft wire up to 2 mm and hard wire up to 1.6 mm. This Pin & Wire Cutter with double action compound hinges is a robustly structured to enhance cutting power with reduced effort, along with a spring

mechanism for improved grip and shock absorption. Gold plated non-slip plier grip handles provides ergonomic grip to the surgeon during long veterinary surgical procedures. Made from high-quality stainless steel, it is sterilizable for reuse.

TC PIN AND WIRE CUTTER FEATURES

- 7 Inch (17.8cm) Overall Length.
- Sharp Angled Blades.
- Maximum Cutting Capacity Up to 2 Mm (5/64") Soft Wire or 1.6 Mm Hard Wire 1.6mm (1/16").
- T/C (Tungsten Carbide Inserts) Cutting Jaws for Effective Cutting Action.
- Double Action Compound Hinges to Enhance Cutting Power

**WI-70001289****T/C Pin & Wire Cutter Double Action 9" Angled Max Cutting Capacity 2.4mm (3/32")****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Pin and Wire Cutter 9" (22.9cm) Double Action CARB-N-SERT Cutting Blades Angled with Spring Action Gold Color Handle Maximum Cutting Capacity 2.4mm (3/32") Premium Quality Stainless Steel Surgical Veterinary Orthopedic Instruments

mechanism for improved grip and shock absorption. Gold plated non-slip plier grip handles provides ergonomic grip to the surgeon during long veterinary surgical procedures. Made from high-quality stainless steel, it is sterilizable for reuse.

TC PIN AND WIRE CUTTER FEATURES

- 9 Inch (22.9cm) Overall Length.
- Sharp Angled Blades.
- Maximum Cutting Capacity Up to 2.4mm (3/32").
- T/C (Tungsten Carbide Inserts) Cutting Jaws for Effective Cutting Action.
- Double Action Compound Hinges to Enhance Cutting Power

**TC PIN AND WIRE CUTTER****WI-70001287****T/C Side Cutter Double Action 9" Max Cutting Capacity 2.4mm 3/32"****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality T/C Side Cutter Double Action 9" (22.8cm) Tungsten Carbide Inserted Pin & Wire Cutter, Maximum Capacity 2.4mm (3/32") Spring Action Gold Color Handle German Grade Stainless Steel Surgical Veterinary Orthopedic Instruments

The self-opening spring action system and non-slip plier grip handles provides better safety an ergonomic grip to the surgeon because of their perfectly balanced handle.

TC SIDE PIN & WIRE CUTTER FEATURES

- Straight Sharp Cutting Edges
- Overall Length 9 Inches.
- Maximum Cutting Capacity 2.4mm (0.094").
- T/c (Tungsten Carbide Inserts) Cutting Jaws for Effective Cutting Action.
- Double Action
- Double Spring Action Plier Grip Handle.
- Reusable.
- Corrosion-resistant.

**T/C SIDE CUTTER DOUBLE ACTION****T/C SIDE CUTTER DOUBLE ACTION**

Side Cutter utilized in fracture repair, bone stabilization, and implant orthopedic surgeries. The strong and sharp jaws used to remove wires and pins from bones. Side cutter double-action tungsten carbide inserted orthopedic plier features a sharp straight cutting edge and scratch-resistant body. This sharp instrument has been assembled with German stainless to become corrosion resistant and robust for many years.

**WI-70001067****T/C Torquing Plier Female/ Outer 5" (12.7cm)****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Torquing Plier Female/Outer 5" (12.7cm) Tungsten Carbide Dental Torque Adjustment Arch Wire Torquing Pliers German Grade Stainless Steel Dentist Orthodontic Instruments

The outer torquing pliers (female torquing Plier) feature an outer tip design and are used in conjunction with male/inner torquing pliers to deliver appropriate torquing to an arch wire to facilitate single tooth movement. Easily and accurately apply torque to arch wires measuring up to .028". Our torquing plier female / outer has tungsten carbide inserts tips for superior grip and longevity. The pliers should be able to apply necessary torque without nicking or fraying the wire.

Our stainless steel Tungsten carbide torquing plier feature a streamlined design, sturdy construction and improved system that you'd expect in an orthodontic torquing tool made for professional orthodontists.

T/C TORQUING PLIER FEMALE FEATURES

- Torquing Plier Female / Outer
- Tungsten Carbide Inserts For Better Grip & Longevity.
- 5" (12.7cm) Overall Length
- Perfect for Applying Torque Without Nicking or Fraying the Wire.
- Durable, Strong, and Corrosion Resistant.
- Fully Autoclavable/Reusable
- Made with German Grade Stainless Steel.

**WI-70003140****T/C Universal Cut & Hold Distal End Cutter 4.5"****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Orthodontic Universal Cut & Hold Distal End Cutter 4.5" (11.4cm) with Tungsten Carbide Inserts Designed to Cut Soft Wire Pins and Ligature Wire from 0.012"/.30mm. Slim profile is a plus to access hard-to-reach areas.

TC UNIVERSAL DISTAL END CUTTER PLIER FEATURES

TC UNIVERSAL CUT & HOLD DISTAL END CUTTER: The Universal Cut And Hold Distal End Cutters are a sturdy handheld orthodontic tool designed to cut soft wire pins and ligature wire from 0.012"/.30mm to 0.021" x 0.025" (0.53mm x 0.64mm)

DURABLE STAINLESS STEEL CONSTRUCTION: The Premium Quality Stainless Steel and Tungsten Carbide composition of the Universal Cut &

Hold Distal End Cutter offers a long life of strong performance. It is built to last and reliable for long term use and will not bend or break on you. Made with anti-corrosion properties + tarnish free.

FULLY REUSABLE + AUTOCLAVABLE:

A hygienic and valuable choice, the TC Distal End Cutters can be fully autoclaved without fear, for long term use without degradation. They are ultrasonic cleaned and have a beautiful high mirror finish.

DESIGNED FOR PRECISE + EASY OPERATION:

The TC Distal End Cutters are designed for easy accessibility, with a comfortable grip and less chance of getting loose at the joints even after long usage. The precision joint with adjustable screw makes it possible to disassemble the plier/cutter and adjust it for wear and tear. These



**WI-70001927****T/C Wire & Pin Cutter 18.5\"/>****\$149.99** per pieceTrade price from **\$137.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Wire & Pin Cutter 18.5\"/>

Heavy Duty Plate Tungsten Carbide Inserts Pin and Wire Cutter with Adjustable Bolt Gold Handle assist surgeons to dissect and remove bone plates, wires and pins during orthopedic surgeries. This instrument has sterilizable material due to high quality stainless steel crafting. It is quite resilient to last for many years and reused after sterilization.

- Premium Grade Stainless Steel For A Reliable Service.
- Overall Length: 18 1/2"
- Max. cutting capacity: 3/16" (4.8 mm)
- Side cut
- Made with High Quality Stainless Steel (German Grade)

PACKAGE INCLUDE:

- T/C Wire & Pin Cutter 18.5\"/>

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Orthopedic Instruments**APPLICATION:** Bone Pin Cutting**CURVATURE:** Straight**CUTTING END:** Upto 4.8mm Diameter**PIN CUTTER FEATURES**

- Large Handles For forcible & controlled grips.
- Tungsten Carbide Inserts Sharp Jaws For effective cutting action.

**WI-70003329****TC Adams Pliers 5\"/>****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

TC Adams Pliers 5\"/>

Adam Pliers tungsten carbide inserts flat pyramid beaks are precision tools designed for right angle bends and appliance adaptation on stainless steel wires up to .040\"/>

wires up to .040\"/>

ADAMS PLIER FEATURES

- 5\"/>
- Tungsten Carbide Inserts Tips.
- 2 Smooth Pyramidal Shaped Beaks.
- Max. Wire Thickness Up to .040\"/>
- Designed for Right Angle Bends and Appliance Adaptation.
- Enables Construction of Adam Clasps, Precise Wire Adjustment, and Efficient Tooth Alignment.



**WI-70003328****TC Bracket Remover Plier 5.5" Straight****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

TC Bracket Remover Plier 5.5" (14cm) Straight Beaks with Gold Handle Direct Bond Braces Removing Premium Quality Stainless Steel Dental Orthodontic Instruments.

The Direct Bond Bracket Removal Plier with tungsten carbide inserts is a handheld tool designed for quick and efficient bracket removal. It is suitable for steel, ceramic, and plastic brackets. The tool features a comfortable grip, smooth maneuvering, and precise operation, making it a valuable addition to any toolkit. The TC tip ensures precision and long-term accuracy in every procedure. The premium quality stainless steel bracket removing pliers are designed for long-lasting performance, ensuring they won't bend

or break, and are also anti-corrosion properties.

TC BRACKET REMOVING PLIER STRAIGHT FEATURES

- 5.5" (14cm) Total Length.
- Straight Tungsten Carbide Inserted Tips.
- Suitable for Stainless Steel, Ceramic, and Plastic Brackets for All Types.
- Gold Plated Smooth Grip Handles.
- Corrosion Resistance.
- Completely Sterilizable.
- Fully Reusable / Autoclavable.
- Crafted With Premium Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc TC Bracket Remover Plier 5.5" (14cm) Straight

WI-70001030**TC Flush Cut Distal End Cutter 4.5" Orthodontic Plier****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Tip Distal End Cutter Flush Cut 4.5" (11.4cm) with Tungsten Carbide Inserts Universal Cutter for All Types of Wires Up to Maximum of 0.021" x 0.025 (.53mm x .64mm) Premium Grade Stainless Steel Dental Orthodontic Instrument.

Flush cut distal end cutter is designed to cut wires as small as 0.30 mm, and as large as 0.53x0.64. The extra thin beaks can cut even the smallest amount of jutting wire. Ensuring no trauma to the patients tissues. The tungsten carbide inserts features gives the beaks maximum sharpness and ensures a long lifetime of operation. The ergonomically designed handles of our instrument, along with the unique L-key screw design, ensure the easiest cutting with the smallest amount of pressure. Manufactured with premium

quality stainless steel to ensure durability. Flush cut distal end cutter is sterilizable and can be reused.

TC FLUSH CUT DISTAL END CUTTER FEATURES

- Universal Cutter with Ultra Thin Beaks Leaving No Trace Of The Wire.
- Maximum Cutting from 0.30mm to 0.53x0.64mm.
- Tungsten Carbide Inserts Cutting Edge For Maximum Efficiency & Durability.
- L-Key Screw Design.
- Ergonomic Grip To Decrease Effort.
- Resuable.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDES



WI-70003326

TC Hard Wire Cutter 15 Degree Curved Orthodontic Pliers

\$12.99 per piece

Trade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

TC Hard Wire Cutter 5" (12.7cm) with 15 Degree Curved Tungsten Carbide Inserts Tip #151 Orthodontic Plier – German Grade Stainless Steel Dental Orthodontic Instruments.

Orthodontic hard wire cutters are designed to precisely cut and trim high-strength archwires and ligatures (up to 0.022" x 0.028"). They effectively trim heavy stainless steel or nickel-titanium wires, clip fine ligature ties or soft wire pins, and assist in creating specialized loops or springs. Featuring tungsten carbide tips for enduring sharpness, these cutters may include angled options for accessing tight spaces and are made from surgical-grade stainless steel to ensure sterilizability through repeated autoclaving.

TC HARD WIRE CUTTER FEATURES:

- 5 Inch (12.7cm) Overall Length.
- Hard Wires Cutting Up to .022" X .028"
- 15-degree Curved Offset Tip Design Improves Access in Difficult Areas.
- Tungsten Carbide Inserts for Enhanced Cutting and Lasting Sharpness.
- Gold Color Plated Handle With Smooth Grip.
- Capable of Easily Trimming Hard Orthodontic Wires and Cutting Ligatures.
- Fully Reusable, Easy to Clean, and Autoclavable.
- Made From German Grade Stainless Steel.



WI-70003325

TC How Crown Plier Curved Dental Orthodontics Surgical

\$30.53 per piece

Trade price from **\$28.09**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

TC How Crown Plier Curved Dental Orthodontics Surgical is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$30.53, falling to \$28.09 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003324

TC How Crown Plier Straight Dental Orthodontics Surgical

\$27.46 per piece

Trade price from **\$25.26**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

TC How Crown Plier Straight Dental Orthodontics Surgical is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$27.46, falling to \$25.26 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003313

TC Pin and Ligature Hardwire Cutter Straight Tip

\$12.99 per piece

Trade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

TC Hard Wire Cutter 5" (12.7cm) with 15 Degree Curved Tungsten Carbide Inserts Tip #151 Orthodontic Plier – German Grade Stainless Steel Dental Orthodontic Instruments.

Orthodontic hard wire pin ligature cutter with TC (Tungsten Carbide Tips) is a dental orthodontic plier for precisely cutting pin wires, metal ligatures, and hard archwires. Featuring a tungsten carbide inserts straight tip design facilitates access to brackets and oral tissues while blunt tip prevents injury to oral tissues. It can trim soft and hard wires up to .022" x .028". Premium quality stainless steel construction ensures strength and resistance to multiple sterilization cycles. The ergonomic, non-slip gold

color handles enhance comfort and reduce user fatigue.

TC PIN AND LIGATURE HARD WIRE CUTTER FEATURES:

- 5 Inch (12.7cm) Overall Length.
- Trims Soft and Hard Wires (Typically Up to .022" X .028").
- Straight Tip Design Allows Precise, Close Access to Brackets and Patient Oral Tissues.
- Tungsten Carbide Inserts for Enhanced Cutting and Lasting Sharpness.
- Gold Color Plated Handle With Smooth Grip.
- Designed to Precisely Cut Pin Wires, Metal Ligatures, and Hard Archwires
- Fully Reusable, Easy to Clean, and Autoclavable.



**WI-70003323****TC Pin Puller 5.5" Orthodontic Veterinary Plier****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

T/C Pin Puller Pliers 5.5" (14cm) with Tungsten Carbide Inserts Serrated Jaws & Gold Plated Handles Premium Quality Stainless Steel Orthodontic Dental Surgical Instruments.

TC pin pulling pliers are designed for cutting and holding wires in orthodontic appliances like braces and wires. This plier feature TC (Tungsten Carbide) inserts serrated jaws for sharpness and durability. They provide high control and safety for practitioners and patients during orthodontic procedures. Pin puller can also be used for wire manipulation, bracket removal, and adhesive removal. They are also suitable for various applications, including wire manipulation, bracket removal, and adhesive removal. TC pin pulling pliers are expertly designed and

crafted from premium quality stainless steel, ensuring durability, strength, and resistance to corrosion.

TC PIN PULLER ORTHODONTIC PLIER FEATURES

- Crafted with 5.5" (14cm) Overall Length.
- 3.0mm Tungsten Carbide Inserts Serrated Jaws.
- Designed for Cutting and Holding Wires in Orthodontic Appliances.
- Provides High Control and Safety During Orthodontic Procedures.
- Can Be Used for Wire Manipulation, Bracket Removal, and Adhesive Removal.
- Highly Durable and Corrosion Resistance.
- Reusable.

**WI-70003322****TC Weingart Pliers 5.5" Utility Dental Orthodontic Instrument****\$12.49** per pieceTrade price from **\$11.49**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Weingart Pliers with Tungsten Carbide Inserts L-Key Dental Wire Bending Orthodontics Braces Placement Stainless Steel Ortho Dental Instrument

Versatile utility pliers with serrated tips hold wires firmly at a convenient working angle. The slim Weingart plier provides easy access to lingual and posterior areas. Tapered beaks fit easily between brackets and hard-to-reach areas.

ADVANTAGES OF WEINGART PLIERS WITH A HEXAGONAL SCREW/L-KEY

- Precise function.
- Angular Design is convenient for the placement and removal of arch wires.
- Longevity of life

- Less chance of getting loose at the joints even after long usage.
- When found loose, can be tightened with the aid of an Allen-key easily.
- The precision joint with a screw makes it adjustable.

TC WEINGART PLIERS FEATURES

PERFECT WEINGART PLIER: TC Weingart utility orthodontic plier is a versatile serrated tool. The tapered angular design is convenient for placement and removal of arch wires. **WEINGART L-KEY PLIERS WITH HEXAGONAL SCREW:** Used to shorten arch wire. Provides an accurate 1mm stop.

EASY TO USE: This Weingart Plier is user-friendly and easy to clean and sterilize.

BEAUTIFULLY FORGED: High mirror

**WI-70003321**

Three Prong Plier Orthodontic Dental Instrument

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Three Prong Plier – Premium Quality Stainless Steel Three Jaw Design Orthodontists Plier for Orthodontics and Prosthodontics Procedures.

Three Prong plier is a specialized instrument for bending, contouring, and adjusting orthodontic wires and clasps. It is designed for precise clasp adjustment, wire bending, and appliance fitting, suitable for both orthodontic and prosthodontic applications. The pliers provide controlled force, accuracy, and reliability in clinical and laboratory settings. They are specifically used for creating sharp angles and precise curves in archwires, adjusting metal clasps on appliances, and forming loops on springs. Made from premium

stainless steel, these pliers are durable and fully sterilizable for repeated use.

THREE PRONG PLIERS FEATURES

- Three-prong Configuration: One Middle Prong and Two Opposing Outer Prongs.
- Suitable for Heavy or Standard Orthodontic Wires Sizes Up to 0.030" (0.76 mm).
- 3 Jaws for Balanced Pressure, Reducing Risks of Wire Damage or Kinking.
- Box-joint Configuration Facilitates Stable Alignment.
- Ergonomic Handles Designed to Minimize Hand Fatigue.
- Fully Autoclavable and Heat-sterilizable.

**WI-70003489**

Toe Nail Cutter 5.5" (14cm) Moon Shaped with Back Lock

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Moon Shape Toe Nail Cutter 5.5" (14cm) with Back Handle Lock Heavy Duty Thick Nails Clipper Nipper Tools Premium Quality Stainless Steel Chiropody Podiatry Instruments.

MOON SHAPE TOE NAIL CUTTER FEATURES

HIGH QUALITY STAINLESS STEEL:

Made with superior quality durable stainless steel. The Premium Quality Stainless Steel composition of these tools offers a long life of strong performance, as they are built to last

FUNCTIONAL & ERGONOMIC DESIGN:

Features an ergonomic design to complete the task and make your life easy. Provides a comfortable experience and offers more control for precise and easy operation. Specially designed zero clearance hinge makes

handles better connected for safe operation

CORROSION RESISTANCE: Highly Polished Finish for Aesthetics, Corrosion Resistance and Anti-Oxidation

CUSTOMER SATISFACTION: We take pride in the quality of our products, each item is inspected under strike quality control system to deliver the finest tools. 100% satisfaction guaranteed.

PACKAGE INCLUDE:


- 1 Piece Nail Nipper Cutter Moon Shaped 5.5" with Handle Lock

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Podiatry Instruments





WI-70000433

Torquing Plier Female/ Outer 5" (12.7cm)

\$10.99 per piece

Trade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Torquing Plier Female/Outer 5" (12.7cm) Dental Torque Adjustment Arch Wire Torquing Pliers German Grade Stainless Steel Dentist Orthodontic Instruments

Female Torquing Plier used by orthodontists to apply precise torque to orthodontic wires, ensuring correct tooth alignment and positioning. Female torquing orthodontic plier is designed to create controlled torque on the wire, enabling adjustments for specific teeth. This torquing plier is commonly utilized for rectangular or square wires in orthodontic treatments. Outer torquing plier feature an outer tip design and are used in conjunction with male/inner torquing pliers to deliver appropriate torquing to an arch wire to facilitate single tooth movement. This high-

grade stainless steel orthodontic plier is designed for durability and sterilization, allowing for precise adjustments without damaging the wire or braces.

FEMALE TORQUING PLIER FEATURES

- 5" (12.7cm) Overall Length
- Torquing Plier Female / Outer.
- Perfect for Applying Torque Without Nicking or Fraying the Wire.
- Durable, Strong, and Corrosion Resistant.
- Fully Autoclavable/Reusable.
- Made with German Grade Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Torquing Plier Female/Outer 5" (12.7cm)



WI-14719250

Ultimate 6-Piece Stainless Steel Pliers Set 8"

\$12.00 per piece

Trade price from **\$8.00**
Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Ultimate 6-Piece Stainless Steel Pliers Set 8" is a pliers & cutter we forge, grind and finish in our own Sialkot workshop. Instruments of this type grip, bend or cut wire, pins and hard material. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$12.00, falling to \$8.00 at 2 sets — about 33% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003327

Universal Cut and Hold Distal End Cutter 4.5\"/>

\$10.99 per piece

Trade price from \$10.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Universal Cut & Hold Distal End Cutter 4.5\"/>

Universal Cut & Hold Distal End Cutter is an orthodontic instrument designed to cut archwire flush with the buccal tube while securely holding the end to prevent ingestion or aspiration of wire segments. It accommodates various wire sizes (up to .021\"/>

Cutter, designed for orthodontic procedures, ensures safety and hygiene. This tool effectively holds the cut-off part, preventing injury to the cheek during use. Nade of high-grade stainless steel for durability and precision in wire-adjustment procedures.

BRACES WIRE ORTHODONTIC DISTAL END CUTTER PLIER FEATURES

- Size: 4.5 Inch (11.4cm) Long.
- Cut wires from 0.012\"/>
- Slim, Tempared Cutting Heads Improve Access.
- Align the Tip With Wire to Avoid Breaking the Tips.



WI-70003320

Universal Pliers Orthodontic Surgical Dental

\$22.78 per piece

Trade price from \$20.96
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Universal Pliers Orthodontic Surgical Dental is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$22.78, falling to \$20.96 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003319****V-Stop Pliers 5.5" with V Notch Tip****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

V-Stop Pliers 5.5" (14cm) Making Precise V Notch on Archwires Up to .022" x .025" Premium Quality Dental Orthodontic Instrument.

V stop pliers are used to place a V bend on archwires up to .022" x .025" (.56 mm x .64 mm) to prevent archwire from traveling, including nickel titanium. V-stop pliers have a "V"-shaped groove at the beak tip, used for wire bending in orthodontic treatment. Orthodontic V-stop plier enables orthodontic professionals to customize archwires to the patient's treatment plan, creating bends and stops for tooth rotation, tipping, and other movements. Orthodontists use the V stop notch wire bending plier to create precise angles in orthodontic wires, crucial for controlling tooth movement and

achieving desired alignment. These pliers made from premium stainless steel, create precise V-bends in braces' arch wires, enhancing daily performance and reshaping orthodontic care for professionals.

V-STOP WIRE FORMING PLIERS FEATURES

- 5.5" (14cm) Overall Length.
- Creates 1mm stop.
- For wires up to .025".
- V Shaped Tip on the Inner Edge.
- Perfect to Create V Bends on Arch Wires.
- Autoclavable / Reusable
- Made with Premium Great Stainless Steel

PACKAGE INCLUDES

- 1 Pc V-Stop Pliers 5.5" (14cm)

**WI-70002447****Weingart Plier 5.5" Wire Bending Orthodontic Instrument****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Weingart Pliers L-Key 5.5" (14cm) Dental Wire Bending Orthodontics Braces Placement Oral Tooth Stainless Steel Ortho Dental Instrument

Versatile utility pliers with serrated tips hold wires firmly at a convenient working angle. Weingart plier provides easy access to lingual and posterior areas. Tapered beaks fit easily between brackets and hard-to-reach areas.

ADVANTAGES OF WEINGART PLIERS WITH A HEXAGONAL SCREW/L-KEY

- Precise function.
- Angular Design is convenient for the placement and removal of arch wires.
- Longevity of life

- Less chance of getting loose at the joints even after long usage.
- When found loose, can be tightened with the aid of an Allen-key easily.
- The precision joint with a screw makes it adjustable.

WEINGART PLIERS FEATURES

PERFECT WEINGART PLIER: Weingart Utility Orthodontic Plier is a versatile serrated tool. The tapered angular design is convenient for placement and removal of arch wires.

WEINGART L-KEY PLIERS WITH HEXAGONAL SCREW: Used to shorten arch wire. Provides an accurate 1mm stop.

EASY TO USE: This Weingart Plier is user-friendly and easy to clean and sterilize.

BEAUTIFULLY FORGED: High mirror





WI-44552379

Wire Mesh Surgical Tools Sterilization Tray Medical

\$65.00 per piece

Trade price from **\$56.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Wire Mesh Surgical Tools Sterilization Tray Medical — one of 4,770 instruments we manufacture ourselves rather than buy in. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship.

Listed at \$65.00, falling to \$56.00 at 10 pieces — about 14% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003317

Young Plier 13.5cm Loop Forming Bending Orthodontic Instrument

\$10.99 per piece

Trade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Young Loop Forming Plier 5.3" (13.5cm) with Three Sectional Beak and Three Grooves for Wires Up to 0.028" Premium Quality Stainless Steel Universal Orthodontic Plier.


Young orthodontic plier, also known as Young's loop-forming plier, is a specialized dental tool designed for bending and shaping orthodontic wires into precise loops (labial, lingual, or canine). Its key uses include loop formation with multi-stepped cylindrical beaks sized at 1.3mm, 3mm, and 4.2mm for accurate loops, wire bending for sharp, controlled bends, and customization of removable or fixed orthodontic appliances. Made of premium quality stainless steel, this orthodontic plier offers long-lasting

performance and being fully sterilizable and reusable.

YOUNG LOOP FORMING PLIER FEATURES

- Overall Length 5.3" (13.5cm)
- For Wires up to 0.028 inches (0.7 mm).
- Three Graduated Cylindrical Sections (1.3mm, 3.0mm, 4.2mm) for Distinct Loop Sizes For Versatility.
- Opposing Contoured or Flat Grooved Jaws Stabilize Hard Wires and Prevent Slipping.
- Ergonomically Designed Handles To Decrease Fatigue.
- Fully Autoclaved
- Made with Anti-corrosion German Grade Stainless Steel

PACKAGE INCLUDE





WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

Request a quotation

Send us product codes from this catalogue, quantities and your country.
We reply the same working day on WhatsApp.

Online shop	shop.waqarintl.com
Email	sales@waqarintl.com
WhatsApp /	
Phone	+92 307 6173627
Works	Roras Road, Opposite Anwar Khawaja Industries, Sialkot 51310, Pakistan

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.