



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

Scalpels & Knives Surgical Instruments

Reusable stainless-steel scalpels & knives forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



10
PCS



10
PCS



207

PRODUCTS

\$0.49

PRICES FROM

OEM

PRIVATE LABEL

Global

WORLDWIDE
EXPORT

Contents

207 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

Product	Code	Price	Page
10 Fahnestock Knife 13cm Small – Double Ended	WI-70001128	\$34.99	7
10 Round Scalpel Handles #3 Fits Blade 10 to 16	WI-70001005	\$39.90	7
10 Scalpel Blades #60 with FREE Scalpel Knife Handle #8	WI-70002151	\$14.99	8
10 Sterile Surgical Blades #10 with Scalpel Knife Handle #3	WI-70002882	\$9.99	8
10 Sterile Surgical Blades #11 with Scalpel Knife Handle #3	WI-70002881	\$9.99	9
10 Sterile Surgical Blades #12 with Scalpel Knife Handle #3	WI-70002880	\$9.99	9
10 Sterile Surgical Blades #15 with Scalpel Knife Handle #3	WI-70002879	\$9.99	10
10 Sterile Surgical Blades #15C with Scalpel Knife Handle #3	WI-70002878	\$9.99	10
10 Sterile Surgical Blades #20 with Scalpel Knife Handle #4	WI-70002887	\$9.99	11
10 Sterile Surgical Blades #21 with Scalpel Knife Handle #4	WI-70002888	\$9.99	11
10 Sterile Surgical Blades #22 with Scalpel Knife Handle #4	WI-70002886	\$9.99	12
10 Sterile Surgical Blades #23 with Scalpel Knife Handle #4	WI-70002885	\$9.99	12
10 Sterile Surgical Blades #24 with Scalpel Knife Handle #4	WI-70002884	\$9.99	13
10 Sterile Surgical Blades #25 with Scalpel Knife Handle #4	WI-70002883	\$9.99	13
100 Scalpel Blades #60 with FREE Scalpel Knife Handle #8	WI-70002187	\$76.99	14
100 Scalpel Blades #70 with Scalpel Knife Handle #8	WI-70001341	\$76.99	14
12 Dental Lab Plaster Knife #12R New	WI-70003825	\$71.68	15
12 Dental Lab Plaster Knife #7R New	WI-70003824	\$69.30	15
2 Molt Mouth Gag 3.5" + 4.5" Pediatric & Child Curved Blades	WI-70003981	\$16.99	16
2 Round Scalpel Handles #3 Fits Blade 10 to 16	WI-70001007	\$13.99	16
20 Pcs Advanced Stainless Steel Medical	WI-01210325	\$4.00	17
20 Scalpel Blades #60 with FREE Scalpel Knife Handle #8	WI-70002150	\$22.99	17
20cm Autopsy Knife Fixed Blade Stainless Steel Post Mortem	WI-36732508	\$18.00	18
20cm Fixed Blade Autopsy & Medical Dissection Knife	WI-36732510	\$18.00	18
20cm Fixed Blade Autopsy Knife	WI-36707577	\$18.00	19
3 Dental Lab Plaster Knife #12R New	WI-70003820	\$22.79	19
3 Dental Lab Plaster Knife #7R New	WI-70003819	\$20.21	20
3 Finger Ring Cutter Replacement Blade – Small	WI-70001317	\$10.99	20
3 McCoy Fiber Optic Flex Tip Laryngoscope Blades #2, #3 & #4	WI-70001329	\$49.99	21
4 Macintosh Conventional Blades #1, #2, #3, #4	WI-70000203	\$35.99	21
4 Macintosh Fiber Optic Blades #1, #2, #3, #4	WI-70001345	\$44.99	22
5 Fiber Optic Miller Blades #0, #1, #2, #3, #4	WI-70001358	\$54.99	22
5 Finger Ring Cutter Replacement Blade – Small	WI-70001316	\$12.99	23
5 Miller Laryngoscope Blades #0, #1, #2, #3, #4 – Conventional	WI-70001359	\$39.99	23
5 Round Scalpel Handles #3 Fits Blade 10 to 16	WI-70001006	\$29.99	24
50 Scalpel Blades #60 with FREE Scalpel Handle #8	WI-70001959	\$43.99	24
6 Dental Lab Plaster Knife #12R New	WI-70003815	\$41.24	25
6 Dental Lab Plaster Knife #7R New	WI-70003814	\$37.64	25
Adjustable Scalpel Handle #3 – 180° Rotation with 7 Adjustable Angles	WI-70001004	\$10.99	26
Bard Parker Scalpel Handle #3 with Ruler	WI-70001016	\$9.99	26
Bard Parker Scalpel Handle #4 with Ruler	WI-70001015	\$9.99	27

Contents (continued)

Product	Code	Price	Page
Black Scalpel Blade Handle #3	WI-70000003	\$9.99	28
Buck Ear Curette 6.5" Straight Sharp Ring Blade Size #0	WI-70000470	\$8.99	28
Buck Ear Curette 6.5" Straight Sharp Ring Blade Size #00	WI-70000471	\$8.99	29
Buck Ear Curette 6.5" Straight Sharp Ring Blade Size #1	WI-70000469	\$8.99	29
Buck Ear Curette 6.5" Straight Sharp Ring Blade Size #2	WI-70000468	\$8.99	30
Buck Ear Curette 6.5" Straight Sharp Ring Blade Size #3	WI-70000467	\$8.99	30
Con Scodellino Knife – Large	WI-70003854	\$9.99	31
Con Scodellino Knife – Small	WI-70003855	\$8.99	31
Dental Lab Plaster Knife #12R New	WI-70003811	\$8.99	32
Dental Lab Plaster Knife #7R New	WI-70003810	\$8.99	32
Design Disposable Surgical Scalpels Steel	WI-28678810	\$1.50	33
Disposable Steel Surgical Scalpels and Blades Instrument	WI-28667834	\$1.50	33
Disposable Sterile Surgical Scalpel Stainless Steel	WI-43380071	\$0.88	34
Dissecting Scalpel Knife 6" – Fixed Blade	WI-70000253	\$8.99	34
Dog Grooming Coat Stripping Knife Set 6	WI-08790844	\$12.25	35
Electrosurgical Pencil 2mm Wide Blade – Single Use	WI-70000035	\$14.99	35
Ergonomic Medical Instrument Set Stainless Steel No. 3/4/7	WI-34506927	\$1.50	36
Fahnestock Wax Knife 5.25" – Small	WI-70003850	\$8.99	36
Fahnestock Wax Knife 7" – Large	WI-70003849	\$8.99	37
Fiber Optic Miller Blade #0 with Laryngoscope Medium Handle	WI-70001355	\$28.99	37
Fiber Optic Miller Blade #1 with Laryngoscope Medium Handle	WI-70001354	\$28.99	38
Fiber Optic Miller Blade #2 with Laryngoscope Medium Handle	WI-70001353	\$29.99	38
Fiber Optic Miller Blade #3 with Laryngoscope Medium Handle	WI-70001352	\$29.99	39
Fiber Optic Miller Blade #4 with Laryngoscope Medium Handle	WI-70001351	\$29.99	39
Finger Ring Cutter 6.5" with Saw Blade Hollow Handle	WI-70001315	\$9.99	40
Finger Ring Cutter Replacement Blade – Large	WI-70000013	\$8.99	40
Finger Ring Cutter Replacement Blade – Small	WI-70001318	\$8.99	41
Foot Callus Shaver Hard Skin Corn Remover with 10 Replacement Blades – Black	WI-70001676	\$8.99	41
Foot Callus Shaver Hard Skin Corn Remover with 10 Replacement Blades – Blue	WI-70001675	\$8.99	42
Foot Callus Shaver Hard Skin Corn Remover with 10 Replacement Blades – White	WI-70001674	\$8.99	42
Forged Stainless Steel Scalpel Surgical Cutting Instrument for Controlled	WI-46940065	\$1.80	43
General Surgery Scalpel Handle Reusable Medical Blade	WI-45279261	\$2.00	43
German Stainless Steel Scalpel Handle 145mm 5.34" Straight	WI-34512907	\$1.50	44
Heavy Duty Finger Ring Cutter Plier 6.5" with 30mm OD Saw Blade	WI-70001068	\$11.99	44
High Carbon Steel No.3 Surgery Scalpel BP Handle	WI-34495866	\$1.50	45
Knife Sharpening Steel 8" Black Handle	WI-70001141	\$14.99	45
Lesmann Knife 13cm – Small	WI-70003842	\$8.99	46
Lesmann Knife 17cm – Large	WI-70003841	\$9.99	46
Lucas Curette #85 Double Ended 2.7mm Spoon Shaped Blades	WI-70001080	\$9.99	47
Mac Blade #0 with Conventional Laryngoscope Medium Handle	WI-70001975	\$20.99	47
Mac Blade #1 with Conventional Laryngoscope Medium Handle	WI-70001974	\$21.99	48
Mac Blade #2 with Conventional Laryngoscope Medium Handle	WI-70001973	\$23.99	48

Contents (continued)

Product	Code	Price	Page
Mac Blade #4 with Conventional Laryngoscope Medium Handle	WI-70001972	\$25.99	49
Mac Blade #5 with Conventional Laryngoscope Medium Handle	WI-70000031	\$29.99	49
Macintosh Blade #0 with Fiber Optic Laryngoscope Medium Handle	WI-70001343	\$28.99	50
Macintosh Blade #1 with Fiber Optic Laryngoscope Medium Handle	WI-70001336	\$28.99	50
Macintosh Blade #2 with Fiber Optic Laryngoscope Medium Handle	WI-70001335	\$29.99	51
Macintosh Blade #3 with Fiber Optic Laryngoscope Medium Handle	WI-70001334	\$29.99	51
Macintosh Blade #4 with Fiber Optic Laryngoscope Medium Handle	WI-70001333	\$29.99	52
Macintosh Conventional Laryngoscope Blade #0	WI-70001980	\$11.99	52
Macintosh Conventional Laryngoscope Blade #1	WI-70001979	\$11.99	53
Macintosh Conventional Laryngoscope Blade #2	WI-70001978	\$12.99	53
Macintosh Conventional Laryngoscope Blade #4	WI-70001977	\$12.99	54
Macintosh Conventional Laryngoscope Blade #5	WI-70000033	\$14.99	54
Macintosh Fiber Optic Laryngoscope Blade #0 Anesthesia Intubation	WI-70001350	\$15.99	55
Macintosh Fiber Optic Laryngoscope Blade #1 Anesthesia Intubation	WI-70001349	\$15.99	55
Macintosh Fiber Optic Laryngoscope Blade #2 Anesthesia Intubation	WI-70001348	\$16.99	56
Macintosh Fiber Optic Laryngoscope Blade #3 Anesthesia Intubation	WI-70001347	\$16.99	56
Macintosh Fiber Optic Laryngoscope Blade #4 Anesthesia Intubation	WI-70001346	\$16.99	57
McCoy Fiber Optic Flex Tip Laryngoscope Blade #2	WI-70001332	\$19.99	57
McCoy Fiber Optic Flex Tip Laryngoscope Blade #3	WI-70001331	\$19.99	58
McCoy Fiber Optic Flex Tip Laryngoscope Blade #4	WI-70001330	\$19.99	58
McCoy Laryngoscope Blade #2 with Fiber Optic Medium Handle	WI-70001328	\$26.99	59
McCoy Laryngoscope Blade #3 with Fiber Optic Medium Handle	WI-70001327	\$26.99	59
McCoy Laryngoscope Blade #4 with Fiber Optic Medium Handle	WI-70001326	\$26.99	60
Medical Sterile Stainless Steel Surgical Scalpel Blades	WI-09286277	\$2.00	60
Micro Finger Ring Cutter 4.75" with Saw Blade	WI-70001322	\$8.99	61
Miller Blade #0 with Conventional Laryngoscope Medium Handle	WI-70001364	\$20.99	61
Miller Blade #1 with Conventional Laryngoscope Medium Handle	WI-70001363	\$20.99	62
Miller Blade #2 with Conventional Laryngoscope Medium Handle	WI-70001362	\$21.99	62
Miller Blade #3 with Conventional Laryngoscope Medium Handle	WI-70001361	\$21.99	63
Miller Blade #4 with Conventional Laryngoscope Medium Handle	WI-70001360	\$21.99	63
Miller Fiber Optic Laryngoscope Blade #0	WI-70001369	\$15.99	64
Miller Fiber Optic Laryngoscope Blade #1	WI-70001368	\$15.99	64
Miller Fiber Optic Laryngoscope Blade #2	WI-70001367	\$16.99	65
Miller Fiber Optic Laryngoscope Blade #3	WI-70001366	\$16.99	65
Miller Fiber Optic Laryngoscope Blade #4	WI-70001365	\$16.99	66
Miller Laryngoscope Blade #0 – Conventional	WI-70001971	\$11.99	66
Miller Laryngoscope Blade #1 Conventional	WI-70001970	\$11.99	67
Miller Laryngoscope Blade #2 Conventional	WI-70001969	\$12.99	67
Miller Laryngoscope Blade #3 Conventional	WI-70001968	\$12.99	68
Miller Laryngoscope Blade #4 Conventional	WI-70001967	\$12.99	68
No. 4 Surgical Scalpel Handle	WI-34516917	\$1.50	69
No-Touch Scalpel Blade Remover Box Holds Up to 100 Blades	WI-70002761	\$9.99	69

Contents (continued)

Product	Code	Price	Page
Pack of 10 Scalpel Blades #60 Autopsy Post Mortem Surgical Blade	WI-70001956	\$12.99	70
Pack of 10 Scalpel Blades No. 10 with Black Scalpel Handle #3	WI-70003148	\$10.99	70
Pack of 100 Scalpel Blades #60	WI-70002186	\$69.99	71
Pack of 100 Scalpel Blades #70	WI-70001342	\$69.99	71
Pack of 100 Surgical Blades #10 with Scalpel Knife Handle #3	WI-70002877	\$11.99	72
Pack of 100 Surgical Blades #11 with Scalpel Knife Handle #3	WI-70002876	\$11.99	72
Pack of 100 Surgical Blades #12 with Scalpel Knife Handle #3	WI-70002875	\$11.99	73
Pack of 100 Surgical Blades #15 with Scalpel Knife Handle #3	WI-70002874	\$11.99	73
Pack of 100 Surgical Blades #15C with Scalpel Knife Handle #3	WI-70002873	\$11.99	74
Pack of 100 Surgical Blades #20 with FREE Scalpel Handle #4	WI-70002872	\$11.99	74
Pack of 100 Surgical Blades #21 with Scalpel Knife Handle #4	WI-70002871	\$11.99	75
Pack of 100 Surgical Blades #22 with Scalpel Knife Handle #4	WI-70002870	\$11.99	75
Pack of 100 Surgical Blades #23 with Scalpel Knife Handle #4	WI-70002869	\$11.99	76
Pack of 100 Surgical Blades #24 with Scalpel Knife Handle #4	WI-70002868	\$11.99	76
Pack of 100 Surgical Blades #25 with Scalpel Knife Handle #4	WI-70002867	\$11.99	77
Pack of 20 Scalpel Blades #60 Autopsy Post Mortem Surgical Blade	WI-70001957	\$17.99	77
Pack of 50 Scalpel Blades #60 Autopsy Post Mortem Surgical Blade	WI-70001958	\$35.99	78
PCS Scalpel Handle No 3 Stainless Steel Derma	WI-28189918	\$1.50	78
Post Mortem Knife 6.5" Smooth Handle	WI-70000396	\$19.99	79
Post Mortem Scalpel Autopsy Cartilage Knives	WI-30864085	\$10.00	79
Reusable Stainless Steel Scalpel Handle Surgical Blade Holder	WI-45266366	\$2.00	80
Round Scalpel Handle #3 Straight Fits Blade 10-16	WI-70001008	\$9.99	80
Scalpel Blade Handle	WI-45143074	\$0.80	81
Scalpel Handle #3	WI-70003147	\$8.99	81
Scalpel Handle #3 & #4 – Double Ended	WI-70000162	\$10.99	82
Scalpel Handle #3L	WI-70003146	\$9.99	82
Scalpel Handle #3LA Angled	WI-70003145	\$9.99	83
Scalpel Handle #4	WI-70003144	\$8.99	83
Scalpel Handle #4L	WI-70003143	\$9.99	84
Scalpel Handle #4LA Angled	WI-70003142	\$9.99	84
Scalpel Handle #5 with Green Plastic Grip	WI-70001018	\$8.99	85
Scalpel Handle #6 with Red Plastic Grip	WI-70001019	\$8.99	85
Scalpel Handle #7 Surgical Scalpel Knife Handle Medical Dental ENT	WI-70003141	\$9.99	86
Scalpel Handle #8 for Scalpel Blades #60 & #70	WI-70002185	\$9.99	86
Scalpel Handle No 4 Surgical Instrument Set / Sialkot	WI-34528863	\$1.50	87
Sickle Knife Surgical Cutting Tool for Fine Tissue Separation	WI-47139992	\$4.80	87
Sickle Knife Surgical Instrument for Precision Tissue Cutting Dissection	WI-47143926	\$5.00	88
Sickle Surgical Knife Precision Cutting Instrument for Controlled Tissue	WI-47157822	\$4.90	88
Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #0	WI-70002410	\$9.99	89
Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #00	WI-70002411	\$11.24	89
Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #1	WI-70002409	\$9.99	90
Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #2	WI-70002408	\$9.99	90

Contents (continued)

Product	Code	Price	Page
Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #3	WI-70002407	\$9.99	91
Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #4	WI-70002406	\$9.99	91
Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #5	WI-70002405	\$9.99	92
Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #6	WI-70002404	\$9.99	92
Single-Use Surgical Scalpel Handle Stainless Steel Blade	WI-28643901	\$1.50	93
Sinus Sickle Knife 7.5" Curved Blunt Tip	WI-70000626	\$8.99	93
Smith Peterson Osteotome 8" Straight 13mm Blade	WI-70001722	\$26.99	94
Smith Peterson Osteotome 8" Straight 16mm Blade	WI-70001721	\$27.99	94
Smith Peterson Osteotome 8" Straight 19mm Blade	WI-70001720	\$28.99	95
Smith Peterson Osteotome 8" Straight 25mm Blade	WI-70001719	\$29.99	95
Smith Peterson Osteotome 8" Straight 32mm Blade	WI-70001718	\$32.99	96
Smith Peterson Osteotome 8" Straight 6.4mm Blade	WI-70003272	\$24.99	96
Smith Peterson Osteotome 8" Straight 9.5mm Blade	WI-70001723	\$25.99	97
Solid Forged Scalpel Stainless Steel Surgical Instrument	WI-46905810	\$2.00	97
Solid Forged Surgical Scalpel Reusable Stainless Steel	WI-46897888	\$1.90	98
Stainless Steel Scalpel Graduated Handle #3	WI-34473973	\$1.50	98
Stainless Steel Scalpel Handle No. 3 Surgical BP	WI-34488920	\$1.50	99
Stainless Steel Scalpel-Low Cost	WI-06301454	\$2.00	99
Stainless Steel Scalpels Dull Polish Finish	WI-06318441	\$1.00	100
Stainless Steel Surgical Scalpel Blade	WI-35052119	\$0.49	100
Stainless Steel Surgical Scalpel Blade Holder	WI-45144024	\$0.50	101
Stainless Steel Surgical Scalpels Instrument Set Customizable	WI-28610950	\$1.50	101
Steel Disposable Surgical Scalpel Handle Instrument	WI-28226866	\$1.50	102
Steel Double Ended Scalpel Handle 16cm/21cm Blade Size	WI-34469968	\$1.50	102
Steel Scalpel and Surgical Knife Handles No.3 & No.4	WI-28292250	\$1.50	103
Steel Surgical Scalpel Precision Instrument Set	WI-28623886	\$1.50	103
Sterile Stainless Steel #10 Scalpel Blades #3 Handle	WI-34501066	\$1.50	104
Surgical Instruments Medical Steel Scalpel Blade	WI-28603915	\$1.50	104
Surgical Scalpel Handle Precision Blade Holder	WI-45257446	\$2.00	105
Surgical Sterile Steel Scalpel Blade Box Compatible Handle	WI-28617921	\$1.50	105
Thomas Uterine Curette 11" (27.9cm) Blunt Blade Size #1	WI-70001112	\$8.99	106
Thomas Uterine Curette 11" (27.9cm) Blunt Blade Size #2	WI-70001111	\$8.99	106
Thomas Uterine Curette 11" (27.9cm) Blunt Blade Size #3	WI-70001110	\$8.99	107
Thomas Uterine Curette 11" (27.9cm) Blunt Blade Size #4	WI-70001109	\$8.99	107
Thomas Uterine Curette 11" (27.9cm) Blunt Blade Size #5	WI-70001108	\$8.99	108
Thomas Uterine Curette 11" (27.9cm) Blunt Blade Size #6	WI-70001107	\$8.99	108
Waterproof Stainless Steel Luxury Battery Powered Curved Blade	WI-40141624	\$2.73	109
Wieder Tongue Depressor Large Blade 36mm Wide Serraed	WI-70003428	\$10.99	109
Wieder Tongue Depressor Small Blade 28mm Wide Serraed	WI-70003473	\$9.99	110

**WI-70001128****10 Fahnestock Knife 13cm Small – Double Ended****\$34.99** per pieceTrade price from **\$32.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Fahnestock Knife 5.25" (13cm) Small Wax & Porcelain Carvers Double Ended with Wooden Handle Fahnestock Spatula Plaster Alginate Lab Mixing Tools Dental Instruments.

Fahnestock Wax and Porcelain Carvers is also known as "Fahnestock Knife" used to place and carve modeling wax to produce dentures. Fahnestock Wax Knives are commonly used in a dental lab to carve wax models, mixing of dental plaster, powders & alginate, impression materials. Sharp edged knife used to cut and mix plasters and other laboratory materials. Fahnestock features a comfortable wooden handle for easy grip and double metal working ends for the mixing of dental, plaster, and other dental materials. You can reuse this dental mixing spatula over

time as they will not bend or tarnish on you. made with high quality stainless steel to meet international quality standards. Fahnestock wax knives can be used in all areas of dentistry.

FAHNESTOCK WAX KNIFE FEATURES:

- Wax & Porcelain Carvers for Dental Lab Technicians.
- Overall Length 5.25" (13cm) – Small.
- Double-ended, One Knife End & Other is Spoon End.
- Smooth Edges Minimize the Risk of Hand Injury.
- Flat Wooden Handle for Easy Grip.
- Reusable.
- High Quality Stainless Steel Blades.

PACKAGE INCLUDE**WI-70001005****10 Round Scalpel Handles #3 Fits Blade 10 to 16****\$39.90** per pieceTrade price from **\$36.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Round Scalpel Handles #3 Straight 6" (15.2cm) Long Premium Quality Stainless Steel compatible with Surgical Blades No. 10 to 16 – Reusable.

Scalpel handle, round and tapered style provides greater rotational movement with fingertip control and enhances dexterity. Handle measures 6" (15cm) long. Compatible with standard scalpel blade #10, #11, #12, #15, #15C, #16. This scalpel handle features an ergonomically designed pen-shaped handle that enhances gripping with knurled style cover that minimizes the risk of slippage and reduces fatigue. Made with high quality stainless steel ensure durability. This handle can be sterilized and reused.

- 10 Pcs Round Scalpel Handle #3 Fits Blades No. 10-16

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Scalpel Handle**RESUABLE:** Yes**STERILEABILITY:** Sterile**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.

**PACKAGE INCLUDES**

**WI-70002151****10 Scalpel Blades #60 with FREE Scalpel Knife Handle #8****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Surgical Blades #60, Carbon Steel Autopsy Post Mortem Scalpel Blades with FREE handle #8

These Post Mortem Autopsy Scalpel Blades are made from surgical carbon steel blades and sharpened to perfection for a super clean cut. The precision shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures.

Exceptionally strong and uniformly sharpened to maximum acuity. These Scalpel Sterile Blades used within the fields of Pathology and Histopathology. Used during Autopsy procedures they are also popular with the Morticians.

PACKAGE INCLUDES:

- 10 Pieces of Scalpel Blades #60
- Free Scalpel Handle #8

SCALPEL BLADE #60 FEATURES:

- Convex cutting edge.
- Thicker blades reduce breakage
- Surgical-quality blades
- Individually sealed
- Peelable packing
- Sterile
- Fits #8 scalpel handle

PACKING:

- All blades packed in Individual pouches

Find our extensive selection of **Scalpel Blades and Handles**

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.

**WI-70002882****10 Sterile Surgical Blades #10 with Scalpel Knife Handle #3****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Sterile Surgical Blades #10 with FREE Scalpel Knife Handle #3 ENT Medical Dental Podiatry, Crafts, Dermaplaning & More.

- The Pack of 10 Scalpel Blades #10 is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #3 that can be used for a variety of different blade sizes.
- The #10 Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for the highest quality + value for the price.
- The sharp and precise blades make these Scalpel Blades a great buy for

Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #10 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.



**WI-70002881****10 Sterile Surgical Blades #11 with Scalpel Knife Handle #3****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Sterile Surgical Blades #11 with FREE Scalpel Knife Handle #3 ENT Medical Dental Podiatry, Crafts, Dermaplaning & More.

- The Pack of 10 Scalpel Blades #11 is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #3 that can be used for a variety of different blade sizes.
- The #11 Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for the highest quality + value for the price.
- The sharp and precise blades make these Scalpel Blades a great buy for

Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #11 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.

**WI-70002880****10 Sterile Surgical Blades #12 with Scalpel Knife Handle #3****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Sterile Surgical Blades #12 with FREE Scalpel Knife Handle #3 ENT Medical Dental Podiatry, Crafts, Dermaplaning & More.

- The Pack of 10 Scalpel Blades #12 is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #3 that can be used for a variety of different blade sizes.
- The #12 Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for the highest quality + value for the price.
- The sharp and precise blades make these Scalpel Blades a great buy for

Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #12 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.



**WI-70002879****10 Sterile Surgical Blades #15 with Scalpel Knife Handle #3****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Sterile Surgical Blades #15 with FREE Scalpel Knife Handle #3 ENT Medical Dental Podiatry, Crafts, Dermaplaning & More.

- The Pack of 10 Scalpel Blades #15 is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #3 that can be used for a variety of different blade sizes.
- The #15 Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for the highest quality + value for the price.
- The sharp and precise blades make these Scalpel Blades a great buy for

Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #15 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.

**WI-70002878****10 Sterile Surgical Blades #15C with Scalpel Knife Handle #3****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Sterile Surgical Blades #15C with FREE Scalpel Knife Handle #3 ENT Medical Dental Podiatry, Crafts, Dermaplaning & More.

- The Pack of 10 Scalpel Blades #15C is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #3 that can be used for a variety of different blade sizes.
- The #15C Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for the highest quality + value for the price.
- The sharp and precise blades make these Scalpel Blades a great buy for

Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #15C scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.



**WI-70002887****10 Sterile Surgical Blades #20 with Scalpel Knife Handle #4****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Sterile Surgical Blades #20 with FREE Scalpel Knife Handle #4 ENT Medical Dental Podiatry, Crafts, Dermaplaning & More.

- The Pack of 10 #20 Scalpel Blades is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #4 that can be used for a variety of different blade sizes.
- The #20 Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for the highest quality + value for the price.
- The sharp and precise blades make these Scalpel Blades a great buy for

Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #20 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.

**WI-70002888****10 Sterile Surgical Blades #21 with Scalpel Knife Handle #4****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


10 Sterile Surgical Blades #21 with FREE Scalpel Knife Handle #4 ENT Medical Dental Podiatry, Crafts, Dermaplaning & More.

- The Pack of 10 #21 Scalpel Blades is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #4 that can be used for a variety of different blade sizes.
- The #21 Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for the highest quality + value for the price.
- The sharp and precise blades make these Scalpel Blades a great buy for

Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #21 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.



**WI-70002886****10 Sterile Surgical Blades #22 with Scalpel Knife Handle #4****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Sterile Surgical Blades #22 with FREE Scalpel Knife Handle #4 ENT Medical Dental Podiatry, Crafts, Dermaplaning & More.

- The Pack of 10 #22 Scalpel Blades is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #4 that can be used for a variety of different blade sizes.
- The #22 Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for the highest quality + value for the price.
- The sharp and precise blades make these Scalpel Blades a great buy for

Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #22 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.

**WI-70002885****10 Sterile Surgical Blades #23 with Scalpel Knife Handle #4****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Sterile Surgical Blades #23 with FREE Scalpel Knife Handle #4 ENT Medical Dental Podiatry, Crafts, Dermaplaning & More.

- The Pack of 10 #23 Scalpel Blades is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #4 that can be used for a variety of different blade sizes.
- The #23 Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for the highest quality + value for the price.
- The sharp and precise blades make these Scalpel Blades a great buy for

Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #23 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.

**WI-70002884****10 Sterile Surgical Blades #24 with Scalpel Knife Handle #4****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Sterile Surgical Blades #24 with FREE Scalpel Knife Handle #4 ENT Medical Dental Podiatry, Crafts, Dermaplaning & More.

- The Pack of 10 #24 Scalpel Blades is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #4 that can be used for a variety of different blade sizes.
- The #24 Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for the highest quality + value for the price.
- The sharp and precise blades make these Scalpel Blades a great buy for

Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #24 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.

**WI-70002883****10 Sterile Surgical Blades #25 with Scalpel Knife Handle #4****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

10 Sterile Surgical Blades #25 with FREE Scalpel Knife Handle #4 ENT Medical Dental Podiatry, Crafts, Dermaplaning & More.

- The Pack of 10 #25 Scalpel Blades is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #4 that can be used for a variety of different blade sizes.
- The #25 Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for the highest quality + value for the price.
- The sharp and precise blades make these Scalpel Blades a great buy for

Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #25 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.



**WI-70002187****100 Scalpel Blades #60 with FREE Scalpel Knife Handle #8****\$76.99** per pieceTrade price from **\$70.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 Scalpel Blades #60 with FREE Scalpel Knife Handle #8 Carbon Steel Autopsy Post Mortem Scalpel Blades Set for Pathology and Histopathology.

These post mortem PM60 autopsy scalpel blades are made from surgical carbon steel sharpened to perfection for a super clean cut. The precision shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Exceptionally strong and uniformly sharpened to maximum acuity. These scalpel sterile blades used within the fields of Pathology and Histopathology. Used during Autopsy procedures they are also popular with the Morticians.

SET INCLUDED:

- 100 Pcs Scalpel Blades #60
- Free Scalpel Handle #8

SCALPEL BLADE #60 FEATURES:

- Convex Cutting Edge.
- Thicker Blades Reduce Breakage.
- Surgical-quality Blades.
- Individually Sealed.
- Peelable Packing.
- Sterile.
- Fits #8 Scalpel Handle.

PACKING:

- All blades packed in Individual pouches

Find our extensive selection of **Scalpel Blades and Handles**

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger

**WI-70001341****100 Scalpel Blades #70 with Scalpel Knife Handle #8****\$76.99** per pieceTrade price from **\$70.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 Surgical Blades #70, Carbon Steel Autopsy Post Mortem Scalpel Blades with Scalpel Knife Handle #8

The number 70 blade is primarily used for postmortem, autopsy and other research purposes. A smaller version of the number 60B blade, this blade features a rounded tip to reduce puncture risk. The number 70 blade will fit the number 8 scalpel handles. Exceptionally strong and uniformly sharpened to maximum acuity. These post mortem autopsy scalpel blades are made from surgical carbon steel blades and sharpened to perfection for a super clean cut. These scalpel sterile blades used within the fields of Pathology and Histopathology. Used

during Autopsy procedures they are also popular with the Morticians.

PACKAGE INCLUDED:

- 100 Pcs Scalpel Blades #70
- 1 Pc Scalpel Handle #8

SCALPEL BLADE #70 FEATURES:

- Rounded tip.
- Thicker blades reduce breakage
- Surgical-quality blades
- Individually sealed
- Peelable packing
- Sterile
- Fits #8 scalpel handle

PACKING:

- All blades packed in Individual pouches

Find our extensive selection of **Scalpel Blades and Handles**





WI-70003825

12 Dental Lab Plaster Knife #12R New

\$71.68 per piece

Trade price from \$65.95
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

12 Dental Lab Plaster Knife #12R New is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$71.68, falling to \$65.95 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003824

12 Dental Lab Plaster Knife #7R New

\$69.30 per piece

Trade price from \$63.76
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

12 Dental Lab Plaster Knife #7R New is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$69.30, falling to \$63.76 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003981****2 Molt Mouth Gag 3.5" + 4.5" Pediatric & Child Curved Blades****\$16.99** per pieceTrade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Molt Mouth Gag Set 3.5" & 4.5" Stainless Steel Surgical Pediatric & Child Mouth Opener for Anesthesia Dental, ENT, Oral & Maxillofacial Instruments.

Molt Mouth Gag used to hold open patient's mouth during extraction or other surgical procedures and reduce TMJ trauma for sedated patients. This is a very useful dentist device when working in the oral cavity. It is fitted between the upper and lower jaws to prevent the mouth from closing during operative procedures of the mouth or throat.

DENTAL MOLT MOUTH GAG FEATURES

UNIQUE DESIGN: The highly polished curved metal ends provide the same

feel as silicone pads without having to worry about the damage that silicone pads face in the autoclave.

LOCKING MECHANISM: The locking mechanism offers more controlled working and the locking mechanism is designed with very small teeth that offer more range in open positions.

HIGH-QUALITY STAINLESS STEEL: Made with medical grade durable stainless steel. The Premium Quality Stainless Steel composition of Dental Molt Mouth Gag 3.5" & 4.5" offers a long life of strong performance, as they are built to last.

FUNCTIONAL & PRECISION: Features a functional design to efficiently complete the task at hand and make your and your subject's life easy. High Degree of Precision and Flexibility while conducting the Clinical Procedure

**2
PCS****WI-70001007****2 Round Scalpel Handles #3 Fits Blade 10 to 16****\$13.99** per pieceTrade price from **\$12.87**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 2 Round Scalpel Handles #3 Straight 6" (15.2cm) Long Premium Quality Stainless Steel compatible with Surgical Blades No. 10 to 16 – Reusable.

Scalpel handle, round and tapered style provides greater rotational movement with fingertip control and enhances dexterity. Handle measures 6" (15cm) long. Compatible with standard scalpel blade #10, #11, #12, #15, #15C, #16. This scalpel handle features an ergonomically designed pen-shaped handle that enhances gripping with knurled style cover that minimizes the risk of slippage and reduces fatigue. Made with high quality stainless steel ensure durability. This handle can be sterilized and reused.

PACKAGE INCLUDES

- 2 Pcs Round Scalpel Handle #3 Fits Blades No. 10-16

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Scalpel Handle

RESUABLE: Yes

STERILEABILITY: Sterile

RUST PREVENTION PROCEDURE:

Passivated

ULTRASONIC CLEANED: Yes

TESTS PERFORMED: Boil Test,

Performance Test, Shape Test

QC PASSED: Yes

PACKING DETAILS: Polybag Packing

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.



WI-01210325

20 Pcs Advanced Stainless Steel Medical

\$4.00 per piece

Minimum order 2 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

20 Pcs Advanced Stainless Steel Medical is a scalpels & knife we forge, grind and finish in our own Sialkot workshop. Instruments of this type make precise incisions and separate tissue planes.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$4.00.

Minimum order 2 sets, and the price moves with quantity. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-70002150

20 Scalpel Blades #60 with FREE Scalpel Knife Handle #8

\$22.99 per piece

Trade price from **\$21.15**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 20 Surgical Blades #60, Carbon Steel Autopsy Post Mortem Scalpel Blades with FREE handle #8

These Post Mortem Autopsy Scalpel Blades are made from surgical carbon steel blades and sharpened to perfection for a super clean cut. The precision shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Exceptionally strong and uniformly sharpened to maximum acuity. These Scalpel Sterile Blades used within the fields of Pathology and Histopathology. Used during Autopsy procedures they are also popular with the Morticians.

PACKAGE INCLUDES:

- 20 Pieces of Scalpel Blades #60
- Free Scalpel Handle #8

SCALPEL BLADE #60 FEATURES:

- Convex cutting edge.
- Thicker blades reduce breakage
- Surgical-quality blades
- Individually sealed
- Peelable packing
- Sterile
- Fits #8 scalpel handle

PACKING:

- All blades packed in Individual pouches

Find our extensive selection of **Scalpel Blades and Handles**

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.





WI-36732508

20cm Autopsy Knife Fixed Blade Stainless Steel Post Mortem

\$18.00 per piece

Trade price from **\$12.00**
Minimum order 20 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture 20cm Autopsy Knife Fixed Blade Stainless Steel Post Mortem for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type make precise incisions and separate tissue planes. The 20cm size shown is our standard; other sizes in the same pattern are made to order.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$18.00, falling to \$12.00 at 20 pieces — about 33% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36732510

20cm Fixed Blade Autopsy & Medical Dissection Knife

\$18.00 per piece

Trade price from **\$12.00**
Minimum order 20 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

20cm Fixed Blade Autopsy & Medical Dissection Knife — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type make precise incisions and separate tissue planes. The 20cm size shown is our standard; other sizes in the same pattern are made to order.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$18.00, falling to \$12.00 at 20 pieces — about 33% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36707577

20cm Fixed Blade Autopsy Knife

\$18.00 per piece

Trade price from \$12.00
Minimum order 20 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of 20cm Fixed Blade Autopsy Knife, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes. The 20cm size shown is our standard; other sizes in the same pattern are made to order.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$18.00, falling to \$12.00 at 20 pieces — about 33% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003820

3 Dental Lab Plaster Knife #12R New

\$22.79 per piece

Trade price from \$20.97
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

3 Dental Lab Plaster Knife #12R New is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$22.79, falling to \$20.97 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003819

3 Dental Lab Plaster Knife #7R New

\$20.21 per piece

Trade price from **\$18.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

3 Dental Lab Plaster Knife #7R New is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$20.21, falling to \$18.59 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001317

3 Finger Ring Cutter Replacement Blade – Small

\$10.99 per piece

Trade price from **\$10.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Finger Ring Cutter Replacement Blade Small 7.8" (22mm) Diameter Emergency Rescue First Aid Circular Saw Blade Stainless Steel Surgical Jewelry Instrument.

The ring cutter blade is a special instrument used to cut and remove rings from swollen fingers with minimum discomfort and effort. The stainless-steel blade can quickly and safely remove rings from fingers.

Finger ring cutter is ideal for cutting the ring or band safely off from finger. Great for doctors, emergency rooms, law enforcement. Used by jewelers, doctors, emergency rooms, law enforcement agencies and more.

Warning: Product (rotary blade) has functional sharp edges. Contact with

blade may result in injury.

FINGER RING CUTTER BLADE FEATURES

- Circular Saw Blade.
- Razor Sharp Micro Cutting.
- Diameter: 7/8" (22mm).
- Reusable
- Made with High Quality Stainless Steel.

PACKAGE INCLUDES

- 3 Pcs Replacement Circular Saw Blade ONLY – Small


TECHNICAL INFORMATION :

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Jewelry Tools

APPLICATION: Ring Cutting Emergency Tool

REUSABLE: Yes





WI-70001329

3 McCoy Fiber Optic Flex Tip Laryngoscope Blades #2, #3 & #4

\$49.99 per piece

Trade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 McCoy Laryngoscope Flexi-Tip Fiber-Optic Led Blade #2, #3 & #4 for Difficult Intubation German Grade Stainless Steel Anesthesiology Medical Instrument

Fibre Optic Laryngoscope Blades (McCoy Flex Tip Profile) feature a low shoulder profile for better visualization and easier intubation of patients.


McCoy Fiber Optic Flex Tip blade have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Moreover, the McCoy blade slimer than the Macintosh blade from the tip and handle sides.

These are reusable blades, extremely useful for better visualization of the vocal cords during difficult intubations. McCoy laryngoscope is used in

patients when difficulty in elevating epiglottis is encountered due to limited mouth opening and activation of blade tip elevates the epiglottis and visualization of vocal cords.

MCCOY FIBER OPTIC FLEX TIP BLADE FEATURES

- Curved Blade With Low Shoulder Profile for Better Visualization
- Mccoy Flex Tip Blades are Slimer Than Macintosh Curved Blade From Tip and Handle Sides.
- Hinged Tip Allows Elevation of the Distal Structures (Epiglottis).
- Extremely Usfull for Difficult Intubation When Difficulty in Elevating Epiglottis.
- Mccoy Blades Used to Manage the Airway of Patients With Limited



WI-70000203

4 Macintosh Conventional Blades #1, #2, #3, #4

\$35.99 per piece

Trade price from **\$33.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 4 Macintosh Conventional Blades Size #1, #2, #3 and #4 – Reusable – Autoclavable German Grade Stainless Steel Laryngoscope Blades EMT Anesthesia Intubation Surgical Instruments

Macintosh Laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.


MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.

- Flange Reduced for Less Trauma.
- Blades Compatible with All Conventional Laryngoscope Systems..
- High Quality Stainless Steel.
- Integrated Light Bulb Directly on the Blad to Providing Light to the Distal Tip.
- Reusable Blade.
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula.
- Available In A Range Of Sizes To Meet Patient Needs.

PACKAGE INCLUDES:

- 1 Pc Conventional Macintosh Blade – Size 1 – Infant
- 1 Pc Conventional Macintosh Blade – Size 2 – Child



**WI-70001345****4 Macintosh Fiber Optic Blades #1, #2, #3, #4****\$44.99** per pieceTrade price from **\$41.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 4 Macintosh Fiber Optic Blades Size #1, #2, #3 and #4. Reusable, Autoclavable German Grade Stainless Steel Laryngoscope Blades EMT Anesthesia Intubation Surgical Instruments

Macintosh Laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.


MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.

- Flange Reduced for Less Trauma.
- Blades Compatible with All Fiber Optic Laryngoscope Systems..
- High Quality Stainless Steel.
- Fiberoptic light transmission provides cool, obstruction free viewing
- Reusable Blade
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula
- Removable light pipe facilitates cleaning and service.

PACKAGE INCLUDES:

- 1 Pc Fiber Optic Laryngoscope Blade – Macintosh; Size 1, Infant; Overall Length: 94 mm
- 1 Pc Fiber Optic Laryngoscope Blade – Macintosh; Size 2, Child; Overall Length: 114 mm

**WI-70001358****5 Fiber Optic Miller Blades #0, #1, #2, #3, #4****\$54.99** per pieceTrade price from **\$50.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Standard Fiber Optic Miller Laryngoscope Blades #0, #1, #2, #3 & #4 Set EMT Anesthesia Intubation Surgical Instruments.

The Miller Blades are the most popular type of the straight blade. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-intubate patients provides a better view of the vocal cords. The straight Miller laryngoscope blade is traditionally recommended for intubation in infants, due to the large size and flexibility of the infant epiglottis

MILLER FIBER OPTIC LARYNGOSCOPE BLADE**FEATURES:**

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.
- Blades Compatible with All Conventional Laryngoscope Systems.
- High Quality German Grade Stainless Steel.
- Fiber optic light transmission provides cool, bright, and focused obstruction-free viewing
- Matt Satin Finish.
- Fully Autoclavable Blade.

SET INCLUDES:

- 1 Miller Fiber Optic Blade Size #0 Neonate; Overall Length: 80 mm



**5
PCS****WI-70001316****5 Finger Ring Cutter Replacement Blade – Small****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Finger Ring Cutter Replacement Blade Small 7.8" (22mm) Diameter Emergency Rescue First Aid Circular Saw Blade Stainless Steel Surgical Jewelry Instrument.

The ring cutter blade is a special instrument used to cut and remove rings from swollen fingers with minimum discomfort and effort. The stainless-steel blade can quickly and safely remove rings from fingers.

Finger ring cutter is ideal for cutting the ring or band safely off from finger. Great for doctors, emergency rooms, law enforcement. Used by jewelers, doctors, emergency rooms, law enforcement agencies and more.

Warning: Product (rotary blade) has functional sharp edges. Contact with

blade may result in injury.

FINGER RING CUTTER BLADE FEATURES

- Circular Saw Blade.
- Razor Sharp Micro Cutting.
- Diameter: 7.8" (22mm).
- Reusable
- Made with High Quality Stainless Steel.

PACKAGE INCLUDES

- 5 Pcs Replacement Circular Saw Blade ONLY – Small

TECHNICAL INFORMATION :

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Jewelry Tools

APPLICATION: Ring Cutting Emergency Tool

REUSABLE: Yes

**WI-70001359****5 Miller Laryngoscope Blades #0, #1, #2, #3, #4 – Conventional****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Standard Conventional Miller Laryngoscope Blades #0, #1, #2, #3 & #4 Set EMT Anesthesia Intubation Surgical Instruments.

The Miller laryngoscope is a straight blade designed to obtain a view of the vocal cords by directly lifting the epiglottis. It has useful application in 'floppy' airways making it popular within paediatric anaesthesia. The unique blade design also makes it ideal for a range of medical approaches. This set allows the physicians to treat the patient without removing them from the MRI suite.

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve

the Lifting of the Epiglottis.

- Round Edges Prevent From Injury Risk.
- Blades Compatible with All Conventional Laryngoscope Systems.
- High Quality Stainless Steel.
- LED Light Bulb Can be Easily Replaced
- Matt Satin Finish.
- Fully Autoclavable Blade.

NOTE: Replacement bulbs sold separately.

PACKAGE INCLUDE:

- 1 Miller Conventional Blade Size #0 Neonate; Overall Length: 80 mm
- 1 Miller Conventional Blade Size #1 Infant; Overall Length: 105 mm



**5
PCS****WI-70001006****5 Round Scalpel Handles #3 Fits Blade 10 to 16****\$29.99** per pieceTrade price from **\$27.59**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Round Scalpel Handles #3 Straight 6" (15.2cm) Long Premium Quality Stainless Steel compatible with Surgical Blades No. 10 to 16 – Reusable.

Scalpel handle, round and tapered style provides greater rotational movement with fingertip control and enhances dexterity. Handle measures 6" (15cm) long. Compatible with standard scalpel blade #10, #11, #12, #15, #15C, #16. This scalpel handle features an ergonomically designed pen-shaped handle that enhances gripping with knurled style cover that minimizes the risk of slippage and reduces fatigue. Made with high quality stainless steel ensure durability. This handle can be sterilized and reused.

PACKAGE INCLUDES

- 5 Pcs Round Scalpel Handle #3 Fits Blades No. 10-16

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Scalpel Handle**RESUABLE:** Yes**STERILEABILITY:** Sterile**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.

**50
PCS****WI-70001959****50 Scalpel Blades #60 with FREE Scalpel Handle #8****\$43.99** per pieceTrade price from **\$40.47**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 50 Surgical Blades #60, Carbon Steel Autopsy Post Mortem Scalpel Blades with FREE handle #8

These Post Mortem Autopsy Scalpel Blades are made from surgical carbon steel sharpened to perfection for a super clean cut. The precision shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Exceptionally strong and uniformly sharpened to maximum acuity. These Scalpel Sterile Blades used within the fields of Pathology and Histopathology. Used during Autopsy procedures they are also popular with the Morticians.

SCALPEL BLADES #60 PACKAGE INCLUDES:

- 50 Pcs Scalpel Blades #60

- Free Scalpel Handle #8

SCALPEL BLADE #60 FEATURES:

- Convex cutting edge.
- Thicker blades reduce breakage
- Surgical-quality blades
- Individually sealed
- Peelable packing
- Sterile
- Fits #8 scalpel handle

PACKING:

- All blades packed in Individual pouches

Find our extensive selection of **Scalpel Blades and Handles**

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger





6 Pieces

WI-70003815

6 Dental Lab Plaster Knife #12R New

\$41.24 per piece

Trade price from \$37.94
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

6 Dental Lab Plaster Knife #12R New is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$41.24, falling to \$37.94 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



6 Pieces

WI-70003814

6 Dental Lab Plaster Knife #7R New

\$37.64 per piece

Trade price from \$34.63
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

6 Dental Lab Plaster Knife #7R New is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$37.64, falling to \$34.63 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001004****Adjustable Scalpel Handle #3 – 180° Rotation with 7 Adjustable Angles****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Adjustable Scalpel Handle #3 Fits Surgical Blades #10, #11, #12, #14, #15, #16 – Round Pattern 7 Degrees Adjustment Locking Points with Rotation Angles 0°, 30°, 60°, 90°, 120°, 150°, 180° Stainless Steel Surgical Knife Handle for Medical, Surgical and DIY Projects – Reusable.

Adjustable Scalpel Handle No. 3 is a multi-functional tool which can be used with a wide range of surgical blades #10, #11, #12, #14, #15, and #16. This is multi-angled scalpel handle with 180 degree range to easily maneuver around corners, hard to reach areas, and tight spaces. With 7 point locking rotations include 0°, 30°, 60°, 90°, 120°, 150°, 180° angles handle position can be adjusted for user-preferred angle resulting in high

degree of precision and flexibility while conducting procedures. No-slip grooved handle provides a strong grip to perform delicate work and/or detailed inspections. Ergonomic design helps reduce hand fatigue during repetitive tasks. Can be easily clean and disinfected for reused purpose.

ADJUSTABLE SCALPEL HANDLE #3 FEATURES

- Scalpel Handle No. 3 Fits Surgical Blades #10, #11, #12, #14, #15, #16.
- 0°-180° Degree Rotation
- 7 Adjustable Point Locking Positions Including 0°, 30°, 60°, 90°, 120°, 150°, 180° Angles.
- Blades Can Slides Securely and Smoothly

**WI-70001016****Bard Parker Scalpel Handle #3 with Ruler****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Bard Parker Scalpel Handle #3 with Ruler Stainless Steel Flat Handle with 0 – 6 Inch (120mm) Scale, Overall Length 5" – Reusable.

The Bard Parker style scalpel handle #3 with metric ruler is used for cutting gingival tissue and making surgical incisions during oral surgeries. Handle is 5" long compatible with standard scalpel blade 10-15 sizes. The standard scalpel handle graduated 1-6 scale with serrated textured grip section which allow precise control during surgical procedures, woodworking, or crafting. Can also used for surgery and art and craft work. B.P scalpel can be reused over time as they will not bend or tarnish on you. Made with high quality stainless steel to meet excellent international quality standard.

PACKAGE INCLUDES

- 1 Pc B.P Scalpel Handle #3 with Ruler

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Scalpel Handle**RESUABLE:** Yes**STERILEABILITY:** Sterile**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.





WI-70001015

Bard Parker Scalpel Handle #4 with Ruler

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Bard Parker Scalpel Handle #4 with Ruler Stainless Steel Flat Handle with 0 – 6 Inch (120mm) Scale, Overall Length 5.5" – Reusable.

The Bard Parker style scalpel handle #4 with metric ruler is used for cutting and making incisions during surgical procedures. Handle is 5.5" long compatible with standard scalpel blade 20- 25 sizes. The standard scalpel handle graduated 1-6 scale with serrated textured grip section which allow precise control during surgical procedures, woodworking, or crafting. Can also used for surgery and art and craft work. B.P scalpel can be reused over time as they will not bend or tarnish on you. Made with high quality stainless steel to meet excellent international quality standard.

PACKAGE INCLUDES

- 1 Pc B.P Scalpel Handle #4 with Ruler

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENT TYPE: Scalpel Handle
RESUABLE: Yes
STERILEABILITY: Sterile
RUST PREVENTION PROCEDURE: Passivated
ULTRASONIC CLEANED: Yes
TESTS PERFORMED: Boil Test, Performance Test, Shape Test
QC PASSED: Yes
PACKING DETAILS: Polybag Packing

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.



WI-70003400

Barraquer Blade Breaker Surgical Supplies

\$18.59 per piece

Trade price from **\$17.10**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Barraquer Blade Breaker Surgical Supplies is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$18.59, falling to \$17.10 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000003****Black Scalpel Blade Handle #3****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Tactical Black Scalpel Handle #3
Competible Scalpel Blades from #10 to #15 Premium Quality Stainless Steel Reusable Matte Black Surgical Knife Handle.

Premium quality, rust proof black scalpel handle no.3 is lightweight and durable, fits for surgical blades No. 10, 11, 12, 15 and 15C. The scalpel handle #3 ergonomically designed for hand precise fingertip control required for hand surgery and other delicate cutting or dissection. Surgical Knife Handle #3 is usually used with smaller blades that have a curved cutting edge, like the #10 or #15 replaceable surgical blades. These handles and blade work well for microdissection or fine surgeries when making precise incisions. With a matte finish for improved grip and visibility in

sterile lab environments, this handle manufactured from surgical grade stainless steel and suitable for autoclaving to promote reusability.

BLACK SCALPEL HANDLE #3 FEATURES

- 5" (12.7cm) Overall Legth.
- Matte Black Finish.
- Knife Handle Fits Blades From Size 10 to 15.
- Textured Gripping Surface for a Firm Hold.
- Blades Slide on Smoothly and Securely.
- Reusable Handle.
- Crafted with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Black Scalpel Blade Handle #3

**WI-70000470****Buck Ear Curette 6.5" Straight Sharp Ring Blade Size #0****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Buck Ear Curette 6.5" (16.5cm) Straight Sharp Ring Blade Size #0 (2mm) with Round Pattern Knurled Grip Handle Stainless Steel ENT Surgical Veterinary Instruments.


Buck Ear Curettes are specialized surgical tools used in ear surgeries to safely remove tumors, clean infections, treat pain and prevent hearing loss through biopsies. This Ear Curette offer surgical options with blunt ring blades for sampling and sharp blades for deep cleaning, available with straight or angled necks.

Buck Ear 6.5" (16.5cm) curettes, made from premium grade stainless steel which making it reusable and autoclavable, are designed to effectively treat deep diseases and have slender shafts for clear working

tips during delicate samples. The veterinary practitioners also use Buck Ear Curettes for ENT procedures like vertical cleaning cealing to remove ear wax and other foreign objects. The long handle to give precise functioning during multiple ENT practices.

BUCK EAR CURETTE FEATURES

- 6.5" (16.5cn) Long.
- Straight Sharp Profile.
- Size #0 (2mm) Tip.
- Slender Shafts For Deep Surgical Access.
- Non-Slip Handle For Controlled Grips.
- Multiple Model and Sizes For Diverse Indications.
- Designed for Ear Canal Cleaning & Removes Cerumen.



**WI-70000471****Buck Ear Curette 6.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Buck Ear Curette 6.5\" (16.5cm) Straight Sharp Ring Blade Size #00 (1.5mm) with Round Pattern Knurled Grip Handle Stainless Steel ENT Surgical Veterinary Instruments.


Buck Ear Currettes are specialized surgical tools used in ear surgeries to safely remove tumors, clean infections, treat pain and prevent hearing loss through biopsies. This Ear Curette offer surgical options with blunt ring blades for sampling and sharp blades for deep cleaning, available with straight or angled necks.

Buck Ear 6.5\" (16.5cm) currettes, made from premium grade stainless steel which making it reusable and autoclavable, are designed to effectively treat deep diseases and have slender shafts for clear working

tips during delicate samples.The veterinary practitioners also use Buck Ear Currettes for ENT procedures like vertical cleaning cealing to remove ear wax and other foreign objects. The long handle to give precise functioning during multiple ENT practices.

BUCK EAR CURETTE FEATURES

- 6.5\" (16.5cn) Long.
- Straight Sharp Profile.
- Size #00 (1.5mm) Tip.
- Slender Shafts For Deep Surgical Access.
- Non-Slip Handle For Controlled Grips.
- Multiple Model and Sizes For Diverse Indications.
- Designed for Ear Canal Cleaning & Removes Cerumen.

**WI-70000469****Buck Ear Curette 6.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Buck Ear Curette 6.5\" (16.5cm) Straight Sharp Ring Blade Size #1 (2.5mm) with Round Pattern Knurled Grip Handle Stainless Steel ENT Surgical Veterinary Instruments.


Buck Ear Currettes are specialized surgical tools used in ear surgeries to safely remove tumors, clean infections, treat pain and prevent hearing loss through biopsies. This Ear Curette offer surgical options with blunt ring blades for sampling and sharp blades for deep cleaning, available with straight or angled necks.

Buck Ear 6.5\" (16.5cm) currettes, made from premium grade stainless steel which making it reusable and autoclavable, are designed to effectively treat deep diseases and have slender shafts for clear working

tips during delicate samples.The veterinary practitioners also use Buck Ear Currettes for ENT procedures like vertical cleaning cealing to remove ear wax and other foreign objects. The long handle to give precise functioning during multiple ENT practices.

BUCK EAR CURETTE FEATURES

- 6.5\" (16.5cn) Long.
- Straight Sharp Profile.
- Size #1 (2.5mm) Tip.
- Slender Shafts For Deep Surgical Access.
- Non-Slip Handle For Controlled Grips.
- Multiple Model and Sizes For Diverse Indications.
- Designed for Ear Canal Cleaning & Removes Cerumen.



**WI-70000468****Buck Ear Curette 6.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Buck Ear Curette 6.5" (16.5cm) Straight Sharp Ring Blade Size #2 (3mm) with Round Pattern Knurled Grip Handle Stainless Steel ENT Surgical Veterinary Instruments.


Buck Ear Currettes are specialized surgical tools used in ear surgeries to safely remove tumors, clean infections, treat pain and prevent hearing loss through biopsies. This Ear Curette offer surgical options with blunt ring blades for sampling and sharp blades for deep cleaning, available with straight or angled necks.

Buck Ear 6.5" (16.5cm) currettes, made from premium grade stainless steel which making it reusable and autoclavable, are designed to effectively treat deep diseases and have slender shafts for clear working

tips during delicate samples. The veterinary practitioners also use Buck Ear Currettes for ENT procedures like vertical cleaning cealing to remove ear wax and other foreign objects. The long handle to give precise functioning during multiple ENT practices.

BUCK EAR CURETTE FEATURES

- 6.5" (16.5cn) Long.
- Straight Sharp Profile.
- Size #2 (3mm) Tip.
- Slender Shafts For Deep Surgical Access.
- Non-Slip Handle For Controlled Grips.
- Multiple Model and Sizes For Diverse Indications.
- Designed for Ear Canal Cleaning & Removes Cerumen.

**WI-70000467****Buck Ear Curette 6.5\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Buck Ear Curette 6.5" (16.5cm) Straight Sharp Ring Blade Size #3 (3.5mm) with Round Pattern Knurled Grip Handle Stainless Steel ENT Surgical Veterinary Instruments.


Buck Ear Currettes are specialized surgical tools used in ear surgeries to safely remove tumors, clean infections, treat pain and prevent hearing loss through biopsies. This Ear Curette offer surgical options with blunt ring blades for sampling and sharp blades for deep cleaning, available with straight or angled necks.

Buck Ear 6.5" (16.5cm) currettes, made from premium grade stainless steel which making it reusable and autoclavable, are designed to effectively treat deep diseases and have slender shafts for clear working

tips during delicate samples. The veterinary practitioners also use Buck Ear Currettes for ENT procedures like vertical cleaning cealing to remove ear wax and other foreign objects. The long handle to give precise functioning during multiple ENT practices.

BUCK EAR CURETTE FEATURES

- 6.5" (16.5cn) Long.
- Straight Sharp Profile.
- Size #3 (3.5mm) Tip.
- Slender Shafts For Deep Surgical Access.
- Non-Slip Handle For Controlled Grips.
- Multiple Model and Sizes For Diverse Indications.
- Designed for Ear Canal Cleaning & Removes Cerumen.



**WI-70003854****Con Scodellino Knife – Large****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Con Scodellino Wax Knife Large 17cm Double Ended Blade with Wooden Handle Premium Quality Stainless Steel Dental Technician Laboratory Instrument.

Con Scodellino (conscodellino) knife is a specialized dental instrument for cutting, mixing, shaping, and carving dental materials such as plaster and wax. Its key uses include plaster work, wax sculpting, and eliminating undercuts on dental casts. The knife features a double-ended design, with a sharp edge for cutting and a scoop for handling materials. Made from surgical-grade stainless steel and equipped with an ergonomic wooden handle for balanced handling and reduced hand fatigue. Con Scodellino wax knife


designed to withstand regular sterilization protocols.

CON SCODELLINO WAX KNIFE FEATURES:

- Wax & Porcelain Carvers for Dental Lab Technicians.
- Overall Length 17cm – Large.
- Double-ended, One Knife End & Other is Spoon End.
- Smooth Edges Minimize the Risk of Hand Injury.
- Lightweight Wooden Handle for Easy Grip.
- Reusable.
- High Quality Stainless Steel Blades.

PACKAGE INCLUDE

- 1 Pc Con Scodellino Knife Large – Double Ended

**WI-70003855****Con Scodellino Knife – Small****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Con Scodellino Wax Knife Small 13cm Double Ended Blade with Wooden Handle Premium Quality Stainless Steel Dental Technician Laboratory Instrument.

Con Scodellino (conscodellino) knife is a specialized dental instrument for cutting, mixing, shaping, and carving dental materials such as plaster and wax. Its key uses include plaster work, wax sculpting, and eliminating undercuts on dental casts. The knife features a double-ended design, with a sharp edge for cutting and a scoop for handling materials. Made from surgical-grade stainless steel and equipped with an ergonomic wooden handle for balanced handling and reduced hand fatigue. Con Scodellino wax knife

designed to withstand regular sterilization protocols.

CON SCODELLINO WAX KNIFE FEATURES:

- Wax & Porcelain Carvers for Dental Lab Technicians.
- Overall Length 5" (13cm) – Large.
- Double-ended, One Knife End & Other is Spoon End.
- Smooth Edges Minimize the Risk of Hand Injury.
- Lightweight Wooden Handle for Easy Grip.
- Reusable.
- High Quality Stainless Steel Blades.

PACKAGE INCLUDE

- 1 Pc Con Scodellino Knife Small – Double Ended





WI-70003811

Dental Lab Plaster Knife #12R New

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Lab Plaster Knife #12R New is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$8.99, falling to \$8.27 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003810

Dental Lab Plaster Knife #7R New

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Dental Lab Plaster Knife #7R New is an instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for action and finish.

Listed at \$8.99, falling to \$8.27 at 21 pieces - about 8% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28678810

Design Disposable Surgical Scalpels Steel

\$1.50 per piece

Trade price from \$0.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Design Disposable Surgical Scalpels Steel is a scalpels & knife we forge, grind and finish in our own Sialkot workshop. Instruments of this type make precise incisions and separate tissue planes.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28667834

Disposable Steel Surgical Scalpels and Blades Instrument

\$1.50 per piece

Trade price from \$0.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Disposable Steel Surgical Scalpels and Blades Instrument is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type make precise incisions and separate tissue planes. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-43380071

Disposable Sterile Surgical Scalpel Stainless Steel

\$0.88 per piece

Trade price from \$0.85
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Disposable Sterile Surgical Scalpel Stainless Steel — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type make precise incisions and separate tissue planes.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$0.88, falling to \$0.85 at 100 pieces — about 3% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000253

Dissecting Scalpel Knife 6" – Fixed Blade

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Dissecting Scalpel Knife 6" (15.2cm) Long with Sharp Fixed Blade Flat Handle Scalpel Knife for Medical Surgical Tissue Dissection Biology Study Lab.

Dissecting scalpel fixed blade knife is a ultra-sharp tool essential for precise tissue cutting in autopsies, forensic pathology, and biological laboratories. These 6-inch dissecting scalpel blade have a rigid, non-detachable blade and flat handle, offering superior control for deep incisions compared to standard scalpels. Dissecting scalpel blade fixed handle knife used in laboratory settings for trimming and sectioning fixed tissues, as well as for educational dissections, detailed work,specimen preparation, gross trimming and anatomical dissection. Made of medical


grade stainless steel, fixed blade scalpel knife are lightweight and designed for continuous heavy-duty use.

DISSECTING SCALPEL FIXED BLADE KNIFE FEATURES

- 6 Inches (15.2 Cm) in Length.
- Sharp, Fixed Blade Designed for Precision.
- Equipped With a Flat Handle for Enhanced Control.
- Ideal for Medical Surgical Tissue Dissection and Biology Study.
- Suitable for Laboratory Use.
- Made with Stainless Steel.

PACKAGE INCLUDED:

- 1 Pc Dissecting Scalpel Knife 6" (15.2cm)





WI-08790844

Dog Grooming Coat Stripping Knife Set 6

\$12.25 per piece

Trade price from \$7.45
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dog Grooming Coat Stripping Knife Set 6, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$12.25, falling to \$7.45 at 10 sets — about 39% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000035

Electrosurgical Pencil 2mm Wide Blade – Single Use

\$14.99 per piece

Trade price from \$13.79
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Electrosurgical Pencil 2mm Wide Blade – Single Use – Diathermy ESU Cautery Pencil Used to Cut and Coagulate Biological Tissues During Electrosurgical Procedures.

Electrosurgical Pencil also known as Cautery Pencil / Monopolar cautery is a disposable, sterile medical device used for cutting and coagulating biological tissues during surgeries. Designed for precision, it features an ergonomic grip, a smooth activation button, and a universal 3-pin connector for compatibility with standard electrosurgical generators. Suitable for various medical fields including general surgery and plastic surgery, these lightweight and individually packed pencils enhance hygiene, efficiency,

and reduce the risk of cross-contamination.

ELECTROSURGICAL PENCIL / CAUTERY PENCIL FEATURES

- 8.25" (21cm) Overall Length.
- 2mm Blade Tip for Cut and Coagulate.
- Hand Switching Push Buttons to Switch Between Power and Precision Settings.
- Three-prong Plug for Compatibility With Standard Electrosurgical Generators.
- Available in 3' Cords for Easy Set Up and Drapeability.
- Single Use, Sterilized and Individually Packed.
- Lightweight and Portable.





WI-34506927

Ergonomic Medical Instrument Set Stainless Steel No. 3/4/7

\$1.50 per piece

Trade price from **\$0.50**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Ergonomic Medical Instrument Set Stainless Steel No. 3/4/7 is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type make precise incisions and separate tissue planes.

The name tells you what you are getting: an ergonomic handle for long lists without fatigue. Supplied as a set — the contents can be adjusted to your

list. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one.

Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price. Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003850

Fahnestock Wax Knife 5.25" – Small

\$8.99 per piece

Trade price from **\$8.27**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Fahnestock Knife 5.25" (13cm) Small Wax & Porcelain Carvers Double Ended with Wooden Handle Fahnestock Spatula Plaster Alginate Lab Mixing Tools Dental Instruments.

Fahnestock Wax and Porcelain Carvers is also known as "Fahnestock Knife" used to place and carve modeling wax to produce dentures. Fahnestock Wax Knives are commonly used in a dental lab to carve wax models, mixing of dental plaster, powders & alginate, impression materials. Sharp edged knife used to cut and mix plasters and other laboratory materials. Fahnestock features a comfortable wooden handle for easy grip and double metal working ends for the mixing of dental, plaster, and other dental materials. You can reuse this dental mixing spatula over

time as they will not bend or tarnish on you. made with high quality stainless steel to meet international quality standards. Fahnestock wax knives can be used in all areas of dentistry.

FAHNESTOCK WAX KNIFE FEATURES:

- Wax & Porcelain Carvers for Dental Lab Technicians.
- Overall Length 5.25" (13cm) – Small.
- Double-ended, One Knife End & Other is Spoon End.
- Smooth Edges Minimize the Risk of Hand Injury.
- Flat Wooden Handle for Easy Grip.
- Reusable.
- High Quality Stainless Steel Blades.

PACKAGE INCLUDE



DOUBLE ENDED

DOUBLE SIZED BLADES MORE CONVENIENCE



WAX & PORCELAIN CARVERS

IDEAL FOR WAX MODELLING MIXING IMPRESSION MATERIAL, POWDERS & ALGinate



WI-70003849

Fahnestock Wax Knife 7" – Large

\$8.99 per piece

Trade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Fahnestock Knife 7" (17.8cm) Large Wax & Porcelain Carvers Double Ended with Wooden Handle Fahnenstock Spatula Plaster Alginate Lab Mixing Tools Dental Instruments.

Fahnestock Wax and Porcelain Carvers is also known as "Fahnenstock Knife" used to place and carve modeling wax to produce dentures. Fahnenstock Wax Knives are commonly used in a dental lab to carve wax models, mixing of dental plaster, powders & alginate, impression materials. Sharp edged knife used to cut and mix plasters and other laboratory materials. Fahnestock features a comfortable wooden handle for easy grip and double metal working ends for the mixing of dental, plaster, and other dental materials. You can reuse this dental mixing spatula over

time as they will not bend or tarnish on you. made with high quality stainless steel to meet international quality standards. Fahnenstock wax knives can be used in all areas of dentistry.

FAHNESTOCK WAX KNIFE FEATURES:

- Wax & Porcelain Carvers for Dental Lab Technicians.
- Overall Length 7" (17.8cm) – Large.
- Double-ended, One Knife End & Other is Spoon End.
- Smooth Edges Minimize the Risk of Hand Injury.
- Flat Wooden Handle for Easy Grip.
- Reusable.
- High Quality Stainless Steel Blades.

PACKAGE INCLUDE



DOUBLE ENDED
DOUBLE SIZED BLADES MORE CONVENIENCE



WAX & PORCELAIN CARVERS
IDEAL FOR WAX MODELING MIXING IMPRESSION MATERIAL, POWDERS & ALGINATE



WI-70001355

Fiber Optic Miller Blade #0 with Laryngoscope Medium Handle

\$28.99 per piece

Trade price from **\$26.67**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Miller Blade #0 with Fiber Optic Laryngoscope Medium Handle 2 Pcs Set EMT Anesthesia Intubation Surgical Instruments.

The Miller laryngoscope is a straight blade designed to obtain a view of the vocal cords by directly lifting the epiglottis. It has useful application in 'floppy' airways making it popular within paediatric anaesthesia. This set allows the physicians to treat the patient without removing them from the MRI suite.

2 PCS MILLER LARYNGOSCOPE SET:

FIBER OPTIC MILLER BLADES: The Miller laryngoscope has a straight blade with a long, curved tip, providing better exposure of the larynx. The

unique blade design also makes it ideal for a range of medical approaches. It comes with different variations such as Miller Laryngoscope Blade, Size 0, 1, 2, 3, and 4 Conventional Style.

CONVENTIONAL LARYNGOSCOPE HANDLE: Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.



**WI-70001354****Fiber Optic Miller Blade #1 with Laryngoscope Medium Handle****\$28.99** per pieceTrade price from **\$26.67**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Miller Blade #1 with Fiber Optic Laryngoscope Medium Handle 2 Pcs Set EMT Anesthesia Intubation Surgical Instruments.

The Miller laryngoscope is a straight blade designed to obtain a view of the vocal cords by directly lifting the epiglottis. It has useful application in 'floppy' airways making it popular within paediatric anaesthesia. This set allows the physicians to treat the patient without removing them from the MRI suite.

2 PCS MILLER LARYNGOSCOPE SET:

FIBER OPTIC MILLER BLADES: The Miller laryngoscope has a straight blade with a long, curved tip, providing better exposure of the larynx. The

unique blade design also makes it ideal for a range of medical approaches. It comes with different variations such as Miller Laryngoscope Blade, Size 0, 1, 2, 3, and 4 Conventional Style.

CONVENTIONAL LARYNGOSCOPE HANDLE: Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.

**WI-70001353****Fiber Optic Miller Blade #2 with Laryngoscope Medium Handle****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Miller Blade #2 with Fiber Optic Laryngoscope Medium Handle 2 Pcs Set EMT Anesthesia Intubation Surgical Instruments.

The Miller laryngoscope is a straight blade designed to obtain a view of the vocal cords by directly lifting the epiglottis. It has useful application in 'floppy' airways making it popular within paediatric anaesthesia. This set allows the physicians to treat the patient without removing them from the MRI suite.

2 PCS MILLER LARYNGOSCOPE SET:

FIBER OPTIC MILLER BLADES: The Miller laryngoscope has a straight blade with a long, curved tip, providing better exposure of the larynx. The

unique blade design also makes it ideal for a range of medical approaches. It comes with different variations such as Miller Laryngoscope Blade, Size 0, 1, 2, 3, and 4 Conventional Style.

CONVENTIONAL LARYNGOSCOPE HANDLE: Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.



**WI-70001352****Fiber Optic Miller Blade #3 with Laryngoscope Medium Handle****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Miller Blade #3 with Fiber Optic Laryngoscope Medium Handle 2 Pcs Set EMT Anesthesia Intubation Surgical Instruments.

The Miller laryngoscope is a straight blade designed to obtain a view of the vocal cords by directly lifting the epiglottis. It has useful application in 'floppy' airways making it popular within paediatric anaesthesia. This set allows the physicians to treat the patient without removing them from the MRI suite.

2 PCS MILLER LARYNGOSCOPE SET:

FIBER OPTIC MILLER BLADES: The Miller laryngoscope has a straight blade with a long, curved tip, providing better exposure of the larynx. The

unique blade design also makes it ideal for a range of medical approaches. It comes with different variations such as Miller Laryngoscope Blade, Size 0, 1, 2, 3, and 4 Conventional Style.

CONVENTIONAL LARYNGOSCOPE HANDLE: Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.

**WI-70001351****Fiber Optic Miller Blade #4 with Laryngoscope Medium Handle****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Miller Blade #4 with Fiber Optic Laryngoscope Medium Handle 2 Pcs Set EMT Anesthesia Intubation Surgical Instruments.

The Miller laryngoscope is a straight blade designed to obtain a view of the vocal cords by directly lifting the epiglottis. It has useful application in 'floppy' airways making it popular within paediatric anaesthesia. This set allows the physicians to treat the patient without removing them from the MRI suite.

2 PCS MILLER LARYNGOSCOPE SET:

FIBER OPTIC MILLER BLADES: The Miller laryngoscope has a straight blade with a long, curved tip, providing better exposure of the larynx. The

unique blade design also makes it ideal for a range of medical approaches. It comes with different variations such as Miller Laryngoscope Blade, Size 0, 1, 2, 3, and 4 Conventional Style.

CONVENTIONAL LARYNGOSCOPE HANDLE: Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.



**WI-70001315**

Finger Ring Cutter 6.5" with Saw Blade Hollow Handle

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Finger Ring Cutter 6.5 Inches (16.5cm) with 7/8" (22mm Diameter) Saw Blade Non-Slip Hollow Grip Handle Stainless Steel Emergency Jewelry Cutting Rescue First Aid Paramedic EMT/EMS Surgical Instrument.

Finger Ring cutter is an emergency tool designed to quickly remove a ring stuck in place due to swelling or injury this ring remover can cut through gold, silver, brass, titanium, and other bands quickly and safely. Ringer cutter made with an extra-large twisting thumb to create more leverage for the user to cut with minimal discomfort and effort.

The small finger guard that slides under the ring to shield it from the ultra-sharp, stainless-steel blade, so you can remove the band while protecting the skin. Our ring cutter can be used by a

wide range of first responders and medical personnel including doctors, paramedics, EMTs, police officers, fire rescue, or even jewelers. This is a ring cutter that saws (utilizes a circular saw blade 7/8"D) ring while still on or off the finger. Made with German surgical stainless steel and can be reused after sterilization.

FINGER RING CUTTER FEATURES

- Overall Approximately 6.5" (16.5cm) in Length.
- Blade Diameter : 7/8" (22mm)
- Circular Razor Sharp Micro Cutting Saw Blade.
- Non-Slip Hollow Grip Handle.
- Easy to Use on Rings on or Off the Finger.

**WI-70000013**

Finger Ring Cutter Replacement Blade – Large

\$8.99 per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Finger Ring Cutter Replacement Blade Large 1-3/16" (30mm Diameter) Emergency Rescue First Aid Circular Saw Blade Stainless Steel Surgical Jewelry Instrument.

The ring cutter blade is a special instrument used to cut and remove rings from swollen fingers with minimum discomfort and effort. The stainless-steel blade can quickly and safely remove rings from fingers.

Finger ring cutter is ideal for cutting the ring or band safely off from finger. Great for doctors, emergency rooms, law enforcement. Used by jewelers, doctors, emergency rooms, law enforcement agencies and more.

Warning: Product (rotary blade) has functional sharp edges. Contact with

blade may result in injury.

FINGER RING CUTTER BLADE FEATURES

- Circular Saw Blade.
- Razor Sharp Micro Cutting.
- Diameter: 1-3/16" (30mm) OD.
- Reusable
- Made with High Quality Stainless Steel.

PACKAGE INCLUDES

- 1 Pc Replacement Circular Saw Blade ONLY – Large


TECHNICAL INFORMATION :

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Jewelry Tools

APPLICATION: Ring Cutting Emergency Tool

AUTOCCLAVABLE/REUSABLE: Yes



**WI-70001318****Finger Ring Cutter Replacement Blade – Small****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Finger Ring Cutter Replacement Blade – Small – Emergency Rescue First Aid Circular Saw Blade Stainless Steel Surgical Jewelry Instrument.

The ring cutter blade is a special instrument used to cut and remove rings from swollen fingers with minimum discomfort and effort. The stainless-steel blade can quickly and safely remove rings from fingers.

Finger ring cutter is ideal for cutting the ring or band safely off from finger. Great for doctors, emergency rooms, law enforcement. Used by jewelers, doctors, emergency rooms, law enforcement agencies and more.


Warning: Product (rotary blade) has functional sharp edges. Contact with blade may result in injury.

FINGER RING CUTTER BLADE FEATURES

- Circular Saw Blade.
- Razor Sharp Micro Cutting.
- Diameter: 7/8" (22mm) OD.
- Reusable
- Made with High Quality Stainless Steel.

PACKAGE INCLUDES

- 1 Pc Replacement Circular Saw Blade ONLY – Small

TECHNICAL INFORMATION :**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Jewelry Tools**APPLICATION:** Ring Cutting Emergency Tool**REUSABLE:** Yes**STERILE:** Non-Sterile**WI-70001676****Foot Callus Shaver Hard Skin Corn Remover with 10 Replacement Blades – Black****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Salon Quality Surgical Grade Stainless Steel Black Pedicure Callous Callus Shaver Foot Hard Tough Skin Corn Remover Tool with 10 Blades. Ideal for removing and reducing Removes Calluses, Corns & Rough Skin. Long easy to use handle comes with Ten (10) replacement blades.

CALLOUS SHAVERS RAZORS FEATURES

MATERIAL: The shaver set is made of good quality steel and plastic. No rust, no deformed foot shaver callus remover, and a long easy to use non-slip handle.

USED FOR: This callus shaver set can remove hard, cracked, dead skin, heel skin calluses, and so on. Perfect gift idea for elders.

EASY TO USE: Removable head design,


very convenient to change, plastic handle touch smooth and comfortable. **PORTABLE:** Our heel shaver for feet is portable and lightweight which is easy to operate, can be used for indoor and outdoor pedicures. the perfect solution to remove callus, helps you save your money and time.

ATTENTION:

Not recommended for people who have thin skin, stay away from children and be cautious while using.

PACKAGE INCLUDE:

- 1 Pc Callous Callus Shaver with Ten (10) Replacement blades – Black

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Pedicure Instruments

**WI-70001675****Foot Callus Shaver Hard Skin Corn Remover with 10 Replacement Blades – Blue****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Salon Quality Surgical Grade Stainless Steel Blue Pedicure Callous Callus Shaver Foot Hard Tough Skin Corn Remover Tool with 10 Blades. Ideal for removing and reducing Removes Calluses, Corns & Rough Skin. Long easy to use handle comes with Ten (10) replacement blades.

very convenient to change, plastic handle touch smooth and comfortable.

PORTABLE: Our heel shaver for feet is portable and lightweight which is easy to operate, can be used for indoor and outdoor pedicures. the perfect solution to remove callus, helps you save your money and time.

CALLOUS SHAVERS RAZORS FEATURES

MATERIAL: The shaver set is made of good quality steel and plastic. No rust, no deformed foot shaver callus remover, and a long easy to use non-slip handle.

USED FOR: This callus shaver set can remove hard, cracked, dead skin, heel skin calluses, and so on. Perfect gift idea for elders.

EASY TO USE: Removable head design,

ATTENTION:

Not recommended for people who have thin skin, stay away from children and be cautious while using.


PACKAGE INCLUDE:

- 1 Pc Callous Callus Shaver with Ten (10) Replacement Blades – Blue Color Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Pedicure

**WI-70001674****Foot Callus Shaver Hard Skin Corn Remover with 10 Replacement Blades – White****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Salon Quality Surgical Grade Stainless Steel White Pedicure Callous Callus Shaver Foot Hard Tough Skin Corn Remover Tool with 10 Blades. Ideal for removing and reducing Removes Calluses, Corns & Rough Skin. Long easy to use handle comes with Ten (10) replacement blades.

very convenient to change, plastic handle touch smooth and comfortable.

PORTABLE: Our heel shaver for feet is portable and lightweight which is easy to operate, can be used for indoor and outdoor pedicures. the perfect solution to remove callus, helps you save your money and time.

CALLOUS SHAVERS RAZORS FEATURES

MATERIAL: The shaver set is made of good quality steel and plastic. No rust, no deformed foot shaver callus remover, and a long easy to use non-slip handle.

USED FOR: This callus shaver set can remove hard, cracked, dead skin, heel skin calluses, and so on. Perfect gift idea for elders.

EASY TO USE: Removable head design,

ATTENTION:

Not recommended for people who have thin skin, stay away from children and be cautious while using.

PACKAGE INCLUDE:


- 1 Pc Callous Callus Shaver with Ten (10) Replacement Blades – White Color Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Pedicure





WI-46940065

Forged Stainless Steel Scalpel Surgical Cutting Instrument for Controlled

\$1.80 per piece

Trade price from **\$0.90**
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Forged Stainless Steel Scalpel Surgical Cutting Instrument for Controlled, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.80, falling to \$0.90 at 50 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-45279261

General Surgery Scalpel Handle Reusable Medical Blade

\$2.00 per piece

Trade price from **\$1.40**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

General Surgery Scalpel Handle Reusable Medical Blade is a scalpels & knife we forge, grind and finish in our own Sialkot workshop. Instruments of this type make precise incisions and separate tissue planes. The name tells you what you are getting: full reusability — clean, autoclave and use again.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$2.00, falling to \$1.40 at 25 pieces — about 30% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34512907

German Stainless Steel Scalpel Handle 145mm 5.34\"/>

\$1.50 per piece

Trade price from \$0.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

German Stainless Steel Scalpel Handle 145mm 5.34\"/>

The name tells you what you are getting: a straight working end for direct, in-line access.

The 145mm size shown is our standard; other sizes in the same pattern are

made to order. Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001068

Heavy Duty Finger Ring Cutter Plier 6.5\"/>

\$11.99 per piece

Trade price from \$11.03
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Heavy Duty Finger Ring Cutter Plier 6.5\"/>

Finger Ringer cutting plier is an emergency tool designed to quickly remove a ring stuck in place due to swelling or injury this ring remover can cut through gold, silver, brass, and other bands quickly and safely. This handheld Ring Cutter is similar to pliers, ring cutter body is made of Stainless Steel. The thin base segment slides between the finger and the band of the ring, acts as protective base plate. Utilizing this tool is easy, slight pressure by squeezing the pliers


handle gently and simple turns by rotating the thumb wheel (Clockwise) operates the cutting action of the blade. After a few turns, the shank will be cut open the ring to be removed.

Our ring cutter can be used by a wide range of first responders and medical personnel including doctors, paramedics, EMTs, police officers, fire rescue, or even jewelers. Made with German surgical stainless steel and can be reused after sterilization.

FINGER RING CUTTING PLIER FEATURES

- Overall Approximately 6.5\"/>
- Saw Blade Diameter: 1-3/16\"/>
- Circular Razor Sharp Micro Cutting Saw Blade.





WI-34495866

High Carbon Steel No.3 Surgery Scalpel BP Handle

\$1.50 per piece

Trade price from **\$0.50**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of High Carbon Steel No.3 Surgery Scalpel BP Handle, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001141

Knife Sharpening Steel 8" Black Handle

\$14.99 per piece

Trade price from **\$13.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Knife Sharpening Steel 8" (20.3cm) Comfortable Black Handle – Kitchen Honing Steel Rod Home & Chiefs Cutlery Blade Sharpening Stick, Carbon Steel, Nickel-Chrome Plating, Polypropylene Handle with Hanging Ring

8" Knife sharpening steel maintains the edge on knives and extends the life of your cutlery with regular use. Steel rod is made of carbon steel and plated with nickel-chrome. Ergonomic, easy-to-grip polypropylene handle with hanging hole. Wipe clean with a damp cloth.

Normally used for dull knives, Best choice for chef home use. Perfect for all types of knives, including kitchen, household, sporting and pocket knives.

KNIFE SHARPENING STEEL FEATURES

- Professional Knife Sharpener 8" (20.3cm).
- Long and Slender Honing Steel Rod.
- Lengthwise Ridges.
- Metal Sharpening Protection Patch.
- Non-Slip Comfortable Polypropylene Black Handle.
- Steel Tip Protection.
- Large Hanging Ring.
- Made With Carbon Steel and Plated With Nickel-chrome.
- Wipe Clean With Damp Cloth.

PACKAGE INCLUDE:

- 1 Pc Knife Sharpening Honing Steel Rod 8" (20.3cm) with Black Handle

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Home & Kitchen

**WI-70003842****Lessmann Knife 13cm – Small****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lessmann Knife Small 13cm Double Ended Blade with Wooden Handle Premium Quality Stainless Steel Dental Technician Laboratory Instrument.

Lessmann knife (Lessmann Knife) is a specialized dental laboratory tool used for carving, shaping, trimming, and modeling dental wax patterns for crowns, bridges, and dentures. Its common applications include sculpting wax-ups for dental work, refining contours of dentures, adding detailing to biting surfaces of model teeth, and using heat to manipulate wax.

Lessmann Knife features a double ended stainless steel blade ensuring sharpness for precise cuts with lightweight wooden handle for balanced handling and reduced hand fatigue. Lessmann wax knife designed

to withstand regular sterilization protocols.

LESSMANN WAX KNIFE FEATURES:

- Wax & Porcelain Carvers for Dental Lab Technicians.
- Overall Length 5.25" (13cm) – Large.
- Double-ended, One Knife End & Other is Spoon End.
- Smooth Edges Minimize the Risk of Hand Injury.
- Lightweight Wooden Handle for Easy Grip.
- Reusable.
- High Quality Stainless Steel Blades.

PACKAGE INCLUDE

- 1 Pc Lessmann Knife 13cm Small – Double Ended

**WI-70003841****Lessmann Knife 17cm – Large****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lessmann Knife Large 17cm Double Ended Blade with Wooden Handle Premium Quality Stainless Steel Dental Technician Laboratory Instrument.

Lessmann knife (Lessmann Knife) is a specialized dental laboratory tool used for carving, shaping, trimming, and modeling dental wax patterns for crowns, bridges, and dentures. Its common applications include sculpting wax-ups for dental work, refining contours of dentures, adding detailing to biting surfaces of model teeth, and using heat to manipulate wax.

Lessmann Knife features a double ended stainless steel blade ensuring sharpness for precise cuts with lightweight wooden handle for balanced handling and reduced hand fatigue. Lessmann wax knife designed

to withstand regular sterilization protocols.

LESSMANN WAX KNIFE FEATURES:

- Wax & Porcelain Carvers for Dental Lab Technicians.
- Overall Length 6.7" (17cm) – Large.
- Double-ended, One Knife End & Other is Spoon End.
- Smooth Edges Minimize the Risk of Hand Injury.
- Lightweight Wooden Handle for Easy Grip.
- Reusable.
- High Quality Stainless Steel Blades.

PACKAGE INCLUDE

- 1 Pc Lessmann Knife 17cm Large – Double Ended



**WI-70001080****Lucas Curette #85 Double Ended 2.7mm Spoon Shaped Blades****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Lucas Curette #85 6.5" (16.5cm) Double-Ended Spoon Shaped 2.7mm Blades with Round Pattern Hollow Handle German Grade Stainless Steel Periodontal Surgical Dental Lab Instrument.

Lucas bone curette used to remove unwanted tissue, cysts or other debris from the tooth socket. These surgical curettes are designed with 2.7mm wide spoon shaped blades that have an elongated radius. Made with premium quality stainless steel ensure maximum durability and dependability, last long and can be reused overtime.

LUCAS CURETTE #85 FEATURES

- 85 Lucas Surgical Curette Overall 6.5" (16.5cm) Length.
- 2.7mm Spoon Shaped Blade Ends.

- Double Ended.
- Non-Slip Round Shaped Hollow Handle.
- Reusable.
- For Cyst Removal, Curettage, and Tooth Socket Debridement.
- Made with Premium Quality Stainless Steel.

PACKAGE INCLUDE

- 1 Pc Lucas Bone Curette #85 Double Ended Hollow Handle

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Dental Lab Instrument**APPLICATION:** Curettage, Cyst Removal Tooth Socket Debridement.**SHAPE:** Double Ended**AUTOCCLAVABLE/REUSABLE:** Yes**WI-70001975****Mac Blade #0 with Conventional Laryngoscope Medium Handle****\$20.99** per pieceTrade price from **\$19.31**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Mac Blade #0 with Conventional Laryngoscope Medium Handle Set Premium Quality EMT Anesthesia Intubation Surgical Instruments.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

2 PCS LARYNGOSCOPE MAC SET:

MAC BLADE #0: Macintosh conventional blade #0 measuring approximately 80mm, designed for

neonates and premature newborns. It helping to visualize the vocal cords and place a breathing tube with difficult airways. Its curved shape fits into the vallecula to indirectly lift the epiglottis. The blade features conventional lighting provided by a standard bulb and is supplied with a backup blade and handle to ensure sterility and immediate availability in emergencies.

CONVENTIONAL LARYNGOSCOPE

HANDLE Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

PACKAGE INCLUDE:

MACINTOSH CONVENTIONAL LARYNGOSCOPE BLADE #0



**WI-70001974****Mac Blade #1 with Conventional Laryngoscope Medium Handle****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Mac Blade #1 with Conventional Laryngoscope Medium Handle Set Premium Quality EMT Anesthesia Intubation Surgical Instruments.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

2 PCS LARYNGOSCOPE MAC SET:

MAC BLADE #1: Macintosh conventional blade #1 measuring approximately 94mm, designed for

infants. It helping to visualize the vocal cords and place a breathing tube with difficult airways. Its curved shape fits into the vallecula to indirectly lift the epiglottis. The blade features conventional lighting provided by a standard bulb and is supplied with a backup blade and handle to ensure sterility and immediate availability in emergencies.

CONVENTIONAL LARYNGOSCOPE HANDLE

Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

PACKAGE INCLUDE:**WI-70001973****Mac Blade #2 with Conventional Laryngoscope Medium Handle****\$23.99** per pieceTrade price from **\$22.07**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Mac Blade #2 with Conventional Laryngoscope Medium Handle Set Premium Quality EMT Anesthesia Intubation Surgical Instruments.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

2 PCS LARYNGOSCOPE MAC SET:

MAC BLADE #2: Macintosh conventional blade #2 measuring approximately 114mm, designed for

pediatric patients, older children, or small-stature adults. It helping to visualize the vocal cords and place a breathing tube with difficult airways. Its curved shape fits into the vallecula to indirectly lift the epiglottis. The blade features conventional lighting provided by a standard bulb and is supplied with a backup blade and handle to ensure sterility and immediate availability in emergencies.

CONVENTIONAL LARYNGOSCOPE HANDLE

Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

PACKAGE INCLUDE:

**WI-70001972****Mac Blade #4 with Conventional Laryngoscope Medium Handle****\$25.99** per pieceTrade price from **\$23.91**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Mac Blade #4 with Conventional Laryngoscope Medium Handle Set Premium Quality EMT Anesthesia Intubation Surgical Instruments.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

2 PCS LARYNGOSCOPE MAC SET:

MAC BLADE #4: Macintosh conventional blade #4 is an large blade, measuring approximately 155mm, is

designed for large, tall individuals, or patients with larger anatomical features, helping to visualize the vocal cords and place a breathing tube with difficult airways. Its curved shape fits into the vallecula to indirectly lift the epiglottis. The blade features conventional lighting provided by a standard bulb and is supplied with a backup blade and handle to ensure sterility and immediate availability in emergencies.

CONVENTIONAL LARYNGOSCOPE HANDLE: Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

**WI-70000031****Mac Blade #5 with Conventional Laryngoscope Medium Handle****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Mac Blade #5 Large Adult with Conventional Laryngoscope Medium Handle Set Premium Quality EMT Anesthesia Intubation Surgical Instruments.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

2 PCS LARYNGOSCOPE MAC SET:

MAC BLADE #5: Macintosh conventional blade #5 is an extra-large blade, measuring approximately 179

mm, is designed for large adults or those with difficult airways. Its curved shape fits into the vallecula to indirectly lift the epiglottis. The blade features conventional lighting provided by a standard bulb and is supplied with a backup blade and handle to ensure sterility and immediate availability in emergencies.

CONVENTIONAL LARYNGOSCOPE HANDLE: Medium size conventional laryngoscope handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

PACKAGE INCLUDE:

**WI-70001343****Macintosh Blade #0 with Fiber Optic Laryngoscope Medium Handle****\$28.99** per pieceTrade price from **\$26.67**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Macintosh Blade #0 with Fiber Optic Laryngoscope Medium Handle Set EMT Anesthesia Intubation Surgical Instruments.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

2 PCS MAC LARYNGOSCOPE SET:

MACINTOSH BLADES: Macintosh laryngoscope blades with a whiter & brighter Fiberoptic light source provides a high level of illumination

allowing instant illumination of the oral cavity. Removable light pipe facilitates cleaning and service. The green color coded to signify compatibility only with fiber optic handle.

CONVENTIONAL LARYNGOSCOPE

HANDLE: Fiber Optic Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.
- Blades Compatible with All Fiber Optic Laryngoscope Systems..

**WI-70001336****Macintosh Blade #1 with Fiber Optic Laryngoscope Medium Handle****\$28.99** per pieceTrade price from **\$26.67**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Macintosh Blade #1 with Fiber Optic Laryngoscope Medium Handle Set EMT Anesthesia Intubation Surgical Instruments.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

2 PCS MAC LARYNGOSCOPE SET:

MACINTOSH BLADES: Macintosh laryngoscope blades with a whiter & brighter Fiberoptic light source provides a high level of illumination

allowing instant illumination of the oral cavity. Removable light pipe facilitates cleaning and service. The green color coded to signify compatibility only with fiber optic handle.

CONVENTIONAL LARYNGOSCOPE

HANDLE: Fiber Optic Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.
- Blades Compatible with All Fiber Optic Laryngoscope Systems..



**WI-70001335****Macintosh Blade #2 with Fiber Optic Laryngoscope Medium Handle****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Macintosh Blade #2 with Fiber Optic Laryngoscope Medium Handle Set EMT Anesthesia Intubation Surgical Instruments.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

2 PCS MAC LARYNGOSCOPE SET:

MACINTOSH BLADES: Macintosh laryngoscope blades with a whiter & brighter Fiberoptic light source provides a high level of illumination

allowing instant illumination of the oral cavity. Removable light pipe facilitates cleaning and service. The green color coded to signify compatibility only with fiber optic handle.

CONVENTIONAL LARYNGOSCOPE

HANDLE: Fiber Optic Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.
- Blades Compatible with All Fiber Optic Laryngoscope Systems..

**WI-70001334****Macintosh Blade #3 with Fiber Optic Laryngoscope Medium Handle****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Macintosh Blade #3 with Fiber Optic Laryngoscope Medium Handle Set EMT Anesthesia Intubation Surgical Instruments.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

2 PCS MAC LARYNGOSCOPE SET:

MACINTOSH BLADES: Macintosh laryngoscope blades with a whiter & brighter Fiberoptic light source provides a high level of illumination

allowing instant illumination of the oral cavity. Removable light pipe facilitates cleaning and service. The green color coded to signify compatibility only with fiber optic handle.

CONVENTIONAL LARYNGOSCOPE

HANDLE: Fiber Optic Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.
- Blades Compatible with All Fiber Optic Laryngoscope Systems..



**WI-70001333****Macintosh Blade #4 with Fiber Optic Laryngoscope Medium Handle****\$29.99** per pieceTrade price from **\$27.59**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Macintosh Blade #4 with Fiber Optic Laryngoscope Medium Handle Set EMT Anesthesia Intubation Surgical Instruments.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

2 PCS MAC LARYNGOSCOPE SET:

MACINTOSH BLADES: Macintosh laryngoscope blades with a whiter & brighter Fiberoptic light source provides a high level of illumination

allowing instant illumination of the oral cavity. Removable light pipe facilitates cleaning and service. The green color coded to signify compatibility only with fiber optic handle.

CONVENTIONAL LARYNGOSCOPE HANDLE: Fiber Optic Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.
- Blades Compatible with All Fiber Optic Laryngoscope Systems..

**WI-70001980****Macintosh Conventional Laryngoscope Blade #0****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Macintosh Conventional Blade #0 EMT Anesthesia Intubation Surgical Instruments

Macintosh laryngoscope blades with a whiter & brighter light source provides a high level of illumination allowing instant illumination of the oral cavity.

The light bulb in the handle can be easily change while the distance from bulb to tip protects keeps patient safe from the accidental burns.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to

facilitate intubation in unparalyzed patients.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.
- Blades Compatible with All Conventional Laryngoscope Systems..
- High Quality Stainless Steel.
- Conventional Light Source.
- Satin Finish.
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula

Replacement Bulbs:

- Small LED Replacement Bulb (for MacIntosh size 1 & 2)



**WI-70001979**

Macintosh Conventional Laryngoscope Blade #1

\$11.99 per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Macintosh Conventional Blade #1 EMT Anesthesia Intubation Surgical Instruments

Macintosh laryngoscope blades with a whiter & brighter light source provides a high level of illumination allowing instant illumination of the oral cavity. The light bulb in the handle can be easily change while the distance from bulb to tip protects keeps patient safe from the accidental burns.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to

facilitate intubation in unparalyzed patients.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.
- Blades Compatible with All Conventional Laryngoscope Systems..
- High Quality Stainless Steel.
- Conventional Light Source.
- Satin Finish.
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula

Replacement Bulbs:

- Small LED Replacement Bulb (for MacIntosh size 1 & 2)

**WI-70001978**

Macintosh Conventional Laryngoscope Blade #2

\$12.99 per pieceTrade price from **\$11.95**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Macintosh Conventional Blade #2 EMT Anesthesia Intubation Surgical Instruments

Macintosh laryngoscope blades with a whiter & brighter light source provides a high level of illumination allowing instant illumination of the oral cavity. The light bulb in the handle can be easily change while the distance from bulb to tip protects keeps patient safe from the accidental burns.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to

facilitate intubation in unparalyzed patients.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.
- Blades Compatible with All Conventional Laryngoscope Systems..
- High Quality Stainless Steel.
- Conventional Light Source.
- Satin Finish.
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula

Replacement Bulbs:

- Small LED Replacement Bulb (for MacIntosh size 1 & 2)



**WI-70001977****Macintosh Conventional Laryngoscope Blade #4****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Macintosh Conventional Blade #4 EMT Anesthesia Intubation Surgical Instruments

Macintosh laryngoscope blades with a whiter & brighter light source provides a high level of illumination allowing instant illumination of the oral cavity. The light bulb in the handle can be easily change while the distance from bulb to tip protects keeps patient safe from the accidental burns.

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to

facilitate intubation in unparalyzed patients.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.
- Blades Compatible with All Conventional Laryngoscope Systems..
- High Quality Stainless Steel.
- Conventional Light Source.
- Satin Finish.
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula

Replacement Bulbs:

- Small LED Replacement Bulb (for MacIntosh size 1 & 2)

**WI-70000033****Macintosh Conventional Laryngoscope Blade #5****\$14.99** per pieceTrade price from **\$13.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Macintosh Conventional Blade #5 Large Adult 7" (179mm) EMT Anesthesia Endotracheal Intubation Surgical Instruments.

Macintosh laryngoscope blades with a whiter & brighter light source provides a high level of illumination allowing instant illumination of the oral cavity. The light bulb in the handle can be easily change while the distance from bulb to tip protects keeps patient safe from the accidental burns.


Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to

facilitate intubation in unparalyzed patients. Mac conventional blade #5 is specifically designed for large adult airways, measuring approximately 179 mm to accommodate individuals with larger anatomies or difficult oversized airways.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.
- Blades Compatible with All Conventional Laryngoscope Systems..
- High Quality Stainless Steel.
- Conventional Light Source.
- Satin Finish.
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula



**WI-70001350**

Macintosh Fiber Optic Laryngoscope Blade #0 Anesthesia Intubation

\$15.99 per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Macintosh Fiber Optic Blade #0
Premie Reusable, Autoclavable
German Grade Stainless Steel
Laryngoscope Blades EMT Anesthesia
Intubation Surgical Instruments

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.

- Blades Compatible with All Fiber Optic Laryngoscope Systems..
- High Quality Stainless Steel.
- Fiberoptic light transmission provides cool, obstruction free viewing
- Reusable Blade
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula
- Removable light pipe facilitates cleaning and service.

PACKAGE INCLUDES:

- Fiber Optic Laryngoscope Blade – MacIntosh; Size 0, Premie; Overall Length: 80 mm

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel**BLADES:** Curved Blade**REUSABLE:** Yes**WI-70001349**

Macintosh Fiber Optic Laryngoscope Blade #1 Anesthesia Intubation

\$15.99 per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Macintosh Fiber Optic Blade #1 Infant
Size Reusable, Autoclavable German
Grade Stainless Steel Laryngoscope
Blades EMT Anesthesia Intubation
Surgical Instruments

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.

- Blades Compatible with All Fiber Optic Laryngoscope Systems..
- High Quality Stainless Steel.
- Fiberoptic light transmission provides cool, obstruction free viewing
- Reusable Blade
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula
- Removable light pipe facilitates cleaning and service.

PACKAGE INCLUDES:

- Fiber Optic Laryngoscope Blade – MacIntosh; Size 1, Infant; Overall Length: 94 mm

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel**BLADES:** Curved Blade**REUSABLE:** Yes

**WI-70001348****Macintosh Fiber Optic Laryngoscope Blade #2
Anesthesia Intubation****\$16.99** per pieceTrade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Macintosh Fiber Optic Blade #2 Child Size Reusable, Autoclavable German Grade Stainless Steel Laryngoscope Blades EMT Anesthesia Intubation Surgical Instruments

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.

- Blades Compatible with All Fiber Optic Laryngoscope Systems..
- High Quality Stainless Steel.
- Fiberoptic light transmission provides cool, obstruction free viewing
- Reusable Blade
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula
- Removable light pipe facilitates cleaning and service.

PACKAGE INCLUDES:

- Fiber Optic Laryngoscope Blade – MacIntosh; Size 2, Child; Overall Length: 114 mm

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**BLADES:** Curved Blade**REUSABLE:** YesMACINTOSH LARYNGOSCOPE
FIBER OPTIC BLADE # 2**WI-70001347****Macintosh Fiber Optic Laryngoscope Blade #3
Anesthesia Intubation****\$16.99** per pieceTrade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Macintosh Fiber Optic Blade #3 Medium Adult Size Reusable, Autoclavable German Grade Stainless Steel Laryngoscope Blades EMT Anesthesia Intubation Surgical Instruments

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.

- Flange Reduced for Less Trauma.
- Blades Compatible with All Fiber Optic Laryngoscope Systems..
- High Quality Stainless Steel.
- Fiberoptic light transmission provides cool, obstruction free viewing
- Reusable Blade
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula
- Removable light pipe facilitates cleaning and service.

PACKAGE INCLUDES:

- Fiber Optic Laryngoscope Blade – MacIntosh; Size 3, Medium Adult; Overall Length: 133 mm

TECHNICAL INFORMATION:**MATERIAL:** Stainless SteelMACINTOSH LARYNGOSCOPE
FIBER OPTIC BLADE # 3

**WI-70001346****Macintosh Fiber Optic Laryngoscope Blade #4
Anesthesia Intubation****\$16.99** per pieceTrade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Macintosh Fiber Optic Blade #4 Large Adult Size Reusable, Autoclavable German Grade Stainless Steel Laryngoscope Blades EMT Anesthesia Intubation Surgical Instruments

Mac laryngoscope has a curved blade which allows exposure of the larynx by positioning the tip in the vallecula, anterior to the epiglottis, lifting it out of view. The Mac laryngoscope is probably the most successful durable in the history of anesthesia to date. Macintosh blade originally designed to facilitate intubation in unparalyzed patients.

MACINTOSH BLADES FEATURES:

- Curved Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Flange Reduced for Less Trauma.

- Blades Compatible with All Fiber Optic Laryngoscope Systems..
- High Quality Stainless Steel.
- Fiberoptic light transmission provides cool, obstruction free viewing
- Reusable Blade
- Fully Autoclavable Blade.
- Smooth and Gentle Curved Spatula
- Removable light pipe facilitates cleaning and service.

PACKAGE INCLUDES:

- Fiber Optic Laryngoscope Blade – Macintosh; Size 4, Large Adult; Overall Length: 155 mm

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**BLADES:** Curved Blade**REUSABLE:** Yes

MACINTOSH LARYNGOSCOPE FIBER OPTIC BLADE # 4

**WI-70001332****McCoy Fiber Optic Flex Tip Laryngoscope
Blade #2****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

McCoy Laryngoscope Flexi-Tip Fiber-Optic Led Blade #2 for Difficult Intubation German Grade Stainless Steel Anesthesiology Medical Instrument

Fibre Optic Laryngoscope Blades (McCoy Flex Tip Profile) feature a low shoulder profile for better visualization and easier intubation of patients. McCoy Fiber Optic Flex Tip blade have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Moreover, the McCoy blade slimmer than the Macintosh blade from the tip and handle sides.

These are reusable blades, extremely useful for better visualization of the vocal cords during difficult intubations. McCoy laryngoscope is used in

patients when difficulty in elevating epiglottis is encountered due to limited mouth opening and activation of blade tip elevates the epiglottis and visualization of vocal cords.

**MCCOY FIBER OPTIC FLEX TIP
BLADE FEATURES**

- McCoy Flex Tip Blade #2
- Curved Blade With Low Shoulder Profile for Better Visualization
- McCoy Flex Tip Blades are Slimer Than Macintosh Curved Blade From Tip and Handle Sides.
- Hinged Tip Allows Elevation of the Distal Structures (Epiglottis).
- Extremely Usfull for Difficult Intubation When Difficulty in Elevating Epiglottis.



MCCOY FIBEROPTIC LARYNGOSCOPE BLADE # 2



**WI-70001331**

McCoy Fiber Optic Flex Tip Laryngoscope Blade #3

\$19.99 per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

McCoy Laryngoscope Flexi-Tip Fiber-Optic Led Blade #3 for Difficult Intubation German Grade Stainless Steel Anesthesiology Medical Instrument

Fibre Optic Laryngoscope Blades (McCoy Flex Tip Profile) feature a low shoulder profile for better visualization and easier intubation of patients. McCoy Fiber Optic Flex Tip blade have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Moreover, the McCoy blade slimer than the Macintosh blade from the tip and handle sides.

These are reusable blades, extremely useful for better visualization of the vocal cords during difficult intubations. McCoy laryngoscope is used in

patients when difficulty in elevating epiglottis is encountered due to limited mouth opening and activation of blade tip elevates the epiglottis and visualization of vocal cords.

MCCOY FIBER OPTIC FLEX TIP BLADE FEATURES

- McCoy Flex Tip Blade #3
- Curved Blade With Low Shoulder Profile for Better Visualization
- Mccoy Flex Tip Blades are Slimer Than Macintosh Curved Blade From Tip and Handle Sides.
- Hinged Tip Allows Elevation of the Distal Structures (Epiglottis).
- Extremely Usfull for Difficult Intubation When Difficulty in Elevating Epiglottis.

WI-70001330

McCoy Fiber Optic Flex Tip Laryngoscope Blade #4

\$19.99 per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

McCoy Laryngoscope Flexi-Tip Fiber-Optic Led Blade #4 for Difficult Intubation German Grade Stainless Steel Anesthesiology Medical Instrument

Fibre Optic Laryngoscope Blades (McCoy Flex Tip Profile) feature a low shoulder profile for better visualization and easier intubation of patients. McCoy Fiber Optic Flex Tip blade have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Moreover, the McCoy blade slimer than the Macintosh blade from the tip and handle sides.

These are reusable blades, extremely useful for better visualization of the vocal cords during difficult intubations. McCoy laryngoscope is used in

patients when difficulty in elevating epiglottis is encountered due to limited mouth opening and activation of blade tip elevates the epiglottis and visualization of vocal cords.

MCCOY FIBER OPTIC FLEX TIP BLADE FEATURES

- McCoy Flex Tip Blade #4
- Curved Blade With Low Shoulder Profile for Better Visualization
- Mccoy Flex Tip Blades are Slimer Than Macintosh Curved Blade From Tip and Handle Sides.
- Hinged Tip Allows Elevation of the Distal Structures (Epiglottis).
- Extremely Usfull for Difficult Intubation When Difficulty in Elevating Epiglottis.



**WI-70001328****McCoy Laryngoscope Blade #2 with Fiber Optic Medium Handle****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pc Fiber Optic McCoy Laryngoscope Set with Flex Tip Blade #2 and Medium Handle for Difficult Intubation German Grade Stainless Steel Anesthesiology Medical Instrument

Fibre Optic Laryngoscope Blades (McCoy Flex Tip Profile) feature a low shoulder profile for better visualization and easier intubation of patients. McCoy Fiber Optic Flex Tip blade have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Moreover, the McCoy blade slimer than the Macintosh blade from the tip and handle sides.

These are reusable blades, extremely useful for better visualization of the vocal cords during difficult intubations. McCoy laryngoscope is used in

patients when difficulty in elevating epiglottis is encountered due to limited mouth opening and activation of blade tip elevates the epiglottis and visualization of vocal cords.

2 PCS MCCOY LARYNGOSCOPE SET:

MCCOY FLEX TIP BLADES: McCoy Flex Tip blades have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Extremely useful for better visualization of the vocal cords during difficult intubation.

CONVENTIONAL LARYNGOSCOPE HANDLE: Fiber Optic Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday

WI-70001327**McCoy Laryngoscope Blade #3 with Fiber Optic Medium Handle****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pc Fiber Optic McCoy Laryngoscope Set with Flex Tip Blade #3 and Medium Handle for Difficult Intubation German Grade Stainless Steel Anesthesiology Medical Instrument

Fibre Optic Laryngoscope Blades (McCoy Flex Tip Profile) feature a low shoulder profile for better visualization and easier intubation of patients. McCoy Fiber Optic Flex Tip blade have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Moreover, the McCoy blade slimer than the Macintosh blade from the tip and handle sides.

These are reusable blades, extremely useful for better visualization of the vocal cords during difficult intubations. McCoy laryngoscope is used in

patients when difficulty in elevating epiglottis is encountered due to limited mouth opening and activation of blade tip elevates the epiglottis and visualization of vocal cords.

2 PCS MCCOY LARYNGOSCOPE SET:

MCCOY FLEX TIP BLADES: McCoy Flex Tip blades have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Extremely useful for better visualization of the vocal cords during difficult intubation.

CONVENTIONAL LARYNGOSCOPE HANDLE: Fiber Optic Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday





WI-70001326

McCoy Laryngoscope Blade #4 with Fiber Optic Medium Handle

\$26.99 per piece

Trade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pc Fiber Optic McCoy Laryngoscope Set with Flex Tip Blade #4 and Medium Handle for Difficult Intubation German Grade Stainless Steel Anesthesiology Medical Instrument

Fibre Optic Laryngoscope Blades (McCoy Flex Tip Profile) feature a low shoulder profile for better visualization and easier intubation of patients. McCoy Fiber Optic Flex Tip blade have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Moreover, the McCoy blade slimer than the Macintosh blade from the tip and handle sides.

These are reusable blades, extremely useful for better visualization of the vocal cords during difficult intubations. McCoy laryngoscope is used in

patients when difficulty in elevating epiglottis is encountered due to limited mouth opening and activation of blade tip elevates the epiglottis and visualization of vocal cords.

2 PCS MCCOY LARYNGOSCOPE SET:

MCCOY FLEX TIP BLADES: McCoy Flex Tip blades have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Extremely useful for better visualization of the vocal cords during difficult intubation.

CONVENTIONAL LARYNGOSCOPE HANDLE: Fiber Optic Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday



WI-09286277

Medical Sterile Stainless Steel Surgical Scalpel Blades

\$2.00 per piece

Trade price from **\$0.80**

Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Medical Sterile Stainless Steel Surgical Scalpel Blades, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$2.00, falling to \$0.80 at 10 pieces — about 60% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001322****Micro Finger Ring Cutter 4.75" with Saw Blade****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Micro Finger Ring Cutter 4.75 Inches with Saw Blade & Non-Slip Grip Handle Jewelry Emergency Rescue First Aid Paramedic EMT EMS German Grade Stainless Steel Surgical Instrument.

Finger ring cutter is widely used to cut off and remove rings that won't come off from fingers. This instrument have twisting thumb screws that allow the rotary blade to cut the circle while keeping the finger from the edge. This is a ring cutter that saws (utilizes a circular saw blade 7/8" OD) ring while still on or off the finger. Ring cutter is styled like the french ring cutting tool, its style permits easy access and a firm grip on the ring.

Finger ring cutter have a feature one black PVC rubber grip handle and the other is non-dlip knurled handle. The

micro finger ring cutter is versatile and durable but small enough to fit into most rescue holsters without being bulky and awkward. Made with German surgical stainless steel and can be reused after sterilization.

MICRO FINGER RING CUTTER FEATURES

- Circular Razor Sharp Micro Cutting Saw Blade 7/8" (22mm) OD.
- 1 Non-slip Rubberized Grip Handle and 1 Knurled Handle.
- Easy to Use on Rings on or Off the Finger.
- Twisting Thumb Screw to Operates the Razor Sharp Serrated Saw.
- Safety Lever Slides Protective Plate Between Finger and Ring.

**WI-70001364****Miller Blade #0 with Conventional Laryngoscope Medium Handle****\$20.99** per pieceTrade price from **\$19.31**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Miller Blade #0 with Conventional Laryngoscope Medium Handle Set EMT Anesthesia Intubation Surgical Instruments.

The Miller laryngoscope is a straight blade designed to obtain a view of the vocal cords by directly lifting the epiglottis. It has useful application in 'floppy' airways making it popular within paediatric anaesthesia. This set allows the physicians to treat the patient without removing them from the MRI suite.

MILLER LARYNGOSCOPE SET:

MILLER BLADES: The Miller laryngoscope has a straight blade with a long, curved tip, providing better exposure of the larynx. The unique blade design also makes it ideal for a

range of medical approaches. It comes with different variations such as Miller Laryngoscope Blade, Size 0, 1, 2, 3, and 4 Conventional Style.

CONVENTIONAL LARYNGOSCOPE HANDLE: Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.



**WI-70001363****Miller Blade #1 with Conventional Laryngoscope Medium Handle****\$20.99** per pieceTrade price from **\$19.31**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Miller Blade #1 with Conventional Laryngoscope Medium Handle Set EMT Anesthesia Intubation Surgical Instruments.

The Miller laryngoscope is a straight blade designed to obtain a view of the vocal cords by directly lifting the epiglottis. It has useful application in 'floppy' airways making it popular within paediatric anaesthesia. This set allows the physicians to treat the patient without removing them from the MRI suite.

MILLER LARYNGOSCOPE SET:

MILLER BLADES: The Miller laryngoscope has a straight blade with a long, curved tip, providing better exposure of the larynx. The unique blade design also makes it ideal for a

range of medical approaches. It comes with different variations such as Miller Laryngoscope Blade, Size 0, 1, 2, 3, and 4 Conventional Style.

CONVENTIONAL LARYNGOSCOPE

HANDLE: Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.

**WI-70001362****Miller Blade #2 with Conventional Laryngoscope Medium Handle****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Miller Blade #2 with Conventional Laryngoscope Medium Handle Set EMT Anesthesia Intubation Surgical Instruments.

The Miller laryngoscope is a straight blade designed to obtain a view of the vocal cords by directly lifting the epiglottis. It has useful application in 'floppy' airways making it popular within paediatric anaesthesia. This set allows the physicians to treat the patient without removing them from the MRI suite.

MILLER LARYNGOSCOPE SET:

MILLER BLADES: The Miller laryngoscope has a straight blade with a long, curved tip, providing better exposure of the larynx. The unique blade design also makes it ideal for a

range of medical approaches. It comes with different variations such as Miller Laryngoscope Blade, Size 0, 1, 2, 3, and 4 Conventional Style.

CONVENTIONAL LARYNGOSCOPE

HANDLE: Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.



**WI-70001361****Miller Blade #3 with Conventional Laryngoscope Medium Handle****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Miller Blade #3 with Conventional Laryngoscope Medium Handle Set EMT Anesthesia Intubation Surgical Instruments.

The Miller laryngoscope is a straight blade designed to obtain a view of the vocal cords by directly lifting the epiglottis. It has useful application in 'floppy' airways making it popular within paediatric anaesthesia. This set allows the physicians to treat the patient without removing them from the MRI suite.

MILLER LARYNGOSCOPE SET:

MILLER BLADES: The Miller laryngoscope has a straight blade with a long, curved tip, providing better exposure of the larynx. The unique blade design also makes it ideal for a

range of medical approaches. It comes with different variations such as Miller Laryngoscope Blade, Size 0, 1, 2, 3, and 4 Conventional Style.

CONVENTIONAL LARYNGOSCOPE

HANDLE: Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.

**WI-70001360****Miller Blade #4 with Conventional Laryngoscope Medium Handle****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

2 Pcs Standard Miller Blade #4 with Conventional Laryngoscope Medium Handle Set EMT Anesthesia Intubation Surgical Instruments.

The Miller laryngoscope is a straight blade designed to obtain a view of the vocal cords by directly lifting the epiglottis. It has useful application in 'floppy' airways making it popular within paediatric anaesthesia. This set allows the physicians to treat the patient without removing them from the MRI suite.

MILLER LARYNGOSCOPE SET:

MILLER BLADES: The Miller laryngoscope has a straight blade with a long, curved tip, providing better exposure of the larynx. The unique blade design also makes it ideal for a

range of medical approaches. It comes with different variations such as Miller Laryngoscope Blade, Size 0, 1, 2, 3, and 4 Conventional Style.

CONVENTIONAL LARYNGOSCOPE

HANDLE: Conventional Laryngoscope Handle are made of Stainless Steel material with a knurled grip to insure strong grip. Ideal for crash carts, emergency situations and everyday O.R. intubation kits. These are reusable laryngoscope handle build to last long.

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.



**WI-70001369****Miller Fiber Optic Laryngoscope Blade #0****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Reusable Fiber Optic Miller Blade #0 Neonate German Grade Stainless Steel EMT Anesthesia Intubation Surgical Instruments. This blade is ideal for neonatal intubations

Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube. Miller Blades originally developed by Dr. Robert Miller for pediatric use, their matte finish makes for easy clean up and use and are available in a range of sizes.

The Fiber Optic Laryngoscope Miller Blades are the most popular type of the straight blade. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater

exposure of the larynx in difficult-to-intubate patients provides a better view of the vocal cords. The straight Miller laryngoscope blade is traditionally recommended for intubation in infants, due to the large size and flexibility of the infant epiglottis

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.
- Blades Compatible with All Conventional Laryngoscope Systems.
- High Quality German Grade Stainless Steel.



MILLER FIBEROPTIC LARYNGOSCOPE BLADE # 0

**WI-70001368****Miller Fiber Optic Laryngoscope Blade #1****\$15.99** per pieceTrade price from **\$14.71**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Reusable Fiber Optic Miller Blade #1 Infant German Grade Stainless Steel EMT Anesthesia Intubation Surgical Instruments.

Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube. Miller Blades originally developed by Dr. Robert Miller for pediatric use, their matte finish makes for easy clean up and use and are available in a range of sizes.

The Fiber Optic Laryngoscope Miller Blades are the most popular type of the straight blade. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-

intubate patients provides a better view of the vocal cords. The straight Miller laryngoscope blade is traditionally recommended for intubation in infants, due to the large size and flexibility of the infant epiglottis

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.
- Blades Compatible with All Conventional Laryngoscope Systems.
- High Quality German Grade Stainless Steel.



MILLER FIBEROPTIC LARYNGOSCOPE BLADE # 1



**WI-70001367****Miller Fiber Optic Laryngoscope Blade #2****\$16.99** per pieceTrade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Reusable Fiber Optic Miller Blade #2 Child German Grade Stainless Steel EMT Anesthesia Intubation Surgical Instruments.

Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube. Miller Blades originally developed by Dr. Robert Miller for pediatric use, their matte finish makes for easy clean up and use and are available in a range of sizes.

The Fiber Optic Laryngoscope Miller Blades are the most popular type of the straight blade. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-

intubate patients provides a better view of the vocal cords. The straight Miller laryngoscope blade is traditionally recommended for intubation in infants, due to the large size and flexibility of the infant epiglottis

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.
- Blades Compatible with All Conventional Laryngoscope Systems.
- High Quality German Grade Stainless Steel.

MILLER FIBEROPTIC LARYNGOSCOPE BLADE # 2

**WI-70001366****Miller Fiber Optic Laryngoscope Blade #3****\$16.99** per pieceTrade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Reusable Fiber Optic Miller Blade #3 Adult German Grade Stainless Steel EMT Anesthesia Intubation Surgical Instruments.

Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube. Miller Blades originally developed by Dr. Robert Miller for pediatric use, their matte finish makes for easy clean up and use and are available in a range of sizes.

The Fiber Optic Laryngoscope Miller Blades are the most popular type of the straight blade. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-

intubate patients provides a better view of the vocal cords. The straight Miller laryngoscope blade is traditionally recommended for intubation in infants, due to the large size and flexibility of the infant epiglottis

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.
- Blades Compatible with All Conventional Laryngoscope Systems.
- High Quality German Grade Stainless Steel.

MILLER FIBEROPTIC LARYNGOSCOPE BLADE # 3



**WI-70001365****Miller Fiber Optic Laryngoscope Blade #4****\$16.99** per pieceTrade price from **\$15.63**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Reusable Fiber Optic Miller Blade #4 Large Adult German Grade Stainless Steel EMT Anesthesia Intubation Surgical Instruments.

Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube. Miller Blades originally developed by Dr. Robert Miller for pediatric use, their matte finish makes for easy clean up and use and are available in a range of sizes.

The Fiber Optic Laryngoscope Miller Blades are the most popular type of the straight blade. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-

intubate patients provides a better view of the vocal cords. The straight Miller laryngoscope blade is traditionally recommended for intubation in infants, due to the large size and flexibility of the infant epiglottis

MILLER LARYNGOSCOPE BLADE FEATURES:

- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.
- Blades Compatible with All Conventional Laryngoscope Systems.
- High Quality German Grade Stainless Steel.

WI-70001971**Miller Laryngoscope Blade #0 – Conventional****\$11.99** per pieceTrade price from **\$11.03**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Reusable Miller Conventional Blade #0 Neonate 80mm Long EMT Anesthesia Intubation Surgical Instruments.

Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube. Miller Blades originally developed by Dr. Robert Miller for pediatric use, their matte finish makes for easy clean up and use and are available in a range of sizes.

Conventional Laryngoscope Miller Blades are the most popular type of the straight blade. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-

intubate patients provides a better view of the vocal cords. The straight Miller laryngoscope blade is traditionally recommended for intubation in infants, due to the large size and flexibility of the infant epiglottis

MILLER LARYNGOSCOPE BLADE FEATURES:

- Miller Laryngoscope Blade Size #0 – Neonate
- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.
- Blades Compatible with All Conventional Laryngoscope Systems.
- High Quality Stainless Steel.

**WI-70001970****Miller Laryngoscope Blade #1 Conventional****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Reusable Miller Conventional Laryngoscope Blade #1 Infant 105mm EMT Anesthesia Intubation Surgical Instruments.

Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube. Miller Blades originally developed by Dr. Robert Miller for pediatric use, their matte finish makes for easy clean up and use and are available in a range of sizes.

Conventional Laryngoscope Miller Blades are the most popular type of the straight blade. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-

intubate patients provides a better view of the vocal cords. The straight Miller laryngoscope blade is traditionally recommended for intubation in infants, due to the large size and flexibility of the infant epiglottis

MILLER LARYNGOSCOPE BLADES FEATURES:

- Miller Laryngoscope Blade Size #1 – Infant
- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.
- Blades Compatible with All Conventional Laryngoscope Systems.
- High Quality Stainless Steel.

**WI-70001969****Miller Laryngoscope Blade #2 Conventional****\$12.99** per pieceTrade price from **\$11.95**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Reusable Miller Conventional Conventional Blade #2 Child 155mm EMT Anesthesia Intubation Surgical Instruments.

Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube. Miller Blades originally developed by Dr. Robert Miller for pediatric use, their matte finish makes for easy clean up and use and are available in a range of sizes.

Conventional Laryngoscope Miller Blades are the most popular type of the straight blade. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-intubate patients provides a better view of the vocal cords. The straight Miller laryngoscope blade is traditionally recommended for intubation in infants, due to the large size and flexibility of the infant epiglottis

**MILLER LARYNGOSCOPE BLADE FEATURES:**

- Miller Laryngoscope Blade Size #2 – Child
- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the

**WI-70001968****Miller Laryngoscope Blade #3 Conventional****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Reusable Miller Conventional Laryngoscope Blade #3 Adult 195mm Long EMT Anesthesia Intubation Surgical Instruments.

Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube. Miller Blades originally developed by Dr. Robert Miller for pediatric use, their matte finish makes for easy clean up and use and are available in a range of sizes.

Conventional Laryngoscope Miller Blades are the most popular type of the straight blade. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-

intubate patients provides a better view of the vocal cords. The straight Miller laryngoscope blade is traditionally recommended for intubation in infants, due to the large size and flexibility of the infant epiglottis

MILLER LARYNGOSCOPE BLADES FEATURES:

- Miller Laryngoscope Blade Size #3 – Adult
- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.
- Blades Compatible with All Conventional Laryngoscope Systems.
- High Quality Stainless Steel.

WI-70001967**Miller Laryngoscope Blade #4 Conventional****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Standard Reusable Miller Conventional Laryngoscope Blade #4 Adult Large 205mm Long EMT Anesthesia Intubation Surgical Instruments

Laryngoscope blades are used as the primary tool for examination of the interior of the larynx and for placement of an endotracheal tube. Miller Laryngoscope Blades originally developed by Dr. Robert Miller for pediatric use, their matte finish makes for easy clean up and use and are available in a range of sizes.

Conventional Laryngoscope Miller Blades are the most popular type of the straight blade. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater

exposure of the larynx in difficult-to-intubate patients provides a better view of the vocal cords. Straight Miller laryngoscope blade is traditionally recommended for intubation in infants, due to the large size and flexibility of the infant epiglottis

MILLER LARYNGOSCOPE BLADES FEATURES:

- Miller Laryngoscope Blade Size #4 – Large Adult
- Straight Blades with Slight Curve at the Beak or Tip Extended to Improve the Lifting of the Epiglottis.
- Round Edges Prevent From Injury Risk.
- Blades Compatible with All Conventional Laryngoscope Systems.



WI-34516917

No. 4 Surgical Scalpel Handle

\$1.50 per piece

Trade price from \$0.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

No. 4 Surgical Scalpel Handle — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type make precise incisions and separate tissue planes. Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price. Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002761

No-Touch Scalpel Blade Remover Box Holds Up to 100 Blades

\$9.99 per piece

Trade price from \$9.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Disposable Blade Remover Handheld No-Touch Blades Cutter Blade Insertion Slot Remover Box Case Holder for Scalpels, Dermaplaning, Razor, ENT, Medical, Dental, Podiatry, Crafts, & More

This device is designed for the safe and easy removal of used disposable blades from any size or type of scalpel handle. It improves safety and hygiene by the elimination of direct hand-blade contact. Depending on the size, it can store 100+ used blades.

DISPOSE OF OLD SCALPEL BLADES

Concerned over how to dispose of old blades in your lock-on scalpel handles?

- Easily Remove Blades from Lock-on Handles.

- Spent Blades Are Stored Inside the Box.
- Safe Disposal with No Contact.
- Simply insert the scalpel into the remover, push a button, and pull out the handle.

DISPOSABLE BLADE REMOVER BENEFITS:

EASE OF USE: This Blade Remover permits the safe and simple removal of used disposable blades from any kind and size of scalpel handle.
SAFE AND HYGIENIC: No direct contact with blades provides the greatest possible protection for medical staff against infection. Sharps are securely and safely retained in the disposable unit.
DISPOSABLE PRODUCT: Once the cost-effective blade remover is filled



**WI-70001956****Pack of 10 Scalpel Blades #60 Autopsy Post Mortem Surgical Blade****\$12.99** per pieceTrade price from **\$11.95**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 10 Scalpel Blades #60, Carbon Steel Autopsy Post Mortem Surgical Blade

These Post Mortem Autopsy Scalpel Blades are made from surgical carbon steel and sharpened to perfection for a super clean cut. The precision shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Exceptionally strong and uniformly sharpened to maximum acuity. These Scalpel Sterile Blades used within the fields of Pathology and Histopathology. Used during Autopsy procedures they are also popular with the Morticians.

PACKAGE INCLUDES:

- 10 Pcs Scalpel Blades #60

SCALPEL BLADE #60 FEATURES:

- Convex cutting edge.
- Thicker blades reduce breakage
- Surgical-quality blades
- Individually sealed
- Peelable packing
- Sterile
- Fits #8 scalpel handle

PACKING:

- All blades packed in Individual pouches

Find our extensive selection of **Scalpel Blades and Handles**

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.

**WI-70003148****Pack of 10 Scalpel Blades No. 10 with Black Scalpel Handle #3****\$10.99** per pieceTrade price from **\$10.11**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack 10 Pcs Sterile Surgical Blades #10 with FREE Black Scalpel Knife Handle #3 – Perfect for Various Surgical Medical Procedures and Lab Work.

Scalpel blades No. 10 with a #3 black handle is designed for precise medical incisions, including C-sections, and is also suitable for beauty tasks like dermaplaning. Additionally, it serves in arts and crafts for cutting wood, leather, fabric, clay, and paper. The sterile blades are individually wrapped for safety and convenience, and the black handle features a reliable slip lock, compatible with various blades for detailed surgical work.

PACKAGE INCLUDE

- 10 Pcs of Surgical Scalpel Blade #10
- 1 Pc Blade Scalpel Handle #3

TECHNICAL INFORMATION:

MATERIAL: Carbon Steel, Stainless Steel

INSTRUMENT TYPE: Scalpel Blade & Handle

APPLICATION: Surgical Medical, Dermaplaning, Hobby & Crafts

RESUABLE: YES

STERILE: Non-Sterile

RUST PREVENTION

PROCEDURE: Passivated

ULTRASONIC CLEANED: Yes

TESTS PERFORMED: Boil Test,

Performance Test, Shape Test

QC PASSED: Yes

PACKING DETAILS: Individually Packed

Discover more **Scalpel Blades and Handles**





WI-70002186

Pack of 100 Scalpel Blades #60

\$69.99 per piece

Trade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 Surgical Blades #60
Carbon Steel Autopsy Post Mortem
Scalpel Blades fully Compatible with
Scalpel Handle #8.

These Post Mortem Autopsy Scalpel
Blades are made from surgical carbon
steel blades and sharpened to
perfection for a super clean cut. The
precision shaped carbon steel blades
are sharp and easy to use for a variety
of different tasks and/or procedures.
Exceptionally strong and uniformly
sharpened to maximum acuity. These
Scalpel Sterile Blades used within the
fields of Pathology and Histopathology.
Used during Autopsy procedures they
are also popular with the Morticians.

SET INCLUDED:

- 100 Pcs Scalpel Blades #60

SCALPEL BLADE #60 FEATURES:

- Convex cutting edge.
- Thicker blades reduce breakage
- Surgical-quality blades
- Individually sealed
- Peelable packing
- Sterile
- Fits #8 scalpel handle

PACKING:

- All blades packed in Individual
pouches

Find our extensive selection of **Scalpel
Blades and Handles**

Check our entire product range by
visiting this website's **Shop** link. Please
do **contact us** if your wish to buy larger
quantities, we can always arrange great
discounts for large order levels.



WI-70001342

Pack of 100 Scalpel Blades #70

\$69.99 per piece

Trade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 Surgical Blades #70
Carbon Steel Autopsy Post Mortem
Scalpel Blades fully Compatible with
Scalpel Handle #8.

The number 70 blade is primarily used
for postmortem, autopsy and other
research purposes. A smaller version of
the number 60B blade, this blade
features a rounded tip to reduce
puncture risk. The number 70 blade will
fit the number 8 scalpel handles.
Exceptionally strong and uniformly
sharpened to maximum acuity. These
post mortem autopsy scalpel blades
are made from surgical carbon steel
blades and sharpened to perfection for
a super clean cut. These scalpel sterile
blades used within the fields of
Pathology and Histopathology. Used

during Autopsy procedures they are
also popular with the Morticians.

PACKAGE INCLUDED:

- 100 Pcs Scalpel Blades #70

SCALPEL BLADE #70 FEATURES:

- Rounded tip.
- Thicker blades reduce breakage
- Surgical-quality blades
- Individually sealed
- Peelable packing
- Sterile
- Fits #8 scalpel handle

PACKING:

- All blades packed in Individual
pouches

Find our extensive selection of **Scalpel
Blades and Handles**




**WI-70002877****Pack of 100 Surgical Blades #10 with Scalpel Knife Handle #3****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 Sterile Surgical Blades #10 with Scalpel Knife Handle #3 Medical ENT Dental.


- The #10 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.
- They come clean and ready to use, 100% brand new and sealed.
- Precisely shaped and balanced for high degree of flexibility while conducting procedures.
- Scalpel handle No 3 can be used with carbon steel or stainless steel blades #10, #10A, #11, #12, #14, #15, #15C and #16.
- Scalpel handle has a graduation scale which is etched with measurements from 0-6cm and division value of 1mm.
- Grooved handle provides a strong grip to perform delicate work and/or detailed inspections.
- Ergonomic design helps reduce hand fatigue during repetitive tasks.
- Non- Reflective smooth finish for easy cleaning and disinfecting the instrument.
- Autoclavable and passivated for prolonging the instrument's life span.
- Premium quality aesthetic finish resulting from superior

**WI-70002876****Pack of 100 Surgical Blades #11 with Scalpel Knife Handle #3****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 Sterile Surgical Blades #11 with Scalpel Knife Handle #3 Medical ENT Dental.

- The #11 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.
- They come clean and ready to use, 100% brand new and sealed.
- Precisely shaped and balanced for high degree of flexibility while conducting procedures.
- Scalpel handle No 3 can be used with carbon steel or stainless steel blades #10, #10A, #11, #12, #14, #15, #15C and #16.
- Scalpel handle has a graduation scale which is etched with measurements from 0-6cm and division value of 1mm.
- Grooved handle provides a strong grip to perform delicate work and/or detailed inspections.
- Ergonomic design helps reduce hand fatigue during repetitive tasks.
- Non- Reflective smooth finish for easy cleaning and disinfecting the instrument.
- Autoclavable and passivated for prolonging the instrument's life span.
- Premium quality aesthetic finish resulting from superior




**WI-70002875****Pack of 100 Surgical Blades #12 with Scalpel Knife Handle #3****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 Sterile Surgical Blades #12 with Scalpel Knife Handle #3 Medical ENT Dental.


- The #12 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.
- They come clean and ready to use, 100% brand new and sealed.
- Precisely shaped and balanced for high degree of flexibility while conducting procedures.
- Scalpel handle No 3 can be used with carbon steel or stainless steel blades #10, #10A, #11, #12, #14, #15, #15C and #16.
- Scalpel handle has a graduation scale which is etched with measurements from 0-6cm and division value of 1mm.
- Grooved handle provides a strong grip to perform delicate work and/or detailed inspections.
- Ergonomic design helps reduce hand fatigue during repetitive tasks.
- Non- Reflective smooth finish for easy cleaning and disinfecting the instrument.
- Autoclavable and passivated for prolonging the instrument's life span.
- Premium quality aesthetic finish resulting from superior

**WI-70002874****Pack of 100 Surgical Blades #15 with Scalpel Knife Handle #3****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 Sterile Surgical Blades #15 with Scalpel Knife Handle #3 Medical ENT Dental.

- The #15 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.
- They come clean and ready to use, 100% brand new and sealed.
- Precisely shaped and balanced for high degree of flexibility while conducting procedures.
- Scalpel handle No 3 can be used with carbon steel or stainless steel blades #10, #10A, #11, #12, #14, #15, #15C and #16.
- Scalpel handle has a graduation scale which is etched with measurements from 0-6cm and division value of 1mm.
- Grooved handle provides a strong grip to perform delicate work and/or detailed inspections.
- Ergonomic design helps reduce hand fatigue during repetitive tasks.
- Non- Reflective smooth finish for easy cleaning and disinfecting the instrument.
- Autoclavable and passivated for prolonging the instrument's life span.
- Premium quality aesthetic finish resulting from superior




**WI-70002873****Pack of 100 Surgical Blades #15C with Scalpel Knife Handle #3****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 Sterile Surgical Blades #15C with Scalpel Knife Handle #3 Medical ENT Dental.

- The #15C scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries.
- Including health and wellness, beauty and skincare, mechanical and tech servicing, processioning laboratories, taxidermy and more.
- The scalpel blades are individually wrapped in foil pouches and are safe, effective and convenient.
- They come clean and ready to use, 100% brand new and sealed.
- Precisely shaped and balanced for high degree of flexibility while conducting procedures.
- Scalpel handle No 3 can be used with carbon steel or stainless steel blades #10, #10A, #11, #12, #14, #15, #15C and #16.
- Scalpel handle has a graduation scale which is etched with measurements from 0-6cm and division value of 1mm.
- Grooved handle provides a strong grip to perform delicate work and/or detailed inspections.
- Ergonomic design helps reduce hand fatigue during repetitive tasks.
- Non- Reflective smooth finish for easy cleaning and disinfecting the instrument.
- Autoclavable and passivated for prolonging the instrument's life span.
- Premium quality aesthetic finish resulting from superior

**WI-70002872****Pack of 100 Surgical Blades #20 with FREE Scalpel Handle #4****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Box of 100 sterile surgical blades #20 with **FREE** scalpel knife handle #4 Medical ENT Dental. All blades packed in Individual pouches.

SURGICAL SCALPEL BLADE NO.20 FEATURES

- Surgical blades size 20 are a larger version of the size 10 surgical blade with a curved cutting edge.
- #20 blade is used for larger incisions through thicker tissues or for dissection of the soft tissues in larger animals.
- Curved cutting edge with a flat back and a "belly" shaped cutting point.
- No. 20 rib-back surgical blades packaged sterile and individually pouched for safe use.
- Made with precision carbon steel.


SET INCLUDED:

- 100 Pieces of #20 Scalpel Blades
- Free Sterile Scalpel Handle #4 Included.

SCALPEL BLADES PERFECT FOR:

- Facial Dermaplaning.
- Podiatry.
- Surgery.
- Dead skin/ Pimple/ Blackhead removal.
- Arts and Crafts.
- Clay Modeling.
- Carving/ Wood Carving.
- Dissecting.
- Fabric cutting

Find more **Scalpel Blades and Handles**Check our entire product range by visiting this website's **Shop** link. Please

**WI-70002871****Pack of 100 Surgical Blades #21 with Scalpel Knife Handle #4****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

The pack of 100 sterile surgical blades #21 with scalpel knife handle #4
Medical ENT Dental

- So, sharp enough or not ? The perception of sharpness for #21 rounded blades can be different because depends on how experienced you are with dermaplaning or dermatology.
- They are very sharp for beginners and maybe not sharp enough for experienced dermatologists. Either way, it works if used correctly.
- Save some money- pack of 100 sterilized scalpel blades #21 help you learn to dermaplaning by yourself, and you can safely use it for a lot of small aesthetic procedures: great as a scalpel for callus, podiatry scalpel, dermatology scalpel pedicure, cyst

removal blade, blackheads removals, acne scalpel, dissection and suture practice, and more.

- While your other disposable scalpel blades #21 may be too dull, not sharp enough, not sterile, flimsy or difficult to use, we brings a high quality carbon steel scalpel blade 21 that is safe and easy to use. This means a younger brighter face for you and lots of other skin problems solved.
- A ton of other uses – mycology, Comedones removal, prototyping electronics, model making, wood and clay sculpture, desktop work, small cuts and trimming, biology projects in school and labs, arts and crafts, trimming, gardening, vinyl cuts, wood projects, medical procedures, dissection, suture and surgery

**WI-70002870****Pack of 100 Surgical Blades #22 with Scalpel Knife Handle #4****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

The pack of 100 sterile surgical blades #22 with scalpel knife handle #4
Medical ENT Dental

- So, sharp enough or not ? The perception of sharpness for #22 rounded blades can be different because depends on how experienced you are with dermaplaning or dermatology.
- They are very sharp for beginners and maybe not sharp enough for experienced dermatologists. Either way, it works if used correctly.
- Save some money- pack of 100 sterilized scalpel blades #22 help you learn to dermaplaning by yourself, and you can safely use it for a lot of small aesthetic procedures: great as a scalpel for callus, podiatry scalpel, dermatology

scalpel pedicure, cyst removal blade, blackheads removals, acne scalpel, dissection and suture practice, and more.

- While your other disposable scalpel blades #22 may be too dull, not sharp enough, not sterile, flimsy or difficult to use, we brings a high quality carbon steel scalpel blade 22 that is safe and easy to use. This means a younger brighter face for you and lots of other skin problems solved.
- A ton of other uses – mycology, Comedones removal, prototyping electronics, model making, wood and clay sculpture, desktop work, small cuts and trimming, biology projects in school and labs, arts and crafts, trimming, gardening, vinyl cuts,

**WI-70002869****Pack of 100 Surgical Blades #23 with Scalpel Knife Handle #4****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

The pack of 100 sterile surgical blades #23 with scalpel knife handle #4
Medical ENT Dental

- So, sharp enough or not ? The perception of sharpness for #23 rounded blades can be different because depends on how experienced you are with dermaplaning or dermatology.
- They are very sharp for beginners and maybe not sharp enough for experienced dermatologists. Either way, it works if used correctly.
- Save some money- pack of 100 sterilized scalpel blades #23 help you learn to dermaplaning by yourself, and you can safely use it for a lot of small aesthetic procedures: great as a scalpel for callus, podiatry scalpel, dermatology

scalpel pedicure, cyst removal blade, blackheads removals, acne scalpel, dissection and suture practice, and more.

- While your other disposable scalpel blades #23 may be too dull, not sharp enough, not sterile, flimsy or difficult to use, we bring a high quality carbon steel scalpel blade 23 that is safe and easy to use. This means a younger brighter face for you and lots of other skin problems solved.
- A ton of other uses – mycology, Comedones removal, prototyping electronics, model making, wood and clay sculpture, desktop work, small cuts and trimming, biology projects in school and labs, arts and crafts, trimming, gardening, vinyl cuts,

**WI-70002868****Pack of 100 Surgical Blades #24 with Scalpel Knife Handle #4****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

The pack of 100 sterile surgical blades #24 with scalpel knife handle #4
Medical ENT Dental

- So, sharp enough or not ? The perception of sharpness for #24 rounded blades can be different because depends on how experienced you are with dermaplaning or dermatology.
- They are very sharp for beginners and maybe not sharp enough for experienced dermatologists. Either way, it works if used correctly.
- Save some money- pack of 100 sterilized scalpel blades #24 help you learn to dermaplaning by yourself, and you can safely use it for a lot of small aesthetic procedures: great as a scalpel for callus, podiatry scalpel, dermatology

scalpel pedicure, cyst removal blade, blackheads removals, acne scalpel, dissection and suture practice, and more.

- While your other disposable scalpel blades #24 may be too dull, not sharp enough, not sterile, flimsy or difficult to use, we bring a high quality carbon steel scalpel blade 24 that is safe and easy to use. This means a younger brighter face for you and lots of other skin problems solved.
- A ton of other uses – mycology, Comedones removal, prototyping electronics, model making, wood and clay sculpture, desktop work, small cuts and trimming, biology projects in school and labs, arts and crafts, trimming, gardening, vinyl cuts,



WI-70002867

Pack of 100 Surgical Blades #25 with Scalpel Knife Handle #4

\$11.99 per piece

Trade price from **\$11.03**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

The pack of 100 sterile surgical blades #25 with scalpel knife handle #4 Medical ENT Dental

- So, sharp enough or not ? The perception of sharpness for #25 rounded blades can be different because depends on how experienced you are with dermaplaning or dermatology.
- They are very sharp for beginners and maybe not sharp enough for experienced dermatologists. Either way, it works if used correctly.
- Save some money- pack of 100 sterilized scalpel blades #25 help you learn to dermaplaning by yourself, and you can safely use it for a lot of small aesthetic procedures: great as a scalpel for callus, podiatry scalpel, dermatology

scalpel pedicure, cyst removal blade, blackheads removals, acne scalpel, dissection and suture practice, and more.

- While your other disposable scalpel blades #25 may be too dull, not sharp enough, not sterile, flimsy or difficult to use, we brings a high quality carbon steel scalpel blade 25 that is safe and easy to use. This means a younger brighter face for you and lots of other skin problems solved.
- A ton of other uses – mycology, Comedones removal, prototyping electronics, model making, wood and clay sculpture, desktop work, small cuts and trimming, biology projects in school and labs, arts and crafts, trimming, gardening, vinyl cuts,



WI-70001957

Pack of 20 Scalpel Blades #60 Autopsy Post Mortem Surgical Blade

\$17.99 per piece

Trade price from **\$16.55**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 20 Scalpel Blades #60, Carbon Steel Autopsy Post Mortem Surgical Blade

These Post Mortem Autopsy Scalpel Blades are made from surgical carbon steel and sharpened to perfection for a super clean cut. The precision shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Exceptionally strong and uniformly sharpened to maximum acuity. These Scalpel Sterile Blades used within the fields of Pathology and Histopathology. Used during Autopsy procedures they are also popular with the Morticians.

PACKAGE INCLUDES:

- 20 Pcs Scalpel Blades #60

SCALPEL BLADE #60 FEATURES:

- Convex cutting edge.
- Thicker blades reduce breakage
- Surgical-quality blades
- Individually sealed
- Peelable packing
- Sterile
- Fits #8 scalpel handle

PACKING:

- All blades packed in Individual pouches

Find our extensive selection of **Scalpel Blades and Handles**

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.





WI-70001958

Pack of 50 Scalpel Blades #60 Autopsy Post Mortem Surgical Blade

\$35.99 per piece

Trade price from \$33.11
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 50 Scalpel Blades #60, Carbon Steel Autopsy Post Mortem Surgical Blade

These Post Mortem Autopsy Scalpel Blades are made from surgical carbon steel and sharpened to perfection for a super clean cut. The precision shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Exceptionally strong and uniformly sharpened to maximum acuity. These Scalpel Sterile Blades used within the fields of Pathology and Histopathology. Used during Autopsy procedures they are also popular with the Morticians.

PACKAGE INCLUDES:

- 50 Pcs Scalpel Blades #60

SCALPEL BLADE #60 FEATURES:

- Convex cutting edge.
- Thicker blades reduce breakage
- Surgical-quality blades
- Individually sealed
- Peelable packing
- Sterile
- Fits #8 scalpel handle

PACKING:

- All blades packed in Individual pouches

Find our extensive selection of **Scalpel Blades and Handles**

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.



WI-28189918

PCS Scalpel Handle No 3 Stainless Steel Derma

\$1.50 per piece

Trade price from \$0.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

PCS Scalpel Handle No 3 Stainless Steel Derma — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type make precise incisions and separate tissue planes.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000396

Post Mortem Knife 6.5" Smooth Handle

\$19.99 per piece

Trade price from \$18.39
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Stainless Steel Autopsy Post Mortem Dissection Knife Blade 6.5" (16.5cm) Solid Smooth Handle for Histopathology and General Post Mortem Procedures.

Post mortem dissection knives, with bull nose tips, are specialized surgical-grade tools used in autopsies and postmortem examinations for precise tissue and organ examination. Post mortem knives are essential for precise incisions in tissues and organs, enabling pathologists to examine them thoroughly during autopsy, determining the cause of death. Specialized blades are available for cutting through cartilage and soft bone structures, as well as deep cuts and scraping. Post Mortem Autopsy scalpel blades are of superior quality, crafted from surgical

steel and expertly sharpened for a precise and clean cut.

POST MORTEM DISSECTION KNIFE BLADE FEATURES

- 6.5" (16.5cm) Overall Length.
- "Bull Nose Tip" Blade.
- Heavy-duty Blade for Postmortem Procedures.
- Helps in Preventing Accidental Sharp Injuries.
- Smooth Solid Handle
- Reusable.
- Made with High-quality Stainless Steel.

PACKAGE INCLUDE:

- 1 Pc Post Mortem Knife 6.5" Smooth Handle

TECHNICAL INFORMATION:



WI-30864085

Post Mortem Scalpel Autopsy Cartilage Knives

\$10.00 per piece

Trade price from \$7.80
Minimum order 3 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


This is our pattern of Post Mortem Scalpel Autopsy Cartilage Knives, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$10.00, falling to \$7.80 at 3 pieces — about 22% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-45266366

Reusable Stainless Steel Scalpel Handle Surgical Blade Holder

\$2.00 per piece

Trade price from **\$1.40**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Reusable Stainless Steel Scalpel Handle Surgical Blade Holder, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$2.00, falling to \$1.40 at 25 pieces — about 30% off the single-unit price. Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001008

Round Scalpel Handle #3 Straight Fits Blade 10-16

\$9.99 per piece

Trade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Round Scalpel Handle #3 Straight 6" (15.2cm) Long Premium Quality Stainless Steel compatible with Surgical Blades No. 10 to 16 – Reusable.

Scalpel handle, round and tapered style provides greater rotational movement with fingertip control and enhances dexterity. Handle measures 6" (15cm) long. Compatible with standard scalpel blade #10, #11, #12, #15, #15C, #16. Can be reused after proper cleaning/sterilization. This scalpel handle features an ergonomically designed pen-shaped handle that enhances gripping with knurled style cover that minimizes the risk of slippage and reduces fatigue. Made with high quality stainless steel ensure durability. This handle can be sterilized and reused.

PACKAGE INCLUDES

- 1 Pc Round Scalpel Handle #3 Fitt Blades No. 10-16

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENT TYPE: Scalpel Handle
RESUABLE: Yes
STERILEABILITY: Sterile
RUST PREVENTION PROCEDURE: Passivated
ULTRASONIC CLEANED: Yes
TESTS PERFORMED: Boil Test, Performance Test, Shape Test
QC PASSED: Yes
PACKING DETAILS: Polybag Packing

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger





WI-45143074

Scalpel Blade Handle

\$0.80 per piece

Trade price from \$0.30
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Scalpel Blade Handle is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type make precise incisions and separate tissue planes. Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$0.80, falling to \$0.30 at 100 pieces — about 62% off the single-unit price. Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003147

Scalpel Handle #3

\$8.99 per piece

Trade price from \$8.27
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Scalpel Handle #3
Stainless Steel Surgical Knife Handle – Reusable.

Premium quality, rust proof stainless steel scalpel knife handle is lightweight and durable, fits for surgical blades No. 10, 11, 12, 15 and 15C. The scalpel handle #3 ergonomically designed for hand precise fingertip control required for hand surgery and other delicate cutting or dissection. The Knife Handle #3 is usually used with smaller blades that have a curved cutting edge, like the #10 or #15 replaceable surgical blades. These handles and blade work well for microdissection or fine surgeries when making precise incisions. These handles are manufactured from surgical grade


stainless steel and suitable for autoclaving or chemical sterilization.


SCALPEL HANDLE #3 FEATURES

- 5" (12.7cm) Overall Legth.
- Knife Handle Fits Blades From Size 10 to 15.
- Textured Gripping Surface for a Firm Hold.
- Blades Slide on Smoothly and Securely.
- Reusable Handle.
- Crafted with Premium Quality Stainless Steel.

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel
INSTRUMENT TYPE: Scalpel Handle
HANDLE: Non Slip
RESUABLE: YES



**WI-70000162****Scalpel Handle #3 & #4 – Double Ended****\$10.99** per pieceTrade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Scalpel Handle #3 & #4 Double Ended Can Holds #10 – #15 Scalpel Blades on One End (#3) and #21 – #25 Scalpel Blades on Other End (#4) Premium Quality Precise Surgery instrument.

Double-ended scalpel blade handle combine #3 and #4 configurations, allowing for versatility in surgical procedures by accommodating different blade sizes on each end. This scalpel handle enable surgeons to switch between a #3 handle for precision dissection and a #4 handle for deeper cuts, streamlining operations and minimizing instrument changes. It features a double-ended solid construction that provides versatility, secure blade attachment, and balanced handling, suitable for a range of surgical and clinical

procedures. This handle facilitate rapid, clean incisions and reduced bleeding, making them suitable for both delicate and larger surgical tasks, and are typically constructed from durable stainless steel for sterilization purposes.

DOUBLE END SCALPEL HANDLE FEATURES

- 6.5 Inch (16.5cm) Overall Length
- Double Ended Design
- Handle #3 End (Small): Uses blades 10, 11, 12, 15 for delicate, precise work.
- Handle #4 End (Large): Uses blades 20, 21, 22, 23, 24, 25 for broader incisions.
- Suitable for Procedures Needing Various Blade Sizes.

**WI-70003146****Scalpel Handle #3L****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Long Scalpel Handle #3L compatible with scalpel blades 10, 11, 12, 15, 15C made with Stainless Steel. 3L scalpel handle ergonomically designed for hand precise fingertip control required for hand surgery and other delicate cutting along curved lines or dissection. These handles and blade work well for microdissection or fine surgeries when making precise incisions. These handles are manufactured from surgical grade stainless steel and suitable for autoclaving or chemical sterilization.

PACKAGE INCLUDE

- 1 Pc Long Scalpel Handle #3L

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Scalpel Handle**HANDLE:** Non Slip**RESUABLE:** YES**STERILEABILITY:** Sterile**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing

Discover more **Scalpel Blades and Handle**

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.

**WI-70003145****Scalpel Handle #3LA Angled****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Scalpel Handle #3LA Long Angled compatible with scalpel blades 10, 11, 12, 15 and 15C made with Stainless Steel. The scalpel handle 3LA ergonomically designed for hand precise fingertip control required for hand surgery and other delicate cutting along curved lines or dissection.

STERILEABILITY: Sterile**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Individually Packed**TECHNICAL INFORMATION:****MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Scalpel Handle**HANDLE:** Non Slip**RESUABLE:** YES

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.

**WI-70003144****Scalpel Handle #4****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Scalpel Handle #4 Stainless Steel Surgical Knife Handle – Reusable.

The Scalpel Handle #4 is a little larger and heavier than the #3. It gives you more control and the ability to apply more force. The handle #4 is often used with larger blades like the #20, #21, #22, #23, #24, and #25 blades. These blades have a similar shapes as the blades used with the #3 handle, but the blades are larger in size. Premium quality, rust proof stainless steel scalpel knife handle is lightweight and durable. The scalpel handle #4 ergonomically designed for hand precise fingertip control required for hand surgery and other delicate cutting or dissection.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Scalpel Handle**HANDLE:** Non Slip**RESUABLE:** YES**STERILEABILITY:** Sterile**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Individually Packed

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.



**WI-70003143****Scalpel Handle #4L****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Scalpel Handle #4L compatible with scalpel blades 20, 21, 22, 23, 24 and 25 made with Stainless Steel. 4L scalpel handle is similar to the #4 scalpel handle, only longer. The scalpel handle #4L ergonomically designed for hand precise fingertip control required for hand surgery and other delicate cutting or dissection. Scalpel handles are versatile tools with slotted tips that securely hold replaceable blades for various cutting tasks. Scalpel Handle No 4L provide reliable support for multiple blades to perform precise incisions.

We have surgical knife handle No.4 in standard or long versions.

SCALPEL HANDLE #4L FEATURES

- Fits Various Blade Profiles.

- Reduces Slippage Through Transverse Ridges.
- Offers Ergonomic Flat Design for Enhanced Fingertip Control.
- Reusable.
- Made From High-quality Stainless Steel.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Scalpel Handle**HANDLE:** Non Slip**RESUABLE:** YES**STERILEABILITY:** Sterile**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**WI-70003142****Scalpel Handle #4LA Angled****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Scalpel Handle #4LA Long Angled Surgical Knife handle is compatible with scalpel blades 20, 21, 22, 23, 24 and 25 made with stainless steel. The scalpel handle 4LA angled ergonomically designed for hand precise fingertip control required for hand surgery and other delicate cutting along curved lines or dissection.

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Scalpel Handle**HANDLE:** Non Slip**RESUABLE:** YES**STERILEABILITY:** Sterile**RUST PREVENTION PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Individually Packed

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger quantities, we can always arrange great discounts for large order levels.

**WI-70001018****Scalpel Handle #5 with Green Plastic Grip****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Scalpel Handle #5 with Green Hard Plastic Grip, Stainless steel Tip Compatible with Surgical Blade No. 10 to 15C – Reusable.

Scalpel handle No. 5 is designed ergonomically for hand surgery and other delicate cutting along curved lines or dissection. Wide plastic handle provides a strong grip to perform delicate and detailed work with ease and reduce hand fatigue during repetitive tasks. This plastic grip scalpel handle is perfect for wide range of applications including dissection experimentation, medical/surgical procedures and for hobby/ arts & crafts. Handle scientifically shaped and blanded and blades slides on smoothly and securely. Can be reused after proper cleaning/sterilization.

PACKAGE INCLUDES

- 1 Pc Scalpel Handle #5 with Green Plastic Grip

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Scalpel Handle**HANDLE:** Hard Plastic Grip**RESUABLE:** Yes**STERILEABILITY:** Sterile**RUST PREVENTION PROCEDURE:**



Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test,

Performance Test, Shape Test

QC PASSED: Yes**PACKING DETAILS:** Individually Packed

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger

**WI-70001019****Scalpel Handle #6 with Red Plastic Grip****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality Scalpel Handle #6 with Red Hard Plastic Grip, Stainless steel Tip Compatible with Surgical Blade No. 20 to 25 – Reusable.

Scalpel handle No. 6 is designed ergonomically for hand surgery and other delicate cutting along curved lines or dissection. Wide plastic handle provides a strong grip to perform delicate and detailed work with ease and reduce hand fatigue during repetitive tasks. This plastic grip scalpel handle is perfect for wide range of applications including dissection experimentation, medical/surgical procedures and for hobby/ arts & crafts. Handle scientifically shaped and blanded and blades slides on smoothly and securely. Can be reused after proper cleaning/sterilization.

PACKAGE INCLUDES

- 1 Pc Scalpel Handle #6 with Red Plastic Grip

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Scalpel Handle**HANDLE:** Hard Plastic Grip**RESUABLE:** Yes**STERILEABILITY:** Sterile**RUST PREVENTION PROCEDURE:**



Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test,

Performance Test, Shape Test

QC PASSED: Yes**PACKING DETAILS:** Individually Packed

Check our entire product range by visiting this website's **Shop** link. Please do **contact us** if your wish to buy larger



**WI-70003141**

Scalpel Handle #7 Surgical Scalpel Knife Handle Medical Dental ENT

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Scalpel Handle # 7, Premium Quality, Rust proof Stainless Steel Scalpel Knife Handle, Lightweight and Durable, Fits Surgical Blades No. 10, 11, 12, 13, 14 and 15

Premium Quality – High Polished thick metal handle durable and comfortable to use. Created with Superior craftsmanship, with high degree of precision and flexibility for smooth and easy use in conducting clinical procedures.

Rust Proof – Scalpel handle No. 7 is designed ergonomically, hand precise fingertip control required for hand surgery and other delicate cutting along curved lines or dissection. The High Quality Stainless Steel assure a rust proof handle.

Reusable Handle – Yes, this can be use several times, as Waqar International scalpel handle #7 is made of stainless steel, so it is durable and long lasting. Can be used by professional, experts and amateurs, alike.

Polished Finish – Highly polished finish for aesthetic and corrosion resistance. High Degree of Precision and Flexibility while conducting the Clinical Procedure.

ADDITIONAL INFORMATION

Material: Stainless Steel

Rust Prevention Procedure: Passivated

Ultrasonic Cleaned: Yes

Tests Performed: Boil Test, Performance Test, Shape Test

Quality Checking Passed: Yes

**WI-70002185**

Scalpel Handle #8 for Scalpel Blades #60 & #70

\$9.99 per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Scalpel Handle No. 8 is designed to hold and manoeuvre blades during a variety of surgical procedures.

SCALPEL HANDLE NO.8 FEATURES

PREMIUM QUALITY: High Polished thick metal handle durable and comfortable to use. It is created with Superior craftsmanship, precision, and flexibility for smooth and easy use in conducting clinical procedures.

RUST PROOF: Scalpel handle No. 8 is designed ergonomically, hand precise fingertip control is required for hand surgery and other delicate cutting along curved lines or dissection. The High-Quality Stainless Steel assure a rust-proof handle.

REUSABLE HANDLE: Yes, this can be used several times, as this scalpel

handle #8 is made of stainless steel, so it is durable and long-lasting. It can be used by professionals, experts and amateurs, alike.

POLISH FINISH: Highly polished finish for aesthetic and corrosion resistance. High Degree of Precision and Flexibility while conducting the Clinical Procedure.

PACKAGE INCLUDE:

- 1 Pc Scalpel Handle #8

TECHNICAL INFORMATION:

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Surgical Blades Handle

APPLICATION: Surgical Procedures
COMPATIBILITY: with Scalpel Blade #60 & #70

FINISH: Dull/Finish





Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34528863

Scalpel Handle No 4 Surgical Instrument Set / Sialkot

\$1.50 per piece

Trade price from **\$0.50**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

This is our pattern of Scalpel Handle No 4 Surgical Instrument Set / Sialkot, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-47139992

Sickle Knife Surgical Cutting Tool for Fine Tissue Separation

\$4.80 per piece

Trade price from **\$1.80**
Minimum order 60 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

We manufacture Sickle Knife Surgical Cutting Tool for Fine Tissue Separation for general and specialty operating theatre use — made from bar stock to finished instrument in house.

Instruments of this type make precise incisions and separate tissue planes. The name tells you what you are getting: a fine pattern for delicate tissue.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$4.80, falling to \$1.80 at 60 pieces — about 62% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



WI-47143926

Sickie Knife Surgical Instrument for Precision Tissue Cutting Dissection

\$5.00 per piece

Trade price from **\$2.00**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


We manufacture Sickie Knife Surgical Instrument for Precision Tissue Cutting Dissection for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type make precise incisions and separate tissue planes. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$5.00, falling to \$2.00 at 25 pieces — about 60% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-47157822

Sickie Surgical Knife Precision Cutting Instrument for Controlled Tissue

\$4.90 per piece

Trade price from **\$1.90**
Minimum order 50 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Sickie Surgical Knife Precision Cutting Instrument for Controlled Tissue — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type make precise incisions and separate tissue planes. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$4.90, falling to \$1.90 at 50 pieces — about 61% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70002410****Sims Uterine Curette 11\"/>****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Sims Uterine Curette Size #0 Overall Length 11\"/>

Sims Uterine Curettes are versatile gynecology surgical instruments designed to remove tough tissue from the lining of the uterus for biopsy or excision. Overall curette's length is 11\"/>

handle with vertical grooves to ensure convenient handling.

SIMS UTERINE CURETTE FEATURES


- Sharp fenestrated oval shape loop tip available in various sizes.
- Working End Sizes 5mm (#00) – 14mm (#6)
- Malleable shaft for manipulation to achieve easy access.
- Autoclaved, Reusable, and Easy to clean.
- Single-Ended OB/GYN Surgical Instrument with Finger Grip Hollow Handle.

PACKAGE INCLUDE:

- 1 Pc Sims Uterine Curette 11\"/>

SIZE 0

SIMS UTERINE CURETTES

**WI-70002411****Sims Uterine Curette 11\"/>****\$11.24** per pieceTrade price from **\$10.34**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Sims Uterine Curette Size #00 Overall Length 11\"/>

Sims Uterine Curettes are versatile gynecology surgical instruments designed to remove tough tissue from the lining of the uterus for biopsy or excision. Overall curette's length is 11\"/>

handle with vertical grooves to ensure convenient handling.

SIMS UTERINE CURETTE FEATURES


- Sharp fenestrated oval shape loop tip available in various sizes.
- Working End Sizes 5mm (#00) – 14mm (#6)
- Malleable shaft for manipulation to achieve easy access.
- Autoclaved, Reusable, and Easy to clean.
- Single-Ended OB/GYN Surgical Instrument with Finger Grip Hollow Handle.

PACKAGE INCLUDE:

- 1 Pc Sims Uterine Curette 11\"/>

SIZE 00

SIMS UTERINE CURETTES



**WI-70002409****Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #1****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Sims Uterine Curette Size #1 Overall Length 11" (27.9cm) Sharp Blades on Malleable Shank Oval Shaped Fenestrated Loop Tip Single-Ended Stainless Steel OB/Gynecology Surgical Instrument

Sims Uterine Curettes are versatile gynecology surgical instruments designed to remove tough tissue from the lining of the uterus for biopsy or excision. Overall curette's length is 11" and features an oval-shaped fenestrated sharp look tip and Malleable shaft. The small sharp loop allows quick pickup of the sample without causing any damage to nearby tissues or surgical areas. Multiple sizes to accommodate various cases. These uterine curettes consist of a finger-grip

handle with vertical grooves to ensure convenient handling.

SIMS UTERINE CURETTE FEATURES



- Sharp fenestrated oval shape loop tip available in various sizes.
- Working End Sizes 5mm (#00) – 14mm (#6)
- Malleable shaft for manipulation to achieve easy access.
- Autoclaved, Reusable, and Easy to clean.
- Single-Ended OB/GYN Surgical Instrument with Finger Grip Hollow Handle.

PACKAGE INCLUDE:

- 1 Pc Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #1 on Malleable Shank

SIZE 1

SIMS UTERINE CURETTES

**WI-70002408****Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #2****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Sims Uterine Curette Size #2 Overall Length 11" (27.9cm) Sharp Blades on Malleable Shank Oval Shaped Fenestrated Loop Tip Single-Ended Stainless Steel OB/Gynecology Surgical Instrument

Sims Uterine Curettes are versatile gynecology surgical instruments designed to remove tough tissue from the lining of the uterus for biopsy or excision. Overall curette's length is 11" and features an oval-shaped fenestrated sharp look tip and Malleable shaft. The small sharp loop allows quick pickup of the sample without causing any damage to nearby tissues or surgical areas. Multiple sizes to accommodate various cases. These uterine curettes consist of a finger-grip

handle with vertical grooves to ensure convenient handling.

SIMS UTERINE CURETTE FEATURES


- Sharp fenestrated oval shape loop tip available in various sizes.
- Working End Sizes 5mm (#00) – 14mm (#6)
- Malleable shaft for manipulation to achieve easy access.
- Autoclaved, Reusable, and Easy to clean.
- Single-Ended OB/GYN Surgical Instrument with Finger Grip Hollow Handle.


PACKAGE INCLUDE:

- 1 Pc Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #2 on Malleable Shank

SIZE 2

SIMS UTERINE CURETTES



**WI-70002407****Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #3****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Sims Uterine Curette Size #3 Overall Length 11" (27.9cm) Sharp Blades on Malleable Shank Oval Shaped Fenestrated Loop Tip Single-Ended Stainless Steel OB/Gynecology Surgical Instrument

Sims Uterine Curettes are versatile gynecology surgical instruments designed to remove tough tissue from the lining of the uterus for biopsy or excision. Overall curette's length is 11" and features an oval-shaped fenestrated sharp look tip and Malleable shaft. The small sharp loop allows quick pickup of the sample without causing any damage to nearby tissues or surgical areas. Multiple sizes to accommodate various cases. These uterine curettes consist of a finger-grip

handle with vertical grooves to ensure convenient handling.

SIMS UTERINE CURETTE FEATURES

- Sharp fenestrated oval shape loop tip available in various sizes.
- Working End Sizes 5mm (#00) – 14mm (#6)
- Malleable shaft for manipulation to achieve easy access.
- Autoclaved, Reusable, and Easy to clean.
- Single-Ended OB/GYN Surgical Instrument with Finger Grip Hollow Handle.

PACKAGE INCLUDE:

- 1 Pc Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #3 on Malleable Shank

SIZE 3

SIMS UTERINE CURETTES

**WI-70002406****Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #4****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Sims Uterine Curette Size #4 Overall Length 11" (27.9cm) Sharp Blades on Malleable Shank Oval Shaped Fenestrated Loop Tip Single-Ended Stainless Steel OB/Gynecology Surgical Instrument

Sims Uterine Curettes are versatile gynecology surgical instruments designed to remove tough tissue from the lining of the uterus for biopsy or excision. Overall curette's length is 11" and features an oval-shaped fenestrated sharp look tip and Malleable shaft. The small sharp loop allows quick pickup of the sample without causing any damage to nearby tissues or surgical areas. Multiple sizes to accommodate various cases. These uterine curettes consist of a finger-grip

handle with vertical grooves to ensure convenient handling.

SIMS UTERINE CURETTE FEATURES

- Sharp fenestrated oval shape loop tip available in various sizes.
- Working End Sizes 5mm (#00) – 14mm (#6)
- Malleable shaft for manipulation to achieve easy access.
- Autoclaved, Reusable, and Easy to clean.
- Single-Ended OB/GYN Surgical Instrument with Finger Grip Hollow Handle.

PACKAGE INCLUDE:

- 1 Pc Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #4 on Malleable Shank

SIZE 4

SIMS UTERINE CURETTES



**WI-70002405****Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #5****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Sims Uterine Curette Size #5 Overall Length 11" (27.9cm) Sharp Blades on Malleable Shank Oval Shaped Fenestrated Loop Tip Single-Ended Stainless Steel OB/Gynecology Surgical Instrument

Sims Uterine Curettes are versatile gynecology surgical instruments designed to remove tough tissue from the lining of the uterus for biopsy or excision. Overall curette's length is 11" and features an oval-shaped fenestrated sharp look tip and Malleable shaft. The small sharp loop allows quick pickup of the sample without causing any damage to nearby tissues or surgical areas. Multiple sizes to accommodate various cases. These uterine curettes consist of a finger-grip

handle with vertical grooves to ensure convenient handling.

SIMS UTERINE CURETTE FEATURES


- Sharp fenestrated oval shape loop tip available in various sizes.
- Working End Sizes 5mm (#00) – 14mm (#6)
- Malleable shaft for manipulation to achieve easy access.
- Autoclaved, Reusable, and Easy to clean.
- Single-Ended OB/GYN Surgical Instrument with Finger Grip Hollow Handle.

PACKAGE INCLUDE:

- 1 Pc Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #5 on Malleable Shank

SIZE 5

SIMS UTERINE CURETTES

**WI-70002404****Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #6****\$9.99** per pieceTrade price from **\$9.19**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Sims Uterine Curette Size #6 Overall Length 11" (27.9cm) Sharp Blades on Malleable Shank Oval Shaped Fenestrated Loop Tip Single-Ended Stainless Steel OB/Gynecology Surgical Instrument

Sims Uterine Curettes are versatile gynecology surgical instruments designed to remove tough tissue from the lining of the uterus for biopsy or excision. Overall curette's length is 11" and features an oval-shaped fenestrated sharp look tip and Malleable shaft. The small sharp loop allows quick pickup of the sample without causing any damage to nearby tissues or surgical areas. Multiple sizes to accommodate various cases. These uterine curettes consist of a finger-grip

handle with vertical grooves to ensure convenient handling.

SIMS UTERINE CURETTE FEATURES

- Sharp fenestrated oval shape loop tip available in various sizes.
- Working End Sizes 5mm (#00) – 14mm (#6)
- Malleable shaft for manipulation to achieve easy access.
- Autoclaved, Reusable, and Easy to clean.
- Single-Ended OB/GYN Surgical Instrument with Finger Grip Hollow Handle.

PACKAGE INCLUDE:

- 1 Pc Sims Uterine Curette 11" (27.9cm) Sharp Blade Size #6 on Malleable Shank

SIZE 6

SIMS UTERINE CURETTES





WI-28643901

Single-Use Surgical Scalpel Handle Stainless Steel Blade

\$1.50 per piece

Trade price from **\$0.50**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Single-Use Surgical Scalpel Handle Stainless Steel Blade, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70000626

Sinus Sickle Knife 7.5" Curved Blunt Tip

\$8.99 per piece

Trade price from **\$8.27**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Sinus Sickle Knife 7.5" (19.1cm) Curved Blunt Tip 10mm Blade with Flat Handle Premium Quality Stainless Steel ENT Surgical Rhinology Instruments.

Sinus sickle knife is an Otolaryngology instrument that is standardly used in nasal interventions and is commonly used in turbinate reductions, uncinectomy, anterior ethmoidectomy and during septoplasty procedures. The Sinus Sickle Knife's principal use is to perform atraumatic dissection of the delicate nasal and middle ear tissues. For this purpose, the knife features a sharp, sickle-shaped blade that fits small cavities and helps in creating precise endaural incisions/cuts. The straight, tapering shaft with a flat handle facilitates easy access into narrow spaces. Crafted from high-

grade stainless steel, our Nasal Sickle Knife guarantees durability and longevity, which are essential qualities for instruments used in surgical settings.

SINUS SICKLE KNIFE FEATURES

- Crafted with 7.5" (19.1cm) Overall Length.
- Blunt Tip Profile.
- 10mm Long Blade.
- Sharp Sickle-Shaped Blade for Highly Precise Incisions.
- Flat Handle with Straight Tapering Shaft.
- Various Pattern Tip Profile.
- Sterilized and Reusable.
- Made With Premium Quality Stainless Steel.

PACKAGE INCLUDES:

**WI-70001722****Smith Peterson Osteotome 8" Straight 13mm Blade****\$26.99** per pieceTrade price from **\$24.83**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Smith-Peterson Osteotome 8" (20.3cm), Straight 1/2" (13mm) wide Blade with Round Handle Stainless Steel Surgical Veterinary Orthopedic Instruments


Smith-Peterson Bone Osteotome is an orthopedic surgical tool that surgeons frequently use to cut, carve and reshape bone structures, in order to perform osteotomy and repair nonunion fractures or collapsed joints. These benefits play a vital role in dividing bone tissues and reshaping bone surfaces.

For these purposes, the instrument features a sharp beveled blade. Moreover, the blade is available in multiple sizes, ranging from 6.4 to 32mm. As a result, surgeons can pick the adequate size for treating small or

large bones, also contemplating the deep access to narrow spaces. In addition, Smith-Peterson Bone Osteotome features a solid, bulky handle supplied with a smooth surface. Also, the handle includes a bottom circular impacting platform. Surgeons can push the blade against the bone and hammer the platform with a mallet to cut cortical bone tissue. Specialty Used in Neurosurgery / Spine, Orthopedic, Post Mortem

SMITH PETERSON OSTEOTOMES FEATURES

- Atraumatic Sharp Beveled Cutting Edge for Smooth Cutting Action.
- From 6.4 to 32mm Different Working End Sizes for Deep Surgical Access.
- Ergonomic Solid Handle for Superior Control.

**WI-70001721****Smith Peterson Osteotome 8" Straight 16mm Blade****\$27.99** per pieceTrade price from **\$25.75**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Smith-Peterson Osteotome 8" (20.3cm), Straight 5/8" (16mm) wide Blade with Round Handle Stainless Steel Surgical Veterinary Orthopedic Instruments



Smith-Peterson Bone Osteotome is an orthopedic surgical tool that surgeons frequently use to cut, carve and reshape bone structures, in order to perform osteotomy and repair nonunion fractures or collapsed joints. These benefits play a vital role in dividing bone tissues and reshaping bone surfaces.


For these purposes, the instrument features a beveled blade. Moreover, the blade is available in multiple sizes, ranging from 6.4 to 32mm. As a result, surgeons can pick the adequate size for treating small or large bones, also

contemplating the deep access to narrow spaces. In addition, Smith-Peterson Bone Osteotome features a solid, bulky handle supplied with a smooth surface. Also, the handle includes a bottom circular impacting platform. Surgeons can push the blade against the bone and hammer the platform with a mallet to cut cortical bone tissue. Specialty Used in Neurosurgery / Spine, Orthopedic, Post Mortem

SMITH PETERSON OSTEOTOMES FEATURES

- Atraumatic Sharp Beveled Cutting Edge for Smooth Cutting Action.
- From 6.4 to 32mm Different Working End Sizes for Deep Surgical Access.
- Ergonomic Solid Handle for Superior Control.



**WI-70001720****Smith Peterson Osteotome 8" Straight 19mm Blade****\$28.99** per pieceTrade price from **\$26.67**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Smith-Peterson Osteotome 8" (20.3cm), Straight 3/4" (19mm) wide Blade with Round Handle Stainless Steel Surgical Veterinary Orthopedic Instruments



Smith-Peterson Bone Osteotome is an orthopedic surgical tool that surgeons frequently use to cut, carve and reshape bone structures, in order to perform osteotomy and repair nonunion fractures or collapsed joints. These benefits play a vital role in dividing bone tissues and reshaping bone surfaces.

For these purposes, the instrument features a sharp beveled blade. Moreover, the blade is available in multiple sizes, ranging from 6.4 to 32mm. As a result, surgeons can pick the adequate size for treating small or

large bones, also contemplating the deep access to narrow spaces. In addition, Smith-Peterson Bone Osteotome features a solid, bulky handle supplied with a smooth surface. Also, the handle includes a bottom circular impacting platform. Surgeons can push the blade against the bone and hammer the platform with a mallet to cut cortical bone tissue. Specialty Used in Neurosurgery / Spine, Orthopedic, Post Mortem

SMITH PETERSON OSTEOTOMES FEATURES

- Atraumatic Sharp Beveled Cutting Edge for Smooth Cutting Action.
- From 6.4mm to 32mm Different Working End Sizes for Deep Surgical Access.

**WI-70001719****Smith Peterson Osteotome 8" Straight 25mm Blade****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Smith-Peterson Osteotome 8" (20.3cm), Straight 1" (25mm) wide Blade with Round Handle Stainless Steel Surgical Veterinary Orthopedic Instruments

Smith-Peterson Bone Osteotome is an orthopedic surgical tool that surgeons frequently use to cut, carve and reshape bone structures, in order to perform osteotomy and repair nonunion fractures or collapsed joints. These benefits play a vital role in dividing bone tissues and reshaping bone surfaces.

For these purposes, the instrument features a sharp beveled blade. Moreover, the blade is available in multiple sizes, ranging from 6.4 to 32mm. As a result, surgeons can pick the adequate size for treating small or

large bones, also contemplating the deep access to narrow spaces. In addition, Smith-Peterson Bone Osteotome features a solid, bulky handle supplied with a smooth surface. Also, the handle includes a bottom circular impacting platform. Surgeons can push the blade against the bone and hammer the platform with a mallet to cut cortical bone tissue. Specialty Used in Neurosurgery / Spine, Orthopedic, Post Mortem

SMITH PETERSON OSTEOTOMES FEATURES

- Atraumatic Sharp Beveled Cutting Edge for Smooth Cutting Action.
- From 6.4mm to 32mm Different Working End Sizes for Deep Surgical Access.



**WI-70001718****Smith Peterson Osteotome 8" Straight 32mm Blade****\$32.99** per pieceTrade price from **\$30.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Smith-Peterson Osteotome 8" (20.3cm), Straight 1 1/4" (32mm) wide Blade with Round Handle Stainless Steel Surgical Veterinary Orthopedic Instruments

Smith-Peterson Bone Osteotome is an orthopedic surgical tool that surgeons frequently use to cut, carve and reshape bone structures, in order to perform osteotomy and repair nonunion fractures or collapsed joints. These benefits play a vital role in dividing bone tissues and reshaping bone surfaces.

For these purposes, the instrument features a sharp beveled blade. Moreover, the blade is available in multiple sizes, ranging from 6.4 to 32mm. As a result, surgeons can pick the adequate size for treating small or

large bones, also contemplating the deep access to narrow spaces. In addition, Smith-Peterson Bone Osteotome features a solid, bulky handle supplied with a smooth surface. Also, the handle includes a bottom circular impacting platform. Surgeons can push the blade against the bone and hammer the platform with a mallet to cut cortical bone tissue. Specialty Used in Neurosurgery / Spine, Orthopedic, Post Mortem

SMITH PETERSON OSTEOTOMES FEATURES

- Atraumatic Sharp Beveled Cutting Edge for Smooth Cutting Action.
- From 6.4mm to 32mm Different Working End Sizes for Deep Surgical Access.

**WI-70003272****Smith Peterson Osteotome 8" Straight 6.4mm Blade****\$24.99** per pieceTrade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Smith-Peterson Osteotome 8" (20.3cm), Straight 1/4" (6.4mm) wide Blade with Round Handle Stainless Steel Surgical Veterinary Orthopedic Instruments.



Smith-Peterson Bone Osteotome is an orthopedic surgical tool that surgeons frequently use to cut, carve and reshape bone structures, in order to perform osteotomy and repair nonunion fractures or collapsed joints. These benefits play a vital role in dividing bone tissues and reshaping bone surfaces.

For these purposes, the instrument features a sharp beveled blade. Moreover, the blade is available in multiple sizes, ranging from 6.4 to 32mm. As a result, surgeons can pick the adequate size for treating small or

large bones, also contemplating the deep access to narrow spaces. In addition, Smith-Peterson Bone Osteotome features a solid, bulky handle supplied with a smooth surface. Also, the handle includes a bottom circular impacting platform. Surgeons can push the blade against the bone and hammer the platform with a mallet to cut cortical bone tissue. Specialty Used in Neurosurgery / Spine, Orthopedic, Post Mortem

SMITH PETERSON OSTEOTOMES FEATURES

- Atraumatic Sharp Beveled Cutting Edge for Smooth Cutting Action.
- From 6.4 to 32mm Different Working End Sizes for Deep Surgical Access.
- Ergonomic Solid Handle for Superior Control.





WI-70001723

Smith Peterson Osteotome 8" Straight 9.5mm Blade

\$25.99 per piece

Trade price from **\$23.91**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Smith-Peterson Osteotome 8" (20.3cm), Straight 3/8" (9.5mm) wide Blade with Round Handle Stainless Steel Surgical Veterinary Orthopedic Instruments


Smith-Peterson Bone Osteotome is an orthopedic surgical tool that surgeons frequently use to cut, carve and reshape bone structures, in order to perform osteotomy and repair nonunion fractures or collapsed joints. These benefits play a vital role in dividing bone tissues and reshaping bone surfaces.

For these purposes, the instrument features a beveled blade. Moreover, the blade is available in multiple sizes, ranging from 6.4 to 32mm. As a result, surgeons can pick the adequate size for treating small or large bones, also

contemplating the deep access to narrow spaces. In addition, Smith-Peterson Bone Osteotome features a solid, bulky handle supplied with a smooth surface. Also, the handle includes a bottom circular impacting platform. Surgeons can push the blade against the bone and hammer the platform with a mallet to cut cortical bone tissue. Specialty Used in Neurosurgery / Spine, Orthopedic, Post Mortem

SMITH PETERSON OSTEOTOMES FEATURES

- Atraumatic Sharp Beveled Cutting Edge for Smooth Cutting Action.
- From 6.4 to 32mm Different Working End Sizes for Deep Surgical Access.
- Ergonomic Solid Handle for Superior Control.



WI-46905810

Solid Forged Scalpel Stainless Steel Surgical Instrument

\$2.00 per piece

Trade price from **\$1.00**

Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Solid Forged Scalpel Stainless Steel Surgical Instrument is a scalpels & knife we forge, grind and finish in our own Sialkot workshop. Instruments of this type make precise incisions and separate tissue planes.


Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.00, falling to \$1.00 at 25 pieces — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46897888

Solid Forged Surgical Scalpel Reusable Stainless Steel

\$1.90 per piece

Trade price from **\$0.90**
Minimum order 30 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Solid Forged Surgical Scalpel Reusable Stainless Steel — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type make precise incisions and separate tissue planes. The name tells you what you are getting: full reusability — clean, autoclave and use again.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$1.90, falling to \$0.90 at 30 pieces — about 53% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34473973

Stainless Steel Scalpel Graduated Handle #3

\$1.50 per piece

Trade price from **\$0.50**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Scalpel Graduated Handle #3 is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type make precise incisions and separate tissue planes. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-34488920

Stainless Steel Scalpel Handle No. 3 Surgical BP

\$1.50 per piece

Trade price from \$0.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Scalpel Handle No. 3 Surgical BP is a scalpels & knife we forge, grind and finish in our own Sialkot workshop. Instruments of this type make precise incisions and separate tissue planes.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-06301454

Stainless Steel Scalpel-Low Cost

\$2.00 per piece

Trade price from \$0.80
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Scalpel-Low Cost for general and specialty operating theatre use — made from bar stock to finished instrument in house. Instruments of this type make precise incisions and separate tissue planes.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.00, falling to \$0.80 at 1 piece — about 60% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



WI-06318441

Stainless Steel Scalpels Dull Polish Finish

\$1.00 per piece

Trade price from \$0.25
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Scalpels Dull Polish Finish is a scalpels & knife we forge, grind and finish in our own Sialkot workshop. Instruments of this type make precise incisions and separate tissue planes.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.00, falling to \$0.25 at 10 pieces — about 75% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-35052119

Stainless Steel Surgical Scalpel Blade

\$0.49 per piece

Trade price from \$0.46
Minimum order 4 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Scalpel Blade is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type make precise incisions and separate tissue planes.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$0.49, falling to \$0.46 at 4 pieces — about 6% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-45144024

Stainless Steel Surgical Scalpel Blade Holder

\$0.50 per piece

Trade price from \$0.30
Minimum order 100 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Scalpel Blade Holder is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type make precise incisions and separate tissue planes.

We work the instrument by hand at the finishing stage — box locks, hinges and

tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$0.50, falling to \$0.30 at 100 pieces — about 40% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28610950

Stainless Steel Surgical Scalpels Instrument Set Customizable

\$1.50 per piece

Trade price from \$0.50
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Scalpels Instrument Set Customizable — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type make precise incisions and separate tissue planes. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28226866

Steel Disposable Surgical Scalpel Handle Instrument

\$1.50 per piece

Trade price from **\$0.50**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Steel Disposable Surgical Scalpel Handle Instrument, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34469968

Steel Double Ended Scalpel Handle 16cm/21cm Blade Size

\$1.50 per piece

Trade price from **\$0.50**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Steel Double Ended Scalpel Handle 16cm/21cm Blade Size, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes. The name tells you what you are getting: two working ends on one handle, halving what you carry.

made to order. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

The 16cm size shown is our standard; other sizes in the same pattern are

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



Detail



WI-28292250

Steel Scalpel and Surgical Knife Handles No.3 & No.4

\$1.50 per piece

Trade price from **\$0.50**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Steel Scalpel and Surgical Knife Handles No.3 & No.4 is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type make precise incisions and separate tissue planes.

Each unit is inspected by hand for alignment and action before packing,

and we can photograph your actual pieces before they ship. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-28623886

Steel Surgical Scalpel Precision Instrument Set

\$1.50 per piece

Trade price from **\$0.50**

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Steel Surgical Scalpel Precision Instrument Set is a scalpels & knife we forge, grind and finish in our own Sialkot workshop. Instruments of this type make precise incisions and separate tissue planes. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34501066

Sterile Stainless Steel #10 Scalpel Blades #3 Handle

\$1.50 per piece

Trade price from \$0.50

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Sterile Stainless Steel #10 Scalpel Blades #3 Handle, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28603915

Surgical Instruments Medical Steel Scalpel Blade

\$1.50 per piece

Trade price from \$0.50

Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Surgical Instruments Medical Steel Scalpel Blade, made for general and specialty operating theatre use. Instruments of this type make precise incisions and separate tissue planes.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-45257446

Surgical Scalpel Handle Precision Blade Holder

\$2.00 per piece

Trade price from **\$1.40**
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Surgical Scalpel Handle Precision Blade Holder is a scalpels & knife we forge, grind and finish in our own Sialkot workshop. Instruments of this type make precise incisions and separate tissue planes.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$2.00, falling to \$1.40 at 25 pieces — about 30% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28617921

Surgical Sterile Steel Scalpel Blade Box Compatible Handle

\$1.50 per piece

Trade price from **\$0.50**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Surgical Sterile Steel Scalpel Blade Box Compatible Handle — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type make precise incisions and separate tissue planes. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$1.50, falling to \$0.50 at 1 piece — about 67% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70001112****Thomas Uterine Curette 11\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thomas Uterine Curette 11\"/>


Thomas uterine curettes are incredible gynecological instruments used to remove the uterine tissues during surgical procedures. They come with multiple variations, such as a curette with 1/2\", Blunt, Size 0, 1, 2, 3, 4, 5, and 6. The long hollow vertical groove handle extends toward the versatile loop that allows easy grasping and removal of the uterine tissues. These instruments are manufactured from high-grade German surgical stainless steel and thus can be reused after sterilization.

THOMAS UTERINE CURETTES FEATURES

- Fenestrated Teardrop Loop Designed Tip
- Blunt Blades Working Tip From #1 to #6
- Malleable Shaft for Manipulation to Achieve Easy Access.
- Vertical Grooves Hollow Handle Easy to Hold and Maneuver During Procedures.
- Autoclaved, Reusable, and Easy to Clean..
- Single-ended OB/GYN Surgical Instrument.
- Made with Top Quality Stainless Steel

PACKAGE INCLUDE:**SIZE 1**

THOMAS UTERINE CURETTES

**WI-70001111****Thomas Uterine Curette 11\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thomas Uterine Curette 11\"/>


Thomas uterine curettes are incredible gynecological instruments used to remove the uterine tissues during surgical procedures. They come with multiple variations, such as a curette with 1/2\", Blunt, Size 0, 1, 2, 3, 4, 5, and 6. The long hollow vertical groove handle extends toward the versatile loop that allows easy grasping and removal of the uterine tissues. These instruments are manufactured from high-grade German surgical stainless steel and thus can be reused after sterilization.


THOMAS UTERINE CURETTES FEATURES

- Fenestrated Teardrop Loop Designed Tip
- Blunt Blades Working Tip From #1 to #6
- Malleable Shaft for Manipulation to Achieve Easy Access.
- Vertical Grooves Hollow Handle Easy to Hold and Maneuver During Procedures.
- Autoclaved, Reusable, and Easy to Clean..
- Single-ended OB/GYN Surgical Instrument.
- Made with Top Quality Stainless Steel

PACKAGE INCLUDE:**SIZE 2**

THOMAS UTERINE CURETTES



**WI-70001110****Thomas Uterine Curette 11\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thomas Uterine Curette 11\"/>


Thomas uterine curettes are incredible gynecological instruments used to remove the uterine tissues during surgical procedures. They come with multiple variations, such as a curette with 1/2\", Blunt, Size 0, 1, 2, 3, 4, 5, and 6. The long hollow vertical groove handle extends toward the versatile loop that allows easy grasping and removal of the uterine tissues. These instruments are manufactured from high-grade German surgical stainless steel and thus can be reused after sterilization.

THOMAS UTERINE CURETTES FEATURES

- Fenestrated Teardrop Loop Designed Tip
- Blunt Blades Working Tip From #1 to #6
- Malleable Shaft for Manipulation to Achieve Easy Access.
- Vertical Grooves Hollow Handle Easy to Hold and Maneuver During Procedures.
- Autoclaved, Reusable, and Easy to Clean..
- Single-ended OB/GYN Surgical Instrument.
- Made with Top Quality Stainless Steel

PACKAGE INCLUDE:**SIZE
3**

THOMAS UTERINE CURETTES

**SIZE
4****WI-70001109****Thomas Uterine Curette 11\"/>****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thomas Uterine Curette 11\"/>


Thomas uterine curettes are incredible gynecological instruments used to remove the uterine tissues during surgical procedures. They come with multiple variations, such as a curette with 1/2\", Blunt, Size 0, 1, 2, 3, 4, 5, and 6. The long hollow vertical groove handle extends toward the versatile loop that allows easy grasping and removal of the uterine tissues. These instruments are manufactured from high-grade German surgical stainless steel and thus can be reused after sterilization.

THOMAS UTERINE CURETTES FEATURES


- Fenestrated Teardrop Loop Designed Tip
- Blunt Blades Working Tip From #1 to #6
- Malleable Shaft for Manipulation to Achieve Easy Access.
- Vertical Grooves Hollow Handle Easy to Hold and Maneuver During Procedures.
- Autoclaved, Reusable, and Easy to Clean..
- Single-ended OB/GYN Surgical Instrument.
- Made with Top Quality Stainless Steel


PACKAGE INCLUDE:

THOMAS UTERINE CURETTES



THOMAS UTERINE CURETTES



**WI-70001108****Thomas Uterine Curette 11" (27.9cm) Blunt Blade Size #5****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thomas Uterine Curette 11" (27.9cm)
Blunt Blade Size #5 on Malleable Shank
Oval Shaped Fenestrated Loop Tip
Single-Ended Stainless Steel
OB/Gynecology Surgical Instrument



Thomas uterine curettes are incredible gynecological instruments used to remove the uterine tissues during surgical procedures. They come with multiple variations, such as a curette with 1/2", Blunt, Size 0, 1, 2, 3, 4, 5, and 6. The long hollow vertical groove handle extends toward the versatile loop that allows easy grasping and removal of the uterine tissues. These instruments are manufactured from high-grade German surgical stainless steel and thus can be reused after sterilization.

THOMAS UTERINE CURETTES FEATURES

- Fenestrated Teardrop Loop Designed Tip
- Blunt Blades Working Tip From #1 to #6
- Malleable Shaft for Manipulation to Achieve Easy Access.
- Vertical Grooves Hollow Handle Easy to Hold and Maneuver During Procedures.
- Autoclaved, Reusable, and Easy to Clean..
- Single-ended OB/GYN Surgical Instrument.
- Made with Top Quality Stainless Steel

PACKAGE INCLUDE:**SIZE 5**

THOMAS UTERINE CURETTES

**WI-70001107****Thomas Uterine Curette 11" (27.9cm) Blunt Blade Size #6****\$8.99** per pieceTrade price from **\$8.27**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Thomas Uterine Curette 11" (27.9cm)
Blunt Blade Size #6 on Malleable Shank
Oval Shaped Fenestrated Loop Tip
Single-Ended Stainless Steel
OB/Gynecology Surgical Instrument


Thomas uterine curettes are incredible gynecological instruments used to remove the uterine tissues during surgical procedures. They come with multiple variations, such as a curette with 1/2", Blunt, Size 0, 1, 2, 3, 4, 5, and 6. The long hollow vertical groove handle extends toward the versatile loop that allows easy grasping and removal of the uterine tissues. These instruments are manufactured from high-grade German surgical stainless steel and thus can be reused after sterilization.

THOMAS UTERINE CURETTES FEATURES

- Fenestrated Teardrop Loop Designed Tip
- Blunt Blades Working Tip From #1 to #6
- Malleable Shaft for Manipulation to Achieve Easy Access.
- Vertical Grooves Hollow Handle Easy to Hold and Maneuver During Procedures.
- Autoclaved, Reusable, and Easy to Clean..
- Single-ended OB/GYN Surgical Instrument.
- Made with Top Quality Stainless Steel

PACKAGE INCLUDE:**SIZE 6**

THOMAS UTERINE CURETTES





WI-40141624

Waterproof Stainless Steel Luxury Battery Powered Curved Blade

\$2.73 per piece

Trade price from **\$1.50**
Minimum order 5 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Waterproof Stainless Steel Luxury Battery Powered Curved Blade is part of our surgical instruments range, produced start to finish at our Sialkot facility. Instruments of this type make precise incisions and separate tissue planes. The name tells you what you are getting: a curved working end for reach around structures.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$2.73, falling to \$1.50 at 5 pieces — about 45% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003428

Wieder Tongue Depressor Large Blade 36mm Wide Serrad

\$10.99 per piece

Trade price from **\$10.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Wieder Tongue Depressor Large Blade 36mm Wide Serrated – Double Ended 9.5" Long Surgical Dental Oral Instrument Ergonomic L- Shaped Curved Handle Stainless Steel Tool for Dentists, Oral Surgeons, ENT and Maxillofacial Surgeons and Dental Hygienists.


surgeons to pass suction tubes. As a result, they can absorb fluids that may block the view of the surgical field. The ergonomically designed has a bottom curve for resting the hand, which avoids slippage. Manufactured in premium grade stainless steel to ensure the highest quality and product performance during surgical procedures.


WIEDER TONGUE DEPRESSOR FEATURES

- Crafted With 9.5" Overall Length.
- Inner Serrations to Prevent Mucosal Damage.
- 28 & 36mm Blade Size to Suit All Applications From Child to Adult.
- 3 Small Fenestrations to Drain Fluids.

Wieder is a specialized device that dentists, ENT and maxillofacial surgeons can use to manipulate, press down and protect the tongue, especially during procedures that require an optimal view of the oral cavity. This tongue depressor features an L-shaped curved blade that accommodates to the jawbone's contour with ease.

Moreover with inner serration, blade has three fenestrations that allow



**WI-70003473****Wieder Tongue Depressor Small Blade 28mm Wide Serrated****\$9.99** per pieceTrade price from **\$9.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Wieder Tongue Depressor Small Blade 28mm Wide Serrated – Double Ended 9.5" Long Surgical Dental Oral Instrument Ergonomic L- Shaped Curved Handle Stainless Steel Tool for Dentists, Oral Surgeons, ENT and Maxillofacial Surgeons and Dental Hygienists.

Wieder is a specialized device that dentists, ENT and maxillofacial surgeons can use to manipulate, press down and protect the tongue, especially during procedures that require an optimal view of the oral cavity. This tongue depressor features an L-shaped curved blade that accommodates to the jawbone's contour with ease.

Moreover with inner serration, blade has three fenestrations that allow

surgeons to pass suction tubes. As a result, they can absorb fluids that may block the view of the surgical field. The ergonomically designed has a bottom curve for resting the hand, which avoids slippage. Manufactured in premium grade stainless steel to ensure the highest quality and product performance during surgical procedures.

WIEDER TONGUE DEPRESSOR FEATURES

- Crafted With 9.5" Overall Length.
- Inner Serrations to Prevent Mucosal Damage.
- 28 & 36mm Blade Size to Suit All Applications From Child to Adult.
- 3 Small Fenestrations to Drain Fluids.



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

Request a quotation

Send us product codes from this catalogue, quantities and your country.
We reply the same working day on WhatsApp.

Online shop	shop.waqarintl.com
Email	sales@waqarintl.com
WhatsApp /	
Phone	+92 307 6173627
Works	Roras Road, Opposite Anwar Khawaja Industries, Sialkot 51310, Pakistan

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.