





WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

PRODUCT CATALOGUE · OCTOBER 2026

Instrument Sets Veterinary Instruments

Reusable stainless-steel instrument sets forged, ground and finished in Sialkot, Pakistan, for distributors, hospitals, clinics and OEM partners worldwide.



114

PRODUCTS

\$4.00

PRICES FROM

OEM

PRIVATE LABEL

Global

WORLDWIDE
EXPORT

Contents

114 products, A to Z. Prices in USD per piece; trade prices for minimum order quantities are shown on each product.

Product	Code	Price	Page
100 Pcs Scalpel Blades #10 for Handle No. 3	WI-70003728	\$11.99	5
100 Pcs Scalpel Blades #20 for Handle No. 4	WI-70003727	\$11.99	5
11 Pcs Advanced Pig Castration Set with Carrying Case	WI-70001379	\$32.99	6
12 Pcs Small Animal Dental Extraction Kit	WI-70001269	\$79.99	6
12 Pcs Small Animal Dental Extraction Set with Sterilization Cassette	WI-70003725	\$109.99	7
12 Pcs Spay Snook Hook Set 8" Curved 5-6mm Wide Blunt Tip Flat Handle	WI-70000587	\$65.99	7
13 Pcs Advance Dental Prophy Kit	WI-70003724	\$79.99	8
15 Pcs Birds Dissection Kit	WI-70004108	\$29.99	8
15 Pcs Birds Dissection Kit Surgical Veterinary Instruments	WI-70000509	\$29.99	9
15 Pcs Feline Spay Instruments Pack	WI-70004001	\$69.99	9
157 Piece Stainless Steel Reusable Veterinary Instrument Student Kit	WI-28886663	\$70.00	10
160 PCS Stainless Steel Anatomy Dissection	WI-25582029	\$4.00	10
17 Pcs Dog Ear Cropping Clamps Set with Carrying Case	WI-70000725	\$79.99	11
17 Pcs Feline Spay Kit Economy Set	WI-70001380	\$49.99	11
17 Pcs Large Animal Deep Surgery Set	WI-70001378	\$149.99	12
18 Pcs Dog Dental Extraction Forceps Elevator Kit Veterinary Set	WI-70003718	\$94.99	12
18 Pcs Dog Dental Extraction Kit	WI-70004002	\$89.99	13
18 Pcs Economy Feline Spay Kit	WI-70001382	\$43.99	13
2 Pcs Calf Puller OB Chain Hook Handle	WI-70001273	\$19.99	14
2 Pcs Double Crush Emasculators Set 9" & 12"	WI-70001749	\$44.99	14
2 Pcs Ratchet Style Calf Puller with Extra 60" OB Chain Set	WI-70001275	\$179.99	15
20 Pcs Standard Veterinary Orthopedic Instrument Set	WI-28956327	\$700.00	15
20-Piece Veterinary Orthopedic Surgical Kit Stainless Steel Bone	WI-30897559	\$185.00	16
22 Pcs Feline Spay Kit with Stainless Steel Box	WI-70001381	\$64.99	16
24 Pcs Feline Spay Pack Economy Grade	WI-70001923	\$49.99	17
24 Pcs T/C Feline Spay Pack	WI-70001924	\$89.99	17
27 Pcs Male Dog Scrotal Castration – Neuter Kit	WI-70003712	\$89.99	18
3 Equine Dental Float Rasps Set Angle Up, Down & Straight Adjustable Wooden Handle	WI-70000457	\$39.99	18
3 Pcs Cat Nail Clipper Kit with Carrying Case	WI-70004003	\$17.99	19
3 Pcs Dog Nail Clipper Kit	WI-70003715	\$17.99	19
3 Pcs Dog Nail Clipper Kit Guillotine Type	WI-70003714	\$17.99	20
3 Pcs Nylon OB Calf Pulling Straps 30", 45" 60" with Hook Handle Set	WI-70000452	\$43.99	20
3 Pcs Ratchet Calf Puller with OB Handle & 60" Chain Cattle Delivery Veterinary Set	WI-70001276	\$194.99	21
3 Pcs Serra Emasculators 35.5cm with Ratchet	WI-70003683	\$164.99	21
3 Pcs Spay Snook Hook Set 8" Curved 5-6mm Wide Blunt Tip Flat Handle	WI-70000589	\$21.99	22
34 Pcs Entomology Dissection Kit	WI-70003634	\$25.99	22
4 Pcs Calf Puller with 60" OB Chain Handle & Carrying Bag Cattle Delivery Set	WI-70001149	\$229.99	23
4 Pcs Canine Mouth Gag Set 3.5", 4.5", 5.5" & 6.5"	WI-70003720	\$39.99	23
4 Pcs Dog Ear Cropping Clamps Set "C", "I", "S", "L" Type with Teeth	WI-70000726	\$49.99	24
4 Pcs Winged Elevators Set 1mm-4mm Tips Standard Handle	WI-70001402	\$29.99	24
47 Pcs Deluxe Veterinary Kit Spay Surgery Set	WI-70001396	\$139.99	25

Contents (continued)

Product	Code	Price	Page
5 Pcs Covault Spay Hook Set 9" with Curved Probe Tip	WI-70000402	\$39.99	26
5 Pcs Dental Scaling Kit For Large Dogs	WI-70001417	\$21.99	26
5 Pcs Dental Scaling Kit For Medium Dogs	WI-70001418	\$21.99	27
5 Pcs Dental Scaling Kit For Small Dogs	WI-70001420	\$21.99	27
5 Pcs Veterinary Otoscope Set Animal Diagnostic Instruments	WI-70001910	\$22.99	28
6 Pcs Basic Dental Prophy Instruments Kit for Small Animals	WI-70004007	\$39.99	28
6 Pcs Basic Dental Prophy Kit	WI-70003723	\$49.99	29
6 Pcs Basic Exodontia Kit	WI-70004006	\$59.99	29
6 Pcs Dog Ear Suture Kit	WI-70003717	\$35.16	30
6 Pcs Dog Tail Docking Kit	WI-70003713	\$35.99	30
6 Pcs Equine Periodontal Instrument Kit	WI-70003631	\$123.74	31
6 Pcs First Aid Kit	WI-28887630	\$40.00	31
6 Pcs Horse First Aid Kit Veterinary Surgical Instruments	WI-28806200	\$30.00	32
6 Pcs Horse Pets First Aid Kit	WI-23883602	\$160.00	32
6 Pcs Pets First Aid Kit	WI-70004000	\$24.99	33
6 Pcs Spay Snook Hook Set 8" Curved 5-6mm Wide Blunt Tip Flat Handle	WI-70000588	\$39.99	33
6 Pcs Stripping and Carding Knives Set – Blue Wooden Handle	WI-70000172	\$32.99	34
6 Pcs Stripping and Carding Knives Set – Garnish Red Handle	WI-70000173	\$32.99	34
6 Pcs Stripping and Carding Knives Set – Multi Color Handle	WI-70000170	\$32.99	35
6 Pcs Stripping and Carding Knives Set – Red Wooden Handle	WI-70000171	\$32.99	35
7 Pcs Canine Periodontal Pack	WI-70003719	\$25.99	36
7 Pcs Veterinarian Umbilical Cord Scissor Kit	WI-70002745	\$35.99	36
72 PC US Minor Surgery Surgical Veterinary	WI-08756437	\$65.00	37
7-Piece Horse Cattle Farrier	WI-25860394	\$30.00	37
7-Piece Reusable Veterinary Horse Care Farrier	WI-25894370	\$30.00	38
9 Pcs Large Animal Eye Kit	WI-70001373	\$54.99	38
Advanced 30-Piece Dissection Kit Medical Students Veterinary Surgical	WI-23513026	\$7.00	39
Advanced 45-Piece Stainless Steel Dissecting Kit	WI-21725912	\$13.00	39
All-In-One 14-Piece Stainless Steel Minor Surgery Veterinary Dental	WI-05867732	\$12.65	40
Animal Orthopedic Surgical Instrument 1.5/2.0/2.4mm Vet LC-DCP Bone	WI-28812440	\$840.00	40
Basic Starter Pets First Aid Kit	WI-28795310	\$30.00	41
Bitch Spay-Neuter Kit 25 Pcs Spay Surgery Instruments	WI-70003722	\$85.99	41
Botany Dissection Kit for Medical Biology Veterinary Students	WI-29411879	\$8.00	42
Canine Feline Spay Pack Veterinary Surgical Instruments	WI-29308398	\$55.00	42
Cat Spay Neuter Surgical Instrument Set	WI-36187616	\$180.00	43
Complete Orthopedic Stainless Steel Suture Kit for Medical	WI-32658802	\$10.00	43
Complete Suture Practice Kit Stainless Steel Instruments Pre-Cut	WI-32067215	\$7.00	44
Deluxe Veterinary Dissection Kit	WI-29399967	\$15.00	44
Dog Dental Extraction Kit 18Pcs Veterinary Surgical Instruments Set	WI-29552734	\$8.00	45
Economy 18-Piece Canine Spay Veterinary Instrument Manikin Kit	WI-24106098	\$30.00	45
Equine Dental McPherson Horse Mouth Speculum & Float Rasp Kit	WI-70002717	\$124.98	46
Extendable Stainless Steel Snake Catching Holder Stick Instrument	WI-27930440	\$20.00	46

Contents (continued)

Product	Code	Price	Page
Farriers Hoof Grooming Tool Kit 6-Piece	WI-28996408	\$20.00	47
Feline Spay Neuter Instrument Set Stainless Steel Veterinary	WI-36133667	\$180.00	47
Feline Spay Neuter Surgical Pack Stainless Steel Veterinary	WI-42807811	\$180.00	48
Feline Sterilization Surgery Instrument Set	WI-46255664	\$180.00	48
Feline Sterilization Surgical Instrument Set Cat Spay Neuter Veterinary	WI-46244729	\$180.00	49
Hoof Farrier Trimming Tools 8 Pcs Set Reusable Breathable	WI-25896343	\$30.00	49
Horse Farrier Hoof Trimmer Tool Kit	WI-25915247	\$8.00	50
Horse Hoof Trimming Farrier Tools	WI-25935362	\$30.00	50
Large Animal Veterinary First Aid Kit 6 Pcs	WI-23955371	\$220.00	51
Pro Eco-Friendly Reusable Stainless Steel Complete Farrier	WI-28999003	\$40.00	51
Pro Horse Hoof Trimming Farrier Tool	WI-25865404	\$30.00	52
Puppy Tail Banding Elastrator with FREE 100 Pcs Bands	WI-70003711	\$18.69	52
Rabbit Rat Rodent Small Animal Cheek Teeth Molar Incisor	WI-29429154	\$15.00	53
Reusable Eco-Friendly Veterinary Orthopedic Instrument Kit Complete System Dog	WI-34073285	\$600.00	53
Reusable Veterinary Orthopedic Instrument Kit Complete System Small Animal	WI-34051245	\$600.00	54
Shorty Hog Holder Instrument Set Veterinary Pig Restraining	WI-27932199	\$18.00	54
Small Animal Veterinary Dental Instrument Set	WI-28939165	\$35.00	55
Stainless Steel Anatomical Model Orthopedics Veterinary	WI-24082237	\$30.00	55
Stainless Steel Feline Spay Neuter Instrument Set for Cat	WI-46260656	\$180.00	56
Stainless Steel Hoof Clincher Instrument Set for Horseshoe	WI-25892288	\$8.00	56
Stainless Steel Orthopedic Surgical Instrument Set Blue Silicone	WI-36641610	\$55.00	57
Stainless Steel Surgical Skin Hooks 15cm Instrument	WI-27567272	\$7.00	57
Suture Practice Kit for Medical Veterinary Students	WI-32054322	\$7.00	58
Suture Practice Kit Medical Veterinary Students	WI-28543634	\$12.00	58
Veterinary Cardiovascular Surgery Anatomical Model Kit-Stainless Steel Medical	WI-24092177	\$30.00	59
Veterinary Dissecting Kit Anatomical Model 14 Pcs	WI-28634030	\$50.00	59
Veterinary Orthopedic Trauma PUSM Elastic Nail Instrument	WI-34691746	\$200.00	60
Veterinary Pig Holder Instrument Set	WI-27906596	\$15.00	60
Veterinary Surgery Horse Spay Pack with Instrument Box	WI-23908573	\$50.00	61
ZebraFish Veterinary Dissecting Kit Surgical Instruments Set	WI-29236743	\$9.00	61

**WI-70003728****100 Pcs Scalpel Blades #10 for Handle No. 3****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Scalpel Blade No. 10 has a curved belly with sharpened edge. This blade is used in making skin incisions. Pack of 100 Pcs Scalpel Blades #3 ideal for ENT Medical Dental Podiatry, Crafts, Dermaplaning & More

SCALPEL BLADE #10 FEATURES

- The Pack of 100 #10 Scalpel Blades is an ideal choice to restock on blades for medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #3 that can be used for a variety of different blade sizes.
- The #10 Scalpel Blades are forged from corrosion-resistant and ultrasonic cleaned carbon steel for

the highest quality + value for the price.

- The sharp and precise blades make these Scalpel Blades a great buy for Designers, Hobbyists, Graphic Artists, Model Makers & more.
- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #10 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries including health and wellness, beauty and skincare, mechanical and tech servicing,

**WI-70003727****100 Pcs Scalpel Blades #20 for Handle No. 4****\$11.99** per pieceTrade price from **\$11.03**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 100 pcs Scalpel Blades #20. The No.20 is a large version of the No.10 blade with a curved cutting edge and a flat, unsharpened back edge. Used for orthopaedic and general surgical procedures it fits handle No. 4, 4L.


SCALPEL BLADE #20 FEATURES

- The Pack of 100 #20 Scalpel Blades is an ideal choice to restock on blades for veterinary, medical, beauty & home uses. The precision-shaped carbon steel blades are sharp and easy to use for a variety of different tasks and/or procedures. Includes a free Handle #4 that can be used for a variety of different blade sizes.
- The #20 Scalpel Blades are forged from corrosion-resistant and

ultrasonic cleaned carbon steel for the highest quality + value for the price.

The sharp and precise blades make these Scalpel Blades a great buy for Designers, Hobbyists, Graphic Artists, Model Makers & more.

- These medical disposables are ready to use and cleaned with gamma radiation. They can be used professionally for dermaplaning to dull layers of skin are removed from the face in order to reveal a brighter complexion.
- The #20 scalpel blades are multi-functional tools which can be used for a wide range of tasks in industries including health and wellness, beauty and skincare, mechanical and tech servicing,



**WI-70001379****11 Pcs Advanced Pig Castration Set with Carrying Case****\$32.99** per pieceTrade price from **\$30.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



Advanced Pig Castration Set – 11 Pcs
Advanced Pig Swine Castration
Neutering Dissection Kit Set German
Tools Veterinary Instruments with Free
Carrying Case.

Surgical castration of male piglets has been a traditional measure for centuries to avoid boar taint, the sex-specific off odor of pork, and to reduce problems due to male-specific behavior.

Castration is essentially universal in pigs; only a few male pigs are chosen as potential breeder boars. One of the primary reasons is behavioral, with age, body size and sexual maturity boars tend to be more aggressive with pen mates and become more difficult to handle than barrows of similar age and weight. The other and perhaps most important reason is that meat from

boars that are nearing sexual maturity (about 200 lbs. body weight) has a high potential for an odor and flavor problem commonly called “boar taint.” that many consumers detect in cooked pork from intact males.

Most of the veterinarians acknowledge that castration helps control aggressive behavior in pigs and recommends that the procedure. Studies have also indicated that pigs castrated at 2 weeks of age displayed fewer indicators of the pain than pigs castrated at 7 weeks of age. Our advanced Pigs Castration kit has all the tools to perform the Castration (De-ballification). All tools are made of surgical grade stainless steel used by professionals.

ADVANCED PIG CASTRATION SET INCLUDES**WI-70001269****12 Pcs Small Animal Dental Extraction Kit****\$79.99** per pieceTrade price from **\$73.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Small Animal Extraction Kit With FREE Scalpel Blades German Grade Stainless Steel Dental Veterinary Instruments Set

Small animals dental extraction is a common procedure in veterinary practice. Extraction can be Non-Surgical (due to some periodontal disease) or Surgical (due to damaged or abnormal teeth extraction). Our Dental Extraction kit contains all the necessary tools for both types of extraction.

Teeth that are extracted in small animal patients may or may not be mobile. This 12 Pcs small animals dental extraction kit has all appropriate instruments to perform selective or full mouth dental extraction. All Tools are

made from Surgical Grade stainless steel and used by professionals.

NOTE: We also make custom design surgery kit to meet with specific needs.

DENTAL EXTRACTION SET INCLUDES:

- 1 TC Mayo Hager Needle Holder 16cm
- 1 T/C Iris Scissors 11.5cm Straight
- 1 Adson Brown Forceps 7×7 Teeth
- 1 Extraction Forceps Miltex DEF E (FIG. #217)
- 1 Scalpel Handle #3
- 1 Periosteal Elevator 2/4
- 1 Pc Winged Elevator 1mm – Standard Handle
- 1 Pc Winged Elevator 2mm – Standard Handle



WI-70003725

12 Pcs Small Animal Dental Extraction Set with Sterilization Cassette

\$109.99 per piece

Trade price from **\$101.19**
Minimum order 21 pieces


- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

sSmall Animal Dental Extraction Kit of 12 Pcs with Sterilization Cassette for tooth removal, Animal Dental Extraction Set Veterinary Instruments.

Waqar International Dental Extraction kit contains all the necessary tools for both types of extraction. Teeth that are extracted in small animal patients may or may not be mobile. This 12 Pcs dental extraction kit has all appropriate instruments to perform selective or full mouth dental extraction. All Tools are made from Surgical Grade stainless steel and used by professionals. We also make custom design surgery kit to meet with specific needs.

DENTAL EXTRACTION SET INCLUDES:

- 1 TC Mayo Hager Needle Holder 16cm
- 1 Iris Scissors 11.5cm TC Straight
- 1 Adson Brown Forceps 7x7
- 1 Pc Extraction Forceps Miltex DEF E (FIG. #217)
- 1 Scalpel Handle #3
- 1 Periosteal Elevator 2/4
- 1 Winged Elevator Stubby Handle 1mm
- 1 Winged Elevator Stubby Handle 2mm
- 1 Winged Elevator Stubby Handle 3mm
- 1 Winged Elevator Stubby Handle 4mm
- 1 Winged Elevator Stubby Handle 5mm



WI-70000587

12 Pcs Spay Snook Hook Set 8" Curved 5-6mm Wide Blunt Tip Flat Handle

\$65.99 per piece

Trade price from **\$60.71**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 12 Pcs Spay Snook Hook Set 8" (20.3cm) Curved 5-6mm Wide Blunt Tip Premium Quality Stainless Steel Ovariectomy Veterinary Surgical Instrument.

Spay Snook Hook is a specialized surgical tool used largely in veterinary surgeries for neutering dogs and cats. Ovarian Snook Hook designed for maximum performance and longevity, assists veterinarians in the safe and precise removal of reproductive organs. Snook hook features a curved blunt tip, long shank and a flat solid handle for gripping. The hook-shaped tip instrument is also used to retrieve the horn of the uterus while performing ovariectomy in felines and canines. This veterinary spay hook made of german-grade stainless steel,


which keeps it rust-free and autoclavable for easy sterilization and can be reused.

SPAY SNOOK HOOK FEATURES

- Crafted with 8" (20.3cm) Overall Length.
- 5-6mm Wide Curved Blunt Tip.
- Hook Shaped Working End.
- Long Straight Shank.
- Flat Smooth Grip Handle.
- Used in Neutering Dogs and Cats During Ovariectomy Procedure.
- Reusable.
- Made with Premium Quality Stainless Steel;

PACKAGE INCLUDE

- 12 Pcs Spay Snook Hook 8" (20.3cm) Curved Blunt Tip



**WI-70003724****13 Pcs Advance Dental Propy Kit****\$79.99** per pieceTrade price from **\$73.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

13 Pcs Advance Dental Propy Kit Premium Quality Stainless Steel Oral Hygien Prophylaxis Procedure Dental Instruments.



Dental prophyl treatment can help the dentist to diagnose any problems at an early stage. Dental prophyl or prophylaxis procedure promotes oral health. There are many dental problems that may not show early visible signs or cause pain. In a dental prophyl procedure, a dentist or a dental hygienist examines the gum tissue, perform a periodontal probing to measure the depth of the gum tissue, examine any decaying teeth that need attention and perform a thorough cleaning. Our 13 Pcs Advance Propy Kit contains all the necessary tools to perform prophylaxis. All instruments

are made from medical-grade stainless steel used by professionals ensure durability and precision performance.

ADVANCE DENTAL PROPHY KIT INCLUDES

Our 13 pcs advance dental prophylaxis set comes with the following instruments;

- 1 Mouth Mirror Handle with Mirrors
- 1 Periodontal Probe CP15/23
- 1 Explorer 11/12 ODU
- 1 Towner-Jacquette Scaler U15/30
- 1 Columbia Curette 13/14
- 1 Sickle Scaler 204S
- 1 Hygienist Scaler H6/H7
- 1 Gracy Curette 1/2
- 1 Gracy Curette 7/8
- 1 Gracy Curette 11/12

**WI-70004108****15 Pcs Birds Dissection Kit****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Bird Dissection Kit – 15 Pcs Bird Anatomy Dissecting Set Biology Medical Veterinary Student Stainless Veterinary Lab Instruments Tools Set

The birds dissection Instruments kit is used for the dissection of birds. This kit has all the necessary instruments to perform the External and Internal anatomy of a bird. All instruments made with high quality stainless steel so students can perform dissection precisely.

External Anatomy of Bird

The dissection kit helpful in below external anatomy procedure;

- To examine an external feature of the bird-like beak, cere, eye, ear, primary wings, hallux, Tarsus.

- To measure the length, from beak to tip of tail feathers, in cm of the bird.
- To examine primary feathers, which are attached to the bones of the "fingers".
- To examine the external Muscular System.

Internal Anatomy of Bird

The dissection kit helpful in below internal anatomy procedure;

- Examine internal organs like esophagus, trachea, crop, liver, gizzard, pancreas, intestine, cloaca, testes, ovary and heart
- Renal and Urinary System
- Examine the Reproductive System. Male birds have two white testes while the female bird has an ovary.



**WI-70000509****15 Pcs Birds Dissection Kit Surgical Veterinary Instruments****\$29.99** per pieceTrade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Birds Dissection 15 Pieces Kit with FREE Carrying Case Top Quality Stainless Steel Birds Anatomy Biology Medical Veterinary Lab Instruments Tools Set

The Birds Dissection kit is a high-grade stainless steel tool used for dissectioning birds, allowing for examination of external features like beak, eye, ear, wings, and tarsus, as well as internal organs like esophagus, trachea, liver, pancreas, intestine, testes, ovary, heart, renal and urinary system, and reproductive system. All instruments made with high quality stainless steel so students can perform dissection precisely.

- 1 Teaser Needle Angled
- 1 Iris Scissors 11.5cm
- 1 Dressing Forceps 16cm
- 1 Scalpel Handle #3
- 1 Adson Forceps Dressing 12cm
- 1 Adson Forceps 1x2 Teeth 12cm
- 1 Micro Dissection Forceps 10cm
- 1 Tweezers with Magnifying Glass
- 2 Scalpel Blades #10
- 2 Scalpel Blades #11
- 1 Mosquito Forceps 12.5cm
- 1 Ruler 6" Stainless Steel

With **FREE** Carrying Case**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Veterinary Dissection Kit**APPLICATION:** Bird Anatomy**15 PCS BIRDS DISSECTION KIT INCLUDES**

- 1 Teaser Needle Straight

**WI-70004001****15 Pcs Feline Spay Instruments Pack****\$69.99** per pieceTrade price from **\$64.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


15 Pcs Feline Spay Pack Neuter Surgery Kit Surgical Veterinary Instruments

Feline Spay/Neuter Pack includes all the surgical instruments that veterinarians require for performing a cat spay. Canine and Feline Spay Neuter Pack Kit is essential during the spay and neutering procedure. This kit is comprised of all instruments to fulfill spaying procedures. All instruments are crafted from 100% premium German grade surgical stainless material and thus can be reused after sterilization

- 1 Snook Hook
- 1 Thumb Dressing Forceps
- 1 Adson Brown Forceps 4.75" (12cm)
- 2 TC Olsen Hager Needle Holder 5.5" (14cm)
- 1 TC Mayo Scissor 6.75" (17.12cm)
- 4 Backus Towel Clamp 3.5" (8.89cm)
- 1 Scalpel Handle #3
- **FREE** Scalpel Blades

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Surgical Veterinary Instruments**APPLICATION:** Spaying Procedure**RUST PREVENTION****PROCEDURE:** Passivated**REUSABLE:** YES**AUTOCCLAVED:** Yes**ULTRASONIC CLEANED:** Yes**FELINE SPAYING PACK INCLUDES**

- 2 Mosquito Forceps 5" (12.7cm) Curved
- 2 Mosquito Forceps 5" (12.7cm) Straight





WI-28886663

157 Piece Stainless Steel Reusable Veterinary Instrument Student Kit

\$70.00 per piece

Trade price from \$25.00
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

157 Piece Stainless Steel Reusable Veterinary Instrument Student Kit is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Supplied as a set — the contents can be adjusted to your list. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$70.00, falling to \$25.00 at 10 sets — about 64% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-25582029

160 PCS Stainless Steel Anatomy Dissection

\$4.00 per piece

Trade price from \$3.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

160 PCS Stainless Steel Anatomy Dissection is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$4.00, falling to \$3.00 at 1 set — about 25% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70000725****17 Pcs Dog Ear Cropping Clamps Set with Carrying Case****\$79.99** per pieceTrade price from **\$73.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

17 Pcs Dog Ear Cropping Clamps Set with FREE Carrying Case- Stainless Steel Professional Veterinarian Ear Trim Clamps Retractor for Various Dog Breeds Livestock Medical Surgical Veterinary Instruments Kit.

Dog ear cropping clamps set ensuring precise cuts for an aesthetically pleasing appearance. The premium-grade clamps provides comfort for the dog and streamlines the process for veterinarians, leading to a smoother recovery. Ear cropping instruments have sharp cutting edges allow clean and accurate incisions, ensuring uniform and precise cuts, resulting in aesthetically pleasing ear shapes that meet breed standards. Ergonomic handles provide a comfortable grip, reducing hand fatigue during prolonged

procedures. Premium quality stainless steel construction ensure lasting performance to meet the needs of veterinary professionals, making the ear-cropping process safer and more efficient. All these instruments are fully autoclavable, rust free and can be reused over time.

DOG EAR CROPPING CLAMPS SET INCLUDES:

- 1 Ear Cropping Clamp 5.12" (13cm) – C Type
- 1 Ear Cropping Clamp 5.5" (14cm) – I Type
- 1 Ear Cropping Clamp 7" (18cm) – S Type
- 1 Ear Cropping Clamp 7" (18cm) – L Type

**WI-70001380****17 Pcs Feline Spay Kit Economy Set****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Feline Spay Economy Set 17 Pieces Kit Stainless Steel Surgical Veterinary Spay Surgery Instruments Set

In female cats, desexing (also known as spaying) eliminates the risk of pyometra (infection of the womb) mammary tumors and eliminates the risk of tumors in the ovaries, uterus, and cervix. It also prevents other medical conditions. Feline Spay/Neuter Pack includes all the surgical instruments that veterinarians require for performing a cat spay.

The Feline Spay Neuter Pack Kit is essential during the spay and neutering procedure. This kit is comprised of all instruments to fulfill spaying procedures and for safe and effective spay surgery.. All instruments are crafted from 100% premium grade

surgical stainless material and thus can be reused after sterilization.

FELINE SPAY KIT INCLUDES

Our 17 Pcs economy feline spay kit contains all the necessary instruments used in spay surgery.

- 1 Mayo Hager Needle Holder 16cm
- 1 Metzenbaum Scissors Straight 15cm
- 2 Kelly Forceps Curved 14cm
- 2 Mosquito Forceps Curved 12.5cm
- 1 Scalpel Handle #3
- 4 Scalpel Blades #10
- 1 Snook Hook
- 3 Backhaus Forceps 11cm
- 1 Adson Brown Forceps 12cm 7x7 Teeth
- 1 Dressing Forceps 16cm





WI-70001378

17 Pcs Large Animal Deep Surgery Set

\$149.99 per piece

Trade price from **\$137.99**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality 17 Pcs Large Animal Deep Surgery Kit with Free Scalpel Blades Surgical Veterinary Instruments Set

Large animals deep soft tissue surgery includes surgery of the abdomen, urogenital system, respiratory system, from a simple lump removal to complex airway or colic surgery. Waqar International's deep tissue surgery kit is designed for large animals. The kit has extra large instruments to reach deep.

All instruments are premium quality manufactured from surgical grade stainless steel.

LARGE ANIMAL DEEP SURGERY KIT INCLUDES

- 1 Mayo Hegar Needle Holder TC 25.5cm

- 1 Allis Tissue Forceps 5×6 Teeth 25.5cm
- 4 Backhaus Towel Forceps 13cm
- 2 Rochester-Carmalt Forceps Straight 20cm
- 2 Rochester-Carmalt Forceps Curved 20cm
- 2 Kelly Forceps 25.5cm Curved
- 1 Metzenbaum Scissors 28cm Curved
- 1 Operating Scissors 18cm Straight S/B
- 1 3L Scalpel Handle Long 21cm
- 1 4L Scalpel Handle Long 21cm
- 1 Dressing Thumb Forceps 25.5cm

with **FREE** Scalpel Blades.

TECHNICAL INFORMATION

MATERIAL: Stainless Steel



WI-70003718

18 Pcs Dog Dental Extraction Forceps Elevator Kit Veterinary Set

\$94.99 per piece

Trade price from **\$87.39**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 18 Pcs Dog Dental Extraction Kit Stainless Steel Extracting Forceps Elevators Surgical Veterinary Instrument Set

Tooth extraction is a common dental procedure performed in most veterinary practices. This 18 Pcs dog dental extraction kit has all appropriate instruments to perform selective or full mouth dental extraction. Dental Extraction can be Non-Surgical (due to some periodontal disease) or Surgical (due to damaged or abnormal teeth extraction). All Tools are made from Surgical Grade stainless steel and used by veterinarian professionals.

DOG DENTAL EXTRACTION KIT INCLUDES

- 1 Pc Scalpel Handle #3

- 1 Pc TC Olsen Hager Needle Holder 15cm
- 1 Pc TC Mayo Scissor 15cm
- 1 Pc LaGrange Scissor 11cm
- 1 Pc Molt Periosteal Elevator 16cm
- 1 Pc Mosquito Forceps Curved 12.5cm
- 1 Pc Mosquito Forceps Straight 12.5cm
- 1 Pc College Cotton Tweezer 15cm
- 1 Pc Adson Tissue Forceps 1×2 Teeth 12cm
- 1 Pc Adson Tissue Forceps 7×7 Teeth 12cm
- 1 Pc Elevator 3mm Straight
- 1 Pc Elevator 3mm Down Angled
- 1 Pc Elevator 5mm Straight
- 1 Pc Elevator 5mm Down Angled



**WI-70004002****18 Pcs Dog Dental Extraction Kit****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Small Animal Dentistry Tools Set – 18 Pcs Dog Dental Extraction Kit Stainless Steel Extracting Forceps Elevators Surgical Veterinary Instrument Set

Tooth extraction is a common dental procedure performed in most veterinary practices. This 18 Pcs dog dental extraction kit has all appropriate instruments to perform selective or full mouth dental extraction. Dental Extraction can be Non-Surgical (due to some periodontal disease) or Surgical (due to damaged or abnormal teeth extraction). All Dental Tools are made from Surgical Grade stainless steel and used by veterinarian professionals.

DOG DENTAL EXTRACTION KIT INCLUDES

- 1 Pc Scalpel Handle #3

- 1 Pc TC Olsen Hager Needle Holder 15cm
- 1 Pc TC Mayo Scissor 15cm
- 1 Pc LaGrange Scissor 11cm
- 1 Pc Molt Periosteal Elevator 16cm
- 1 Pc Mosquito Forceps Curved 12.5cm
- 1 Pc Mosquito Forceps Straight 12.5cm
- 1 Pc College Cotton Tweezer 15cm
- 1 Pc Adson Tissue Forceps 1x2 Teeth 12cm
- 1 Pc Adson Tissue Forceps 7x7 Teeth 12cm
- 1 Pc Elevator 3mm Straight
- 1 Pc Elevator 3mm Down Angled
- 1 Pc Elevator 5mm Straight
- 1 Pc Elevator 5mm Down Angled

**WI-70001382****18 Pcs Economy Feline Spay Kit****\$43.99** per pieceTrade price from **\$40.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Economy Feline Spay Kit 18 Pieces Surgical Veterinary Set Stainless Steel Feline Spay Surgery Instruments Pack with FREE Blades.


Feline Spay/Neuter Pack includes all the surgical instruments that veterinarians require for performing a cat spay. Canine and Feline Spay Neuter Pack Kit is essential during the spay and neutering procedure. This kit is comprised of all instruments to fulfill spaying procedures. All instruments are crafted from 100% premium quality surgical stainless material and thus can be reused after sterilization

ECONOMY FELINE SPAY KIT INCLUDES

- 1 Mayo Scissors Straight 6.75"
- 1 Metzenbaum Scissors Curved 5.5"

- 1 Olsen-Hegar Needle Holder 5.70"
- 1 Allis Tissue Forceps 4x5 5.90"
- 2 Mosquito Forceps Curved 5"
- 2 Kelly Forceps Curved 5.5"
- 2 Rochester Carmalt Forceps Curved 6.25"
- 1 Snook Hook
- 1 Tissue Forceps 6.29"
- 1 Adson-Brown Forceps 4.75" 7x7 Teeth
- 1 Scalpel Handle #3
- 4 Backhaus Towel Clamps 4.3"

with **FREE** Scalpel Blade**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Surgical Veterinary Instruments**APPLICATION:** Spaying Procedure

**WI-70001273****2 Pcs Calf Puller OB Chain Hook Handle****\$19.99** per pieceTrade price from **\$18.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Highly Durable 2 Pcs Calf Puller OB Hook Handle Set with Blue PVC Grip for Cow, Goat, and Cattle Delivery Instruments on Dairy Farms.

Obstetrical Chain (OB) Hook Handle is a specialized instrument used in veterinary obstetrical procedures. This handle features a hook that perfectly holds the OB chain. OB Chain hook handle is an ideal instrument for birth assistance in calves and other cattle animals. This hook handle is crafted from German grade stainless and features blue pvc grip to minimize the risk of slippage during use. You can use these hooks with 30", 45" & 60" long OB Apparatus. Robust and strongly constructed, these OB chains hook handles will not rust, providing

years of service. Use with the OB chain handle for ease and grip while calving.

OB CHAIN HANDLE FEATURES

- 6" Long and 5" Wide.
- Anti-slip Pvc Grip.
- Compatible With 30", 45" & 60" Long Obstetrical Apparatus
- Ideal for Obstetrical Use
- Resuable
- High Quality Comfortable Fit for Any Delivery.
- Constructed With Stainless Steel Cutting Heads

PACKAGE INCLUDE

- 2 Pcs Calf Puller OB Chain Hook Handle

TECHNICAL INFORMATION**IDEALSIZE**
FITS ANY HAND**DURABLE POLYCOATED**
STAINLESS STEEL HANDLE**WI-70001749****2 Pcs Double Crush Emasculators Set 9" & 12"****\$44.99** per pieceTrade price from **\$41.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Set of 2 Double Crush Emasculators 9" and 12" – White's Style – Livestock, Sheep, Dog Tails Docking Castration Plier Stainless Steel Surgical Veterinary Instruments.

The double crush emasculators crush the cord and blood vessels followed by a sharp cut, reducing post-op bleeding and infection. The Removable wing nut simplifies disassembly for cleaning, and serves as adjustment for tightening or approximating the blades. Made with top quality stainless steel.

DOUBLE CRUSH EMASCULATOR FEATURES:

- White Emasculators Double Crush Jaw
- Less Trauma
- Serrated Edge Blades

- Tension Adjuster
- Smooth Grip Plier Handle
- Crushes the Arteries Severing them Cleanly.
- Minimizes the Risk of Serious Bleeding.
- No Blood Loss Because Blood Vessels are Crushed.
- Made With High Quality Stainless Steel

PACKAGE INCLUDES:

- 1 Pc Double Crush Emasculator 9" (22.9cm) – White's Style
- 1 Pc Double Crush Emasculator 12" (30.5cm) – White's Style

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Emasculator – Veterinary Instruments**DOUBLE CRUSH EMASCULATOR**



WI-70001275

2 Pcs Ratchet Style Calf Puller with Extra 60" OB Chain Set

\$179.99 per piece

Trade price from \$165.59
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
 - ✓ OEM / private-label marking
 - ✓ Reusable, autoclavable
- Set of 2 Ratchet Style Calf Puller Heavy Duty Calving Cattle Birthing Extractor 65" Standard Rod with 60 Inches Long OB Chain Veterinary Surgical Cattle Farming Instruments
- High Quality Stainless Steel Jack Handle
 - Bree chain Handle Length: 25"
 - Calf Puller Total Length: 65" (165 cm)


RATCHET STYLE CALF PULLER BENEFIT

2 PCS CALF PULLER CATTLE DELIVERY SET INCLUDED

- 1 Ratchet Style Calf Puller 65" Set
- 1 Pc Calf Puller OB Chain 60" Long

CALF PULLER FEATURES:

- Aero Space Heat Treated Steel Alloy Jack
- State-of-the-Art Handle Release
- Hardened Steel Cams
- Precision Two-Section Rod Assembly
- Rotating Jack and Bree chain
- Corrosion resistant, lightweight, durable for longer life.
- One person operation permitting a smooth release on both sets of cams.
- Guards against slipping and provides uniform tension and minimal rod damage.
- Smooth for easy cleaning and disinfecting – Convenient to carry and assemble.
- For proper positioning of the Breechain below pin bones



WI-28956327

20 Pcs Standard Veterinary Orthopedic Instrument Set

\$700.00 per piece

Trade price from \$480.00
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ OEM / private-label marking
- ✓ Reusable, autoclavable

ABOUT THIS INSTRUMENT

20 Pcs Standard Veterinary Orthopedic Instrument Set — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 20-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. We work the instrument by

hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$700.00, falling to \$480.00 at 10 sets — about 31% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-30897559

20-Piece Veterinary Orthopedic Surgical Kit
Stainless Steel Bone

\$185.00 per piece

Trade price from **\$130.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

20-Piece Veterinary Orthopedic Surgical Kit Stainless Steel Bone is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$185.00, falling to \$130.00 at 1 set — about 30% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001381

22 Pcs Feline Spay Kit with Stainless Steel Box

\$64.99 per piece

Trade price from **\$59.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality 22 Feline Spay Kit Including Scalpel Blades and Stainless Steel Box German Grade Feline Spay Surgery Instruments Pack

Spaying female pets is recommended by most vets. It reduces the number of unwanted puppies and kittens and also prevents health problems such as pyometra, mammary tumors, and ovarian, uterine, and cervical cancer. Feline Spay/Neuter Pack includes all the surgical instruments that veterinarians require for performing a cat spay.


The Feline Spay Neuter Pack Kit is essential during the spay and neutering procedure. This kit is comprised of all instruments to fulfill spaying procedures and for safe and effective spay surgery.. All instruments are

crafted from 100% premium German grade surgical stainless material and thus can be reused after sterilization.

22 PCS FELINE SPAY KIT
INCLUDES

Our Feline Spay Kit contains all the necessary instruments used for spay surgery including stainless steel instruments box for easy storage.

- 1 TC Olsen Hager Needle Holder 14cm
- 1 Metzenbaum Scissors 17cm Curved
- 2 Kelly Forceps Straight 14cm
- 2 Mosquito Forceps Curved 12.5cm
- 2 Carmalt Forceps Straight 16cm
- 2 Carmalt Forceps Curved 16cm
- 1 Adson Brown Forceps 12cm 7:7 Teeth
- 1 Scalpel Handle #3



**WI-70001923****24 Pcs Feline Spay Pack Economy Grade****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking



24 Pcs Canine and Feline Spay Neuter Pack Surgery Kit Surgical Veterinary Instruments

Feline Spay/Neuter Pack includes all the surgical instruments that veterinarians require for performing a cat spay. Canine and Feline Spay Neuter Pack Kit is essential during the spay and neutering procedure. This kit is comprised of all instruments to fulfill spaying procedures. All instruments are crafted from 100% premium grade surgical stainless material and thus can be reused after sterilization

FELINE SPAY KIT INCLUDES

- 2 Kelly Hemostatic Forceps 5.5"(14cm) Curved
- 2 Kelly Hemostatic Forceps 5.5"(14cm) Straight

- 2 Mosquito Forceps 5" (12.7cm) Curved
- 2 Mosquito Forceps 5" (12.7cm) Straight
- 4 Backhaus Towel Clamp 3.5" (8.89cm)
- 2 Thumb Dressing Forceps 5.5"(14cm)
- 1 Adson-Brown Forceps, 4.75" (12cm)
- 1 Mayo Scissors, 5.75" (15cm)
- 1 Mayo-Hegar Needle Holder 6.25" (15.87cm)
- 1 Olsen-Hegar Needle Holder 5.5"(14cm)
- 1 Snook Hook
- 1 Scalpel Handle #3
- 4 Scalpel Blades #10

TECHNICAL INFORMATION**WI-70001924****24 Pcs T/C Feline Spay Pack****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


24 Pcs T/C Canine and Feline Spay Neuter Pack Surgery Kit Surgical Veterinary Instruments

Feline Spay/Neuter Pack includes all the surgical instruments that veterinarians require for performing a cat spay. Canine and Feline Spay Neuter Pack Kit is essential during the spay and neutering procedure. This kit is comprised of all instruments to fulfill spaying procedures. All instruments are crafted from 100% premium German grade surgical stainless material and thus can be reused after sterilization

T/C FELINE SPAYING PACK INCLUDES

- 2 Kelly Hemostatic Forceps 5.5"(14cm) Curved

- 2 Kelly Hemostatic Forceps 5.5"(14cm) Straight
- 2 Mosquito Forceps 5" (12.7cm) Curved
- 2 Mosquito Forceps 5" (12.7cm) Straight
- 4 Backhaus Towel Clamp 3.5" (8.89cm)
- 2 Thumb Dressing Forceps 5.5"(14cm)
- 1 Adson-Brown Forceps, 4.75" (12cm)
- 1 T/C Mayo Scissors, 5.75" (15cm)
- 1 T/C Mayo-Hegar Needle Holder 6.25" (15.87cm)
- 1 T/C Olsen-Hegar Needle Holder 5.5"(14cm)
- 1 Snook Hook
- 1 Scalpel Handle #3



**WI-70003712****27 Pcs Male Dog Scrotal Castration – Neuter Kit****\$89.99** per pieceTrade price from **\$82.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Male Dog Scrotal Castration – 27 Pieces Neuter Kit German Grade Stainless Steel Surgical Veterinary Instruments Set

Male dog castration or neutering is the surgical removal of the testicles. Most of the veterinarians agree that the best age for castration is between 6-12 months of age. At this age, the procedure is very low risk and easy to perform. Early castration is the responsible thing to do. It could save the dog's life and make him a much more desirable housemate.

Our professional quality 27 Pcs male dog castration Neutering kit contains all the necessary tools to perform a complete procedure. All tools are made from surgical grade stainless steel and used by professionals.

MALE DOG SCROTAL CASTRATION INCLUDES

Dog Neutering Kit consists of;

- 1 T/C Olsen Hagar Needle Holder 16cm
- 1 T/C Mayo Scissors 15cm Straight
- 2 Kelly Forceps 14cm Straight
- 2 Kelly Forceps 14cm Curved
- 2 Mosquito Forceps 12.5cm Straight
- 2 Mosquito Forceps 12.5cm Curved
- 2 Dressing Forceps 16cm
- 2 Tissue Forceps 16cm 1:2 Teeth
- 1 Scalpel Handle #3
- 4 Scalpel Handle Blades
- 4 Backhaus Towel Forceps 11cm
- 1 Adson Forceps 12cm
- 1 Adson Forceps 12cm 1:2 Teeth
- 1 Adson-Brown Forceps 12cm 7:7 Teeth

**WI-70000457****3 Equine Dental Float Rasps Set Angle Up, Down & Straight Adjustable Wooden Handle****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Equine Dental Float Rasps Set 24" (61cm) Angle Up, Down & Straight with Adjustable Wooden Handle Premium Quality Stainless Steel Horse Teeth Veterinary Instruments.

Float rasp used to file overgrowths on the teeth or reducing sharp points in order to allow the horse to chew more comfortably. Horse dental float rasp has a high degree of accuracy and flexibility, improving efficiency when using it. This equine dental float rasp set is made of stainless steel, rust proof, wear resistant, resistant, and can be used for a long time. High flexibility fine polishing having smoother surfaces equine dental kit for animal farm and professional veterinarian use. Float rasps are a crucial dental tool for horses, providing efficient tooth filing

and maintenance, enhancing the veterinary toolkit for better dental outcomes.

EQUINE DENTAL FLOAT RASP FEATURES

- Crafted 24" (61cm) Overall Length.
- Available in Up, Down and Straight Blade Profile.
- Adjustable Wooden Handle.
- Ideal for Animal Farm and Professional Veterinarian Use.
- Reusable.
- Made of Stainless Steel

PACKAGE INCLUDE:

- 1 Pc Equine Dental Float Rasp – Up – Adjustable Wooden Handle
- 1 Pc Equine Dental Float Rasp – Down – Adjustable Wooden Handle



**WI-70004003****3 Pcs Cat Nail Clipper Kit with Carrying Case****\$17.99** per pieceTrade price from **\$16.55**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Stainless Steel 3 Pcs Cat Nail Clipper with Nail File Kit – Professional Pet Nail Clipper Suitable for Cat, Dog & Small Pets to Avoid Over Cutting.

Trimming cats nail is challenging without practice and proper tools. Hold the clipper on the right hand and a paw in other hand and press a toe pad gently to extend the claw. Catchy pet nail clipper with nail file set for dogs and cats has a unique design to avoid overcutting. This pet clippers set suitable for most small pets such as kittens, cats, rabbits, puppy, pig and so on, practical tools for animal trainers and pet keepers.

Keep your pet nails clean and trim regularly with this nail clippers set. Cats nail clipper kit is designed for cats and small animals. Easily remove dirt hiding

in nails, and then file their nails to make nails more neat, give your pet a clean look.

CAT NAIL CLIPPER FEATURES

- Fast & Easy Trimming
- Premium Quality Pet Grooming Tools
- Ensures Swift & Clean Cut
- Suitable for Cat, Dogs and Small Pets
- Made With High Quality Stainless Steel

PACKAGE INCLUDE

- 1 Cat Nail Clipper
- 1 Nail File
- 1 Carrying Case

TECHNICAL INFORMATION**MATERIAL:** Plastic, Stainless Steel**INSTRUMENTS TYPE:** Veterinary Instruments**WI-70003715****3 Pcs Dog Nail Clipper Kit****\$17.99** per pieceTrade price from **\$16.55**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Stainless Steel 3 Pcs Dog Nail Clipper with Nail File Kit – Professional Pet Nail Clipper Suitable for Cat, Dog & Small Pets to Avoid Over Cutting.

Trimming dog nail is challenging without practice and proper tools. Hold the clipper on the right hand and a paw on the other hand and press a toe pad gently to extend the claw. Catchy pet nail clipper with nail file set for dogs and cats has a unique design to avoid overcutting. This pet clippers set suitable for most small pets such as kittens, cats, rabbits, puppy, pig and so on, practical tools for animal trainers and pet keepers.

Keep your pet nails clean and trim regularly with this nail clippers set. This nail clipper kit is designed for dog and small animals. Easily remove dirt hiding

in nails, and then file their nails to make nails more neat, give your pet a clean look. This kit include excellent quality tools used by professionals.

DOG NAIL CLIPPER FEATURES

- Fast & Easy Trimming
- Smooth Action
- Premium Quality Pet Grooming Tools
- Ensures Swift & Clean Cut
- Suitable for Cat, Dogs and Small Pets
- Made With High Quality Stainless Steel

PACKAGE INCLUDE

- 1 Dog Nail Clipper
- 1 Pet Nail File
- 1 Carrying Case

TECHNICAL INFORMATION**MATERIAL:** Plastic, Stainless Steel

**WI-70003714****3 Pcs Dog Nail Clipper Kit Guillotine Type****\$17.99** per pieceTrade price from **\$16.55**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Guillotine Dog Nail Clipper with Nail File and Carrying Case Set Stainless Steel Groom Professional and Pet Keeper Nail Clipper Kit is Best Fitted for Smaller Breeds with Lighter Nails.

This dog nail clipper kit is perfect for the size for medium to large-sized dogs. Guillotine type nail clipper/trimmer with spring action is easy to use, The nail should be cut from top to bottom, not side to side. The tip of the nail is placed in the hole in the trimmer and the handles are pressed to advance the cutting blade through the nail. They have a comfortable, anti-slip handle and are safe and accurate to use designed for dogs and small animals.

GUILLOTINE DOG NAIL CLIPPER KIT FEATURES

- Sharp & Long-lasting Blades
- Double-bladed Guillotine Style Dog Nail Trimmer
- Owner-friendly Nail Clipper With Anti-slip Handle
- Suitable for Cat, Dog and Smaller Breeds
- Safe, Accurate and Easy to Use
- Made With High Quality Stainless Steel

PACKAGE INCLUDE

- 1 Guillotine Nail Clipper
- 1 Nail File
- 1 Carrying Case

TECHNICAL INFORMATION**MATERIAL:** Plastic, Stainless Steel**WI-70000452****3 Pcs Nylon OB Calf Pulling Straps 30", 45" 60" with Hook Handle Set****\$43.99** per pieceTrade price from **\$40.47**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Nylon OB Calf Pulling Straps 30", 45" 60" with OB Hook Handle Set for Aid in Difficult Births Cow, Goat, and Cattle Delivery on Dairy Farms Highly Durable Veterinary Calving Supplies.

Nylon OB Calf Strap & OB Hook Handle used in the cow-calf industry for difficult births especially for use with goats and cattle delivery on dairy farms. OB Calf Pulling Strap with stainless-steel D rings made with soft yet sturdy, washable nylon. This obstetrical strap is 1" wide and spreads pulling tension over a wider area of the calf's legs to help prevent bone fractures and soft-tissue injury.

Calf strap has a hand loop that allows the user to apply pressure with Hook Handle, an ideal instrument for birth

assistance in calves and other cattle animals. The robust and strongly constructed calf strap with handle provides years of service with an easy and strong grip while calving.

OB CALF STRAP WITH HANDLE SET FEATURES

- **Calving Strap with Hook Set:** Nylon Calf Pulling Strap with OB Hook Handle are ideal cattle delivery supplies for aid in difficult births.
- **Spreads Stress:** Wide strap distributes pulling tension over calf's legs, reducing risk of injury.
- **Attaches to Calf Puller:** Metal loops allow attachment to calf puller or traction device.
- **Hand Loop:** Allows pressure application.



WI-70001276

3 Pcs Ratchet Calf Puller with OB Handle & 60" Chain Cattle Delivery Veterinary Set

\$194.99 per piece

Trade price from **\$179.39**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Set of 3 Ratchet Style Calf Puller Heavy Duty Calving Cattle Birthing Extractor 65" Standard Rod with 60 Inches Long OB Chain & Handle Veterinary Surgical Cattle Farming Instruments

- High Quality Stainless Steel Jack Handle
- Bree chain Handle Length: 25"
- Calf Puller Total Length: 65" (165 cm)

RATCHET STYLE CALF PULLER BENEFIT

- Corrosion resistant, lightweight, durable for longer life.
- One person operation permitting a smooth release on both sets of cams.
- Guards against slipping and provides uniform tension and minimal rod damage.
- Smooth for easy cleaning and disinfecting – Convenient to carry and assemble.
- For proper positioning of the Breechain below pin bones

3 PCS CALF PULLER CATTLE DELIVERY SET INCLUDED

- 1 Ratchet Style Calf Puller 65" Set
- 1 Pc Calf Puller OB Chain 60" Long
- 1 Pc Calf Puller OB Chain Handle

CALF PULLER FEATURES:

- Aero Space Heat Treated Steel Alloy Jack
- State-of-the-Art Handle Release
- Hardened Steel Cams
- Precision Two-Section Rod Assembly
- Rotating Jack and Bree chain



WI-70003683

3 Pcs Serra Emasculators 35.5cm with Ratchet

\$164.99 per piece

Trade price from **\$151.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Heavy Duty Curved Handle Serra Emasculator 14" (35.5cm) Triple Crush Castration Clamp Surgical Veterinary Instrument

Serra Emasculator has three actions: clamping, crushing and cutting comes with curved handles ratchet lock. It draws the cord concentrically into the jaws, thus preventing slipping of the cord, as the jaws are closed on the spermatic cord it is gripped by the clamping part. Ratchets hold emasculator locked and Back-lock retains the clamping part closed unless unlocked.


The stainless steel castration tool is designed for castration minimal blood loss as blood vessels are crushed immediately after they are being cut. It reduces trauma with minimal blood loss.

SERRA EMASCULATOR FEATURES:

- No Blood Loss Because Blood Vessels are Crushed.
- Triple Crush Jaws with Clamping, Crushing & Cutting Action
- Less Trauma
- Crushes the Arteries Severing them Cleanly.
- Minimizes the Risk of Serious Bleeding.
- Work Efficiency is High, Have More Work Done in Less Time.
- Back-lock Retains the Crushing Portion Locked Unless Unlocked.
- Made With High Quality Stainless Steel

PACKAGE INCLUDES:

- 3 Pcs Serra Emasculator Castration Clamp 14" (35.5cm) Triple Crush –



**WI-70000589****3 Pcs Spay Snook Hook Set 8\"/>****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 3 Pcs Spay Snook Hook Set 8\"/>

Spay Snook Hook is a specialized surgical tool used largely in veterinary surgeries for neutering dogs and cats. Ovarian Snook Hook designed for maximum performance and longevity, assists veterinarians in the safe and precise removal of reproductive organs. Snook hook features a curved blunt tip, long shank and a flat solid handle for gripping. The hook-shaped tip instrument is also used to retrieve the horn of the uterus while performing ovariectomy in felines and canines. This veterinary spay hook made of german-grade stainless steel,


which keeps it rust-free and autoclavable for easy sterilization and can be reused.

SPAY SNOOK HOOK FEATURES

- Crafted with 8\"/>
- 5-6mm Wide Curved Blunt Tip.
- Hook Shaped Working End.
- Long Straight Shank.
- Flat Smooth Grip Handle.
- Used in Neutering Dogs and Cats During Ovariectomy Procedure.
- Reusable.
- Made with Premium Quality Stainless Steel;

PACKAGE INCLUDE

- 3 Pcs Spay Snook Hook 8\"/>

**WI-70003634****34 Pcs Entomology Dissection Kit****\$25.99** per pieceTrade price from **\$23.91**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 34 Pcs Entomology Dissecting Kit Premium Quality Stainless Steel Insects Dissection Instruments Set with Free Carrying Case.

Entomology Dissecting Kit is used to study the anatomy, morphology, and physiology of butterflies, flies, bees, and other insects. Entomology Dissection Instruments Kit is the perfect tool set for all your entomology needs. Our kit includes a variety of high-quality tools designed specifically for the dissection of insects.


The kit comes in a sturdy travel case ensuring safe transport and storage. Ultrasonically cleaned and Passivated to enhance the durability. This Kit allows you to perform precise and delicate cuts, grips, and manipulations without damaging the specimens.

Whether you are a student, a teacher, a researcher, or a hobbyist, this entomology dissection kit is an essential tool for your insect studies. The ergonomically designed tools are made of premium grade stainless steel material and can be reused.

NOTE: We can make custom packs according to your requirements. Contact us for details.

ENTOMOLOGY DISSECTION KIT INCLUDES

- 2 Teaser Needles Straight
- 2 Teaser Needle Angled
- 1 Iris Scissors 11.5cm
- 1 Scalpel Handle #3
- 1 Adson Micro Forceps 12cm
- 1 Dropper



**WI-70001149****4 Pcs Calf Puller with 60" OB Chain Handle & Carrying Bag Cattle Delivery Set****\$229.99** per pieceTrade price from **\$211.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

4 Pcs Calf Puller Veterinary Cattle Delivery Set – Heavy Duty Ratchet Style Calf Puller Cattle Birthing Calving Extractor 65" Standard Rod with 60 Inches Long OB Chain & Hook Handle with Zipper Carrying Case Veterinary Surgical Cattle Farming Instruments

- Precision Two-Section Rod Assembly
- Rotating Jack and Bree chain
- High Quality Stainless Steel Jack Handle
- OB Chain Length: 39.37" (100cm)
- Bree chain Handle Length: 25"
- Calf Puller Total Length: 65" (165 cm)
- Water Resistant Black Zipper Carrying Case

4 PCS CALF PULLER CATTLE DELIVERY SET INCLUDED

- 1 Ratchet Style Calf Puller 65" Set
- 1 Pc Calf Puller OB Chain 60" Long
- 1 Pc Calf Puller OB Chain Handle
- 1 Pc Zipper Carrying Case

CALF PULLER FEATURES:

- Aero Space Heat Treated Steel Alloy Jack
- State-of-the-Art Handle Release
- Hardened Steel Cams

RATCHET STYLE CALF PULLER BENEFIT

- Corrosion resistant, lightweight, durable for longer life.
- One person operation permitting a smooth release on both sets of cams.
- Guards against slipping and provides uniform tension and minimal rod

**WI-70003720****4 Pcs Canine Mouth Gag Set 3.5", 4.5", 5.5" & 6.5"****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


Pack of 4 Canine Mouth Gag Set 3.5", 4.5", 5.5" and 6.5" Premium Quality Stainless Steel Dog Oral Dental Examination Surgical Veterinary Instruments.

Canine Veterinary Mouth Gag is a spring-loaded stainless steel tool used to facilitate oral examinations and procedures in dogs by keeping their mouths open. It provides stable access to the oral cavity, prevents the jaws from closing, and aids in intubation while minimizing trauma through its nylon inserts. These mouth gags are easy to sterilize and autoclavable, making it reliable for veterinary practitioners.

- X-large Size Measuring 6.5 Inches (16.5 cm).
- Spring-loaded Design With Nylon Inserts to Minimize Trauma.
- Prevents the Jaws From Closing and Aids in Intubation.
- Designed for Dog's Oral Dental Examinations and Surgeries.
- Available Various Sizes to Fit Different Breeds and Dog Sizes.
- Easy to Sterilize and Fully Autoclavable.
- Made With Corrosion Resistant High-quality Stainless Steel.

PACKAGE INCLUDES:

- 1 Pc Canine Mouth Gag Small Size 3.5"
- 1 Pc Canine Mouth Gag Medium Size 4.5"

CANINE MOUTH GAG FEATURES:



WI-70000726

4 Pcs Dog Ear Cropping Clamps Set “C”, “I”, “S”, “L” Type with Teeth

\$49.99 per piece

Trade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 4 Dog Ear Cropping Clamps Guide C, I, S & L Type Veterinart Retractor with Teeth – Stainles Steel Professional Veterinarian Ear Trim Clamps Retractor for Various Dog Breeds Livestock Medical Surgical Veterinary Instruments.

Ear cropping clamps are specially designed for the individual shape of each breed and this allows for better cosmetic and more consistent ear trims. The practice of cropping a dog's ears involves surgically altering the shape of ears for aesthetic reasons. It is mainly done for certain dog breeds to conform to breed standards or traditional preferences. Moreover, cropped ears may reduce the risk of ear infections and improve overall ear hygiene by allowing better airflow. This

ear cropping clamps set consists of 4 sizes to accommodate all breeds.

Our ear cropping guid retractor set constructed from premium-grade materials. Equipped with meticulously crafted tiny serrations. It confidently grips the ear, ensuring precise cuts resulting in aesthetically pleasing ear shapes that meet breed standards. Not only does the clamp streamline the process for veterinarians, but it also minimizes discomfort leading to a smoother and quicker recovery. Highly resistant to corrosion and rust. Precisely crafted for accurate use.

NOTE: Generally Ear Cropping Guide Clamps are best used on 7 to 8 week pups. Must follow professional veterinarian guide before cropping.



WI-70001402

4 Pcs Winged Elevators Set 1mm-4mm Tips Standard Handle

\$29.99 per piece

Trade price from **\$27.59**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 4 Veterinary Winged Elevator Set 1mm, 2mm, 3mm & 4mm Tips with Standard Handle – Luxating, Dental Extraction, Small Animals Veterinary Specific Elevators

Winged Elevators are especially designed for small breed animal dental work. The curved edges are designed to improve control and follows the contour of the tooth itself allowing for better leverage when compared to other elevators. The short winged blades with lateral edges which curve upward farther than traditional elevators. The winged shape of the elevator tip matches the curve of the tooth and the extra leverage aids in extraction procedure. The shorter narrow neck enhances visibility of the root, handling and control. Our winged


elevators with increased surface area help to sever the periodontal ligament easily. They have a surgical grade German stainless body and are available 6 inches with various tip variations.


SET INCLUDES

- 1 Pc Winged Elevator 1mm – Standard Handle
- 1 Pc Winged Elevator 2mm – Standard Handle
- 1 Pc Winged Elevator 3mm – Standard Handle
- 1 Pc Winged Elevator 4mm – Standard Handle

WINGED ELEVATORS FEATURES

- Short Winged Blades.
- Thin and Sharp Working Ends.





WI-70001396

47 Pcs Deluxe Veterinary Kit Spay Surgery Set

\$139.99 per piece

Trade price from **\$128.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Set of 47 Pcs Deluxe Veterinary Dissection Kit. Excellent Quality Surgery Pack for Routine Spay and General Surgery for Veterinarian. Our spay and general surgery veterinary instruments kit is packed in a stainless steel box. All instruments are made of German-grade stainless steel can be reused after sterilization.

DELUXE VETERINARY DISSECTION KIT INCLUDE

This comprehensive veterinary surgical pack compiled for vet clinics and professionals

- 2 Mayo Hegar Needle Holder 7"
- 1 Crile Wood Needle Holder 6"
- 1 Mayo Scissors 6.75" (straight)
- 2 Bandage Scissors 7.25"
- 3 Metzenbaum Scissors 7" (curved)

- 4 Carmalt Forceps 6.25" (straight)
- 1 Babcock Forceps 6.25"
- 2 Tissue Forceps 5" 1x2 teeth
- 2 Doyen Intestinal Forceps 9" (straight)
- 2 Adson Brown Forceps 4.75" 7X7 teeth
- 3 Mosquito Forceps 5" (curved)
- 3 Mosquito Forceps 5" (straight)
- 1 U.S. Army-Navy Retractor
- 12 Backhaus Towel Clamp 3.25"
- 2 Gelpi Retractor Sharp Points 7"
- 1 Periosteal Elevator
- 1 Spay Hook Snook Style
- 1 Kelly Forceps 5.25" (curved)
- 1 Grooved Director 6"
- 1 Instrument Tray with Cover
- 1 Medicin Bowl 50ml

TECHNICAL INFORMATION



WI-70001274

5 Pcs Calf Puller, Polansky Speculum, OB Chain & Handle Veterinary Instruments Set

\$274.99 per piece

Trade price from **\$252.99**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Set of 5 Ratchet Style Calf Puller Calving Extractor 65" Standard Rod with 60 Inches Long OB Chain, 2 Handles and Polansky Speculum Veterinary Surgical Cattle Farming Instruments

5 PCS CALF PULLER VETERINARY SET INCLUDED

- 1 Ratchet Style Calf Puller 65" Set
- 1 Pc Calf Puller OB Chain 60" Long
- 2 Pc Calf Puller OB Chain Handle
- 1 Pc Polansky Vaginal Speculum

CALF PULLER FEATURES:

- Aero Space Heat Treated Steel Alloy Jack
- State-of-the-Art Handle Release
- Hardened Steel Cams

- Precision Two-Section Rod Assembly
- Rotating Jack and Bree chain
- High Quality Stainless Steel Jack Handle
- Bree chain Handle Length: 25"
- Calf Puller Total Length: 65" (165 cm)

RATCHET STYLE CALF PULLER BENEFIT

- Corrosion resistant, lightweight, durable for longer life.
- One person operation permitting a smooth release on both sets of cams.
- Guards against slipping and provides uniform tension and minimal rod damage.
- Smooth for easy cleaning and disinfecting – Convenient to carry



**5
PCS****WI-70000402****5 Pcs Covault Spay Hook Set 9" with Curved Probe Tip****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Covault Spay Hook Set 9" (22.9cm) Octagon Handle with Curved Probe Tip Premium Quality Stainless Steel Veterinary General Surgery Instruments.

Covault spay hook is a surgical tool used in ovariectomy for dogs and cats to retrieve the uterus horn during spay procedures. A long, thin hook, designed with a curved end and buttoned tip, is inserted into the incision to grasp the uterine horn, enabling the surgeon to exteriorize and remove the uterus. The Covault Ovariectomy Hook, made from high-quality stainless steel, ensures safe and effective manipulation of the uterine horn, making it a reliable tool for long-term use in veterinary or medical settings. This veterinary spay hook

made of german-grade stainless steel, which keeps it rust-free and autoclavable for easy sterilization and can be reused.

COVAULT OVARIECTOMY HOOK FEATURES

- Crafted with 9" (22.9cm) Overall Length.
- Long, Thin Hook With a Curved End and Buttoned Tip.
- Enables Surgeon to Exteriorize and Remove the Uterus.
- Used in Spay Procedures for Dogs and Cats.
- Rust-free and Autoclavable for Reuse Purpose.
- Made From High-quality Stainless Steel for Safe Manipulation.

PACKAGE INCLUDE

9 INCH

WI-70001417**5 Pcs Dental Scaling Kit For Large Dogs****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Dental Scaling Kit For Large Dogs Pet High Quality Stainless Steel Tartar Scrapping Dental Tools Kit Surgical Veterinary Instruments

Tooth and gum problems are the most common in pets. Bad breath in dogs and cats go hand in hand with other health problems it is important to treat this problem even if the pet's breath is not objectionable.

Placid temperament is not difficult with our 5 Pcs Tartar Scrapping Kit. Many pet professionals perform excellent tooth cleaning at home. Kit Scalars remove the tartar and other unnecessary foreign bodies from teeth, dental mirror help to view deep teeth while dental tweezers remove an unnecessary object from teeth. This dental scaling kit includes a specially

designed spring action Mouth Gag opens the mouth with reduced trauma.

DENTAL SCALING KIT INCLUDES

Get your pet's teeth cleaned professionally with our 5 Pcs Dental Scaling Kit.

- 1 Pc Dental Tweezers
- 1 Pc Dental Scalars
- 1 Pc Dog Scalars
- 1 Pc Dental Explorer
- 1 Pc Mouth Gag 16.5cm

Note: Large dogs mouth gag is designed for smaller breed having height upto 25 inches at shoulders.

TECHNICAL INFORMATION

MATERIAL: Stainless Steel
INSTRUMENTS TYPE: Surgical Veterinary Instruments

**LARGE DOG
DENTAL SCALING**

**WI-70001418****5 Pcs Dental Scaling Kit For Medium Dogs****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Dental Scaling Kit For Medium Dogs Pet High Quality Stainless Steel Tartar Scrapping Dental Tools Kit Surgical Veterinary Instruments

Tooth and gum problems are the most common in pets. Bad breath in dogs and cats go hand in hand with other health problems it is important to treat this problem even if the pet's breath is not objectionable.

Placid temperament is not difficult with our 5 Pcs Tartar Scrapping Kit. Many pet professionals perform excellent tooth cleaning at home. Kit Scalars remove the tartar and other unnecessary foreign bodies from teeth, dental mirror help to view deep teeth while dental tweezers remove an unnecessary object from teeth. This

dental scaling kit includes a specially designed spring action Mouth Gag opens the mouth with reduced trauma.

DENTAL SCALING KIT INCLUDES

Get your pet's teeth cleaned professionally with our 5 Pcs Dental Scaling Kit.

- 1 Pc Dental Tweezers
- 1 Pc Dental Scalars
- 1 Pc Dog Scalars
- 1 Pc Dental Explorer
- 1 Pc Mouth Gag 14cm

Note: Medium dogs mouth gag is designed for smaller breed having height upto 18-25 inches at shoulders.

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

WI-70001420**5 Pcs Dental Scaling Kit For Small Dogs****\$21.99** per pieceTrade price from **\$20.23**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 5 Pcs Dental Scaling Kit For Small Dogs Pet High Quality Stainless Steel Tartar Scrapping Tools Kit Surgical Veterinary Instruments

Tooth and gum problems are the most common in pets. Bad breath in dogs and cats go hand in hand with other health problems it is important to treat this problem even if the pet's breath is not objectionable.

Placid temperament is not difficult with our 5 Pcs Tartar Scrapping Kit. Many pet professionals perform excellent tooth cleaning at home. Kit Scalars remove the tartar and other unnecessary foreign bodies from teeth, dental mirror help to view deep teeth while dental tweezers remove an unnecessary object from teeth. This dental scaling kit includes a specially

designed spring action Mouth Gag opens the mouth with reduced trauma.

DENTAL SCALING KIT INCLUDES

Get your pet's teeth cleaned professionally with our 5 Pcs Dental Scaling Kit.


- 1 Pc Dental Tweezers
- 1 Pc Dental Scalars
- 1 Pc Dog Scalars
- 1 Pc Dental Explorer
- 1 Pc Mouth Gag 11cm

Note: Small dogs mouth gag is designed for smaller breed having height upto 16 inches.

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Surgical Veterinary Instruments



**WI-70001910****5 Pcs Veterinary Otoscope Set Animal Diagnostic Instruments****\$22.99** per pieceTrade price from **\$21.15**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

3.2V Bright White LED Veterinary Otoscope Set Includes 3 Sizes Of Reusable Specula with FREE Zipper Carrying Case ENT Medical Surgical Diagnostic Instruments.

Ear infections are common in animals, especially those with long ears. Ear infections can cause extreme discomfort and further health problems if left untreated. A thorough clinical examination allows veterinarians to determine the cause of the infection. This brand new professional veterinary otoscope set comes with a handle, an otoscope head, 3 different sizes of specula and a protective zippered case. The disposable specula adaptor in sizes 4mm, 5mm, and 7mm – disposable specula works with our otoscope reducing the risk of



spreading infection. You will not be disappointed with the look, feel, and function of this otoscope.

5 PCS VETENIARY OTOSCOPE DIAGNOSTIC SET INCLUDES

- 1 Otoscope
- 1 Pc 4mm Specula
- 1 Pc 5mm Specula
- 1 Pc 7mm Specula
- 1 Pc 3.2 volts LED bulb

FREE Zippered Carrying Case**PROFESSIONAL OTOSCOPE SET FEATURES**

A COMPLETE DIAGNOSTIC LED OTOSCOPE KIT: Superior quality brand new professional diagnostic veterinary otoscope kit with 3 reusable specula and zippered carrying case. Also

**WI-70004007****6 Pcs Basic Dental Prophy Instruments Kit for Small Animals****\$39.99** per pieceTrade price from **\$36.79**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

6 Pieces Basic Dental Prophy Kit Stainless Steel Surgical Veterinary Instruments Set for Dental Prophylaxis Procedure



Periodontal disease is common in dogs and cats. The only way to prevent periodontal disease is to start before periodontal disease exists. The Prophy Kit is a set of dental surgical instruments used to perform dental prophylaxis on small animals. Our prophy kits feature an ergonomic design and robust structure that assist veterinary dental professionals in cleaning calculus, bacteria, and plaque from the teeth of small animals. Our Basic Dental Prophy Pack is designed to provide complete instrumentation for providing proper dental prophy and exam. All instruments are made from

medical grade stainless steel used by professionals.

BASIC DENTAL PROPHY KIT INCLUDES

Our basic dental prophylaxis set comes with the following instruments;

- 1 #5 Mouth Mirror Handle with Mirrors
- 1 Periodontal Probe CP15
- 1 Explorer 17/23
- 1 Towner-Jacquette Scaler U15/30
- 1 Columbia Curette 13/14
- 1 Columbia Curette 4R/4L

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Veterinary Dental Instruments**APPLICATION:** Dental Prophylaxis

**WI-70003723****6 Pcs Basic Dental Propy Kit****\$49.99** per pieceTrade price from **\$45.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

6 Pieces Basic Dental Propy Kit
Stainless Steel Surgical Veterinary
Instruments Set for Dental Prophylaxis
Procedure

Periodontal disease is common in adult dogs and cats. The only way to prevent periodontal disease is to start before periodontal disease exists. The Propy Kit is a set of dental surgical instruments used to perform dental prophylaxis on small animals. Our propy kits feature an ergonomic design and robust structure that assist veterinary dental professionals in cleaning calculus, bacteria, and plaque from the teeth of small animals. Our Basic Dental Propy Pack is designed to provide complete instrumentation for providing proper dental propy and exam. All instruments are made from

medical grade stainless steel used by professionals.

BASIC DENTAL PROPHY KIT INCLUDES

Our basic dental prophylaxis set comes with the following instruments;

- 1 #5 Mouth Mirror Handle with Mirrors
- 1 Periodontal Probe CP15
- 1 Explorer 17/23
- 1 Towner-Jacquette Scaler U15/30
- 1 Columbia Curette 13/14
- 1 Columbia Curette 4R/4L

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Veterinary Dental Instruments**APPLICATION:** Dental Prophylaxis**WI-70004006****6 Pcs Basic Exodontia Kit****\$59.99** per pieceTrade price from **\$55.19**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Basic Exodontia 6 Pieces Small Animal
Dental Extraction Kit Stainless Steel
Surgical Veterinary Instruments Set

Exodontia is the removal of a tooth from its alveolar bone socket. Pets with advanced periodontal disease may experience mobile teeth, pulp infection, and bone loss. As a result, veterinarians must extract severely damaged teeth. We offer a Basic Exodontia Kit to perform tooth extractions.

Tooth extraction is a common dental procedure performed in most veterinary practices. Our Basic Exodontia Kit or Basic Dental Extraction Kit contains necessary tools to perform extraction in cat, dog, and other small animals. All instruments are made from

medical grade stainless steel used by professionals.

BASIC EXODONTIA KIT INCLUDES

- 1 Apexo Root Elevator 301
- 1 Seldin Root Elevator 34
- 1 Heidbrink Root Elevator Fig. 1
- 1 Molt Periosteal Elevator Fig. 9
- 1 Hemingway Bone Curette Fig. 2
- 1 Parly Extracting Forceps 32A

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Dental Veterinary Instruments**APPLICATION:** Small Animal Dental Extraction Procedure**RUST PREVENTION****PROCEDURE:** Passivated

**WI-70003717****6 Pcs Dog Ear Suture Kit****\$35.16** per pieceTrade price from **\$32.35**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

6 Pcs Dog Ear Suture Kit with FREE Scalpel Blades Stainless Steel Surgical Veterinary Instruments.

Suture Kit assists veterinarians to suture ear after cropping. Our advanced Surgical Suture kit comprises of all the necessary instruments to suture ears after cropping. This dog suture kit can also be used for minor surgery. All instruments are made of surgical grade stainless steel to resist rust and stand up during regular daily usage. They are easy to clean, fully autoclavable and reusable.

DOG EAR SUTURE KIT FEATURES

This kit has professional quality instruments used by the veterinarian.

- Mayo-Hager Needle Holder 14cm
- Adson Forceps 12cm

- Stitch Removal Scissors 11.4cm
- Kelly Forceps Straight 14cm
- Operating Scissors 14cm Sharp/Blunt
- Scalpel Handle #3
- Disposable Scalpel Blades

TECHNICAL INFORMATION:**MATERIAL:** Stainless Steel**INSTRUMENT TYPE:** Surgical Veterinary Instrument**APPLICATION:** Ear Cropping Suturing Procedure**AUTOCCLAVABLE/REUSABLE:** Yes**RUST PREVENTION****PROCEDURE:** Passivated**ULTRASONIC CLEANED:** Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag Packing**DRESSING THUMB**
Stainless Steel, Straight, Coated**SUTURE STITCH SCISSOR**
Durable, 100% High Grade Stainless Steel**WI-70003713****6 Pcs Dog Tail Docking Kit****\$35.99** per pieceTrade price from **\$33.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

6 Pcs Dog Tail Docking Kit with Free Scalpel Blades Stainless Steel Surgical Veterinary Instruments Set for Veterinarian

Tail Docking Kit is used to dock the tail of puppies. Most veterinarians and senior breeders recommend that this be done between birth and up to five days of age so that puppies experience the least amount of pain or trauma. Our dog tail docking set consist of all nesseary instruments to performa a safe and effective tail docking.

TAIL DOCKING KIT INCLUDES

- 1 Tail Docking Scissors Heavy Pattern
- 1 Adson Forceps 12cm
- 1 Spencer Scissors 12cm
- 1 Mayo Hagar Needle Holder 16cm

- 1 Dressing Forceps 16cm
- 1 Scalpel Handle #4

with **FREE** Scalpel Blades**TECHNICAL INFORMATION****MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Veterinary Instrument**APPLICATION:** Dog Tail Docking**REUSABLE:** Yes**AUTOCCLAVED:** Yes**RUST PREVENTION PROCEDURE:**

Passivated

ULTRASONIC CLEANED: Yes**TESTS PERFORMED:** Boil Test, Performance Test, Shape Test**QC PASSED:** Yes**PACKING DETAILS:** Polybag PackingFind more **Dogs Veterinary Instruments****TAIL DOCKING SCISSORS**



WI-70003631

6 Pcs Equine Periodontal Instrument Kit

\$123.74 per piece

Trade price from **\$113.84**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

6 Pcs Equine Dental Preiodontal Instrument Kit with Zipper Carrying Case Premium Quality Horse Dentistry Tools for Professional Veterinary Instruments.

Regular oral examinations in horses are crucial for early detection of dental problems, reducing the risk of severe health issues. This equine periodontal examination tool set used to identify and address potential oral health issues in horses, enhancing their overall periodontal health. Periodontal disease, primarily caused by calculus, is a leading cause of tooth loss. Our equine dental periodontal kit includes surgical grade stainless steel tools for examination and calculus treatment in horses. Our horse dentistry equipment


set is crafted from high-quality stainless steel with precision.

EQUINE DENTAL PERIODONTAL KIT FEATURES

- 1 Equine Dental Banded Periodontal Probe 14.75" (37.4cm)
- 1 Equine Dental Explorer Fine Point 15.2" (38.6cm)
- 1 Equine Dental Sickle Scaler Large Hook 14.75" (37.4cm)
- 1 Equine Dental Sickle Scaler Small Hook 14.75" (37.4cm)
- 1 Equine Needle Pick, 15" (40cm)
- 1 Equine Dental Mouth Mirror 18" Head 50mm Dia. with T Handle

with **FREE** Zipper Carrying Case

TECHNICAL INFORMATION



WI-28887630

6 Pcs First Aid Kit

\$40.00 per piece

Trade price from **\$25.00**
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of 6 Pcs First Aid Kit, made for veterinary practice and animal surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$40.00, falling to \$25.00 at 10 sets — about 38% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28806200

6 Pcs Horse First Aid Kit Veterinary Surgical Instruments

\$30.00 per piece

Trade price from \$18.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of 6 Pcs Horse First Aid Kit Veterinary Surgical Instruments, made for veterinary practice and animal surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$30.00, falling to \$18.00 at 10 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-23883602

6 Pcs Horse Pets First Aid Kit

\$160.00 per piece

Trade price from \$100.00
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

6 Pcs Horse Pets First Aid Kit — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$160.00, falling to \$100.00 at 10 sets — about 38% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70004000****6 Pcs Pets First Aid Kit****\$24.99** per pieceTrade price from **\$22.99**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

6 Pcs Pets First Aid Instruments Kit with FREE Carrying Case Stainless Steel Small Animals Surical Veterinary Instruments Set

Our Pet First Aid kit is the must-have pet care product for everyday and emergency use. Pet's emergency can strike anytime; anyone who has pets must keep a pet first aid instrumetns kit in easy reach. Our Pets first aid kit is specifically designed for small animals; perfect for cats, pigs, dogs, and many other small animals.

The kit is designed to be used in the home as well as camping outside with pets. Besides an emergency, some instruments in the kit can be used regularly. This instrumetn set contain all necessary tools for first aid. All Tools

are made from Surgical Grade stainless steel and used by professionals.

NOTE: We also make custom design surgery kit to meet with specific needs.

PETS FIRST AID INSTRUMENTS KIT INCLUDES


- 1 Kelly Forceps 14cm Curved.
- 1 Operating Scissors 14cm Sharp/Blunt
- 1 Dressing Forceps 16cm
- 1 Dropper
- 1 Pc EMT Truama Scissors 14cm
- 1 Pet Toe Nail Clipper/ Trimmer

with **FREE** Carrying Case

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENTS TYPE: Veterinary Instrument

**WI-70000588****6 Pcs Spay Snook Hook Set 8" Curved 5-6mm Wide Blunt Tip Flat Handle****\$39.99** per pieceTrade price from **\$36.79**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Spay Snook Hook Set 8" (20.3cm) Curved 5-6mm Wide Blunt Tip Premium Quality Stainless Steel Ovariectomy Veterinary Surgical Instrument.

Spay Snook Hook is a specialized surgical tool used largely in veterinary surgeries for neutering dogs and cats. Ovarian Snook Hook designed for maximum performance and longevity, assists veterinarians in the safe and precise removal of reproductive organs. Snook hook features a curved blunt tip, long shank and a flat solid handle for gripping. The hook-shaped tip instrument is also used to retrieve the horn of the uterus while performing ovariectomy in felines and canines. This veterinary spay hook made of german-grade stainless steel,

which keeps it rust-free and autoclavable for easy sterilization and can be reused.

SPAY SNOOK HOOK FEATURES

- Crafted with 8" (20.3cm) Overall Length.
- 5-6mm Wide Curved Blunt Tip.
- Hook Shaped Working End.
- Long Straight Shank.
- Flat Smooth Grip Handle.
- Used in Neutering Dogs and Cats During Ovariectomy Procedure.
- Reusable.
- Made with Premium Quality Stainless Steel;

PACKAGE INCLUDE

- 6 Pcs Spay Snook Hook 8" (20.3cm) Curved Blunt Tip



**WI-70000172****6 Pcs Stripping and Carding Knives Set – Blue Wooden Handle****\$32.99** per pieceTrade price from **\$30.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Stripping & Carding Knives Set Premium Quality Stainless Steel Blades & Blue Wooden Handle. Efficient Undercoat Removal, Hand Stripping Tools for Dog Grooming.

Professional quality right-handed stripping and carding knife set is a specialized grooming tool for hand-stripping wire-coated dogs (like terriers and schnauzers). It effectively removes dead topcoat hairs without cutting, featuring knives of various tooth sizes to manage coat density and maintain breed-specific texture. This 6-piece set, made from high-quality German-grade steel with durable, ergonomic wooden handles, is suitable for both professionals and home use, providing efficient grooming for various coat types and ensuring comfort and control

during use. High quality German-grade steel blades ensure sharpness and rust resistance, enhancing grooming efficiency for breeders, groomers, and pet owners. Order this perfect pet grooming set today and keep your pet healthy and happy. Available in different colors

STRIPPING AND CARDING KNIVES SET FEATURES

- Multiple Knives Overall Lengths from 5" to 7".
- Knives With Coarse, Medium, Fine, and Extra Fine Teeth Sizes.
- Perfect for Flat and Fine Work, Rolling, Detailing Hand-stripping.
- Right-handed Model.
- Blades are Made From a High-quality Steel Ensure Sharpness and Rust

**WI-70000173****6 Pcs Stripping and Carding Knives Set – Garnish Red Handle****\$32.99** per pieceTrade price from **\$30.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Stripping & Carding Knives Set Premium Quality Stainless Steel Blades & Garnish Red Wooden Handle. Efficient Undercoat Removal, Hand Stripping Tools for Dog Grooming.

Professional quality right-handed stripping and carding knife set is a specialized grooming tool for hand-stripping wire-coated dogs (like terriers and schnauzers). It effectively removes dead topcoat hairs without cutting, featuring knives of various tooth sizes to manage coat density and maintain breed-specific texture. This 6-piece set, made from high-quality German-grade steel with durable, ergonomic wooden handles, is suitable for both professionals and home use, providing efficient grooming for various coat

types and ensuring comfort and control during use. High-quality German-grade steel blades ensure sharpness and rust resistance, enhancing grooming efficiency for breeders, groomers, and pet owners. Order this perfect pet grooming set today and keep your pet healthy and happy. Available in different colors

STRIPPING AND CARDING KNIVES SET FEATURES

- Multiple Knives Overall Lengths from 5" to 7".
- Knives With Coarse, Medium, Fine, and Extra Fine Teeth Sizes.
- Perfect for Flat and Fine Work, Rolling, Detailing Hand-stripping.
- Right-handed Model.



**WI-70000170****6 Pcs Stripping and Carding Knives Set – Multi Color Handle****\$32.99** per pieceTrade price from **\$30.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Stripping & Carding Knives Set Premium Quality Stainless Steel Blades & Multi Color Wooden Handle. Efficient Undercoat Removal, Hand Stripping Tools for Dog Grooming.

Professional quality right-handed stripping and carding knife set is a specialized grooming tool for hand-stripping wire-coated dogs (like terriers and schnauzers). It effectively removes dead topcoat hairs without cutting, featuring knives of various tooth sizes to manage coat density and maintain breed-specific texture. This 6-piece set, made from high-quality German-grade steel with durable, ergonomic wooden handles, is suitable for both professionals and home use, providing efficient grooming for various coat

types and ensuring comfort and control during use. High quality German-grade steel blades ensure sharpness and rust resistance, enhancing grooming efficiency for breeders, groomers, and pet owners. Order this perfect pet grooming set today and keep your pet healthy and happy. Available in different colors

STRIPPING AND CARDING KNIVES SET FEATURES

- Multiple Knives Overall Lengths from 5" to 7".
- Knives With Coarse, Medium, Fine, and Extra Fine Teeth Sizes.
- Perfect for Flat and Fine Work, Rolling, Detailing Hand-stripping.
- Right-handed Model.

**WI-70000171****6 Pcs Stripping and Carding Knives Set – Red Wooden Handle****\$32.99** per pieceTrade price from **\$30.35**

Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Pack of 6 Pcs Stripping & Carding Knives Set Premium Quality Stainless Steel Blades & Red Wooden Handle. Efficient Undercoat Removal, Hand Stripping Tools for Dog Grooming.

Professional quality right-handed stripping and carding knife set is a specialized grooming tool for hand-stripping wire-coated dogs (like terriers and schnauzers). It effectively removes dead topcoat hairs without cutting, featuring knives of various tooth sizes to manage coat density and maintain breed-specific texture. This 6-piece set, made from high-quality German-grade steel with durable, ergonomic wooden handles, is suitable for both professionals and home use, providing efficient grooming for various coat types and ensuring comfort and control

during use. High quality German-grade steel blades ensure sharpness and rust resistance, enhancing grooming efficiency for breeders, groomers, and pet owners. Order this perfect pet grooming set today and keep your pet healthy and happy. Available in different colors

STRIPPING AND CARDING KNIVES SET FEATURES

- Multiple Knives Overall Lengths from 5" to 7".
- Knives With Coarse, Medium, Fine, and Extra Fine Teeth Sizes.
- Perfect for Flat and Fine Work, Rolling, Detailing Hand-stripping.
- Right-handed Model.
- Blades are Made From a High-quality Steel Ensure Sharpness and Rust



**WI-70003719****7 Pcs Canine Periodontal Pack****\$25.99** per pieceTrade price from **\$23.91**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

7 Pcs Canine Periodontal Instrument Pack Small Animals Stainless Steel Surgical Dental Veterinary Instruments Kit



Periodontal diseases affect one or more of the periodontal tissues alveolar bone, periodontal ligament, cementum and gingiva. The treatment of periodontal disease begins with the removal of sub-gingival calculus (tartar). Our canine periodontal instrument kit has all the necessary tools to perform root planning and scaling.

The periodontal pack is ideally crafted to help veterinary dental professionals during root planing and scaling for small animals. The set includes seven periodontal instruments, including curettes, scalers, probes, explorers,

and mirror handles. All these instruments are long-lasting and reusable.

CANINE PERIODONTAL INSTRUMENT KIT INCLUDES

- Mouth Mirror Handle
- Mouth Mirror #4
- Explorer 17/23
- Columbia Curette 13/14
- Columbia Curette 4L/4R
- Gracy Curette 1/2
- Towner Jacquette U15/30

TECHNICAL INFORMATION**MATERIAL:** Stainless Steel**INSTRUMENTS TYPE:** Veterinary Dental Instrument**APPLICATION:** Small Animals Periodontal Procedure**REUSABLE:** Yes**WI-70002745****7 Pcs Veterinarian Umbilical Cord Scissor Kit****\$35.99** per pieceTrade price from **\$33.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Set of 7 Pcs Veterinarian Umbilical Cord Scissors Kit Including Carry Case High-Grade Stainless Steel Instruments Perfect for Anatomy, Biology, and Veterinary Students.


7-piece veterinarian umbilical cord scissor kit is a stainless steel instrument set for safe cutting and management of umbilical cords during animal births, such as puppies and kittens. It is used for whelping, anatomy studies, and veterinary education, helping to prevent injury when cutting the umbilical cord of newborn small animals. The stainless steel tools are typically autoclavable, making them suitable for sterilization.

VETERINARIAN UMBILICAL CORD SCISSOR KIT FEATURES

DESIGN FOR STUDENTS: This is a veterinary student's umbilical cord scissors set, perfect for dogs and puppies. The two included scissors are suitable for large, medium, and small breeds of animals.

A KIT YOU'LL BE PROUD TO USE: With 6 reusable tools and 1 premium carrying case, this kit is the perfect balance of functionality, value, and quality. A complete kit for all your needs: our kit is suitable for teachers, students, and hobbyists, and suitable for studies in animal anatomy, zoology, veterinary and biology, and much more. This is a student-level medical instrument set and not a clinical grade.

OUR QUALITY: Our tools are made from high-quality German stainless steel and will not corrode, tarnish, or rust. All conventional sterilization





WI-08756437

72 PC US Minor Surgery Surgical Veterinary

\$65.00 per piece

Trade price from \$48.00
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

72 PC US Minor Surgery Surgical Veterinary is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$65.00, falling to \$48.00 at 10 sets — about 26% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25860394

7-Piece Horse Cattle Farrier

\$30.00 per piece

Trade price from \$15.00
Minimum order 10 bags

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture 7-Piece Horse Cattle Farrier for veterinary practice and animal surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$30.00, falling to \$15.00 at 10 bags — about 50% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25894370

7-Piece Reusable Veterinary Horse Care Farrier

\$30.00 per piece

Trade price from **\$15.00**
Minimum order 10 bags

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

7-Piece Reusable Veterinary Horse Care Farrier is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$30.00, falling to \$15.00 at 10 bags — about 50% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70001373

9 Pcs Large Animal Eye Kit

\$54.99 per piece

Trade price from **\$50.59**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Large Animal Eye Kit 9 Pieces German Grade Stainless Steel Animal Ophthalmology Veterinary Instruments Set

Large Animal Eye Kit is a high-quality and versatile kit for veterinary ophthalmic surgery. This kit contains 9 instruments that are specially designed for eye procedures on large animals such as horses, cattle, and zoo animals. Whether you need to perform diagnostics, therapeutics, or surgery, this animal ophthalmology kit has everything you need to ensure a successful outcome.

LARGE ANIMAL EYE KIT INCLUDES

Our Large Animal Ophthalmology kit provides a range of instruments used

for veterinary ophthalmic surgery. The set consists of;

- 1 Bishop Harmon Forceps 8.5cm
- 1 Mosquito Forceps 12.5cm Straight
- 1 Mosquito Forceps 12.5cm Curved
- 1 Colibri Forceps 7.3cm
- 1 T/C Castroviejo Forceps 14cm
- 1 Eyelid Speculum Large
- 1 Micro Blade Handle
- 1 T/C Iris Scissors 11cm Curved
- 1 Castroviejo Corneal Scissors 11cm Right

TECHNICAL INFORMATION

MATERIAL: Stainless Steel

INSTRUMENT TYPE: Veterinary Ophthalmology Instrument

APPLICATION: Large Animals Eye Procedures

REUSABLE: Yes



TC IRIS SCISSORS
CURVED



MOSQUITO
FORCEPS CURVED



WI-23513026

Advanced 30-Piece Dissection Kit Medical Students Veterinary Surgical

\$7.00 per piece

Trade price from **\$4.50**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Advanced 30-Piece Dissection Kit Medical Students Veterinary Surgical is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$7.00, falling to \$4.50 at 1 set — about 36% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-21725912

Advanced 45-Piece Stainless Steel Dissecting Kit

\$13.00 per piece

Trade price from **\$6.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Advanced 45-Piece Stainless Steel Dissecting Kit — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$13.00, falling to \$6.00 at 1 set — about 54% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-05867732

All-In-One 14-Piece Stainless Steel Minor Surgery Veterinary Dental

\$12.65 per piece

Trade price from \$7.86
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of All-In-One 14-Piece Stainless Steel Minor Surgery Veterinary Dental, made for veterinary practice and animal surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$12.65, falling to \$7.86 at 10 sets — about 38% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28812440

Animal Orthopedic Surgical Instrument 1.5/2.0/2.4mm Vet LC-DCP Bone

\$840.00 per piece

Trade price from \$760.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Animal Orthopedic Surgical Instrument 1.5/2.0/2.4mm Vet LC-DCP Bone — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The 2.4mm size shown is our standard; other sizes in the same pattern are made to order.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$840.00, falling to \$760.00 at 1 set — about 10% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-28795310****Basic Starter Pets First Aid Kit****\$30.00** per pieceTrade price from **\$18.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Basic Starter Pets First Aid Kit is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$30.00, falling to \$18.00 at 10 pieces — about 40% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-70003722****Bitch Spay-Neuter Kit 25 Pcs Spay Surgery Instruments****\$85.99** per pieceTrade price from **\$79.11**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Premium Quality 25 Pcs Bitch Spay-Neuter Kit Stainless Steel Surgical Veterinary Spay Surgery Neuter Pack – Gold


This kit provides all the necessary instruments for veterinary surgeons to perform a bitch spay.

Ovariohysterectomy is one of the most common procedures performed in veterinary clinics. This procedure helps prevent mammary tumours, pyometra and unwanted pregnancy. This Bitch Spay-Neuter Kit having comprehensive set of instruments for bitch desexing (spaying) All instruments are made of surgical grade stainless steel to resist rust and stand up during regular daily usage. They are easy to clean, fully autoclavable and reusable.

BITCH SPAY-NEUTER KIT INCLUDES

The Bitch Spay Pack (Bitch Spay-Neuter Kit) have professional quality instruments used by the veterinarian.

- 2 Kelly Hemostatic Forceps 5.5" Curved
- 2 Kelly Hemostatic Forceps 5.5" Straight
- 2 Mosquito Forceps 5" Curved
- 2 Mosquito Forceps 5" Straight
- 4 Backhaus Towel Clamp 3.5"
- 1 Dressing Forceps 5.5"
- 1 Scalpel Handle #3
- 5 Scalpel Blades
- 1 Adson-Brown Forceps 4.75"
- 1 T/C Mayo Scissors 6.75" Straight
- 1 Mayo-Hegar Needle Holder 6"





WI-29411879

Botany Dissection Kit for Medical Biology Veterinary Students

\$8.00 per piece

Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Botany Dissection Kit for Medical Biology Veterinary Students is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$8.00.

Minimum order 10 sets, and the price moves with quantity. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29308398

Canine Feline Spay Pack Veterinary Surgical Instruments

\$55.00 per piece

Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Canine Feline Spay Pack Veterinary Surgical Instruments is a veterinary instrument we forge, grind and finish in our own Sialkot workshop. The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish.

We state the grade we used on your quotation. Listed at \$55.00.

Minimum order 10 sets, and the price moves with quantity. Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36187616

Cat Spay Neuter Surgical Instrument Set

\$180.00 per piece

Trade price from \$164.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Cat Spay Neuter Surgical Instrument Set — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$180.00, falling to \$164.00 at 1 set — about 9% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32658802

Complete Orthopedic Stainless Steel Suture Kit for Medical

\$10.00 per piece

Trade price from \$7.00
Minimum order 25 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Complete Orthopedic Stainless Steel Suture Kit for Medical is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$10.00, falling to \$7.00 at 25 pieces — about 30% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32067215

Complete Suture Practice Kit Stainless Steel Instruments Pre-Cut

\$7.00 per piece

Trade price from \$4.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Complete Suture Practice Kit Stainless Steel Instruments Pre-Cut is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$7.00, falling to \$4.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29399967

Deluxe Veterinary Dissection Kit

\$15.00 per piece

Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Deluxe Veterinary Dissection Kit for veterinary practice and animal surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$15.00.

Minimum order 10 sets, and the price moves with quantity. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29552734

Dog Dental Extraction Kit 18Pcs Veterinary Surgical Instruments Set

\$8.00 per piece

Trade price from **\$5.00**
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Dog Dental Extraction Kit 18Pcs Veterinary Surgical Instruments Set, made for veterinary practice and animal surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a 18-piece set.

We can add, drop or substitute items so the tray matches your protocol rather than ours. The steel is heat treated and

passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation.

Listed at \$8.00, falling to \$5.00 at 10 sets — about 38% off the single-unit price. Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-24106098

Economy 18-Piece Canine Spay Veterinary Instrument Manikin Kit

\$30.00 per piece

Trade price from **\$19.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Economy 18-Piece Canine Spay Veterinary Instrument Manikin Kit is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$30.00, falling to \$19.00 at 1 set — about 37% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70002717

Equine Dental McPherson Horse Mouth Speculum & Float Rasp Kit

\$124.98 per piece

Trade price from **\$114.98**
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Equine Dental Full Mouth Speculum
4000 Horse Mouth Gag with Float
Rasps Excellent Quality Forged
Ratchets with Leather Straps Stainless
Steel McPherson (Hausmann) Style
Speculum with Spare Bite Plate

performance. Regular oral and dental checkups can help identify and correct any issues that may have developed over time.

ERGONOMIC COMFORTABLE DESIGN: Professional Equine Dental McPherson Speculum is designed for full mouth examination. It has forged and deep ratchet to ensure comfortable fitting. Its made from high-grade stainless steel with "A" grade leather backings. Spring-loaded ratchets on both sides of the speculum keep the mouth open at the desired level with an extra secure non-slippery locking mechanism.

FLOAT RASP: Equine Horse Dental float rasp adjustable straight angled doubled with Wood handle stainless steel veterinary horse instruments. Highly polished finish for aesthetic and corrosion resistant.

MCPHERSON HORSE MOUTH SPECULUM FEATURES

PREMIUM GRADED STAINLESS STEEL: This stainless steel McPherson Horse Full Mouth Speculum is a necessary piece of equipment for the practice of equine dentistry.

INDISPENSABLE EQUINE DENTISTRY TOOL: In Horses, regular oral examination and care are critical to support overall health and avoid any potential problems. Any oral irregularity may lead to weight loss and poor



WI-27930440

Extendable Stainless Steel Snake Catching Holder Stick Instrument

\$20.00 per piece

Trade price from **\$12.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Extendable Stainless Steel Snake Catching Holder Stick Instrument — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$20.00, falling to \$12.00 at 1 piece — about 40% off the single-unit price.

Larger quantities are priced on request. Order it on its own or as part of a mixed shipment; we do not insist on a full container.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28996408

Farriers Hoof Grooming Tool Kit 6-Piece

\$20.00 per piece

Trade price from **\$18.00**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Farriers Hoof Grooming Tool Kit 6-Piece — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$20.00, falling to \$18.00 at 10 pieces — about 10% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36133667

Feline Spay Neuter Instrument Set Stainless Steel Veterinary

\$180.00 per piece

Trade price from **\$164.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Feline Spay Neuter Instrument Set Stainless Steel Veterinary is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$180.00, falling to \$164.00 at 1 set — about 9% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-42807811

Feline Spay Neuter Surgical Pack Stainless Steel Veterinary

\$180.00 per piece

Trade price from **\$170.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Feline Spay Neuter Surgical Pack Stainless Steel Veterinary — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$180.00, falling to \$170.00 at 1 set — about 6% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46255664

Feline Sterilization Surgery Instrument Set

\$180.00 per piece

Trade price from **\$164.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Feline Sterilization Surgery Instrument Set — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$180.00, falling to \$164.00 at 1 set — about 9% off the single-unit price. Larger quantities are priced on request.

Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46244729

Feline Sterilization Surgical Instrument Set Cat Spay Neuter Veterinary

\$180.00 per piece

Trade price from **\$164.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Feline Sterilization Surgical Instrument Set Cat Spay Neuter Veterinary for veterinary practice and animal surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$180.00, falling to \$164.00 at 1 set — about 9% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25896343

Hoof Farrier Trimming Tools 8 Pcs Set Reusable Breathable

\$30.00 per piece

Trade price from **\$15.00**
Minimum order 10 bags

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

This is our pattern of Hoof Farrier Trimming Tools 8 Pcs Set Reusable Breathable, made for veterinary practice and animal surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: full reusability — clean, autoclave and use again.

matches your protocol rather than ours. Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge.

Listed at \$30.00, falling to \$15.00 at 10 bags — about 50% off the single-unit price. Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Supplied as a 8-piece set. We can add, drop or substitute items so the tray

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25915247

Horse Farrier Hoof Trimmer Tool Kit

\$8.00 per piece

Trade price from \$4.50
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Horse Farrier Hoof Trimmer Tool Kit — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$8.00, falling to \$4.50 at 10 pieces — about 44% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25935362

Horse Hoof Trimming Farrier Tools

\$30.00 per piece

Trade price from \$15.00
Minimum order 10 bags

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Horse Hoof Trimming Farrier Tools is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$30.00, falling to \$15.00 at 10 bags — about 50% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-23955371

Large Animal Veterinary First Aid Kit 6 Pcs

\$220.00 per piece

Trade price from \$150.00
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Large Animal Veterinary First Aid Kit 6 Pcs is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$220.00, falling to \$150.00 at 10 sets — about 32% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28999003

Pro Eco-Friendly Reusable Stainless Steel Complete Farrier

\$40.00 per piece

Trade price from \$25.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Pro Eco-Friendly Reusable Stainless Steel Complete Farrier for veterinary practice and animal surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$40.00, falling to \$25.00 at 10 pieces — about 38% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25865404

Pro Horse Hoof Trimming Farrier Tool

\$30.00 per piece

Trade price from \$15.00
Minimum order 10 bags

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Pro Horse Hoof Trimming Farrier Tool — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Forging, milling, heat treatment, hand finishing and passivation all happen

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$30.00, falling to \$15.00 at 10 bags — about 50% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-70003711

Puppy Tail Banding Elastrator with FREE 100 Pcs Bands

\$18.69 per piece

Trade price from \$17.19
Minimum order 21 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

Elastrator Plastic Castration Tail Docking Dogs Livestock + 100 FREE Bands Veterinary instrument. Ideal for Cattle, Calves, Cows, Dogs, Goats, Lambs, Sheep, and other Livestock. Multi-Function Snap Pliers Easy to operate and convenient to handle.

PUPPY TAIL BANDING ELASTRATOR FEATURES

CASTRATION BANDER WITH STRIPS:

You will get 1 pcs castration bander pliers, 100 pcs castrator rings in red. A large number of castrator rings with castration bander pliers are perfect for young animals, meeting your demands.

HIGH ELASTICITY RUBBER:

Compared with other castration rings on the market, our castration rings are made of premium rubber with good

elasticity and strong flexibility, avoiding breaking. Safer to dock the tail.

PREMIUM CASTRATION SNAP PLIERS:

Grappling hook in 4 directions designed for pulling the castrator rings to the max with ease, the adjustable range: 0.5-5cm. Ergonomics handle for comfortable to grip, adopted to lever principle, one-hand operation, save labor to operate.

SIMPLE TO USE:

Firstly put the rubber band on the pliers and open the pliers; secondly, put the lamb's tail in and put the rubber band between the 3 and 4 vertebrae of the tail, safe to use (you can feel the joint seams by pinching the root of the lamb's tail). Finally, operate the pliers, work done!

GOOD HELPER FOR FARMER:





WI-29429154

Rabbit Rat Rodent Small Animal Cheek Teeth Molar Incisor

\$15.00 per piece

Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Rabbit Rat Rodent Small Animal Cheek Teeth Molar Incisor — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: interlocking teeth for a positive hold on tough tissue.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$15.00. Minimum order 10 sets, and the price moves with quantity.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34073285

Reusable Eco-Friendly Veterinary Orthopedic Instrument Kit Complete System Dog

\$600.00 per piece

Trade price from \$597.00
Minimum order 5 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

Reusable Eco-Friendly Veterinary Orthopedic Instrument Kit Complete System Dog — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Supplied as a set — the contents can be adjusted to your list. We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$600.00, falling to \$597.00 at 5 sets — about 1% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-34051245****Reusable Veterinary Orthopedic Instrument Kit Complete System Small Animal****\$600.00** per pieceTrade price from **\$597.00**
Minimum order 5 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Reusable Veterinary Orthopedic Instrument Kit Complete System Small Animal is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The name tells you what you are getting: full reusability — clean, autoclave and use again.

Supplied as a set — the contents can be adjusted to your list. Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$600.00, falling to \$597.00 at 5 sets — about 1% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

**WI-27932199****Shorty Hog Holder Instrument Set Veterinary Pig Restraining****\$18.00** per pieceTrade price from **\$12.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Shorty Hog Holder Instrument Set Veterinary Pig Restraining — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$18.00, falling to \$12.00 at 1 piece — about 33% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28939165

Small Animal Veterinary Dental Instrument Set

\$35.00 per piece

Trade price from \$33.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Small Animal Veterinary Dental Instrument Set — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$35.00, falling to \$33.00 at 1 set — about 6% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-24082237

Stainless Steel Anatomical Model Orthopedics Veterinary

\$30.00 per piece

Trade price from \$19.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Anatomical Model Orthopedics Veterinary is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$30.00, falling to \$19.00 at 1 set — about 37% off the single-unit price.

Larger quantities are priced on request. Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.

Forging, milling, heat treatment, hand finishing and passivation all happen

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-46260656

Stainless Steel Feline Spay Neuter Instrument Set for Cat

\$180.00 per piece

Trade price from **\$164.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Feline Spay Neuter Instrument Set for Cat — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

The steel is heat treated and passivated for corrosion resistance, then each

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$180.00, falling to \$164.00 at 1 set — about 9% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-25892288

Stainless Steel Hoof Clincher Instrument Set for Horseshoe

\$8.00 per piece

Trade price from **\$4.50**
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Hoof Clincher Instrument Set for Horseshoe is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$8.00, falling to \$4.50 at 10 pieces — about 44% off the single-unit price. Larger quantities are priced on request.

We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-36641610

Stainless Steel Orthopedic Surgical Instrument Set Blue Silicone

\$55.00 per piece

Trade price from \$45.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Stainless Steel Orthopedic Surgical Instrument Set Blue Silicone for veterinary practice and animal surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$55.00, falling to \$45.00 at 1 piece — about 18% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-27567272

Stainless Steel Surgical Skin Hooks 15cm Instrument

\$7.00 per piece

Trade price from \$4.00
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Stainless Steel Surgical Skin Hooks 15cm Instrument is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. The 15cm size shown is our standard; other sizes in the same pattern are made to order.

Each unit is inspected by hand for alignment and action before packing, and we can photograph your actual pieces before they ship. Listed at \$7.00, falling to \$4.00 at 1 piece — about 43% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-32054322

Suture Practice Kit for Medical Veterinary Students

\$7.00 per piece

Trade price from \$4.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking


ABOUT THIS INSTRUMENT

This is our pattern of Suture Practice Kit for Medical Veterinary Students, made for veterinary practice and animal surgery. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$7.00, falling to \$4.00 at 10 pieces — about 43% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28543634

Suture Practice Kit Medical Veterinary Students

\$12.00 per piece

Trade price from \$10.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT


Suture Practice Kit Medical Veterinary Students — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$12.00, falling to \$10.00 at 1 set — about 17% off the single-unit price.

Larger quantities are priced on request. We ship door to door worldwide, and can quote DDP if you would rather not handle the import yourself.

The steel is heat treated and passivated for corrosion resistance, then each

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-24092177

Veterinary Cardiovascular Surgery Anatomical Model Kit-Stainless Steel Medical

\$30.00 per piece

Trade price from \$19.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Veterinary Cardiovascular Surgery Anatomical Model Kit-Stainless Steel Medical is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

Forging, milling, heat treatment, hand finishing and passivation all happen under one roof, which is why we can change a dimension or a finish without a tooling charge. Listed at \$30.00, falling to \$19.00 at 1 set — about 37% off the single-unit price. Larger quantities are priced on request.

Order it on its own or as part of a mixed shipment; we do not insist on a full container.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-28634030

Veterinary Dissecting Kit Anatomical Model 14 Pcs

\$50.00 per piece

Trade price from \$35.00
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Veterinary Dissecting Kit Anatomical Model 14 Pcs is an instrument set we forge, grind and finish in our own Sialkot workshop. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

instrument is checked by hand for jaw meet, smooth action and finish. We state the grade we used on your quotation. Listed at \$50.00, falling to \$35.00 at 1 set — about 30% off the single-unit price.

Larger quantities are priced on request. Small trial orders are welcome — most buyers start with a sample or a single carton before committing.

The steel is heat treated and passivated for corrosion resistance, then each

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-34691746

Veterinary Orthopedic Trauma PUSM Elastic Nail Instrument

\$200.00 per piece

Trade price from **\$120.00**
Minimum order 1 set

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Veterinary Orthopedic Trauma PUSM Elastic Nail Instrument for veterinary practice and animal surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Every piece is hand-finished and inspected before it leaves the bench: joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$200.00, falling to \$120.00 at 1 set — about 40% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking

WI-27906596

Veterinary Pig Holder Instrument Set

\$15.00 per piece

Trade price from **\$8.00**
Minimum order 1 piece

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

We manufacture Veterinary Pig Holder Instrument Set for veterinary practice and animal surgery — made from bar stock to finished instrument in house. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$15.00, falling to \$8.00 at 1 piece — about 47% off the single-unit price. Larger quantities are priced on request.

Tell us the quantity you actually need and we will price it, including quantities below the listed minimum.



Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-23908573

Veterinary Surgery Horse Spay Pack with Instrument Box

\$50.00 per piece

Trade price from \$35.00
Minimum order 10 sets

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

Veterinary Surgery Horse Spay Pack with Instrument Box — one of 4,770 instruments we manufacture ourselves rather than buy in. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table.

Every piece is hand-finished and inspected before it leaves the bench:

joint action, tip alignment and surface finish are checked on each unit, not on a sample. Listed at \$50.00, falling to \$35.00 at 10 sets — about 30% off the single-unit price.

Larger quantities are priced on request. Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WI-29236743

ZebraFish Veterinary Dissecting Kit Surgical Instruments Set

\$9.00 per piece

Trade price from \$5.00
Minimum order 10 pieces

- ✓ Surgical-grade stainless steel
- ✓ Reusable, autoclavable
- ✓ OEM / private-label marking

ABOUT THIS INSTRUMENT

ZebraFish Veterinary Dissecting Kit Surgical Instruments Set is part of our veterinary instruments range, produced start to finish at our Sialkot facility. Instruments of this type cover a complete procedure in one tray, so nothing is missing at the table. Supplied as a set — the contents can be adjusted to your list.

We work the instrument by hand at the finishing stage — box locks, hinges and tips are set individually, which is what separates a good instrument from a cheap one. Listed at \$9.00, falling to \$5.00 at 10 pieces — about 44% off the single-unit price. Larger quantities are priced on request.

Samples go out on request so you can check the steel and the finish before you place a real order.

Material	Stainless steel — we state the grade on your quotation
Reusable	Yes, autoclavable
Finish	Satin, mirror or black — your choice at order
Marking	Your logo or part number can be applied — ask us what is possible for this item
Lead time	Typically 15–30 days depending on quantity and marking



WAQAR INTERNATIONAL

PRECISION SURGICAL INSTRUMENTS

Request a quotation

Send us product codes from this catalogue, quantities and your country.
We reply the same working day on WhatsApp.

Online shop	shop.waqarintl.com
Email	sales@waqarintl.com
WhatsApp /	
Phone	+92 307 6173627
Works	Roras Road, Opposite Anwar Khawaja Industries, Sialkot 51310, Pakistan

Prices in USD, ex-works Sialkot, subject to change. OEM marking, private label and custom packaging on request.